# T-CLASSIC T-FLAP

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**ENG** T-Classic / T-Flat, is a Mercury Free thermometer, useful for clinical measurement of body temperature when positioned in axillary, oral or rectal position.

## Use and Storage conditions

Measurement principle:	Thermal expansion of a calibrated liquid alloy
Data Display:	Displayed mechanically on a calibrated graduated scale.
Measurement range:	35 - 42°C
Minimum scale interval:	0.10°C
Measurement accuracy (at 37°C and 41°C):	+0.10 to -0.15°C
Environmental operational limits:	+20°C - +30°C
Shipping/storage environmental limits:	0°C - 42°C

- The ideal environment temperature for correct measurement should be between 20 and 30°C. If the thermometer has been kept in an environment with a temperature below 20°C, before using it, it is necessary to keep it for 1-2 hours in an environment with a temperature that is between 20 and 30°C or kept for approximately ten minutes in a room with a temperature between 30 and 35°C.

- If environment temperature exceeds 35°C, immerse the tip of the thermometer in cold water for approximately 10-20 seconds and then shake it in order to lower the reference column and then position it immediately in the site used for measurement.

- Temperature must be below 42 °C when stored and transported in order to avoid irreparable damage.

- After use, always immediately lower the liquid metal column, then correctly sanitise/clean it and place it back in the protective plastic cover, in order to avoid any damage or breakage that may place it in a dangerous situation.

**IMPORTANT:** To lower the liquid metal column in the T-Classic model, carefully tighten the final part of the thermometer between the thumb and index finger and, paying attention not to crush it and not hitting it against anything, shake the thermometer downwards multiple times, until the liquid column does not go below 36°C. The same procedure can also be used with the T-Flat model that is however equipped with a plastic package that was specifically designed to safely and quickly lower the liquid metal column. If you have purchased T-Flat, it is sufficient to insert the device into its plastic container, tighten the cover well and then tightly hold the plastic bag located at its end between thumb and index finger, allowing the thermometer to swing quickly as if it was a pendulum, until the temperature shown falls below 36°C.

## Measurement method

It is important to remember that the body temperature reading depends on where it is measured.

## Axillary (positioning it under the arm)

If the thermometer is in axillary position, a surface temperature is obtained that can oscillate approximately between 0.5°C and 1.5°C compared to rectal temperatures in adults. Once the thermometer has been positioned, the patient must keep his arm tight against his body, without applying pressure but ensuring contact between the bulb and the axilla. The average measurement time for this method is approximately 5 minutes. It is worthwhile to note however that it is not possible to obtain an exact reading if the axilla are cooled down. In this case, it is recommended to extend measurement time by approximately 2 minutes in order to obtain the most precise as possible reading that corresponds as much as possible to body temperature.

## Oral (positioned in the mouth)

In the mouth there are various hot areas. In general, oral temperature is 0,3° C - 0,8° C lower than rectal temperature. In order to guarantee a reading as precise as possible, position the tip of the thermometer on the left or right of the root of the tongue. The tip of the thermometer must be in constant contact with tissue while reading and must be positioned under the tongue. Keep the mouth closed while reading and breath regularly using nose. Do not eat or drink anything before measurement. Average measurement time is approximately 3 minutes

## Rectal (positioning it in the rectum)

This method is the most precise from a medical point of view, since it is the closest to body temperature. The tip of the thermometer is carefully inserted into the rectum a maximum of 0.5 cm. The average measurement time is approximately 3 minutes.

#### Before proceeding with a measurement, always remember that:

1. The liquid metal column must always be brought under 36°C before starting a measurement.
2. Before each use it is necessary to correctly clean / disinfect/sterilise the thermometer (depending on the type of application, the environment where it is used and single or promiscuous use). Sterilise in medical use alcohol or in plastic sterilisation tubes.
3. In case of rectal measurement, it is recommended to apply a specific lubricant before insertion.
4. It is best to avoid physical activity, as well as taking a bath, eating or smoking at least 30 minutes before carrying out the measurement.
5. Avoid axillary measurement after applying refreshing deodorants or deodorants with particles of silver or other. Their application may interfere in part with a correct measurement.

#### IMPORTANT PRECAUTIONS

1. The use of clinical thermometers for oral temperature measurement is extremely inadvisable in children that are under 10 years of age;
2. Patients to whom the thermometer is applied orally, must pay careful attention to avoid chewing, talking, playing, running or moving until the thermometer is removed from the mouth;
3. Rectal use of the thermometer on children must be applied in a silent environment, under constant control and sustaining the thermometer with a hand in order to keep it from penetrating too deep or sliding outdoors.
4. Always verify, at least once a year, structural integrity of the device.

#### Safety precautions

- This thermometer must be used exclusively to read bodily temperature (fever).
- This thermometer was not designed for special applications (for example thermometers for premature babies, ovulation thermometers, etc.);
- The device must be used following the information provided in the attached documents.
- This device must only be used for its destined use and according to instructions supplied in this manual. A different use must be considered incorrect and dangerous.
- Manufacturer and Distributor cannot be considered liable for damages due to an improper use or storage.
- In case of malfunctioning, only contact the distributor's after-sales service or an after-sales service centre authorised by the Manufacturer.
- The device must not come into contact with hot water.
- Do not expose it to high temperatures or direct sunlight.
- Do not let the thermometer fall. It is not impact resistant.
- Do not clean it using solvents, petrol or benzene. Only clean using water or disinfectant.
- The thermometer body is made of glass, a very fragile material that is potentially dangerous in case of breakage. In the same way it is also dangerous if swallowed, the metal liquid contained inside may come out if the glass body were to break. For this reason, never leave the thermometer in children's hands.
- Persistent fever, especially in children, must be treated by a physician. In presence of fever, always consult your physician.

 USING THIS THERMOMETER DOES NOT REPLACE YOUR PHYSICIAN'S ADVICE

#### DISPOSAL

The body of this device is glass and inside it contains an alloy of liquid metals made of gallium (68% - 0.544g), indium (20% - 0.16g) and tin (12% - 0.096g) and considering its fragility, in order to avoid dispersing the materials, this clinical thermometer must always be handled with extreme care. In particular, one must pay attention to never let the instrument fall and never subject it to treatments or force that may cause damage or breakage. In case of breakage, always proceed by picking up glass and liquid metal that may have come out, using a soft cloth soaked with water and soap. Carefully pick up any residue and throw away even the cloth used for this operation.

#### Technical Specification

Model References :	TR.100200 – T-CLASSIC TR.100300 – T-FLAP
Size:	Bulb length= 14 ± 3mm Product length = 125±5mm H/W = 12.1/8.8 mm
Weight:	approx 7.2 g without plastic cover
Degree of protection against liquids: :	Waterproof

#### Legal requirements and guidelines


This product conforms to the European Directive on Medical Devices 93/42/EEC. The device conforms to the requirements foreseen by standard EN12470-1: 2000 + A1: 2009, "Clinical Thermometers - Part 1: Metallic liquid in glass thermometer with maximum degree".

#### Calibration control

This thermometer does not require any additional periodical or long term calibration.

#### Warranty

This product is covered by a 2 year warranty starting from the date the purchase document was issued. Damages caused by an incorrect use are not covered by warranty. The warranty is void in case of tampering, if the defect is caused by an improper use or if the damage is not a result of manufacturing defects (accidental fall, careless transport, etc.).

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CA-MI S.r.l. - Via Ugo La Malfa nr.13 - Frazione Pilastrò