# STORZ

## KARL STORZ — ENDOSKOPE

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TL400 Kaltlicht-Fontäne power LED rubina

### INSTRUCTION MANUAL

TL400 Cold Light Fountain power LED rubina

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the cold light fountain power LED rubina. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG

Utilization and reproduction by third parties require express authorization by KARL STORZ SE & Co. KG.

All rights reserved.

**Rx only** Federal (USA) law restricts this device to sale by or on the order of a physician.

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

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

**2. Controls, displays,  
connectors, and  
their uses**

- ① Standby button
- ② Touch screen
- ③ Light outlet, can be unscrewed
- ④ Footswitch connector
- ⑤ Potential equalization connector
- ⑥ External sync interface
- ⑦ Power cord socket
- ⑧ Line fuses (fuse holder)
- ⑨ USB connector
- ⑩ KS HIVE connection (only in combination with KARL STORZ devices)

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ICG = Indocyaningrün  
IR = Infrarot  
NIR = Nahinfrarot

**Abbreviations used**

ICG = IndoCyanineGreen  
IR = Infrared  
NIR = near infrared

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
**3. Symbols employed**

**3.1 Symbols on the device**

	Follow instructions for use
	Standby/ON/OFF
	Potential equalization connector
	Type CF applied part according to IEC 60601-1
	Alternating current
	Fuse
	Ethernet
	USB
	Connection socket for footswitch
	Electronic information product pollution control (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Hersteller
	Hot surface

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**3.2 Label and packaging symbols**

	Manufacturer
	Catalogue number
	Serial number
	Number of products in the product packaging
	Keep dry
	Fragile, handle with care
	Consult instructions for use
	CE mark

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


	<b>White light</b>
①	Mode display for white light and NIR applications
②	Footswitch (only shown if the footswitch is connected)
③	Menu (only shown if the light source is stand-by)
④	ON/OFF (white = Off) ON/OFF (orange = On)
⑤	Manual light intensity regulation
⑥	Light intensity (only shown if manual light intensity regulation is active)
⑦	Automatic light intensity regulation (can only be used if the light source is connected with the IMAGE1 S via KS HIVE).
⑧	Lower light intensity manually
⑨	Increase light intensity manually
⑩	Slide control for manual light intensity regulation

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## 4. General information

### 4.1 Device description

The cold light fountain power LED rubina is a high-performance cold light source for diagnostic and therapeutic interventions in endoscopy, microscopy and open surgery.

The device offers numerous functions and features. These include:

- White light and NIR: High and constant light intensity with low current consumption thanks to LED technology.
- NIR: Fluorescent modes to aid the user by showing only the fluorescent signal or by overlapping the white light display with the fluorescent information in real time.
- Comfortable working environment: Innovative cooling concept for minimal noise emission.
- Durable: No costs or effort required to change the lamp if used properly.
- Clear: Direct operation via a color screen with touch-sensitive surface (touchscreen)
- Integration: The KS HIVE function enables the light source to be controlled remotely and the system status information to be retrieved, while providing simple configuration during initialization.
- Up-to-date: Software updates can be installed easily and conveniently via an interface.

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## 5. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 5.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**NOTE:** A Note indicates special information about operating the product, or clarifies important information.

**NOTE:** Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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**WARNING:** Read this instruction manual thoroughly before putting the device into operation. Read the section on safety instructions carefully to avoid putting your patients, personnel, or yourself at risk.

**WARNING:** The electrical installation of the operating room in which the device is connected and operated must comply with the applicable IEC standards.

**WARNING:** Keep device out of reach of patients.

**WARNING:** For safety reasons, do not simultaneously touch the device output sockets and the patient.

**WARNING:** The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**WARNING:** Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. Clause 16 of IEC 60601-1).

**WARNING:** Additional devices which are connected to electrical medical devices must comply with the relevant IEC or ISO standards (e.g., IEC 60950 for data processing devices). Furthermore, all configurations must comply with the requirements for medical electrical systems (see Clause 16 of the 3rd Ed. of the IEC 60601-1). Anybody connecting additional devices to electrical medical devices is configuring a medical system, and is therefore responsible for ensuring that the system complies with the standard requirements for medical electrical systems. It should be noted that the local legal regulations have priority over the above-mentioned requirements. If in doubt, contact your representative on site or the technical customer service.

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**WARNING:** Grounding reliability can only be achieved when the unit is connected to a properly installed 'Hospital Only' or 'Hospital Grade' receptacle (i.e. approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**WARNING:** To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

**WARNING:** At the end of its useful operating life, dispose of the unit as electronic scrap.

**WARNING:** Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

**WARNING:** Ensure sufficient air circulation for the device.

**WARNING:** It is not permitted to press more than one touch point simultaneously.

**WARNING:** Only operate the device with the voltage stated on the manufacturer's identification plate.

**WARNING:** DO NOT use the device in the presence of flammable anesthetics.

**WARNING:** Servicing may only be performed by specialists authorized by KARL STORZ.

**WARNING:** Test this device prior to each surgical procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

**WARNING:** The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

**WARNING:** Do not use the device if the input controls or the displays are defective.

**WARNING:** The device must not be used in oxygenated environments.

**WARNING:** The device may be subject to sudden failures. Keep a spare device at hand or alternate procedures in mind.

**WARNING:** Applied parts of devices used in combination with this device must be of type BF or CF.

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**WARNING:** Patient leakage currents may be additive when energized endoscopes are used with energized endotherapy devices. This is particularly important if a CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

**WARNING:** The compatibility of this device with any accessories and/or energized endotherapy devices should be checked in the individual instructions for use.

**WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

**WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

**WARNING:** To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

**WARNING:** Always use the lowest possible settings for your intervention.

**WARNING:** Danger of glare! Never look into the light outlet when the light system is switched on.

**WARNING:** The light source is transmitting high energy light leading to high temperatures at the light window of the endoscope. Tissue damage is possible.

**WARNING:** Insufficient ventilation of the unit can cause an internal build-up of heat which leads to a safety shut-off.

**WARNING:** The device may only be operated within specified environmental conditions.

**WARNING:** The correct configuration of the device must be checked after each software update.

**WARNING:** Certain devices or accessories can cause danger in lower power settings. For example, the risk of gas embolism in argon assisted coagulation rises if the HF power is insufficient for the fast creation of an impenetrable eschar layer on the target tissue.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

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## 5.2 Intended use

It is only allowed to use the power LED rubina light source in combination with approved accessories from KARL STORZ or with equipment that is approved as safety-related uncritical.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

**WARNING:** No modification of this equipment is allowed.

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### 5.2.1 Indication

NIR/ICG fluorescent light sources are intended for generating light in endoscopic and microscopic diagnostic examinations and in surgical procedures.

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#### Indication (applies only to USA)

Light sources are used for illumination during endoscopic, open or microscopic procedures. Additionally they are used for excitation during near infrared fluorescence imaging.

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### 5.2.2 Contraindication

When using NIR/ICG fluorescent light sources, the medical contraindications that are listed for the suitable marker substance that may be used apply, and use is contraindicated if endoscopic or microscopic examination techniques contradict the latest state of scientific knowledge.

NIR/ICG fluorescent light sources are contraindicated for ophthalmologic interventions. For other interventions, the responsible physician must decide whether the prescribed application is admissible based on the general condition of the patient.

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### 5.3 Intended purpose

NIR/ICG fluorescent light sources are used for white light illumination in diagnostic and therapeutic interventions in endoscopy, microscopy and open surgery.

Additionally, they are used for stimulation during fluorescence imaging in the near infrared spectral range.

### 5.4 User qualification

Only qualified doctors and medical assistant staffs are allowed to use the cold light fountain power LED rubina light source after training. Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ sales representative is available to provide training on your premises and to give you information about further training options.

### 5.5 Patient profile

The use of the power LED rubina is not limited to a certain patient profile (sex, age, weight, etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

### 5.6 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means

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## 5.7 Intended conditions of use

### 5.7.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operates properly before using it. During treatment using the power LED rubina cold light fountain, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

**WARNING:** Any treatment may only be performed if visual observation of the action of the device is ensured.

**WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

### 5.7.2 Other intended conditions

Frequency of use: one or more times a day

Length of use: from a few minutes to several hours a day.

Place of installation: Positioning on a level, vibration-free surface

Mobility: can be moved if positioned on a cart

Combination: can be used on the patient at the same time as other equipment required for the operations

Control: can be controlled via KS HIVE

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**5.8 User position**

For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$  at a distance of approx. 30 – 70 cm to the front panel.

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**5.9 Safety precautions at the site of installation**

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones.

**WARNING:** The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

**WARNING:** The device should be placed on a flat surface or a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.

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**5.10 Safety features**

The cold light fountain power LED rubina is equipped with the following safety features:

**5.10.1 Self-test**

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

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## 6. Installation and operating instructions

### 6.1 Unpacking the equipment

Carefully unpack the cold light fountain power LED rubina and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

### 6.2 Basic equipment

KARL STORZ Set No.: TL400

- 1 power LED rubina cold light fountain
- 1 400A power cord, length 300 cm
- 1 Instruction manual 96286010D
- 1 Patch cable CAT6a, 2.0 m, UL, WO10275
- 1 Sync connecting cable, 100 cm, TL006
- 1 OPAL1 NIR/ICG demonstration card 96240726

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
### 6.3 Installation

**i NOTE:** The cold light fountain power LED rubina including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

**! WARNING:** The device is not intended for use in potentially explosive zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

**! WARNING:** Risk of injury! The touchscreen is made of glass. This glass can break if the device is dropped on a hard surface or receives a substantial impact. Place the device away from the patient (1 meter)! This minimizes the risk of fragments coming into contact with the patient if the touchscreen is damaged during use. Make sure that no pieces of glass are left inside the patient!

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If the glass breaks, do not touch the glass parts of the device or attempt to remove the broken glass from the device. Stop using the device until the glass is replaced by qualified service personnel.

Set the device on a horizontal, flat surface (Fig. 1, Page 12).

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**6.3.1 Connecting the ground line**

**! WARNUNG:** The device is equipped with a connector for attaching a ground line. Carry out potential equalization as per the applicable national regulations.

The device is equipped with a connector ⑤ for attaching a ground line.

The device's ground line should be installed by a qualified electrician (Fig. 2).

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**6.3.2 Connecting the power cord**

**! WARNUNG:** Only operate the device with the voltage stated on the manufacturer's identification plate.

**! WARNUNG:** Only connect/disconnect the power plug to/from the power supply outside of areas subject to explosion hazards.

**! WARNUNG:** The device must be set up in such a way that the power cord can be simply pulled out at any time.

Connect power cord and insert the power plug fully into the power cord receptacle ⑦ (Fig. 3).

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**6.3.3 KS HIVE**

**i NOTE:** To prevent the KS HIVE connecting cable from being pulled out accidentally, the KS HIVE connector possesses a protection device.

Insert the connector into the KS HIVE socket ⑩ of the power LED rubina (Fig. 4) until you hear a click—the connector is now correctly mechanically engaged. Connect the other end of the cable with a KS HIVE control device.

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6

**6.3.4 Connecting sync connecting cable**

**i NOTE:** To ensure full functionality of the TIPCAM1 Rubina with the power LED rubina, an IMAGE1 S 4U-LINK (TC304) is required.

Insert the connector into the Ext Sync socket ⑥ of the power LED rubina (Fig. 5). Insert the other end of the cable into the Ext Sync socket of the IMAGE1 S 4U-LINK (TC304) (Fig. 6).

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**i NOTE:** Application parts for this device are light cables and endoscopes (see the following Sections 6.3.5 and 6.3.6).

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7

**6.3.5 Connecting light cable**

Plug the light cable into the socket ③ as far as it will go (Fig. 6). Only handle the light cable using the handle piece. Never pull on the cable.

**i NOTE:** The light outlet ③ is provided with an antiglare flap, which does not allow any direct outlet of light. Light is only emitted when the light cable is inserted.

**! WARNING:** Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up (can exceed 41 °C). This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

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**6.3.6 Connecting the light cable to the endoscope**

Connect the light cable to the endoscope (Fig. 8).

**NOTE:** The use of original KARL STORZ light cables only is recommended. Light cables from other manufacturers may not be optimally adjusted for light transmission.

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**6.3.7 Connecting the footswitch**

Connect the cable from the footswitch to its mating connector ④.

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**General applications**

The power LED rubina is suitable for use in a wide range of applications. The correct light intensity depends on the particular application, the instrument being used, and the connected light cable. If the automatic controller for the light source (see section 6.4.2) isn't used, the intensity should always be kept as low as possible.

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**6.4 Commissioning the power LED rubina light source**

Switch the device on at the Standby button ① (Fig. 9).

**NOTE:** Do not touch the touchscreen when starting the device.

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Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo (Fig. 10).



The device is now running a self test. Once the test has successfully completed, an acoustic availability signal sounds and the light source is ready for use. If the self test fails for critical components, the device goes into 'Malfunction' mode (Safe State) and an acoustic information signal sounds (Fig. 11). Verify that the device is properly connected and that the preparations have been made according to the instructions.

For further information on troubleshooting, please refer to the section 'Technical description – Troubleshooting'.

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After switching on, the power LED rubina is in manual white light mode (standby mode with the last set light intensity) (Fig. 12).


**NOTE:** The icon for calling up the setting menu is only shown if the light is switched off (the symbol on the ON/OFF button is white).

Press button **3** to switch to the Settings menu. Further details can be found in the section 'Settings'.

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
#### 6.4.1 Switching the light on and off

Press the ON/OFF button **4** to switch the light on and off. The color of the symbol on the ON/OFF button shows the current status of the light (orange = on, white = off) (Fig. 13).

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
**6.4.2 Using the automatic light source control**

If the power LED rubina is connected to the camera control unit IMAGE1 S via KS HIVE cable, automatic light intensity regulation can be used.

Prerequisite: IMAGE1 S with software release version 4.0 as a minimum

Press the button **5** and select 'auto' **7** (Fig. 14). Press the ON/OFF button **4** to set the light (the symbol on the ON/OFF button **4** is orange – Fig. 15).

The light intensity is then controlled automatically via the IMAGE1 S, ensuring optimal illumination in every situation.

**i NOTE:** The automatic light regulation is deactivated as soon as 'manual' **5** is pressed during operation. For this press 'auto' and then select 'manual' which appears to the right of the 'auto' field (opposite of Fig. 14). The automatic light intensity regulation can be switched on again via the camera control unit. Further information can be found in the IMAGE1 S instruction manual.

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
**6.4.3 Using manual light intensity regulation**

The power LED rubina is in manual mode after being switched on, and the light is switched off (standby mode).

The brightness can be modified using the buttons **8** and **9** both when the light is switched on and off (Fig. 16).

Press the button **9** to increase the brightness.

Press the button **8** to reduce the brightness.

The light intensity can also be reduced (push the slide control to the left) or increased (push slide control to the right) directly with the slide control **10**. The light intensity can be set between 5% and 100%.

The current set brightness is shown in % on the display **6**.

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#### 6.4.4 Changing the light mode

For decommissioning, switch off the unit at the power switch. The light mode (white light and NIR applications) can only be changed using the buttons of a compatible camera head connected to the IMAGE1 S, via the IMAGE1 S (see IMAGE1 S instruction manual) or with a footswitch connected to the light source.

The selected mode is shown in the mode display for white light and NIR applications ❶ (see, for example, Fig. 16).

**❶ NOTE:** It is not intended to change the light mode via the power LED rubina screen.

The following white light and NIR applications are available:

- **White light:** Only the white light information is shown on the IMAGE1 S.
- **Monochromatic:** Only the NIR information is shown on the IMAGE1 S.
- **Overlay:** A white light display is overlaid with the fluorescent information and shown on the IMAGE1 S.
- **Intensity Map:** A white light display with the fluorescent information (fluorescent information in false colors depending on the fluorescent intensity) is shown overlaid on the IMAGE1 S.

**❶ NOTE:** When using a H3-Z FI camera head (TH102) or with a missing KS HIVE connection, 'Fluorescence' is shown as the NIR application in the mode display for white light and NIR applications ❶.

#### 6.4.5 Switching off the device (decommissioning)

Switch off the device at the Standby button.

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### 6.5 Settings


Click on the menu icon (3) to call up the Settings menu (Fig. 17).

**NOTE:** The setting menu symbol is only visible and the setting menu can only be called up if the light source is in stand-by (symbol on the ON/OFF button (4) is white).

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
In the Settings menu, you can select various setting and management fields (Fig. 18).

- Service (only authorized service technicians)
- Language
- Event protocol
- Device Info
- 'Audio' (on/off)

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#### 6.5.1 Service


The 'Service' window is opened by pressing the 'Service' button in the Settings menu (Fig. 18).

**NOTE:** This area is password protected (Fig. 19).

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#### 6.5.2 Language

The 'Language' window is opened by pressing the 'Language' button in the Settings menu (Fig. 18).


The following screen languages are available: Danish, English, Finnish, French, German, Italian, Polish, Portuguese, Russian, Spanish, Korean, Japanese, Chinese, Swedish, Dutch and Norwegian (displayed in the appropriate language). Select a new language if required, and then confirm.

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### 6.5.3 Event protocol

The corresponding window is opened by pressing the 'Event protocol' button in the Settings menu (Fig. 18).

The Event protocol contains the saved information and alarms.

Each entry occupies one row and consists of date, time and alarm/information ID.

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
### 6.5.4 Device Info


The 'Device Info' window is opened with device-specific data by pressing the 'Device Info' button in the Settings menu (Fig. 18).

In this window, you will find information about your device (e.g., the serial number, software version, and operating hours) (Fig. 22).

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### 6.5.5 Audio settings

The acoustic signals can be switched on or off (Audio on/Audio off) directly in the Settings menu (Fig. 23).

When the crossed-out bell is orange, the acoustic signals are switched off.

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**7. Maintenance**

**7.1 Fuse replacement**

**!** **WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Switch off the device and remove the power plug from the electrical outlet (Fig. 24).

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25

Remove the fuse holder ⑧ with a screwdriver or other suitable tool (Fig. 25).

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26

**!** **WARNING:** Only use fuses of the correct rating. Replace defective fuses and screw the line fuse holder back on.

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27

Refit the fuse holder ⑧.  
Reconnect the power cord.  
Test the device for proper operation (Fig. 27).

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**7.2 Reprocessing**

**7.2.1 General warnings for medical devices**

**WARNING:** Risk of infection: Incorrectly prepared medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**WARNING:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**WARNING:** Country-specific laws and regulations must be observed.

**7.2.2 Wipe-down disinfection of device**

**WARNING:** Always unplug the unit and disconnect the endoscope before carrying out any maintenance or cleaning.

**CAUTION:** Make absolutely sure that no liquid can get into the device housing.

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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#### 7.3.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

#### 7.3.2 **Safety checks/repeat inspections as per IEC 62353**

**WARNING:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

#### **Visual inspection**

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

#### **Electric measurements**

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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**Test for proper functioning**

Carry out the functional check as per the instruction manual (see section 6.4).

**Documentation**

The safety check/repeat inspection and results must be documented.

**WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**NOTE:** For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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**7.4 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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**7.5 Disposal**

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

**WARNING:** At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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### **7.6 Repairing devices**

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-Off Gewerbepark 83  
78579 Neuhausen, Germany

Service hotline: +49 7461 708-980  
E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

### **7.7 Infection prevention**

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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**7.8 Limitation of liability**

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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**7.9 Manufacturer's warranty**

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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**8.1 Information signals**

**8.1.1 Configuration**

All of the information signals created, with the exception of those which display the cause of an inoperable device ('Malfunction', see e.g. page 15 Self test), are not self-retaining. This means that the signal is only issued whilst the signal conditions are present.

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration, with which the described information signals also comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by a signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If numerous signal conditions with the same priority occur, only the condition most recently detected is signaled.

**8.1.2 Visual signals**

Information signals are displayed without flashing (Fig. 28).

'Light blue field'	11 Text display of an information signal
	12 Acknowledge information signal
	13 Switch off acoustic signals

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 10 seconds.

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
**8.1.3 Acoustic signals**


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. But at least one complete signal sequence is played.

The acoustic signal can be switched on or off in the Settings menu.

But the acoustic signal can only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

**8.1.4 Information signals (general)**

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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


**8.1.5 Availability signal**

On completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

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Complete failure of the device.

**Possible cause**

1. Power supply failure.
2. Defective power fuse.
3. Power cord connector is not properly connected to device socket.

**Recommended remedy**

1. Check that there is electricity to the wall outlet.
2. Change fuses as described in the instruction manual. Use only fuses of the correct rating.
3. Push power cord connector firmly into device socket.

**Issue**

No light is emitted.

**Possible cause**

1. Electronics are defective.
2. Overheating of the power LED rubina (LED has switched off due to overheating of the unit caused by covered air vent).

**Recommended remedy**

1. Notify Service.
2. Remove objects from in front of the air vents. Switch off the power LED rubina, let it cool down (approx. 10–15 minutes) and switch it back on.  
Ensure that it is sufficiently ventilated.

**Issue**

No light emitted, Standby button lights up (on).

**Possible cause**

1. Strong impacts or vibrations during operation.
2. Power supply or LED is defective.

**Recommended remedy**

1. Turn the unit off and then on again.
2. Have the unit repaired.

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<p><b>Issue</b> Too little light.</p> <p><b>Possible cause</b></p> <ol style="list-style-type: none"> <li>1. Ends of light cable and/or endoscope are soiled.</li> <li>2. Light cable or endoscope is defective.</li> <li>3. LED service life has been exceeded.</li> </ol> <p><b>Recommended remedy</b></p> <ol style="list-style-type: none"> <li>1. Clean ends of light cables, and light entrances and exits on the endoscope.</li> <li>2. Replace light cable or endoscope.</li> <li>3. Have the unit repaired.</li> </ol>
<p><b>Issue</b> Light intensity regulation shows 'errors'.</p> <p><b>Possible cause</b> Internal error in input processing.</p> <p><b>Recommended remedy</b> Press the Standby button once, in order to restore 'normal performance'.</p>
<p><b>Problem</b> Touchscreen does not react.</p> <p><b>Possible cause</b></p> <ol style="list-style-type: none"> <li>1. Touchscreen was actuated during the start process.</li> <li>2. The touchscreen is very wet.</li> </ol> <p><b>Recommended remedy</b></p> <ol style="list-style-type: none"> <li>1. Turning the device off and on again.</li> <li>2. Remove any moisture on the front of the device with a clean, dry cloth.</li> </ol>
<p><b>Issue</b> Device certificate expired.</p> <p><b>Possible cause</b> The device certificate is no longer up-to-date or invalid.</p> <p><b>Recommended remedy</b> Contact KARL STORZ Customer Service.</p>

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**8.3 Technical data**

power LED rubina	TL400
Line voltage	100...240 V~
Line frequency	50/60 Hz
Power consumption	220 VA
Line fuse	2 x T4.0 AH 250V
Dimensions (w x h x d)	305 mm x 120 mm x 370 mm
Weight	7.4 kg
Operating mode	Continuous operation
Light emitted (dependent on the mode)	Wavelength
White light	400...700 nm
Near-infrared light (NIR)	700...800 nm
Operating conditions	
Temperature	0 °C...40 °C (32 °F...104 °F)
Relative humidity (non-condensing)	5%...85%
Max. altitude	3,000 m
Storage/transport conditions	
Temperature	-10 °C...60 °C (14 °F...140 °F)
Relative humidity (non-condensing)	5%...95%
Atmospheric pressure	+500 hPa...+1080 hPa

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shock: Applied part type CF
- Type of protection against moisture: Drip-water protection as per IPX0

**According to IEC 60601-1-2:**

Observe the information on electromagnetic compatibility in chapter 10 (p. 37-53).

**Directive compliance** (for TL400)

**According to Medical Device Directive (MDD):**

This medical device belongs to Class IIa  
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

**NOTE:** The code number after the CE mark indicates the responsible notified body.

**Conformidad con la norma** (para TL400)  
**Según CEI 60601-1, IEC 60601-1-6, ANSI/AAMI ES60601-1, CAN/CSA-C22.2 No. 60601-1:**

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#### 8.4 Technical documentation

On request, the manufacturer will provide available parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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### 8.5 Software ownership and licensing

With respect to KARL STORZ products containing software components, Customer has a non-exclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ Products System (Software) internally, only in the form in which delivered to Customer and for the sole purpose of operating in accordance with KARL STORZ's written instructions for such products sold to Customer (and for no other product). The Software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the property of KARL STORZ. Customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the Software or in any way modify the Software. Customer has no right to, and may not, create derivatives of the software, and Customer may not attempt to copy, create or re-create the source code of the Software. Any and all such modifications or enhancements to the Software by Customer, in contravention of this license, will immediately become the sole property of KARL STORZ.

Customer hereby acknowledges and agrees (I) that the purchase of KARL STORZ products does not constitute a sale of the Software, (II) that the Software is the property of KARL STORZ, (III) that Customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (IV) KARL STORZ retains all right, title, and interest in and to the Software, at all times, regardless of the form or media in or on which the original or other copies of the Software may exist. In the event of a failure of Customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the License will, without any further action by KARL STORZ or any other party, immediately terminate.

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**Attribution notice**

You can view a list of the open source software included in the product, and the open source software license terms, on the power LED rubina device itself. Select the tab 'Device info' and then 'Open source licenses' in the Settings menu.

To obtain a copy of the source code being made publicly available by KARL STORZ related to open source software used in this KARL STORZ product, you should send your request, along with the device and serial number, in writing by email to [opensource@karlstorz.com](mailto:opensource@karlstorz.com) or by postal mail to:

KARL STORZ Imaging, Inc.  
Attention: Open Source Software Inquiries  
1 S. Los Carneros  
Goleta, CA 93117  
USA

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Item	Cat. no.
<b>All KARL STORZ light cables are suitable for use</b>	(see catalog; corresponding specialities)
<b>Power cord</b> (Schuko)	400A
<b>Power cord</b> 'Hospital Grade' (USA)	400B
<b>Triple Adapter</b>	TL005
<b>One-pedal footswitch</b> , one-stage	UF101
<b>Patch cable</b> , CAT6a, length 2.0 m, UL-listed	W010275
<b>Sync connecting cable</b> , dia. 3.5 mm, length 100 cm	TL006*

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\* Not all accessories are available in all regions.

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**10. Electromagnetic Compatibility (EMC) Information**

**i NOTE:** In this chapter, the term device is used instead of 'cold light fountain power LED rubina TL400'.

The device is suitable for use in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

**! WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and cables used for EMC Compliance				
Type	Shielded	Length [m]	Ferrite	Used for
PA (Potential equalization)	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

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## Electromagnetic Compatibility (EMC) Information

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**NOTE:** The emission characteristics of this device make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This device offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting the device.

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<b>Table 1 – Compliance level for Immunity tests</b>			
<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.			
<b>Immunity tests</b>	<b>EN/IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the device further from sources of electromagnetic fields or to install magnetic shielding. Before installing the device, the electromagnetic field should be measured to ensure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	6 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	

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**Table 3 – Test levels for Radiated and conducted Immunity Tests**  
Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	6 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> . $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the device is used exceeds the above compliance levels, the device should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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<b>Table 4 – Emission class and group</b> <b>Guidelines and manufacturer’s declaration – electromagnetic emissions</b>		
The device is intended for use in such an environment as specified below. The customer or user of the device should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The device is suitable for use in all establishments including domestic ones and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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<b>Table 5</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the device</b>			
The device is intended for use in an electromagnetic environment in which RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
<b>Rated maximum output power of transmitter [W]</b>	<b>Separation distance d [m] as a function of the transmitter frequency</b>		
	<b>150 kHz to 80 MHz</b> $d = 1.2 \sqrt{P}$	<b>80 MHz to 800 MHz</b> $d = 1.2 \sqrt{P}$	<b>800 MHz to 2.5 GHz</b> $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.			
Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

**HF surgical instrument compatibility**

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga, Ontario  
L5N 3R3 Canada  
Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8525  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Av. Ejercito Nacional No. 453 Piso 2,  
Colonia Granada, Alcaaldia Miguel Hidalgo  
C.P. 11520 Ciudad de México  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6° Piso, B1638 CAA – Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Stamveien1  
1483 Hagan, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
127 39 Skårholmen, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)96824775  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representative Office  
Zalgirio St. 94, LT9300  
Vilnius, Lithuania  
Phone: +370 5 272 0448  
Mobile: +370 685 67 000  
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk  
KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopiae Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com  
KARL STORZ Endoscopiae France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info-hrv@karlstorz.com

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.  
ul. Bójkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoszkóp Magyarország Kft.  
Toberek utca 2. fsz. 17/b  
HU-1112 Budapest, Hungary  
Phone: +36 195 096 31, Fax: +36 195 096 31  
E-Mail: info-hu@karlstorz.com  
KARL STORZ Endoscopia Romania srl  
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4  
041393 Bukarest, Romania  
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801  
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.\*  
Patriarhou Grigoriou E' 34  
54248 Thessaloniki, Greece  
Phone: +30 2310 304868, Fax: +30 2310 304862  
E-Mail: info-gr@karlstorz.com

\*Repair & Service Subsidiary

KARL STORZ Industrial\*\*  
Gedik Is Merkezi B Blok  
Kat 5, D 38-39, Bagdat Cad. No: 162  
Maltepe Istanbul, Turkey  
Phone: +90 216 442 9500, Fax: +90 216 442 9030

\*\*Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK  
Derbenyevskaya nab. 7, building 4  
115114 Moscow, Russia  
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41  
E-Mail: info-ru@karlstorz.com  
TOV LLC KARL STORZ Ukraine  
Avenue Geroyiv Stalingrada Str. 2D, office 717  
Kyiv, 04210/Ukraine  
Phone: +38 095 000-895-0, +38-097-000-895-0,  
+38 073 000-895-0  
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG  
Representation Office  
Sabit Orudschow 1184, apt. 23  
1025 Baku, Azerbaijan  
Phone: +99 450 613 30 60  
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –  
East Mediterranean and Gulf (Offshore) S.A.L.  
Spark Tower 1<sup>st</sup> floor  
Charles Helou St., Horch Tabet – Sin El Fil  
Beirut, Lebanon  
Phone: +961 1 501105, Fax: +961 1 501950  
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.  
P.O. 6061, Roggebaai 8012  
Cape Town, South Africa  
Phone: +27 21 417 2600, Fax: +27 21 421 5103  
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan  
Saryarka, 6, BC "Arman", off. 910  
010000 Astana, Republic of Kazakhstan  
Phone: +7 7172 552-549, 552-788, Fax: -444  
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE  
East Mediterranean & Gulf (branch)  
Building West Side 7A – Unit 7WA – 3008  
Dubai Airport Free Zone, P.O. Box 54983  
Dubai – United Arab Emirates  
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282  
Service Hotline: +971 (0)4 3415882  
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG  
Interchange 21 Tower, Level 33  
399 Sukhumvit Road,  
North Klongtoey, Wattana,  
10110 Bangkok, Thailand  
Phone: +66 2 660 3669  
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch  
Room 701A Building 5 & Room 501 Building 7,  
No. 3000 Longdong Avenue, Pilot Free Trade Zone,  
201203, Shanghai, People's Republic of China  
Phone: +86 21 60339888, Fax: +86 21 60339808  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch  
Room 803-805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
610041, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: sales@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch  
Room 02B & 03 & 04A, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China  
Phone: +86 20 87321281, Fax: +86 20 87321286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Indonesia  
Sinarmas MSIG Tower Level 37  
Jl. Jend. Sudirman No. Kav. 21  
Jakarta Selatan  
DKI Jakarta 12920  
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ Endoscopy Philippines Inc.  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.  
31 Morningside Drive Mt Albert  
Auckland, 1025, New Zealand  
PO Box 56 511, Dominion Rd  
Auckland, 1446, New Zealand  
Phone: +64 9 846 6044, Fax: +64 9 846 6808  
Toll free: +64 508 84 84 84 (New Zealand only)  
E-Mail: sales@karlstorz.co.nz

KARL STORZ Endoscopy Australia Pty. Ltd.  
68 Waterloo Road, Macquarie Park NSW 2113, Australia  
PO Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: karlstorz@karlstorz.com.au

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 7461 708-0  
Telefax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)