AUML

AUML-0427	R 6 x 50 mL + Std 1 x 5 mL
AUML-0497	R 1 x 100 mL + Std 1 x 5 mL
AUML-0507	R 6 x 100 mL + Std 1 x 5 mL
AUML-0707	R 4 x 250 mL + Std 1 x 5 mL
AUML-0420	R 6 x 50 mL
AUML-0500	R 6 x 100 mL
AUML-0250	R 12 x 20 mL
AUML-5405	R 50 mL
AUML-5505	R 100 mL
AUML-5710	R 250 mL
AUML-5220	R 20 mL

FTRO-AUML-v28 (12/2024)_ (PIT-AUML-4-v28)

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BIBLIOGRAFIE

1. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 66.
2. Dufour D. R., The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
3. Johnson, A. M., Amino Acids and Proteins. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 286.
4. Doumas, B. T., et al., Albumin standards and the measurement of serum albumin with bromocresol green, Clin Chim Acta (1971), **31**, 87.
5. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
6. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
7. Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Dufour, D.R., *The Liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
3. *German Society for Clinical Chemistry, Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
4. Wu, A.H.B., *Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 78.
5. Rosalki, S.B., *et al.*, *Clin. Chem.*, (1993), **39**, 648.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 64.
3. Dufour, R., *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
4. Schumann, G., et al., *Clin Chem Lab Med.* (2002), **40**, 718.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Panteghini, M. & Biais, R., *Enzymes*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Dufour, D.R., *The Pancreas: Function and Chemical Pathology*, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 651 and appendix.
3. Winn-Deen, E.S., et al, Clin. Chem., (1988), **34**, 2005
4. Guder W.G, et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, (2002). WHO/DIL/LAB/99.1 Rev.2.
5. Schumann, G., et al. Clin. Chem. Lab Med., (2006), **44**, 1146.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACCC Press, (1997).
8. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACCC Press, (1995).

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 64.
3. Dufour, R., *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
4. Schumann, G., et al., *Clin Chem Lab Med.*, (2002), **40**, 718.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Higgins, T., *et al.*, *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R. Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
2. Dufour, D.R. *The liver: Function and Chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
3. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
4. Wu, A.H.B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 172.
5. Howard, C.R., *Breastmilk & Breastfeeding Jaundice*, Pediatric Clinical Advisor: Instant Diagnosis and Treatment, 2nd Ed., Garfunkel, L.C., Kaczorowski, J.M. & Christy, C., (Mosby Inc. eds.), (2007), 82.
6. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
8. Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid metabolism*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.
3. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
4. Allain, C.C., *et al.*, Clin. Chem., (1974), **20**, 470.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
8. Young D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Wu, A. H. B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 306
2. Sanhai, W.R. et al., Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis. Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.
3. Apple, F. S. & Jaffe, A. S., Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.
4. Boot, E., de Haan, L., Psychopharmacology (Berl), (2000), 150, 347.
5. Schumann, G., et al., Clin. Chem. Lab. Med., (2002), 40, 635.
6. Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), 59, 263.
7. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
8. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995)

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the quality control of the performances of reagents, used for CK-MB quantitative determination.

This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum and human CK-MB.
- Sodium azide < 0.1 % (w/w)
- CK-MB concentrations are lot-specific.

MATERIALS REQUIRED

BUT NOT PROVIDED

- VitalScientific reagent, used for CK-MB quantitative determination.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated on the instructions for use of the reagent used in combination with the recommended calibration factor.

PRECAUTIONS FOR USE AND WARNINGS

- Harmful to aquatic life with long lasting effects.
- Avoid release into the environment.
- Consult Safety Data Sheet (SDS) for a proper handling.
- Each unit of human blood used in the manufacture of this product was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA-approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.

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PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.
- Add exactly 3 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

PRODUCT DETERIORATION

- The product after reconstitution should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

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After reconstitution :

- These products are stable 5 days at 2-8 °C.
- These products should be immediately and tightly capped to prevent contamination and evaporation.

PROCEDURE

To use CK-MB CONTROL, follow the procedure described in the instructions for use of the reagent used.

LIMITATIONS

Using CK-MB CONTROL has been validated with the VitalScientific systems (Analyzers and reagents used for CK-MB quantitative determination).

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Immuno-inhibition/IFCC method - Kinetic.

CK-MB SL/CK-MB reagent contains an antibody inhibiting specifically CK-M subunits (i.e. 100% of CK-MM and 50% of CK-MB). The remaining activity, corresponding to CK-B fraction activity, is quantified according to IFCC reference method for total CK measurement. CK-MB activity is then obtained by multiplying by 2 the remaining activity.

CK-B

Creatine Phosphate + ADP $\xrightarrow{\text{CK-B}}$ Creatine + ATP

ATP + D-Glucose $\xrightarrow{\text{Hexokinase}}$ D-Glucose-6-Phosphate + ADP

G-6-P + NADP⁺ $\xrightarrow{\text{G-6-PDH}}$ D-Gluconate-6-Phosphate + NADPH + H⁺

G-6-P : D-Glucose-6-Phosphate

G-6-PDH : Glucose-6-Phosphate Dehydrogenase

The rate of increase of NADPH concentration measured at 340 nm is directly proportional to the enzymatic activity.

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- Stabilitatea la bord

Refer to § PREPARATION AND STABILITY OF WORKING REAGENT.

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BIBLIOGRAFIE

1. Wu, A. H. B., Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds.), (2006), 306
2. Sanhai, W.R. *et al.*, Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.
3. Apple, F. S. & Jaffe, A. S., Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.
4. Schumann, G., *et al.*, Clin Chem Lab Med., (2002), **40**, 635.
5. Neumeier, D., *et al.*, Clin.Chim.Acta., (1976), **73**, 445.
6. Klein, G., *et al.*, Clin. Chem., (2001), **47**, Suppl. A30.
7. Berth, M. & Delanghe, J., Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
9. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Kaplan, J.M., First, M.R., *Renal Function. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds), (2010), 567 and appendix.
2. Delaney, M.P., Price, C.P., Lamb, E.J., *Kidney Function and Disease. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 631.
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 316.
4. Vasiliades, J., *Clin. Chem.*, (1976), **22**, 1664.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Ceriotti, F. *et al.*, *Reference Intervals for Serum Creatinine Concentrations: Assessment of Available Data for Global Application*, *Clin. Chem.*, (2008), **54**, 559.
7. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press (1995).
8. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
9. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press (1997).

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English - EN

INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of quantitative tests listed in the value sheet.

These *in vitro* diagnostic devices are for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.
- Concentrations for each analyte to test are lot-specific.

MATERIALS REQUIRED BUT NOT PROVIDED

- VitalScientific reagents listed in the value sheet.
- General Laboratory equipment.

TRACEABILITY

The traceability is indicated on the instructions for use of the reagent used in combination with the recommended calibrator.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.
- Carefully open the vial, avoiding the loss of lyophilizate.
- Add exactly 5 mL of distilled or deionized water.
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

Important: To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.
- Do not use after the expiry date stated on the vial label.

After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.
- Stability of the components:
Between 2 and 8 °C : 5 days
Between -25 and -15 °C : 4 weeks (when frozen once)

Exceptions:

- Stability of total bilirubin (when stored protected from light) :
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light) :
Between 2 and 8 °C : 8 hours
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of alkaline phosphatase:
Between 2 and 8 °C : 24 hours
Between -25 and -15 °C : 4 weeks (when frozen once)

Note:

- Store the vials tightly capped and protected from light after reconstitution.
- Mix content thoroughly prior to use.
- Thawing frozen specimen at 37°C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for the use of the reagent used.

LIMITATIONS

Using ELITROL I & II has been validated with the VitalScientific systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.

Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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BIBLIOGRAFIE

1. Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycemia*, (2006).
5. Trinder, P., *Ann. Clin. Biochem.*, (1969), 6, 24.
6. Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Dufour D. R., *The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.

2. Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 916.

3. Doumas, B.T., *et al.*, *Clin. Chem.*, (1981), **27**, 1642.

4. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.

5. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed, AACC Press (1997).

7. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed, AACC Press (1995).

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BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.
3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*, JAMA, (2001), **285**, 2486.
4. Fossati, P. & Prencipe, L., Clin. Chem., (1982), **28**, 2077.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Langlois, M.R., et al. for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
7. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press (1997).
9. Young D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press (1995).

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