MANAGEMENT SYSTEM CERTIFICATE

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10.02.2022

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STATEMENT

We, XEMA Co., Ltd. Having a registered office at 48, 9th Parkovaya st., 104264 Moscow, Russia, assign Sanmedico Srl. Having a registered office at srt. A. Corobceanu 7A, apt. 9, Chişinău MD 2012, Moldova, as authorized representative in correspondence with the conditions of directive 93/42/EEC, 98/79/EEC and 90/385/EEC.

We declare that the company mentioned above is authorized to register, notify, renew or modify the registration of medical devices on the territory of the Republic of Moldova.

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«Cytomegalovirus IgG-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF IgG ANTIBODIES TO CYTOMEGALOVIRUS
IN HUMAN SERUM OR PLASMA**

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
CONTENT

1. INTENDED USE	11
2. SUMMARY AND EXPLANATION	11
3. PRINCIPLE OF THE TEST	11
4. WARNINGS AND PRECAUTIONS	12
5. KIT COMPONENTS	13
6. SPECIMEN COLLECTION AND STORAGE	14
7. TEST PROCEDURE	14
8. QUALITY CONTROL	16
9. CALCULATION OF RESULTS	16
10. EXPECTED VALUES	16
11. LITERATURE	16

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF IGG ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of IgG antibodies to Cytomegalovirus in blood serum or plasma.

This kit is designed for measurement of IgG antibodies to Cytomegalovirus in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 45 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Cytomegalovirus (CMV) belongs to herpesviruses and often causes clinically asymptomatic or mild infection, mostly in young children. It can be transmitted via stool, saliva, and breast milk. CMV can also be transmitted via the placenta and cause severe fetal malformations. Specific IgG-antibodies to CMV are evaluated in women before or during pregnancy to assess and manage the risk of transplacental fetus involvement.

In immunocompromised hosts, CMV reactivation or primary infection may have serious and even life-threatening consequences. Therefore, the absence of specific IgG-antibodies to CMV (seronegativity) in organ transplant recipients requires the seronegativity of the donor.

Specific IgG-antibodies to CMV do not protect from virus reactivation, and usually raise in titer during reactivation caused by decrease of immune system capacity to control the virus replication.

3. PRINCIPLE OF THE TEST

This test is based on indirect enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	cytomegalovirus IgG EIA strips, 8x12 wells	1	pcs		until exp. date
2	Calibrator set, 1.5 ml each. The set contains 5 calibrators: 0; 0.5; 1.5; 3; 6 U/ml	5	pcs	blue(C1 - colourless)	2 months
3	Control serum (1.5 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	red	until exp. date
5	EIA buffer, 50 ml	1	pcs	blue	until exp. date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp. date
7	Washing solution concentrate 26x, 22 ml	1	pcs	colourless	Concentrate - until exp date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp. date
9	Plate sealing tape	2	pcs		N/A
10	Instruction Cytomegalovirus IgG EIA	1	pcs		N/A
11	QC data sheet Cytomegalovirus IgG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0;
- Dry thermostat for 37 °C ±0.1 °C.

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

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Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at – 20 °C or lower. Repeated freezing/thawing should be avoided.

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- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

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7.4. Alternative units:

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1 - 5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 101 fold. (5 µl of sample + 500 µl of diluent) . Do not dilute control sample and calibrators.
3	Pipet 100 µl of calibrators, control sample and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at +37 °C.
5	Prepare washing solution by 26x dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10-20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.
14	Apply point-by-point method for data reduction.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	Sample into the well, µl	Calculation factor
blood serum or plasma	grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing	5 µl of sample + 500 µl of diluent	100	1

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The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of *GLP* (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS

Some laboratories, based on their population studies, set up a second cutoff, which stands between anamnestic ('normal') IgG antibody level and 'high' IgG antibody level characteristic of reactivation or late period of primary infection. Recommended values for this second cut-off for two age groups (8 months – 3 year, > 3 years) are presented in the table below.

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CMV IgG. Based on data obtained by XEMA, the following normal range is recommended (see below).

Sex, age	Units, U/ml	
	Lower limit	Upper limit
Seronegative	<0.1	0.5
Seropositive > 3 years	0.6	2.7
newborn*	<0.1	0.9
under 8 months*	<0.1	1.9
8 months – 3 years	<0.1	3.1

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Russian Diagnostic
Manufacturers Association

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Russian Association
of Medical Laboratory
Diagnosticians

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«Cytomegalovirus IgM-ИФА»

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF IgM ANTIBODIES TO CYTOMEGALOVIRUS
IN HUMAN SERUM OR PLASMA**

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
CONTENT

1. INTENDED USE	10
2. SUMMARY AND EXPLANATION	10
3. PRINCIPLE OF THE TEST	10
4. WARNINGS AND PRECAUTIONS	11
5. KIT COMPONENTS	12
6. SPECIMEN COLLECTION AND STORAGE	13
7. TEST PROCEDURE	13
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	15
11. PERFORMANCE CHARACTERISTICS	16
12. LITERATURE	16

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IGM ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to Cytomegalovirus in serum or plasma.

This kit is designed for measurement of IgM antibodies to Cytomegalovirus in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Cytomegalovirus (CMV) infection is a human infectious disease which may progress both with and without clinical signs and cause local or multiorgan damage. Due to a life-long persistence of CMV in the infected organism, the disease may have recurrences. Most cases of CMV infection are asymptomatic, clinical signs being evident only in immunocompromised patients (HIV infection, chemotherapy and other immunosuppressive treatment). CMV infection in pregnant women may cause intra-uterine infection of the fetus and birth defects. That is why in vitro diagnostics of CMV infection plays an important role - especially, before organ transplantation and during pregnancy.

In case of a primary infection, first anti-CMV antibodies of IgM isotype are found several weeks after infection. They persist in circulation for several weeks and then, from month 2-3, began to gradually decrease to reach minimal concentration by the month 8. CMV-specific IgG antibodies appear ca. 1 week following appearance of IgM and persist in circulation life-long. Recurrence of the infection induces raise of CMV-specific IgG. Sometimes it is accompanied by elevation of CMV-specific IgM level; however, this level never reaches titers seen during primary infection.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated with murine monoclonal antibodies to human IgM. Antibodies from the specimen bind coated murine monoclonal antibodies to IgM on the microwell surface. Unbound material is removed by washing procedure. CMV antigen, labelled with peroxidase enzyme, is then added into the microwells and binds to anti-CMV IgM if present in total IgM fixed. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations). Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Cytomegalovirus IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2 CONTROL -, CONTROL +	polystyrene microwells coated with CMV antigen pp150 (CONTROL-) dilution of preselected human serum, not containing IgM antibodies to Cytomegalovirus with preservative - 0.01% Bronidox L, 0.01% 2-Methyl-4-isothiazolin-3-one-hydrochloride, colourless; (CONTROL+) dilution of preselected human serum with high content of human IgM antibodies to Cytomegalovirus with preservative - 0.01% Bronidox L, 0.01 % 2-Methyl-4-isothiazolin-3-one-hydrochloride; also contains red dye	2	pcs	colour-less ; red	2 months
3 CONJ HRP	aqueous solution of murine monoclonal to IgM coupled with horseradish peroxidase diluted on phosphate buffered solution with casein from bovine milk and detergent (Tween-20), contains 0.1 % phenol as preservative and green dye	1	pcs	green	until exp.date
4 DIL	Saline; contains blue dye	1	pcs	blue	until exp.date
5 SUBS TMB	ready-to-use single-component tetramethylbenzidine (TMB) solution	1	pcs	colour-less	until exp.date
6 BUF WASH 26X	aqueous solution of sodium chloride and detergent (Tween 20), contains proClin300 as a preservative	1	pcs	colour-less	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
7 STOP	5.0% vol/vol solution of sulphuric acid	1	pcs	colour-less	until exp.date
8 N003	Plate sealing tape	2	pcs		N/A
9 K103MI	Instruction Cytomegalovirus IgM EIA	1	pcs		N/A
10 K103MQ	QC data sheet Cytomegalovirus IgM EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat or thermostat shaker for 37 °C ±2 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.

(Continuation see page 14)

K103MI

5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C and continuous shaking at 500-600 rpm.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

Alternative incubation:

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at +37 °C.
5	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at +37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10–20 minutes at +18...+25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

For the assay to be valid, the following requirements should be met:

1. OD450 for CONTROL+ should be ≥ 1 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.
3. OD450 for any CONTROL- replicate should be within 50%-150% of the mean OD450 value for CONTROL-. If any value lies outside this range (although meets requirement #2), it should be discarded and not used for calculation of the mean OD450 value for CONTROL-

9. CALCULATION OF RESULTS

1. Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.
2. Calculate the cut-off value: (mean OD450 for CONTROL-) + 0.3.
3. Calculate Positivity Index (PI) for each sample: $PI = \text{mean OD450}(\text{sample}) / \text{Cut-off}$

10. EXPECTED VALUES AND ASSAY LIMITATIONS

If PI value is greater than 1.1, the result is POSITIVE.

If PI value is less than 0.9, the result is NEGATIVE.

If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative.

Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures. NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity

Specificity of the test was evaluated on 55 serum specimens found negative in DiaSorin (Italy) and NovaTec (Germany). All tested specimens were found negative. Based on these data, specificity of the test is 100%.

11.2. Analytical sensitivity

The kit Cytomegalovirus IgM EIA was evaluated using panel of naturally occurring serum and plasma specimens from Boston Biomedica, Inc (Anti-CMV Mixed Titer Performance Panel PTC202(M)).

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