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Model: EPK-i7010 + EC-3890LZi +EC34-i10M+ ED34-i10T2;

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To whom it may concern,

STATEMENT

We, HOYA Corporation (6-10-1 Nishi-shinjuku Shinjuku-ku, Tokyo 160-0023 Japan), hereby certify that attached document listed in above is the copy of original document.

Seiya Raiju
Vice President
Global Quality Assurance & Regulatory Affairs
PENTAX Lifecare Division
HOYA Corporation

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Product Service

EC Certificate

Production Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex V

(Devices in class I in sterile conditions, sterilised systems or procedure packs)

No. G2S 068357 0029 Rev. 00

Manufacturer

HOYA Corporation
6-10-1 Nishi-shinjuku
Shinjuku-ku
Tokyo
160-0023 JAPAN

EC-Representative:

PENTAX Europe GmbH
Julius-Vosseler-Str. 104,
22527 Hamburg,
GERMANY

Product Category(ies):

**Balloon for Ultrasound Endoscopes,
Sterile Endoscope Channel Seal for the Prevention
of Reflux of Body Fluids from Endoscopes,
Sterile Irrigation Tap for Use
in Biopsy Endoscopic Procedures**

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture in accordance with MDD Annex V. This quality assurance system covers those aspects of manufacture concerned with securing and maintaining sterile conditions of the respective devices / device categories and conforms to the requirements of this Directive. It is subject to periodical surveillance. See also notes overleaf.

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Product Service

EC Certificate

Production Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex V

(Devices in class I in sterile conditions, sterilised systems or procedure packs)

No. G2S 068357 0029 Rev. 00

Facility(ies):

HOYA Corporation
6-10-1 Nishi-shinjuku, Shinjuku-ku, Tokyo,
160-0023 JAPAN

HOYA Corporation PENTAX Lifecare Division
Showanomori Technology Center
1-1-110 Tsutsujigaoka, Akishima-shi, Tokyo,
196-0012 JAPAN

HOYA Corporation PENTAX Lifecare Division
Production Technology Center / Ogawa Factory
395 Oaza-kakuyama, Ogawa-machi, Hiki-gun,
Saitama, 355-0316 JAPAN

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EC Certificate

Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)
(Devices in Class IIa, IIb or III)

No. G1 068357 0028 Rev. 00

Manufacturer:

HOYA Corporation
6-10-1 Nishi-shinjuku
Shinjuku-ku
Tokyo
160-0023 JAPAN

EC-Representative:

PENTAX Europe GmbH
Julius-Vosseler-Str. 104,
22527 Hamburg,
GERMANY

Product Category(ies): Endoscopes, Ultrasound Endoscopes,
their Related Equipment and Accessories (IIa, IIb)

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

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
1. Country: JAPAN

This public document

2. has been signed by **Teruo Takayanagi**

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of **Teruo Takayanagi**, Notary

Certified

5. at Tokyo

6. JAN. 9, 2019

7. by the Ministry of Foreign Affairs

8. 19- **N9 022801**

9. Seal/stamp:

10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

NOTARIAL CERTIFICATE

This is to certify that Lee Jinwoo, an agent of Seiya Raiju, Vice President of HOYA Corporation PENTAX Lifecare Division Global Quality Assurance & Regulatory Affairs has stated in my very presence that said Seiya Raiju acknowledged himself to have signed and sealed to the attached document.

Dated this 9th day of January, 2019.

Teruo Takayanagi
Teruo Takayanagi

NOTARY

9 - 21, 3 - CHOME, SHIBAZAKICHO

TACHIKAWA, TOKYO, JAPAN

TOKYO LEGAL AFFAIRS BUREAU

EC Declaration of Conformity
Dossier No. CE-SC-001

PRODUCT DESCRIPTION

1. Product Name: Radiance 19, 24, 26, 27, 32, 42 and 55" inch Flat-Panel Displays
2. Model Number: SC-WUxx-wxyzz, where w may be A, L, R, S, M or N; x may be 0-9, A or B; y may be 1-9 or A, and zz may be 00-23.
3. Classification: Class I per Annex IX, Rule 12 of the Council Directive 93/42/EEC include amending directive 2007/47/EC
4. Safety Standards: The Radiance series Flat-Panel Display is tested to meet the medical safety requirements of IEC 60601-1:2005 3rd edition, ANSI/AAMI ES60601-1, CAN/CSA C22.2 No. 60601-1, and the Medical Device Directive, 93/42/EEC and 2007/47/EC.
5. Intended use: The Radiance display is a Class I medical device intended for use in a medical environment to display high quality video and graphic images.

NDS QUALITY SYSTEM ASSESSMENT

6. Certification Body: Intertek Certification Number: 0019394-00 and, 0019396-00
7. Conformity Route: NDS Surgical Imaging, LLC, has established a full quality system in accordance with ISO 9001:2008, ISO 13485:2003 and the Medical Device Directive, Annex VII. FDA Registration #2954921

The management of NDS Surgical Imaging, LLC, declares that the devices described above are in compliance with the following directives, standards and normative documents.

- Council Directive 93/42/EEC and 2004/108/EEC include directive 2007/47/EC
- EN 60601-1-2:2001 +A1:2006
- Council Directive 2011/65/EU (also known as "RoHS Recast"). F-21-005, 75B0101

Under 93/42/EEC

- In accordance with EN 60601-1-2:2007
 - CISPR 11:1998 +A1:1999 + A2:2002 Radiated & Conducted Emissions
 - EN 61000-3-2: Harmonics
 - EN 61000-3-3: Flicker
 - EN 61000-4-2: Electrostatic Discharge
 - EN 61000-4-3: Radiated Immunity
 - EN 61000-4-4: Transients
 - EN 61000-4-5: Surges Line to Line
 - EN 61000-4-6: Conducted Immunity
 - EN 61000-4-8: Magnetic Immunity
 - EN 61000-4-11: Voltage Variations & Interruptions

Technical documentation for conformity assessment of the aforementioned device is maintained on file at the following location: NDS Surgical Imaging, 5750 Hellyer Ave., San Jose, CA 95138, U.S.A.

Signed: Chad Chang
Name: Chad Chang

Function: Chief Technology Officer
Date: 2015-03-9 Page 1 of 1

DEC™ Video Duodenoscope ED34-i10T2

Lock in. Elevate your care.

PENTAX
MEDICAL

Raising the standard of health

PENTAX Medical is delivering a unique and innovative advancement in cleaning capabilities for infection prevention and control: the Video Duodenoscope ED34-i10T2, with the single-patient use, sterile, disposable elevator cap (DEC™).

This pioneering solution for diagnostic and therapeutic endoscopic retrograde cholangio-pancreatography (ERCP) procedures combines improved medical device hygiene clinical control with HD image quality (HD+) and procedural performance.

DEC™ with its disposable elevator responds to the need for enhanced patient safety in endoscopy. This need comes as a result of increasing incidences of Carbapenem-resistant Enterobacteriaceae (CRE) and other infections, which may be linked also to improper cleaning and/or disinfection of a duodenoscope's elevator mechanism.

The PENTAX Medical Video Duodenoscope ED34-i10T2 combines a sterile disposable elevator cap (DEC™) for single-patient use and simple disposal that advances cleaning capability of the duodenoscope to help reducing risk of cross contamination and improve patient care with High-Definition image quality for detailed endoscopic visualization during ERCP procedures.

1

Maximize infection prevention, reduce risk

The DEC™ allows simplified reprocessing and increased cleaning capability, and thus helps **reducing the risk for cross-contamination**.

2

Reliable therapeutic performance

PENTAX Medical's first duodenoscope with a disposable elevator that offers **HD image quality** for detailed visualization of the papilla for **high precision** in ERCP procedures.

3


Reprocessing operational efficiency

35% reduction in distal end reprocessing* due to better access for cleaning and disinfection as well as disposability of the elevator.

*In comparison to the standard duodenoscopes of the major manufacturers. Source: PENTAX Medical internal benchmarking.

Maximize infection prevention, reduce risk

Risks factors associated with endoscope infection outbreaks

Bacterial contamination of duodenoscope and accessories in U.S. Report dated from 2012 to 2017.
The data might vary in the future according to new reports.
Source: FDA Maude Report.

“ERCP is associated with the highest risk of infectious complications among all gastrointestinal endoscopic procedures.”

Verfaillie et al, Endoscopy 2015; 47: 493–502

The PENTAX Medical solution

The innovative design of the PENTAX Medical DEC™ Duodenoscope introduces a unique solution in the market and simplifies the overall reprocessing process.

The elevator cap is designed as a disposable component

The elevator area of the duodenoscope is one of the most difficult components to clean. With the DEC™ Duodenoscope design, this challenge is overcome through the elevator integration into the disposable distal end cap.

“Manufacturers may consider enhanced design to ensure the adequate reprocessing of ERCP endoscopes and other complex reusable instrumentation, to further minimize the risk of disease transmission during GI endoscopy.”

Muscarella et al. World J Gastrointest Endosc 2014 October 16; 6(10): 457-474 ISSN 1948-5190

Reliable therapeutic performance

Confidence when delivering therapy

The ED34-i10T2 is a duodenoscope providing High-Definition endoscopic images for routine and therapeutic ERCP procedures. HD+ combined with i-scan allows for better visualization of the lesion and facilitates easier cannulation. Combined with i-scan and i-scan OE, detailed visualization of the mucosal structures is possible to support detection and characterization of lesions and help locate the minor papillae.

Precise control

Material couplings and the disposable design allow the use of a new elevator every time a procedure is performed, with no risk of wear-out-effects.

The improved elevator shape and ergonomic design of the duodenoscope facilitate easy cannulation as well as precise and controlled positioning of devices during ERCP procedures.

Optimized visualization

HD+ imaging allows better visualization and anatomical interpretation of papillary structures, for an easier cannulation of the ducts. i-scan can be used to highlight and assess lesions.


i-scan facilitates the visualization of the orifice of the papilla

i-scan OE increases contrast and visualizes vascular patterns

i-scan can be used to highlight adenomas pit pattern structures

i-scan OE can help to precisely visualize the area to be ablate

Reprocessing operational efficiency

The re-engineered and innovative distal end of the PENTAX Medical DEC™ Duodenoscope allows direct access and exposure to critical surfaces for brushing, cleaning and disinfection.

The disposable elevator cap: better benefit through disposability

Thanks to the distal end cap's single-patient-use design with its integrated elevator there is no need for elevator reprocessing.

Distal end with open access

The DEC™ Duodenoscope distal end, a crucial element for the prevention of patient infection, has been designed to offer easier access for brushing.

35%
reduction in steps
for distal end
reprocessing*

*In comparison to the standard duodenoscopes of the major manufacturers. Source: PENTAX Medical internal benchmarking.

PROFILE single-use cleaning brushes were made to meet the high-quality standards of PENTAX Medical engineering. These brushes have been specifically designed for PENTAX Medical and were tested to reduce residual soil compared to legacy brushes.

PENTAX Medical recommends to reprocess the DEC™ Duodenoscope by using PROFILE single-use cleaning brushes.

Find out more about PENTAX Medical hygiene solutions here:

www.the-hygiene-solution-that-fits.com

Training for Excellence

PENTAX Medical Training for Excellence is a unique platform with the aim of offering professional education and supporting the exchange of best practices. This training platform forms part of PENTAX Medical's drive to improve quality in endoscopy.

www.training-for-excellence.eu

my PENTAX Medical: The exclusive online customer portal

The "my PENTAX Medical" customer portal provides you with 24/7 access to all online services, including manual downloads, repair requests and tracking, as well as training and events registration.

- Convenient access to all PENTAX Medical online services via one platform
- Allows 24/7 access from desktop and mobile devices
- Always updated with the latest information
- Easy access to the contact details of all your PENTAX Medical partners

How to join

www.mypentaxmedical.eu

1

As our customer, you receive your personal set-up credentials together with your personal membership card.

2

You can create a personalized account on 'my PENTAX Medical' using the initial set-up credentials.

3

Please ensure you record your personal password on the membership card or in a safe place.

4

For using your account in the most convenient way, you can organize and set up your departments and equipment as you need.

EC	REP
----	-----

**EMEA Headquarter
Germany**
PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

United Kingdom
PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical.uk@pentaxmedical.com

France
PENTAX France Life Care S.A.S.
116 quai de Bezons
B.P. 204
95106 Argenteuil Cedex
Tel: +33 1 / 3025 7575
Fax: +33 1 / 3025 7445
E-mail: clients.fr@pentaxmedical.com

Netherlands
PENTAX Nederland B.V.
Edisonring 6
6669NB Dodewaard
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: info.nl@pentaxmedical.com

Italy
PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.it@pentaxmedical.com

Spain
SIMMEDICA – Sistemas Integrales
de Medicina, S.A.
Avenida del Sistema Solar 25
28830 San Fernando de Henares · Madrid
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

Russia
Moscow Representative office
of PENTAX Europe GmbH (Germany)
13, 4-y Lesnoy pereulok
125047 Moscow
Tel.: +7 495 114 52 31
Fax: +7 495 114 52 31

Turkey
PENTAX Turkey
Veko Giz Plaza, Meydan Sokak No:3/43
343396 Maslak – Istanbul
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

Japan
 HOYA Corporation
6-10-1 Nishi-shinjuku
Shinjuku-ku
Tokyo 160-00231

NEW

The **MagniView** endoscopes.

A step forward in characterisation.

The PENTAX Medical MagniView endoscopes join the outstanding HD+ image series with an optical zoom function. In combination with i-scan, the MagniView endoscopes provides a new standard for daily practice, supporting enhanced characterisation in a wide range of procedures, with its 136 times magnified image.

Advanced ergonomics and graduated flexibility of the insertion tube, increases both patient and physician comfort.

MagniView: HD+ and Optical Zoom technology combined.

HD+, Optical Zoom technology and i-scan - the perfect combination for improved characterisation.

Powerful magnification

- The optical zoom technology magnifies the image 136 times, allowing for detailed surface structure observation and supporting enhanced characterisation.

Excellent imaging capabilities


- The combination of i-scan, HD+ and magnification in a standard endoscope provides state-of-the-art technology in detection and characterisation in daily clinical routine.

Fully functional endoscope

- Designed to have the same working channel and equivalent outer diameter as a standard diagnostic endoscope, with optical magnification capabilities. Easy to position, easy control of magnification and advanced ergonomic design.

Magnified endoscopic image

Magnified endoscopic image

MagniView VIDEO ENDOSCOPES

	Product	Insertion tube [ø mm]	Distal end [ø mm]	Instrument channel [ø mm]	Working length [mm]	Field of view [mm]	Tip deflection up/down [°]	Tip deflection right/left [°]	Remarks
Gastroscope	EG-2990Zi	9.8	10.6	2.8	1,050	140	210/120	120/120	HD+, MagniView, Water Jet System
Colonoscopes	EC-3890LZi	13.2	13.0	3.8	1,700	140	180/180	160/160	HD+, MagniView, GDF, Water Jet System
	EC-3890FZi	13.2	13.0	3.8	1,500	140	180/180	160/160	HD+, MagniView, GDF, Water Jet System
	EC-3890MZi	13.2	13.0	3.8	1,300	140	180/180	160/160	HD+, MagniView, GDF, Water Jet System

*Distal Rubber Hood: OE-A58 compatible with EG-2990Zi, OE-A59 compatible with EC-3890Zi series

PENTAX doesn't take any liability for the accuracy of the content as laid out in this brochure. In the interest of technical progress specifications are subject to change without notice. Details concerning the compatibility with other medical devices are to be found within the corresponding Instructions for Use.

PENTAX Europe GmbH

Julius-Vosseler-Straße 104 · 22527 Hamburg · Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13 · www.pentaxmedical.com

PENTAX
MEDICAL

Leveraging clinical evidence and expertise

Clinical leaflet

PENTAX
MEDICAL

A quarantine process for the resolution of duodenoscope-associated transmission of multidrug-resistant *Escherichia coli*

Ross et al. *Gastrointestinal Endoscopy*, Volume 82, No. 3: 2015

Background:

Because of their complex design, duodenoscopes have been long recognized to be difficult to fully disinfect and may play a role in transmission of bacteria between patients. Recent reports of duodenoscope associated carbapenem-resistant enterobacteriaceae (ERCP) transmission have confirmed these suspicions. An outbreak of a multidrug-resistant strain of *Escherichia coli* was recently reported at our institution. Herein we report the results of our investigation and the process improvements that we deployed in an effort to contain the outbreak.

Methods:

A full investigation into the environment, endoscopists, infection control practices, high-level disinfection process, and endoscopes was undertaken in conjunction with the local county health authority and the Centers for Disease Control and Prevention. Duodenoscopes were cultured and quarantined for 48 hours until negative cultures were obtained. Ergonomic changes were made to the endoscope reprocessing area, duodenoscopes were returned for routine maintenance, and surveillance cultures were obtained from all patients undergoing ERCP.

Results:

Between November 2012 and August 2013, 32 patients were found to harbor one of two clonal strains of multidrug-resistant *E. coli*, all of whom had undergone ERCP or duodenoscopy. A total of 1149 ERCPs were performed during this time period. Seven patients died within 31 days of the organism being identified in culture, 16 patients died overall by March 2015. The exact contribution of *E. coli* to death is unclear because most patients had underlying late-stage malignancy or other severe medical comorbidities. No breach in high-level disinfection protocol or infection control practices was identified. The clonal strain of *E. coli* was identified in culture on four of eight duodenoscopes, three of which required critical repairs despite lack of obvious malfunction. The defect rate in high-level disinfection of duodenoscopes was 2 % over a one year period. The implemented quality improvements, subsequent to which 1625 ERCPs have been performed, were successful in halting the outbreak.

Conclusions:

The existing manufacturer-recommended high-level disinfection protocols for duodenoscopes are inadequate. Although the ultimate solution may be a design change to the instrument, the timeline for such a change appears long and potentially difficult to exact. In the interim, a reliable method to ensure that bacterial pathogens are not present on the duodenoscope after high-level disinfection is needed (*Gastrointest Endosc* 2015;82:477-83).

This publication reports the results of the investigation following an outbreak of a multidrug-resistant strain of *E. coli* in this institution. Authors took in consideration environment, protocols and equipment in the hospital as possible root causes of the outbreak.

Key facts

- 1149 ERCPs performed (November 2012 and August 2013)
- 32 patients infected with multidrug-resistant *E. coli*
- 7 patients died within 31 days
- 2 % defect rate in high-level disinfection of duodenoscopes

Additional considerations

32 patients accounts for the 2,7 % of the procedures. This can be considered a remarkable number if we correlate it with the early mortality associated with these patients, which in this publication is reported as more than 20 %, most likely due to the additional comorbidities characterizing this patient population.

Two percent defect rate in high-level disinfection of duodenoscopes is a rate that has concerned the authors and led them to question policy and process in the institution. Of note is the fact that out of the eight duodenoscopes that were returned to the manufacturer, four required critical repairs despite the devices' lack of functional defects.

Author remarks

"Because of design complexities, the duodenoscope has been implicated as a vector for bacterial transmission since its introduction into clinical practice."

"[...]solution may be a design change to the instrument."

"Potential solutions include scheduled maintenance after a certain number of uses."

EC	REP
----	-----

**EMEA Headquarter
Germany**
PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com

United Kingdom
PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical.uk@pentaxmedical.com

France
PENTAX France Life Care S.A.S.
116 quai de Bezons
B.P. 204
95106 Argenteuil Cedex
Tel.: +33 1 / 3025 7575
Fax: +33 1 / 3025 7445
E-mail: clients.fr@pentaxmedical.com

Netherlands
PENTAX Nederland B.V.
Edisonring 6
6669NB Dodewaard
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: info.nl@pentaxmedical.com

Italy
PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.it@pentaxmedical.com

Spain
SIMMEDICA – Sistemas Integrales
de Medicina, S.A.
Avenida del Sistema Solar 25
28830 San Fernando de Henares · Madrid
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

Russia
Moscow Representative office
of PENTAX Europe GmbH (Germany)
13, 4-y Lesnoy pereulok
125047 Moscow
Tel.: +7 495 114 52 31
Fax: +7 495 114 52 31

Turkey
PENTAX Turkey
Veko Giz Plaza, Meydan Sokak No:3/43
343396 Maslak – Istanbul
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

Japan
 HOYA Corporation
6-10-1 Nishi-shinjuku
Shinjuku-ku
Tokyo 160-00231

DEC™ Video Duodenoscope ED34-i10T2

Lock in. Elevate your care.

The PENTAX Medical Video Duodenoscope ED34-i10T2 combines a disposable elevator cap (DEC™) that advances cleaning capabilities with High-Definition image quality for detailed endoscopic visualization during ERCP procedures.

Setting the new standard in ERCP

1

Maximize infection prevention, reduce risk

The DEC™ allows simplified reprocessing and increased cleaning capability, and thus helps **reduce the risk of cross-contamination**.

2

Reliable therapeutic performance

PENTAX Medical's first duodenoscope with a disposable elevator that offers **HD+ image quality** for detailed visualization of the papilla for **high precision** in ERCP procedures.

3

Reprocessing operational efficiency

35 % reduction in distal end reprocessing* due to better access for cleaning and disinfection as well as disposability of the elevator.

*In comparison to the standard duodenoscopes of the major manufacturers. Source: PENTAX Medical internal benchmarking.

DEC™ Video Duodenoscope ED34-i10T2

Type	ED34-i10T2
Direction of view [°]	100 (10 retro)
Field of view [°]	100
Depth of field [mm]	4–60
Tip angulation up - down [°]	120–90
Tip deflection right - left [°]	105–90
Distal end width [mm]	13.6
Insertion tube width [mm]	11.6
Minimum instrument channel width [mm]	4.2
Insertion tube working length [mm]	1,250

EC REP EMEA Headquarter Germany
PENTAX Europe GmbH
Julius-Vosseler-Straße 104 . 22527 Hamburg
Tel.: +49 40 / 5 61 92 - 0 . Fax: +49 40 / 5 60 42 13
E-mail: info.emea@pentaxmedical.com . www.pentaxmedical.com

TÜV Süd CE0123 • Medical device class: IIa • This product must be used only by healthcare professionals. Before use and for detailed product specifications, please refer to the instructions for use. In the interest of technical process, specifications may change without notice.

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PENTAX Medical
The i10 endoscope series.

High-Definition and high performance at its best.


The high-end i10 endoscope series is the natural choice to meet your need for extraordinarily High-Definition image (HD+) and high endoscopic performance. The i10 series offers impressive visualisation and excellent therapeutic options. These endoscopes allow you to fully exploit their wide range of features to undertake the most sophisticated procedures. The graduated flexibility allows optimal insertion control of the endoscope, while the advanced ergonomics deliver increased maneuverability, control and comfort.

The i10 endoscope series at a glance.

High performance i10 series HD+ endoscopes

Next generation megapixel CCD for improved HD+ image quality

- Crisp, clear and bright endoscopic image combined with exceptional field of view for a superior visualisation of the mucosa. Gastrosopes, colonoscopes and, the first-to-market outstanding HD+ duodenoscope.

Exclusive Close Focus

- State-of-the-art endoscopy combined with Close Focus allows detailed examination and close visualisation of suspicious mucosal details.

Expanded Therapeutic Option

- Excellent Therapeutic Option in all standard and slim endoscopes – For your daily challenges, the best ratio of external diameter & working channel combined with Water Jet System.

Improved Graduated Flexibility (GDF)

- Improved flexibility of the insertion tube to increase patient comfort and examination performance during every colonoscopy.

Optimized control body

- Lighter weight and improved size to serve you during challenging procedures.

Product	Insertion tube [ø mm]	Instrument channel [ø mm]	Working length [mm]	Total length [mm]	Angle of view [mm]	Focal range [mm]	Tip deflection up/down [°]	Tip deflection right/left [°]
Gastrosopes i10								
EG16-K10	5.4	2.0	1,100	1,416	140	4-100	210/120	120/120
EG27-i10	9.0	2.8	1,050	1,366	140	2-100	210/120	120/120
EG29-i10	9.8	3.2	1,050	1,366	140	2-100	210/120	120/120
Duodenoscope i10								
ED34-i10T	11.6	4.2	1,250	1,566	100	4-60	120/90	105/90
Colonoscopes i10								
EC34-i10L	11.6	3.8	1,700	2,016	140	2-100	180/180	160/160
EC34-i10F	11.6	3.8	1,500	1,816	140	2-100	180/180	160/160
EC34-i10M	11.6	3.8	1,300	1,616	140	2-100	180/180	160/160
EC38-i10L	13.2	3.8	1,700	2,016	140	4-100	180/180	160/160
EC38-i10F	13.2	3.8	1,500	1,816	140	4-100	180/180	160/160
EC38-i10F2	13.2	3.8	1,500	1,816	140	4-100	180/180	160/160
EC38-i10M	13.2	3.8	1,300	1,616	140	4-100	180/180	160/160
EC38-i10M2	13.2	3.8	1,300	1,616	140	4-100	180/180	160/160

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Certification Mark:

Scope of Certificate:

Design and Development, Production,
Distribution and Servicing of
Medical Endoscopes and Related Products

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.:

JAQ235030361

Valid from:

2018-02-01

Valid until:

2021-01-31

Date, 2018-01-26

Stefan Preiß

Page 1 of 2

Product Service

CERTIFICATE**No. Q1N 17 08 68357 024****Applied Standard(s):**

EN ISO 13485:2012 + AC:2012
Medical devices - Quality management systems -
Requirements for regulatory purposes
(ISO 13485:2003 + Cor. 1:2009)
DIN EN ISO 13485:2012
Upgrade required until 2019-03-31

Facility(ies):

HOYA Corporation
6-10-1 Nishi-shinjuku, Shinjuku-ku, Tokyo, 160-0023
JAPAN

HOYA Corporation PENTAX Lifecare Division
Showanomori Technology Center
1-1-110 Tsutsujigaoka, Akishima-shi, Tokyo, 196-0012
JAPAN

HOYA Corporation PENTAX Lifecare Division Medical
Instrument SBU Japan Sales Headquarters
3-35-1 Shimo-renjaku, Mitaka-shi, Tokyo, 181-0013
JAPAN

HOYA Corporation PENTAX Lifecare Division
Production Technology Center / Ogawa Factory
395 Oaza-kakuyama, Ogawa-machi, Hiki-gun, Saitama,
355-0316 JAPAN

HOYA Corporation PENTAX Miyagi Factory
30-2 Okada, Aza-Shimomiyano, Tsukidate, Kurihara-shi,
Miyagi, 987-2203 JAPAN

HOYA Corporation PENTAX Yamagata Factory
4-1 Hinode-cho, Nagai-shi, Yamagata, 993-0012 JAPAN

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To whom it may concern,

STATEMENT

We, HOYA Corporation (6-10-1 Nishi-shinjuku Shinjuku-ku, Tokyo 160-0023 Japan),
hereby certify that attached document listed in above is the original document.

Seiya Raiju
Vice President
Global Quality Assurance & Regulatory Affairs
PENTAX Lifecare Division
HOYA Corporation

MANUFACTURER's CERTIFICATE

September 9, 2019

We, HOYA CORPORATION (6-10-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 160-0023, Japan) being the manufacturer of PENTAX/PENTAX Medical endoscopic equipment, instruments and accessories, herewith certify that our affiliate:

PENTAX EUROPE GmbH
(Julius-Vosseler-Str.104, 22527, Hamburg, Germany)

is the exclusive distributor and supplier of PENTAX/ PENTAX Medical endoscopic equipment, instruments and accessories to the territory of the Russian Federation and Republics of Ukraine, Belarus, Moldova, Georgia, Armenia, Kazakhstan, Tajikistan, Uzbekistan, Turkmenia and Kirgizia. PENTAX Europe GmbH's exclusivity also extends to connected sales, registration, service and training activities.

Hereby we grant our warranty for our products offered by PENTAX Europe GmbH according to our warranty conditions.

HOYA Corporation also certifies that PENTAX Europe GmbH sets pricelist for the sales of all PENTAX/PENTAX Medical products by itself.

The Authorization Letter is valid until September 9, 2022.

Seiya Raiju
Vice President
Global Quality Assurance & Regulatory Affairs
PENTAX Lifecare Division
HOYA Corporation

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
APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by **Shigeru Kaneda**
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of **Shigeru Kaneda**, Notary
Certified
5. at Tokyo
6. SEP. 12, 2019
7. by the Ministry of Foreign Affairs
8. 19- **Nº 060561**
9. Seal/stamp:
10. Signature

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NOTARIAL CERTIFICATE

This is to certify that Hayato Hasegawa, an agent of Seiya Raiju, Vice President of HOYA Corporation PENTAX Lifecare Division Global Quality Assurance & Regulatory Affairs has stated in my very presence that said Seiya Raiju acknowledged himself to have signed and sealed to the attached document.

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CERTIFICATE OF REGISTRATION

This is to certify that the quality management system of:

NDS Surgical Imaging LLC

Main Site: 5750 Hellyer Avenue

San Jose, California 95138, United States

(refer to Appendix for additional sites)

has been assessed by Intertek as conforming to the requirements of:

ISO 13485:2016

The quality management system is applicable to:

The Design and Development, Manufacture, and Servicing of Medical Grade Visualization, Wireless and Integration Products.

(An appendix is also part of this certificate of registration)

Certificate Number:

0019395-06

Initial Certification Date:

17 January 2015

Certificate Issue Date:

13 March 2019

Certificate Expiry Date:

16 January 2021

A handwritten signature in black ink, appearing to read "Calin Moldovean", written over a horizontal line.

Calin Moldovean

President

Intertek Testing Services NA Ltd.,
1829, 32nd avenue, Lachine, QC, H8T 3J1,
Canada

APPENDIX TO CERTIFICATE OF REGISTRATION

This appendix identifies the locations by the management system of

NDS Surgical Imaging LLC

This appendix is linked to the Main Certificate #0019395-06 and cannot be shown nor reproduced without it.

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Scope: The design and development of medical grade visualization, wireless and integration products

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Scope: The servicing of medical grade visualization, wireless and integration products

NDS Surgical Imaging LLC

Novanta Europe GmbH

Parking 57-59, 85748 Garching, Germany

Scope: The order management of medical grade visualization, wireless and integration products.

NDS Surgical Imaging LLC

4531 36th Street, Orlando, Florida, 32811-6527, United States

Scope: The servicing of medical grade visualization, wireless and integration products

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Scope: The manufacture of medical grade visualization, wireless and integration products.

Instructions for Use

PENTAX Video Upper G.I. Scope

PENTAX Video Colonoscope

Operation

EG-2990Zi

EC-3890LZi/FZi/MZi

For cleaning, high-level disinfection, and sterilization of the product after use, refer to the separate Instructions for Use (Reprocessing) with the model name of the endoscope.

Instructions for Use

This Instructions for Use (hereinafter referred to as “IFU”) contains essential information, such as operating procedures and handling precautions, on using this endoscope safely and effectively. Before use, fully understand the contents of, and properly follow, this IFU and the instruction manuals of all equipment that are going to be used in combination. Do not use this endoscope for any purpose other than its intended use.

In addition, review and fully understand the contents of the separate IFU for reprocessing (hereinafter referred to as “IFU (Reprocessing)”). Inappropriate use of the product may result in damage to the equipment or injuries, including, but not limited to, burns, electric shock, perforation, infection, and bleeding.

This IFU does not describe specific endoscopic procedures. The specific procedures should be determined according to the discretion of a medical professional.

If you have any questions or concerns about any information in this IFU, contact your local PENTAX Medical service facility.

The content of the IFU may be changed without prior notice.

Unauthorized reproduction of any part of this IFU is prohibited.

Keep this IFU and all related instruction manuals in a safe, accessible location.

Signal words and symbols

Signal words

The following signal words are used throughout this IFU.

Warning

Indicates a situation that could result in death or serious injury if not avoided.

Caution

Indicates a potentially hazardous situation that could result in minor or moderate injury or damage to equipment if not avoided.

Note

Indicates supplementary or useful information regarding use.

Symbols

The meaning(s) of the symbol(s) on the endoscope, accessories, and/or on their packaging are as follows:

Symbol	Description
	Caution
	Date of Manufacture
	Type BF applied part
	Follow the Instructions for Use
	Manufacturer
	Authorized representative in the European Community
	This product complies with the applicable standards harmonised under the Directive 93/42/EEC and Directive 2011/65EU.

Contents

Instructions for Use	3
Signal words and symbols	3
Important information: Please read before use	6
Product summary	6
Intended use	6
Application	6
Classification	6
Specifications	7
Compatible products	7
Reprocessing before the initial use, reprocessing, and storage after use	8
General warnings and cautions	9
Maintenance management	10
1 Package contents	12
1-1. Package contents	12
2 Nomenclature and functions	14
2-1. Control body, insertion portion	14
2-2. Connector	16
3 Preparation and inspection	17
3-1. Preparation of the equipment	18
3-2. Inspection of the endoscope	20
3-3. Inspection of accessories and attachment to the endoscope	28
3-4. Inspection and connection of ancillary equipment to the endoscope	42
3-5. Inspection of the endoscopic system	45
4 Directions for use	56
4-1. Preparation immediately before insertion of the endoscope	58
4-2. Insertion and observation	59
4-3. Using an endoscopic device	65
4-4. Using a non-flammable gas	68
4-5. Laser cauterization	70
4-6. Electrosurgery	71
4-7. Withdrawal of the endoscope	72
4-8. Care after use	

5 Troubleshooting	74
5-1. Withdrawal of an endoscope with an abnormality	74
5-2. Returning the endoscope for repair	75
Disposal	76
Electromagnetic compatibility (EMC)	77
Endoscope specifications	80
System chart	81

Important information: Please read before use

Product summary

This endoscope visualizes subjects under illumination transmitted from a dedicated video processor with a solid-state image sensor located at the distal end of the endoscope and provides images for observation of the target anatomy through the images reproduced on the video monitor via the video processor.

It can be used with endoscopic devices which are introduced from the instrument channel inlet of the control body.

The endoscope also allows for angulation operation of the bending sections via operation of the angulation control knob; air/water feeding from the distal end of the endoscope via operation of the air/water feeding valve; and suction through the channel at the end of the endoscope via operation of the suction control valve.

Intended use

■ EG-2990Zi:

This PENTAX Video Upper G.I. scope is intended to provide optical visualization of (via a video monitor), and therapeutic access to, the upper gastrointestinal tract. This anatomy includes, but is not restricted to, the organs; tissues; and subsystems: esophagus, stomach, and duodenum.

This endoscope is introduced via the mouth or the nose, as decided by the physician, when indications consistent with the need for the procedure are observed in adult and pediatric patient populations.

■ EC-3890Zi series:

These PENTAX Video Colonoscopes are intended to provide optical visualization of (via a video monitor), and therapeutic access to, the lower gastrointestinal tract. This anatomy includes, but is not restricted to, the organs, tissues, and subsystems: large bowel to the cecum.

These endoscopes are introduced via the rectum when indications consistent with the need for the procedure are observed in adult and pediatric patient populations.

Application

Medical purposes	Provide images for observation, diagnosis, visualization, and treatment.
Patient population	Patients who are considered suitable for the application of this endoscope by the physicians (pediatric to adult patients).
Intended anatomical area	EG-2990Zi; Upper gastrointestinal tract (esophagus, stomach) duodenum EC-3890Zi series; Large intestine and terminal ileum
User qualifications	Physicians (Experts who have been approved by the endoscopic medical safety administrator at each medical facility. If the eligibility requirements are defined by an official body, such as a government entity and/or an academic society, follow such requirements). Specific training to use this endoscope is not required.
Location of use	A medical facility

Classification

Degree of protection against electric shock for the applied parts	TYPE BF applied part (when connected to a compatible PENTAX Medical video processor)
Degree of protection against water	IPX7 (with the soaking cap attached)
Mode of operation	Continuous operation

Specifications

■ Environment

Operating environment	Ambient temperature	10 to 40 °C
	Relative humidity	30 to 85 % RH
	Air pressure	700 to 1,060 hPa

Storage/transportation environment	Ambient temperature	–20 to 60 °C
	Relative humidity	0 to 85 % RH
	Air pressure	700 to 1,060 hPa

■ Software version

For the software version for each model, refer to the back cover.

■ Endoscope specifications

For details, refer to “Endoscope specifications” (p. 80).

Compatible products

This section describes the equipment that can be used in combination with this endoscope. For more details, refer to “System chart” (p. 81).

For the equipment used in combination during cleaning/high-level disinfection/sterilization, refer to the separate IFU (Reprocessing) of this endoscope.

The combinations of equipment and accessories that can be used with this product are listed below. Prior to use, the product must be prepared and inspected according to its IFU.

Warning

PENTAX Medical does NOT warrant compatibility with unlisted products.

If products are NOT listed, contact the manufacturer of the equipment or accessory to confirm the compatibility and instructions for use with PENTAX Medical products.

Note

When this endoscope is used in combination with other equipment, depending on how it is connected, it may result in malfunction and/or unforeseen events to patients and/or medical professionals. Pre-use operation checks and risk management associated with such changes are recommended, particularly when the equipment used in combination is changed, added, or upgraded.

■ Video processor

Video processor models that can be connected with these endoscopes are shown below. For instructions on video processor operation, refer to the IFU of the respective video processor.

Model Name	Brand Name
EPK-i7000 series	PENTAX Medical
EPK-i5000 series	
EPK-i	

Reprocessing before the initial use, reprocessing, and storage after use

■ Reprocessing before the initial use

The endoscope identified in this IFU is a reusable semi-critical device. Since it is packaged non-sterile, it must be cleaned and high level disinfected, or cleaned and sterilized, or cleaned and disinfected and additionally can be sterilized (if applicable) according to the separate IFU (Reprocessing) of this product before initial use. Insufficient reprocessing may increase the risk of cross contamination.

Note

The wording "high-level disinfection" in this IFU defines the disinfection of the endoscope and the accessories with a completely virucidal disinfectant.

■ Reprocessing

After use, the endoscope must be appropriately cleaned, high-level disinfected and/or sterilized, and stored. Insufficient and/or incomplete cleaning, high-level disinfection, and/or sterilization of this endoscope may result in its non-optimal function of and/or damage to the endoscope and may pose a risk of infection to the patient and/or users.

Warning

When using an endoscope and its accessories on patients with Creutzfeldt-Jakob disease (CJD) or variant Creutzfeldt-Jakob disease (vCJD), use only dedicated instruments and equipment. The instruments and equipment used on these patients must be discarded so that they cannot be used again on another patient. The pathogenic agents that cause this disease, which are called "prions", cannot be destroyed or inactivated using the cleaning, disinfection, and sterilization methods presented in this IFU. Please consult the guidelines that apply to your country or region for more detailed information regarding the handling of prion-contaminated instruments.

■ Storage after use

Warning

Observe the following guidelines. Failure to do so may result in contamination of the endoscope with bacteria or pose a risk of infection to patients and/or users.

- Ensure that all removable accessories, such as air/water feeding valve, suction control valve, inlet seal, cleaning adapter, and soaking caps are removed from the endoscope when storing.
- Do NOT store the endoscope in areas of high humidity or high temperature.
- Do NOT store the endoscope, its components, and accessories in the carrying case.
- Ensure that the endoscope, its components, and accessories are completely moisture-free before storage.
- Before the next use, the endoscope, its components, and accessories that have been stored inappropriately or for a prolonged period of time must be subjected to appropriate cleaning, high-level disinfection, and/or sterilization processes according to the separate IFU (Reprocessing).

Caution

Observe the following precautions when storing the endoscopes, its accessories, or device. Failure to do so may result in damage to property.

- Endoscope insertion portion, umbilical cord, and endoscopic devices should be kept as straight as possible during storage.
- Keep away from chemicals, direct sunlight, or ultraviolet rays.
- Maintain adequate distances between the endoscope, its accessories, and devices, so that they do NOT hit against each other.

Note

It is recommended to store the endoscope hanging down straight in a well-ventilated room or cabinet dedicated for endoscope storage.

For storage after use, also refer to the separate IFU (Reprocessing) of this endoscope.

General warnings and cautions

Warning

- The medical facility should determine restrictions or non-use of the endoscope in patients suspected of having lowered immunity.
- Users as well as the assisting personnel should always wear protective equipment (e.g., gloves, goggles, masks, medical gowns, etc.) to minimize the risk of cross contamination, as patient's body fluids may be dispersed from instrument components such as the instrument channel inlet and the suction control valve.
- Do NOT use this endoscope for any purpose other than its intended use. Doing so may result in patient injury.
- Do NOT use this endoscope with equipment other than those that have been specified for combined use. Endoscope operation in the freeze or magnification mode may result in damage to the endoscope and patient injury.
- Do NOT drop this endoscope or apply a strong shock to it. Doing so may result in damage to the endoscope. In particular, do NOT apply a strong shock to the lens surface at the distal end. Visual abnormalities may occur, which may result in unforeseen events.
- Ensure to attach/connect an appropriate device to the connectors of the PVE connector such as suction nipple, air/water port, venting connector, or feedback terminal according to the IFU.
- Always check the endoscopic image during endoscope angulation, air/water feeding, and suctioning, use of endoscopic devices, and endoscope insertion and withdrawal. Ensure that these operations are performed in the normal (non-frozen, non-magnified) mode. Endoscope operation in the freeze or magnification mode may result in damage to the endoscope and/or patient injury.
- Do NOT forcefully insert and withdraw the endoscope. Doing so may result in patient injuries, including bleeding and perforation.
- Do NOT perform reverse observations inside a narrow lumen. Doing so may cause patient injury or make it impossible to withdraw the endoscope.

Warning

After using operational/cleaning accessories (e.g., forceps, needles, snares, brushes etc.) with the endoscope, carefully check that all accessories are intact and that no parts have fallen off and become lodged within the endoscope's instrument/suction channel. Furthermore, ensure that any therapeutic devices (e.g., clips, stents, etc.) passed through the channel are accounted for after use. If the instrument/suction channel becomes blocked or clogged due to the accumulation of debris, an accessory that cannot be removed, or other cause, do NOT attempt to correct the blockage or continue to use the endoscope. In such a case, contact your local PENTAX Medical service facility to have the endoscope repaired. The use of an endoscope with a blocked internal channel may result in ineffective reprocessing and/or the introduction of debris and/or device components into a patient during a subsequent procedure, posing a risk of cross contamination.

Caution

- Do NOT excessively twist, rotate, or bend any of the insertion portions, strain relief boots, or umbilical cord. Doing so may damage the endoscope.
- Do NOT hit the remote buttons with hard objects or pull or twist them. Doing so may cause internal damage to the endoscope that may lead to water leaks.
- Do NOT attach or remove the PVE connector of the endoscope while the power of the video processor is turned on. Doing so may damage the endoscope.
- Electromagnetic interference may occur with equipment labeled with the following symbol or near mobile RF communication equipment such as mobile phones. If electromagnetic interference occurs, reorient or relocate the endoscope or shield the location of use.

Maintenance management

The service life of this endoscope is 6 years after date of shipment with the following conditions.

- Perform inspection before use, care after use, storage, and replacement of consumables according to this IFU.
- Have a specialist specified by PENTAX Medical perform repairs and at least annual periodic inspections.

1 Package contents

1-1 . Package contents

Check the package contents according to the separate Standard Accessories List provided with this product. For detail picture of the contents/accessories, refer to Figure 1.1 and 1.2 in the IFU.

If there are any damaged or missing components, do not use the endoscope; immediately contact your local PENTAX Medical service facility. (Optional depending on the model.).


Figure 1.1


Figure 1.2









■ Accessories		
		
Cleaning Adapter (OF-B153)	Irrigation Tube (OF-B113)	Ventilation Cap (OF-C5)
		
O Ring Set (OF-B192*) * For Air/Water Feeding Valve (OF-B188)	O Ring Set (OF-B127**) ** For Suction Control Valve (OF-B120)	Check Valve Sets (OE-C15)
		
Distal Hood (For EG-2990Zi, OE-A58) (For EC-3890Zi series, OE-A59)	Bite Block (OF-Z5) (For EG-2990Zi)	Air/Water Feeding Valve (OF-B188) (Installed)
		
Suction Control Valve (OF-B120) (Installed)	Water Jet Check Valve Adapter (OE-C12) (Installed)	Water Jet Connector Cap (OF-B118) (Installed)
■ Endoscopic device*** *** This is an optional device depending on the sales region.  Biopsy Forceps (EC-3890LZi/3890FZi: KW2422R), (EC-3890MZi: KW2418R), (EG-2990Zi: KW2415R)		
■ Others		
		
IFU (Operation; this document)	IFU (Reprocessing)	Standard Accessories List

Figure 1.2

2

Nomenclature and functions

2

2-1 . Control body, insertion portion

Nomenclature and functions


Figure 2.1

1. Remote Buttons 1-3
Functions assigned to each button can be remotely controlled by pressing each of the remote buttons. Functions of the remote buttons 1-3 are assigned from the video processor.
Refer to the IFU of the respective video processor for assignment of functions to each remote button.
2. Magnification control lever
The optical zooming function is available through the magnification control lever operation.
3. Up/Down Angulation Control Knob
By turning in the "▲ U" direction, the bending section moves upwards.
By turning in the "▲ D" direction, the bending section moves downwards.
4. Up/Down Angulation Lock Lever
By turning counterclockwise, upward/downward bending of the bending section is locked.
By turning in the "F ►" direction, the bending lock is released.
5. Right/Left Angulation Control Knob
By turning in the "▲ R" direction, the bending section moves to the right.
By turning in the "▲ L" direction, the bending section moves to the left.
6. Right/Left Angulation Lock Knob
By turning counterclockwise, right/left bending of the bending section is locked.
By turning in the "F ►" direction, the bending lock is released.
7. Suction Control Valve (OF-B120)
Attach to the suction cylinder. Depress it to suction fluids or air through the instrument channel of the endoscope.
8. Air/Water Feeding Valve (OF-B188)
Attach to the air/water feeding cylinder. Covering the hole on the top of the valve feeds air from the air nozzle at the distal end of the endoscope. Depressing the valve feeds water from the water nozzle.
9. Strain Relief Boot
The strain relief boot protects the connecting parts.
10. Inlet Seal (OF-B190)
The inlet seal is attached to the instrument channel inlet to avoid fluid/air leakage.
11. Model Name Label
The model name label shows the model name, minimum instrument channel width, and other related information. (Figure 2.2)


Figure 2.2

2-2. Connector


Figure 2.3

12. Light Guide Plug
The light guide plug transmits the light received from the light source to the distal end of the endoscope..
13. Suction Nipple
Connect the suction tube on the suction source to the suction nipple.
14. Air/Water Port
Connect the air/water feeding hose on the water bottle assembly to the air/water port.
15. Water Jet Port
Attach a water jet check valve adapter (OE-C12).
16. Water Jet Check Valve Adapter (OE-C12)
Use it by attaching to the water jet port.
Connect the irrigation tube (OF-B113) to send sterile water from a syringe or irrigation pump to the water jet nozzle at the distal end of the endoscope.
When an irrigation tube is not connected, close it with the water jet connector cap (OF-B118).
17. Venting Connector
Attach the ventilation cap (OF-C5) or the endoscope connector of the leakage tester here.
18. Feedback Terminal
When using an electrosurgical device, connect the electrosurgical unit's endoscope feedback cord (S-cord) or the condenser earth cable (OL-Z4/OL-Z3: check the package contents according to the separate standard accessories list provided with this product.).

3

Preparation and inspection

Before use, the endoscope, accessories, video processor, and other components must be prepared and carefully inspected according to the IFU. Any equipment used in combination with the endoscope must also be prepared and inspected according to the respective instruction manuals.

Always perform pre-use inspection before each use.

If any abnormality is suspected during inspection, do not use the endoscope; send it for repair according to "5-2. Returning the endoscope for repair" (p. 75).

Warning

Always perform pre-use inspection before each use. NEVER use an endoscope with a suspected abnormality. Doing so may result in malfunction, endoscope damage, and/or injury to the patient and/or user.


Caution

Ensure that another endoscope is also prepared to avoid interruption of the procedure due to endoscope failure or unforeseen events.

3-1 . Preparation of the equipment

Prepare the endoscope, accessories, ancillary equipment, and protective equipment. Refer to the "Compatible products" to prepare the ancillary equipment as necessary and to the IFU provided with the video processor for its inspection.

■ Endoscope


■ Accessories

Air/Water Feeding Valve (OF-B188)

Suction Control Valve (OF-B120)

Inlet Seal (OF-B190)


Silicone Oil Lubricant (OF-Z11)

Water Jet Check Valve Adapter (OE-C12)

Irrigation Tube (OF-B113)


Water Jet Connector Cap
(OF-B118)

Bite Block (OF-Z5)
(For EG-2990Zi)

Distal Hood
(For EG-2990Zi, OE-A58)
(For EC-3890Zi series, OE-A59)

3-2. Inspection of the endoscope

Prepare an endoscope that has been reprocessed according to the procedure specified in the separate IFU (Reprocessing) of this endoscope.

Warning

- NEVER disassemble or modify the endoscope. Doing so may impair its original functionality and possibly result in serious injury to the patient and/or user.
- NEVER use an endoscope with any abnormality. Doing so may result in endoscope damage, detachment of parts into the patient's body cavity, malfunction during use, and/or injury to the patient and/or user.
- Use only sterile water for inspection. Failure to do so may result in contamination of the endoscope with waterborne bacteria and other microorganisms. Do NOT use water that has been left uncovered for a prolonged period of time.

Caution

- Do NOT excessively twist, bend, or rotate any of the strain relief boots on the instrument (See Figure 3.2 (1) and 3.2 (2) to identify the strain relief boots). Doing so may result in instrument damage. Pay special attention to the careful handling of the strain relief boot of the insertion portion (See Figure 3.2 (1)) of the endoscope, because it has a small diameter and is more likely to suffer damage due to mishandling.
- When carrying the endoscope, do NOT grasp or carry it only by its umbilical cord or insertion portion. Moreover, do NOT squeeze or forcefully bend the bending section. (Figure 3.3) Doing so may result in equipment damage.


Figure 3.2


Figure 3.3

**Note**

In case the endoscope is hot/cold immediately after cleaning, high-level disinfection, and/or sterilization, wait until it returns to room temperature before using it. Observation may be hindered by fogging of lens and other effects because of the difference in the temperature of the endoscope and the room.

■ Carrying the endoscope by hand

When carrying the endoscope by hand, loosely loop the umbilical cord and insertion portion, hold the control body and insertion portion (near the bending section) in one hand, and hold the PVE connector in the other hand as shown in Figure 3.4.


Figure 3.4


Inspection of the entire endoscope

Warning

NEVER use the endoscope with any abnormality in function. Doing so may result in endoscope damage, detachment of parts into the patient's body cavity, endoscope malfunction during use, and/or injury to the patient and/or user.


1. Check the entire surface of the endoscope for any visible adhered material.
2. Check the entire surface of the insertion portion for abnormalities such as wrinkles, scars, sharp edges, clouding of the surface, dents, catching, protrusions, attachment of foreign materials, detachment of parts, etc.
3. Check the surface of the adhesive bands on both ends of the bending section for abnormalities such as scratches, clouding, and peeling. With clean gauze, lightly wipe the surface of the adhesive bands to ensure that there is no catching and/or attachment of the adhesive on the gauze.

(1) Adhesive Bands

Figure 3.5

4. Check the case of the distal end of the endoscope (especially around the periphery of the instrument channel) for any abnormalities such as deformation or chipping.

- (1) Objective Lens
- (2) Light Guide
- (3) Water Jet Nozzle
- (4) Instrument Channel
- (5) Air Nozzle
- (6) Water Nozzle
- (7) Case

Figure 3.6

5. Check the objective lens at the distal end of the endoscope and the light guides for any abnormalities such as attachment of foreign material, scratches, or chipping, and ensure that there is no gap on the periphery of the lens.
6. Ensure that there are no scratches, clouding, or peeling on the surface of the adhesive glue around the objective lens at the distal end of the endoscope and it looks glossy.
7. Gently clean the objective lens and light guides with clean gauze or a cotton-tip applicator moistened with 70%–90% medical grade ethyl or isopropyl alcohol. Check that there is no attachment of the adhesive on the gauze.

Note

Clear images cannot be obtained when there are any attachment of foreign material or residuals on the objective lens or light guides. When an endoscope with any attachment of foreign material or residuals on the objective lens or light guides is used, vapor may be generated from the water contents of the attachment of foreign material or residuals being heated by the light.

8. Check the air nozzle and water nozzle at the distal end of the endoscope for any abnormalities such as clogging, dents, deformations, chipping, etc.
9. Form an arch with the insertion tube as shown in Figure 3.7 using both hands. By sliding the insertion tube in the direction of the arrow in Figure 3.7, check that the entire insertion tube can be bent smoothly and easily to form an arch.


Figure 3.7

10. Check the entire surface of the umbilical cord for abnormalities such as wrinkles, scars, sharp edges, clouding of the surface, catching, protrusions, attachment of foreign materials, detachment of parts, etc.

11. Check the control body, PVE connector, and electrical contacts for abnormalities such as scratches, deformities, loose parts, etc. Pay special attention when checking the parts shown in the Figure 3.8. Take clean, lint-free gauze to gently hold these parts, and move them in various directions to ensure that there are no abnormalities such as looseness.


Figure 3.8

Inspection of the angulation mechanism

Ensure that there is nothing near the bending section that would hinder its operation, and inspect the angulation mechanism while the insertion portion is kept straight.

■ Inspection of bending function

Warning

Do NOT use the endoscope with angulation abnormalities such as non-smooth operation, excessive play in the angulation control lever, or excessive angulation attenuation, as the inside of the endoscope may be damaged. Using the endoscope in this state may result in worsening of endoscope damage, malfunction during use, and patient injury.

1. Turn the up/down angulation lock lever and right/left angulation lock knob in the "F ►" direction until they stop to release the lock of the angulation control knobs.


Figure 3.9

- Turn the up/down and right/left angulation control knobs slowly in each direction until they stop, and return them to their original position. Check that the angulation control knobs operate smoothly with no roughness or catching.


Figure 3.10

- Check that the bending section angulates in the direction in which the angulation control knobs are turned and that the maximum angulation can be achieved.


Figure 3.11

- Turn the angulation control knobs back to the neutral position. Check that the bending section returns to an approximately straight condition.


Figure 3.12

■ Inspection of the up/down bending lock mechanism

1. Turn the up/down angulation lock lever counterclockwise until it stops.


Figure 3.13

2. Turn the up/down angulation control knob slowly in the "▲ U" or "▲ D" direction until it stops.


Figure 3.14

3. Check that the curved form of the bending section is fixed when releasing the angulation control knob.
4. Turn the up/down angulation lock lever in the "F ►" direction until it stops to release the lock. Check that the bending section returns to an approximately straight condition.


Figure 3.15

■ Inspection of the right/left bending lock mechanism

1. Turn the right/left angulation lock knob counterclockwise until it stops.


Figure 3.16

2. Turn the right/left angulation control knob slowly in the "▲ R" or "▲ L" direction until it stops.


Figure 3.17

3. Check that the curved form of the bending section is fixed when releasing the angulation control knob.
4. Turn the right/left angulation lock knob in the "F ►" direction until it stops to release the lock. Check that the bending section returns to an approximately straight condition.


Figure 3.18

3-3. Inspection of accessories and attachment to the endoscope

When using reusable accessories, ensure that they have been cleaned, high-level disinfected, and/or sterilized according to the separate IFU (Reprocessing) of this endoscope.

Warning

NEVER disassemble or modify the accessories and endoscopic devices. Doing so may impair their original functionality and possibly result in serious injury to the patient and/or user.

3

Preparation and inspection

Inspection of the air/water feeding valve (OF-B188)


Warning

- If any abnormality is suspected with the check-valve, replace the air/water feeding valve with a new one. Continue use of the air/water feeding valve with abnormalities could cause unintended continuous air feeding and pose a risk of pain or perforation to the patient. It could also pose a risk of infection to the user as result of reflux or dispersal of patient's body fluids from the air/water feeding valve.
- The O-ring of the air/water feeding valve is a consumable. If any abnormality is suspected with the O-ring, stop use immediately and replace it with a new one. Use the compatible O-ring set for replacement. Using an O-ring with abnormalities or non-compatible O-ring could lower the function of air/water feeding, cause unintended continuous air feeding, and pose a risk of pain or perforation to the patient. It could also pose a risk of infection to the user as result of reflux or dispersal of patient's body fluids from the air/water feeding valve.
- The replacement O-ring is NOT sterilized or disinfected before shipment. Perform cleaning and high-level disinfection, and/or sterilization of the air/water feeding valve after O-ring replacement.

Note

Use the O ring set (OF-B192) for air/water feeding valve (OF-B188) for replacement. For details on the O-ring replacement method, refer to the IFU provided with the O ring set (OF-B192).

- (1) O-ring
- (2) Check Valve
- (3) Hole

Air/Water Feeding Valve (OF-B188)

Figure 3.19

1. Check the air/water feeding valve (OF-B188) for any abnormalities such as attachment of foreign materials, deformation, cracks, or hole blockage.

2. Check that the O-ring is properly attached and that there is no chipping, breaks, or peeling in the O-ring or check-valve.

Inspection of the suction control valve (OF-B120)


Warning

- If any abnormalities are suspected in the rubber seal, replace the suction control valve with a new one. Using a suction control valve with any abnormality could result in continuous weak aspiration, which may hinder the procedure. It could also result in potential reflux or dispersal of patient's body fluids, posing a risk of infection.
- The O-ring of the suction control valve is a consumable. If any abnormality is suspected with the O-ring, stop use immediately and replace it with a new one. Use the compatible O-ring set for replacement. Using an O-ring with abnormalities or non-compatible O-ring could result in unintended continuous suction and may hinder the examination. It could also pose a risk of infection to the user as a result of reflux or dispersal of patient's body fluids from the suction control valve.
- The replacement O-ring is NOT sterilized or disinfected before shipment. Perform cleaning and high-level disinfection, or sterilization of the suction control valve after O-ring replacement.

Note

Use the O ring set (OF-B127) for the suction control valve (OF-B120) for replacement.

- (1) Rubber Seal
- (2) O-ring
- (3) Hole

Suction Control Valve (OF-B120)

Figure 3.20

1. Check the suction control valve (OF-B120) for any abnormalities such as attachment of foreign materials, deformation, cracks, or hole blockage.
2. Check that the O-ring is properly attached and that there is no chipping, breaks, or peeling in the O-ring or sealing rubber part.

Inspection of the inlet seal (OF-B190)

Warning

NEVER use an inlet seal (OF-B190) that has any abnormality. Replace it with a new one. Inlet seals are consumables. Using a damaged and/or worn inlet seal may result in lowered suction function and potential reflux or dispersal of patient's body fluids, posing a risk of infection.

1. Check the slit in the cap of the inlet seal (OF-B190) and the hole of the body of the inlet seal for any abnormalities such as cracks, wear, chipping, and attachment or presence of foreign materials. Check that the light does not shine through the slit of the cap.


Figure 3.21

2. Put the cap to the body of the inlet seal and check that the cap is correctly attached.


Figure 3.22

Inspection of the water jet check valve adapter (OE-C12)


Warning

- NEVER use a water jet check valve adapter (OE-C12) that has any abnormality. Replace it with a new one. Check valve adapters are consumables. Using a damaged check valve adapter may result in potential reflux or dispersal of patient's body fluids, posing a risk of infection.
- Replacement check valve sets (OE-C15) are NOT sterilized or disinfected before shipment. Perform cleaning and high-level disinfection or sterilization of the water jet check valve adapter after check valve replacement.

Note

Use the check valve set (OE-C15) for replacement.

- (1) O-ring
- (2) Hole
- (3) Check Valve (OE-C14)

Figure 3.23

1. Check the water jet check valve adapter for any abnormalities such as attachment of foreign materials, deformation or cracks, or hole blockage.
2. Ensure that the check valve is attached correctly to the water jet check valve adapter without any gaps or pinching.


Figure 3.24

3. Check the O-ring and check valve for any abnormalities such as cracks, breaks, and peeling.

Inspection of the irrigation tube (OF-B113)

Warning

NEVER use the irrigation tube (OF-B113) when abnormality is suspected in inspection. Replace it with a new one. Using the OF-B113 with abnormality in the process of cleaning, high level disinfection, or sterilization may cause leaking of detergent from the connection part and detachment of the OF-B113. The cleaning, high-level disinfection or sterilization may NOT be effective due to the insufficient processing.


Figure 3.25

1. Check the entire surface of the irrigation tube (OF-B113) for abnormalities such as bending/breakage/looseness of the connector, cut/chip of the O-ring, buckling/deterioration/hardening of the tube, and/or broken luer.
2. Attach the syringe filled with the sterile water to the luer connector of the irrigation tube (OF-B113) and flush sterile water.
3. Check that sterile water flows in a steady stream from the connector of the irrigation tube (OF-B113).

Inspection of the distal hood (OE-A58/OE-A59)

Warning

- The distal hoods (OE-A58 / OE-A59) are provided non-sterile and non-disinfected. Prior to the initial use, it must be subjected to appropriate cleaning and high-level disinfection, or sterilization according to the respective IFU for the distal hood (OE-A58 / OE-A59). Insufficient processing prior to use may increase the risk of cross contamination.
- NEVER use a distal hood that has any abnormality. Replace it with a new one. Using a distal hood with an abnormality may result in endoscope damage and patient injury.

Note

Distal hoods are consumables. Even while it is still within its service life, avoid reuse whenever possible, instead replace it with a new one on every use.

Check the distal hood for any abnormalities such as attachment of foreign materials, cracks, deformity, chipping, and discoloration.

Figure 3.26

Inspection of the bite block (OF-Z5) (EG-2990Zi only)

Caution

NEVER use a bite block with any abnormality. Replace it with a new one.

Using a bite block with an abnormality may result in endoscope damage and injury to the oral cavity of patients.

Check the bite block for any abnormalities such as attachment of foreign materials, cracks, deformity, chipping, and discoloration.


Figure 3.27

Inspection of the endoscopic devices

For details on the inspection of each endoscopic device, refer to the instruction manual provided with the specific endoscopic device. For reusable endoscopic devices, prepare ones that have been cleaned and sterilized by following the instruction manual for the respective endoscopic device.

Warning

- NEVER use an endoscopic device with signs of damage and/or operational abnormality. Doing so may result in malfunction during use, endoscope damage, and/or patient injury.
- All reusable endoscopic devices must be cleaned and sterilized before initial use as well as before every subsequent use.
- Use endoscopic devices specified by PENTAX Medical whose compatibility has been confirmed. Using endoscopic devices whose compatibility has NOT been confirmed may result in endoscope damage and/or patient injury caused by failure during use.

This section describes the use of a biopsy forceps.

1. Check the entire surface of the forceps for any visible adhered material.
2. Check the insertion portion and control body of the biopsy forceps for abnormalities such as wrinkles, scars, sharp edges, clouding of the surface, dents, catching, protrusions, attachment of foreign materials, falling of parts, etc.


Figure 3.28

3. Check that the cups of the biopsy forceps open/close smoothly by operating its handle.


Figure 3.29

4. Form a loop with a diameter of 20 - 30 cm with the flexible shaft at approximately 20 - 30 cm from the tip of the insertion portion of the biopsy forceps. Check that the cups of the biopsy forceps open/close smoothly by operating its handle.


Figure 3.30

5. Check that the cups align with each other when closed.


Figure 3.31

Attachment of accessories

Warning

- Attach the accessories properly to the endoscope. Failure to do so may result in lowered function and potential reflux or dispersal of patient's body fluids, posing a risk of infection.

■ Attachment of the air/water feeding valve (OF-B188) and suction control valve (OF-B120)

Caution

- Ensure to apply silicone oil lubricant (OF-Z11) onto the O-ring of each valve and the rubber seal of the suction control valve (OF-B120). Using the valves without applying the oil or applying a silicone oil other than the specified one could deteriorate the functions and may result in damage to the endoscope and/or patient injury.
- Attach the air/water feeding valve (OF-B188) and suction control valve straight to the cylinder of the endoscope. Pressing them at an angle may result in damage to the O-ring and check valve.


Figure 3.32

1. Apply a small amount of silicone oil lubricant (OF-Z11) onto the O-rings of the air/water feeding valve (OF-B188) and suction control valve (OF-B120), and the rubber seal. Place a small droplet of oil onto a sterile gloved forefinger, gently swirl the oil between the thumb and the forefinger, and apply it onto the necessary parts. Wipe off the excess lubricant with soft gauze.
2. Attach the air/water feeding valve to the air/water feeding cylinder of the endoscope.


Figure 3.33

3. Ensure that the air/water feeding valve (OF-B188) is firmly attached. Press down the air/water feeding valve a few times to ensure that it moves smoothly.



Figure 3.34

4. Align the metal tab on the shaft of the suction control valve (OE-B120) with the notch on the suction cylinder of the endoscope.

- (1) Notch
(2) Metal Tab

Figure 3.35

5. Attach the suction control valve to the suction cylinder of the endoscope.


Figure 3.36

6. Check that the suction control valve is firmly attached. Press down the suction control valve a few times to ensure that it moves smoothly.


Figure 3.37

■ Attachment of the inlet seal (OF-B190)

1. Attach the inlet seal (OF-B190) to the instrument channel inlet.


Figure 3.38

2. Ensure that the inlet seal is tightly attached to the instrument channel inlet without gaps.


Figure 3.39

■ Attachment of the water jet check valve adapter (OE-C12) and water jet connector cap (OF-B118)

1. Attach the water jet connector cap (OF-B118) and water jet check valve adapter (OE-C12) to the water jet port of the endoscope.


Figure 3.40

2. Ensure that the water jet check valve adapter is firmly attached to the water jet port without gaps. (Close the lid of the water jet connector cap.)


Figure 3.41

■ Attachment of the distal hood (OE-A58/OE-A59)

Warning

The distal hood (OE-A58/OE-A59) is provided non-sterile and non-disinfected. Prior to use, it must be subjected to appropriate cleaning, high-level disinfection, and/or sterilization processes according to the respective IFUs.

Caution

Lightly hold the distal end of the endoscope when attaching/removing the distal hood (OE-A58/OE-A59).

Do NOT apply excessive force to hold the distal end of the endoscope as it may deform the bending section and damage the endoscope.

Note

- The distal hood OE-A58 can be attached to EG-2990Zi, and the distal hood OE-A59 can be attached to EC-3890Zi series. By attaching the distal hood, the distance with the observed section is kept constant, making it easier for magnified observation.
- The distal hood may appear in the monitor in the non-magnified state.

1. Check the distal end of the endoscope for presence of lubricant or water.
2. Hold the distal end gently and push the distal hood into the flange portion inside.


Figure 3.42

3. Ensure that there is no gap between the flange portion of the distal hood and the distal end of the endoscope.

3-4. Inspection and connection of ancillary equipment to the endoscope

Inspect the ancillary equipment prepared in “3-1. Preparation of the equipment”, such as the video processor, monitor, and suction source, according to their respective manuals.

Video processor
Monitor
Water bottle assembly
Endoscopic device
Suction source, etc.

Inspection of the video processor

Only use compatible PENTAX Medical video processors.

For compatible video processors, refer to “Compatible products” (p. 7) or “System chart” (p. 81).

For details on the preparation and inspection of the video processor, refer to the IFU of the respective video processor.

Connection of the endoscope and ancillary equipment

■ Connection to the video processor

1. Ensure that all ancillary equipments are turned off.
2. Ensure that the endoscope locking lever is placed in the “OPEN” position. Securely insert the endoscope electrical contacts and light guide plug into the processor connector and receptacle.


Figure 3.43

3. Turn the endoscope locking lever to the “LOCK” position.


Figure 3.44

4. Lightly pull the PVE connector to ensure that it is connected securely.


Figure 3.45

■ Connection of the water bottle assembly, suction tube, and irrigation tube (OF-B113)

Warning

- Use only sterile water in the water bottle assembly. Failure to do so may pose a risk of infection.
- Connect the suction tube of the suction source firmly to the suction nipple. Failure to do so may result in disconnection of the suction tube during use and pose a risk of cross contamination to the user as a result of reflux or dispersal of patient's body fluids.

Caution

Do NOT use defoaming agents in the water bottle assembly. Such agents cling to the inside channel of the endoscope. A blocked channel can lower the air/water feeding function and cause damage to the endoscope.

Note

Turn off the air/water feeding pump of the video processor beforehand.

1. Attach the water bottle assembly correctly according to the IFU of the video processor.
2. Insert the air/water connector of the water bottle assembly into the air/water port of the endoscope until it clicks.


Figure 3.46

Note

Failure to connect the water bottle assembly correctly not only lowers the air/water feeding function, but may also cause insufficient cleaning of the objective lens.

3. Connect the suction tube of the suction source to the suction nipple of the endoscope.


Figure 3.47

4. Remove the lid of the water jet connector cap (OF-B118) and push the irrigation tube (OF-B113) into the water jet check valve adapter (OE-C12) until it clicks.


Figure 3.48

Caution

Have the irrigation tube (OF-B113) exactly parallel to the water jet check valve adapter (OE-C12) when attaching or removing it. Attaching or removing it diagonally may break the irrigation tube.

Note

Do not use the irrigation tube (OF-B113) if you have any difficulty in attaching or you could not feel a click feeling when attaching it to the endoscope. The use of damaged luer connector may result in water leakage from the connected part or tube disconnection.

3-5. Inspection of the endoscopic system

Inspection of the endoscopic image

Caution

- Do NOT look directly at the light emitted from the distal end of the endoscope or the video processor unit. The intense light may cause eye injuries. Turn off the lamp when looking directly at the distal end of the endoscope.
- Do NOT attach or remove the PVE connector while the video processor power is turned on. Doing so may damage the endoscope.

Note

- If the video processor connected has scope eject function, the PVE connector can be removed when the power of the video processor is turned on by using the scope eject function of the video processor. For details, refer to the IFU of the respective video processor.
- The video processor shown in the Figure is just an example. For details on the operation of each video processor, refer to the IFU of the respective video processor.

1. Turn on the video processor by pushing the Power Switch.


Figure 3.49

2. Turn on the lamp button of the video processor control panel.
3. Ensure that the lamp is lit and light is being emitted from the distal end of the connected endoscope. After pushing the lamp button, it takes several seconds for the lamp to light up.


Figure 3.50

4. Check that the endoscopic image is clear and is displayed normally.

Note

In case the image is not clear, gently clean the endoscope objective lens with clean gauze moistened with 70%–90% medical grade ethyl or isopropyl alcohol.

5. On the control panel of the video processor, check that the exposure control is set to [AUTO].
6. While checking the image displayed on the monitor and following the IFU of the respective video processor, adjust the brightness level as appropriate.
7. Adjust the white balance by following the IFU of the respective video processor.
8. While observing the palm of your hand, check that the brightness adjustment is made normally by moving your palm closer to and away from the distal end of the endoscope.


Figure 3.51

Caution

Do NOT directly touch the distal end of the endoscope (particularly the light guide) for a prolonged period of time when the light is being emitted. Doing so may result in burn injury.

9. Operate the angulation control knobs of the endoscope to move the bending section, and check if the image of the intended direction is displayed by corresponding to the angulation of the bending section. Also check for abnormalities such as appearance of noise in the endoscopic image or disappearance of the image.


Figure 3.52

Inspection of the remote buttons

Warning

Always inspect the remote buttons even if they are NOT expected to be used. During a procedure, the endoscopic image may freeze or other abnormalities may occur, which may result in patient injury.


Figure 3.53

1. Press each remote button.
2. Check that the function assigned to each remote button is operating normally.

Inspection of the optical magnification function


Figure 3.54

1. When the endoscope is connected to the processor, the lens position initialization starts automatically. When the initialization is finished, the standard (non magnified) image will appear on the monitor.
2. Check that turning the magnification control lever on the control body clockwise zooms the image, and turning it counter-clockwise returns the display to the standard image size.
3. Check that the magnification changes one by one with each operation of the magnification control lever, and that the image keeps the same state when the finger is removed from the lever.

Caution

When the monitor image stays zoomed after connecting with the video processor, remove the endoscope from the processor once, make sure that the bending section and insertion section are in the straight state, and then connect it again. The initialization may NOT have finished correctly. If the image remains zoomed, stop using this product and contact your local PENTAX Medical service facility.

Note

The optical magnification function and electrical magnification function can be used together in this product, however, the magnification control lever of the 90Zi series can be used only for the optical magnification function. The magnification control lever cannot be used for the electrical magnification function.

Inspection of the air/water feeding function

Warning

Use sterile water for inspection of the air/water feeding function. Failure to do so may pose a risk of infection.

Note

Refer to the separate IFU of water bottle assembly for details of the operating procedure.

1. Set the A/W-drain lever of the water bottle assembly at the "A/W" position.


Figure 3.55

2. Press the pump button on the control panel of the video processor.
3. Set the pump level to "5" using the pump level setting button on the control panel of the video processor.
4. Insert the distal end of the endoscope into a container filled with sterile water, and check that air bubbles are not continuously discharged from the air nozzle at the distal end of the endoscope. (The example shown in the figure is of the EC-3890Zi series)


Figure 3.56

Warning

If air bubbles are continuously discharged from the air/water nozzle at the distal end of the endoscope when the hole on the top of the air/water feeding valve is NOT closed, stop use immediately and replace the air/water feeding valve with a new one. Continuous use of an air/water feeding valve with abnormalities could cause unintended continuous air feeding and pose a risk of pain or perforation to the patient.

5. Block the hole in the top of the air/water feeding valve. Check whether air bubbles come out vigorously from the air/water nozzle of the distal end of the endoscope.


Figure 3.57

6. Check that the discharge of air bubbles stops when releasing the finger from the air/water feeding valve.


Figure 3.58

7. Pull the endoscope out of the container, and depress the air/water feeding valve. Check that a certain amount of water flows out from the air/water nozzle. (It takes a few seconds until water comes out the first time.)


Figure 3.59

8. Remove the finger from the air/water feeding valve. Check that the air/water feeding valve returns to the original position smoothly and that the water stops at the same time as the finger is removed from the valve.


Figure 3.60

Caution

Do NOT attempt to clear the air or water nozzles with a needle or any other sharp object if nozzle blockage is suspected. It may result in lowered performance as well as damage to the endoscope.

Note

Do not attempt to correct the blockage or continue to use the endoscope, if air/water cannot be fed smoothly and blockage in the nozzle or channel of the endoscope is suspected. In such a case, contact your local PENTAX Medical service facility to have the endoscope repaired.

Inspection of the irrigation function

Warning

Use sterile water for inspection of the irrigation function. Failure to do so may pose a risk of infection.

1. Fill a syringe with sterile water.
2. Put the distal end of the endoscope to the clean container, attach the syringe filled with sterile water, and insert it into the inlet seal (OF-B190).


Figure 3.61

3. Check that sterile water comes out of the opening of the instrument channel when the syringe is pressed to flush the channel. Moreover, check that no foreign materials come out.
4. Remove the syringe from the inlet seal (OF-B190).
5. Fill the syringe with air and insert it into the inlet seal.
6. Flush the sterile water remaining inside the channel by pressing the syringe.
7. Remove the syringe from the inlet seal.

Inspection of the suction function

Warning

Use sterile water for inspection of the suction function. Failure to do so may pose a risk of infection.

Note

Before inspecting the suction function, attach the cap to the inlet seal. Failure to do so may cause lowered suction function.

1. Turn on the suction source and adjust to an appropriate pressure setting.
2. Insert the distal end of the endoscope into a container filled with sterile water and press the suction control valve (OF-B120). Check that water is being suctioned up. (The example shown in the figure is of the EC-3890Zi series)


Figure 3.62

3. Check that when the suction control valve is released, it smoothly returns to the initial position and the suctioning stops.


Figure 3.63

4. Repeat steps 2 and 3 several times to check that there is no water leakage from the suction control valve or the inlet seal.
5. Pull the distal end of the endoscope out of the container. Press the suction control valve, and suction air in order to remove the water remaining inside the instrument channel.

Inspection of the instrument channel

Use a biopsy forceps for inspection of the instrument channel.

Prepare a biopsy forceps which has been cleaned and sterilized according to the manual provided with that product and ensure to perform a pre-use inspection.

Warning

Do NOT use the endoscope if you feel a significant resistance when inserting a biopsy forceps. The inside of the channel may be damaged and it may result in unforeseen events to patients and/or medical professionals.

Caution

- Slowly and gently insert and withdraw the forceps from the inlet seal (OF-B190). Applying strong force may cause endoscope damage.
- Keep the endoscope bending section as straight as possible when inserting the forceps. When the endoscope bending section is greatly angulated, it may NOT be possible to insert the forceps.

1. Close the biopsy forceps cups by operating its handle.


Figure 3.64

Note

Do not close the biopsy forceps cups tightly. Doing so may make its insertion into the instrument channel difficult.

2. Insert the biopsy forceps into the inlet seal (OF-B190). When the cups are first passed through the inlet seal, temporary resistance will be encountered. Hold the shaft at approximately 5cm from the cups and push the biopsy forceps through.


Figure 3.65

3. Slowly advance the biopsy forceps and check that its tip appears from the distal end of the endoscope. Moreover, check that no foreign materials come out.
4. Check that the biopsy forceps can be smoothly withdrawn from the inlet seal.

Inspection of the water jet feeding function

Warning

Use sterile water for inspection of the water jet feeding function. Failure to do so may pose a risk of infection.


1. To use the irrigation pump, prepare to feed sterile water by following the instructions for use for the irrigation pump.
2. Check the entire surface of the irrigation tube (OF-B113) for abnormalities such as bending/breakage/looseness of the connector, cut/chip of the O-ring, buckling/deterioration/hardening of the tube, and/or broken luer.
3. Open the water jet connector cap (OF-B118) and push tube for water pump irrigation or the irrigation tube (OF-B113) in until the water jet check valve adapter (OE-C12) clicks.

Note

Do not use the irrigation tube (OF-B113) if you have any difficulty in attaching or you could not feel a click feeling when attaching it to the endoscope. The use of damaged luer connector may result in water leakage from the connected part or tube disconnection.

4. Feed sterile water using the irrigation pump or syringe filled with sterile water attached to the luer connector of the irrigation tube (OF-B113).

- (1) Water Jet Connector Cap (OF-B118)
- (2) Irrigation Tube (OF-B113)
- (3) Water Jet Check Valve Adapter (OE-C12)

Figure 3.66

Note

When the irrigation tube (OF-B113) is connected with locking type of luer connector, ensure that the luer connectors are properly locked. Do not use the irrigation tube (OF-B113) if there is damage on the luer connector of the irrigation tube (OF-B113) and/or luer connection is not properly locked.

5. Check that a certain amount of water flows out forward from the water jet nozzle at the distal end of the endoscope. (It takes a few seconds until water comes out the first time.) (The example shown in the figure is of the EC-3890Zi series)


Figure 3.67

6. Check that there is no water leakage from the connection between the water jet port of the endoscope and the water jet check valve adapter (OE-C12), or from the connection between the water jet check valve adapter and the irrigation tube (OF-B113).


Figure 3.68

4

Directions for use

This endoscope should only be used by a physician authorized by the medical safety administrator at each medical facility to perform endoscopy.

The device should never be used by individuals who are not licensed medical professionals or used at facilities other than medical facilities.

This section describes the essential information, such as operating procedures and handling precautions, on using this endoscope safely and effectively. This IFU does not describe specific endoscopic procedures. The specific procedures should be determined according to the discretion of a medical professional.

Warning

- Users as well as the assisting personnel should always wear protective equipment (e.g., gloves, goggles, masks, medical gowns, etc.) to minimize the risk of infection, as the patient's body fluids may be dispersed from endoscope components such as the instrument channel inlet and the suction control valve.
- Immediately stop the endoscopic procedure and slowly and cautiously withdraw the endoscope when any abnormality, such as roughness, is felt during angulation. NEVER forcefully turn the angulation control lever. Continuing to use the endoscope with any abnormality may result in endoscope damage and/or patient injuries, including bleeding and perforation.
- NEVER withdraw the endoscope while the bending section is angulated. Doing so may result in patient injury.
- Always check the endoscopic image during endoscope angulation, air/water feeding, and suctioning, use of endoscopic devices, and endoscope insertion and withdrawal. Ensure that these operations are performed in the normal (non-frozen, non-magnified) mode. Endoscope operation in the freeze or magnification mode may result in damage to the endoscope and patient injury.
- Do NOT forcefully insert and withdraw the endoscope. Doing so may result in patient injury.
- When inserting the endoscope into narrow lumina, Do NOT rotate the endoscopic image using the video processor. Doing so may result in patient injury or make the withdrawal of the endoscope impossible.
- Immediately stop the endoscopic procedure if the endoscopic image disappears unexpectedly because of blackout and/or damage to the lamp, video processor, and/or endoscope. Slowly withdraw the endoscope following the instructions in "5-1. Withdrawal of an endoscope with an abnormality" (p. 74). Continuing to use the endoscope may result in patient injury.
- Set the brightness to the minimum necessary. Maintain an appropriate distance between the distal end of the endoscope and the mucosa in order to avoid prolonged illumination of the mucosa. The temperature at the distal end of the endoscope may exceed 41°C and even reach 50°C due to the light emitted from it. This may result in mucosal injury to the patient.

Caution

- Do NOT look directly at the light emitted from the endoscope or direct it at the eyes of other individuals as the intense light may cause eye injuries.
- Do NOT use the endoscope with suspected adherence of patient's body fluid, blood, etc., on the light guide, as this causes the observation image to become dark. The temperature at the distal end of the endoscope may increase and result in mucosal injury to the patients.
- Use the minimum pressure necessary for suctioning. Do NOT suction from the mucosa for a prolonged period of time. Doing so may result in patient injury.
- Do NOT use a water supply device that can exert 30kPa or greater of water pressure to the suction channel (suction valve) during endoscopic examination.

Note

- Prior to a procedure, remove any debris or secretions from the observation area as much as possible to obtain a clear image.
- The objective lens may be cleaned during a procedure by performing air/water feeding and suctioning simultaneously or alternately.

4-1 . Preparation immediately before insertion of the endoscope

Perform appropriate patient preparation for endoscopy as necessary.

Caution


Do NOT spray or wipe the surface of the endoscope insertion portion with an anesthetic (particularly anesthetic spray containing alcohol) or non-medical lubricant (such as petroleum jelly). Doing so could cause cracking or peeling of the external surface of the insertion portion and may result in endoscope damage.

1. Apply a medical grade lubricant to the insertion portion, as necessary.
2. Place a bite block (OF-Z5) into the patient's mouth (only for EG-2990Zi).

Note

- Do not apply lubricants to the objective lens for getting clear observation images.
- When using lens cleaner, ensure to follow the instructions of that product.

(1) Lubricant

Figure 4.1

4-2 . Insertion and observation

Insertion of the endoscope

Caution

Do NOT forcefully bend the strain relief boot as shown (Figure 4.2). Doing so may result in endoscope damage.


Figure 4.2

Slowly and cautiously insert the endoscope.

Adjust the brightness as appropriate for observation with the video processor.

Note

Clear images cannot be obtained if any foreign material is attached to the objective lens or the light guide. Continued use of the light guide with any foreign material attached to it might cause visible steam-like vaporization associated with water vaporization of the organic material heated by the light.

If this vapor is observed, stop the procedure immediately and withdraw the endoscope from the patient. Using clean gauze, clean off any foreign material that has attached and then resume endoscopy.

Angulation operation

Warning

Immediately stop the endoscopic procedure and slowly and cautiously withdraw the endoscope when an abnormality, such as roughness, is felt during angulation. NEVER forcefully turn the angulation control knob as it may result in endoscope damage and/or patient injuries, including bleeding and perforation.

1. Slowly and cautiously operate the angulation control knobs according to the position of the endoscope.
2. Turn the up/down angulation lock lever and right/left angulation lock knob to hold the bending angle of the distal end of the endoscope, as necessary.

Air/Water feeding

Warning

Be careful NOT to feed too much air and to properly control air insufflation into the body cavity. Excessive air insufflation into the patient's body cavity may pose a risk of pain, embolism, and perforation to the patient.

1. Set the appropriate pump level using the pump level setting button on the control panel of the video processor.
2. Cover the hole on top of the air/water feeding valve with a finger to feed air through the air/water nozzle at the distal end of the endoscope.
3. Press in the air/water feeding valve to feed water from the air/water nozzle onto the objective lens.

Suction

Warning

- Do NOT aspirate solid materials as it may cause a clogging in the suction control valve and/or the suction channel.
- If the instrument/suction channel becomes blocked or clogged due to the accumulation of debris, an accessory that cannot be removed, or other cause, do NOT attempt to correct the blockage or continue to use the endoscope. In such a case, contact your local PENTAX Medical service facility to have the endoscope repaired. The use of an endoscope with a blocked internal channel may result in ineffective reprocessing and/or the introduction of debris and/or device components into a patient during a subsequent procedure, causing patient injury and/or posing a risk of cross contamination.
- The use of an endoscope with a blocked internal channel may result in ineffective reprocessing and/or the introduction of debris and/or device components into a patient during a subsequent procedure, causing patient injury and/or posing a risk of cross-contamination.
- Securely attach the cap to the inlet seal. Failure to do so may result in lowered suction function as well as potential reflux or dispersal of patient's body fluids, posing a risk of infection.

Caution

- Do NOT use a cleaning brush or biopsy forceps to remove a foreign object that has occluded the suction channel. This may result in damage to the channel.
- Observe these precautions when suctioning. Failure to do so may result in mucosal injury to the patient.
 - Do NOT apply excessive suction pressure.
 - Maintain distance between the distal end of the endoscope and the mucosa to ensure that the instrument channel opening of the distal end of the endoscope does NOT suction the mucosa.
 - Immediately stop suctioning if the mucosa is suctioned. Do NOT suction mucosa for a prolonged period of time.
 - Stop use immediately when any abnormality is suspected of the suction control.

Suction fluid from inside the body cavity through the instrument channel by pressing the suction control valve.

Water jet feeding

Warning

Use sterile water for water jet feeding. Failure to do so may pose a risk of infection.

Caution

Use minimum pressure for water feeding while observing the condition of the patient's mucosa. Water feeding with the excessive pressure may result in mucosal injury to the patient.

Use the irrigation pump by following its operation manual for the irrigation pump or by attaching a syringe to the luer connector of the irrigation tube (OF-B113) and delivering water into it.

Note

When the irrigation tube (OF-B113) is connected with locking type of luer connector, ensure that the luer connectors are properly locked. Do not use the irrigation tube (OF-B113) if there is damage of the luer connector of the irrigation tube (OF-B113) and/or luer connection is not properly locked.

4

Directions for use

Remote control

Caution

Do NOT apply strong force to the remote button from its side or in an oblique direction, as the button may get stuck and become inoperable.

Operate the remote button for image capture, hardcopy, VCR recording, etc., as necessary.

Note

Leaving a finger on the remote button may result in unintentional pressing of the remote button, causing it to operate.

Optical Magnifying System

Warning

Ensure to insert the endoscope after turning the magnification control lever counterclockwise to return the standard size image. Under the magnified vision the field of view is limited, which worsens insertion operability and could result in a health hazard to the patient.

Note

The optical magnification system gives change in the focal distance and the depth of field. Adjust the distance between the distal end of the endoscope and the observation area to obtain the optimum image. With the maximum magnification rate, the optimum image can be obtained with an observation distance of 2 to 3 mm from the observation area.

4

Directions for use


Figure 4.3

1. Control the optical magnification system by the magnification control lever on the control body. Turning it clockwise magnifies the image, and turning it counterclockwise returns the image to the normal rate.
2. The operation speed of the magnification system can be adjusted in three stages: Low, Medium, and High. For how to change the operation speed, refer to the IFU of the respective video processor.
3. The figure below shows the monitor screen when the optical magnification system is in use. The information for the magnification is displayed in the lower right of the monitor screen. The "+" position in the magnification information shows the approximate magnification rate. On the normal rate of the image the mark "+" is positioned in the left end, and moves to the right with the magnified operation. With the maximum magnification rate, the mark "+" is displayed in the right end. When using EPK-i7000, the mark "▼" moves from left to right.


Figure 4.4

Note

The optical magnification system is also available using the foot switch. Each pressing of the foot switch can switch the magnification rate in four steps, from the normal rate (or the largest magnification rate) to the largest magnification rate (or the normal rate). For how to allocate the magnification system to the foot switch, refer to the IFU of the respective video processor.

The magnification system allocable to the foot switch can be allocated to the remote buttons 1 to 3 of the endoscope. For how to allocate the magnification system to the remote button of the endoscope, refer to the IFU of the respective video processor.

Electrical magnifying system

The electrical magnification system is available when this magnification system is allocated to the control panel, keyboard, or foot switch of the video processor, as well as to the remote buttons of the endoscope. For how to allocate the magnification system to the remote buttons of the endoscope, refer to the IFU of the respective video processor.

4-3. Using an endoscopic device

4

Directions for use

Warning

- NEVER use an endoscopic device showing signs of damage and/or operational abnormality. Doing so may result in endoscope malfunction or damage and/or patient injury.
- All reusable endoscopic devices must be cleaned and sterilized before initial use as well as before every subsequent use.
- Before using the endoscopic device, check its compatibility with the endoscope and read and understand the respective IFU of the endoscopic device. Incorrect use of the endoscopic device may result in damage to the endoscopic device and patient injury.
- When inserting or withdrawing the endoscopic device, ensure that its distal tip is closed or retracted within the sheath. Straighten the endoscopic device and slowly withdraw it. Failure to do so may result in inlet seal damage and/or falling of the broken inlet seal particle(s) into the patient's body cavity.
- Ensure that the distal tip of the endoscopic device is adequately projecting from the distal end of the endoscope before operating it. Failure to do so may result in damage to the instrument channel and/or falling of the broken instrument channel particle(s) inside the patient's body cavity.
- After the endoscopic device is inserted into the inlet seal, NEVER let it hang down. Ensure that the endoscopic device is supported with a hand and no load is applied to the inlet seal. Failure to do so may result in lowered suction function as well as potential reflux or dispersal of patient's body fluids, posing a risk of infection.
- Use only compatible endoscopic devices specified by PENTAX Medical. Using non-compatible endoscopic devices NOT specified by PENTAX Medical may result in clogging and/or damage to the instrument channel and/or endoscopic device. If a liquid such as sterile water or physiological saline is injected with a syringe from the instrument channel inlet while the instrument channel is clogged, the suction control valve may detach and result in potential reflux or dispersal of patient fluids, posing a risk of infection.
- Immediately stop the endoscopic procedure if the endoscopic device cannot be withdrawn from the endoscope. Do NOT attempt to forcefully withdraw the endoscopic device. Slowly and cautiously withdraw the endoscope in which the endoscopic device is inserted. Failure to do so may result in damage to the endoscopic device and/or instrument channel as well as potential reflux or dispersal of patient's body fluids, posing a risk of infection.

Caution

- Do NOT forcefully insert the endoscopic device when the instrument channel is clogged, as this may result in damage to the endoscope.
- Constantly check the endoscopic image while cautiously inserting and withdrawing the endoscopic device.
- Keep the endoscope bending section as straight as possible when inserting and withdrawing the endoscopic device. Forcefully inserting and withdrawing the endoscopic device may result in damage to the instrument channel and endoscopic device and/or patient injury.

Note

The minimum instrument channel width is found on the model name label.

Insertion and operation of the endoscopic device

1. Ensure that the distal tip of the endoscopic device is closed or retracted into the sheath. In case of biopsy forceps, operate the forceps to fully close the cups at the tip. There is a certain amount of resistance when inserting for the first time. Hold the shaft at approximately 5 cm away from its end, and press it in.


Figure 4.5

2. Slowly insert the endoscopic device, and check that the distal tip of the endoscopic device is within the field of view.


Figure 4.6

3. Operate the endoscopic device according to the IFU provided with it.

Withdrawal of the endoscopic device

Warning

- Do NOT forcefully withdraw the endoscopic device or in an oblique direction. Doing so may result in lowered suction function caused by inlet seal damage, falling of the broken inlet seal particle(s) into the patient's body cavity, and potential reflux or dispersal of patient's body fluids, posing a risk of infection. When withdrawing the endoscopic device, prevent the dispersal of patient's body fluids by covering the inlet seal with clean gauze, and withdraw the device slowly in a straight direction against the inlet seal.
- Immediately stop the therapeutic procedure if significant resistance is encountered when withdrawing the endoscopic device or if the endoscopic device cannot be withdrawn from the endoscope. Do NOT attempt to forcefully withdraw the endoscopic device. Failure to do so may result in equipment damage. Close or retract the distal tip of the endoscopic device and slowly withdraw the endoscope into which the endoscopic device is inserted.

1. Ensure that the distal tip of the endoscopic device is closed or retracted into the sheath.
2. Slowly withdraw the endoscopic device in a straight direction against the inlet seal.


Figure 4.7

4-4. Using a non-flammable gas

If there is a possibility of an inflammable gas being present within a body cavity, convert the gas to a non-flammable gas using carbon dioxide prior to laser cautery or electrosurgery.

Warning

- Do NOT use non-flammable gas cylinders whose pressure and flow settings cannot be controlled. Set the gas pressure to 49 kPa or less, and the flow to 4 L/min or less. Using a gas cylinder whose settings cannot be controlled or whose settings are uncertain may result in damage to the endoscope and excessive insufflation of gas into the patient's body cavity.
- Be careful NOT to deliver too much gas and to properly control gas delivery into the channel. Excessive insufflation of gas into the patient's body cavity may pose a risk of a pain, embolism, and perforation to the patient.

Caution

- Perform adequate ventilation when using a non-flammable gas in a small room prolong period of time. An elevated CO₂ concentration in the room may pose a risk of negative effects to the physical condition of the patient and/or user.
- Turn off the air/water feeding pump of the video processor before opening/closing the gas cylinder. If the gas cylinder is opened without turning it off, excessive load is applied to the air/water feeding pump of the video processor, which may result in damage to the air/water feeding pump.

Note

Use of the optionally available gas/water feeding valve (OF-B194) is recommended instead of the air/water feeding valve to prevent gas leakage into the room. Use the gas/water feeding valve according to the IFU provided with it.

1. Prepare a gas cylinder and the optionally available gas adapter (OF-G11). Ensure that the gas cylinder valve is closed. Turn off the air/water feeding pump of the video processor.
2. Remove the air/water connector of the water bottle assembly from the air/water port of the endoscope, and connect the gas adapter instead.


Figure 4.8

3. Connect the gas cylinder to the gas adapter (OF-G11).


Figure 4.9

4. Connect the air/water connector of the water bottle assembly to the gas adapter.


Figure 4.10

5. Ensure that all the devices are securely connected before opening the gas cylinder valve.

4-5. Laser cauterization

Warning

- Laser equipment should be used only by experts who have thorough knowledge of the laser equipment and endoscopic laser treatment.
- Before using the laser equipment, thoroughly read the manual supplied with it, and always perform pre-use inspection. Ensure that the laser equipment is ready for use by performing the safety checks specified in the manual.
- Use Nd:YAG laser (wavelength 1064 nm) or a laser with a wavelength of 800–1000 nm only.
- When using laser equipment, both users and the assisting personnel should wear goggles. Failure to do so may result in eye injuries.
- Do NOT use the laser equipment in flammable surroundings, such as an oxygen-rich environment. If there is a possibility of a flammable gas being present within a body cavity, convert the gas to a non-flammable gas prior to laser cauterization. Using the laser equipment in flammable surroundings may result in combustion or an explosion.
- Set the laser output to the minimum level necessary.
 - If the laser is continuously emitted at a high level, the endoscopic image may become white (whiteout). Do NOT perform laser cautery during whiteout, as it may result in patient injury.
 - Continuously emitting the laser at a high level may damage the instrument.
- Maintain an adequate distance between the distal end of the endoscope and the patient's body cavity wall. Before activation of the laser, ensure that the distal tip of the laser probe emerges from the distal end of the endoscope. Failure to do so may result in instrument damage and patient injury.

4

Directions for use

1. Insert the laser probe into the inlet seal (OF-B190) as described in "4-3. Using an endoscopic device".
2. Operate the laser probe according to the manual provided with it.
3. When the procedure is complete, withdraw the laser probe via the inlet seal as described in "4-3. Using an endoscopic device".

Note

- It is normal for the guide beam to appear white in the video endoscopic image.
- If the distal end of the endoscope is within 20 mm of the irradiated target surface, the guide beam may create a smear in the image. If such a smear affects visibility of the target, decrease the intensity of the light.
- When operating the laser at a high power and/or if the distal end of the endoscope is moved within 10 mm of the irradiated target, a flare may appear at the corners of the image.


Figure 4.11

4-6. Electrosurgery

4

Directions for use

Warning

- Thoroughly read the manual provided with the high frequency generator and device before using them, and always perform a pre-use inspection. Ensure that the high frequency generator and device are ready for use by performing the safety checks specified in the manual. Use in combination with the electrosurgical device could result in increased leakage of current to the patient.
- Ensure that the released energy from the high-frequency does not affect the peripheral device such as pacemaker and to use the minimum necessary output level of high-frequency when using it near the heart. It may stimulate the heart.
- Do NOT use the high frequency generator in flammable surroundings, such as an oxygen-rich environment. If there is a possibility of a flammable gas being present within a body cavity, convert the gas to a non-flammable gas prior to electrosurgery. Using the high frequency generator in flammable surroundings may result in combustion or an explosion.
- Ensure that the active portion of the endoscopic device does NOT come in contact with the peripheral tissues, as it may result in patient injuries.
- Set the high-frequency output level and waveform mode appropriately according to usage. Minimize the energization time, as it can result in patient injury.
- Before using an electrosurgical device, check the entire surface of the endoscope for any abnormalities such as cracks and exposure of internal metals. Failure to do so may result in burns from high-frequency current.

Caution

- Users as well as the assisting personnel should always wear insulated gloves. Failure to do so may result in burns from high-frequency current.
 - High frequency generators may be of the floating (Type BF or Type CF) or non-floating (Type B) types. Only use floating-type high frequency generators to avoid patient and user burns.
 - There are two types of floating-type high frequency generators: those with an endoscope feedback cord (S-cord) and those without.
 - 1) High frequency generator with an S-cord: Connect the S-cord to the feedback terminal located on the PVE connector of the endoscope.
 - 2) High frequency generator without an S-cord: Connect the condenser earth cable (OL-Z4/OL-Z3) between the endoscope feedback terminal and the video processor potential equalization terminal.
- Follow the directions for each type described above, as failure to do so may result in burns from high-frequency current.
- Only use insulated devices. Failure to do so may result in burns from high-frequency current.
 - Ensure that the distal tip of the electrosurgical device is adequately projecting from the distal end of the endoscope before operating it. Failure to do so may result in endoscope damage.
 - During use, follow the precautions below, as failure to do so may result in endoscope damage, burns, and/or mucosal injury.
 - Maintain an adequate distance between the distal end of the endoscope and the insulated tip and active portion of the endoscopic device. Ensure that the distal tip of the endoscopic device is adequately projecting from the distal end of the endoscope before operating it.
 - Users and the assisting personnel should NOT touch the patient during device use.
 - Turn on the high frequency generator just before the procedure and turn it off immediately after the procedure.

1. Insert the electrosurgical device into the inlet seal as described in “4-3. Using an endoscopic device”.
2. Operate the electrosurgical device according to the IFU provided with it.
3. When the procedure is complete, withdraw the electrosurgical device via the inlet seal as described in “4-3. Using an endoscopic device”.

4-7 . Withdrawal of the endoscope

Warning

- When withdrawing the endoscope, prevent dispersal of patient’s body fluids by holding clean gauze along the insertion portion. Failure to do so may pose a risk of infection.
- Before withdrawing the endoscope, do NOT remove the water bottle assembly from the video processor while the water bottle assembly is connected to the endoscope. Doing so may result in reflux of patient’s body fluids into the water bottle assembly.
- Do NOT withdraw the endoscope while the bending section is angulated. Doing so may result in patient injury.

1. Operate the suction control valve to suction any fluid remaining inside the patient’s body cavity.
2. If the electrical magnifying function used, set it back to standard image size.
3. Turn the up/down angulation lock lever and right/left angulation lock knob in the “F►” direction until they stop to release the lock of the angulation control knobs.
4. While checking the endoscopic image, slowly and cautiously withdraw the endoscope.
5. Remove the bite block from the patient’s mouth. (Only for EG-2990Zi)
6. Turn the video processor lamp off.

4

Directions for use

4-8. Care after use

Caution

Do NOT touch the light guide plug and electrical contacts after use. It could result in burns.

■ Endoscope:

Perform cleaning, high-level disinfection, and/or sterilization according to the procedure specified in the separate IFU (Reprocessing) of this endoscope.

■ Accessories:

Air/water feeding valve (OF-B188), suction control valve (OF-B120), inlet seal (OF-B190), bite block (OF-Z5) (Only for EG-2990Zi), water jet check valve adapter (OE-C12), water jet connector cap (OF-B118), irrigation tube (OF-B113), distal hoods (OE-A58/OE-A59), and other optional equipment: Perform cleaning, high-level disinfection, and/or sterilization according to the procedure specified in the respective IFU provided with them.

■ Endoscopic devices:

Reusable endoscopic devices:

All reusable devices must be cleaned and sterilized according to the respective IFUs provided with them.

Single use endoscopic devices:

Follow the national or local laws/guidelines to appropriately dispose of single use endoscopic devices.

■ Video Processors/irrigation pump:

Follow the IFU provided with it for its care after use.

■ Water bottle assembly:

For cleaning and disinfection and/or sterilization of the water bottle assembly, refer to IFU provided with the water bottle assembly.

Disconnecting the endoscope from the video processor

Caution

Do NOT attach or remove the PVE connector while the video processor power is powered on. Doing so may damage the endoscope.

1. Immediately after use, perform pre-cleaning according to the separate IFU (Reprocessing) of this endoscope.
2. After completion of cleaning in the examination room, turn off power of the video processor.


Figure 4.12

3. Turn the endoscope locking lever to the "OPEN" position; then, hold the endoscope PVE connector and remove the endoscope electrical contacts and light guide plug from the processor connector and receptacle.


Figure 4.13

Note

- If the video processor connected has a scope eject function, the PVE connector can be removed when the power of the video processor is turned on by using the scope eject function of the video processor. For details, refer to the IFU of the respective video processor.

5 Troubleshooting

After inspecting the endoscope according to “Chapter 3 Preparation and inspection”, if any abnormality is suspected, do not use the endoscope; send the endoscope for repair according to “5-2. Returning the endoscope for repair” (p. 75).

Warning

Do NOT use an endoscope with any apparent abnormality. Continuing to use an endoscope with an abnormality may result in endoscope damage, malfunction, and/or injury to the patient and/or user.

5-1 . Withdrawal of an endoscope with an abnormality

Immediately stop the endoscopic procedure and slowly and cautiously withdraw the endoscope when any abnormality occurs.

When the endoscopic image is displayed

1. When using the endoscopic device, close the distal tip or retract it within the sheath. Then, slowly withdraw the endoscopic device from the endoscope.
2. Operate the suction control valve to suction any fluid remaining inside the patient's body cavity.
3. If the electrical magnifying function used, set it back to standard image size.
4. Turn the up/down angulation lock lever and right/left angulation lock knob in the “F ►” direction until they stop to release the lock of the angulation control knobs.
5. While checking the endoscopic image, slowly and cautiously withdraw the endoscope.

When the endoscopic image is not displayed

1. When using the endoscopic device, close the distal tip or retract it within the sheath. Then, slowly withdraw the endoscopic device from the endoscope.
2. Turn the up/down angulation lock lever and right/left angulation lock knob in the “F ►” direction until they stop to release the lock of the angulation control knobs.
3. Remove your hand from the up/down and right/left angulation control knobs.
4. Slowly and cautiously withdraw the endoscope.

5-2 . Returning the endoscope for repair

When returning the endoscope for repair, follow the instructions below. For more details, contact your local PENTAX Medical service facility. Always subject the endoscope to cleaning and high-level disinfection before returning it for repair.

Warning

Only qualified personnel from PENTAX Medical are authorized to repair this endoscope. PENTAX Medical is NOT liable for any damage or injury that occurs as a result of repairs attempted by non-PENTAX Medical personnel. It must be recognized that PENTAX Medical does NOT evaluate non-PENTAX Medical parts, components, materials and/or servicing methods and therefore questions regarding material compatibility and/or functionality of PENTAX Medical endoscopes built with these unauthorized, untested and unapproved items, materials, repair/assembly methods must be referred to the third party service organization and/or device remanufacturer.

1. Place this endoscope in the dedicated carrying case. Ensure that the PVE soaking cap is included for water leakage test.
2. When transporting by air, ensure that the ventilation cap is attached to prevent damage to the endoscope.
3. Include any PENTAX Medical accessory that is suspected to be associated with the damage.
4. Contact your local PENTAX Medical service facility for shipping address and inform us the description of failures that need repair, model name, serial number, and name/phone number/address of the contact person.

Disposal

Warning

Follow the national or local laws/guidelines to appropriately dispose of the consumables. Failure to do so may create a risk of cross contamination or infection.

Contact your local PENTAX Medical service facility when disposing the endoscope(s).

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Guidance and manufacturer's declaration – electromagnetic emissions

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Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. This product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycle 70 % U_T (30 % dip in U_T) for 25 cycle <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycle 70 % U_T (30 % dip in U_T) for 25 cycle <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued operation during power mains interruptions, it is recommended that this product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	It is recommended that this product be used apart from other devices operated with large current.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	The recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	The recommended separation distance: $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. d is the recommended separation distance in metres (m).			

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Rated maximum output power of transmitter (W)	Recommended distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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
Note

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Endoscope specifications

Model Name		EG-2990Zi	EC-3890FZi	EC-3890MZi	EC-3890LZi
Direction of view		Forward (0°)			
Field of view		Normal observation (WIDE): 140°, Observation with the maximum zooming (TELE): 49°			
Depth of field		Normal (WIDE): 5-100mm	Normal (WIDE): 4-100mm		
Magnification method		Optical/electrical			
Tip angulation	Up-Down	210°-120°	180°-180°		
	Right-Left	120°-120°	160°-160°		
Rigid distal width		Ø10.6mm [Ø12.6mm]	Ø13.0mm [Ø15.0mm]		
Distal end width		Ø10.2mm [Ø12.0mm]	Ø12.0mm [Ø14.4mm]		
Insertion tube width		Ø9.8mm	Ø13.2mm		
Maximum insertion portion width *1		Ø11.6mm [Ø12.9mm]	Ø14.65mm [Ø15.3mm]		
Minimum instrument channel width *2		Ø2.8mm	Ø3.8mm		
Endoscopic device view on the endoscopic image					
Insertion Tube Working Length *1		1,050mm [1,053mm]	1,500mm [1,503mm]	1,300mm [1,303mm]	1,700 mm [1,703mm]
Total Length		1,376mm [1,379mm]	1,825mm [1,828mm]	1,625mm [1,628mm]	2,025 mm [2,028mm]
Laser cauterization		Available			
Electrosurgery treatment		Available			
Water jet feeding function		Available			



[]: Values when a distal hood is attached

Specifications are subject to change without prior notice and without any obligation on the part of the manufacturer.

*1 There is no guarantee that equipment selected solely using the maximum insertion portion width and insertion portion working length will be compatible when used in combination.

*2 There is no guarantee that equipment selected solely using this minimum instrument channel width will be compatible when used in combination.

Distal End

- (1) Objective Lens
- (2) Light Guide
- (3) Water Jet Nozzle
- (4) Instrument Channel
- (5) Air Nozzle
- (6) Water Nozzle
- (7) Case

System chart

This section shows the system chart (configuration) for this endoscope and the ancillary equipment.

Warning


Use this product in combination only with compatible products shown in “Compatible products” (p. 7) and the “System chart”. Failure to do so may result in lowered function and patient/user injury or damage to the equipment.

Note

When this endoscope is used in combination with other equipment, depending on how it is connected, it may result in malfunction and/or unforeseen events to patients and/or medical professionals.

Pre-use operation check and risk management associated with changes are recommended, particularly when the equipment(s) used in combination is changed, added, or upgraded.

1. Suction Control Valve (OF-B120)
2. Air/Water Feeding Valve (OF-B188)
3. Gas/Water Feeding Valve (OF-B194*) *Optional Equipment
4. Inlet Seal (OF-B190)
5. Water Jet Check Valve Adapter (OE-C12)
6. Irrigation Tube (OF-B113)
7. Gas Adapter (OF-G11*) *Optional Equipment
8. Distal Hood OE-A58 (for EG-2990Zi) OE-A59 (for EC-3890Zi series)
9. Bite Block (OF-Z5) (Only for EG-2990Zi)

Software Version

EG-2990Zi	00B3C-4
EC-3890Zi series	

Contacts

Manufacturer

 **HOYA Corporation**
6-10-1 Nishi-shinjuku,
Shinjuku-ku, Tokyo
160-0023 Japan

Distributors

PENTAX Europe GmbH

Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel: +49 40 561 92-0
Fax: +49 40 560 42 13

PENTAX Medical A Division of PENTAX of America, Inc.

3 Paragon Drive
Montvale, NJ 07645-1782
USA
Tel: +1 201 571 2300
Toll Free: +1 800 431 5880
Fax: +1 201 391 4189

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80566

LCPM: 04/2018/06/35010302 2018. 03 6217001 S069 R01

In the interest of technical progress, specifications may change without notice.

PENTAX
MEDICAL
Excellence in Focus

INSTRUCTIONS FOR USE

EGA-500P
PENTAX Medical
Irrigation Pump

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PENTAX® is a registered trademark of Hoya Corporation.
ENDOGATOR™ is a trademark of Medivators Inc.

IFU-66 Revision N

© 2020 Medivators Inc.

CONTENT

English (EN)	5 - 20
Croatian (HR)	21 - 37
Czech (CS).....	38 - 54
Danish (DA).....	55 - 71
Dutch (NL)	72 - 88
Estonian (EE)	89 - 105
Finnish (FI)	106 - 122
French (FR)	123 - 139
German (DE).....	140 - 156
Hungarian (HU).....	157 - 173
Italian (IT)	174 - 190
Kazakh (KK)	191 - 207
Lithuanian (LT)	208 - 224
Norwegian (NO)	225 - 241
Polish (PL).....	242 - 258
Portuguese (PT).....	259 - 275
Romanian (RO)	276 - 292
Russian (RU).....	293 - 309
Serbian (SR).....	310 - 326
Slovakian (SK)	327 - 343
Slovenian (SL).....	344 - 360
Spanish (ES)	361 - 377
Swedish (SV)	378 - 394
Ukrainian (UA)	395 - 411

Table Of Contents

I. Unpacking And Inspection:	3
II. Introduction:.....	3
Definitions:.....	3
Indications For Use:.....	3
Contra-indications:.....	3
III. Warnings And Cautions	4
IV. Product Features	6
Theory Of Operation	6
V. Technical Specifications	6
Electrical Specifications.....	6
Mechanical Specifications	7
Flow Rates	7
Environmental requirements.....	7
Accessory Items.....	8
VI. Front Panel Controls/Connections	9
VII. Rear Panel Connections	10
VIII. Setup And Operation.....	11
Initial Setup.....	11
Pre-Procedure setup	12
Water Pre-Warming.....	13
Operation	13
Shutdown.....	13
IX. Care and Maintenance	14
Cleaning	14
Maintenance	14
X. Electromagnetic Compatibility (EMC) Information.....	15

Instruction Manual

I. Unpacking And Inspection:

Upon receiving the PENTAX® Medical IRRIGATION PUMP, ensure that the following items are contained in the shipping box:

- Base tray and water bottle holder
- Hardware to assemble base tray and water bottle holder
- Operator's manual
- Power cord
- Foot pedal
- Pump Head

READ THIS MANUAL THOROUGHLY BEFORE PROCEEDING WITH THE OPERATION OF THIS EQUIPMENT.

These instructions should be stored and used for reference as needed. If there are any questions, please contact your local PENTAX service facility.

II. Introduction:

Definitions:

- Throughout this document the PENTAX Medical IRRIGATION PUMP may be referred to as "unit" or "device"
- mL/min – milliliters per minute (flow)
- VAC – Volts Alternating Current (electrical potential)
- kPa – kilo Pascal (pressure)
- W – Watts
- VA – Volt Amps
- Hz – Hertz (frequency)

Indications For Use:

Federal (USA) law restricts this device to sale by or on the order of physician.

The PENTAX Medical IRRIGATION PUMP is indicated for GI endoscopic irrigation with room temperature or warm water for use with washing catheters, integral GI endoscope water jet channels and GI endoscope working channels.

The device contains a heated water bottle cradle that is intended to maintain a sterile water bottle for GI endoscopic irrigation at a maximum temperature set point of 37°C with a tolerance of $\pm 3^{\circ}\text{C}$.

Contra-indications:

The device should only be used for GI endoscopic procedures requiring irrigation of the gastrointestinal tract and should not be used for any other treatments or procedures.

The device is not intended for use with or exposure to magnetic resonance imaging (MRI) systems. Do not use the device in an environment where it could be exposed to MRI systems.

III. Warnings And Cautions

SIGNAL WORDS:

WARNING:

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION:

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.

Safety Symbols:

Attention

Refer To Operating Instructions

Warning: Dangerous Voltage

Equipotentiality

DO NOT Allow Fingers To Contact Moving Parts.

Hot Surface

MR Unsafe

WARNINGS

- A. To reduce risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.
- B. To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth ground.
- C. The PENTAX® Medical IRRIGATION PUMP is not suitable for use in the presence of a flammable anesthetic mixture with oxygen.
- D. The PENTAX Medical IRRIGATION PUMP shall be sold only by prescription for use by physicians/clinicians who are trained regarding the amount of water to use for irrigation and infusion.
- E. Always watch for abnormal flow rates and adjust the flow rate as necessary. If an unusual flow rate remains after adjustments have been made, then stop the procedure and investigate.
- F. Only qualified medical personnel in an acceptable medical facility should operate the PENTAX Medical IRRIGATION PUMP.
- G. The PENTAX Medical IRRIGATION PUMP should be connected to a properly grounded 120V receptacle marked "Hospital Grade" or "Hospital Only," otherwise grounding reliability cannot be achieved.
- H. Extreme precaution must be taken when handling liquids around electrical equipment. DO NOT operate the PENTAX Medical IRRIGATION PUMP if liquid has been spilled on the unit.
- I. The PENTAX Medical IRRIGATION PUMP should not be used adjacent to or stacked with other equipment other than another PENTAX Medical device. Electromagnetic or other interference may occur between the PENTAX Medical IRRIGATION PUMP and other electronic devices. If stacked or adjacent use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.
- J. The PENTAX Medical IRRIGATION PUMP should only be used in conjunction with other equipment whose safety against leakage currents has been established.
- K. The instructions contained in the operating manuals of any equipment to be used in conjunction with the PENTAX Medical IRRIGATION PUMP must be followed to avoid any possible hazard from incompatibility.
- L. The instructions for use described in this manual MUST be followed. Otherwise, compromised safety, malfunction, injury to the operator and/or patient, or costly damage to the unit and other equipment may occur.
- M. The PENTAX Medical IRRIGATION PUMP must be connected to an appropriate power source when loss of power source would result in an unacceptable risk.

CAUTION

- A. If emergency or abnormal function occurs, immediately turn off the power to the unit.
- B. Field-servicing of the PENTAX Medical IRRIGATION PUMP is limited to the replacement of the power cords, water bottle holder and heater assemblies, pump heads, foot pedals and fuses.
- C. Remove power from the device before initiating any field servicing of the replacement parts.
- D. There are no user serviceable parts inside the PENTAX Medical IRRIGATION PUMP unit. Repairs to the PENTAX Medical IRRIGATION PUMP should only be performed by qualified service personnel.
- E. Do not use the device if the enclosure is damaged or enclosure integrity has been compromised.
- F. Do not attempt to operate the device before reading and understanding all sections of this manual.
- G. Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in Section X.
- H. Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment. Do not expose the device to sources of electromagnetic interference such as CT equipment, diathermy equipment, cellular phones, RFID tags and metal detectors.

IV. Product Features

- Built in water bottle warmer maintains water at up to 37°C with an allowable tolerance of $\pm 3^{\circ}\text{C}$
- “Prime” feature runs pump for 20 seconds with automatic timer
- Includes side water bottle clip for traditional water bottle holder
- Safety feature prevents pump from running when pump head is open
- Robust motor design for consistent performance
- Water bottle holder is angled for easier access to water bottle and tubing

Theory Of Operation

The PENTAX® Medical IRRIGATION PUMP works by turning a peristaltic-roller-type pump head to move liquid through a tube set and into a GI endoscope system. The pump head will not operate if the pump head is open, and will cease to operate if the pump head is opened while the motor is activated.

The unit incorporates a “prime” feature which allows the unit to operate for a pre-set period of time. This feature can be canceled while in progress by pressing the foot pedal.

The water bottle warmer is intended to maintain the warmth of the water in a sterile water bottle that has been pre-warmed. It is controlled by dual temperature sensors for redundancy and safety. There is a heating element integrated into the water bottle holder which then radiates heat through the bottle and into the water. The temperature of the water bottle holder is maintained at a maximum of $37^{\circ}\text{C} \pm 3^{\circ}\text{C}$.

V. Technical Specifications

Electrical Specifications

Input Voltage:	100-240 VAC
Input Frequency:	50-60 Hz
Power Consumption:	110 VA
Fuse Rating:	M10AL250V
	Medium acting, 10 amp, low breaking capacity, 250 volt

Replace fuses only with those of same type and rating

Certifications:	IEC-60601-1, IEC-60601-1-2, IEC-60601-2-18
Classification:	Class 1 Type B
Water Bottle Warmer:	20W resistive element with dual temperature sensors.
IP Rating (Ingress Protection):	IP24

WARNING:

Grounding reliability can only be achieved when connected to a receptacle marked “Hospital Grade”.

Mechanical Specifications

Physical Dimensions:

Height	4 ¾" inches	121 mm
Width	7 ¾" inches	197 mm
Depth	13 ¾" inches	349 mm
Weight	10.5 pounds	4.8 kg

Flow Rates

Auxiliary Water Channel	0-300 ml/min*
Biopsy Channel	0-650 ml/min*

*These values are approximate and are based on an average flow obtained with three different manufacturer's GI endoscopes with biopsy channel diameters ranging from 2.8 mm - 4.2 mm and auxiliary water channel diameters ranging from 0.8mm - 0.1 mm. User results may vary depending on scope used, channel size and working channel length.

CAUTION:

The PENTAX® Medical IRRIGATION PUMP is intended for use only by physicians/clinicians who are trained regarding the amount of water to use for irrigation and infusion during GI endoscopic procedures.

Environmental requirements

Operating Temperature:	16° to 24°C (+61° to 75°F)
Operating Relative Humidity:	30% to 75% non-condensing
Operating Pressure:	70 kPa-106 kPa_(10.2 PSI -15.4 PSI)

Accessory Items

The PENTAX® Medical IRRIGATION PUMP, EGA-500P, is intended to be used only with the specific models of compatible tubing sets, connectors and accessories identified in the table below and/or officially recommended by PENTAX. Use of the device with accessories not identified in the table below or officially recommended by PENTAX may result in incompatibility and/or the risk of cross-contamination and infection transmission.

<u>Order Number</u>	<u>Description</u>
100150CO2P, 100150CO2UP (24 hr Use)	PENTAX Medical DispoCap CO ₂
100160P, 100160UP (24 hr Use)	PENTAX Medical DispoCap Air
100130P, 100130UP (24 hr Use)	PENTAX Medical Irrigation Tubing For EGA-500P
200230P, 200230UP (24 hr Use)	PENTAX Medical Irrigation Tubing For EGP-100P
100242P (Single use only)	PENTAX Medical Single Use Waterjet Connector
100116P (24 hr Use)	PENTAX Medical 24 hr Use Waterjet Connector
100551P	PENTAX Medical CO ₂ Source Tubing

Replacement Parts

EGA-500P

<u>Item</u>	<u>Order Number</u>
Power cord	EGA-7014
Water Bottle holder And Warmer Assembly	EGA-7010
Pump Head	EGA-7016
Foot Pedal	47049-942

VI. Front Panel Controls/Connections

NOTE:

All status indicator lights will be green during normal operation.

1. **On/Off Push Button** Turns main power on or off to the unit
2. **Water Bottle Warmer On/Off** Push to turn water warmer on or off
3. **Pump Flow Control** The flow control knob is a continuous variable flow control that has settings between 0 and 100% flow. The lowest knob setting corresponds to no flow while the max flow control setting represents 100% flow. It contains bracketed visual iterations to represent flow magnitude in increasing increments of approximately 10%, however the flow control knob response is not entirely linear. Incremental flow control knob increases between 50-100% have a larger impact on flow rate response as compared to incremental changes between 0-50%.

Flow Rates:

Auxiliary Water Channel	0-300 ml/min*
Biopsy Channel	0-650 ml/min*


*These values are approximate and are based on an average flow obtained with three different manufacturer's GI endoscopes with biopsy channel diameters ranging from 2.8 mm-4.2 mm and auxiliary water channel diameters ranging from 0.8 mm-1.0mm. User results may vary depending on scope used, channel size and working channel length.

CAUTION:

The PENTAX® Medical IRRIGATION PUMP is intended for use only by physicians/clinicians who are trained regarding the amount of water to use for irrigation and infusion during GI endoscopic procedures.

4. **Foot Pedal Connection**
5. **Prime Button**
When pressed, pump will run for a period of 20 seconds. The cycle can be canceled at any time by either pressing the foot pedal or prime button again. See Section VIII for full explanation.
6. **Pump Head**
Refer to Section VIII for assembly instructions

(The appearance of your device may differ slightly from picture shown above)

VII. Rear Panel Connections

1. **Water Bottle Warmer Connection**
Power connection for water bottle warmer
2. **Equipotentiality (Earth Ground)**
3. **AC Power Connection**
AC power input 100-240 VAC, 50-60HZ

WARNING:

Use only the hospital grade power line cord supplied with this unit. Connect only to a power receptacle marked as hospital grade.

4. **Fuse Receptacle**


WARNING:

Replace only with type and rating marked M10AL250V

5. **Water Bottle Holder And Warmer**

The water bottle holder contains the heating element and an integral hold-down bracket to secure the water bottle in place.


(The appearance of your device may differ slightly from picture shown above)

VIII. Setup And Operation

Initial Setup

1. Before beginning, inspect the PENTAX® Medical IRRIGATION PUMP for any signs of damage.
2. Install the pump head onto the device by aligning the flat shaft on the front of the device with the slotted hole in the pump head and rotate the pump head until it locks in place.
3. Assemble the water bottle holder and warmer assembly in the following manner:
 - A. The following supplies are necessary and can be found in a bag attached to the water bottle holder:
 - a.) Short Screws (2)
 - b.) Long Screws (2)
 - c.) Large Hex Key
 - d.) Small Hex Key

- B. Attach the water bottle holder to the base tray using the (2) small screws and small hex key.

- C. Place the PENTAX Medical IRRIGATION PUMP on the base tray and install the (2) long screws using the large hex key.

- D. Adjust the water bottle hold-down bracket to fit the water bottles used by your facility. This is done by loosening the nut on the back and sliding the bracket up or down so it contacts the water bottle and then re-tightening the nut.

4. Plug the water warmer cable into the receptacle marked "water warmer" on the back of the unit.
5. Place the PENTAX® Medical IRRIGATION PUMP on a flat surface such as an endoscope accessory cart or other suitable working surface.
6. Before connecting the power cord to the hospital grade wall receptacle, make sure the power is off to the device, and no accessories are connected.
7. Connect the power cord to the back of the device.

Pre-Procedure setup

1. Place the foot pedal on the floor and insert the foot pedal cord into the designated front panel receptacle.
2. Insert a water bottle containing sterile water into the water bottle holder

CAUTION:

Use of irrigation fluid other than sterile water can pose an infection risk to patients.

Adjust the hold-down bracket by loosening the knob, lowering the bracket until it contacts the top of the bottle and then tightening the adjusting screw. After initial setup, if the same size and type of bottle is used, it should require no further adjustment.

3. Open a sterile pouch of the 100130P PENTAX Medical irrigation tubing. Read and follow its instructions for use.

CAUTION:

The PENTAX Medical IRRIGATION PUMP EGA-500P is intended to be used only with specific models of compatible tubing sets, connectors and accessories identified in the table in Section V or officially recommended by PENTAX. Use of the device with accessories not identified in the accessories table or officially recommended by PENTAX may result in incompatibility and/or the risk of cross contamination and infection transmission.

4. Connect the tubing to the water bottle.
5. Open the pump head and place the tubing into the pump head rollers so that the water flow will be from right to left facing the pump.

Note:

The device is equipped with a safety feature so that the pump will not operate if the pump head is open.

6. Turn the flow rate indicator to its lowest speed and turn the unit on.
7. Perform the following procedure to ensure the tubing is fully filled with sterile water each time before connecting to a GI endoscope:
 - A. Direct the end of the tube into a receptacle capable of catching expelled water.
 - B. Depress the foot switch and turn up the flow rate adjustment to the desired level based on the irrigation connection to the GI endoscope.

PRIME FUNCTION: The prime function allows the user to perform this feature with the assistance of an automatic timer in lieu of pressing the foot pedal. Pressing the prime button will run the pump for a period of 20 seconds. If at any time the prime cycle needs to be stopped the button can be pressed again, or the foot pedal can be pressed.

The prime function is not intended to be a definite means of fully priming the tubing. It is the responsibility of the end user to ensure the tubing and GI endoscope channel are fully filled with water before proceeding.

8. When priming is complete, attach the irrigation tube set to the GI endoscope using the proper adapter for the endoscope model and endoscope channel. (See Section V for a complete list of irrigation tubing and adaptors for various GI endoscopes.)
9. The PENTAX Medical IRRIGATION PUMP is now ready to use.

CAUTION:

NEVER allow the sterile water bottle to empty completely to avoid air entering the irrigation tubing.

Water Pre-Warming

1. The PENTAX® Medical irrigation pump has an integral water warming system that can maintain the temperature of the water at up to $37 \pm 3^{\circ}\text{C}$.

NOTE:

The system is not intended to raise the temperature of the water; therefore pre-warming of the water bottle is required.

2. Place sterile water bottle in a fluid warming oven set to 37°C for at least 2 hours until the water has been warmed to at least $37 \pm 3^{\circ}\text{C}$. Do not exceed this temperature. If warm water is not desired, turn the water warmer on the pump off and place a bottle of room temperature sterile water on the water bottle holder.

NOTE:

If the water warmer has been enabled and has exceeded the upper temperature set point limit, the water warmer's switch indicator light will turn blue and the warmer will automatically shut off. If the switch indicator light is flashing blue, ensure the water warmer cable is plugged into the back of the unit. If it is connected and the light is flashing, this may indicate a problem with the unit. Contact customer service at 800-444-4729.

Read and become familiar with all manufacturer's instructions on warming ovens and water bottles regarding maximum fluid temperatures. Never exceed 40°C during water pre-warming.

CAUTION:

Never use a microwave oven to warm a water bottle as this could heat the water to a dangerously high temperature or could result in uneven warming.

CAUTION:

The surface of the water bottle warmer may be hot to the touch.

Operation

1. Water flow is initiated by pressing the foot pedal and releasing when no longer desired.
2. If more water flow is required, increase the flow control setting on the front panel. Refer to Section VI for flow rates.

Note:

The device is equipped with a safety feature so that the pump will not operate if the pump head is open.

The physician must assess the condition of the patient and use clinical judgment to set the flow rate from the pump to a suitable level to avoid patient trauma. The flow control should always be set at the lowest setting at the start of the procedure and be increased incrementally to a level commensurate with the clinical condition of the patient and the degree of irrigation required.

Shutdown

1. Press the power button to turn the unit off. The light should extinguish.
2. If the warming element has been activated, press the power button for the warmer to turn off. The blue light should extinguish.
3. After the procedures are complete for the day, remove the tubing from the pump head and discard the tubing and water bottle.

Note:

ENDOGATOR™ Irrigation tubing may be used for 24 hours.

IX. Care and Maintenance

Cleaning

- Before cleaning ensure the power to the unit is turned off and the electrical cord unplugged.
- The outside surface of the PENTAX® Medical IRRIGATION PUMP can be cleaned with a damp cloth, 70% isopropyl alcohol solution, or a 10% bleach with water solution as often as deemed necessary.
- Do not use abrasive or harsh cleaners.
- Do not allow liquid to enter the unit.
- To disinfect the outside surface, use a mild disinfectant according to the manufacturer's instructions.
- Do not sterilize the unit by any means.

Maintenance

Pump Head Maintenance

If liquids are spilled onto the pump head, it should be removed and cleaned with a mild detergent. Remove any tubing from the pump head prior to removal for cleaning.

The same cleaning procedure should be used to limit the build-up of dust (which can become electrostatically charged and/or heated by friction).

Servicing

The unit is not user serviceable. Please contact your local PENTAX service facility.

Limited Warranty

Materials and workmanship of this product, PENTAX Medical IRRIGATION PUMP will be warranted for a period of one (1) year from the original date of purchase. If this product should become inoperable due to a defect in materials or workmanship during this one (1) year warranty period, the product will be repaired or replaced. This limited warranty does NOT include replacement or service due to an accident, natural disaster, use of the product other than specified, external electrical error, improper installation, negligence, modification, unauthorized service or normal wear.

Disposal

No special care is needed when disposing of the device.

X. Electromagnetic Compatibility (EMC) Information

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this section.

Portable and mobile RF communications equipment can affect medical electrical equipment. Do not expose the device to sources of electromagnetic interference such as CT equipment, diathermy equipment cellular phones, RFID tags and metal detectors.

The use of accessories, transducers and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**Table 1 – Guidance And MANUFACTURER'S Declaration
ELECTROMAGNETIC EMISSIONS
For All ME EQUIPMENT And ME SYSTEMS**

Guidance and Manufacturer's Declaration – Emissions		
The EGA-500P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-500P should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF Emissions CISPR 11	Group 1	The EGA-500P uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The EGA-500P is suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

**Table 2 – Guidance And MANUFACTURER'S Declaration
Electromagnetic IMMUNITY
For All ME EQUIPMENT And ME SYSTEMS**

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The EGA-500P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-500P should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	±6kV Contact ±8kV Air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%.
EFT EN/IEC 61000-4-4	±2kV Mains ±1kV I/Os	±2kV Mains ±1kV I/Os	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout EN/IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EGA-500P requires continued operation during power mains interruptions, it is recommended that the EGA-500P be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field EN/IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

**Table 3 – Guidance And MANUFACTURER'S Declaration
Electromagnetic IMMUNITY
For ME EQUIPMENT And ME SYSTEMS That Are Not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – Electromagnetic Immunity			
The EGA-500P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-500P should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	(V1)Vrms = 3 Vrms	Portable and mobile communications equipment should be separated from the EGA-500P by no less than the distances calculated/listed below: $D = (3.5/V1)(\text{Sqrt } P)$ $D = (3.5/E1)(\text{Sqrt } P)$ 80 to 800 MHz $D = (7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	(E1)V/m = 3 V/m	Where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

**Table 4 – Recommended Separation Distances Between Portable And
Mobile RF Communications Equipment And
The ME EQUIPMENT Or ME SYSTEM
For ME EQUIPMENT And ME SYSTEMS That Are Not LIFE-SUPPORTING**

Recommended Separations Distances for the EGA-500P			
The EGA-500P is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the EGA-500P can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the EGA-500P as recommended below, according to the maximum output power of the communications equipment.			
Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D=(3.5/\sqrt{P})(\sqrt{P})$	Separation (m) 80 to 800MHz $D=(3.5/\sqrt{E1})(\sqrt{P})$	Separation (m) 800MHz to 2.5GHz $D=(7/\sqrt{E1})(\sqrt{P})$
0.01	.1166	.1166	.2333
0.1	.3689	.3689	.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.333

Approved accessories		
3.7m	Manufacturer/Model	Maximum length
Power cable	Interpower/86610610	3.7m

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MEDICAL

INSTRUCTIONS FOR USE

EGA-501P
PENTAX Medical
CO₂ Insufflator

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PENTAX® is a registered trademark of Hoya Corporation.

IFU-65 Revision L

© 2017 Medivators Inc.

CONTENT

English (EN)	5 - 21
Croatian (HR)	22 - 38
Czech (CS)	39 - 55
Danish (DA)	56 - 72
Dutch (NL)	73 - 89
Estonian (EE)	90 - 106
Finnish (FI)	107 - 123
French (FR)	124 - 140
German (DE)	141 - 157
Hungarian (HU)	158 - 174
Italian (IT)	175 - 191
Kazakh (KK)	192 - 208
Norwegian (NO)	209 - 225
Polish (PL)	226 - 242
Portuguese (PT)	243 - 259
Romanian (RO)	260 - 276
Russian (RU)	277 - 293
Serbian (SR)	294 - 310
Slovakian (SK)	311 - 327
Slovenian (SL)	328 - 344
Spanish (ES)	345 - 361
Swedish (SV)	362 - 378
Ukrainian (UA)	379 - 395

Table of Contents

I. Unpacking and Inspection	6
II. Introduction	6
Definitions	6
Indications for Use	6
Contraindications	6
III. Warnings and Cautions	7
Safety Symbols	7
IV. Product Features	9
Safety Features	9
Theory of Operation	9
V. Technical Specifications	9
Electrical Specifications	9
Mechanical Specifications	10
CO ₂ Specifications	10
Environmental requirements	11
Accessory Items	11
Replacement parts	11
VI. Front panel controls/connections	12
VII. Rear Panel Connections	13
VIII. Setup And Operation	14
Initial Setup	14
Air hose attachment	15
Pre-Procedure setup	15
Water Pre-Warming	16
Operation	16
Shutdown	16
IX. Care and Maintenance	17
Cleaning	17
Maintenance	17
X. EMC Tables	18

Instruction Manual

I. Unpacking and Inspection

Upon receiving the PENTAX® Medical CO₂ INSUFFLATOR, ensure that the following items are included in the shipping box:

- Base tray and water bottle holder
- Hardware for assembling water bottle holder
- Operator's manual
- Power cord

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II. Introduction

Definitions

- Throughout this document the PENTAX Medical CO₂ INSUFFLATOR may be referred to as “unit” or “device”
- L/min – Liters per minute (flow)
- kPa – kilo Pascal (pressure)
- VAC – Volts Alternating Current (electrical potential)
- W – Watts
- Hz – Hertz (frequency)
- PSI – Pounds Per Square Inch (pressure)

Indications for Use

The PENTAX Medical CO₂ INSUFFLATOR is designed to use CO₂ as a distention media in the gastrointestinal tract when used in conjunction with a gastrointestinal endoscope.

The PENTAX Medical CO₂ INSUFFLATOR contains a CO₂ warming feature that allows the clinician to modulate and control delivery of CO₂ to the patient that has been warmed to a target temperature 37°C.

Contraindications

The PENTAX Medical CO₂ INSUFFLATOR should only be used by or under the direct guidance of a physician experienced in the standard practices of gastrointestinal endoscopy procedures. It should only be used for procedures where insufflation of the gastrointestinal tract is necessary to support navigation of the GI endoscope and performance of any necessary evaluation procedures. The device should not be used in any other method or for any other treatments or procedures.

The device is contraindicated for laparoscopic or hysteroscopic insufflation. It must not be used for intrauterine distension.

The device is contraindicated for CT colonography. Also, the device is not intended for use with or exposure to magnetic resonance imaging (MRI) systems. Do not use the device in an environment where it could be exposed to MRI systems.

III. Warnings and Cautions

SIGNAL WORDS

WARNING:

Indicates a potentially hazardous situation which if not avoided, could result in death or serious injury.

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MR Unsafe

WARNINGS

- A. To reduce risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.
- B. To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- C. The PENTAX® Medical CO₂ INSUFFLATOR is not suitable for use in the presence of a flammable anesthetic mixture with oxygen.
- D. The PENTAX Medical CO₂ INSUFFLATOR shall be sold only by prescription for use by physicians/clinicians who are trained and qualified.
- E. Only qualified medical personnel in an acceptable medical facility should operate the PENTAX Medical CO₂ INSUFFLATOR.
- F. The PENTAX Medical CO₂ INSUFFLATOR should be connected to a properly grounded receptacle marked "Hospital Grade" or "Hospital Only," otherwise grounding reliability cannot be achieved.
- G. Extreme precaution must be taken when handling liquids around electrical equipment. DO NOT operate the PENTAX Medical CO₂ INSUFFLATOR if liquid has been spilled on the unit.
- H. Never place or stack the PENTAX Medical CO₂ INSUFFLATOR on other electrical equipment other than the PENTAX Medical irrigation pump. Electromagnetic or other interference may occur between the PENTAX Medical CO₂ INSUFFLATOR and other electronic devices. The equipment or system should be observed to verify normal operation in the configuration in which it will be used.
- I. The PENTAX Medical CO₂ INSUFFLATOR should only be used in conjunction with other equipment whose safety against leakage currents has been established.
- J. The instructions contained in the operating manuals of any equipment to be used in conjunction with the PENTAX Medical CO₂ INSUFFLATOR must be followed to avoid any possible hazard from incompatibility.
- K. The instructions for use described in this manual MUST be followed. Otherwise, compromised safety, malfunction, injury to the operator and/or patient, or costly damage to the unit and other equipment may occur.
- L. The PENTAX Medical CO₂ INSUFFLATOR must be connected to an appropriate power source when loss of power source would result in an unacceptable risk.

CAUTION

- A. If emergency or abnormal function occurs, immediately turn off the power to the unit.
- B. Use only medical USP medical grade CO₂ from a "D" or "E" size CO₂ tank.
- C. Read and understand all warnings that come with your commercially available "D" or "E" sized CO₂ tanks.
- D. Field-servicing of the PENTAX Medical CO₂ INSUFFLATOR is limited to the replacement of the power cords, high pressure hoses, yoke adapters, water bottle holder and heater assemblies, and fuses.
- E. Remove power from the device before initiating any field servicing of the replacement parts.
- F. There are no user serviceable parts inside this unit. Repairs to the PENTAX Medical CO₂ INSUFFLATOR should only be performed by qualified service personnel.
- G. When the low pressure warning indicator illuminates, exchange the CO₂ tank as soon as possible to avoid loss of function.
- H. Always keep a spare full tank of CO₂ nearby for quick access.
- I. Always keep the CO₂ tank in an upright position to avoid fluid entering the unit.
- J. Do not use the device if the enclosure is damaged or enclosure integrity has been compromised.
- K. Do not attempt to operate the device before reading and understanding all sections of this manual.
- L. Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in section X.
- M. Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment. Do not expose the device to sources of electromagnetic interference such as CT equipment, diathermy equipment, cellular phones, RFID tags and metal detectors.

IV. Product Features

- Internal CO₂ warmer eliminates cold CO₂
- Includes control for optional water bottle warmer to maintain water at close to body temperature
- Stackable with the PENTAX® Medical IRRIGATION PUMP
- Standard male luer lock output
- Uses standard “D” or “E” size CO₂ cylinders
- Accepts input from wall CO₂ with standard DISS type connection
- Usable at CO₂ input pressures lower than most comparable units, making efficient use of tank supply
- Works with disposable PENTAX Medical DispoCap CO₂ (100150CO₂P) and PENTAX Medical CO₂ Source Tubing (100551P) for safety and compliance

Safety Features

- Dual inline pressure regulators ensure a constant pressure
- Additional mechanical pressure relief at ≤12 PSI prevents overpressure
- Low CO₂ input pressure is indicated by a yellow light on the front panel
- On initial power up, gas flow is off until the flow knob is turned on
- CO₂ warmer will not operate unless flow is initiated

Theory of Operation

The PENTAX Medical CO₂ INSUFFLATOR works by regulating CO₂ from a tank source to a preset pressure and then controlling the output of CO₂ to an GI endoscope system for ultimate delivery as a distention media in the gastrointestinal tract. The clinician will use the air/water valve on the GI endoscope and visual feedback of the GI endoscopic system to manually distend the gastrointestinal tract with CO₂.

The PENTAX Medical CO₂ INSUFFLATOR contains a CO₂ warming feature that allows the clinician to modulate and control delivery of CO₂ to the patient that has been warmed to a target temperature of 37°C with an upper tolerance of +3°C. Warming is achieved as CO₂ flows across a heating element which can be switched on or off by the user, depending on if warming of the CO₂ is desired. The heating of the CO₂ is accomplished through an internal 20W heating element with dual redundant temperature sensors.

The flow rate can be adjusted to low, medium or high settings depending on the desired rate of insufflation, with flow rates of 1.4 L/min, 2.4 L/min and 3.5 L/min maximum, respectively.

When the unit is turned on, the CO₂ flow and heating function are off by default. To initiate gas flow the user must turn the flow setting knob from the “off” position. The unit will detect gas flow and will not allow the heater to operate unless gas flow is detected. In addition, there is a low gas pressure sensor and the unit will indicate when the input CO₂ pressure reaches 25PSI by illuminating a yellow light on the front panel.

The unit also has an output to power the water bottle warmer. The water bottle heating element maintains the temperature of the water in a sterile water bottle at 37±3°C. It is controlled by dual temperature sensors for redundancy and safety.

V. Technical Specifications

Electrical Specifications

Input Voltage:	100-240 VAC
Input Frequency:	50-60 Hz
Power Consumption:	82 VA
Fuse Rating:	M10AL250V
	Medium acting, 10 amp, low breaking capacity, 250 volt

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Certifications:	IEC-60601-1, IEC-60601-1-2, IEC-60601-2-18
Classification:	Class 1 Type B
IP Rating (Ingress Protection):	IP24

**WARNING:**

Grounding reliability can only be achieved when connected to a receptacle marked "Hospital Grade".

Mechanical Specifications

Physical Dimensions:

Height	4 ¾" inches	121 mm
Width	7 ¾" inches	197 mm
Depth	13 ¾" inches	349 mm
Weight	10.5 pounds	4.8 kg

CO₂ Specifications

Free Flow rates (maximum unrestricted output to GI endoscope system):**Front panel setting**

Low: 1.4 L/min

Medium: 2.4 L/min

High: 3.5 L/min

GI Endoscopic Flow rates (approximate flow rate at GI endoscope output):**Front panel setting**

Low: 0.6-0.9 L/min*

Medium: 1.2-1.4 L/min*

High: 1.5-1.8 L/min*

** These values are approximate and are based on testing with respective models of GI endoscopes. User results may vary depending on GI endoscope used, channel diameter, and channel length.*

Output Connection: Male luer lock

Input Connection: ¼" male flare fitting.

CO₂ Heater: 20W, internal, flow through, with dual redundant temperature sensors.

CO₂ Output temperature: 37°C (+3°C)

Operating Pressure:

Input Pressure 1900 PSI Maximum
25 PSI Minimum*

** The low pressure warning light will illuminate at 25 PSI, although the unit may still flow. It is recommended to replenish the CO₂ supply when the light illuminates.*

Pressure regulator set point 8 PSI Maximum

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Accessory Items

The PENTAX® Medical CO₂ INSUFFLATOR is intended to be used only with the specific models of compatible tubing sets, connectors and accessories identified in the table below and/ or officially recommended by PENTAX. Use of the device with accessories not identified in the table below or officially recommended by PENTAX may result in incompatibility and/ or the risk of cross contamination and infection transmission.

Order Number	Description
100150CO ₂ P, 100150CO ₂ UP (24 hr Use)	PENTAX Medical DispoCap CO ₂
100160P, 100150CO ₂ UP (24 hr Use)	PENTAX Medical DispoCap Air
100130P, 100130UP (24 hr Use)	PENTAX Medical Irrigation Tubing For EGA-500P
200230P, 200230UP (24 hr Use)	PENTAX Medical Irrigation Tubing For EGP-100P
100242P (Single use only)	PENTAX Medical Single Use Waterjet Connector
100116P (24 hr Use)	PENTAX Medical 24hr Use Waterjet Connector
100551P	PENTAX Medical CO ₂ Source Tubing

Replacement parts

Order Number	Description
EGA-7014	Power cord (E.U.)
EGP-042	Power cord (U.K.)
EGA-7011	High pressure hose
EGA-7012	Yoke adapter
EGA-7010	Water bottle holder and warmer assembly

VI. Front panel controls/connections

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
- 1. On/Off push button**
Turns main power on or off to the unit
- 2. Water warmer On/Off**
Push to turn water warmer on or off
- 3. CO₂ Flow control**
Controls flow rate of CO₂ output. The three levels are represented by three green indicators of "L", "M", and "H" which are illuminated to indicate the flow rate setting. The flow rates are as follows:
Low – 1.4L/min
Medium – 2.4L/min
High – 3.5L/min

**Flow rates shown are maximum free flow measurements at output of unit*

Low gas pressure is indicated by the words "LOW PRESSURE" flashing to notify the user that the gas input pressure has fallen below 25PSI.

- 4. CO₂ gas output connector**
Connect to tubing and GI endoscope.
Standard male luer lock connector.
- 5. CO₂ gas heater On/Off**
Push to toggle CO₂ gas heater On or Off

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VII. Rear Panel Connections

1. Water bottle warmer connection

Power connection for water bottle warmer
Refer to section VIII for warmer explanation

2. CO₂ input connection

For connecting a "D" or "E" sized CO₂ tank. Use only pre-filtered medical grade CO₂ gas.

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5. Fuse receptacle


WARNING:

Replace only with type and rating marked M10AL250V.

6. Water bottle holder

The water bottle holder contains the warming system for the water bottle and connects to the main PENTAX® Medical CO₂ INSUFFLATOR by means of an electrical cable which plugs into the receptacle marked "water warmer".

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
VIII. Setup And Operation

Initial Setup

1. Place the PENTAX® Medical CO₂ INSUFFLATOR on a flat surface such as an accessory cart or other suitable working surface.
2. Before beginning, inspect the PENTAX Medical CO₂ INSUFFLATOR for any signs of damage.
3. Connect the power cord to the back of the device.
4. Before connecting the power cord to the hospital grade wall receptacle, make sure the power is off, and no accessories are connected. Then plug the power cord into the wall receptacle.
5. Assemble the water bottle holder and warmer assembly in the following manner:
 - A. The following supplies are necessary and can be found in a bag attached to the water bottle holder:
 - a.) Short Screws (2)
 - b.) Long Screws (2)
 - c.) Large Hex Key
 - d.) Small Hex Key

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- C. Place the PENTAX Medical CO₂ INSUFFLATOR on the base tray and install the (2) long screws using the large hex key.

- D. Adjust the water bottle hold-down bracket to fit the water bottles used by your facility. This is done by loosening the nut on the back and sliding the bracket up or down so it contacts the water bottle and then re-tightening the nut.
6. Plug the water warmer cable into the receptacle marked "water warmer" on the back of the unit.

Air hose attachment

Tools required: (not included)

- 9/16" open-end wrench
- 3/4" open-end wrench

Accessories required:

- "D" or "E" size CO₂ tank filled with medical grade CO₂ (not included)
- High pressure hose
- Yoke adapter

1. Assemble the high pressure hose and the yoke by tightening the hose with the 3/4" wrench while holding the yoke with the 9/16" open-end wrench.
2. Unscrew the T-handle on the yoke to allow it to slip over the top post of the CO₂ tank. The correct orientation is so that the two pins on the yoke line up with the two alignment holes on the tank. (CO₂ tank is not included with system, purchase separately)
3. After aligning the two pins, tighten the T-handle firmly.
4. Tighten the other end of the high pressure hose and yoke assembly to the back of the PENTAX® Medical CO₂ INSUFFLATOR. Hold the brass fitting with wrench and use a second wrench to tighten the hose nut. DO NOT USE TEFLON TAPE OR THREAD SEALANT ON THE FITTINGS.

Pre-Procedure setup

1. Open the valve on the CO₂ tank approximately 1 turn.
2. Turn the power switch to "On"
3. Test the setup by turning the flow control knob position clockwise. Verify that there is CO₂ flowing out the front output of the unit. If there is no CO₂ observed out the front of the unit verify that the CO₂ tank is full and check the connection outlined in air hose attachment above.
4. Turn off the flow control knob.
5. Connect the tubing to the front output of the device. Use only PENTAX Medical tubing manufactured by MEDIVATORS.
6. When using PENTAX Medical CO₂ source tubing, attach PENTAX Medical CO₂ source tubing to PENTAX Medical tubing. Attach PENTAX Medical tubing to the GI endoscope according to the manufacturer's instructions.

Water Pre-Warming

1. The PENTAX® Medical CO₂ INSUFFLATOR can be used to control the water bottle warmer system. The warmer system can maintain the temperature of bottles of sterile water at 37°C (±3°C). However, the system is not intended to raise the temperature of the water; therefore, pre-warming of the water bottle is required.
2. The water bottle must be pre-warmed with a method such as a fluid warmer if warm water is desired. Water bottle pre-warming should be performed according to the following instructions:

Place sterile water bottle in a fluid warming oven set to 37°±3°C. Do not exceed this temperature. If warm water is not desired, turn the water warmer on the pump off and place a bottle of room temperature sterile water on the water bottle holder.

NOTE:

If the water warmer has been enabled and has exceeded the upper temperature set point limit, the water warmer's switch indicator light will turn blue and the warmer will automatically shut off. If the switch indicator light is flashing blue, ensure the water warmer cable is plugged into the back of the unit. If it is connected and the light is flashing, this may indicate a problem with the unit. Contact your local service facility.

Read and become familiar with all manufacturer's instructions on warming ovens and water bottles regarding maximum fluid temperatures. Never exceed 40°C during water pre-warming. Never use a microwave oven to warm a water bottle, as this could heat the water to a dangerously high temperature or could result in uneven warming.

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Operation

1. After the pre-procedure setup has been completed as outlined in section VIII, flow can be initiated.
2. Turn the flow control knob to the first position. "L" will illuminate for "LOW". If a higher flow rate is desired, turn the knob to the next position for "MEDIUM" and one more position for "HIGH". "M" will illuminate for "MEDIUM" and "H" will illuminate for "HIGH". The knob will turn freely in each direction, and will increase or decrease the flow rate accordingly. See section V for flow rates.
3. If warming of the CO₂ is desired, press the button for the CO₂ heater. The button will illuminate green.
4. If the device is attached to a water bottle warmer, and water warming is desired, press the button for water bottle warmer. The button will illuminate green.
5. Use the appropriate air/water valve on the GI endoscope and visual feedback of the GI endoscopic system to manually distend the gastrointestinal tract with CO₂.

Shutdown

1. When the procedure is over, turn the knob counter clockwise until all lights are off, turn off the water bottle warmer (if attached) and turn off the CO₂ heater.
2. Turn off the main power switch and disconnect the tubing from the front of the unit.
3. Discard the tubing according to the manufacturer's instructions for proper disposal.

IX. Care and Maintenance

Cleaning

- Before cleaning ensure the power to the unit is turned off and the electrical cord unplugged.
- The outside surface of the unit can be cleaned with a damp cloth, 70% isopropyl alcohol solution, or a 10% bleach with water solution as often as deemed necessary.
- To disinfect the outside surface, use a mild disinfectant as needed according to the manufacturer's instructions.
- Do not use abrasive or harsh cleaners.
- Do not allow liquid to enter the unit.
- Do not sterilize the unit by any means.

Maintenance

Servicing

The device is not user serviceable (excluding the replacement parts listed in section V). Please contact your local PENTAX® service facility.

Limited Warranty

Materials and workmanship of this product, PENTAX Medical CO₂ INSUFFLATOR will be warranted for a period of one (1) year from the original date of purchase. If this product should become inoperable due to a defect in materials or workmanship during this one (1) year warranty period, the product will be repaired or replaced. This limited warranty does NOT include replacement or service due to an accident, natural disaster, use of the product other than specified, external electrical error, improper installation, negligence, modification, unauthorized service or normal wear.

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Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The EGA-501P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-501P should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF Emissions CISPR 11	Group 1	The EGA-501P uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The EGA-501P is suitable for use in all establishments, other than domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

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Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The EGA-501P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-501P should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	±6kV Contact ±8kV Air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%
EFT EN/IEC 61000-4-4	±2kV Mains ±1kV I/Os	±2kV Mains ±1kV I/Os	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout EN/IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	100% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles Note 1	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EGA-501P requires continued operation during power mains interruptions, it is recommended that the EGA-501P be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field EN/IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Note 1: The EUT shuts off and must be manually restarted by the operator during a 5 second loss of AC Mains power.

**Table 4 – Guidance And MANUFACTURER'S Declaration
Electromagnetic IMMUNITY
For ME EQUIPMENT And ME SYSTEMS That Are Not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity			
The EGA-501P is intended for use in the electromagnetic environment specified below. The customer or user of the EGA-501P should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	(3)Vrms	Portable and mobile communications equipment should be separated from the EGA-501P by no less than the distances calculated/listed below: $D=(3.5/V1)(\text{Sqrt } P)$ $D=(3.5/E1)(\text{Sqrt } P)$ 80 to 800 MHz $D=(7/E1)(\text{Sqrt } P)$ 800 MHz to 2.5 GHz
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	(3)V/m	Where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter..

Table 6 – Recommended Separation Distances Between Portable And Mobile RF Communications Equipment And The EGA-501P EQUIPMENT and SYSTEMS That Are Not LIFE-SUPPORTING

Recommended Separations Distances for the EGA-501P			
The EGA-501P is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the EGA-501P can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the EGA-501P as recommended below, according to the maximum output power of the communications equipment.			
Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D=(3.5/\sqrt{P})$	Separation (m) 80 to 800MHz $D=(3.5/\sqrt{P})$	Separation (m) 800MHz to 2.5GHz $D=(7/\sqrt{P})$
0.01	.1166	.1166	.2333
0.1	.3689	.3689	.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.333

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Fax: +86-21-61701655

PENTAX Medical Singapore Pte. Ltd.

438A Alexandra Road, #08-06 Alexandra Technopark Singapore 119967

Tel: +65-6271-1669

Fax: +65-6271-1691

**PENTAX Medical Company
A Division of PENTAX of America, inc.**

3 Paragon Drive Montvale, New Jersey 07645-1782 USA

Tel: +1-201-571-2300 Toll Free: +1-800-431-5880

Fax: +1-201-391-4189

PENTAX Europe GmbH

Julius Vosseler Strasse 104, 22527 Hamburg, Germany

Tel: +49-40-561-920

Fax: +49-40-560-4213

Manufactured in the USA by:

CE  Medivators Inc.
14605 28th Avenue North
Minneapolis, MN 55447 USA
Toll Free: +1.800.444.4729
0050 www.medivators.com

Medivators BV
Sourethweg 11
6422PC Heerlen
The Netherlands
Tel: +31.45.5.471.471

printed in USA

PENTAX

MEDICAL

PENTAX
MEDICAL

INSTRUCTIONS FOR USE

EPIK-i7010

PENTAX Medical VIDEO PROCESSOR

INDICATION FOR USE:

The PENTAX Medical EPK-i7010 Video Processor is intended to be used with the PENTAX camera heads, endoscopes, light sources, monitors and other ancillary equipments for endoscopic diagnosis, treatment and video observation.

The PENTAX Medical EPK-i7010 includes a digital post-processing imaging enhancement technology (PENTAX i-Scan™) and an optical imaging enhancement technology (OE). These imaging enhancement technologies are intended to be used as an optional adjunct following traditional white light endoscopy and is not intended to replace histopathological sampling. i-Scan and OE are compatible with PENTAX video endoscopes.

DEVICE DESCRIPTION:

The PENTAX Medical EPK-i7010 video processor consists of a video system, integrated light source, monitor, and ancillary equipment. This processor is intended for endoscopic diagnostic, treatment and video observation.

The PENTAX Medical EPK-i7010 video processor contains two types of contrast enhancement techniques: PENTAX i-Scan technology, and optical enhancement (OE) technology. PENTAX i-Scan technology is a digital filter-based contrast enhancement technique, while OE technology combines the band limited light and the digital processing.

PENTAX i-Scan technology has three modes, i-Scan 1, 2 and 3. i-Scan 1 enhances image topography and edges and i-Scan 2 and 3 enhance the color tone of the image by dissecting and recombining the individual red, green and blue (RGB) components of a white light image.

PENTAX i-Scan modes 1, 2, and 3, are intended to give the user an enhanced view of the texture of the mucosal surface and blood vessels. i-Scan 1 provides the user with a view that sharpens surface vessels and enhances surface texture of the mucosa. i-Scan 2 provides the user with increased visibility of blood vessels while also providing the same enhancements to the mucosa achieved in i-Scan 1. i-Scan 3 provides the user with increased visibility of blood vessels including dimly illuminated far-field regions while also providing the same enhancement to the mucosa achieved in i-Scan 1. The user can select either white light image or i-Scan modes by pressing a pre-programmed button on the scope, by using a pre-programmed foot pedal or by pressing a keyboard button. i-Scan is intended to be used as an optional adjunct following traditional white light endoscopy and is not intended to replace histopathological sampling.

The PENTAX Medical EPK-i7010 video processor is equipped with two optical filters placed inside Xenon lamp light path to provide optical enhancement (OE), which combines the band limited light and the digital processing. PENTAX OE technology has two modes, OE Mode1, 2.

In addition to i-Scan, OE Mode1 and 2 are intended to provide alternative methods to improve blood vessel visibility (emphasizes mucosal microvasculature and fine mucosal structures) on the mucosal surface by combination of band limited light illumination source. OE Mode1 uses a spectral filter that transmits band limited green and blue light and is intended to provide the user with enhanced image of blood vessels and fine structure of mucosa. OE Mode2 uses a different spectral filter that transmits band limited red, green and blue light and is intended to provide the user with enhanced image of blood vessels and fine structure of mucosa in an image closer to white light image.

The user can select either white light image or OE modes by pressing a pre-programmed button on the scope, by using a pre-programmed foot pedal or by pressing a keyboard button. OE is intended to be used as an optional adjunct following traditional white light endoscopy and is not intended to replace histopathological sampling.

White light is captured from a 300 Watt xenon lamp housed in the PENTAX Medical EPK-i7010 video processor. All visualization is done with the white light mode first. White light (BGR) illuminates the tissue and transfers the captured light through the video scope or a charged coupled device (CCD). Note that the white light visualization mode is always used first by the physician. For i-Scan image enhancement, the modification of the combination of RGB components for each pixel occurs when the i-Scan function is turned on in the PENTAX Medical EPK-i7010 video processor. The resulting i-Scan image is then displayed on the observation monitor. For OE image enhancement, one of the two optical filters corresponding to Mode1 and Mode2 are inserted into illumination light path when the OE function is turned on in the PENTAX Medical EPK-i7010 video processor. The resulting OE image is then displayed on the observation monitor.

FUNCTIONS USED FREQUENTLY:

The functions used frequently on this equipment are the still image display function, image capture function, white balance adjustment function, and patient information registration/display function.

APPLICATION:

- ① Intended medical use: Observation of body cavities (Lighting and air/water supply through a connected endoscope and monitor output of observed images from the endoscope itself)
- ② Patient populations: Patients ranging from pediatric to adult (Patients who are considered suitable for the application of equipment by the physician)
- ③ Intended anatomical area (tissues/organs):
Pharynx; Larynx; Esophagus; Stomach; Duodenum; Small Intestine; Large Intestine; Pancreatic and Biliary Duct; Trachea; Bronchial Tree; Ear; Nasal Passage; Kidney; Renal Cavity; Renal Calyces; Urethra; Endometrial Canal (all require endoscope connection)
- ④ User: An operator with experience, or an operator working under the supervision of someone with experience
- ⑤ Usage environment: Hospitals, Ambulatory Surgery Centers, and Medical Clinics

IMPORTANT:

This Instructions for Use (IFU) describes the recommended procedures for inspecting and preparing this video processor prior to its use and the care and maintenance after its use. It does not describe how an actual procedure is to be performed, nor does it attempt to teach the beginner the proper technique or any medical aspects regarding the use of the equipment.

Failure to follow the instructions in this IFU may result in damage to and/or malfunction of the equipment. If you have any questions regarding any of the information in this IFU or concerns pertaining to the safety and/or use of this equipment, please contact your local PENTAX Medical service facility.

CAUTION

After the product has been delivered, promptly install, set up, and start using the equipment. Avoid storing the equipment over an extended period of time without using.

TRADEMARKS:

- PENTAX is either registered trademark or trademark of HOYA Corporation.
- PENTAX Medical is the trademark of HOYA Corporation.
- Windows 7, Windows 7 Ultimate, Windows Media Player are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.
- All other product names mentioned in this IFU provided with this processor are registered trademarks or trademarks of their respective manufacturers.

Contents

1	Safety Precautions - Important	1
2	Nomenclature, Controls, and Functions	5
	Main Unit	5
	Touch Panel	6
	Rear Panel.....	7
	PENTAX Medical Water Bottle Assembly (OS-H5)	8
3	Preparation	9
	System Configuration	9
	Installation	10
	Connecting a Video Endoscope	17
4	Pre-use Inspections	19
	Power-On and Lamp Inspection	19
	On-screen Image Check.....	20
	Automatic Brightness Control Inspection.....	20
	Endoscope Control Button.....	20
	Exposure and Measurement Inspection	21
	Color Adjustment	21
	Suction Inspection	21
	Air/Water Supply Inspection	22
5	Operating Procedure	23
	Turning On/Off the Processor.....	23
	Monitor Screens	23
	Touch Panel Operation.....	25
	i-scan	30
	Optical Enhancement (OE).....	32
	Image Processing Settings.....	33
	Setup Menu Operations.....	35
	[i-scan] Tab	35
	[Image] Tab.....	39
	[Customize] Tab.....	44
	[System] Tab.....	58
	[Peripheral] Tab	71
	Keyboard Operations.....	72
	Operations in the Keyboard Menu Screen	75
	Operations in the [Patient list preset] Tab.....	75
	Operations in the [User list preset] Tab	77
	Settings in the [Setup] Tab	78
6	Maintenance	80
	After Use.....	80
	Cleaning and Storage of the Water Bottle Assembly	81
	Storage	81
	Replacing the Lamp.....	82
	Replacing the Fuses	84
	Repair	85

7	Troubleshooting	86
8	Specifications	91
9	Electromagnetic Compatibility	92

The following precautions should always be exercised with the use of all electro-medical equipment to ensure safety to all involved parties - user(s), patient(s), etc.

Please carefully read and follow this IFU.

1-1. TRAINING

1. This equipment should only be used under the supervision of a trained physician in a medical facility. Do not use in other locations or for any other purposes than the intended application. Failure to observe this precaution may result in personal injury.

1-2. INSTALLATION

1. This equipment should NEVER be installed or used in areas where the unit could get wet or be exposed to any environmental conditions such as high temperature, humidity, direct sunlight, dust, salt, etc., which could adversely affect the equipment.
2. This equipment should NEVER be installed or used in the presence of flammable or explosive gases or chemicals.
3. This equipment should NEVER be installed, used or transported in an inclined position nor should it be subjected to impact or vibration.
4. For safety reasons, this equipment must be properly grounded.
5. Ensure that all power requirements are met and conform to those specified on the rating plate located on the rear panel.
6. Do not block the air intake vent grid of this equipment.
7. Do not allow the power cord to become twisted, crushed or pulled taut.
8. When using an isolation transformer for any ancillary equipment, ensure the power requirements of the devices do not exceed the capacity of the isolation transformer. For further information, contact your local PENTAX Medical service facility.

1-3. PRIOR TO USE

1. Confirm that this equipment functions properly and check the operation of all switches, indicators, etc.
2. To prevent electrical shock when used with endoscopes, this equipment is insulated (type BF electro-medical equipment). Do not allow it to be grounded to other electrical devices being used on the patient. Rubber gloves should always be worn to prevent grounding through user(s).
3. Confirm that other devices used in conjunction with this equipment function properly and that these other devices will not adversely affect the operation or safety of this equipment. If any component of the endoscopic system is not properly functioning, the procedure should not be performed.
4. Check and confirm that all cords or cables are connected correctly and securely.
5. The lamp life when used in this equipment is 500 hours. Prior to use, check the lamp life indicator on the operation panel to ensure the indicator is lit green. After 500 hours of use, the indicator turns red and the image quality will deteriorate. The lamp life could be affected by frequency of use. In which case, the lamp life might become shorter than 500 hours.
6. Right after turning on the power of the processor, wait till the backlight of the touch panel comes on. Using the operation panel or the keyboard before the backlight comes on would result in system malfunction.

1-4. DURING USE

1. To prevent electric shock, the endoscope and/or any other ancillary device should NEVER be applied directly to the heart.
2. Make sure that no contact is made between the patient and this equipment.

3. To avoid damage to the luminous display and flat membrane switches, do not press any keys with any sharp or pointed objects.
4. The light emitted by the Xenon lamp is extremely intense. Avoid looking directly at the light exiting the endoscope and/or this equipment.
5. To protect the users eyes and avoid risk of thermal injury during an endoscopic examination, use only the minimum amount of brightness required.
6. During clinical procedures, avoid unnecessary prolonged use which could compromise patient/user safety.
7. Continually monitor this equipment and the patient for any signs of irregularities.
8. In the event that some type of irregularity is noted to the patient or this equipment, take the appropriate action to ensure patient safety.
9. If the operation of any of the components of the endoscopic system fails during the procedure and the visualization of the procedure is lost or compromised, place the endoscope in the neutral position and slowly withdraw the endoscope.
10. This equipment should only be used according to the instruction and operating conditions described in this IFU. Failure to do so could result in compromised safety, equipment malfunction or instrument damage.
11. The equipment may become hot during use. Take special care when touching the equipment.
12. If an error message is displayed on the monitor, follow the message instruction. Failure to do so could result in compromised safety for patient and/or user or equipment malfunction.
13. When the observation object moves quickly, the after image might appear on the monitor.
14. The processor's performance could be negatively affected by electromagnetic waves emitted from cellular phones or other portable wireless equipments. Recommended separation distance between the processor and the portable/mobile equipment is shown on a table near the end of this IFU. Follow these recommendation to avoid electromagnetic interference by the portable/mobile equipments.

1-5. AFTER USE

1. The equipment may become hot immediately after use. Take special care when touching the equipment.
2. Refer to the operating instructions provided with all the components of the endoscopic system to establish the right order in which components should be turned off. Some peripheral devices may have to be turned off first to avoid compromising their operation.
3. Wipe all surfaces clean with gauze slightly dampened with alcohol.
4. Be sure connector interfaces and air intake vents are not allowed to become wet or splashed with liquids.
5. Clean this equipment regularly. Remove substances on movable units as well as dirt on the surface. Substances on endoscope locking lever and/or endoscope electrical connector may cause malfunction.

1-6. STORAGE

1. This equipment should NEVER be stored in areas where the unit could get wet or be exposed to any environmental conditions such as high temperature, humidity, direct sunlight, dust, salt, etc., which could adversely affect the equipment.
2. This equipment should NEVER be stored in the presence of flammable or explosive gases or chemicals.
3. This equipment should NEVER be stored or transported in an inclined position, nor should it be subjected to impact or vibration.
4. Cords, accessories, etc., should be cleaned by wiping the surface with alcohol-dampened gauze and neatly stored.
5. This equipment should be maintained in a clean condition during storage and be ready for subsequent use.

1-7. SERVICE

1. Alterations/modifications to the equipment should NEVER be made. Repairs should only be performed by an authorized PENTAX Medical service facility.
2. When replacing the lamp, use only the lamp recommended by PENTAX Medical and follow all PENTAX Medical instructions provided.

1-8. MAINTENANCE

1. Periodically this equipment and any applicable accessories should be inspected for operation and safety.

1-9. DISPOSAL

1. The equipment should be returned for disposal to PENTAX Medical.
Contact your local PENTAX Medical representative or service facility.

An information on Disposal for users in the European Union

This product is a medical device. In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted waste, but should be collected separately. Contact your local PENTAX Medical service facility for correct disposal and recycling.

By disposing of this product correctly you will help ensure that the waste undergoes the necessary treatment, recovery and recycling and thus prevent potential negative effects on the environment and human health which could otherwise arise due to inappropriate waste handling.

POWER REQUIREMENTS

Check the standard power plug configurations that are used in your country. If the appropriate power cord is not included in your product, notify your local PENTAX Medical service facility.

SOFTWARE VERSION

The software version of the products corresponding to this IFU is described on its back cover.

SYMBOLS ON MARKING

Manufacturer

Date of manufacture

Alternating current

Type BF applied part (Safety degree specified by IEC 60601-1)

Switches the power on and off.

Warning

Warning, electricity

Caution

Warning; Hot surface

Protective earth (ground)

Equipotentiality

Instructions for Use

Follow the Instructions for Use

Endoscope

Illumination lamp model

Authorized representative in the European Community

CONVENTIONS

The following conventions have been established in the text of this IFU to aid in the identification of potential hazards of operation:

: Could result in death or serious injury.


: May result in minor or moderate injury or property-damage.

: May result in property-damage. Also, advises owner/operator about important information on the use of this equipment.

Main Unit

NOTE

- Always insert an endoscope prior to turning on the processor. Also, always turn off the processor before removing an endoscope.
- Before turning on the processor, make sure the air intake vents are not blocked.
- The lamp of the processor should be turned off except during pre-use inspections and clinical use.
- After connecting an endoscope to the processor, check that the endoscope is firmly secured and that the endoscope locking lever is in the LOCK position.

CAUTION

- Immediately after use, the metal light guide prong and the electrical contacts/pins of the endoscope may be hot. To avoid burns, do not touch these parts immediately after use. For safe handling after a procedure, grasp the PVE connector housing of the endoscope.
- The light emitted by the Xenon lamp is extremely intense. Avoid looking directly at the light exiting this equipment without Light Guide Receptacle being attached.

Touch Panel

2

COLOR BALANCE Button

Use this button to display the color balance adjustment screen.

 "Color Balance Button"(P.28)

Customize Buttons

Each press of a button changes the setting of the indicated function. You can set the functions you use very frequently for these four buttons. Depending on the function, the setting may be displayed on the top right of a button.

 "Customize Buttons"(P.29)

XLUM ON/OFF Button

Use this button to maximize the light intensity of the endoscope.

 "XLUM ON/OFF Button"(P.29)

Brightness Setting Buttons

Use  to adjust the brightness level.

 "Brightness Setting Buttons"(P.26)

Pump Level Setting Buttons

Use  to adjust the strength of the pump.

 "Pump Level Setting Buttons"(P.27)

LAMP Button

Use this button to turn the lamp on or off.

 "Lamp Button"(P.26)

PUMP Button

Use this button to operate or stop the air/water pump.

 "Pump Button"(P.26)

WHITE BALANCE Button

Use this button to adjust the white balance.

 "White Balance Button"(P.27)

EXPOSURE CONTROL Button

Use this button to change the method of adjusting the screen brightness.

AUTO: Adjust the screen brightness automatically.

MANUAL: Adjust the screen brightness manually.

 "Exposure Control Button"(P.28)

SETUP Button

Use this button to display the setting menus.

 "Setup Button"(P.29)

Lamp Life Indicator [LIFE]

The lamp usage time is indicated by this indicator.

3 green circles: Less than 400 hours

2 green circles: 400 hours to 449 hours

1 yellow circle: 450 hours to 499 hours

1 red circle: 500 hours or more (the lamp needs to be replaced)

 "Replacing the Lamp"(P.82)

MODE Button

Use this button to select the light measuring method.

AVERAGE:

Adjust the [BRIGHTNESS] level in relation to the average value of the brightness of the video signal.

PEAK: Adjust the [BRIGHTNESS] level in relation to the maximum value of the brightness of the screen.


 "Mode Button"(P.28)

NOTE

- The lamp life is 500 hours. Check the LIFE indicator (lamp life indicator) on the touch panel before using the processor. Replace the lamp cartridge if a red circle is displayed on the LIFE indicator. The lamp life may not be as long as 500 hours depending on the method of use and operating environment.
 "Replacing the Lamp"(P.82)
- When pressing a button on the touch panel, make sure to press the center of the button. If the edge of the button is pressed, the system might activate the button that is next to the button you intended to press.
- When pressing a button on the touch panel, make sure that the setting/configuration change that was intended to be made in the previous press has been correctly applied. Pressing a button in succession may cause a system error. If a system error occurs, turn the power off and then on again.

Rear Panel

2


WARNING

- To avoid the risk of electric shock, this equipment must only be connected to a power supply with protective earth.
- This equipment must be connected to an appropriate power supply.

CAUTION

Make sure to firmly connect each connector to the specified position. Failing to do so could cause certain function to be unavailable.

PENTAX Medical Water Bottle Assembly (OS-H5)


2

NOTE

When cleaning and sterilizing the water bottle assembly (OS-H5), follow the instruction given in the Instructions for Use for OS-H5.

System Configuration

The following shows examples of system configurations for using the processor.

※

NOTE

Recording devices specified in the diagram have been tested for compatibility with the EPK-i7010. Equivalent devices can be used but may not be compatible with the EPK-i7010.

NOTE

- For peripheral devices, use the recommended devices shown above or equivalent products.
- PENTAX Medical recommends the devices shown above based on the independent test results.
- The processor or the peripheral devices may not operate normally by the combination of their connections including the recommended devices.
- When any connected peripheral devices are added, use the processor and the peripheral devices after confirming their normal operations by checking their connected combinations.
- The data may not be recorded properly if non-compatible external hard drive or non-compatible USB flash memory is used, or if the USB flash memory is connected to the USB port of this product using USB extension cable.

Installation

■ Installing the Processor in a Cart



Install the processor and peripheral devices in a cart while referring to the installation examples in "System Configuration" (P.9). Place the cart on a stable and level surface.

WARNING

- The processor is an electro-medical device incorporating precision components and sophisticated circuitry so do not install it in a location where it will be subject to harsh conditions, including excessive vibrations and/or severe impact. Never drop the processor or subject it to a severe impact as doing so could compromise the functionality and/or safety of the unit. Should this equipment be mishandled or dropped, do not use it. Return it to an authorized PENTAX Medical service facility for inspection and repair.
- Because of the adverse effect that one equipment could cause to another equipment, it is possible that the correct operation of this processor could be compromised if another equipment is placed right nearby or stacked on top of the processor. If such situations cannot be avoided for extenuating reasons, be sure to check the correct operation of each equipment before using them in clinical use.
- Do NOT block the ventilation grids on the sides and at the back of the processor. In particular make sure to separate the ventilation grids on the processor's left side at least 15 cm (6 inches) from the wall. Blocking ventilation may cause equipment temperature to rise.

CAUTION

- When you raise this processor, hold the portion from the both sides shown in the figure.
- Make sure to firmly connect each connector to the specified position. Failing to do so could cause certain function to be unavailable.
- Install the processor on a stable and level surface (cart, counter, stand, etc.).
 - Avoid places where the processor may be splashed with liquid.
 - Never use the processor where explosive or flammable gases are present.
 - Do not install the processor where it will be exposed to high temperature and humidity or direct sunlight.
 - Do NOT use this equipment by installing it to stand vertically. Doing so will hinder the air to vent through the ventilation grid, and eventually it will become the cause of product damage.
- Never install, operate, or store electro-medical equipment in a dusty place. The accumulation of dust within the equipment may cause a malfunction, smoke, or ignition.
- Do not block the air intake vents on the back and the sides of the processor.
- When moving the processor, do not hold the endoscope locking lever.
- Connect a water bottle assembly (OS-H5) before using the processor.
- When loading the processor on a shelf of the cart, provide sufficient space for attaching or detaching the water bottle assembly (OS-H5) easily.

■ Connecting a Power Supply

Connect the power cord to a suitable power outlet with a protective earth terminal that meets the power rating indicated on the rating plate or a medical grade isolation transformer. When peripheral devices are used, connect the devices to the medical grade isolation transformer of the appropriate power rating in order to reduce the risk of electric shock.


⚠ WARNING

- To reduce the risk of electric shock, connect the power cord of the equipment into a power outlet with protective earth.
- When peripheral devices such as a monitor are used, connect the devices to the medical grade isolation transformer of the appropriate power rating. Do not connect the devices which are not used to the isolation transformer.

⚠ CAUTION

- Check that the total power consumption of all connected devices does not exceed the isolation transformer's power rating.
- Make sure that output is in compliance with IEC 60601-1.
- Make sure that the power cord is connected to a suitable power outlet with a protective earth terminal.

3

- 1 **Connect the AC power cord to the power input socket of the processor.**
- 2 **Connect the power cord to a suitable power outlet with a protective earth terminal that meets the power rating indicated on the rating plate or a medical grade isolation transformer.**
 - Make sure the processor is turned off beforehand. If the switch has “looseness”, the processor is turned on. If the switch has “looseness”, push the switch again and confirm there is no “looseness”.
 - When connecting peripheral devices to an isolation transformer, make sure the peripheral devices are turned off beforehand. For details, refer to the IFU for the peripheral devices.

■ Connecting Peripheral Devices

Connect the peripheral devices to the processor. Make sure the processor and peripheral devices are turned off beforehand.

! WARNING


When connecting peripheral equipments to the processor, use only the cables specified in this IFU. If different cables are used, the equipments' resistance to the electromagnetic waves will change, and the equipments could cause adverse effect to each other. If you use different cables, please contact your local PENTAX Medical service facility.

! CAUTION

- Before connecting peripheral devices, check that the devices operate correctly.
- Before connecting or removing peripheral devices, make sure the processor and peripheral devices are turned off.
- To avoid malfunction caused by changing the connecting configuration, use after confirming the normal operation of the processor and the peripheral devices by checking adequately.
- Since the processor may not operate normally if the peripheral devices are removed while in use, be sure to remove them after turning off the power.
- To avoid malfunction caused by updating or upgrading of the peripheral devices, use after confirming the normal operation of the processor and the peripheral devices by checking adequately.
- When used near a radio or TV receiver in clinical or residential areas, the processor may be subjected to radio interference.
- To reduce electromagnetic interference, do not keep turning on the main power switch of the processor while an endoscope is connected but not ready for use.
- To avoid and resolve adverse electromagnetic effects, do not operate the processor near RF (Radio Frequency) generating equipment.

● Connecting a Monitor or Recorder

Connect a monitor or recorder to the processor. For details, refer to the IFU for the monitor and recorder.

! CAUTION

- The range of display may not match with the following connection methods because the image output specifications differ.

Connecting a monitor or recorder with an HD-SDI input connector (BNC) (recommended)

Use the PENTAX Medical HD-SDI cable (OS-A85) to connect an HD-SDI connector to an HD-SDI input connector on the monitor(26"Radiance G2 HB).

Connecting a monitor or recorder with a DVI or VGA input connector

Use the PENTAX Medical DVI cable (OS-A78) to connect the DVI connector to a DVI input connector of the recorder(HVO-1000MD), or use the PENTAX VGA cable (OS-A74) to connect the DVI connector to a VGA input connector on the monitor(Radiance19).

Connecting a monitor or recorder with a composite video input connector

Use the PENTAX Medical BNC video cable (OS-A17) to connect the VIDEO OUT/SYNC OUT connector to a video input connector on the monitor or recorder.

Connecting a monitor or recorder with an RGB input connector (BNC)

Use the PENTAX Medical RGB cable (OS-A25) to connect the RGB connector to the RGB input connector on the monitor.

Connecting a monitor or recorder with a 4-pin female mini-DIN connector

Use the PENTAX Medical Y/C cable (OS-A24) to connect a Y/C connector to the 4-pin female mini-DIN connector on the monitor or recorder.


NOTE

Depending on the setting of the monitor screen, endoscopic image or text information might spill out of the screen. Adjust the display position by performing overscan/underscan adjustment or vertical/horizontal adjustment and make sure that the image as well as text is correctly displayed on the screen.

● Connecting a Printer, USB Flash Memory or External Hard Disk Drive

Connect a printer or "Image Recording Media" (P.91) to the processor.

For details on a printer, refer to the IFU for the printer.

3

Using RGB output method

Use an RGB cable (OS-A25) to connect the RGB connector to the connector on the printer.

Using Y/C output method

Use a Y/C cable (OS-A24) to connect the Y/C connector to the connector on the printer(UP-55MD).

Connecting a printer with an video input connector

Connect the VIDEO OUT/SYNC OUT terminal of the processor with video input terminal of the printer (UP-25MD) using a BNC cable (OS-A17).

Connecting a USB printer

Use a USB cable(2.0m) to connect a USB port on the front or rear of the processor to the port on the USB printer(UP-D25MD).

⚠ CAUTION

- When using a printer, be sure to adjust the settings of the printer so that the colors in images appear the same on the monitor and in printouts.
- Compatible USB printer with the processor is UP-D25MD (Sony Business Solutions Corporation).

Connecting USB flash memory

Connect the USB flash memory(TS32GJF600) to a USB port on the front or rear of the processor.

Connecting an external hard disk drive


Connect the external hard disk drive (WDBBGB0030HBK) to a USB port on the front or rear of the processor.

NOTE

- Some Image Recording Media may not be recognized by the processor.
- It is possible that the correct data recording could be compromised because of using the USB3.0 compatible cable.

● Connecting Other Peripheral Devices

For details, refer to the corresponding IFU.

Connecting a keyboard

Connect the PENTAX Medical keyboard (OS-A83) to the KEYBOARD connector.

Connecting foot switches

Connect the PENTAX foot switches (OS-A61) to the FTSW connector.

Connecting a peripheral devices supporting remote (trigger) input

When you want to use a peripheral device that includes support for remote (trigger) input, use a control cable (OS-A58) to connect the REMOTE connector to the remote connector on the peripheral device. For details on a peripheral device, refer to the IFU for the peripheral device.

Connecting an external strobe device

Connect the VIDEO OUT/SYNC OUT connector to the signal connector on an external strobe device recommended by PENTAX Medical. For details on connecting a compatible device, refer to the IFU for the compatible device.

Connecting a microphone

If you want to record video images with audio, connect a microphone to the AUDIO IN connector. Use a microphone (with the following specifications) recommended by PENTAX Medical.

Type: Stereo mini-plug (three pins)

Support for plug-in-power

Microphone characteristics: Sensitivity of -40 to -50 dB(0dB = 1V/1Pa, 1kHz)

Connecting a computer with an RS-232C connector

Use an RS-232C cable (2.0m) to connect the RS-232C connector to the connector on the computer with Endonet installed. For details on connecting an Endonet equipment, refer to the Instructions for Use for the equipment.

Connecting a computer with Endoimage2 installed

Connect one end of a LAN cable (5.0m) to the RJ45 port of a computer with Endoimage2 installed and connect the other end to the processor. For details on connecting compatible devices, refer to IFU for each compatible device.

Connecting a compatible device recommended by PENTAX Medical

Connect the RJ45 connector and a compatible device recommended by PENTAX Medical. For details on connecting a compatible device, refer to the IFU for the compatible device.

CAUTION

- When the processor is connected to an existing network environment through the RJ45 connector, it is necessary to use a router equipped with security functions between the RJ45 connector of the processor and network.
- Do not use a router to connect to other than the processor and the existing network.
- For details on the specification of the router, contact your PENTAX Medical service facility.
- This equipment is a Class B Medical Equipment (specified CISPR 11) and is intended for hospitals, ambulatory surgery centers, and medical clinics.
- Use the connection cable of length specified below:
HD-SDI cable (2.0 m), DVI-A cable (2.0 m), Composite video cable (1.5 m), Y/C video cable (1.5 m), Control Cable (1.5 m), RS-232C cable (2.0 m), Audio cable (1.5 m), and Keyboard cable (1.65 m)
- Refer to [Peripheral] Tab(P.71) for a connection setup of the peripheral devices. Also refer to [Network setup](P.78) for a setup of the network connected with RJ45 terminal.
- This processor transmits an observation image, patient information, etc. to Endoimage2, and receives a still image capture command. To endoPRO, it transmits an observation image and receives patient information, a still picture capture command, etc..

PENTAX Medical Policy on Antivirus Software

Connecting the processor to a network may result in the processor being exposed to computer viruses and security vulnerabilities. If the processor is connected to a security-compromised network or virus-infected computers and USB flash drives, this could cause the processor to be vulnerable to computer virus infection. For these reasons, you should consult with the IT coordinators in your hospital or clinic before connecting the processor to a network or USB flash drives via the LAN port or USB port. It is the responsibility of the customer to establish a secure network connection and ensure appropriate virus protection software is regularly updated and the network is configured with security. Due to the increasing major antivirus solutions and licensing issues, it is not possible for PENTAX Medical to predict and prepare appropriate antivirus solutions for each customer. Therefore, PENTAX Medical customers typically have to install virus protection software at each medical facility according to the policies of the facility.

Connecting a Video Endoscope

Connect an endoscope to the processor. Connect the water bottle assembly and suction unit to the endoscope before use. For details on an endoscope, refer to the IFU for the endoscope.

- 1 Check to ensure the endoscope locking lever is open.
- 2 Insert the endoscope firmly into connector ① of the processor. Close the endoscope locking lever in the direction of ② until it clicks into position. (Refer to Figure 3.1.)

NOTE

- After connecting the endoscope to the processor, make sure that the endoscope is firmly secured and that the endoscope locking lever is in the LOCK position.
- On rare occasions, noise may be temporarily displayed on endoscopic images immediately after the start-up video processor or endoscope connection.


Figure 3.1

- 3 Fill the bottle approximately 2/3 full (to the "200" marking on the bottle scale) with sterile water, and then screw the cap onto the bottle and securely tighten it.

WARNING

Check whether there are any defects such as cracks on the water bottle assembly exterior. If there is a defect on the water bottle assembly, do not use it.

NOTE

Do not overtighten the water bottle cap.

- 4 Align the A/W DRAIN lever with the A/W (air/water supply) position. (Refer to Figure 3.2.)
- 5 Insert the air pipe of the water bottle assembly into the water bottle socket of the processor until it clicks into position.

NOTE

Do not use excessive force when pushing the water bottle towards the processor. Rough handling may cause water to leak from the bottle onto/into the video processor.


Figure 3.2

- Insert the air/water connector into the holder on the water bottle cap until the endoscope is connected. (Refer to Figure 3.3.)

NOTE

- Remove the water bottle assembly before moving the processor to a place other than the usual place of use.
- Remove the water bottle assembly before packing the processor for transportation.


Figure 3.3

- 6 Connect the air/water connector to the air/water socket on the PVE connector of the endoscope. (Refer to Figure 3.4.)
- 7 Connect the suction tube of the suction unit to the suction nipple on the PVE connector of the endoscope.


Figure 3.4

■ When Using an Endoscopic Electrosurgery Device

WARNING

When performing endoscopic electrosurgery, use the PENTAX Medical condenser earth cable OL-Z4 as shown below to reduce noise appearing on the screen.

- 1 Connect the condenser earth cable OL-Z4 to the potential equalization terminal on the processor. (Refer to Figure 3.5.)


Figure 3.5

- 2 Connect the other end of condenser earth cable to the feedback terminal on the PVE connector of the endoscope. (Refer to Figure 3.6.)


Figure 3.6

NOTE

- Some of the 90i series endoscopes are necessary to update the scope's software. If a 90i series endoscope with outdated software is connected to this video processor, the error message is displayed on the monitor. In this case, contact your local PENTAX Medical service facility to have them inspect the endoscope and ask them to update the endoscope's software.
- Depending on the endoscope model that is being used, there is a possibility that slight oscillating noise may be visible in the top and bottom sides of the endoscopic image.

4 Pre-use Inspections

To ensure the equipment is in good working condition for use with patients, check that the endoscope, processor, and accessories are clean and operate properly before use. Also, provide spare equipment in case use of the equipment must be stopped as result of a malfunction or other unforeseen problem.

Be sure to perform the following inspection steps before using the processor. If you discover any abnormality, stop using the processor immediately and contact a PENTAX Medical service facility.

Check the following before you begin the inspection.

- The processor is turned off.
- The processor is installed on a stable and level surface.
- The water bottle assembly is prepared properly and connected correctly.
- The power cord is connected properly.
- The endoscope is connected properly.
- The keyboard is connected properly.
- The USB flash memory or other storage media is connected properly.

Before saving any images, check that

- (a) the format for the data to be saved to the USB flash memory is the desired format, and
- (b) the USB flash memory or similar is selected as the storage location.

 "Select drive (Still image)" (P.60), "Select drive (Video image)" (P.63)

WARNING

Be sure to perform the following inspections before you use the processor. In the event that any function or device in the video endoscope system does not work properly, do not perform the endoscopic examination. Contact a PENTAX Medical service facility before using the equipment.

Power-On and Lamp Inspection

1 Turn on the monitor and peripheral devices.

2 Press the power switch to turn on the processor, and check that the LED illuminates.

- The touch panel turns on.
- Check that the LIFE indicator (lamp life indicator) on the touch panel is lit green.
- If pressing  does not turn on the processor, check the fuses.

 "Replacing the Fuses" (P.84)

CAUTION

- The lamp life is 500 hours. Check the LIFE indicator (lamp life indicator) on the touch panel before using the processor. Replace the lamp cartridge if one red circle on the LIFE indicator is lit red and a message prompting you to replace the lamp is displayed on the monitor. The lamp life may not be as long as 500 hours depending on the method of use and operating environment.
 "Replacing the Lamp" (P.82)
- If the touch panel does not lit up after more than one and a half minutes from the system startup, restart the system.

3 Press on the touch panel.

-  illuminates green.

4 Check that the main lamp of the processor illuminates, and the distal end of the endoscope emits light.

WARNING

- Do not look directly at the light emitted from the distal end of the endoscope. Doing so can cause eye injury.
- Do not open the lamp housing cover to avoid electric shock.

NOTE

- If the lamp fails to illuminate, stop the inspection and contact a PENTAX Medical service facility.
- If a PENTAX Medical video endoscope or video module that is not compatible with the processor is connected, endoscopic images will not be displayed.
- The auxiliary lamp is for the purpose of ensuring that you can see what you are doing while removing the endoscope in the event of the main lamp fails during a procedure. If this occurs, release the angulation of the endoscope, remove the endoscope from the patient's body cavity, and stop the examination. Do not perform a procedure using the auxiliary lamp.

On-screen Image Check

- 1 Check the image displayed on the monitor.

WARNING

Check the live image is properly output on the monitor in the right direction.

Automatic Brightness Control Inspection

Check the video images from the endoscope on the monitor while the lamp is lit.

- 1 Exercise the endoscope's automatic iris. Bring the tip of the endoscope within 1 cm of the palm of your hand and move it to about 5 cm away from your palm. Watch the image displayed on the monitor to ensure the brightness at both distances is similar. Lift the distal end of the endoscope close to the room lights, and check that the light being emitted from the distal end of the endoscope lowers significantly (dependent on the ambient light levels in the room). Return the distal end of the endoscope to a point on your palm and ensure that the light is being emitted from the distal end of the endoscope. (Refer to Figure 4.1.)


Figure 4.1

4

CAUTION

- Wear rubber gloves when performing the inspection.
- If items are brought up close to the distal end of the endoscope and left there for a long time, it could cause heat to be generated by the light being emitted from the distal end of the endoscope.
- If dark-colored gloves are brought close to the distal end of the endoscope, it may cause heat to be generated, which may result in a burn injury.

NOTE

- It may take more than 10 seconds for the endoscopic image to appear on the monitor.
- If the endoscopic image does not appear on the monitor, turn the power off and then on again.

Endoscope Control Button

- 1 Press each button on the control body of the endoscope and check that a confirmation tone sounds. (Refer to Figure 4.2.)
 - For details on each of the buttons on the control body of the endoscope, refer to the IFU for the endoscope.


Figure 4.2

Exposure and Measurement Inspection

1 Press **EXPOSURE CONTROL** on the touch panel, and select [MANUAL].

- The selected item illuminates green. Each press switches the item between [AUTO] and [MANUAL].

2 Press **BRIGHTNESS** of [BRIGHTNESS] and check that the level of the indicator changes.

3 Press **EXPOSURE CONTROL**, and select [AUTO].

- The selected item illuminates green. Each press switches the item between [AUTO] and [MANUAL].
- Selection is not possible when an endoscope is not connected.

4 Press **MODE**, and select [AVERAGE] or [PEAK].

- The selected item illuminates green. Each press switches the item between [AVERAGE] and [PEAK].
- Selection is not possible when an endoscope is not connected.

Color Adjustment

1 Adjust the white balance before adjusting the color balance.

☞ "White Balance Button" (P.27)

2 Press **COLOR BALANCE** on the touch panel.

- The color balance adjustment screen appears.

3 Press **BLUE** or **RED** and check the level value. Gently place your hand around the distal end of the endoscope, and check the changes in the image of your hand on the monitor.

Suction Inspection

1 Place the distal end of the endoscope in a basin of sterile water, and press the suction control valve. (Refer to Figure 4.3.)

- The sterile water should be rapidly aspirated into the suction system.

NOTE

Be sure to use freshly distilled or sterilized water.


Figure 4.3

Air/Water Supply Inspection

- 1 Set the A/W DRAIN lever to A/W position. (Refer to Figure 4.4)
Insert the air pipe of the water bottle assembly into the water bottle socket of the video processor and press until the water bottle assembly 'clicks' into position.


Figure 4.4

- 2 Press  (ON/OFF) on the touch panel so that  illuminates green and the operating sound of the air pump can be heard.
- 3 Press  of [AIR FLOW] to select the strength level.

- 4 Air/water is supplied through the endoscope. If you insert the distal end of the endoscope into water and gently cover the air hole on the air/water feeding valve, air will come out of the distal end of the endoscope. You can confirm the air supply by checking whether air bubbles are generated. (Refer to Figure 4.5.)

NOTE

Be sure to use freshly distilled or sterilized water.


Figure 4.5

- 5 Fully press the air/water feeding valve and check that water comes out of the distal end of the endoscope. (Refer to Figure 4.6.)


Figure 4.6

NOTE

For details about air and water delivery, refer to the Instructions for Use for the endoscope.

CAUTION

If all of the inspection items are working properly, then the endoscopic examination can be performed. If there is a problem with any of the functions, do not perform the endoscopic examination. Contact a PENTAX Medical service facility.

5 Operating Procedure

Turning On/Off the Processor

1 Press ① on the front of the processor. (Refer to Figure 5.1.)

- The processor turns on, the touch panel illuminates, and the screen for normal mode appears on the monitor.
-  "Normal Screen" (P.23)
- The patient information and user information displayed immediately prior to turning off the processor last time appear on the monitor.

2 Press ① again to turn off the processor.


Figure 5.1

CAUTION

- The processor can be shut down the power supply by turning off a power switch.
- After the processor has been turned off, wait at least 30 seconds before turning it back on again. A fuse may be blown if you turn the power on and off repeatedly within a short period of time.  "Replacing the Fuses" (P.84)
- Do not use a sharp object such as the tip of a pen to press buttons on the touch panel.
- On rare occasions, the processor may not start up when it is turned on. If this happens, turn the processor off and then on again.
- Do not use the end of the light guide of the endoscope to operate buttons or switches.
- Data may not be saved correctly to an external device if power is turned off to the processor while the data is being transferred.

NOTE

On rare occasions, noise may be temporarily displayed on endoscopic images immediately after the start-up video processor or endoscope connection.


5

Monitor Screens

Normal Screen


The normal screen shown below appears on the monitor when you turn on the processor.

HD Monitor Display

- ① Date
- ② Time (Hours/Minutes/Seconds)
- ③ Patient's Name (Up to 24 Characters)
- ④ ID (Up to 12 Characters)
- ⑤ Age (Up to 3 Characters)
- ⑥ Sex (1 Character)

SD Monitor Display

- ⑦ Comment (Up to 40 Characters)
- ⑧ Doctor's Name (Up to 24 Characters)
- ⑨ Facility (Up to 16 Characters)
- ⑩ Endoscope Name
- ⑪ Endoscope Serial No.
- ⑫ Film counter

NOTE

- You can change the display for the HD monitor and the SD monitor in [Display setting].  "Display setting" (P.66)
- Depending on the setting of the monitor screen, endoscopic image or text information might spill out of the screen. Adjust the display position by performing overscan/underscan adjustment or vertical/horizontal adjustment and make sure that the image as well as text is correctly displayed on the screen.
- Displayed/undisplayed of each item can be set up by [Character display].  "Settings in the [User data edit No.*] Screen" (P.78)

Freeze Screen P-in-P Mode Screen

When the freeze function is executed, the monitor display is split into a main screen and sub screen, and still images are displayed in the main screen and endoscope video images are displayed in the sub screen.

When images are input from an external device, you can continue to display the endoscope video images in the main screen and display the images input from the external device in the sub screen. (P-in-P mode)

HD Monitor Display

SD Monitor Display

CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time before attempting to use the equipment or after changing the image mode from freeze screen to normal screen or twin mode screen.

NOTE

If you set the [Sub screen position] to [Upper Right] on an SD monitor, [Date] and [Time] will be hidden.

Twin Mode Screen

This screen displays two endoscope video images together. It enables you to confirm an image for which image processing has been performed and an unprocessed image at the same time.

HD Monitor Display

SD Monitor Display

Touch Panel Operation

⚠ CAUTION

- Do not use a sharp object such as the tip of a pen to press the buttons on the touch panel.
- Do not use the end of the light guide of the endoscope to operate the touch panel.

Lamp Button

1 Press on the touch panel.

-  illuminates green. The processor's main lamp illuminates, and the distal end of the endoscope connected to the processor emits light.
- If the lamp does not illuminate, the auxiliary lamp illuminates and  flashes green.

CAUTION

- To protect your eyes, avoid looking directly at the light emitted from the endoscope and/or processor.
- The lamp life used by the processor is 500 hours. Prior to use, check the lamp life indicator on the touch panel and confirm that it is lit green or yellow. After 500 hours of use, the indicator will turn red and image quality will deteriorate. Excessive use of the lamp beyond its rated 500 hours (approaching a thousand hours of use or more) could cause the lamp to explode resulting in damage to the video processor. The lamp life may not be as long as 500 hours depending on the method of use and operating environment.
 "Replacing the Lamp" (P.82)

NOTE

- If the main lamp does not illuminate even if you press , check the lamp life indicator on the touch panel. Replace the lamp if the indicator is red ( "Replacing the Lamp" (P.82)). If the indicator is green or yellow, the lamp system or lamp cartridge is broken.
 "Troubleshooting" (P.86)
- The lamp of the processor should be turned off except during pre-use inspection and clinical use.
- The auxiliary lamp is for the purpose of ensuring that you can see what you are doing while removing the endoscope in the event of the main lamp failing during a procedure. Release the angulation of the endoscope, remove the endoscope from the patient's body cavity, and stop the examination. Do not perform a procedure using the auxiliary lamp.
- If the auxiliary lamp is lit, all image processing is disabled.

5

Brightness Setting Buttons

You can manually adjust the brightness of images appearing on the monitor.

1 Press or of [BRIGHTNESS].

- Each press of  increases the brightness by one step and each press of  decreases the brightness by one step.
- When you adjust the brightness, the level is indicated by the indicator and the value is displayed on the monitor. Note that the standard value is [0] when you make adjustments.

NOTE

The brightness setting is stored in memory even after you turn off the processor.

Pump Button

1 Press on the touch panel.

-  illuminates green, and the air/water pump operates. Press the button again to stop operation.
- Press  or  of [AIR FLOW] to set the strength of the pump air flow.
 "Pump Level Setting Buttons" (P.27)

NOTE

When pump button becomes inoperable, stop the procedure. Release the finger tip from the air/water feeding valve of the endoscope to stop feeding of the air and water into the patient's body. Then, carefully withdraw the endoscope from the patient.

Pump Level Setting Buttons

Set the strength of the air and water supply of the pump.

1 Press or of [AIR FLOW].

- When you set the pump level, the level is indicated by the indicator and the value is displayed on the monitor.

CAUTION

- Regardless of pump level setting selected, avoid delivering too much air to minimize the potential for perforation or other medical complications.
- Monitor the patient continuously, and ensure that an air embolism does not occur due to excessive air supply.

NOTE

- If the objective lens is difficult to clean, the level of the pump setting on the processor can be temporarily raised but simultaneously press the air/water feeding valve and suction control valve on the endoscope while the level is raised to minimize the potential of perforation. After removing the dirt from the lens, return the pump setting to a normal level.
- When the pump level setting is [1] to [3], the setting is stored in memory even after you turn off the processor. When the setting is [4] or [5], the setting is automatically changed to [3] when you turn off the processor.

White Balance Button

Adjust the white balance. If the message "Check white balance." appears when the processor is turned on or an endoscope is connected to the processor, adjust the white balance. You can also adjust it whenever you want.

1 Connect the endoscope to the processor, and turn on the processor.

2 Set the exposure control mode to [AUTO].

3 Press on the touch panel to emit light from the distal end of the endoscope.

4 Insert the distal end of the endoscope into the supplied white balance adjuster. Move the adjuster so that the circle at the bottom can be viewed completely on the monitor screen, and verify that the inside of the adjuster can be seen clearly. (Refer to Figure 5.2.)

CAUTION

- Be careful not to hit the distal end of the endoscope against the bottom of the white balance adjuster.
- Be sure to use the PENTAX Medical white balance adjuster OS-A43H when adjusting the white balance. Do not attempt to adjust the white balance under any other condition.


Figure 5.2


5 Hold the distal end of the endoscope in place and press on the touch panel for at least 2 seconds.

- The white balance is adjusted. After the white balance is set, the message “White balance OK!” appears on the monitor for about 5 seconds.
- If adjustment of the white balance fails, the message “White balance failed! Please retry.” appears on the monitor. Adjust the white balance again.

! CAUTION

Be sure to press  and adjust white balance before the procedure.

Normal Screen

Circle at the
Bottom of the
Adjuster

**EXPOSURE
CONTROL**

Exposure Control Button

You can select from two exposure control modes for adjusting the brightness of images appearing on the monitor.

[AUTO]: The brightness set with the brightness setting buttons is automatically adjusted so that it is uniform. Selection is possible when an endoscope is connected.

[MANUAL]: The brightness set with the brightness setting buttons is applied without any automatic adjustments.

1 Press on the touch panel and select a setting.

- Each press switches the item between [AUTO] and [MANUAL].
- The selected item illuminates green.
- When changing the brightness, press  or  of [BRIGHTNESS] to set the level.

MODE

Mode Button

1 When the exposure control method is set [AUTO], select the light measuring method from [AVERAGE] and [PEAK]. The indicator of the selected method illuminates.

[AVERAGE]: Adjust the [BRIGHTNESS] level in relation to the average value of the brightness of the video signal.

[PEAK]: Adjust the [BRIGHTNESS] level in relation to the maximum value of the brightness of the screen.

NOTE

The light measuring method setting is stored in memory even after you turn off the processor.

**COLOR
BALANCE**

Color Balance Button

You can manually adjust the color tones of images appearing on the monitor.

1 Press on the touch panel.

- The color balance adjustment screen appears.

2 Press or of [BLUE] or [RED].

- Each press of  or  changes the color tone.
- When you adjust a color tone, the level is indicated by the indicator and the value is displayed on the monitor. Note that the standard value is [0] when you make adjustments.

NOTE

The color balance setting is stored in memory even after you turn off the processor.

XLUM ON/OFF Button

Use this button when you want to know the position of the distal end of the endoscope from outside the body. When this is set to on, the light intensity of the endoscope is set to maximum and you can confirm the position of the distal end of the endoscope from outside the body.

1 Press on the touch panel.

- When this is turned on, [ON] of the button illuminates green and the [BRIGHTNESS] indicator indicates the Maximum value.
- Press the button again to turn it off.

CAUTION

- Do not use the XLUM ON function during routine use.
- When XLUM is used, Brightness cannot be adjusted.
- XLUM is the function that enables the physician to locate the distal end of the endoscope within the patient's body by intensifying the light emitted by the endoscope to its maximum level. However, the effectiveness of XLUM function can be limited by such factors as patient's body size and so on.

Customize Buttons

Pressing a customize button activates the function assigned to it. Depending on the registered function, you may be able to change, for example, the setting value by just pressing the button repeatedly or you may be able to display the setting screen for that function by just pressing the button once.

You can change the functions set for the customize buttons in [Customize button settings] on the [Customize] tab of the Setup menu.

 "Customize button settings" (P.44)

1 Press the desired customize button.

- The setting value changes or the screen for setting the function appears.

Setup Button

Pressing the Setup button displays the menu for configuring the settings of each of [i-scan], [Image], [Customize], [System], and [Peripheral].

1 Press on the touch panel.

- The Setup menu appears.
 "Setup Menu Operations" (P.35)

i-scan

This video processor includes three, optional preprogrammed modes of i-scan.

i-scan is an optional digital, post-processing imaging enhancement technology, which is intended to be used as an adjunct to traditional white light endoscopy.

i-scan 1, 2, and 3 are preprogrammed default settings. Only one default setting can be activated at a time. The user is able to either choose the default settings or change the volumetric options of the settings for customizability.

The i-scan modes 1, 2 and 3 can be accessed by using a pre-programmed button on the endoscope, by using a pre-programmed foot switch or by pressing a keyboard button.

Surface Enhancement (SE), Contrast Enhancement (CE) and Tone Enhancement (TE) are software based image enhancement functions that work in combination to produce the enhancements of i-scan 1, 2 and 3.

■ Surface Enhancement (SE)

The software based SE mode allows for enhancement of detailed patterns and micro-unevenness (texture) of the surface structures of the image.

■ Contrast Enhancement (CE)

The software based CE mode allows for enhancement of micro-unevenness (texture) by identifying locally dark areas in the images and further enhancing them with a slight blue tone.

NOTE

- If [Enhancement], [SE], or [CE] is used, noise may become visible in images on the monitor screen.
- [Enhancement], [CE], [SE] settings are stored in memory even if the processor is turned off.

5

■ Tone Enhancement (TE)

The software based TE mode allows for the adjustment of color contrast in the image to enhance contrast between mucosa and blood vessels. Each TE setting displays a unique combination of red, green and blue contrast.

! WARNING

- If any newly identified irregularities in the endoscopic image are observed during any of the image enhancement observational modes, immediately return to the normal white light observation mode.
- Do not rely on the image enhancement observational modes alone for primary detection of lesions or for decisions regarding any potential diagnostic or therapeutic intervention.

! CAUTION

The tone and the brightness of the endoscopic image of the i-scan observation are different from the white light observation. Please use the i-scan observation after a thorough white light examination and reviewing the property of each mode. i-scan observation is an adjunct to white light endoscopy, and is not a substitute for either a thorough traditional examination or for histopathological results.

i-scan Setting

No.	Setting Item	Setting Value	Explanation
1	Enhancement (Enh)	Off, +1, +2, +3, +4, +5, +6, +1/a1, +2/a2, +3/a3, +4/a4, +5/a5, +6/a6	Enhances the sharpness of edges and emphasizes the mucosal structures (the enhancement level is increased with the setting value).
2	Surface Enhancement (SE)	Off, +1, +2, +3, +4, +5, +6	Enhances detailed patterns and micro-unevenness (texture) of the surface structures of the image.
3	Contrast Enhancement (CE)	Off, +1, +2, +3, +4, +5, +6	Enhances micro-unevenness (texture) by identifying locally dark areas in the images and further enhancing them with a slight blue tone.
4	Brightness (B)	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5	Allows for adjustment of brightness of image (a larger value increases brightness).
5	Color Balance	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5	Allows for adjustment of red and blue tone (a larger value increases red or blue color intensity).
6	Tone Enhancement (TE)	Off, p, v, b, e, g, c	Allows for adjustment of color contrast in image to enhance contrast between mucosa and blood vessels. Each TE setting displays a unique combination of red, green and blue contrast.
7	Light Measuring Method (LMM)	Ave, Peak	Allows for adjustment of brightness level with respect to an averaging of the brightness of the video signal or with respect to the peak brightness of the screen.
8	Noise Reduction (NR)	Off, Low, Medium, High	Allows for correction of noise appearing by image processing.

Optical Enhancement (OE)

This video processor is equipped with two optical filters placed inside Xenon lamp light path to provide optical enhancement (OE), which combines the band limited light and the digital processing.

PENTAX OE technology has two modes, OE Mode1, 2.

OE Mode1 uses a spectral filter that transmits band limited green and blue light and is intended to provide the user with enhanced image of blood vessels and fine structure of mucosa.

OE Mode2 uses a different spectral filter that transmits band limited red, green and blue light and is intended to provide the user with enhanced image of blood vessels and fine structure of mucosa in an image closer to white light image.

NOTE

OE mode cannot be processed simultaneously with [SE], [CE], [TE].

WARNING

- If any newly identified irregularities in the endoscopic image are observed during any of the image enhancement observational modes, immediately return to the normal white light observation mode.
- Do not rely on the image enhancement observational modes alone for primary detection of lesions or for decisions regarding any potential diagnostic or therapeutic intervention.

Optical Enhancement (OE) setting

Setting Item	Setting Value	Explanation
OE (Optical Enhancement)	Off, Mode1, Mode2	Inserts one of two optical filters corresponding to Mode1 and Mode2 into the illumination light path to provide a band-limited spectrum. The band-limited spectrum enhances blood vessels and fine structure of mucosa.
Enhancement (Enh)	Off, +1, +2, +3, +4, +5, +6, +1/a1, +2/a2, +3/a3, +4/a4, +5/a5, +6/a6	Enhances the sharpness of edges and emphasizes the mucosal structures (the enhancement level is increased with the setting value).
Brightness (B)	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5	Allows for adjustment of brightness of image (a larger value increases brightness).
Ave / Peak (Light Measuring Method)	Average, Peak	Allows for adjustment of brightness level with respect to an averaging of the brightness of the video signal or with respect to the peak brightness of the screen.
Red / Blue (Color Balance)	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5	Allows for adjustment of red and blue tone (a larger value increases red or blue color intensity).
Noise Reduction (NR)	Off, Low, Medium, High	Allows for correction of noise appearing by image processing.

Image Processing Settings

CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time after changes have been made to image processing settings.

■ Using SE

1 Press the [SE] key on the keyboard / customize button on the touch panel or configure setting on [i-scan] Tab.

- Each press changes the level of SE.
 "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)

CAUTION

SE is an observation mode that employs software image enhancement of regional concavity and convexity after the image has been collected by the endoscope.
All mucosal areas are to be viewed using traditional white light. SE should not be used as a substitute for a thorough traditional white light examination of the mucosa.

■ Using CE

1 Press the [CE] key on the keyboard / customize button on the touch panel or configure setting on [i-scan] Tab.

- Each press changes the level of CE.
 "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)

NOTE

- If [Enhance level], [SE], or [CE] is set to high level, noise may become visible in images on the monitor screen.
- [Enhance level], [CE], [SE] settings are stored in memory even if the processor is turned off.
- K series with CE: Equivalent effectiveness to i series with CE mode.

■ Using TE

1 Press the [TE] key on the keyboard / customize button on the touch panel or configure setting on [i-scan] Tab.

- Each press changes the mode of TE.
- Select from [p], [v], [b], [e], [g], [c], and [off].
 "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)

NOTE

K series with TE [e], [g], [c]: Equivalent effectiveness to i series with TE [e], [g], [c] mode.

■ Using OE

OE modes can be activated during i-scan or white light observations.

Switching from White light image to OE image:

- 1 If i-scan functions are activated: Press [i-scan1], [i-scan2], [i-scan3] key on the keyboard, or customize button on the touch panel, or configure setting on [i-scan] Tab.**
 - Select OFF of i-scan1, 2 or 3.
☞ "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)
- 2 Press the [OE] key on the keyboard / customize button on the touch panel or configure setting on [i-scan] Tab.**
 - Each press changes the mode of OE.
 - Select from [Mode1], [Mode2] and [off].
☞ "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)

Switching from i-scan1, i-scan2, or i-scan3 image to OE image:

- 1 Press the [OE] key on the keyboard / customize button on the touch panel or configure setting on [i-scan] Tab with i-scan1, 2 or 3 function are been activating.**
 - Each press changes the mode of OE.
 - Select from [Mode1], [Mode2] and [off].
☞ "Function Keys" (P.73), "Customize button settings" (P.44), "[i-scan] Tab" (P.35)

NOTE

- 70K and 80K series endoscopes are not compatible with OE, therefore OE function cannot be used.
- 90K, K10 and J10 series endoscopes are not compatible with Twin mode viewing in OE mode image observation.
- Some of 90K series endoscope require software update to be used with OE. Please contact your local PENTAX Medical service facility.


■ Using Enhancement

Configure the edge enhancement setting.

- 1 Press the [Enhance level] key on the keyboard / customize button on the touch panel.**
 - Each press changes the intensity of enhancement.
☞ "Function Keys" (P.73), "Customize button settings" (P.44)

Setup Menu Operations

Pressing  on the touch panel displays the Setup menu.



The Setup menu contains the following categories.

- i-scan** Configure settings for image processing related to i-scan mode.
 "[i-scan] Tab" (P.35)
- Image** Configure settings for image processing performed on endoscopic images, and settings related to twin mode.
 "[Image] Tab" (P.39)
- Customize** Configure settings for functions assigned to customize buttons, endoscope buttons, and foot switches, settings related to monitor display, and other settings.
 "[Customize] Tab" (P.44)
- System** Configure settings related to saving endoscopic images, and settings related to the processor unit. Be sure to check the setting options before using the processor for the first time.
 "[System] Tab" (P.58)
- Peripheral** Configure settings related to the peripheral device.
 "[Peripheral] Tab" (P.71)

[i-scan] Tab

Pressing the [i-scan] tab displays the following screen. The screen has two pages.

SE

Set the level for SE (surface enhancement).

1 Press [Off] or one of the six levels.

- A higher numerical value sets a higher level.

CE

Set the level for CE (contrast enhancement).

1 Press [Off] or one of the six levels.

- A higher numerical value sets a higher level.

NOTE

- If [SE] or [CE] is used, noise may become visible in images on the monitor screen.
- The [SE] and [CE] settings are stored in memory even after you turn off the processor.

TE

Set the mode for TE (tone enhancement).

1 Select [Off], [p], [v], [b], [e], [g], or [c].

CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time after changes have been made to [SE], [CE], and [TE].

OE

Set the mode for OE (optical enhancement).

1 Select [Off], [Mode1], or [Mode2].

- It cannot be processed simultaneously with [SE] [CE] [TE].
- If an endoscope is not connected to the processor, the message "Connect endoscope." appears.

NOTE

- [OE] setting remains even when an endoscope is taken out and inserted by [Scope eject] with a processor's main body turned on.
- When changing the mode of OE, an image may be disturbed until its completion.

i-scan profile

Set the functions of [i-scan1], [i-scan2], and [i-scan3].

1 Press [>].

- The [i-scan profile] screen appears.
- The setting value in the main screen and other screens is displayed for each item.

2 Select [i-scan1], [i-scan2], or [i-scan3].

- If one of the i-scan modes is already selected, the currently enabled button will be lit.

3 Press the parameter box for which you want to change the setting.

- When a parameter box is selected, it turns blue and the buttons for changing the selected parameter appear at the bottom of the screen.

4 Change the setting with the buttons at the bottom of the screen.


- [SE] [CE] [TE] and [OE] cannot be processed simultaneously.

5 Press [Save].

- The changes are saved.
- If [Load] was pressed when the setting values were changed, the saved setting values are loaded.

Select the i-scan mode to set.

Load the setting values of a saved i-scan mode.

Save the changes.

Load the recommended setting values.

Change the parameter settings.

Press the parameter box for which you want to change the setting.

5

List for assignable i-scan setting

Setting Item	Setting Value
Enhancement (Enh)	Off, +1, +2, +3, +4, +5, +6, +1/a1, +2/a2, +3/a3, +4/a4, +5/a5, +6/a6
SE (Surface Enhancement)	Off, +1, +2, +3, +4, +5, +6
CE (Contrast Enhancement)	Off, +1, +2, +3, +4, +5, +6
TE (Tone Enhancement)	Off, p, v, b, e, g, c
OE (Optical Enhancement)	Off, Mode1, Mode2
Brightness (B)	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5
Ave / Peak (Light Measuring Method)	Average, Peak
Red / Blue (Color Balance)	-5, -4, -3, -2, -1, 0, +1, +2, +3, +4, +5
Noise Reduction (NR)	Off, Low, Medium, High

NOTE

- With [Brightness] and [Ave/Peak], [No change] can be selected in addition to the normal levels. [No change] enables the level for when the [i-scan] setting is [Off] to be used.
- If you attempt to perform a keyboard operation other than pressing the [Lamp] or [Pump] key, endoscope button operation, or foot switch operation in the [i-scan profile] screen, the message "Close profile menu." appears on the monitor and the operation cannot be performed.

i-scan1-3

Select the i-scan mode.

1 Select [Off], [i-scan1], [i-scan2], or [i-scan3].

- The selected button illuminates, and the settings registered in [i-scan profile] are loaded.
- If one of the i-scan modes is already selected, the currently enabled button will be lit.
- If an endoscope is not connected to the processor, the message “Connect endoscope.” appears.




NOTE

- Configure each setting of [i-scan1], [i-scan2], and [i-scan3] in [i-scan profile].
- The [i-scan1-3] settings are retained even if the endoscope is removed and inserted using [Scope eject] while the processor is on.

[Image] Tab

Pressing the [Image] tab displays the following screen. The screen has three pages.

Display the setup menu for each of the categories.

5

Enhancement


Set the level for edge enhancement/structure enhancement display.

1 Press [>].

- The [Enhancement] screen appears.

When [Scope enhancement] on the [Customize] tab is set to [Enable] (Only 70K series)

2 Select [Off], [Low], [Medium], or [High].

When [Scope enhancement] on the [Customize] tab is set to [Disable]

3 Press [Off] or one of two types × six levels.

- A higher numerical value sets a higher level for edge enhancement / structure enhancement display.

[+1 to +6] Edge enhancement mode

[+1/a1 to +6/a6] Edge enhancement / Structure enhancement mode : allows for enhancing the mucous membrane structure and contour.

When connecting the 70K, 80K, 90K, K10 and J10 series endoscope, the edge enhancement mode [+1 to +6] is effective, and when the i series endoscope is connected, the structure enhancement mode [a1 to a6] becomes effective.

Structure enhancement mode and [SE] cannot be processed simultaneously.

! CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time after changes have been made to [Enhancement].

Noise reduction

Set the level for noise reduction.

1 Select [Off], [Low], [Medium], or [High].

! CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time after changes have been made to [Noise reduction].

NOTE

If [SE], [CE], or [Enhancement] is used, noise may become visible in the image on the monitor screen.

Contrast

Select the mode for contrast processing.

1 Select [Normal], [Mode1], or [Mode2].

! CAUTION

Be sure to check that the endoscopic images on the monitor are not static and that the images are displayed with the correct orientation each time after changes have been made to [Contrast].

NOTE

The [Enhancement], [Noise reduction], and [Contrast] settings are stored in memory even after you turn off the processor.

Digital zoom

Set the magnification ratio of endoscopic images.

1 Select [Off], [x1.2], [x1.5], or [x2.0].


NOTE

As the optical or electronic expanded observation display has a narrow view, please secure a view enough when performing an operation.

Image rotate

Set the display orientation of the endoscopic image.

1 Press [Normal] or [Rotate].

Twin mode

Set or cancel twin mode.

"Twin Mode Screen" (P.24)

1 Press [Off] or [On].

- If an endoscope is not connected to the processor, the message "Connect endoscope." appears.


NOTE

- If [SE] or other image processing is performed in twin mode, the processing is only reflected on the image on the right side of the screen. The image processing parameters are inherited by the normal screen when twin mode is cancelled.
- When the Twin mode setup is made after freezing, an image processing is not reflected in the left side screen. In the case where the USB capture is performed as is, an image-processing parametric information is added to the still picture even if not checking [Add image processing parameter information.]of [Option (Still image)]. And the image-processing parametric information of the left side screen displays "-". Image processing of the left side screen will become effective after freezing release.
- 90K, K10 and J10 series endoscopes do not correspond to the Twin mode of white light and OE.
- (When executed an image-processing)In case that the display is too dark to obtain a proper observation result, we recommend you to return to a normal observation display or change one display into white light.

Twin mode display

Set the sizes of the image on the left side of the screen and on the right side of the screen in twin mode.

1 Select [1:1], [3:2], or [2:3].

NOTE

The [Twin mode display] setting is stored in memory even after you turn off the processor.

Twin mode profile

Set the image processing for twin mode.

1 Press [➤].

- The [Twin mode profile] screen appears.
- The setting value in the main screen and the other screens is displayed for each item.

2 Select [Off] or [On] for [Profile setting].

- [Off]: The [Twin mode profile] screen is not displayed in twin mode. Image processing parameters other than [Enhancement] are not reflected on the left side of the screen.
- [On]: The [Twin mode profile] screen is always displayed on the left side of screen in twin mode. "TM profile" is displayed in the upper left of the left side of the screen.

3 Press the parameter box for which you want to change the setting.

- When a parameter box is selected, it turns blue and buttons for changing the selected parameter appear at the bottom of the screen.

4 Change the setting with the buttons at the bottom of the screen.


- If the parameter box does not turn blue, the buttons will not function properly.

5 Press [Save].

- The changes are saved.
- If [Load] was pressed when the setting values were changed, the saved setting values are loaded.

Set whether to display the [Twin mode profile] screen.

Load the saved twin mode setting values.

Save the changes.

Load the recommended setting values.

Change the setting of the parameter.

Press the parameter box for which you want to change the setting.

NOTE

- If you attempt to perform a keyboard operation other than pressing the [Lamp] or [Pump] key, endoscope button operation, or foot switch operation in the [Twin mode profile] screen, the message "Close profile menu." appears on the monitor and the operation cannot be performed.
- When the Twin mode setup is made after freezing, an image processing is not reflected in the left side screen. In the case where the USB capture is performed as is, an image-processing parametric information is added to the still picture even if not checking [Add image processing parameter information.] of [Option (Still image)]. And the image-processing parametric information of the left side screen displays "-". Image processing of the left side screen will become effective after freezing release.

■ ND mode

Set the level for ND mode.

1 Select [Off], [Low], [Medium], or [High].

⚠ CAUTION

- If an option other than [Off] is set, a normal screen will not be displayed on the monitor when a PDT or other strong light is not utilized.
- When [Low], [Medium], or [High] is set, the current user cannot be changed.  "Operations in the [User list preset] Tab" (P.77)
- When the setting screen is displayed by pressing the [Menu] key, [ND mode] operation is not possible.
- If [ND mode] is set to an option other than [Off], [i-scan1], [i-scan2], and [i-scan3] cannot be used.
- When the [ND mode] is set to anything other than [Off], the settings for [Shutter mode], [OE] will be canceled.
- Depending on the endoscope model being used, the time it takes to resume from the ND mode varies. It takes maximum of 10 seconds.
- There are certain models among 90K series endoscopes that cannot utilize ND mode. In order to use the ND mode, it is necessary to update the integrated software of the endoscope. Contact your local PENTAX Medical service facility to arrange the software update and inspection of the endoscope.

■ D-range expansion

Select the mode for D-range expansion processing.

1 Select [Off], [Low], [Medium], or [High].

NOTE

The [D-range expansion] settings are stored in memory even after you turn off the processor.

[Customize] Tab

Pressing the [Customize] tab displays the following screen. The screen has four pages.

Display the setup menu for each of the categories.




Return to the previous screen.

Return to the main screen.

Change the page.

Display the advanced settings screen.

Press any of the buttons to configure the setting.

5

Customize button settings

Set the functions for the customize buttons (four buttons at the top of the main screen).

- 1 Press [>].
 - The [Customize button settings] screen appears.
- 2 Press the customize button for which you want to set the function.
 - The selected customize button illuminates blue.


The function assigned to each customize button is displayed.

3 Press the button of the function to set.

- You can press ◀▶ to display additional functions.
- The name of the assigned function is displayed below the customize button.

Assignable Functions

Display additional functions.

The following functions can be assigned.

Function Name	Function Description	Reference Page
SE	Switch the SE level to a level available in [SE setting]. It cannot be processed simultaneously with [OE] or the structure enhancement mode [a1 to a6].	P.55
CE	Switch the CE level to a level available in [CE setting]. It cannot be processed simultaneously with [OE].	P.56
TE	Switch the TE mode to a mode available in [TE setting]. It cannot be processed simultaneously with [OE].	P.56
OE*	Switch the OE mode to a mode available in [OE setting]. It cannot be processed simultaneously with [SE] [CE] [TE].	P.56
i-scan1*	Use the settings of [i-scan1]. When the function is configured, a check mark is displayed in the check box of the corresponding button.	P.37
i-scan2*	Use the settings of [i-scan2]. When the function is configured, a check mark is displayed in the check box of the corresponding button.	P.37
i-scan3*	Use the settings of [i-scan3]. When the function is configured, a check mark is displayed in the check box of the corresponding button.	P.37
i-scan Off-1-2-3*	Switch the i-scan mode between [Off], [i-scan1], [i-scan2], and [i-scan3].	P.38
i-scan profile	Open the i-scan profile menu.	P.37
Remote1*	Output a signal from the REMOTE1 connector on the back of the processor.	—
Remote2*	Output a signal from the REMOTE2 connector on the back of the processor.	—
Remote3*	Output a signal from the REMOTE3 connector on the back of the processor.	—
Freeze*	Switch between video mode and freeze mode.	—

Function Name	Function Description	Reference Page
USB capture*	Save a still image to Image Recording Media in accordance with the setting of [File storage setting (Still image)].	P.60
USB printer*	Send a still image to the connected USB printer.	P.71
Video Rec/Stop*	Start and stop recording video. The video files are saved to Image Recording Media in accordance with the setting of [File storage setting (Video image)].	P.62
Enhancement*	Switch the edge enhancement/structure enhancement level available in [Enhancement setting].	P.39
Noise reduction	Switch the level for noise reduction between [Off], [Low], [Medium], and [High].	P.40
Contrast	Switch the mode for contrast processing between [Normal], [Mode1], and [Mode2].	P.40
Twin mode*	Set or cancel twin mode.	P.41
Twin mode display	Switch the size of the images on the left and right sides of the screen in twin mode between [1:1], [3:2], and [2:3].	P.41
Digital zoom*	Switch the magnification ratio for the endoscopic image in accordance with the setting of [Digital zoom setting].	P.54
Image rotate	Switch the display orientation for the endoscopic image between [Normal] and [Rotate].	P.41
ND mode*	Switch the level for ND mode between [Off], [Low], [Medium], and [High].	P.43
Film counter reset	Reset the film counter to 1.	P.68
Optical zoom speed	Switch the zoom speed for an optical zoom endoscope between [Low], [Medium], and [High].	P.54
P-in-P display	Switch P-in-P mode display between [Off] and [On].	P.54
External video input	Display images input from an external device on the touch panel.	P.54
P-in-P Main/Side change	Switch between the endoscopic image and external input video screens in P-in-P mode.	P.55
Scope eject*	Press and hold the button to remove the endoscope while the processor is on.	P.66
Image size	Switch the display size of the endoscopic image between [Full] and [Medium].	P.66

Function Name	Function Description	Reference Page
Stopwatch	Start or stop the stopwatch. Press and hold the button to reset the stopwatch.	P.66
System information	Display the [System information] screen.	P.68
Shutter mode*	Switch the shutter mode for the endoscope between [Off] and [On] (K series endoscope) or between [Off], [Low], [Medium], and [High] (i series endoscope).	P.52
Copy*	Capture a still image and output it to the Endonet system or endoPRO.	—
Video*	Start/stop recording video to the Endonet system or endoPRO.	—
PC capture*	Output the image to PC via Endoimage2.	—
D-range expansion	Switch the level for D-range expansion between [Off], [Low], [Medium], and [High].	P.43
Manual brightness level	Set the brightness adjustment level by endoscope type while [MANUAL] is set for [EXPOSURE CONTROL].	P.67
iDoc button	Send the current settings for i-Doc management to endoPRO.	—
Cecum marker	Send a command for the Cecum marker to endoPRO.	—

*: The button does not function when an endoscope is not connected to the processor.

Assign scope buttons

Set the functions for the endoscope buttons.

1 Press [>].

- The [Assign scope buttons] screen appears.

The function assigned to each endoscope button is displayed.

2 Select the endoscope buttons ([1], [2], [3], or [4]) for which you want to assign a function.

- The selected endoscope button illuminates blue.

3 Press the button of the function you wish to configure.


- You can press ◀▶ to display additional functions.
- The name of the assigned function is displayed below the endoscope button.

Assignable Functions

Display additional functions.

- Pressing [Peripheral] displays the [Peripheral] screen. On this screen, each external output can be switched on or off. Press the button of the function to output so that the button illuminates blue. Simultaneous output of multiple functions is possible.
- Pressing the button again switches off the output.

The following functions can be assigned.

Function Name	Function Description	Reference Page
No operation	A function is not assigned to the button.	—
Freeze	Switch between video mode and freeze mode.	—

Function Name		Function Description	Reference Page
Peripheral	Remote1	Output a signal from the REMOTE1 connector on the back of the processor.	—
	Remote2	Output a signal from the REMOTE2 connector on the back of the processor.	—
	Remote3	Output a signal from the REMOTE3 connector on the back of the processor.	—
	Copy	Capture a still image and output it to the Endonet system, or endoPRO. This cannot be configured simultaneously with [USB memory], [USB printer], [Endoimage2], and [Video].	—
	USB memory	Save a still image to Image Recording Media in accordance with the setting of [File storage setting (Still image)]. This cannot be configured simultaneously with [Copy], [USB printer], [Endoimage2], and [Video].	P.60
	USB printer	Send a still image to the connected USB printer. This cannot be configured simultaneously with [Copy], [USB memory], [Endoimage2], and [Video].	P.71
	Endoimage2	Output the image to a PC via Endoimage2. This cannot be configured simultaneously with [Copy], [USB memory], [USB printer], and [Video].	—
	Video	Start/stop recording video to the Endonet system or endoPRO. This cannot be configured simultaneously with [Copy], [USB memory], [USB printer], and [Endoimage2].	—
Video Rec/Stop		Start and stop recording video. The video files are saved to Image Recording Media in accordance with the setting of [File storage setting (Video image)].	P.62
Auto/Manual		Switch between [AUTO] and [MANUAL] for exposure control.	P.28
Average/Peak		Switch between [AVERAGE] and [PEAK] for the light measuring mode.	P.28
Enhancement		Switch the edge enhancement/structure enhancement level to a level available in [Enhancement setting].	P.39
Noise reduction		Switch the level for noise reduction between [Off], [Low], [Medium], and [High].	P.40
SE		Switch the SE level to a level available in [SE setting]. It cannot be processed simultaneously with [OE] or the structure enhancement mode [a1 to a6].	P.55
CE		Switch the CE level to a level available in [CE setting]. It cannot be processed simultaneously with [OE].	P.56


Function Name	Function Description	Reference Page
TE	Switch the TE mode to a mode available in [TE setting]. It cannot be processed simultaneously with [OE].	P.56
OE	Switch the OE mode to a mode available in [OE setting]. It cannot be processed simultaneously with [SE] [CE] [TE].	P.56
i-scan1	Use the settings of [i-scan1].	P.37
i-scan2	Use the settings of [i-scan2].	P.37
i-scan3	Use the settings of [i-scan3].	P.37
i-scan Off-1-2-3	Switch the i-scan mode between [Off], [i-scan1], [i-scan2], and [i-scan3].	P.38
Digital zoom	Switch the magnification ratio for the endoscopic image in accordance with the setting of [Digital zoom setting].	P.54
Stopwatch	Start or stop the stopwatch.	P.66
Stopwatch reset	Reset the stopwatch.	P.66
ND mode	Switch the level for ND mode between [Off], [Low], [Medium], and [High].	P.43
P-in-P Main/Side change	Switch between the endoscopic image and external input video screens in P-in-P mode.	P.55
Optical zoom (T)	Enlarge the display of the endoscopic image for an optical zoom endoscope.	—
Optical zoom (W)	Reduce the display of the endoscopic image for an optical zoom endoscope.	—
Twin mode	Set or cancel twin mode.	P.41
Twin mode display	Switch the size of the images on the left and right sides of the screen in twin mode between [1:1], [3:2], and [2:3].	P.41
Shutter mode	Switch the shutter mode for the endoscope between [Off] and [On] (K series endoscope) or between [Off], [Low], [Medium], and [High] (i series endoscope).	P.52
D-range expansion	Switch the level for D-range expansion between [Off], [Low], [Medium], and [High].	P.43
iDoc button	Send the current settings for i-Doc management to endoPRO.	—
Cecum marker	Send a command for the Cecum marker to endoPRO.	—

Assign foot switch

Set the functions for the foot switches.


- 1 Press [➤].
 - The [Assign foot switch] screen appears.
- 2 Select the foot switch for which you want to set a function from [L] and [R].
 - The selected foot switch button illuminates blue.

The function assigned to each foot switch is displayed.

- 3 Press the button of the function you wish to configure.
 - You can press ◀▶ to display additional functions.
 - The name of the assigned function is displayed below the foot switch.
 - The functions that can be assigned are the same as "Assign scope buttons" (P.48).

Assignable Functions

Display additional functions.

5

Scope enhancement

Set whether or not to enable enhancement processing in the endoscope.

[Enable] Enable endoscope enhancement processing in the endoscope. Enhancement processing in the processor is disabled.

[Disable] Disable endoscope enhancement processing in the endoscope. Enhancement processing in the processor is enabled.

- 1 Press [Enable] or [Disable].

CAUTION

When a 90i/i10/90K/K10/J10 series endoscope is connected, only enhancement processing in the processor is used, regardless of the setting of [Scope enhancement].

NOTE

For use in combination with a different model color video processor, we recommend setting [Enable]. When an endoscope that was used with this set to [Disable] is connected to a different model color video processor, the [Enhancement] setting will be set to [Off].

■ Character display

Set character display for the monitor.

[Off] All characters are not displayed.

[Date/Time] Only the date and time are displayed.

[On] The information set to [ON] in [Character display] of the keyboard menu is displayed.

1 Select [Off], [Date/Time], or [On].

■ Shutter setting

Configure the shutter speed settings of the endoscope.

1 Press [>].

- The [Shutter setting] screen appears.

2 Configure the settings of [K series scope] or [i series scope].

● K series scope

Set whether or not to automatically change the shutter speed of a K series endoscope when capturing normal images and close-up images.

“K series” consists of 70K, 80K, 90K, K10 and J10 series.

[Off] The shutter speed is not changed. Select [Off] when using an external light source.

[On] The shutter speed is changed. When capturing close-up images, this enables clear images to be obtained if sufficient amount of light is emitted from the lamp.

1 Press [Off] or [On].

● i series scope

Set whether or not to automatically change the shutter speed of an i series endoscope when capturing normal images and close-up images.

[Off] The shutter speed is not changed. Select [Off] when using an external light source.

[Low] The shutter speed is changed in low speed mode. Select [Low] if limited amount of light is emitted from the lamp when capturing close-up images.

[Medium] The shutter speed is changed in medium speed mode. When capturing close-up images, this enables clear images to be obtained if sufficient amount of light is emitted from the lamp.

[High] The shutter speed is changed in high speed mode. When capturing close-up images, this enables the clearest images to be obtained if sufficient amount of light is emitted from the lamp.

1 Select [Off], [Low], [Medium], or [High].

NOTE

When turned off the lamp, the shutter speed is not changed.

Press any of the buttons to configure the setting.

■ Sub screen setting

Configure the settings for the sub screen displayed when, for example, the freeze function is executed.

1 Press [>].

- The [Sub screen setting] screen appears.

2 Configure the [Sub screen delay], [Sub screen position], and [Sub screen position (ERCP)] settings.

● Sub screen delay

Set the display of the sub screen to be delayed.

[Off] The sub screen is displayed immediately after the freeze function is executed.

[On] The sub screen is displayed shortly after the freeze function is executed.

1 Press [Off] or [On].

● Sub screen position

Set the position to display the sub screen.

[Upper Left] Display the sub screen on the upper left of the monitor.

[Upper Right] Display the sub screen on the upper right of the monitor.

1 Press [Upper left] or [Upper right].

● Sub screen position (ERCP)

When using a side viewing endoscope (ED series), set the position to display the sub screen.


[Upper Left] Display the sub screen on the upper left of the monitor.

[Upper Right] Display the sub screen on the upper right of the monitor.

1 Press [Upper left] or [Upper right].

NOTE

The [Sub screen position] and [Sub screen position (ERCP)] settings are enabled only on SD monitor displays.

Press any of the buttons to configure the setting.

Digital zoom setting


Configure the zoom settings for when [Digital zoom] is assigned to a customize button on the touch panel, endoscope button, or foot switch, and when zoom operations are performed with the keyboard.

- 1 Press [➤].
 - The [Digital zoom setting] screen appears.
- 2 Select [Off - x1.2], [Off - x1.5], [Off - x2.0], [Off - x1.2 - x1.5], [Off - x1.5 - x2.0], [Off - x1.2 - x2.0], or [Off - x1.2 - x1.5 - x2.0] to set the zoom ratio.
- 3 Select [↗↗] or [↘↘] to set the up/down order of the zoom.

- [↗↗] When the zoom ratio has reached the maximum, pressing the customize button will return the zoom ratio to actual size.
- [↘↘] When the zoom ratio has reached the maximum, pressing the customize button will incrementally reduce the zoom ratio until it returns to the actual size.

An animation is played at the set zoom ratio and up/down order.

Indicates the set zoom ratio and up/down order.

Set the up/down order of the zoom.

Set the zoom ratio.

Optical zoom speed

You can adjust the zoom speed of an optical zoom endoscope.

- 1 Select [Low], [Medium], or [High].

NOTE

As the optical or electronic expanded observation display has a narrow view, please secure a view enough when performing an operation.

5

Film counter type

Set the display format of the film counter.

- 1 Select [1-99], [1/2], [1/4], [1/8], or [1/16].
 - If [1/2], [1/4], [1/8], or [1/16] is selected, the counter is reset to 1 when the number in the numerator becomes equal to the number in the denominator.

P-in-P Display

Configure the display settings for P-in-P mode.

[Off] The images input from an external device are not displayed on the monitor.

[On] The images input from an external device are displayed on the monitor.

- 1 Press [Off] or [On].

External video input

Display images input from an external device on the touch panel. Only images that are input from an external device to the VIDEO IN connector can be displayed.

- 1 Press [Switch].

- The images input from an external device are displayed on the touch panel. Touching the touch panel displays the previous screen.

■ P-in-P Main/Side change

Switch between the endoscopic image and external input video screens in P-in-P mode.

- [Main] The images input from an external device are displayed on the main screen of the monitor. The endoscopic images from the endoscope are displayed in the sub screen.
- [Side] The endoscopic images from the endoscope are displayed on the main screen of the monitor. The images input from an external device are displayed in the sub screen.

1 Press [Main] or [Side].

! CAUTION

Confirm that the message "Main:External video" is displayed on the monitor while [Side] is selected.

NOTE

During the twin mode, the setting for [P-in-P Main/Side change] is fixed to [Main].

■ Enhancement setting

Specify the edge enhancement/structure enhancement display levels that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [>].

- The [Enhancement setting] screen appears.
- The level that is currently specified is indicated by green characters.

2 Press the buttons of the levels you want to use.

- The selected level buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.

The level that is currently specified is indicated by green characters.

Select the levels you want to use.

Confirm the selection.

■ SE setting

Specify the SE display levels that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [>].

- The [SE setting] screen appears.
- The level that is currently specified is indicated by green characters.

2 Press the buttons of the levels you want to use.

- The selected level buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.

The level that is currently specified is indicated by green characters.

Select the levels you want to use.

Confirm the selection.

CE setting

Specify the CE display levels that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [>].

- The [CE setting] screen appears.
- The level that is currently specified is indicated by green characters at the top of the screen.

2 Press the buttons of the levels you want to use.

- The selected level buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.

The level that is currently specified is indicated by green characters.

Select the levels you want to use.

Confirm the selection.

TE setting

Specify the TE display modes that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [>].

- The [TE setting] screen appears.
- The mode that is currently specified is indicated by green characters at the top of the screen.

2 Press the buttons of the modes you want to use.

- The selected mode buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.

The mode that is currently specified is indicated by green characters.

Select the modes you want to use.

Confirm the selection.

OE setting

Specify the OE modes that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [>]

- The [OE setting] screen appears.
- The mode that is currently specified is indicated by green characters at the top of the screen.

2 Press the buttons of the modes you want to use.

- The selected mode buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.

The mode that is currently specified is indicated by green characters.

Select the modes you want to use.

Confirm the selection.

■ i-scan setting

Specify the i-scan modes that can be switched using a customize button on the touch panel, endoscope button, foot switch, or keyboard.

1 Press [➤].

- The [i-scan setting] screen appears.
- The mode that is currently specified is indicated by green characters at the top of the screen.

2 Press the buttons of the modes you want to use.

- The selected mode buttons illuminate blue.
- Pressing a button again cancels the selection.

3 Press [OK].

- If you switch to another screen without pressing [OK], the settings will not be applied.

The mode that is currently specified is indicated by green characters.

Select the modes you want to use.

Confirm the selection.

■ Freeze scan

Set whether or not to use the freeze scan function.

[Off]: The image is displayed as is when the freeze function is executed.

[Short][Middle][Long]:

The image with the least blur is selected automatically from the images up to a few seconds before the freeze function was executed.

The amount of time for image selection increases in the following order:[Short],[Middle], and [Long].

1 Press [Off],[Short],[Middle] or [Long].

NOTE




When you use the freeze scan function, wait for the endoscopic image on the LCD monitor to show the effect of desired image processing before capturing the image.

[System] Tab

Pressing the [System] tab displays the following screen. The screen has four pages.

Display the setup menu for each of the categories.

Display the advanced settings screen.


5

Date & Time

Set the date and time to be displayed on the monitor.

- Press [>]**
 - The [Date & Time] screen appears.
- Set the time zone setting.**
 - When pressing [>] in the [Time zone setting], the [Time zone setting] will be displayed. "Time zone setting" (P.59)
- Press the box of the item for which you want to change the setting from the year, month, day, hour, minute, and second items.**
 - The selected box turns blue, and ▲▼ appear above and below the box.
- Press ▲▼ to change the setting.**
- Repeat steps 3 and 4 to set all of the items you want to change.**
 - If you press [Date style], you can change the display format of the date. Each press of the button switches the order between "year/month/day," "month/day/year," and "day/month/year".
- Press [OK].**
 - If you switch to another screen without pressing the [OK] button, the settings will not be applied.

Press the box of the item for which you want to change the setting.

Confirm the settings.

● Time zone setting

Set the time zone setting.


1 Press [➤] of [Time zone setting] in the [Date & Time] screen.

- The [Time zone setting] screen appears.

2 ▲▼ to change the setting.

- Adding a check mark in the [Daylight saving time] check box adds 1 hour to the current time. Please add a check mark into the [Daylight saving time] check box manually during Daylight saving time period.

Press a check box to add or clear a check mark.

Change the setting.

Confirm the settings.

The time zone can be selected from the followings.

Time zone No.[01]	(GMT -12:00)	IDL (west side)
Time zone No.[02]	(GMT -11:00)	UTC-11
Time zone No.[03]	(GMT -10:00)	Hawaii
Time zone No.[04]	(GMT -09:00)	Alaska
Time zone No.[05]	(GMT -08:00)	PST (in the United States and Canada)
Time zone No.[06]	(GMT -07:00)	Arizona, Chihuahua, MST (in the United States and Canada)
Time zone No.[07]	(GMT -06:00)	Guadalajara, Mexico City, Monterrey, Saskatchewan, Central America, CST (in the United States and Canada)
Time zone No.[08]	(GMT -05:00)	Eastern Indiana, Bogota, Lima, Quito, EST(in the United States and Canada)
Time zone No.[09]	(GMT -04:30)	Caracas
Time zone No.[10]	(GMT -04:00)	Asuncion, Cuiaba, Santiago, AST (in Canada)
Time zone No.[11]	(GMT -03:30)	Newfoundland
Time zone No.[12]	(GMT -03:00)	Greenland, Salvador, Buenos Aires, Brasilia, Montevideo
Time zone No.[13]	(GMT -02:00)	UTC-2, Central Atlantic Ocean
Time zone No.[14]	(GMT -01:00)	Azores, Cape Verde
Time zone No.[15]	(GMT)	UTC, GMT, Dublin, Edinburgh, Lisbon, London, Monrovia, Reykjavik
Time zone No.[16]	(GMT +01:00)	Amsterdam, Berlin, Bern, Rome, Stockholm, Sarajevo, Windhoek, Warsaw, Zagreb, Brussels, Paris, Copenhagen, Madrid, Belgrade, Budapest, Ljubljana, West and Central Africa
Time zone No.[17]	(GMT +02:00)	Athens, Bucharest, Amman, Istanbul, Jerusalem, Cairo, Helsinki, Kiev, Sofia, Riga, Tallin
Time zone No.[18]	(GMT +03:00)	Kuwait, Riyadh, Nairobi, Baghdad, Kaliningrad, Minsk
Time zone No.[19]	(GMT +03:30)	Tehran
Time zone No.[20]	(GMT +04:00)	Abu Dhabi, Muscat, Yerevan, Caucasus Standard Time, Moscow, Baku, Saint Petersburg
Time zone No.[21]	(GMT +04:30)	Kabul
Time zone No.[22]	(GMT +05:00)	Islamabad, Karachi, Tashkent
Time zone No.[23]	(GMT +05:30)	Sri Jayewardenepura Kotte
Time zone No.[24]	(GMT +05:45)	Kathmandu
Time zone No.[25]	(GMT +06:00)	Astana, Ekaterinburg, Dhaka
Time zone No.[26]	(GMT +06:30)	Yangon
Time zone No.[27]	(GMT +07:00)	Novosibirsk, Bangkok, Hanoi, Jakarta
Time zone No.[28]	(GMT +08:00)	Ulan Bator, Kuala Lumpur, Singapore, Krasnoyarsk, Perth, Taipei, Beijing, Chongqing, Hong Kong, Urumqi
Time zone No.[29]	(GMT +09:00)	Irkutsk, Seoul, Osaka, Sapporo, Tokyo
Time zone No.[30]	(GMT +09:30)	Darwin
Time zone No.[31]	(GMT +10:00)	Canberra, Melbourne, Sydney, Guam, Port Moresby

Time zone No.[32]	(GMT +11:00)	Vladivostok, Solomon Islands, New Caledonia
Time zone No.[33]	(GMT +12:00)	Auckland, Wellington, Fiji, UTC+12
Time zone No.[34]	(GMT +13:00)	Samoa, Nuku'alofa

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.
- Make sure to set the date correctly.

■ File storage settings (Still image)

Set the method to save still images.

1 Press [>].

- The [File storage settings (Still image)] screen appears.

2 Configure the settings of [Select drive], [Folder and file name], [Disk space warning], and [Option].

Display the advanced settings screen.

● Select drive (Still image)

Set the "Image Recording Media" (P.91) as the destination for saving still images.


1 Press [>] of [Select drive] in the [File storage settings (Still image)] screen.

- The [Select drive (Still image)] screen appears. The Image Recording Media currently connected and the total size and free space of each Image Recording Media is indicated on the screen.
- If you connect Image Recording Media after displaying the [Select drive (Still image)] screen, press [Refresh] to perform recognition of the Image Recording Media.

2 Press the button of the Image Recording Media to use as the save destination.

- The selected Image Recording Media button illuminates blue.

Icon indicating that the Image Recording Media is selected as the save destination for video images.

Indicates the Image Recording Media currently connected and the total size and free space of each Image Recording Media.

Perform Image Recording Media recognition again.

NOTE

- Only one Image Recording Media can be set to save images at one time.
- Up to three Image Recording Media can be connected at one time. However, do not connect Image Recording Media via a USB hub.
- After replacing Image Recording Media, make sure you recheck that the [Select drive] setting is configured properly.
- High capacity Image Recording Media tends to take longer to be recognized by the system. Wait till the system recognizes the external module before start operating.

● Folder and file name (Still image)

Set the names of the folders and files generated in the PENTAX folder of the "Image Recording Media" (P.91) when still images are saved.

1 Press [➤] of [Folder and file name] in the [File storage settings (Still image)] screen.


- The [Folder and file name (Still image)] screen appears.

2 Select [Name], [ID], or [Date] of [Current folder] to set the name.

3 Press [Name] or [ID] of [File name] to switch it on or off.

- Each press of a button switches the name on or off.
- In addition to [Name] or [ID], a consecutive number is automatically added to a file name. A four-digit number is added when either [Name] or [ID] is on, and an eight-digit number when both of them are off.

The current folder and file name settings are displayed in a folder tree.

Set the folder and file names.

● Disk space warning (Still image)

Set the threshold for displaying the "Image Recording Media" (P.91) free space warning message.

1 Press [➤] of [Disk space warning] in the [File storage settings (Still image)] screen.

- The [Disk space warning (Still image)] screen appears.


2 Set the threshold value for displaying a warning message by pressing the number buttons.

- You can set a value within the range of 10 to 65,535.
- Pressing [Clear] resets the value to 0.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.
- If the free space of the Image Recording Media reaches the set value, the message "Check memory space (Still image)." appears on the monitor.

Set the value for which to display a warning message.

Reset the setting value.

Confirm the setting.

● Option (Still image)

Set options related to saving still images.


1 Press [➤] of [Option] in the [File storage settings (Still image)] screen.

- The [Option (Still image)] screen appears.

2 Set each of the items.

- [File format]
Select the setting for the format and compression rate for saving files from [BMP], [JPEG High quality], and [JPEG Standard quality].
- [Add external.] check box
Add a check mark to save endoscopic images and external input images as one file while external input images are displayed in P-in-P mode.
- [Add patient information.] check box
Add a check mark to save the patient's name, ID, and other patient information together with still images.
- [Add image processing parameter information.] check box
Add a check mark to save the SE, CE, and other image processing parameter information together with still images.

Press any of the buttons to configure the setting.

Press a check box to add or clear a check mark.

NOTE

- Since lossy compression is used for JPEG images, it may result in the loss of details in images.
- Please make sure to check the JPEG image quality prior to use if it is stored in JPEG format.
- When the setting allows for simultaneous saving of endoscopic image and the external input image in a single file, the file containing both image (endoscopic image and external input image) will be saved in peripheral device also. (The peripheral device must be specified under "Add external." menu)
- While still image data remains in the processor during the process of capturing to a USB printer or saving to Image Recording Media, a white "○" is displayed on the bottom right of the monitor. Do not disconnect the USB printer or Image Recording Media or turn off the processor until the "○" disappears.
- Please try not to execute the next capture until the "○" disappears.
- When trying to capture an image while system is processing the image at the same time, wait for the endoscopic image on the LCD monitor to be switched before proceeding to capture the next image.
- When the Twin mode setup is made after freezing, an image processing is not reflected in the left side screen. In the case where the USB capture is performed as is, an image-processing parametric information is added to the still picture even if not checking [Add image processing parameter information.] of [Option (Still image)]. And the image-processing parametric information of the left side screen displays "-". Image processing of the left side screen will become effective after freezing release.

■ File storage settings (Video image)

Set the method to save video images.

1 Press [➤].

- The [File storage setting (Video image)] screen appears.

2 Configure the settings of [Select drive], [Folder and file name], [Disk space warning], and [Option].

Display the advanced settings screen.


● Select drive (Video image)

Set the "Image Recording Media" (P.91) as the destination to save video images.

1 Press [>] of [Select drive] in the [File storage settings (Video image)] screen.

- The [Select drive (Video image)] screen appears. The Image Recording Media currently connected and the total size and free space of each Image Recording Media is indicated on the screen.
- If you connect Image Recording Media after displaying the [Select drive (Video image)] screen, press [Refresh] to perform recognition of the Image Recording Media.

Icon indicating that the Image Recording Media is selected as the save destination for still images.

Indicates the Image Recording Media currently connected and the total size and free space of each Image Recording Media.

Perform Image Recording Media recognition again.

2 Press the button of the Image Recording Media to use as the save destination.

- The selected Image Recording Media button illuminates blue.

NOTE

- Only one Image Recording Media can be set to save images at one time.
- Up to three Image Recording Media can be connected at one time. However, do not connect Image Recording Media via a USB hub.
- After replacing Image Recording Media, make sure you recheck that the [Select drive] setting is configured properly.
- High capacity Image Recording Media tends to take longer to be recognized by the system. Wait till the system recognizes the external module before start operating.

5

● Folder and file name (Video image)

Set the names of the folders and files generated in the PENTAX folder of the "Image Recording Media" (P.91) when video images are saved.

1 Press [>] of [Folder and file name] in the [File storage settings (Video image)] screen.


- The [Folder and file name (Video image)] screen appears.

2 Select [Name], [ID], or [Date] of [Current folder] to set the name.

3 Press [Name] or [ID] of [File name] to switch it on or off.

- Each press of a button switches the name on or off.
- In addition to [Name] or [ID], a consecutive number is automatically added to a file name. A four-digit number is added when either [Name] or [ID] is on, and an eight-digit number when both of them are off.

The current folder and file name settings are displayed in a folder tree.

Set the folder and file names.

● Disk space warning (Video image)

Set the threshold for displaying the "Image Recording Media" (P.91) free space warning message.

1 Press [➤] of [Disk space warning] in the [File storage setting (Video image)] screen.

- The [Disk space warning (Video image)] screen appears.


2 Set the threshold value for displaying a warning message by pressing number buttons.

- You can set a value within the range of 100 to 65,535.
- Pressing [Clear] resets the value to 0.

3 Press [OK].

- If you switch to another screen without pressing the [OK] button, the settings will not be applied.
- If the free space of the Image Recording Media reaches the set value, the message "Check memory space (Video image)." appears on the monitor.

Set the value for which to display a warning message.

Reset the setting value.

Confirm the setting.

● Option (Video image)

Set options related to saving video images.


1 Press [➤] of [Option] in the [File storage settings (Video image)] screen.

- The [Option (Video image)] screen appears.

2 Set each of the items.

- [Compressibility]
 - Select the setting for the compression rate for saving files from [Best quality], [High quality], and [Standard quality].
- [Add character display.] check box
 - Add a check mark to save the character information together with the video images.

Press any of the buttons to configure the setting.

Press a check box to add or clear a check mark.

NOTE

- "○ Rec" or the recording time is displayed on the monitor during the recording of video images. If the Image Recording Media is removed and inserted or the processor is turned off while "○ Rec" is displayed, the video images may not be played normally or other problems may occur.
- The maximum recording time varies depending on the setting of [Compressibility]. The maximum recording time is displayed in parenthesis on the monitor during recording.
The maximum recording times for each [Compressibility] setting are as follows:
[Best quality]: 15 minutes / [High quality]: 30 minutes / [Standard quality]: 50 minutes
- Minimum of 2 GB or more of availability is required in a Image Recording Media to record videos in each compression rate for maximum time period.
- The video images are automatically saved to Image Recording Media at the point in time that the maximum recording time is reached.
- Do not remove the Image Recording Media or turn off the processor during the message "Now saving." when displayed on the monitor after stopping the video image recording.
- Once the recording is started, do not stop recording at least for 10 seconds. If the recording is terminated without recording certain amount of time, abnormal values may be displayed on the computer monitor. When the video recording time is short, the length of video recording time indicated in Windows 7 might differ from the actual recorded time.

● Operation checked Software

- 1 The video images (MPEG-2 TS format) generated by this video processor are checked their operation with the following software in “Operation check list”. But this check only guarantees against its performance in our testing environment (PC) but not for your usage environment.

Test environment	
Item	Contents
CPU	Core2 Quad @2.66GHz
Memory	4GB
OS	Windows 7 Ultimate 32bit
Monitor	1920 x 1200 resolution

- 2 For some software, the operation check is done after installing Service Pack.
- 3 As for the device drivers like Graphics Card, etc., replay and edit in the latest condition.
- 4 As to the each operation method of software, contact to the customer service of each software manufacturer.
- 5 In following “Operation check list”, “○” means that check has been done, and “-” means inoperative.
- 6 “○” of “Sound replay” means the sound replayable case during the movie replay.
- 7 “○” is marked in case that all processes like loading, trimming, saving and format conversion are possible against the movie files (MPEG-2 TS type) generated by this product or in the case replayable the edited movie/audio file by each software.

5

Operation check list					
Name of software	Version Edition	Manufacturer	Operation check		
			Movie replay	Sound replay	Edit
Windows Media Player	11	Microsoft Corp.	○	-	-
PowerDVD	Ver.11, 12 Ultra	CyberLink Corp.	○	○	-
Adobe Premiere	Ver.10 Elements	Adobe System Inc.	○	○	○

! CAUTION

Don't make a diagnosis based on the digital file of still image or video image generated by this video processor.

NOTE

- The usage of the monitor of this video processor's system configuration example is recommendable when replaying the still images and video images generated by this video processor by PC. If not, please be in mind that its reproduction characteristic may be inadequate and its color reproducibility and gradation may be variable.
- Use DVI-D output for outputting from PC when replaying the still images and video images. Or its color reproducibility and gradation may become inadequate.
- On rare occasions, the seek bar or the seek time does not be displayed appropriately.

■ Display setting

Set the display format of the monitor.

- [16:9 mode1] Setting for HD monitor display (The character information is displayed on the right side of the endoscopic image.)(The DVI-D signal is output from the DVI connector.)
- [16:9 mode2] Setting for HD monitor display (The character information is displayed at the top and bottom of the endoscopic image.)(The DVI-D signal is output from the DVI connector.)
- [4:3] Setting for SD monitor display(The VGA signal is output from the DVI connector.)

1 Select [16:9 mode1], [16:9 mode2], or [4:3].

NOTE

If [4:3] is selected, images cannot be output from the HD-SDI connector.

■ Scope eject

Use this button when you want to remove and insert an endoscope while the processor is on.

1 Press and hold [Eject] for 2 seconds or longer.

- The message “OK to eject.” appears on the touch panel and the monitor.
- If an endoscope is not connected, the button will not function.

2 Remove the endoscope from the processor.

NOTE

- Reconnect the endoscope when you want to use it again.
- If the “Eject” button does not change from blue to black even when the endoscope is reconnected, please try reconnecting the endoscope.

■ Stopwatch

Use the stopwatch function.

1 Press [Start].

- The stopwatch “S 00:00:00” indication is displayed below the time indication on the monitor.
- To pause the stopwatch, press [Stop]. To cancel the pause, press [Restart].
- Pressing [Reset] resets the counter to “00:00:00.” Pressing [Reset] when the counter is “00:00:00” ends the stopwatch function.

■ Image size

Change the display size of endoscopic images.

- [Full] Set the display size of the endoscopic image to normal size.
- [Medium] Set the display size of the endoscopic image to smaller than normal size.

1 Press [Full] or [Medium].

NOTE

- In twin mode, the display size does not change even if you change the [Image size] setting. The [Image size] setting will be applied when you set [Twin mode] to [Off].
- During analog output(exclude VGA output) of the endoscopic image, the image is always displayed in [Full] size even if the [Image size] setting is changed.

■ Sync out

Set the output from the VIDEO OUT/SYNC OUT connector to a synchronous signal for an external strobe device.

- [Off] Video signals are output.
- [On] Synchronous signals for a strobe device are output.

1 Press [Off] or [On].

■ Freeze release

Set the method for freeze release.

- [Off] The frozen image is maintained even after the still image is saved. Release the frozen image manually.
[On] Freeze is released automatically when the still image is saved.

1 Press [Off] or [On].

■ Manual brightness level

Set the brightness adjustment level by endoscope type for when the exposure control method is set to [MANUAL].

- [Normal] Use this mode for the observation with normal brightness adjustment level.
[Boost] Use this mode for the observation with higher brightness adjustment level. This mode is intended for the use with an otolaryngology endoscope with small diameter.

1 Press [Normal] or [Boost].

⚠ CAUTION

- Use [Normal] mode at first for any observation.
- If the [Boost] mode is used when a gastrointestinal endoscope with large diameter is connected, do not perform observation for a long time with brightness level +4 or higher. Failure to do so could result in a burn injury.

NOTE

If an electronic endoscope is connected, [Boost] cannot be selected.

5

■ Lamp data

Confirm the lamp usage time, configure the lamp ignition method setting, etc.

1 Press [>].

- The [Lamp data] screen appears.

2 Configure the [Lamp ignition] setting or reset the lamp life counter with [Lamp life reset].

● Lamp ignition

Set the lamp ignition method.

- [Manual] The lamp can be turned on manually from the touch panel or keyboard.
[Power on] The lamp turns on automatically when the processor is turned on.

1 Press [Manual] or [Power on].

● Lamp life reset

Reset the lamp life counter.

1 Press [Reset].


- The confirmation screen appears.

2 Press [Yes].

- The lamp life counter is reset to "00:00:00."
- Pressing [No] returns to the previous screen without resetting the lamp life counter.

The lamp life counter is displayed.

The number of times the lamp life counter has been reset is displayed.


Press any of the buttons to configure the setting.

System information

Displays the system information of the processor.

1 Press [>].

- The [System information] screen appears.
- The "Serial No." "Software Ver." "System life" "Lamp life" "Scope name" "Scope serial No." "Scope software Ver." and "Scope connection count" information appear on the screen.

Factory default

Reset all settings to their default values.

1 Press [Default].

- The confirmation screen appears.

2 Press [Yes].

- Pressing [No] returns to the previous screen without resetting the settings.

Film counter reset

Reset the film counter.

1 Press [Reset].

- The film counter indication on the monitor is reset to "1."

Procedure history

Save the patient's name, recording date and time, and other procedure history to "Image Recording Media" (P.91).

1 Press [Output].

- The [Procedure history] screen appears. The Image Recording Media currently connected and the total size and free space of each Image Recording Media is indicated on the screen.
- If you connect Image Recording Media after displaying the screen, press [Refresh] to perform recognition of the Image Recording Media.


2 Press the button of the Image Recording Media to be used as the save to destination.

3 Press [Output].

- The procedure history is saved to Image Recording Media.

Indicates the Image Recording Media currently connected and the total size and free space of each Image Recording Media.

Perform Image Recording Media recognition again.

Output the procedure history.

NOTE

- The procedure history for up to 1,000 procedures is saved to a file in CSV format.
- If the procedure history is saved properly, the message "Procedure history saved." appears on the monitor. If it is not saved properly, the message "Copy failed." appears on the monitor.

■ Configuration copy

Copy the system configuration data of the processor to "Image Recording Media" (P.91).

1 Press [Output].

- The [Configuration copy] screen appears. The Image Recording Media currently connected and the total size and free space of each Image Recording Media is indicated on the screen.
- If you connect Image Recording Media after displaying the screen, press [Refresh] to perform recognition of the Image Recording Media.


2 Press the button of the Image Recording Media to be used as the copy to destination.

3 Press [Output].

- The system configuration data is copied to Image Recording Media.

Indicates the Image Recording Media currently connected and the total size and free space of each Image Recording Media.

Perform Image Recording Media recognition again.

Copy the system configuration.

NOTE

If the configuration data is copied properly, the message "Configuration saved." appears on the monitor. If it is not copied properly, the message "Copy failed." appears on the monitor.

■ Configuration load

Load the system configuration data saved in "Image Recording Media" (P.91) to the processor.

1 Press [Load].

- [Configuration load] screen appears. The Image Recording Media currently connected and the total size and free space of each Image Recording Media is indicated on the screen.
- If you connect Image Recording Media after displaying the screen, press [Refresh] to perform recognition of the Image Recording Media.


2 Press the buttons of the Image Recording Media with the system configuration data to load.

3 Press [Load].

- The system configuration data is loaded to the processor.

Indicates the Image Recording Media currently connected and the total size and free space of each Image Recording Media.

Perform Image Recording Media recognition again.

Load the system configuration data.

NOTE

- The system configuration data is loaded from the file "EPK-i7010ConfigurationDataEx.dat" in Image Recording Media.
The following files cannot be read.
"EPK-i7000ConfigurationData.dat", "EPK-i7010ConfigurationData.dat"
- The files saved other model's system settings cannot be read.
- If there is no data to load in Image Recording Media, the message "No data." appears.
- Configuration load is not available while the endoscope is being connected.

■ Language

Change the display language of the touch panel menus.

1 Press [>].

- The [Language] screen appears.

2 Select [English], [Français], [Deutsch], [Español], [Italiano], [Magyar], [Polski], [Русский], or [中文].

- [English] The menus are displayed in English.
- [Français] The menus are displayed in French.
- [Deutsch] The menus are displayed in German.
- [Español] The menus are displayed in Spanish.
- [Italiano] The menus are displayed in Italian.
- [Magyar] The menus are displayed in Hungarian.
- [Polski] The menus are displayed in Polish.
- [Русский] The menus are displayed in Russian.
- [中文] The menus are displayed in Chinese.
(Only products for China display Chinese.)

The language that is currently specified is indicated by green characters.

NOTE

In case that any incomprehensible language is displayed, user(s) should refer to the English IFU to manipulate.

■ Input method

Set the input method of the keyboard.

[Alphabet input] Input English alphabet letters.

[French input] Input French alphabet letters.


[German input] Input German alphabet letters.

1 Select [Alphabet input], [French input], or [German input].

[Peripheral] Tab

Pressing the [Peripheral] tab displays the following screen. The screen has one page.

Display the setup menu for each of the categories.

Return to the previous screen.

Return to the main screen.

Press any of the buttons to configure the setting.

USB printer control

Capture and output endoscopic images to the USB printer or print captured images.

1 Press [>].

- The [USB printer control] screen appears.

2 Select [1-up], [2-up], or [4-up] to set [Multi picture mode].

- The set split method is displayed in a preview.

3 Press [Capture].


- The endoscopic images are captured and output to the USB printer. "○" is displayed on the monitor during the capture.
- The boxes in the preview display illuminate blue. When [2-up] or [4-up] is set for [Multi picture mode], the boxes for the number of captured images illuminate blue in numerical order, and the box number in which the next image will be captured flashes.

4 Press [Print].

- The captured endoscopic images are printed.
- An animation indicates the print progress.

Split method preview display

Press any of the buttons to configure the setting.

Print a captured image.

Capture an endoscopic image.

CAUTION

- Do not connect more than one USB printer or connect a printer via a USB hub.
- A white "○" is displayed on the monitor until printing has finished. Do not turn the processor off while "○" is displayed.

NOTE

If a USB printer is not connected, [USB printer control] is not available.

RS-232C

Set the device connected to the RS-232C connector.

1 Press [Off] or [Endonet].

Network

Set the network connected to via the RJ45 connector.

1 Press [Off] or [Endoimage2].

Remote

Set the active period for output from the REMOTE connector.

1 Select [200ms], [300ms], or [400ms].

Keyboard Operations


! WARNING

Do not use the keyboard as well as the keyboard menu while the message “Network/File system loading” is being displayed on the bottom right corner of the monitor screen.

! CAUTION

For the EPK-i7010 processor, use the following PENTAX Medical keyboard OS-A83.

Key Assignment

5

! CAUTION

The cursor does not appear when typing the text. In case you need to correct the text, delete the text and re-type it again.




NOTE

- Do not use a sharp object such as the tip of a pen to press the keys on the keyboard.
- Do not use the end of the light guide of the endoscope to press keys.
- Only the keyboard designed specifically for use with the processor can be used (OS-A83).
- The above illustration may not be true to the real keyboard.

Function Keys

The following list shows the functions assigned to the function keys of the keyboard designed specifically for use with the processor.

Function Keys	Key Name	Function Outline
Keyboard OS-A83		
	[Patient]	Display the patient list.
	[New patient]	Register new patient data.
	[Clear]	Clear patient data or user data.
	[User]	Display the user list.
	[SE]	Switch the SE level to a level available in [SE setting]. (It cannot be processed simultaneously with [OE] or the structure enhancement mode [a1 to a6].)
	[CE]	Switch the CE level to a level available in [CE setting]. (It cannot be processed simultaneously with [OE].)
	[TE]	Switch the TE mode to a mode available in [TE setting]. (It cannot be processed simultaneously with [OE].)
	[OE]	Switch the OE mode to a mode available in [OE setting]. (It cannot be processed simultaneously with [SE] [CE] [TE].)
	[Enhance level]	Switch the edge enhancement/structure enhancement level to a level available in [Enhancement setting].
	[Menu]	Display the Setup screen.
	[Twin mode]	Switch twin mode on or off.
	[Freeze]	Switch between still image mode and video mode.
	[Remote out]	Output a signal from the REMOTE connector.
	[Video Rec/Stop]	Start/stop recording video.
	[PC capture]	Output the image to a PC via Endoimage2.
	[USB capture]	Output images to Image Recording Media.
	[i-scan1]	Switch the [i-scan1] i-scan mode on or off.
	[i-scan2]	Switch the [i-scan2] i-scan mode on or off.

Function Keys	Key Name	Function Outline
Keyboard OS-A83		
	[i-scan3]	Switch the [i-scan3] i-scan mode on or off.
	[Noise reduction]	Change the noise reduction level.
	[Image size]	Change the endoscopic image display size.
	[Image rotate]	Change the endoscopic image display orientation.
	[Digital zoom]	Change the endoscopic image display magnification.
	[Counter reset]	Reset the film counter.
	[Character display]	Change the character display on the monitor.
	[Shutter mode]	Change the shutter mode.
	[White balance]	Display the [White balance] screen.
	[Pump]	Display the [Pump] screen.
	[Lamp]	Display the [Lamp] screen.

CAUTION


Do not press multiple keys in combinations not described above.

Operations in the Keyboard Menu Screen

Pressing the [Patient], [User], [Menu], [White balance], [Pump], or [Lamp] key displays the keyboard menu screen on the monitor.

The keyboard menu screen consists of the [Patient list preset], [User list preset], and [Setup] tabs.

- 1 Press the [←] or [→] key to select a tab.
- 2 Press the [↑] or [↓] key to select an item and press the [→] key.
 - Some settings require pressing the [←] or [→] key repeatedly to change the value and some settings require advanced settings to be configured in a sub screen.
- 3 Press the [Enter] or [Esc] key.
 - Pressing the [Enter] key returns to the normal screen where any changes to the settings are reflected.
 - Pressing the [Esc] key returns to the normal screen where any changes to the settings are not reflected.

5 Operations in the [Patient list preset] Tab

■ Registering New and Editing Existing Patient Data

Display the patient list. Then, register new patient data or edit registered patient data. You can register the following patient data items: [Name], [ID], [Age], [Sex], and [Comment].

- 1 Press the [←] or [→] key in the keyboard menu screen to select the [Patient list preset] tab.
 - The [Patient list preset] screen (registration numbers 1 to 10) appears.
- 2 Press the [↑] or [↓] key to select a patient. To register a new patient, select a number with no patient registered. To edit existing patient data, select the corresponding patient.
 - To change the page of the patient list, press the [↑] or [↓] key to align the cursor with [Page */*] and then press [←] or [→].
 - To delete the selected patient data, press the [Backspace] or [Delete] key. Also, to delete all patient data, press the [Clear] key. Pressing each key displays the deletion confirmation message. Select [Yes] or [No] and press the [Enter] key. If you select [Yes], the patient data is deleted from the patient list.

3 Press the [→] key.

- The [Patient data edit No.*] screen for a new patient or the selected patient appears.
- Follow the procedure below to enter or change data in the [Patient data edit No.*] screen. Be sure to enter a name in the [Name] field.

① Press the [↑] or [↓] key to select an item.

② Press the alphanumeric keys to enter or change data.

Press the [Backspace] key in the [Patient data edit No.*] screen to delete one character on the left of the cursor. Press the [Delete] key to delete one character on the right of the cursor.

③ Press the [Enter] or [Esc] key to return to the [Patient list preset] screen.

Pressing the [Enter] key reflects any data entered and changes made and pressing the [Esc] key discards any data entered and changes made.

Pressing the [Esc] key after changing data displays the change cancellation confirmation message. Select [Yes] to return to the [Patient list preset] screen. Select [No] to return to the [Patient data edit No.*] screen.

[Patient data edit No.*] Screen

4 To enter or change data for another patient, press the [↑] or [↓] key to select a patient and repeat Step 3.

5 Press the [Enter] or [Esc] key.

- Press the [Enter] key to return to the normal screen where the information of the patient selected in the [Patient list preset] screen is displayed.
- Press the [Esc] key to return to the normal screen displayed prior to the [Patient list preset] screen.

CAUTION

- If the capture function is enabled, do not use any special symbols such as the ones below.
Examples of symbols not to use: / . * ? " < > : \ |

5

Deleting Patient Data

Press the [Backspace] or [Delete] key while the [Patient list preset] is displayed to delete the data for that patient.

1 Press the [Backspace] or [Delete] key.

- A confirmation message for deleting the data appears.
- Select [Yes] to delete the displayed patient data. Select [No] to keep the registered data.

2 Press the [←] or [→] key to select [Yes] or [No] and press the [Enter] key.

- The [Patient list preset] screen reappears.
- Press the [Esc] key to return to the [Patient list preset] screen without deleting patient data.

[Patient list preset] Screen

■ Registering New and Editing Existing User Data

Display the user preset list. Then, register new user data or edit registered user data. You can register the following user data items: [User name], [Character display], [Assign scope buttons & foot SW], [Picture data], and [Picture data initialize].

1 Press the [←] or [→] key in the keyboard menu screen to select the [User list preset] tab.

- The [User list preset] screen (registration numbers 1 to 10) appears.

[User list preset] Screen

2 Press the [↑] or [↓] key to select a user. To register a new user, select [New] or a number with no user registered. To edit existing user data, select the corresponding user.

- To change the page of the user list, press the [↑] or [↓] key to align the cursor with [Page */*] and then press [←] or [→].
- To delete the selected user data, Press the [Backspace] or [Delete] key. Pressing each key display the deletion confirmation message. Select [Yes] or [No] and press the [Enter] key. If you select [Yes], the selected user data is deleted from the user list.
- To delete all user data, press the [Clear] key. Pressing the key displays the deletion confirmation message. Select [Yes] or [No] and press the [Enter] key. The deletion confirmation message appears again. If you select [Yes] again, all user data is deleted from the user list.

3 Press the [→] key.

- The [User data edit No.*] screen appears.
- Follow the procedure below to enter or change data in the [User data edit No.*] screen.
 - Press the [↑] or [↓] key to select an item.
 - Press the [→] key. The corresponding setup screen appears.
 - Press the alphanumeric keys to enter or change data. Press the [↑], [↓], [←], and [→] keys to select items and configure the settings.
 - Press the [Enter] or [Esc] key. The [User list preset] screen reappears.

Pressing the [Enter] key reflects any data entered and changes made and pressing the [Esc] key discards any data entered and changes made.

[User data edit No.*] Screen

4 To enter or change other data, press the [↑] or [↓] key to select a user and repeat Step 3.

5 Press the [Enter] or [Esc] key.

- The normal screen reappears.
- Press the [Enter] key to return to the normal screen where the information of the user selected in the [User list preset] screen is displayed.
- Press the [Esc] key to return to the normal screen displayed prior to the [User list preset] screen.

■ Settings in the [User data edit No.*] Screen

You can configure the following settings for the user selected in the [User list preset] screen.

[User name]

Enter the name of the user.

During user input, press the [Backspace] key to delete one character on the left of the cursor. Press the [Delete] key to delete one character on the right of the cursor.

[Character display]

Set the items to display on the monitor. Select [On] for each of the [Patient name], [ID], [Age], [Sex], [User name], [Facility name], [Comment], [Date], and [Film counter] items you want to display, and [Off] for each of the items you do not want to display.

[Assign scope buttons & foot SW]

Assign a function to each endoscope buttons 1, 2, 3 and 4 and foot switches L and R.

☞ "Assign scope buttons" (P.48)

☞ "Assign foot switch" (P.51)

[Picture data]

Configure the image processing settings for endoscopic images.

- [SE] Select [Off] or one of six levels for the SE (surface enhancement) display setting.
- [CE] Select [Off] or one of six levels for the CE (contrast enhancement) display setting.
- [Enhancement] Select [Off] or one of two types × six levels for the edge enhancement/structure enhancement display setting.
- [Average/Peak] Select [Average] or [Peak] for the light measuring method.
- [Noise reduction] Select [Off], [Low], [Medium], or [High] for the noise reduction level.
- [Sub screen setting]
 - [Sub screen delay] Set the display of the sub screen to be delayed when the freeze function is executed.
 - [Off] The sub screen is displayed immediately after the freeze function is executed.
 - [On] The sub screen is displayed shortly after the freeze function is executed.
 - [Sub screen position] Select [Upper Right] or [Upper Left] for the position of the sub screen displayed on the monitor when using an endoscope (except when using a side viewing endoscope).
 - [Sub screen position (ERCP)] Select [Upper Right] or [Upper Left] for the position of the sub screen displayed on the monitor when using a side viewing endoscope (ED series).
- [Manual brightness level] Select one of 11 levels for the brightness setting when exposure control is set to [MANUAL].
- [Auto brightness level] Select one of 11 levels for the brightness setting when exposure control is set to [AUTO].

[Picture data initialize]

Restore the endoscopic image settings to their initial values. Select [Yes] and then press the [Enter] key in the confirmation screen that appears. Select [No] or press the [Esc] key to return to the [User data edit No.*] screen without initializing the settings.

Settings in the [Setup] Tab

Configure the settings related to the processor unit and peripheral devices.

[Facility name]

Enter the hospital name.

[Save new patient data]

Set the registration method for registering patient information with the [New patient] key.

- [List+Temporary] The information is saved to the patient list.
- [Temporary] The registered information is only reflected in the monitor screen as temporary data.

[Network setup]

Configure the settings for the network to connect to via the RJ45 connector. The items that can be set in [Network setup] vary depending on the setting of [Network] in the [Peripheral] tab on the touch panel.

- [Endoimage2 setup] (when [Network] is set to [Endoimage2])
 - [Processor setup]
 - [DHCP mode]

[Setup] Screen

Select [Enable] or [Disable]. If you select [Enable], the settings for [IP address], [Subnet mask], and [Default gateway] become invalid.

[IP address] / [Subnet mask] / [Default gateway]

Follow the procedure below to enter or change settings.

- ① Press the [←] or [→] key to select an item when entering the four items separated by “.”
- ② Press the alphanumeric keys or the keys on the numeric keypad to enter or change the numbers. Each press of the [Backspace] key deletes one character on the left of the cursor. Each press of the [Delete] key deletes one character on the right of the cursor.

[Password]

Press the alphanumeric keys to enter a password. Only alphanumeric characters can be used (a to z, A to Z, and 0 to 9) and up to 12 characters can be entered. The password is case sensitive. Leave the field blank if you do not wish to set a password.

[Authorization]

Configure the network settings for the connection destination. Use the numeric keypad to configure of the [Connect ID] settings. Use alphanumeric keys to enter the Connect ID. The number of characters that can be used for the connect ID is between 6 and 12. Alphabet use in the connect ID is case sensitive.

- [iPro setup] (when [Network] is set to [TCP/IP])

[Processor setup]

[DHCP mode]

Select [Enable] or [Disable]. If you select [Enable], the settings for [IP address], [Subnet mask], and [Default gateway] become invalid.

[IP address] / [Subnet mask] / [Default gateway]

Follow the procedure below to enter or change settings.

- ① Press the [←] or [→] key to select an item when entering the four items separated by “.”
- ② Press the alphanumeric keys or the keys on the numeric keypad to enter or change the numbers. Each press of the [Backspace] key deletes one character on the left of the cursor. Each press of the [Delete] key deletes one character on the right of the cursor.

[Password]

Press the alphanumeric keys to enter a password. Only alphanumeric characters can be used (a to z, A to Z, and 0 to 9) and up to 12 characters can be entered. The password is case sensitive. Leave the field blank if you do not wish to set a password.

[iGate setup]

[iGate address] / [Gateway port] / [Active port] / [Maintenance port]

Follow the procedure below to enter or change settings.

- ① Press the [←] or [→] key to select an item when entering the four items separated by “.”
- ② Press the alphanumeric keys or the keys on the numeric keypad to enter or change the numbers. Each press of the [Backspace] key deletes one character on the left of the cursor. Each press of the [Delete] key deletes one character on the right of the cursor.

CAUTION

- Do not operate the processor when there is an IP address conflict.
- Once an IP address conflict occurs, set the correct IP address, and turn the processor off and then on again.

[Lamp]

Turn the lamp on or off.

- [On] Turn on the lamp.
- [Off] Turn off the lamp.

[Pump]

Set the pump level or operate/stop the pump.

- [Pump level] Press [←] or [→] key to set the pump to one of five levels.
- [On] Operate the pump.
- [Off] Turn off the pump.

[White balance]

Set the white balance of the endoscope.

- [White balance] Adjust the white balance of the connected endoscope.
- [R gain] Display the red gain value of the endoscope (K series only). Press the [←] or [→] key to adjust the value. Pressing the [Alt] key and [←] or [→] key simultaneously increases or decreases the value in increments of 10.
- [B gain] Display the blue gain value of the endoscope (K series only). Press the [←] or [→] key to adjust the value. Pressing the [Alt] key and [←] or [→] key simultaneously increases or decreases the value in increments of 10.
- [Reload] Restore the state before the white balance was adjusted.(K series only)

After Use

Follow the procedure below after use.

! CAUTION

Make sure no water spills on the processor when you are performing the after care procedure. In particular, getting the connector parts and ventilation grids wet may cause a malfunction.

1 Turn off the peripheral devices.

! CAUTION

Turning off the processor prior to turning off any connected peripheral devices may cause the peripheral devices to malfunction.

NOTE

When turning off the peripherals of the processor, follow the instructions in the IFU for each peripheral device.

2 Press ① to turn off the processor.

3 Disconnect the power plug, endoscope, and water bottle assembly. (Refer to Figure 6.1.)

NOTE

Always turn OFF the processor BEFORE disconnecting an endoscope.

! WARNING

Recommend irrigation devices that are reprocessed after every patient use or recommend disposable irrigation devices.


Figure 6.1

4 Wipe the surfaces of the processor with a gauze dampened with 70-90% ethyl or isopropyl alcohol.

- The procedure below should be followed to decontaminate the surface of this video processor that is suspected of being contaminated.
- ① Wipe the surface with gauze lightly dampened with a detergent (neutral pH, low foaming, with or without enzymes, mild/compatible with PENTAX Medical product, etc.)
- ② Wipe the same areas with gauze lightly dampened with clean water (to remove residual detergent).
- ③ Wipe the processor with dry gauze or lint-free cloth to remove remaining water/fluids.
- ④ Wipe all of the surfaces with gauze dampened with 70-90% ethyl or isopropyl alcohol.

NOTE

- Do not allow liquid to splash onto/into the processor. In particular, make absolutely sure that the connector interfaces and air intake vents do not get wet. Do not use harsh chemicals, cleaning agents, or the like when wiping the front panel of the processor because doing so may result in a malfunction of the processor. Be sure to use ethanol for disinfection when wiping the surfaces.
- Do not use spray-type medical agents such as rubbing alcohol directly to the video processor. Medical agents might enter the video processor through the air intake vents, etc. and may result in property-damage.

CAUTION

- If the processor is turned off prior to turning off any connected peripheral devices, it may result in malfunction of the peripheral devices.
- Make sure that no water spills on the processor when you perform the aftercare procedure. In particular, getting water into the connector parts and air intake vents may result in malfunction in the processor.

Observe the following points when storing the processor.

- Be sure to turn off the processor and disconnect the power plug prior to storage.
- Do not store the processor in a location with high temperature and humidity or a location where the processor could get wet or be exposed to direct sunlight.
- Store the processor in a location where dust does not enter inside the processor. Wipe off any dust on the processor. Excessive amount of dust accumulation inside the processor may result in malfunction in the processor, smoke or fire from the processor.
- Be careful not to drop the processor or subject it to an excessive mechanical shock. The processor may become unsafe to use or damaged. If the processor is subjected to an excessive mechanical shock, stop using it and contact your PENTAX Medical service facility.

Cleaning and Storage of the Water Bottle Assembly

When washing, disinfecting, and sterilizing the water bottle assembly, be sure to follow the instructions in the IFU for OS-H5.

Storage

Observe the following points on storage

- Be sure to turn off the processor and disconnect the power plug prior to storage.
- Do not store the processor in a very hot and humid location, or where it is likely to be exposed to water or direct sunlight.

WARNING

- Store the processor in a location where dust will not enter inside the unit. Wipe off any dust on the processor. Furthermore, when storing the processor for a long term, take precautions to reduce dust build up within the processor. Excessive amount of dust accumulation inside the unit may cause the processor to malfunction, emit smoke, or catch fire.
- Do NOT store the unit in direct sunlight or where temperature and humidity are high or where it can be exposed to liquids.

CAUTION

Be careful not to drop the processor or subject it to an excessive mechanical shock. The processor may become unsafe to use or damaged. If the processor is subjected to an excessive mechanical shock, stop using it and contact our service personnel.

Replacing the Lamp

Check the LIFE indicator (lamp life indicator) on the touch panel before using the processor. Follow the procedure below to replace the lamp cartridge (OL-X29) if one circle on the LIFE indicator illuminates red and a message prompting you to replace the lamp appears on the monitor.

! CAUTION

- The lamp life is 500 hours.
- When the total usage time of the lamp becomes 500 hours or longer, one of the circles on the LIFE indicator illuminates red and the message "Please replace the lamp." appears on the monitor.
- The lamp life may not be as long as 500 hours depending on the frequency of use.
- Please confirm lighting when you replace the lamp cartridge.  "Power-On and Lamp Inspection" (P.19)
- Replacement lamp should not be stored over an extended period of time.

■ Removing the Lamp Cartridge

- 1 Press  to turn off the processor and disconnect the plug from the power outlet.
- 2 Turn the screw on the lamp protective cover counterclockwise, open the cover, and check the lamp cartridge. (Refer to Figure 6.2.)

! CAUTION

Immediately after use, the metal lamp protective cover and the lamp bulb may be hot. To avoid burns, do not touch these parts immediately after use.


Figure 6.2

- 3 Turn the two screws at the top of the lamp cartridge counterclockwise until the fixed lamp cartridge loosens. (Refer to Figure 6.3.)
- 4 Pull the lamp cartridge out towards you while making sure it does not bang against the lamp cartridge housing.
- 5 Follow local regulations for discarding the old lamp. If unsure of appropriate procedures for lamp disposal, return the lamp cartridge/module to PENTAX Medical after packing it appropriately to avoid any damage during shipping.


Figure 6.3

! WARNING

When the lamp needs to be replaced, we strongly recommends replacing both the lamp and lamp cartridge as a set. Use the PENTAX replacement lamp cartridge OL-X29, which has been developed to ensure optimal brightness and safety. Using a lamp not authorized by PENTAX Medical may result in excessive brightness and/or heat, and it is not known what kind of adverse effect that may have on the patients.

! CAUTION

- To avoid electric shock, do not touch both parts on the processor and the patient at the same time when performing this operation.
- During these procedures, always wear rubber gloves.

■ Installing a Lamp Cartridge

⚠ CAUTION

When performing the following procedure, always wear rubber gloves to prevent oil from your hands from getting on the glass surfaces. Do not directly touch the glass surfaces of the new lamp with your fingers.

- 1 **Hold the lamp cartridge so that the two screws on the lamp cartridge are at the front and fully insert it into the housing while keeping it straight. (Refer to Figure 6.4.)**
- 2 **Turn the screws of the lamp cartridge clockwise to secure the lamp cartridge into place.**
- 3 **Close the lamp protective cover and turn the screw on the right side of the cover clockwise.**

⚠ CAUTION

Be sure to turn screw ① on the lamp cartridge first to secure the lamp cartridge in place.


Figure 6.4

■ Resetting the Lamp Life Indicator

After replacing the lamp cartridge, reset the lamp life indicator on the front panel of the processor.

- 1 **Plug the power cord into the power outlet and turn on the processor.**
- 2 **Press  on the touch panel.**
 - The setup screen appears.
- 3 **Press the [System] tab at the top of the screen.**
 - The system setup screen appears.
- 4 **Press the page [3] tab on the right side of the screen.**
- 5 **Press [>] of [Lamp data].**
 - The [Lamp data] screen appears.
- 6 **Press [Reset] of [Lamp life reset].**
 - The [Lamp life reset] confirmation screen appears.
- 7 **Press [Yes].**
- 8 **Press the ① power switch to turn off the processor and then press ① again to turn the processor back on. Confirm that the lamp life indicator illuminates green.**

Replacing the Fuses

If the power does not come on when you press the power switch, turn off the power switch and then check the following.

- The power cord is securely plugged in.
- The power breaker is activated. When activated, turn on the power switch.

If the power still fails to come on, follow the procedure below to replace the fuses.

NOTE

A flat blade screwdriver is required to replace the fuses.

- 1 **Unplug the power cord.**
- 2 **Insert a screwdriver into the notch on either or both sides of the fuse box immediately above the power inlet socket. (Refer to Figure 6.5.)**
- 3 **Move the inserted screwdriver slightly in the inward direction.**
- 4 **When the fuse box pops out a little, pull it out by hand.**


Figure 6.5

- 5 **Inspect the fuses, and remove and replace any fuse that has blown with a supplied spare fuse. (Refer to Figure 6.6.)**

WARNING

- Always replace a fuse with a supplied spare fuse (Ø5 x 20, T3.15AH, 250V).
- Never bypass a fuse.
- If the supplied spare fuse specified on the rating plate on the rear panel is not available, contact your local PENTAX Medical service facility.
- To avoid electric shock, do not touch both parts on the processor and the patient at the same time when performing this operation.


Figure 6.6

- 6 **Reinsert the fuse box, and push it with a finger until it clicks into place.**
- 7 **Connect the power cord, press the ① power switch on the processor to turn on the power, and check that the LED of ① illuminates.**

NOTE

If the processor cannot be turned on even after the fuses have been replaced with new ones, immediately turn off the power, disconnect the power cord, and contact a PENTAX Medical service facility.

The service life of this product is six years.

Follow the instruction given in the IFU on care after use, maintenance, and usage of this product. If the processor is maintained and used improperly, it could shorten the service life of a product. Some consumable components of this product have the service life of less than six years. Contact your nearest PENTAX Medical service facility every three years for inspection and maintenance of all consumable components except lamp and fuse. Replacement of consumable components may be necessary. If the product is kept being used without periodic maintenance, it might cause a problem in switching between modes from Exposure control mode to Observation mode.

Repair

CAUTION

When the abnormalities of the processor are not improved after it is performed the maintenance according to the Instructions for Use or exchanged a lamp or a fuse, please be sure to request repair to your local PENTAX Medical service facility.

When returning any instrument for repair to PENTAX Medical, please be aware of the following points. Check with your local PENTAX Medical distributor for more details.

- (1) All equipment requiring repair should be shipped in a carrying box with appropriate packing together with details on the instrument damage or problem.
- (2) Write down the repair order number and your contact name, phone number, and shipping address on a piece of paper and include it in the package.

WARNING

- Be sure to submit the instrument to a PENTAX Medical service facility for repairs. Note that PENTAX Medical is in no way liable for any injuries to patients or users, damage or malfunction of the processor, and inadequate disinfection or sterilization as a result of repairs being undertaken by an unauthorized personnel or agency. It must be recognized that PENTAX Medical does not evaluate non-PENTAX Medical parts, components, materials and/or servicing methods. Therefore, questions regarding material compatibility and/or functionality of PENTAX Medical instruments built with the unauthorized, untested and unapproved items, materials, repair/assembly methods must be referred to the third party service organization and/or device remanufacturer. It is unknown to PENTAX Medical if serviced or remanufactured instruments (performed by unauthorized PENTAX Medical entities) which still bear a PENTAX Medical label are within PENTAX Medical device specifications and/or if unauthorized activities have significantly changed the instrument's performance, intended use, safety and/or effectiveness. Ultimately, owners of these medical devices are responsible for selecting an appropriate service facility or vendor whose activities render an instrument to the same expectations and quality of a finished device supplied by the endoscope OEM.
- Be sure to delete any user and patient data from the processor before shipping in order to prevent the disclosure of personal information.
- Technical information related to repairs is listed in the service manual.
- For details about repairs, contact your local PENTAX Medical service facility.

Symptom	Possible Cause	Possible Solution	Reference Page
Cannot turn on the processor.	The power cord is not connected properly.	Securely connect the power cord to the power input socket and isolation transformer.	P.11
	A fuse has blown.	Replace the fuse.	P.84
	Lamp protective cover is not closed.	Close the lamp protective cover and turn the screw knob on the cover clockwise	—
The lamp does not illuminate.	The lamp life is over the limit.	Replace the lamp cartridge if the lamp life indicator on the touch panel is lit red and a message prompting you to replace the lamp is displayed on the monitor.	P.82
No image appears on the monitor.	The monitor settings are not configured properly.	Turn on the monitor.	P.23
		Set the proper video input signal.	P.12
		Reset the adjustment values for the monitor.	*1
	A cable is not connected properly.	Connect all cables properly.	P.12
	The endoscope is not connected properly.	Align the endoscope locking lever to the LOCK position.	P.17
	The lamp is not lit.	Press  on the touch panel to turn on the lamp.	P.26
	The endoscope is not compatible with the processor.	Contact a PENTAX Medical service facility.	—
	Abnormality in the connectors.	Contact a PENTAX Medical service facility.	—
	The monitor is connected to VIDEO OUT/SYNC OUT connector, and [Sync out] is set to [On].	Set [Sync out] to [Off].	P.66
Images are dark.	The lamp is not lit.	Press  on the touch panel to turn on the lamp.	P.26
Text cannot be typed to screen.	The keyboard is not connected properly.	Connect the keyboard to the processor properly.	P.15
The peripheral devices cannot be operated.	The peripheral devices are not connected properly.	Connect the peripheral devices to the processor properly.	P.12 -16
A long beep tone sounds and the [Copy] and [Video] functions do not work.	The settings of the customize buttons, endoscope buttons, or foot switches are not configured properly, or the peripheral devices are not connected properly.	① Change the settings of the [Copy] and [Video] functions in the [Customize button settings], [Assign scope buttons], and [Assign foot switch] screens. If the device that is the output destination is not connected or turned off, cancel the settings. ② Connect the cables of the peripheral devices properly and then turn on the peripheral devices.	P.44 P.48 P.51 P.15
The endoscopic image stopped being displayed on the LCD monitor. Or observation is difficult due to insufficient illumination.	Power outage or there may be a damage in equipment such as lamp, processor, or LCD monitor.	Release the endoscope locking lever and the endoscope's distal tip should be straightened to its neutral position. Then, carefully withdraw the endoscope from the patient.	*2
Software version cannot be displayed correctly.	Attempting to display software version or failed to display software version.	Exit the software version display menu, and enter another menu. Then, reenter the software version display menu. Or, restart the processor.	P.68

Symptom	Possible Cause	Possible Solution	Reference Page
The endoscope buttons do not work.	Cannot communicate properly with the endoscope.	Press [Scope eject] button and remove the endoscope from the processor by aligning the endoscope locking lever to the OPEN position, and then connect the endoscope again. If the function is still not available, release the endoscope locking lever and the endoscope's distal tip should be straightened to its neutral position. Then, carefully withdraw the endoscope from the patient.	P.66
Cannot release freeze.	System malfunction caused by static generation.	Release the endoscope locking lever and the endoscope's distal tip should be straightened to its neutral position. Then, carefully withdraw the endoscope from the patient.	—
The time stamp of the saved still image file/video image file does not match the time when the image was captured.	This processor's "Time zone setting (P.59)" or "Date & Time (P.58)" differs from PC settings.	Set up this processor's "Time zone setting (P.59)" and "Date & Time (P.58)" to become the same as PC settings.	P.58 P.59
		When replayed the still image file/the video image file captured by PC having been set a different "Date & Time" or "Time zone setting", the time stamp on the still image file/video image file may be different from the date and time captured. Therefore, confirm the date and time displayed on the PC monitor.	—
Seek bar or seek time of video image file is not displayed appropriately by reproduction software.	As time property of video image file can't be read out appropriately by reproduction software.	Use other reproduction software to replay.	P.65
Image processing is not reflected in the left side screen after the Twin mode setup.	Setup of Twin mode was made after freezing.	Release freezing and check Twin mode being set, then execute freezing again.	P.41
It cannot be copied from Endonet during a "VCU ready." display.	RS-232C cable is not connected.	Connect RS-232C cable correctly and try to copy again.	—
Image is deteriorated when the saved recording file is played with PC.	This is generated originating in a compression codec processing while setting up the compression ratio Best quality.	Use setting it up with the compression rate for High quality or Standard quality.	P.64
"No Signal" is displayed in a still image.	Checked [Add external.] box of [Option (Still image)] while a still image was captured in a state of no inputting an external input image in P-in-P mode.	Capture a still image again after check that an external input image is inputted correctly. Or when an external input is not used, cancel P-in-P mode and capture a still image again.	—

*1: For details, refer to the Instructions for Use for the monitor.

*2: For details, refer to the Instructions for Use for the endoscope.

NOTE

- If there is no response even if you press any key on the keyboard, turn off the processor and contact a PENTAX Medical service facility.
- If the problem persists even after taking the measures described in the corresponding solutions above, contact a PENTAX Medical service facility.

Error Messages

When an error message is displayed on the monitor, try the following solutions. If there is a problem that cannot be resolved with the following solutions, contact your PENTAX Medical service facility.

Message Displayed on Monitor	Meaning	Solution	Reference Page
Please be prepared to replace the lamp soon.	It is almost time to replace the lamp.	Prepare to replace the lamp.	P.82
Please replace the lamp.	It is time to replace the lamp cartridge.	Replace the lamp.	P.82
Scope is not connected.	An endoscope is not connected.	Connect an endoscope.	P.17
The scope is not supported.	The endoscope is not compatible to this processor.	Connect a compatible endoscope.	P.91
Please contact your PENTAX service facility.	90i series endoscope with outdated software is connected to this video processor.	Contact a PENTAX Medical service facility.	—
Check scope connection(No**).	There is a problem communicating with the endoscope.	Reconnect the endoscope.	P.17
White balance failed! Please retry.	The white balance is abnormal.	Set the white balance again.	P.27
Change exposure control to Auto.	Set the Exposure control to [MANUAL] and adjusted the white balance of the i series endoscope.	Change the method of Exposure control to Auto.	—
Change OE mode to Off.	OE is set up except "Off" and the white balance was performed.	Set up OE "Off" and perform white balance.	—
Check USB memory.	Image Recording Media is not connected or there is insufficient space.	Connect Image Recording Media with sufficient space to the processor.	P.14
Capture failed.	Saving of the still image failed.	Try saving the still image again. If problem is not resolved after the above procedure, contact a PENTAX Medical service facility.	—
Check memory space(Still image).	Available space in the Image Recording Media for saving still images has become as low as the value set for the disk space warning.	Connect Image Recording Media with sufficient space to the processor.	P.14
Check memory space(Video image).	Available space in the Image Recording Media for saving video images has become as low as the value set for the disk space warning.	Connect Image Recording Media with sufficient space to the processor.	P.14
Check video printer.	USB printer is not connected.	Connect USB printer to the processor.	P.14
No response from VCU.	Communication with Endonet is abnormal.	Connect Endonet properly.	—
No VCU connection.	No Endonet is connected.	Connect Endonet properly.	—
AUX lamp ON. Check main lamp.	The auxiliary lamp is lit.	Check main lamp.	P.26
Lamp cooler fan stopped. Please contact your PENTAX service facility.	There is a problem with lamp cooler fan rotation.	Contact a PENTAX Medical service facility.	—
Close profile menu.	The [i-scan profile] or [Twin mode profile] screen is open.	Close the [i-scan profile] or [Twin mode profile] screen.	P.37, P.42
Retry - C function. Freeze release - F function or Frz-key on KBD.	Capturing an endoscope still image failed.	- To try capturing again: Press the Copy button. - To cancel capturing: Press the Freeze button. - When the error message does not disappear even if the copy button and the freeze button are pressed, please try reconnecting the endoscope.	—

Message Displayed on Monitor	Meaning	Solution	Reference Page
System error. Please restart the processor(No**-**).	There was a problem with the system startup.	Restart the processor. When the same message is displayed repeatedly, contact a PENTAX Medical service facility.	P.23
System error. Please contact your PENTAX service facility(No**-**).	There was a problem with the system startup.	Contact a PENTAX Medical service facility.	—
AUX lamp ON. Please contact your PENTAX service facility.	Either the lamp was turned on when an illumination function error occurred, or an illumination function error occurred when turning on the lamp.	Contact a PENTAX Medical service facility.	—
Copy failed.	Copying the procedure history or system configuration to Image Recording Media failed.	Try copying the data again.	—
Data transfer failed.	Failed to transfer data to peripheral devices.	Retry to transfer data after confirming network connection.	—
PC:Not ready.	Capture is done while unconnected to Endoimage2.	Confirm the connection with a PC activating Endoimage2.	—
PC:HDD alert.	There is no sufficient space at PC in which installed Endoimage2.	Connect the processor to a PC with sufficient memory.	—
PC:HDD warning.	The amount of memory in PC installed Endoimage2 decreased to the alert level of insufficient space.	Connect the processor to a PC with sufficient memory.	—
PC:Transfer error.	Failed to transfer data to Endoimage2.	Contact a PENTAX Medical service facility.	—
IP CONFLICT	Setting value of IP address conflicts.	Set the IP address not to conflict and reboot the processor.	—
Video rec error.	Failed to save a movie.	Confirm the correct connection with an image record media or reboot the processor. If problem is not resolved after the above procedure, contact a PENTAX Medical service facility.	—
		When image processing level is set up highly, an error may occur to the compression codec technology. In that case, use setting it up with the compression ratio for High quality or Standard quality.	P.64
No data.	There is no data file to load in an image record media in spite of [Configuration load].	Confirm the existing file to load in an image record media.	—
Changing OE mode failed! Please retry.	When OE was changed into Mode1 or Mode2, the mode change failed.	Perform the mode change again.	—
Changing OE mode Off failed! Check scope connection.(No.**)	Because of the abnormalities occurred in communication with an endoscope when changed OE into Off, the mode change of OE failed.	Reconnect the endoscope.	—
Out of internal memory space.	It has reached the continuous-capture ceiling of the memory capacity.	Try capture again after a while.	—

Other Messages

Message Displayed on Monitor	Meaning	Reference Page
Now printing.	USB printer is printing now.	—
Now saving.	Saving to Image Recording Media.	—
Connect to VCU.	Connected to Endonet.	—
VCU ready.	Ready to communicate with Endonet.	—
Connected to Net	Connected to the network.	—
Capture OK.	Successfully captured the endoscope still image.	—
Check white balance.	Adjust the white balance.	P.27
White balance OK!	Successfully adjusted the white balance.	P.27
Not available.	The auxiliary lamp is lit. (Some functions such as image processing cannot be used.)	—
Procedure history saved.	Saving of the procedure history is complete.	P.68
Configuration saved.	Saving of system configuration is complete.	P.69
Network/File system loading...	Loading the network and file systems required to save the endoscopic still images.  WARNING Do not use the keyboard as well as the keyboard menu while this message is being displayed on the monitor screen.	—
Data transfer OK.	Successfully data transferred to peripheral devices.	—
Waiting for the mode change.	Light control function is switching.	—
PC:Ready.	Ready to start transporting to Endoimage2.	—
PC:Capture OK.	Successfully captured endoscope static image to Endoimage2.	—
PC:Buffer full.	File is being saved to Endoimage2.	—
PC:Command from PC.	Capture required from Endoimage2.	—
OK to eject.	Possible to eject the endoscope.	—
This scope is not supported OE.	The connected K series endoscope is non-OE-compliant.	—
Main:External video	The images input from an external device are displayed on the main screen of the monitor.	P.54

8 Specifications

Model Name			EPK-i7010
Power Requirements	Voltage		230 VAC
	Frequency		50 - 60 Hz
	Input current		2.4 A
	Voltage fluctuation		±10%
Mode of Operation			Continuous Operation
Operating Environment	Ambient temperature		10 - 40°C
	Relative humidity		30 - 85%
	Air pressure		700 - 1060 hPa
Storage and Transport Environment	Ambient temperature		-20 - 60°C
	Relative humidity		0 - 85%
	Air pressure		700 - 1060 hPa
Illumination	Lamp		Xenon lamp (300 W) Model: OL-X29
	Color temperature		6,000 K
	Brightness control		Selection: Automatic/Manual
	Auxiliary lamp		3W White LED
Endoscope Compatibility	Colonoscope Family Gastroscope Family Duodenoscope Family Bronchoscope Family Naso-Pharyngo-Laryngoscope Family		K series video endoscopes, i series video endoscopes, PENTAX Fiberscopes (via video module) ("K series" consists of 70K, 80K, 90K, K10 and J10 series. "i series" consists of 90i and i10 series.)
Video Module Compatibility			PENTAX video module PVK-1070, PVK-1070Z
Air Supply Function	Air pump system		DC diaphragm
	Air supply pressure (at flow rate of 0)		45 - 70 kPa (6.5 - 10.2 PSI)
	Air supply rate (at inlet of water bottle assembly)		Level 1: 2.0 - 2.8 L/min Level 2: 2.9 - 3.4 L/min Level 3: 3.5 - 4.0 L/min Level 4: 4.1 - 4.5 L/min Level 5: 4.6 - 7.2 L/min
Water Supply Function	Water bottle assembly		Capacity = 200 ml
Brightness Control System	Auto		Selection: Average/Peak ±5 step adjustment
	Manual		±5 step adjustment
Color System	Color correction		Red/Blue ±5 step adjustment
Freeze Function			Live video images displayed while freeze mode activated.
Cooling			Forced air cooling
Video Signals	Digital output	HD-SDI	2 sets
		DVI	1 set (combined with VGA)
	Analog output	RGB	1 set
		Y/C	2 sets
		VIDEO OUT	1 set (combined with SYNC OUT)
		VGA	1 set (combined with DVI)
	Analog input	VIDEO IN	1 set
Synchronous Signal	SYNC OUT		1 set (combined with VIDEO OUT)
Voice Signal	Microphone		1 set
Control Signals	RS-232C		1 set
	RJ45		1 set
	Remote		3 sets
	Keyboard		1 set
	Foot switches		1 set
Image Recording Media	USB flash memory, external hard disk drive		
Classification as Electro Medical Equipment	Type of protection against electric shock		Class I equipment, 3 prong plug
	Degree of protection against electric shock		Type BF (Body Floating), using insulated endoscope. Use on a heart is prohibited.
	Degree of explosion proofing		Do not use in potentially flammable surroundings.* * The processor is not suitable for use in a mixture of air and flammable anesthetic gas or a mixture of oxygen/nitrous oxide and flammable anesthetic gas.
Compliance	Designed in accordance with		IEC 60601-1:2005+A1:2012 IEC 60601-2-18:2009
Size (Excluding Protrusions)			400 mm (W) 205 mm (H) 520 mm (D)
Weight			21.5 kg
Electro Magnetic Compatibility			IEC 60601-1-2:2007
IP Classification			IPX0

This product conforms to IEC60601-1-2: 2007: Medical electrical equipment, EMC standard.

Guidance and manufacturer's declaration-electromagnetic emissions

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. This product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration-electromagnetic immunity

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued operation during power mains interruptions, it is recommended that this product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	It is recommended that this product be used apart from other devices operated with large current.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	The recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	The recommended separation distance: $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. d is the recommended separation distance in metres (m).			

Note:

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- Interference may occur in the vicinity of equipment marked with the following symbol:

- Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^{a)} should be less than the compliance level in each frequency range.^{b)}
 - a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this product is used exceeds the applicable RF compliance level above, this product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating this product.
 - b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this product

This product is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note:

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

Contacts

Manufacturer

HOYA Corporation
6-10-1 Nishi-shinjuku,
Shinjuku-ku, Tokyo
160-0023 Japan

Distributors

PENTAX Europe GmbH

Julius-Vosseler-Straße 104
22527 Hamburg, Germany
Tel: +49 40 561 92-0
Fax: +49 40 560 42 13

PENTAX Medical
A Division of PENTAX of America, Inc.

3 Paragon Drive
Montvale, NJ 07645-1782
USA
Tel: +1 201 571 2300
Toll Free: +1 800 431 5880
Fax: +1 201 391 4189

PENTAX Medical Shanghai Co., Ltd.

Room 701, 291 Fumin Road, Shanghai
200031 P. R. China
Tel: +86 21 6170 1555
Fax: +86 21 6170 1655

PENTAX Medical Singapore Pte. Ltd.

438A Alexandra Road, #08-06
Alexandra Technopark, 119967 Singapore
Tel: +65 6507 9266
Fax: +65 6271 1691
Customer Service Toll Free:
400 619 6570 (within China)
1800 2005 968 (within India)
1300 PENTAX (within Australia)

81240

2018. 07 6217001 P227 R03

LCPM 02/2018/09/35011802

Specifications are subject to change without notice and without any obligation on the part of the manufacturer.

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Seite 2

OPERATING MANUAL

for mobile equipment carts uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart and endo-cart with and without isolating transformer

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This is a Class I medical product in accordance with the European Medical Device Directive (MDD) 93/42/EEC, Appendix IX.

The manufacturer declares the conformity of this product with the basic requirements according to MDD 93/42/EEC, Appendix VII, and documents this with the CE designation.

Legal manufacturer :

ITD GmbH
Sportplatzstr. 3
84381 Johanniskirchen, Germany
Tel: + 49 89 61 44 25- 0
Fax: +49 89 61 44 25- 200
Email: sales@itd-cart.com

English

Sales and service:

North America

ITD Corporation
50 RADO DR UNIT A
Naugatuck, CT 06770-2211, USA
Tel: +1 203 714 67-00
Email: sales.usa@itd-cart.com

Europe

ITD GmbH - Office Unterhaching
Gruenwalder Weg 13a
82008 Unterhaching, Germany
Tel: +49 89 61 44 25 0
Email: sales@itd-cart.com

China

ITD Medical Technology Products
(Shanghai) Co., Ltd.
Section A, Room 101, 1st Floor
11 South RiYing Road
Pilot Free Trade Zoon
Shanghai 200131, China
Tel: +86 21 6442 5704
Email: s.pan@itd-cart.com

Australia

ITD Australia Pty Ltd
Aviation Place 3
3043 Tullamarine, Victoria, Australia

Tel: +61 3 9330 50 50
Email: s.oconnor@itd-cart.com

We work constantly to further develop our products. Please understand that we must reserve the right to make changes to the delivery package in terms of form, equipment and technology at any time.

Reproduction, multiplication or translation, in part or in entirety, is not permitted without the written permission of ITD GmbH!

All rights according to copyright law are expressly reserved for ITD GmbH.

Version © 11/2018

Print date 12/2018

These instructions for usage are valid for the following products:

Part number	Description
CC.02xx.xxx / 03xx.xxx	System components and accessories compact-cart
CC.45xx.xxx / 50xx.xxx / 65xx.xxx	Equipment cart compact-cart, 30 U – 40 U
CC.45xx.xxx / 50xx.xxx / 65xx.xxx	System components and accessories compact-cart
DC.42xx.xxx / 53xx.xxx / 64xx.xxx	Equipment cart duo-cart, 21 U – 30 U
DC.40xx.xxx / 50xx.xxx / 60xx.xxx	System components and accessories duo-cart
EB.45xx.xxx / 50xx.xxx / 65xx.xxx	Equipment cart compact-cart Economy
EC.04xx.xxx	Equipment cart endo-cart, 30 U
GW.52xx.xxx / 53xx.xxx / 54xx.xxx	Equipment cart classic-cart, 21 U – 40 U
GW.62xx.xxx / 63xx.xxx / 64xx.xxx	
GW.01xx.xxx / 02xx.xxx / 03xx.xxx	System components and accessories classic-cart
GW.04xx.xxx / 05xx.xxx / 06xx.xxx	
GW.07xx.xxx / 08xx.xxx	
GW.40xx.xxx / 41xx.xxx / 45xx.xxx	
GW.50xx.xxx / 51xx.xxx / 55xx.xxx	
GW.60xx.xxx / 61xx.xxx / 65xx.xxx	
GW.80xx.xxx / 98xx.xxx / 99xx.xxx	
HA.1xx.xxx / 2xx.xxx	System components and accessories flexion-port
HA.45xx.xxx / 5xx.xxx	
HA.60xx.xxx / 65xx.xxx	
NT.20xx.xxx / 40xx.xxx / 41xx.xxx	System components and accessories vexio-cart
NT.50xx.xxx	Equipment cart vexio-cart, 21 U – 50 U
PT.20xx.xxx / 40xx.xxx / 41xx.xxx	System components and accessories pro-cart
PT.50xx.xxx / PT.90xx.xxx	Equipment cart pro-cart, 21 U – 50 U
RS.41xx.xxx / 48xx.xxx / 49xx.xxx	Equipment cart uni-cart, 21 U – 50 U
RS.00xx.xxx / 01xx.xxx / 02xx.xxx	System components and accessories uni-cart
RS.4xx.xxx / 5xx.xxx	
TS.03xx.xxx / 05xx.xxx	System components and accessories flexion-port

These instructions for usage are valid for the following products:

English

Part number	Description
VS.53xx.xxx / 54xx.xxx / 63xx.xxx	Video rack cart classic-cart, 30 U – 40 U
VS.63xx.xxx	
VT.43xx.xxx	Video cart with isolating transformer pro-cart, 30 U
VT.45xx.xxx	Video rack cart with isolating transformer compact-cart, 40 U
VT.54xx.xxx / 64xx.xxx	Video rack cart with isolating transformer classic-cart, 40 U
VW.54xx.xxx / 64xx.xxx	Video cart classic-cart, 40 U
ZV.94xx.xxx / 96xx.xxx / 97xx.xxx	System components and accessories universal
ZV.98xx.xxx / 99xx.xxx	

KD.7xxx.xxx / 8xxx.xxx / 9xxx.xxx	Customer-specific, mobile equipment carrier, series uni-cart, pro-cart, duo-cart, classic-cart, compact-cart and endo-cart
KU.0xxx.xxx – KU.9xxx.xxx	Customer-specific, mobile equipment carrier, series uni-cart, vexio-cart, pro-cart, duo-cart, classic-cart, compact-cart and endo-cart
OC.0xxx.xxx – OC.9xxx.xxx OC.1 – OC.999	
OM.0xxx.xxx – OM.9xxx.xxx	

Content

1	Important information	
1.1	Correct usage	30
1.2	General explanation of the symbols.....	31
1.3	Safety instructions.....	33
2	Assembly	
2.1	Completeness.....	34
2.2	Loading	34
2.3	Loading sequence.....	34
2.4	Risk through mechanical instability.....	35
2.5	Castors.....	35
2.6	Load capacity	35
2.7	Assembly / handling	35
2.8	Additional mounting of system components.....	36

3	Electrical safety	
3.1	Positioning of the electrical equipment	36
3.2	Power columns (compact-cart, classic-cart, endo-cart) and vertical columns (uni-cart, vexio-cart and pro-cart)	36
3.3	Gases	37
3.4	Equipotential bonding (POAG)	37
3.5	Isolating transformer - leakage current	37
3.6	Earth-Leakage-Guard (ELG)	37
3.7	Cable connection	39
3.8	Combination of equipment	39
3.9	EMC	40
3.10	Excluded from the final electrical inspection of system components and accessories	40
3.11	Minimum safety	41
4	Transportation	
4.1	Safe transportation on castors	41
4.2	Safe transportation by carrying	41
5	Mechanical and electrical height adjustment	42
6	Support arms	
6.1	Wiring	43
6.2	Horizontal swivel	43
6.3	Tilting / rotating equipment	43
6.4	Variable height support arms (flexion-port)	44
6.5	Dismantling and repositioning system components and accessories	44
6.6	Intended use when manoeuvring the mobile equipment carts	44
6.7	Operation of the height adjustable support arms (flexion-port) on the monitor shelf, 2-fold	45
7	Miscellaneous	
7.1	Cleaning and disinfection	45
7.2	Service / Repair	46
7.3	Environmental conditions	46
7.4	Disposal	46
7.5	Spare parts	46
8	Accessories	46
9	Maintenance	47
10	Technical data	
10.1	Load capacity uni-cart	48
10.2	Load capacity vexio-cart	48
10.3	Load capacity pro-cart	48
10.4	Load capacity duo-cart	48
10.5	Load capacity compact-cart	48
10.6	Load capacity classic-cart / endo-cart	48
10.7	Load capacity modul-port	49
10.8	Load capacity flexion-port	49

1 Important information

All products from ITD GmbH are manufactured for a long and trouble-free service life. Development, construction, sales and production are certified at ITD GmbH according to ISO 13485.

This is a basis for:

- highest quality and a long service life
- easy, safe and ergonomic operation
- functional design
- optimisation for the planned usage

The products satisfy the requirements of the European Medical Device Directive (EMDD) and bear the CE designation.

- Carefully read these instructions for usage from the beginning in order to become familiar with the functions step-by-step.
- Please be sure to address all questions or concerns to the manufacturer.
- The mobile equipment carriers are only intended for the use as described.
- These instructions are to be kept for the service life of the product.

The system configurator is to make the instructions for usage of the overall configuration available to the final customer.

We expressly note here that the system configurator is responsible for the observance of IEC 60601-1 and the EMC norm IEC 60601-1-2 in the valid version!

1.1 Correct usage

The functions of the mobile equipment carriers of ITD GmbH are:

- mounting medical and IEC-tested devices according to the permitted load information in compliance with the requirements of IEC 60601-1 in the currently valid edition.
- Connection and distribution of mains voltage from the local electricity outlet as well as from data lines.
- Consolidation of original ITD system components and accessories.

Using the mobile equipment carrier, the medical equipment can be transported inside the building or arbitrarily positioned in the area before and after the application. Therefore, a flexible, economical efficiency of all equipment is possible. In addition, it is easier to clean the floor area.

1.2 General explanation of the symbols

"ON" (voltage)
lights up green

"OFF" (voltage)

"ON" lights up
green/"OFF"
(pressure actuated)

Equipotential bonding: Identifies equipotential bonding terminals on the housing of the isolating transformer; equipotential bonding ensures that resistance between all conductive materials is sufficiently low.

Connection to protective conductor:
Connects conductors, equipment units, conducting parts, main earth terminals and earth

Conductive castors:
Conductive castors are indicated by a flash or a yellow point

Move only with arm folded up


Use the handle to push


Follow the instruction manual


Alternating current


Total output:
Total power that can be supplied by the individual plug in locations must not exceed the total power.

Total load rating (basic frame)
max. total load rating (= total of all load ratings of mounted system components)
Please refer to the label for appropriate load.

Load rating (system components)
Please refer to the label for appropriate load.

Humidity limit

IMPORTANT INFORMATION

English

General warning sign:

This symbol is used at the socket strips. The overall rating given on the name plate must not be exceeded.

Only suitable for the interior

Heavy object:

It is recommended that two people lift the uni-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart to prevent injury.

Distributed by


Adjustment of the clamping force (tilt and swivel unit)

Set load:

Describes the load range and also the direction of rotation for setting them.


Risk of tilting:

It is imperative to observe the sequence of loading and unloading mobile equipment carts.


Brake:

Information about positioning of the support arm in the event of device acceptance and indicates the direction for release and blocking the locking function.

Operating position:

This symbol describes the permitted operating positions (right/left) and indicates a prevailing danger of tilting when changing over.

Locking function:

This describes the closing / opening direction for components.

Do not push:

It is forbidden to push the equipment cart above the handle because of a risk of tipping over.

1.3 Safety instructions

General:

- Only those mobile equipment carriers may be operated whose main voltage equipment has been tested and approved by appropriate, qualified personnel!
- Ensure the isolating transformer is only connected to a power supply with a functioning protective earth connection that complies with the provisions of the IEC 60364-7-710 2012-10 "Electrical installations in buildings, part 7-710: Requirements for special installations or locations – Medical locations". If in doubt, contact a specialist in the electrical trade or an authorised employee of the hospital's engineering team.
- Personnel (hospital and service personnel) working directly or indirectly with a mobile equipment carrier must be instructed!
- Setting adjustments may only be carried out by qualified personnel.
- Repairs may only be carried out by ITD GmbH.

Safe working with the equipment cart:

- The device is only disconnected from the mains supply when the power plug is removed from the socket.

Operation:

- With every change of location it must be ensured that no one is injured or objects damaged!

Connections:

- In the USA and Canada, use a hospital-grade connection cable when connecting isolating transformers in the 115 V mode and in Japan use the optional Japanese connection cable.
- When connecting to the multiple-socket strip, secure the plug using an optional available plug pull-out lock.
- Only devices that fulfil the requirements of IEC 60601-1 or are IEC-tested may be connected to the sockets/connecting lines.
- Additional medical equipment with connecting bolts for equipotential bonding are to be connected with the green-yellow cable to the optional equipotential bonding connecting bolts!

- **Caution:** The overall rating given on the nameplate must not be exceeded. Please observe that no further multiple sockets must be connected to an existing multiple socket.

Load capacity:

- The total weight of the equipment and the accessories on the mobile equipment carrier may not exceed the permitted payload weight (see load capacity sticker on the base frame).
- The surface load imprinted on the system components may not be exceeded!
- The load shown on the fittings (e.g. infusion stand, jointed arms) may not be exceeded!

Infection protection:

- Hygiene regulations are to be observed when cleaning!
- Give only cleaned and disinfected equipment and fittings to a service technician for maintenance and repair work!

Environmental protection:

- Dispose of all cleaning and disinfection agent residue in a manner not harmful to the environment!

2 Assembly

2.1 Completeness

Unpack the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system and check whether all the parts you have ordered have been included in the consignment.

2.2 Loading

Place the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system on a level, horizontal surface. Place the appliances in the cart taking into consideration the loading sequences. Loading and assembly may only be effected if the cart is separated from current. We recommend the optional accessory (e. g. tension belt) to be applied for the safety of each appliance.


The maximum loading is specified in section 10. If the tilt stability is not adequate, the cart must be clearly labelled in conformity with the 5° labelling specification.

2.3 Loading sequence

Make sure that all the equipment placed on the uni-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system is stowed securely to prevent slipping, tipping, falling down, etc. by taking appropriate measures (especially also during transport). It is advisable to position all heavy components on the uni-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system using two people. You should take into account the fact that the centre of gravity changes with the loading.


The cart should be loaded in the following sequence:

- shelves and drawers from the bottom to top
- support arm system (rigid, traversable, tiltable, height adjustable, single or multiple) should be burdened last.

The cart should be unloaded in the following sequence:

- support arm system (rigid, traversable, tiltable, height adjustable, single or multiple) should be relieved first.
- shelves and drawers should be unloaded from top to bottom.


For transportation please refer to section 4 of this manual.

2.4 Risk of mechanical instability

The overall system must be in accordance with IEC 60601-1.

2.5 Castors

The uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system has twin swivel castors with brakes. You should check to see that the wheel locks are working properly before using the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system. When parking the cart or stopping during transport, all roller brakes on the equipment cart must be pressed (roller locking device). Accordingly, all roller brakes must be released before moving and transport.

The castor must be tested every 12 months for safety as well as the castor attachment pin. If these should be lost, please contact your supplier immediately.

2.6 Load capacity

It is not permitted to exceed the load rating of the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system.

You should observe the maximum load capacity of the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system (see section 10 of these operating instructions).

2.7 Assembly / Handling

2.7.1 Shelves

Shelves can be removed or installed in another position. Loosen the screws, re-position the shelf and then screw the shelf tightly in position again. Check the safety-earth resistance.

2.7.2 Drawers

Drawer unit blocks (pro-cart) are provided with a latch mechanism.

To release the lock on the pro-cart, the operating handle on the front cover must be pulled upwards.

The drawers can be lifted out in the extended position.

A label strip can be affixed to the front trim (except pro-cart).

Insert drawer during transport.

2.7.3 The isolating transformer

Follow the instruction manual for the isolating transformer. The isolating transformer is mounted in a housing under the base of the basic frame. This is assembled at the factory.

2.8 Additional mounting of system components


Retrospective installation of ITD system components must only be carried out by specialists in accordance with the specifications of the installation instructions supplied. The modified overall system must be rechecked in accordance with the IEC 60601-1.

3 Electrical safety

3.1 Positioning of the electrical equipment

Make sure that electrical equipment on the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system does not get wet.

Under no condition should you position the products, which may lead to loss of their fluidness, over the electrical equipment or extension lead in which fluidity may be permeate.

3.2 Energy column (compact-cart, classic-cart, endo-cart) and vertical column (uni-cart, vexio-cart, pro-cart)

The socket strips or the cabling system are located in the left or right power columns of classic-cart, compact-cart or endo-cart, with uni-cart and pro-cart it is in the vertical column. The vexio-cart's socket strips and cabling system are located beneath the base (cabling can be routed through optional cable ducts, which are attached laterally to the vertical column).

Power columns are located on the left and right behind the vertical extrusion and permit optimum accommodation of the power equipment cables available. Do not drill holes in the power columns / vertical columns because power cables may be located inside the columns.

3.3 Gases

Electrical equipment should not be operated in the vicinity of gases, e.g. flammable gas used in anaesthesia or similar gases. The user is responsible for maintaining this requirement and for compliance with IEC 60601-1-2 standard.

3.4 Equipotential bonding (POAG)

Equipotential bonding should be carried out for the uni-cart, pro-cart, classic-cart, compact-cart or endo-cart carrier system with an isolating transformer. vexio-cart offers an optional POAG set without isolating transformer. The equipotential bonding cable should be first connected to the base of the uni-cart, vexio-cart, pro-cart, classic-cart, compact-cart or endo-cart carrier system and then be connected to the equipotential bonding plug in the room. Next, connect the POAG supply cables to the POAG pins of the multiple sockets and the appliance.

English

3.5 Isolating transformer – leakage current

The uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system is intended to provide a practical and mobile workstation for medical electrical equipment. Total earth leakage currents must not exceed a maximum value of 0.5 mA so that the entire medical electrical system can be compliant with IEC 60601-1. If the sum of the ground leakage currents would exceed the tolerance range, the system must be supplied through a safety transformer.

If an isolating transformer is not mounted, the multiple socket strip / auxiliary socket strip system of the uni-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system must not be used to connect up equipment that is not in conformity with requirements for earth leakage currents under IEC 60601-1.

If an isolating transformer is mounted, the total power uptake of all the connected units must remain within the rated value for the isolating transformer.


3.6 Earth-leakage-guard (ELG)

The following descriptions are only applicable for the carrier systems in which the item "isolating transformer with insulating monitor" is installed.

3.6.1 Control panel of the earth-leakage-guard (ELG)

At the devices with integrated-earth-leakage guard the evaluation unit is in the transformer case, the separate control panel is placed in one of the shelves. Both components are connected with an interface cord, which lays in one of the profiles.

Control panel of the earth-leakage-guard (optional):

- 1 Power on control (green)
- 2 Signals exceeding dielectric resistance (yellow)
- 3 Signals exceeding temperature (yellow)
- 4 Reset-button for fault release
- 5 Test-button starting self-test function manually

3.6.2 Application requirements

The earth-leakage-guard monitors the dielectric resistance for equipment (where safety isolation is necessary) which are connected to the isolation transformers. Simultaneously, the temperature of the transformer is evaluated by a built-in thermal switch which responds when triggered. The evaluation procedure is processor controlled.

3.6.3 Application requirements

Turn the main power switch of the isolating transformer to "ON". Within 5 s the self-test of the ELG is automatically activated. After the self-test of the ELG is concluded, the green LED (Power on control) glows permanently. During operation this test is repeated (automatically) in a time cycle of 8 hours. This self-test may also be activated manually through the test-button.

The control procedure for the manual self-test is as follows:

- An insulation error is simulated, the yellow LED (signals exceeding dielectric resistance) glows, a warn signal (2.4 kHz) is activated. Both signals stop after approx. 5 s.
- Following this, a temperature error is simulated. The yellow LED (signals exceeding temperature) glows, a warn signal (2.4 kHz pulse tone) is activated. Both signals stop after approx. 5 s.

An error situation is recognized as follows:

- In case of an insulation error, the yellow LED (signals exceeding dielectric resistance) glows permanently, the warn signal with 2.4 kHz is activated permanently. The warn signal can be deactivated through the reset-button for fault release. The LED glows until the error is eliminated.
- If the earth-leakage-guard is turned "OFF" in the meantime and the error is not eliminated, when turning the ELG back "ON" the above described procedure begins anew. With the occurrence of an isolation error the acoustic and optical alarm stays activated (even if the isolation error is corrected on its own) until deactivated through the fault release:

First activation of the fault release: acoustic alarm is deactivated

Second activation of the fault release: optical alarm is deactivated

- In case of a temperature error, the yellow LED (signals exceeding temperature) glows permanently and the warn signal with 2.4 kHz is activated intermittent. The warn signal can be deactivated through the reset-button for fault release. The LED glows until the error is eliminated. If the earth-leakage-guard is turned "OFF" in the meantime and the error is not eliminated, when turning the earth-leakage guard back ON the above described procedure begins anew.
- Should an insulation error and a temperature error occur simultaneously, the dielectric resistance fault has priority (acoustic).

Self test for the Earth-Leakage Guard

In addition to the manual self-test, during operation the earth-leakage-guard repeats the self-test in a time cycle of 8 hours and each time the unit is turned "ON".

From outside, you will not notice the self-test, the time period for this process is approx. 5 seconds. In case of fault, the green LED (power on control) blinks in a frequency of 0.5 Hz and with the same frequency an acoustic alarm is activated. The fault signal can not be deactivated through the reset-button.

3.6.4 Trouble shooting

Do not attempt to repair on your own. Warranty will be invalid if unauthorized repair has been carried out. Due to safety reasons repair and maintenance may only be done by the manufacturer.

Please note:

For further technical data and information, please refer to the separate and comprehensive instructions for use, enclosed with the product, isolating transformer and insulating monitor.

3.7 Cable connection

Users of the uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier system without an isolating transformer must ensure that the detachment of the cable plug connection between the cart's multiple power outlet strip and devices can only occur using tools. Our line of accessories includes a selection of covers for multiple power outlet strips.

3.8 Combination of equipment

The following should be observed for combination of equipment on the carrier system:

- Auxiliary equipment connected to analogue and digital interfaces of the equipment must be certified in compliance with the relevant EN specifications (e. g. IEC 60950 for data-processing equipment and IEC 60601-1 for medical electrical equipment).
- Furthermore, all configurations must comply with the valid version of the standard IEC 60601-1. Anyone connecting additional equipment to the signal input or signal output is a system configurer and is therefore responsible for ensuring compliance with the valid version of the standard IEC 60601-1.

If you have any questions, you should contact your local dealer or technical services.

Please note:


Make sure that this is also the case for adaptation of equipment in the power supply circuit (e. g. multiple socket strip).

3.9 EMC


Electromagnetic compatibility of the medical electrical equipment positioned on the cart has to be checked by the overall system configurer. If you use a different equipment combination, you should check electromagnetic compatibility between individual items of equipment before using in a medical setting.

Customer-specific equipment trollies used within the nuclear spin environment must be tested by the customer for suitability for use due to the ferromagnetic materials they contain.

ITD GmbH excludes any liability in this respect.

Using the carrier system with an isolating transformer within an environment subject to explosive hazard is not permissible.

3.10 Excluded from the final electrical inspection of system components and accessories

ITD GmbH exclude the following system components and accessories from the final electrical inspection:

- Multiple socket strips without additional protective conductors that are not wired in the mounting
- ME cables and appliance cables included
- POAG plates and cables included

- Non-electrified equipment carts and support systems
- Height adjustments and attachment parts to height adjustments
- Handles, mouse pads, drawers, drawer bodies and attachments (bottle mounting brackets, baskets, camera mounting brackets, infusion tripods, ...)
- Isolation transformers, which are not mounted but instead leave the ITD as a single part.
- Keyboard extensions and extendable shelves
- Computer mounts at the top and bottom
- Conductive castors
- Support arms installed and monitor mounting brackets
- Secondary power circuit with insulation monitors are only excluded from the dielectric strength inspection!

3.11 Minimum safety

ITD GmbH is not aware of any item of equipment or accessory that reduces the minimum safety of the system. Only equipment not presenting a hazard may be used.
If necessary, this should be clarified by means of a risk analysis (ISO 14971).

4 Transportation

4.1 Safe transportation on castors

Please ensure that before the mobile use of the equipment cart:

- All equipment / products in transit are secured against detachment.
- All swivel arms must be swivelled in and secured.
- The feeder from the local outlet must be disconnected.
- The brakes of the castors are released.

The equipment cart should only be moved at a maximum walking pace of 6 km/h. The velocity is significantly reduced when knocked down by bedplates, cables as well as hoses. It is safer to deliver by driving on ramps having a max. gradient of 10° so that the equipment cart can be brought to a standstill at any time.

Although all precautionary measures have been observed to guarantee maximum stability for this product, it is important to pay attention to unevenness in floors, lift door frames, cables, etc. in order to prevent accidents.

Fundamentally, the requirements of the IEC 60601-1 apply.

4.2 Safe transportation by carrying

Handles are not provided for lifting the mobile equipment carts, but exclusively have a function for pushing. Lifting and carrying the mobile device carrier must only be carried out by two persons using the extensions of the base.

Fundamentally, the requirements of the IEC 60601-1 apply.

5 Mechanical and electrical height adjustment

Specific safety regulations must be observed for the „mechanical height adjustment“ using the „gas pressure“, as well as for the electro-mechanical height adjustment using „Linear drive“, in accordance with IEC 60601-1 „Mechanical Risk in Conjunction with Moving Parts“. Hereby:

- Take into consideration and adhere to the permitted distance between moving parts, in accordance with IEC 60601-1 in Table 20 (ISO 13857:2008).
- Products with height adjustment are manufactured and supplied ex works conforming to the standard, under consideration of the permitted safety distances. Due to the equipping of or replacement with ME devices and / or component, these distances change. This can result in a mechanical risk. The respective person who configures the system is responsible for adhering to the minimum distances required.
- The overall weight of the devices and accessories installed must not exceed the specified maximum overall payload of the height adjustment. Overloading results in damage to the height adjustment and loss of the warranty.
- Stored energy is released with the mechanical height adjustment using gas pressure. Thereby, for unloaded systems sudden, unbraked activation of the height adjustment can result in injuries and damage.
 - o In order to prevent injuries and damage, before installation and removal of the devices, place the height adjustment at the uppermost position („energy-free“).
 - o The height adjustable supporting arm system „flexion-port“ must also be fixed and secured at the uppermost position („energy-free“) with the aid of the clamping lever [refer to the separate user instructions for the „flexion-port“, as well as the risk notes on the supporting arm system].
- Unintended activation of the electro-mechanical height adjustment using the manual button can also result in injuries and damage.
 - o In order to prevent injuries and damage, before installation and removal of the devices, disconnect the height adjustment from the power supply.
 - o Servicing and maintenance tasks in the „interior area“ of the height adjustment, i.e. in the covered area within the supporting column not accessible from outside, must only be carried out by specialists.
 - o Caution: If operating the height adjustment using a remote control, make sure that no persons are in the hazard area.

6 Support arms

6.1 Wiring


Please observe the following instructions:

- So that there is no damage to the cable or failure of the device when swivelling, the cable must be sufficiently dimensioned.
- Possible sagging cables must under no circumstances be used as a handle.
- Please make sure that the enclosed assembly material is correctly applied, in accordance with the assembly instructions.
- When swivelling the arms, pay attention to possible cable loops present.

6.2 Horizontal swivel


Make sure that the swivel range of the system components matches the dimensions of the equipment and the ambient conditions in the working environment.

When system components and attached equipment are swivelled horizontally, they must not collide with other equipment, with other system components or with the wall. Any collision may result in damage to equipment and injuries to persons.


6.3 Tilting / rotating equipment

If system components are tiltable or rotatable, it is important to ascertain whether the clamping force is appropriate for the unit being fixed in place. If the force applied is incorrectly adjusted, the equipment is at risk of tilting over. Adjustment must therefore be carried out to ensure that the unit can be slightly tilted or rotated, while the unit remains stable in any desired position.

6.4 Variable height support arms (flexion-port)

When loading the system components that can be height adjusted, unconditionally observe the minimum and maximum permitted total weight. Also, due to safety reasons, please make sure that the space below the height-adjustable support arm (flexion-port) remains clear. In order to be able to set the support arm to the load, it must be placed into the horizontal position.

6.5 Dismantling and repositioning system components and accessories

When system components and accessories are dismantled or repositioned, it is important to remove any units mounted on these elements before any changes are made. If this concerns the disassembly / assembly of the height-adjustable support arms flexion-port, these must first be placed in the uppermost position and the clamping (brake) must be determined (refer to decal).

6.6 Intended use when manoeuvring the mobile equipment carts

When manoeuvring the mobile equipment carts it is imperative to make sure that the support arms attached (swivelling, height adjustable) are placed into the respective park position and, as necessary, locked. Otherwise, the stability cannot be guaranteed (refer to Point 2.4).

6.7 Operation of the height adjustable, support arm system (flexion-port) on the monitor shelf, 2-fold

When operating the height adjustable support arm system (flexion-port) on the monitor shelf, 2-fold, always make sure that the swivelling arm that provides the changeover, is always placed parallel to the front side of the cart (engaged position). Thereby, please observe that only two operating positions of the support arm (right / left) are permitted (refer to „Operating position“ in Chapter „1.2 General symbol explanation“, Page 31).

By a changeover of the operating position from left to right, or vice versa, the locking elements must be pushed downwards and the height adjustable support arm swivelled to the other side. Hereby, for flexion-port with swivelling arm make sure that it is folded when changing over. Otherwise, the tilting stability cannot be guaranteed (refer to Point 2.4).

7 Miscellaneous

7.1 Cleaning and disinfection

Caution: Disconnect from power before cleaning and disinfection!

The equipment cart uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart must be cleaned using commercial all-purpose cleaning agents (neutral cleaning agents).

For disinfecting, commercial disinfectants approved for disinfecting surfaces or wipe disinfection can be used. The disinfection agents must be used solely as disinfection for wiping, in accordance with the manufacturer's specification.

For example, ITD have carried out tests using the following disinfectants:

Product	Manufacturer
Bacillol plus	Bode
Cleanisept Wipes	Dr. Schumacher
Mikrobac Tissues	Bode
Mikrozid Sensitive Wipes	Schülke
Terralin Protect	Schülke
Incidin PLUS	Ecolab
Incidin Foam	Ecolab

If complete disinfection is required, assemblies can be disassembled by a specialist and wipe disinfected when disassembled.

7.2 Service / Repair

The cart should be always be cleaned and disinfected with a suitable cleaning agent before any service operations are undertaken, and before the cart is returned for purposes of repair. Repairs to the cart should only be effected by professional personnel. We recommend consulting ITD GmbH on all matters relating to service activities.

7.3 Environmental conditions

English

The uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart carrier systems are designed for standard operation in hospitals and medical practices.

Operation:

Ambient temperature:	10° C to 40° C
Air humidity:	30 % to 75 %
Air pressure:	700 hPa to 1060 hPa
Protection class:	IP20

Transport/storage:

Ambient temperature:	-25°C to 70°C
Air humidity:	10% to 95%
Air pressure:	500 hPa to 1200 hPa

7.4 Disposal

Separate Collection for Electrical and Electronic Equipment in compliance with Waste Electrical and Electronic Equipment Directive WEEE (registration number for Germany: DE35464575). All electrical and electronic equipment provided with systems released after 13 August 2005 is marked with a Separate Collection for Electrical and Electronic Equipment symbol, indicating that this equipment must undergo separate collection for disposal, in countries where EU directive 2002/96/EC is in effect.

7.5 Spare parts

Only spare parts authorized by ITD may be used. A sticker with an order number is attached to the base of your cart. All order numbers and the associated spare parts are archived at ITD GmbH. Spare parts may be obtained from ITD GmbH.

8 Accessories

A comprehensive range of accessories is provided in our catalogues or under www.itd-cart.com (information for dealers).

9 Maintenance

The uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart or endo-cart has been developed and constructed for many years of trouble-free use.

Check the functional capability of the following parts every 12 months in order to guarantee safety.

Monitor shelf:

- Swivels and tilts easily without too much play.

Shelves:

- Check whether the mounting screws have been tightened and whether the shelf is stable and flat.

Castors:

- Check that the castors run freely and the brakes operate correctly.
- Check that the 4 bolts holding the castors to the bottom side of the base and the castors themselves are seated firmly in their holders.
- The running surfaces of castors must be free from contamination, so that they can function correctly.

Socket strips:

- Check the main cable for damage and firm seating.

Auxiliary sockets:

- Check the cable for damage and firm seating.

Variable height support arms (flexion-port):

- The height adjustment functions freely, the raising force is adjusted to the weight of the device.

Support arms:

- Rotating and tilting functions smoothly, without excess play.

Isolating transformers:

- Safety-relevant checks of isolating transformers.

Serial number:

- Compare the serial number of the cart with the data of the equipment log book.

Fuses:

- Check whether the correct fuses have been installed.

If you encounter any problems during these checks you should contact your supplier immediately.

TECHNICAL DATA

10 Technical data

10.1 Load capacity uni-cart

- Basic frame, total payload from 50 kg / 110 lbs
- Shelf 10 kg / 22 lbs
- Rack shelf 20 kg / 44 lbs
- Drawer 3 kg / 6.6 lbs (lockable: 10 kg / 22 lbs)
- Monitor holder 14 kg / 30.8 lbs
- Mouse pad 3 kg / 6.6 lbs

10.2 Load capacity vexio-cart

- Basic frame, total payload 65 kg / 143 lbs
- Shelf 15 kg / 33 lbs
- Drawer unit 3 kg / 6.6 lbs
- Monitor holder 14 kg / 30.8 lbs
- Mouse pad 3 kg / 6.6 lbs

10.3 Load capacity pro-cart

- Basic frame, total payload 80 kg / 176 lbs
- Shelf 20 kg / 44 lbs
- Drawer unit 15 kg + 3 kg / 33 lbs + 6.6 lbs
- Monitor holder 14 kg / 30.8 lbs
- Mouse pad 3 kg / 6.6 lbs

10.4 Load capacity duo-cart

- Basic frame, total payload 80 kg / 176 lbs
- Shelf 50 kg / 110 lbs (pull-out: 20 kg / 44 lbs)
- Drawer unit 3 kg / 6.6 lbs (lockable: 20 kg / 44 lbs)
- Mouse pad 3 kg / 6.6 lbs

10.5 Load capacity compact-cart

- Basic frame „Profi“, total payload 180 kg / 396 lbs
- Basic frame „Economy“, total payload 150 kg / 330 lbs
- Shelf 50 kg / 110 lbs
- Drawer unit 3 kg / 6.6 lbs
- Monitor shelf max. 35 kg / 77 lbs (depends on type)
- Mouse pad 3 kg / 6.6 lbs

10.6 Load capacity classic-cart / endo-cart

- Basic frame, total payload 150 kg / 330 lbs
- Shelf 50 kg / 110 lbs (pull-out: 20 kg / 44 lbs)
- Drawer unit 3 kg / 6.6 lbs
- Monitor shelf max. 35 kg / 77 lbs (depends on type)
- Mouse pad 3 kg / 6.6 lbs

10.7 Load capacity modul-port (stationary carrier systems)

- Support profile, total added load 25-150 kg / 55-330 lbs
depends on length
- Support arm up to 23 kg / 50.6 lbs
- Swivel arm, 1-fold up to 23 kg / 50.6 lbs
- Swivel arm, 2-fold up to 18 kg / 39.6 lbs
- Monitor mount with VESA 75/100 adapter up to 18 kg / 39.6 lbs
- Monitor mount with universal adapter up to 14 kg / 30.8 lbs
- Monitor mount with Table Top Mount Adapter up to 14 kg / 30.8 lbs
- Shelf 10 kg / 22 lbs
- Drawer unit 3 kg / 6.6 lbs
- Keyboard holder 5 kg / 11 lbs
- Mouse pad 3 kg / 6.6 lbs

10.8 Load capacity flexion-port (variable height support arms)

- flexion-port (depends on modell) 3-10 kg / 6.6-22 lbs or
8-14 kg / 17.6-30.8 lbs
11-20 kg / 22.4-44 lbs
- Tilt and swivel unit bis 14 kg / up to 30.8 lbs
- Post ("Down-Post") 10 kg / 22 lbs
- Mouse pad 3 kg / 6.6 lbs

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China

ITD Medical Technology Products
(Shanghai) Co., Ltd.
Section A, Room 101, 1st Floor
11 South RiYing Road
Pilot Free Trade Zoon
Shanghai 200131, China
Tel: +86 21 6442 5704
Email: s.pan@itd-cart.com

Australia

ITD Australia Pty Ltd
Aviation Place 3
3043 Tullamarine, Victoria, Australia

Tel: +61 3 9330 50 50
Email: s.oconnor@itd-cart.com

Responsible for content: ITD GmbH.

Operating manual uni-cart, vexio-cart, pro-cart, duo-cart, compact-cart, classic-cart and endo-cart 11/2018

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LEAKAGE TESTER (SHA-P5)

STANDARD SET LIST

This package includes the following items:

- Leakage Tester (Main Body and Tube)
- Instructions for Use

⚠ CAUTION:

Never drop the Leakage Tester or subject it to severe impact, as this could compromise the functionality and/or safety of the unit. Should this equipment be mishandled or dropped, do not use it. Return it to an authorized PENTAX service facility for inspection or repair.


Symbols

Symbol	Description	Symbol	Description
	Date of manufacture		Batch Code
	The CE mark confirms the compliance to applicable European(EU) requirements.		
	Medical device in the EU		

Symbol	Description
	*If this bar code symbol is displayed on the product, see the following explanation. This is a UDI (Unique Device Identification) code required by Unique Device Identification System designed to adequately identify devices through distribution and use. The following information is coded in 2D bar code (GS1 Data Matrix). - (01) GS1 Commodity code (Global Trade Item Number) - (11) Production date - (10) Lot Number

Intended use


This product is intended to be used to perform leakage testing on PENTAX Medical flexible endoscopes.

- The PENTAX Leakage Tester allows for simple two-stage testing of the watertight integrity of the PENTAX Endoscopes. Air pressure is delivered by means of a hand operated bulb, eliminating the need for any electro-mechanical devices. The main body and the tube are packed separately. Before using, connect the tube to the air connector of the main body by turning the knurled screw of the tube clockwise.

STAGE I - DRY LEAKAGE TEST

BEFORE IMMERSION, PENTAX Endoscopes should be tested for any major loss of integrity in their watertight construction (example: major tear in the instrument channel). If applicable, all valves, inlet seals, distal rubber tips and other removable components should be detached from the endoscope prior to leakage test.


- 1) Secure the endoscope connector to the venting connector on the endoscope light guide. If applicable, the soaking cap should be securely attached to/ over the electrical contacts. The endoscope connector and the venting connector on the endoscope MUST be dry before connecting. Proper connection will require alignment of the venting connector pin and clockwise rotation of the Leakage Tester connector.
- 2) Depress the upper half of the toggle switch to retain pressure within the endoscope.
- 3) Pressurize the endoscope by pumping the hand bulb until the indicator on the gauge is in the GREEN zone. Do NOT pressurize beyond the GREEN zone into the RED zone, as it may cause serious damage to the endoscope.

⚠ CAUTION:

Do NOT squeeze the hand bulb too abruptly or forcefully. Doing so will cause the Leakage Tester to be pressurized to the point where the indicator enters into the red zone. Be careful not to let the indicator enter into the red zone.

NOTE:


During the leakage test procedure, the insertion tube of the endoscope should be flexed in various positions, and the distal bending section should be fully angulated in all directions to confirm the absence of a leak.


- 4) Observe the gauge pressure to determine if the indicator is standing still in the GREEN zone. If the indicator drops from the GREEN zone rapidly, a major leak may be indicated.

NOTE:

Avoid depressing the lower half of the toggle switch during the leakage test, as this will cause the endoscope to deflate and the leak tester gauge needle to return to its baseline zero position. If this happens, you must repeat the leakage test. NEVER IMMERSE the entire endoscope if the gauge indicator does not remain in the GREEN zone. Instead, contact your PENTAX Service Facility.

STAGE II - WET LEAKAGE TEST


- After determining the absence of any major leak in Stage I testing, PENTAX Endoscopes should be immersed to test for loss of integrity in their watertight construction.

- With the Leakage Tester securely attached to the endoscope, the soaking cap (if applicable) securely attached to/over the endoscope, ALL valves and inlet seals removed endoscope, and the endoscope pressurized so that the gauge indicator is in the GREEN zone, immerse the entire endoscope in clean water.

NOTE:

Only the endoscope connector and the distal portion of its rubber tubing should be immersed. NEVER immerse the entire Leakage Tester.

- Observe the instrument carefully while fully angulating the distal tip of the endoscope. A few bubbles may occur initially released from recessed areas of the endoscope. This is normal. If a continuous stream of bubbles is observed from the same spot, a leak is indicated. Immediately remove the endoscope from the water. Do NOT use the endoscope. Thoroughly dry the endoscope and contact your PENTAX Service Facility.
- For a product equipped with a forceps elevator, inspect the forceps elevator while raising and lowering it to confirm if air bubbles are present. If a continuous stream of air bubbles is observed from the same spot, stop using the endoscope and send it to PENTAX Medical for repair.

- After removing the endoscope from water, release the air pressure within the endoscope by pushing the toggle switch downward. After the gauge indicates 'zero', disconnect the Leakage Tester from the endoscope.

You may proceed with cleaning and either disinfection or sterilization of the endoscope as described in the endoscope Instructions for Use (Reprocessing).

NOTE:

NEVER connect or disconnect the Leakage Tester under water. This will cause leakage of water into the endoscope and Leakage Tester.

CARE AFTER USE

- Wipe the rubber tubing with a gauze lightly dampened with a cleaning detergent (neutral pH, low foaming, mild/compatible with PENTAX products, etc.).
- Wipe the rubber tubing with a gauze lightly dampened with potable water (to remove residual detergent). Do NOT allow any fluids to enter the distal endoscope connector (green) or the tubing.
- Wipe the rubber tubing with a dry gauze or lint-free cloth to remove remaining water/fluids.
- Wipe the entire Leakage Tester with a gauze dampened with 70-90% medical grade ethyl alcohol.


⚠ CAUTION:

Keep the entire hand-held section of the Leakage Tester (gauge, bulb, etc) from getting wet. NEVER completely immerse the entire Leakage Tester.

- After cleaning and drying as per the steps above, store the Leakage Tester in a clean, dry, and well-ventilated area. NEVER store it in an area where it could get wet or be exposed to any adverse environmental conditions such as high temperature, humidity, direct sunlight, dust, salt, etc.

⚠ WARNING:

Follow the national or local laws/guidelines to appropriately dispose of the consumables. Failure to do so may create a risk to the environment.


⚠ CAUTION:

When storing the Leakage Tester, raise the toggle switch. Storing the Leakage Tester with the toggle switch in the locked position for extended periods can damage the device.

⚠ CAUTION:

Never attempt to unscrew/rotate the gray ring above the hand bulb. This can damage the Leakage Tester.

NOTE:

As per the Regulation (EU) 2017/745 (EU-MDR) users should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user is established.

Manufacturer

HOYA Corporation

6-10-1 Nishi-shinjuku, Shinjuku-ku, Tokyo 160-0023 Japan

Distributors

PENTAX Medical (Shanghai) Co. Ltd.

Room 701, 291 Fumin Road, Shanghai, China 200031

Tel: +86-21-61701555 Fax: +86-21-61701655

PENTAX Medical Singapore Pte. Ltd.

438A Alexandra Road, #08-06 Alexandra Technopark, Singapore 119967

Tel: +65-6271-1669 Fax: +65-6271-1691

PENTAX Medical Company

A Division of PENTAX of America, Inc.

3 Paragon Drive, Montvale, New Jersey 07645-1782 USA

Toll Free: +1-800-431-5880

Tel: +1-201-571-2300 Fax: +1-201-391-4189

PENTAX Canada Inc.

6715 Millcreek Drive, Unit1 Mississauga, Ontario L5N 5V2 Canada

Tel: +1-905-286-5570 Fax: +1-905-286-5571

PENTAX Europe GmbH

Julius Vosseler Strasse 104, 22527 Hamburg, Germany

Tel: +49-40-561-920 Fax: +49-40-560-4213

PENTAX U.K. Ltd.

PENTAX House, Heron Drive, Langley SLOUGH SL3 8PN, Great Britain

Tel: +44-1753-792792 Fax: +44-1753-792794

Authorized representative
in the European Community

88228

PENTAX™

LEAKAGE TESTER (SHA-P5) 2/2

LCPM: 02/2016/09/35102702 2020.02 6217001 Z935 R05

PENTAX Medical
CO₂ Insufflator and Irrigation Pump

PENTAX Medical is committed to improving the quality of patient care and productivity of healthcare delivery around the globe by providing the best endoscopic products and services, with a focus on clinically relevant innovation, simplicity, and quality.

The CO₂ Solution for GI Endoscopy

The PENTAX Medical CO₂ Insufflator is a versatile system that connects to either a wall or tank CO₂ source. Combined with the full line of the disposable PENTAX Medical DispoCap CO₂, the PENTAX Medical CO₂ Insufflator allows clinicians to offer the optimal procedure for the patient.

The PENTAX Medical CO₂ Insufflator offers a complete range of CO₂ solutions to suit individual facility needs for PENTAX Medical GI endoscopes. With direct connection to sterile water and any CO₂ tank or wall source – PENTAX Medical CO₂ Insufflator has the comprehensive solution for your facility and patients.

Benefits of Using CO₂

- CO₂ is easier on the patient both during and after the procedure due to rapid absorption by surrounding tissue.¹
- Minimises patient discomfort associated with bloating and cramping from room air.¹
- Less sedation required during the procedure.²
- Reduces the number of patients that need to be seen by the physician in recovery.²

Daily Use PENTAX Medical DispoCap CO₂ Full disposable system minimises infection control risks.

Warm Water Technology


Deliver irrigation for GI endoscopy procedures with a warm water option by utilising the PENTAX Medical Irrigation Pump. Combine the pump with PENTAX Medical GI endoscopes and irrigate through the auxiliary water jet or biopsy channel to improve procedural visibility.

- Greater visibility for improved diagnosis and polyp detection
- Auxiliary water jet or biopsy channel irrigation capability
- Warm water technology
- Compatible with PENTAX Medical Irrigation Tubing and Connectors
- Easy to use water bottle cradle
- Automatic prime button provides instant irrigation upon foot pedal depression


Auxiliary Water Jet Connection

CO₂ Insufflator Features

Irrigation Pump Features

Specifications

PENTAX Medical CO₂ Insufflation

REF#	EGA-501P
Dimensions [mm]	197x 121x 349
Weight [Kg]	4.8
Power Requirements [V]	100-240

PENTAX Medical Irrigation Pump

REF#	EGA-500P
Dimensions [mm]	197x 121x 349
Weight [Kg]	4.8
Power Requirements [V]	100-240
Max Flow Rate:	
Auxiliary Water Channel [ml/min]	0-300
Biopsy Channel [ml/min]	0-600

Easy Set Up and Installation


1. Choose a CO₂ Source

2. Connect Hose and Tubing to Insufflator

3. Ready for Procedure

Compatible Accessories

Tubing

REF# 100150CO2P	PENTAX Medical DispoCap CO ₂
REF # 100160P	PENTAX Medical DispoCap Air
REF# 100130P	PENTAX Medical Irrigation Tubing for EGA-500P
REF# 100551P	PENTAX Medical CO ₂ Source Tubing

Connectors

REF# 100116P	PENTAX Medical 24hr Use Waterjet Connector
REF# 100242P	PENTAX Medical Single Use Waterjet Connector

Source: 1. Dellon, E.S. (2009). The use of carbon dioxide for insufflation during GI endoscopy: a systematic review. Gastrointestinal Endoscopy, 69(4), 843-849.
2. Author: Mirosław Szura, MD PhD, Rafał Solecki, MD PhD, Krzysztof Bucki, MD, Andrzej Matyja, MD PhD. MEDICINA Specialist Diagnostic & Therapeutic Centre
Title: Helpfulness of Carbon Dioxide Insufflation During Colonoscopy - SAGES 2011 Program Number P211

PENTAX Europe GmbH
Julius-Vosseler-Straße 104
22527 Hamburg
Germany
Tel.: +49 40 / 5 61 92 - 0
Fax: +49 40 / 5 60 42 13
E-mail: medical@pentax.de

PENTAX U.K. Limited
PENTAX House
Heron Drive, Langley
Slough SL3 8PN
United Kingdom
Tel.: +44 17 53 / 79 27 33
Fax: +44 17 53 / 79 27 94
E-mail: medical@pentax.co.uk

PENTAX France Life Care S.A.S.
112 quai de Bezons
B.P. 204
95106 ARGENTEUIL CEDEX
France
Tel.: +33 1 / 30 25 94 78
Fax: +33 1 / 30 25 74 45
E-mail: contact.medical@pentax.fr

PENTAX Nederland B.V.
Amsterdamseweg 29
1420 CA Uithoorn
Netherlands
Tel.: +31 88 / 5 30 30 30
Fax: +31 88 / 5 30 30 40
E-mail: medical@pentax.nl

PENTAX Italia S.r.l.
Via Dione Cassio, 15
20138 Milano
Italy
Tel.: +39 / 02 50 99 58 1
Fax: +39 / 02 50 99 58 60
E-mail: marketing.lifecare@pentaxitalia.it

**SIMMEDICA – Sistemas Integrales
de Medicina, S.A.**
Avenida del Sistema Solar 25
28830 San Fernando de Henares - Madrid
Spain
Tel.: +34 91 / 301 62 40
Fax: +34 91 / 751 31 15
E-mail: sim@simmedica.com

PENTAX Europe GmbH
Representative office in Moscow
Sadovnicheskaya str, 82, build 2, entrance 6
Regus, Business center „Aurora”
Office 2012, 2013
Moscow, 115035
Russian Federation
Tel.: +7 495 / 792 52 00
Fax: +7 495 / 792 35 66

PENTAX Turkey
Veko Giz Plaza, Meydan Sokak
No:3/43 343396
Maslak – Istanbul
Turkey
Tel.: +90 212 / 705 05 26
Fax: +90 212 / 705 05 00

HOYA Corporation
PENTAX Life Care Division
1-1-110, Tsutsujigaoka
Akishima-shi
196-0012 Tokyo
Japan
Tel.: +81 3 / 39 60 51 55
Fax: +81 3 / 53 92 67 24

Manufactured in the U. S. A. by: MEDIVATORS, 3150 Pollok Drive Conroe, TX 77303
U. S. A. Phone: +1-973-890-7220
PENTAX® is a registered trademark of HOYA Corporation.

© 2013 ©HOYA Corporation
PENTAX Lifecare Div