

## Detro® Activ

Pre - cleaner and Disinfectant for Medical Instrument and Endoscopes

## Detro® Cid Activ

Pre - cleaner and Disinfectant for Medical Instrument and Endoscopes

### USAGE AREAS:

- Used for high level disinfection of rigid and flexible endoscopes, surgical instruments, anesthesia and dental surgery instruments, other medical instruments made of stainless steel, plastic, silicon, glass, acrylate and rubber.

### PRODUCT PROPERTIES:



- Concentrated product in powder form.
  - Generates peracetic acid by reaction of sodium percarbonate and TAED (tetraacetylenediamine) with water.
  - Cleaning and disinfection in one step.
  - Effective cleaning power.
  - Extends instrument life thanks to corrosion inhibitors.
  - Aldehyde and phenol free.
- 93/42 EEC Medical device directive class IIb.

### MICROBIOLOGICAL ACTIVITY:

- Bactericidal and fungicidal in 5 minutes at %0,5 concentration,
- Tuberculocidal in 5 minutes at %2 concentration.
- Virucidal and sporicidal in 10 minutes at %2 concentration.

### USAGE:

- The solution is prepared by solving the product in water in a clean and closed container at the suitable concentration (5 gr product into 1 lt water for 0,5% concentration, 20 gr product into 1 lt water for 2% concentration).
- The instruments are immersed thoroughly so that all the surfaces and cavities are covered by the solution.
- The container is closed.
- The contact time is completed.
- Then the instruments are rinsed with sterile or distilled water.
- The ready to use solution in the container can be used daily.
- It is recommended to change the solution in case of intensive contamination.
- Solution activity can be controlled by test strips.
- Do not leave instruments in solution longer than the contact time.
- Do not use on aluminium, nickel and chromium instruments.



### COMPOSITION:

- Peracetic acid generated by reaction of sodium percarbonate and TAED (tetraacetylenediamine) with water, enzymes (amylase, lipase, protease).
- Generates >1000 ppm peracetic acid at %2 concentration.

### PACKAGE:

QUANTITY	PACKING TYPE	BARCODE
2,5 KG	4 x BUCKET	8680097100344
750 GR	4 x BUCKET	8680097103277

### COMPOSITION:

- Peracetic acid generated by reaction of sodium percarbonate and TAED (tetraacetylenediamine) with water.
- Generates >1000 ppm peracetic acid at %2 concentration.

### PACKAGE:

QUANTITY	PACKING TYPE	BARCODE
2,5 KG	4 x BUCKET	8680097100603
750 GR	4 x BUCKET	8680097103215