# AUTOCAL

## Calibration Serum

**REF** 13160 4 x for 5 ml Calibrator  
**IVD**

### Intended use

AUTOCAL is intended for calibration of HUMAN methods on clinical-chemistry analyzers of AUTOHUMALYZER and HumaStar series.

### Summary

AUTOCAL is a lyophilised calibrator based on human serum.

The concentrations and activities have been selected to ensure optimum calibration.

### Reagents and Contents

Lyophilisate for 5.0 ml per vial AUTOCAL. Reactive components in the lyophilisate: Human serum with chemical additives and extracts of animal origin.

The concentrations / activities of the components are lot-specific. The exact calibration values are listed in the attached table.

### Reagent Preparation

Carefully open a bottle, avoiding the loss of lyophilisate, and add exactly 5.0 ml of distilled/deionised water. Carefully close the bottle with the rubber stopper and dissolve the contents completely by occasional gentle swirling within 30 minutes at room temperature. Avoid the formation of foam.

### Storage and Stability

Store at 2...8°C. Stability of the lyophilised calibrator serum at 2...8°C: up to the stated expiration date.

Stability of the components in the reconstituted calibrator:

at	25°C	at least 8 hours
at	2...8°C	at least 7 days
at	-20°C	at least 1 month (when frozen once)

### Limitations:

Stability of bilirubin in the reconstituted calibrator (when stored protected from light): at 2...8°C: 8 hours.

For acid phosphatase determinations, the reconstituted material should be stabilised by adding one drop (25 - 30 µl) of acetic acid (0.7 mol/l) to 1 ml of the calibrator.

At 2...8°C stabilised acid phosphatase is stable for 2 days.

Alkaline phosphatase levels will rise over the stability period after reconstitution. The reconstituted serum should be allowed to stand for 2 hours at 25 °C before measurement.

Store calibrators tightly capped when not in use.

### Calibration Values

The calibration values were determined using the analytical procedure described in the "Method" section of the respective HUMAN test kit (see assigned values list).

### Precautions and Warnings

This product has been prepared exclusively from the blood of donors tested individually and shown by approved methods to be non-reactive for HBsAg and antibodies to HCV and HIV 1 and 2. However, since no test method can rule out the potential risk of infection with absolute certainty, the material should be treated just as carefully as a patient specimen. In the event of exposure the directives of the responsible health authorities should be followed.<sup>1,2</sup>

### Performance Characteristics

Typical performance data can be found in the Verification Report, accessible via

[www.human.de/data/gb/vr/ca-auto.pdf](http://www.human.de/data/gb/vr/ca-auto.pdf) or  
[www.human-de.com/data/gb/vr/ca-auto.pdf](http://www.human-de.com/data/gb/vr/ca-auto.pdf)

### References

1. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. Fed. Register July 1, 6, 267-280 (1998)
2. Council Directive (90/679/EEC). Official Journal of the European Communities. No. L374 from Dec. 31, 1-12 (1990)

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