


CE-IMMUNDIAGNOSTIKA GmbH Karl-Landsteiner-Straße 6 D-69151 Neckargemünd Phone: +49 6223 80094 00 Fax: +49 6223 80094 99 www.ce-immundiagnostika.com	
Instruction for use Use by professionals only Tests per ml: max. 25 with drop size 40µl when using separate volumetric pipettes.	 0483
Revision:	29/07-2019
Product-Name:	Product-Code:
AHG polyspez.-rabbit-	AHG-green-Ka
Humanglobulin - Coombsserum	
All described test methods are only valid for manual applications as recommended in this instruction. The user must determine their suitability for use in other techniques (automates, gel-cards, others) according to recognized techniques in individual responsibility. Only for in-vitro diagnostic laboratory use. Store at + 2 - 8 °C when not in use.	

Product description:	The reagent is produced from rabbit. The rabbit serum pool is absorbed to ensure that no antibodies against other human antigens are present. These sera are adjusted with bovine albumin in isotonic saline.
Indirect anti-human globulin test:	This method is a two-step test for detection of incomplete antibodies and demonstrations of Antigens with sera which can only be used in coombs test. 1. Place 1 drop of the patient serum to be tested or one drop of the specific blood typing test serum into a tube. 2. Add 1 drop of a 2-3% erythrocytes suspension in 0,9% isotonic saline solution or 22% BSA. (test erythrocytes or donor erythrocytes for antigen testing) 3. Incubate 15-60 minutes at 37° C and test for agglutination. 4. Wash erythrocytes 3 times in isotonic saline and discard last supernatant. 5. Add to sediment 2 drops of anti human globulin serum and centrifuge for one minute at 400 xg (1500 UpM) or at a suitable other rpm for an appropriate time. 6. Shake gently and check for agglutination. 7. The addition of Liss or Bovine Albumin allows reduction of the incubation time and the reaction may be enhanced.
Direct anti-human globulin test:	This is a one step test for demonstration of erythrocyte-bound antibodies (for example in haemolytic syndrome in the new born, haemolytic anaemia) 1. Wash erythrocytes to be tested three times in isotonic saline and discard the supplement. 2. Prepare a 2-3% erythrocyte suspension in isotonic saline solution. 3. Place 1 drop of erythrocyte suspension into a tube. 4. Add 1-2 drops of coombs serum and mix. 5. Centrifuge for one minute at 400 xg (1500 UpM) or other time with adequate rpm. 6. Shake gently and check for agglutination.
Cross match:	The cross match is an antibody diagnostic test to reveal confusion and detect irregular antibodies. The cross match is used in three forms: Mayor test: Recipient serum + donor erythrocytes Minor test: Donor serum + recipients erythrocytes Controls: Recipient serum + recipient erythrocytes All three are performed as outlined below in three-step tube test:
Tube Centrifuge test – Three step-test:	First step: Saline milieu at room temperature 1. For each vial of an antibody detection test mark one test tube and additional one for control of autoantibodies. 2. Place 1 drop of the serum under examination into each test tube. 3. Add 1 drop of test-erythrocytes suspension to the appropriate labelled test tubes and 1 drop of the one washed patients cells (3-5%) into the control tube and mix well. 4. Incubate for five minutes, centrifuge for one minute at 400 xg or other time with adequate rpm, shake tubes gently and read. Second step: Albumin milieu at 37° C 1. Pipette 1 drop of 30% bovine albumin in each tube and mix well. 2. Incubate the tests for 30 - 60 minutes at 37° C 3. Centrifuge for 1 minute at 400 xg or other time with adequate rpm, shake tubes gently and read. Third step: Coombstest 1. Wash the erythrocytes 2 or 3 times in saline. 2. Add 1 drop of Anti-Human globulin serum and mix well. 3. Centrifuge for 1 minute at 400 xg or other time with adequate rpm, shake tubes gently and read. If the results of the Coombstest are doubtful or negative run a check with a coombs-control.
Performance data:	The reagent shows equal or better performance criteria as comparable reagents on the market. It fulfills the requirements of the specifications of CE-Immundiagnostika GmbH.
Other Tets methods:	As alternative method it is recommended to run with Albu-Liss (modified Liss) If there is suspicion of cold antibodies, it is recommended to incubate in a saline milieu at 4° C for 30 minutes. Conservation: Na-Azid (< 0,1%) Shelf-life: to be used until expiry date.