MACOPHARMA SSP/SSP+ (FR)

DESCRIPTION DU DISPOSITIF

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MACOPHARMA SSP/SSP+ (EN)

DESCRIPTION OF THE DEVICE

The devices include a bag designed with one or two connection modes, according to the references. Each bag contains a certain volume of platelet additive solution SSP or SSP+ according to the following composition :

- **Composition SSP** : Na₂-citraté 2H₂O : 2.94g, Na-acetate 3H₂O : 4.08g, NaCl : 6.75g, Water For Injection 1000 ml, pH 7.2
- **Composition SSP+** : Na₂-citraté 2H₂O : 3.18g, Na-acetate 3H₂O : 4.42g, NaH₂PO₄ 2H₂O : 1.05g, Na₂HPO₄ : 3.05g, KCl : 0.37g, MgCl₂ 6H₂O : 0.30g, NaCl : 4.05g, Water For Injection 1000 ml, pH 7.2

All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male fetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia. All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free

INDICATIONS

The platelet additive solutions SSP and SSP+ are designed to partially replace plasma in the preparation and storage of buffy-coat derived platelet concentrates or apheresis platelet units. The recommended ratios are as follows: up to 70% SSP / 30% plasma, up to 80% SSP+ / 20% plasma. These solutions enable platelets to be stored at 22°C ± 2°C, under gentle agitation, for up to 5 days (SSP) or 7 days (SSP+) following collection and according to local regulations.

SIDE EFFECTS

There are no known side effects with the use of SSP or SSP+.

WARNING

- The device is intended to be used by authorised staff in transfusion centers and hospitals.
- Avoid contact with pointed or sharp-edged objects.
- Check the integrity of the device and of the packaging (secondary and tertiary) before use.
- The secondary packaging is a protective packaging and not the sterile barrier. The solution bag is the sterile barrier.
- Do not use if there are any visible signs of deterioration of the tertiary or secondary packaging and/ or of the device or if they have been accidentally opened.
- Open the tertiary packaging by locating one of the two notches below the two heat sealings and pulling on the notch to tear open the tertiary packaging.
- Do not remove the tertiary packaging from the cardboard box.
- Moisture or micro-droplets may be present inside the secondary and tertiary packaging.
- Open the secondary packaging by pulling on the tab at the corner.
- Use the device immediately after opening the secondary packaging.
- Do not use the device.
- Do not use if the solution is cloudy and/ or any particles can be seen in the solution.
- Do not use if any labels are missing.
- Do not directly infuse the solution
- Do not re-use the device as there could be a risk of contamination. And additionally, the remaining quantity of solution might be insufficient.
- Take out the luer's cap.
- Use the solution at a temperature of between 18°C and 24°C.
- Store the device at a temperature of between 0°C and 35°C.
- Following Carelide's situation and the conditions of safety of use and maintenance of performance which we are committed, we are producing under derogation with the agreement of ANSM.

INSTRUCTIONS FOR USE

Manifold pooling (use of 4 to 6 Buffy-Coats)

- Make a sterile connection between each buffy-coat tubing and each end-line tubing of the manifold device.
- Connect by sterile connection the SSP or SSP+ bag tubing to the remaining manifold set end-tubing.
- Hang all the buffy-coat bags at the same level and the SSP or SSP+ bag at a higher level. Check the integrity of all sterile connections.
- Allow the buffy-coats to drain into the pooling bag.
- Close the clamp above the pooling bag and then allow the SSP or SSP+ solution to drain into each buffy-coat bag.
- Agitate in order to recover a maximum of buffy-coat & open the clamp, allow the SSP/SSP+ solution to drain into the pooling bag.
- Remove the manifold set and then follow standard procedures to obtain a pooled buffy-coat derived platelet concentrate

Cascade pooling (4-6 Buffy Coats)

- Connect by sterile connection the SSP or SSP+ bag tubing to the bleed line of the first buffy-coat bag in the train.
- Make a sterile connection between the buffy-coat plasma line (tube equipped with a breakaway cannula) of the first buffy-coat bag and the bleed line of the second buffy-coat bag.
- Then make a sterile connection between the second buffy-coat and the third buffy-coat. Repeat the operation with each buffy-coat until a maximum of 6 buffy-coats attached.
- (Option: connect the tubing of the plasma line of the last buffy-coat bag to a platelet reco-very set).
- Hang the chain vertically by suspending the SSP or SSP+ bag with the chain of buffy-coats underneath. Check the integrity of all the connections and allow the buffy-coats to drain into the pooling bag.
- Use the SSP or SSP+ solution to drain one bag after each other.
- Remove the set and then follow standard procedures to obtain a pooled buffy-coat derived platelet concentrate.

Apheresis

- Connected to the apheresis set, the SSP or SSP+ solution can be used as a substitution-solution for a part of the plasma in platelet concentrates.
- Connect by screwing the luer of the device to the apheresis kit. Check the integrity of the connection
- Break the breakable connector in the Additive Solution line by moving the connector back and forth 90° until an audible "snap" sound is heard in both directions (see drawing (A)).
- Bend tubing 180° to force the arrow piece of the frangible upward into the soft tubing closer to the bag opening and away from the base (B). This prevents this piece from resealing. To ensure proper flow, check to see if a gap exists between the arrow and base of the frangible (C). If not repeat step B.
- Hang the Additive Solution bag on the IV hook with the platelet product bag supported (not hanging below).
- Follow the valid apheresis protocol of your institution. Allow the appropriate amount of storage solution to flow into the platelet product bag and mix gently.

DISCARD WARNINGS

Discard according to the local regulations' on the potentially dangerous biological products.

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