

STRUCTURAL HEART THERAPY

AMPLATZER™ SEPTAL OCCLUDERS

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THE PROVEN STANDARD

FOR TRANSCATHETER ATRIAL SEPTAL DEFECT CLOSURE^{1, 2}

Amplatzer™ Septal Occluders are the standard of care for minimally invasive atrial septal defect (ASD) closure. These double-disc occluders are comprised of Nitinol mesh with polyester fabric. They are designed to securely appose the septal wall on each side of the defect and create a platform for tissue in-growth after implantation.

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WHEN SIMPLICITY MATTERS

The primary treatment option for closure of secundum atrial septal defects is with transcatheter devices.² This is when trust and simplicity matter. The Amplatzer Septal Occluder is the proven standard of care in transcatheter ASD closure.^{1,2}

SAFETY IN NUMBERS

98.5% Closure Rate: With no significant residual shunt (less than 2 mm) at 1-year follow-up¹

Low Major and Minor Complication Rates:

Postprocedure and long-term based on 5-year follow-up¹


SPECIFICALLY DESIGNED FOR ASD CLOSURE

Wide Waist: Centers device and fills the ASD³

Shape-memory Nitinol Mesh: Designed to securely appose both sides of the septal wall³

Polyester Material: Promotes occlusion and tissue in-growth³

Precise Placement: Device can be easily recaptured and redeployed³

The waist of the Amplatzer™ Septal Occluder fills the defect for optimal occlusion.

WHEN FLEXIBILITY MATTERS

Not all septal defects are the same. They come in a variety of shapes, sizes and often there can be multiple communications between the left and right atria. These types of atrial septal defects require a special device designed with the flexibility to meet multiple needs.

The addition of the Amplatzer Multifenestrated Septal Occluder - “Cribriform” to the Amplatzer family of occluders enables transcatheter closure for the majority of atrial septal defects.


SPECIFICALLY DESIGNED FOR MULTI-FENESTRATED ASD CLOSURE

Narrow Waist: Allows for placement through a central defect³

Matched Disc Diameters: Maximizes coverage of multiple fenestrations³

Polyester Material: Promotes occlusion and tissue in-growth³

Precise Placement: Device can be easily recaptured and redeployed³

The Amplatzer™ Cribriform Multi-Fenestrated Septal Occluder enables occlusion of the defect by covering the fenestrations with a single device.

Stocking both **Amplatzer™ Septal and Cribriform Occluders** provides you the ability to choose the appropriate device for a wide size-range of defects.

ORDERING INFORMATION

AMPLATZER™ SEPTAL OCCLUDER

Model / Reorder Number	Waist Diameter (mm)	Waist Width (mm)	Right Atrial Disc Diameter (mm)	Left Atrial Disc Diameter (mm)
9-ASD-004	4	3	12	16
9-ASD-005	5	3	13	17
9-ASD-006	6	3	14	18
9-ASD-007	7	3	15	19
9-ASD-008	8	3	16	20
9-ASD-009	9	3	17	21
9-ASD-010	10	3	18	22
9-ASD-011	11	4	21	25
9-ASD-012	12	4	22	26
9-ASD-013	13	4	23	27
9-ASD-014	14	4	24	28
9-ASD-015	15	4	25	29
9-ASD-016	16	4	26	30
9-ASD-017	17	4	27	31
9-ASD-018	18	4	28	32
9-ASD-019	19	4	29	33
9-ASD-020	20	4	30	34
9-ASD-022	22	4	32	36
9-ASD-024	24	4	34	38
9-ASD-026	26	4	36	40
9-ASD-028	28	4	38	42
9-ASD-030	30	4	40	44
9-ASD-032	32	4	42	46
9-ASD-034	34	4	44	50
9-ASD-036	36	4	46	52
9-ASD-038	38	4	48	54
9-ASD-040	40	4	50	56

AMPLATZER™ CRIBRIFORM MULTI-FENESTRATED SEPTAL OCCLUDER

Model / Reorder Number	Right & Left Atrium Disc Diameter (mm)	Waist Width (mm)
9-ASD-MF-018	18	3
9-ASD-MF-025	25	3
9-ASD-MF-030	30	3
9-ASD-MF-035	35	3
9-ASD-MF-040	40	3

REFERENCE:

1. Amplatzer Septal Occluder Instructions for Use.
2. Kashour TS, Latroche B, Elhoury ME, et al. Successful Percutaneous Closure of a Secundum Atrial Septal Defect through Femoral Approach in a Patient with Interrupted Inferior Vena Cava. Congenital Heart Disease. 2010;5(6):620-623.
3. Test(s) performed by and data on file at Abbott.

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Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

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© 2020 Abbott. All Rights Reserved. MAT-2004327 v1.0 | Item approved for Global OUS use only.

AMPLATZER™

Septal Occluder

Instructions for Use

bg: Септално устройство за затваряне Инструкции за употреба	hr: Septalni okluder Upute za uporabu	ro: Dispozitiv de închidere de sept Instrucțiuni pentru utilizare
bs: Septalni okluder Uputstvo za upotrebu	hu: Szeptumtapasz Használati utasítás	ru: окклюдер для перегородок Инструкции по применению
cs: Septální okludér Návod k použití	it: Occlusore settale Istruzioni per l'uso	sk: Septálny okluzor ventrikulárnych septálnych defektov Návod na použitie
da: Septal okklusionsanordning Brugsanvisning	ko: 중격 폐색기 사용 설명서	sl: Septalno zapiralo Navodila za uporabo
de: Septum-Okluder Gebrauchsanweisung	lt: Pertvaros okluderis Naudojimo instrukcija	sq: Okludues septal Udhëzime për përdorim
el: Συσκευή απόφραξης διαφράγματος Οδηγίες χρήσης	lv: Starpsienas Nosprostotājs Lietošanas instrukcija	sr: Septalni okluder Uputstvo za upotrebu
es: Ocluser septal Instrucciones de uso	nl: Septale occluder Gebruiksaanwijzing	sv: Anordning för slutning av septum Bruksanvisning
et: Vaheseina sulgur Kasutusjuhend	no: Septal lukkeanordning Bruksanvisning	tr: Septal Oklüder Kullanma Talimatı
fi: Väliseinäaukon sulkija Käyttöohjeet	pl: Korek do zamykania ubytków przegrody Instrukcja użytkowania	zh: 房间隔缺损封堵器 使用说明
fr: Dispositif d'occlusion septal Mode d'emploi	pt: Dispositivo de oclusão septal Instruções de Utilização	ar: سداة الحاجز تعليمات الاستخدام

STERILE EO


CE 0086

R ONLY

LATEX
Does not contain
natural rubber
latex components

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MORE FLEXIBILITY,
MORE CONTROL

DISCOVER THE UNIQUE
AMPLATZER™ TREVISIO™

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

RELIABLE PRECISION WHEN IT MATTERS MOST


The Amplatzer™ Trevisio™ Intravascular Delivery System is an ultra-flexible delivery system enabling interventional cardiologists to perform their work with complete confidence. It leverages the one-piece cable design utilized by the Amplatzer™ TorqVue™ Delivery System, also known as the Classic Amplatzer™ Delivery System*. Trevisio is designed for no compromises on torque strength, sheath diameter and pushability.

ULTRA-FLEXIBLE TIP

- Improves assessment of device position before cable release.
- Reduces bias on the device.

FLEXIBLE TRANSITION SECTION

Maintains sheath position during deployment of the device.

CLASSIC AMPLATZER™ DELIVERY SYSTEM*


AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM


ULTRA-
FLEXIBLE TIP

FLEXIBLE
TRANSITION
SECTION

STIFF
PROXIMAL
SECTION

Flexible tip reduces the bias on the device and improves the assessment of device position prior to cable release.

STIFF PROXIMAL SECTION

Maintains pushability of the
delivery system.

DELIVERY SYSTEM DIMENSIONS

DELIVERY SYSTEM (SHEATH SIZE)	INNER DIAMETER OF SHEATH	OUTER DIAMETER OF SHEATH	MODEL NUMBER/DELIVERY SYSTEM SIZE (MM)
6 Fr	2.11 mm (0.08 in)	2.79 mm (0.11 in)	9-ATV06F45/60
7 Fr	2.44 mm (0.10 in)	3.18 mm (0.13 in)	9-ATV07F45/60
7 Fr	2.44 mm (0.10 in)	3.18 mm (0.13 in)	9-ATV07F45/80
8 Fr	2.69 mm (0.11 in)	3.45 mm (0.14 in)	9-ATV08F45/60
8 Fr	2.69 mm (0.11 in)	3.45 mm (0.14 in)	9-ATV08F45/80
9 Fr	3.00 mm (0.12 in)	3.81 mm (0.15 in)	9-ATV09F45/80
10 Fr	3.30 mm (0.13 in)	4.14 mm (0.16 in)	9-ATV10F45/80
12 Fr	3.99 mm (0.16 in)	4.80 mm (0.19 in)	9-ATV12F45/80
13 Fr	4.32 mm (0.17 in)	5.13 mm (0.20 in)	9-ATV13F45/80

*Amplatzer™ TorqVue™ Delivery System

DOWNLOAD NOW

Download the free Amplatzer™ Portfolio App from your favorite app store.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Abbott

Park Lake, Culliganlaan 2b, 1831 Diegem, Belgium
www.cardiovascular.abbott

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2020 Abbott. All Rights Reserved. MAT-2001046 v1.0 | Item approved for OUS use only.

AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM


DEVICE DESCRIPTION

The Amplatzer™ Trevisio™ Intravascular Delivery System is an ultra-flexible delivery system enabling interventional cardiologists to perform their work with complete confidence. It leverages the one-piece cable design utilized by the Amplatzer™ TorqVue™ Delivery System, also known as the Classic Amplatzer™ Delivery System. Trevisio is designed for no compromises on torque strength, sheath diameter and pushability.

The Amplatzer™ Trevisio™ Intravascular Delivery System comes with a 45° sheath curve in 60 or 80 cm lengths from 6Fr to 13Fr.

INDICATIONS AND USAGE


The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices. The Amplatzer™ Trevisio™ Intravascular Delivery System is indicated for use by interventional cardiologists to place an Amplatzer™ Occluder device.

DEVICE SPECIFICATIONS

Model Number/ Delivery System Size mm	Delivery System (Sheath Size) Fr	Inner Diameter of Sheath mm (inch)	Outer Diameter of Sheath mm (inch)	Length cm
9-ATV06F45/60	6	2.11 (0.08)	2.79 (0.11)	60
9-ATV07F45/60	7	2.44 (0.10)	3.18 (0.13)	60
9-ATV07F45/80	7	2.44 (0.10)	3.18 (0.13)	80
9-ATV08F45/60	8	2.69 (0.11)	3.45 (0.14)	60
9-ATV08F45/80	8	2.69 (0.11)	3.45 (0.14)	80
9-ATV09F45/80	9	3.00 (0.12)	3.81 (0.15)	80
9-ATV10F45/80	10	3.30 (0.13)	4.14 (0.16)	80
9-ATV12F45/80	12	3.99 (0.16)	4.80 (0.19)	80
9-ATV13F45/80	13	4.32 (0.17)	5.13 (0.20)	80

PACKAGING INFORMATION

- A. Loader – Introduces an Amplatzer™ device into the sheath
- B. Hemostasis valve with extension tube and stopcock – Allows flushing of the delivery system and controls blood backflow
- C. Sheath – Provides a pathway through which an Amplatzer™ device is delivered
- D. Dilator – Eases penetration of tissue and minimizes vessel trauma
- E. Delivery cable – Attaches to the device to control its movement through the sheath
- F. Plastic vise – Attaches to the delivery cable and serves as a handle for disconnecting (unscrewing) the delivery cable from a device

Note that the device and the delivery system are being sold separately.

ASSOCIATED AMPLATZER™ PRODUCTS

Amplatzer™ Trevisio™ Intravascular Delivery System Sizes							
	6 Fr	7 Fr	8 Fr	9 Fr	10 Fr	12 Fr	13 Fr
Amplatzer™ Septal (ASD) Occluder	9-ASD-004	9-ASD-011	9-ASD-018	9-ASD-020	9-ASD-026	9-ASD-032	9-ASD-032
	9-ASD-005	9-ASD-012	9-ASD-019	9-ASD-022	9-ASD-028	9-ASD-034	9-ASD-034
	9-ASD-006	9-ASD-013	–	9-ASD-024	9-ASD-030	9-ASD-036	9-ASD-036
	9-ASD-007	9-ASD-014	–	–	–	9-ASD-038	9-ASD-038
	9-ASD-008	9-ASD-015	–	–	–	9-ASD-040	9-ASD-040
	9-ASD-009	9-ASD-016	–	–	–	–	–
	9-ASD-010	9-ASD-017	–	–	–	–	–
Amplatzer™ Multi-fenestrated ASD (Cribriform) Occluder	–	–	9-ASD-MF-018	–	–	–	–
	–	–	9-ASD-MF-025	9-ASD-MF-035	9-ASD-MF-040	–	–
	–	–	9-ASD-MF-030	–	–	–	–
Amplatzer™ PFO Occluder	–	–	9-PFO-018	–	–	–	–
	–	–	9-PFO-025	9-PFO-035	–	–	–
	–	–	9-PFO-030	–	–	–	–
Amplatzer™ Muscular VSD Occluder	9-VSD-MUSC-004	–	–	–	–	–	–
	9-VSD-MUSC-006	9-VSD-MUSC-012	9-VSD-MUSC-014	9-VSD-MUSC-018	–	–	–
	9-VSD-MUSC-008	–	9-VSD-MUSC-016	–	–	–	–
	9-VSD-MUSC-010	–	–	–	–	–	–
Amplatzer™ Muscular PI VSD Occluder	–	–	–	9-VSDMUSCPI-016	9-VSDMUSCPI-020	–	–
	–	–	–	9-VSDMUSCPI-018	9-VSDMUSCPI-022	–	–
	–	–	–	–	9-VSDMUSCPI-024	–	–

LATEX-FREE INFORMATION

These Amplatzer™ products do not contain latex.

Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, B-1831 Diegem, Belgium, Tel: +32 2 714 14 11

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott. Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

For more information, visit our web site at www.abbott.com

© 2020 Abbott. All rights reserved. 9-EH-11038-01 03-2020

AMPLATZER™

SIZING BALLOON II

DEFECT VISUALIZATION FOR ACCURATE DEVICE SELECTION

Featuring a radiopaque marker system for better visualization and positioning, the Amplatzer™ Sizing Balloon II provides a simple, precise method for visualization of cardiovascular defects in patients undergoing a transcatheter implant procedure. It features compliant balloon material to accurately visualize cardiovascular defects of various shapes and tissue properties.

DESIGNED FOR ACCURATE DEFECT MEASUREMENT

Radiopaque Marker Bands

- Are strategically placed to assist with alignment of imaging and support proper sizing calibration (Figure 1)

Soft Distal Tip and Flexible Shaft

- Offer smooth access for atraumatic balloon positioning

Compliant Balloon Material

- Allows precise visualization of a variety of defect shapes and tissue properties when using a stop-flow technique
- Supports low-pressure inflation

The balloon shall be available in 3 sizes, 18mm, 24mm and 34mm.

- A minimum of 20mm when inflated with a max inflation volume of 15cc of saline
- A minimum of 25mm when inflated with a max inflation volume of 30cc of saline
- A minimum of 40mm when inflated with a max inflation volume of 90 cc of saline


Figure 1: The most proximal marker band is for sizing calibration. The two central marker bands provide visual alignment for catheter placement in the cardiovascular defect.


ORDERING INFORMATION

Amplatzer™ Sizing Balloon II Specifications

Model/Reorder Number	Nominal Dimension (mm)	Nominal Inflation Volume (cc)	Max Defect Size ^a (mm)	Useable Length (cm)	Guidewire (inch)	Max Balloon Inflation Volume (cc)
9-SB-018	18	10	20	70	0.035	15
9-SB-024	24	25	25	70	0.035	30
9-SB-034	34	55	40	70	0.035	90

a. Maximum defect size refers to diameter of native defect as measured by echocardiography.
Note: Do not insert balloon catheter through an introducer sheath.

Figure 2

SUPPLEMENTAL INFORMATION

1: Catheter distal tip – Eases penetration of tissue and minimizes vessel trauma.

2/3: Catheter – The catheter shaft is composed of a tri-lumen tube. Two lumens provide a pathway for inflating and deflating the balloon. There are several inflation holes (3) along the distal end of the shaft to facilitate uniform inflation and deflation of the balloon with inflation medium. The third lumen accepts a guidewire which facilitates advancement of the balloon through the vasculature and positioning of the balloon within the defect.

4: Balloon – The balloon is filled with contrast medium within the defect and visualized via imaging for defect assessment.

5: Marker bands – Provide visual markers for orienting imaging equipment.

6/7/8: Bifurcation hub – Molded component (6) that provides a pathway for a guidewire lumen (7) and a balloon inflation/deflation lumen (8) for balloon inflation/deflation.

9/10/11: Inflation/deflation sidearm (9) – Attaches to the bifurcation hub balloon inflation/deflation lumen. A stopcock (10) on the inflation/deflation sidearm allows for control of balloon inflation volume. A balloon inflation/deflation port (11) allows for a luer lock syringe to be attached for inflation/deflation of the balloon.

12: Protective balloon cover – protective cover to prevent balloon damage during handling prior to the procedure.

NOTE: Figure 2 shows the balloon both uninflated (inside the protective balloon cover) and inflated.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.
Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.
Photo(s) on file at Abbott.

Abbott

3200 Lakeside Dr., Santa Clara, CA. 95054 USA

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

www.structuralheart.abbott

© 2024 Abbott. All Rights Reserved. MAT-2314177 v1.0 | Item approved for OUS use only. (MDR Version)

AMPLATZER™ GUIDEWIRES

Ordering Information

Contents: 1 guidewire

Model/ Reorder No	Diam (in)	Body	Floppy Tip Lgth (cms)	Tip Description	Usable Lgth (cm)
9-GW-001	0.035	Super Stiff	6	7.5 mm, Modified J-tip	260
9-GW-002	0.035	Super Stiff	5	1.5 mm, Modified J-tip	260
9-GW-003	0.035	Super Stiff	20	6 mm, J-tip	300
9-GW-004*	0.035	Soft Tip, Fixed Core	NA	6 mm, J-tip	300

*9-GW-004, also referred to as “Noodlewire,” is a soft tip, flexible guidewire recommended for establishing an arterial-venous loop, facilitating closure of ventricular septal defects.

EU MDR Declaration of Conformity
Amplatzer™ Septal Occluder and
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”

This declaration of conformity is issued under the sole responsibility of the manufacturer.

I, the undersigned, for and on behalf of Abbott Medical, hereby declare that the medical device(s) specified above conform(s) with the applicable *General Safety & Performance Requirements* listed in Annex I and all relevant provisions of Regulation (EU) 2017/745.

Directive 2006/42/EC on Machinery and Directive 89/686/EEC (and the superseding Regulation (EU) 2016/425) on Personal Protective Equipment do not apply.

Common Specifications Utilized:

No applicable common specifications

Notified Body:

BSI Group The Netherlands B.V.
NB #: 2797
Say Building
John M. Keynesplein 9
1066 EP Amsterdam
The Netherlands

Conformity Assessment Procedure:

Conformity assessment for the Amplatzer™ Septal Occluder and Amplatzer™ Multifenestrated Septal Occluder – “Cribiform” is based on a quality management system and assessment of the technical documentation as per Annex IX.

Supporting Certificate(s):

- EU Quality Management System Certificate: MDR 750915
- Technical Documentation Assessment Certificate Number: MDR 777149

EU MDR Declaration of Conformity
Amplatzer™ Septal Occluder and
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”

This signature section is applicable to all declarations of conformity, including to other European Union legislation, if applicable:

Prepared By: Nicolette Pedersen Date: 05-FEB-2024
Nicolette Pedersen, Regulatory Affairs Specialist II

Authorized Signatory: Christopher Gallivan Date: 05-Feb-2024
Christopher Gallivan, DVP, Quality Person Responsible for Regulatory Compliance (PRRC)

Place of Issue: Santa Clara, California, USA Issue Date: 05-Feb-2024

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity
Amplatzer™ Septal Occluder and
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”

Declaration of Conformity Product List

Model Numbers	Product Trade Name	Basic UDI-DI	Original MDR CE Mark Date
9-ASD-004	Amplatzer™ Septal Occluder	5415067AMP1100DV	07-Dec-2023
9-ASD-005			
9-ASD-006			
9-ASD-007			
9-ASD-008			
9-ASD-009			
9-ASD-010			
9-ASD-011			
9-ASD-012			
9-ASD-013			
9-ASD-014			
9-ASD-015			
9-ASD-016			
9-ASD-017			
9-ASD-018			
9-ASD-019			
9-ASD-020			
9-ASD-022			
9-ASD-024			
9-ASD-026			
9-ASD-028			
9-ASD-030			
9-ASD-032			
9-ASD-034			
9-ASD-036			
9-ASD-038			
9-ASD-040			

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity
Amplatzer™ Septal Occluder and
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”

Model Numbers	Product Trade Name	Basic UDI-DI	Original MDR CE Mark Date
9-ASD-MF-018	Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”	5415067AMP1100DV	15-Nov-2023
9-ASD-MF-025			
9-ASD-MF-030			
9-ASD-MF-035			
9-ASD-MF-040			

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity

Amplatzer™ Delivery Systems

Manufacturer's Name:	Abbott Medical
Manufacturer's Single Registration Number (SRN):	US-MF-000018613
Manufacturer's Address:	177 County Road B East St. Paul, MN 55117 USA
Authorized Representative's Name, Address, and Single Registration Number (SRN):	Abbott Medical The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium SRN: BE-AR-000008744
Product Trade Name(s):	• Amplatzer™ TorqVue™ Delivery System • Amplatzer™ TorqVue™ Exchange System • Amplatzer™ Trevisio™ Intravascular Delivery System • Amplatzer™ TorqVue™ 2 Delivery Sheath • Amplatzer™ Amulet™ Delivery Sheath • Amplatzer™ TorqVue™ LP Delivery System
Model Number(s):	See attached Product List
Intended Purpose:	See attached Product List
Risk Classification and Rule:	See attached Product List
EMDN Code(s):	See attached Product List
GMDN Code(s):	45419 – Cardiac Occluder Delivery Kit
Basic UDI-DI:	See attached Product List

EU MDR Declaration of Conformity
Amplatzer™ Delivery Systems

This declaration of conformity is issued under the sole responsibility of the manufacturer.

I, the undersigned, for and on behalf of Abbott Medical, hereby declare that the medical device(s) specified above conform(s) with the applicable *General Safety & Performance Requirements* listed in Annex I and all relevant provisions of Regulation (EU) 2017/745.

Directive 2006/42/EC on Machinery and Directive 89/686/EEC (and the superseding Regulation (EU) 2016/425) on Personal Protective Equipment do not apply.

Common Specifications Utilized:

- No applicable common specifications

Notified Body:

BSI Group The Netherlands B.V.
NB #: 2797
Say Building
John M. Keynesplein 9
1066 EP Amsterdam
The Netherlands

Conformity Assessment Procedure:

Conformity assessment for the Amplatzer™ Delivery Systems is based on a quality management system and assessment of the technical documentation as per Annex IX.

Supporting Certificate(s):

- EU Quality Management System Certificate: MDR 750915
- Technical Documentation Assessment Certificate Number: MDR 750953

This signature section is applicable to all declarations of conformity, including to other European Union legislation, if applicable:

Prepared By: Ashley Starr Date: 10 July 2024
Ashley Starr, Associate Director, Regulatory Affairs

Authorized Signatory: C. Gallivan Date: 10 July 2024
Christopher Gallivan, DVP, Quality Person Responsible for Regulatory Compliance (PRRC)

Place of Issue: St. Paul, Minnesota, USA

Issue Date: 10 July 2024

EU MDR Declaration of Conformity

Amplatzer™ Delivery Systems

Declaration of Conformity Product List

Model Numbers	Product Name	Original MDR CE Mark Date	Intended Purpose	Risk Classification and Rule	EMDN Code	Basic UDI-DI
9-ITV06F45/60	Amplatzer™ TorqVue™ Delivery System	21-Feb-2023	The Amplatzer™ TorqVue™ Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of occluders.	Class III, Rule 7, MDR per Annex VIII	C0502 - Cardiovascular Introducing Sheaths, Valved	5415067AMP1500EH
9-ITV07F45/60						
9-ITV07F45/80						
9-ITV08F45/60						
9-ITV08F45/80						
9-ITV09F45/80						
9-ITV10F45/80						
9-ITV12F45/80						
9-ITV13F45/80						
9-ITV05F180/60						
9-ITV06F180/60						
9-ITV06F180/80						
9-ITV07F180/80						
9-ITV08F180/80						
9-ITV09F180/80						
9-EITV09F45/80	Amplatzer™ TorqVue™ Exchange System	21-Feb-2023	The Amplatzer™ TorqVue™ Exchange System is intended to facilitate the removal and exchange of a delivery sheath with one of equal or larger diameter to enable the recapture and replacement of an occluder.	Class III, Rule 7, MDR per Annex VIII	C0502 - Cardiovascular Introducing Sheaths, Valved	5415067AMP1500EH
9-EITV12F45/80						
9-EITV06F180/80						
9-EITV08F180/80						

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity

Amplatzer™ Delivery Systems

Model Numbers	Product Name	Original MDR CE Mark Date	Intended Purpose	Risk Classification and Rule	EMDN Code	Basic UDI-DI
9-ATV06F45/60	Amplatzer™ Trevisio™ Intravascular Delivery System	21-Feb-2023	The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of occluders.	Class III, Rule 7, MDR per Annex VIII	C0502 - Cardiovascular Introducing Sheaths, Valved	5415067AMP1500EH
9-ATV07F45/60						
9-ATV07F45/80						
9-ATV08F45/60						
9-ATV08F45/80						
9-ATV09F45/80						
9-ATV10F45/80						
9-ATV12F45/80						
9-ATV13F45/80						
9-TV2-05F120	Amplatzer™ TorqVue™ 2 Delivery Sheath	20-Mar-2023	The Amplatzer™ TorqVue™ 2 Delivery Sheath is intended to provide a pathway through which a device is introduced within the chambers and coronary vasculature of the heart.	Class III, Rule 7, MDR per Annex VIII	C0501 - Cardiovascular Introducing Sheaths, Non-Valved	5415067AMP1500EH
9-TV2-06F120						
9-TV2-07F120						

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity

Amplatzer™ Delivery Systems

Model Numbers	Product Name	Original MDR CE Mark Date	Intended Purpose	Risk Classification and Rule	EMDN Code	Basic UDI-DI
DS-TV45X45-12F-080	Amplatzer™ Amulet™ Delivery Sheath	2023-May-17	The Amplatzer™ Amulet™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.	Class III, Rule 6, MDR per Annex VIII	C0501 - Cardiovascular Introducing Sheaths, Non-Valved	5415067AMP1001DS
DS-TV45X45-14F-080						
9-TVLP4F90/060	Amplatzer™ TorqVue™ LP Delivery System	2023-Dec-07	The Amplatzer™ TorqVue™ LP Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of occluders.	Class III, Rule 7, MDR per Annex VIII	C0501 - Cardiovascular Introducing Sheaths, Non-Valved	5415067AMP1401EE
9-TVLP4F90/080						
9-TVLP5F90/060						
9-TVLP5F90/080						

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity
Amplatzer™ Septal Occluder and
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”

Manufacturer's Name:	Abbott Medical
Manufacturer's Single Registration Number (SRN):	US-MF-000018613
Manufacturer's Address:	177 County Road B East St. Paul, MN 55117 USA
Authorized Representative's Name, Address, and Single Registration Number (SRN):	Abbott Medical The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium SRN: BE-AR-000008744
Product Trade Name(s):	See attached product list
Model Number(s):	See attached product list
Intended Purpose:	The Amplatzer™ Septal Occluder is a device intended for the nonsurgical occlusion of ostium secundum atrial septal defects or fenestrations. The Amplatzer™ Multifenestrated Septal Occluder – “Cribiform” is a device intended for the occlusion of a multifenestrated atrial septal defect.
Risk Classification and Rule:	Class III, Rule 8, MDR per Annex VIII
EMDN Code(s):	P07040302
GMDN Code(s):	45418
Basic UDI-DI:	See attached product list

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity
Amplatzer™ Sizing Balloon II

Manufacturer's Name:	Abbott Medical
Manufacturer's Single Registration Number (SRN):	US-MF-000018613
Manufacturer's Address:	177 County Road B East St. Paul, MN 55117 USA
Authorized Representative's Name, Address, and Single Registration Number (SRN):	Abbott Medical The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium SRN: BE-AR-000008744
Product Trade Name:	Amplatzer™ Sizing Balloon II
Model Number:	9-SB-018 9-SB-024 9-SB-034
Intended Purpose:	The Amplatzer™ Sizing Balloon II is intended to facilitate sizing of a cardiovascular defect.
Risk Classification and Rule:	Class III, Rule 7, MDR per Annex VIII
EMDN Code:	C0104020103
GMDN Code:	42417
Basic UDI-DI:	5415067AMP1101DX

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU MDR Declaration of Conformity

Amplatzer™ Sizing Balloon II

This declaration of conformity is issued under the sole responsibility of the manufacturer.

I, the undersigned, for and on behalf of Abbott Medical, hereby declare that the medical device(s) specified above conform(s) with the applicable *General Safety & Performance Requirements* listed in Annex I and all relevant provisions of Regulation (EU) 2017/745.

Directive 2006/42/EC on Machinery and Directive 89/686/EEC (and the superseding Regulation (EU) 2016/425) on Personal Protective Equipment do not apply.

Common Specifications Utilized:

No applicable common specifications

Notified Body:

BSI Group The Netherlands B.V.
NB #: 2797
Say Building
John M. Keynesplein 9
1066 EP Amsterdam
The Netherlands

Conformity Assessment Procedure:

Conformity assessment for the Amplatzer™ Sizing Balloon II is based on a quality management system and assessment of the technical documentation as per Annex IX.

Supporting Certificate(s):

- EU Quality Management System Certificate: MDR 750915
- Technical Documentation Assessment Certificate Number: MDR 777443

This signature section is applicable to all declarations of conformity, including to other European Union legislation, if applicable:

Prepared By: Nicolette Pedersen Date: 28-May-2024
Nicolette Pedersen, Regulatory Affairs Specialist II

Authorized Signatory: Christopher Gallivan Date: 30-May-2024
Christopher Gallivan, DVP, Quality Person Responsible for Regulatory Compliance (PRRC)

Place of Issue: Santa Clara, California, USA

Issue Date: 30-May-2024

This confidential document is the property of Abbott and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777149 R000

Manufacturer: Abbott Medical

Address:

177 County Road B East
St. Paul
Minnesota
55117
USA

Single Registration Number: US-MF-000018613

EU Authorised Representative: Abbott Medical


Address:

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium

Scope: See attached **Device Schedule**

On the basis of our assessment of the technical documentation in accordance with Regulation (EU) 2017/745, Annex IX Chapter II, the technical documentation meets the requirements of the Regulation. For the placing on the market of these devices an additional Annex IX Chapter I and III certificate is required.

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Graeme Tunbridge, Senior Vice President Medical Devices

First Issue Date: **2023-11-15**

Current Issue Date: **2023-12-07**

Starting Validity Date: **2023-12-07**

Expiry Date: **2028-11-14**

...making excellence a habit.™

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777149 R000

Device Schedule:

Intended Purpose as per the Instructions for Use:

The Amplatzer™ Multifenestrated Septal Occluder – “Cribiform” is a device intended for the occlusion of a multifenestrated atrial septal defect.

Type (Codes as per (EU) 2017/2185): MDN 1101

Risk Classification: Class III Implantable

Basic UDI-DI: 5415067AMP1100DV

Device Name	Model	Left Atrial Disc Diameter mm	Right Atrial Disc Diameter mm	Waist Length mm
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”	9-ASD-MF-018	18	18	3
	9-ASD-MF-025	25	25	3
	9-ASD-MF-030	30	30	3
	9-ASD-MF-035	35	35	3
	9-ASD-MF-040	40	40	3

First Issue Date: **2023-11-15**

Current Issue Date: **2023-12-07**

Starting Validity Date: **2023-12-07**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 2 of 5

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777149 R000

Intended Purpose as per the Instructions for Use:

The Amplatzer™ Septal Occluder is a device intended for the nonsurgical occlusion of ostium secundum atrial septal defects or fenestrations.

Type (Codes as per (EU) 2017/2185): MDN 1101

Risk Classification: Class III Implantable

Basic UDI-DI: 5415067AMP1100DV

Device Name	Model	Waist Diameter mm	Left Atrial Disc Diameter mm	Right Atrial Disc Diameter mm	Waist Length mm
Amplatzer™ Septal Occluder	9-ASD-004	4	16	12	3
	9-ASD-005	5	17	13	3
	9-ASD-006	6	18	14	3
	9-ASD-007	7	19	15	3
	9-ASD-008	8	20	16	3
	9-ASD-009	9	21	17	3
	9-ASD-010	10	22	18	3
	9-ASD-011	11	25	21	4
	9-ASD-012	12	26	22	4
	9-ASD-013	13	27	23	4
	9-ASD-014	14	28	24	4
	9-ASD-015	15	29	25	4
	9-ASD-016	16	30	26	4
	9-ASD-017	17	31	27	4
	9-ASD-018	18	32	28	4
	9-ASD-019	19	33	29	4
	9-ASD-020	20	34	30	4

First Issue Date: **2023-11-15**

Current Issue Date: **2023-12-07**

Starting Validity Date: **2023-12-07**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 3 of 5

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777149 R000

Device Name	Model	Waist Diameter mm	Left Atrial Disc Diameter mm	Right Atrial Disc Diameter mm	Waist Length mm
Amplatzer™ Septal Occluder	9-ASD-022	22	36	32	4
	9-ASD-024	24	38	34	4
	9-ASD-026	26	40	36	4
	9-ASD-028	28	42	38	4
	9-ASD-030	30	44	40	4
	9-ASD-032	32	46	42	4
	9-ASD-034	34	50	44	4
	9-ASD-036	36	52	46	4
	9-ASD-038	38	54	48	4
	9-ASD-040	40	56	50	4

First Issue Date: **2023-11-15**

Current Issue Date: **2023-12-07**

Starting Validity Date: **2023-12-07**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 4 of 5

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777149 R000

Certificate History

(References to applicable Common Specifications, Harmonized Standards complied with, and the relevant test and audit reports that support any of the below certificate changes may be requested from Certificate.Verification@bsigroup.com)

Date	Reference Number	Action
2023-11-15	3746939	Issued
Current	3746937	Supplemented – Addition of Amplatzer™ Septal Occluder devices.

First Issue Date: **2023-11-15**

Current Issue Date: **2023-12-07**

Starting Validity Date: **2023-12-07**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 5 of 5

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Manufacturer: Abbott Medical

Address:


177 County Road B East
St. Paul
Minnesota
55117
USA

Single Registration Number: US-MF-000018613

Scope: See attached **Device Schedule**

On the basis of our examination of the quality system in accordance with Regulation (EU) 2017/745, Annex IX Chapter I and III, the quality system meets the requirements of the Regulation. For the placing on the market of Class III devices, and Class IIb implantable devices that are not considered well-established technologies as specified in Article 52(4) an additional Annex IX Chapter II certificate is required.

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Graeme Tunbridge, Senior Vice President Global Regulatory & Quality

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

EU Authorised Representative: Abbott Medical

Address:

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium

EU Authorised Representative: Abbott Vascular International BVBA

Address:

Park Lane
Culliganlaan, 2B
1831 Diegem
Belgium

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 2 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Device Schedule: Class III and Class IIb devices

Class III, Implantable	Intended purpose
Amplatzer™ Talisman™ PFO Occlusion System	See MDR 751010
Navitor™ Transcatheter Aortic Valve	See MDR 751017
Navitor Titan™ Transcatheter Aortic Valve	
Amplatzer™ Amulet™ Left Atrial Appendage Occluder	See MDR 751008
Portico™ Transcatheter Aortic Valve	See MDR 751019
MitraClip G4 System	See MDR 751009
Amplatzer™ Duct Occluder	See MDR 777151
Amplatzer™ Duct Occluder II	See MDR 777152
Amplatzer Piccolo™ Occluder	
Masters Series HP™ Valved Graft with Gelweave Valsalva™ technology (VAVGJ)	See MDR 772712
Masters Series™ Aortic Valved Graft with Hemashield™ Graft Technology (CAVGJ)	See MDR 772713
Tendyne Transcatheter Mitral Valve System	See MDR 772768
TriClip™ G4 System	See MDR 777146
Master Series™ Mechanical Heart Valves	See MDR 775965
Amplatzer™ Multifenestrated Septal Occluder – “Cribiform”	See MDR 777149
Amplatzer™ Septal Occluder	
Epic™ Max Stented Tissue Valves	See MDR 781775
Epic™ Plus Stented Tissue Valves	See MDR 773424
Séguin™ Annuloplasty Ring	See MDR 777154
SJM Tailor™ Annuloplasty Ring	
SJM Tailor™ Annuloplasty Band	
Rigid Saddle Ring	
Amplatzer™ Valvular Plug III	See MDR 777143
Amplatzer™ Muscular VSD Occluder	See MDR 777441
Amplatzer™ P.I. Muscular VSD Occluder	
Regent™ Mechanical Heart Valve	See MDR 775967
Regent™ Mechanical Heart Valve with FlexCuff™	

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 3 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Class III	Intended purpose
FlexNav™ Delivery System	See MDR 751005
Amplatzer™ Amulet™ Delivery Sheath	See MDR 750953
Amplatzer™ Trevisio™ Intravascular Delivery System	
Amplatzer™ TorqVue™ Delivery System	
Amplatzer™ TorqVue™ Exchange System	
Amplatzer™ TorqVue™ 2 Delivery Sheath	
Amplatzer™ TorqVue™ LP Catheter and Delivery System	
Amulet™ Steerable Delivery Sheath	
Amplatzer™ Sizing Balloon II	See MDR 777443
Epic™ Max Heart Valve Sizer Set and Holder Handle	See MDR 781775
Epic™ Plus Heart Valve Sizer Set and Holder Handles	See MDR 773424
Mechanical Heart Valve Accessories – Holder Handles, Leaflet Tester, Replacement Holder/Rotators and Sizer Set	See MDR 775966
Séguin™ Annuloplasty Ring Sizer Set, Holder Handle and Extension Handle	See MDR 777154
SJM Tailor™ Annuloplasty Ring Sizer Set	
Rigid Saddle Ring Sizer Set	
Regent™ Mechanical Heart Valve Accessories - Replacement Holder/Rotators and Sizer Set	See MDR 775967
Class IIb, Implantable	Intended purpose
Amplatzer™ Vascular Plug, Amplatzer™ Vascular Plug II, Amplatzer™ Vascular Plug 4	See MDR 767903

Device Schedule: Class IIa, Custom-made and other devices

Device(s)	Risk Classification
Transcatheter Aortic Heart Valve Loading System	Class Is
Tendyne Loading System	Class Is
Tendyne Stand Components	Class Is
For Class Is devices, the Notified Body conformity assessment is limited to the aspects relating to establishing, securing and maintaining sterile conditions.	

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 4 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Certificate History

(References to applicable Common Specifications, Harmonized Standards complied with, and the relevant test and audit reports that support any of the below certificate changes may be requested from Certificate.Verification@bsigroup.com)

Date	Reference Number	Action
2022-04-21	3447260	Issued
2022-08-16	3682747	Amended – Addition of subcontractors. Addition of service for subcontractor: Microbiology service. Removal of services for subcontractor: Labelling and packaging. Supplemented - Addition of Navitor™ Transcatheter Aortic Valve, FlexNav™ Delivery System, Amplatzer™ Amulet™ Left Atrial Appendage Occluder, and Transcatheter Heart Valve Loading System.
2022-09-29	3677048	Amended – Addition of subcontractor.
2022-11-11	3795938	Supplemented – Addition of Portico™ Transcatheter Aortic Valve. Amended – Administrative update to previous entries Reference Number 3682747 and Reference Number 3677048.
2022-11-29	3766365	Supplemented – Addition of MitraClip G4 System, Amplatzer™ Steerable Delivery Sheath, and Amplatzer™ Vascular Plug, Amplatzer™ Vascular Plug II, Amplatzer™ Vascular Plug 4. Amended – Addition of subcontractor.
2023-02-21	3854349	Supplemented – Addition of Amplatzer™ TorqVue™ Delivery System, Amplatzer™ TorqVue™ Exchange System, Amplatzer™ Trevisio™ Intravascular Delivery System
2023-03-20	3873040	Supplemented – Addition of Amplatzer™ TorqVue™ 2 Delivery Sheath.
2023-05-17	3894086	Supplemented – Addition of Amplatzer™ Amulet™ Delivery Sheath, Amplatzer™ Duct Occluder, Amplatzer™ Duct Occluder II, and Amplatzer™ Piccolo™ Occluder.
2023-06-14	3869869	Supplemented – Addition of Masters HP™ Valved Graft with Gelweave Valsalva™ technology (VAVGJ) and Masters Series™ Aortic Valved Graft with Hemashield™ Technology (CAVGJ).

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 5 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Date	Reference Number	Action
2023-08-25	30002038	Amended – Minor update to clarify the device name from “Transcatheter Heart Valve Loading System” to “Transcatheter Aortic Heart Valve Loading System”. Supplemented – Addition of TriClip™ G4 System, Tendyne Transcatheter Mitral Valve System, Master Series™ Mechanical Heart Valves, Tendyne Loading System and Tendyne Stand Components.
2023-11-15	30007555	Supplemented – Addition of Amplatzer™ Multifenestrated Septal Occluder – “Cribiform” and Amplatzer™ Sizing Balloon II to device schedule. Amended – Administrative correction of device names in device schedule for Masters Series HP™ Valved Graft with Gelweave Valsalva™ Technology (VAVGJ) and Masters Series™ Aortic Valved Graft with Hemashield™ Graft Technology (CAVGJ)
2023-12-07	30057065	Amended – Minor formatting update to consolidate devices with same MDR TDA certificate number into the same column of the device schedule. Supplemented – Addition of Amplatzer™ Septal Occluder and Amplatzer™ TorqVue™ LP Catheter and Delivery System.
2023-12-22	30072385	Supplemented – Addition of Epic™ Max Stented Tissue Valves, Epic™ Plus Stented Tissue Valves, Epic™ Max Heart Valve Sizer Set and Holder Handle, and Epic™ Plus Heart Valve Sizer Set and Holder Handles.
2024-03-08	30092058	Supplemented – Addition of Séguin™ Annuloplasty Ring, Séguin™ Annuloplasty Ring Sizer Set, Holder Handle and Extension Handle, Mechanical Heart Valve Accessories – Holder Handles, Leaflet Tester, Replacement Holder/Rotators and Sizer Set, Amplatzer™ Valvular Plug III and Navitor Titan™ Transcatheter Aortic Valve.

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 6 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Quality Management System Certificate

Regulation (EU) 2017/745, Annex IX Chapter I and III

MDR 750915 R000

Date	Reference Number	Action
2024-04-26	30106034	Supplemented – Addition of Amplatzer™ Muscular VSD Occluder and Amplatzer™ P.I. Muscular VSD Occluder, Regent™ Mechanical Heart Valve and Regent™ Mechanical Heart Valve with FlexCuff™, and Regent™ Mechanical Heart Valve Accessories – Replacement Holder/Rotators and Sizer Set.
2024-06-06	30167623	Supplemented – Addition of SJM Tailor™ Annuloplasty Ring, SJM Tailor™ Annuloplasty Band, Rigid Saddle Ring, SJM Tailor™ Annuloplasty Ring Sizer Set, Rigid Saddle Ring Sizer Set
Current	30169950	Supplemented – Addition of Amulet™ Steerable Delivery Sheath Restricted – Removal of Amplatzer™ Steerable Delivery Sheath, no longer placed on the market. Amended – Correction of typo for Piccolo device to remove the trademark (™) symbol from the word Amplatzer

First Issue Date: **2022-04-21**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

Page 7 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Manufacturer: Abbott Medical

Address:

177 County Road B East
St. Paul
Minnesota
55117
USA

Single Registration Number: US-MF-000018613

EU Authorised Representative: Abbott Medical


Address:

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium

Scope: See attached **Device Schedule**

On the basis of our assessment of the technical documentation in accordance with Regulation (EU) 2017/745, Annex IX Chapter II, the technical documentation meets the requirements of the Regulation. For the placing on the market of these devices an additional Annex IX Chapter I and III certificate is required.

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Graeme Tunbridge, Senior Vice President Global Regulatory & Quality

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Device Schedule:

Risk Classification: Class III

Type (Codes as per (EU) 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1001DS

Device Name	Model	Sheath Size (Fr)	Usable Length (cm)	Intended purpose (as per the IFU)
Amulet™ Steerable Delivery Sheath	STRBL-14F-075	14	75	The Amulet™ Steerable Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 2 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Risk Classification: Class III

Type (Codes as per (EU) 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1500EH

Device Name	Model	Sheath Size (Fr)	Curve Radius	Usable Length (cm)	Intended purpose (as per the IFU)
Amplatzer™ TorqVue™ Delivery System	9-ITV06F45/60	6	45°	60	The Amplatzer™ TorqVue™ Delivery System is intended to facilitate the attachment, loading, delivery and deployment of occluders.
	9-ITV07F45/60	7	45°	60	
	9-ITV07F45/80	7	45°	80	
	9-ITV08F45/60	8	45°	60	
	9-ITV08F45/80	8	45°	80	
	9-ITV09F45/80	9	45°	80	
	9-ITV10F45/80	10	45°	80	
	9-ITV12F45/80	12	45°	80	
	9-ITV13F45/80	13	45°	80	
	9-ITV05F180/60	5	180°	60	
	9-ITV06F180/60	6	180°	60	
	9-ITV06F180/80	6	180°	80	
	9-ITV07F180/80	7	180°	80	
	9-ITV08F180/80	8	180°	80	
	9-ITV09F180/80	9	180°	80	

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 3 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Risk Classification: Class III

Type (Codes as per (EU) 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1500EH

Device Name	Model	Sheath Size (Fr)	Curve Radius	Usable Length (cm)	Intended purpose (as per the IFU)
Amplatzer™ TorqVue™ Exchange System	9-EITV09F45/80	9	45°	80	The Amplatzer™ TorqVue™ Exchange System is intended to facilitate the removal and exchange of a delivery sheath with one of equal or larger diameter to enable the recapture and replacement of an occluder.
	9-EITV12F45/80	12	45°	80	
	9-EITV06F180/80	6	180°	80	
	9-EITV08F180/80	8	180°	80	

Risk Classification: Class III

Type (Codes as per (EU) 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1500EH

Device Name	Model	Sheath Size (Fr)	Curve Radius	Usable Length (cm)	Intended Purpose (as per the IFU)
Amplatzer™ Trevisio™ Intravascular Delivery System	9-ATV06F45/60	6	45°	60	The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of occluders.
	9-ATV07F45/60	7	45°	60	
	9-ATV07F45/80	7	45°	80	
	9-ATV08F45/60	8	45°	60	
	9-ATV08F45/80	8	45°	80	
	9-ATV09F45/80	9	45°	80	
	9-ATV10F45/80	10	45°	80	
	9-ATV12F45/80	12	45°	80	
	9-ATV13F45/80	13	45°	80	

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 4 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Risk Classification: Class III

Type (per EU 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1500EH

Device Name	Model	Sheath Size (Fr)	Usable Length (cm)	Intended purpose (as per the IFU)
Amplatzer™ TorqVue™ 2 Delivery Sheath	9-TV2-05F120	5	120	The Amplatzer™ TorqVue™ 2 Delivery Sheath is intended to provide a pathway through which a device is introduced within the chambers and coronary vasculature of the heart.
	9-TV2-06F120	6	120	
	9-TV2-07F120	7	120	

Risk Classification: Class III

Type (Codes as per (EU) 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1001DS

Device Name	Model	Sheath Size (Fr)	Usable Length (cm)	Intended purpose (as per the IFU)
Amplatzer™ Amulet™ Delivery Sheath	DS-TV45X45-12F-080	12	80	The Amplatzer™ Amulet™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.
	DS-TV45X45-14F-080	14	80	

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 5 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Risk Classification: Class III

Type (per EU 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1401EE

Device Name	Model	French Size	Length Dimension (cm)	Delivery Wire Length (cm)	Intended purpose (as per the IFU)
Amplatzer™ TorqVue™ LP Delivery System	9-TVLP4F90/060	4F	60	160	The Amplatzer™ TorqVue™ LP Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of occluders.
	9-TVLP4F90/080	4F	80	195	
	9-TVLP5F90/060	5F	60	160	
	9-TVLP5F90/080	5F	80	195	

Risk Classification: Class III

Type (per EU 2017/2185): MDN 1203

Basic UDI-DI: 5415067AMP1401EE

Device Name	Model	French Size	Length Dimension (cm)	Intended purpose (as per the IFU)
Amplatzer™ TorqVue™ LP Catheter	9-TVLPC4F90/080	4F	80	The Amplatzer™ TorqVue™ LP Catheter is intended to facilitate the loading, delivery, and deployment of occluders.

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 6 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 750953 R000

Certificate History

(References to applicable Common Specifications, Harmonized Standards complied with, and the relevant test and audit reports that support any of the below certificate changes may be requested from Certificate.Verification@bsigroup.com)

Date	Reference Number	Action
2022-11-29	3586410	Issued
2023-02-21	3447531	Supplemented – Addition of Amplatzer™ TorqVue™ Delivery System and Amplatzer™ TorqVue™ Exchange System.
	3447788	Supplemented – Addition of Amplatzer™ Trevisio™ Intravascular Delivery System.
2023-03-20	3447784	Supplemented – Addition of Amplatzer™ TorqVue™ 2 Delivery Sheath.
2023-05-17	3447786	Supplemented – Addition of Amplatzer™ Amulet™ Delivery Sheath.
2023-12-07	3746933	Supplemented – Addition of Amplatzer™ TorqVue™ LP Delivery System, and Amplatzer™ TorqVue™ LP Catheter
Current	30003596	Supplemented – Addition of Amulet™ Steerable Delivery Sheath (STRBL-14F-075) Restricted – Removal of Amplatzer™ Steerable Delivery Sheath (ASDS-14F-075) no longer placed on the market.

First Issue Date: **2022-11-29**

Current Issue Date: **2024-06-21**

Starting Validity Date: **2024-06-21**

Expiry Date: **2027-11-28**

...making excellence a habit.™

Page 7 of 7

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777443 R000

Manufacturer: Abbott Medical

Address:

177 County Road B East
St. Paul
Minnesota
55117
USA

Single Registration Number: US-MF-000018613

EU Authorised Representative: Abbott Medical


Address:

The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium

Scope: See attached **Device Schedule**

On the basis of our assessment of the technical documentation in accordance with Regulation (EU) 2017/745, Annex IX Chapter II, the technical documentation meets the requirements of the Regulation. For the placing on the market of these devices an additional Annex IX Chapter I and III certificate is required.

For and on behalf of BSI, a Notified Body for the above Regulation (Notified Body Number 2797):

Graeme Tunbridge, Senior Vice President Medical Devices

First Issue Date: **2023-11-15**

Current Issue Date: **2023-11-15**

Starting Validity Date: **2023-11-15**

Expiry Date: **2028-11-14**

...making excellence a habit.™

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777443 R000

Device Schedule:

Device Name: Amplatzer™ Sizing Balloon II

Type (Codes as per (EU) 2017/2185): MDN 1203

Intended Purpose as per the Instructions for Use: The Amplatzer™ Sizing Balloon II is intended to facilitate sizing of a cardiovascular defect.

Risk Classification: Class III

Basic UDI-DI: 5415067AMP1101DX

Model Number	Nominal Outer Diameter (mm)	Nominal Inflation Volume (cc)	Usable Length (cm)	Minimum Outer Diameter at Maximum Inflation Volume (mm)	Maximum Inflation Volume (cc)
9-SB-018	18	10	70	20	15
9-SB-024	24	25	70	25	30
9-SB-034	34	55	70	40	90

First Issue Date: **2023-11-15**

Current Issue Date: **2023-11-15**

Starting Validity Date: **2023-11-15**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 2 of 3

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

EU Technical Documentation Assessment Certificate

Regulation (EU) 2017/745, Annex IX Chapter II

MDR 777443 R000

Certificate History

(References to applicable Common Specifications, Harmonized Standards complied with, and the relevant test and audit reports that support any of the below certificate changes may be requested from Certificate.Verification@bsigroup.com)

Date	Reference Number	Action
Current	3748975	Issued

First Issue Date: **2023-11-15**

Current Issue Date: **2023-11-15**

Starting Validity Date: **2023-11-15**

Expiry Date: **2028-11-14**

...making excellence a habit.™

Page 3 of 3

Validity of this certificate is conditional on the Manufacturer's quality system being maintained to the requirements of the Regulation as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

NB Contact: BSI Group The Netherlands B.V., Say Building, John M. Keynesplein 9, 1066 EP, Amsterdam, Netherlands. Tel: + 31 (0) 20 346 07 80
Corporate Contact: BSI Group Assurance Limited, registered in England under number 05435540 at 389 Chiswick High Road, London, W4 4AL, UK.
A Member of the BSI Group of Companies.

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THE GOLD STANDARD FOR OVER 20 YEARS

OFTEN IMITATED, NEVER MATCHED

Unique design features
set the Amplatzer™ Talisman™
PFO Occluder apart from the rest.

INTAGLIO™
WIRE TREATMENT

ASYMMETRIC
DOUBLE DISC
DESIGN*

>94%
CLOSURE RATE†
at 6 months in
RESPECT trial¹

DOWNLOAD NOW

Download the free
Amplatzer™ Portfolio App
from your favorite app store.

DURABLE NITINOL
WIRE MESH WITH
POLYESTER
FABRIC THREAD

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† Effective closure

* Note: Only the 18mm size is not asymmetric.

¹ Saver JL, Carroll JD, Thaler DE, et al. Long-term outcomes of patent foramen ovale closure or medical therapy after stroke. *N Engl J Med* 2017; 377:1022-32

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AMPLATZER™ TALISMAN™ PFO OCCLUDER

INDICATION FOR USE

The Amplatzer™ Talisman™ PFO Occluder is indicated for percutaneous transcatheter closure of a patent foramen ovale (PFO) to reduce the risk of recurrent ischemic stroke in patients who have had an ischemic stroke due to a presumed paradoxical embolism.

PRODUCT HIGHLIGHTS

MR Conditional Information

Non-clinical testing has demonstrated the Amplatzer™ Talisman™ PFO Occluder is MR Conditional. A patient with the Amplatzer™ Talisman™ PFO Occluder can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm).
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)


Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 16 mm from the Amplatzer™ Talisman™ PFO Occluder when imaged with a gradient echo pulse sequence in a 3.0T MR system.

SIZING AND DEVICE SELECTION

Amplatzer™ Talisman™ PFO Occluder

Model/Reorder Number	Right Atrial Disc Diameter A	Left Atrial Disc Diameter B	Recommended Sheath Size
9-PFO-1818	18 mm	18 mm	8 F; 45° curve
9-PFO-2518	25 mm	18 mm	8 F; 45° curve
9-PFO-3025	30 mm	25 mm	9 F; 45° curve
9-PFO-3525	35 mm	25 mm	9 F; 45° curve

DELIVERY SHEATH

Amplatzer™ Talisman™ Delivery Sheath

Model/Reorder Number	Sheath Size	Tip Angle	Sheath Inner Diameter	Sheath Outer Diameter	Usable Length
9-TDS-08F45-80	8F	45°	2.69 mm	3.45 mm	80 cm
9-TDS-09F45-80	9F	45°	3.00 mm	3.81 mm	80 cm

ANCILLARY PRODUCTS

Amplatzer™ Guidewire

Model/Reorder Number	Diameter	Body	Tip Description	Usable Length
9-GW-002	0.035 inch	Super Stiff	1.5 mm, Modified J-tip	260 cm

DEVICE SIZING GUIDELINES

PFO Morphology	Example Anatomical Characteristics	Suggested Amplatzer™ Talisman™ PFO Occluder Size (mm)
Simple PFO or PFO with a non-prominent ASA PFO where a secure device position and effective PFO closure can be achieved when using the 25 mm device size	1. Absence of ASA, long tunnel, and thickened septum secundum 2. Non-prominent ASA (<20 mm total excursion) without a long tunnel (≥ 10 mm length) and without a thickened septum secundum (≥ 10 mm thickness)	25
Complex PFO PFO with one or more anatomical characteristics that may complicate the ability to achieve a secure device position and effective PFO closure when using the 25 mm device size	1. ASA (≥ 10 mm excursion) with long tunnel (≥ 10 mm length) 2. ASA (≥ 10 mm excursion) with thickened septum secundum (≥ 10 mm thickness) 3. Prominent ASA with excessive mobility (≥ 20 mm total excursion) 4. Lipomatous hypertrophy of septum secundum (≥ 15 mm thickness)	30 or 35
PFO with small anatomy Anatomy not suitable for 25 mm device size secondary to interference with adjacent cardiac structures	1. Septal primum length < 20 mm	18

Note: Evaluate the position of the device after deployment, but before detachment. Use echocardiography to ensure that the device does not impinge on the free atrial wall or aortic root. If the device interferes with an adjacent cardiac structure (such as free atrial wall or aortic root), recapture the device and redeploy. If device position remains unsatisfactory, recapture the device and either replace with a smaller device (18 mm or 25 mm) or refer the patient for alternative treatment.

The presence of an ASA alone does not necessarily prevent successful PFO closure with a 25 mm device size. In RESPECT¹, 180 patients (36%) in the device closure group had an ASA. The 25 mm device size was used in the majority of patients with an ASA (77%) to close the PFO, and at 6-months post-implant, effective closure was achieved in 95% of these patients. There were no cases of device embolization in any patient in the study.

¹. Saver JL, Carroll JD, Thaler DE, et al. Long-term outcomes of patent foramen ovale closure or medical therapy after stroke. *N Engl J Med*. 2017;377:1022-32.

PFO CLOSURE: FROM RARE TO ROUTINE

*“...closing the PFO may reduce the risk of having another stroke better than medication alone.”*¹ —Steven R. Messé, M.D.


Updated guidance from physician societies around the world supports percutaneous closure of a patent foramen ovale (PFO) to prevent stroke recurrence in select patients.

An expanded body of evidence prompted the American Academy of Neurology (AAN) and American Heart Association/American Stroke Association (AHA/ASA) to update their guidances regarding percutaneous PFO closure and risk of recurrent stroke.^{2,3}

A multidisciplinary team of neurologists, internists and cardiologists was convened for the AAN practice advisory to review the results from a meta-analysis of eight clinical studies.

In addition to the AAN findings, the updated AHA/ASA guideline also favors PFO closure, recommending it for patients with a PFO that have high-risk anatomical features, such as atrial septal aneurysm and/or large right-to-left shunt.

AAN SUMMARY OF EFFECTS OF PERCUTANEOUS PFO CLOSURE²

Relative risk reduction for
recurrent stroke compared to
medical management:
59%

- ▶ Absolute risk reduction of stroke at 5 years:
3.4%
- ▶ Periprocedural complication risk:
3.9%
- ▶ Increased absolute rate of non-periprocedural atrial fibrillation:
0.33% Per Year

AAN CLINICAL IMPLICATIONS²

PFO closure may be recommended for people < 60 years of age:

- When a PFO is present and no other mechanism of stroke has been identified after a thorough diagnostic evaluation.
- After discussing the potential benefits and risks.

PFO closure may be offered for people 60-65 years of age:

- After a thorough evaluation, including prolonged monitoring for atrial fibrillation.
- With very limited degree of traditional vascular risk factors (hypertension, diabetes, hyperlipidemia, smoking).
- In whom no other mechanism of stroke has been detected.

AHA/ASA CLINICAL IMPLICATIONS³

PFO closure may be recommended in patients 18-60 years of age:

- When stroke is thought to be caused by a PFO with high-risk anatomic features, such as atrial septal aneurysm and/or large shunt.
- If PFO is considered low risk anatomically and RoPE Score has been factored into clinical decision.

For more information about the Amplatzer™ Talisman™ PFO Occluder or updated guidelines, contact your Abbott sales representative.

RECOMMENDED BY MULTIPLE PHYSICIAN SOCIETIES FROM AROUND THE WORLD

Beyond the U.S., an increasing number of physician societies recommend PFO closure for well selected patients with a PFO-associated stroke, including:

- Japanese Stroke Society⁴ - 2021
- Australian Stroke Council⁵ - 2019
- Chinese Stroke Association⁶ - 2019
- German Society of Neurology/DSG German Stroke Society/ German Cardiology Society⁷ - 2018

EFFECTIVE PFO CLOSURE MADE EASIER⁸, WITH THE AMPLATZER™ TALISMAN™ PFO OCCLUDER

EXTENSIVE EXPERIENCE


OVER 180,000 PATIENTS TREATED GLOBALLY⁹

#1 device selected for PFO closure

An unmatched track record with over two decades of experience

Wide range of device sizes including a re-engineered 30mm with a smaller left atrial disc

CONFIDENCE IN CLOSURE

>94% CLOSURE RATE⁺
at 6 months in RESPECT trial¹⁰

SIMPLIFIED PREPARATION AND ENHANCED EASE OF USE

- The Amplatzer™ Talisman™ PFO Occluder comes assembled and ready for use, reducing preparation time.
- Fully recapturable and repositionable, ensuring optimal anatomical placement of the device for effective closure.

LONG-TERM PATIENT FOLLOW-UP¹⁰

5,810
patient-years of data

5.9
years median patient follow up

EXCELLENT SAFETY¹⁰

ZERO

- Device erosions
- Device thrombus
- Device embolization events
- Wire frame fractures

< 1% AF*
Low risk of serious atrial fibrillation

For more information about the Amplatzer™ Talisman™ PFO Occluder or updated guidelines, contact your Abbott sales representative.

THE LANDMARK *RESPECT* TRIAL

Studied in the landmark *RESPECT* trial, the AMPLATZER™ PFO Occluder is a double-disc device comprised of Nitinol mesh and polyester fabric.

RIGOROUS DESIGN DEMONSTRATING SIGNIFICANT STROKE REDUCTION¹

- MOST EXTENSIVE patient follow-up (13+ years of data)
- LARGEST trial, with broad inclusion criteria (including patients on anticoagulation therapy)
- EXCELLENT safety profile
- SIGNIFICANTLY reduced risk of recurrent ischemic stroke

THE RESULTS OF THE LANDMARK *RESPECT* TRIAL WERE PUBLISHED IN THE SEPTEMBER 14, 2017 ISSUE OF THE *NEW ENGLAND JOURNAL OF MEDICINE*.

THE MOST EXTENSIVE PATIENT FOLLOW-UP

The *RESPECT* trial enrolled patients at 69 centers across the U.S. and Canada, collecting a total of 5,810 patient-years of data.

~2
TIMES

More Patient-Years

The most extensive patient follow-up of any PFO closure trial ever conducted.

5.9
YEARS

Longer Median Follow-up

Patients in the *RESPECT* trial were followed significantly longer (nearly 2x more follow-up time) than patients in any other PFO trial.

THE LANDMARK PFO TRIAL

THE LARGEST TRIAL, WITH BROAD INCLUSION CRITERIA

The trial enrolled 980 patients with a PFO who had previously suffered a cryptogenic ischemic stroke. Patients on anticoagulation therapy were not excluded from the trial.

980

Total Patients Enrolled

The landmark trial included more patients than any other PFO trial conducted.

1:1 Randomization

Patients were randomly assigned to receive either the AMPLATZER™ PFO Occluder or medical therapy.


20%

Anticoagulation Therapy

Trial participants included approximately 20% of patients on anticoagulation therapy (a patient population thought to have a higher risk of venous thromboembolism).

EXCELLENT PROCEDURAL RESULTS

The device demonstrated significant technical and procedural success, and achieved excellent effective closure.

Technical Success

Designed for optimized ease of use, the device showed highly successful delivery and release.

Procedural Success

An excellent rate of successful implantation without in-hospital SAEs was demonstrated.

Effective Closure

(n≤9 Bubbles at 6 months)

Effective closure was achieved according to highly stringent criteria (2x more stringent than the criteria utilized in other PFO device trials).

DEMONSTRATING INVALUABLE PATIENT BENEFITS

AN EXCELLENT SAFETY PROFILE

Low rates of Serious Adverse Events (SAEs) were reported, including any SAE (0.064 per pt-yr rate), device related SAE (0.004 per pt-yr rate) and death (0.002 per pt-yr rate). Additionally, the trial showed no device embolization, no aortic erosion/dissection, no thrombus formation and low risk of atrial fibrillation, consistent with medical therapy.

0%

DEVICE EMBOLIZATION
AORTIC EROSION/DISSECTION
DEVICE THROMBUS


LOW RISK OF SERIOUS ATRIAL FIBRILLATION (AF)

There was no statistical difference in serious AF events between the closure and medical therapy groups (0.48 vs. 0.34 per 100 patient-years, $P = 0.36$).


SIGNIFICANTLY REDUCED RISK OF RECURRENT STROKE

The RESPECT trial demonstrated a clear benefit of closure with AMPLATZER™ PFO Occluder for reducing the risk of recurrent stroke and helping patients live healthy through significant stroke reduction.

The RESPECT trial showed there was a **45% RELATIVE RISK REDUCTION** for any recurrent ischemic stroke over nearly six years of follow-up with the use of the AMPLATZER™ PFO Occluder.

The AMPLATZER™ PFO Occluder also demonstrated a **62% RELATIVE RISK REDUCTION** for recurrent ischemic stroke of unknown mechanism.

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
Abbott Medical
177 County Road B East
St. Paul, MN 55117, USA
Tel: +1 855 478 5833

0072181 Rev. B

MDR Declaration of Conformity

Manufacturer:	Abbott Medical
Manufacturer SRN:	US-MF-000018613
Address:	177 County Road B East St. Paul, MN 55117 USA
Manufacturing Site(s):	Abbott Medical 5050 Nathan Lane North Plymouth, Minnesota 55442 USA
European Authorized Representative:	Abbott Medical The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium
European Authorized Representative SRN:	BE-AR-000008744

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Signature:  Teleshia Taylor De Sanchez Sr. Director Quality, Abbott Medical Structural Heart	<u>26 April 2022</u> Issue Date On behalf of Abbott Medical, signed at 177 County Road B East, St. Paul, Minnesota 55117
--	---

MDR Declaration of Conformity


Intended Purpose:	<u>Amplatzer™ Talisman™ PFO Occluder:</u> The Amplatzer™ Talisman™ PFO Occluder is a percutaneous, transcatheter occlusion device intended to close a patent foramen ovale. <u>Amplatzer™ Talisman™ Delivery Sheath:</u> The Amplatzer™ Talisman™ Delivery Sheath is intended to facilitate the delivery and deployment of an Amplatzer™ Talisman PFO Occluder.
Risk Classification:	<u>Amplatzer™ Talisman™ PFO Occluder:</u> Class III, per Medical Device Regulation 2017/745, Annex VIII, Chapter III, Section 5.4 <u>Amplatzer™ Talisman™ Delivery Sheath:</u> Class III, per Medical Device Regulation 2017/745, Annex VIII, Chapter III, Section 5.2
Classification Rationale:	<u>Amplatzer™ Talisman™ PFO Occluder:</u> Rule 8 - All implantable devices and long-term surgically invasive devices are classified as class IIb unless they are intended to be used in direct contact with the heart, the central circulatory system or the central nervous system, in which case they are classified as class III. <u>Amplatzer™ Talisman™ Delivery Sheath:</u> Rule 6 - All surgically invasive devices intended for transient use are classified as class IIa unless they are intended to be used in direct contact with the heart, the central circulatory system or the central nervous system, in which case they are classified as class III.
EMDN Code(s):	P07040303 - Cardiac Occluders, Patent Foramen Ovale And Patent Ductus Arteriosus C05 - Cardiovascular Introducing Sheaths
Basic UDI-DI:	Structural Intervention, Talisman System: 5415067AMP2000DX Structural Intervention, Talisman System (implant): 5415067AMP2001DZ

MDR Declaration of Conformity

	Structural Intervention, Talisman Sheath (DS): 5415067AMP2002E3
--	--

The products described in this declaration are in conformity with all applicable EU harmonized legislation, including:

- Regulation (EU) 2017/745, and the applicable *General Safety & Performance Requirements* in Annex 1 Directive 2006/42/EC on Machinery and Directive 89/686/EEC (and the superseding Regulation (EU) 2016/425) on Personal Protective Equipment do not apply.

Common Specifications Applied:	No applicable common specifications
STED #	0000044760
Notified Body:	BSI Group The Netherlands B.V. NB #: 2797 Say Building John M. Keynesplein 9 1066 EP Amsterdam The Netherlands
Supporting Certificate(s):	EC Certificate No: MDR 751010 R000 Expiration Date: 2027-04-20 QMS Certificate: MDR 750915 R000 Expiration: 2027-04-20
Original CE Mark Date:	2022-04-21
Conformity Assessment:	Annex IX
Device Photograph:	 Figure 1. Talisman Occluder  Figure 2. Talisman Delivery Sheath

Abbott Medical
177 County Road B East
St. Paul, MN 55117, USA
Tel: +1 855 478 5833

0072181 Rev. B

MDR Declaration of Conformity

The products in the attached Declaration of Conformity Product List are approved under EC Certificate MDR 751010.

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The signature is applied on page 1
88136 MDR Declaration of Conformity Template Rev F

Page 4 of 4

This confidential document is the property of Abbott Medical and shall not be reproduced, distributed, disclosed or used without the express written consent of Abbott Medical.

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MDR 750915 R000

Date	Reference Number	Action
Current	30106034	Supplemented – Addition of Amplatzer™ Muscular VSD Occluder and Amplatzer™ P.I. Muscular VSD Occluder, Regent™ Mechanical Heart Valve and Regent™ Mechanical Heart Valve with FlexCuff™, and Regent™ Mechanical Heart Valve Accessories – Replacement Holder/Rotators and Sizer Set.

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
The Amplatzer™ Talisman™ PFO Occluder is a percutaneous, transcatheter occlusion device intended to close a patent foramen ovale.

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The Amplatzer™ Talisman™ Delivery Sheath is intended to facilitate the delivery and deployment of an Amplatzer™ Talisman™ PFO Occluder.

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Date	Reference Number	Action
Current	3447814	Issued

First Issued: **2022-04-21**

Date: **2022-04-21**

Expiry Date: **2027-04-20**

...making excellence a habit.™

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© 2020 Abbott. All Rights Reserved. MAT-2003588 v1.0 | Item approved for Global OUS use only
© 2022 Abbott. All rights reserved. 9-EH-2-13625-01 09-2022