Artisse™ Detachment Device

fr	Dispositif de détachement Artisse™	no	Artisse™-frigjøringsenhet
de	Vorrichtung zur Ablösung des Implantats Artisse™	sk	Odpájacie zariadenie Artisse™
it	Dispositivo di distacco Artisse™	ro	Dispozitiv de detaşare Artisse™
es	Dispositivo de desprendimiento Artisse™	bg	Устройство за освобождаване Artisse™
sv	Artisse™ fråkopplingsenhet	et	Lahutusseade Artisse™
nl	Artisse™-loskoppelingsapparaat	hr	Proizvod za odvajanje Artisse™
pt	Dispositivo de destacamento Artisse™	lt	„Artisse™“ atskyrimo prietaisas
fi	Artisse™-irrotuslaite	lv	Artisse™ atvienošanas ierīce
da	Artisse™-frigørelsesenhet	sl	Pripomoček za odklop Artisse™
el	Συσκευή απόσπασης Artisse™	sr	Artisse™ uređaj za odvajanje
cs	Odpojovací nástroj Artisse™	pt-br	Dispositivo de desprendimento Artisse™
hu	Artisse™ leválasztó eszköz	ko	Artisse™ 분리 장치
ru	Устройство для отделения имплантата Artisse™	vi	Thiết bị tháo Artisse™
pl	Urządzenie do odłączania Artisse™	kk	Artisse™ ажырату құрылғысы
tr	Artisse™ Ayırma Cihazı	ar	Artisse™ جهاز الفصل

TABLE OF CONTENTS

Artisse™ Detachment Device

English.....	3
Français.....	5
Deutsch.....	8
Italiano.....	11
Español.....	14
Svenska.....	17
Nederlands.....	19
Português.....	22
Suomi.....	25
Dansk.....	28
Ελληνικά.....	30
Česky.....	33
Magyar.....	36
Русский.....	39
Polski.....	42
Türkçe.....	45
Norsk.....	47
Slovenčina.....	50
Română.....	52
Български.....	55
Eesti.....	58
Hrvatski.....	61
Lietuvių.....	63
Latviski.....	66
Slovenščina.....	68
Srpski.....	71
Português (Brasil).....	73
한국어.....	76
Tiếng Việt.....	79
Қазақша.....	82
العربية.....	87
Symbol Glossary.....	88

English

Instructions for Use

en

Artisse™ Detachment Device

CAUTIONS

- It is important to read the instructions for use with careful attention to cautions, notes and warnings prior to using this device.
- This device should be used only in the hospital operating room environment by physicians with a thorough understanding of angiographic and percutaneous neuro-interventional procedures.
- Do not use the device after the expiration date (Use-by Date) printed on the product label.
- Always confirm there are several backup Artisse™ Detachment Devices, in their unopened sterile package, as well as extra Sterile needles (20 G or 22 G) on the shelf before procedure. Verify that all Artisse™ Detachment Devices are within their indicated shelf life.
- Never open the main housing. Opening the main housing affects proper device functioning.
- This device is designed to minimize the effects of uncontrolled electromagnetic interference and other types of interference from external sources. Avoid use of other equipment that may cause erratic operation or degradation in the performance of this device.
- The device must be kept dry. It must not be submerged or be used in wet environment.
- The device should be used in the temperature range 10°C to 34°C.
- The device should be used in the Humidity range 30% to 85%.
- Do not use if sterile packaging is compromised or damaged.
- The device is intended for SINGLE PATIENT USE ONLY. DO NOT RESTERILIZE AND/OR REUSE. Re-sterilization may corrode the device, resulting in device malfunction.

ESSENTIAL PERFORMANCE

Artisse™ Detachment Device shall not allow an output of unintended current >2.1 mA.

DESCRIPTION


The Artisse™ Detachment Device is a hand held, battery operated device designed to initiate and control the detachment of the Artisse™ Implant. The Artisse™ Detachment Device is designed to apply a constant current through the device and to detect detachment.

The device maintains constant current by:

- Sensing the amount of resistance to current flow through the device.
- Adjusting the voltage required to maintain a constant current. Once detachment is completed, the device indicates detachment with an audible and visual "Detach" signal and stops the flow of current to the device.

WARNINGS

- No modification of this device is allowed. Device could be damaged.
- Use of this device adjacent to or stacked with other devices should be avoided because it could result in improper operation. If such use is necessary, this device and the other devices should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of this device. Otherwise, degradation of the performance of this device could result.

INTENDED PURPOSE / INDICATIONS FOR USE

The Artisse™ Detachment Device is intended for use only with the Artisse™ Intrascapular Device.

General Technical Data	
Device Dimensions:	130 mm x 40 mm x 30 mm (L x W x D)
Device Weight:	Approximately 100 grams
Climatic Class:	Operation at 10°C to 34°C
Atmospheric pressure:	68 kPa to 110 kPa
Humidity (relative, non –condensing):	30% to 85%
Maximum Operating Voltage:	28 Volts
Maximum Direct Current Output:	2.1 mA
Identification Plate:	Artisse™ Detachment Device

Artisse™ Detachment Device Visual and Audible Indication Chart			
Visual Indication		Audible Indication	Device Status
Battery Indicator	Detachment Indicator		
N/A	N/A	Single acoustic tone	- At power ON. - On single push button press to start detachment.
GREEN	OFF	N/A	Device ready for use.
GREEN	Continuously glowing and dimming BLUE	N/A	Detachment process is in progress. Note: Detachment process will last for 60 seconds.
GREEN	Interrupted flashing BLUE	Interrupted acoustic tone for 20 seconds	SUCCESSFUL DETACHMENT Precaution: Successful detachment must be verified by fluoroscopy that Artisse™ Implant has detached. On successful detachment follow instructions in Note 1 below.
GREEN	AMBER	Continuous acoustic tone for 20 seconds	Device encountered error in detachment or timed out after 60 seconds of continuous detachment process. If error in detachment is encountered follow instructions in Note 2 below.
AMBER	AMBER	N/A	Device has encountered system error. Discontinue use of the device.
AMBER	OFF	N/A	Device battery level is insufficient. Discontinue use of the device.

Note 1: User may acknowledge the indication by pressing the Push Button for at least 1/10 Sec but not more than 3 Sec.

Note 2:

- Ensure that the cables are properly connected.
 - If cables are not connected properly, connect the cables properly.
 - If cables are connected properly, check for the detachment using fluoroscopy. If not detached, restart detachment process.
- Acknowledge the error by pressing the push button for at least 1/10 Sec but not more than 3 Sec.
- After acknowledgement, Battery Indicator will light in GREEN and the device is ready for use for next detachment.

Artisse™ Detachment Device Push Button Function Chart	
To power ON device	Press and hold Push Button for minimum 3 Sec.
To power OFF device	Press and hold Push Button for minimum 4.5 Sec.
To start detachment	Press Push Button for at least 1/10 Sec but not more than 3 Sec.
To acknowledge detachment success	Press Push Button for at least 1/10 Sec but not more than 3 Sec.
To acknowledge detachment Error	Press Push Button for at least 1/10 Sec but not more than 3 Sec.
Note: The device shall power OFF, If kept idle for 15 minutes.	

CONTRAINDICATIONS

There are no known contraindications for this device.

DIRECTIONS FOR USE

Patient Connection and Cable Sets

Cables connection shall only be made after the Artisse™ Implant is placed at the desired location and ready for detachment.

- Insert a sterile needle (20 G or 22 G = 0.7 or 0.9 mm, respectively) at the shoulder.

- b. Press the end of the connector of the sterile black cable so that the J-hook protrudes from the distal end and clip it on the needle.
- c. Press the end of the connector of the sterile red cable so that the J-hook protrudes from the distal end and clip to the proximal end of the pusher wire of the Artisse™ Implant.

Note: If the polarity is incorrect, detachment will not be possible.

Switching "ON" the Artisse™ Detachment Device and Testing:

Press and hold the Push button for approximately 3 seconds to switch ON the device. Battery Indicator and Detachment Indicator shall flash. Upon successful self-diagnosis of the device, the Battery Indicator shall turn GREEN to confirm the device is ready for use. If the Battery and Detachment Indicators are AMBER, the device is in System Error state and shall not be used.

After powering ON, the Artisse™ Detachment Device shall remain ready until it receives user input, maximum idle time of 15min is reached, or insufficient battery capacity reached.

STARTING THE DETACHMENT

- a. Once the cable connection to the patient has been made, the detachment may start. This is performed by pressing the Push button on the device for at least 1/10 Sec but not more than 3 Sec.
- b. The Detachment Indicator shall continuously glow and dim in BLUE without any audible tone throughout the detachment cycle.

PRECAUTION

- Before starting the detachment process, verify proper connection of the Anode ("red") and Cathode ("black") Clips from the Artisse™ Detachment Device to the respective Pusher wire & Needle.

STOP THE DETACHMENT

In case of emergency, Artisse™ Detachment Device detachment process may be stopped by pressing and holding the Push Button for minimum 4.5 Sec.

- "Starting a new detachment" – Artisse™ Detachment Device to be powered ON by pressing push button for minimum of 3 Sec.

SUCCESSFUL DETACHMENT

The device shall indicate the successful detachment to user with interrupted flashing BLUE Detachment Indicator and interrupted acoustic tone. User may acknowledge the indication by pressing the Push Button for at least 1/10 Sec but not more than 3 Sec. After acknowledgment, the device goes back into Ready state and a new detachment process can be started.

CAUTION

- DO NOT remove cathode/anode cable connector clips from needle/pusher wire during detachment process.**
- Accidental connector Clip removal/disconnection from needle/pusher wire during the detachment process may elicit a false "Successful Detachment" indication. Always visually confirm implant non detachment using fluoroscopy, before starting a new detachment.

PRECAUTION

- Before pulling out Pusher wire, visually confirm detachment status using fluoroscopy.
- Successful detachment of the Artisse™ Implant must be verified by fluoroscopy to ensure that the Artisse™ implant has detached.
- Once Artisse™ implant detachment has been detected and visually confirmed using fluoroscopy, slowly remove the Clips from the needle and pusher wire.

INTERRUPTION OF THE ELECTRICAL CONNECTION TO THE PATIENT

If the connection to the patient is interrupted in the first 10 seconds during the detachment process, the device will recognize this and indicate with a continuous AMBER on Detachment Indicator and a continuous acoustic tone for 20 seconds approximately. User may acknowledge the indication by pressing the Push Button for at least 1/10 Sec but not more than 3 Sec. Ensure that the cable clips are properly connected prior to starting a new detachment process.

MAXIMUM DETACHMENT TIME

The maximum detachment time has been set to 60 seconds. After maximum detachment time has elapsed, the device will sound a continuous acoustic tone with a continuous AMBER light ON visual indication.

POWERING OFF THE Artisse™ Detachment Device

The Artisse™ Detachment Device may be powered OFF by pressing and holding the Push Button for approximately 5 seconds. Confirm Device has Powered OFF by visually observing Battery Indicator light turning OFF.

REQUIRED ADDITIONAL ITEMS AND ACCESSORIES

Always have an additional Artisse™ Detachment Device (in sterile pouch) and extra Sterile needles (20 G or 22 G) available as backup.

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to:

- Access site complications such as infection
- Complications of radiation exposure such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia
- Device complications such as premature detachment or non-detachment
- Hematologic complications such as intracranial hemorrhage, hemolysis
- Neurological deficit/dysfunctions such as cerebral infarction, stroke
- Systemic complications such as hypersensitivity, toxicity, shock
- Vascular complications such as embolism, vasoconstriction, inflammation, necrosis, stenosis, thrombosis

*Consult Instructions for Use for other required therapy devices and medications for additional potential complication information. If a serious incident related to the device occurs, contact your Medtronic representative and the competent authority in your respective country/region.

WARNINGS

- The Artisse™ Detachment Device must not be operated in the proximity of or in connection with equipment that emits short waves or microwave radiation nor with HF surgical equipment.
- The Artisse™ Detachment Device is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Visually inspect all sterile barrier systems, that are labeled as sterile, immediately prior to use. Do not use the device if breaches in sterile barrier system integrity are evident.
- For additional Materials of Concerns information such as REACH, CA Prop 65 or other product stewardship programs, go to www.medtronic.com/productstewardship

PRECAUTIONS

The Artisse™ Detachment Device has been tested and found to comply with the limits for medical devices to the **IEC 60601-1-2:2014/02/01 Ed: 4.0**. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference

to other devices in the vicinity. However there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to other devices, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the devices.
- Consult Medtronic.

TECHNICAL SPECIFICATIONS

Declaration about the Electromagnetic Compatibility of Artisse™ Detachment Device according to EN60601-1-2 Ed: 4.0

Electromagnetic Emission

The Artisse™ Detachment Device is intended to be used in the electromagnetic environment specified below. The customer or user of the device should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
RF-Emissions CISPR 11	Group 1	The device uses RF Energy only for its internal function. Therefore the RF Emission is very low and it is more unlikely that nearby equipment is affected.
RF-Emissions CISPR 11	Class [A]	The Artisse™ Detachment Device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic Immunity

The Artisse™ detachment Device is intended to be used in the Electromagnetic environment specified below.

The customer or user of the device should ensure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic Environment Guidance
ESD 61000-4-2	± 8 kV Contact ±15 kV Air	±2/4/6/8 kV Contact ±2/4/8/15 kV Air	Floors should be wood, concrete or Ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not Applicable	Not Applicable
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not Applicable	Not Applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Not Applicable	Not Applicable
Power Frequency 50/60Hz magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Suggested distances between mobile RF-Communication Devices and Artisse™ Detachment Device

The Artisse™ Detachment Device is intended for use in the electromagnetic environment in which radiated Disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Distance depending on the Sending Frequency:

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150kHz to 80MHz $d=1.17 \sqrt{P}$	80MHz to 800Mhz $d=1.17 \sqrt{P}$	800MHz to 2.5 GHz $d=2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150kHz to 80MHz $d=1.17 \sqrt{P}$	80MHz to 800MHz $d=1.17 \sqrt{P}$	800MHz to 2.5 GHz $d=2.33 \sqrt{P}$
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Electromagnetic Immunity

The Artisse™ Detachment Device is intended to be used in the electromagnetic Environment specified below. The customer or user of the device should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test level	Compliance Level
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms (V1)
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m (E1)

Electromagnetic Environment Guidance

Portable and mobile RF communications equipment should be used no closer to any part of the Artisse™ Detachment Device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

$d=1.17 \sqrt{P}$ for 150 kHz to 80 MHz

$d=1.17 \sqrt{P}$ for 80 MHz to 800 MHz

$d=2.33 \sqrt{P}$ for 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^a.

Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Artisse™ Detachment Device is used exceeds the applicable RF compliance level above, the Artisse™ Detachment Device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Artisse™ Detachment Device.
- b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

HOW SUPPLIED

- This device is supplied STERILE using ethylene oxide. This device is non-pyrogenic.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Micro Therapeutics, Inc.
Medtronic B.V.