

VANCOGENX[®] - SPACE KNEE


TECRES[®]
ADVANCING HIGH TECHNOLOGY

Vancogen[®]-Space Knee

Ref	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	Gentamicin base (g)	Vancomycin base (g)	Overall weight (g)
SPK0030	54	40	60	36	16	0,9	0,9	47,71
SPK0130	64	47	70	42	17	1,3	1,3	72,38
SPK0230	74	54	80	48	18	1,8	1,8	94,94
SPK0330	84	61	90	54	19	2,7	2,7	140,65

Vancogen[®]-Space Knee Ats

Ref	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	Gentamicin base (g)	Vancomycin base (g)	Overall weight (g)
SPK0430	60	36	7	32	11	0,3	0,3	17,02
SPK0530	60	36	12	32	11	0,5	0,5	26,77
SPK0630	80	48	7	40	11	0,5	0,5	27,06
SPK0730	80	48	12	40	11	0,8	0,8	43,78

VANCOGENX[®]-SPACE KNEE

EN

INSTRUCTIONS FOR USE

To the Medical Staff's attention

TEMPORARY SPACER FOR KNEE & TEMPORARY AUGMENTED TIBIAL STEM (ATS) SPACER WITH GENTAMICIN and VANCOMYCIN

These IFU refer generically to "Vancogenx[®]-Space Knee" when they consider all three components (femoral, tibial and ATS); on the contrary they refer to femoral and tibial components as "Vancogenx[®]-Space Knee (F&T)" and to ATS component as "Vancogenx[®]-Space Knee ATS" when they refer to these specific device parts.

PRESENTATION

Vancogenx[®]-Space Knee is a temporary implantable spacer device for the knee joint.

Vancogenx[®]-Space Knee includes and releases two medicinal substances: gentamicin sulfate and vancomycin hydrochloride.

Vancogenx[®]-Space Knee is a sterile and single use medical device. The device is sterilized with gamma irradiation.

Vancogenx[®]-Space Knee devices are available in the following variants:

- Vancogenx[®]-Space Knee (F&T) in 4 sizes (S M, L, XL).
- Vancogenx[®]-Space Knee ATS in two sizes with two different thicknesses.

The use of Vancogenx[®]-Space Knee ATS is optional and it is recommended in case of a large tibial defect. This component is available in different sizes and shall be coupled, using Vancogenx[®] bone cement, with the tibial component of the following Vancogenx[®]-Space Knee (F&T):

Vancogenx [®] -Space Knee ATS	Combination with Vancogenx [®] -Space Knee (F&T)
Vancogenx [®] -Space Knee ATS 60/07 (ref. SPK0430)	Vancogenx [®] -Space Knee 6054 (ref. SPK0030)
Vancogenx [®] -Space Knee ATS 60/12 (ref. SPK0530)	or Vancogenx [®] -Space Knee 7064 (ref. SPK0130)
Vancogenx [®] -Space Knee ATS 80/07 (ref. SPK0630)	Vancogenx [®] -Space Knee 8074 (ref. SPK0230)
Vancogenx [®] -Space Knee ATS 80/12 (ref. SPK0730)	or Vancogenx [®] -Space Knee 8094 (ref. SPK0330)

DESCRIPTION

Vancogenx[®]-Space Knee has a three-piece design, which resembles an ultra-congruent condylar knee prosthesis.

- Vancogenx[®]-Space Knee (F&T): two articulating independent elements, a tibial component made of a flat base upon which a femoral component articulates.
- Vancogenx[®]-Space Knee ATS: an additional optional component, which can be used in case of a large tibial defect.

Vancogenx[®]-Space Knee is composed of fully formed polymethylmethacrylate (PMMA) radiopaque bone cement loaded with gentamicin sulphate and vancomycin hydrochloride. The PMMA bone cement includes barium sulphate as radiopaque agent and gentamicin sulphate and vancomycin hydrochloride as antibiotics.

The formulation of the fully formed resin in the finished product is:

	Formulation [% w/w]
PMMA*	84,90
BaSO ₄	10,00
Gentamicin sulphate	3,20
Vancomycin hydrochloride	1,90

* Polymethyl Methacrylate, N,N-dimethyl-p-toluidine, Hydroquinone, Benzoyl peroxide, Methyl methacrylate

The design allows for limited motion of the leg for basic mobility (sitting, standing, walking) under limited weight bearing conditions (e.g. protected loading with use of crutches, walkers, canes etc.).

Dimensions, gentamicin content and vancomycin content of each size and variant of Vancogenx[®]-Space Knee devices are reported in the tables placed at the beginning of this leaflet.

INTENDED USE

Vancogenx[®]-Space Knee is intended for temporary use (maximum 180 days) as an adjunct to total knee replacement (TKR) in skeletally mature patients undergoing a two-stage procedure due to a septic process.

Vancogenx[®]-Space Knee is not intended for use for more than 180 days, at which time it must be explanted and a permanent device implanted or another appropriate treatment performed (e.g. resection arthroplasty, fusion etc.).

It is particularly intended for operations which present risks of or existing infections caused by organisms sensitive to Gentamicin and/or Vancomycin.

Vancogenx[®]-Space Knee ATS is intended for temporary use (maximum 180 days) in a total knee replacement (TKR) in patients undergoing a two-stage procedure due to a septic process and where a large tibial defect is present.

The device must be combined with the tibial component of Vancogenx[®]-Space Knee.

Vancogenx[®]-Space Knee ATS is not intended for use for more than 180 days, at which time it must be explanted and a permanent device implanted or another appropriate treatment performed (e.g. resection arthroplasty, fusion etc.).

It is particularly intended for operations which present risks of or existing infections caused by organisms sensitive to Gentamicin and/or Vancomycin.

INDICATIONS FOR USE

Peri-prosthetic joint infection (PJI) of the knee joint.

CONTRAINDICATIONS

- The patient's condition is such that a two-stage arthroplasty procedure is contraindicated due to decreased immune response or other relevant systemic clinical conditions.
- Bone loss precluding adequate support of the prosthesis.
- Lack of adequate competence (anatomical and functional) of peripheral ligamentous apparatus and extensor mechanism.
- The procedure is unjustified due to deficiencies in the patient's muscular, nervous or vascular systems.
- Poor bone quality (as in osteoporosis) could cause the prosthesis to migrate or to fracture host bone.
- Infection of the TKR cannot be confirmed.
- The infected TKR devices cannot be removed.
- The infecting pathogens are resistant to gentamicin and vancomycin.
- The patient is sensitive (allergic) to gentamicin, aminoglycosides, vancomycin, glycopeptides or PMMA bone cement.
- A systemic or secondary remote infection is suspected or confirmed.
- The patient does not have a TKR and the infection is secondary to trauma, septic arthritis or other surgical procedures.
- The patient does not have sufficient bone stock to allow insertion and fixation of the prosthesis
- The patient has neuromuscular disorders that do not allow control of the knee joint.
- The patient's weight, age or activity level would cause the surgeon to expect early failure of the system.
- The patient is unwilling or unable to use protected weight bearing mobility throughout the implantation period (e.g. crutches, canes, walkers etc).
- Myasthenia gravis.
- Severe renal impairment.
- Pre-existing calcium metabolism disorder.

INTENDED USERS

The device is intended to be used by suitably qualified orthopaedic surgeons, experienced in two-stage knee arthroplasty procedures, in an operating theatre environment.

PATIENT TARGET GROUP(S)

The device is intended to be used in skeletally mature patients suffering from a peri-prosthetic joint infection (PJI) of the Knee joint.

CLINICAL BENEFITS & PERFORMANCES

Clinical performance

Main

- Maintain the correct joint space between femur and tibia following prosthesis removal.
- Allow partial weight-bearing in relation to device stability and anatomical condition.
- Restore a (partial) joint functionality which allows movement and

possibility to perform physical therapy, providing a higher quality of life.

- Facilitate the second stage revision as muscle retraction and joint stiffness are avoided.
- Reduce functional recovery times after implantation of the revision prosthesis thanks to the maintenance of a certain muscle tonicity

Ancillary

- Help to oppose the present infection with the controlled and effective release of the antibiotics directly in the infection site all over the implantation period.

Clinical benefits

- The patient can maintain a certain joint mobility, walk with partial weight-bearing (according to surgeon indications based on device stability and clinical/anatomical condition) and perform some physical therapy (active and passive) and muscle strengthening exercises: this will make easier the second-stage and quicker the functional recovery after the second stage since muscle retraction and joint stiffness are avoided.
- The patient can perform basic movements, generally without the need of help: self-independence is a great psychological help; moreover, this

may limit the need of an external help (caregiver).

- The first operation (1st stage) is shorter since the device is ready-made, and the second operation is generally facilitated and shorter as well (2nd stage) thanks to the muscle tonicity. This means that also the anaesthesiological risks and any surgery-related risk are reduced.
- Earlier hospital dismissal is generally allowed to the patient.
- The presence of the antibiotics protects from bacterial adhesion and the local release of the antibiotics participates to the treatment of the infection.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCES

The Summary of Safety and Clinical Performances (SSCP) is available on the European Database on medical devices (EUDAMED) at the link <https://ec.europa.eu/tools/eudamed>, within six months after the date of publication of the notice of its full functionality.

ADVERSE EVENTS

The following adverse events may arise when using a spacer. Some are typical of any surgery, others are specific of arthroplasty operations. Most are not directly attributable to the spacer as such; however, surgeons must be aware of these complications and be ready to treat them if they occur.

General Surgical Risks

- Arrhythmias
- Myocardial infarction
- Pulmonary embolism
- Sudden death
- Transitory hypotension
- Venous thrombosis

Knee replacement Surgery Risks

- Damage to femur or tibia
- Arthrofibrosis
- Damage to blood vessels
- Extensive blood loss
- Hematoma
- Limb length discrepancy
- Nerve damage, bone bed damage
- Phlebitis, thrombophlebitis
- Wound healing problems

Vancogen[®]-Space Knee Device Risks

- Adverse soft tissue reaction
- Allergic pyrexia
- Debris release
- Device breakage
- Device dislocation
- Device sub-luxation

- Device loosening
- Difficult device removal
- Gentamicin-associated toxicity (nephro- / ototoxicity)
- Recurrent infection
- Vancomycin-associated toxicity (nephro- / ototoxicity)

Other adverse reactions associated with systemic use of Gentamicin therapy include hypersensitivity, anaphylactic reactions, nausea, vomiting, urticaria, reversible granulocytopenia, anemia, blood dyscrasia, convulsions, central nervous toxicity, abnormal hepatic function, hypomagnesaemia, stomatitis, purpura, allergic contact sensitization and neuromuscular blockade.

Adverse reactions associated with systemic use of Vancomycin include hypersensitivity, anaphylactic reactions, rash, urticarial, pruritic, toxic epidermal necrolysis.

The postoperative follow-up should be closely monitored by the surgeon for patients reporting a previous experience of sensitivity or allergic reaction to other antibiotics, in particular of the aminoglycosides' family.

Attention: since the device must be fixed with antibiotic loaded bone cement, the surgeon must be aware of its negative effects as it may increase the potential of toxic drug reactions (gentamicin/vancomycin-associated nephro-/oto- toxicity). Please refer to section WARNINGS - PRECAUTIONS FOR USE.

INFORMATION FOR USE

The present instructions for use leaflet must be read carefully.

Familiarity with and attention to Vancogen[®]-Space Knee surgical technique and familiarity with proper use of the Vancogen[®]-Space Knee is essential for successful use of the device.

Only surgeons who have read Vancogen[®]-Space Knee surgical technique and are aware of the limitations of its application should utilize the device.

WARNINGS - PRECAUTIONS FOR USE

Do not use after the expiration date since the effectiveness of the device may be compromised.

Altering the structural identity of the device can cause damages to the same (creation of cracks or fragments and risk of sterility loss leading to infection risks for the patient).

Surgeons must make sure that the device to be implanted is of the right size for the site of implantation.

Remove the spacer not later than six months after its implantation. The device was tested to be safely used for not more than 6 months. If this period is extended for too long this can lead to wearing, development of debris and eventually to breakage that can cause pain, inflammation and bone re-absorption.

Aggressive assembly methods are not required for proper implantation of the device.

During the application do not subject the device to excessive forces (e.g. hammer strikes) that could cause damage.

Any damage to the device may affect the fatigue strength and lead to failure under load.

Using the device under conditions other than the indicated use is unlikely to provide benefit to the patient and increases the risk of the development of drug-resistant bacteria.

Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient.

Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the device's morphology, the efficiency of the antibiotics and of the mechanical features, causing a malfunction of the same with serious risks for the patient's health.

Re-using the device after its extraction must by all means be avoided since it can cause contamination and worsening of the patient's infection.

During the extraction, the spacer can also be damaged and residues of the cement can remain in the operatory site: in this case all residues of cement shall be fully removed.

The device is foreseen to be used in conjunction with systemic antimicrobial antibiotic therapy (standard treatment approach to an infection) as defined by physician and/or surgeon.

Vancogen[®]-Space Knee should be used with caution:

- in conjunction with other nephrotoxic or ototoxic drugs (e.g. other aminoglycosides, amphotericin B, cephalosporins, ciclosporins, cisplatin, clindamycin, colistin, polymyxin B, foscarnet, vancomycin, viomycin), since they can increase gentamicin toxicity (the low serum level of gentamicin makes this interaction very unlikely to occur);
- in patients who are predisposed to or who have pre-existing clinical conditions that would put them at risk for gentamicin and vancomycin toxicity (e.g. renal dysfunction, hearing difficulties, dehydration, advanced age, taking drugs which may affect kidneys, undergoing general anaesthesia, etc).

In such cases, patients should be monitored for toxic blood levels of gentamicin and vancomycin (mainly in the first days following implantation), as well as for kidney function.

Simultaneous use of gentamicin and/or vancomycin and strong diuretics, such as ethacrynic acid or furosemide, may intensify the ototoxic effect of gentamicin and vancomycin by altering the antibiotic concentrations in serum and tissue.

In cases of significant obesity gentamicin serum concentrations should be closely monitored and a reduction in dose should be considered.

Co-administration with the following agents should be avoided:

- Neuromuscular blocking agents such as succinylcholine, botulinum toxin risk of toxicity due to enhanced neuromuscular block
- Anticoagulants such as warfarin and phenindione as Gentamicin has been known to potentiate them
- Cholinergics: antagonism of effect of neostigmine and pyridostigmine.
- Bisphosphonates: increased risk of hypocalcaemia.
- Sulfites in susceptible people, especially those with a history of asthma or allergy: risk of allergic-type reactions including anaphylactic symptoms and bronchospasm.
- Piperacilline and tazobactam as they can potentiate vancomycin nephrotoxicity.

The use of Vancogen[®]-Space Knee may result in overgrowth of non-susceptible organisms. If overgrowth of non-susceptible organism occurs, appropriate therapy should be initiated.

The use of antibiotic-loaded bone cement for fixation of the device may increase the potential of toxic drug reactions.

In vitro elution studies (microbiological method) has shown that the daily release of gentamicin and vancomycin never exceeds the recommended systemic adult dose for gentamicin (5-7 mg/kg/day, adult with normal renal function) and vancomycin (30-45 mg/Kg/day) according to Goodman and Gilman's recommendations¹. It is therefore unlikely that the amount of gentamicin and vancomycin absorbed locally from Vancogenx®-Space Knee will result in serum levels in the toxic range.

Owing to the presence of gentamicin and vancomycin, adverse events of these substances may occur, even if occurrence is unlikely because of the very low serum levels reached.

Gentamicin and vancomycin are potentially nephrotoxic and/or ototoxic.

Despite the risk of accumulation is low, care should be taken, and the level of gentamicin and vancomycin should be monitored in patients with impaired renal function, and in patients with previous neuromuscular disorders (e. g. myasthenia gravis, Parkinson's disease), in particular if they also have a history of renal insufficiency.

In case of administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified (using calcium salts it may be reversed). However, this is quite unlikely in patients with healthy kidneys, whereas in patients with impaired renal function, the probability of interaction increases in proportion to serum levels.

Plasma trough concentrations of gentamicin higher than 2µg/ml for periods longer than 10 days have been associated to toxicity.

Ototoxicity is associated to an excessively high concentration of vancomycin in plasma (between 60 and 100 µg/ml). Nephrotoxicity induced by vancomycin is unusual, when normal doses are used. It depends on kidney function and measurement of antibiotic concentration in blood.

Use of Vancogenx®-Space Knee ATS in combination with Vancogenx®-Space Knee (F&T)

The use of Vancogenx®-Space Knee ATS is expected to lead to a higher systemic exposure of gentamicin and vancomycin compared to using the Vancogenx®-Space Knee alone.

USE DURING PREGNANCY, BREAST-FEEDING AND IN CHILDREN

There are no available tests which demonstrate the safe use of the device during pregnancy or breast-feeding, and in children.

During pregnancy and breast-feeding the surgeon should weight the benefit for the mother against the potential risk to the child before using the device. It is advisable not to implant the device during pregnancy and breast-feeding. In these patients, the device should only be used in critical, life endangering situations.

In the event that the device is used during pregnancy or if the patient becomes pregnant while the device is in use, she should be notified of the possible dangers to the foetus. Women of childbearing potential should consider the benefits and dangers associated with the use of the product. The device is indicated for use in children (skeletally young patients) only when the surgeon deems that saving the joint through other forms of intervention is not possible.

STERILE PACKAGING INFORMATION

The device is contained inside a double blister packaging configuration and it is sterilised by gamma radiation.

Please discard the device if the blisters are damaged or unintentionally opened before use because the sterility of the device may be compromised. Do not re-sterilize.

STORAGE & TRANSPORT CONDITIONS

Store and transport at a temperature below 25°C and in a cool and dry environment.

However, if the temperature during the transport cannot be maintained under 25 °C with a relative humidity < 90%, the medical device can tolerate a temperature between 0 °C and 40 °C with a relative humidity less than 90% for maximum 5 days.

The device must be shipped and transported in a double wall corrugated box with ECT ≥ 44 lb/in (or ECT ≥ 7.71 kN/m).

UTILIZATION, IMPLANTATION and EXPLANATION

It is important to maintain strictly aseptic surgical techniques.

The surgeon shall become thoroughly familiar with the technique of implantation of the prostheses by: (1) reading appropriate literature, (2) training in the operative skills and techniques required for Knee arthroplasty revision surgery, and (3) reviewing information regarding use of instrumentation for sizing, implantation and explanation of the components.

In order to prevent dislocation/sub-luxation, the same measures taken for a permanent total Knee replacement (TKR) are advised, plus other specifics such as:

- 1) choice of the correct size (refer to Section "Vancogenx®-Space Knee Size Choice");
- 2) proper cement fixation (Vancogenx® bone cement) of the components;
- 3) insertion with appropriate joint tension of the soft tissues around the knee joint;
- 4) application, in cases at risk, of an orthopaedic abduction brace (possibly articulated) to assist flexion without dislocation/sub-luxation.

No specific instrumentation is required to successfully implant the device; however, it is recommended that the following instruments and devices be available for the procedure:

- Trial and Templates,
- Vancogenx® bone cement for fixation purposes,
- Equipment for knee arthroplasty surgery (such as osteotomies, rasps, reamers, impactor without metallic tips, etc.).

For the application of the device the front median or the medial para-patellar route may be used.

The infected prosthesis components and the residual bone cement (if any) must be removed from the femur, tibia and patella before the device implantation.

Clean the host site by aggressive debridement, removing all residues.

Do not insert the device if the joint hosts other osteosynthesis devices that can interfere mechanically with the device.

Vancogenx®-Space Knee Size Choice

Selection of the proper size of the device depends on the judgement of the surgeon with relationship to the requirements of the patient.

¹ *The Pharmacological Basis of Therapeutics 13th Edition (Brunton L, Hilal-Dandan R, Knollmann B) - McGraw Hill, New York. Chapter 58: Aminoglycosides (McDougall C); Chapter 59: Protein Synthesis Inhibitors and Miscellaneous Antibacterial Agents (McDougall C)*

For the implant size selection, templates and trial devices are available.

Size appropriateness can be evaluated preoperatively with the superimposition of Vancogen[®]-Space Knee templates on patient x-ray imaging.

Size appropriateness should be confirmed through trial reduction, achieving a compromise between stability and joint mobility.

Size selection for Vancogen[®]-Space Knee device can be selected by considering:

- Dimensions of the removed implant
- Remaining bone stock
- State of ligamentous apparatus
- Flexion and extension spaces

Trial Placement

Position the Trials into the joint space and reduce the knee. The knee should not be too tight as it will tighten further upon cement fixation of components.

Tightness may be relieved through downsizing and/or recontouring the femoral bone to achieve a satisfactory fit.

Once the trial has been used, remove it from the joint, and use it as a reference to select the implantable spacer device size.

Caution: Trials must not be implanted. The thickness generated by the cement is not reflected in the trial.

Vancogen[®]-Space Knee Implantation

Vancogen[®]-Space Knee is applied using the front median or the medial para-patellar route.

After the removal of the pre-existent prosthesis, the host site must be irrigated with Ringer or physiological solution and a thorough debridement must be performed before inserting the Vancogen[®]-Space Knee device. In particular, eliminate, wherever possible, all residuals of cement left by the previous implant.

Avoid washing with aqueous solutions the spacer before or after implantation in order to maintain optimal levels of antibiotic release.

Care should be taken to keep the wound dry before the spacer is inserted to avoid loss of antibiotics at the surface of the implant, as the spacer antibiotics are activated in an aqueous environment.

The most suitable size to use is selected on the basis of:

- the size of the middle lateral and antero-posterior femoral and tibial prosthetic components removed;
- the type of the bone defect (in case of lack of a significant part of the tibial plateau, the use of Vancogen[®]-Space Knee ATS is recommended);
- the ligamentous apparatus state;
- the flex and extension spaces.

For the selection of the size transparent radiograph overlays and trial devices are available.

Note that the thickness generated by the cement is not reflected in these templates.

Among the sizes, the size to be implanted is that which:

- is nearest to the size of the removed implant;
- achieves the best compromise between stability and joint mobility: during the operation, it is important to achieve full extension and a 90° flex, with a flex area sufficiently close to avoid anteroposterior movement of the flexed knee.

Vancogen[®]-Space Knee components must be tightly fixed with bone cement antibiotic loaded Vancogen[®].

The entire contact surface of the components must be cemented to create continuity between spacer and bone. Missing or insufficient cement - especially in the posterior portion of the condyles - may weaken the structure of the device.

Following removal of the existing knee implant and debridement, apply the femoral part first on the femoral condyles, and wait for the polymerization of the cement, then proceed with the application of the tibial part on the tibial plateau, avoiding in both cases that excess cement adheres to the joint surfaces.

In case of large tibial defect, connect **Vancogen[®]-Space Knee ATS** to the tibial component of **Vancogen[®]-Space Knee** applying a layer of Vancogen[®] bone cement between the two components.

Apply a layer of bone cement to the non-articulating surface of combined tibial and ATS components as well as the proximal surface of the tibia. Then, manually position the combined components on the proximal tibia, taking care to remove all excess of bone cement.

Reduce eventually the joint before the cement of the tibial part has cured and perform some flex-extension movements to allow the self-centering of the tibial component in relation to the femoral component.

Clean the area from any debris or other with pulsing washing systems.

With the suture and the reconstruction of the extensor apparatus completed, the knee must be stable, but not too tight, and have a joint excursion ranging from 0° and 90°.

POST-OPERATIVE TREATMENT

As a general rule, post-operative treatment is superimposable with a primary knee prosthesis, with the difference that the weight-bearing can be only partial (use of crutches).

Partial weight bearing must be assessed on an individual basis in relation to the anatomic conditions of the femur and tibia, the bone growth and the clinical conditions of the patient. In particular, one should avoid the risk that excessive weight-bearing or forced mobilisation cause the structure of the spacer to damage the biological structure.

The articular recovery objective is 90° flex, or at least the flex obtained during operation when the suture is completed.

If the surgeon deems it necessary, on the basis of the stability achieved and, especially, the condition of the extensor apparatus, temporary use of an articulated post-operative orthosis may be prescribed.

The physician shall make the patient aware of all the warnings, precautions and measures to be taken and shall provide patient with the filled implant card delivered with the device.

The follow-up may include x-rays and/or removal of fluid from the knee joint to see if the infection is healing. The follow-up may also include drug monitoring tests to make sure the level of drugs in the system is correct. After the infection is healed, a second surgery to remove the device and implant a permanent knee joint implant will be needed.

MRI SAFETY INFORMATION

Vancogen[®]-Space Knee is considered as "MR Safe" in accordance with the definition of the standard ASTM F2503, therefore it does not pose any risk in the Magnetic Resonance Environment.

IMPLANT CARD

The device is provided with an implant card to be filled by the Healthcare Institution/provider with patient data, and to be supplied to the patient. This implant card reports the information allowing the identification of the implanted device. Further information to be supplied to the patient about warnings, precautions and measures to be taken is reported on the patient leaflet and at the following website: <https://ifu.tecres.com>. The patient leaflet included in the device package must be supplied to the patient.

Instruction for Implant Card completion

The following data shall be filled by the HealthCare Institution/provider on the implant card front side:

1. Name of the patient or patient ID.
2. Date of implantation (day/month/year).
3. Name and address of the Healthcare Institution/provider.

The image shows a rectangular implant card. On the left side, there are three rows of input fields. Each row starts with a small icon: a person with a question mark, a calendar, and a person with a plus sign. To the right of each icon is a horizontal line for text. Above each line is a small circle containing a number: 1, 2, and 3. At the bottom left, there is a small icon of a computer monitor followed by the URL <https://ifu.tecres.com>. On the top left is the Tecres logo (a cube) and the text 'TECRES' and 'IMPLANT CARD'. On the right side, there is a large, stylized black and white geometric logo consisting of several overlapping triangles.

Once filled the Implant Card shall be detached from its support following the dotted lines and provided to the patient.

Explanation of overall symbols present on the implant card (front and back sides) are reported under section "SYMBOLS on Implant Card" of the present leaflet.

EXPLANATION

The device is not intended for permanent implantation and must be removed within 180 days of implantation.

Osteotomes, mallets and other revision instruments may be used to aid in the explantation procedure.

Care should be taken to assure that the wound site is thoroughly cleaned of all bone cement debris prior to implantation of a definitive prosthesis or performing an alternative surgical procedure (e.g. resection arthroplasty, fusion etc.). Failure to remove cement and/or bone debris may shorten the survival of the revision implant.

NOTE: It may happen that the device will be left indefinitely in the patient when the patient becomes unfit for further surgery or when the surgeon agrees with the patient under their responsibility that this is the best treatment option in a benefit/risk evaluation.

DISPOSAL

Vancogenx[®]-Space Knee once explanted shall be disposed as biological waste according to local waste regulations, since the device could have been contaminated with potentially infectious substances of human origin.

In case of expired shelf-life, damaged or unintentionally opened packaging, please dispose Vancogenx[®]-Space Knee according to local waste regulations for unused or expired drugs since the device contain antibiotics.

SERIOUS INCIDENT REPORTING

Any serious incident potentially attributable to Vancogenx[®]-Space Knee shall be reported to the manufacturer, Tecres S.p.A., and to the competent authority in which the patient is established.

The manufacturer can be reached at the following contacts:

- email: info@tecres.it
- tel.: +39 045 9217311

The contact information of the competent authority located in the European Member State in which the patient is established is reported on European Commission website, under section "medical devices" at the following link: https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en

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The image shows a rectangular form titled "TECRES IMPLANT CARD". On the left side, there is a logo for TECRES (ADVANCED BONE TECHNOLOGY) and three icons: a person with a question mark, a calendar showing the number 31, and a person with a plus sign. To the right of these icons are three horizontal lines for text entry, each with a circled number (1, 2, 3) above it. Below these lines is a small computer icon and the URL <https://ifu.tecres.com>. On the right side of the form, there is a large, stylized graphic of a knee joint.

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The image shows a rectangular 'IMPLANT CARD' form. On the left side, there are three rows of input fields, each preceded by a small icon: a person icon with a question mark, a calendar icon, and a person icon with a plus sign. To the right of these fields are three numbered circles (1, 2, 3) corresponding to the instructions. Below the third field is a QR code icon and the URL 'https://ifu.tecres.com'. On the right side of the card, there is a large, stylized graphic of a knee joint. The top left corner features the Tecres logo and the text 'TECRES IMPLANT CARD'.

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TECRES S.p.A.