***H. pylori* IgG**

**For the Detection of IgG
Antibodies to *H. pylori*
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 *H. pylori* IgG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the qualitative detection of IgG antibodies to *Helicobacter pylori* in human serum from symptomatic adults, as an aid in the diagnosis of *Helicobacter pylori* infection.

Catalog Number: **L2KHPG2** (200 tests),
L2KHPG6 (600 tests)

Test Code: **HGP** Color: **Dark Green**

Summary and Explanation

In 1983, Warren and Marshall described a curved bacillus associated with gastric mucosa in cases of chronic gastritis, tentatively drawing comparisons between that organism and the genus *Campylobacter*. More recently, investigators have shown a correlation between colonization with this organism (now termed *Helicobacter pylori*) and gastric and duodenal ulcers and chronic gastritis. *H. pylori* colonization is chronic in nature and appears to cause histological inflammation of the gastric mucosa. When *H. pylori* is eliminated from the gastric mucosa, the inflammation tends to decrease. If the organism recolonizes, the inflammation increases in severity and coincides with the appearance of gastrointestinal symptoms. A lack of gastrointestinal symptoms in the presence of *H. pylori* indicates colonization rather than infection. The most recent clinical studies have confirmed *H. pylori* to be the causative agent for most cases of chronic gastritis and ulcers. Evidence has been presented that *H. pylori* is also associated with gastric carcinoma.

Current procedures for detecting *H. pylori* rely on isolation of the organism from tissue obtained from endoscopic biopsy. The tissue is also tested by histology, gram staining, and/or testing for the enzyme urease, which is produced by the organism in abundance. Each of the methods, including histology and gram staining, require biopsy material taken from multiple sites. In addition to the

necessity for performing an invasive gastroscopy, each of these methods suffers from other inadequacies.⁶ The presence of *H. pylori* has also been detected with a urea breath test (using radioactive or nonradioactive isotopes) and by serological methods.¹⁴⁻¹⁶

A positive serological response to *H. pylori* antigens has been determined in individuals with duodenitis, chronic gastritis, and gastric or duodenal ulcer. Further, many people without clinical symptoms are seropositive for *H. pylori* antibodies, with the prevalence increasing with age. Therefore, while *H. pylori* serology is a sensitive method to determine colonization, differentiation between colonization and active disease is not possible.

Principle of the Procedure

IMMULITE 2000 *H. pylori* IgG is a solid-phase, chemiluminescent immunometric assay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 *H. pylori* IgG has not been tested with all possible variations of tube types. Consult the section on

Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum.

Automatic Predilution Factor: 20

Dispense required volume of diluent (L2HPGZ2) to a suitable test tube with barcode label applied.

Storage: 3 days at 2–8°C, or 6 months at –20°C.¹⁷

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²¹⁻²³

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

The bead is coated with *inactivated* *H. pylori* antigen. However, caution is advised because of the possible presence of residual organism when working with, or disposing of, the materials supplied.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

***H. pylori* IgG Bead Pack (L2HPG12)**

With barcodes. 200 beads, coated with inactivated, partially purified *H. pylori* antigen. Stable at 2–8°C until expiration date.

L2KHPG2: 1 pack. **L2KHPG6:** 3 packs.

***H. pylori* IgG Reagent Wedge (L2HPGA2)**

With barcodes. 2 reagents. 11.5 mL buffer solution. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in buffer. Stable at 2–8°C until expiration date.

L2KHPG2: 1 wedge.

L2KHPG6: 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

***H. pylori* IgG Adjustors (LHPGL, LHPGH)**

Two vials (Low and High), 4 mL each, of human serum with IgG reactive to *H. pylori*, in a buffer, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHPG2: 1 set. **L2KHPG6:** 2 sets.

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

***H. pylori* IgG Controls (LHPGC1, LHPGC2, LHPGC3)**

Three vials, (Negative, Low Positive and Positive), 2 mL each. **LHPGC1 (Negative Control):** One vial containing human

serum with IgG nonreactive to *H. pylori*, with preservative. **LHPGC2, LHPGC3 (Low Positive Control, Positive Control):** Two vials containing human serum with IgG reactive to *H. pylori*, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHGP2: 1 set. **L2KHGP6:** 2 sets.

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

For the current control ranges, please refer to the Control insert.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

***H. pylori* IgG Sample Diluent (L2HPGZ2)**

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) buffer solution, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHGP2: 1 vial. **L2KHGP6:** 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2HPGZ2: 5 labels.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LHPGCM: Tri-level *H. pylori* IgG Control Module

Also required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine

maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval: 1 week.

Quality Control Samples: Use controls supplied with kit.

Interpretation of Results

Positive: A result greater than or equal to 1.1 U/mL is considered to be "positive", and indicates that *H. pylori* IgG antibodies were detected in the sample.

Indeterminate: A result greater than or equal to 0.9 U/mL and less than 1.1 U/mL, is considered to be "indeterminate". Indeterminate samples should be retested. Samples still yielding "indeterminate" values should be examined by an alternate method, or a second sample should be collected, if possible, within a reasonable time period (e.g., one week).

Negative: A result less than 0.9 U/mL is considered to be "negative", and indicates that *H. pylori* IgG antibodies were not detected in the sample. Negative results by this test do not preclude recent primary infection.

The presence of IgG antibodies to *H. pylori* is an indication of previous exposure to the organism. A single specimen can only be used to determine the immune status of the individual.

The anti-*H. pylori* IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 *H. pylori* IgG EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

The majority of individuals exposed to *H. pylori* possess IgG antibodies to the organism. In addition, the presence of *H. pylori* antibodies is a function of age,

race, geography and clinical condition.¹⁶ Age-specific rates for *H. pylori* infection are similar for males and females.

The prevalence of *H. pylori* infection is between 30 and 40% in the US and Canada, about 20% in Australia, 70% in Europe, and between 70 and 90% in Africa, South America and Asia.¹⁸

H. pylori has also been identified as a risk factor for gastric cancer. Nearly 100% of patients with duodenal ulcer, 85% of those with gastric ulcer,¹⁹ and 50 – 100% of patients with gastric cancer are infected with *H. pylori*.²⁰

A relatively large proportion of patients who have positive levels of antibody are asymptomatic, even though they are colonized with the *H. pylori* organism. Therefore, antibody levels do not necessarily correlate with the severity of clinical symptoms.¹⁶

Limitations

For the determination of seroconversion from *nonreactive* to *reactive*, two serum samples should be drawn three to four weeks apart, during the acute and convalescent stages of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Individuals with an acute *H. pylori* infection may not exhibit any detectable IgG antibodies at the early stage of infection.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from cord blood, neonates, pediatric patients, or pretransplant patients.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an

anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in U/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: 0.4 to 8.0 U/mL.

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 0.4 U/mL.

Precision: Six samples were assayed in triplicate over the course of 6 days, two runs per day, for a total of 12 runs and 36 replicates. (See "Precision" table.)

Specificity: The assay is highly specific for *H. pylori* IgG antibodies, with no observed crossreactivity to *Campylobacter jejuni*, *Campylobacter fetus* and *Campylobacter coli* microorganisms. This conclusion was drawn from a study using the IMMULITE *H. pylori* IgG assay.

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 809 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST® vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of *H. pylori* IgG, to obtain values throughout the calibration range of the assay, and

then assayed by the IMMULITE 2000 *H. pylori* IgG procedure.

(Heparin) = 1.06 (Serum) – 0.18 U/mL
r = 0.997

(EDTA) = 1.00 (Serum) – 0.09 U/mL
r = 0.993

(SST) = 1.00 (Plain Tubes) – 0.04 U/mL
r = 0.999

Means:

2.5 U/mL (Serum)

2.5 U/mL (Heparin)

2.4 U/mL (EDTA)

2.4 U/mL (SST)

The results show that heparinized and EDTA plasma yield virtually the same results as serum in the IMMULITE 2000 *H. pylori* IgG procedure.

Method Comparisons: The assay was compared to clinical diagnosis by biopsy on 155 randomly selected retrospective specimens from patients suspected of *H. pylori* infection. A specimen was considered clinically positive if the culture test, or both CLO and histology tests, were positive. A specimen was considered clinically negative if none of the biopsy tests were positive.

IMMULITE 2000

Clinical Diagnosis	Pos	Ind	Neg	Sensitivity	Specificity
Pos	65	2	2	97.0%	100.0%
Neg	0	2	84		

Agreement: 98.7%

Central 95% Confidence Limits for Sensitivity and Specificity, respectively: 89.6% – 99.6% and 95.7% – 100%.

The assay was also compared to IMMULITE *H. pylori* IgG (LKHEQ) procedure on 197 samples.
(Concentration range: approximately 0.4 to 7.0 U/mL. See graph.) By linear regression:

(IML 2000) = 0.96 (IML) + 0.05 U/mL
r = 0.992

Means:

2.19 U/mL (IMMULITE 2000)

2.23 U/mL (IMMULITE)

References

1) Antonescu CG, Marshall BJ. *Helicobacter pylori*: A potentially curable form of peptic ulcer disease. *Gastroenterol J Club* 1990;1:3-12. 2) Barthel JS, Everett ED. Diagnosis of *campylobacter pylori* infections: the "Gold

Standard" and the alternatives. *Rev Infec Dis* 1990;12 (suppl.1):107-14. 3) DeCross AJ, Peura DA. Role of *H. pylori* in peptic ulcer disease. *Contemp Gastroenterol* 1992;5:18-28. 4) George LL, Borody TJ, Andrews P, Devine M, Moore-Jones D et al. Cure of duodenal ulcer after eradication of *H. pylori*. *Med J Australia* 1990;153. 5) Goodwin CS, Armstrong JA, Marshall BJ. *Campylobacter pylori*, gastritis, and peptic ulceration. *J Clin Pathol* 1986;39:353-65. 6) Graham DY, Klein PD. What you should know about the methods, problems, interpretations, and uses of urea breath tests. *Amer J Gastroenterol* 1991;86 (9):1118-22. 7) Kosunen TU, Sepälä K, Sarna S, Sipponen P. Diagnostic value of decreasing IgG, IgA, and IgM antibody titers after eradication of *Helicobacter pylori*. *Lancet* 1992;339:893-5. 8) Parsonnet J, Friedman GD, Vandersteen DP, Chang Y, Vogelman JH et al. *Helicobacter pylori* infection and the risk of gastric carcinoma. *New Engl J Med* 1991;325:1127-31. 9) Rauws EA, Tytgat GNJ. Cure of duodenal ulcer associated with eradication of *Helicobacter pylori*. *Lancet* 1990;335:1233-5. 10) Talley NJ, Newell DG, Ormand JE, Carpenter HA, Wilson WR et al. Serodiagnosis of *Helicobacter pylori*: Comparison of enzyme-linked immunosorbent assays. *J Clin Microbiol* 1991;29:1635-9. 11) Jones DM et al. Antibody to the gastric *Campylobacter*-like organism (*Campylobacter pylori*) — clinical correlations and distribution in the normal population. *J Med Microbiol* 1986;22:57-62. 12) Thomas JE, Whatmore AM, Barer MR, Eastham EJ, Kehoe MA. Serodiagnosis of *Helicobacter pylori* infection in childhood. *J Clin Microbiol* 1990;28:2641-6. 13) Warren JR, Marshall BJ. Unidentified curved bacillus on gastric epithelium in active chronic gastritis. *Lancet* 1983;1:1273-5. 14) Fennerty MB. A review of tests for the diagnosis of *Helicobacter pylori* infection. *Lab Med* 1998;29:561-6. 15) Lichtenstein DR. Peptic ulcer disease: clinical options for detecting and eradicating *H. pylori* infection. *Mod Med* 1994;62. 16) Dunn BE, et al. *Helicobacter pylori*. *Clin Microbiol Rev* 1997;7:20-41. 17) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:358. 18) Bonney RC. The developing market for *Helicobacter pylori* testing. *PJB Publications Ltd*, 1997:43-4. 19) Marshall BJ. *Helicobacter pylori*: the etiological agent for peptic ulcer. *Journal of the American Medical Association* 1995; 274:1064-1067. 20) Vanagunas A. The link between *Helicobacter pylori* and gastric cancer. *Infect Med* 1998;15:644-640. 21) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377-82, 387-8. 22) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards

Institute; 2005. NCCLS Document M29-A3.
23) Federal Occupational Safety and Health
Administration, Bloodborne Pathogens
Standard, 29 CFR 1910.1030.

Technical Assistance

In the United States, contact Siemens
Healthcare Diagnostics Technical
Services department. Tel: 877.229.3711.
Outside the United States, contact your
National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare
Diagnostics Products Ltd. is certified to ISO
13485.

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
*Results were consistently below 0.4 U/mL. ⁶

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
© 2009 Siemens Healthcare Diagnostics. All rights reserved.