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
Federal Budget Institute of  
Science "Central Research  
Institute for Epidemiology"  
3A Novogireevskaya Street  
Moscow 111123 Russia

**TABLE OF CONTENTS**

- 1. INTENDED USE ..... 3
- 2. PRINCIPLE OF PCR DETECTION ..... 3
- 3. CONTENT ..... 3
- 4. ADDITIONAL REQUIREMENTS. .... 4
- 5. GENERAL PRECAUTIONS..... 4
- 6. SAMPLING AND HANDLING ..... 5
- 7. WORKING CONDITIONS..... 7
- 8. PROTOCOL ..... 7
- 9. DATA ANALYSIS ..... 9
- 10. TROUBLESHOOTING..... 11
- 11. TRANSPORTATION..... 11
- 12. STABILITY AND STORAGE..... 11
- 13. SPECIFICATIONS..... 12
- 14. REFERENCES ..... 12
- 15. QUALITY CONTROL..... 13
- 16. KEY TO SYMBOLS USED ..... 13

## 1. INTENDED USE

**AmpliSens® *Leptospira*-FRT** PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of 16S RNA of pathogenic *Leptospira* genospecies in the clinical material (blood and cerebrospinal fluid), autopsy material (brain, kidney, liver, lung tissue, and mesenteric lymph nodes) and biological material (material obtained from died animals (lung, brain, and kidney tissue) and animals suffering from acute leptospirosis (blood) or *Leptospira* persisting in kidneys (urine)) using real-time hybridization-fluorescence detection of amplified products.

The results of PCR analysis are taken into account in complex diagnostics of disease.

## 2. PRINCIPLE OF PCR DETECTION

Detection of 16S RNA of pathogenic *Leptospira* genospecies by the polymerase chain reaction (PCR) is based on the amplification of the pathogen genome specific region using specific primers. In the real-time PCR, the amplified product is detected with the use of fluorescent dyes. These dyes are linked to oligonucleotide probes, which bind specifically to the amplified product during thermocycling. The real-time monitoring of the fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

**AmpliSens® *Leptospira*-FRT** PCR kit is a qualitative test that contains the Internal Control (Internal Control STI-87-rec (IC)). It must be used in the isolation procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

**AmpliSens® *Leptospira*-FRT** PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by separation of nucleotides and Taq-polymerase using chemically modified polymerase (TaqF). Chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

## 3. CONTENT

**AmpliSens® *Leptospira*-FRT** PCR kit is produced in 1 form:

**AmpliSens® *Leptospira*-FRT** PCR kit variant FRT **REF** R-B49(RG,iQ)-CE.

AmpliSens® *Leptospira*-FRT PCR kit variant FRT includes:

<i>Reagent</i>	<i>Description</i>	<i>Volume, ml</i>	<i>Quantity</i>
RT-G-mix-2	colorless clear liquid	0.01	2 tubes
RT-PCR-mix-1-FRT <i>Leptospira</i>	colorless clear liquid	0.6	1 tube
RT-PCR-mix-2-FEP/FRT	colorless clear liquid	0.3	1 tube
Polymerase (TaqF)	colorless clear liquid	0.03	1 tube
TM-Revertase (MMIv)	colorless clear liquid	0.015	1 tube
Positive Control cDNA <i>Leptospira</i> (C+ <i>Leptospira</i> )	colorless clear liquid	0.1	1 tube
RNA-eluent	colorless clear liquid	0.07	2 tubes
Negative Control (C-)*	colorless clear liquid	1.6	2 tubes
Positive Control <i>Leptospira</i> -rec	colorless clear liquid	0.03	5 tubes
Internal Control STI-87-rec (IC) **	colorless clear liquid	0.12	5 tubes

\* must be used in the extraction procedure as Negative Control of Extraction.

\*\* add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture (see RIBO-sorb **REF** K2-1-Et-50-CE, RIBO-zol-C **REF** K2-13-50-CE, RIBO-prep **REF** K2-9-Et-50-CE protocols).

AmpliSens® *Leptospira*-FRT PCR kit is intended for 60 reactions, including controls.

#### 4. ADDITIONAL REQUIREMENTS

- RNA extraction kit.
- Disposable powder-free gloves and a laboratory coat.
- Pipettes (adjustable).
- Sterile RNase-free pipette tips with aerosol filters (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with a rotor for 2-ml reaction tubes.
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia); iCycler iQ5 (Bio-Rad, USA), or equivalent).
- Disposable polypropylene PCR tubes (0.2-ml).
- Refrigerator for 2–8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.

- Reservoir for used tips.

## 5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use the PCR kit if the internal packaging was damaged or its appearance was changed.
- Do not use the PCR kit if the transportation and storage conditions according to the Instruction Manual were not observed.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid samples and reagents contact with the skin, eyes, and mucous membranes. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice immediately.
- Safety Data Sheets (SDS) are available on request.
- The PCR kit is intended for single use for PCR analysis of specified number of samples (see the section “Content”).
- The PCR kit is ready for use in accordance with the Instruction Manual. Use the PCR kit strictly for intended purpose.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment

and reagents in the area where the previous step was performed.

Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

## 6. SAMPLING AND HANDLING

Obtaining samples of biological materials for PCR-analysis, transportation and storage are described in manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

**AmpliSens<sup>®</sup> Leptospira-FRT** PCR kit is intended for analysis of RNA extracted with RNA extraction kits from the following material:

### Human material

- clinical material (blood, cerebrospinal fluid);
- autopsy material (brain, kidney, liver, and lung tissue and mesenteric lymph nodes).

### Animal material (biological material)

- urine;
- blood;
- brain, kidney, and lung tissue.

The material can be stored at 2-8 °C for 1 day. The autopsy material can be stored at ≤ -16 °C for 1 week, at ≤ -68 °C for a long time.

### Sampling and pretreatment

6.1. ~~Bsp~~ . Whole blood is taken in the morning after overnight fasting into the tube with 6 % EDTA solution in the ratio 1 : 20. The closed tube with whole peripheral blood should be rotated several times. The tube with blood should be centrifuged at 1,000 g for 10 min to obtain blood plasma (if the blood was stored at 2–8 °C more than 1 hour, it should be mixed carefully by inverting the tube). Transfer 1 ml of plasma into two tubes. Two tubes with 1.0 ml of plasma should be centrifuged at 13,000 rpm for 10 min to concentrate bacterial cells. Then, 900 µl of the supernatant plasma should be removed with a filter tip into the container with disinfectant. The pellet and 100 µl of the supernatant are tested for the presence of ~~b~~ 16S RNA. The second pellet prepared in the same way should be stored at ≤ -16 °C for repeated extraction (if any technological procedure is performed incorrectly). The pellet obtained from blood plasma can be stored at ≤ -16 °C for 1 week or at ≤ -68 °C for a long time.

When cerebrospinal fluid is analyzed, the pellet is obtained by the same procedure by centrifugation at 13,000 rpm for 10 min. The pellet and 100 µl of the supernatant are analyzed.

6.2. ~~U~~ for analyses is taken into a sterile container. If there is no chance to test material within 24 h after sampling, urine is transferred to a centrifuge tube or an Eppendorf tube. The contents of the tube is mixed with glycerol (~10 % v/v) and frozen. It can be stored at  $\leq -16$  °C for 1 week or at  $\leq -68$  °C for a long time.

If a cooling centrifuge (4 °C) with a speed of 9,000–10,000 g intended for 30-ml tubes is available, the following sample preparation procedure is used. The sample is centrifuged at 9,000–10,000 g for 10 min, the supernatant is transferred to a container with disinfectant. Leave ~1 ml of the supernatant over the pellet in the tube and resuspend it. Transfer the suspension to a new tube and concentrate it by centrifugation at 13 000,rpm for 10 min. Then, 900  $\mu$ l of the supernatant is transferred to the container with disinfectant, and the pellet and 100  $\mu$ l of the supernatant is used for RNA isolation. In case of large quantities of salts and mucus, 100  $\mu$ l of the supernatant and the upper layer of cells should be carefully taken from the salt pellet and transferred into a new tube for RNA isolation.

If you have no centrifuge for 30-ml tubes and a speed of 9,000–10,000 g, bacteria are concentrated from 1 ml of urine as described above using 1.5-ml tubes and a microcentrifuge for Eppendorf tubes. The remaining urine should be decontaminated in a disinfectant.

6.3. ~~U~~ ~~U~~ ~~U~~ ~~U~~ is to be homogenized in sterile porcelain mortars with pestles. Then, 10 % suspension in sterile saline or phosphate buffer is prepared; 30  $\mu$ l of the suspension is taken for RNA extraction.

## 7. WORKING CONDITIONS

AmpliSens® *Leptospira*-FRT PCR kit should be used at 18–25 °C.

## 8. PROTOCOL

### 8.1. RNA extraction

It is recommended to use the following nucleic acid extraction kits:

- RIBO-prep, **REF** K2-9-Et-50-CE.
- RIBO-zol-C **REF** K2-13-50-CE and RIBO-sorb **REF** K2-1-Et-50-CE.

Extract the RNA according to the manufacturer's protocol.

**If using the RIBO-prep kit** extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- If RNA is extracted from the tissues homogenates. Add **10  $\mu$ l** of **Internal Control STI-87-rec (IC)** to each tube and then add **300  $\mu$ l** of **Solution for Lysis**. Label test tubes. Then, add **50  $\mu$ l** of test suspensions to the tubes

with **Solution for Lysis** and **Internal Control STI-87-rec (IC)**.

- If RNA is extracted from the cell pellets after **β** concentration, **300 µl** of **Solution for Lysis** and **10 µl** of **Internal Control STI-87-rec (IC)** should be added directly into the tubes with the pellets.
- Add only **10 µl** of **Internal Control STI-87-rec (IC)** and **300 µl** of **Solution for Lysis** into the tube for the Negative Control of extraction (C–).
- Add **300 µl** of **Solution for Lysis**, and then add **90 µl** of **Negative Control (C–)**, **10 µl** of **Internal Control STI-87-rec (IC)** and **10 µl** of **Positive Control *Leptospira-rec*** into the tube for the Positive Control of extraction (PCE).
- **40 µl** of **RNA-buffer** is to be used.

**If using the RIBO-zol-C and RIBO-sorb kit** extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- Add **300 µl** of **Solution D** and then **10 µl** of **Internal Control STI-87-rec (IC)** into the tubes with the test pellets of blood, cerebrospinal fluid (when extracting from cerebrospinal fluid add 50 µl of Negative control (C–) into the tube after lysis) or urine, or with tissues homogenate.
- In case of analyzing the suspensions of viscera tissues, at the 1<sup>st</sup> stage **300 µl** of **Solution D** and **10 µl** of **Internal Control STI-87-rec (IC)** can be added first into empty tubes and then **30 µl** of test sample can be added.
- Add **300 µl** of **Solution D**, **10 µl** of **Internal Control STI-87-rec (IC)** and **100 µl** of **Negative control (C–)** into the tube for the Negative Control of extraction (C–).
- Add **300 µl** of **Solution D**, and then add **90 µl** of **Negative Control (C–)**, **10 µl** of **Internal Control STI-87-rec (IC)** and **10 µl** of **Positive Control *Leptospira-rec*** into the tube for the Positive Control of extraction (PCE).
- After the first stage of RNA extraction with **RIBO-zol-C** kit the top phase of obtained samples (**400 µl** when extracting RNA from cerebrospinal fluid, urine pellets and tissues homogenates, **200 µl** – from blood pellets, 450 µl – from C– and PCE samples) should be transferred into new 1.5-ml tubes with **400 µl** of **Lysis Solution** (from RIBO-sorb kit).
- **40 µl** of **RNA-buffer** is to be used.

## 8.2. Preparing reverse transcription and PCR

The total reaction volume is **25 µl**, the volume of DNA sample is **10 µl**.

Only RNase-free, DNase-free disposable plastic consumables must be used when working with RNA.

### 8.2.1 Preparing tubes for reverse transcription and PCR

1. Prepare the required number of tubes for amplification of cDNA obtained from clinical and control samples. The type of tubes depends on the PCR instrument used for analysis. For carrying out N reactions with 2 controls, N+2 tubes are required.
2. Prepare the **reaction mixture**, calculating per **one** reaction:

- **10 µl** of RT-PCR-mix-1-FRT *Leptospira*
- **5 µl** of RT-PCR-mix-2-FEP/FRT
- **0.5 µl** of polymerase (TaqF)
- **0.25 µl** of TM-Revertase (MMIv)
- **0.25 µl** of RT-G-mix-2

3. Transfer **15 µl** of the prepared mixture to each tube. Discard the unused mixture.
4. Using filter tips add **10 µl** of **RNA** samples obtained at the RNA extraction stage into prepared tubes.

Avoid transferring of sorbent together with the RNA samples into the reaction mixture for RT-PCR.

5. Carry out the control amplification reactions:

- NCA** – Add **10 µl** of **RNA-eluent** to the tube labeled NCA (Negative Control of Amplification).
- C+** – Add **10 µl** of **Positive Control cDNA *Leptospira* (C+*Leptospira*)** to the tube labeled C+ (Positive Control of Amplification).
- C–** – Add **10 µl** of **the sample extracted from the Negative Control of Extraction sample** to the tube labeled C– (Negative control of Extraction).
- PCE** – Add **10 µl** of **the sample extracted from the Positive control of Extraction sample** to the tube labeled PCE (Positive control of Extraction).

Amplification is to be carried immediately after mixing the reaction mixture, RNA-sample and controls. The time period between addition of RNA-samples into the reaction mixture and amplification starting is to be not more than 10-15 min.

### 8.2.2. Amplification

1. Create a temperature profile on your instrument as follows:

Table 1

**Amplification program for *Leptospira* cDNA**

Step	Rotor-type instruments <sup>1</sup>			Plate-type instruments <sup>2</sup>		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	50 °C	30 min	1	50 °C	30 min	1
2	95 °C	15 min	1	95 °C	15 min	1
3	95 °C	20 s	10	95 °C	20 s	10
	65 °C	50 s		65 °C	50 s	
	72 °C	20 s				
4	95 °C	20 s	38	95 °C	20 s	40
	61 °C	50 s Fluorescence acquiring		61 °C	50 s Fluorescence acquiring	
	72 °C	20 s		65 °C	20 s	

<sup>1</sup> For example, Rotor-Gene 3000/Rotor-Gene 6000 (Corbett Research, Australia).

<sup>2</sup> For example, iCycler iQ5 (Bio-Rad, USA).

Fluorescent signal is detected in the channels for the FAM and JOE fluorophores.

2. Adjust the fluorescence channel sensitivity according to the Guidelines [2].
3. Insert tubes into the reaction module of the device.
4. Run the amplification program with fluorescence detection.
5. Analyze results after the amplification program is completed.

## 9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in two channels:

- The signal of the IC cDNA amplification product is detected in the channel for the FAM fluorophore.
- The signal of the  $\beta$  cDNA amplification product is detected in the channel for the JOE fluorophore.

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a  $C$  value of the cDNA sample in the corresponding column of the results grid.

Principle of interpretation is the following:

- $\beta$  cDNA is **detected** if the  $C$  value determined in the results grid in the channel for the JOE fluorophore is less than the boundary  $C$  value specified in the Guidelines.
- $\beta$  cDNA is **not detected** in a sample if the  $C$  value is not determined (absent) in the channel for the JOE fluorophore, whereas the  $C$  value determined in the channel for the FAM fluorophore is less than the boundary  $C$  value specified in the Guidelines.
- The result is **invalid** if the  $C$  value is not determined (absent) in the channel for the JOE fluorophores whereas the  $C$  value in the channel for the FAM fluorophore is greater than the specified boundary  $C$  value. In such cases, the PCR analysis of this sample should be repeated starting from the RNA extraction stage.
- The result is **equivocal** if the  $C$  value determined in the channel for the JOE fluorophore is greater than the boundary  $C$  value specified in the Guidelines, whereas the  $C$  value determined in the channel for the FAM fluorophore is less than the boundary  $C$  value specified in the Guidelines. In such cases, the PCR analysis of this sample should be repeated two times starting from the RNA extraction stage.

Boundary  $C$  values are specified in the Guidelines [2]

**The result of the analysis is considered reliable only if the results obtained for the**

Positive and Negative Controls of amplification and extraction are correct (see Table 2).

Table 2

**Results for controls**

Control	Stage for control	Ct value in the channel for the fluorophore	
		FAM	JOE
PCE	DNA extraction	Present	< boundary value
C-	DNA extraction	Present	Absent
NCA	PCR	Absent	Absent
C+	PCR	Absent	< boundary value

**10. TROUBLESHOOTING**

Results of analysis are not being registered in the following cases:

1. If the C value is determined for the Negative Control of extraction (C-) in the channel for the JOE fluorophore and/or for the Negative Control of amplification (NCA) in the channels for the FAM and JOE fluorophores in the results grid, it indicates contamination of reagents or samples. In such cases, the results of analysis are considered to be irrelevant. Analysis should be repeated and measures to detect and eliminate the source of contamination should be taken.
2. If no signal is detected for the Negative Control of extraction (C-) in the channel for the FAM fluorophore and/or for the Positive Control of extraction (PCE) in the channels for the FAM and JOE fluorophores, the results of analysis are considered invalid. Analysis of all samples should be repeated starting from the extraction stage.
3. If no signal is detected for Positive Control of amplification (C+) in the channel for the JOE fluorophore, the results of analysis are considered invalid. Analysis of all samples should be repeated starting from the RT-PCR stage.

If you have any further questions or if you encounter problems, please contact our Authorized representative in the European Community.

**11. TRANSPORTATION**

**AmpliSens® Leptospira-FRT** PCR kit should be transported at 2–8 °C for no longer than 5 days.

**12. STABILITY AND STORAGE**

All components of the **AmpliSens® Leptospira-FRT** PCR kit are to be stored at 2–8 °C when not in use (except for RT-G-mix-2, RT-PCR-mix-1-FRT ~~2~~, RT-PCR-mix-2-FEP/FRT, Polymerase (TaqF), and TM-Revertase (MMIv)). All components of the

**AmpliSens® Leptospira-FRT** PCR kit are stable until the expiry date stated on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

RT-G-mix-2, RT-PCR-mix-1-FRT ~~Ⓢ~~, RT-PCR-mix-2-FEP/FRT, Polymerase (TaqF), and TM-Revertase (MMIv) are to be stored at the temperature from minus 24 to minus 16 °C.

RT-PCR-mix-1-FRT ~~Ⓢ~~ is to be stored away from light

## 13. SPECIFICATIONS

### 13.1. Sensitivity

Analytical Sensitivity of **AmpliSens® Leptospira-FRT** PCR kit is not less than  $5 \times 10^3$  copies per 1 ml of sample (copies/ml).

The claimed analytical features of **AmpliSens® Leptospira-FRT** PCR kit are guaranteed only when additional reagents kits **RIBO-sorb** and **RIBO-zol-C** or **RIBO-prep** (manufactured by Federal Budget Institute of Science “Central Research Institute for Epidemiology”) are used.

### 13.2. Specificity

The analytical specificity of **AmpliSens® Leptospira-FRT** PCR kit is ensured by the selection of specific primers and probes as well as stringent reaction conditions. The primers and probes have been checked for possible homologies to all sequences published in gene banks by sequence comparison analysis.

The clinical specificity of **AmpliSens® Leptospira-FRT** PCR kit was confirmed in laboratory clinical trials.

## 14. REFERENCES

1. Handbook “Sampling, Transportation, and Storage of Clinical Material for PCR diagnostics”, developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology” of Federal Service for Surveillance on Consumers’ Rights Protection and Human Well-Being, Moscow, 2010.
2. Guidelines to the **AmpliSens® Leptospira-FRT** PCR kit for qualitative detection of 16S RNA of pathogenic ~~Ⓢ~~ genospecies in the clinical material, autopsy material and biological material by the polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science “Central Research Institute for Epidemiology”.

## 15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science “Central Research Institute for Epidemiology” ISO 13485-Certified Quality Management System, each lot of the **AmpliSens® *Leptospira-FRT*** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

## 16. KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Sufficient for
	<i>Itb</i> diagnostic medical device		Expiration Date
	Version		Consult instructions for use
	Temperature limitation		Keep away from sunlight
	Manufacturer	<b>NCA</b>	Negative control of amplification
	Date of manufacture	<b>C-</b>	Negative control of extraction
	Authorised representative in the European Community	<b>C+</b>	Positive control of amplification
		<b>PCE</b>	Positive control of extraction
		<b>IC</b>	Internal control

### List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
13.12.10	Cover page	The phrase “For Professional Use Only” was added
	Intended use	The phrase “The results of PCR analysis are taken into account in complex diagnostics of disease” was added
	Content	New sections “Working Conditions” and “Transportation” were added
		The “Explanation of Symbols” section was renamed to “Key to Symbols Used”
	Stability and Storage	The information about the shelf life of reagents before and after the first use was added
		Information that RT-PCR-mix-1-FRT <del>β</del> is kept away from light was added
	Key to Symbols Used	The explanation of symbols was corrected
	Footer	Reference number was changed from R-B49(RG)-CE to R-B49(RG,iQ)-CE
Sampling and Handling	Duration of the urine sample storage was changed	
	Information about storage and disposal of urine samples was added	
01.07.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science “Central Research Institute for Epidemiology”
16.10.15 ChA	3. Content	Quantity of Negative Control (C-) tubes was changed from 1 to 2
	8.1. RNA Isolation	For RIBO-prep the procedure of positive control of extraction preparation was changed: 90 µl of Negative Control (C-) is to be added additionally
17.01.17 ME	Text	Corrections according to the template
	1. Intended use	Types of biological material was specified
	8.1. RNA extraction	The sections were rewritten
	8.2.2. Amplification	
	9. Data analysis	
14. References	The reference to Guidelines was added	

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# CERTIFICATE

No. 18 0028 SJ

We confirm on the basis of a performed audit that company

**Federal Budget Institute of Science “Central  
Research Institute for Epidemiology“**

3a, Novogireevskaya str., 111123 Moscow, Russian Federation

has implemented and documented a functional quality management system  
in compliance with the requirements of the standard

**EN ISO 13485:2016**

Covering the following activities:

Design and development, manufacturing and final control of *in vitro* diagnostic  
medical devices

The Certificate is issued on the basis of the results mentioned in Audit Report No. 233404641/2018.  
The Certificate validity is conditioned by positive results of surveillance audits, which the certified  
company committed to undergo.

During use of the Certificate the Certificate Holder undertakes to follow the Rules of Use of the Certificate. This  
document is publicly available on [www.itczlin.cz](http://www.itczlin.cz)

Date of Issue: 07. 05. 2018

Valid until: 07. 05. 2021

Date of the first certification awarding: 20. 05. 2015

  
Ing. Pavel Vaněk

Head of Certification Body

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

## EC Certificate - Full Quality Assurance System No. 11 0040 QS/NB

The quality system of manufacturer

### **Federal Budget Institute of Science “Central Research Institute for Epidemiology”**

3A Novogireevskaya Street, Moscow 111123, Russia

has been certified as meeting the requirements of

### **Directive 98/79/EC**

**on in vitro diagnostic medical devices, Annex IV excluding (4, 6)**

for the following product category(ies):

### **AmpliSens® PCR kits**

The Notified Body No. 1023 declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with IVDD Annex IV. This quality assurance system conforms to the requirements of this Directive and is subjected to periodical surveillance. For placing on the market of List A devices covered by this certificate, an EC Design-Examination Certificate according to Annex IV (Section 4) is required.

**Valid from:** 2018-10-31  
**Valid until:** 2021-06-16  
**First Issued:** 2011-01-24  
**Revision:** h

Date: 2018-10-31

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

## Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**

**3A Novogireevskaya Street, Moscow 111123, Russia**

### Product(s):

**Name:** **AmpliSens® Rubella virus-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT-50 F

**Classification:** List B

**GMDN:** 30793

**Name:** **AmpliSens® Toxoplasma gondii-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT-50 F

**Classification:** List B

**GMDN:** 52428

**Name:** **AmpliSens® CMV-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5-ml tubes), variant FEP (0.2-ml tubes)

**Classification:** List B

**GMDN:** 30798

**Name:** **AmpliSens® CMV-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 30798

Date: 2018-10-31

Revision: h

Mgr. Jiří Heš

Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name: AmpliSens® HSV / CMV-MULTIPRIME-FEP  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5-ml tubes), variant FEP (0.2-ml tubes)

**Classification:** List B

**GMDN:** 61348

**Name: AmpliSens® HSV / CMV-MULTIPRIME-FRT  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 61348

**Name: AmpliSens® CMV-screen/monitor-FRT PCR  
kit**

**Trade name(s):** -

**Model(s):** variant FRT-100 F

**Classification:** List B

**GMDN:** 30798

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name: AmpliSens® EBV / CMV / HHV6-screen-FRT  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT-100 F

**Classification:** List B

**GMDN:** 61348

**Name: AmpliSens® Chlamydia trachomatis-FEP  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5-ml tubes), variant FEP (0.2-ml tubes)

**Classification:** List B

**GMDN:** 30677


**Name: AmpliSens® Chlamydia trachomatis-FRT  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 30677

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name: AmpliSens® C.trachomatis / Ureaplasma-  
MULTIPRIME-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 50409

**Name: AmpliSens® C.trachomatis / M.genitalium  
MULTIPRIME-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 50409


**Name: AmpliSens® C.trachomatis / Ureaplasma /  
M.genitalium-MULTIPRIME-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 50409

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute  
for Epidemiology”**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name:** **AmpliSens® C.trachomatis / Ureaplasma /  
M.genitalium-MULTIPRIME-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 50409

**Name:** **AmpliSens® C.trachomatis / Ureaplasma / M.  
hominis-MULTIPRIME-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 50409


**Name:** **AmpliSens® C.trachomatis / Ureaplasma /  
M.hominis-MULTIPRIME-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 50409

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name:** **AmpliSens® C.trachomatis / Ureaplasma /  
M.genitalium / M.hominis-MULTIPRIME-FRT  
PCR kit**

**Trade name(s):** -  
**Model(s):** variant FRT, variant FRT-100 F  
**Classification:** List B  
**GMDN:** 50409

**Name:** **AmpliSens® N.gonorrhoeae / C.trachomatis /  
M.genitalium / T.vaginalis-MULTIPRIME-FRT  
PCR kit**

**Trade name(s):** -  
**Model(s):** variant FRT, variant FRT-100 F  
**Classification:** List B  
**GMDN:** 50409

**Name:** **AmpliSens® N.gonorrhoeae / C.trachomatis /  
M.genitalium-MULTIPRIME-FEP PCR kit**

**Trade name(s):** -  
**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)  
**Classification:** List B  
**GMDN:** 50409

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name:** **AmpliSens® N.gonorrhoeae / C.trachomatis /  
M.genitalium-MULTIPRIME-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 50409

**Name:** **AmpliSens® Genoscreen HLA B\*5701-FRT  
PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT

**Classification:** List B

**GMDN:** 56403

**Name:** **AmpliSens® Mycoplasma pneumoniae /  
Chlamydomphila pneumoniae-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 58957

Date: 2018-10-31

Revision: h

Mgr. Jiří Heš

Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Name:** **AmpliSens<sup>®</sup> Mycoplasma pneumoniae /  
Chlamydophila pneumoniae-FRT PCR kit**

**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 58957

**Name:** **AmpliSens<sup>®</sup> T.vaginalis / N.gonorrhoeae /  
C.trachomatis-MULTIPRIME-FEP PCR kit**

**Trade name(s):** -

**Model(s):** variant FEP (0.5 ml tubes), variant FEP (0.2 ml tubes)

**Classification:** List B

**GMDN:** 61144

**Name:** **AmpliSens<sup>®</sup> T.vaginalis / N.gonorrhoeae /  
C.trachomatis-MULTIPRIME-FRT PCR kit**


**Trade name(s):** -

**Model(s):** variant FRT, variant FRT-100 F

**Classification:** List B

**GMDN:** 61144

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

Notified Body 1023  
**INSTITUTE FOR TESTING AND CERTIFICATION, Inc.**,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute  
for Epidemiology”**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Facility(ies):**

Federal Budget Institute of Science “Central Research Institute for Epidemiology”  
3A Novogireevskaya Street, Moscow 111123, Russia

Date: 2018-10-31  
Revision: h

**Mgr. Jiří Heš**  
Representative of the Notified Body No. 1023

Notified Body 1023  
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,  
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

**Annex to EC Certificate No. 11 0040 QS/NB**  
issued for manufacturer:

**Federal Budget Institute of Science "Central Research Institute  
for Epidemiology"**  
**3A Novogireevskaya Street, Moscow 111123, Russia**

**Certificate History:**

Revision	Date	Reference Number	Action
	2011-01-24	813600111	Certification
a	2011-07-21	813600161	Change of manufacturer name
b	2012-02-13	343601304	Product scope extension
c	2014-05-13	343602568	Product scope extension
d	2016-01-15	813600504a	Prolongation of certificate validity
e	2016-06-17	813600504	Re-certification process
f	2016-08-29	343603690	Change of manufacturer facility address
g	2017-11-30	343603888	Changes of product compositions, packaging and quality system documentation
h	2018-10-31	813600754	Change of product labelling, shelf life extension and quality system documentation

Date: 2018-10-31  
Revision: h

Mgr. Jiří Heš  
Representative of the Notified Body No. 1023

FEDERAL SERVICE FOR SUPERVISION OF CONSUMER RIGHTS PROTECTION AND HUMAN WELFARE

FEDERAL BUDGET INSTITUTE OF SCIENCE  
«CENTRAL RESEARCH INSTITUTE FOR EPIDEMIOLOGY»

111123, Moscow, 3A Novogireevskaya street, Tel.: +7 495 974 96 42, Fax: +7 495 305 54 23,  
e-mail: obtk@pcr.ru

**EC DECLARATION OF CONFORMITY**

Directive 98/79/EC of the European Parliament and of the Council of 27<sup>th</sup> of October 1998 on  
*In Vitro* Diagnostic Medical Devices

Federal Budget Institute of Science “Central Research Institute for Epidemiology” hereby under own responsibility declares that the products covered by the declaration conform with Essential Requirements listed in Annex I of EC Directive 98/79/EC (IVD Directive). Supporting documentation is retained under the premises of the manufacturer.

The quality management system meets the requirements of the standard EN ISO 13485 “Medical devices – Quality management systems – Requirements for regulatory purposes” and is certified by Institute for testing and certification, Inc. (certificate No. 15 0125 SJ/a, valid until 2018.05.20).

<b>Manufacturer:</b>	Federal Budget Institute of Science “Central Research Institute for Epidemiology”
<b>Authorized Representative:</b>	Ecoli s.r.o. Studenohorska 12 841 03 Bratislava 47 Slovak Republic tel.: (+421) 02/64 789 336 fax: (+421) 02/64 789 040 email: ecoli@ecoli.sk
<b>Product Name:</b>	Annex for this Declaration
<b>Description:</b>	Reagent kits for qualitative detection and quantification of DNA (RNA) of different infectious agents or HLA B*5701 DNA in human specimens
<b>Classification:</b>	Article 9, paragraph 3 of EC Council Directive 98/79/EC on <i>in Vitro</i> Diagnostic Devices Annex II List B IVDs (According to EC Declaration of Conformity List)
<b>Conformity Assessment Route:</b>	Annex IV (IVDD) Full QA System
<b>Notified Body:</b>	Institute for testing and certification, Inc. třída Tomáše Bati 299 Louky, 763 02 Zlin, Czech Republic e-mail: itc@itczlin.cz Notified Body No. 1023
<b>EC Certificate:</b>	No. 11 0040 QS/NB revision g, valid until 2021.06.16
<b>Place, Date of Issue:</b>	Zlin, Czech Republic, 2017.11.30

Signed \_\_\_\_\_

Valid from 2018.01.26  
Valid until 2021.06.16

Full name: Vasiliy G. Akimkin  
Title: Interim Director

№№	Description	Product Code
1.	AmpliSens® <i>Rubella virus</i> -FRT PCR kit	R-V24-S(RG,iQ,Mx)-CE
2.	AmpliSens® <i>Toxoplasma gondii</i> -FRT PCR kit	R-P1(RG,iQ,Mx)-CE
3.	AmpliSens® CMV-FEP PCR kit	V7-100-R0,2-FEP-CE V7-100-R0,5-FEP-CE
4.	AmpliSens® CMV-FRT PCR kit	R-V7(RG)-CE R-V7-F(RG,iQ)-CE
5.	AmpliSens® HSV / CMV-MULTIPRIME-FEP PCR kit	V60-100-R0,5-FEP-CE V60-100-R0,2-FEP-CE
6.	AmpliSens® HSV / CMV-MULTIPRIME-FRT PCR kit	R-V60(RG)-CE R-V60-F(RG,iQ)-CE
7.	AmpliSens® CMV-screen/monitor-FRT PCR kit	R-V7-100-S(RG,iQ,Mx)-CE
8.	AmpliSens® EBV / CMV / HHV6-screen-FRT PCR kit	R-V48(RG,iQ,Mx)-CE
9.	AmpliSens® <i>Chlamydia trachomatis</i> -FEP PCR kit	B1-100-R0,2-FEP-CE B1-100-R0,5-FEP-CE
10.	AmpliSens® <i>Chlamydia trachomatis</i> -FRT PCR kit	R-B1(RG)-CE R-B1-F(RG,iQ)-CE
11.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> -MULTIPRIME-FEP PCR kit	B47-100-R0,2-FEP-CE B47-100-R0,5-FEP-CE
12.	AmpliSens® <i>C.trachomatis</i> / <i>M.genitalium</i> -MULTIPRIME-FEP PCR kit	B66-100-R0,5-FEP-CE B66-100-R0,2-FEP-CE
13.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> / <i>M.genitalium</i> -MULTIPRIME-FEP PCR kit	B46-100-R0,2-FEP-CE B46-100-R0,5-FEP-CE
14.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> / <i>M.genitalium</i> -MULTIPRIME-FRT PCR kit	R-B46(RG)-CE R-B46-F(RG,iQ)-CE
15.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> / <i>M.hominis</i> -MULTIPRIME-FEP PCR kit	B43-100-R0,2-FEP-CE B43-100-R0,5-FEP-CE
16.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> / <i>M.hominis</i> -MULTIPRIME-FRT PCR kit	R-B43(RG)-CE R-B43-F(RG,iQ)-CE
17.	AmpliSens® <i>C.trachomatis</i> / <i>Ureaplasma</i> / <i>M.genitalium</i> / <i>M.hominis</i> -MULTIPRIME-FRT PCR kit	R-B60(RG)-CE R-B60-F(RG)-CE
18.	AmpliSens® <i>N.gonorrhoeae</i> / <i>C.trachomatis</i> / <i>M.genitalium</i> / <i>T.vaginalis</i> -MULTIPRIME-FRT PCR kit	R-B61(RG)-CE R-B61-F(RG)-CE
19.	AmpliSens® <i>N.gonorrhoeae</i> / <i>C.trachomatis</i> / <i>M.genitalium</i> -MULTIPRIME-FEP PCR kit	B67-100-R0,5-FEP-CE B67-100-R0,2-FEP-CE
20.	AmpliSens® <i>N.gonorrhoeae</i> / <i>C.trachomatis</i> / <i>M.genitalium</i> -MULTIPRIME-FRT PCR kit	R-B67(RG)-CE R-B67-F(RG,iQ)-CE
21.	AmpliSens® Genoscreen HLA B*5701-FRT PCR kit	R-O2(RG,iQ)-CE
22.	AmpliSens® <i>Mycoplasma pneumoniae</i> / <i>Chlamydophila pneumoniae</i> -FEP PCR kit	B42-50-R0,5-FEP-CE B42-50-R0,2-FEP-CE B42-Mod-50-R0,2-FEP-CE; B42-Mod-50-R0,5-FEP-CE
23.	AmpliSens® <i>Mycoplasma pneumoniae</i> / <i>Chlamydophila pneumoniae</i> -FRT PCR kit	R-B42-4x(RG)-CE R-B42-100-F-CE
24.	AmpliSens® <i>T.vaginalis</i> / <i>N.gonorrhoeae</i> / <i>C.trachomatis</i> -MULTIPRIME-FEP PCR kit	B83-100-R0,5-FEP-CE B83-100-R0,2-FEP-CE
25.	AmpliSens® <i>T.vaginalis</i> / <i>N.gonorrhoeae</i> / <i>C.trachomatis</i> -MULTIPRIME-FRT PCR kit	R-B83(RG)-CE R-B83-F(RG,iQ)-CE

### EC DECLARATION OF CONFORMITY

Directive 98/79/EC of the European Parliament and of the Council of 27<sup>th</sup> of October 1998 on  
In Vitro Diagnostic Medical Devices

Federal Budget Institute of Science "Central Research Institute for Epidemiology" hereby under own responsibility declares that the products covered by the declaration conform with Essential Requirements listed in Annex I of EC Directive 98/79/EC (IVD Directive). Supporting documentation is retained under the premises of the manufacturer.

The quality management system meets the requirements of the standard EN ISO 13485 "Medical devices – Quality management systems – Requirements for regulatory purposes" and is certified by Institute for testing and certification, Inc. (certificate No. 18 0028 SJ, valid until 2021.05.07).

<b>Manufacturer:</b>	Federal Budget Institute of Science "Central Research Institute for Epidemiology"
<b>Authorized Representative:</b>	Ecoli s.r.o. Studenohorska 12 841 03 Bratislava 47 Slovak Republic Tel.: +421 264 789 336 Fax: +421 264 789 040
<b>Product Name:</b>	Annex for this Declaration
<b>Description:</b>	Reagent kits for qualitative detection and quantification of DNA (RNA) of different infectious agents
<b>Classification:</b>	Article 9, paragraph 1 of EC Council Directive 98/79/EC on <i>in Vitro</i> Diagnostic Devices
<b>Conformity Assessment Route:</b>	Annex III (IVDD)

Signed \_\_\_\_\_

Valid from 2018.10.05

Full name: Vasily G. Akimkin  
Title: Director

№№	Description	Product Code
1.	<b>Hemolytic</b>	137-CE
2.	<b>Mucolysin</b> reagent for sputum preliminary treatment	180-CE
3.	<b>RNA-medium</b>	981-CE
4.	<b>Transport Medium for Storage and Transportation of Respiratory Swabs</b>	957-CE 958-CE 959-CE
5.	<b>Transport Medium with Mucolytic Agent</b>	952-CE 953-CE
6.	<b>Transport medium TM-EDEM</b>	1533-CE
7.	<b>Citolizin</b> nucleic acid extraction kit	K1-3-100-CE
8.	<b>AmpliSens® DNA-sorb-D</b> nucleic acid extraction kit	K8-2331-100-CE
9.	<b>AmpliSens® MAGNO-sorb-URO</b> nucleic acid extraction kit	K4-2181-100-CE
10.	<b>DNA-sorb-AM</b> nucleic acid extraction kit	K1-11-50-CE K1-11-100-CE K1-12-100-CE
11.	<b>DNA-sorb-B</b> nucleic acid extraction kit	K1-2-50-CE K1-2-100-CE
12.	<b>DNA-sorb-C</b> nucleic acid extraction kit	K1-6-50-CE
13.	<b>EDEM</b> reagents kit for extraction of DNA by express method	K2-17-100-CE
14.	<b>MAGNO-sorb</b> Nucleic Acid Extraction Kit	K2-16-200-CE K2-16-1000-CE
15.	<b>RIBO-prep</b> nucleic acid extraction kit	K2-9-Et-50-CE K2-9-Et-100-CE
16.	<b>RIBO-sorb</b> nucleic acid extraction kit	K2-1-Et-50-CE K2-1-Et-100-CE
17.	<b>RIBO-zol-B</b> nucleic acid extraction kit	K2-3-50-CE K2-3-100-CE
18.	<b>RIBO-zol-C</b> nucleic acid extraction kit	K2-13-50-CE K2-13-100-CE
19.	<b>REVERTA-L</b> RT reagents kit	K3-4-50-CE K3-4-100-CE
20.	<b>EPh</b> detection agarose kit	K5-200-CE K5-300-CE K6-200-CE K6-300-CE
21.	<b>AmpliSens® Adenovirus-EPh</b> PCR kit	V23-50-R0,2-CE V23-50-R0,5-CE
22.	<b>AmpliSens® All-screen-FEP</b> PCR kit	B45-FEP-CE
23.	<b>AmpliSens® All-screen-FRT</b> PCR kit	R-B45(RG,iQ)-CE
24.	<b>AmpliSens® ARVI-screen-FRT</b> PCR kit	R-V57-100-F(RG,iQ,Dt)-CE
25.	<b>AmpliSens® Astrovirus-EPh</b> PCR kit	V19-50-R0,2-CE V19-50-R0,5-CE
26.	<b>AmpliSens® Bacillus anthracis-FRT</b> PCR kit	R-B41(RG)-CE
27.	<b>AmpliSens® Bordetella multi-FRT</b> PCR kit	R-B84-100-F(RG,iQ,Dt)-CE
28.	<b>AmpliSens® Borrelia burgdorferi sensu lato-EPh</b> PCR kit	B37-50-R0,2-CE B37-50-R0,5-CE

29.	<b>AmpliSens® <i>Borrelia burgdorferi sensu lato</i>-FEP PCR kit</b>	B37-FEP-CE
30.	<b>AmpliSens® <i>Borrelia burgdorferi sensu lato</i>-FRT PCR kit</b>	R-B37(RG)-CE
31.	<b>AmpliSens® <i>Borrelia miyamotoi</i>-FRT PCR kit</b>	H-2791-1-CE H-2792-1-4-CE
32.	<b>AmpliSens® <i>Brucella</i> spp.-FEP PCR kit</b>	B10-50-R0,2-FEP-CE B10-50-R0,5-FEP-CE
33.	<b>AmpliSens® <i>Brucella</i> spp.-FRT PCR kit</b>	R-B10(RG)-CE
34.	<b>AmpliSens® <i>C.albicans</i> / <i>C.glabrata</i> / <i>C.krusei</i>-MULTIPRIME-FRT PCR kit</b>	R-F3-F(RG,iQ)-CE R-F3(RG)-CE
35.	<b>AmpliSens® <i>Campylobacter</i> spp.-EPh PCR kit</b>	B35-50-R0,2-CE B35-50-R0,5-CE
36.	<b>AmpliSens® <i>Candida albicans</i>-EPh PCR kit</b>	F1-100-R0,2-CE F1-100-R0,5-CE
37.	<b>AmpliSens® <i>Candida albicans</i>-FEP PCR kit</b>	F1-100-R0,2-FEP-CE F1-100-R0,5-FEP-CE
38.	<b>AmpliSens® <i>Candida albicans</i>-FRT PCR kit</b>	R-F1-F(RG,iQ)-CE R-F1(RG)-CE
39.	<b>AmpliSens® CCHFV-FRT PCR Kit</b>	R-V22-50-F(RG,iQ,Mx,Dt)-CE
40.	<b>AmpliSens® <i>Corynebacterium diphtheriae</i>-EPh PCR kit</b>	B8-50-R0,2-CE B8-50-R0,5-CE
41.	<b>AmpliSens® <i>Corynebacterium diphtheriae</i> / <i>tox-genes</i>-FRT PCR kit</b>	H-2842-1-CE H-2843-1-4-CE
42.	<b>AmpliSens® Cov-Bat-FRT PCR kit</b>	R-V65-F-CE
43.	<b>AmpliSens® <i>Coxiella burnetii</i>-FRT PCR kit</b>	R-B85-50-F(RG,iQ,Mx,Dt)-CE
44.	<b>AmpliSens® <i>Cronobacter sakazakii</i>-FEP PCR kit</b>	B58-FEP-CE
45.	<b>AmpliSens® <i>Cronobacter sakazakii</i>-FRT PCR kit</b>	R-B58(RG,iQ)-CE
46.	<b>AmpliSens® <i>Cryptococcus neoformans</i>-FRT PCR kit</b>	R-F4-F(RG,iQ)-CE
47.	<b>AmpliSens® <i>Dengue virus</i> type-FRT PCR kit</b>	R-V63(RG,CFX)-CE
48.	<b>AmpliSens® <i>Dengue virus</i>-FRT PCR kit</b>	H-2391-1-CE H-2392-1-4-CE
49.	<b>AmpliSens® EBOV Zaire-FRT PCR kit</b>	R-V69-50-F-CE
50.	<b>AmpliSens® EBV-EPh PCR kit</b>	V9-100-R0,2-CE V9-100-R0,5-CE
51.	<b>AmpliSens® EBV-screen/monitor-FRT PCR kit</b>	R-V9-100-S(RG,iQ,Mx)-CE
52.	<b>AmpliSens® EHEC-FEP PCR kit</b>	B59-FEP-CE
53.	<b>AmpliSens® EHEC-FRT PCR kit</b>	R-B59(RG,iQ)-CE
54.	<b>AmpliSens® Enterovirus 71-FRT PCR kit</b>	R-V64-F-CE
55.	<b>AmpliSens® Enterovirus-EPh PCR kit</b>	V16-50-CE V16-50-R0,2-CE V16-50-R0,5-CE
56.	<b>AmpliSens® Enterovirus-FEP PCR kit</b>	V16-50-R0,2-FEP-CE V16-50-R0,5-FEP-CE
57.	<b>AmpliSens® Enterovirus-FRT PCR kit</b>	R-V16(RG)-CE R-V16-F(RG)-CE R-V16-F-CE
58.	<b>AmpliSens® Escherichioses-FEP PCR kit</b>	B62-FEP-CE
59.	<b>AmpliSens® Escherichioses-FRT PCR kit</b>	R-B62(RG,iQ)-CE

60.	<b>AmpliSens® FiloA-screen-FRT PCR kit</b>	H-2781-1-4-CE
61.	<b>AmpliSens® Florocenosis / Aerobes-FRT PCR kit</b>	R-B88-100-FT-CE
62.	<b>AmpliSens® Florocenosis / Bacterial vaginosis-FRT PCR kit</b>	R-B74-100-FT(RG)-CE
63.	<b>AmpliSens® Florocenosis / Candida-FRT PCR kit</b>	R-F5-100-FT(RG,CFX)-CE
64.	<b>AmpliSens® Florocenosis / Mycoplasma-FRT PCR kit</b>	R-B75-100-FT(RG,iQ,Mx)-CE
65.	<b>AmpliSens® Gardnerella vaginalis-EPh PCR kit</b>	B7-100-R0,2-CE B7-100-R0,5-CE
66.	<b>AmpliSens® Gardnerella vaginalis-FEP PCR kit</b>	B7-100-R0,2-FEP-CE B7-100-R0,5-FEP-CE
67.	<b>AmpliSens® Gardnerella vaginalis-FRT PCR kit</b>	R-B7-F(RG,iQ)-CE R-B7(RG)-CE
68.	<b>AmpliSens® Genoscreen-IL28B-FRT PCR kit</b>	R-O5-100-F(RG,iQ,Dt,CFX)-CE
69.	<b>AmpliSens® Hantavirus-EPh PCR kit</b>	V39-50-R0,2-CE V39-50-R0,5-CE
70.	<b>AmpliSens® HAV-EPh PCR kit</b>	V4-50-R0,2-CE V4-50-R0,5-CE
71.	<b>AmpliSens® HAV-FEP PCR kit</b>	V4-FEP-CE
72.	<b>AmpliSens® HAV-FRT PCR kit</b>	R-V4(RG,iQ)-CE
73.	<b>AmpliSens® Helicobacter pylori-EPh PCR kit</b>	B9-50-R0,2-CE B9-50-R0,5-CE
74.	<b>AmpliSens® Helicobacter pylori-FEP PCR kit</b>	B9-FEP-CE
75.	<b>AmpliSens® Helicobacter pylori-FRT PCR kit</b>	R-B9(RG,iQ)-CE
76.	<b>AmpliSens® HGV-EPh PCR kit</b>	V2-100-R0,2-CE V2-100-R0,5-CE
77.	<b>AmpliSens® HGV-FRT PCR kit</b>	R-V2-50-F(RG,iQ,Mx,Dt)-CE
78.	<b>AmpliSens® HHV6-EPh PCR kit</b>	V10-100-R0,2-CE V10-100-R0,5-CE
79.	<b>AmpliSens® HHV6-screen-titre-FRT PCR kit</b>	R-V10-T(RG,iQ,Mx)-CE
80.	<b>AmpliSens® HHV7-screen/monitor-FRT PCR kit</b>	H-2431-1-1-CE
81.	<b>AmpliSens® HPV 16/18-EPh PCR kit</b>	V12-100-R0,2-CE V12-100-R0,5-CE
82.	<b>AmpliSens® HPV 16/18-FEP PCR kit</b>	V12-FEP-CE V12-100-R0,5-FEP-CE V12-100-R0,2-FEP-CE V12-Mod-FEP-CE
83.	<b>AmpliSens® HPV 16/18-FRT PCR kit</b>	R-V12(RG,iQ,Mx)-CE R-V12-100(RG)-CE R-V12-100(iQ,Mx,Dt)-CE R-V12-Mod(RG,iQ,Mx)-CE
84.	<b>AmpliSens® HPV 31/33-EPh PCR kit</b>	V13-100-R0,2-CE V13-100-R0,5-CE
85.	<b>AmpliSens® HPV 35/45-EPh PCR kit</b>	V14-100-R0,2-CE V14-100-R0,5-CE
86.	<b>AmpliSens® HPV 6/11-EPh PCR kit</b>	V11-50F-CE
87.	<b>AmpliSens® HPV 6/11-FEP PCR kit</b>	V11-FEP-CE V11-100-R0,5-FEP-CE V11-100-R0,2-FEP-CE V11-Mod-FEP-CE

88.	<b>AmpliSens<sup>®</sup> HPV 6/11-FRT PCR kit</b>	R-V11(RG,iQ,Mx)-CE R-V11-100(RG)-CE R-V11-100(iQ,Mx,Dt)-CE R-V11-Mod(RG,iQ,Mx)-CE
89.	<b>AmpliSens<sup>®</sup> HPV HCR genotype-EPh PCR kit</b>	V25-50F-CE
90.	<b>AmpliSens<sup>®</sup> HPV HCR genotype-FRT PCR kit</b>	R-V25(RG,iQ,Mx)-CE
91.	<b>AmpliSens<sup>®</sup> HPV HCR genotype-titre-FRT PCR kit</b>	R-V67-F-CE R-V71-F-CE
92.	<b>AmpliSens<sup>®</sup> HPV HCR screen-EPh PCR kit</b>	V31-100F-CE
93.	<b>AmpliSens<sup>®</sup> HPV HCR screen-FEP PCR kit</b>	V31-3x-FEP-CE V31-FEP-CE
94.	<b>AmpliSens<sup>®</sup> HPV HCR screen-titre-14-FRT PCR kit</b>	H-2311-1-13-CE
95.	<b>AmpliSens<sup>®</sup> HPV HCR screen-titre-FRT PCR kit</b>	R-V31-T-2x(RG,iQ,SC)-CE R-V31-T-4x(RG,iQ,Mx)-CE R-V31-F-CE
96.	<b>AmpliSens<sup>®</sup> HSV I, II-EPh PCR kit</b>	V8-100-R0,2-CE V8-100-R0,5-CE
97.	<b>AmpliSens<sup>®</sup> HSV I, II-EPh-M PCR kit</b>	V8-100-R0,2-Mod-CE V8-100-R0,5-Mod-CE
98.	<b>AmpliSens<sup>®</sup> HSV I, II-FRT PCR kit</b>	R-V8-F(RG,iQ)-CE R-V8(RG)-CE
99.	<b>AmpliSens<sup>®</sup> HSV-typing-FEP PCR kit</b>	V38-100-R0,2-FEP-CE V38-100-R0,5-FEP-CE
100.	<b>AmpliSens<sup>®</sup> HSV-typing-FRT PCR kit</b>	R-V38-F(RG,iQ)-CE R-V38(RG)-CE
101.	<b>AmpliSens<sup>®</sup> Human enterovirus-FEP PCR kit</b>	H-2771-2-2-CE H-2772-2-5-CE H-2773-2-CE
102.	<b>AmpliSens<sup>®</sup> Human enterovirus-FRT PCR kit</b>	H-2771-1-2-CE H-2773-1-CE
103.	<b>AmpliSens<sup>®</sup> Influenza virus A/B-EPh PCR kit</b>	V36-50-R0,2-CE V36-50-R0,5-CE
104.	<b>AmpliSens<sup>®</sup> Influenza virus A/H1-swine-FEP PCR kit</b>	V55-50-R0,5-FEP-CE V55-50-R0,2-FEP-CE
105.	<b>AmpliSens<sup>®</sup> Influenza virus A/H1-swine-FRT PCR kit</b>	R-V55(RG)-CE R-V55-F(RG,iQ,Dt,SC)-CE R-V55-F(SC)-CE
106.	<b>AmpliSens<sup>®</sup> Influenza virus A-type-FEP PCR kit</b>	V54-50-R0,5-FEP-CE V54-50-R0,2-FEP-CE
107.	<b>AmpliSens<sup>®</sup> Influenza virus A-type-FRT PCR kit</b>	R-V54(RG)-CE R-V54(iQ,Dt)-CE R-V54-100-F(RG,iQ,Dt,SC)-CE
108.	<b>AmpliSens<sup>®</sup> Influenza virus A-type-H5, H7, H9-FRT PCR kit</b>	R-V66-F-CE
109.	<b>AmpliSens<sup>®</sup> Influenza virus A H5N1-FEP PCR kit</b>	V33-50-R0,2-FEP-CE V33-50-R0,5-FEP-CE
110.	<b>AmpliSens<sup>®</sup> Influenza virus A H5N1-FRT PCR kit</b>	R-V33(RG)-CE R-V33(SC)-CE
111.	<b>AmpliSens<sup>®</sup> Influenza virus A/B-FEP PCR kit</b>	V36-Mod-50-R0,2-FEP-CE V36-Mod-50-R0,5-FEP-CE
112.	<b>AmpliSens<sup>®</sup> Influenza virus A/B-FRT PCR kit</b>	R-V36-50-Mod-CE R-V36-100-F-Mod(RG,iQ,Dt,CFX,SC)-CE

113.	<b>AmpliSens® JCV-BKV screen-monitor-FRT</b> PCR kit	H-2441-1-1-CE
114.	<b>AmpliSens® Legionella pneumophila-FEP</b> PCR kit	B50-R0,2-FEP-CE B50-50-R0,5-FEP-CE
115.	<b>AmpliSens® Legionella pneumophila-FRT</b> PCR kit	R-B50(RG)-CE
116.	<b>AmpliSens® Leptospira-FRT</b> PCR kit	R-B49(RG,iQ)-CE
117.	<b>AmpliSens® Leucosis Quantum M-bcr-FRT</b> PCR kit	TR-O1(RG,iQ,Mx,A)-CE
118.	<b>AmpliSens® Leucosis Quantum M-bcr- FRT-M</b> PCR kit	SK9-3651-1-1-CE
119.	<b>AmpliSens® Listeria monocytogenes-EPh</b> PCR kit	B14-50-R0,2-CE B14-50-R0,5-CE
120.	<b>AmpliSens® Listeria monocytogenes- screen-titre-FRT</b> PCR kit	R-B14-100-FT(RG,iQ)-CE
121.	<b>AmpliSens® Listeria monocytogenes- screen/monitor-FRT</b> PCR kit	H-2161-1-1-CE
122.	<b>AmpliSens® M.genitalium-screen-titre-FRT</b> PCR kit	R-B4-100-FT(RG,iQ,Mx)-CE
123.	<b>AmpliSens® M.hominis / G.vaginalis- MULTIPRIME-FEP</b> PCR kit	B48-100-R0,2-FEP-CE B48-100-R0,5-FEP-CE
124.	<b>AmpliSens® M.hominis-screen-titre-FRT</b> PCR kit	R-B3-100-FT(RG,iQ,Mx)-CE
125.	<b>AmpliSens® MBT-EPh</b> PCR kit	B15-100-R0,2-CE B15-100-R0,5-CE
126.	<b>AmpliSens® MDR KPC/OXA-48-FRT</b> PCR kit	R-C2(RG,CFX)-CE
127.	<b>AmpliSens® MDR MBL-FRT</b> PCR kit	R-C1(RG,CFX)-CE
128.	<b>AmpliSens® MRSA-screen-titre-FRT</b> PCR kit	R-B78-100-FT(RG,iQ)-CE
129.	<b>AmpliSens® MTC-diff-FRT</b> PCR kit	R-B80(RG,iQ,Dt,SC)-CE
130.	<b>AmpliSens® MTC-FEP</b> PCR kit	B57-FEP-CE
131.	<b>AmpliSens® MTC-FRT</b> PCR kit	R-B57(RG,iQ,SC,Dt)-CE
132.	<b>AmpliSens® Mycoplasma genitalium-EPh</b> PCR kit	B4-100-R0,2-CE B4-100-R0,5-CE
133.	<b>AmpliSens® Mycoplasma genitalium-FEP</b> PCR kit	B4-100-R0,2-FEP-CE B4-100-R0,5-FEP-CE
134.	<b>AmpliSens® Mycoplasma genitalium-FRT</b> PCR kit	R-B4-F(RG,iQ)-CE R-B4(RG)-CE
135.	<b>AmpliSens® Mycoplasma hominis-EPh</b> PCR kit	B3-100-R0,2-CE B3-100-R0,5-CE
136.	<b>AmpliSens® Mycoplasma hominis-FEP</b> PCR kit	B3-100-R0,2-FEP-CE B3-100-R0,5-FEP-CE
137.	<b>AmpliSens® Mycoplasma hominis-FRT</b> PCR kit	R-B3-F(RG,iQ)-CE R-B3(RG)-CE
138.	<b>AmpliSens® Mycoplasma pneumoniae-EPh</b> PCR kit	B40-50-R0,2-CE B40-50-R0,5-CE
139.	<b>AmpliSens® N.meningitidis / H.influenzae / S.pneumonia-FEP</b> PCR kit	B25-FEP-CE
140.	<b>AmpliSens® N.meningitidis / H.influenzae / S.pneumonia-FRT</b> PCR kit	R-B25(RG,iQ)-CE
141.	<b>AmpliSens® Neisseria gonorrhoeae-EPh</b> PCR kit	B5-100-R0,2-CE B5-100-R0,5-CE B24-100-R0,2-CE B24-100-R0,5-CE
142.	<b>AmpliSens® Neisseria gonorrhoeae- screen-FEP</b> PCR kit	B51-100-R0,2-FEP-CE B51-100-R0,5-FEP-CE

143.	<b>AmpliSens® <i>Neisseria gonorrhoeae</i>-screen-FRT PCR kit</b>	R-B51-F(RG,iQ)-CE R-B51(RG)-CE
144.	<b>AmpliSens® <i>Neisseria gonorrhoeae</i>-test-FEP PCR kit</b>	B56-100-R0,2-FEP-CE B56-100-R0,5-FEP-CE
145.	<b>AmpliSens® <i>Neisseria gonorrhoeae</i>-test-FRT PCR kit</b>	R-B56(RG)-CE R-B56-F(RG,iQ)-CE
146.	<b>AmpliSens® <i>Neisseria meningitidis</i> A, B, C-EPh PCR kit</b>	B26-50-R0,2-CE B26-50-R0,5-CE
147.	<b>AmpliSens® <i>Neisseria</i> spp., <i>Haemophilus</i> spp., <i>Streptococcus</i> spp.-EPh PCR kit</b>	B25-50-R0,2-CE B25-50-R0,5-CE
148.	<b>AmpliSens® <i>Norovirus</i> GI / GII-FRT PCR kit</b>	H-2751-1-3-CE
149.	<b>AmpliSens® <i>Norovirus</i> genotypes 1, 2-EPh PCR kit</b>	V27-50-R0,2-CE V27-50-R0,5-CE
150.	<b>AmpliSens® <i>Parvovirus</i> B19-FRT PCR kit</b>	R-V49(RG,iQ,Mx)-CE
151.	<b>AmpliSens® <i>Pneumocystis jirovecii</i> (<i>carinii</i>)-FRT PCR kit</b>	R-F2(RG,iQ,Mx)-CE R-F2-Mod(RG,iQ,Mx)-CE
152.	<b>AmpliSens® <i>Poliovirus</i>-FEP PCR kit</b>	V58-FEP-CE
153.	<b>AmpliSens® <i>Poliovirus</i>-FRT PCR kit</b>	R-V58(RG,iQ)-CE
154.	<b>AmpliSens® <i>Pseudomonas aeruginosa</i>-screen-titre-FRT PCR kit</b>	R-B76-50-FT(RG,iQ)-CE
155.	<b>AmpliSens® <i>Rickettsia conorii</i>-FRT PCR kit</b>	H-2741-1-CE H-2742-1-4-CE
156.	<b>AmpliSens® <i>Rotavirus</i> / <i>Norovirus</i> / <i>Astrovirus</i>-FEP PCR kit</b>	V40-FEP-CE
157.	<b>AmpliSens® <i>Rotavirus</i> / <i>Norovirus</i> / <i>Astrovirus</i>-FRT PCR kit</b>	R-V40(RG,iQ)-CE
158.	<b>AmpliSens® <i>Rotavirus</i>-EPh PCR kit</b>	V15-50-R0,2-CE V15-50-R0,5-CE
159.	<b>AmpliSens® <i>Salmonella</i> spp.-EPh PCR kit</b>	B11-50-R0,2-CE B11-50-R0,5-CE
160.	<b>AmpliSens® <i>Salmonella</i> spp.-FEP PCR kit</b>	B11-FEP-CE
161.	<b>AmpliSens® <i>Salmonella</i> spp.-FRT PCR kit</b>	R-B11(RG,iQ)-CE
162.	<b>AmpliSens® <i>Salmonella typhi</i>-FEP PCR kit</b>	B63-FEP-CE
163.	<b>AmpliSens® <i>Salmonella typhi</i>-FRT PCR kit</b>	R-B63(RG,iQ)-CE
164.	<b>AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i> / <i>Salmonella</i> spp. / <i>Campylobacter</i> spp.-FRT PCR kit</b>	R-B44(RG,iQ)-CE
165.	<b>AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i> / <i>Salmonella</i> spp. / <i>Campylobacter</i> spp.-FEP PCR kit</b>	B44-FEP-CE
166.	<b>AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i>-EPh PCR kit</b>	B12-50-R0,2-CE B12-50-R0,5-CE
167.	<b>AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i>-FEP PCR kit</b>	B12-FEP-CE
168.	<b>AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i>-FRT PCR kit</b>	R-B12(RG,iQ)-CE
169.	<b>AmpliSens® <i>Streptococcus agalactiae</i>-screen-titre-FRT PCR kit</b>	R-B77-100-FT(RG,iQ)-CE
170.	<b>AmpliSens® <i>Streptococcus pyogenes</i>-screen-titre-FRT PCR kit</b>	R-B82-100-FT(RG,iQ)-CE
171.	<b>AmpliSens® <i>Streptococcus pyogenes</i>-screen/monitor-FRT PCR kit</b>	H-2171-1-1-CE H-2172-1-14-CE
172.	<b>AmpliSens® <i>Streptococcus</i> spp.-EPh PCR kit</b>	B18-50-R0,2-CE B18-50-R0,5-CE

173.	<b>AmpliSens® <i>T.vaginalis</i> / <i>N.gonorrhoeae</i>-MULTIPRIME-FEP</b> PCR kit	B65-100-R0,5-FEP-CE B65-100-R0,2-FEP-CE
174.	<b>AmpliSens® <i>T.vaginalis</i> / <i>N.gonorrhoeae</i>-MULTIPRIME-FRT</b> PCR kit	R-B65(RG)-CE R-B65-F(RG,iQ)-CE
175.	<b>AmpliSens® TBE-FRT</b> PCR kit	R-V52(RG)-CE
176.	<b>AmpliSens® TBEV, <i>B.burgdorferi</i> sl, <i>A.phagocytophilum</i>, <i>E.chaffeensis</i> / <i>E.muris</i>-FRT</b> PCR kit	R-V59(RG,iQ,Mx,Dt)-CE R-V59-50-F(RG,iQ,Mx,Dt)-CE
177.	<b>AmpliSens® <i>Treponema pallidum</i>-FEP</b> PCR kit	B20-100-R0,2-FEP-CE B20-100-R0,5-FEP-CE
178.	<b>AmpliSens® <i>Treponema pallidum</i>-FRT</b> PCR kit	R-B20-F(RG,iQ)-CE R-B20(RG)-CE
179.	<b>AmpliSens® <i>Trichomonas vaginalis</i>-EPh</b> PCR kit	B6-100-R0,2-CE B6-100-R0,5-CE
180.	<b>AmpliSens® <i>Trichomonas vaginalis</i>-FEP</b> PCR kit	B6-100-R0,2-FEP-CE B6-100-R0,5-FEP-CE
181.	<b>AmpliSens® <i>Trichomonas vaginalis</i>-FRT</b> PCR kit	R-B6-F(RG,iQ)-CE R-B6(RG)-CE
182.	<b>AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i>-EPh</b> PCR kit	B19-100-R0,2-CE B19-100-R0,5-CE
183.	<b>AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i>-FEP</b> PCR kit	B19-100-R0,2-FEP-CE B19-100-R0,5-FEP-CE
184.	<b>AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i>-FRT</b> PCR kit	R-B19-F(RG,iQ)-CE R-B19(RG)-CE
185.	<b>AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i>-screen-titre-FRT</b> PCR kit	R-B19-100-FT(RG,iQ,Mx)-CE
186.	<b>AmpliSens® <i>Ureaplasma</i> spp.-EPh</b> PCR kit	B2-100-R0,2-CE B2-100-R0,5-CE
187.	<b>AmpliSens® <i>Ureaplasma</i> spp.-FEP</b> PCR kit	B2-100-R0,2-FEP-CE B2-100-R0,5-FEP-CE
188.	<b>AmpliSens® <i>Ureaplasma</i> spp.-FRT</b> PCR kit	R-B2-F(RG,iQ)-CE R-B2(RG)-CE
189.	<b>AmpliSens® <i>Ureaplasma</i> spp.-screen-titre-FRT</b> PCR kit	R-B2-100-FT(RG,iQ,Mx)-CE
190.	<b>AmpliSens® <i>Vibrio cholerae</i>-FRT</b> PCR kit	R-B53(RG)-CE
191.	<b>AmpliSens® VZV-FRT</b> PCR kit	R-V61-50-F(RG)-CE
192.	<b>AmpliSens® WNV-FRT</b> PCR kit	R-V53(RG,iQ,Mx)-CE
193.	<b>AmpliSens® Yellow fever virus-FRT</b> PCR kit	H-2461-1-CE H-2462-1-4-CE
194.	<b>AmpliSens® <i>Yersinia enterocolitica</i> / <i>Y.pseudotuberculosis</i>-FEP</b> PCR kit	B64-FEP-CE
195.	<b>AmpliSens® <i>Yersinia enterocolitica</i> / <i>Y.pseudotuberculosis</i>-FRT</b> PCR kit	R-B64(RG,iQ)-CE
196.	<b>AmpliSens® <i>Yersinia pestis</i>-FRT</b> PCR kit	R-B79(RG,iQ,Dt)-CE
197.	<b>AmpliSens® Zika virus-FRT</b> PCR kit	H-2411-1-CE
198.	<b>BRCA-screen</b> kit	S-1619-6-CE
199.	<b>PEERO-prep</b> reagent kit for sample preparation	K15-1611-40-CE