www.mobiakcare.com

Our Company

First Aid Products - Stretchers


Folding Aluminium Stretcher

Ref: **0804591**

Technical Specifications

- Frame: Aluminium
- Material: Canvas
- Includes: 2 Safety Belts
- Packing Dimensions (LxWxH) : 104 x 18 x 11cm
- Dimensions Open (LxWxH) : 208 x 55 x 14cm
- Product Weight: 4.78Kg
- Weight Tolerance: 159Kg

Reinforced Safety Belt

Ref: **0806482** *One Piece*

- With Car-type Safety Lock
- Suitable for All Wheelchair Types & Stretcher
- Secured around the Seat or the Backrest

Length **(A)**: 162cm
Width **(B)**: 5cm

Soft Stretcher

Ref: **0804592**

Technical Specifications

- Frame: Canvas, Waterproof, Washable
- Includes: 2 Safety Belts
- Dimensions when deployed (LxW) : 205 x 70 cm (measurement without the handles)
- Dimensions when closed (LxW) : 74 x 35 cm
- Product Weight: 1.72Kg
- Weight Tolerance: 159Kg

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INTENDED USE:

The product is used for the transportation of patients in emergency or in situations demanding a quick and safe evacuation of the patient. It has to be used exclusively by professional and qualified personnel, duly trained for emergency care. Can be used in ambulances and in hospital and sports environments.

CAUTION (!):

- In case of damage or malfunction of your product DOT NOT use the product and please contact the authorized and trained representative.
- Do not use the product for a purpose not indicated in this manual.
- MOBIAK S.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration on the frame of the product.
- MOBIAK S.A reserves the right to change the information contained in this document without previous notice
- There is a ± 3%. Deviation for the Dimensions of the Products.
- It is always recommended to follow the procedures designed for emergency response.

GENERAL WARNINGS:

- If you don't read this user manual, it is preferable not to use this device or another available part.
- If you don't understand the usage or the precautions, please contact with the dealer or the proper technical person before you use the device because may be caused damage.
- READ carefully the entire manual before using the product.
- Save this document for future reference.
- The product should be used only by qualified personnel.
- Keep away the product from heat sources.
- Be careful when children are nearby and do not allow the children to play or operate with the device.
- Special care is always required where there are moving parts that can cause limb entrapment or injury.

BEFORE EVERY USE:

- Check if the frame is damaged in order to guarantee a safe use of the product (there are no cracks or fractures in the frame).
- In case of damage do not use the product and contact your dealer for further instructions.
- Check, if the product is properly assembled and all screws are secure and well screwed on.
- Always check the condition of wear for the mechanical parts to ensure the absolute safe use of the product for people and objects.
- Periodically check its sturdiness and the stability.

USE OF PRODUCT:

Open the stretcher.

Successfully move the patient on the stretcher and secure with the safety belts.

The presence of handles allows you to lift the patient by grabbing the stretcher from several sides.

CLEANING AND DISINFECTION:

- Use only a wet cloth and mild soap on surfaces where they need cleaning. Then dry thoroughly before you use the product.
- Dry away from direct sunlight.
- DO NOT leave the product wet.
- Do not use chemical cleaners to clean the frame. This could cause damage to the surface of the product.
- Never use acids, alkalis or solvents such as acetone
- If you need to disinfect the product, use a common, mild disinfectant.
- Protect the item from scratches, cuts and punctures.

DECLARATION OF CONFORMITY:

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 745/2017 and where appropriate, the standards and legislation referred to.

ENVIRONMENTAL PROTECTION:

If one day you find that your product needs to be replaced or it is no longer working for you, consider protecting the environment:

- Do not dispose your product along with the rest of the public waste (this is also the meaning of the shown recycling sign).
- Contact your Public Authorities and they will instruct you of the Recycling centers to which your product must be disposed.
- Correct disposal of your product helps the protection of the environment as well as the recycling of the product's components.

INCIDENT REPORTING:

Any serious incident that occurs in relation to use of the device must be reported to the manufacturer and a competent authority of the Member State in which the user and / or patient resides.

WARRANTY:

The product is guaranteed for a period of two (2) years from the purchasing date. The warranty apply on factory defects & DOES NOT apply on damages resulting from misuse, poor maintenance, modification, overuse, or non-compliance with the product's operating instructions. Also DOES NOT apply on the parts, that can wear out during use or over time. Damages or defects caused by natural disasters, unauthorized maintenance or repair, power supply problems (where applicable) & carriers are not covered by the warranty. No service or spare part is applied unless it has been inspected by an authorized service department or the manufacturer of the product. Shipping costs for goods of warranty are to be paid by the end user- customer or by dealer.

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We hereby certify that the under mentioned organization has established and maintains for the listed activities a quality management system that complies with the requirements of the standard.

Validity of the certificate is based upon the successful completion of surveillance audits.

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Address: **KATHIANA AKROTIRIOU, CHANIA CRETE GREECE.**

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Activities: **• MANUFACTURE AND TECHNICAL SUPPORT OF MEDICAL DEVICES
(as listed in annex).**
**• DISTRIBUTION OF MEDICAL DEVICES AND THEIR CONSUMABLES.
(as listed in annex).**

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ACTIVITIES:

• MANUFACTURE AND TECHNICAL SUPPORT OF MEDICAL DEVICES:

- RESPIRATORY EQUIPMENT FOR DIAGNOSIS, TREATMENT AND HOME CARE.
- ORTHOPEDICS FOR EXTERNAL USE AND PHYSIOTHERAPY.
- BEDS FOR HOME CARE.
- ANTI DECUBITUS AIDS.
- MOBILITY AIDS.
- WALKING AIDS.
- LIFTING AND BATHROOM AIDS.
- WHEELCHAIRS.
- MEDICAL EQUIPMENT.
- HOSPITAL HOTEL MEDICAL EQUIPMENT.

• DISTRIBUTION OF MEDICAL DEVICES:

- RESPIRATORY EQUIPMENT FOR DIAGNOSIS, TREATMENT AND HOME CARE.
- ORTHOPEDICS FOR EXTERNAL USE AND PHYSIOTHERAPY.
- BEDS FOR HOME CARE.
- ANTI DECUBITUS AIDS.
- MOBILITY AIDS.
- WALKING AIDS.
- LIFTING AND BATHROOM AIDS.
- WHEELCHAIRS.
- MEDICAL EQUIPMENT.
- HOSPITAL HOTEL MEDICAL EQUIPMENT.
- CONSUMABLES FOR GENERAL MEDICAL USE.
- PORTABLE DEFIBRILLATORS.

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DECLARATION OF CONFORMITY

We are solely responsible for declaring that the Medical Devices mentioned in this statement are of Low-Risk Class (Class I) and comply with the requirements of the European Regulation 2017/745 and where appropriate, the standards and legislation referred to.

MANUFACTURER:	MOBIAK S.A
SRN:	GR-MF-000016256
SEAT ADDRESS:	KATHIANA AKROTIRIOU-CHANIA-CRETE-GREECE
E-MAIL:	mobiakcarequality@mobiak.com
TELEPHONE:	+30 2821063222
WEB SITE:	www.mobiakcare.com

COMPETENT AUTHORITY:	National Organization for Medicines
CONTACT DETAILS:	Address: Mesogeion 284, PS 15562, Cholargos Phone : 2132040000 Website: https://www.eof.gr/

LIST OF PRODUCTS COVERED BY THIS DECLARATION

PRODUCT	CODE	UDI-DI	INTENDED USE	RULE
PATIENT TRANSFER STRETCHER	0806458	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
SOFT STRETCHER	0806459	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
FOLDING STRETCHER, SINGLE FOLD	0806460	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
FOLDING STRETCHER, SINGLE FOLD	0804591	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
STAIR CHAIR STRETCHER	0804593	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
SOFT STRETCHER	0804592	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1

BASKET STRETCHER	0805906	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
SPINE BOARD	0805907	521300690STRETCHERS184Z	IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1
SCOOP STRETCHER	0805910		IT IS SUITABLE FOR TRANSPORTATION OF PATIENTS IN EMERGENCY CONDITIONS OR IN CONDITIONS THAT REQUIRE QUICK AND SAFE REMOVAL OF THE PATIENT.	1

CONFORMITY ASSESSMENT PROCEDURE

According to Annexes II & III of Regulation (EU) 2017/745

APPLIED STANDARDS & LEGAL REQUIREMENTS

ISO 13485:2016, ISO 9001:2015, EN ISO 15223-1: 2016, EN ISO 14971:2019, EN ISO 10993-1:2018, EN ISO 10993-1:2009, ISO 10993-10:2021 , (EU) 2017/745

FOR APPROVAL	
NAME:	SVOURAKI MARIA
POSITION:	CEO
PLACE:	CHANIA
DATE:	20/12/2023
SIGN:	