



(>200 mg/L).

hs-CRP Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1003

Getein1600: Cat.# IF2003

User Manual

INTENDED USE

hs-CRP Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of C-reactive protein (CRP) in serum, plasma whole blood, or fingertip blood. Measurement of CRP is useful for the detection and evaluation of infection, tissue injury and inflammatory disorders. Measurement of high sensitivity CRP (hs-CRP), when used in conjunction with traditional clinical laboratory evaluation of acute coronary syndromes (ACS), may be useful as an independent marker of prognosis for recurrent events in patients with stable coronary disease or ACS.

SUMMARY

C-reactive protein is an acute-phase reactant that precipitated with Pneumococcal C-polysaccharide, and is a non-specific immune response component. CRP has wide distribution in our body, and is an acute-phase protein produced in the liver in response to microbial infection or tissue injury, and the hs-CRP can be used to detect lower concentrations of CRP in serum or plasma. Studies revealed hs-CRP levels seem to be correlated with Atherosclerosis and Acute Myocardial Infarction. And the hs-CRP is an inflammation "marker" for ACS patient and is helpful for primary prevention and risk assessment of cardiovascular disease. Its combination with the ratio of total cholesterol to HDL-C is more accurate than other risk factor in predicting cardiovascular disease.

The American Heart Association and US Centers for Disease Control and Prevention have advocated hs-CRP as a predictor of cardiovascular disease (CVD) to define risk groups: less than 1.0 mg/L indicates low risk, 1.0 to 3.0 mg/L means moderate risk, and the amount above 3.0 mg/L (lower than 10 mg/L) strongly suggests a high risk of CVD. Moreover, higher CRP levels are found in late pregnant women, mild inflammation and viral infections (10-40 mg/L), active inflammation, bacterial infection (40-200 mg/L), severe bacterial infections and burns

with another anti-human hs-CRP monoclonal antibody and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months.

Store the sample diluent/whole blood buffer at 2~8°C for better results.

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch or the cartridge is damaged.
5. Do not open pouches or the cartridge until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for *serum, plasma, whole blood and fingertip blood samples*. *Heparin, sodium citrate and EDTA* can be used as the anticoagulant for plasma, whole blood and fingertip blood. Samples should be free of hemolysis.

PRINCIPLE

The test uses an anti-human hs-CRP monoclonal antibody conjugated with fluorescence latex and another anti-human hs-CRP monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human hs-CRP monoclonal antibody binds with the hs-CRP in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human hs-CRP monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of hs-CRP in sample.

Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of hs-CRP in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100 contains:
 - Getein hs-CRP test card in a sealed pouch with desiccant 25
 - Disposable pipet 25
 - Sample diluent 25
 - SD card 1
 - User manual 1
2. A kit for Getein1600 contains:
 - Sealed cartridge with 24/48 Getein hs-CRP test cards .. 2
 - User manual 1
 - Package specifications:
 - 2×24 tests/kit, 2×48 tests/kit
 - Materials required for Getein1600:
 - Sample diluent 1
 - Box with pipette tips 1
 - Mixing plate 1
3. Sample diluent composition:
 - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
4. A test card consists of:
 - A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human hs-CRP monoclonal antibody, the test line is coated

- Suggest using serum or plasma for better results.
- If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2–8°C or stored at -20°C for 6 months before testing (whole blood sample may be stored up to 3 days at 2–8°C).
- Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples.
- SAMPLE VOLUME (for Getein1100): 10 µl.

TEST PROCEDURE

- Collect specimens according to user manual.
- Test card, sample and reagent should be brought to room temperature before testing.

For Getein1100:

- Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
- On the main interface of Getein1100, press "ENT" button to enter testing interface.
- Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
- Put the test card on a clean table, horizontally placed.
- Using sample transfer pipette, deliver 10 µl of sample into one tube of sample diluent, mix gently and thoroughly. Then drop 100 µl of sample mixture (or 3–4 drops of sample mixture when using disposable pipet) into the sample port on the test card.
- Reaction time: 3 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1600:

- Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
- Place samples in the designated area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

- It is required to perform "SD Card Calib" calibration when using a new batch of kits.
- It is suggested to calibrate once for one batch of kits for Getein1100.
- Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

EXPECTED VALUE

hs-CRP: The expected normal value for hs-CRP was determined by testing samples from 500 apparently healthy individuals. The 95th percentile of the concentration for hs-CRP is 3 mg/L. (The probability that hs-CRP value of a normal person below 3 mg/L is 95%.)

CRP: The expected normal value for CRP was determined by testing samples from 500 apparently healthy individuals. The 95th percentile of the concentration for CRP is 10 mg/L. (The probability that CRP value of a normal person below 10 mg/L is 95%.)

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.5~200 mg/L
Lower Detection Limit	≤0.5 mg/L
Within-Run Precision	≤10%
Between-Run Precision	≤15%

Method Comparison:

The assay was compared with HITACHI 7600/OLYMPUS AU5400 and its matching hs-CRP test kits with 200 serum samples (61 positive samples and 139 negative samples). The correlation coefficient (r) for hs-CRP is 0.941.

LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Samples containing interferences may influence the results. The table below listed the maximum allowance of these potential interferences.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	10 g/L	10 g/L	0.2 g/L

REFERENCES

- Danesh J, Whincup P, Wslker M, et al. Low grade inflammation

and coronary heart disease: prospective study and updated meta-analysis. *BJM* 2000; 321:199–204.

- Rifai N, Ridker PM. Proposed cardiovascular risk assessment algorithm using high-sensitivity C-reactive protein and lipid screening. *Clin Chem* 2001; 47:28–30.
- EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
- EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on hs-CRP Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use	LOT	Batch code
	Temperature limitation	IVD	<i>In vitro</i> diagnostic medical device
	Sufficient for	EC REP	Authorized representative in the European Community
CE	CE mark		Do not use if package is damaged

Thank you for purchasing hs-CRP Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF04-S-02



Getein Biotech, Inc.
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.bio-GP.com.cn



Total IgE Fast Test Kit (Immunofluorescence Assay)

IF1069 for Getein1100
IF5069 for Getein1160
IF3069 for Getein1180
IF4069 for Getein1200
IF2069 for Getein1600

REF

User Manual

INTENDED USE

Total IgE Fast Test Kit (Immunofluorescence Assay) is intended for in vitro quantitative determination of circulating total IgE antibodies in human serum, plasma and whole blood samples. This test can be used as an aid in the clinical diagnosis of IgE mediated allergic disorders in conjunction with other clinical findings.

SUMMARY

IgE is an immunoglobulin with a molecular weight of approximately 190,000 daltons. Produced by plasma cells, IgE has a significant role in atopic diseases such as allergic rhinitis, allergic asthma, and atopic dermatitis. IgE has a high affinity for receptors on mast cells and basophils, mediating the binding of allergens to these cells. The subsequent release of vasoactive amines, such as histamine, produce the clinical manifestations associated with atopic disease. Measurement of IgE serum levels can be important in the diagnosis and treatment of these disorders.

In most nonatopic patients, IgE serum levels are relatively low. However, certain parasitic or helminth infections have been associated with elevated IgE levels due to IgE sensitization of macrophages, eosinophils, and other inflammatory cells. The IgE concentration in a patient is dependent on both the extent of the allergic reaction and the number of different allergens to which the patient is sensitized. Nonallergic normal individuals have IgE concentrations that vary widely and increase steadily during childhood, reaching their highest levels at age 15 to

20, and thereafter remaining constant until about age 60 when they slowly decline.

PRINCIPLE

The test uses an anti-human IgE monoclonal antibody I conjugated with fluorescence latex coated on the sample pad, and another anti-human IgE monoclonal antibody II coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human IgE monoclonal antibody I binds with IgE in sample and forms a marked antigen-antibody complex. This complex moves to the test detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human IgE monoclonal antibody II. The fluorescence intensity of test line increases in proportion to the amount of IgE concentration in sample.

Then insert test card into Getein1100/Getein1160/Getein1180 Immunofluorescence Quantitative Analyzer/ Automatically inserted by Getein1200/Getein1600 Immunofluorescence Quantitative Analyzer (hereafter referred to as Getein1100, Getein1160, Getein1180, Getein1200 and Getein1600), the concentration of IgE in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100/Getein1160/Getein1180 contains:
Package specifications: 25 tests/kit, 10 tests/kit

- 1) Getein Total IgE test card in a sealed pouch with desiccant
- 2) Disposable pipet
- 3) Sample diluent
- 4) User manual: 1 piece/kit
- 5) SD card: 1 piece/kit

2. A kit for Getein1200/Getein1600 contains:
Package specifications: 2x24 tests/kit, 2x48 tests/kit
Sealed cartridge with 24/48 Getein Total IgE test cards

User manual: 1 piece/kit

Materials required for Getein1200/Getein1600:

- 1) Sample diluent: 1 bottle/kit
- 2) Box with pipette tips: 96 tips/kit
- 3) Mixing plate: 1 piece/kit

3. Sample diluent composition :

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

4. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad (coated with fluorescence latex-labeled anti-human IgE monoclonal antibody I), nitrocellulose membrane (the test line is coated with another anti-human IgE monoclonal antibody II, and the control line is coated with goat anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1180 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer
Getein1160 Immunofluorescence Quantitative Analyzer
Getein1200 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test kit at 4~30°C with a valid period of 24 months.

Use the test card for Getein1100/Getein1160/Getein1180 within 1 hour once the foil pouch is opened.

For test card of Getein1200/Getein1600: If the cartridge is opened, it could be stable within 24 hours once exposure to air. Please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

PRECAUTIONS

1. For in vitro diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.

4. Do not use the test card if the foil pouch or the cartridge is damaged.
5. Do not open pouches or the cartridge until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for **serum, plasma and whole blood samples**. Heparin, sodium citrate and EDTA can be used as the anti-coagulant for plasma. Samples should be free of hemolysis.
2. Suggest using serum and plasma for better results. The test should be performed within 4 hours after blood collection.
3. If testing is delayed, serum and plasma samples may be stored up to 5 days at 2~8°C or stored at -20°C for 6 months before testing (whole blood samples may be stored up to 3 days at 2~8°C).
4. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
5. Do not use heat-inactivated or hemolysis samples.
6. **SAMPLE VOLUME (for Getein1100/Getein1160/Getein1180): 100 µL**

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should reach to room temperature before test.

For Getein1100:

1. Confirm SD card lot No. in accordance with test kit lot No. Perform "SD card" calibration when necessary.
2. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
3. Put the test card on a clean table, horizontally

placed.
 4. Using sample transfer pipette, deliver **100 µL** of sample into one tube of sample diluent and mix thoroughly. Then drop **100 µL** of sample mixture into sample well on the test card.

5. **Reaction time: 15 minutes.** Insert the test card into Getein1100 and click on “Start” icon after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1160/Getein1180:

1. Confirm SD card lot No. in accordance with test kit lot No.. Perform “SD card” calibration when necessary.
 2. Enter testing interface of Getein1160/Getein1180.
 3. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.

4. Put the test card on a clean table, horizontally placed.

5. Using sample transfer pipette, deliver **100 µL** of sample into one tube of sample diluent and mix thoroughly. Then drop **100 µL** of sample mixture into sample well on the test card.

6. Insert the test card into Getein 1160/Getein 1180 immediately after sample loading. The analyzer will count down the reaction time (15 minutes) and automatically test the card after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1200/Getein1600:

1. Each cartridge for Getein1200/Getein1600 contains a specific RFID card which can calibrate automatically.

2. Put the sample diluent at the correct position in Getein1200/Getein1600.

3. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1200/Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform “SD card” calibration when using a new batch of kits.
 2. It is suggested to calibrate once for one batch of kits for Getein1100/Getein1160/Getein1180.

3. Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1160/Getein1180/Getein1200/Getein1600.

Others: Measuring range of the Total IgE test kit is 1.00-2000.00 IU/mL. Dilute the sample which concentration is higher than the upper limit with negative samples, and the dilution ratio should be less than 3 times.

EXPECTED VALUE

The expected normal value for total IgE was determined by testing samples from 240 apparently healthy individuals. The reference range of total IgE is 1.00 IU/mL~165.00 IU/mL calculated by using normal distribution methods(95% confidence interval). It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	1.00-2000.00 IU/mL
Lower Detection Limit	≤1.00 IU/mL
Within-Run Precision	≤10%
Between-Run Precision	≤15%

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
 2. Interferents in samples may influence the results. The table below listed the maximum allowance of these potential interferent.

Interferent	Hemoglobin	Triglyceride
Concentration (Max)	50 g/L	0.2 g/L

REFERENCES

- Jones HE, Inouye JC, Mcgerity JL, et al. Atopic disease and serum immunoglobulin-E [J]. Br J Dermatol, 2010, 92 (1): 17-25.
- Stone KD, Prussin C, Metcalfe DD. IgE, Mast Cells, Basophils, and Eosinophils [J]. J Allergy Clin Immunol, 2010, 125 (2): S73-80.
- Johansson SGO, Berglund A, Kjellman NIM. Comparison of IgE values as determined by different solid phase radioimmunoassay methods [J]. Clin Allergy, 2010, 6 (1): 91-98.
- Burney PG, Malmberg E, Chinn S, et al. The distribution of total and specific serum IgE in the European Community Respiratory Health Survey [J]. J Allergy Clin Immunol, 1997, 99 (3):314-322.
- Sampson, Hugh A. Utility of food-specific IgE concentrations in predicting symptomatic food allergy [J]. J Allergy Clin Immunol, 2001, 107 (5): 891-896.
- Ceska M, Lundkvist U. A new and simple radioimmunoassay method for the determination of IgE [J]. Immunochemistry, 1972, 9 (10): 1021-1030.
- Barbee A, Halonen M, Lebowitz M, et al. Distribution of IgE in a community population sample: correlations with age, sex, and allergen skin test reactivity [J]. J Allergy Clin Immunol, 1981, 68 (2): 106-111.
- Wickman M, Lilja G, Soderstrom L, et al. Quantitative analysis of IgE antibodies to food and inhalant allergens in 4-year-old children reflects their likelihood of allergic disease [J]. Allergy, 2015, 60 (5): 650-657.
- Asero R, Ballmer-Weber BK, Beyer K, et al. IgE-mediated food allergy diagnosis: Current status and new perspectives [J]. Mol Nutr Food Res, 2010, 51 (1): 135-147.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on Total

IgE Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more detail in the European Standard EN ISO 15223-1:2021.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		In vitro diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/ European Union
	CE mark		Do not use if package is damaged and consult instructions for use
	Catalogue number		

Thank you for purchasing Total IgE Fast Test Kit (Immunofluorescence Assay).

Please read this user manual carefully before operating to ensure proper use.

Version: WIF94-S-10

Getein Biotech, Inc.
 Add.: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail:tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.getein.com

CMC Medical Devices & Drugs S.L.
 Add.: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain
 Tel: +34951214054



LH

Fast Test Kit

(Immunofluorescence Assay)

IF1055 for Getein 1100
 IF5055 for Getein 1160
 IF3055 for Getein 1180
 IF2055 for Getein 1600
 IF4055 for Getein 1200



Instructions for Use

INTENDED USE

LH Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of LH in human serum and plasma. This test is used to determine menopause, pinpoint ovulation and monitor endocrine therapy. For professional and laboratory use only.

SUMMARY

Luteinizing hormone (LH) is produced in both men and women from the anterior pituitary gland in response to luteinizing hormone-releasing hormone (LH-RH or Gn-RH), which is released by the hypothalamus.

LH is a glycoprotein hormone having two subunits. The alpha subunit is similar to FSH, HCG and TSH. The beta subunit is different from those of the other glycoprotein hormones and confers its biochemical specificity.

In males, LH is also called interstitial cell-stimulating hormone (ICSH). In both males and females, LH secretion is regulated by a balance of positive and negative feedback mechanisms involving the hypothalamic-pituitary axis, the reproductive organs and the pituitary and sex steroid hormones. LH and the other pituitary gonadotropin, FSH, play a critical role in maintaining the normal function of the male and female reproductive systems.

Patients suffering from hypogonadism show increased concentration of serum LH. A decrease in steroid hormone production in females is a result of immature ovaries, primary ovarian failure, polycystic ovary disease or menopause. In these cases, LH secretion is not regulated. A similar loss of regulatory hormones occurs in males when the testes develop abnormally or anorchia exists. Increased concentrations of LH may be found in primary

testicular failure, Klinefelter syndrome, renal failure, cirrhosis, hyperthyroidism and severe starvation. LH is a useful marker in determining the homeostasis fertility regulation via the hypothalamic-pituitary-gonadal axis.

PRINCIPLE

LH Fast Test Kit (Immunofluorescence Assay) is a lateral flow immunoassay designed in a sandwich format. After the sample is applied to the test strip, the fluorescence-labelled LH monoclonal antibody binds with the LH in the sample to form a marked antigen-antibody complex. This complex moves to the detection zone on the test card by capillary action and is captured by another LH monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of LH in the sample. The fluorescent signal intensity can then be analyzed by an appropriate device to quantitatively detect the LH in the sample.

APPLICABLE DEVICE

Getein 1100 Immunofluorescence Quantitative Analyzer
 Getein 1160 Immunofluorescence Quantitative Analyzer
 Getein 1180 Immunofluorescence Quantitative Analyzer
 Getein 1200 Immunofluorescence Quantitative Analyzer
 Getein 1600 Immunofluorescence Quantitative Analyzer

CONTENTS

Materials provided	Getein 1100/Getein 1160/ Getein 1180		Getein 1200/ Getein 1600	
	10 T/kit	25 T/kit	2*24 T/kit	2*48 T/kit
LH test card*	10 pcs	25 pcs	24 test cards in 1 cartridge, and 2 cartridges in 1 box	48 test cards in 1 cartridge, and 2 cartridges in 1 box
Disposable pipet	10 pcs	25 pcs	/	/
Instructions for use	1 pc	1 pc	1 pc	1 pc
SD card	1 pc	1 pc	1 pc in each cartridge	1 pc in each cartridge

* LH test card

A test card mainly consists of: Fluorescence-labelled LH monoclonal antibody, LH monoclonal antibody.

Consumables for Getein 1200/ Getein 1600

- Box with pipette tips (96 tips/box)

- Fixing plate (1 piece/box)

Note:

- The standard curve data can be written to RFID card in the kit. According to the function of RFID card, we define it as "Standard Curve Data Card", short for "SD Card".
- Do not mix or interchange different batches of kits.

STORAGE AND STABILITY

Realtime stability:

Store the kit at 4~30°C with a valid period of 24 months. The test kits are stable until the expiry date printed on the labels.

In-use stability:

-For the test card of Getein 1100/Getein 1160/Getein 1180: Use the test card within 1 hour once the foil pouch is opened.

-For test card of Getein 1200/Getein 1600: If the cartridge is opened, it could be stable within 24 hours once exposure to air. The valid period after opening is 7 days, it is recommended to put the cartridge back to the foil bag and re-seal along the entire edge of zip-seal if not used up.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card or pipet.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- Carefully read and follow instructions for use to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- This test can be used for **serum and plasma samples**.
- Heparin, EDTA and sodium citrate can be used as the anticoagulant for plasma samples.
- It is recommended to test the sample within 4 hours after collection. If testing is delayed, serum and plasma samples are stable for 7 days when stored at 2~8°C and

6 months when stored at -20°C.

- Refrigerated or frozen sample should reach room temperature before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples or hemolysis samples.
- Sample volume (**Getein 1100/Getein 1160/Getein 1180**): 100 µL.

TEST PROCEDURE

- User must carefully read and operate in strict accordance with the instructions for use before testing, otherwise reliable results cannot be guaranteed.
- Test kit and sample should be brought to room temperature before testing.

For Getein 1100:

- Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
- Select the corresponding "Sample" on the analyzer according to the sample type (see the user manual of analyzer for details).
- Remove the test card from the sealed pouch before use and put the test card on a clean table, horizontally placed.
- Use disposable pipet or pipette to drop 100 µL of sample into the sample well on the test card.
- Reaction time: **15 minutes**. Insert the test card into Getein 1100 and click on "Start" icon after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein 1160/Getein 1180:

- Confirm SD card lot No. in accordance with test kit lot No. Perform "SD card" calibration when necessary.
- Select the corresponding "Sample" on the analyzer according to the sample type (see the user manual of analyzer for details).
- Remove the test card from the sealed pouch immediately before use and put the test card on a clean table, horizontally placed.
- Use disposable pipet or pipette to drop 100 µL of sample into the sample well on the test card.
- Insert the test card into Getein 1160/Getein 1180 **immediately** after sample loading. The analyzer will count down the reaction time (15 minutes) and automat-

ically test the card after reaction time is elapsed. The result will be shown on the screen and displayed automatically.

For Getein 1200/Getein 1600:

- 1) Each cartridge for Getein 1200/Getein 1600 contains a specific RFID card which can calibrate automatically.
- 2) Put the consumables at the correct position in Getein 1200/Getein 1600.
- 3) Place samples in the designed area of the sample holder, insert the holder, set parameters (more operational details refer to the user manual of analyzer) and run the instrument, Getein 1200/ Getein 1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform “SD card” calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein 1100/Getein 1160/Getein 1180.
3. Make sure the insertion of test card and the sample are correct and complete.

TEST RESULTS

Getein 1100/Getein 1160/Getein 1180/Getein 1200/Getein 1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein 1100/Getein 1160/Getein 1180/Getein 1200/Getein 1600.

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Some substances in blood as listed below may interfere with the test and cause erroneous results. The maximum allowance concentration of each is as follows:

Interferent	Concentration (Max)
Triglyceride	20.32 mmol/L
Bilirubin	171.0 umol/L
Hemoglobin	300 mg/dL
Human serum albumin	3 g/dL

EXPECTED VALUE

The expected normal value for LH was determined by testing blood samples from apparently healthy individuals. Reference range of LH:

Group	No.	Reference Range (mIU/mL)
Male	221	1.21-8.69
Female	Mid-follicle	210 2.03-11.21
	Mid-cycle peak	188 19.41-103.58
	Mid-luteal phase	211 1.28-12.82
	Postmenopausal	195 11.01-58.99

It is recommended that each laboratory determine the applicability of the reference ranges through experimentation and establish their own laboratory-specific reference ranges if necessary.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.20~150.00 mIU/mL
Detection of Limit	≤0.20 mIU/mL
Within-Run Precision	≤10%
Between-Lot Precision	≤15%

REFERENCES

1. Beastall GH, Ferguson KM, O'Reilly DSJ, et al. Assays for follicle stimulating hormone and luteinizing hormone: Guidelines for the provision of a clinical biochemistry service. Ann Clin Biochem 1987, 24:246-262.
2. Collip JB. William Henry Welch lectures: Some recent advances in physiology of anterior pituitary. J.MA Sinai Hosp 1934; 1:28-71.
3. Mia VG, Julia R, et al. Persistent organic pollutants as predictors of increased FSH: LH ratio in naturally cycling, reproductive age women. Environmental Research 2018, 164.
4. Samipoor Forough, Paksereshd Sedigheh, Rezasoltani Parvaneh, et al. The association between hypogonadism symptoms with serum testosterone, FSH and LH in men. The Aging Male, 2018, 21(1):1-8.
5. Mansour Alizadeh, Ali Nasebakht, Rohollah Valizadeh, et al. A preliminary evaluation of serum level of testosterone, LH, and FSH in patients with varicocele after varicocelectomy as a kidney-related disease.

- Therapeutics and Clinical Risk Management. 2018, 2018(default):1585-1590.
6. Yonggang Huang, Xiaosheng Lu, Zhaoxia Huang, et al. Effects of human chorionic gonadotropin combined with clomiphene on Serum E₂, FSH, LH and PRL levels in patients with polycystic ovarian syndrome. Saudi Journal of Biological Sciences. 2017, 24(2):241-245.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on LH Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2021.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		In vitro diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/ European Union
	CE mark		Do not use if package is damaged and consult instructions for use
	Catalogue number		Caution

Thank you for using LH Fast Test Kit (Immunofluorescence Assay). Please read the instructions for use carefully before operating to ensure proper use. Please report any product problems or adverse events to the below manufacture or authorized representative in the European Community in time.

Version: WIF53-S-09

Getein Biotech, Inc.
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.getein.com

CMC Medical Devices & Drugs S.L.
 Add: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain
 Tel: +34951214054



mAlb Fast Test Kit

(Immunofluorescence Assay)

User Manual

Cat.# IF1009

INTENDED USE

mAlb Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of microalbuminuria (mAlb) in urine. An elevated mAlb concentration below the proteinuric level has long been recognized as a marker of kidney disease and increased cardiovascular risk in diabetic nephropathy.

SUMMARY

Albumin is one of the major plasma proteins. In normal circumstances, albumin molecules are too large to cross the glomerular basement membrane. Therefore, albumin is usually present in very low concentration in urine. Damage to the glomerular basement membrane can alter its permeability. Albumin is then able to enter the urine. Sustained elevation of urinary albumin concentration is called microalbuminuria (mAlb). mAlb arises from increased leakage of glomerular basement membrane. So, mAlb is recognized as a marker of kidney damage. The epidemiology of microalbuminuria reveals a close association between systemic endothelial dysfunction and vascular disease, also implicating glomerular endothelial dysfunction in microalbuminuria.

Recent years, determination of mAlb is linked with increased risk for cardiovascular events rather than progression to end-stage kidney diseases. It is a valuable tool for the detection of cardiovascular risk in diabetic nephropathy. Early detection of microalbuminuria in diabetics is critical because immediate intervention can slow the progression of disease.

PRINCIPLE

The test is based on the competition immune-detection method and uses an anti-human mAlb monoclonal antibody conjugated with fluorescence latex and recombinant mAlb antigen coated on the test line. After the sample has been applied to the test strip, mAlb in the sample will compete with recombinant mAlb antigen on nitrocellulose matrix for fluorescence latex-labelled mAlb monoclonal antibody. As a result, the concentration of mAlb antigen in specimen shows inverse proportionally with the fluorescence intensity of mAlb.

Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100), the concentration of mAlb in sample will be measured and displayed on the screen. The value will be stored in Getein1100 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

A kit contains:

- 1. Getein mAlb test card in a sealed pouch with desiccant 25
- 2. Disposable pipet 25
- 3. User manual 1
- 4. SD card 1

A test card consists of:

A plastic shell and a reagent strip which is composed of sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human mAlb monoclonal antibody, the test line is coated with mAlb recombinant antigen, and the control line is coated with rabbit anti-goat IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card within 1 hour once the foil pouch is opened.

PRECAUTIONS

- 1. For *in vitro* diagnostic use only.
- 2. For professional use only.
- 3. Do not use the kit beyond the expiration date.
- 4. Do not use the test card if the foil pouch is damaged.
- 5. Do not open pouches until ready to perform the test.
- 6. Do not reuse the test card.
- 7. Do not reuse the pipet.
- 8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- 9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- 1. This test can be used for *urine sample*.
- 2. *Urine sample* can be preserved at room temperature for 4 hours, please test it as soon as possible. If testing will be delayed, urine sample may be stored up to 3 days at 2~8°C before testing .
- 3. Do not use frozen urine sample.
- 4. Samples should be brought to room temperature before testing.
- 5. Do not use heat-inactivated samples.
- 6. SAMPLE VOLUME: **100 µl**.

TEST PROCEDURE

- 1. Collect specimens according to user manual.
- 2. Test card, sample should be brought to room temperature before testing.
- 3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
- 4. On the main interface of Getein1100, press "ENT" button to

enter testing interface.

- Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
- Put the test card on a clean table, horizontally placed.
- Using sample transfer pipette, deliver **100 µl** of sample (or 3~4 drops of sample when using disposable pipet) into the sample port on the test card.
- Reaction time: 3 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

Notes:

- It is required to perform "SD Card Calib" calibration when using a new batch of kits.
- It is suggested to calibrate once for one batch of kits.
- Make sure the test card insertion is correct and complete.

TEST RESULTS

Getein1100 can scan the test card automatically and display the result on the screen. Please follow the procedure in user manual of Getein1100 for result printing. For additional information, please refer to the user manual of Getein1100.

EXPECTED VALUE

The expected normal value for mAlb was determined by testing samples from 500 apparently healthy individuals. The 95th percentile of the concentration for mAlb is 20.0 mg/L. (The probability that value of a normal person below 20.0 mg/L is 95%.) It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	10.0~200.0 mg/L
Lower Detection Limit	≤10 mg/L
Within-Run Precision	≤10%
Between-Run Precision	≤15%

Method Comparison:

The assay was compared with OLYMPUS AU5400 analyzer

and its matching Randox mAlb test kits with 200 urine samples (62 positive samples and 138 negative samples). The correlation coefficient (r) is 0.984.

LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Samples containing interferents may influence the results. The table below listed the maximum allowance of these potential interferents.













Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	10 g/L	10 g/L	100 g/L

REFERENCES

- Cöl M, Ocaktan E, Ozdemir O, et al. Microalbuminuria: prevalence in hypertensives and diabetics. *Acta Med Austriaca*. 2004, 31(1):23-29.
- McTaggart MP, Price CP, Pinnock RG, et al. The diagnostic accuracy of a urine albumin - creatinine ratio point-of-care test for detection of albuminuria in primary care. *Am J Kidney Dis*. 2012, 60(5):787-794.
- Denis Sviridov, Glen L. Hortin. Urine albumin measurement: Effects of urine matrix constituents. *Clinica Chimica Acta*. 2009, 404(2):140-143.
- Reboldi G, Gentile G, Angeli F, et al. Microalbuminuria and hypertension. *Minerva Med*. 2005, 96(4):261-75.
- EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
- EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on mAlb Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing mAlb Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF10-S-01



Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508

Fax: +86-25-68568500

E-mail: tech@getein.com.cn

overseas@getein.com.cn

Website: www.bio-GP.com.cn



NT-proBNP Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1002

Getein1600: Cat.# IF2002

User Manual

INTENDED USE

NT-proBNP Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of N-terminal B-type natriuretic peptide precursor (NT-proBNP) in serum, plasma or whole blood. This test is used as an aid in the clinical diagnosis, prognosis and evaluation of Heart Failure (HF).

SUMMARY

N-terminal B-type natriuretic peptide precursor (NT-proBNP) is secreted from the left cardiac ventricle in response to volume and pressure overload. It's an inactive N-terminal fragment that split from BNP prohormone. NT-proBNP can be used to evaluate heart contractile, diastolic dysfunction, and ventricular segmental wall motion coordination. Besides, it has high sensitivity and negative predictive value (>97%). As a gold standard recommended by the European Society of Cardiology, American Heart Association, and American College of Cardiology for the diagnosis and prognosis of heart failure, NT-proBNP is used to indicate heart failure patient at the early stage, determine HF risk levels, monitor medical efficiency of HF drug, evaluate prognosis of HF patient and to distinguish dyspnea that caused by HF from other diseases. Furthermore, NT-proBNP is a risk assessment indicator for Acute Coronary Syndrome.

PRINCIPLE

The test uses an anti-human NT-proBNP monoclonal antibody conjugated with fluorescence latex and an anti-human NT-proBNP polyclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human NT-proBNP monoclonal antibody binds with the NT-proBNP in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human NT-proBNP polyclonal antibody. The fluorescence

intensity of the test line increases in proportion to the amount of NT-proBNP in sample.

Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of NT-proBNP in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

- A kit for Getein1100 contains:
 - Getein NT-proBNP test card in a sealed pouch with desiccant 25
 - Disposable pipet 25
 - Whole blood buffer 1
 - SD card 1
 - User manual 1
- A kit for Getein1600 contains:
 - Sealed cartridge with 24/48 Getein NT-proBNP test cards 2
 - User manual 1
 - Package specifications:
 - 2×24 tests/kit, 2×48 tests/kit
 - Materials required for Getein1600:
 - Sample diluent 1
 - Box with pipette tips 1
 - Mixing plate 1
- Sample diluent/Whole blood buffer composition:
 - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
- A test card consists of:
 - A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human NT-proBNP monoclonal antibody, the test line is coated with another anti-human NT-proBNP polyclonal antibody and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Components from different batches must not be interchanged.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months.

Store the sample diluent/whole blood buffer at 2~8°C for better results.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card.
- Do not reuse the pipet.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- Carefully read and follow the manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- This test can be used for *serum, plasma and whole blood samples*. *Heparin and sodium citrate* should be used as the anticoagulant for plasma and whole blood. Samples should be free of hemolysis.
- Suggest using serum or plasma for better results.
- Serum or plasma can be used directly. For whole blood sample, one drop of whole blood buffer must be added before testing.
- If testing will be delayed, serum and plasma samples may be stored up to 1 day at 2~8°C or stored at -20°C for 3 months before testing (whole blood sample may be stored up to 3 days at 2~8°C).
- Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples.
- SAMPLE VOLUME (for Getein1100): 100 µl.**

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.

For Getein1100:

3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
4. On the main interface of Getein1100, press "ENT" button to enter testing interface.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver 100 µl of sample (or 3~4 drops of sample when using disposable pipet) into the sample port on the test card (for whole blood sample, one drop of whole blood buffer must be added after loading 100 µl sample on the test card).
8. **Reaction time: 10 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1600:

9. Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD Card Calib" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein1100.
3. Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

EXPECTED VALUE

The expected normal value for NT-proBNP was determined by testing samples from 2,500 apparently healthy individuals. The 95th percentile of the concentration for NT-proBNP is 185 pg/ml and the 97.5th percentile of the concentration for NT-proBNP is 300 pg/ml. Because of the apparent difference of the concentration

of NT-proBNP among different age groups, the reference values of the NT-proBNP are reported in groups. Details refer to Table 1. Clinical diagnosis value: refer to Roche criterion, details see Table 2.

Table 1 NT-proBNP reference value

Age Percentile	≤44	45-54	55-64	65-74	≥75	Statistic analysis
95	98.5	130	215	290	530	185
97.5	116	170	270	350	740	300

Table 2 Standard of excluding/diagnosing HF by NT-proBNP

Age	<50	50-75	≥75	Diagnosis of HF
NT-proBNP (pg/ml)	≥450	≥900	≥1800	High probability of HF
	300-450	300-900	300-1800	Low probability of HF, need to combine with other clinical evaluation
	<300	<300	<300	Exclude HF

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	100~35000 pg/ml
Lower Detection Limit	≤100 pg/ml
Within-Run Precision	≤10%
Between-Run Precision	≤15%

Method Comparison:

The assay was compared with Roche MODULAR ANALYTICS E170 and its matching NT-proBNP test kits with 200 serum samples (63 positive samples and 137 negative samples). The correlation coefficient (r) for NT-proBNP is 0.959.

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferents may influence the results. The table below listed the maximum allowance of these potential interferents.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	10 g/L	0.2 g/L













REFERENCES

1. de Lemos JA, McGuire DK, Drazner MH. B-type natriuretic peptide in cardiovascular disease. *Lancet* 2003; 362:316~322.

2. Pfister R, Scholz M, Wielckens K, Erdmann E, Schneider CA. The value of natriuretic peptides NT-pro-BNP and BNP for the assessment of left-ventricular volume and function. A prospective study of 150 patients. *Deutsche medizinische Wochenschrift* (1946) 2002; 127(49):2605.
3. EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
4. EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on NT-proBNP Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing NT-proBNP Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF03-S-02



Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn
overseas@getein.com.cn
Website: www.bio-GP.com.cn



PCT Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1007
Getein1600: Cat.# IF2007

User Manual

INTENDED USE

PCT Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of Procalcitonin (PCT) in serum, plasma or whole blood. The test is used as an aid in the assessment and evaluation of patients suspected of bacterial infection, trauma or shock.

SUMMARY

PCT is a peptide precursor of the hormone calcitonin, the latter being involved with calcium homeostasis. It is composed of 116 amino acids and is produced by parafollicular cells (C cells) of the thyroid and by the neuroendocrine cells of the lung and the intestine.

Measurement of PCT can be used as a marker of severe sepsis and generally grades well with the degree of sepsis, although levels of PCT in the blood are very low. PCT has the greatest sensitivity and specificity for differentiating patients with systemic inflammatory response syndrome (SIRS) from those with sepsis.

PCT levels may be useful to distinguish bacterial infections from nonbacterial infections. It has shown that PCT may help guide therapy and reduce antibiotic use, which can help save on cost of antibiotic prescriptions and drug resistance.

PRINCIPLE

The test uses an anti-human PCT monoclonal antibody conjugated with fluorescence latex. For PCT product, test line 1 was coated with anti-human PCT polyclonal antibody and test line 2 was coated with another anti-human PCT monoclonal antibody. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human PCT monoclonal antibody binds with the PCT in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen

antibody complex is captured on the test line by the other anti-human PCT monoclonal antibody or the polyclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of PCT in sample. Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of PCT in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100 contains:
 - Getein PCT test card in a sealed pouch with desiccant 25
 - Disposable pipet 25
 - Whole blood buffer 1
 - SD card 1
 - User manual 1
2. A kit for Getein1600 contains:
 - Sealed cartridge with 24/48 Getein PCT test cards 2
 - User manual 1
 - Package specifications:
 - 2×24 tests/kit, 2×48 tests/kit
 - Materials required for Getein1600:
 - Sample diluent 1
 - Box with pipette tips 1
 - Mixing plate 1
3. Sample diluent/Whole blood buffer composition:
 - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
4. A test card consists of:
 - A plastic shell and a reagent strip which is composed of a sample pad, fluorescence latex pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human PCT monoclonal antibody, the test line are coated with another anti-human PCT monoclonal antibody and polyclonal antibody, and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days. Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months. Store the sample diluent/whole blood buffer at 2~8°C for better results.

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch or the cartridge is damaged.
5. Do not open pouches or the cartridge until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow the manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for *serum, plasma and whole blood samples*. *Heparin and sodium citrate* should be used as the anticoagulant for plasma and whole blood. Samples should be free of hemolysis.
2. Suggest using serum or plasma for better results.
3. Serum or plasma can be used directly. For whole blood sample, one drop of whole blood buffer must be added before testing.
4. If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2~8°C or stored at -20°C for 6 months before testing (whole blood sample may be stored up to 3 days at 2~8°C).
5. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
6. Do not use heat-inactivated samples.
7. **SAMPLE VOLUME (for Getein1100): 100 µl.**

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.

For Getein1100:

3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD Card Calib" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
 4. On the main interface of Getein1100, press "ENT" button to enter testing interface.
 5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
 6. Put the test card on a clean table, horizontally placed.
 7. Using sample transfer pipette, deliver 100 µl of sample (or 3-4 drops of sample when using disposable pipet) into the sample port on the test card (for whole blood sample, one drop of whole blood buffer must be added after loading 100 µl sample on the test card).
 8. **Reaction time: 15 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.
- ### For Getein1600:
9. Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
 10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD Card Calib" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein1100.
3. Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

EXPECTED VALUE

The expected normal value for PCT was determined by testing samples from 500 apparently healthy individuals. The 99th percentile of the concentration for PCT is 0.1 ng/ml. (The probability that value of a normal person below 0.1 ng/ml is 99%.) The table below comes from the research of ACCP/SCCM (American College of Chest Physicians/Society of Critical Care

Medicine), showing the PCT value and its clinical meaning^[4]:

PCT concentration	Clinical significance
< 0.5 ng/ml	Local bacterial infection is possible, systemic infection (sepsis) is not likely.
≥ 0.5 and < 2.0 ng/ml	Systemic infection (sepsis) is possible, a moderate risk of severe sepsis and/or septic shock.
≥ 2.0 ng/ml	Systemic infection (sepsis) is likely, a high risk of severe sepsis and/or septic shock.

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.1~50.0 ng/ml
Lower Detection Limit	≤0.1 ng/ml
Within-Run Precision	≤10%
Between-Run Precision	≤15%

Method Comparison:

The assay was compared with Roche MODULAR ANALYTICS E170 automatic immunoassay system and its matching PCT test kits with 200 serum samples (68 positive samples and 132 negative samples). The correlation coefficient (r) for PCT is 0.983.

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferent may influence the results. The table below lists the maximum allowance of these potential interferents.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	10 g/L	0.2 g/L













REFERENCES

1. Balci C, Sungurtekin H, Gürses E, Sungurtekin U, Kaptanoglu B. Usefulness of procalcitonin for diagnosis of sepsis in the intensive care unit, Crit Care. 2003 February 7 (1):85-90.
2. Schuetz P, Christ-Crain M, Thomann R, et al. Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: the ProHOSP randomized controlled trial. JAMA. Sep 9 2009; 302(10):1059-66.
3. Briel M, Schuetz P, Mueller B, et al. Procalcitonin-guided antibiotic use vs a standard approach for acute respiratory

- tract infections in primary care. Arch Intern Med. Oct 13 2008; 168(18):2000-7; discussion 2007-8.
4. Meisner M. Procalcitonin (PCT) - A New innovative infection parameter. Biochemical and clinical aspects. Thieme Stuttgart, New York 2000, ISBN: 3-13-105503-0.
5. EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
6. EN ISO 18113-2:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on PCT Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing PCT Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF06-S-02



Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn
overseas@getein.com.cn
Website: www.bio-GP.com.cn



PRL

Fast Test Kit

(Immunofluorescence Assay)

IF1048 for Getein1100
 IF5048 for Getein1160
 IF3048 for Getein1180
 IF4048 for Getein1200
 IF2048 for Getein1600

REF

User Manual

INTENDED USE

PRL Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of prolactin (PRL) in serum or plasma samples. This test can be used as an aid in the diagnosis of male and female infertility and pituitary dysfunction, monitoring of male and female gonadal disorders and management of amenorrhea and galactorrhea.

SUMMARY

Prolactin is a single chain polypeptide of 199 amino acids and a molecular weight of approximately 23,000 Daltons. It is produced by the anterior pituitary and its secretion is regulated physiologically by inhibitory and releasing factors of the hypothalamus.

The major physiologic action of PRL is the initiation and maintenance of lactation in women. In adults, basal circulating prolactin is present in concentrations up to 30 µg/L. During pregnancy and postpartum lactation, serum prolactin can increase 10 to 20 times. Exercise, stress and sleep also cause transient increases in prolactin levels.

Hyperprolactinemia has been established as a common cause of infertility and gonadal disorders in men and women. Prolactin has been shown to inhibit the secretion of ovarian steroids and to interfere with follicle maturation and the secretion of LH and FSH in the human female.

PRINCIPLE

This test uses an anti-human PRL monoclonal antibody I conjugated with fluorescence latex and another anti-human PRL monoclonal antibody II coated on the test line. After

sample is applied to the test strip, the fluorescence latex-labelled anti-human PRL monoclonal antibody I binds with the PRL antigen in sample and forms marked antibody-antigen complex. The complex moves to the detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human PRL monoclonal antibody II. The fluorescence intensity of the test line increases in proportion to the amount of PRL antigen in the sample.

Insert test card into Getein1100/Getein1160/Getein1180 Immunofluorescence Quantitative Analyzer/Automatically inserted by Getein1200/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100, Getein1160, Getein1180, Getein1200 and Getein1600), the concentration of PRL antigen in the sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 and available for downloading. Results can be easily transmitted to a laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100/Getein1160/Getein1180 contains:

Package specifications: 25 tests/kit, 10 tests/kit

- 1) Getein PRL test card in a sealed pouch with desiccant
- 2) Disposable pipet
- 3) User manual: 1 piece/kit
- 4) SD card: 1 piece/kit

2. A kit for Getein1200/Getein1600 contains:

Package specifications: 2×24 tests/kit, 2×48 tests/kit

- 1) Sealed cartridge with 24/48 Getein PRL test cards
- 2) User manual: 1 piece/kit

Materials required for Getein1200/Getein1600:

- 1) Box with pipette tips: 96 tips/kit
- 2) Mixing plate: 1 piece/kit

3. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with fluorescence latex-labelled PRL antibody I, the test line is coated with another PRL antibody II and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
 Getein1180 Immunofluorescence Quantitative Analyzer
 Getein1600 Immunofluorescence Quantitative Analyzer
 Getein1160 Immunofluorescence Quantitative Analyzer
 Getein1200 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test kit at 4~30°C with a valid period of 24 months. Use the test card for Getein1100/Getein1160/Getein1180 within one hour once the foil pouch is opened. For test card of Getein1200/Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Do not use the kit beyond the expiration date.
3. Do not use the test card if the foil pouch is damaged.
4. Do not open pouches until ready to perform the test.
5. Do not reuse the test card.
6. Do not reuse the pipet.
7. Handle all specimens as potentially infectious. Proper handling and disposal methods should follow by local regulations.
8. Carefully read and follow the manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for **serum and plasma samples**. **Heparin, sodium citrate and EDTA** can be used as the anticoagulant for plasma. Samples should be free of hemolysis.
2. The test should performed within 4 hours after blood collection.
3. If testing is delayed, serum and plasma samples may be stored up to 7 days at 2~8°C or stored at -20°C for 6 months before testing.

4. Refrigerated or frozen sample should be reached to room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
5. Do not use heat inactivated or hemolysis samples.
6. Sample volume: (**for Getein1100/Getein1160/Getein1180**): **100 µL**

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.
For Getein1100:
 1. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
 2. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
 3. Put the test card on a clean table, horizontally placed.
 4. Using sample transfer pipette, deliver **100 µL** of sample into the sample well on the test card.
 5. Reaction time: **15 minutes**. Insert the test card into Getein1100 and press "ENT" button or click on "Start" icon (for Android Getein1100) after reaction time is elapsed. The result will be shown on the screen and printed automatically.
For Getein1160/Getein1180 :
 1. Confirm SD card lot No.in accordance with test kit lot No.. Perform "SD card"calibration when necessary.
 2. Enter testing interface of Getein1160/Getein1180 .
 3. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
 4. Put the test card on a clean table, horizontally placed.
 5. Using sample transfer pipette, deliver 100 µL of sample into the sample well on the test card.
 6. Insert the test card into Getein1160/Getein1180 immediately after sample loading. The analyzer will count down the reaction time (15 minutes) and automatically test the card after reaction time is elapsed. The result will be shown on the screen and printed automatically.
For Getein1200/Getein1600:
 1. Each cartridge for Getein1200/Getein1600 contains a

- specific RFID card which can calibrate automatically.
- Place the sample diluent at the correct position in Getein 1200/Getein1600.
 - Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1200/Getein1600 will do the testing and print the result automatically.

Notes:

- It is required to perform “SD card” calibration when using a new batch of kits.
- It is suggested to calibrate once for one batch of kits for Getein1100/Getein1160/Getein1180.
- Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

- Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1160/Getein1180/Getein1200/Getein1600.
- Due to different methodologies or antibody specificity, there may be deviations between the test results of different manufacturers, so they can't be compared directly.

EXPECTED VALUE

The expected reference range for PRL was determined by testing serum samples from apparently 355 healthy adult males and females.

PRL expected range:

Reference Group		N	Median	Range(ng/mL)
Male		146	5.53	2.64-13.13
Female	Premenopausal (< 50 years)	121	8.28	3.34-26.72
	Postmenopausal (> 50 years)	88	6.20	2.74-19.64

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.50 ng/mL~200.00 ng/mL
Lower Detection Limit	≤ 0.50 ng/mL

Within-Run Precision	≤ 10%
Between-Run Precision	≤ 15%

LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Interferents in samples may influence the results. The table below listed the maximum allowance of these potential interferents.














Interferent	Hemoglobin	Triglyceride	Bilirubin	Cholesterol
Concentration (Max)	500 mg/dL	1800 mg/dL	10 mg/dL	400 mg/dL

REFERENCES

- Zhang SY, Li X, Chen XH, et al. Significant associations between prolactin gene polymorphisms and growth traits in the channel catfish (*Ictalurus punctatus* Rafinesque, 1818) core breeding population. *Meta Gene*, 2018.
- Bo-Young Yun, Chunghee Cho, Byung-Nam Cho. Differential activity of 16K rat prolactin in different organic systems [J]. *Animal Cells and Systems*, 2019, 23 (2).
- Cacho ND, Butjosa A, Cuadras D, et al. Prolactin levels in Drug-Naïve First Episode nonaffective Psychosis Patients compared with healthy controls. Sex differences [J]. *Psychiatry Research*, 2019.
- Mohammad JA, Friedrich P. Prolactin and Prolactin-inducible protein (PIP) in the pathogenesis of primary acquired nasolacrimal duct obstruction (PANDO) [J]. *Medical Hypotheses*, 2019, 125.
- Satoshi F, Motoshi K, Kishio T, et al. Orexin A modulates prolactin production by regulating BMP-4 activity in rat pituitary lactotrope cells [J]. *Peptides*, 2019, 113.
- Siew H, Papillon G, David J, et al. Elevated Prolactin during Pregnancy Drives a Phenotypic Switch in Mouse Hypothalamic Dopaminergic Neurons [J]. *Cell Reports*, 2019, 26 (7).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on PRL Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more detail in the European Standard EN ISO 15223-1: 2021.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		<i>In vitro</i> diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/ European Union
	CE mark		Do not use if package is damaged and consult <i>instructions for use</i>
	Catalogue number		

Thank you for purchasing PRL Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF55-S-07



Getein Biotech, Inc.
 Add.: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.getein.com



CMC Medical Devices & Drugs S.L.
 Add.: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain
 Tel: +34951214054



Prog Fast Test Kit (Immunofluorescence Assay)

IF1071 for Getein 1100
IF5071 for Getein 1160
IF5071 for Getein 1180
IF2071 for Getein 1600
IF4071 for Getein 1200

REF

Instructions for Use

INTENDED USE

Prog Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of progesterone in human serum and plasma samples. The test is used clinically as an aid in the diagnosis of diseases related to abnormal progesterone levels, such as threatened miscarriage and luteal phase defect. For professional and laboratory use only.

SUMMARY

Progesterone (P) is a steroid hormone secreted by the corpus luteum of the ovary, and its level varies with the physiological cycle, which is important for ovulation detection, luteal function, normal pregnancy, and placental function. Progesterone is secreted by the ovaries and adrenal glands in non-pregnant women, and the level of progesterone is very low in the follicular phase. After ovulation, the corpus luteum of the ovary produces a large amount of progesterone, which peaks 8-9 days after ovulation and then declines rapidly to follicular phase levels about 4 days before the next menstrual period. If the level of progesterone in the luteal phase is significantly lower than the physiologic value, it suggests possible luteal insufficiency. In early pregnancy, progesterone is secreted by the corpus luteum of pregnancy and serum progesterone levels continue to rise. If the corpus luteum secretes low levels of progesterone during this period, preeclampsia is likely to occur. After 8-10 weeks of pregnancy, the placenta replaces the corpus luteum as the main secretory gland of progesterone, and at this time, the level of progesterone is stable and gradually increases, if the serum progesterone level drops significantly during this period, suggesting that the placenta may be functionally impaired. Progesterone levels will be higher in late pregnancy.

During a woman's pregnancy, progesterone provides support and security for the early growth and development of the fetus. Progesterone levels have a close relationship with the activity of the pregnancy, and low levels of progesterone often indicate an abnormal pregnancy, and can even lead to miscarriage. Determination of progesterone in the blood is clinically important for regulating the female menstrual cycle, detecting luteal function and maintaining pregnancy.

PRINCIPLE

Prog Fast Test Kit (Immunofluorescence Assay) is a lateral flow immunoassay in a competitive design. After the sample has been applied to the test strip, the fluorescence-labelled progesterone monoclonal antibody binds with the progesterone in sample and forms a marked antigen-antibody complex. Meanwhile, the uncombined fluorescence-labelled progesterone monoclonal antibody binds with the progesterone antigen on the test line. The fluorescence intensity of the test line decreases in proportion to the amount of progesterone in sample. Fluorescent signals intensity can be analyzed by applicable device thus the progesterone in sample be detected quantitatively.

APPLICABLE DEVICE

Getein 1100 Immunofluorescence Quantitative Analyzer
Getein 1160 Immunofluorescence Quantitative Analyzer
Getein 1180 Immunofluorescence Quantitative Analyzer
Getein 1200 Immunofluorescence Quantitative Analyzer
Getein 1600 Immunofluorescence Quantitative Analyzer

CONTENTS

Materials provided	Getein 1100/Getein 1160/ Getein 1180		Getein 1200/ Getein 1600	
	10 T/kit	25 T/kit	2*24 T/kit	2*48 T/kit
Prog test card*	10 pcs	25 pcs	24 test cards in 1 cartridge, and 2 cartridges in 1 box	48 test cards in 1 cartridge, and 2 cartridges in 1 box
Disposable pipet	10 pcs	25 pcs	/	/
Sample diluent**	10 pcs	25 pcs	1 box	1 box
Instructions for use	1 pc	1 pc	1 pc	1 pc
SD card	1 pc	1 pc	1 pc in each cartridge	1 pc in each cartridge

*Prog test card

A test card main consists of: Fluorescence-labelled prog monoclonal antibody and prog antigen.

**Sample diluent

(1) Sample diluent for Getein 1100/Getein 1160/Getein 1180 is 0.2 mL contained in each tube consists of:

- Sample diluent main contains phosphate buffer (20 mmol/L), Na₂S₂O₃ (<0.1%).

(2) Sample diluent for Getein 1200/Getein 1600 is an independent packing box main consists of:

- Phosphate buffer (20 mmol/L), Na₂S₂O₃ (<0.1%) (25 mL/bottle for Getein 1200, 40 mL/bottle for Getein 1600),

- Box with pipette tips (96 tips/box),

- Mixing plate (1 piece/box)

Note:

- The standard curve data can be written to RFID card in the kit. According to the function of RFID card, we define it as "Standard Curve Data Card", short for "SD Card".
- Do not mix or interchange different batches of kits.

STORAGE AND STABILITY

Realtime stability:

Store the kit at 4~30°C with a valid period of 24 months. The test kits are stable until the expiry date printed on the labels.

In-use stability:

-For the test card of Getein 1100/Getein 1160/Getein 1180: Use the test card within 1 hour once the foil pouch is opened.

-For test card of Getein 1200/Getein 1600: If the cartridge is opened, it could be stable within 24 hours once exposure to air. The valid period after opening is 7 days, it is recommended to put the cartridge back to the foil bag and reseal along the entire edge of zip-seal if not used up.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card or pipet.
- Handle all specimens as potentially infectious. Proper

handling and disposal methods should be followed in accordance with local regulations.

- Carefully read and follow instructions for use to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- This test can be used for **serum and plasma samples**.
- Heparin, EDTA and sodium citrate can be used as the anticoagulant for plasma samples.
- It is recommended to test the sample within 4 hours after collection. If testing is delayed, serum and plasma samples are stable for 5 days when stored at 2~8°C and 6 months when stored at -20°C.
- Refrigerated or frozen sample should reach room temperature before testing. Avoid multiple freeze-thaw cycles.
- Do not use heat-inactivated samples or hemolysis samples.
- Sample volume (**Getein 1100/Getein 1160/Getein 1180**): 100 µL.

TEST PROCEDURE

- User must carefully read and operate in strict accordance with the instructions for use before testing, otherwise reliable results cannot be guaranteed.
- Test kit and sample should be brought to room temperature before testing.

For Getein 1100:

- Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
- Select the corresponding "Sample" on the analyzer according to the sample type (see the user manual of analyzer for details).
- Remove the test card from the sealed pouch before use and put the test card on a clean table, horizontally placed.
- Use disposable pipet or pipette to drop 100 µL of sample into one tube of sample diluent, mix gently and thoroughly, then add 100 µL of sample mixture into the sample well on the test card.
- Reaction time: **15 minutes**. Insert the test card into Getein 1100 and click on "Start" icon after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein 1160/Getein 1180:

1. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
2. Select the corresponding "Sample" on the analyzer according to the sample type (see the user manual of analyzer for details).
3. Remove the test card from the sealed pouch immediately before use and put the test card on a clean table, horizontally placed.
4. Use disposable pipet or pipette to drop 100 µL of sample into one tube of sample diluent, mix gently and thoroughly, then add **100 µL** mixture into the sample well on the test card.
5. Insert the test card into Getein 1160/Getein 1180 immediately after sample loading. The analyzer will count down the reaction time (15 minutes) and automatically test the card after reaction time is elapsed. The result will be shown on the screen and displayed automatically.

For Getein 1200/Getein 1600:

1. Each cartridge for Getein 1200/Getein 1600 contains a specific RFID card which can calibrate automatically.
2. Put the sample diluent at the correct position in Getein 1200/Getein 1600.
3. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein 1200/Getein 1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein 1100/Getein 1160/Getein 1180.
3. Make sure the insertion of test card and the sample are correct and complete.

RESULTS

Getein 1100/Getein 1160/Getein 1180/Getein 1200/Getein 1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein 1100/Getein 1160/ Getein 1180/Getein 1200/Getein 1600.

Prog Fast Test Kit (Immunofluorescence Assay) results are provided in ng/mL.

LIMITATIONS

1. The test results of this reagent are for clinical reference only, and cannot be used as the basis for diagnosis or exclusion of cases alone additional tests should be performed accordingly.
2. Some substances in blood as listed below may interfere with the test and cause erroneous results. The maximum allowance concentration of each is as follows:

Interferent	Concentration (Max)
Triglyceride	25 g/L
Bilirubin	0.1 g/L

EXPECTED VALUE

The expected normal value for Progesterone was determined by testing blood samples from apparently healthy individuals. Reference range of Progesterone:

Group		n	95% Reference range (ng/mL)
Healthy men		200	0.15-1.97
Healthy women	Follicular phase	121	0.34-1.52
	Luteal phase	58	5.20-18.76
	Post menopause	156	0.12-0.76
Healthy pregnant women	1 st trimester	135	4.73-50.21
	2 nd trimester	49	19.27-45.34

It is recommended that each laboratory determine the applicability of the reference ranges through experimentation and establish their own laboratory-specific reference ranges if necessary.

PERFORMANCE CHARACTERISTICS

1. Measuring Range 0.10–40.00 ng/mL
2. Limit of Detection ≤0.10 ng/mL
3. Within-Run Precision ≤10%
4. Between-lot Precision ≤15%















REFERENCES

1. Julia S, R J W, Pawel B, et al. Role of progesterone and

- progesterin therapy in threatened abortion and preterm labour[J]. *Frontiers in bioscience: a journal and virtual library*, 2008, 13:1981-90.
2. Julia S. Progesterone orchestrates foetal-maternal interactions in pregnancy[J]. *Journal of Reproductive Immunology*, 2023, 159.
 3. Clinical and Laboratory Standards Institute. Protocols for determination of limits of quantitation; approved guideline-second edition, EP17-A, CLSI,2004.
 4. Clinical and Laboratory Standards Institute. Evaluation of precision performance of quantitative measurement method; approved guideline-second edition, EP17-A, CLSI,2004.
 5. National Committee for Clinical Laboratory. Method comparison and bias estimation using patient samples; approved guideline. EP9-A2, NCCLS, 2002.


DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on Prog Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2021.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		In vitro diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/ European Union
	CE mark		Do not use if package is damaged and consult instructions for use
	Catalogue number		Caution

Thank you for using Prog Fast Test Kit (Immunofluorescence Assay). Please read the instructions for use carefully before operating to ensure proper use. Please report any product problems or adverse events to the below manufacture or authorized representative in the European Community in time.

 Getein Biotech, Inc.
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.getein.com

 CMC Medical Devices & Drugs S.L.
 Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain
 Tel: +34951214054



RF Fast Test Kit (Immunofluorescence Assay)

IF1075 for Getein1100
IF5075 for Getein1160
IF3075 for Getein1180
IF4075 for Getein1200
IF2075 for Getein1600

REF

User Manual

INTENDED USE

RF Fast Test Kit (Immunofluorescence Assay) is intended for in vitro quantitative determination of rheumatoid factor (RF) in serum, plasma or whole blood samples. The test is designed to aid in the diagnosis of autoimmune diseases, such as rheumatoid arthritis (RA).

SUMMARY

Rheumatoid factor is an autoantibody targeting the Fc fragment of human or animal denatured IgG molecule. RF mainly includes four types of IgM, IgG, IgA and IgE, and IgM is the main type of RF. Under the direct stimulation of denatured IgG or Epstein-Barr virus, B cells in patients with rheumatoid arthritis will synthesize RF in large quantities. On the contrary, in healthy people, there are few clones of B cells that produce RF, and the soluble factors secreted by monocytes can inhibit the production of RF, which is generally difficult to be detected. RF is mainly used in the clinical diagnosis of RA. RF has a positive detection rate of 80% in RA patients. Positive RF is one of the criteria for RA classification by the American College of Rheumatology, but positive RF is not the sole basis for the diagnosis of RA.

PRINCIPLE

The test kit adopts a double-antigen sandwich method to quantitatively detect the concentration of RF in human serum, plasma and whole blood samples. After the sample has been applied to the test card, the fluorescence latex-labelled RF antigen binds with RF in

sample and forms a marked antigen-antibody complex. The complex moves to the detection area by capillary action, then it is captured by RF antigen coated on the detection area of nitrocellulose membrane, forming a double- antigen complex. The complex generates a fluorescent signal and the intensity increases in proportion to the amount of RF in sample. Then insert test card into Getein1100/Getein1160/Getein1180 Immunofluorescence Quantitative Analyzer/ automatically inserted by Getein1200/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100, Getein1160, Getein1180, Getein1200 and Getein1600), the concentration of RF in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100/Getein1160/Getein1180 contains:

Package specifications: 25 tests/kit, 10 tests/kit

- 1) Getein RF test card in a sealed pouch with desiccant
- 2) Capillary pipet
- 3) Sample diluent
- 4) User manual: 1 piece/kit
- 5) SD card: 1 piece/kit

2. A kit for Getein1200/Getein1600 contains:

Package specifications: 2x24 tests/kit, 2x48 tests/kit

- 1) Sealed cartridge with 24/48 Getein RF test cards
- 2) User manual: 1 piece/kit

Materials required for Getein1200/Getein1600:

- 1) Sample diluent: 1 bottle/kit
- 2) Box with pipette tips: 96 tips/kit
- 3) Mixing plate: 1 piece/kit

3. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad ,nitrocellulose membrane (one end of the membrane is coated with fluorescence latex-labelled RF antigen, the test line is coated with RF antigen and the control line is coated with polyclonal mouse anti human IgG antibody),absorbent paper and liner.

4. Sample diluent composition:

Phosphate buffered saline, proteins, detergent,

preservative, stabilier.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1160 Immunofluorescence Quantitative Analyzer
Getein1180 Immunofluorescence Quantitative Analyzer
Getein1200 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test kit at 4~30°C with a valid period of 24 months.

Use the test card for Getein1100/Getein1160/Getein1180 within one hour once the foil pouch is opened. For test card of Getein1200/Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

PRECAUTIONS

1. For in vitro diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch or the cartridge is damaged.
5. Do not open pouches or the cartridge until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. **Serum, plasma or whole blood** samples can be used for the test. Other body fluids and samples may not give accurate results. Samples should be

free of hemolysis.

2. Venous blood should be collected under aseptic conditions; serum or plasma is preferred for testing.
3. Heparin, sodium citrate or EDTA can be used as the anticoagulant for plasma and whole blood samples.
4. The test should be performed at room temperature (15~30°C) within 4 hours after sample collection.
5. If testing is delayed, serum and plasma samples may be stored up to 7 days at 2~8°C and 3 months at -20°C before testing. Whole blood samples should not be frozen and can be stored at 2~8°C for 3 days. Do not heat inactivated samples or use hemolyzed blood samples.
6. Refrigerated or frozen sample should be reached to room temperature (15~30°C) before testing. Frozen samples must be completely thawed, rewarmed and evenly mixed. Serum and plasma samples can freeze and thaw twice at most. Avoid multiple freeze-thaw cycles.
7. **SAMPLE VOLUME (for Getein1100/Getein1160/Getein1180): 10µL.**

TEST PROCEDURE

1. Collect specimens according to user manual.
 2. Test card, sample and reagent should reach to room temperature before test.
- For Getein1100:
1. Confirm SD card lot No. in accordance with test kit lot No. Perform "SD card" calibration when necessary.
 2. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
 3. Put the test card on a clean table, horizontally placed.
 4. Using sample transfer pipette, deliver **10 µL** of sample into one tube of sample diluent, mix gently and thoroughly. Then drop **100µL** of sample mixture into the sample well on the test card.
 5. **Reaction time: 10 minutes.** Insert the test card into Getein1100 and click on "Start" icon (for Android Getein1100) after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1160/Getein1180:

1. Confirm SD card lot No.in accordance with test kit lot No..Perform "SD card" calibration when necessary.
2. Enter testing interface of Getein1160/Getein1180.
3. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
4. Put the test card on a clean table, horizontally placed.
5. Using sample transfer pipette, deliver **10 µL** of sample into one tube of sample diluent, mix gently and thoroughly. Then drop **100 µL** of sample mixture into the sample well on the test card.
6. Insert the test card into Getein1160/Getein1180 immediately after sample loading. The analyzer will count down the reaction time (10 minutes) and automatically test the card after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1200/Getein1600:

1. Each cartridge for Getein1200/Getein1600 contains a specific RFID card which can calibrate automatically.
2. Put the sample diluent at the correct position of Getein1200/Getein1600.
3. Place samples in the designed area of the sample holder,insert the holder and select the right test item, Getein1200/Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein1100/Getein1160/Getein1180.
3. Make sure the test card and the sample insertion are correct and complete.

TEST RESULTS

Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1160/Getein1180/Getein1200/Getein1600.

Others: Samples whose concentration exceeds the upper limit should be diluted no more than 4 times.

EXPECTED VALUE

The expected normal value for RF is determined by testing samples from 282 apparently healthy individuals. The upper 97.5th percentile value is 15.9 IU/mL.

It is recommended that each laboratory determine the applicability of the reference value through experiments and establish its own reference ranges if necessary.

PERFORMANCE CHARACTERISTICS

Measuring Range	10.0-640.0 IU/mL
Lower Detection Limit	≤10.0 IU/mL
Within-run Precision	≤10%
Between-run Precision	≤15%

LIMITATIONS

1. Bilirubin and triglyceride in the sample may interfere with the test results, and the maximum allowable concentrations are 0.2 g/L and 10 g/L respectively.
2. The test results of this kit are for clinical reference only, and should not be used as the sole criteria for clinical diagnosis. It is recommended to conduct a comprehensive analysis on the condition in combination with symptoms/signs, history and other laboratory tests.

REFERENCES














1. Eberhardt KB, Truedson L, Pettersson H, et al. Disease activity and joint damage progression in early rheumatoid arthritis: relation to IgG, IgA and IgM rheumatoid factor [J]. Ann Rheum Dis, 1990, 49(11):906.
2. Kolarz B, Podgorska D, Podgorski R. Insights of rheumatoid arthritis biomarkers.[J]. Biomarkers : biochemical indicators of exposure, response, and susceptibility to chemicals,2020: 1-34.
3. Nithya L, Jeremy S, Lauren J, et. al. Combination of anticitrullinated protein antibodies and rheumatoid factor is associated with increased systemic inflammatory mediators and more rapid progression from preclinical to clinical rheumatoid arthritis[J]. Clinical Immunology,2018,195:

119-126.

4. Veerle I, Xavier B, Daniël B, Ellen DL. Prevalence and clinical correlates of rheumatoid factor and anticitrullinated protein antibodies in patients with idiopathic inflammatory myopathy[J]. RMD Open,2018,4(2).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on the test kit are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2021.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		<i>In vitro</i> diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/ European Union
	CE mark		Do not use if package is damaged and consult instructions for use
	Catalogue number		

Thank you for purchasing RF Fast Test Kit (Immunofluorescence Assay).

Please read this user manual carefully before operating to ensure proper use.

Version: WIF82-S-05



Getein Biotech, Inc.
 Add.: No.9 Bofu Road, Luhe District, Nanjing,
 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.getein.com



CMC Medical Devices & Drugs S.L.
 Add.: C/ Horacio Lengo N° 18, CP 29006,
 Málaga, Spain



T3

Fast Test Kit

(Immunofluorescence Assay)

User Manual

REF IF1022 for Getein1100
IF2022 for Getein1600

INTENDED USE

T3 (Triiodothyronine) Fast Test Kit (Immunofluorescence Assay) is intended for in vitro quantitative determination of T3 in serum and plasma. It can be used in the monitoring of hyperthyroidism and hypothyroidism, and also used as an aid in the functional diagnosis of thyroidea.

SUMMARY

The thyroid gland exerts powerful and essential regulatory influences on growth, differentiation, cellular metabolism and general hormonal balance, as well as on the maintenance of metabolic activity and the development of the skeletal and organ system. The hormones thyroxine (T4) and triiodothyronine (T3) circulate in the blood stream, mostly bound to the plasma protein, thyroxine binding globulin (TBG). The concentration of T3 is much less than that of T4, but its metabolic potency is much greater.

T3 is produced by the thyroid and secreted in response to TSH. T3 determination is an important factor in the diagnosis of thyroid disease. Its measurement has uncovered a variant of hyperthyroidism in thyrotoxic patients with elevated T3 levels and normal T4 levels. An increase in T3 without an increase in T4 is frequently a forerunner of recurrent thyrotoxicosis in previously treated patients. In other patients, euthyroidism attributable to normal T3, although their T4 values are subnormal.

In women, T3 levels are elevated during pregnancy, during estrogen treatment, and contraceptive hormone therapy. When T3 levels parallel TBG increases in a manner analogous to T4 levels, these changes are not reflection of altered thyroid status.

PRINCIPLE

The test uses an anti-human T3 monoclonal antibody

conjugated with fluorescence latex coated on fluorescent pad and another antibody conjugated with fluorescein isothiocyanate (FITC) coated on the test line, also a T3 antigen conjugated with FITC is coated on the junction of nitrocellulose membrane and sample pad. After the sample has been applied to the test strip, the T3 in sample and the T3 antigen conjugated with FITC will compete the fluorescence latex-labelled T3 antibody and form marked T3-FITC-T3 antibody and T3-T3 antibody, respectively. These complexes move to the test card detection zone by capillary action. Then marked T3-FITC-T3 antibody complexes will be captured on the test line by antibodies conjugated with FITC. The fluorescence intensity of test line increases in proportion to the amount of T3 in sample.

Then Insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Automatically inserted by Getein1600 Immunofluorescence Quantitative Analyzer (hereafter referred to as Getein1100 and Getein1600), the concentrations of T3 in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100 contains:

Package specifications: 25 tests/box, 10 tests/box

- 1) Getein T3 test card in a sealed pouch with desiccant
- 2) Disposable pipet
- 3) Sample treatment solution A and B
- 4) User manual: 1 piece/box
- 5) SD card/RFID card: 1 piece/box

2. A kit for Getein1600 contains:

Package specifications: 2x24 tests/kit, 2x48 tests/kit

- 1) Sealed cartridge with 24/48 Getein T3 test cards
- 2) User manual: 1 piece/box

Materials required for Getein1600:

- 1) Sample diluent: 1 bottle/box
- 2) Box with pipette tips: 96 tips/box
- 3) Mixing plate: 1 piece/box

3. Sample treatment solution/sample diluent composition:

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

4. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the

membrane is coated with fluorescence latex-labelled anti-human T3 monoclonal antibody and T3 antigen, the test lines is coated with FITC antibody, and the control line C is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months.

Use the test card within 1 hour once the foil pouch is opened. For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample treatment solution/sample diluent at 0~30°C with a valid period of 24 months.

Store the sample treatment solution/sample diluent at 2~8°C for better results.

PRECAUTIONS

1. For in vitro diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch is damaged.
5. Do not open pouches until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for **serum and plasma samples**. **Heparin, EDTA and sodium citrate** can be used as the anticoagulant for plasma. Samples should be free of hemolysis.
2. The test should be performed within 4 hours after whole blood collection.
3. If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2-8 °C or stored at -20 °C for 6 months before testing.
4. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
5. Do not use heat-inactivated samples or hemolysis samples.
6. **SAMPLE VOLUME (for Getein1100): 100 µL**.

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.
For Getein1100:
3. Confirm SD card or RFID card lot No. in accordance with test kit lot No.. Perform "SD card or RFID card" calibration when necessary.
4. Enter testing interface of Getein1100.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver **10 µL** sample treatment B into sample treatment A, and then deliver **100 µL** of sample into sample treatment A. Mixing for 1 minute and then drop **100 µL** of the sample mixture to the test card.
8. **Reaction time: 15 minutes**. Insert the test card into Getein1100 and start test after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1600:

9. Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card or RFID card" calibration

when using a new batch of kits.

2. It is suggested to calibrate once for one batch of kits for Getein1100.
3. Make sure the test card insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

Others: Measuring range of the T3 test kit is 0.30 nmol/L~10.00 nmol/L. Dilute the sample which concentration is higher than the upper limit with calf serum or negative samples, and the dilution ratio should be less than 4 times.

EXPECTED VALUE

The expected normal value for T3 and was determined by testing samples from 391 apparently healthy individuals. The reference range of T3 is 1.30 nmol/L~3.10 nmol/L calculated by using normal distribution methods (95% confidence interval).

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.30~10.00 nmol/L
Lower Detection Limit	≤0.30 nmol/L
Within-Run Precision	≤10%
Between-Run Precision	≤15%

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferent may influence the results. The table below listed the maximum allowance of these potential interferent.













Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25 g/L	0.1 g/L

REFERENCES

1. Spencer C A, LoPresti J S, Patel A, et al. Applications of a new chemiluminometric thyrotropin assay to subnormal measurement. *J Clin Endocrinol Metab.* 1990, 70(2):453-460.
2. Spencer C A, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyrolobulin assays. *Clin Chem.* 1996,42(1):164-173.
3. Koidl J, Hödl H, Schmid MG et al. Enantio-recognition of triiodothyronine and thyroxine enantiomers using different chiral selectors by HPLC and micro-HPLC. *J. Biochem Biophys Methods.* 2008, 70(6):1254-1260.
4. Keffer JH. Preanalytical considerations in testing thyroid function. *Clin. Chem.* 1996, 42(1):125-134.
5. Lindstedt G, Berg G, Jansson S, et al. Clinical use of laboratory thyroid tests and investigations. *J Int Fed Clin Chem.* 1994, 6(4):136-141.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on T3 Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2016

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		In vitro diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing T3 Fast Test Kit (Immunofluorescence Assay).

Please read this user manual carefully before operating to ensure proper use.

 Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508 Fax: +86-25-68568500
E-mail: tech@getein.com.cn overseas@getein.com.cn
Website: www.bio-GP.com.cn



T4

Fast Test Kit

(Immunofluorescence Assay)

User Manual

REF IF1023 for Getein1100
IF2023 for Getein1600

INTENDED USE

T4 (Thyroxine) Fast Test Kit (Immunofluorescence Assay) is intended for in vitro quantitative determination of T4 in serum and plasma. It can be used in the monitoring of hyperthyroidism and hypothyroidism, and also used as an aid in the functional diagnosis of thyroidea.

SUMMARY

The thyroid gland exerts powerful and essential regulatory influences on growth, differentiation, cellular metabolism and general hormonal balance, as well as on the maintenance of metabolic activity and the development of the skeletal and organ system.

T4 is the most commonly measured substance in the assessment of thyroid function. Most of the thyroxine secreted into the bloodstream is bound to a transport protein, thyroxine binding globulin (TBG), and to albumin and pre-albumin. Only less than 1% of thyroxine remains unbound as free T4 in blood. Elevated total thyroxine levels have been associated with hyperthyroidism, a condition with an excess amount of circulating thyroid hormone and decreased total thyroxine levels have been associated with hypothyroidism, a condition with insufficient levels of thyroxine concentration. Primary malfunction of the thyroid gland or any diseases affecting the thyroid-pituitary-hypothalamus system may result in the abnormal thyroxine concentration in blood. Measurement of thyroxine has been one of the most widely used method for evaluation of an individual's thyroid status.

PRINCIPLE

The test uses an anti-human T4 monoclonal antibody conjugated with fluorescence latex coated on fluorescent pad

and another antibody conjugated with fluorescein isothiocyanate (FITC) coated on the test line, also a T4 antigen conjugated with FITC is coated on the junction of nitrocellulose membrane and sample pad. After the sample has been applied to the test strip, the T4 in sample and the T4 antigen conjugated with FITC will compete the fluorescence latex-labelled T4 antibody and form marked T4-FITC-T4 antibody and T4-T4 antibody, respectively. These complexes move to the test card detection zone by capillary action. Then marked T4-FITC-T4 antibody complexes will be captured on the test line by antibody conjugated with FITC. The fluorescence intensity of test line increases in proportion to the amount of T4 in sample.

Then Insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Automatically inserted by Getein1600 Immunofluorescence Quantitative Analyzer (hereafter referred to as Getein1100 and Getein1600), the concentrations of T4 in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100 contains:

Package specifications: 25 tests/box, 10 tests/box

- 1) Getein T4 test card in a sealed pouch with desiccant
- 2) Disposable pipet
- 3) Sample treatment solution A and B
- 4) User manual: 1 piece/box
- 5) SD card/RFID card: 1 piece/box

2. A kit for Getein1600 contains:

Package specifications: 2x24 tests/kit, 2x48 tests/kit

- 1) Sealed cartridge with 24/48 Getein T4 test cards
- 2) User manual: 1 piece/box

Materials required for Getein1600:

- 1) Sample diluent: 1 bottle/box
- 2) Box with pipette tips: 96 tips/box
- 3) Mixing plate: 1 piece/box

3. Sample treatment solution/sample diluent composition:

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

4. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with fluorescence latex-labelled

anti-human T4 monoclonal antibody and T4 antigen, the test lines is coated with FITC antibody, and the control line C is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months.

Use the test card within 1 hour once the foil pouch is opened. For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.

Store the sample treatment solution/sample diluent at 0~30°C with a valid period of 24 months.

Store the sample treatment solution/sample diluent at 2~8°C for better results.

PRECAUTIONS

1. For in vitro diagnostic use only.
2. For professional use only.
3. Do not use the kit beyond the expiration date.
4. Do not use the test card if the foil pouch is damaged.
5. Do not open pouches until ready to perform the test.
6. Do not reuse the test card.
7. Do not reuse the pipet.
8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
9. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for *serum and plasma samples*. *Heparin, EDTA and sodium citrate* can be used as the anticoagulant for plasma. Samples should be free of hemolysis.
2. The test should be performed within 4 hours after whole blood collection.
3. If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2~8℃ or stored at -20℃ for 6 months before testing.
4. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
5. Do not use heat-inactivated samples or hemolysis samples.
6. SAMPLE VOLUME (for *Getein1100*): 25 µl.

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.

For *Getein1100*:

3. Confirm SD card or RFID card lot No. in accordance with test kit lot No.. Perform "SD card or RFID card" calibration when necessary.
4. Enter testing interface of *Getein1100*.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver 10 µl sample treatment B into sample treatment A, and then deliver 25 µl of sample into sample treatment A. Mixing for 1 minute and then drop 100 µl of the sample mixture to the test card.
8. **Reaction time: 15 minutes.** Insert the test card into *Getein1100* and start test after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For *Getein1600*:

9. Each cartridge for *Getein1600* contains a specific RFID card which can calibrate automatically.
10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, *Getein1600* will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card or RFID card" calibration

- when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for *Getein1100*.
3. Make sure the test card insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of *Getein1100/Getein1600*.

Others: Measuring range of the T4 test kit is 5.4 nmol/L~320.0 nmol/L. Dilute the sample which concentration is higher than the upper limit with calf serum or negative samples, and the dilution ratio should be less than 4 times.

EXPECTED VALUE

The expected normal value for T4 and was determined by testing samples from 391 apparently healthy individuals. The reference range of T4 is 59.0 nmol/L~154.0 nmol/L calculated by using normal distribution methods (95% confidence interval).

It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	5.4~320.0 nmol/L
Lower Detection Limit	≤5.4 nmol/L
Within-Run Precision	≤10%
Between-Run Precision	≤15%

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferer may influence the results. The table below listed the maximum allowance of these potential interferer.













Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25 g/L	0.1 g/L

REFERENCES

1. Spencer C A, LoPresti J S, Patel A, et al. Applications of a new chemiluminometric thyrotropin assay to subnormal measurement. *J Clin Endocrinol Metab.* 1990, 70(2):453-460.
2. Spencer C A, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyrolobulin assays. *Clin Chem.* 1996,42(1):164-173.
3. Koidl J, Hödl H, Schmid MG et al. Enantio recognition of triiodothyronine and thyroxine enantiomers by different chiral selectors by HPLC and micro-HPLC. *J Biochem Biophys Methods.* 2008, 70(6):1254-1260.
4. Keffer JH. Preatalytical considerations in testing thyroid function. *Clin. Chem.* 1996, 42(1):125-134.
5. Lindstedt G, Berg G, Jansson S, et al. Clinical use of laboratory thyroid tests and investigations. *J Int Fed Clin Chem.* 1994, 6(4):136-141.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on T4 Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2016.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		In vitro diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing T4 Fast Test Kit (Immunofluorescence Assay).

Please read this user manual carefully before operating to ensure proper use.



Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508 Fax: +86-25-68568500
E-mail: tech@getein.com.cn overseas@getein.com.cn
Website: www.bio-GP.com.cn



Testosterone Fast Test Kit (Immunofluorescence Assay)

REF IF1073 for Getein1100
IF5073 for Getein1160
IF3073 for Getein1180

Instructions for Use

INTENDED USE

Testosterone Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of Testosterone in human serum and plasma. Testosterone tests are used to assist in the diagnosis of female polycystic ovary syndrome and male testosterone insufficiency. For professional and laboratory use.

SUMMARY

Testosterone, one of the main sex hormones in the human body, is a steroid hormone with 19 carbon atoms synthesized from cholesterol. About 90% of male testosterone comes from leydig cells, and the rest is produced in adrenal cortex and other tissues; Testosterone in women synthesized by ovarian stromal cells, portal cells and adrenal cortex reticular zone, account for about 50% of total. The remaining 50% is mainly converted from androstenedione in liver, fat, skin and other tissues. It is the precursor of estradiol synthesis in the ovary, and has a certain role in maintaining the female gonadal function. Testosterone exists in blood in two forms: bound and free. The binding state is mainly bound with sex hormone binding protein, albumin and other proteins, accounting for 97%~99%; Free testosterone accounts for 1%~3%, which is the main form of biological activity *in vivo*. Detection of female testosterone can assist in the diagnosis of polycystic ovary syndrome (PCOS), masculinized tumors in women, congenital adrenal hyperplasia and adrenal tumors. Detection of male testosterone content can assist in the diagnosis of diseases with insufficient testosterone production, such as congenital testicular dysplasia, chromosomal abnormalities, and in addition, it can assist in the diagnosis of hypopituitarism, cirrhosis.

PRINCIPLE

Testosterone Fast Test Kit (Immunofluorescence Assay) is a lateral flow immunoassay in a competitive design. After the

sample has been applied to the test strip, the fluorescence latex-labelled testosterone monoclonal antibody binds with the testosterone in sample and forms a marked antigen-antibody complex. The uncombined fluorescence latex-labelled testosterone monoclonal antibody binds with the testosterone on the test line. The fluorescence intensity of the test line decreases in proportion to the amount of testosterone in sample. Fluorescent signals intensity can be analyzed by applicable device thus the testosterone in sample be detected quantitatively.

APPLICABLE DEVICE

Getein 1100 Immunofluorescence Quantitative Analyzer
Getein 1160 Immunofluorescence Quantitative Analyzer
Getein 1180 Immunofluorescence Quantitative Analyzer

CONTENTS

Materials provided	Getein 1100/ Getein 1160/ Getein 1180	
	10 T/kit	25 T/kit
Testosterone test card	10 pcs	25 pcs
Disposable pipet	10 pcs	25 pcs
Sample diluent 3	10 tube	25 tube
Instructions for use	1 pc	1 pc
SD card	1 pc	1 pc

Sample diluent 3 for Getein 1100/ Getein 1160/ Getein 1180 consists of:

Sample diluent 3 contains Carbonate buffer (50 mmol/L), ProClin™ 300 (0.1%).

A test card consists of:

Fluorescence latex-labelled testosterone monoclonal antibody, Fluorescence latex-labelled Chicken immunoglobulin Y natural protein, testosterone antigen and Goat anti chicken immunoglobulin Y polyclonal antibody.

Note:

- The standard curve data can be written to RFID card in the kit. According to the function of RFID card, we define it as "Standard Curve Data Card", short for "SD Card".
- Do not mix or interchange different batches of kits.

STORAGE AND STABILITY

Realtime stability:

Store the kit at 4~30°C with a valid period of 24 months. The test kits are stable until the expiry date printed on the labels.

In-use stability:

For the test card of Getein 1100/Getein 1160/Getein 1180: Use the test card within 1 hour once the foil pouch is opened.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- Handle all specimens as potentially infectious. The foil bag is nondegradable. Proper handling and disposal methods should be followed in accordance with local regulations.

SPECIMEN COLLECTION AND PREPARATION

- This test can be used for **plasma and serum samples**.
- Heparin, sodium citrate or EDTA can be used as the anticoagulant for plasma samples.
- It is recommended to test the sample within 8 hours after collection. Stable in plasma for 7 days when stored at 2~8°C and 6 months when stored at -20°C.
- Refrigerated or frozen sample should reach room temperature before testing. Avoid multiple freeze-thaw cycles.

CALIBRATION AND TRACEABILITY


- Calibration: The regression equation fitting the concentration value of the working calibrator with the reaction signal value is written into the SD card in advance. Before detection, the SD card is written into the instrument, which can automatically read the calibration curve information in the SD card. During detection, the content of analyte can be calculated by substituting the obtained signal value into the regression equation.

Calibration Frequency: A new calibration is required when using a new reagent lot or a new instrument.

TEST PROCEDURE

- Before use, you must carefully read the instructions for use and operate in strict accordance with the instructions, otherwise reliable results cannot be guaranteed.
- Test kit and sample should be brought to room temperature before testing.
- Confirm SD card lot No. in accordance with test kit lot No. Perform calibration using the SD card when necessary.
- Select the corresponding model on the analyzer according to

- the sample type (see the instructions of analyzer for details).
- Remove the test card from the sealed pouch before use. Horizontally place the test card.
- Deliver 100 μ L of sample into one tube of sample diluent 3 using **disposable pipet**, mix gently and thoroughly. (Samples must be added using the disposable pipet in the kit to avoid incorrect results).




Note 1: Press the top of the disposable pipet to the bottom with your finger during sampling. Ensure that the slit is fully submerged in the sample.

Note 2: Thoroughly press the disposable pipette only once to take a sample, not repeatedly.

Note 3: Insert the disposable pipet into the sample diluent 3 tube to mix the sample by pushing the cap at the top of the disposable pipet for 4-6 times and wait 5~10 minutes.

Note 4: It is recommended to wait for 5-10 minutes after mixing the samples, and the results of early or overtime testing are inaccurate.



Insert the disposable pipet to the bottom to prevent bubbles

8. Deliver the sample mixture by pushing the cap at the top of the disposable pipet and dispense the sample mixture into the sample port "S" on the test card.

For Getein 1100:

Reaction time: **15 minutes**. Insert the test card into Getein 1100 and press "ENT" button (click on "Start" icon for Android Getein 1100) after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein 1160/Getein 1180:

Insert the test card into Getein 1160/Getein 1180 immediately after sample loading. The analyzer will count down the reaction time (15 minutes) and automatically test the card after reaction time is elapsed. The result will be shown on the screen and printed automatically.

RESULTS

Getein1100/Getein1160/Getein1180 can scan the test card automatically and display the result on the screen. For additional information, please refer to the instructions for use of Getein1100/Getein1160/ Getein1180.

Testosterone Fast Test Kit (Immunofluorescence Assay) results are provided in ng/mL.

Results in ng/mL may be converted to nmol/mL as shown with an example below.

Testosterone Fast Test Kit (Immunofluorescence Assay) result as reported by the system (example) 1.00 ng/mL
The reported example result equals: 3.47 nmol/L

Others: Measuring range of the Testosterone Fast Test Kit is 0.10 ~ 16.00 ng/mL. Samples initially outside the measuring range may be diluted with 1% bovine serum albumin, upper limit can be up to 48.00 ng/mL through dilution.

PERFORMANCE CHARACTERISTICS

- 1. Measuring Range 0.10 ~ 16.00 ng/mL
- 2. Limit of Detection ≤0.10 ng/mL
- 3. Within-Run Precision ≤10%
- 4. Between-Run Precision ≤15%

LIMITATIONS

- 1. The test results of this kit are only for clinical reference and cannot be used as the basis for confirming or excluding cases alone. In order to achieve the purpose of diagnosis, this test result should be used in combination with clinical examination, medical history and other examination results.
- 2. Do not use the test card if the foil pouch or the cartridge is damaged.

- 3. Do not open pouches until performing the test.
- 4. Patient samples may contain heterophilic antibodies (e.g. human anti-mouse antibodies (HAMA) and rheumatoid factors) that could react in immunoassays to give a falsely elevated or depressed result. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed.
- 5. Triglyceride and bilirubin in the sample may interfere with the test results, and the maximum allowable concentrations are 18 g/L and 0.1 g/L respectively.

EXPECTED VALUE

The expected normal value for Testosterone was determined by testing samples from apparently healthy individuals.

Reference range of Testosterone:

Group	n	95% Reference range
		(ng/mL)
Healthy men	260	1.75-7.81
Healthy women	280	0.10-0.75

Expected values may vary with age, sex, sample type, diet and geographical location. Each laboratory should verify the transferability of the expected values to its own population, and if necessary, determine its own expected values according to good laboratory practice.














REFERENCES

- 1. Payne AH, Hales DB. Overview of steroidogenic enzymes in the pathway from cholesterol to active steroid hormones. *Endocr Rev.* 2004 Dec;25(6):947-70. doi: 10.1210/er.2003-0030. PMID: 15583024.
- 2. RICHARD, Sharpe M. Intratesticular Factors Controlling Testicular Function[J]. *Biol Reprod.* 1984,30 (1) : 29-49
- 3. EN ISO 18113-1:2011 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
- 4. EN ISO 18113-2:2011 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on Testosterone Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European

Standard EN ISO 15223-1:2021.


Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use or consult electronic instructions for use		Batch code
	Temperature limit		In vitro diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community/European Union
	CE mark		Do not use if package is damaged and consult instructions for use
	Catalogue number		

Thank you for using Testosterone Fast Test Kit (Immunofluorescence Assay). Please read this instructions for use carefully before operating to ensure proper use. Please report any product problems or adverse events to the below manufacture or authorized representative in the European Community in time.

Document no.: WIF80-S-07

Effective date: 2023.08.29

 Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn
overseas@getein.com.cn
Website: www.getein.com

 CMC Medical Devices & Drugs S.L.
Add: C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain
Tel: +34951214054



tPSA Fast Test Kit (Immunofluorescence Assay)

IF2053 for Getein1600
IF1053 for Getein1100
IF3053 for Getein1180
IF4053 for Getein1200
IF5053 for Getein1160



User Manual

INTENDED USE

tPSA (total Prostate Specific Antigen) Fast Test Kit (Immunofluorescence Assay) is intended for *in vitro* quantitative determination of tPSA in human serum and plasma samples. It can be used as an aid in the diagnosis and management of patients with prostate cancer.

SUMMARY

Prostate-specific antigen (PSA) is a single-chain glycoprotein with molecular weight of 34 kilodaltons. As a serine protease with chymotrypsin-like activity, PSA belongs to the kallikrein family. PSA exists as a free or complex form with protease inhibitors such as α -1-antichymotrypsin (ACT) in blood. Total PSA represents the sum of both free and complex forms. Elevated PSA in serum or plasma is found in patients with prostate cancer, benign prostatic hypertrophy, or inflammatory tissues. PSA is uniquely associated with prostate tissues from normal, inflamed or cancerous states.

PSA has been found in normal, benign hyperplastic, malignant prostatic tissue, metastatic prostatic carcinoma and also in prostatic fluid as well as in seminal fluid. PSA is not found in any other tissues in men, and it is not produced by cancers originating in the lung, colon, rectum, stomach, pancreas or thyroid. PSA measurement is an essential tool in assessing the status of disease in patients with prostate cancer when serial samples are measured over time. The clinical value realized by monitoring tPSA concentration in patients with prostate cancer regardless of the treatment regimen is well known.

PRINCIPLE

The test uses an anti-human PSA monoclonal antibody conjugated with fluorescence latex coated on the junction of nitrocellulose membrane and sample pad and another anti-human PSA monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled antihuman PSA antibody binds with the PSA in sample and forms marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by antihuman PSA antibody. The

fluorescence intensity of test lines increases in proportion to the amount of tPSA in sample. Then insert test card into Getein1100/Getein1160/Getein1180 Immunofluorescence Quantitative Analyzer/automatically inserted by Getein1200/Getein1600 Immunofluorescence Quantitative Analyzer (hereafter referred to as Getein1100, Getein1160, Getein1180, Getein1200 and Getein1600), the concentration of tPSA in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100/Getein1160/Getein1180 contains:

- Package specifications: 25 tests/box, 10 tests/box
- 1) Getein tPSA test card in a sealed pouch with desiccant
 - 2) Disposable pipet
 - 3) User manual: 1 piece/box
 - 4) SD card: 1 piece/box

2. A kit for Getein1200/Getein1600 contains:

- Package specifications: 2×24 tests/box, 2×48 tests/box
Sealed cartridge with 24/48 Getein tPSA test cards
User manual: 1 piece/box
Materials required for Getein1200/Getein1600:

- 1) Sample diluent: 1 bottle/box
- 2) Box with pipette tips: 96 tips/box
- 3) Mixing plate: 1 piece/box

3. Sample diluent composition:

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

4. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad (the junction of sample pad and nitrocellulose membrane is coated with fluorescence latex-labelled anti-human PSA monoclonal antibody), nitrocellulose membrane (test line is coated with another fluorescence latex-labelled anti-human PSA monoclonal antibody and the control line C is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1180 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer
Getein1160 Immunofluorescence Quantitative Analyzer
Getein1200 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100/Getein1160/Getein1180 within 1 hour once the foil pouch is opened.

For test card of Getein1200/Getein1600: if the cartridge is

opened, it could be stable within 24 hours once exposed to air. If the test cards can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days.
Store the sample diluent at 0~30°C with a valid period of 24 months.

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Do not use the kit beyond the expiration date.
3. Do not use the test card if the foil pouch is damaged.
4. Do not open pouches until ready to perform the test.
5. Do not reuse the test card.
6. Do not reuse the pipet.
7. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
8. Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

1. This test can be used for **serum and plasma samples**. Heparin, EDTA and sodium citrate can be used as the anticoagulant for plasma. Samples should be free of hemolysis.
2. The test should be performed within 4 hours after whole blood collection.
3. If testing is delayed, serum and plasma samples may be stored up to 7 days at 2~8°C or stored at -20°C for 6 months before testing.
4. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze-thaw cycles.
5. Do not use heat-inactivated samples or hemolysis samples.
6. **SAMPLE VOLUME (for Getein1100/Getein1160/Getein1180): 100 μ l.**

TEST PROCEDURE

1. Collect specimens according to user manual.
2. Test card, sample and reagent should be brought to room temperature before testing.
For Getein1100:
3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
4. Enter testing interface of Getein1100.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver **100 μ l** of sample into the sample port on the test card.
8. **Reaction time: 15 minutes.** Insert the test card into

Getein1100/Getein1160/Getein1180/Getein1200/Getein1600. The result will be shown on the screen and printed automatically.

For Getein1160/Getein1180:

9. Confirm SD card lot No. in accordance with test kit lot No.. Perform "SD card" calibration when necessary.
10. Enter testing interface of Getein1160/Getein1180.
11. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
12. Put the test card on a clean table, horizontally placed.
13. Using sample transfer pipette, deliver **100 μ l** of sample into the sample port on the test card.
14. **Reaction time: 15 minutes.** Insert the test card into Getein1160/Getein1180 immediately after sample loading. The analyzer will count down the reaction time and automatically test the card after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1200/Getein1600:

15. Each cartridge for Getein1200/Getein1600 contains a specific RFID card which can calibrate automatically.
16. Place the sample diluent at the correct position in Getein1200/Getein1600.
17. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1200/Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein1100/Getein1160/Getein1180.
3. Make sure the test card insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1160/Getein1180/Getein1200/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1160/Getein1180/Getein1200/Getein1600.

Others:

Measuring range of the tPSA test kit is 0.50 ng/mL~100.00 ng/mL. Dilute the sample which concentration is higher than the upper limit with negative samples, and the dilution ratio should be less than 4 times.

EXPECTED VALUE

The expected normal value for tPSA was determined by testing samples from 1000 apparently healthy individuals. The reference range of tPSA is 4.00 ng/mL calculated by using normal distribution methods (95% confidence interval). It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.50~100.00 ng/mL
Lower Detection Limit	≤ 0.50 ng/mL
Within-Run Precision	≤ 10%
Between-Run Precision	≤ 15%
With-run Precision: Test tPSA with same batch for 10 times using tPSA control 1 (3.20~4.80 ng/mL) and tPSA control 2 (24.00~36.00 ng/mL) respectively, then calculate within-run precision which should not greater than 10%.	
Between-run Precision: Randomly select 3 consecutive batches of tPSA products, and take 10 strips for each batch to test the quality control (24.00~36.00 ng/mL), calculate between-run precision which should not greater than 15%.	

LIMITATIONS

- As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
- Samples containing interferent may influence the results. The table below listed the maximum allowance of these potential interferent.













Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25g/L	0.1 g/L

REFERENCES

- Mc Jimpsey EL. Molecular Form Differences Between Prostate-Specific Antigen (PSA) Standards Create Quantitative Discordances in PSA ELISA Measurements. Scientific Reports. 2016, 6: 22050.
- Jun Seok Kim, Je-Guk Ryu, Jin Woong Kim, et al. Prostate-Specific Antigen fluctuation: what does it mean in diagnosis of prostate cancer? Br J Cancer. Int Braz J Urol. 2015, 41(2): 258-264.
- Yasuhide Kitagawa, Mikio Namiki. Prostate-specific antigen-based population screening for prostate cancer: current status in Japan and future perspective in Asia. Asian J Androl. 2015, 17(3): 475-480.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on tPSA Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2016.

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limit		<i>In vitro</i> diagnostic medical device
	Contains sufficient for <n> tests		Do not use if package is damaged
	Catalogue number		Authorized representative in the European Community

Thank you for purchasing tPSA Fast Test Kit (Immunofluorescence Assay).
Please read this user manual carefully before operating to ensure proper use.

Version: WIF48-S-07



Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn
overseas@getein.com.cn
Website: www.getein.com



Lotus NL B.V.
Add: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
E-mail: peter@lotusnl.com
Tel: +31644168999



TSH Fast Test Kit

(Immunofluorescence Assay)

Getein1100: Cat.# IF1024

Getein1600: Cat.# IF2024

User Manual

INTENDED USE

TSH Fast Test Kit (Immunofluorescence Assay) is intended for in vitro quantitative determination of thyroid-stimulating hormone (TSH) in serum, plasma or whole blood. This test is used in the screening, clinical diagnosis, prognosis and therapeutic effect evaluation of thyroid diseases.

SUMMARY

Thyroid-stimulating hormone (TSH) is the main regulator of thyroid cell growth, thyroid hormone synthesis and secretion. TSH(MW 30 kDa) is synthesized and secreted by tsh cells of pituitary gland, it has negative feedback to the synthesis and secretion process. The fluctuation of TSH is faster and more significant than thyroid hormones when thyroid function was changed, it is a sensitive biomarker of hypothalamic-pituitary-thyroid function.

PRINCIPLE

The test uses an anti-human TSH monoclonal antibody conjugated with fluorescence latex and another anti-human TSH monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human TSH monoclonal antibody binds with the TSH in sample and forms marked antigen-antibody complex. The complex moves to the detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human TSH monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of TSH in sample. Insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Automatically inserted by Getein1600 Immunofluorescence Quantitative Analyzer (hereafter referred

to as Getein1100 and Getein1600), the concentration of TSH in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

1. A kit for Getein1100 contains:

Getein TSH test card in a sealed pouch with desiccant	25
Disposable pipet	25
User manual	1
SD card/RFID card	1
Whole blood buffer	1

2. A kit for Getein1600 contains:

Package specifications: 2×24 tests/kit, 2×48 tests/kit	
Sealed cartridge with 24/48 Getein TSH test cards	2
User manual	1
Sample diluent	1
Box with pipette tips	1
Mixing plate	1

3. Sample diluent / Whole blood buffer composition:

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

4. A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human TSH monoclonal antibody, the test line is coated with another anti-human TSH monoclonal antibody and the control line is coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

For test card of Getein1600: if the cartridge is opened, it could be stable within 24 hours once exposed to air. If the test cards

can't be used up at a time, please put the cartridge back to the foil pouch and reseal along the entire edge of zip-seal. The remaining test cards should be used up within 7 days. Store the sample diluent/whole blood buffer at 0~30°C with a valid period of 24 months. Store the sample diluent/whole blood buffer at 2~8°C for better results.

PRECAUTIONS

- 1.For in vitro diagnostic use only.
- 2.For professional use only.
- 3.Do not use the kit beyond the expiration date.
- 4.Do not use the test card if the foil pouch is damaged.
- 5.Do not open pouches until ready to perform the test.
- 6.Do not reuse the test card.
- 7.Do not reuse the pipet.
- 8.Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- 9.Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- 1.This test can be used for *serum, plasma and whole blood. Heparin, EDTA and sodium citrate* can be used as the anticoagulant for plasma and whole blood samples. Samples should be free of hemolysis.
2. Suggest using serum or plasma for better results.
3. Serum or plasma can be used directly. For whole blood sample, one drop of whole blood buffer must be added before testing.
4. If testing will be delayed, serum and plasma samples may be stored up to 7 days at 2~8°C or stored at -20°C for 6 months before testing (whole blood sample may be stored up to 3 days at 2~8°C).
5. Refrigerated or frozen sample should reach room temperature and be homogeneous before testing. Avoid multiple freeze- thaw cycles.
6. Do not use heat-inactivated samples.
7. SAMPLE VOLUME (for Getein1100): 100 µl.

TEST PROCEDURE

- 1.Collect specimens according to user manual.
- 2.Test card, sample and reagent should be brought to room temperature before testing.

For Getein1100:

3. Confirm SD card or RFID card lot No. in accordance with test kit lot No.. Perform "SD card or RFID card Calib" calibration when necessary.
4. Enter testing interface of Getein1100.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver 100 µl of sample (or 3~4 drops of sample when using disposable pipet) into the sample port on the test card (for whole blood sample, one drop of whole blood buffer must be added after loading 100 µl sample on the test card).

8. **Reaction time: 15 minutes.** Insert the test card into Getein1100 and start test after reaction time is elapsed. The result will be shown on the screen and printed automatically.

For Getein1600:

9. Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
10. Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

1. It is required to perform "SD card or RFID card Calib" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits for Getein1100.
3. Make sure the test card insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

Others: Measuring range of the test kit is 0.10 µIU/mL~50.0 µIU/mL, dilute the sample which concentration is higher than the upper limit, the dilution ratio should be less than 4 times.

EXPECTED VALUE

The expected normal value for TSH was determined by testing samples from serum of 391 apparently healthy individuals. The reference range of TSH is 0.27 µIU/mL~ 4.20 µIU/mL

calculated by using normal distribution methods (95% confidence interval). It is recommended that each laboratory should establish its expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	0.10 µIU/mL~ 50.00 µIU/mL
Lower Detection Limit	0.10 µIU/mL
Within-Run Precision	≤10%
Between-Run Precision	≤15%

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferences may influence the results. The table below listed the maximum allowance of these potential interferences.

Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25 g/L	0.1 g/L













REFERENCES

1. Spencer C A, LoPresti J S, Patel A, et al. Applications of a new chemiluminometric Thyrotropin assay to subnormal assessment. Clin Endocrinol Metab. 1990, 70(2):453-460.
2. Sakai H, Fukuda G, Suzuki N, et al. Falsely Elevated Thyroid-Stimulating Hormone (TSH) Level Due to Macro-TSH. Endocr J. 2009, 56(3):435-440.
3. Abalovich M, Amino N, Barbour LA, et al. Management of thyroid dysfunction during pregnancy and postpartum: an Endocrine Society Clinical Practice Guideline. J Clin Endocrinol Metab. 2007, 92(8):1-47.
4. Spencer C A, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyrotropin assays. Clin Chem. 1996, 42(1):164-173.
5. EN ISO 18113-1:2011 In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
6. EN ISO 18113-2:2011 In vitro diagnostic medical devices -

Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use (ISO 18113-2:2011).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on TSH Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2012.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		In vitro diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing TSH Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF26-S-02



Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508 Fax: +86-25-68568500

E-mail: tech@getein.com.cn overseas@getein.com.cn

Website: www.bio-GP.com.cn