

DICOM Conformance Statement

iQ-VIEW 4.1

A component of iQ-SYSTEM PACS v1

MD **CE** 0482

Table of Contents

1. Conformance statement overview.....	7
2. Introduction	12
2.1 Audience	12
2.2 Remarks	12
2.3 Terms and definitions	12
2.4 Abbreviations	15
2.5 References.....	16
3. Networking	17
3.1 Implementation model	17
3.1.1 Application data flow	17
3.1.2 Functional definition of AEs	18
3.1.2.1 Functional definition of "STORAGE-SCU"	18
3.1.2.2 Functional definition of "WORKLIST-SCU"	18
3.1.2.3 Functional definition of "PRINT-SCU"	18
3.1.2.4 Functional definition of "FIND-SCU"	18
3.1.2.5 Functional definition of "MOVE-SCU"	18
3.1.2.6 Functional definition of "STORAGE-SCP"	18
3.1.2.7 Functional definition of "ECHO-SCP"	19
3.1.3 Sequencing of real world activities	19
3.2 AE specifications	19
3.2.1 STORAGE-SCU.....	19
3.2.1.1 SOP classes.....	19
3.2.1.2 Association policies.....	23
3.2.1.2.1 General	23
3.2.1.2.2 Number of associations.....	24
3.2.1.2.3 Asynchronous nature	24
3.2.1.2.4 Implementation identifying information	24
3.2.1.3 Association initiation policy.....	24
3.2.1.3.1 Activity – Send storage request	24
3.2.1.3.1.1 Description and sequencing of activities.....	24
3.2.1.3.1.2 Proposed presentation contexts	24
3.2.1.3.1.2.1 Extended negotiation.....	25
3.2.1.3.1.3 SOP specific conformance for SOP class(es).....	26
3.2.1.3.1.3.1 SOP specific conformance to storage SOP classes.....	26
3.2.1.3.1.3.2 Presentation context acceptance criterion	26
3.2.1.3.1.3.3 Transfer syntax selection policies	26
3.2.1.3.1.3.4 Response status.....	26
3.2.1.4 Association acceptance policy	26
3.2.2 WORKLIST-SCU	27
3.2.2.1 SOP classes.....	27
3.2.2.2 Association policies.....	27
3.2.2.2.1 General.....	27
3.2.2.2.2 Number of associations	27
3.2.2.2.3 Asynchronous nature.....	27
3.2.2.2.4 Implementation identifying information	27
3.2.2.3 Association initiation policy.....	27

3.2.2.3.1	Send worklist request.....	28
3.2.2.3.1.1	Description and sequencing of activities	28
3.2.2.3.1.2	Proposed presentation contexts	28
3.2.2.3.1.2.1	Extended negotiation	28
3.2.2.3.1.3	SOP specific conformance for SOP class(es)	28
3.2.2.3.1.3.1	SOP specific conformance to workflow management SOP classes	28
3.2.2.3.1.3.2	Presentation context acceptance criterion	30
3.2.2.3.1.3.3	Transfer syntax selection policies.....	30
3.2.2.3.1.3.4	Response status.....	30
3.2.2.4	Association acceptance policy	31
3.2.3	PRINT-SCU.....	31
3.2.3.1	SOP classes.....	31
3.2.3.2	Association policies.....	31
3.2.3.2.1	General.....	31
3.2.3.2.2	Number of associations.....	31
3.2.3.2.3	Asynchronous nature.....	31
3.2.3.2.4	Implementation identifying information	32
3.2.3.3	Association initiation policy.....	32
3.2.3.3.1	Send print request.....	32
3.2.3.3.1.1	Description and sequencing of activities	32
3.2.3.3.1.2	Proposed presentation contexts	32
3.2.3.3.1.2.1	Extended Negotiation.....	32
3.2.3.3.1.3	SOP specific conformance for SOP class(es)	33
3.2.3.3.1.3.1	SOP specific conformance for the printer SOP class	33
3.2.3.3.1.3.2	Printer SOP class operations (N-GET)	33
3.2.3.3.1.3.3	SOP specific conformance for the film session SOP class.....	34
3.2.3.3.1.3.4	Film session SOP class operations (N-CREATE)	34
3.2.3.3.1.3.5	Film session SOP class operations (N-DELETE)	35
3.2.3.3.1.3.6	SOP specific conformance for the film box SOP class.....	35
3.2.3.3.1.3.6.1	Film box SOP class operations (N-CREATE)	35
3.2.3.3.1.3.7	Film box SOP class operations (N-ACTION)	37
3.2.3.3.1.3.8	SOP specific conformance for the image box SOP class	38
3.2.3.3.1.3.9	Image box SOP class operations (N-SET).....	38
3.2.3.4	Association acceptance policy	39
3.2.4	FIND-SCU.....	39
3.2.4.1	SOP classes.....	39
3.2.4.2	Association policies.....	39
3.2.4.2.1	General.....	39
3.2.4.2.2	Number of associations.....	39
3.2.4.2.3	Asynchronous nature	40
3.2.4.2.4	Implementation identifying information	40
3.2.4.3	Association initiation policy.....	40
3.2.4.3.1	Send query request.....	40
3.2.4.3.1.1	Description and sequencing of activities	40
3.2.4.3.1.2	Proposed presentation contexts	40
3.2.4.3.1.2.1	Extended negotiation	41
3.2.4.3.1.3	SOP specific conformance for SOP class(es)	41
3.2.4.3.1.3.1	SOP specific conformance to query SOP classes.....	41
3.2.4.3.1.3.2	Presentation context acceptance criterion	42
3.2.4.3.1.3.3	Transfer syntax selection policies	42
3.2.4.4	Association acceptance policy	43

3.2.5	MOVE-SCU.....	43
3.2.5.1	SOP classes.....	43
3.2.5.2	Association policies.....	44
3.2.5.2.1	General.....	44
3.2.5.2.2	Number of associations.....	44
3.2.5.2.3	Asynchronous nature.....	44
3.2.5.2.4	Implementation identifying information.....	44
3.2.5.3	Association initiation policy.....	44
3.2.5.3.1	Send retrieve request.....	45
3.2.5.3.1.1	Description and sequencing of activities.....	45
3.2.5.3.1.2	Proposed presentation contexts.....	45
3.2.5.3.1.2.1	Extended negotiation.....	45
3.2.5.3.1.3	SOP specific conformance for SOP class(es).....	45
3.2.5.3.1.3.1	SOP specific conformance to retrieve SOP classes.....	45
3.2.5.3.1.3.2	Presentation context acceptance criterion.....	46
3.2.5.3.1.3.3	Transfer syntax selection policies.....	46
3.2.5.3.1.3.4	Response status.....	46
3.2.5.3.1.3.5	Sub-operation dependent behavior.....	47
3.2.5.4	Association acceptance policy.....	47
3.2.6	STORAGE-SCP.....	47
3.2.6.1	SOP classes.....	47
3.2.6.2	Association policies.....	52
3.2.6.2.1	General.....	52
3.2.6.2.2	Number of associations.....	52
3.2.6.2.3	Asynchronous nature.....	53
3.2.6.2.4	Implementation identifying information.....	53
3.2.6.3	Association initiation policy.....	53
3.2.6.4	Association acceptance policy.....	53
3.2.6.4.1	Receive storage request.....	53
3.2.6.4.1.1	Description and sequencing of activities.....	53
3.2.6.4.1.2	Accepted presentation contexts.....	53
3.2.6.4.1.2.1	Extended negotiation.....	54
3.2.6.4.1.3	SOP specific conformance for SOP class(es).....	54
3.2.6.4.1.3.1	SOP specific conformance to storage SOP classes.....	54
3.2.6.4.1.3.2	Presentation context acceptance criterion.....	54
3.2.6.4.1.3.3	Transfer syntax selection policies.....	54
3.2.6.4.1.3.4	Response status.....	55
3.2.7	ECHO-SCP.....	55
3.2.7.1	SOP classes.....	55
3.2.7.2	Association policies.....	55
3.2.7.2.1	General.....	55
3.2.7.2.2	Number of associations.....	56
3.2.7.2.3	Asynchronous nature.....	56
3.2.7.2.4	Implementation identifying information.....	56
3.2.7.3	Association initiation policy.....	56
3.2.7.4	Association acceptance policy.....	56
3.2.7.4.1	Receive echo request.....	56
3.2.7.4.1.1	Description and sequencing of activities.....	56
3.2.7.4.1.2	Accepted presentation contexts.....	56
3.2.7.4.1.2.1	Extended negotiation.....	57
3.2.7.4.1.3	SOP specific conformance for SOP class(es).....	57

3.2.7.4.1.3.1	SOP specific conformance to verification SOP classes	57
3.2.7.4.1.3.2	Presentation context acceptance criterion	57
3.2.7.4.1.3.3	Transfer syntax selection policies.....	57
3.3	Network interfaces	57
3.3.1	Physical network interface	57
3.3.2	Additional protocols	58
3.3.3	IPv4 and IPv6 support.....	58
3.4	Configuration	58
3.4.1	AE title / presentation address mapping	58
3.4.2	Parameters	58
4.	Media interchange.....	59
4.1	Implementation model	59
4.1.1	Application data flow diagram	59
4.1.2	Functional definition of AEs	59
4.1.2.1	Functional definition of "MEDIA-FSR"	59
4.1.2.2	Functional definition of "MEDIA-FSC"	59
4.1.3	Sequencing of real world activities	60
4.1.3.1	Activity "import from medium"	60
4.1.3.2	Activity "export to medium"	60
4.1.4	File meta information for implementation class and version	60
4.2	AE specifications	60
4.2.1	MEDIA-FSR application entity specification	60
4.2.1.1	File meta information for the application entity	60
4.2.1.2	Real-world activities.....	61
4.2.1.2.1	Activity – Import from medium.....	61
4.2.2	MEDIA-FSC application entity specification	61
4.2.2.1	File meta information for the application entity	61
4.2.2.2	Real-world activities	61
4.2.2.2.1	Activity – Export to medium.....	61
4.3	Augmented and private application profiles	61
4.4	Media configuration.....	62
5.	Transformation of DICOM to CDA	63
6.	Support of character sets	64
6.1	Overview	64
6.2	Character sets.....	64
6.3	Character set configuration.....	65
7.	Security	66
7.1	Security profiles.....	66
7.2	Association level security	66
7.3	Application level security	66
8.	Annexes.....	67
8.1	IOD contents	67
8.1.1	Created SOP instance(s)	67
8.1.1.1	Secondary capture image storage	67
8.1.1.2	Grayscale softcopy presentation state storage	68
8.1.1.3	Basic text SR storage	68
8.1.1.4	Encapsulated PDF storage.....	69

8.1.1.5	Common modules	69
8.1.1.6	Secondary capture image modules	71
8.1.1.7	Grayscale softcopy presentation state modules.....	73
8.1.1.8	Basic text structured report modules	77
8.1.1.9	Encapsulated PDF modules	80
8.1.2	Usage of attributes from received IODs.....	81
8.1.3	Attribute mapping.....	81
8.1.4	Coerced / modified fields	81
8.2	Data dictionary of private attributes.....	81
8.3	Coded terminology and templates	82
8.4	Grayscale image consistency	82
8.5	Standard extended / specialized / private SOP classes.....	82
8.6	Private transfer syntaxes.....	82

1. Conformance statement overview

This software and version implements the necessary DICOM services to transfer DICOM objects via network transmission, to query other DICOM systems for workflow and data information and send print requests to DICOM Print aware systems (printer).

Table 1-1 provides an overview of the network services supported by this software and version:

Table 1-1 Network services

Network services		
SOP Classes	User of service (SCU)	Provider of service (SCP)
Transfer		
(RETIRED) Stored Print Storage	(Yes)*	Yes
(RETIRED) Hardcopy Grayscale Image Storage	Yes	Yes
(RETIRED) Hardcopy Color Image Storage	Yes	Yes
Computed Radiography Image Storage	Yes	Yes
Digital X-Ray Image Storage For Presentation	Yes	Yes
Digital X-Ray Image Storage For Processing	Yes	Yes
Digital Mammography X-Ray Image Storage For Presentation	Yes	Yes
Digital Mammography X-Ray Image Storage For Processing	Yes	Yes
Digital Intra Oral X-Ray Image Storage For Presentation	Yes	Yes
Digital Intra Oral X-Ray Image Storage For Processing	Yes	Yes
CT Image Storage	Yes	Yes
Enhanced CT Image Storage	Yes	Yes
Legacy Converted Enhanced CT Image Storage	(Yes)*	Yes
(RETIRED) Ultrasound Multiframe Image Storage	Yes	Yes
Ultrasound Multiframe Image Storage	Yes	Yes
MR Image Storage	Yes	Yes
Enhanced MR Image Storage	Yes	Yes
MR Spectroscopy Storage	(Yes)*	Yes
Enhanced MR Color Image Storage	Yes	Yes
Legacy Converted Enhanced MR Image Storage	(Yes)*	Yes
(RETIRED) Nuclear Medicine Image Storage	Yes	Yes
(RETIRED) Ultrasound Image Storage	Yes	Yes
Ultrasound Image Storage	Yes	Yes
Enhanced US Volume Storage	(Yes)*	Yes
Secondary Capture Image Storage	Yes	Yes
Multiframe Single Bit Secondary Capture Image Storage	Yes	Yes
Multiframe Grayscale Byte Secondary Capture Image Storage	Yes	Yes

Network services

Multiframe Grayscale Word Secondary Capture Image Storage	Yes	Yes
Multiframe True Color Secondary Capture Image Storage	Yes	Yes
(RETIRED) Standalone Overlay Storage	(Yes)*	Yes
(RETIRED) Standalone Curve Storage	(Yes)*	Yes
Twelve Lead ECG Waveform Storage	(Yes)*	Yes
General ECG Waveform Storage	(Yes)*	Yes
Ambulatory ECG Waveform Storage	(Yes)*	Yes
Hemodynamic Waveform Storage	(Yes)*	Yes
Cardiac Electrophysiology Waveform Storage	(Yes)*	Yes
Basic Voice Audio Waveform Storage	(Yes)*	Yes
General Audio Waveform Storage	(Yes)*	Yes
Arterial Pulse Waveform Storage	(Yes)*	Yes
Respiratory Waveform Storage	(Yes)*	Yes
(RETIRED) Standalone Modality LUT Storage	(Yes)*	Yes
(RETIRED) Standalone VOI LUT Storage	(Yes)*	Yes
Grayscale Softcopy Presentation State Storage	Yes**	Yes
Color Softcopy Presentation State Storage	(Yes)*	Yes
Pseudo Color Softcopy Presentation State Storage	(Yes)*	Yes
Blending Softcopy Presentation State Storage	(Yes)*	Yes
XA XRF Grayscale Softcopy Presentation State Storage	(Yes)*	Yes
X-Ray Angiographic Image Storage	Yes	Yes
Enhanced XA Image Storage	Yes	Yes
X-Ray Radiofluoroscopic Image Storage	Yes	Yes
Enhanced XRF Image Storage	(Yes)*	Yes
(RETIRED) X-Ray Angiographic BiPlane Image Storage	Yes	Yes
X-Ray 3D Angiographic Image Storage	(Yes)*	Yes
X-Ray 3D Craniofacial Image Storage	(Yes)*	Yes
Breast Tomosynthesis Image Storage	Yes	Yes
Intravascular Optical Coherence Tomography Image Storage For Presentation	(Yes)*	Yes
Intravascular Optical Coherence Tomography Image Storage For Processing	(Yes)*	Yes
Nuclear Medicine Image Storage	Yes	Yes
Raw Data Storage	(Yes)*	Yes
Spatial Registration Storage	(Yes)*	Yes
Spatial Fiducials Storage	(Yes)*	Yes
Deformable Spatial Registration Storage	(Yes)*	Yes
Segmentation Storage	(Yes)*	Yes
Surface Segmentation Storage	(Yes)*	Yes
Real World Value Mapping Storage	(Yes)*	Yes
Surface Scan Mesh Storage	(Yes)*	Yes
Surface Scan Point Cloud Storage	(Yes)*	Yes

Network services		
(RETIRED) VL Image Storage	Yes	Yes
VL Endoscopic Image Storage	Yes	Yes
Video Endoscopic Image Storage	Yes	Yes
VL Microscopic Image Storage	Yes	Yes
Video Microscopic Image Storage	Yes	Yes
VL Slide Coordinates Microscopic Image Storage	(Yes)*	Yes
VL Photographic Image Storage	Yes	Yes
Video Photographic Image Storage	Yes	Yes
Ophthalmic Photography 8Bit Image Storage	Yes	Yes
Ophthalmic Photography 16Bit Image Storage	Yes	Yes
Stereometric Relationship Storage	(Yes)*	Yes
Ophthalmic Tomography Image Storage	(Yes)*	Yes
VL Whole Slide Microscopy Image Storage	(Yes)*	Yes
(RETIRED) VL Multi Frame Image Storage	Yes	Yes
Lensometry Measurements Storage	(Yes)*	Yes
Autorefracton Measurements Storage	(Yes)*	Yes
Keratometry Measurements Storage	(Yes)*	Yes
Subjective Refraction Measurements Storage	(Yes)*	Yes
Visual Acuity Measurements Storage	(Yes)*	Yes
Spectacle Prescription Report Storage	(Yes)*	Yes
Ophthalmic Axial Measurements Storage	(Yes)*	Yes
Intraocular Lens Calculations Storage	(Yes)*	Yes
Macular Grid Thickness And Volume Report Storage	(Yes)*	Yes
Ophthalmic Visual Field Static Perimetry Measurements Storage	(Yes)*	Yes
Ophthalmic Thickness Map Storage	(Yes)*	Yes
Corneal Topography Map Storage	(Yes)*	Yes
Basic Text SR Storage	Yes	Yes
Enhanced SR Storage	Yes	Yes
Comprehensive SR Storage	Yes	Yes
Comprehensive 3D SR Storage	(Yes)*	Yes
Procedure Log Storage	(Yes)*	Yes
Mammography CAD SR Storage	Yes	Yes
Key Object Selection Document Storage	Yes	Yes
Chest CAD SR Storage	(Yes)*	Yes
X-Ray Radiation Dose SR Storage	Yes	Yes
Colon CAD SR Storage	(Yes)*	Yes
Implantation Plan SR Document Storage	(Yes)*	Yes
Encapsulated PDF Storage	Yes	Yes
Encapsulated CDA Storage	(Yes)*	Yes
Positron Emission Tomography Image Storage	Yes	Yes
Legacy Converted Enhanced PET Image Storage	(Yes)*	Yes

Network services		
(RETIRED) Standalone PET Curve Storage	(Yes)*	Yes
Enhanced PET Image Storage	(Yes)*	Yes
Basic Structured Display Storage	(Yes)*	Yes
RT Image Storage	(Yes)*	Yes
RT Dose Storage	(Yes)*	Yes
RT Structure Set Storage	(Yes)*	Yes
RT Beams Treatment Record Storage	(Yes)*	Yes
RT Plan Storage	(Yes)*	Yes
RT Brachy Treatment Record Storage	(Yes)*	Yes
RT Treatment Summary Record Storage	(Yes)*	Yes
RT Ion Plan Storage	(Yes)*	Yes
RT Ion Beams Treatment Record Storage	(Yes)*	Yes
Generic Implant Template Storage	(Yes)*	Yes
Implant Assembly Template Storage	(Yes)*	Yes
Implant Template Group Storage	(Yes)*	Yes
RT Beams Delivery Instruction Storage	(Yes)*	Yes
Color Palette Storage	(Yes)*	Yes
(DRAFT) SR Text Storage	(Yes)*	Yes
(DRAFT) SR Audio Storage	(Yes)*	Yes
(DRAFT) SR Detail Storage	(Yes)*	Yes
(DRAFT) SR Comprehensive Storage	(Yes)*	Yes
(DRAFT) Waveform Storage	(Yes)*	Yes
(DRAFT) RT Beams Delivery Instruction Storage	(Yes)*	Yes
Query/Retrieve		
Study Root Query/Retrieve Information Model - FIND	Yes	No
Study Root Query/Retrieve Information Model - MOVE	Yes	No
Workflow Management		
Modality Worklist Information Model FIND	Yes	No
Print Management		
Basic Grayscale Print Management Meta	Yes	No

* These SOP classes are not supported by the viewer application, but network transfer (C-STORE) will be possible.

** The viewer application only supports Grayscale Softcopy Presentation State Objects created by iQ-VIEW itself. Other GSPS objects will not be stored or displayed.

Table 1-2 provides an overview of the Media Storage Application Profiles supported by iQ-VIEW:

Table 1-2 Media services

Media services		
Media Storage Application Profile	Write Files (FSC)	Read Files (FSR)

Media services		
Compact Disk - Recordable		
General Purpose CD-R Interchange	Yes	Yes
DVD		
General Purpose Interchange on DVD-RAM Media	Yes	Yes
BD		
General Purpose Interchange on BD Media	Yes	Yes

2. Introduction

2.1 Audience

This document is written for the people that need to understand how this software and version will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

2.2 Remarks

The scope of this DICOM Conformance Statement is to facilitate integration between this software and version and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard [DICOM]. DICOM by itself does not guarantee interoperability. The Conformance Statement does, however, facilitate a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

2.3 Terms and definitions

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Term	Definition
Abstract Syntax	The information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.
Application Entity (AE)	An end point of a DICOM information exchange, including the DICOM network or media interface software; i.e., the software that sends or receives DICOM information objects

Term	Definition
	or messages. A single device may have multiple Application Entities.
Application Entity Title (AET)	The externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network.
Application Context	The specification of the type of communication used between Application Entities. Example: DICOM network protocol.
Association	A network communication channel set up between Application Entities.
Attribute	A unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).
Information Object Definition (IOD)	The specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.
Joint Photographic Experts Group (JPEG)	A set of standardized image compression techniques, available for use by DICOM applications.
Media Application Profile	The specification of DICOM information objects and encoding exchanged on removable media(e.g., CDs)
Module	A set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.
Negotiation	First phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded.
Presentation Context	The set of DICOM network services used over an Association, as negotiated between Application Entities; includes Abstract Syntaxes and Transfer Syntaxes.
Protocol Data Unit (PDU)	A packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.
Security Profile	A set of mechanisms, such as encryption, user authentication, or digital signatures, used by an Application Entity to ensure confidentiality, integrity, and/or availability of exchanged DICOM data

Term	Definition
Service Class Provider (SCP)	<p>Role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User).</p> <p>Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).</p>
Service Class User (SCU)	<p>Role of an Application Entity that uses a DICOM network service; typically, a client.</p> <p>Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)</p>
Service/Object Pair Class (SOP Class)	<p>The specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification.</p> <p>Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.</p>
Service/Object Pair Instance (SOP Instance)	<p>An information object; a specific occurrence of information exchanged in a SOP Class.</p> <p>Examples: a specific x-ray image.</p>
Tag	<p>A 32-bit identifier for a data element, represented as a pair of four digit hexadecimal numbers, the "group" and the "element". If the "group" number is odd, the tag is for a private (manufacturer-specific) data element.</p> <p>Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210)[private data element]</p>
Transfer Syntax	<p>The encoding used for exchange of DICOM information objects and messages.</p> <p>Examples: JPEG compressed (images), little endian explicit value representation.</p>
Unique Identifier (UID)	<p>A globally unique "dotted decimal" string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier.</p> <p>Examples: Study Instance UID, SOP Class UID, SOP Instance UID.</p>
Value Representation (VR)	<p>The format type of an individual DICOM data element, such as text, an integer, a person's name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.</p>

2.4 Abbreviations

Abbreviations	Meaning
AE	Application Entity
AET	Application Entity Title
CAD	Computer Aided Detection
CDA	Clinical Document Architecture
CD-R	Compact Disk Recordable
CSE	Customer Service Engineer
CR	Computed Radiography
CT	Computed Tomography
DHCP	Dynamic Host Configuration Protocol
DICOM	Digital Imaging and Communications in Medicine
DIT	Directory Information Tree (LDAP)
DN	Distinguished Name (LDAP)
DNS	Domain Name System
DX	Digital X-ray
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
GSDF	Grayscale Standard Display Function
GSPS	Grayscale Softcopy Presentation State
HIS	Hospital Information System
HL7	Health Level 7 Standard
IHE	Integrating the Healthcare Enterprise
IOD	Information Object Definition
IPv4	Internet Protocol version 4
IPv6	Internet Protocol version 6
ISO	International Organization for Standards
IO	Intra-oral X-ray
JPEG	Joint Photographic Experts Group
LDAP	Lightweight Directory Access Protocol
LDIF	LDAP Data Interchange Format
LUT	Look-up Table
MAR	Medication Administration Record
MPEG	Moving Picture Experts Group
MG	Mammography (X-ray)
MPPS	Modality Performed Procedure Step
MR	Magnetic Resonance Imaging
MSPS	Modality Scheduled Procedure Step
MTU	Maximum Transmission Unit (IP)
MWL	Modality Worklist
NM	Nuclear Medicine

Abbreviations	Meaning
NTP	Network Time Protocol
O	Optional (Key Attribute)
OP	Ophthalmic Photography
OSI	Open Systems Interconnection
PACS	Picture Archiving and Communication System
PET	Positron Emission Tomography
PDU	Protocol Data Unit
R	Required (Key Attribute)
RDN	Relative Distinguished Name (LDAP)
RF	Radiofluoroscopy
RIS	Radiology Information System.
RT	Radiotherapy
SC	Secondary Capture
SCP	Service Class Provider
SCU	Service Class User
SOP	Service-Object Pair
SPS	Scheduled Procedure Step
SR	Structured Reporting
TCP/IP	Transmission Control Protocol/Internet Protocol
U	Unique (Key Attribute)
UL	Upper Layer
US	Ultrasound
VL	Visible Light
VR	Value Representation
XA	X-ray Angiography


2.5 References

- [DICOM]: Digital Imaging and Communications in Medicine (DICOM) Standard, <http://dicom.nema.org>
- [IHE]: IHE Radiology Technical Framework, <http://www.ihe.net>
- [HL7]: Health Level Seven (HL7) Standard, <http://www.hl7.org>

3. Networking

3.1 Implementation model

3.1.1 Application data flow



The application provides a user interface, internal database and network listener that spawns additional threads as necessary to handle incoming connections.

Conceptually the network services may be modeled as the following separate AEs, though in fact all the AEs share a single (configurable) AE Title:

STORAGE-SCU, which sends outbound instances
WORKLIST-SCU, which queries remote AEs for worklist items
PRINT-SCU, which sends print requests to remote AEs
FIND-SCU, which queries remote AEs for lists of studies, series and instances
MOVE-SCU, which retrieves Hanging Protocol Instances
STORAGE-SCP, which receives incoming instances
ECHO-SCP, which responds to verification requests

3.1.2 Functional definition of AEs

3.1.2.1 Functional definition of "STORAGE-SCU"

STORAGE-SCU is activated through the user interface when a user selects studies, series or instances from the local database and requests that they be sent to a remote AE (selected from a pre-configured list).

3.1.2.2 Functional definition of "WORKLIST-SCU"

WORKLIST-SCU is activated through the user interface when a user requests the patient and study information being updated from a list of worklist items queried from a remote AE (from a pre-configured list).

3.1.2.3 Functional definition of "PRINT-SCU"

PRINT-SCU is activated through the user interface when a user requests a selected study, series or set of images to be printed on a remote AE (from a pre-configured list).

3.1.2.4 Functional definition of "FIND-SCU"

FIND-SCU is activated through the user interface when a user selects a remote AE to query (from a pre-configured list), then initiates a query. Queries are performed on the study level first and continue recursively through the series and instance levels when requested by the user.

3.1.2.5 Functional definition of "MOVE-SCU"

MOVE-SCU is activated through the user interface when a user selects a study, series or instance for retrieval. A connection to the remote AE is established to initiate and monitor the retrieval and the STORAGE-SCP AE receives the retrieved instances.

3.1.2.6 Functional definition of "STORAGE-SCP"

STORAGE-SCP waits in the background for connections, will accept associations with Presentation Contexts for SOP Classes of the Storage Service Class, and will store the

received instances to the local database where they may subsequently be listed and viewed through the user interface.

3.1.2.7 Functional definition of "ECHO-SCP"

ECHO-SCP waits in the background for connections, will accept associations with Presentation Contexts for SOP Class of the Verification Service Class, and will respond successfully to echo requests.

3.1.3 Sequencing of real world activities

All SCP activities are performed asynchronously in the background and not dependent on any sequencing.

All SCU activities are sequentially initiated in the user interface, and another activity may not be initiated until the prior activity has completed.

3.2 AE specifications

3.2.1 STORAGE-SCU

3.2.1.1 SOP classes

STORAGE-SCU provide Standard Conformance to the following SOP Classes:

Table 3-1 SOP Classes Supported By STORAGE-SCU

SOP Class Name	SOP Class UID	SCU	SCP
(RETIRED) Stored Print Storage	1.2.840.10008.5.1.1.27	Yes	No
(RETIRED) Hardcopy Grayscale Image Storage	1.2.840.10008.5.1.1.29	Yes	No
(RETIRED) Hardcopy Color Image Storage	1.2.840.10008.5.1.1.30	Yes	No
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	Yes	No
Digital X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.1	Yes	No
Digital X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	Yes	No
Digital Mammography X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.2	Yes	No
Digital Mammography X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	Yes	No
Digital Intra Oral X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.3	Yes	No

SOP Class Name	SOP Class UID	SCU	SCP
Digital Intra Oral X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.3.1	Yes	No
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	Yes	No
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1	Yes	No
Legacy Converted Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.2	Yes	No
(RETIRED) Ultrasound Multiframe Image Storage	1.2.840.10008.5.1.4.1.1.3	Yes	No
Ultrasound Multiframe Image Storage	1.2.840.10008.5.1.4.1.1.3.1	Yes	No
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	Yes	No
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1	Yes	No
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2	Yes	No
Enhanced MR Color Image Storage	1.2.840.10008.5.1.4.1.1.4.3	Yes	No
Legacy Converted Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.4	Yes	No
(RETIRED) Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.5	Yes	No
(RETIRED) Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6	Yes	No
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	Yes	No
Enhanced US Volume Storage	1.2.840.10008.5.1.4.1.1.6.2	Yes	No
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Yes	No
Multiframe Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1	Yes	No
Multiframe Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	Yes	No
Multiframe Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3	Yes	No
Multiframe True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	Yes	No
(RETIRED) Standalone Overlay Storage	1.2.840.10008.5.1.4.1.1.8	Yes	No
(RETIRED) Standalone Curve Storage	1.2.840.10008.5.1.4.1.1.9	Yes	No
Twelve Lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	Yes	No
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	Yes	No
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	Yes	No
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	Yes	No
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	Yes	No
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	Yes	No
General Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.2	Yes	No
Arterial Pulse Waveform Storage	1.2.840.10008.5.1.4.1.1.9.5.1	Yes	No
Respiratory Waveform Storage	1.2.840.10008.5.1.4.1.1.9.6.1	Yes	No

SOP Class Name	SOP Class UID	SCU	SCP
(RETIRED) Standalone Modality LUT Storage	1.2.840.10008.5.1.4.1.1.10	Yes	No
(RETIRED) Standalone VOI LUT Storage	1.2.840.10008.5.1.4.1.1.11	Yes	No
Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.1	Yes	No
Color Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.2	Yes	No
Pseudo Color Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.3	Yes	No
Blending Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.4	Yes	No
XA XRF Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.5	Yes	No
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	Yes	No
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1	Yes	No
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	Yes	No
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1	Yes	No
(RETIRED) X-Ray Angiographic BiPlane Image Storage	1.2.840.10008.5.1.4.1.1.12.3	Yes	No
X-Ray 3D Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.13.1.1	Yes	No
X-Ray 3D Craniofacial Image Storage	1.2.840.10008.5.1.4.1.1.13.1.2	Yes	No
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3	Yes	No
Intravascular Optical Coherence Tomography Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.14.1	Yes	No
Intravascular Optical Coherence Tomography Image Storage For Processing	1.2.840.10008.5.1.4.1.1.14.2	Yes	No
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	Yes	No
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66	Yes	No
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1	Yes	No
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2	Yes	No
Deformable Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.3	Yes	No
Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.4	Yes	No
Surface Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.5	Yes	No
Real World Value Mapping Storage	1.2.840.10008.5.1.4.1.1.67	Yes	No
Surface Scan Mesh Storage	1.2.840.10008.5.1.4.1.1.68.1	Yes	No
Surface Scan Point Cloud Storage	1.2.840.10008.5.1.4.1.1.68.2	Yes	No
(RETIRED) VL Image Storage	1.2.840.10008.5.1.4.1.1.77.1	Yes	No
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	Yes	No

SOP Class Name	SOP Class UID	SCU	SCP
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1	Yes	No
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	Yes	No
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1	Yes	No
VL Slide Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	Yes	No
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	Yes	No
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1	Yes	No
Ophthalmic Photography 8Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	Yes	No
Ophthalmic Photography 16Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2	Yes	No
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3	Yes	No
Ophthalmic Tomography Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.4	Yes	No
VL Whole Slide Microscopy Image Storage	1.2.840.10008.5.1.4.1.1.77.1.6	Yes	No
(RETIRED) VL Multi Frame Image Storage	1.2.840.10008.5.1.4.1.1.77.2	Yes	No
Lensometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.1	Yes	No
Autorefractometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.2	Yes	No
Keratometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.3	Yes	No
Subjective Refraction Measurements Storage	1.2.840.10008.5.1.4.1.1.78.4	Yes	No
Visual Acuity Measurements Storage	1.2.840.10008.5.1.4.1.1.78.5	Yes	No
Spectacle Prescription Report Storage	1.2.840.10008.5.1.4.1.1.78.6	Yes	No
Ophthalmic Axial Measurements Storage	1.2.840.10008.5.1.4.1.1.78.7	Yes	No
Intraocular Lens Calculations Storage	1.2.840.10008.5.1.4.1.1.78.8	Yes	No
Macular Grid Thickness And Volume Report Storage	1.2.840.10008.5.1.4.1.1.79.1	Yes	No
Ophthalmic Visual Field Static Perimetry Measurements Storage	1.2.840.10008.5.1.4.1.1.80.1	Yes	No
Ophthalmic Thickness Map Storage	1.2.840.10008.5.1.4.1.1.81.1	Yes	No
Corneal Topography Map Storage	1.2.840.10008.5.1.4.1.1.82.1	Yes	No
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	Yes	No
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	Yes	No
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	Yes	No
Comprehensive 3D SR Storage	1.2.840.10008.5.1.4.1.1.88.34	Yes	No
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40	Yes	No
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	Yes	No
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	Yes	No
Chest CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.65	Yes	No

SOP Class Name	SOP Class UID	SCU	SCP
X-Ray Radiation Dose SR Storage	1.2.840.10008.5.1.4.1.1.88.67	Yes	No
Colon CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.69	Yes	No
Implantation Plan SR Document Storage	1.2.840.10008.5.1.4.1.1.88.70	Yes	No
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Yes	No
Encapsulated CDA Storage	1.2.840.10008.5.1.4.1.1.104.2	Yes	No
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	Yes	No
Legacy Converted Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.128.1	Yes	No
(RETIRED) Standalone PET Curve Storage	1.2.840.10008.5.1.4.1.1.129	Yes	No
Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.130	Yes	No
Basic Structured Display Storage	1.2.840.10008.5.1.4.1.1.131	Yes	No
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	Yes	No
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	Yes	No
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	Yes	No
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4	Yes	No
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	Yes	No
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6	Yes	No
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7	Yes	No
RT Ion Plan Storage	1.2.840.10008.5.1.4.1.1.481.8	Yes	No
RT Ion Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.9	Yes	No
Generic Implant Template Storage	1.2.840.10008.5.1.4.43.1	Yes	No
Implant Assembly Template Storage	1.2.840.10008.5.1.4.44.1	Yes	No
Implant Template Group Storage	1.2.840.10008.5.1.4.45.1	Yes	No
RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.7	Yes	No
Color Palette Storage	1.2.840.10008.5.1.4.39.1	Yes	No
(DRAFT) SR Text Storage	1.2.840.10008.5.1.4.1.1.88.1	Yes	No
(DRAFT) SR Audio Storage	1.2.840.10008.5.1.4.1.1.88.2	Yes	No
(DRAFT) SR Detail Storage	1.2.840.10008.5.1.4.1.1.88.3	Yes	No
(DRAFT) SR Comprehensive Storage	1.2.840.10008.5.1.4.1.1.88.4	Yes	No
(DRAFT) Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1	Yes	No
(DRAFT) RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.1	Yes	No

3.2.1.2 Association policies

3.2.1.2.1 General

STORAGE-SCU initiates but never accepts associations.

Table 3-2 Maximum PDU Size Received as SCP for STORAGE-SCU

Maximum PDU size received as SCP	
Maximum PDU size received	Unlimited

3.2.1.2.2 Number of associations

Table 3-3 Number of Associations as SCP for STORAGE-SCU

Number of associations as SCP	
Maximum number of simultaneous associations	1

3.2.1.2.3 Asynchronous nature

STORAGE-SCU will only allow a single outstanding operation on an Association. Therefore, STORAGE-SCU will not perform asynchronous operations window negotiation.

3.2.1.2.4 Implementation identifying information

Table 3-4 DICOM Implementation Class and Version for STORAGE-SCU

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

3.2.1.3 Association initiation policy

STORAGE-SCU attempts to initiate a new association for each series it attempts to transfer.

3.2.1.3.1 Activity – Send storage request

3.2.1.3.1.1 Description and sequencing of activities

For each series selected from the user interface to be transferred, a single attempt will be made to transmit it to the selected remote AE. If the send fails, for whatever reason, no retry will be performed, and an attempt will be made to send the next series.

3.2.1.3.1.2 Proposed presentation contexts

Table 3-5 Proposed Presentation Contexts for STORAGE-SCU

Presentation context table

Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-1	See Table 4-1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		JPEG Baseline - Process 1 (default 8-bit lossy JPEG compression)	1.2.840.10008.1.2.4.50		
		JPEG Extended - Process 2 & 4 (default 12-bit lossy JPEG compression)	1.2.840.10008.1.2.4.51		
		JPEG Lossless, Non-Hierarchical - Process 14	1.2.840.10008.1.2.4.57		
		JPEG Lossless, Non-Hierarchical, First-Order Prediction - Process 14 [Selection Value 1] (default lossless JPEG compression)	1.2.840.10008.1.2.4.70		
		JPEG-LS Lossless	1.2.840.10008.1.2.4.80		
		JPEG-LS Lossy (Near-Lossless)	1.2.840.10008.1.2.4.81		
		JPEG 2000 (Lossless Only)	1.2.840.10008.1.2.4.90		
		JPEG 2000 (Lossless or Lossy)	1.2.840.10008.1.2.4.91		
		RLE Lossless	1.2.840.10008.1.2.5		

STORAGE-SCU will propose Presentation Contexts only for the SOP Class of the series that is to be transferred.

For that SOP Class, STORAGE-SCU will propose multiple Presentation Contexts, one for the pre-configured "preferred syntax" of the remote AE and two more for each little endian syntax.

For Abstract Syntaxes that do not contain pixel data, in particular that do not contain "PixelData" elements, only Presentation Contexts for uncompressed transfer syntaxes will be proposed.

3.2.1.3.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.1.3.1.3 SOP specific conformance for SOP class(es)

3.2.1.3.1.3.1 SOP specific conformance to storage SOP classes

STORAGE-SCU provides standard conformance to the Storage Service Class.

3.2.1.3.1.3.2 Presentation context acceptance criterion

STORAGE-SCU does not accept associations.

3.2.1.3.1.3.3 Transfer syntax selection policies

STORAGE-SCU prefers the pre-configured "preferred" Transfer Syntax for a remote AE. If offered a choice of Transfer Syntaxes in the accepted Presentation Contexts, it will apply the following priority to the choice of Presentation Context to use for the C-STORE operation:

- The "preferred" Transfer Syntax
- Transfer Syntax of first accepted presentation context
- First encountered explicit Transfer Syntax
- Default Transfer Syntax (Little Endian Implicit)

3.2.1.3.1.3.4 Response status

STORAGE-SCU will behave as described in the Table below in response to the status returned in the C-STORE response command message.

Table 3-6 Response Status for STORAGE-SCU

Service status	Further meaning	Status code(s)	Behavior
Refused	Out of Resources	A7xx	Ignored
Error	Data Set does not match SOP Class	A9xx	Ignored
	Cannot understand	Cxxx	Ignored
Warning	Coercion of Data Elements	B000	Ignored
	Data Set does not match SOP Class	B007	Ignored
	Elements Discarded	B006	Ignored
Success		0000	Ignored

3.2.1.4 Association acceptance policy

STORAGE-SCU does not accept associations.

3.2.2 WORKLIST-SCU

3.2.2.1 SOP classes

WORKLIST-SCU provide Standard Conformance to the following SOP Classes:

Table 3-7 SOP Classes Supported by WORKLIST-SCU

SOP Class Name	SOP Class UID	SCU	SCP
Modality Worklist Information Model - FIND	1.2.840.10008.5.1.4.31	Yes	No

3.2.2.2 Association policies

3.2.2.2.1 General

WORKLIST-SCU initiates but never accepts associations.

Table 3-8 Maximum PDU Size Received as SCP for WORKLIST-SCU

Maximum PDU size received as SCP	
Maximum PDU size received	Unlimited

3.2.2.2.2 Number of associations

Table 3-9 Number of Associations as SCP for WORKLIST-SCU

Number of associations as SCP	
Maximum number of simultaneous associations	1

3.2.2.2.3 Asynchronous nature

WORKLIST-SCU will only allow a single outstanding operation on an Association. Therefore, WORKLIST-SCU will not perform asynchronous operations window negotiation.

3.2.2.2.4 Implementation identifying information

Table 3-10 DICOM Implementation Class and Version for WORKLIST-SCU

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

3.2.2.3 Association initiation policy

WORKLIST-SCU attempts to initiate a new association when the user performs the worklist query action from the user interface.

3.2.2.3.1 Send worklist request

3.2.2.3.1.1 Description and sequencing of activities

The request for a Worklist is initiated by user interaction, i.e. pressing the "DICOM Modality Worklist" button in the "Import" or "Modify" dialog. A dialog to enter search criteria is opened and an interactive query can be performed. Within the interactive queries a single attempt will be made to query the remote AE. If the query fails, for whatever reason, no retry will be performed.

3.2.2.3.1.2 Proposed presentation contexts

Table 3-11 Proposed Presentation Contexts for WORKLIST-SCU

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-7	See Table 4-7	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

WORKLIST-SCU will propose a single Presentation Context for Abstract Syntax with all of the supported Transfer Syntaxes.

3.2.2.3.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.2.3.1.3 SOP specific conformance for SOP class(es)

3.2.2.3.1.3.1 SOP specific conformance to workflow management SOP classes

WORKLIST-SCU provides standard conformance to the C-FIND Service Class.

No CANCEL requests are ever issued.

Unexpected attributes returned in a C-FIND response (those not requested) are ignored. Requested return attributes not returned by the SCP are ignored. Non-matching responses returned by the SCP due to unsupported (hopefully optional) matching keys are not filtered locally by the WORKLIST-SCU and thus will still be presented in the Worklist query result. No attempt is made to filter out duplicate responses.

The pre-configured Specific Character Set will always be included in each query. If present in the response, the Specific Character Set will be used to identify character sets other than the default character set for display of strings and the application of the information to the corresponding instances.

Table 3-12 Worklist Request Identifier

Module name Attribute name	Tag	VR	M	R	Q	D
SOP Common						
Specific Character Set	(0008,0005)	CS		X		
Scheduled Procedure Step						
Scheduled Procedure Step Sequence	(0040,0100)	SQ				
> Modality	(0008,0060)	CS	S	X	X	X
> Scheduled Station AET	(0040,0001)	AE				
> Scheduled Procedure Step Start Date	(0040,0002)	DA	S,R	X	X	X
> Scheduled Procedure Step Start Time	(0040,0003)	TM		X		X
> Scheduled Procedure Step Description	(0040,0007)	LO		X		X
> Scheduled Station Name	(0040,0010)	SH	(S)	X		
Requested Procedure						
Study Instance UID	(0020,000D)	UI		X		
Imaging service Request						
Accession Number	(0008,0050)	SH	S,*	X	X	X
Patient Identification						
Patient Name	(0010,0010)	PN	S,*	X	X	X
Patient ID	(0010,0020)	LO	S,*	X	X	X
Patient Demographic						
Patient's Birth Date	(0010,0030)	DA		X		X
Patient's Sex	(0010,0040)	CS		X		X

The above table should be read as follows:

- Module Name** The name of the associated module for supported worklist attributes.
- Attribute Name** Attributes supported to build a WORKLIST-SCU Worklist Request Identifier.
- Tag** DICOM tag for this attribute.
- VR** DICOM VR for this attribute.
- M** Matching key. An "S" will indicate that WORKLIST-SCU will supply an attribute value for Single Value Matching, an "R" will indicate Range Matching and a "*" will denote wild card matching. It can be configured if "Scheduled Station Name" is additionally supplied "(S)".
- R** Return keys. An "x" will indicate that WORKLIST-SCU will supply this attribute as Return Key with zero length for Universal Matching.
- Q** Interactive Query Key. An "x" will indicate that WORKLIST-SCU will supply this attribute as matching key, if entered in the Worklist Query dialog. For

example, the Patient Name can be entered thereby restricting Worklist responses to the patient.

- D** Displayed keys. An "x" indicates that this worklist attribute is displayed to the user in the Worklist Query dialog.

3.2.2.3.1.3.2 Presentation context acceptance criterion

WORKLIST-SCU does not accept associations.

3.2.2.3.1.3.3 Transfer syntax selection policies

WORKLIST-SCU prefers explicit Transfer Syntaxes. Only a single Presentation Context with all supported Transfer Syntaxes is offered and therefore, if accepted, WORKLIST-SCU will use the explicit syntax.

3.2.2.3.1.3.4 Response status

WORKLIST-SCU will behave as described in the Table below in response to the status returned in the C-FIND response command message.

Table 3-13 Response Status for WORKLIST-SCU

Service status	Further meaning	Status code(s)	Behavior
Refused	Out of Resources	A700	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
Error	Data Set does not match SOP Class	A900	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
	Unable to process	Cxxx	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
Success	Matching is complete - No final Identifier is supplied	0000	Current query is terminated.
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier used to populate query result browser.
	Matches are continuing - Warning that one or more Optional Keys were not	FF01	Identifier used to populate query result browser.

Service status	Further meaning	Status code(s)	Behavior
	supported for existence and/or matching for this Identifier		

3.2.2.4 Association acceptance policy

WORKLIST-SCU does not accept associations.

3.2.3 PRINT-SCU

3.2.3.1 SOP classes

PRINT-SCU provide Standard Conformance to the following SOP Classes:

Table 3-14 SOP Classes Supported by PRINT-SCU

SOP Class Name	SOP Class UID	SCU	SCP
Basic Grayscale Print Management Meta	1.2.840.10008.5.1.1.9	Yes	No

3.2.3.2 Association policies

3.2.3.2.1 General

PRINT-SCU initiates but never accepts associations.

Table 3-15 Maximum PDU Size Received as SCP for PRINT-SCU

Maximum PDU size received as SCP	
Maximum PDU size received	Unlimited

3.2.3.2.2 Number of associations

Table 3-16 Number of Associations as SCP for PRINT-SCU

Number of associations as SCP	
Maximum number of simultaneous associations	1

3.2.3.2.3 Asynchronous nature

PRINT-SCU will only allow a single outstanding operation on an Association. Therefore, PRINT-SCU will not perform asynchronous operations window negotiation.

3.2.3.2.4 Implementation identifying information

Table 3-17 DICOM Implementation Class and Version for PRINT-SCU

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

3.2.3.3 Association initiation policy

PRINT-SCU attempts to initiate a new association for each page when the user performs the print action from the user interface.

3.2.3.3.1 Send print request

3.2.3.3.1.1 Description and sequencing of activities

A user composes images onto pages and requests them to be sent to a specific hardcopy device (remote AE). The user can select the desired film format, orientation and whether to print additional information like patient and study information. Each page is internally processed and converted to a "STANDARD/1,1" page.

The PRINT-SCU AE is invoked by the application for each page. Only a single page is handled per association. If no association to the printer can be established, the print-job of the individual page is switched to a failed state and the user informed.

3.2.3.3.1.2 Proposed presentation contexts

Table 3-18 Proposed Presentation Contexts for WORKLIST-SCU

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-14	See Table 4-14	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

PRINT-SCU will propose a single Presentation Context for the Abstract Syntax with all of the supported Transfer Syntaxes.

3.2.3.3.1.2.1 Extended Negotiation

No extended negotiation is performed.

3.2.3.3.1.3 SOP specific conformance for SOP class(es)

3.2.3.3.1.3.1 SOP specific conformance for the printer SOP class

PRINT-SCU supports the following DIMSE operations and notifications for the Printer SOP Class:

- N-GET

Details of the supported attributes and status handling behavior are described in the following subsections.

3.2.3.3.1.3.2 Printer SOP class operations (N-GET)

PRINT-SCU uses the Printer SOP Class N-GET operation to obtain information about the current printer status. The attributes obtained via N-GET are listed in the Table below:

Table 3-19 Printer SOP Class N-GET Request Attributes

Attribute name	Tag	VR	Value	Presence of value	Source
Printer Status	(2110,0010)	CS	Provided by Printer	ALWAYS	Printer
Printer Status Info	(2110,0020)	CS	Provided by Printer	ALWAYS	Printer

The Printer Status information is evaluated as follows:

- If Printer status (2110,0010) is NORMAL, the print-job continues to be printed.
- If Printer status (2110,0010) is FAILURE, the print-job is marked as failed. The contents of Printer Status Info (2110,0020) is logged.
- If Printer status (2110,0010) is WARNING, the print-job continues to be printed. The contents of Printer Status Info (2110,0020) is logged.

The behavior of PRINT-SCU when encountering status codes in a N-GET response is summarized in the Table below:

Table 3-20 Printer SOP Class N-GET Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The request to get printer status information was success.
Warning	Warning	Any status WARNING	The request to get printer status information was successful. Warning information are ignored.

Service status	Further meaning	Status code(s)	Behavior
*	*	Any other status code	The Association is aborted using A-RELEASE and the print-job is marked as failed. The status meaning is logged.

3.2.3.3.1.3.3 SOP specific conformance for the film session SOP class

PRINT-SCU supports the following DIMSE operations and notifications for the Film Session SOP Class:

- N-CREATE
- N-DELETE

Details of the supported attributes and status handling behavior are described in the following subsections.

3.2.3.3.1.3.4 Film session SOP class operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below:

Table 3-21 Film Session SOP Class N-CREATE Request Attributes

Attribute name	Tag	VR	Value	Presence of value	Source
Number of Copies	(2000,0010)	IS	1	ALWAYS	Configuration*
Print Priority	(2000,0020)	CS	HIGH MED (default) LOW	ALWAYS	Configuration*
Medium Type	(2000,0030)	CS	PAPER CLEAR FILM BLUE FILM (default) MAMMO CLEAR FILM MAMMO BLUE FILM	ALWAYS	Configuration*
Film Destination	(2000,0040)	CS	PROCESSOR (default) MAGAZINE	ALWAYS	Configuration*

* The values can be changed in the application configuration by the administrator. The user cannot modify the configuration in the user interface.

The behavior of PRINT-SCU when encountering status codes in an N-CREATE response is summarized in the Table below:

Table 3-22 Film Session SOP Class N-CREATE Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The SCP has completed the operation successfully. The print job continues.
Warning	Warning	Any status WARNING	The SCP has completed the operation successfully. Warning information are ignored. The print job continues.
*	*	Any other status code	The Association is aborted using A-RELEASE and the print-job is marked as failed. The status meaning is logged.

3.2.3.3.1.3.5 Film session SOP class operations (N-DELETE)

The behavior of PRINT-SCU when encountering status codes in an N-DELETE response is summarized in the Table below:

Table 3-23 Film Session SOP Class N-DELETE Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The SCP has completed the operation successfully.
Warning	Warning	Any status WARNING	The SCP has completed the operation successfully. Warning information are ignored.
*	*	Any other status code	The SCP failed completing the operation. The error is ignored.

3.2.3.3.1.3.6 SOP specific conformance for the film box SOP class

PRINT-SCU supports the following DIMSE operations and notifications for the Film Box SOP Class:

- N-CREATE
- N-ACTION

Details of the supported attributes and status handling behavior are described in the following subsections.

3.2.3.3.1.3.6.1 Film box SOP class operations (N-CREATE)

The attributes supplied in an N-CREATE Request are listed in the Table below:

Table 3-24 Film Box SOP Class N-CREATE Request Attributes

Attribute name	Tag	VR	Value	Presence of value	Source
Image Display Format	(2010,0010)	CS	STANDARD\1,1	ALWAYS	Auto
Referenced Film Session Sequence	(2010,0500)	SQ		ALWAYS	Auto
> Referenced SOP Class UID	(0008,1150)	UI	1.2.840.10008.5.1.1.1	ALWAYS	Auto
> Referenced SOP Instance UID	(0008,1155)	UI	From created Film Session SOP Instance	ALWAYS	Auto
Film Orientation	(2010,0040)	CS	PORTRAIT LANDSCAPE	ALWAYS	User
Film Size ID	(2010,0050)	CS	8INX10IN 8_5INX11IN 8_5INX35IN 8_5INX47IN 10INX12IN 10INX14IN 11INX14IN 11INX17IN 12INX35IN 12INX47IN 13INX18IN 13INX35IN 13INX47IN 14INX14IN 14INX17IN 24CMX24CM 24CMX30CM A4 A3 A3 NOBI	ALWAYS	User
Magnification Type	(2010,0060)	CS	REPLICATE BILINEAR CUBIC (BICUBIC*) NONE		User**
Border Density	(2010,0100)	CS	150 20 BLACK WHITE		User**
Empty Image Density	(2010,0110)	CS	0		Configuration** *
Minimum Density	(2010,0120)	US	FFFF		Fixed

Attribute name	Tag	VR	Value	Presence of value	Source
Maximum Density	(2010,0130)	US	FFFF		Fixed
Trim	(2010,0140)	CS	NO		Configuration** *

* For compatibility reasons with old SCP implementations the user can also configure BICUBIC which is not a DICOM defined value for the Magnification Type.

** Values are pre-configured per remote AE.

*** The values can be changed in the application configuration by the administrator. The user cannot modify the configuration in the user interface.

The behavior of PRINT-SCU when encountering status codes in an N-CREATE response is summarized in the Table below:

Table 3-25 Film Box SOP Class N-CREATE Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The SCP has completed the operation successfully. The print job continues.
Warning	Warning	Any WARNING status	The SCP has completed the operation successfully. Warning information are ignored. The print job continues.
*	*	Any other status code	The Association is aborted using A-RELEASE and the print-job is marked as failed. The status meaning is logged.

3.2.3.3.1.3.7 Film box SOP class operations (N-ACTION)

The behavior of PRINT-SCU when encountering status codes in an N-DELETE response is summarized in the Table below:

Table 3-26 Film Box SOP Class N-ACTION Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The SCP has completed the operation successfully. The page has been accepted for printing.
Warning	Warning	Any WARNING status	The SCP has completed the operation successfully. Warning information are ignored.
*	*	Any other status code	The SCP failed completing the operation. The error is logged.

3.2.3.3.1.3.8 SOP specific conformance for the image box SOP class

PRINT-SCU supports the following DIMSE operations and notifications for the Image Box SOP Class:

- N-SET

Details of the supported attributes and status handling behavior are described in the following subsections.

3.2.3.3.1.3.9 Image box SOP class operations (N-SET)

The attributes supplied in an N-SET request are listed in the Table below:

Table 3-27 Image Box SOP Class N-CREATE Request Attributes

Attribute name	Tag	VR	Value	Presence of value	Source
Image Box Position	(2010,0010)	US		ALWAYS	Auto
Requested Image Size	(2020,0030)	DS	Width in mm	Only when "Real-size" printing is selected	Auto
Requested Decimate Crop Behavior	(2020,0040)	CS	CROP	Only when "Real-size" printing is selected	Auto
Basic Grayscale Image Sequence	(2020,0110)	SQ		ALWAYS	Auto
> Samples Per Pixel	(0028,0002)	US	1	ALWAYS	Auto
> Photometric Interpretation	(0028,0004)	CS	MONOCHROME2	ALWAYS	Auto
> Rows	(0028,0010)	US	Depends on image	ALWAYS	Auto
> Columns	(0028,0011)	US	Depends on image	ALWAYS	Auto
> Bits Allocated	(0028,0100)	US	8	ALWAYS	Auto
> Bits Stored	(0028,0101)	US	8	ALWAYS	Auto
> High Bit	(0028,0102)	US	7	ALWAYS	Auto
> Pixel Representation	(0028,0103)	US	0	ALWAYS	Auto
> Pixel Data	(7FE0,0010)	OB	Pixels of rendered page	ALWAYS	Auto

The behavior of PRINT-SCU when encountering status codes in an N-SET response is summarized in the Table below:

Table 3-28 Image Box SOP Class N-CREATE Response Status Behavior

Service status	Further meaning	Status code(s)	Behavior
Success	Success	0000	The SCP has completed the operation successfully. Image successfully stored in the Image Box.
Warning	Warning	Any status WARNING	The SCP has completed the operation successfully. Image successfully stored in the Image Box. Warning information are ignored.
*	*	Any other status code	The Association is aborted using A-RELEASE and the print-job is marked as failed. The status meaning is logged.

3.2.3.4 Association acceptance policy

PRINT-SCU does not accept associations.

3.2.4 FIND-SCU

3.2.4.1 SOP classes

FIND-SCU provide Standard Conformance to the following SOP Classes:

Table 3-29 SOP Classes Supported by FIND-SCU

SOP Class Name	SOP Class UID	SCU	SCP
Study Root Query/Retrieve Information Model -FIND	1.2.840.10008.5.1.4.1.2.2.1	Yes	No

3.2.4.2 Association policies

3.2.4.2.1 General

FIND-SCU initiates but never accepts associations.

Table 3-30 Maximum PDU Size Received as SCP for FIND-SCU

Maximum PDU size received as SCP	
Maximum PDU size received	Unlimited

3.2.4.2.2 Number of associations

Table 3-31 Number of Associations as SCP for FIND-SCU

Number of associations as SCP

Maximum number of simultaneous associations	1
---	---

3.2.4.2.3 Asynchronous nature

FIND-SCU will only allow a single outstanding operation on an Association. Therefore, FIND-SCU will not perform asynchronous operations window negotiation.

3.2.4.2.4 Implementation identifying information

Table 3-32 DICOM Implementation Class and Version for FIND-SCU

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

3.2.4.3 Association initiation policy

FIND-SCU attempts to initiate a new association when the user performs the query action from the user interface. Recursive queries for lower query levels in the hierarchy will be performed on a new association.

3.2.4.3.1 Send query request

3.2.4.3.1.1 Description and sequencing of activities

A single attempt will be made to query the remote AE. If the query fails, for whatever reason, no retry will be performed.

3.2.4.3.1.2 Proposed presentation contexts

Table 3-33 Proposed Presentation Contexts for FIND-SCU

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-31	See Table 4-31	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

FIND-SCU will propose a single Presentation Context for the Abstract Syntax with all of the supported Transfer Syntaxes.

3.2.4.3.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.4.3.1.3 SOP specific conformance for SOP class(es)

3.2.4.3.1.3.1 SOP specific conformance to query SOP classes

FIND-SCU provides standard conformance to the Query Service Class.

Only a single information model, Study Root, is supported.

All queries are initiated at the highest level of the information model (the STUDY level). For each study unfolded by the user new query at the next lower levels (the SERIES and then IMAGE levels) is initiated in order to completely elucidate the "tree" of instances available on the remote AE (from which the user may subsequently request a retrieval at any level).

No CANCEL requests are ever issued.

Unexpected attributes returned in a C-FIND response (those not requested) are ignored. Requested return attributes not returned by the SCP are ignored. Non-matching responses returned by the SCP due to unsupported (hopefully optional) matching keys are not filtered locally by the FIND-SCU and thus will still be presented in the browser. No attempt is made to filter out duplicate responses.

Specific Character Set will always be included at every query level. If present in the response, Specific Character Set will be used to identify character sets other than the default character set for display of strings in the browser.

Table 3-34 Study Root Identifier for FIND-SCU

Name	Tag	Type of matching
STUDY level		
Patient Name	(0010,0010)	S,*,U
Patient ID	(0010,0020)	S,*,U
Patient Date of Birth	(0010,0030)	S,*,U,R
Patient Sex	(0010,0040)	S,*,U
Study ID	(0020,0010)	S,*,U
Study Description	(0008,1030)	S,*,U
Modalities in Study	(0008,0061)	S,*,U
Study Date	(0008,0020)	S,*,U,R
Study Time	(0008,0030)	S,*,U,R
Referring Physician's Name	(0008,0090)	S,*,U
Accession Number	(0008,0050)	S,*,U
Study Instance UID	(0020,000D)	UNIQUE
SERIES level		

Name	Tag	Type of matching
Series Number	(0020,0011)	S,*,U
Series Description	(0008,103E)	S,*,U
Modality	(0008,0060)	S,*,U
Series Date	(0008,0021)	S,*,U,R
Series Time	(0008,0031)	S,*,U,R
Body Part Examined	(0018,0015)	S,*,U
Series Instance UID	(0020,000E)	UNIQUE
IMAGE level		
Instance Number	(0020,0013)	S,*,U
Image Type	(0008,0008)	S,*,U
Instance Creation Date	(0008,0012)	S,*,U,R
Instance Creation Time	(0008,0013)	S,*,U,R
Acquisition Date	(0008,0022)	S,*,U,R
Acquisition Time	(0008,0023)	S,*,U,R
Acquisition Number	(0020,0012)	S,*,U
Image Position Patient	(0020,0032)	S,*,U
Image Orientation Patient	(0020,0037)	S,*,U
Slice Location	(0020,1041)	S,*,U
Number of Frames	(0028,0008)	S,*,U
SOP Instance UID	(0008,0018)	UNIQUE
SOP Class UID	(0008,0016)	NONE
All levels		
Specific Character Set	(0008,0005)	NONE
Query Retrieve Level	(0008,0052)	S

Types of Matching: the types of Matching supported by the C-FIND SCU.

- S** Single Value Matching
- R** Range Matching
- *** Wildcard matching
- U** Universal Matching
- NONE** No matching supported
- UNIQUE** Unique Key for that query level (Universal Matching or Single Value Matching depending on the query level)

3.2.4.3.1.3.2 Presentation context acceptance criterion

FIND-SCU does not accept associations.

3.2.4.3.1.3.3 Transfer syntax selection policies

FIND-SCU prefers explicit Transfer Syntaxes. Only a single Presentation Context with all supported Transfer Syntaxes is offered and therefore, if accepted, FIND-SCU will use the explicit syntax.

Response status

FIND-SCU will behave as described in the Table below in response to the status returned in the C-FIND response command message.

Table 3-35 Response Status for WORKLIST-SCU

Service status	Further meaning	Status code(s)	Behavior
Refused	Out of Resources	A700	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
Error	Data Set does not match SOP Class	A900	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
	Unable to process	Cxxx	Current query is terminated. The status meaning is logged. Any additional error information in the Response will be logged.
Success	Matching is complete - No final Identifier is supplied	0000	Current query is terminated.
Pending	Matches are continuing - Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys	FF00	Identifier used to populate query result browser.
	Matches are continuing - Warning that one or more Optional Keys were not supported for existence and/or matching for this Identifier	FF01	Identifier used to populate query result browser.

3.2.4.4 Association acceptance policy

FIND-SCU does not accept associations.

3.2.5 MOVE-SCU

3.2.5.1 SOP classes

MOVE-SCU provide Standard Conformance to the following SOP Classes:

Table 3-36 SOP Classes Supported by MOVE-SCU

SOP Class Name	SOP Class UID	SCU	SCP
Study Root Query/Retrieve Information Model -MOVE	1.2.840.10008.5.1.4.1.2.2.2	Yes	No

3.2.5.2 Association policies

3.2.5.2.1 General

MOVE-SCU initiates but never accepts associations.

Table 3-37 Maximum PDU Size Received as SCP for MOVE-SCU

Maximum PDU size received as SCP	
Maximum PDU size received	Unlimited

3.2.5.2.2 Number of associations

Table 3-38 Number of Associations as SCP for MOVE-SCU

Number of associations as SCP	
Maximum number of simultaneous associations	1

3.2.5.2.3 Asynchronous nature

MOVE-SCU will only allow a single outstanding operation on an Association. Therefore, MOVE-SCU will not perform asynchronous operations window negotiation.

3.2.5.2.4 Implementation identifying information

Table 3-39 DICOM Implementation Class and Version for MOVE-SCU

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

3.2.5.3 Association initiation policy

MOVE-SCU attempts to initiate a new association when the user performs the retrieve action from the user interface.

3.2.5.3.1 Send retrieve request

3.2.5.3.1.1 Description and sequencing of activities

For the entity (study or series) selected from the user interface to be retrieved, a single attempt will be made to retrieve it from the selected remote AE. If the retrieve fails, for whatever reason, no retry will be performed.

3.2.5.3.1.2 Proposed presentation contexts

Table 3-40 Proposed Presentation Contexts for MOVE-SCU

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-36	See Table 4-36	Implicit VR Little Endian	1.2.840.10008.1.2	SCU	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

MOVE-SCU will propose a single Presentation Context for the Abstract Syntax with all of the supported Transfer Syntaxes.

3.2.5.3.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.5.3.1.3 SOP specific conformance for SOP class(es)

3.2.5.3.1.3.1 SOP specific conformance to retrieve SOP classes

MOVE-SCU provides standard conformance to the Retrieve Service Class.

Only a single information model, Study Root, is supported.

A retrieval will be performed at the STUDY or SERIES level depending on what level of entity has been selected by the user in the browser.

No CANCEL requests are ever issued.

The instances are retrieved to the current application's local database by specifying the destination as the AE Title of the STORE-SCP AE of the local application. This implies that the remote C-MOVE SCP must be preconfigured to determine the presentation address corresponding to the STORE-SCP AE. The STORE-SCP AE will accept storage requests addressed to it from anywhere, so no pre-configuration of the local application to accept from the remote AE is necessary (except in so far as it was necessary to configure FIND-SCU).

Table 3-41 Study Root Identifier for MOVE-SCU

Name	Tag	Type of matching
STUDY level		
Study Instance UID	(0020,000D)	UNIQUE
SERIES level		
Series Instance UID	(0020,000E)	UNIQUE
IMAGE level		
SOP Instance UID	(0008,0018)	UNIQUE
All levels		
Query Retrieve Level	(0008,0052)	S

Types of Matching: the types of Matching supported by the C-FIND SCU.

- S** Single Value Matching
- UNIQUE** Unique Key for that query level (Universal Matching or Single Value Matching depending on the query level)

3.2.5.3.1.3.2 Presentation context acceptance criterion

MOVE-SCU does not accept associations.

3.2.5.3.1.3.3 Transfer syntax selection policies

MOVE-SCU prefers explicit Transfer Syntaxes. Only a single Presentation Context with all supported Transfer Syntaxes is offered and therefore, if accepted, MOVE-SCU will use the explicit syntax.

3.2.5.3.1.3.4 Response status

MOVE-SCU will behave as described in the Table below in response to the status returned in the C-MOVE response command message.

Table 3-42 Response Status for WORKLIST-SCU

Service status	Further meaning	Status code(s)	Behavior
Refused	Out of Resources	A70x	Retrieval is terminated
Error	Data Set does not match SOP Class	A900	Retrieval is terminated
	Unable to process	Cxxx	Retrieval is terminated
Success	Sub-operations Complete – No Failures	0000	Retrieval is terminated

Service status	Further meaning	Status code(s)	Behavior
Cancel Warning	Sub-operations terminated due to Cancel Indication	FE00	Retrieval is terminated (should never occur, since cancel is never issued)
	Sub-operations Complete - One or more Failures	B000	Retrieval is terminated
Pending	Sub-operations are continuing	FF00	Retrieval continues

3.2.5.3.1.3.5 Sub-operation dependent behavior

Since the C-MOVE operation is dependent on completion of C-STORE sub-operations that are occurring on a separate association, the question of failure of operations on the other association(s) must be considered.

MOVE-SCU completely ignores whatever activities are taking place in relation to the STORAGE-SCP AE that is receiving the retrieved instances. Once the C-MOVE has been initiated it runs to completion (or failure) as described in the C-MOVE response command message(s). There is no attempt by MOVE-SCU to confirm that instances have actually been successfully received or locally stored.

Whether or not completely or partially successful retrievals are made available in the local database to the user is purely dependent on the success or failure of the C-STORE sub-operations, not on any explicit action by MOVE-SCU.

Whether or not the remote AE attempts to retry any failed C-STORE sub-operations is beyond the control of MOVE-SCU.

If the association on which the C-MOVE was issued is aborted for any reason, whether or not the C-STORE sub-operations continue is dependent on the remote AE; the local STORAGE-SCP will continue to accept associations and storage operations regardless.

3.2.5.4 Association acceptance policy

MOVE-SCU does not accept associations.

3.2.6 STORAGE-SCP

3.2.6.1 SOP classes

STORAGE-SCP provides Standard Conformance to the following SOP Classes:

SOP class name	SOP class UID	SCU	SCP
(RETIRED) Stored Print Storage	1.2.840.10008.5.1.1.27	No	(Yes)*
(RETIRED) Hardcopy Grayscale Image Storage	1.2.840.10008.5.1.1.29	No	Yes

SOP class name	SOP class UID	SCU	SCP
(RETIRED) Hardcopy Color Image Storage	1.2.840.10008.5.1.1.30	No	Yes
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1	No	Yes
Digital X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.1	No	Yes
Digital X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.1.1	No	Yes
Digital Mammography X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.2	No	Yes
Digital Mammography X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.2.1	No	Yes
Digital Intra Oral X-Ray Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.1.3	No	Yes
Digital Intra Oral X-Ray Image Storage For Processing	1.2.840.10008.5.1.4.1.1.1.3.1	No	Yes
CT Image Storage	1.2.840.10008.5.1.4.1.1.2	No	Yes
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1	No	Yes
Legacy Converted Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.2	No	(Yes)*
(RETIRED) Ultrasound Multiframe Image Storage	1.2.840.10008.5.1.4.1.1.3	No	Yes
Ultrasound Multiframe Image Storage	1.2.840.10008.5.1.4.1.1.3.1	No	Yes
MR Image Storage	1.2.840.10008.5.1.4.1.1.4	No	Yes
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1	No	Yes
MR Spectroscopy Storage	1.2.840.10008.5.1.4.1.1.4.2	No	(Yes)*
Enhanced MR Color Image Storage	1.2.840.10008.5.1.4.1.1.4.3	No	Yes
Legacy Converted Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.4	No	(Yes)*
(RETIRED) Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.5	No	Yes
(RETIRED) Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6	No	Yes
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1	No	Yes
Enhanced US Volume Storage	1.2.840.10008.5.1.4.1.1.6.2	No	(Yes)*
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	No	Yes
Multiframe Single Bit Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.1	No	Yes
Multiframe Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2	No	Yes
Multiframe Grayscale Word Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.3	No	Yes
Multiframe True Color Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.4	No	Yes

SOP class name	SOP class UID	SCU	SCP
(RETIRED) Standalone Overlay Storage	1.2.840.10008.5.1.4.1.1.8	No	(Yes)*
(RETIRED) Standalone Curve Storage	1.2.840.10008.5.1.4.1.1.9	No	(Yes)*
Twelve Lead ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.1	No	(Yes)*
General ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.2	No	(Yes)*
Ambulatory ECG Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1.3	No	(Yes)*
Hemodynamic Waveform Storage	1.2.840.10008.5.1.4.1.1.9.2.1	No	(Yes)*
Cardiac Electrophysiology Waveform Storage	1.2.840.10008.5.1.4.1.1.9.3.1	No	(Yes)*
Basic Voice Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.1	No	(Yes)*
General Audio Waveform Storage	1.2.840.10008.5.1.4.1.1.9.4.2	No	(Yes)*
Arterial Pulse Waveform Storage	1.2.840.10008.5.1.4.1.1.9.5.1	No	(Yes)*
Respiratory Waveform Storage	1.2.840.10008.5.1.4.1.1.9.6.1	No	(Yes)*
(RETIRED) Standalone Modality LUT Storage	1.2.840.10008.5.1.4.1.1.10	No	(Yes)*
(RETIRED) Standalone VOI LUT Storage	1.2.840.10008.5.1.4.1.1.11	No	(Yes)*
Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.1	No	Yes**
Color Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.2	No	(Yes)*
Pseudo Color Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.3	No	(Yes)*
Blending Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.4	No	(Yes)*
XA XRF Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.5	No	(Yes)*
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1	No	Yes
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1	No	Yes
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2	No	Yes
Enhanced XRF Image Storage	1.2.840.10008.5.1.4.1.1.12.2.1	No	(Yes)*
(RETIRED) X-Ray Angiographic BiPlane Image Storage	1.2.840.10008.5.1.4.1.1.12.3	No	Yes
X-Ray 3D Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.13.1.1	No	(Yes)*
X-Ray 3D Craniofacial Image Storage	1.2.840.10008.5.1.4.1.1.13.1.2	No	(Yes)*
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3	No	Yes

SOP class name	SOP class UID	SCU	SCP
Intravascular Optical Coherence Tomography Image Storage For Presentation	1.2.840.10008.5.1.4.1.1.14.1	No	(Yes)*
Intravascular Optical Coherence Tomography Image Storage For Processing	1.2.840.10008.5.1.4.1.1.14.2	No	(Yes)*
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20	No	Yes
Raw Data Storage	1.2.840.10008.5.1.4.1.1.66	No	(Yes)*
Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.1	No	(Yes)*
Spatial Fiducials Storage	1.2.840.10008.5.1.4.1.1.66.2	No	(Yes)*
Deformable Spatial Registration Storage	1.2.840.10008.5.1.4.1.1.66.3	No	(Yes)*
Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.4	No	(Yes)*
Surface Segmentation Storage	1.2.840.10008.5.1.4.1.1.66.5	No	(Yes)*
Real World Value Mapping Storage	1.2.840.10008.5.1.4.1.1.67	No	(Yes)*
Surface Scan Mesh Storage	1.2.840.10008.5.1.4.1.1.68.1	No	(Yes)*
Surface Scan Point Cloud Storage	1.2.840.10008.5.1.4.1.1.68.2	No	(Yes)*
(RETIRED) VL Image Storage	1.2.840.10008.5.1.4.1.1.77.1	No	Yes
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1	No	Yes
Video Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1.1	No	Yes
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2	No	Yes
Video Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2.1	No	Yes
VL Slide Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3	No	(Yes)*
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4	No	Yes
Video Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4.1	No	Yes
Ophthalmic Photography 8Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1	No	Yes
Ophthalmic Photography 16Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2	No	Yes
Stereometric Relationship Storage	1.2.840.10008.5.1.4.1.1.77.1.5.3	No	(Yes)*
Ophthalmic Tomography Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.4	No	(Yes)*
VL Whole Slide Microscopy Image Storage	1.2.840.10008.5.1.4.1.1.77.1.6	No	(Yes)*
(RETIRED) VL Multi Frame Image Storage	1.2.840.10008.5.1.4.1.1.77.2	No	Yes
Lensometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.1	No	(Yes)*
Autorefractometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.2	No	(Yes)*
Keratometry Measurements Storage	1.2.840.10008.5.1.4.1.1.78.3	No	(Yes)*
Subjective Refraction Measurements Storage	1.2.840.10008.5.1.4.1.1.78.4	No	(Yes)*

SOP class name	SOP class UID	SCU	SCP
Visual Acuity Measurements Storage	1.2.840.10008.5.1.4.1.1.78.5	No	(Yes)*
Spectacle Prescription Report Storage	1.2.840.10008.5.1.4.1.1.78.6	No	(Yes)*
Ophthalmic Axial Measurements Storage	1.2.840.10008.5.1.4.1.1.78.7	No	(Yes)*
Intraocular Lens Calculations Storage	1.2.840.10008.5.1.4.1.1.78.8	No	(Yes)*
Macular Grid Thickness And Volume Report Storage	1.2.840.10008.5.1.4.1.1.79.1	No	(Yes)*
Ophthalmic Visual Field Static Perimetry Measurements Storage	1.2.840.10008.5.1.4.1.1.80.1	No	(Yes)*
Ophthalmic Thickness Map Storage	1.2.840.10008.5.1.4.1.1.81.1	No	(Yes)*
Corneal Topography Map Storage	1.2.840.10008.5.1.4.1.1.82.1	No	(Yes)*
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	No	Yes
Enhanced SR Storage	1.2.840.10008.5.1.4.1.1.88.22	No	Yes
Comprehensive SR Storage	1.2.840.10008.5.1.4.1.1.88.33	No	Yes
Comprehensive 3D SR Storage	1.2.840.10008.5.1.4.1.1.88.34	No	(Yes)*
Procedure Log Storage	1.2.840.10008.5.1.4.1.1.88.40	No	(Yes)*
Mammography CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.50	No	Yes
Key Object Selection Document Storage	1.2.840.10008.5.1.4.1.1.88.59	No	Yes
Chest CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.65	No	(Yes)*
X-Ray Radiation Dose SR Storage	1.2.840.10008.5.1.4.1.1.88.67	No	Yes
Colon CAD SR Storage	1.2.840.10008.5.1.4.1.1.88.69	No	(Yes)*
Implantation Plan SR Document Storage	1.2.840.10008.5.1.4.1.1.88.70	No	(Yes)*
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	No	Yes
Encapsulated CDA Storage	1.2.840.10008.5.1.4.1.1.104.2	No	(Yes)*
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128	No	Yes
Legacy Converted Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.128.1	No	(Yes)*
(RETIRED) Standalone PET Curve Storage	1.2.840.10008.5.1.4.1.1.129	No	(Yes)*
Enhanced PET Image Storage	1.2.840.10008.5.1.4.1.1.130	No	(Yes)*
Basic Structured Display Storage	1.2.840.10008.5.1.4.1.1.131	No	(Yes)*
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1	No	(Yes)*
RT Dose Storage	1.2.840.10008.5.1.4.1.1.481.2	No	(Yes)*
RT Structure Set Storage	1.2.840.10008.5.1.4.1.1.481.3	No	(Yes)*
RT Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.4	No	(Yes)*
RT Plan Storage	1.2.840.10008.5.1.4.1.1.481.5	No	(Yes)*

SOP class name	SOP class UID	SCU	SCP
RT Brachy Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.6	No	(Yes)*
RT Treatment Summary Record Storage	1.2.840.10008.5.1.4.1.1.481.7	No	(Yes)*
RT Ion Plan Storage	1.2.840.10008.5.1.4.1.1.481.8	No	(Yes)*
RT Ion Beams Treatment Record Storage	1.2.840.10008.5.1.4.1.1.481.9	No	(Yes)*
Generic Implant Template Storage	1.2.840.10008.5.1.4.43.1	No	(Yes)*
Implant Assembly Template Storage	1.2.840.10008.5.1.4.44.1	No	(Yes)*
Implant Template Group Storage	1.2.840.10008.5.1.4.45.1	No	(Yes)*
RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.7	No	(Yes)*
Color Palette Storage	1.2.840.10008.5.1.4.39.1	No	(Yes)*
(DRAFT) SR Text Storage	1.2.840.10008.5.1.4.1.1.88.1	No	(Yes)*
(DRAFT) SR Audio Storage	1.2.840.10008.5.1.4.1.1.88.2	No	(Yes)*
(DRAFT) SR Detail Storage	1.2.840.10008.5.1.4.1.1.88.3	No	(Yes)*
(DRAFT) SR Comprehensive Storage	1.2.840.10008.5.1.4.1.1.88.4	No	(Yes)*
(DRAFT) Waveform Storage	1.2.840.10008.5.1.4.1.1.9.1	No	(Yes)*
(DRAFT) RT Beams Delivery Instruction Storage	1.2.840.10008.5.1.4.34.1	No	(Yes)*

* These SOP classes are not supported by the viewer application, but network transfer (C-STORE) and storage to the local file system will be possible. The acceptance of these classes is disabled by default but can be activated by the administrator if acceptance is desired.

3.2.6.2 Association policies

3.2.6.2.1 General

STORAGE-SCP accepts but never initiates associations.

Table 3-43 Maximum PDU Size Received as SCP for STORAGE-SCP

Maximum PDU size received as SCP	
Maximum PDU size received	128 KB

3.2.6.2.2 Number of associations

Table 3-44 Number of Associations as SCP for STORAGE-SCP

Number of associations as SCP	
Maximum number of simultaneous associations	Configurable

3.2.6.2.3 Asynchronous nature

STORAGE-SCP will only allow a single outstanding operation on an Association. Therefore, STORAGE-SCP will not perform asynchronous operations window negotiation.

3.2.6.2.4 Implementation identifying information

Table 3-45 DICOM Implementation Class and Version for STORAGE-SCP

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0
Implementation Version Name	IIS_368

3.2.6.3 Association initiation policy

STORAGE-SCP does not initiate associations.

3.2.6.4 Association acceptance policy

When STORAGE-SCP accepts an association, it will respond to storage requests.

3.2.6.4.1 Receive storage request

3.2.6.4.1.1 Description and sequencing of activities

As instances are received they are copied to the local file system and a record inserted into the local database. If the received instance is a duplicate of a previously received instance, the old file and database record will not be overwritten and the new one will be discarded.

3.2.6.4.1.2 Accepted presentation contexts

Table 3-46 Acceptable Presentation Contexts for STORAGE-SCP and Receive Storage Request

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-42	See Table 4-42	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		
		JPEG Baseline - Process 1	1.2.840.10008.1.2.4.50		

Presentation context table

	(default 8-bit lossy JPEG compression)	
	JPEG Extended - Process 2 & 4 (default 12-bit lossy JPEG compression)	1.2.840.10008.1.2.4.51
	JPEG Lossless, Non-Hierarchical - Process 14	1.2.840.10008.1.2.4.57
	JPEG Lossless, Non-Hierarchical, First-Order Prediction - Process 14 [Selection Value 1] (default lossless JPEG compression)	1.2.840.10008.1.2.4.70
	JPEG-LS Lossless	1.2.840.10008.1.2.4.80
	JPEG-LS Lossy (Near-Lossless)	1.2.840.10008.1.2.4.81
	JPEG 2000 (Lossless Only)	1.2.840.10008.1.2.4.90
	JPEG 2000 (Lossless or Lossy)	1.2.840.10008.1.2.4.91
	RLE Lossless	1.2.840.10008.1.2.5

3.2.6.4.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.6.4.1.3 SOP specific conformance for SOP class(es)

3.2.6.4.1.3.1 SOP specific conformance to storage SOP classes

STORAGE-SCP provides standard conformance to the Storage Service Class.

3.2.6.4.1.3.2 Presentation context acceptance criterion

STORAGE-SCP will always accept any Presentation Context for the supported SOP Classes with the supported Transfer Syntaxes. More than one proposed Presentation Context will be accepted for the same Abstract Syntax if the Transfer Syntax is supported, whether or not it is the same as another Presentation Context.

3.2.6.4.1.3.3 Transfer syntax selection policies

ECHO-SCP, if offered a choice of Transfer Syntaxes in a Presentation Context, will apply the following priority to the choice of Transfer Syntax:

- first encountered Transfer Syntax,
- default Transfer Syntax

STORAGE-SCP will accept duplicate Presentation Contexts, that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same priority for selecting a Transfer Syntax for each.

3.2.6.4.1.3.4 Response status

STORAGE-SCP will behave as described in the Table below when generating the C-STORE response command message.

Table 3-47 Response Status for STORAGE-SCP

Service status	Further meaning	Status code(s)	Reason
Refused	Out of Resources	A7xx	Internal errors or failed to storage in file system
Error	Data Set does not match SOP Class	A9xx	Never sent
	Cannot understand	Cxxx	Missing or invalid SOP Class/ Instance UIDs
Warning	Coercion of Data Elements	B000	Never sent - no coercion is ever performed
	Data Set does not match SOP Class	B007	Mismatch between C-STORE request and dataset
	Elements Discarded	B006	Never sent
Success		0000	

3.2.7 ECHO-SCP

3.2.7.1 SOP classes

ECHO-SCP provides Standard Conformance to the following SOP Classes:

Table 3-48 SOP Classes Supported By ECHO-SCP

SOP class name	SOP class UID	SCU	SCP
Verification SOP Class	1.2.840.10008.1.1	No	Yes

3.2.7.2 Association policies

3.2.7.2.1 General

ECHO-SCP accepts but never initiates associations.

Table 3-49 Maximum PDU Size Received as SCP for ECHO-SCP

Maximum PDU size received as SCP	
Maximum PDU size received	128 KB

3.2.7.2.2 Number of associations

Table 3-50 Number of Associations as SCP for ECHO-SCP

Number of associations as SCP	
Maximum number of simultaneous associations	Configurable

3.2.7.2.3 Asynchronous nature

ECHO-SCP will only allow a single outstanding operation on an Association. Therefore, ECHO-SCP will not perform asynchronous operations window negotiation.

3.2.7.2.4 Implementation identifying information

Table 3-51 DICOM Implementation Class and Version for ECHO-SCP

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0
Implementation Version Name	IIS_368

3.2.7.3 Association initiation policy

ECHO-SCP does not initiate associations.

3.2.7.4 Association acceptance policy

When ECHO-SCP accepts an association, it will respond to echo requests.

3.2.7.4.1 Receive echo request

3.2.7.4.1.1 Description and sequencing of activities

Each echo request received over an accepted association will result in an echo response.

3.2.7.4.1.2 Accepted presentation contexts

Table 3-52 Acceptable Presentation Contexts for ECHO-SCP and Receive Echo Request

Presentation context table					
Abstract syntax		Transfer syntax		Role	Extended negotiation
Name	UID	Name List	UID list		
See Table 4-48	See Table 4-48	Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		
		Explicit VR Big Endian	1.2.840.10008.1.2.2		

3.2.7.4.1.2.1 Extended negotiation

No extended negotiation is performed.

3.2.7.4.1.3 SOP specific conformance for SOP class(es)

3.2.7.4.1.3.1 SOP specific conformance to verification SOP classes

ECHO-SCP provides standard conformance to the Verification Service Class.

3.2.7.4.1.3.2 Presentation context acceptance criterion

ECHO-SCP will always accept any Presentation Context for the supported SOP Classes with the supported Transfer Syntaxes. More than one proposed Presentation Context will be accepted for the same Abstract Syntax if the Transfer Syntax is supported, whether or not it is the same as another Presentation Context.

3.2.7.4.1.3.3 Transfer syntax selection policies

ECHO-SCP prefers explicit Transfer Syntaxes. If offered a choice of Transfer Syntaxes in a Presentation Context, it will apply the following priority to the choice of Transfer Syntax:

- first encountered explicit Transfer Syntax,
- default Transfer Syntax.

ECHO-SCP will accept duplicate Presentation Contexts, that is, if it is offered multiple Presentation Contexts, each of which offers acceptable Transfer Syntaxes, it will accept all Presentation Contexts, applying the same priority for selecting a Transfer Syntax for each.

3.3 Network interfaces

3.3.1 Physical network interface

The application is indifferent to the physical medium over which TCP/IP executes, which is dependent on the underlying operating system and hardware.

3.3.2 Additional protocols

When host names rather than IP addresses are used in the configuration properties to specify presentation addresses for remote AEs, the application is dependent on the name resolution mechanism of the underlying operating system.

3.3.3 IPv4 and IPv6 support

This product supports IPv4 only.

3.4 Configuration

The configuration can either be changed in the Server Administration application or in the DICOM Settings within the main application.

3.4.1 AE title / presentation address mapping

The Calling AE Title of the local application is configurable in the Server Administration application, and is shared by all of the AEs. The mapping of the logical name by which remote AEs are described in the user interface to Called AE Titles as well as presentation address (hostname or IP address and port number) is configurable in the DICOM Settings of the application.


3.4.2 Parameters

The network listen port for SCP AEs can also be changed in the Server Administration application. Only a single port will be used by all SCP AEs.

4. Media interchange

4.1 Implementation model

4.1.1 Application data flow diagram



The application is a single application that provides a user interface, network support and media support as a File Set Reader (FSR) and File Set Creator (FSC).

Conceptually it may be modeled as the following AEs:

- MEDIA-FSR, which imports a user-selected file or directory into the local database, either from the local file system or from a PS3.12 compliant media according to one of the General Purpose Media Application Profiles of PS3.11 (CD-R, DVD-RAM or BD). It is associated with the real-world activity "Import from Medium".
- MEDIA-FSC, which exports a user-selected set of patients, studies, series, either to the local file system or to a writeable storage medium according to one of the General Purpose Media Application Profiles of PS3.11 (CD-R, DVD-RAM or BD). It is associated with the real-world activity "Export to Medium".

4.1.2 Functional definition of AEs

4.1.2.1 Functional definition of "MEDIA-FSR"

MEDIA-FSR is activated through the user interface by scanning selected directories in the local file system and import compliant files into the local database. If a DICOMDIR file is found it will be used to import referenced instances.

4.1.2.2 Functional definition of "MEDIA-FSC"

MEDIA-FSC is activated through the user interface by exporting selected patients, studies or series to a medium or directory. The SOP Instances associated with the selection will be

collected and written to a single available media. Depending on the media type selected for export a corresponding profile is selected (CD-R, DVD-RAM or BD).

4.1.3 Sequencing of real world activities

4.1.3.1 Activity "import from medium"

All import activities are sequentially initiated in the user interface, and another activity may not be initiated until the prior activity has completed.

4.1.3.2 Activity "export to medium"

At least one patient, study or series must exist and be selected before the MEDIA-FSC Application Entity can be invoked. The operator can insert a new CD-R, DVD-RAM or BD media at any time before or after invocation of the MEDIA-FSC Application Entity. The MEDIA-FSC Application Entity will wait indefinitely for a media to be inserted before starting to write to the media device. If no media is available the export job can be canceled.

4.1.4 File meta information for implementation class and version

Table 4-1 DICOM Implementation Class and Version for Media Storage

Implementation class and version	
Implementation Class UID	1.2.826.0.1.3680043.2.360.0.3.6.8
Implementation Version Name	OFFIS_DCMTK_368

4.2 AE specifications

4.2.1 MEDIA-FSR application entity specification

The MEDIA-FSR Application Entity provides standard conformance to the Media Storage Service Class. The Application Profiles and roles are listed below:

Table 4-2 Application Profiles, Activities and Roles for MEDIA-FSR

Application profiles supported	Real-world activity	Role
STD-GEN-CD	Import from Medium	FSR
STD-GEN-DVD-RAM	Import from Medium	FSR
STD-GEN-BD	Import from Medium	FSR

4.2.1.1 File meta information for the application entity

Not applicable, since MEDIA-FSR is not an FSC or FSU.

4.2.1.2 Real-world activities

4.2.1.2.1 Activity – Import from medium

MEDIA-FSR is activated through the user interface when a user selects the “Filesystem – Scan folder” operation.

A browser will be displayed, from which studies or series may be selected and imported into the local database.

4.2.2 MEDIA-FSC application entity specification

The MEDIA-FSC Application Entity provides standard conformance to the Media Storage Service Class. The Application Profiles and roles are listed below:

Table 4-3 Application Profiles, Activities and Roles for MEDIA-FSC

Application profiles supported	Real-world activity	Role
STD-GEN-CD	Export to Medium	FSC
STD-GEN-DVD-RAM	Export to Medium	FSC
STD-GEN-BD	Export to Medium	FSC

4.2.2.1 File meta information for the application entity

No additional values pertain in the File Meta Information of the exported instances.

4.2.2.2 Real-world activities

4.2.2.2.1 Activity – Export to medium

MEDIA-FSC is activated through the user interface when a user selects the “Export – Export to medium” operation.

A dialog will be presented allowing the user to modify the suggested media configuration and provides control over the available media capacity. Additional content can be selected, e.g. if the iQ-LITE viewer application will be included or if JPEG versions of images are created additionally.

The dialog will wait in case a medium is not present. The content of the export job will be written together with a corresponding DICOMDIR to the medium in a single session. Writing in multi-session mode is not supported. The user can cancel an export job.

4.3 Augmented and private application profiles

The iQ-VIEW application does not support any augmented or private application profiles.

4.4 Media configuration

The iQ-VIEW application does not provide any media specific configuration.

5. Transformation of DICOM to CDA

The iQ-VIEW application does not support the export of CDA documents.

6. Support of character sets

6.1 Overview

The application supports all in table 7-1 listed Specific Character Sets, including single-byte and multi-byte character sets as well as code extension techniques using ISO 2022 escapes.

Support extends to correctly decoding and displaying the correct symbol for all names and strings found in the DICOMDIR, in storage instances from media and received over the network, and in the local database.

No specific support for sorting of strings other than in the default character set is provided in the browsers.

For the creation of DICOM instances only character sets without code extension are supported.

6.2 Character sets

In addition to the default character repertoire, the following Specific Character Sets and their Defined Terms are supported:

Table 6-1 Supported Specific Character Set Defined Terms

Character set description	Defined term
Single-byte character sets	
Latin alphabet No. 1	ISO_IR 100
Latin alphabet No. 2	ISO_IR 101
Latin alphabet No. 3	ISO_IR 109
Latin alphabet No. 4	ISO_IR 110
Cyrillic	ISO_IR 144
Arabic	ISO_IR 127
Greek	ISO_IR 126
Hebrew	ISO_IR 138
Latin alphabet No. 5	ISO_IR 148
Japanese	ISO_IR 13
Thai	ISO_IR 166
Single-byte character sets with code extensions*	
Default repertoire	ISO 2022 IR 6
Latin alphabet No. 1	ISO 2022 IR 100
Latin alphabet No. 2	ISO 2022 IR 101
Latin alphabet No. 3	ISO 2022 IR 109
Latin alphabet No. 4	ISO 2022 IR 110
Cyrillic	ISO 2022 IR 144
Arabic	ISO 2022 IR 127

Character set description	Defined term
Greek	ISO 2022 IR 126
Hebrew	ISO 2022 IR 138
Latin alphabet No. 5	ISO 2022 IR 148
Japanese	ISO 2022 IR 13
Thai	ISO 2022 IR 166
Multi-byte character sets with code extensions*	
Japanese	ISO 2022 IR 87
Japanese	ISO 2022 IR 159
Korean	ISO 2022 IR 149
Multi-byte character sets without code extensions	
Unicode in UTF-8	ISO_IR 192
GB18030	GB18030

* These Specific Character Sets are only supported for decoding and display.

6.3 Character set configuration

Whether or not characters are displayed correctly depends on the presence of font support in the underlying operating system.

7. Security

7.1 Security profiles

None supported.

7.2 Association level security

None supported.

Any Calling AE Titles and/or IP addresses may open an Association.

7.3 Application level security

None supported.

8. Annexes

8.1 IOD contents

8.1.1 Created SOP instance(s)

The iQ-VIEW application supports the creation of the following SOP Instances:

Table 8-1 Created SOP Instances

SOP Class Name	SOP Class UID	Table
Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7	Table 8-3
Grayscale Softcopy Presentation State Storage	1.2.840.10008.5.1.4.1.1.11.1	Table 8-4
Basic Text SR Storage	1.2.840.10008.5.1.4.1.1.88.11	Table 8-5
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Table 8-6

The following tables use a number of abbreviations. The abbreviations used in the "Presence of ..." column are:

Table 8-2 Abbreviations

Abbreviation	Meaning
"Presence of ..." column	
VNAP	Value Not Always Present (attribute sent zero length if no value is present)
ANAP	Attribute Not Always Present
ALWAYS	Always Present
EMPTY	Attribute is sent without a value
"Source" column	
MWL	Attribute value source is Modality Worklist
USER	Attribute value source is from User input
AUTO	Attribute value source is generated automatically
COPY	Attribute value source is from source instances
MPPS	Attribute value source is the same as that use for Modality Performed Procedure Step
CONFIG	Attribute value source is a configurable parameter

8.1.1.1 Secondary capture image storage

Table 8-3 IOD of created SC SOP Instances

IE	Module	Reference	Presence of module
Patient	Patient	Table 8-7	ALWAYS
Study	General Study	Table 8-8	ALWAYS

IE	Module	Reference	Presence of module
Series	General Series	Table 8-10	ALWAYS
Equipment	General Equipment	Table 8-11	ALWAYS
	SC Equipment	Table 8-13	ALWAYS
Image	General Image	Table 8-14	ALWAYS
	Image Pixel	Table 8-15	ALWAYS
	SC Image	Table 8-16	ALWAYS
	Modality LUT	Table 8-12	ALWAYS
	VOI LUT	Table 8-17	Only for 24 bit
	SOP Common	Table 8-18	ALWAYS

8.1.1.2 Grayscale softcopy presentation state storage

Table 8-4 IOD of created Presentation State SOP Instances

IE	Module	Reference	Presence of module
Patient	Patient	Table 8-7	ALWAYS
Study	General Study	Table 8-8	ALWAYS
Series	General Series	Table 8-10	ALWAYS
	Presentation Series	Table 8-19	ALWAYS
Equipment	General Equipment	Table 8-11	ALWAYS
Image	Presentation State Identification	Table 8-20	ALWAYS
	Presentation State Relationship	Table 8-21	ALWAYS
	Presentation State Shutter	Table 8-22	Only with shutter
	Display Shutter	Table 8-23	Only with shutter
	Displayed Area	Table 8-24	ALWAYS
	Graphic Annotation	Table 8-25	Only with annotation
	Spatial Transformation	Table 8-26	ALWAYS
	Graphic Layer	Table 8-27	Only with annotation
	Modality LUT	Table 8-12	ALWAYS
	Softcopy VOI LUT	Table 8-28	ALWAYS
	Softcopy Presentation LUT	Table 8-29	ALWAYS
	SOP Common	Table 8-30	ALWAYS

8.1.1.3 Basic text SR storage

Table 8-5 IOD of created Basic Text SR SOP Instances

IE	Module	Reference	Presence of module
Patient	Patient	Table 8-7	ALWAYS
Study	General Study	Table 8-8	ALWAYS
Series	SR Document Series	Table 8-31	ALWAYS

IE	Module	Reference	Presence of module
Equipment	General Equipment	Table 8-11	ALWAYS
Image	SR Document General	Table 8-32	ALWAYS
	SR Document Content	Table 8-33	ALWAYS
	SOP Common	Table 8-34	ALWAYS

8.1.1.4 Encapsulated PDF storage

Table 8-6 IOD of created Encapsulated PDF SOP Instances

IE	Module	Reference	Presence of module
Patient	Patient	Table 8-7	ALWAYS
Study	General Study	Table 8-8	ALWAYS
	Patient Study	Table 8-9	ALWAYS
Series	Encapsulated Document Series	Table 8-35	ALWAYS
Equipment	General Equipment	Table 8-11	ALWAYS
	SC Equipment	Table 8-36	ALWAYS
Image	Encapsulated Document	Table 8-37	ALWAYS
	SOP Common	Table 8-38	ALWAYS

8.1.1.5 Common modules

Table 8-7 Patient Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Patient's Name	(0010,0010)	PN	From Modality Worklist, user input or source instances. Values supplied via Modality Worklist will be entered as received. Values supplied via user input will contain all 5 components (some possibly empty). Maximum 64 characters per component.	VNAP	MWL/ USER/ COPY
Patient ID	(0010,0020)	LO	From Modality Worklist, user input or source instances. Maximum 64 characters.	VNAP	MWL/ USER/ COPY
Patient's Birth Date	(0010,0030)	DA	From Modality Worklist, user input or source instances.	VNAP	MWL/ USER/ COPY

Attribute name	Tag	VR	Value	Presence of value	Source
Patient's Sex	(0010,0040)	CS	From Modality Worklist, user input or source instances.	VNAP	MWL/ USER/ COPY

Table 8-8 General Study Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Study Instance UID	(0020,000D)	UI	From Modality Worklist, user input, source instances or generated by device.	ALWAYS	MWL/ USER/ COPY/ AUTO
Study Date	(0008,0020)	DA	From Modality Worklist, user input, source instances or generated by device. <yyyymmdd>	ALWAYS	MWL/ USER/ COPY/ AUTO
Study Time	(0008,0030)	TM	From Modality Worklist, user input, source instances or generated by device. <hhmmss>	ALWAYS	MWL/ USER/ COPY/ AUTO
Referring Physician's Name	(0008,0090)	PN	From Modality Worklist, user input or source instances.	VNAP	MWL/ USER/ COPY
Study ID	(0020,0010)	SH	From Modality Worklist or source instances.	VNAP	MWL/ COPY
Accession Number	(0008,0050)	SH	From Modality Worklist, user input or source instances.	VNAP	MWL/ USER/ COPY
Study Description	(0008,1030)	LO	From Modality Worklist, user input or source instances.	VNAP	MWL/ USER/ COPY

Table 8-9 Patient Study Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Patient's Age	(0010,1010)	AS	Calculated from Patient's Birth Date on base of actual Date	ALWAYS	AUTO

Table 8-10 General Series Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Series Instance UID	(0020,000E)	UI	Generated.	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated.	ALWAYS	AUTO
Series Description	(0008,103E)	LO	"Secondary Capture Sequence" or "IIS Presentation State".	ALWAYS	AUTO

Table 8-11 General Equipment Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Manufacturer	(0008,0070)	LO	IMAGE Information Systems	ALWAYS	AUTO
Institution Name	(0008,0080)	LO	From Configuration	VNAP	CONFIG

Table 8-12 Modality LUT Module of created SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Rescale Intercept	(0028,1052)	DS	From source image or 0.	ALWAYS	AUTO
Rescale Slope	(0028,1053)	DS	From source image or 1.	ALWAYS	AUTO
Rescale Type	(0028,1054)	LO	From source image or US.	ALWAYS	AUTO

8.1.1.6 Secondary capture image modules

Table 8-13 SC Equipment Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Conversion Type	(0008,0064)	CS	WSD	ALWAYS	AUTO
Modality	(0008,0060)	CS	OT	ALWAYS	AUTO

Table 8-14 General Image Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Instance Number	(0020,0013)	IS	Generated	ALWAYS	AUTO
Patient Orientation	(0020,0020)	CS	Zero length	EMPTY	AUTO
Image Type	(0008,0008)	CS	ORIGINAL\PRIMARY\OTHER	ALWAYS	AUTO

Table 8-15 Image Pixel Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Samples per Pixel	(0028,0002)	US	1 or 3. Images with burned-in annotations have always value 3.	ALWAYS	AUTO
Photometric Interpretation	(0028,0004)	CS	MONOCHROME1 (for Samples per Pixel = 1) or RGB (for Samples per Pixel = 3).	ALWAYS	AUTO
Rows	(0028,0010)	US	Generated	ALWAYS	AUTO
Columns	(0028,0011)	US	Generated	ALWAYS	AUTO
Bits Allocated	(0028,0100)	US	8	ALWAYS	AUTO
Bits Stored	(0028,0101)	US	8	ALWAYS	AUTO
High Bit	(0028,0102)	US	7	ALWAYS	AUTO
Pixel Representation	(0028,0103)	US	0000H	ALWAYS	AUTO
Planar Configuration	(0028,0006)	US	0	Only for Samples per Pixel Value = 3	AUTO
Pixel Data	(7FE0,0010)	OB	The Pixel Data itself may contain burned-in annotation.	ALWAYS	

Table 8-16 SC Image Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Date of Secondary Capture	(0018,1012)	DA	Generated	ANAP	AUTO
Time of Secondary Capture	(0018,1014)	TM	Generated	ANAP	AUTO

Table 8-17 VOI LUT Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Window Center	(0028,1050)	DS	127	ALWAYS	AUTO
Window Width	(0028,1051)	DS	256	ALWAYS	AUTO

Table 8-18 SOP Common Module of created SC SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Specific Character Set	(0008,0005)	CS	From source instances or configuration	ALWAYS	AUTO/CONFIG
SOP Class UID	(0008,0016)	UI	1.2.840.10008.5.1.4.1.1.7	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated	ALWAYS	AUTO

8.1.1.7 Grayscale softcopy presentation state modules

Table 8-19 Presentation Series Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Modality	(0008,0060)	CS	PR	ALWAYS	AUTO

Table 8-20 Presentation State Identification Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Presentation Creation Date	(0070,0082)	DA	Generated by device	ALWAYS	AUTO
Presentation Creation Time	(0070,0083)	TM	Generated by device	ALWAYS	AUTO

Table 8-21 Presentation State Relationship Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Referenced Series Sequence	(0008,1115)	SQ	One or more items	ALWAYS	AUTO
>Series Instance UID	(0020,000E)	UI	From referenced image	ALWAYS	AUTO
> Referenced Image Sequence	(0008,1140)	SQ	From referenced image	ALWAYS	AUTO
>>Referenced SOP Class UID	(0008,1150)	UI	From referenced image	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1155)	UI	From referenced image	ALWAYS	AUTO
>>Referenced Frame Number	(0008,1160)	IS	Frame number if referenced image is a multi-frame image, or 1 for single-frame image	ALWAYS	AUTO
>> Private Creator	(0009,0010)	LO	IIS_IQ-VIEW	ALWAYS	AUTO

Attribute name	Tag	VR	Value	Presence of value	Source
>> Private Referenced Filename	(0009,1010)	UT	Filename of referenced file in local database	ALWAYS	AUTO

Table 8-22 Presentation State Shutter Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Shutter Presentation Value	(0018,1622)	US	0 if shutter is present.	VNAP	AUTO

Table 8-23 Display Shutter Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Shutter Shape	(0018,1600)	CS	RECTANGULAR if shutter is applied.	VNAP	AUTO
Shutter Left Vertical Edge	(0018,1602)	IS	If shutter is applied.	VNAP	AUTO
Shutter Right Vertical Edge	(0018,1604)	IS	If shutter is applied.	VNAP	AUTO
Shutter Upper Horizontal Edge	(0018,1606)	IS	If shutter is applied.	VNAP	AUTO
Shutter Lower Horizontal Edge	(0018,1608)	IS	If shutter is applied.	VNAP	AUTO

Table 8-24 Displayed Area Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Displayed Area Selection Sequence	(0070,005A)	SQ	One or more items	ALWAYS	AUTO
> Displayed Area Top Left Hand Corner	(0070,0052)	SL	From current display setting	ALWAYS	AUTO
> Displayed Area Bottom Right Hand Corner	(0070,0053)	SL	From current display setting	ALWAYS	AUTO
> Presentation Size Mode	(0070,0100)	CS	From current display setting. SCALE TO FIT/ TRUE SIZE/ MAGNIFY	ALWAYS	AUTO

Attribute name	Tag	VR	Value	Presence of value	Source
>Presentation Pixel Spacing	(0070,0101)	DS	From current display setting.	ANAP	AUTO
> Presentation Pixel Aspect Ratio	(0070,0102)	IS	1\1	ANAP	AUTO
> Presentation Pixel Magnification Ratio	(0070,0103)	FL	From current display setting.	ANAP	AUTO

Table 8-25 Graphic Annotation Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Graphic Annotation Sequence	(0070,0001)	SQ	One or more items	ANAP	AUTO
> Referenced Image Sequence	(0008,1140)	SQ	One or more items	ALWAYS	AUTO
>> Referenced SOP Class UID	(0008,1155)	UI	From referenced image	ALWAYS	AUTO
>> Referenced SOP Instance UID	(0008,1160)	UI	From referenced image	ALWAYS	AUTO
>>Referenced Frame Number	(0008,1160)	IS	Frame number if referenced image is a multi-frame image, or 1 for single-frame image	ALWAYS	AUTO
> Graphic Layer	(0070,0002)	CS	LAYER	ALWAYS	AUTO
> Text Object Sequence	(0070,0008)	SQ	One or more items if text annotation present.	ANAP	AUTO
>> Specific Character Set	(0008,0005)	CS	ISO_IR 192	ALWAYS	AUTO
>> Bounding Box Annotation Units	(0070,0003)	CS	PIXEL	ANAP	AUTO
>> Anchor Point Annotation Units	(0070,0004)	CS	PIXEL	ANAP	AUTO
>> Unformatted Text Value	(0070,0006)	ST	From user input.	ALWAYS	USER
>> Bounding Box Top Left Hand Corner	(0070,0010)	FL	From user input.	ANAP	USER
>> Bounding Box Bottom Right Hand Corner	(0070,0011)	FL	From user input.	ANAP	USER
>> Bounding Box Text Horizontal Justification	(0070,0012)	CS	LEFT	ANAP	AUTO
>> Anchor Point	(0070,0014)	FL	From user input.	ANAP	USER

Attribute name	Tag	VR	Value	Presence of value	Source
>> Anchor Point Visibility	(0070,0015)	CS	From user input.	ANAP	USER
> Graphic Object Sequence	(0070,0009)	SQ	One or more items if graphic annotation present	ANAP	AUTO
>> Graphic Annotation Units	(0070,0005)	CS	PIXEL	ALWAYS	AUTO
>> Graphic Dimensions	(0070,0020)	US	2	ALWAYS	AUTO
>> Number of Graphic Points	(0070,0021)	US	From user input.	ALWAYS	USER
>> Graphic Data	(0070,0022)	FL	From user input.	ALWAYS	USER
>> Graphic Type	(0070,0023)	CS	POLYLINE or ELLIPSE	ALWAYS	USER
>> Graphic Filled	(0070,0024)	CS	N	ALWAYS	AUTO
>> Private Creator	(0071,0010)	LO	IIS_IQ-VIEW	ALWAYS	AUTO
>> Private Measurement Descriptor	(0071,0110)	UT	Filename of referenced file in local database	ALWAYS	AUTO

Table 8-26 Spatial Transformation Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Image Rotation	(0070,0042)	US	From current display setting	ALWAYS	AUTO
Image Horizontal Flip	(0070,0041)	CS	From current display setting	ALWAYS	AUTO

Table 8-27 Graphic Layer Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Graphic Layer Sequence	(0070,0060)	SQ	One or more items	ANAP	AUTO
> Graphic Layer	(0070,0002)	CS	LAYER	ALWAYS	AUTO
> Graphic Layer Order	(0070,0062)	IS	1	ALWAYS	AUTO
> Graphic Layer Recommended Display Grayscale Value	(0070,0066)	US	56535 (0xFFFF - white)	ALWAYS	AUTO

Table 8-28 Softcopy VOI LUT Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Softcopy VOI LUT Sequence	(0028,3110)	SQ	One or more items	ALWAYS	AUTO
> Window Center	(0028,1050)	DS	From current display setting	ALWAYS	AUTO
> Window Width	(0028,1051)	DS	From current display setting	ALWAYS	AUTO

Table 8-29 Softcopy Presentation LUT Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Presentation LUT Shape	(2050,0020)	CS	IDENTITY or INVERSE	ALWAYS	AUTO

Table 8-30 SOP Common Module of created GSPS SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Specific Character Set	(0008,0005)	CS	ISO_IR 100	ALWAYS	AUTO
SOP Class UID	(0008,0016)	UI	1.2.840.10008.5.1.4.1.1.11.1	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated	ALWAYS	AUTO

8.1.1.8 Basic text structured report modules

Table 8-31 SR Document Series Module of created SR SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Modality	(0008,0060)	CS	SR	ALWAYS	AUTO
Series Description	(0008,103E)	LO	IIS Structured Report	ALWAYS	AUTO
Referenced Performed Procedure Step Sequence	(0008,1111)	SQ	Sequence without item	EMPTY	AUTO
Series Instance UID	(0020,000E)	UI	Generated	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated	ALWAYS	AUTO

Table 8-32 SR Document General Module of created SR SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Content Date	(0008,0023)	DA	Generated	ALWAYS	AUTO
Content Time	(0008,0033)	TM	Generated	ALWAYS	AUTO
Referring Physician's Address	(0008,0092)	ST	From SR source study	VNAP	AUTO
Referring Physician's Telephone Numbers	(0008,0094)	SH	From SR source study	VNAP	AUTO
Patient's Address	(0010,1040)	LO	From SR source study	VNAP	AUTO
Admitting Diagnoses Description	(0008,1080)	LO	From SR source study	VNAP	AUTO
Patient's Telephone Number	(0010,2154)	SH	From SR source study	VNAP	AUTO
Instance Number	(0020,0013)	IS	1	ALWAYS	AUTO
Images in Acquisiton	(0020,1002)	IS	Generated	ALWAYS	AUTO
Requesting Physician	(0032,1032)	PN	From SR source study	VNAP	AUTO
Completion Flag	(0040,A491)	CS	PARTIAL or COMPLETE From user input.	ALWAYS	USER
Verification Flag	(0040,A493)	CS	UNVERIFIED or VERIFIED From user input.	ALWAYS	USER
Verifying Observer Sequence	(0040,A073)	SQ	From user input. If Veri- fication Flag is VERIFIED.	ANAP	AUTO
> Verifying Organization	(0040,A027)	LO	From user input	VNAP	USER
> Verification DateTime	(0040,A030)	DT	Generated	ALWAYS	AUTO
> Verifying Observer Name	(0040,A075)	PN	From user input	VNAP	USER
> Verifying Observer Identification Code Sequence	(0040,A088)	SQ	Sequence without item	EMPTY	AUTO
Performed Procedure Code Sequence	(0040,A372)	SQ	Sequence without item	EMPTY	AUTO

Table 8-33 SR Document Content Module of created SR SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Value Type	(0040,A040)	CS	CONTAINER	ALWAYS	AUTO
Concept Name Code Sequence	(0040,A043)	SQ	Contains single item	ALWAYS	AUTO
> Code Value	(0008,0100)	SH	1111	ALWAYS	AUTO
> Coding Scheme Designator	(0008,0102)	SH	Unknown	ALWAYS	AUTO
> Code Meaning	(0008,0104)	LO	Report	ALWAYS	AUTO
Continuity of Content	(0040,A050)	CS	SEPARATE	ALWAYS	AUTO
Content Sequence	(0040,A730)	SQ	One or more items	ALWAYS	AUTO
> Relationship Type	(0040,A010)	CS	Either HAS OBS CONTEXT or CONTAINS	ALWAYS	AUTO
> Value Type	(0040,A040)	CS	One of COMPOSITE, TEXT or PNAME	ALWAYS	AUTO
> Concept Name Code Sequence	(0040,A043)	SQ	If Value Type is COMPOSITE	ANAP	AUTO
> Text Value	(0040,A160)	UT	If Value Type is TEXT	ANAP	AUTO
> Person Name	(0040,A123)	PN	If Value Type is PNAME	ANAP	USER
Private Creator	(0071,0010)	LO	IIS_IQ-VIEW	ALWAYS	AUTO
IIS Extended SR Private Creator	(0077,0010)	LO	IIS_IQ-VIEW	ALWAYS	AUTO
Private Sequence	(0077,1010)	SQ	One item per series in SR source study	VNAP	AUTO
> Modality	(0008,0060)	CS	From source series	VNAP	AUTO
> Manufacturer	(0008,0070)	LO	From source series	VNAP	AUTO
> Institution Name	(0008,0080)	LO	From source series	VNAP	AUTO
> Series Description	(0008,103E)	LO	From source series	VNAP	AUTO
> Institutional Department Name	(0008,1040)	LO	From source series	VNAP	AUTO
> Performing Physician's Name	(0008,1050)	PN	From source series	VNAP	AUTO
> Manufacturer's Model Name	(0008,1090)	LO	From source series	VNAP	AUTO
> Series Instance UID	(0020,000E)	UI	From source series	VNAP	AUTO

Table 8-34 SOP Common Module of created SR SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Specific Character Set	(0008,0005)	CS	From configuration	ALWAYS	CONFIG
SOP Class UID	(0008,0016)	UI	1.2.840.10008.5.1.4.1.1.88.11	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated	ALWAYS	AUTO

8.1.1.9 Encapsulated PDF modules

Table 8-35 Encapsulated Document Series Module of created PDF SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Modality	(0008,0060)	CS	DOC	ALWAYS	AUTO
Series Instance UID	(0020,000E)	UI	Generated	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated	ALWAYS	AUTO
Series Description	(0008,103E)	LO	IIS PDF Storage	ALWAYS	AUTO

Table 8-36 SC Equipment Module of created PDF SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Conversion Type	(0008,0064)	CS	WSD	ALWAYS	AUTO

Table 8-37 Encapsulated Document Module of created PDF SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Instance Number	(0020,0013)	IS	Generated	ALWAYS	AUTO
Burned In Annotation	(0028,0301)	CS	YES or NO	ALWAYS	AUTO
Concept Name Code Sequence	(0040,A043)	SQ	Sequence without item	EMPTY	AUTO
Performed Procedure Code Sequence	(0040,A372)	SQ	Sequence without item	EMPTY	AUTO
Document Title	(0042,0010)	ST	IIS PDF IMPORT	ALWAYS	AUTO
Encapsulated Document	(0042,0011)	OB	Document data.	ALWAYS	USER

Table 8-38 SOP Common Module of created PDF SOP Instances

Attribute name	Tag	VR	Value	Presence of value	Source
Specific Character Set	(0008,0005)	CS	From configuration	ALWAYS	CONFIG
SOP Class UID	(0008,0016)	UI	1.2.840.10008.5.1.4.1.1.104.1	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated	ALWAYS	AUTO

8.1.2 Usage of attributes from received IODs

No SOP Class specific fields for images are required.

8.1.3 Attribute mapping

Not applicable.

8.1.4 Coerced / modified fields

No coercion is performed.

8.2 Data dictionary of private attributes

Table 8-39 Data Dictionary of Private Attributes

Attribute name	Tag	VR	VM	Description
Private Creator	(0009,0010)	US	1	Private Creator
Private Filter Sharpen	(0009,0100)	US	1	Presentation State Filter
Private Filter Soften	(0009,0101)	US	1	Presentation State Filter
Private Filter Edge Enhance	(0009,0102)	US	1	Presentation State Filter
Private Filter Blur	(0009,0103)	US	1	Presentation State Filter
Private Filter Low Pass	(0009,0104)	US	1	Presentation State Filter
Private Filter High Pass	(0009,0105)	US	1	Presentation State Filter
Private Color Scheme	(0009,0200)	US	1	Color Scheme
Private Referenced Filename	(0009,1010)	UT	1	Presentation State Referenced Image Filename
Private PR Sequence	(0009,3000)	SQ		Presentation State Sequence
Version Name	(0071,0010)	LO	1	iQ-VIEW Version Name
Version Number	(0071,0020)	LO	1	iQ-VIEW Version Number
Private PR Filename	(0071,0030)	UT	1	Presentation Filename
No Series Scope	(0071,0090)	US	1	No Series Scope

Attribute name	Tag	VR	VM	Description
Measurement Descriptor	(0071,0110)	LO	1	Presentation State Measurement
Private Measurement Data	(0071,0120)	UT	1	Presentation State Measurement Data
Specific Scope	(0071,0130)	US	1	Presentation State Specific Scope

8.3 Coded terminology and templates

The iQ-VIEW application is not using any Codes (SNOMED) or Controlled Terminology, such as the use of the DICOM Content Mapping Resource (DCMR).

8.4 Grayscale image consistency

The high resolution display monitor attached to the product can be calibrated according to the Grayscale Standard Display Function (GSDF). External tools can be used together with a luminance meter to measure the Characteristic Curve of the display system and the current ambient light. The result of the calibration procedure is an ICC profile that will be applied to systems graphic subsystem.

8.5 Standard extended / specialized / private SOP classes

None.

8.6 Private transfer syntaxes

None.

 **IMAGE Information Systems Europe GmbH**

Lange Str. 16 | 18055 Rostock | Germany
+49 381 496 58 20 | www.iq-image.com