Medtronic

Endurant™ II/IIs

Stent Graft System

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Contents

1 Device description	15
1.1 Stent graft	15
1.2 Delivery system	18
2 Intended use	20
2.1 Intended users	20
2.2 Indications for use	20
2.3 Contraindications	21
2.4 Clinical benefits	21
2.5 Intended patient population	21
2.6 Device performance characteristics	21
3 Warnings and precautions	22
3.1 General	22
3.2 Patient selection	22
3.3 Before the implant procedure	24
3.4 During the implant procedure	24
3.5 Treatment and follow-up	25
4 Adverse events	26
4.1 Potential adverse events	26
5 Patient selection and treatment	27
5.1 Individualization of treatment	27
6 Patient counseling information	27
7 How supplied	28
7.1 Sterility	28
7.2 Contents	28
8 Device disposal	28
9 Clinical use information	28
9.1 Physician training requirements	28
9.2 Recommended device sizing	29
9.3 Device inspection	30
9.4 Required materials (not included in the stent graft system packaging)	30
9.5 Recommended materials (not included in the stent graft system packaging)	30
9.6 MRI safety information	31
10 Implant instructions	31
10.1 Vascular access and device preparation	31
10.2 Delivery procedure	32
11 Bail-out techniques	47
11.1 Screw gear handle disassembly	47
11.2 Ballooning	48
11.3 Back-end handle disassembly	48
11.4 Snare the tapered tip	48
12 Follow-up imaging recommendations	48
12.1 General	48
12.2 X-ray	49
12.3 CT with contrast	49
12.4 Noncontrast CT	49
12.5 Duplex ultrasound	49
12.6 MRI or MRA	49
12.7 Imaging tests	49

12.8 Supplemental imaging 49

13 Additional surveillance and treatment 50

14 Disclaimer of warranty 50

1 Device description

The Endurant™ II/Endurant™ IIs stent graft system (hereinafter referred to as the Endurant II/IIs stent graft system) is designed for the endovascular repair of aneurysms. When placed within the target lesion, the stent graft provides an alternative conduit for blood flow within the patient's vasculature by excluding the lesion from blood flow and pressure.

The stent graft system is comprised of 2 main components: the implantable stent graft and the disposable delivery system. The stent graft is preloaded into the delivery system and advanced to the aneurysm using fluoroscopic guidance. Upon deployment, the stent graft self-expands to conform to the shape and size of the seal zones above and below the aneurysm.

The Endurant II/IIs stent graft can also be used with the Heli-FX EndoAnchor system (available separately). The Heli-FX EndoAnchor system is designed to provide fixation and augment sealing between the Endurant II/IIs stent graft and the native artery. The system consists of an EndoAnchor implant that is delivered using the Heli-FX applier through the steerable Heli-FX guide.

No accessories are provided with this device. Please see *Section 9.4* for additional required equipment and *Section 9.4* for additional recommended equipment.


1.1 Stent graft

The Endurant II/IIs stent graft (*Figure 1*) has 2 basic configurations: a bifurcated configuration and a limb configuration. Additional configurations include iliac extension, aortic extension, abdominal tube, and aorto-uni-iliac (AUI). After placement of the bifurcated or AUI device, limbs and additional stent grafts are introduced separately into the vessel and mated with the implanted configuration.

All Endurant II/IIs stent graft configurations are composed of nitinol stents sewn to a fabric graft with nonresorbable sutures. Radiopaque markers are sewn onto the stent graft to aid in visualization and to facilitate accurate placement. The nitinol stents may also be visible under fluoroscopy.

Stent graft components should be oversized to be larger than the measured vessel inner diameter (*Section 9.2*). *Table 1* contains a summary of the stent graft materials.

Figure 1. Stent graft configurations and locations of RO markers

- | | |
|---|---|
| 1 Radiopaque marker | 6 Endurant IIs bifurcated configuration |
| 2 'e' marker | 7 Endurant II iliac extension configuration |
| 3 Radiopaque gate marker | 8 Endurant II limb configuration |
| 4 Endurant II aortic extension/abdominal tube configuration | 9 Endurant II aorto-uni-iliac configuration |
| 5 Endurant II bifurcated configuration | 10 Overlap marker |

Note: This and all other product graphics appearing in this manual are not drawn to scale.

Table 1. Stent graft materials

Component	Material	Endurant II bifur (g)	Endurant IIs bifur (g)	AUI (g)	Aortic extension (g)	Limb (g)	Iliac extension (g)
Stents	Nickel-titanium (nitinol) alloy	1.64	1.98	1.26	1.13	1.73	0.84
Button radiopaque markers	Platinum-iridium alloy	0.17	0.15	0.12	0.10	0.10	0.10
Contralateral gate marker							
"e" radiopaque marker	Platinum	0.01	0.01	0.01	0.01	N/A	N/A
Graft material	Polyester	1.37	1.83	0.94	0.93	1.27	0.77
Suture	Polyethylene	0.02	0.02	0.02	0.02	0.01	0.01
Suture	Silicone-coated polyester	0.41	0.54	0.28	0.28	0.36	0.20

The Endurant II/IIs stent graft system does not contain natural rubber latex. However, during the manufacturing/assembly process, it may have incidental contact with latex-containing products.

No substances have been identified in this stent graft that are known to be carcinogenic, mutagenic, or toxic for reproduction (CMR). See *Section 1.2* for information about a CMR substance (cobalt) in the delivery system.

1.1.1 Bifurcated configuration

The bifurcated stent graft is available in 2 configurations: the Endurant II bifurcated configuration and the Endurant IIs bifurcated configuration. The Endurant II bifurcated configuration is an aortic bifurcated configuration that is available in 3 lengths. The Endurant IIs bifurcated configuration is an aortic configuration available in a single, shorter length (*Figure 1*). The proximal end of both bifurcated configurations deploy into the proximal neck and upper section of the aneurysm. The proximal end of the bifurcated configuration is composed of nitinol stents sewn to a fabric graft. The suprarenal portion of the proximal end is not covered with graft fabric (*Figure 1*). The suprarenal stent also contains anchor pins to fix the stent graft in place inside the aorta.

The aortic section bifurcates distally into 2 smaller tubes: an ipsilateral leg and a shorter contralateral leg. In the Endurant II bifurcated configuration, all stents on the ipsilateral leg are sewn to the outside of the fabric creating a smooth inner lumen. In the Endurant IIs bifurcated configuration, the 4 distal stents are sewn to the inside of the ipsilateral leg graft fabric. For all sizes, the stents on the contralateral leg are sewn to the inside of the graft fabric (*Figure 1*).

1.1.2 Limb configuration

The proximal end of the limb configuration deploys within the legs of the bifurcated configuration, while the distal end deploys into the iliac artery. The proximal end of the limb configuration has an open web configuration (*Figure 1*), which contains no graft material in its stent valleys.

Note: A limb device is implanted on both the ipsilateral and contralateral legs of an Endurant IIs bifurcated configuration. See *Section 10.2.12*.

1.1.3 Iliac extension configuration

An iliac extension configuration is available if an additional distal stent graft length is needed. It has an open web configuration on its proximal end (*Figure 1*).

Note: An appropriately sized limb configuration can be used as an iliac extension configuration.

1.1.4 Aortic extension and abdominal tube configurations

Aortic extension and abdominal tube configurations are available if additional proximal stent graft length is needed. The aortic extension and abdominal tube stent grafts have a bare proximal suprarenal stent with anchor pins (*Figure 1*).

1.1.5 Aorto-uni-iliac (AUI) configuration

The proximal end of the aorto-uni-iliac (AUI) configuration is deployed into the proximal neck and upper section of the aneurysm. All stents on the proximal aortic end of the AUI configuration are sewn to the outside of the graft fabric. The proximal stent (suprarenal) of the aortic section is not covered with fabric. As such, this bare stent design allows the AUI stent graft to be fixed above the renal arteries without obstructing them with graft fabric. Refer to *Figure 1* for a pictorial representation of the proximal configuration. The suprarenal stent includes anchor pins to help fix the AUI device in place. The suprarenal stent is joined to the proximal edge of the graft by ultra high molecular weight polyethylene suture.

Distally, the aortic section tapers down to a smaller diameter tube. In the distal end of the tapered AUI device, the stents are sewn to the inside of the graft fabric.

Note: The Talent occluder system (available separately) can be used to assist in stent graft implantation.

1.1.6 Balloon expandable covered stent configuration

The Endurant II/IIIs stent graft can also be used in a parallel graft technique in combination with 6 mm to 7 mm diameter balloon-expandable covered stents that are approved for use in the renal arteries (available separately). When used, these devices are deployed adjacent to the proximal body of the bifurcated or AUI stent graft configurations to perfuse the renal arteries. They are used at the discretion of the physician, per the instructions detailed in the corresponding Instructions for Use (IFU). For additional information on using balloon expandable covered stents in conjunction with the Endurant II/IIIs stent graft system, see *Section 10.2.9*. For additional information about the balloon expandable covered stent, refer to the IFU provided with the device.

Note: The clinical data used to evaluate the Endurant II/IIIs stent graft system in a parallel graft technique consisted primarily of cases in which the device was used with the Atrium Advanta™ V12 balloon expandable covered stent.

1.1.7 Heli-FX EndoAnchor system

The Endurant II/IIIs stent graft can also be used with the Heli-FX EndoAnchor system (available separately). The system consists of an EndoAnchor implant that is delivered using the Heli-FX applier through the steerable Heli-FX guide. The EndoAnchor implant is designed to provide fixation and augment sealing between the Endurant II/IIIs stent graft and the native artery. The EndoAnchor implant may be implanted at the time of the initial stent graft implantation or during a secondary (for example, repair) procedure.

Note: For infrarenal EVAR procedures using the Heli-FX EndoAnchor system, the EndoAnchor implant should be implanted only into areas of aortic tissue that are free of or have insignificant calcified plaque or thrombus, or where such pathology is diffuse and less than 2 mm in thickness. It should not be used in patients with bleeding diathesis or with a known hypersensitivity or allergy to the device materials.

For additional information about using this system with the Endurant II/IIIs stent graft, refer to *Section 10.2.18*. For additional information about the Heli-FX EndoAnchor system, refer to the instructions for use provided with the system.

1.1.8 Talent occluder system

The Talent occluder system is available to be used with Endurant II/IIIs stent graft system. The Talent occluder system is typically used in conjunction with the aorto-uni-iliac component and a femoral-femoral bypass in the event that:

- It is not possible to place a bifurcated component due to patient's anatomy, tortuous anatomy, a significant extension of aneurysmal disease, inability to access the contralateral iliac.
- It is needed to stop retrograde blood flow into the aneurysm sac.
- It is needed for repair.


The Talent occluder system is a tubular stent graft segment sealed at both ends. The Talent occluder system ensures that the blood flowing through the contralateral iliac artery does not “backflow” into the aneurysm sac, but instead flows from the femoral-femoral bypass into the patient's contralateral leg artery. The Talent occluder delivery system consists of a dilator used for delivery system insertion, a cartridge preloaded with the Talent occluder, a sheath for delivering the stent graft to the target location, and a pushrod to advance and help deploy the stent graft. Additional details are available in the Talent occluder system instructions for use.

1.2 Delivery system

The Endurant II delivery system that delivers all stent graft configurations consists of a single-use, disposable catheter with an integrated handle to provide controlled deployment. It is available in 14, 16, 18, and 20 Fr graft cover diameters with a working length of 57 cm ± 2 cm. The catheter assembly is flexible and compatible with a 0.035 in (0.89 mm) guidewire. There are 2 types of delivery systems: the aortic (*Figure 2*) and the iliac (*Figure 3*) delivery systems. The aortic delivery system delivers the bifurcated components, aortic extension, AUI, and abdominal tube stent graft configurations. The iliac delivery system delivers the limb and iliac extension stent graft configurations. The aortic delivery system features a tip capture mechanism, which is not present in the iliac delivery system.

Note: This delivery system contains cobalt (CAS No. 7440-48-4), present at a concentration >0.1% weight/weight. Cobalt is classified as an endocrine disruptor (ED) or as a carcinogen, mutagen, or substance toxic for reproduction (CMR) of Category 1A or 1B.


Figure 2. Aortic delivery system

- 1 Rear handle
- 2 Back-end wheel
- 3 Screw gear
- 4 External slider
- 5 Trigger
- 6 Front grip

- 7 Graft cover
- 8 Markerband
- 9 Spindle
- 10 Sleeve
- 11 Tapered tip

Figure 3. Iliac delivery system

- 1 Rear handle
- 2 Screw gear
- 3 External slider
- 4 Trigger

- 5 Front grip
- 6 Graft cover
- 7 Markerband
- 8 Tapered tip

2 Intended use

The intended use of the Endurant II/IIIs stent graft system is to repair infrarenal abdominal aortic or aortoiliac aneurysms, or juxtarenal abdominal aortic or aortoiliac aneurysms (with a parallel graft technique), in an endovascular approach by providing an alternative conduit for blood flow within the patient's vessel (aorta and iliac arteries) by excluding the aneurysm (target lesion).

2.1 Intended users

The Endurant II/IIIs stent graft system, including when used in conjunction with the Heli-FX EndoAnchor system for short neck or with a balloon-expandable covered stent in a parallel graft technique, should only be used by physicians and teams trained in vascular interventional techniques, including training in the use of the device. Specific training expectations are described in *Section 9.1*.

This device should only be used in a clinical setting in a sterile environment.

2.2 Indications for use

The Endurant II/IIIs stent graft system is indicated for the endovascular treatment of infrarenal or juxtarenal abdominal aortic or aortoiliac aneurysms in patients with the following characteristics:

- Adequate iliac or femoral access vessel morphology that is compatible with vascular access techniques, devices, or accessories
- Iliac distal fixation length of ≥ 15 mm
- Iliac diameters with a range of 8 to 25 mm
- Morphology suitable for aneurysm repair
- One of the following:
 - Aneurysm diameter > 5 cm
 - Aneurysm diameter of 4 to 5 cm, which has also increased in size by 0.5 cm in the last 6 months
 - Aneurysm that is at least 1.5 times the diameter of the normal infrarenal aorta

In addition, for treatment of infrarenal abdominal aortic or aortoiliac aneurysms, the following patient characteristics apply:

- Aortic neck diameters with a range of 19 to 32 mm

- Proximal neck length of ≥ 10 mm, or ≥ 4 mm and < 10 mm when used in conjunction with the Heli-FX EndoAnchor system, with insignificant calcification, or insignificant thrombus with $\leq 60^\circ$ infrarenal and $\leq 45^\circ$ suprarenal angulation and a vessel diameter approximately 10% to 20% smaller than the labeled Endurant II/IIIs stent graft diameter
- Proximal neck length of ≥ 15 mm with insignificant calcification, or insignificant thrombus with $\leq 75^\circ$ infrarenal and $\leq 60^\circ$ suprarenal angulation and a vessel diameter approximately 10% to 20% smaller than the labeled Endurant II/IIIs stent graft diameter

In addition, for treatment of juxtarenal abdominal aortic or aortoiliac aneurysms with the parallel graft technique, the following patient characteristics apply:

- Aortic neck diameters with a range of 19 to 30 mm
- Infrarenal proximal neck length of ≥ 2 mm and total available proximal seal length > 15 mm with a balloon-expandable covered stent deployed in a parallel graft technique in one or two renal arteries and with insignificant calcification, or insignificant thrombus with $\leq 60^\circ$ infrarenal, $\leq 45^\circ$ suprarenal, and $\leq 45^\circ$ suprasuperior mesenteric artery angulation (supraSMA angulation) and an aortic diameter approximately 20% to 30% smaller than the labeled Endurant II/IIIs stent graft diameter
- Adequate brachial or axillary access that is compatible with vascular access techniques, devices, or accessories
- Sufficient renal fixation length
- Renal take-off angulation $< 90^\circ$ relative to the aortic centerline

Note: The angle between the aortic and renal artery blood lumen centerlines is the renal take-off angulation. An angle of $< 90^\circ$ is associated with a downward-pointed renal artery relative to the trajectory of the aortic centerline.

2.3 Contraindications

The Endurant II/IIIs stent graft system is contraindicated for patients with the following conditions:

- Patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation procedures (*Chapter 9 to Chapter 12*).
- Patients who have a condition that threatens to infect the graft
- Patients who are sensitive to or have allergies to the device materials listed in *Section 3.2*

When used with the Heli-FX EndoAnchor system, the Endurant II/IIIs stent graft system is also contraindicated in patients with known sensitivities to the EndoAnchor implant materials.

Also consider patient selection information (*Section 3.2*). For contraindications regarding ancillary devices used with the Endurant II/IIIs stent graft system, refer to the instructions for use provided with each device.

2.4 Clinical benefits

Note: After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00001977U.

The Endurant II/IIIs stent graft system is intended to repair infrarenal abdominal aortic or aortoiliac aneurysms, or juxtarenal abdominal aortic or aortoiliac aneurysms (with a parallel graft technique). This minimally invasive endovascular approach provides an alternative conduit for blood flow within the patient's vessel (aorta and iliac arteries) by excluding the aneurysm (target lesion), which should prevent rupture and aneurysm related mortality (ARM). Therefore, the primary clinical benefits include prevention of rupture and ARM.

2.5 Intended patient population

The Endurant II/IIIs stent graft system is intended to be used in patients with infrarenal or juxtarenal abdominal aortic or aortoiliac aneurysm with anatomic characteristics and vessel morphology compatible with the device. For detailed patient characteristics, please refer to *Section 2.2*.

2.6 Device performance characteristics

The Endurant II/IIIs stent graft system consists of various stent graft configurations and sizes as well as a catheter delivery system. The stent graft configurations can be implanted in combination into the patient's abdominal aorta and iliac arteries. When implanted, the system provides an alternative conduit for blood flow in the patient's aorta and iliac arteries. The Endurant II/IIIs stent graft system provides high flexibility and conformability to adapt to a straight or tortuous anatomy. It combines accurate and controlled positioning with simple deployment, providing adequate stent graft placement in challenging aortic necks.

The Endurant II/IIIs stent graft can also be used with the Heli-FX EndoAnchor system (available separately). The system consists of an EndoAnchor implant that is delivered using the Heli-FX applier through the steerable Heli-FX guide. The EndoAnchor implant is designed to provide fixation and to augment sealing between the Endurant II/IIIs stent graft and the native artery. The EndoAnchor implant may be implanted at the time of the initial stent graft implantation or during a secondary (for example, repair) procedure.

Note: For additional information about using Heli-FX EndoAnchor system with the Endurant II/IIIs stent graft, refer to *Section 10.2.18*. For further information about the Heli-FX EndoAnchor system, consult the instructions for use provided with the system.

The Endurant II/IIIs stent graft can also be used in a parallel graft technique in combination with 6 mm to 7 mm diameter balloon-expandable covered stents that are approved for use in the renal arteries (available separately).

Note: For additional information on using balloon-expandable covered stents in conjunction with the Endurant II/IIIs stent graft system, see *Section 10.2.9*. For further information about the balloon expandable covered stent, consult the instructions for use provided with the device.

There is robust clinical data to support the safety and effectiveness of the Endurant II/IIIs stent graft system, deriving from both Medtronic sponsored studies and/or from published literature. Data demonstrated that the Endurant II/IIIs stent graft system, when used in accordance with the instructions for use, is safe and effective for the endovascular treatment of infrarenal abdominal aortic and aortoiliac aneurysms. This includes the treatment of short necks (when used in conjunction with the Heli-FX EndoAnchor system) and juxtarenal AAA (utilizing parallel graft technique).

3 Warnings and precautions

Caution: Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

3.1 General

- The Endurant II/IIIs stent graft system should only be used by physicians and teams trained in vascular interventional techniques, including training in the use of the device. Specific training expectations are described in *Section 9.1*.
- When the Endurant II/IIIs stent graft system is used with the Heli-FX EndoAnchor system, it should only be used by physicians and teams who are also trained in the use of the EndoAnchor system.
- The Heli-FX EndoAnchor system should be used when augmented radial fixation and/or sealing is required; in particular, in the treatment of abdominal aortic aneurysms with short (≥ 4 mm and < 10 mm) infrarenal necks.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

3.2 Patient selection

- Do not use the Endurant II/IIIs stent graft system in patients with sensitivities or allergies to the stent graft system materials. These materials include (* Denotes materials that are a part of the delivery system but do not contact the patient):
 - Polyether block amide (PEBA)
 - Polyether block amide (PEBA) with barium-sulfate filler
 - Polyamide 12
 - Acrylonitrile-butadiene-styrene (ABS) copolymer*
 - Stainless steel
 - Platinum-iridium alloy
 - Nickel-cobalt-chromium-molybdenum alloy
 - Nickel-titanium (Nitinol) alloy
 - Ethylene propylene rubber*
 - Silicone
 - Polycarbonate
 - Photo-polyacrylamide copolymer
 - Photoreactive quaternary salt
 - Photo-polyvinylpyrrolidone copolymer
 - Polyvinylpyrrolidone
 - Platinum
 - Polyester
 - Polyethylene
 - D&C blue no. 6
 - Tungsten
- Medtronic recommends that the physician disclose to the patient, in written form, all risks associated with treatment using the Endurant II/IIIs stent graft system. Details regarding risks occurring during and after implantation of the device are provided in *Adverse Events (Chapter 4)*.

- Inappropriate patient selection may result in poor device performance or device performance not otherwise in accordance with the specifications.
- The Endurant II/IIIs stent graft system is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging.
- The Endurant II/IIIs stent graft system is not recommended in patients exceeding weight or size limits necessary to meet imaging requirements.
- Key anatomic elements that may affect successful exclusion of the aneurysm include:
 - Severe proximal neck angulation with short proximal aortic neck ($>75^\circ$ with neck length <15 mm or $>60^\circ$ with neck length <10 mm);
 - Thrombus or calcium formation at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface;
 - Narrowing of the aorta at the bifurcation point;
 - Irregular calcification or plaque may compromise the fixation and sealing of the implantation sites.

Necks exhibiting these key anatomic elements may be more conducive to graft migration. The use of a bifurcated device in a patient with a narrowing of the distal aorta may result in reduced flow through the limbs or limb occlusion. The AUI stent graft is recommended for patients with significant distal aorta narrowing. For patients with a short aortic neck length, either a parallel graft technique or an infrarenal EVAR procedure using the Heli-FX EndoAnchor system should be considered.
- When a parallel graft technique is considered, patients should have the following anatomic characteristics:
 - Renal-aortic centerline angulation <90 degrees
 - Adequate renal artery seal zone length, without extensive calcification or thrombus that, in the operator's opinion, would preclude adequate sealing
 - Infrarenal neck ≥ 2 mm
- For infrarenal EVAR procedures using the Heli-FX EndoAnchor system:
 - The EndoAnchor implant should be implanted only into areas of aortic tissue that are free of or have insignificant calcified plaque or thrombus, or where such pathology is diffuse and less than 2 mm in thickness.
 - The EndoAnchor implant should not be used in patients with bleeding diathesis or with a known hypersensitivity or allergy to the device materials.
 - Refer to the Heli-FX EndoAnchor system instructions for use for additional information.
 - Methodologies for measurement of neck length vary. Where neck length is defined as the length over which the aortic diameter remains within 10% of the infrarenal diameter, other definitions may result in an estimation of neck length that is longer than that obtained using this definition.
- Iliac conduits may be used to ensure the safe insertion of the delivery system if the patient's access vessels, as determined by treating physician, preclude safe insertion of the delivery system.
- Deploying the stent graft in an area of vessel calcification may lead to abrasion of the stent graft on calcified plaque, potentially causing development of holes or tears in the graft.
- The long-term safety and effectiveness of the Endurant II/IIIs stent graft system has not been established.
- The safety and effectiveness of the Endurant II/IIIs stent graft system has not been evaluated in patients who:
 - Are less than 18 years old
 - Are pregnant or lactating
 - Have an aneurysm that is:
 - Suprarenal or pararenal
 - Isolated iliofemoral
 - Mycotic
 - Inflammatory
 - Pseudoaneurysm
 - Have a dominant patent inferior mesenteric artery and an occluded or stenotic celiac or superior mesenteric artery
 - Have an untreated thoracic aneurysm >4.5 cm in diameter
 - Requires emergent aneurysm treatment, for example, trauma or rupture
 - Have a history of bleeding diathesis or coagulopathy
 - Have had a myocardial infarction (MI) or cerebrovascular accident (CVA) within 3 months prior to implantation
 - Have a reversed conical neck, which is defined as a >4 mm distal diameter increase over a 10 mm length (infrarenal EVAR only)
 - Have a known hypersensitivity or contraindication to anticoagulants, antiplatelets, or contrast media, which is not amenable to pretreatment

- Have significant (typically >25% of vessel circumference of aortic neck and iliac artery, or >50% of the length of the iliac artery) aortic mural thrombus at either the proximal or distal attachment location that would compromise bilateral fixation and seal of the device
- Have ectatic iliac arteries requiring bilateral exclusion of hypogastric blood flow
- Have iliofemoral arterial access site that is not expected to accommodate the diameter of the device (14 to 20 Fr) due to size or tortuosity
- Have axillary arterial access site that is not expected to accommodate the diameter of sheaths compatible with the selected balloon expandable covered stent.
- Have active infection at the time of the index procedure documented by pain, fever, drainage, positive culture, or leukocytosis (WBC >11,000/mm³) that is treated with antimicrobial agents (nonprophylactic)
- Have congenital degenerative collagen disease
- Have a creatinine >2.0 mg/dL
- Are on dialysis
- Have connective tissue disorder
- The safety and effectiveness of the Endurant II/IIIs stent graft system with the Heli-FX EndoAnchor system has additionally not been evaluated in patients who have an infrarenal aortic neck with significant thrombus or calcium that precludes adequate EndoAnchor implant penetration of the aortic wall.
- All patients should be advised that endovascular treatment requires lifelong, regular follow-up care to assess their health and the performance of their endovascular graft. Patients with specific clinical findings (for example, endoleaks, enlarging aneurysms, stenosis of balloon-expandable covered stents, changes in the structure or position of the endovascular graft), or less than the recommended number of EndoAnchor implants when used in short proximal necks (≥4 mm and <10 mm), should receive enhanced follow-up care. Specific follow-up guidelines are described in *Chapter 12*.
- Patients experiencing reduced blood flow through the graft limbs, balloon-expandable covered stents, or leaks may be required to undergo secondary interventions or surgical procedures.
- Consider intervention or conversion to standard open surgical repair following the initial endovascular repair for patients experiencing enlarging aneurysms or endoleak. An increase in aneurysm size or persistent endoleak may lead to aneurysm rupture.

3.3 Before the implant procedure

- Perform preoperative planning for access and placement before opening the device packaging.
- Refer to the balloon expandable covered stent IFU for device preparation requirements (parallel graft technique only)
- For infrarenal EVAR procedures using the Heli-FX EndoAnchor system, refer to the Heli-FX EndoAnchor system instructions for use for device preparation requirements.
- Carefully inspect the Endurant II/IIIs stent graft system packaging and device for damage or defects before use. Do not use the product if any sign of damage or breach of the sterile barrier is observed. Do not attempt to resterilize the delivery system or stent graft.
- Do not bend, kink, or otherwise alter the Endurant II/IIIs delivery system prior to implantation because it may cause deployment difficulties.
- To reduce the risk of thrombotic complications, administer an additional bolus of IV heparin before inserting the device.

3.4 During the implant procedure

- Exercise care in handling and delivery technique to help prevent vessel rupture.
- Studies indicate that the danger of micro-embolization increases with increased procedure duration.
- Renal complications may occur:
 - from an excess use of contrast agents
 - as a result of embolic or misplaced stent graft or balloon expandable covered stent
- Exercise care and utilize suitable imaging techniques when deploying the aortic endograft into a short proximal neck, to ensure accurate placement. Inaccurate placement could result in unsuccessful implantation of EndoAnchor implants, the need to place a proximal extension, or unintentional artery coverage.
- Do not deploy the stent grafts in a location that will cause an endoleak or occlude arteries necessary to supply blood flow to organs or extremities. This could necessitate surgical removal of the device.
- Use fluoroscopic guidance to advance the delivery system, detect kinking, and assess alignment problems with the stent graft devices. Do not use excessive force to advance or withdraw the delivery system when resistance is encountered. If the delivery system kinks during insertion, do not attempt to deploy the stent graft. Remove the device and insert a new system.
- Do not continue to torque the delivery system if the tip is not rotating with the delivery system.

- Exercise particular care in difficult areas, such as areas of stenosis, intravascular thrombosis, or in calcified or tortuous vessels. Balloon angioplasty at the site of a narrowed or stenotic vessel may be considered prior to attempting to gently reintroduce the catheter delivery system.
- An inadequate seal zone may result in increased risk of leakage into the aneurysm or migration of the stent graft or balloon expandable covered stent.
- Systemic anticoagulation should be used during the implantation procedure based on hospital or physician protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- A stent graft that has been deployed, even if only partially, cannot be repositioned or drawn back into the delivery system.
- If the graft cover is accidentally withdrawn, the device will prematurely deploy and may be incorrectly positioned.
- For the limb stent graft overlap criteria with Endurant IIs ipsilateral leg only, refer to *Table 8*. As noted, do not overlap more than 3 stents for the limb stent graft configurations that have an overlap criteria of 3 stents only.
- When deploying the stent graft, be sure to hold the front grip of the delivery system stationary.
- If a compliant balloon catheter is used, do not overinflate or inflate outside the graft material. Follow all manufacturer instructions regarding catheter operation.
- High-pressure injections of contrast media made at the edges of the stent graft immediately after implantation can cause endoleaks.

3.5 Treatment and follow-up

- Physicians should review the instructions for use for applicable information to be shared with the patient.
- A patient implant card is included in the device package and contains identifying information about the implanted device. After device implant, complete the patient implant card and provide it to the patient before they are discharged. Patients can access information about their device on the website that is listed on their patient implant card. Physicians should communicate the following instructions to their patients:
 - Always carry their implant card with them.
 - Access additional information about their device on the website that is listed on their patient implant card.

Note: If the patient is unable to access the website, the physician must provide the information from the website to the patient.

 - Always inform any healthcare personnel that they have an implanted device before any procedure has begun.
 - Contact their healthcare provider if they notice any new or changing symptoms.
- Any endoleak left untreated during the implantation procedure must be carefully monitored after implantation.
- All patients with endovascular aneurysm repair should undergo periodic imaging to evaluate the stent graft, aneurysm size, and occlusion of vessels in the treatment area. Significant aneurysm enlargement (>5 mm), the appearance of a new endoleak, evidence of perigraft flow, change in aneurysm pulsatility, or migration resulting in an inadequate seal zone should prompt further investigation and may indicate the need for additional intervention or surgical conversion.
- Additional treatment, including endovascular treatment or surgical conversion, should be strongly considered in the following cases:
 - Aneurysm growth >5 mm, with or without endoleak, since last follow-up
 - Change in aneurysm pulsatility, with or without growth or endoleak
 - Persistent endoleak, with or without aneurysm growth
 - Stent graft migration resulting in an inadequate seal zone
 - Stenosis within a balloon expandable covered stent or at its distal attachment site
 - Decrease in renal function due to renal artery occlusion (migration or poor placement)
- Following endovascular aneurysm repair (EVAR), spinal cord ischemia (SCI) may result in a rare complication of paraplegia or paraparesis. Cerebrospinal fluid (CSF) drain is advised if spinal cord ischemia is suspected.
- For infrarenal EVAR procedures using the Heli-FX EndoAnchor system:
 - Always use fluoroscopy for guidance, delivery, and observation of any Heli-FX system components within the vasculature.
 - Medtronic recommends that the implantation of EndoAnchor implants be done after the aortic endograft has been placed and any balloon remodeling of the infrarenal seal zone of the stent graft system has been completed. Exercise care in balloon remodeling of the stent graft system to avoid moving the main body endograft from its intended implant location.
 - EndoAnchor implant locations should be based upon a detailed examination of the preoperative CT imaging in cases involving irregular or eccentric plaque in the intended sealing zone(s). The EndoAnchor implant should be implanted only into areas of aortic tissue that are free of or have insignificant calcified plaque or thrombus, or where such pathology is diffuse and less than 2 mm in thickness. Attempting to place EndoAnchor implants into more severe plaque or thrombus may be associated with implantation difficulty and suboptimal endograft fixation and/or sealing.
 - The recommended number of EndoAnchor implants for a bifurcated endograft is based on native vessel diameter and is independent of the amount of endograft oversizing.

- Stability of the stent graft in short (≥ 4 mm and < 10 mm) infrarenal necks is augmented by the EndoAnchor implants. Ensure successful deployment of the recommended minimum number of EndoAnchor implants. Where the number of successfully deployed EndoAnchor implants is below the minimum recommended, there may be greater risk of postoperative Type 1a endoleak or migration.

4 Adverse events

4.1 Potential adverse events

Adverse events or complications associated with the use of this stent graft system that may occur or require intervention include, but are not limited to, the following:

- Access failure
- Allergic/hypersensitivity reaction (to contrast, antiplatelet therapy, stent graft material)
- Anesthetic complications and subsequent problems (for example, aspiration)
- Aneurysm enlargement
- Aneurysm rupture
- Aortic damage, including perforation, dissection, bleeding, rupture
- Aortoenteric fistula
- Arterial or venous pseudoaneurysm
- Arterial or venous thrombosis
- Bleeding, hematoma, or coagulopathy
- Bowel complications (for example, ileus, transient ischemia, infarction, necrosis, obstruction)
- Cardiac complications and subsequent problems (for example, arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension, angina, cardiac tamponade, pericarditis)
- Catheter breakage
- Claudication (for example, buttock, lower limb)
- Contrast toxicity
- Death
- Deployment difficulties/failures
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleaks
- Excessive or inappropriate radiation exposure
- Fever and inflammation
- Genitourinary complications and subsequent problems (for example, ischemia, erosion, fistula, incontinence, hematuria, infection)
- Hepatic failure
- Infection of the aneurysm, device access site, including abscess formation, transient fever, and pain, possibly leading to sepsis
- Insertion or removal difficulty
- Local or systemic neurological complications and subsequent problems (for example, confusion, stroke, transient ischemic attack, paraplegia, paraparesis, paralysis, nerve injury, neuropathy)
- Lymphatic complications and subsequent problems (for example, lymph fistula)
- Pulmonary complications and subsequent problems (for example, pulmonary embolism, pleural effusion, pneumonia, pneumothorax, pulmonary edema, respiratory depression or failure)
- Renal complications and subsequent problems (for example, artery occlusion, insufficiency, failure)
- Sexual dysfunction
- Stenosis
- Stent graft migration
- Stent graft misplacement
- Stent graft occlusion
- Stent graft rupture or fracture
- Stent graft suture break
- Stent graft twisting or kinking
- Surgical conversion to open repair

- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula, stenosis, dissection
- Vascular spasm or vascular trauma, or vessel damage (for example, extrusion/erosion, iliofemoral vessel dissection, bleeding, rupture, death)
- Vessel occlusion
- Wound complications and subsequent problems (for example, dehiscence, infection, hematoma, seroma, cellulitis)

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

5 Patient selection and treatment

5.1 Individualization of treatment

Each Endurant II/IIIs stent graft system must be ordered in the appropriate size to fit the patient's anatomy. Proper sizing of the device is the responsibility of the physician. The stent graft should be oversized to be larger than the vessel inner diameter. Aortic configurations are recommended to be oversized 10% to 20% (for infrarenal/aortoiliac aneurysm treatment) or 20% to 30% when used with a balloon expandable covered stent (for juxtarenal aneurysm treatment); iliac configurations are recommended to be oversized 10% to 25%. Refer to *Section 9.2* for further details. The stent graft configurations cover aortic diameters ranging from 19 mm to 32 mm or 19 to 30 mm when used with parallel grafts, and iliac diameters from 8 to 25 mm. The recommended overall length of the stent graft, including multiple deployed devices, should extend from the lowest renal artery to just above the internal iliac or hypogastric artery. In the case of use with a balloon expandable covered stent, the aortic configuration can be positioned above the renal arteries, requiring additional treatment length. Refer to the Instructions for Use provided with the balloon expandable covered stent for recommended sizing guidelines. All lengths and diameters of the stent graft devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters or lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

Medtronic may consult with physicians to determine proper stent graft dimensions based on the physician's assessment of the patient's anatomical measurements. The benefits and risks previously described should be carefully considered for each patient before use of the stent graft system.

Patient selection factors to be assessed should include:

- Patient's age and life expectancy
- Comorbidities (for example, cardiac, pulmonary, or renal insufficiency prior to surgery)
- Patient morphologic suitability for endovascular repair
- The risk of aneurysm rupture balanced against the risk of treatment with the stent graft system

Note: Due to the nature of the design and the flexibility of the Endurant II/IIIs stent graft system, the overall length of each stent graft may be shorter when deployed.

Caution: Vessel overdilation and damage, or partial stent graft infolding, may be caused by excessive oversizing of the stent graft in relation to the diameter of the blood vessel.

For infrarenal EVAR procedures using the Heli-FX EndoAnchor system:

- The EndoAnchor implant should be implanted only into areas of aortic tissue that are free of or have insignificant calcified plaque or thrombus, or where such pathology is diffuse and less than 2 mm in thickness.
- The EndoAnchor implant should not be used in patients with bleeding diathesis or with a known hypersensitivity or allergy to the device materials.
- Refer to the Heli-FX EndoAnchor system instructions for use for additional information.

Note: Clinical experience from Medtronic's ANCHOR Registry supports the use of Endurant II/IIIs stent graft system in conjunction with Heli-FX EndoAnchor system in the treatment of abdominal aortic aneurysms with short proximal neck (≥ 4 mm and < 10 mm). For the purposes of evaluating outcomes of short neck subjects in the registry, neck length was defined by the core lab as that length over which the aortic diameter remains within 10% of the infrarenal diameter.

6 Patient counseling information

The physician should review the following risks and benefits when counseling the patient about this endovascular device and procedure:

- Patient age and life expectancy
- Differences between endovascular repair and surgical repair.
- Risks and benefits related to open surgical repair
- Risks and benefits related to endovascular repair
- Risks related to noninterventional treatment or medical management

- Risks of aneurysm rupture compared to endovascular repair
- Possibility that subsequent endovascular or open surgical repair of the aneurysm may be required
- The long-term safety and effectiveness of the Endurant II/IIIs stent graft system or the Endurant II/IIIs stent graft system used in conjunction with the Heli-FX EndoAnchor system has not been established
- Long-term, regular follow-up care is needed to assess patient health status and stent graft performance
- Patients with specific clinical findings (for example, endoleaks, enlarging aneurysms) should be monitored closely
- Symptoms of aneurysm rupture

Medtronic recommends that the physician disclose to the patient, in written form, all risks associated with treatment using the Endurant II/IIIs stent graft system. Details regarding risks occurring during and after implantation of the device are provided in Adverse Events (*Chapter 4*).

7 How supplied

7.1 Sterility

Each stent graft configuration (bifurcated, AUI, limb, aortic and iliac extension, and abdominal tube) is individually contained within a delivery system. It is sterilized using electron beam and is supplied sterile for single use only.

- Do not reuse or attempt to resterilize.
- If the device is damaged or the integrity of the sterile barrier has been compromised, do not use the product and contact your Medtronic representative for return information.

7.2 Contents

- One Endurant II/IIIs stent graft system
- One web reference leaflet

For information on the supply of the Heli-FX EndoAnchor system, please see the instructions for use supplied with the device.

8 Device disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

9 Clinical use information

9.1 Physician training requirements

All physicians should be trained in the use of the Endurant II/IIIs stent graft system before using it.

Caution: The Endurant II/IIIs stent graft system should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device.

The following are the knowledge and skill requirements for physicians using the Endurant II/IIIs stent graft system:

- Natural history of abdominal aortic aneurysms (AAA), aortoiliac aneurysms, and comorbidities associated with AAA repair
- Use of balloon expandable covered stents used for parallel graft technique
- Use of EndoAnchor implants for infrarenal EVAR using the Heli-FX EndoAnchor system
- Radiographic, fluoroscopic, and angiographic image interpretation
- Appropriate use of radiographic contrast material
- Arterial cutdown, arteriotomy, and repair or percutaneous access and closure techniques
- Nonselective and selective guidewire and catheter techniques
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Techniques to minimize radiation exposure
- Device selection and sizing

Caution: The Heli-FX EndoAnchor system should only be used by physicians and teams trained in the use of this device. For physician training recommendations for the Heli-FX EndoAnchor system, please see the instructions for use provided with the device.

9.2 Recommended device sizing

The Endurant II/IIIs stent graft systems are available in the sizes described in *Table 2* through *Table 7*. For questions about Endurant II/IIIs sizing, refer to contact information in the back of this Instructions for Use. For the balloon expandable covered stent sizing guidelines, refer to the Instructions for Use provided with the device.

Table 2. Sizing chart — Endurant II bifurcated configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)	
			Infrarenal EVAR	Parallel graft technique
20	36x20	145, 166	29-32	27-30
	36x16			
	32x20	124, 145, 166	26-28	24-26
	32x16			
28x20				
28x16				
18	28x13	124, 145, 166	23-25	21-23
	25x16			
	25x13			
	23x16		21-22	19-20
	23x13			

Table 3. Sizing chart — Endurant IIs bifurcated configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)	
			Infrarenal EVAR	Parallel graft technique
20	36x14	103	29-32	27-30
	32x14		26-28	24-26
18	28x14		23-25	21-23
	25x14		21-22	19-20
	23x14		19-20	N/A

Table 4. Sizing chart — limb configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)
16	16x28	82, 93, 124, 146, 156, 199	23-25
	16x24		19-22
	16x20		15-18
	16x16	146, 156, 199	12-14
	16x13		10-11
	16x10		8-9
14	16x16	82, 93, 124	12-14
	16x13		10-11
	16x10		8-9

Table 5. Sizing chart — iliac extension configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)
18	28x28	82	23-25
16	24x24		19-22
	20x20		15-18
14	13x13		10-11
	10x10		8-9

Table 6. Sizing chart — aortic extension and abdominal tube configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)	
			Infrarenal EVAR	Parallel graft technique
20	36x36	49, 70	29-32	N/A
	32x32		26-28	N/A
18	28x28		23-25	N/A
	25x25		21-22	N/A
	23x23		19-20	N/A

Table 7. Sizing chart — aorto-uni-iliac (AUI) configuration

OD (Fr)	Proximal x distal diameter (mm x mm)	Covered length (mm)	Vessel inner diameter (mm)	
			Infrarenal EVAR	Parallel graft technique
20	36x14	102	29-32	27-30
	32x14		26-28	24-26
18	28x14		23-25	21-23
	25x14		21-22	19-20
	23x14		19-20	N/A

9.3 Device inspection

Inspect the Endurant II/IIIs device and packaging to verify that no damage or defects exist. If the “Use by” date has elapsed, the device is damaged, or the sterile barrier has been compromised, do not use the device and contact a Medtronic representative for return or replacement information.

9.4 Required materials (not included in the stent graft system packaging)

At the time of surgery have available:

- Additional Endurant II/IIIs stent graft systems (bifurcated, contralateral limb, aortic extension, abdominal tube, and iliac extension components) of various lengths and diameters to customize the implant to fit the anatomy of the individual patient.
- Balloon expandable covered stents in various lengths and diameters (parallel graft technique only)
- For infrarenal EVAR procedures, additional Heli-FX EndoAnchor systems and compatible introducer sheaths.
- Fluoroscope with digital angiographic capabilities (C-arm or fixed unit). Fluoroscopic imaging and the ability to record and recall all imaging.
- Assorted guidewires of adequate length. In addition to guidewires used for accessing the vessel, 0.035 in (0.89 mm) diameter guidewires or equivalents must be used to maximally support the Endurant II/IIIs delivery system into the aortic vasculature.
- Heparinized saline solution

9.5 Recommended materials (not included in the stent graft system packaging)

At the time of surgery have available:

- An appropriately sized introducer sheath to provide an adequate conduit for the Endurant II/IIIs delivery system to be used.
- Sterile introducer sheaths for introduction into femoral arteries during road mapping or further diagnostic imaging
- Power injector for angiographic contrast studies
- Radiopaque ruler with centimeter increments
- Compliant balloon catheters for stent graft expansion following implantation
- Assorted diagnostic catheters
- Assorted balloon catheters to potentially dilate blood vessels prior to insertion of the Endurant II/IIIs delivery system
- Radiopaque contrast media
- Sterile silicone lubricant or sterile mineral oil
- Interventional snare devices
- Endovascular coils, vascular plugs, and vascular occluding devices (for example, Talent occluder — see *Section 1.1.8* and *Section 10.2.16*)

9.6 MRI safety information

Nonclinical testing has demonstrated that the Endurant II/IIIs stent graft is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 or 3.0 T only
- Maximum spatial gradient magnetic field of 2500 gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Mode)

Under the scan conditions defined above, the Endurant II/IIIs stent graft is expected to produce:

- Maximum temperature rise of 1.00°C after 15 minutes of continuous scanning in a 1.5 tesla scanner
- Maximum temperature rise of 3.27°C after 15 minutes of a continuous scanning in a 3.0 tesla scanner

The image artifact extends approximately 5 mm and 8 mm from the device, both inside and outside the device lumen when scanned in nonclinical testing using the sequence: spin echo and gradient echo, respectively, in a 3.0 T Siemens TrioTim (VB 13 Software) MR system with a whole-body coil.

If multiple devices are present, the most restrictive MRI exposure requirements must be used. Contact the appropriate device manufacturers if you have questions.

10 Implant instructions

10.1 Vascular access and device preparation

Correct sizing of the aorta and iliac vessels must be determined before implantation of the aortic and iliac stent graft configurations using contrast-enhanced computer-aided tomography (CT), as well as angiograms of both the iliac arteries and aorta. 3D imaging may also be beneficial. Refer to *Section 9.2*. These images should be available for review during the procedure. Vascular instruments and other surgical supplies necessary to gain access to the artery should also be available.

To reduce the risk of thromboembolism, it is recommended that the patient be heparinized for the duration of the procedure.

Caution: Do not retract the graft cover of the delivery system until it is accurately placed within the vasculature and ready for deployment.

Caution: Never advance or retract equipment from the vasculature without the use of fluoroscopy.

10.1.1 Vascular access

1. Following aseptic procedure, perform a vascular access at the femoral arteries.
2. Place a guidewire in the ipsilateral femoral artery and advance it above the renal arteries.
3. From the contralateral side femoral artery, place a second guidewire directed to the abdominal aorta.
4. Over this guidewire, place an angiography catheter above the renal arteries.
5. Take an angiogram.

Note: An additional incision might be necessary to access the common iliac artery.

10.1.2 Renal cannulation through axillary or brachial access (parallel graft technique only)

1. Following aseptic procedure, perform a vascular access at the left axillary or brachial artery.
2. Cannulate target renal artery or arteries with standard interventional technique
3. Using standard technique, exchange for an appropriately stiff wire with a curved atraumatic tip
4. Following cannulation of the renal artery or arteries, advance an appropriately sized introducer sheath into the renal artery to maintain access.
5. Prepare balloon expandable covered stent delivery system according to the device Instruction for Use
6. Advance the balloon expandable covered stent or stents into the target renal artery or arteries.

Note: Monitor anticoagulation during procedures using activated clotting time or an equivalent method, especially in the case of multiple chimney graft placement or prolonged operation time.


Note: Use superselective renal angiography to verify cannulation and exclude possible injury to the renal artery.

Note: A molding balloon can be helpful to stabilize catheters and sheaths during the cannulation process.

Caution: Brachial or axillary access through a right sided approach is associated with a higher rate of cerebral events.

Caution: Calcium and thrombus may interfere with cannulation of the renal artery or arteries.

Figure 4. Renal cannulation (parallel graft technique only)

10.1.3 Device preparation

1. View the delivery system under fluoroscopy to visualize the radiopaque markers on the stent graft. The radiopaque markers indicate the position of the proximal and distal edges of the graft material.
2. Turn the graft cover to align the radiopaque gate marker on the short leg of the bifurcated configuration with the patient's contralateral iliac artery.
3. Flush the guidewire lumen with heparinized saline solution.
4. Prior to insertion into the vessel, activate the hydrophilic coating by wiping the outer surface of the graft cover with a sterile gauze, saturated in saline, until the graft cover is slippery to touch.

10.2 Delivery procedure

Medtronic recommends using an appropriate caliber introducer sheath to perform diagnostic tests. No sheath is necessary for the introduction of the delivery system or deployment of the stent graft. For infrarenal EVAR procedures using the Heli-FX EndoAnchor system, the access vessel diameter and morphology should be compatible for use with the device and should accommodate a 16 Fr introducer sheath. An infrarenal EVAR procedure using the Heli-FX EndoAnchor system is performed at the discretion of the implanting physician. Medtronic recommends that the implantation of EndoAnchor implants be done after the aortic endograft has been placed and any balloon remodeling of the infrarenal seal zone of the stent graft system has been completed. See *Section 10.2.18*.

Caution: Do not remove the guidewire while the delivery system is in the patient.

Warning: To prevent thrombotic problems, a second bolus of IV heparin is recommended before inserting the device.

10.2.1 Introduction of bifurcated configuration

Warning: Do not advance the delivery system without placing a guidewire.

1. Slowly insert the delivery system.
2. Advance over the guidewire so that the most proximal stents and the radiopaque markers are visualized in the target proximal aortic neck (*Figure 5*).

3. Inject contrast media through an angiographic (pigtail) catheter into the abdominal aorta and mark the position of the target location, either on the imaging screen or on the patient's body.
4. Adjust the position of the bifurcated stent graft configuration so that the top edge of the graft fabric is below the lowest untreated visceral vessel.

Note: The edge of the graft fabric is 0.5 mm to 1.0 mm above the top edge proximal radiopaque markers.

Note: If the top edge of the graft fabric is to be placed close to an untreated visceral vessel, contrast media may be injected to identify the location of the lowest untreated visceral vessel and verify the position before full deployment.


Note: For a parallel graft technique, carefully monitor the position of the stent graft during balloon remodeling as proximal migration of the stent graft may occur.

Caution: Once proximal position has been identified, do not move the patient or imaging equipment, as it may compromise accuracy of stent graft placement.

Caution: The angiographic catheter can be removed prior to deployment. However, if the angiographic catheter is not removed until after deployment, ensure that the tip is straightened (such as with a pigtail catheter) with a guidewire before removal so that the stent graft is not pulled down.

Caution: When aligning the position of stent graft, be sure the fluoroscope is angled perpendicular to the center line of the infrarenal aorta to avoid parallax or other sources of visualization error. Some cranial caudal angulation of the image intensifier (I-I) tube may be necessary, especially if there is anterior angulation of the aneurysm neck.

Figure 5. Introduce the aortic delivery system

1 Infrarenal EVAR

2 Parallel graft technique

10.2.2 Confirm position

1. Ensure that the distal portion of the contralateral stub leg is above the aortic bifurcation and within the aneurysmal sac, and not within the iliac vessel.
2. Rotate the handle until the radiopaque marker on the distal stent of the contralateral stub leg is aligned with the contralateral iliac artery.

Note: When attempting to rotate the system, if the tip does not rotate with the handle, pull back the system and reposition until the intended position is achieved.

10.2.3 Deploy proximal end of bifurcated configuration

1. With 1 hand on the front grip, hold the delivery system stationary.
2. With the other hand, slowly withdraw the graft cover by rotating the external slider counterclockwise (in the direction of the slider arrow), until the constrained suprarenal stent is exposed and 2 to 3 of the Endurant II or Endurant IIs body stents have been fully deployed (*Figure 6*).
3. Use angiography to verify position of the bifurcated configuration in relation to the lowest untreated visceral vessel.
4. If needed, gently push the entire delivery system proximally or pull distally until the proximal end of the graft material is distal to the lowest untreated visceral vessel.


Note: In the unlikely event of delivery system failure that results in partial stent graft deployment due to graft cover severance, the “handle disassembly” technique may permit successful deployment of the stent graft. Refer to *Chapter 11*.

Caution: Do not rotate the graft cover during deployment as this may torque the device and cause it to rotate during deployment.

Caution: If the graft cover is accidentally withdrawn, the stent graft will prematurely deploy and may be incorrectly positioned.

Warning: Failure to properly align the radiopaque markers may result in improper deployment of the stent graft.

Figure 6. Deploy the proximal end of the bifurcated configuration


1 Infrarenal EVAR

2 Parallel graft technique

10.2.4 Deploy contralateral leg of bifurcated configuration

Continue holding the front grip of the delivery system stationary and then rotate the slider handle counterclockwise, stopping immediately after the contralateral leg is released from the graft cover or delivery sheath (*Figure 7*).

Figure 7. Deploy the contralateral leg

- 1 Infrarenal EVAR
- 2 Parallel graft technique


10.2.5 Release proximal end of suprarenal stent

1. Use angiography to verify the position of the bifurcated configuration in relation to the lowest untreated visceral vessel.
2. Continue to hold the delivery system stationary with 1 hand on the front grip.
3. With the other hand, rotate the back-end wheel clockwise, moving the tapered tip forward to release the proximal end of the suprarenal stent (*Figure 8*).
4. Observe the release of the suprarenal stent under fluoroscopy and continue turning the back-end wheel until it is completely clear of the delivery system spindle.

Note: In the unlikely event that the proximal end of the suprarenal stent cannot be released, refer to *Chapter 11*.

Caution: In the unlikely event that the back-end wheel separates during wheel rotation, remove the wheel. Manually advance the exposed tabs forward on the screw gear until all of the suprarenal stents release from the spindle. Refer to *Chapter 11*.

Figure 8. Release the proximal end of the suprarenal stent

- 1 Infrarenal EVAR
- 2 Parallel graft technique


10.2.6 Deploy distal end of bifurcated configuration

Either continue to rotate the external slider counterclockwise or while holding the front grip of the delivery system stationary, use thumb to pull the trigger on the external slider and pull it back until the bifurcated stent graft is completely deployed.

Note: Retract the graft cover past the flexible stent stop tip (approximately 10 mm) to ensure that the graft cover edge does not disturb the graft position during forward advancement of the catheter for tip recapture.

Caution: When using the trigger to rapidly deploy the stent graft, the delivery system must remain stationary. Do not rotate the delivery system during stent graft deployment.

Figure 9. Deploy the distal end of the bifurcated configuration

- 1 Endurant II bifurcated configuration
- 2 Endurant IIs bifurcated configuration

10.2.7 Recapture spindle in tapered tip

Note: For the Endurant II bifurcated configuration, leave the delivery system in situ while deploying the limb stent graft into the contralateral leg.

1. Continue to hold the delivery system stationary with 1 hand on the front grip.
2. Confirm that the spindle has fully separated from the suprarenal stent; gently torque the delivery system if it has not fully separated.
3. Gently rotate the delivery system while pushing the entire delivery system approximately 3 cm proximally so that the tapered tip and spindle are completely clear of the suprarenal stent.
4. With the other hand, rotate the back-end wheel counterclockwise to recapture the spindle in the tapered tip (*Figure 10*).
5. Observe the recapture of the spindle within the sleeve of the tapered tip under fluoroscopy.
6. Continue turning the back-end wheel counterclockwise until the spindle has been completely recaptured and the back-end wheel is at the bottom (*Figure 10*).

Note: When pushing the delivery system forward, be careful not to displace the distal end of the ipsilateral limb.


Note: Ensure that the suprarenal stent is fully disengaged from the spindle before pushing the delivery system forward.

Note: If the spindle catches on the suprarenal stent during advancement, completely advance the back-end wheel clockwise. Using a gentle in-and-out motion with the delivery system, rotate the delivery system until the spindle slips past the suprarenal stent. Then continue with the withdrawal process.

Caution: Stop rotating the back-end wheel when the bottom of the back-end screw gear is reached.

Warning: Failure to adequately advance the delivery system to recapture the spindle can result in the trapping of a suprarenal apex within the tapered tip sleeve. This will alter the proximal landing zone during delivery system withdrawal.

Figure 10. Recapture the spindle in the tapered tip

10.2.8 Remove delivery system (Endurant II bifurcated configuration only)


1. Continue to hold the delivery system stationary with 1 hand on the front grip and the other hand on the external slider.
2. Gently torque and withdraw the delivery system until the spindle is retracted into the fabric portion of the stent graft.
3. Pull back the external slider trigger and hold it stationary while bringing the front grip to the slider (*Figure 11*).

4. Use continual fluoroscopy and watch the top of the stent graft while slowly pulling back the tapered tip into the graft cover of the delivery system.
5. Gently remove the delivery system. Use fluoroscopy to ensure that the bifurcated configuration does not move during withdrawal.

Note: Maintain vessel access and wire placement until all stent graft configurations are in place.

6. For parallel graft technique procedures, deploy the balloon expandable covered stent. See *Section 10.2.9*.

Figure 11. Remove delivery system


10.2.9 Balloon expandable covered stent deployment (parallel graft technique only)

1. Partially withdraw the main delivery system (Endurant IIs only) until the spindle is distal to the proximal fabric edge of the Endurant IIs stent graft.
2. Verify that the balloon expandable covered stent sizing is correct, such that once deployed, the proximal edge is midway between the aortic bifurcate proximal graft edge and the suprarenal anchor pins and that there is sufficient distal seal length within the renal artery.
3. Immediately prior to deployment of the balloon expandable covered stent, withdraw the introducer sheath.
Note: To complete deployment or release the balloon expandable covered stent, the distal tip of the introducer sheath should be positioned below the suprarenal anchor pins.
Caution: If the introducer sheath is withdrawn up to or beyond the suprarenal anchor pins, rupture of the balloon expandable stent catheter is possible.
4. Deploy the parallel stent graft according to the balloon expandable covered stent Instructions for Use (*Figure 12*).

Note: Following deployment of the balloon expandable covered stent, you may insert the introducer sheath into the balloon expandable covered stent to maintain access.


Figure 12. Deploy balloon expandable covered stent or stents (parallel graft technique only)

10.2.10 Deploy limb stent graft into contralateral leg

1. Prepare the iliac stent graft system as described in *Section 10.1.3*.
2. On the patient's contralateral side, insert a guidewire through the contralateral leg and aortic neck of the previously placed bifurcated configuration.
3. Place the delivery system over the guidewire and into the contralateral leg of the bifurcated stent graft.
4. Insert the limb stent graft into the contralateral leg of the bifurcated stent graft. The proximal radiopaque marker of the limb should be aligned to the radiopaque marker at the bifurcation of the bifurcated stent graft configuration (*Figure 13*).

Figure 13. Introduce the iliac delivery system


- 1 Endurant II bifurcated configuration (Infrarenal EVAR)
- 2 Endurant IIs bifurcated configuration (Parallel graft technique)

5. Ensure there is a 3-stent overlap (*Figure 14*).
6. With 1 hand on the front grip, hold the delivery system stationary.
7. With the other hand, slowly withdraw the graft cover by rotating the external slider counterclockwise.
8. At any point, pull the slider trigger and pull the external slider all the way back to finish deploying the limb configuration.
9. Remove the delivery system.

Note: In the unlikely event of delivery system failure that results in partial stent graft deployment, the “handle disassembly” technique may permit the successful deployment of the stent graft configuration. Refer to *Chapter 11*.

Caution: Do not torque the iliac delivery system while it is in the patient.

Figure 14. Deploy the limb configuration


- 1 Endurant II bifurcated configuration
- 2 Endurant IIs bifurcated configuration

10.2.11 Remove delivery system (Endurant IIs bifurcated configuration only)

1. Continue to hold the delivery system stationary with 1 hand on the front grip and the other hand on the external slider.
2. Gently torque and withdraw the delivery system until the spindle is retracted into the fabric portion of the stent graft.
3. Pull back the external slider trigger and hold it stationary while bringing the front grip to the slider (*Figure 15*).
4. Use continual fluoroscopy and watch the top of the bifurcated stent graft while slowly pulling back the tapered tip into the graft cover of the delivery system.
5. Gently remove the delivery system. Use fluoroscopy to ensure that the bifurcated configuration does not move during withdrawal.

Note: Maintain vessel access and wire placement until all stent graft configurations are in place.

Figure 15. Remove the delivery system

10.2.12 Deploy the limb stent graft into the ipsilateral leg (Endurant IIs bifurcated configuration only)

1. Prepare the iliac stent graft system as described in *Section 10.1.3*.
2. On the patient's ipsilateral side, track the delivery system over the guidewire and into the ipsilateral leg of the previously placed bifurcated configuration. The limb stent graft has 2 markers on its proximal edge, 2 markers on its distal edge, and 1 overlap marker approximately 25 mm distal of the proximal markers.
3. Position the device. Overlap criteria between the limb stent graft and the ipsilateral leg of the bifurcated stent graft is dependent on limb selection. See *Table 8* for recommended device overlap.

Table 8. Recommended device overlap — limb stent graft and ipsilateral leg of Endurant IIs bifurcated stent graft

Proximal diameter	Distal diameter	Length	Overlap	Deployment reference
16	10	82	3 stents only	Refer to step a
		93		
	13	82		
		93		
	20	82		
		93		
	24	82		
		93		
28	82			
	93			


Table 8. Recommended device overlap — limb stent graft and ipsilateral leg of Endurant IIs bifurcated stent graft (continued)

Proximal diameter	Distal diameter	Length	Overlap	Deployment reference
16	10	124	3 to 5 stents	Refer to step b
		146		
		156		
		199		
	13	124		
		146		
		156		
		199		
	16	82		
		93		
		124		
		146		
		156		
		199		
	20	124		
		146		
		156		
		199		
	24	124		
		146		
156				
199				
28	124			
	146			
	156			
	199			

a. 3 stent overlap: Align the overlap marker on the limb stent graft with the 2 markers on the distal end of the ipsilateral leg of the Endurant IIs bifurcated stent graft (*Figure 16*).


Warning: For the limb stent graft overlap criteria with Endurant IIs ipsilateral leg only, refer to *Table 8*. As noted in *Table 8*, do not overlap more than 3 stents for the limb stent graft configurations that have an overlap criteria of only 3 stents.

Figure 16. Introduce the iliac delivery system — 3 stent overlap

- b. 3 to 5 stent overlap: For the limb stent grafts that can have an overlap between 3 to 5 stents, minimum 3 stent overlap can be achieved by following the instruction in step a (*Figure 16*) or the maximum 5 stent overlap can be achieved by aligning the proximal edge markers on the limb stent graft to the bifurcation marker on the Endurant IIs bifurcated stent graft (*Figure 17*).

Figure 17. Introduce the iliac delivery system — 5 stent overlap

4. With 1 hand on the front grip, hold the delivery system stationary.
5. With the other hand, slowly withdraw the graft cover by rotating the external slide counterclockwise.
6. At any point, pull the slider trigger and pull the external slider all the way back to finish deploying the limb configuration.
7. Remove the delivery system (*Section 10.2.11*).

Note: In the unlikely event of delivery system failure that results in partial stent graft deployment, the “handle disassembly” technique may permit the successful deployment of the stent graft configuration. Refer to *Chapter 11*.

Caution: Do not torque the iliac delivery system while in the patient.

10.2.13 Iliac or aortic extension stent graft configurations

1. If an aortic extension stent graft configuration is needed for an infrarenal EVAR procedure, ensure that there is a minimum 3-stent overlap between the aortic extension stent graft and the bifurcated stent graft configuration.
2. Follow the bifurcated stent graft deployment process, except rotate the handle to open the extension component entirely before releasing the proximal end of the suprarenal stent of the aortic configuration.

3. If an iliac extension stent graft configuration is needed for an infrarenal EVAR or parallel graft procedure, ensure that there is a minimum 3-stent overlap between the iliac extension and the configuration it is being mated with. This is achieved by aligning the limb overlap marker with the most distal marker on the configuration with which the extension is being mated.
4. Follow the limb stent graft configuration deployment process (*Section 10.2.10*).

10.2.14 Abdominal tube stent graft configuration

If an abdominal tube stent graft configuration is needed, follow the bifurcated stent graft deployment process, with the following exception: Rotate the handle to open the abdominal tube stent graft entirely before releasing the proximal end of the suprarenal stent.

10.2.15 AUI stent graft configuration

Follow the bifurcated stent graft configuration deployment steps, except rotate the handle to open the AUI stent graft covered portion completely before releasing the proximal end of the suprarenal stent of the AUI stent graft. Follow the steps described in *Section 10.2.7* to *Section 10.2.8* to remove the delivery system.

In the event that an additional distal stent graft configuration is required, use the limb stent graft as the AUI distal extension. Place the delivery system over the existing guidewire and follow the limb stent graft deployment steps in *Section 10.2.10*. To ensure proper docking between the 2 stent grafts, align the overlap radiopaque marker on the limb stent graft with the distal radiopaque markers on the AUI stent graft to ensure a 3 stent overlap.

An occluder may be used to block flow through the contralateral iliac artery. Refer to *Section 10.2.16* and the Talent endoluminal occluder system Instructions for Use.

10.2.16 Occluder deployment

The Talent occluder system (packaged separately) is available to be used with the Endurant II/IIIs stent graft system and is typically used in conjunction with the AUI stent graft component. The Talent occluder system is closed at both ends to stop retrograde blood flow into the aneurysm sac.

For details on use and implantation of the Talent occluder system, refer to the appropriate sections of the Talent endoluminal occluder system Instructions for Use.

10.2.17 Smoothing stent graft fabric and modeling stent graft

The Reliant stent graft balloon catheter (packaged separately) can be used to assist in stent graft implantation by modeling the covered portion of the stent graft and removing wrinkles and folds from the graft material, as needed. Use the balloon catheter to model the proximal and distal seal zones as well as any overlapping connection (or junction) areas between the stent graft components. Suboptimal expansion of the self-expanding stent graft components may also be improved by use of the balloon catheter. Refer to the Reliant stent graft balloon catheter Instructions for Use for specific instructions.

Note: The Reliant stent graft balloon catheter is recommended for use with the Endurant II/IIIs stent graft system. Data is not available for the use of other balloon catheters in remodeling stent grafts.

Note: Use a kissing balloon technique to remodel the proximal end of the Endurant II or Endurant IIIs bifurcated or AUI stent graft (Parallel Graft technique only). Before remodeling with the Reliant stent graft balloon, deploy an appropriately sized balloon into each balloon expandable covered stent for dilation and support. Maintain support of the balloon expandable covered stents until after the Reliant balloon has been deflated.

Note: For a parallel graft technique, carefully monitor the position of the stent graft during balloon remodeling as proximal migration of the stent graft may occur.

Note: When using the Heli-FX EndoAnchor system, Medtronic recommends that the implantation of EndoAnchor implants be done after the aortic endograft has been placed and any balloon remodeling of the infrarenal seal zone of the stent graft system has been completed. Exercise care in balloon remodeling of the stent graft system to avoid moving the main body endograft from its intended implant location.

Caution: Over inflation of the balloon can cause graft tears or vessel dissection or rupture.

Warning: When expanding a vascular prosthesis, there is an increased risk of vessel injury or rupture, and possible patient death, if the balloon's proximal and distal radiopaque markers are not completely within the covered (graft fabric) portion of the prosthesis.

Warning: Do not use the Reliant stent graft balloon catheter in the treatment of dissections.

10.2.18 Additional fixation and sealing with the Heli-FX EndoAnchor system (Infrarenal neck \geq 4 mm and $<$ 10 mm in length)

An infrarenal EVAR procedure using the Heli-FX EndoAnchor system is performed at the discretion of the implanting physician. Medtronic recommends that the implantation of EndoAnchor implants be done after the aortic endograft has been placed and any balloon remodeling of the infrarenal seal zone of the stent graft system has been completed.

The Heli-FX EndoAnchor system is required when treating short infrarenal necks (≥ 4 mm and < 10 mm in length). For bifurcated endografts, the recommended minimum number of EndoAnchor implants is based on native vessel diameter and is independent of the amount of endograft oversizing. For aortic neck diameters ≤ 29 mm, the recommended minimum number of EndoAnchor implants is 4. For aortic neck diameters 30 - 32 mm, the recommended minimum number of EndoAnchor implants is 6.

1. Prepare the system, per the Heli-FX EndoAnchor system instructions for use.

Note: The access vessel diameter and morphology should be compatible for use with the Heli-FX EndoAnchor system.

2. Implant the EndoAnchor implant, per the Heli-FX EndoAnchor system instructions for use.

Note: Stability of the stent graft in short (≥ 4 mm and < 10 mm) infrarenal necks is augmented by the EndoAnchor implants. Ensure successful deployment of the recommended minimum number of EndoAnchor implants. Where the number of successfully deployed EndoAnchor implants is below the minimum recommended, there may be greater risk of postoperative Type 1a endoleak or migration.

Note: Always use fluoroscopy for guidance, delivery, and observation of any Heli-FX system components within the vasculature.

Note: The EndoAnchor implant should be implanted only into areas of aortic tissue that are free of or have insignificant calcified plaque or thrombus, or where such pathology is diffuse and less than 2 mm in thickness. It should not be used in patients with bleeding diathesis or with a known hypersensitivity or allergy to the device materials.

3. Remove the system, per the Heli-FX EndoAnchor system instructions for use.

10.2.19 Verify placement and seal

1. At the completion of the procedure, perform angiography to assess the stent graft for proximal and distal endoleaks, to verify position of the implanted stent graft in relation to the aneurysm and renal arteries, to evaluate implantation of EndoAnchor implants (if EndoAnchor implants are used), and if applicable, to assess patency of femoral-femoral bypass graft.

Note: When removing ancillary devices after a parallel graft technique, withdraw the balloons positioned in the renal artery into the sheath prior to withdrawing the sheath proximal to the suprarenal anchor pins.

2. Leaks at the attachment or connection sites should be treated using the balloon to remodel the stent graft against the vessel wall.

Note: For a parallel graft technique, a kissing balloon technique may be required to obtain adequate proximal sealing.

Note: For an infrarenal EVAR procedure using the Heli-FX EndoAnchor system, additional EndoAnchor implants may be used to augment the seal.

3. Major leaks that cannot be corrected by reballoning or adding EndoAnchor implants may be treated by adding aortic or iliac extension components to the previously placed stent graft components. See *Section 10.2.13*.

Caution: Any leak left untreated during the implantation procedure must be carefully monitored after implantation.

Note: For a parallel graft technique, small Type 1a endoleaks can be visualized intraprocedurally; many of these Type 1a endoleaks resolve spontaneously without intervention.

10.2.20 Close entry sites

1. Remove ancillary devices prior to repairing the entry site.
2. Repair the entry site with standard closure technique.

11 Bail-out techniques

In the unlikely event of delivery system failure, the following bail-out techniques may be used.

11.1 Screw gear handle disassembly

If partial stent graft deployment due to graft cover severance occurs, the screw gear handle disassembly technique may permit successful deployment of the stent graft.

1. Pull back the trigger and fully retract the slider.
2. Stabilize the delivery system.
3. Insert the tips of a pair of hemostats into each of the screw gear handle disassembly ports on the front grip.
4. Disengage the front grip from the screw gear by pressing the tips of the hemostats into the handle disassembly ports, simultaneously advancing the front grip away from the screw gear.
5. Advance the front grip until it fully clears the screw gear.
6. Separate the screw gear halves in order to identify the location of graft cover severance.
7. Manually retract the graft cover with your fingers or with hemostats until the stent graft is fully deployed.
8. Follow the instructions for tip capture deployment and delivery system removal.

11.2 Ballooning

If the captured proximal tip of the suprarenal stent cannot be deployed and the back-end wheel section still works, the ballooning technique may permit successful deployment of the suprarenal stent.

1. Use a compliant or semicompliant balloon (Reliant balloon recommended).
2. Insert the balloon and move it to the bifurcated configuration's aortic section.
3. Inflate the balloon inside the stent graft to vessel size to stabilize stent graft.
4. Follow the instructions for tip capture deployment and delivery system removal.

11.3 Back-end handle disassembly

If no or partial deployment of the proximal end of the suprarenal stent occurs due to back-end wheel failure, the back-end handle disassembly technique may permit the successful deployment of the suprarenal stent.

1. Use hemostats to depress the exposed tabs to disassemble the back-end wheel.
2. Insert the tips of hemostats into each of the rear handle disassembly ports.
3. Disengage the rear handle by pressing the tips of the hemostats into the handle disassembly ports simultaneously retracting the rear handle from the delivery system.
4. Stabilize the delivery system.
5. Manually push the exposed tabs of the back-end T-tube to release the suprarenal stent from the spindle.
6. Manually pull back the exposed tabs of the back-end T-tube to recapture the tapered tip after deployment.
7. Follow the instructions for delivery system removal.
8. Hold the exposed tabs of the back-end T-tube so that it remains retracted and the tapered tip remains recaptured during delivery system removal.

11.4 Snare the tapered tip

If the back-end handle disassembly technique is unsuccessful due to an excessively high deployment force, the snare the tapered tip technique may permit successful deployment of the suprarenal stent.

1. Use a snare device.
2. Advance the snare device to the delivery system tapered tip section through upper torso access (such as, brachial).
3. Utilize fluoroscopy to snare the edge of the delivery system tapered tip.
4. Stabilize the delivery system, especially the back-end section.
5. Pull the snare device to separate the suprarenal stent from the tip capture.
6. Manually pull back the back-end T-tube to recapture the tapered tip after deployment.
7. Follow the instructions for delivery system removal.
8. Ensure that the back-end T-tube remains retracted and the tapered tip recaptured during delivery system removal.

12 Follow-up imaging recommendations

12.1 General

Current imaging of stent graft patients includes abdominal X-ray and CT, with and without contrast medium. Alternative imaging modalities such as magnetic resonance imaging should be used in patients with impaired renal function or intolerance to contrast media.

Imaging should be decided based upon the physician's clinical assessment of the patient preimplantation and postimplantation of the stent graft. After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth, or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including 1) abdominal radiographs to examine device integrity (stent fracture, separation between bifurcated device and proximal cuffs or limb extensions, if applicable), and 2) contrast and noncontrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity, and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

Medtronic recommends consideration of an enhanced follow-up schedule for patients with short infrarenal necks (≥ 4 mm and < 10 mm in length), where the minimum recommended number of EndoAnchor implants were not used.

Additional imaging guidelines for Heli-FX EndoAnchor implants can be found in the instructions for use provided with the device.

12.2 X-ray

Abdominal X-rays should be used to assess the presence of stent graft fracture. Four-view kidney, ureter, bladder (KUB) X-rays should be taken. Posterior/anterior (PA) and lateral images are recommended for visualization of the stent graft. Ensure the entire device is captured on images for device assessment.

12.3 CT with contrast

Contrast-enhanced CT should be used to assess stent graft fixation, deformation, apposition to the vessel wall at proximal and distal fixation sites, stent graft migration, stent graft patency, AAA size, occlusion of branch vessels, and endoleak (including source and type if present).

A precontrast scan of 5 mm thick slices is suggested to determine if there are calcifications or areas where metal artifacts may be misinterpreted as endoleak. An arterial phase with <3 mm slice thickness and overlapping images with coverage from the celiac artery to the external iliac is recommended. In aneurysms that are not shrinking and have no apparent endoleak or fixation problems, a late venous phase scan may be performed. The venous phase scan may also be performed with thicker collimation (5 mm). It is recommended that the source data set be archived in case specialized evaluation is needed later (volume measurements, 3-dimensional reconstruction, or computer-aided measurement software). If the aneurysm is not shrinking by more than 5 mm within the first year, volume measurements may be obtained as a more sensitive indicator of AAA size using 3-dimensional software. Patients who are allergic to contrast should be premedicated 12 hours to 24 hours prior to receiving it.

12.4 Noncontrast CT

For patients with impaired renal function or those who are allergic to contrast medium, a spiral CT without contrast may be considered to assess stent graft fixation, deformation, apposition to the vessel wall at proximal and distal fixation sites, stent graft migration, occlusion of vessels, and size of the AAA diameter and volume measurements.

12.5 Duplex ultrasound

For patients with impaired renal function or those who are allergic to contrast medium, a color-duplex ultrasound may be considered to assess size of AAA diameter, endoleaks, and stent graft occlusion and stenosis.

12.6 MRI or MRA

Patients with impaired renal function, such as renal insufficiency, may also be considered for magnetic resonance imaging or angiography (MRI, MRA) in facilities that have expertise in this area. Artifact may occur related to the stent, and care should be used to insure adequate imaging of the outer aneurysm wall to assess AAA size. Volume measurement may be helpful if the aneurysm is not clearly shrinking. If there are concerns regarding imaging of calcified areas, fixation sites, or the outer wall of the aneurysm sac, adjunctive CT without contrast may be needed.

12.7 Imaging tests

Refer to *Table 9* for the recommended follow-up imaging schedule after stent graft implant.

Table 9. Imaging recommendations

Imaging test	12-month and annual follow-up for up to 5 years
CT ^a or MRA ^b with contrast	X
Abdominal X-rays (4-View, KUB)	X

^a CT evaluation may include "3-phase technique," volume studies, 3-D reconstruction, or computer-aided measurements.

^b MRA may be used in patients with impaired renal failure or contrast intolerance.

12.8 Supplemental imaging

Note: Additional radiological imaging may be necessary to further evaluate the stent graft in situ based on findings revealed by one of the surveillance programs. The following recommendations may be considered.

- If there is evidence of poor or irregular position of the stent graft, severe angulation, kinking, or migration of the stent graft on abdominal X-rays, a spiral CT should be performed to assess aneurysm size and the presence or absence of an endoleak.
- If a new endoleak or increase in AAA size is observed by spiral CT, adjunctive studies such as 3-D reconstruction or angiographic assessment of the stent graft and native vasculature may be helpful in further evaluating any changes of the stent graft or aneurysm.
- Spiral CT without contrast, MRI or MRA may be considered in select patients who cannot tolerate contrast media or who have renal function impairment. For centers with appropriate expertise, gadolinium or CO₂ angiography may be considered in patients with renal function impairment requiring angiographic assessment.

13 Additional surveillance and treatment

Additional endovascular repair or open surgical aneurysm repair should be considered for patients with evidence of suboptimal stent graft fixation, proximal endoleak, distal endoleak, junction endoleak, unknown origin of persistent perigraft flow, or increase in AAA size > 5 mm.

14 Disclaimer of warranty

Note: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.