












For *in Vitro* Diagnostic Use


## DNA-Sorb-B

### DNA extraction kit

#### Key to symbols used

	List Number		Store at 2-25°C
	For <i>in Vitro</i> Diagnostic Use		Caution!
	Lot Number		Version
	Expiration Date		Consult instructions for use
	Contains reagents		Manufacturer

#### NAME

DNA-Sorb-B

#### INTENDED USE

The **DNA-Sorb-B** nucleic acid extraction kit is intended for the isolation and purification of DNA from plasma, serum, whole blood, liquor, amniotic liquid, tissue, urine, feces, bronco aspirates and other biological materials.

#### MATERIALS PROVIDED


- Lysis Solution, 15 ml;
- Washing Solution 1, 15 ml;
- Washing Solution 2, 50 ml;
- Sorbent, 1,25 ml;
- DNA-eluent, 5 ml.

Contains reagents for 50 tests.

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Biological cabinet
- Desktop microcentrifuge for “eppendorf” type tubes (RCF max. 16,000 x g); Eppendorf 5415D or equivalent
- 60°C ± 2°C dry heat block
- Vortex mixer
- Pipettors (capacity 5-40 µl; 40-200 µl; 200-1000 µl) with aerosol barrier
- 1,5 ml polypropylene sterile tubes (Sarstedt, QSP, Eppendorf)
- Disposable gloves, powderless
- Tube racks
- Freezer
- Refrigerator

#### WARNINGS AND PRECAUTIONS

-  Lysis Solution contains guanidine thiocyanate. Guanidine thiocyanate is harmful if inhaled, or comes in contact with skin or if swallowed. Contact with acid releases toxic gas. (Xn; R: 20/21/22-36/37/38; S: 36/37/39).
- Wear disposable gloves, laboratory coats and eye protection when handling specimens and reagents. Thoroughly wash hands afterward.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Specimens should be considered potentially infectious and handled in biological cabinet in accordance with Biosafety Level 2 or other appropriate biosafety practices.
- Clean and disinfect all spills of specimens or reagents using a disinfectant such as 0,5% sodium hypochlorite, or other suitable disinfectant.
- Avoid contact of specimens and reagents with the skin, eyes and mucous membranes. If these solutions come into contact, rinse immediately with water and seek medical advice immediately.
- Material Safety Data Sheets (MSDS) are available on request.
- Use of this product should be limited to personnel trained in the techniques of DNA amplification.
- Workflow in the laboratory must proceed in a uni-directional manner, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where you performed previous step.

## STORAGE AND SHIPPING

1. **DNA-Sorb-B** can be stored at 2-25°C storage temperature. Reagent will crystallize upon storage at 2-8°C.
2. **DNA-Sorb-B** can be shipped at room temperature.

## STABILITY

**DNA-Sorb-B** is stable up to the expiration date indicated on the kit label.

## SPECIMEN COLLECTION AND CONSERVATION

**DNA-Sorb-B** Kit can isolate DNA from:

- *whole blood* collected in either ACD or EDTA tubes;
- *serum* collected blood in Serum Separator tubes;
- *plasma* collected blood in ACD or EDTA tubes;
- *bone marrow aspirate* collected in EDTA tube;
- *liquor* stored in "Eppendorf" tube;
- *lacrimal liquid* stored in "Eppendorf" tube;
- *amniotic liquid* stored in "Eppendorf" tube;
- *sinovial liquid* stored in "Eppendorf" tube;
- *peritoneal and pleuric versament* stored in "Eppendorf" tube;
- *tissue* homogenized with mechanical homogenizer and dissolved in PBS sterile;
- *urine (sediment)*;
- *prostatic liquid* stored in "Eppendorf" tube;
- *seminal liquid*: transfer about 30 µl of seminal liquid to a polypropylene tube (1,5 ml) and add 70 µl of sterile saline solution;
- *sputum*: add 1 volume of sputum to 5 volumes of "Mucolisin" (reagent not provided). Vortex well and incubate at room temperature for 30 min, mix batchly. Transfer 1,0 ml of clinical material to a sterile polypropylene tube (1,5 ml) and centrifuge at 10000g/min for 10 min. Discard the supernatant and leave about 100 µl of solution for DNA extraction;
- *bronco aspirate*: transfer 1,0 ml to a polypropylene tube (1,5 ml) and centrifuge at 10000g/min for 10 min. Discard the supernatant and leave about 100 µl of solution for DNA extraction;
- *feces*:
  - ❖ *Bactery*:
    - Prepare required quantity of 1,5 ml polypropylene tubes with 1,0 ml of Saline Solution. Add to each tube 0,1 g of feces. Vortex to get a homogeneous suspension and centrifuge for 5 min to 7000-12000g and using a micropipette with a plugged aerosol barrier tip transfer in a new tube 0,1 ml of the bacterial fraction (white-yellowish line between the sediment and the supernatant). Add 0,8 ml of sterile Saline Solution.
    - Vortex vigorously and centrifuge for 5 min at 7000-12000g. Remove and discard the supernatant.
    - Resuspend the pellet in 0,3 ml of Saline Solution
  - ❖ *Virus*:
    - prepare required quantity of 1,5 ml polypropylene tubes with 1,0 ml of Saline Solution. Add to each tube 0,1 g of feces. Vortex vigorously to get the homogeneous suspension. Centrifuge for 5 min at 7-12000g and use the supernatant for the extraction of the DNA.

It is recommended to process samples immediately after collection. Store samples at 2-8 °C for no longer than 24 hours, or freeze at -20/80°C.

## SPECIMEN AND REAGENT PREPARATION

1. **Lysis Solution** and **Washing Solution** (in case of their storage at +2-8°C) should be warmed up to 60-65°C until disappearance of ice crystals. Prepare required quantity of 1.5 ml polypropylene tubes including one tube for **Negative Control of Extraction**.
2. Add to each tube **10 µl of Internal Control** (if provided with the amplification kit) and **300 µl of Lysis Solution**.
3. Add **100 µl of Samples** to the appropriate tube.
4. Prepare Controls as follows:
  - add **100 µl of C-** (**Neg Control** provided with the amplification kit) to the tube labeled *Cneg*.
5. Vortex the tubes and incubate for 5 min at 65°C. Centrifuge for 7-10 sec. If the sample is not completely dissolved it is recommended to re-centrifuge the tube for 5 min at a maximum speed (12000-16000 g.) and transfer the supernatant into a new tube for DNA extraction.
6. Vortex vigorously **Sorbent** and add **20 µl** to each tube.
7. Vortex for 5-7 sec and incubate all tubes for 3 min at room temperature. Repeat this step.
8. Centrifuge all tubes for 30 sec at 5000g and using a micropipette with a plugged aerosol barrier tip, carefully remove and discard supernatant from each tube without disturbing the pellet. Change tips between the tubes.
9. Add **300 µl of Washing Solution 1** to each tube. Vortex vigorously and centrifuge for 30 sec at 8000g. Remove and discard supernatant from each tube.
10. Add **500 µl of Washing Solution 2** to each tube. Vortex vigorously and centrifuge for 30 sec at 8000g. Remove and discard supernatant from each tube.
11. Repeat step 10 and incubate all tubes with open cap for 5 min at 65°C.
12. Resuspend the pellet in **50 µl of DNA-eluent**. Incubate for 5 min at 65°C and vortex periodically.
13. Centrifuge the tubes for 1 min at 12000g.
14. The supernatant contains DNA ready for amplification. If amplification is not performed the same day of extraction, the processed samples can be stored at 2-8°C for at maximum period of 5 days or frozen at -20°/-80°C.



Sacace Biotechnologies Srl  
44 Scalabrini str., 22100 Como, Italy