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Serial N°:

8. MAINTENANCE

8.1 General recommendations on maintenance

The machine was built only and exclusively for drying endoscopes that have already been handled/sterilized. For this reason it is necessary to provide some useful instructions for the operators who will be performing maintenance on it.

The maintenance technicians, in normal operating conditions, are not subject to risks if they work safely using suitable means of protection.

In order to work safely the maintenance technician must:

- Carefully comply with the instructions set forth in this manual.
- Use safety devices appropriately and with care, as well as group and individual safety gear provided in the workplace.
- Be particularly careful and aware in case of maintenance for the replacement of potentially contaminated bacteriological filters.

Maintenance operations for the machine described in this manual can be divided into "**Routine Maintenance**" and "**Special Maintenance**".

GENERAL GUIDELINES:

MACHINE STATUS

The machine must not be powered electrically and the dedicated safety device must be in the OFF position. The person performing the task must ensure that there is no-one around the machine during this operation.

SAFETY SYSTEMS TO BE ADOPTED

The operation must be carried out in compliance with behavioral rules, in compliance with standards regarding contact with parts of the machine which may be contaminated by pathogenic materials and with use of suitable PPEs.

8.2 Procedure for routine maintenance work

Routine maintenance includes all operations aimed at keeping various parts of the machine clean and functional. They must be performed on a regular basis or when considered necessary due to incorrect performance of drying cycle.

Since these are simple cleaning operations, they are normally performed by the machine operator on his own liability.

8.3 Table of routine maintenance tasks

The following table shows the various routine maintenance tasks, their frequency, who is to perform them and the reference to the specific intervention form.

Each single task is more fully explained in the single reference forms.

Even if the water supply is relatively soft, the high temperature can cause the formation of residues which may create problems with the heating element, compromising the correct wash cycle and the reaching of the disinfection temperature.

For these reasons it is advisable to carry out regular cleaning as described below.

TABLE OF ROUTINE MAINTENANCE TASKS

Steelco		ED 200																REFERENCE
		Programmed maintenance scheme																
		Components	Step	months												Activity		
3	6			9	12	15	18	21	24									
		make every														Replace.	M2	
	Drying F5 filter	make every 4500 hours														Replace.	M2	
	HEPA filter	make every 4500 hours																
	Temperature probes	make every					x								x	During periodic validation, check the sensor status.	M1	
	Safety thermostat	make every					x								x	Verify the sensor.	M1	
	Channels connection pipe	make every					x								x	Replace.		
	Instrument and machine connectors	make every		x			x				x				x	Check the state and if necessary replace.		
	Blower	every week														Check for free rotation.		
	Door gasket	make every		x			x				x				x	Check the state of the seal and replace it every 4 years.		
	Compressor	make every		x			x				x				x	Check for air leaks from the gasket.		
	Compressor air filter unit 0.2 µm cod. ED100085	make every		x			x				x				x	Check the state of the compressor air filter unit. Replace it whenever necessary and at least once a year.		
	Drying heating element	none														Operation is checked by the control system.		
	Pressure switches	make every					x								x	Operation is checked by the control system.		
	Flowmeter	make every					x								x	Check the state of the pipes, the seal as well as the internal filters.		
	White cylindrical component of the nozzle	make every											x			Replace.		
	2 nozzle O-ring, fixed to the white cylindrical component	make every		x			x								x	Replace.		
	UV-C light	make every					x								x	Check the UV-C pipes. The operation is monitored by the control system.		

N.B.:

The time frames for execution of the maintenance programme may vary by +/- 15 days from the period indicated in the table.

N.B.:

Routine maintenance tasks must be performed at the intervals set forth in the table. It is however advisable to carry out single cleaning tasks anytime you feel they may be necessary.

Miele Group
Member

User manual

ENDOSCOPES WASHER DISINFECTOR

EW 2/1

EW 2/2

Serial N°:

CE 0051

7. PROGRAMMED MAINTENANCE

7.1 Operator maintenance

	ENDOSCOPES WASHER EW2								TIME	
	Programmed maintenance scheme									
Components	Step	months							Activity	
	make every.....	3	6	9	12	15	18	24		
Chamber filters	make every day								Take filters and clean.	10'
Leak test filter	make every.....		x						Replace.	1'
Chemical compartment	make every day								Check every day the cleaning of chemical compartment. In case of leakage call the service.	5'
Pre filter dryer F5 (optional)	make every 100 hours								Replace.	2'
Dryer HEPA filter (optional)	make every 300 hours								Replace.	1'
Water filters 0,45 micron	make every.....		x						Replace.	10'
Water filters 0,1 micron	make every.....		x						Replace.	10'
Channel filters 40 micron	make every.....		x						Replace.	10'
Temperature probes	make every.....				x				During periodic validation, check the sensor status.	5'
Safety thermostat	make every.....				x				Verify the sensor.	1'
Chemical/water flowmeter	make every.....		x						Check the impeller, the calibration and the presence of leakages.	30'

Chemical tank level sensor	make every.....		x						Check the functioning and visual control (calcium problems)	5'
Dosing pump connection pipe	make every.....		x						Check of any crashing, leakages or hardening.	2'
Washing arms or nozzles	every week								Check for free rotation. Open the cleaning caps and wash inside; check and in case clean the nozzles.	2'
Door gasket	make every...		x						Verify the gasket and replace after one year.	20'
Washing pumps	make every...				x				Check for water leakage from the arm seal and correct flow (50 to 60 impuls)	5'
Water heating element	make every...				x				Check for water leakage from the gasket and check ampere consumption (10 A per unit)	1'
Dryer heating elements (optional)	none								Operation is checked by the control system and check ampere consumption	5'
Macchine nozzles of basket connection	make every...		x						Check the status and replace the oring.	2'
Drain pump	make every...				x				Check for any leakage, if necessary, remove and clean the membrane seat.	1'
General test	make every...								Operation is checked by the control system. Check for any leakage and the electrical connections.	15'
Unloading water pipe	make every...				x				Check the situation of the pipe and the seal.	5'
Loading water pipes	make every...				x				Check the situation of the pipe and the seal.	5'
Teflon shim inside the chamber	every month								Check the correct position. Fix the screws.	10'

Compressed air	make every...		x						Check the calibration of compressed air for channels purging and pneumatic valves.	2'
Basket/instrument connection pipes	make every day								Check the status of the silicon pipe for the connection of basket/instrument. If necessary replace them.	1'
Machine leakage	make every...		x						Check the leakage during functioning.	5'
Check the security on the doors	make every		x						Control the correct work of the motors and switch	5'
Silicon pipe basket/instrument	make every...				x				Replace.	20'
Luer connectors male/female basket	make every...						x		Replace.	30'
Endoscope oring connectors	make every day								Check the status of the oring. If necessary replace them.	1'
Basket	make every		X						Check the status of screws and if necessary, fix it.	5'