CERTIFIED COMPANY UNI EN ISO 9001 & UNI EN ISO 13485

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BOWIE-DICK TEST

Bowie-Dick type test is a method to control periodically the correct autoclave function. This test is used to verify air removal in a saturated steam autoclave with fractional vacuum during cycles for porous loads at 134°C for 3.5 minutes and 121°C for 15 minutes. Following the sterilization cycle, the color change of the stripes indicates that the residual air has been sufficiently removed to allow the complete penetration of the steam inside the indicator sheet and the effectiveness of the sterilization process.

Laminated indicator sheet with pink to black colour tone. Each package is contained in a practical box to be inserted directly into the autoclave.

The Bowie & Dick test pack is classified as Type 2 in accordance with ISO 11140-1 and ISO 11140-4.

Equivalent to the test pack indicated in EN 285.

The product is latex-free and does not contain toxic substances or heavy metals.

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Bowie-Dick test is used to verify air removal in a saturated steam autoclave with fractional vacuum during cycles for porous loads at 134°C for 3.5 minutes and 121°C for 15 minutes

For use in professional test laboratory only

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The indicator sheet will change from pink to black when exposed to steam sterilization cycle with acceptable vacuum cycles.

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Steam sterilization cycle functioning at 134° C for 3.5 minutes and 121°C for 15 minutes.

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Bowie & Dick Test Pack must be used every day, before the first sterilization cycle.

- *Run a shortened cycle to bring the sterilization chamber to use temperature.*
- *Place the Bowie & Dick Test Pack inside the empty sterilization chamber and on the bottom shelf above the drain. If the package is positioned correctly the writing "THIS SIDE UP" must be visible.*
- *Perform a steam sterilization cycle at 134°C for 3.5 minutes or at 121°C for 15 minutes*

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The indicator sheet is not changed in the center or is not changed completely, indicating that the autoclave is not working correctly and therefore technical assistance is required.

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The indicator sheet has gray stains or rings and sometimes a lifting of the protective film, highlighting the presence of damp steam. The test is considered successful.

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The Bowie & Dick Test could be hot, we recommend handling it carefully to avoid burns.

The indicator does not release toxic substances that could cause health risks or compromise the properties of the product.

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Store the product at temperatures 10 - 30 °C and relative humidity 30 - 60%.

Keep in a cool, dry place, away from heat sources and direct sunlight.

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Classified as Type 2 in accordance with EN ISO 1140-1 and EN ISO 11140-4.

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The product must be used within 4 years from date of manufacture.

IMBALLO / PACKING

Quantità (pz): Quantity (pcs):	20	Confezione interna (pz): Internal packing (pcs):	20 pacchi test 20 test box	QUANTITÀ MINIMA VENDIBILE / MINIMUM SALEABLE QUANTITY
Misura esterna scatola (cm): External box dimensions (cm):	29 x 17,5 x 17	Peso (Kg): Weight (Kg):	2,8	Volume (m ³): Volume (m ³):
				0,009

APTACA S.p.A.

Regione Monforte, 30 - Canelli (Asti) ITALY
TEL. +39 0141/83.50.75 - FAX +39 0141/83.52.92
E-mail: info@aptaca.com - http://www.apataca.com

Certified Company
UNI EN ISO 9001 & UNI CEI EN ISO 13485

Certificate of conformity and quality assurance

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Bowie-Dick test Dim. 130 x 130 x 20 mm

Test di Bowie-Dick Dim. 130 x 130 x 20 mm

We hereby declare that the Device in object is realized and tested in compliance with Quality procedures by controls effected with specific methodologies. Aptaca S.p.A. have Quality System certified in accordance with UNI EN ISO 9001 and UNI CEI EN ISO 13485 Standards. The Device meets the basic requirements of security and conformity provided by current Regulations and Directives.

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The Device, in its integral package and under appropriate storage conditions, has a expiry date of 2 years from manufacturing date.

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Data del certificato

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23/07/2025

Duilio BUONO
Quality & Regulatory Affairs Manager