INSTRUCTIONS FOR USE

Occlusion Balloon System

TABLE OF CONTENTS

Occlusion Balloon System

English	3
Français.	6
Deutsch.	10
Italiano	13
Español	17
Svenska.	20
Nederlands.	24
Português	27
Suomi.	31
Dansk	34
Ελληνικά	38
Česky	41
Magyar	45
Русский.	48
Polski	52
Türkçe.	55
Norsk	59
Slovenčina	62
Română.	66
Български	69
한국어.	73
عربي	79
Symbol Glossary	80

English

Instructions for Use

EN

Occlusion Balloon System

CAUTION

- Federal (USA) law restricts this device to use by or on the order of a physician.
- This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional radiologic procedures.
- Do not use if pouch is open or damaged.

CONTENTS

One (1) Occlusion Balloon Catheter and one (1) 0.010" Hydrophilic Guidewire.

DESCRIPTION

The Occlusion Balloon System is a single lumen balloon catheter that requires the insertion of the 0.010" guidewire to occlude the central lumen to allow inflation of the balloon. When the distal 10 cm platinum coil tip of the guidewire is advanced to or past the catheter tip, it occludes the inflation holes allowing the balloon to inflate through catheter sideholes.

INDICATIONS FOR USE

The Occlusion Balloon Catheters are indicated for use in blood vessels of the peripheral and neuro vasculature where temporary occlusion is desired. These catheters offer a vessel selective technique of temporary vascular occlusion, which is useful in selectively stopping or controlling blood flow; the Occlusion Balloon Catheters may also be used in balloon-assisted embolization of intracranial aneurysms and treatment of vasospasm.

CONTRAINDICATIONS

- Not intended for embolectomy and angioplasty procedures
- Not intended for use in coronary vessels
- Not intended for use as an infusion catheter
- Not intended for pediatric or neonatal use

WARNINGS

- Do not have a continuous flush in place while guidewire is advanced to the catheter tip in order to prevent inadvertent balloon inflation.
- Do not exceed the maximum recommended inflation volume as balloon rupture may occur.
- Verify the size of vessel under fluoroscopy, ensuring that when inflated, the balloon does not exceed size of vessel.
- Ensure accurate vessel dimensions are obtained prior to any Balloon-Assisted coiling (BAC) procedures.
- Do not maintain inflation for prolonged periods of time, as long-term occlusion can lead to increased complications.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Do not withdraw the guidewire to deflate the balloon, instead pull negative on the syringe to deflate.
- Do not steam shape the catheter as it may compromise the integrity of the balloon material.
- Do not use the catheter for subselective angiography as inadvertent balloon inflation could result.
- Do not advance the catheter against resistance until the source of the resistance is identified under fluoroscopy.
- Use only the recommended balloon inflation medium of 50:50 solution by volume of 60% contrast medium and normal, sterile saline.
- Do not use air or any gaseous medium to inflate the balloon.
- Avoid touching the balloon material as contrast residue can stick to gloved hands.
- Ensure the guidewire can be securely tightened in RHV. Inadvertent proximal guidewire movement can cause blood to enter the balloon lumen which may inhibit balloon deflation.

STORAGE

This device should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

COMPLICATIONS

Possible complications include, but are not limited, to the following:

- Vessel Rupture
- Death
- Stroke or Infarction
- Ischemia
- Perforation of Vessel or Arterial Wall
- Infection
- Vasospasm
- Thrombus Formation and Release, Increasing with Prolonged Insertion

COMPATIBILITY


The maximum outer diameter of the Occlusion Balloon System are shown in the table below, use appropriately sized guiding or angiography catheters.

Model Number	Maximum Outer Diameter
104-4112	0.040"
104-4113	0.040"
104-4127	0.040"
104-4132	0.040"
104-4153	0.040"
104-4310	0.040"
104-4315	0.040"
104-4370	0.040"
104-4415	0.040"
104-4420	0.040"
104-4470	0.040"
104-4515	0.040"
104-4520	0.040"
104-4530	0.040"
104-4715	0.046"
104-4770	0.046"

The Balloon System is dependent upon the 0.010" guidewire for balloon inflation.

PREPARATION OF SYSTEM


Figure 1: Proper Set up of System

- 1. 3-Way Stopcock
- 2. 3 cc/ml Syringe
- 3. RHV
- 4. Catheter Hub
- 5. 1 cc/ml Syringe

- Always hydrate wire first, this will allow proper hydration time (30 seconds) while prepping the balloon.
- Attach RHV and syringes to balloon. Use Cadence syringe with 3 way stopcock and 3cc syringe.
- Purge RHV assembly and balloon lumen with 50:50 contrast-saline solution (remove as much air as possible).
- Using wire introducer, advance guide wire through RHV almost to tip of balloon catheter.
- Flush balloon lumen again to ensure air has been removed.
- Advance guidewire past tip.
- Inflate balloon (will stick on first inflation).
- Check for bubbles.
- If bubbles exist, withdraw wire inside catheter and repeat Steps 4 through 7.
- Keep balloon tip in sterile water or saline until used with patient. This will keep wire hydrated and moving free.

Figure 2: Guidewire Position for Balloon Inflation






1. Balloon Inflation Seal



DIRECTIONS FOR USE

- Advance the System through the proximal valve of the introducer sheath or the RHV attached to the hub of a guiding catheter and into the vasculature.
- Advance a portion of the distal radiopaque guidewire past the distal catheter radiopaque marker and track to selected vessel.
- To inflate the balloon, use the 1 cc syringe to gently infuse the desired volume of recommended contrast solution. Balloon cannot inflate unless a portion of the distal 10 cm of the guidewire tip is occluding the distal inflation holes of the catheter. **Do not exceed the maximum recommended inflation volume as balloon rupture may occur.**
- Deflate the balloon by gently retracting the syringe plunger. Do not withdraw the guidewire from the tip of the catheter to deflate.
- Following inflation, if flush of the central lumen is desired, withdraw the guidewire tip proximal to the proximal balloon marker, and gently flush with recommended contrast solution. Watch the balloon for potential inflation during the flush procedure.
- Discontinue flush. Advance the guidewire beyond the catheter tip and track to the selected vessel for additional inflations.
- After multiple inflations, the catheter and guidewire should be removed and inspected. Verify the performance of the balloon and inspect for presence of clot at the tip or inflation holes.

Balloon Inflation Compliance Chart for Models

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TECHNICAL SUGGESTIONS

Guidewire Exchange Technique

If the tip of the guidewire needs reshaping, a guidewire exchange can be accomplished using the following technique:

1. Attach a 3 cc syringe filled with the recommended contrast solution for flush to the stopcock. (As shown in Figure 1).
2. Open the stopcock to the RHV and the catheter; loosen the RHV knob. While slowly withdrawing the guidewire, gently infuse the recommended contrast solution with the 3 cc syringe while observing the balloon for potential inflation.
3. Close off the RHV knob upon complete removal of the guidewire.
4. To insert a new guidewire, open the RHV knob and insert a guidewire introducer into the hub. Advance the guidewire into the proximal shaft leaving the guidewire tip proximal to the proximal balloon marker, remove the wire introducer, and gently flush the RHV.
5. Close off the RHV knob and gently flush the catheter lumen with the 3 cc syringe prior to the guidewire passing the distal marker band.
6. Advance the guidewire tip out of the catheter tip and proceed as directed.

Balloon-Assisted Coiling (BAC) Technique

Literature has shown that Balloon-Assisted Coiling (BAC) can provide the following benefits in the endovascular treatment of aneurysms:

- Enable the placement of initial coils within wide-neck aneurysms;
- Increase the packing density of coils within an aneurysm

If Balloon-Assisted Coiling (BAC) is determined to be an appropriate option, the following procedural technique is recommended:

1. Guide the deflated balloon portion of the balloon catheter across the neck of the aneurysm.
2. Simultaneously access the aneurysm sac via the microcatheter for coil delivery.
3. Inflate the balloon.
4. Place the coils into the aneurysm sac via the microcatheter.
5. Deflate balloon.
6. Evaluate for potential migration of coils and desired packing density.
7. If repositioning or placement of additional coils is required, repeat Steps 3-6 until desired result is obtained.

BALLOON VISUALIZATION

- If the balloon is not visible during an inflation procedure, blood may have entered the catheter lumen. This may occur during extensive guidewire manipulation, or if the guidewire is extended beyond the catheter tip and a strong vacuum was pulled on the syringe to deflate the balloon.
- To clear the lumen of blood, deflate the balloon, withdraw the guidewire proximal to the proximal balloon marker band and perform a gentle flush procedure with the recommended contrast solution.

BALLOON DEFLATION

- The balloon should always be deflated with the guidewire remaining beyond the catheter tip in order to minimize aspiration of blood into the lumen.
- Do not withdraw the guidewire tip to deflate the balloon.
- Never aspirate when the guidewire is withdrawn proximal to the catheter tip.

HYPER BALLOON DEFLATION EXPECTATIONS

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Micro Therapeutics, Inc.