REPUBLICA

MOLDOVA

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Catalog Number 4485261

Pub. No. 4486415 Rev. A

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Product information

Name, intended use, and principle of the procedure

The Avian Influenza Virus RNA Test Kit is a highly sensitive, qualitative, one-step, real-time reverse transcription PCR (real-time RT-PCR) assay to detect Avian Influenza Virus (AIV) RNA isolated from poultry (chicken/turkey) oropharyngeal/tracheal swab samples.

AIV is an enveloped, negative-sense RNA virus of the genus *Influenzavirus A* and family *Orthomyxoviridae*. AIV subtypes are defined by the surface glycoproteins hemagglutinin and neuraminidase. Low pathogenic avian influenza (LPAI) strains exist in avian reservoir hosts and can be transmitted to poultry. H5 and H7 subtypes of LPAI strains are unique in their capability to undergo host adaptation and to evolve into highly pathogenic avian influenza (HPAI). HPAI viral infections arise *de novo* in poultry infected with LPAI H5 and H7 subtypes and become rapidly fatal due to an overwhelming systemic collapse. The VetMAX™-Gold AIV Detection Kit enables diagnosis of AIV in poultry populations.

The assay consists of a single-well/tube, real-time RT-PCR assay in which RNA is reverse-transcribed into cDNA; two viral matrix targets and one nucleoprotein target are amplified and detected in real time using fluorescent TaqMan® probes (hydrolysis probe chemistry). The assay detects sequences that are common to all AIV subtypes. The kit includes:

- Influenza Virus-Xeno™ RNA Control Mix: serves as a positive control for the real-time RT-PCR components and it is also used to set the cycle threshold (C_T) for evaluating test results.
- Xeno™ RNA Control: serves as an internal positive control for the RNA purification process and monitors for the presence of PCR inhibitors.
- Influenza Virus Primer Probe Mix for optimized multiplex real-time RT-PCR amplification of Xeno™ RNA Control and AIV RNA targets.

Limitations

- The kit is not intended for differentiating AIV subtypes.
- Handle samples as recommended in Table 1 to prevent degradation of any AIV RNA that is present.
- RNA extraction methods should yield RNA free of RT-PCR inhibitors, which can prevent amplification of target RNA.
- Follow "Good laboratory practices for PCR and RT-PCR" on page 5 to prevent false positive amplifications due to contamination of test samples with PCR products.

For Veterinary Use Only.

Kit contents and storage conditions

Reagents for 100 25- μ L real-time RT-PCR tests are supplied.

Component	Volume	Storage
2X Multiplex RT-PCR Buffer	1375 μ L	-30°C to -10°C
Multiplex RT-PCR Enzyme Mix	280 μ L	-30°C to -10°C
Influenza Virus Primer Probe Mix	110 μ L	-30°C to -10°C
Xeno™ RNA Control (10,000 copies/ μ L)	250 μ L	-30°C to -10°C
Influenza Virus-Xeno™ RNA Control Mix (1000 copies/ μ L)	80 μ L	-30°C to -10°C
Nuclease-free Water	1.75 mL	-30°C to +25°C

Required materials not supplied

Item	Source
Plates or tubes appropriate for the Applied Biosystems® 7500 Fast Real-Time PCR System (96-well)	<ul style="list-style-type: none">• MicroAmp® Optical 8-Cap Strip (Cat. no. 4323032), or equivalent• MicroAmp® Fast Optical 96-Well Reaction Plate with Barcode, 0.1-mL (Cat. nos. 4366932, 4346906), or equivalent.• MicroAmp® Optical Adhesive Film (Cat. nos. 4311971, 4360954), or equivalent• MicroAmp® Fast 8-Tube Strip, 0.1-mL (Cat. no. 4358293), or equivalent• Precision Plate Holder for 0.1-mL tube strips (Cat. no. 4403809), or equivalent
Nuclease-free pipettes and filtered pipette tips	Major laboratory supplier (MLS)
Nuclease-free reagent tubes for preparing master mixes	MLS
Real-time PCR thermal cycler	Applied Biosystems® 7500 Fast Real-Time PCR System (96-well), running SDS software v1.4
1X Phosphate Buffered Saline (PBS), pH 7.4	MLS
Viral Transport Media	MLS
2 ice buckets:	MLS
<ul style="list-style-type: none">• One for the PCR setup area where the master mix is prepared• One for the area where RNA may be present	

Isolate RNA from samples

Table 1 Sample handling recommendations

Step or process	Recommendation
Transport/storage of samples	Transport oropharyngeal/tracheal swab samples between 4°C and 25°C or in accordance with manufacturer's instructions.
Preparation of swab samples	<ol style="list-style-type: none"> Place one oropharyngeal/tracheal swab sample into a 1.5-mL tube or deep-well 96-well plate, then add 0.75 mL of Viral Transport Media. Vortex vigorously for 3 minutes, then pulse-spin to remove debris from the tube cap. Remove 50 µL of supernatant for RNA isolation.
Preparation of mock-purified samples (for use in extraction control PCRs)	Prepare duplicate mock-purified samples, using 1X PBS as the starting material. Process with the same RNA isolation method that is used for test samples.
Proposed RNA isolation method	MagMAX™-96 Viral RNA Isolation Kit (Cat. nos. AM1836, AMB1836-5) or an equivalent RNA isolation method.
Required modifications to the RNA isolation method	<ul style="list-style-type: none"> Add 2 µL of undiluted Xeno™ RNA Control per isolation to the lysis solution used for RNA isolation. Add carrier RNA to the lysis solution according to the manufacturer recommendations. Carrier RNA is provided in the MagMAX™-96 Viral RNA Isolation Kit (Cat. nos. AM1836, AMB1836-5)

Perform real-time RT-PCR

- Determine the quantity of reactions and thaw the reagents
 - On each plate, include the following control reactions (for step 4 of this procedure):
 - Positive control (prepare duplicate reactions); use 8 µL of the Influenza Virus-Xeno™ RNA Control Mix (1000 copies/µL).
 - No-template control (NTC) (prepare duplicate reactions); use Nuclease-free Water in place of sample RNA.
 - Plan the plate layout so that the wells containing NTCs are located as far as possible from positive controls and test samples to prevent accidental cross-contamination.
 - Thaw RT-PCR master mix reagents in one ice bucket and controls and samples in a separate ice bucket, gently vortex each tube to mix the contents thoroughly, then briefly centrifuge to collect the solution at the bottom of the tube. Keep the reagents on ice.

- Prepare the RT-PCR master mix on ice

Combine the following components for the number of reactions required plus 10% overage.

Component	Volume per reaction
2X Multiplex RT-PCR Buffer	12.5 µL
Multiplex RT-PCR Enzyme Mix	2.5 µL
Influenza Virus Primer Probe Mix	1.0 µL
Nuclease-free Water	1.0 µL
Total volume of RT-PCR master mix	17.0 µL

- Set up the RT-PCR reactions

- Dispense 17 µL of RT-PCR master mix to the appropriate wells of a PCR plate or PCR tubes on ice.
- Add the appropriate component for the reaction type, according to the following table:

Reaction type	Component	Volume per reaction
Test sample	Sample RNA	8.0 µL
NTC	Nuclease-free Water	8.0 µL
Positive control	Influenza Virus-Xeno™-RNA Control Mix (1000 copies/µL)	8.0 µL
Extraction control	Mock-purified sample	8.0 µL

- Seal each reaction vessel, mix, then centrifuge briefly to bring the contents to the bottom.

4 Set up and run the real-time PCR instrument

- a. Following the manufacturer's instructions, set up the run using the following parameters:
- Run mode: Standard 7500
 - Reaction volume: 25 μ L
 - ROX™ passive reference dye: Included in the RT-PCR Buffer
 - TaqMan® probe reporter dyes and quenchers:

Target	Reporter	Quencher
AIV RNA	FAM™ dye ^[1]	Eclipse® Q
Xeno™ RNA Control	VIC® dye ^[2]	Eclipse® Q

^[1] Absorbance maximum of 495 nm; emission maximum of 520 nm.

^[2] Absorbance maximum of 540 nm; emission maximum of 552 nm.

- b. Run the thermal cycler program and collect real-time amplification data during stage 3. Use the following thermal cycler settings:

Stage		Reps.	Temp.	Time
Reverse transcription	1	1	48°C	10 minutes
RT inactivation/initial denaturation	2	1	95°C	10 minutes
Amplification	3	40	95°C	15 seconds
			60°C	45 seconds

Data analysis

Refer to your real-time PCR instrument user guide for instructions on how to analyze your data, using the following method.

Table 2 Data analysis

Method	Details
Use the Control-Based Threshold setting for data analysis.	<ol style="list-style-type: none"> 1. Select Manual C_T. 2. Export ΔRn values for the positive control samples (Influenza Virus-Xeno™ RNA Control, 1000 copies/μL). 3. Average the FAM™ and VIC® values (separately) for the ΔRn at cycle 40 for all replicates of the positive control reaction. 4. Set the threshold for the AIV RNA reactions at 5% of the average maximum fluorescence value of the AIV RNA amplification signal in the positive control reactions. Example: If the average maximum fluorescence value for the AIV RNA target in the positive control reactions is 3.0, set the AIV RNA threshold at 0.15. 5. Repeat step 4 for the Xeno™ RNA Control target using a 5% threshold. Example: If the average maximum fluorescence value for the Xeno™ RNA target in the positive control reactions is 2.0, set the Xeno™ RNA threshold at 0.1.
Check the raw fluorescence data.	Verify that increased fluorescence seen in the normalized data is also evident without mathematical data processing.

Interpretation of test results

Table 3 Criteria for a valid real-time RT-PCR run

Reaction type	C _T value for AIV RNA	C _T value for Xeno™ RNA Control
Positive control	25–29	25–29
NTC	40 (undetermined) ^[1]	40 (undetermined) ^[1]
Extraction control	40 (undetermined) ^[1]	27.5–34

^[1] The run is invalid if the C_T values for either AIV or Xeno™ RNA Control targets in the NTC are <40, see "Troubleshooting" on page 4.

Table 4 Interpretation of sample test results

C _T value for AIV RNA	C _T value for Xeno™ RNA Control	Interpretation
<38	27.5–34 ^[1]	AIV-positive sample
40 [undetermined]	27.5–34	AIV-negative sample
38–40	27.5–34 ^[2]	Suspect result

^[1] If the C_T value falls outside of this range, see “Troubleshooting” on page 4.

^[2] If the C_T value falls outside of this range, see Table 5.

Table 5 Assessment of suspect results

Result	Action	
Suspect result: The sample AIV C _T value is 38–40.	Analyze suspect RNA samples for the presence/absence of RT-PCR inhibitors by calculating the Xeno™ RNA C _T shift: Xeno™ RNA C_T Shift = SS – XEC , where: SS = C _T of Xeno™ RNA in the suspect sample XEC = Average C _T of Xeno™ RNA in the extraction controls	
	Workflow A Xeno™ RNA C _T shift ≥1.5	Workflow B Xeno™ RNA C _T shift <1.5
	<ol style="list-style-type: none"> Repeat the real-time RT-PCR with 2 µL of the suspect RNA sample. [RT-PCR inhibitors may be present in the RNA.] If the AIV C_T value is: <ul style="list-style-type: none"> <38—The sample is AIV positive. No further testing is required. ≥38—Continue with steps 2 through 5 of this procedure. Dilute the original diagnostic sample 1:4. Repeat the RNA purification on triplicate aliquots of the diluted sample. Repeat the real-time RT-PCR with 8 µL of purified RNA from step 3. Determine the number of samples with a AIV C_T value <40: <ul style="list-style-type: none"> 0 of 3: AIV negative 1 of 3: Presumptive positive; confirm with secondary test method ≥2 of 3: AIV positive 	<ol style="list-style-type: none"> Repeat the RNA purification on triplicate aliquots of the original diagnostic sample. Repeat the real-time RT-PCR with 8 µL of purified RNA from step 1. Determine the number of samples with a AIV C_T value <40: <ul style="list-style-type: none"> 0 of 3: AIV negative 1 of 3: Presumptive positive; confirm with secondary test method ≥2 of 3: AIV positive

Confirmatory testing

All samples generating a positive or presumptive positive test result with the VetMAX™-Gold AIV Detection Kit should be submitted to a national laboratory (or equivalent authorized laboratory) for confirmatory testing.

Troubleshooting

Observation	Possible cause	Recommended action
Positive control reaction Influenza Virus-Xeno™Control RNA—no signal Xeno™ RNA Control—no signal	The Influenza Virus-Xeno™ RNA Control Mix was improperly handled, resulting in RNA degradation.	Use appropriate precautions against RNase contamination when handling the control RNAs. For example, wear clean gloves and use nuclease-free barrier pipette tips.
	The Multiplex RT-PCR Enzyme Mix was stored or handled improperly, and it lost activity.	Repeat the RT-PCR with fresh reagents.
	The thermal cycler was not properly set up.	Check the thermal cycler settings. See “Set up and run the real-time PCR instrument” on page 3.
	The RT-PCR master mix was prepared incorrectly.	Repeat the test with correctly prepared RT-PCR master mix.

Observation	Possible cause	Recommended action
NTC or extraction control reaction C _T value is <40	There was contamination during the RNA extraction or PCR.	<ul style="list-style-type: none"> Repeat the RNA isolation or real-time RT-PCR with fresh reagents and freshly decontaminated pipettes. Set up the real-time RT-PCR in an area separate from areas used for RNA isolation and PCR product analysis.
Test samples Xeno™ RNA Control—no or low signal AIV RNA—high signal	The Xeno™ RNA Control primers and probe are at limiting concentrations in the RT-PCR. High levels of AIV RNA in a sample can reduce amplification of Xeno™ RNA Control.	No or low signal from Xeno™ RNA Control is expected in a reaction that has a strong signal for AIV RNA.
Test samples Xeno™ RNA Control—no signal AIV RNA—no signal or inconclusive-range signal	Poor RNA recovery.	Check the C _T values of Xeno™ RNA Control in the mock-purified samples. C _T ≥ 38: indicates that Xeno™ RNA Control was omitted or that RNA recovery was poor. Repeat the RNA purification of the original diagnostic sample.
	The RNA samples contain inhibitors of RT-PCR.	See Table 5.

Explanation of symbols

The symbols present on the product label are explained in the following table.

	MANUFACTURER		USE BY
	CATALOG NUMBER		CONSULT INSTRUCTIONS FOR USE
	BATCH CODE		CAUTION, CONSULT ACCOMPANYING DOCUMENTS
	SERIAL NUMBER		UPPER AND LOWER LIMITS OF TEMPERATURE

Good laboratory practices for PCR and RT-PCR

When preparing samples for PCR or RT-PCR amplification:

- Wear clean gloves and a clean lab coat (not previously worn while handling amplified products or used during sample preparation).
- Change gloves whenever you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation and reaction setup.
 - Amplification and analysis of products.
- Do not bring amplified products into the reaction setup area.
- Open and close all sample tubes carefully. Avoid splashing or spraying samples.
- Keep reactions and components capped as much as possible.
- Use a positive-displacement pipettor or aerosol-resistant barrier pipette tips.
- Clean lab benches and equipment periodically with 10% bleach solution or DNAzap™ Solutions (Cat. no. AM9890).

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Veterinary Use Only.

Manufacturer's address: Life Technologies Corporation | 2130 Woodward Street | Austin, TX 78744

US Veterinary License No. 432.

The information in this guide is subject to change without notice.

DISCLAIMER

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

LIMITED USE LABEL LICENSE No. 460: PCR Veterinary Diagnostics; Notice to Purchaser: For veterinary diagnostic services including reporting the results of such services for a fee and for research purposes only. Human diagnostic uses require a separate license from Roche.

©2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Eclipse is a registered trademark of ELITech Holding B.V.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com

25 March 2014

VetMAX™ NDV Kit

TaqMan® real-time RT-PCR for the detection of paramyxovirus type 1

Catalog Number NDV

Doc. Part No. 100020401 Pub. No. MAN0008826 Rev. B.0

Technology	Species	Samples	Test type
Real-time RT-PCR (RNA) • Duplex assay • Xeno™ RNA Control	Poultry	Tracheal, cloacal swabs Feces Tissues	Individual

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

WARNING! POTENTIAL BIOHAZARD. Read the biological hazard safety information at this product's page at thermofisher.com. Wear appropriate protective eyewear, clothing, and gloves.

Product information

Name, intended use, and principle of the procedure

The Applied Biosystems™ VetMAX™ NDV Kit (Cat. No. NDV) enables real-time reverse transcription PCR (RT-PCR) detection of paramyxovirus type 1 extracted from swabs, feces, or tissues.

NDV is an enveloped, single-stranded RNA type 1 paramyxovirus responsible for Newcastle Disease, an infectious and contagious disease which affects domestic and wild birds.

The VetMAX™ NDV Kit provides assays and reagents required for single-well, real-time RT-PCR in which RNA is transcribed to cDNA, and NDV and Xeno™ RNA targets are amplified and detected using fluorescent TaqMan® probes (hydrolysis probe chemistry). It contains:

- 25X NDV Primer Probe Mix: primers and TaqMan® probes for optimized duplex real-time PCR amplification of NDV target and Xeno™ RNA target.
- 2X qRT-PCR Buffer: RT-PCR buffer.
- 25X qRT-PCR Enzyme Mix: RT-PCR enzyme.
- 25X NDV Control RNA: serves as an external positive control for the real-time RT-PCR components and it is also used to set the cycle threshold (C_t) for evaluating test results.
- Xeno™ RNA Control: the internal positive control is added to each sample and control at the lysis step of the RNA extraction procedure. It serves as an internal positive control for the RNA purification process and monitors for the presence of RT-PCR inhibitors.

Kit contents and storage conditions

Reagents for 100 25-µL real-time RT-PCR tests are supplied.

Component	Volume	Storage
25X NDV Primer Probe Mix	110 µL	-30°C to -10°C
2X qRT-PCR Buffer	1375 µL	-30°C to -10°C
25X qRT-PCR Enzyme Mix	110 µL	-30°C to -10°C
Nuclease-free Water	1750 µL	-30°C to +25°C
25X NDV Control RNA	15 µL	-25°C to -15°C
Xeno™ RNA Control (10,000 copies/µL)	110 µL	-25°C to -15°C
Nucleic Acid Dilution Solution	500 µL	-30°C to -10°C

Required materials not supplied

Unless otherwise indicated, all materials are available through thermofisher.com. MLS: Fisher Scientific (fisherscientific.com) or other major laboratory supplier.

Item	Source
Real-time PCR thermal cycler capable of detecting: • FAM™ (maximum emission: λ515 nm) • CAL Fluor™ Orange 560 (maximum emission: λ561 nm)	MLS
96-well plate, PCR strips (8- or 12-wells), microtubes or capillaries compatible with thermal cycler used	MLS
Suitable plate covers or caps	MLS
Nuclease-free pipettes and filtered pipette tips	MLS
Nuclease-free reagent tubes for preparing master mix	MLS
2 ice buckets or refrigerated racks: • One for the PCR setup area where the RT-PCR master mix is prepared • One for the area where RNA may be present	MLS
1X TE buffer	MLS
DNase/RNase-free water	MLS

Isolate RNA from samples

RNA extraction from initial samples is required prior to real-time PCR analysis.

Step or process	Recommendation
Preparation of mock-purified samples for use as extraction control	Prepare mock-purified samples, using nuclease-free water as the starting material. Process with the same RNA isolation method that is used for test samples.
Proposed RNA isolation methods	MagMAX™-96 Viral RNA Isolation Kit (Cat. Nos. AM1836, AMB18365) or an equivalent RNA isolation method.
Required modifications to the RNA isolation method	<ul style="list-style-type: none"> • Add 1 µL of undiluted Xeno™ RNA Control per isolation to the lysis solution used for RNA isolation. • Add carrier RNA to the lysis solution according to the manufacturer recommendations. Carrier RNA is provided in the MagMAX™-96 Viral RNA Isolation Kit (Cat. Nos. AM1836, AMB18365)

Perform real-time RT-PCR

- 1 Determine the quantity of reactions and thaw the reagents.
- For each real-time RT-PCR run, include the following reactions:
 - Positive control; use 10.5 µL of diluted NDV Positive Control.
 - Extraction control; use 10.5 µL of mock-purified samples.
 - No-template control; use nuclease-free water in place of sample RNA.
 - Thaw PCR master mix reagents in one ice bucket or refrigerated rack and controls and samples in a separate ice bucket or refrigerated rack, gently vortex each tube to mix the contents thoroughly, then briefly centrifuge to collect the solution at the bottom of the tube. Keep the reagents on ice or on a refrigerated rack.

- 2 Prepare the diluted NDV Positive Control. Combine the following components and place the mixture on ice between 2°C and 8°C for immediate use, or below -16°C for later use.

Component	Volume
Nucleic Acid Dilution Solution	9.5 µL
25X NDV Control RNA	1.0 µL
Total volume of diluted NDV Positive Control	10.5 µL

- 3 Prepare the real-time RT-PCR master mix on ice. Combine the following components for the number of reactions required plus 10% overage.

Component	Volume per reaction
2X qRT-PCR Buffer	12.5 µL
25X qRT-PCR Enzyme Mix	1.0 µL
25X NDV Primer Probe Mix	1.0 µL
Total volume of real-time RT-PCR master mix	14.5 µL

- 4 Set up the real-time RT-PCR reactions.
- Dispense 14.5 µL of real-time RT-PCR master mix to the appropriate PCR plate well, PCR strips, or capillaries.
 - Add the appropriate component for the sample type, according to the following table:

Sample type	Component	Volume per reaction
Test sample	Sample RNA	10.5 µL
Positive control	Diluted NDV Positive Control	10.5 µL
Extraction control	Mock-purified sample	10.5 µL
No-template control	Nuclease-free water	10.5 µL

- c. Seal each reaction vessel.

- 5 Set up and run the real-time RT-PCR instrument.
- Following the manufacturer's instructions, set up the following parameters for the real-time RT-PCR run.
 - Reaction volume: 25 µL
 - ROX™ passive reference dye: Included in 2X qRT-PCR Buffer (required if the thermal cycler is capable of ROX™ detection)

- b. Set up and assign TaqMan® probe reporter dyes and quenchers for each well, tube, or capillary used in the analysis:

Target	Reporter	Quencher
NDV	FAM™ dye	BHQ™-1 dye
Xeno™ RNA Control	CAL Fluor™ Orange 560 dye	BHQ™-1 dye

Note: Use CAL Fluor™ Orange 560 dye to calibrate the thermal cycler, if possible. Otherwise use equivalent dye detectors such as VIC™.

- c. Run the thermal cycler program and collect real-time amplification data during the elongation step (45 seconds at 60°C). Use the following thermal cycler settings:

Stage	Repetitions	Temperature	Time
1	1×	48°C	10 minutes
2	1×	95°C	10 minutes
3	40×	95°C	15 seconds
		60°C	45 seconds

Data analysis

Refer to the recommendations of the thermal cycler manufacturer for raw data analysis.

1. Set thresholds separately for each real-time RT-PCR target.
2. Interpret the results for each control and samples according to the obtained C_t values as indicated in the following section.

Validation

The run is validated if the following criteria are met.

Control reaction	C_t value for NDV RNA	C_t Value for Xeno™ RNA Control
Positive control	24–27	>40 (no signal detected)
Extraction control	>40 (no signal detected)	29–34
No-template control	>40 (no signal detected)	>40 (no signal detected)

Interpretation

C_t value for NDV RNA	C_t Value for Xeno™ RNA Control	Interpretation
<40	<40	Sample positive for NDV
>40 (no signal detected)	< C_t Extraction control \pm 3 C_t	Sample negative for NDV
>40 (no signal detected)	>40 (no signal detected)	Not validated ⁽¹⁾

⁽¹⁾See "Troubleshooting".

Troubleshooting

Observation	Possible cause	Recommended action
Test samples Xeno™ RNA Control—no signal NDV RNA—high signal	The Xeno™ RNA Control primers and probes are present at limiting concentrations in the RT-PCR reactions. High levels of NDV RNA in a sample can reduce amplification of Xeno™ RNA Control.	No or low signal from the Xeno™ RNA Control is expected in a reaction that has a strong signal for NDV RNA.
Test samples NDV RNA—no signal Xeno™ RNA Control—no signal	Xeno™ RNA Control was omitted	Verify that Xeno™ RNA Control was added to the lysis solution during RNA extraction.
	Poor RNA recovery	Check the recovery of RNA carrier used in RNA extraction.
	The RNA samples contain inhibitors of RT-PCR	Reduce the quantity of sample added to the RT-PCR reactions. For example: <ul style="list-style-type: none"> • Add 1–2 μL of sample (add nuclease-free water to bring the reaction to the proper volume) • Dilute the RNA sample 1:10 in Nucleic Acid Dilution Solution before adding it to the reaction. • The C_t values for Xeno™ RNA Control and NDV amplifications are expected to decrease proportionally to the decrease in sample quantity (\sim2–3C_t).
	Poor RNA recovery Or The RNA samples contain inhibitors of RT-PCR	Compare the results of amplifying sample RNA using RT-PCR master mix with and without 1 μ L of Xeno™ RNA Control: <ul style="list-style-type: none"> • If the reactions amplified using the RT-PCR master mix with addition of Xeno™ RNA Control return a signal, but the reactions amplified using the RT-PCR master mix without addition of Xeno™ RNA Control return no signal, Xeno™ RNA Control was not recovered. • If Xeno™ RNA Control signal is not detected in either sample, the RNA sample contains inhibitors of real-time RT-PCR.

Observation	Possible cause	Recommended action
Positive control reaction NDV Positive Control—no signal Xeno™ RNA Control—no signal	The NDV Positive Control was improperly handled, resulting in RNA degradation.	Use appropriate precautions against RNase contamination when handling the control RNAs. For example, wear clean gloves and use nuclease-free barrier pipette tips.
	The 25X qRT-PCR Enzyme Mix was stored or handled improperly, and it lost activity.	Repeat the RT-PCR with fresh reagents.
	The thermal cycler was not properly set up.	Check the thermal cycler settings. See "Set up and run the real-time PCR instrument." on page 5.
	The RT-PCR master mix was prepared incorrectly.	Repeat the test with correctly prepared RT-PCR master mix.
No-template control reaction NDV RNA—signal detected or Xeno™ RNA Control—signal detected	Contamination during the PCR.	<ul style="list-style-type: none"> Repeat the real-time RT-PCR with fresh reagents and freshly decontaminated pipettes. Set up the real-time RT-PCR in an area separate from areas used for RNA isolation and PCR product analysis.
Extraction control reaction NDV RNA—signal detected	Contamination during RNA extraction or the PCR.	<ul style="list-style-type: none"> Repeat the RNA isolation or real-time RT-PCR with fresh reagents and freshly decontaminated pipettes. Set up the real-time RT-PCR in an area separate from areas used for RNA isolation and PCR product analysis.

Documentation and support

Customer and technical support

Technical support: visit thermofisher.com/askaquestion
 Visit thermofisher.com/support for the latest in services and support, including:

- Worldwide contact telephone numbers
 - Order and web support
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)
- NOTE:** For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at thermofisher.com/support.

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The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history of Pub. No. MAN0008826 (English)

Revision	Date	Description
B.0	9 May 2017	Updated to the current document template, with associated updates to the warranty, trademarks, and logos.
A.0	10 April 2014	Baseline for revision history

Limited Use Label License No. 460: PCR Veterinary Diagnostics. Notice to Purchaser: For veterinary diagnostic services including the reporting of such services for a fee and for research purposes only. Human diagnostic uses require a separate license from Roche.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche, used under permission and license. CAL Fluor and BHQ are trademarks of BioSearch Technologies, Inc.