

# Medtronic

# Stapling Evidence

Enter

SECTION 1  
Powered Stapling  
VS  
Manual stapling

Select

SECTION 2  
Buttress  
Colorectal

Select

SECTION 3  
Environmental  
Impact

Select

# Acronyms

AB	Anastomotic bleeding
AL	Anastomotic leak
CS	Circular staplers
DCS	double-row circular stapler
LVRS	lung volume reduction
MCS	Two-row manual circular stapler
MUS	Multi-Use Surgical Staplers
NR	Not Reported
NSCLC	non-small cell lung cancer
PCS	Powered circular stapler
PGA	polyglycolic acid
PSW	Propensity score weighting
pts	Patients
RCT	Randomized controlled trial
SUS	Single-Use Stapler
TCS	triple-row circular stapler
TRCS	Three-row manual circular stapler
VATS	video-assisted thoracoscopic surgery





Home

## SECTION 1

# Clinical advantages of powered stapling vs manual stapling



Acronyms

Articles (conclusions)	Tissue	Staplers
<a href="#">Martín-Arévalo J. et al. 2024</a> Powered circular stapler presented a significantly lower risk of leakage and anastomotic bleeding while Three-row manual circular stapler only demonstrated a risk reduction in Anastomotic leak (AL). Risk difference of AL was superior in the Powered circular stapler than in Three-row manual circular stapler.	Colorectal	Two-row manual circular stapler Three-row manual circular stapler Powered circular stapler
<a href="#">McKechnie T. et al. 2024</a> Triple-row staple technology may reduce the risk of anastomotic leak in left-sided colorectal anastomoses.	Colorectal	Tri-staple Technology with EEA circular staplers (Medtronic) Tri-staple Technology with the Endo GIA articulating reload (Medtronic)  2-row circular staplers Medtronic Endo GIA Universal Straight Reload (Medtronic)
<a href="#">Vanstraelen S. et al. 2023</a> The introduction of a powered circular stapler contributed to an 80% risk reduction of anastomotic leakage after esophageal resection for cancer with intrathoracic anastomoses in a propensity score matched monocentric observational series. This consequently, resulted in a reduced length of stay.	Esophagogastric	Manual: EEA™ Circular Stapler (Medtronic) Powered: Echelon Circular™ Powered Stapler (Ethicon)
<a href="#">Wang T. et al. 2023</a> The risk of AL observed with manual two-row circular staplers was similar to that seen with three-row devices. This study affirms the safety of manual two-row circular staplers in colorectal anastomosis.	Colorectal	Ethicon manual circular staplers Medtronic EEA™ circular stapler with DST™ technology Medtronic EEA™ circular stapler with Tri-Staple™ technology
<a href="#">Gan C. et al. 2022</a> Gripping Surface Technology (GST) system was associated with better intraoperative outcomes in clinical practice in China	Lung	Echelon Flex GST powered staplers  manual: NR



PubMed link



PubMed link



PubMed link



PubMed link



PubMed link



Home

SECTION 1

# Clinical advantages of powered stapling vs manual stapling



Acronyms

Articles (conclusions)	Tissue	Staplers
<p><a href="#">Quero G. et al. 2022</a>            Three-row CSs significantly decrease the rate of AL after rectal resection. Further multicenter controlled trials are still needed to confirm the advantages of three-row CSs on anastomotic complications.</p>	Colorectal	Two-row circular stapler: ETHICON™ Circular Stapler (Ethicon) - DST Series EEA (Medtronic)  Three-row circular stapler: Tri-staple™ Technology (Medtronic)
<p><a href="#">McKinnon C. et al. 2022</a>            The Medtronic I-Drive and the Ethicon Echelon are comparable in reload time, stapler misfires, leak test rates, and cost.</p>	Gastrointestinal	Medtronic I-Drive powered stapler with EndoGIA reinforcement arms  Ethicon Echelon powered stapler with Seamguard
<p><a href="#">Shigeeda W. et al. 2021</a>            Using a powered stapler to divide the incomplete interlobular fissure decreased the need for additional intraoperative management using fibrin glue and reduced postoperative air leakage in radical pulmonary resection</p>	Lung	Powered stapler: Signia™ stapling system (Medtronic)  Manual stapler: ENDO GIA™ Ultra Universal Stapler (Medtronic)  Reinforcement: Endo GIA™ Reinforced Reload with Tri-Staple™ technology (Medtronic)
<p><a href="#">Akil A. et al. 2019</a>            The iDrive powered stapling system offers one-handed, push-button operation, which eliminates the manual firing force and possibly enables more precise resection. In the current study, the novel system led to comparable results with the conventional mechanical stapler without any disadvantages in patients undergoing bilateral VATS-LVRS.</p>	Lung	Powered: iDrive powered stapling system (Medtronic) Manual: EndoGIA (Medtronic)



PubMed link



PubMed link



PubMed link



PubMed link



Home

SECTION 1

# Clinical advantages of powered stapling vs manual stapling



Acronyms

Articles (conclusions)	Tissue	Stapler
<a href="#">Mazaki J. et al. 2022</a> Propensity score-matched analysis of clinical data suggested that the use of TCS contributes to a reduced risk of AL after colorectal anastomosis with the double stapling technique	Colorectal	EEA™ circular stapler with Tri-Staple™ technology Double-row circular stapler
<a href="#">Nakanishi R. et al. 2021</a> Anastomoses created using a triple-layered circular stapler had high bursting pressure, which might contribute to a lower incidence of anastomotic leakage after rectal cancer surgery.	Colorectal	Double- or triple-layered circular stapler



PubMed link



PubMed link



Home

## SECTION 1

# Clinical advantages of powered stapling vs manual stapling



Acronyms

### Other suggested readings

	Study type
<p><i>EE145 Cost-Effectiveness Analysis of the Novel Powered Circular Stapler in Colorectal Anastomosis Surgery</i> González de Julián S. et al. 2022 Value in Health (2022) 25:12 Supplement (S81)</p>	Cost-effectiveness study
<p>Indirect comparison between powered and manual circular staplers for left-sided colorectal anastomoses: clinical and economic outcomes in China Bai J. et al. 2022 Cost Effectiveness and Resource Allocation (2022) 20:1 Article Number: 45</p>	Decision analysis model
<p>Economic analysis of leak complications in anastomoses performed with powered versus manual circular staplers in left-sided colorectal resections: A us-based cost analysis Pollack E. et al. 2021 ClinicoEconomics and Outcomes Research (2021) 13 (531-540).</p>	Budget impact model
<p>A Novel Predictive Model for Anastomotic Leakage in Colorectal Cancer Using Auto-artificial Intelligence Mazaki J et al 2021 Anticancer Res. 2021;41(11):5821-5825</p>	Efficacy using AI model
<p>Comparison of pressure resistance of double-rows and triple-rows circular stapler in rectal double stapling technique: In vitro study. Mazaki J et al 2022 Medicine (Baltimore). 2022;101(28):e29600</p>	Pre-clinical study
<p>Evaluating Performance of Circular Staplers Using Comparative Test Methods for Evidence-Based Surgery Strassner H. et al. 2023 Surg Innov. 2023;30(5):576-585</p>	Pre-clinical study



PubMed link



PubMed link



PubMed link



PubMed link



PubMed link



PubMed link

Powered stapling  
vs  
manual stapling

# Two-row, three-row or powered circular stapler, which to choose when performing colorectal anastomosis?

## A systematic review and meta-analysis

Martín-Arévalo J. et al.

*Int J Colorectal Dis.* 2024;39(1):51



The *objective* is to provide the existing evidence on which of these circular staplers would have a lower risk of presenting a leak and/or anastomotic bleeding.



### Results

• *Anastomotic bleeding*  
138/2200 pts (6.27%)

PCS: 5 (0.96%)  
MCS: 133 (7.91%)

Common-effect model:  
OR: 0.19 (95% CI: 0.07-0.048)

Random-effect model:  
0.2 (95% CI: 0.08-0.52) I<sup>2</sup> 0.0%, Q:  
1.29, p = 0.526  
results suggest a protective effect  
of PCS

Risk: - 0.03 (95% CI: - 0.07 to -  
0.01) I<sup>2</sup> 87%  
PCS needed to avoid bleeding  
would be 34 cases

• *Anastomotic leakage/air leakage incidence*  
A) 208/2947 pts  
PCS: 35 (4.1%)  
MCS: 173 (8.23%)  
(p < 0.001)

Common-effect model:  
OR: 0.3695 ( 95%-CI: 0.247-0.553)

Random-effect model:  
OR: 0.402 (0.266-0.608) I<sup>2</sup> 0.0%, H<sup>2</sup> 0.91  
results suggest a protective effect of PCS

Risk: - 0.06 (95% CI: - 0.07 to - 0.04) I<sup>2</sup> 0%  
significantly decreased risk of PCS use  
The number of patients needed to treat with PCS to avoid  
leakage would be 17.

B) 1631/16.490 pts  
TRCS: 123 (6.26%, 1.8-8.19%)  
MCS: 1508 (9.14%, 6.12-11.6%)

Random-effect model:  
0.446 (95%-CI: 0.217 to 0.916) I<sup>2</sup> 75.1%, Q = 12.05, p = 0.007  
results suggest a protective effect of PCS

Risk: - 0.04 (95%-CI: - 0.04 to - 0.08)  
suggests a significantly decreased risk of TRCS  
The number of patients needed to treat with TRCS to avoid  
leakage would be 24.



### Insight

- **Powered circular stapler** could have a **significantly lower risk of leakage and anastomotic bleeding** than a two-row manual circular stapler.
- The **three-row circular stapler** may also have a **reduced risk of anastomotic leakage** compared with the two-row manual circular stapler.
- The reduction of the risk of anastomotic leakage could be greater in the powered circular stapler group than in the three-row circular stapler group.

Meta-analysis  
7 articles

A)  
Two-row manual circular  
stapler (MCS): 2103 cases  
Powered circular stapler  
(PCS): 844 cases

B)  
Three-row manual circular  
stapler (TRCS): 1965 cases  
(10.65%)  
Two-row manual circular  
stapler (MCS): 16.490 cases  
(89.35%)

Indication: colorectal

Acronyms

PubMed link



Powered stapling  
vs  
manual stapling

# Double-row staple technology versus triple-row staple technology for colorectal surgery: A systematic review and meta-analysis

Tyler McKechnie et al.

*Surgery* 176 (2024) 633e644



The *objective* is to compare double- and triple-row staple technology for colorectal anastomoses



## Results

- *Anastomotic bleeding*  
Triple-row stapler vs double-row stapler  
1.2% vs 3.9%, RR 0.47, 95% CI 0.15-1.49, P: .20, I2: 57%
- *Anastomotic leakage/air leakage incidence*  
Pooled analysis  
Significantly reduced with triple-row staple technology (6.3% vs 7.5%, RR: 0.54, 95% CI 0.31-0.94, P: .03, I2: 75%)  
  
Sensitivity analysis  
Significant reduction with triple-row staple technology (3.4% vs 5.9%, RR 0.55, 95% CI 0.36 - 0.85, P: .01, I2: 4%)\*  
\*exclusion of articles with conflict of interest and non-circular staplers  
  
Findings no longer statistically significant (6.4% vs 7.4%, RR 0.65, 95% CI 0.39 -1.09, P: .10, I2: 64%)\*\*  
\*\*exclusion of non-circular staplers  
  
Significant reduction with triple-row staple technology (3.7% vs 6.6%, RR 0.47, 95% CI 0.31 - 0.70, P < .01, I2: 17%\*\*\*  
\*\*\*exclusion of articles with conflict of interest
- *Length of hospital stay*  
Mean: 0.80 days less in the triple-row  
(95% CI -3.59 to 1.99 days, P: .57, I2: 98%)
- *Operation time*  
(95% CI -25.80 to 39.61 minutes, P: .68, I2: 76%)  
Mean difference of 6.91 minutes favoring double-row staple



## Insight

- **Triple-row staple technology** was associated with a **significant reduction** in anastomotic leak risk
- Mortality and reoperation showed no significant difference

Meta-analysis  
6 articles

Double-row staple: 19.372  
Triple-row staple: 2.298

Indication: colorectal  
surgery

Tri-staple Technology with  
EEA circular staplers (tri-  
EEA) (Medtronic)  
Tri-staple Technology with  
the Endo GIA articulating  
reload (Medtronic)

2-row circular staplers  
Medtronic Endo GIA  
Universal Straight Reload  
(Medtronic)

Acronyms

 PubMed link



Powered stapling  
vs  
manual stapling

# Real-life introduction of powered circular stapler for esophagogastric anastomosis: cohort and propensity matched score study

Vanstraelen S. et al.

*Dis Esophagus. 2023;36(5):doac073*



The *objective* is to evaluate the impact of a powered circular stapler on complications after esophagectomy with intrathoracic anastomosis for esophageal cancer.



## Results

- *Anastomotic leakage/air leakage incidence*  
Manual: 13/66 (19.7%)  
Powered: 2/62 (3.2%)  
p=0.004

Relative risk reduction: 83%  
OR: 7.3 (95% confidence interval [CI]: 1.58-33.74).

Propensity score matched  
Manual: 10 (20.4%)  
Powered: 2 (4.1%)  
p: 0.013

- *Length of hospital stay*  
Mean  
Manual: 18.3  
Powered: 11.9  
p: 0.026

Manual: 18.7  
Powered: 11.1  
p: 0.022

- *Adverse Events*  
90-day mortality  
Manual: 4 (6.1%)  
Powered: 1 (1.6%)  
p: 0.366

Comprehensive Complication  
Manual: 33.8  
Powered: 23.0  
P: 0.005



## Insight

- **Powered circular stapler** contributed to an **80% risk reduction of anastomotic leakage** after esophageal resection for cancer with intrathoracic anastomoses in a propensity score matched analysis.
- This consequently, **resulted in a reduced length of stay**

Retrospective study

128 pts  
Powered: 62  
Manual: 66

Propensity score matched cohort  
Powered: 49  
Manual: 49

Indication Esophagectomy

Manual: EEA™ Circular Stapler (Medtronic)  
Powered: Echelon Circular™ Powered Stapler (Ethicon)

Powered stapling  
vs  
manual stapling

## Risk of anastomotic leakage with two-row versus three-row manual circular staplers in colorectal anastomosis: a U.S. cohort study

Wang T. et al.

*Int J Colorectal Dis.* 2023;38(1):264



The *objective* is to compare the risk of anastomotic leak (AL) between Ethicon manual circular staplers (two-row) versus Medtronic EEA™ circular stapler with Tri-Staple™ technology (three-row) and between Medtronic EEA™ circular stapler with DST™ Series technology (two-row) versus Tri-Staple™ technology



### Results

- *Anastomotic leakage/air leakage incidence*

Ethicon manual circular staplers vs. Medtronic EEA™ circular stapler with Tri-Staple™ technology  
RR: 0.72; 95% CI, 0.52-1.01

Medtronic EEA™ circular stapler with DST™ technology vs. Medtronic EEA™ circular stapler with Tri-Staple™ technology  
RR, 0.75; 95% CI, 0.53-1.06



### Insight

- The **risk of AL** observed with **manual two-row circular staplers** was **similar** to that seen **with three-row devices**, after PSW.
- Similar results were found for patients who did not have a diverting stoma prior to or on the same day as a left-sided colorectal surgery during the index admission.
- **Two-row manual circular staplers** are **safe** for **colorectal anastomosis**

### Retrospective study

Ethicon manual circular staplers: 8337  
Medtronic EEA™ circular stapler with DST™ technology: 7928  
Medtronic EEA™ circular stapler with Tri-Staple™ technology: 1306

Indication: Left-sided colorectal surgery

Powered stapling  
vs  
manual stapling

## Powered stapling system with gripping surface technology for pulmonary resection of lung cancer: real-world clinical effectiveness

Gan C. et al.

*Cost Eff Resour Alloc.* 2022;20(1):72



The *objective* is to examine the real-world effectiveness of Gripping Surface Technology system on intraoperative and postoperative outcomes of pulmonary resection



### Results

- *Anastomotic leakage/air leakage incidence*

No bubble

Powered: 48 (88.9%)

Manual: 46 (85.2%)

p 0.556

Countable bubble

Powered: 5 (9.3%)

Manual: 8 (14.8%)

- *Length of hospital stay*

Powered: 7.78 ± 2.72

Manual: 7.71 ± 2.38

p 0.877

- *Operation time*

Powered: 118.04 ± 37.39

Manual: 127.06 ± 46.22

p 0.281

- *Intraoperative bleeding*

Powered: 22.8%

Manual: 51.9%

p = 0.003

- *Use of NEOVEIL (reinforcement material to prevent air leakage from the staple line) intraoperatively*

Powered: 24.1%

Manual: 50%

p = 0.01



### Insight

- **GST system** is associated with better intraoperative outcomes, such as **reduced** bleeding and intervention rates.

Retrospective study  
108 pts  
Powered: 54  
Manual: 54

Indication: Surgical lung  
resection

Echelon Flex GST powered  
staplers

Manual: NR

Powered stapling  
vs  
manual stapling

# Preliminary evaluation of two-row versus three-row circular staplers for colorectal anastomosis after rectal resection: a single-center retrospective analysis

Quero G. et al.

*Int J Colorectal Dis.* 2022;37(12):2501-2510



The *objective* is to evaluate potential advantages of three-row circular staplers (Three-CS) on anastomotic leakage (AL), stenosis (AS), and hemorrhage (AH) rates after rectal resection as compared to two-row circular staplers (Two-CS).



## Results

- *Anastomotic bleeding*  
two-CS: 3 (1.5)  
three-CS: 2 (1.1)  
p 0.73

- *Anastomotic leakage/air leakage incidence*  
two-row CS: 19 (9.6)  
three-row CS: 6 (3.4)  
p 0.01

AL onset was significantly associated to:  
BMI  $\geq 25$  (p = 0.05)  
tumor lesions located in the middle/low rectum (p < 0.0001)  
Two-CS (p = 0.01)

Multivariate analysis  
Tumor location in the middle/low rectum OR: 4 [95% CI: 1.64-9.87] p = 0.002  
Two-CS OR: 2.63 [1.07-6.46] p = 0.03  
Independent prognostic factors

Subanalysis on patients with low rectal cancer  
46 patients  
Two-CS: 24 (52.2%)  
Three-CS: 22 (47.8%)  
p = 0.95

- *Length of hospital stay*  
two-row CS: 8 (6-11)  
three-row CS: 5 (4-6)  
p < 0.0001

- *Surgical approach minimally invasive procedures*  
Three-CS group: 170-95.5%  
Two-CS: 141-71.5%  
p < 0.0001

- *Anastomotic stenosis*  
two-row CS: 1 (0.5)  
three-row CS: 1 (0.6)  
p 0.94

- *in-hospital mortality*  
Two-row CS: 1 (0.5)  
Three-row-CS: 0  
p 0.34



## Insight

- **Less** anastomotic leakage incidence when the **Three-CS** was employed
- The **type of circular stapler did not show any influence** on anastomotic leakage severity
- Furthermore, the **Two-CS** and **Three-CS** were **comparable** for anastomotic bleeding and anastomotic stenosis rate

Retrospective study  
375 pts  
two-row circular: 197  
(52.5%)  
three-row circular: 178  
(47.5%)

Indication: Rectal resection

Two-row circular stapler:  
ETHICON™ Circular Stapler  
(Ethicon) - DST Series EEA  
(Medtronic)  
Three-row circular stapler: Tri-  
staple™ Technology  
(Medtronic)

Powered stapling  
vs  
manual stapling

## Medtronic I-Drive vs Ethicon Echelon: A Head-to-head Randomized Controlled Trial

McKinnon C. et al.

World J Lap Surg 2022; 15 (1):65-68



The *objective* is to compare Echelon with Seamguard (ESG) or I-Drive with EndoGIA reinforcement arms (GIA-R) in sleeve gastrectomies.



### Results

- *Anastomotic bleeding*  
GIA-R system: 15.63  
Echelon: 34.38%  
p = 0.44

- *Operation time*  
Sleeve creation time (minutes)  
Echelon: 12.14  
GIA-R system: 14.26  
P=0.04  
Time to reload (seconds)  
Echelon: 41.77  
GIA-R system: 39.78  
P=0.42

- *Economic outcomes*  
Stapler cost (Government rate)  
Echelon: \$2,449.44  
GIA-R system: \$2,097.66  
P=0.0002  
Stapler cost (Commercial rate)  
Echelon: \$1,982.17  
GIA-R system: \$2,037.25  
P=0.48

- *Stapler loads used*  
Echelon: mean of 5.38  
GIA-R arm: mean of 4.92  
p = 0.052



### Insight

- **Echelon with Seamguard (ESG) or I-Drive with EndoGIA reinforcement arms (GIA-R) were found to be **comparable** regarding reload time, misfires, leak rates, and overall cost when considering commercial pricing.**

Randomized controlled trial (RCT)  
63 pts  
i-Drive: 25  
Ethicon: 26

Indications: laparoscopic sleeve gastrectomies

Medtronic I-Drive powered stapler with EndoGIA reinforcement arms  
Ethicon Echelon powered stapler with Seamguard

## Powered stapling vs manual stapling

# Utility of the powered stapler for radical pulmonary resection: a propensity score-matched analysis

Shigeeda W. et al.

*Surgery Today* (2021) 51:4 (582-588)



The *objective* is to compare the utility of a powered stapler with that of a manual stapler, using the same cartridge, for pulmonary resection.



### Results

- Anastomotic leakage/air leakage incidence*

Powered: 2.8%

Manual: 8.6%

$p = 0.014$

Propensity score analysis

Powered: 2.8%

Manual: 10.7%

$p = 0.003$

- Length of hospital stay*

Powered:  $8.5 \pm 4.0$  days

Manual:  $11.6 \pm 17.9$  days

$p = 0.023$

Propensity score analysis

Powered:  $8.5 \pm 4.0$  days

Manual:  $12.2 \pm 21.27$  days

$p = 0.026$

- Operation time*

no significant differences

Duration of chest tube placement

Powered:  $4.5 \pm 1.7$  days

Manual:  $5.1 \pm 3.1$  days

$p = 0.011$

Propensity score analysis

Duration of chest tube placement

Powered:  $4.5 \pm 1.7$  days

Manual:  $5.3 \pm 3.6$  days

$p = 0.006$

- Economic outcomes*

Respective average costs, those for the powered-stapler group were significantly less than those for the manual-stapler group (JPY  $2,102,090 \pm 21,894.5$  vs.  $2,260,660 \pm 76,199.1$ ,  $p = 0.046$ ).



### Insight

- Powered staplers** are associated with **decreased** [fibrin glue usage](#) and **reduced** [air leakage](#), as well as **shorter** [chest tube placement](#) and [hospital stays](#).
- Powered staplers** offer [advantages in pulmonary procedures](#), potentially leading to **improved** [patient outcomes](#) and [cost savings](#)

Retrospective study  
478 pts

Powered: 179  
Manual: 269

Propensity score analysis  
Powered: 98.9% (177/179)  
Manual: 65.8% (177/269)

Indication: Radical  
pulmonary resection

Powered stapler: Signia™  
stapling system (Medtronic)

Manual stapler: ENDO GIA™  
Ultra Universal Stapler  
(Medtronic)

Reinforcement: Endo GIA™  
Reinforced Reload with Tri-  
Staple™ technology  
(Medtronic)

Powered stapling  
vs  
manual stapling

## Use of a Powered Stapling System for Minimally Invasive Lung Volume Reduction Surgery: Results of a Prospective Double-Blind Single-Center Randomized Trial

Akil A. et al.

*Thoracic and Cardiovascular Surgeon (2019) 67:3 (216-221)*



The *objective* is to test iDrive powered stapling system during video-assisted thoracoscopic surgery (VATS) for lung volume reduction (LVRS) and compare with a non-electronic conventional device



### Results

- *Anastomotic leakage/air leakage incidence*  
iDrive-group: 230 +- 81 (0-2,300)  
EndoGIA-group: 321 +- 116 (0-3,400)  
p = 0.6

Day1  
iDrive-group: 308 +- 94 (0-2,200)  
EndoGIA-group: 255 +- 95 (0-2,300)  
p = 0.7

- *Operation time*  
iDrive-group: 52 +- 2.5 (29-93)  
EndoGIA-group: 54 +- 3.8 (21-105)  
p = 0.5



### Insight

- **iDrive system** performed comparably to the conventional stapler, showing **no significant differences** in surgical duration, air leakage, or length of chest tube therapy.

Prospective, randomized trial  
40 pts

Inpatient control

Indication: advanced lung  
emphysema

Powered: iDrive powered  
stapling system (Medtronic)  
Manual: EndoGIA (Medtronic)

Powered stapling  
vs  
manual stapling

## Effectiveness of a new triple-row circular stapler in reducing the risk of colorectal anastomotic leakage: A historical control and propensity score-matched study

Junichi Mazaki et al.

Medicine (Baltimore). 2022;101(27):e29325



The *objective* is to evaluate the effectiveness of the EEA™ circular stapler, a new triple-row circular stapler (TCS), relative to a conventional, double-row circular stapler (DCS).



### Results

- *Anastomotic leakage/air leakage incidence*  
17/285 (6.0%)

DCS: 15 (7.4%)  
TCS: 2 (2.4%)  
p: 0.17

Case matching  
DCS: 13 (11.6%)  
TCS: 1 (1.8%)  
p: 0.04

- *Diverting stoma*  
AL after case matching, there were 6 patients in the DCS group (46%) and no patients in the TCS group (0%)



### Insight

- **EEA™ circular stapler with Tri-Staple™ technology** is associated with a **significantly lower risk** of anastomotic leak compared to the DCS

Retrospective study  
285 pts

Double-row CS: 202  
Triple-row CS: 83

Propensity score matching  
1:2  
112:56

Indication: Colorectal  
cancer

EEA™ circular stapler with  
Tri-Staple™ technology

Double-row circular stapler

Powered stapling  
vs  
manual stapling

# Clinical impact of the triple-layered circular stapler for reducing the anastomotic leakage in rectal cancer surgery: Porcine model and Multicenter retrospective cohort analysis

Nakanishi R. et al.

*Ann Gastroenterol Surg.* 2021;6(2):256-264



The *objective* is to investigate the impact of the triple-layered circular stapler compared with the double-layered circular stapler on anastomotic leakage after rectal cancer surgery.



## Results

- *Anastomotic leakage/air leakage incidence*  
triple-layered stapler: 0.0%  
double-layered stapler: 5.8%  
P = .0362
- *Length of hospital stay*  
EEA group: 14 days  
tri-EEA: 4 days  
P = .4757
- *Operation time*  
EEA: 227 min  
tri-EEA: 278 min  
p=0.0042
- *Multivariate analysis factors independently associated with a lower incidence of anastomotic leakage*  
triple-layered vs stapler usage  
OR 0.00, 95% confidence interval: 0.00-0.96, P = .0465
- *Postoperative complications*  
(CD grade  $\geq 3$ )  
EEA group: 8%  
tri-EEA: 2%  
P = .1731



## Insight

- Anastomoses created using a **triple-layered circular stapler** had high bursting pressure, which **might contribute to a lower incidence of anastomotic leakage** after rectal cancer surgery.

Retrospective study  
194 pts

Double-layered: 153  
Triple-layered: 41

Indication: rectal cancer

EEA™ Circular Stapler with  
Tri-Staple™ Technology  
(tri-EEA)  
EEA™ Circular Stapler with  
DST Series™ Technology  
(EEA)



Home

SECTION 2

# Buttress reload use for colorectal



Acronyms

Article (conclusions)	Tissue	Stapler
<p><a href="#">Gerbasi I.S. et al. 2024</a>            Tri-Staple™ reduces the risk of leakage in colorectal anastomosis. However, this study provides only insights, and further research is warranted to confirm these findings.</p>	Colorectal tract	DST™ Series technology (Medtronic) circular stapler  EEA Tri-Staple™ (Medtronic) circular stapler
<p><a href="#">Mori S. et al. 2024</a>            The simple technique for constructing a sutureless overlapped anastomosis using a 60 mm linear stapler with a reinforced bioabsorbable material in robotic right colectomy with intracorporeal anastomosis appears to be safe and feasible.</p>	Colorectal tract	Signia Stapling Technology; Endo GIA 60 mm Articulating Medium Thick Reload with Tri-Staple™ Technology (Medtronic)
<p><a href="#">Shibutani M. et al. 2023</a>            Double-stapling technique anastomosis with PGA sheet, which is easy to perform, contributes to the reduction of anastomotic leakage rate by increasing the strength of the anastomotic site.</p>	Colorectal tract	PGA: Neoveil
<p><a href="#">Hokonohara K. et al. 2022</a>            The overlapped delta-shaped anastomosis technique using linear staplers with reinforced bioabsorbable material appears to be both safe and feasible.</p>	Colorectal tract	Endo GIA 60 mm Articulating Medium Thick Reload with Tri-Staple Technology (Medtronic)
<p><a href="#">Shibutani M. et al. 2021</a>            A PGA sheet may be effective for preventing anastomotic leakage in double-stapling technique anastomosis for rectal surgery.</p>	Colorectal tract	PGA: Neoveil
<p><a href="#">Naito M. et al. 2017</a>            The automatic PGA-felt suture device is safe and effective for colorectal resections and anastomoses</p>	Colorectal tract	PGA: Neoveil  Endo GIA Reinforced Reload with Tri-Staple Technology  DST Series EEA Stapler (Medtronic)



PubMed link



PubMed link



PubMed link



PubMed link



PubMed link



PubMed link



Home

SECTION 2

# Buttress reload use



Acronyms

Article (conclusions)	Tissue	Stapler
<p><a href="#">Deguchi H et al. 2020</a> Using the stapler with PGA to divide the incomplete interlobular fissure in pulmonary lobectomy reduced postoperative air leakage and decreased the need for additional intraoperative management using fibrin glue.</p>	Lung	PGA: Neoveil and Gunze
<p><a href="#">Khandhar S.J. et al. 2024</a> The two-row stapler reloads demonstrated favorable safety and efficacy profiles when fired across hilar vessels in the thoracic space with a 99% hemostatic rate, independent of surgical access and stapler handle preference.</p>	Lung	Signia™ Small Diameter Reloads



PubMed link



PubMed link

Buttress reload  
use for  
colorectal

## A sutureless overlapped anastomosis technique using linear staplers with reinforced bioabsorbable material in robotic right colectomy with intracorporeal anastomosis

Mori S. et al.

*Colorectal Disease* (2024) 26:4 (754-759)



The *objective* is to present a simple method for constructing a sutureless overlapped anastomosis using a 60 mm linear stapler with a reinforced bioabsorbable material in RRC with intracorporeal anastomosis.



### Results

- *Length of hospital stay*  
10 days (range 8-15 days)
- *Operation time*  
Median  
281 min (range 228-459 min)
- *Anastomosis construction time*  
12 min (range 11-17 min)
- *Intraoperative adverse events*  
No



### Insight

- **Sutureless technique** appears safe and feasible, offering a potentially faster alternative to traditional handsewn methods for creating anastomosis during RRC with IA

Technical study  
10 pts

Indication: Robotic right  
colectomy (RRC) tumours of  
the caecum

Signia Stapling Technology;  
Endo GIA 60 mm  
Articulating Medium Thick  
Reload with Tri-Staple™  
Technology (Medtronic)

Acronyms

  
PubMed link



Buttress reload  
use for  
colorectal

## Efficacy of the polyglycolic acid sheet for preventing anastomotic leakage in double-stapling technique anastomosis for left-sided colon or rectal cancer surgery: a propensity score-matched study

Shibutani M. et al.

*BMC surgery* (2023) 23:1 (135)



The *objective* is to evaluate the effect of the PGA sheet on preventing anastomotic leakage in patients with left-sided colorectal cancer who underwent DST anastomosis



### Results

- *Anastomotic bleeding*  
PGA: None  
Control: 5 (1.6%)  
No statistical difference

Propensity score matching  
PGA: 0 (0%)  
Control: 0 (0%)  
p NA

- *Anastomotic leakage/air leakage incidence*  
PGA: 1 (2.3%)  
Control: 42 (13.4%)  
p = 0.042

Propensity score matching  
PGA: 1/43 (2.3%)  
Control: 8/43 (18.6%)  
p = 0.030

- *Re-operation*  
PGA: None  
Control: 2.6%  
p = 0.603

Propensity score matching  
PGA: 0/43 (0%)  
Control: 1/43 (2.3%)  
p > 0.999



### Insight

- **PGA sheet** is an **easy-to-perform** and **cost-effective** method to strengthen the anastomosis site and reduce the risk of this common and serious surgical complication

Retrospective study  
356  
PGA: 43 pts  
Control: 313 pts

Indication: Left-sided  
colorectal cancer

PGA: Neoveil

Buttress reload  
use for  
colorectal

## A new overlapped delta-shaped anastomosis technique using linear staplers with reinforced bioabsorbable material for intracorporeal anastomosis during laparoscopic colectomy

Hokonohara K. et al.

*Colorectal Disease* (2022) 24:11 (1427-1429)



The *objective* is to evaluate a construction technique for an overlapped delta-shaped anastomosis using a linear stapler with a reinforced bioabsorbable material.



### Results

- *Length of hospital stay*  
12 days (range 9-15 days)
- *Operation time*  
Median 329 min (range 285-682 min)
- *Anastomosis construction time*  
19 min (range 14-29 min)
- *Perioperative complications*  
No



### Insight

- The aim of this technique is to minimize surgical trauma and wound complications compared to traditional methods.
- The results show the **feasibility** and **safety** of the new **approach** in a small group of patients, suggesting it could lead to quicker recovery and better outcomes.

Technical study  
5

Indication: laparoscopic  
colectomy with  
intracorporeal anastomosis  
Colon tumor

Endo GIA 60 mm  
Articulating Medium Thick  
Reload with Tri-Staple  
Technology (Medtronic)

Acronyms

  
PubMed link



Buttress reload  
use for  
colorectal

## Prevention of anastomotic leakage using a polyglycolic acid sheet in double-stapling technique anastomosis for rectal surgery

Shibutani M. et al.

*Ann Med Surg (Lond)*. 2021;72:103117



The *objective* is to examine whether or not a polyglycolic acid (PGA) sheet is effective for reinforcing rectal anastomosis.



### Results

- *Anastomotic leakage/air leakage incidence*  
None
- *Complications*  
None
- *Anastomotic stenosis*  
None



### Insight

- **PGA sheets** may be an **effective** strategy for preventing **anastomotic leakage** in double-stapling technique anastomosis for rectal surgery.
- **PGA sheets** may **reduce** the risk of this dangerous post-operative issue

Retrospective study  
15 pts

Indication: Rectal surgery

PGA: Neoveil

Buttress reload  
use for  
colorectal

## Safety and efficacy of a novel linear staple device with bioabsorbable polyglycolic acid felt in laparoscopic colorectal surgery

Naito M. et al.

*Asian journal of endoscopic surgery (2017) 10:1 (35-39)*



The *objective* is to assess an automatic suture device with PGA felt to reinforce staple lines during gastrointestinal anastomosis, aiming to reduce postoperative anastomotic complications like leakage.



### Results

- *Anastomotic bleeding*  
None
- *Operation time*  
Median  
159min (range, 113-322min)
- *Post-op Complications*  
None



### Insight

- The **automatic PGA-felt suture device** is **safe** and **effective for colorectal resections and anastomoses**, showing no postoperative bleeding or significant complications.

Retrospective study  
17 pts

Indication: Colorectal  
cancer

PGA: Neoveil

Endo GIA Reinforced Reload with  
Tri-Staple Technology

DST Series  
EEA Stapler (Medtronic)

Buttress reload  
use for  
colorectal

# Does Tri-Staple™ technology reduce the risk of anastomotic leakage in colorectal surgery? A propensity score matching analysis

Gerbasi LS et al.

*J Surg Oncol.* 2024 Sep;130(4):830-837



The *objective* is to evaluate the Tri-Staple™ technology in colorectal anastomosis



## Results

- *Anastomotic leakage/air leakage incidence*

EEA: 11%

Tri-EEA: 4%

p = 0.023

Propensity score matching

EEA: 12%

Tri-EEA: 4%

p = 0.04

- *Length of hospital stay*

EEA: 6.67 (0.66)

Tri-EEA: 5.3 (0.44)

p = 0.085

Propensity score matching

EEA: 6.68 (1.22)

Tri-EEA: 5.3 (0.44)

p = 0.110

- *Re-operations*

EEA: 12%

Tri-EEA: 4%

p = 0.016

Propensity score matching

EEA: 14%

Tri-EEA: 4%

p = 0.014

- *Severe complications*

EEA: 14%

Tri-EEA: 6%

p = 0.026

Propensity score matching

EEA: 15%

Tri-EEA: 6%

p = 0.041



## Insight

- **Tri-Staple™** is associated with a **lower incidence of leakage, reoperations, and severe complications**, even after propensity score matching

336 pts  
EEA: 228  
Tri-EEA: 108

*Propensity score matching*  
DST™ Series  
technology (Medtronic)  
circular stapler (EEA): 108  
EEA Tri-Staple™  
(Medtronic) circular stapler  
(Tri-EEA): 108

Acronyms

  
PubMed link



Buttress reload  
use

# Reduction of air leakage using linear staple device with bioabsorbable polyglycolic acid felt for pulmonary lobectomy

Deguchi H et al.

*Gen Thorac Cardiovasc Surg.* 2020 Mar;68(3):266-272



The *objective* is to investigate the effectiveness of using a linear staple device with bioabsorbable polyglycolic acid (PGA) in reducing air leakage after pulmonary lobectomy for non-small cell lung cancer (NSCLC)



## Results

- *Post-operative air leakage*  
Stapler: 28/125 (22.4%)  
Stapler-PGA: 12/125 (9.6%)  
P= 0.006
- *Length of stay (days)*  
Stapler: 10.7 ± 7.2  
Stapler-PGA: 11.5 ± 12.4  
P= 0.327
- *Operation time (min)*  
Stapler: 248.9 ± 66.4  
Stapler-PGA: 257.3 ± 53.8  
P= 0.055
- *Post-operative complications*  
Stapler: 35/125 (28.0%)  
Stapler-PGA: 31/125 (24.8%)  
P= 0.566



## Insight

- **Stapler with PGA significantly reduced** postoperative air leakage.
- **PGA-reinforced stapler** is a valuable technique for **minimizing complications** and **decreasing the need for additional interventions** like fibrin glue.

Retrospective study  
463 pts

Stapler group: 135 pts  
Stapler-PGA group: 328 pts

*Propensity score matching*  
Stapler group: 125 pts  
Stapler-PGA group: 125 pts

Indication: non-small cell  
lung cancer (NSCLC)

PGA: Neoveil and Gunze

## Buttress reload use

# A real-world study evaluating the safety and utility of a two-row stapler reload on pulmonary vasculature

Sandeep J. Khandhar et al.

*J Thorac Dis.* 2024;16(6):3753-3763.



The *objective* is to evaluate the safety and effectiveness of a two-row stapler reload (Signia™ Small Diameter Reload or SDR) when used on pulmonary vasculature during thoracic surgical procedures.



## Results

- *Intraoperative hemostatic intervention rate*  
0.99% (3/302)

- *Numbers of fires*  
302; 48% of firings on arteries, 52% on veins  
Reinforcement: applied to 14/302 (4.63%)

None of these procedures required conversion to an open procedure, a blood transfusion, suture placement, or the need for more extensive pulmonary parenchymal resection.

Likert scale  
Integrous firings: 99%  
Acceptable: 299  
no intervention: 283

Manual handle: 60%  
Powered handle: 40%

- *Adverse events*  
No device-related injuries  
Possibly device-related: 3/44 pts (6.8%)



## Insight

- **Signia™ Small Diameter Reloads** demonstrates favorable **safety** and **efficacy**, with a low rate of hemostatic intervention

Prospective study  
120 pts

Intervention: pulmonary  
vasculature

Signia™ Small Diameter  
Reloads



Home

SECTION 3

# Environmental impact



Acronyms

Article	Conclusions
<p><a href="#">Bischofberger et al. 2023</a> Assessing the environmental impact of an anastomotic leak care pathway</p>	<p>The findings suggest that efforts to reduce the incidence and severity of surgical complications like anastomotic leak may lead to improvements in the overall environmental sustainability of surgical departments</p>
<p><a href="#">Meissner et al. 2023</a> Evaluating the Environmental Impact of Single-Use (SUS) and Multi-Use Surgical Staplers (MUS) with Staple Line Buttressing in Laparoscopic Bariatric Surgery</p>	<p>The survey indicates that using MUS with pre-attached buttressing significantly reduces waste, resource consumption, and greenhouse gas emissions compared to using SUS with separate buttressing in bariatric surgery</p>
<p><a href="#">Meissner et al. 2021</a> Evaluating the Waste Prevention Potential of a Multi-versus Single-Use Surgical Stapler.</p>	<p>The study suggests that switching from single-use to multi-use surgical stapling systems can substantially reduce waste and resource consumption in surgical procedures, provided that the multi-use components are adequately reused</p>



PubMed link



PubMed link



PubMed link

## Environmental impact

# Assessing the environmental impact of an anastomotic leak care pathway

Bischofberger S. et al.

*Surg Open Sci.* 2023;14:81-86



The *objective* is to assess the environmental impact associated with the care pathway of anastomotic leak (AL)



## Insight

- **Anastomotic leak (AL) patients** have an **increased length of hospital stay**, spending an **additional 8.36 days** in general inpatient care and **1.50 days** in intensive care, on average, compared to those who do not experience leaks.
- The **treatment of AL** is **associated with** a substantial **environmental impact**. The average environmental impact per patient is **1303 kg CO<sub>2</sub>-eq** for [climate impact](#), **1803 m<sup>3</sup>** of water for [water impact](#), and **123 kg** of waste for [waste impact](#).
- **Stoma home management** is the **largest contributor** to the **climate impact** of **AL patients**, accounting for **46.6%** (607 kg CO<sub>2</sub>-eq) of the total, and **47.3%** (58.1 kg) of the total waste impact. The majority of the climate and water impact of stoma home management is associated with the manufacture of ostomy bag material.
- **In-patient days** make the **largest contribution** to the total **water impact**, accounting for **46.7%** (841.1 m<sup>3</sup>) of the total water footprint of AL patients. The environmental impact of in-patient care is largely associated with the use of energy and consumables.
- The findings suggest that **efforts to reduce the incidence and severity of surgical complications like AL** may lead to **improvements in the overall environmental sustainability** of surgical departments.




Fig. 1. Care pathway for the treatment of patients with AL. Based on 100 patients experiencing AL.





Environmental  
impact

# Evaluating the Environmental Impact of Single-Use and Multi-Use Surgical Staplers with Staple Line

## Buttressing in Laparoscopic Bariatric Surgery

Meissner M. et al.

*Risk Manag Healthc Policy. 2023;16:1423-1433*



The *objective* is to assess the environmental impact of single-use versus multi-use surgical staplers in bariatric surgery, focusing on staple line buttressing.

Survey of bariatric surgeons



### Results

- **Environmental Awareness**  
Concerned: 12 (67%)  
Hospital promotes environmentally friendly purchasing: 22%
- **Priorities in Purchasing Decisions**  
Never made changes to their work practice based on environmental concerns: 50%  
Prioritize ease of use and price over environmental impact: 61%
- **Surgeon Practices**  
Staple line reinforcement: 12/18 (67%)  
Powered staplers over manual staplers: 88%
- **Cartridge Usage**  
Sleeve gastrectomy (SG): 0.6 45mm cartridges and 5.0 60mm cartridges.  
Roux-en-Y gastric bypass (RYGB): 0.2 30mm cartridges, 1.0 45mm cartridges, and 5.0 60mm cartridges.  
Bariatric surgery: 60mm cartridges
- **Waste Generation**  
MUS reduce product waste by 40% in both sleeve gastrectomy (SG) and Roux-en-Y gastric bypass (RYGB) procedures.  
MUS reduce packaging waste by 60% in SG and 57% in RYGB.  
Overall, MUS reduce total waste by 50% for SG and 48% for RYGB
- **Total Material Requirement (TMR)**  
The TMR is reduced by more than 90% when using MUS compared to SUS for both SG and RYGB procedures
- **Lithium Content and Greenhouse Gas Emissions**  
Lithium  
MUS: 4.96%  
SUS: 6.72%  
The greenhouse gas impact  
MUS: 0.018 g CO<sub>2</sub>eq  
SUS: 5.904 g CO<sub>2</sub>eq
- **Sensitivity Analysis**  
If buttressing is applied to all staple lines, the waste reduction offered by the MUS increases to 1,022.13 g (54.7%) compared to 649.88 g (45.6%) if no staple line reinforcement is used.



### Insight

- The survey highlights that while many surgeons are aware of the environmental issues associated with surgical practices, environmental concerns are not a primary driver in their purchasing decisions.
- Multi-use staplers (MUS) significantly reduce waste compared to single-use staplers (SUS).
- Preloaded buttressing in MUS produces less waste than separate buttressing used with SUS per stapler firing.
- Although MUS have a higher TMR per cartridge than SUS, the reusability of MUS components leads to a significant overall reduction in TMR.
- MUS contain less lithium than SUS
- Sensitivity analysis, shows that MUS generate less waste and have lower TMR than SUS, regardless of the rate of buttressing applied.

Acronyms

  
PubMed link



## Environmental impact

# Evaluating the Waste Prevention Potential of a Multi- versus Single-Use Surgical Stapler

Meissner M. et al.

*Risk Manag Healthc Policy. 2021;14:3911-3921*



The *objective* is to evaluate the waste prevention potential and extended resource use of multi-use surgical stapling systems (MUS) versus single-use surgical stapling systems (SUS).



## Results

- *Waste Generation*  
Sleeve gastrectomy: reduced by 0.29 kg  
Gastric bypass: reduced by 0.96 kg  
VATS lobectomy: reduced by 0.86 kg  
Waste when using the MUS compared to the SUS.
- *Total Material Requirement (TMR)*  
60 mm power handle and adapter (PSEE60A) of the SUS is the most resource-intensive component, requiring 328 kg of resources  
Power handle (SIGPHANDLE) of the MUS is the most resource-intensive component, requiring 560 kg of resources.  
  
Sleeve gastrectomy: decreased from 329 kg to 27 kg  
Gastric bypass: decreased from 633 kg to 25 kg  
VATS lobectomy: decreased from 633 kg to 34 kg  
when switching from a SUS to a MUS
- *Sensitivity Analysis*  
Single-use MUS scenario, waste generation increases by 1.13 kg for laparoscopic sleeve gastrectomy, 0.46 kg for laparoscopic gastric bypass, and 0.57 kg for VATS lobectomy when compared to SUS.  
To achieve a lower TMR with the MUS than the SUS, the multi-use MUS components must be used at least three times.



## Insight

- Use of the MUS leads to a reduction in the total amount of system waste accumulated per surgical procedure, compared to SUS. 40% for laparoscopic sleeve gastrectomy, 70% for laparoscopic gastric bypass, and 62% for VATS lobectomy procedures.
- TMR decreased significantly when switching from SUS to MUS. TMR reductions were higher, at 92% (sleeve gastrectomy), 96% (gastric bypass), and 95% (VATS lobectomy).
- If the multi-use components of the MUS are used more than once, the reduction in waste and TMR associated with a switch from SUSs to MUSs is maintained for a range of reuses.

Ethicon's SUS:  
ECHELON FLEX™

Medtronic's MUS: Signia™  
Stapling Technology

Mechanically deconstruction