# STORZ

## KARL STORZ — ENDOSKOPE

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HAMOU® ENDOMAT®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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**RxOnly** Federal (USA) law restricts this device to sale by or on the order of a physician.

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**It is recommended to check the suitability of the product for the intended procedure prior to use.**

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




**2. Controls, displays,  
connectors, and  
their uses**

- ① Power switch
- ② TFT touch screen
- ③ Display of field of application
  - ④ Irrigation [ml/min] (**laparoscopy**)
    - ⑤ Reduce the irrigation flow setpoint
    - ⑥ Bar display shows actual irrigation flow rate when the pump is activated (orange 'triangles' indicate setpoint)
    - ⑦ Irrigation flow setpoint
    - ⑧ Increase the irrigation flow setpoint
    - ⑨ Display of pressure limit [mmHg]
    - ⑩ Set pressure limit (100->300->500->100...)
    - ⑪ Irrigation Start/Stop button
  - ⑫ Irrigation [mmHg] (**hysteroscopy**)
    - ⑬ Reduce the irrigation pressure setpoint
    - ⑭ Bar display shows actual irrigation pressure rate when the pump is activated (orange 'triangles' indicate setpoint)
    - ⑮ Irrigation pressure setpoint
    - ⑯ Increase the irrigation pressure setpoint
    - ⑰ Display of flow limit [ml/min]
    - ⑱ Set flow limit (200->400->600->200...)
    - ⑲ Irrigation Start/Stop button
  - ⑲ Suction [bar]
    - ⑳ Reduce the suction pressure setpoint
    - ㉑ Bar display shows actual suction pressure rate when the pump is activated (orange 'triangles' indicate setpoint)

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## Controls, displays, connectors, and their uses

- ②② Increase the suction pressure setpoint
- ②③ Suction pressure setpoint
- ②④ Suction Start/Stop button
- ②⑤ Pump rollers
- ②⑥ Light barrier LEDs
- ②⑦ Pump lever (cartridge locking lever)
- ②⑧ Suction port
- ②⑨ SCB connectors\*
- ③⑩ Potential equalization connector
- ③① Power cord receptacle
- ③② Line fuses
- ③③ Disposable tubing cartridge
- ③④ Suction bottle (optional)
- ③⑤ Suction connection to the device (with bacterial filter)
- ③⑥ Suction connection to the instrument

\* The KARL STORZ SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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	Follow instructions for use
	ON
○	OFF
	Potential equalization connector
	Applied part of type BF
~	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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**Symbols on label and packaging of the device**

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316 DF. This can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).

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6.10	Presetting irrigation parameters .....	23
6.11	Preselecting suction pressure .....	24
6.12	Test for proper operation and air-relief of the tubing system.....	25
6.13	Initiating irrigation.....	26
6.14	Initiating suction .....	26
6.15	Settings .....	27
6.15.1	Audio.....	27
6.15.2	Language .....	27
<b>7.</b>	<b>Maintenance.....</b>	<b>29</b>
7.1	Fuse replacement .....	29
7.2	Reprocessing.....	30
7.2.1	Wipe-down disinfection of device .....	30
7.2.2	General warnings for reusable medical devices.....	30
7.2.3	Accessories for preparation .....	31
7.2.4	Disassembly .....	31
7.2.5	Preparation for cleaning and disinfection .....	32
7.2.6	Manual precleaning.....	32
7.2.7	Manual cleaning.....	32
7.2.8	Manual disinfection .....	33
7.2.9	Machine cleaning and disinfection .....	34
7.2.10	Connecting .....	34
7.2.11	Assembly, inspection and care .....	35
7.2.12	Packaging systems.....	35
7.2.13	Sterilization .....	35
7.2.14	Storage.....	36
7.3	Maintenance and safety check .....	37
7.3.1	Maintenance.....	37
7.3.2	Safety check/repeat inspection according to IEC 62353.....	37
7.4	Servicing and repair .....	38
7.5	Disposal.....	38
7.6	Repair program.....	39
7.7	Limitation of liability .....	40
7.8	Manufacturer's warranty .....	40
<b>8.</b>	<b>Technical description .....</b>	<b>41</b>
8.1	Information signals .....	41
8.1.1	Configuration .....	41
8.1.2	Visual signals .....	41
8.1.3	Acoustic signals.....	42

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## 4. General information

### 4.1 Device description

The HAMOU® ENDOMAT® is a combination irrigation/suction pump for use in gynecological and laparoscopic surgery. Use of the most advanced technologies allows its application in both hysteroscopic and laparoscopic surgery. The device automatically adapts to suit the type of operation being performed by providing optimal operating parameters when the correct tubing cartridge for hysteroscopy or laparoscopy is inserted. The newly developed tubing cartridge makes attachment to the pump considerably easier.

The device is operated and monitored via a touch display. A display for the setpoint and actual values enables optimal monitoring of the current operating status. An electronic safety circuit interrupts irrigation or suction if device performance persistently deviates from the setpoint, and emits an audible signal to inform the user. The high safety standard of the ENDOMAT® is further enhanced by an electronic auto-check system, which tests the system components each time the device is switched on, and notifies the operator of any failures detected.

### 4.2 Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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## 5. Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

### 5.1 Explanation of warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**NOTE:** A Note indicates special information about operating the product, or clarifies important information.

**WARNING:** Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the device, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**NOTE:** Any damage to the device resulting from incorrect operation is not covered by the manufacturer's warranty.

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## Safety instructions

### Warnings and cautions

**WARNING:** The electrical installations in the operating room in which the device is connected and operated must comply with the applicable IEC standards.

**WARNING:** Keep out of reach of patients.

**WARNING:** For safety reasons, when the system is in use the output sockets of the device and the patient must not be touched at the same time.

**WARNING:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**WARNING:** Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this (cf. Clause 16 of the 3rd Ed. of IEC 60601-1).

**NOTE:** The compliance of additional equipment connected to medical electrical equipment with the respective IEC or ISO standards must be verifiable. Furthermore all configurations shall comply with the standard requirements for medical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements.

**WARNING:** Grounding reliability can only be achieved when the equipment is connected to a "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

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**Safety instructions****Warnings and cautions**

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**WARNING:** Always unplug the device before cleaning.

**WARNING:** Always unplug the device before performing any maintenance work on it.

**WARNING:** Always use the lowest possible pressure for sufficient distension, irrigation or suction.

**WARNING:** Treatment may only be performed if visual observation of the action of the device is ensured.

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**WARNING:** The irrigation tubing cartridge should always be unpressurized when it is inserted, as insertion of a pressurized cartridge (membrane extended) may cause the pressure measurement to be falsified.

**WARNING:** Avoid allowing fluids to penetrate the device. Do not store liquids on or directly above the device.

If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

**WARNING:** Before you start the intervention, make sure that the tubing set is not leaking or kinked and that the ends of the tubes are properly secured

**WARNING:** Close the pump lever completely before starting the pump to ensure correct positioning of the cartridge.

**WARNING:** Connect the potential equalization cable to the potential equalization connector. The device's ground line should be installed by a qualified electrician.

**WARNING:** Do not reuse disposable tubing sets. Discard after use.

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**WARNING:** Do not touch the pressure sensor membrane.

**WARNING:** Do not touch the pump reels when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

**WARNING:** Ensure sufficient air circulation for the device.

**WARNING:** It is not permitted to press more than one touch point simultaneously.

**WARNING:** Make sure that you use the cartridges only for applications specified on their packaging.

**WARNING:** Observe the applicable local regulations when disposing of disposable tubing.

**WARNING:** Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

**WARNING:** Only operate the device with the voltage stated on the manufacturer's identification plate.

**WARNING:** Only use fuses of the correct rating. Replace defective fuses and plug the line fuse holder back in.

**WARNING:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, some checks (see chapter 7.3) must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

**WARNING:** Risk of explosion if used in the presence of flammable anesthetics.

**WARNING:** Servicing may only be performed by specialists authorized by KARL STORZ.

**WARNING:** Test this device prior to each surgical procedure to ensure that it is functioning correctly. It must not be used if any damage is evident.

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**WARNING:** The connection of the irrigation tubing set (sheath) on the patient side must remain in the sterile area.

**WARNING:** Try to keep the height difference between device and patient as low as possible. The device controls the pressure measured at the outlet of the irrigation pump. However, the intracavitary pressure depends, among other things, on the height difference between the device and the cavity.

**WARNING:** The device may not be operated in the demarcated hazard zone wherever highly inflammable inhalation anesthetics and their mixtures are used (see 5.10 Safety precautions at the site of installation).

**WARNING:** The device may only be operated under the environmental conditions specified.

**WARNING:** The device may only be operated with the power cable delivered by KARL STORZ or a similar power cable which has a national inspection seal.

**WARNING:** The device must not be operated if the touch screen or the display is defective.

**WARNING:** The device must not be used in oxygenated environments.

**WARNING:** The fluid level in the irrigation bag must be monitored carefully and be exchanged if appropriate. If the irrigation bag runs empty it is possible for air to be pumped into the cavity.

**WARNING:** The fluid level in the suction bottle must be monitored carefully. The suction bottle must be emptied or exchanged in due time to ensure that the vacuum pump functions correctly.

**WARNING:** The device may be subject to sudden failure during use. Keep a replacement device at hand or alternative procedures in mind if the application requires permanent fluid transportation.

**WARNING:** The intracavitary pressure is always equal to or lower than the pressure measured at the outlet of the device, which is used to regulate the device.

**WARNING:** Position the irrigation bag to the side of and above the device. The irrigation pump is not self-priming.

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**WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

**WARNING:** The one-day tubing cartridge for irrigation can remain in the reel pump for a maximum of one surgery day (12h), during which it can be used for several operations. Afterwards, it must be disposed of.

**WARNING:** Check the packaging of the tubing set for damage. Damaged packaging may impair the sterility of the tubing set.

**WARNING:** Check the tubing set for damage etc. before using it.

**WARNING:** The patient tube of a one-day tubing cartridge must be disposed of after each intervention, while the cartridge may be reused.

**WARNING:** The pump tube (inside the cartridge) is a wear part and must be replaced after five interventions.

**WARNING:** The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

**WARNING:** The tubing set must be filled with irrigation fluid before use on a patient.

**WARNING:** The user interface includes a periodically changing signal, which indicates that the device is working properly. If the signal stops changing, the device must not be operated further.

**NOTE:** A visual information message may cover the signal.

**WARNING:** The patient tube connector must be screwed securely into the pump tube coupling to avoid a leak.

**WARNING:** To avoid the risk of electrical shock, do not open the device! Refer servicing to the manufacturer or to personnel authorized by the manufacturer. Removal of covers by unauthorized personnel will void the equipment's warranty.

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**WARNING:** To guarantee the sterility of the **connected** tubing cartridge during patient changeover, the patient tubes are equipped with a separately packed, sterile plug. This must be screwed into the tube coupling immediately after removal of the patient tube from the cartridge and remain there until the next patient tube is connected.

**WARNING:** Use only the approved sterilization methods for the reusable tubing set.

**WARNING:** Valves and tube clamps should be either completely opened or completely closed.

**WARNING:** When cleaning the device surface, the area around the pressure sensor should be wiped carefully due to the sensor's sensitive surface.

**WARNING:** When cleaning the device, always observe the instructions in the instruction manual (section 'Cleaning and Care')

**WARNING: When used in combination with RF surgical devices:** Only RF operating modes with a maximum output voltage of 5.5 kVpp may be selected.

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Before using the device on the patient, it is imperative that you be acquainted with how the device operates and is controlled.

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## 5.2 Intended use

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs and operating fields as well as the suctioning off of irrigation and bodily fluids, secretions and tissue during endoscopic interventions.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

### Contraindication

Suction/irrigation pumps without a pressure monitoring system must not be used for the dilation of hollow organs and joints. Serious injuries to the patient may result if this use exclusion is not observed.

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The HAMOU® ENDOMAT® may only be operated with tubing cartridges or tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

No responsibility for the safe operation of the instrument can be accepted if any tube systems other than those specified are used. The following serious defects may arise in particular:

- Inadequate grip on the tube connectors means a loss of sterility and may lead to the device being damaged by ingress of liquid.
- Premature breakage of the pump tube segment and consequent loss of sterility and damage to the device caused by ingress of liquid.
- Display of an incorrect actual irrigation rate due to deviations in tube diameter.

**WARNING:** *Unauthorized conversions or modifications to the device are not allowed for safety reasons.*

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**5.3 User qualification**

The HAMOU® ENDOMAT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the device.

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**5.4 Training in the operation and function of the device**

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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**5.5 Patient profile**
**5.5.1 Laparoscopy**

Gender: irrelevant  
Age: irrelevant  
Weight: irrelevant  
Nationality: irrelevant  
State of health: suitable for endoscopic interventions, in the physician's judgment

**5.5.2 Gynecology**

Gender: female  
Age: irrelevant  
Weight: irrelevant  
Nationality: irrelevant  
State of health: suitable for endoscopic interventions, in the physician's judgment

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**5.6 Intended region of treatment on the body**

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**5.7 User profile of physician and assistants**

Training: recognized medical skills in the field of application (specialist physician, qualified medical staff)  
Intelligence: adequate powers of comprehension to rationally assess the situation in hand  
Ability to act: familiarity with the use of technical apparatus  
Language: Suitable command of one of the languages used by the device and in the instruction manual.  
Experience: at least one comprehensive session of training in the use of the device.  
Impairments: no physical impairment which makes it difficult or impossible to perceive acoustic or visual signals.

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**5.8 Intended conditions of use**

**5.8.1 Use**

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual.

**5.8.2 Other intended conditions**

Frequency of use: one or more times a day  
Length of use: from a few minutes to several hours a day.  
Place of installation: Positioning on a level, vibration-free surface  
Mobility: can be moved if positioned on a cart  
Combination: can be used on the patient at the same time as other equipment required for the operations  
Control: can be controlled via the KARL STORZ SCB.

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
**5.9 User position**

For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of  $\pm 45^\circ$  at a distance of approx. 30 – 70 cm to the front panel.

For observation of the actual values during the application, a visual distance from the device of 2 m is assumed (the tube length is 2 m).

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**5.10 Safety precautions at the site of installation**

The device may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations. It is not intended for use in hazardous zones.

This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the device must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The device is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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**5.11 Safety precautions when operating the device**

It is the user's responsibility to make sure the device is safe and operates properly before using it.

During treatment with the HAMOU® ENDOMAT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic effects.

Any treatment may only be performed if visual observation of the device is ensured.

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**5.12 Safety features**

The HAMOU® ENDOMAT® is equipped with the following safety features:

**5.12.1 Self-test**

A self-test is performed each time the device is switched on. If any of the tests is negative, a corresponding error message is displayed.

**WARNING:** Before starting an operation, check whether the next self-test is pending, in order to avoid an inoperable pump, and restart the device if necessary. The device must perform a self-test every 24 hours to ensure safe operation.

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**5.12.2 Testing the touch display**

After switching on the supply voltage, a KARL STORZ logo with the device name is displayed for a few seconds, depending on the loading time of the operating system.

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**5.12.3 Monitoring the irrigation and suction pressure**

An electronic safety circuit interrupts irrigation or suction if device performance persistently deviates from the setpoint.

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## 6. Installation and operating instructions

### 6.1 Unpacking the equipment

Carefully unpack the HAMOU® ENDOMAT® and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

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
### 6.2 Basic equipment

KARL STORZ Set No.: 263311 01-1	
263311 20-1	HAMOU® ENDOMAT® SCB incl. Mains cord
20 0901 70	SCB Connecting Cable, 100 cm
031523-03	Cartridge tubing set, PC, for single use
031524-03	Cartridge tubing set, FC, for single use
Instruction manual with wallet	

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### 6.3 Installation

**NOTE:** The HAMOU® ENDOMAT® including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

**! WARNING:** The device is not intended for use in hazardous zones. Do not operate the device within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Set the device on a horizontal, flat surface.

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The device is equipped with a connector ⑳ for attaching a ground line. The device's ground line should be installed by a qualified electrician.

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Connect power cord and insert the power plug fully into the power cord receptacle ㉑.

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-  **WARNING:** Operate the device only with the system voltage stated on the rating plate.
-  **WARNING:** Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.
-  **WARNING:** Ensure that access to the power plug is ensured at the site of installation.

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#### 6.4 KARL STORZ SCB

**i NOTE:** To avoid mistakenly pulling out the SCB connecting cable, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets 29. Connect the other end of the cable to the KARL STORZ SCB (KARL STORZ Communication Bus) control device or other SCB devices (see KARL STORZ SCB control NEO System Instruction Manual).

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#### 6.5 Field of application

The HAMOU® ENDOMAT® device is suitable for both hysteroscopic and laparoscopic applications. The required field of application is automatically selected when the correct tubing cartridge for the respective field of application is inserted. The device determines the intended field of application from the different 'encodings' of the tubing cartridge, and then automatically selects the appropriate operating parameters.

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#### 6.6 Operating the device

Switch on the device at its power switch 1 (green indicator light illuminates). The device will then conduct a self-test, checking all system components for proper operation.

**i NOTE:** Once the self-test at power-up has been successfully completed, an acoustic availability signal is given.

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**6.7 Inserting the tubing cartridge**

**CAUTION:** Only use the tubing cartridges or tubing sets designed for the HAMOU® ENDOMAT®.

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Insertion of the tubing cartridge is explained in an animation. Once all of the conditions have been fulfilled, i.e. a tubing cartridge with a valid coding (HYS or LAP) has been inserted and the pump lever has been closed, the device jumps to the appropriate working level (hysteroscopy or laparoscopy).

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**Disposable or one-day tubing cartridge**

Insert the tubing cartridge and turn the lever counterclockwise.

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**i NOTE:** Reusable tubing cartridges are also available. The instruction manual supplied with the reusable tubing cartridges must be observed.

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
Connect the correct end of the length of tubing to the irrigation-liquid bag (puncture needle) or to the irrigation port on the instrument (LUER lock).

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The field of application preselected by means of the tubing cartridge (in this case 'Laparoscopy' ③) is indicated on the touch display.

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### 6.8 Removing the tubing cartridge

Turn the lever clockwise to the 9 o'clock position and remove the tubing cartridge.

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**6.9 Connecting the tubing set for suction (optional)**

Attach the bacterial filter (see Recommended Accessories) to the vacuum port ③⑤ on the suction bottle and then secure the suction tube to the filter housing.


**CAUTION:** To avoid contaminating the device, a bacterial filter must be fitted between the suction bottle and the device.

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Connect the other end of the suction tube to the device port ②⑧.

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
Use a silicone suction tube to connect the suction port of the suction bottle ③⑥ and the suction port on the instrument.

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### 6.10 Presetting irrigation parameters

The setpoint (digital display – element ⑦\*/⑮\*\*) and actual value (bar display – element ⑥) for the irrigation flow rate\* and ⑭ for the irrigation pressure\*\* are displayed in the operating menu. In addition, the pressure limit\* ⑨ and irrigation flow rate limit\*\* ⑰ are displayed.

The roller pump is activated and deactivated with the associated Start/Stop button (button ⑪).


- \* Field of application – Laparoscopy
- \*\* Field of application – Hysteroscopy

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The setpoint can be adjusted in 100 ml/min steps (field of application – laparoscopy) using the ± buttons ⑤ and ⑧ and in 10 mmHg steps using the ± buttons ⑬ and ⑯ (field of application – hysteroscopy).

This is not affected by the pump status.

Using the ⑩ button, the pressure limit (field of application – laparoscopy) can be set within the range of 100 mmHg --> 300 mmHg --> 500 mmHg --> 100 mmHg...


Using the ⑰ button, the flow limit (field of application – hysteroscopy) can be set within the range of 200 ml/min --> 400 ml/min --> 600 ml/min --> 200 ml/min...

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
When the pump has been activated by means of the Start/Stop button ⑪ the bar indicates the actual irrigation flow rate\* ⑥ or the actual irrigation pressure ⑭. The current setpoint value is indicated by two orange triangles in the bar display. The bar display reflects the entire setting range (Hysteroscopy: 20...200 mmHg and Laparoscopy: 100...1300 ml/min).

**NOTE:** The settings are retained in the device's memory after the device is switched off.

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### 6.11 Preselecting suction pressure


The setpoint (digital display – element 23) and actual value (bar display – element 21) for the negative pressure for the suction are also displayed.

The suction pump is activated and deactivated with the associated Start/Stop button 24.

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Using the ± buttons 20 and 22, the setpoint is adjusted in 0.1 bar steps. This is not affected by the pump status.

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When the suction pump has been activated (by means of the Start/Stop button 24), the actual value for the negative pressure is indicated by a white bar 21. The current setpoint value is indicated by two orange triangles in the bar display. The bar display reflects the entire setting range (Laparoscopy/Hysteroscopy: -0.1... -0.8 bar).

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
**6.12 Test for proper operation  
and air-relief of the tubing  
system**

**WARNING:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

**NOTE:** Have a collecting container ready to catch irrigation liquid as it runs out.

Press the Start/Stop button (11) and let the irrigation pump run until all of the air has been released from the tubing system (i.e. no bubbles).

**Laparoscopy:**

The irrigation flow rate indicated on the bar display (6) must agree with the flow-rate setpoint (see 'triangles'). As soon as the inflow valve on the instrument is closed, the roller pump must stop irrigation.

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**Hysteroscopy:**

Slowly close the inflow valve on the instrument. As soon as the pressure indicated on the bar display (14) rises above the setpoint (see 'triangles'), the roller pump must stop irrigation.

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
Press the suction Start/Stop button (24). Close off the suction tube (instrument tube) with your finger and check whether a vacuum forms. The suction pressure indicated on the bar display (21) must agree with the setpoint (see 'triangles').

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### 6.13 Initiating irrigation

Press the Start/Stop button (11) to initiate irrigation.

**Laparoscopy:** The actual irrigation flow rate is shown continuously on the bar display (6).

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
**Hysteroscopy:** The actual irrigation pressure is shown continuously on the bar display (14).

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### 6.14 Initiating suction

Press the Start/Stop button (15) to initiate suction. The actual suction pressure is shown continuously on the bar display (21).

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
## 6.15 Settings

Press the 'Setting wrench' to call up the 'Settings' window.

In the 'Settings' window, you can select various setting and management fields.

- Audio
- Language
- Device Info
- Event Log
- Service (only authorized service technicians)

### 6.15.1 Audio

To switch off the acoustic signals of the information messages, the tick must be removed.

**NOTE:** When the device is powered up, the acoustic signal function is always active.

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
### 6.15.2 Language


Press the 'Language' button in the 'Settings' window to go to the 'Language' window. In this window, you can select the screen language. The following languages are available via individual buttons: German, English, Spanish, French, Italian, Portuguese, Polish and Russian.

Select a new language, if required, and confirm with Enter. The 'Settings' window will appear again. If you do not choose a new language, only the 'Cancel' button is available.

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
Press the 'Device Info' button in the 'Settings' window to go to the 'About' window.

In this window, you will find information about your device (e.g. the serial number, software version, and operating hours).

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
Press the 'License Info' button to go to the 'Copyrights' window.

In this window, the 'Third Party Software' is listed.

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Press the 'Event Log' button in the 'Settings' window to go to the 'Event Log' window.

Here you will find the saved information messages. The maximum possible number of entries is 1,000. In the scroll bar, the current page number out of the maximum pages present is given (e.g. 1/59). Each entry occupies one line in the list and consists of the date, time, and info ID. The most recent entry is in the top line on page 1.

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Press the 'Service' button in the 'Settings' window to go to the 'Service Password' window.

**NOTE:** The service field is password-protected for security reasons.

Enter the four-digit pin using the number buttons. Entry of the password is concealed, i.e. only asterisks appear in the entry field in place of the numbers entered. After entering the (valid) four-digit pin and confirming with the 'Enter' button, the 'Service' window will appear.

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**NOTE:** The 'Delete' button (see arrow) deletes all of the numbers entered.

'Cancel' calls up the 'Settings' window again. The settings that are possible in the Service Menu are described in the Service Manual.

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**7. Maintenance**

**7.1 Fuse replacement**

**!** **WARNING:** The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Switch off the device and remove the power plug from the electrical outlet.

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Remove the fuse holder ② with a screwdriver or other suitable tool.

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**!** **WARNING:** Only use fuses of the correct rating.

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Refit the fuse holder ②.  
Reconnect the power cord.  
Test the device for proper operation.

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**7.2 Reprocessing**

**7.2.1 Wipe-down disinfection of device**

**WARNING:** Always disconnect from the power supply before cleaning!

**CAUTION:** Make absolutely sure that no liquid can get into the device housing.

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

**7.2.2 General warnings for reusable medical devices**

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

**WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

**WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

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**CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions paying close attention to proper concentration, exposure time, and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**CAUTION:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from [hygiene@karlstorz.com](mailto:hygiene@karlstorz.com).

**CAUTION:** National laws and regulations must be observed.

**NOTE:** The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

**Reusable tubing**

**WARNING:** Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

**7.2.3 Accessories for preparation**

Accessories required for carrying out preparation:

Brushes: Art. 27652 / 27651AL

Miscellaneous: Art. 27660

**7.2.4 Disassembly**

The medical device must be disassembled into its individual components as far as possible or opened.

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**7.2.5 Preparation for cleaning and disinfection**

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

**7.2.6 Manual precleaning****Brushing the surfaces**

Any heavy soiling of the surfaces must be completely removed by cleaning them mechanically under cold running water with the aid of a brush or sponge.

**Brushing the lumina**

In order to remove all visible contamination, the working channel, lumina and hollow spaces of the medical device must be precleaned under cold running water with the aid of a suitable brush.

**Rinsing out with a pressure water gun**

The working channel, lumina and hollow spaces of the medical device must be rinsed using a pressure water gun and a suitable rinsing attachment. Use 3 x 5 second blasts from a pulsed jet of water with a pressure of at least 2.5 bar for rinsing.

**7.2.7 Manual cleaning**

The medical device must be completely immersed in a cleaning solution. Through disassembling and opening the instrument, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

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**7.2.8 Manual disinfection**

The medical device must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, it must be ensured that even surfaces with restricted access are covered and that no air bubbles are present. Following the exposure time, the medical device must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces, channels, and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. No. 27660) is ideal for this purpose.

**Suction bottle and cap**

Carefully pull the tubes from the suction lances and take the bottle out of the holder.

**i NOTE:** *Prompt transportation to the preparation station and a subsequent decontamination process must be ensured.*

**Preparation station / CSSD:** Carefully remove the cap from the suction bottle, empty the contents into the sink and rinse out, repeat rinsing if necessary.

Rinse out the suction bottle with a disinfectant solution and then fill it with the remaining disinfectant solution. Immerse the cap in a plastic basin with disinfectant solution making sure that all air bubbles have been removed and, if necessary, rinse out/fill the attachments/connectors with disinfectant solution using an injection aid.

For further procedures see section 7.2.8

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**7.2.9 Machine cleaning and disinfection**

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments":

- Machine cleaning
- Thermal disinfection

Thermal disinfection is preferable. The relevant national requirements and the A<sub>0</sub> value must be taken into account when using this method.

A suitable slide-in tray or instrument holder, to ensure that the medical device is thoroughly rinsed out or through, must be selected following consultation with the manufacturer of the device.

**i NOTE:** If necessary, the instrument must be dried off afterwards by hand.

**7.2.10 Connecting**

In order to ensure effective machine cleaning and disinfection, the instrument must be connected up to the washer and disinfecter. It must be ensured that the instrument is rinsed out by means of appropriate rinsing lances

**Suction Bottle**

**i NOTE:** The bottles must be rinsed out once and precleaned before they can be loaded into the washer and disinfecter.

Carefully place the suction bottle in a suitable tray holder or connect up (to a rinsing jet). Contact between the glass and other medical devices which could result in damage to the material must be prevented during the decontamination program. It is recommended that the suction lances are connected to the rinsing tubes of the washer and disinfecter.

Start the program (93 °C device-specific, -alkaline cleaner pH value ≤ 11.5).

Perform a visual inspection of the suction bottle after machine preparation.

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**7.2.11 Assembly, inspection and care**

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness.

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

**i NOTE:** During care procedures, use items from the catalog "Care, Sterilization and Storage Techniques".

**Reusable tubing**

**! WARNUNG:** Check the entire tubing set for leaks after preparation. This makes it possible to detect a leak early, through which liquid could escape and enter the device.

To check for leaks use leakage tester 13242 XL, for example. To do this, connect one tube end to the leak tester and block off the other end, with your finger, for example. If there is a leak in the tubing set, a pressure drop will be observed on the pressure gauge of the leakage tester.

**7.2.12 Packaging systems**

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, DIN EN ISO 11607 Parts 1 + 2, DIN 58953).

**7.2.13 Sterilization**

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments". The procedure must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

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**KARL STORZ has validated and approved steam sterilization using the fractionated prevacuum procedure for this medical device**

**i NOTE:** Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C - 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

**i NOTE:** This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

- or 134°C - 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

**i NOTE for suction bottle and cap:** Only cleaned, disinfected, dry and functional suction bottles in suitable packaging may be sterilized. The suction bottle and cap can be steam sterilized, hold time at least 4 minutes, fractionated prevacuum method. To improve sterility efficiency, sterilize the suction bottle in its disassembled state. We recommend sterilizing the suction bottle separately, i.e. packed in sterile goods non-woven material/paper and separated from the other medical devices to be sterilized.

**Repreparation limits**

**WARNING:** For safety reasons, reusable tube cassettes must no longer be used if the type identification is no longer legible or if there are noticeable inconsistencies between the available procedures and the cassette type.

The end of the product's lifespan is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

**7.2.14 Storage**

Store the suction bottles and caps in suitable sterile goods packaging in a closed cupboard to protect them from dust, moisture and temperature fluctuations.

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### 7.3 **Maintenance and safety check**

#### 7.3.1 **Maintenance**

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

#### 7.3.2 **Safety check/repeat inspection according to IEC 62353**

**WARNING:** *Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/ repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.*

#### **Visual inspection**

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

#### **Electric measurements**

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure earth leakage current in accordance with IEC 62353: For limit values, please refer to the current standard.
- Measure touch current according to IEC 62353: For limit values, please refer to the current standard.
- Measure patient leakage current according to IEC 62353: For limit values, please refer to the current standard.

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**Test for proper functioning**

A functional inspection must be carried out in accordance with the instruction manual (see section 6.12).

**Documentation**

The safety check/repeat inspection and results must be documented.

**WARNING:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**NOTE:** For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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**7.4 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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**7.5 Disposal**

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

**WARNING:** At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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**7.6 Repair program**

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In Germany you can refer repairs directly to

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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**Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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**7.7 Limitation of liability**

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

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**7.8 Manufacturer's warranty**

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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### 8.1 Information signals

#### 8.1.1 Configuration

All of the information signals created, with the exception of those which display the cause of an inoperable device, ('Safe State' entry in the 'Reaction' column of the table on page 46 onwards) are not self-retaining. This means that the signal is only issued whilst the signal conditions are present (unless something different is defined under 'Reaction').

If the signal conditions are only very brief, there are rules as to the minimum visual/acoustic signal duration, with which the described information signals also comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by a signal with a higher priority, or a low priority information signal is suppressed for as long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized first will be signalized. Should signals deviate from this rule, this is noted in the corresponding entry form page 46 onwards.

#### 8.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, the minimum display duration is set at 5s.

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**8.1.3 Acoustic signals**


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. But at least one complete signal sequence is played.

The acoustic signal can be switched on or off using a button in the footer of the device's display.

But the acoustic signal can only be switched off temporarily by the user. When the line voltage is switched on, the setting is always 'Audio On'.

**8.1.4 Information signals (general)**

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated. In the case of messages which are related to the Safe State, the signal sequence is repeated every 15s, but otherwise it is just played once.

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


**8.1.5 Availability signal**


On completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by  $\pm 2$  Hz with a modulation frequency of 1.5 Hz. Five different harmonics are generated.

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**8.2 Information messages**

<b>Brief description (1 line with approx. 40 characters possible)</b>	<b>Reaction</b>
<b>Possible cause</b>	
42: Restored default settings Settings/date were reset due to a fault (one or more memory sectors)	None
44: Software version conflict When the power was switched on incompatible software versions for the user interface and process control were found	Safe state
48: Touchscreen short circuit Maximum touchscreen contact time of 25 s exceeded	Deactivation of pump(s), message remains until reactivated
4A: Defective drive path When the power supply was switched on a defective component was found in the irrigation drive path	Safe state
50: Pressure out of tolerance The value of the exposed sensor is outside the tolerance zone	Safe state
52: Irrigation pressure critical Value of the output pressure generated by the pumps exceeds 240 mmHg in HYS or 600 mmHg in LAP, while pump is still in rotation	Safe state
53: Irrigation pressure critical A pressure is measured on the sensor which could lead to its destruction	Safe state, selftest and recalibration of sensor is required
54: Pressure sensor defect Protection system detected a pressure difference between the two sensors of more than 30 mmHg	Safe state
60: Int. temperature high Temperature inside the unit is higher than 60°C	None, message remains until the temperature falls below the limit value or is superseded by message 61
61: Int. temperature critical Temperature inside the unit is outside the 10°C – 70°C range	Deactivation of pump(s), message remains until temperature falls below the limit value
62: Irrigation pump hot Temperature >70 °C	None, message remains until the temperature falls below the limit value or is superseded by message 63
63: Irrigation pump overheat Temperature >80 °C or disconnection of a sensor line	Deactivation of pump(s), message remains until temperature falls below the limit value

<b>Brief description (1 line with approx. 40 characters possible)</b>	<b>Reaction</b>
<b>Possible cause</b>	
64: Suction pump hot Temperature >70°C	None, message remains until the temperature falls below the limit value or is superseded by message 65
65: Suction pump overheat Temperature >80°C or disconnection of a sensor line	Deactivation of pump(s), message remains until temperature is falls below the limit value
70: Sensor test due in <4h The device has been in continuous operation for more than 20 hours since the last sensor test with offset measurement	None, message remains until device starts without cartridge
71: Sensor test overdue Since the last sensor test with offset measurement was conducted more than 24 hours ago, a new sensor test is necessary while at least one pump is activated	Message remains
72: Sensor test overdue Since the last sensor test with offset measurement was conducted more than 24 hours ago, a new sensor test is necessary. If the time is overdue, an inserted cartridge may have lead to this message.	Safe State
82: Supply voltage critical Supply voltages of process control are out of tolerance	Safe State
FE: Internal system error III An unexpected system status was detected	Safe State

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**8.3 Troubleshooting**

**WARNING:** Always unplug the device  
before all maintenance work.

*Symptom:*

- Complete failure of the device.

*Possible causes:*

- Power supply failure.
- Defective power fuse.
- Power cord connector is not properly connected to device socket.

*Remedy:*

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord connector firmly into device socket.

*Symptom:*

- Poor suction.

*Possible causes:*

- Leakage occurring in suction system.

*Remedy:*

- Check all suction lines, and replace any exhibiting evidence of leakage.
- Check the seating of the sealing cap.

*Symptom:*

- No suction.

*Possible causes:*

- Float blocking suction inlet.
- Bacterial filter on suction bottle is moist and therefore impermeable.

*Remedy:*

- Check the fluid level in the vessel and empty, if required.
- Clean float and check for freedom of movement.
- Change the bacterial filter.

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- No irrigation pressure build-up.

*Possible causes:*

- Leakage in tubing or tubing incorrectly connected.
- Failure in device's control circuitry.

*Remedy:*

- Check all tubing, in particular the connections and replace if necessary.
- Send device in for repair.

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
**Standard compliance (for 263311 20-1)**

**According to IEC 60601-1, UL 60601-1,  
CAN/CSA C22.2 No. 601.1-M90:**

- Type of protection against electric shocks:  
Protection Class I
- Degree of protection against electric shocks:  
Applied part type BF
- Type of protection against moisture:  
Drip-water protection as per IPX 1

**According to IEC 60601-1-2:**

Observe the information on electromagnetic compatibility in chapter 10 (p. 55-71).

**Directive compliance (for 263311 20-1)**

**According to Medical Device Directive (MDD):**

This medical device belongs to Class IIb

This medical device bears the CE mark in accordance with MDD 93/42/EEC.

**NOTE:** The code number after the CE mark indicates the responsible notified body.

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**8.5 Technical documentation**

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the device that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

*We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.*

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**9. Spare parts, recommended accessories**

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Item	Cat. no.
SCB connecting cord, length 100 cm	20 0901 70
Line fuses T1.6AH 250V, pack of 10	1973290
Power cord (grounded)	400 A
Power cord 'Hospital Grade' (USA)	400 B
Instruction manual	96116041 D

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Item	Cat. no.
HYS: Tubing Set, Irrigation, PC, sterile, for single use, package of 10	031523-10*
LAP: Tubing Set, Irrigation, FC, sterile, for single use, package of 10	031524-10*
Leakage tester	13242 XL
10 disposable bacterial filters	031124-10*
Suction bottle, 5 l, sterilizable	20 3000 50
Cap for suction bottle 20 3000 50, 5 l, sterilizable	20 3000 34
Bottle stand for suction bottle, 5 l	20 3000 32
Bottle stand holder	20 3000 33
Silicone suction tubing kit, sterilizable	26 3311 42
Tube connectors for suction bottle lid	20 3001 80
Plastic sphere in the nonreturn valve of the suction bottle	2911590
Plastic sleeve for nonreturn valve	20 3005 80

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**10. Electromagnetic Compatibility (EMC) Information**

The HAMOU® ENDOMAT® 26 3311 20-1 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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Accessories and cables used for EMC Compliance				
Type	Shielded	Length [m]	Ferrite	Used for
PE	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

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## Electromagnetic Compatibility (EMC) Information

**WARNING:** The use of an accessory or cable with the HAMOU® ENDOMAT® other than those specified in this manual may result in increased emissions or decreased immunity of the HAMOU® ENDOMAT®. When using an accessory or cable with HAMOU® ENDOMAT® other than those specified in this manual, it becomes the responsibility of the user of the HAMOU® ENDOMAT® to determine compliance with IEC 60601-1-2 when using this item.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the HAMOU® ENDOMAT®, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**NOTE:** The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests			
Guidance and manufacturer’s declaration – electromagnetic immunity			
The HAMOU® ENDOMAT® is intended for use in the electromagnetic environment specified below. The user of HAMOU® ENDOMAT® should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HAMOU® ENDOMAT® require continued operation during power mains interruptions, it is recommended that the HAMOU® ENDOMAT® be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	100 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the HAMOU® ENDOMAT® further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	10 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	10 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	

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**Table 3 – Test levels for Radiated and conducted Immunity Tests**  
Guidance and manufacturer’s declaration – electromagnetic immunity

The HAMOU® ENDOMAT® is intended for use in the electromagnetic environment specified below. The user of HAMOU® ENDOMAT® should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	10 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the HAMOU® ENDOMAT®, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HAMOU® ENDOMAT® is used exceeds the applicable RF compliance level above, the HAMOU® ENDOMAT® should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the HAMOU® ENDOMAT®.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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<b>Table 4 – Emission class and group</b> <b>Guidelines and manufacturer's declaration – electromagnetic emissions</b>		
<p>The HAMOU® ENDOMAT® is intended for use in the electromagnetic environment specified below. The customer or user of the HAMOU® ENDOMAT® should ensure that it is used in such an environment.</p>		
RF emissions as per CISPR 11	Group 1	The HAMOU® ENDOMAT® uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class B	The HAMOU® ENDOMAT® is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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<b>Table 5</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the HAMOU® ENDOMAT®</b>			
The HAMOU® ENDOMAT® is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the HAMOU® ENDOMAT® can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the HAMOU® ENDOMAT® as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga, Ontario  
L5N 3R3 Canada  
Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8525  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Edificio Atlantic, Oficina 3G  
Calle D e/ 1ra y 3ra  
10400 Vedado, Havana, Cuba  
Phone: +537 836 95 06, Fax: +537 836 97 76  
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.  
Lago Constanza No 326,  
Col. Granada, Del. Miguel Hidalgo  
C.P. 11520 Ciudad de México  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar - Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Stamveien1  
1483 Hagan, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
127 39 Skårholmen, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)968247755  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representation Office  
Kęstucio st. 59 / Lenktoji st. 27  
08124 Vilnius, Lithuania  
Phone: +370 5 272 0448  
Mobile: +370 685 67 000  
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraßer Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich - Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & CO. KG  
Interchange 21 Tower, Level 33  
399 Sukhumvit Road,  
North Klontoe, Wattana,  
10110 Bangkok, Thailand  
Phone: +84 28 3823 8000  
Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd., Hong Kong  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch,  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 56381188, Fax: +86 10 56381199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shanghai Branch  
36F, Tower 1 Grand Gateway  
No. 1 Hongqiao Road, Xuhui District,  
200030, Shanghai, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang Intenational Building  
No. 1 West Linjiang Road, Wuhou District,  
6100414, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd .  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)