EIAgen

HCV Ab (v.4) Kit

REF 071067


 96

REF 071064

 192

REF 071068

 480

IVD

CE 0459

This package insert must be read carefully before product use.

Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

Manufacturer:

Adaltis S.r.l

Via Durini, 27

20122 Milano (Italy)

Tel. +39-0774-5791 - Fax +39-0774-353085

www.adaltis.net

en

SYMBOLS USED ON LABELS

English EN							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
				MICROPLATE	CONTROL+	CONTROL-	CAL
	Manufacturer	Keep away from Sunlight	Date of Manufacture	Microplate	Positive Control	Negative Control	Calibrator
	CONJ	DILSPE	SUBSTMB	SOLNSTOP	WASHBUF20X	DILAS	RCNSx mL
	Conjugate	Sample Diluent	Substrate TMB	Stop Solution (H ₂ SO ₄ 0,3 M)	Wash Buffer Conc. (20x)	Assay Diluent	Reconstitute with x mL
	Biological Risk	Danger	Warning				

Attention:

Stop Solution classified as: Skin Corr. 1A

- Signal word:**
Danger
- Hazard-determining components of labelling:**
Sulphuric acid
- Hazard statements:**
H314 Causes severe skin burns and eye damage.
- Precautionary statements:**
P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Attention:

Negative Control, Positive Control, Calibrator, Conjugate, Sample Diluent, Assay Diluent and Wash Buffer Concentrate 20X classified as: Skin Sens. 1

- Signal word:**
Warning
- Hazard-determining components of labelling:**
Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)
- Hazard statements:**
H317 May cause an allergic skin reaction.
- Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P321 Specific treatment (see on this label).
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P302+P352 IF ON SKIN: Wash with plenty of water.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Refer to www.adaltis.net for the Safety Data Sheets.

A. INTENDED USE

Fourth generation Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis C Virus in human plasma (EDTA, Heparin and Citrate) and sera. The kit may be used for the screening of blood units of HCV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) define Hepatitis C infection as follows:

"Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. HCV is spread primarily by direct contact with human blood. The major causes of HCV infection worldwide are use of unscreened blood transfusions, and re-use of needles and syringes that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is too costly for most persons in developing countries to afford. Thus, from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviours (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the flaviviridae family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease. An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high propensity (80%) of inducing chronic infection. HCV is clustered into several distinct genotypes which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Cirrhosis develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 5% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer who do

not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, through the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug-users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioural practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils, or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used, the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or re-use unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (when prevalence data are available) is high compared to some countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme immunoassay (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of active infection and a potential for spread of the infection and/or the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin, can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with ribavirin is effective in about 30% to 50% of patients. Ribavirin does not appear to be effective when used alone.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any protective immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HCV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation of plasma derived products; (c) implementation and maintenance of infection control practices in health care settings, including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce overuse of injections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices. “

The genome encodes for structural components, a nucleocapsid protein and two envelope glycoproteins, and functional constituents involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

C. PRINCIPLE OF THE TEST

Microplates are coated with HCV-specific antigens derived from “core” and “ns” regions encoding for conservative and immunodominant antigenic determinants (Core peptide, recombinant NS3, NS4 and NS5 peptides).

The solid phase is first treated with the diluted sample and HCV Ab are captured, if present, by the antigens. After washing out all the other components of the sample, in the 2nd incubation bound HCV antibodies, IgG and IgM as well, are detected by the addition of polyclonal specific anti hlgG&M antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate TMB mixture, generates an optical signal that is proportional to the amount of anti HCV antibodies present in the sample. A cut-off value let optical densities be interpreted into HCV antibody negative and positive results.

D. COMPONENTS

The kit contains reagents for 96 tests (code 071067), 192 tests (code 071064), or 480 tests (code 071068).

Microplate	1
Negative Control	1x4 mL/vial
Positive Control	1x2 mL/vial
Calibrator	2 vials
Wash Buffer Concentrate 20x	1x50 mL/vial
Conjugate	1x16 mL/vial
Sample Diluent	1x50 mL/vial
Substrate TMB	1x16 mL/vial
Stop Solution	1x15 mL/vial
Assay Diluent	1x8 mL/vial
Plate sealing foils	2
Number of tests	96
Code	071067

Microplate	2
Negative Control	2x4 mL/vial
Positive Control	1x4 mL/vial
Calibrator	3 vials
Wash Buffer Concentrate 20x	2x50 mL/vials
Conjugate	2x16 mL/vials
Sample Diluent	2x50 mL/vials
Substrate TMB	2x16 mL/vials
Stop Solution	2x15 mL/vial
Assay Diluent	2x8 mL/vial
Plate sealing foils	4
Number of tests	192
Code	071064

Microplate	5
Negative Control	1x20 mL/vial
Positive Control	1x10 mL/vial
Calibrator	7 vials
Wash Buffer Concentrate 20x	5x50 mL/vials
Conjugate	2x40 mL/vial
Sample Diluent	5x50 mL/vials
Substrate TMB	2x40 mL/vials
Stop Solution	2x40 mL/vial
Assay Diluent	1x40 mL/vial
Plate sealing foils	10
Number of tests	480
Code	071068

1. Microplate

12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed into a aluminium pouch with desiccant. Bring the microplate to room temperature (18...24°C) before opening the bag. Unused strips have to be returned into the pouch and the pouch has to be sealed and stored back to 2...8°C, in presence of the desiccant.

2. Negative Control

Ready to use control. It contains 10 mM Na-citrate buffer pH 6.0 ± 0.1, 2% casein as protein base and 0.1% Proclin 150 as preservative. The negative control is olive green colour coded.

3. Positive Control

Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HCV, 10 mM Na-citrate buffer pH 6.0 ± 0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Proclin 150 as preservatives. The Positive Control is dark green colour coded.

Important Note: The absence of viable pathogens in the Positive Control can not be fully ensured, and therefore,

the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.

4. Calibrator

Lyophilized calibrator. To be dissolved with the volume of EIA grade water reported on the label. It contains foetal bovine serum proteins, human antibodies to HCV whose content is calibrated on the NIBSC Working Standard code 06/188-006, 10 mM Na-citrate buffer pH 6.0 ± 0.1 , 0.3 mg/mL gentamicin sulfate and 0.1% Proclin 150 as preservatives.

Important Note: *The absence of viable pathogens in the Calibrator can not be fully ensured, and therefore, the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.*

Note: *The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.*

5. Wash Buffer Concentrate 20x

20x concentrated solution. Once diluted, the wash solution (wash buffer diluted) contains 10 mM phosphate buffer pH 7.0 ± 0.2 , 0.05% Tween 20 and 0.05% Proclin 150.

Once diluted, the wash solution is stable for 1 week at $2...8^{\circ}\text{C}$.

6. Conjugate

Ready to use and red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Citrate buffer pH 6.4 ± 0.1 , 0.1% Proclin and 0.05% Tween 20 as preservatives.

7. Substrate TMB

Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H_2O_2 . Mix gently before use.

Note: *To be stored protected from light as sensitive to strong illumination.*

8. Assay diluent

Ready-to-use component. It contains goat serum, 10 mM tris buffered solution pH 8.0 ± 0.1 containing 0.1% Proclin 150 and 0.09% Na azide for the pre-treatment of samples and controls in the plate, blocking interference.

9. Stop Solution

Ready-to-use component.

It contains 0.3 M H_2SO_4 solution. Mix gently before use.

10. Sample Diluent

Ready-to-use component and dark green colour coded. It contains 1% casein, 10 mM Na-citrate buffer pH 6.0 ± 0.1 and 0.1% Proclin 150 as preservative.

To be used to dilute the sample.

Note: *The diluent changes colour from olive green to dark bluish green in the presence of sample.*

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (200 μL and 10 μL) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of $+37^{\circ}\text{C}$.
6. Calibrated ELISA microwell reader with 450nm (reading) and possibly with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
This package insert must be read carefully before product use.
2. Read carefully the Safety Data Sheet (SDS) before product use.
3. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
4. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
5. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
6. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
7. Upon receipt, store the kit at $2...8^{\circ}\text{C}$ into a temperature controlled refrigerator or cold room.
8. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
9. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
10. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
11. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

12. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
13. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
14. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
15. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
16. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
17. The Stop Solution contains 0.3M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
18. Dispose of reagent solutions containing sodium azide as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.
19. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
20. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
21. Do not pipette by mouth.

G. SPECIMEN: PREPARATION AND RECCOMANDATIONS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.

4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 1 re-use of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the pouch. Check that the desiccant has not turned to dark green, indicating a defect in manufacturing. In this case, call Adaltis's customer service.

Unused strips have to be placed back inside the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2...8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

4. Calibrator:

Dissolve carefully the content of the lyophilised vial with the volume of EIA grade water reported on its label. Mix well on vortex before use.

Handle this component as potentially infectious.

Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.

5. Wash Buffer Concentrate 20x (vial of 50 mL):

The whole content of the 20x concentrated solution has to be diluted with EIA grade water up to 1000 mL (the volume is reported on the label) and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at +2...8°C.

6. Conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component

has to be transferred use only plastic, possibly sterile disposable containers.

7. Substrate TMB:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

8. Assay Diluent:

Ready to use. Mix well on vortex before use.

9. Stop Solution:

Ready to use. Mix well on vortex before use.

10. Sample Diluent:

Ready to use. Mix well on vortex before use.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized. Usually 4-5 washing cycles (aspiration + dispensation of $350\text{ }\mu\text{L}$ /well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well-characterized negative and positive reference samples, and check to match the values reported below in sections O "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of $\pm 5\%$.
 - ✓ Short Incubation Method (for 1st / 2nd incubation tolerance between 43 min to 47 min; for 3rd incubation tolerance between 14 and 16 min).
 - ✓ Standard Incubation Method (for 1st incubation tolerance between 57 min to 63 min; for 2nd and 3rd incubation tolerance between 29 and 31 min).
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second

filter (620-630nm) for blanking purposes. Its standard performances should be (a) bandwidth $\leq 10\text{ nm}$; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in sections O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Substrate TMB is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Buffer as described above.
4. Dissolve the Calibrator as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in section I.3.
7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

The Assay Procedure can be done with two times incubation procedures. Choose the one that is required by your regulation:

1. Standard Incubation (1st incubation 60 minutes, 2nd and 3rd incubation 30 minutes)
2. Short Incubation (1st and 2nd incubation 45 minutes, 3rd incubation 15 minutes)

1. Standard Incubation - Manual assay:

1. Place the required number of microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **60 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted wash solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **30 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the microplate at room temperature (18...24°C) for **30 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

11. Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will

turn the positive control and positive samples from blue to yellow.

12. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
3. Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

2. Short Incubation - Manual assay:

1. Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **45 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted washing solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **45 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the

microplate at room temperature (18...24°C) for **15 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
- Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

N. ASSAY SCHEME

Method	Operations (Standard Incubation)	Operations (Short Incubation)
Controls & Calibrator	200 µL	200 µL
Sample Diluent & Sample	200 µL diluent.+ 10 µL sample	200 µL diluent.+ 10 µL sample
Assay Diluent	50 µL	50 µL
1st incubation	60 min (± 3)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Enzyme Conjugate	100 µL	100 µL
2nd incubation	30 min (± 1)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Substrate TMB	100 µL	100 µL
3rd incubation	30 min (± 1)	15 min (± 1)
Temperature	Room temperature (18...24°C)	Room temperature (18...24°C)
Stop Solution	100 µL	100 µL
Reading OD	450/620nm	450/620nm

An example of dispensation scheme is reported below (valid for both incubation time procedures):

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	S1	S9										

Legenda: BLK = Blank NC = Negative Control
CAL = Calibrator PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450/620nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD 450/620nm value
Negative Control (NC)	< 0.050 mean OD450/620nm after blanking
Calibrator	S/Co >1.1
Positive Control	>1.000 OD450/620nm value

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problems	Check
Blank well > 0.100 OD450nm	1. that the Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co<1.1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control <1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150,too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. RESULTS

The tests results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off (Co)} = \text{NC mean} + 0.350$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused. Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.

A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgement errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies **earlier** than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result ! The sample has to be anyway submitted to a confirmation test.
4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV antibodies - lacking anti hIgM conjugate in the formulation of the enzyme tracer and therefore missing IgM reactivity - may be present. The real positivity of the sample for antibodies to HCV should be then confirmed by examining also IgM reactivity, important for the diagnosis of HCV infection.
5. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
6. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user.

Negative Control: 0.019 – 0.020 – 0.021 OD450nm

Mean Value: 0.020 OD450nm Lower than 0.050 – Accepted

Positive Control: 2.189 OD450nm

Higher than 1.000 – Accepted

Cut-Off = 0.020+0.350 = 0.370

Calibrator: 0.550 - 0.530 OD450nm

Mean value: 0.540 OD450nm S/Co = 1.4

S/Co higher than 1.1 – Accepted

Sample 1: 0.070 OD450nm

Sample 2: 1.690 OD450nm

Sample 1 S/Co < 0.9 = negative

Sample 2 S/Co > 1.1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC) and with both incubation procedures (standard and short).

1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated using the short incubation assay by means of the British Working Standard for anti-HCV, NIBSC code 06/188-006. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	3,50	4,00	4,30
2 X	2,10	2,60	2,60
4 X	1,3	1,40	1,30
Negative Plasma	0,25	0,20	0,20

In addition the sample coded Accurun 1 - series 3000 - supplied by Boston Biomedica Inc., USA, has been evaluated "in toto" showing the results below:

Accurun 1 series	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	2,90	3,04	3,40

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

2.1 Diagnostic specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. A total of more 5000 unselected donors, including 1st time donors, were examined.

The diagnostic specificity was assessed against a kit US FDA approved.

5043 blood donors were tested providing a specificity of 99.5%. 210 hospitalized patients were tested for HCV Ab; a diagnostic specificity of 99.5% was found.

Moreover, diagnostic specificity was assessed by testing 162 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.). A value of specificity of 100% was assessed.

No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage. No interference was observed.

2.2 Diagnostic Sensitivity

It defined as the probability of the assay of scoring positive in the presence of specific analyte. The diagnostic sensitivity has been assessed externally on a total number of 348 specimens; a diagnostic sensitivity of 100% was found. Internally more than other 50 positive samples were tested, providing a value of diagnostic sensitivity of again 100%. Positive samples from infections carried out by different genotypes of HCV were tested as well.

Furthermore, most of seroconversion panels available from Boston Biomedica Inc., USA, (PHV) and Zeptometrix, USA, (HCV) have been studied.

Results are reported below for some of them.

Panel	N° samples	Adaltis ¹	Ortho ^{1,2}
PHV 901	11	9	9
PHV 904	7	2	4
PHV 905	9	3	4
PHV 906	7	7	7
PHV 907	7	3	2
PHV 908	13	10	8
PHV 909	3	2	2
PHV 910	5	3	3
PHV 911	5	3	3
PHV 912	3	1	1
PHV 913	4	2	2
PHV 914	9	5	5
PHV 915	4	3	0
PHV 916	8	4	3
PHV 917	10	6	6
PHV 918	8	2	0
PHV 919	7	3	3
PHV 920	10	6	6
HCV 10039	5	2	0
HCV 6212	9	6	7
HCV 10165	9	5	4

Note:

1. Positive samples detected
2. HCV v.3.0

Finally the Product has been tested on the panel EFS Ac HCV, lot n° 06.140817, supplied by the Etablissement Francais Du Sang (EFS), France, with the following results:

EFS Panel Ac HCV

sample	Lot#1 S/Co	Lot#2 S/Co	Lot#3 S/Co	Average value
HCV 1	0,53	0,52	0,55	Negative
HCV 2	3,28	5,91	3,04	Positive
HCV 3	2,17	3,18	2,56	Positive
HCV 4	2,26	2,23	2,35	Positive
HCV 5	6,10	7,06	6,90	Positive
HCV 6	1,66	1,77	1,67	Positive

3. PRECISION

It has been calculated on five samples, one negative and four positives examined in 4 replicates each in six separate runs. The precision, obtained for both protocols, is equivalent.

Results are reported as follows:

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (short incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.03	6.66	10.56	12.48
Positives	1.20	8.52	8.49	12.03
	1.51	7.69	12.22	14.44
	3.57	7.43	11.82	13.97
	11.87	3.42	9.32	9.92

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (long incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.04	4.67	12.34	13.19
Positives	1.47	9.62	11.40	14.92
	1.82	8.92	12.77	15.58
	4.31	4.59	12.88	13.67
	13.78	2.42	8.96	9.26

Inter-lot results: ElAgen HCV Ab (v.4) Kit - 1st, 2nd and 3rd Lot (short incubation assay)

		Precision - %CV		
Sample	Lot 1	Lot 2	Lot 3	
Negative	8,65	8,29	6,13	
Calibrator	4,98	4,44	5,38	
Positive	4,11	3,11	1,37	

The variability shown in the tables above did not result in sample misclassification.

S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ($\pm 50\%$) from OD reported on QC	- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) - error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	- not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) - error in washing or in reading (suggestion: check operating or settings of respective instruments) - contamination of Substrate (suggestion: use only disposable and clean plastic containers) - pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colorimetric reaction after addition of substrate	- some reagent not pipetted - strong contamination of Conjugate or Substrate - errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	- incubation time too short, incubation temperature too low - incorrect conjugate dilution
too high reaction (too high ODs)	- incorrect conjugate dilution - incubation time too long, incubation temperature too high - water quality for wash buffer insufficient (low grade of deionization) - insufficient washing (conjugates not properly removed)
unexplainable outliers	- contamination of pipettes, tips or containers - inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	- reagents and/or strips not pre-warmed to Room Temperature prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	- incubation conditions not constant (time, temperature) - controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) - person-related variation

T. AUTOMATION

The procedures identified in this Instruction for Use are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

U. LIMITATIONS

Repeatable false positive results, not confirmed by RIBA or similar confirmation techniques, were assessed as less than 0.1% of the normal population.

Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

BIBLIOGRAPHY

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon-α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
BIBLIOGRAFIA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhuth MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretsch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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

REFERENCES

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998 ;44 :98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993 ;167 :66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993 ;269 :392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994 ;169 :90-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995 ;171 :768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993 ;153 :1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhuth MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretsch, D. Diagnostic tests for hepatitis C. A cikk megtalálható: <http://www.hepnet.com/nih/gretsch.html>. A Nemzeti Egészségügyi Hivatal által szervezett, C hepatitisz konferenciára készült, amit 1997 március 24-26.-án tartottak Bethesda Maryland-ban
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
BIBLIOGRAFIA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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LITERATURA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, wyd. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turyn, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, wyd. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Waszyngton, DC: US Government Printing Office, 1994; Nr publikacji NIH 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN i in. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH i in. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ i in. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW i in. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB i in. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ i in. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A i in. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK i in. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK i in. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL i in. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW i in. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ i in. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland i A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre i P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. Artykuł dostępny pod adresem: <http://www.hepnet.com/nih/gretch.html>. Napisany w ramach konferencji Narodowego Instytutu Zdrowia dot wirusowego zapalenia wątroby typu C w dniach 24-26 marca 1997 roku w Bethesda w stanie Maryland (USA)
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli i E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams i J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada i T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata i T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIOGRAFIE


1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretsch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JL, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon-α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIOGRAFIJA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhoo MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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BIBLIYOGRAPHI

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. MMWR 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. Hepatology 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. Viral Hepatitis and Liver Disease, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. Digestive diseases in the United States: epidemiology and impact. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. JAMA 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. Ann Intern Med 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. Lancet 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. Am J Med Sci 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. N Engl J Med 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. Science 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. Science 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. N Engl J Med 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. N Engl J Med 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. N Engl J Med 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. Semin Liver Dis 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. Am J Epidemiol 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. West J Med 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. Arch Intern Med 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. Am J Pub Health 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. Blood 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. Blood 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. J Med Virology 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. ASAIO Journal 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. J Infect Dis 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. JAMA 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. J Infect Dis 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. J Infect 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. J Infect Dis 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. Arch Intern Med 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. Infect Control Hosp Epidemiol 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. J. Clin. Microbiol. 40:1656-1659.
32. Dusheiko, G., H. Schilmovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. Hepatology 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. J. Clin. Microbiol. 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. J. Clin. Microbiol. 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. J. Hepatol. 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. Hepatology 16:293-299.

***H. pylori* Antigen ELISA Test Kit**

An enzyme immunoassay (ELISA) for the qualitative and quantitative detection of *Helicobacter pylori* (*H. pylori*)

Antigen in human stool

IVD For *in vitro* diagnostic and professional use only

2°C  8°C Store at (2° to 8°C)

 96 Tests

INTENDED USE

The *H. pylori* Antigen ELISA Test Kit is an enzyme immunoassay for the qualitative and quantitative detection of *H. pylori* antigen in human stool. It is intended as an aid in the diagnosis of possible *H. pylori* infection and in the follow-up of patients undergoing antimicrobial therapy.

INTRODUCTION

Helicobacter pylori are Gram-negative spiral-shaped bacteria that have adapted to living in the harsh acidic conditions of the stomach. These bacteria can alter their surrounding micro-environment by reducing its acidity so they can survive. Their spherical shape facilitates penetration of the epithelial lining, where the bacteria are protected by mucus against cells of the immune system.

Infections with *H. pylori*, though harmless during childhood, manifest as peptic ulcers of the stomach, duodenum and of small intestine, active and chronic gastritis, as well as non-ulcer dyspepsia, in about 60% of the global adult population. The mechanism of bacterial transmission is still unknown, but is thought to be oral and/or fecal borne.

The *H. pylori* Antigen ELISA Test Kit is an immunoassay for the qualitative and quantitative detection of *H. pylori* Antigen in human stool. The test utilizes antibodies to *H. pylori* to selectively detect *H. pylori* Antigen in human stool.

PRINCIPLE OF THE TEST

The *H. pylori* Antigen ELISA Test Kit is a solid phase enzyme immunoassay based on sandwich principle for the qualitative and quantitative detection of *H. pylori* antigen in human stool. The microwell plate is coated with anti-*H. pylori* antibodies. During testing, the antigens are extracted from the specimen with extraction solution and added onto the antibodies coated microwell plate along with the enzyme- conjugated antibodies to *H. pylori*, and then incubated. If specimens contain *H. pylori* antigens, it will bind to the antibodies coated on the microwell plate and simultaneously bind to the conjugate to form immobilized antibody-*H. pylori* antigen-conjugate complexes. If specimens do not contain *H. pylori* antigens, the complexes will not be formed. After initial incubation, the microwell plate is washed to remove unbound materials. Substrate A and substrate B are added and then incubated to produce a blue color indicating the amount of *H. pylori* antigens present in the specimens. Sulfuric acid solution is added to the microwell plate to stop the reaction producing a color change from blue to yellow. The color intensity, which corresponds to the amount of *H. pylori* antigens present in the specimens, is measured with a microplate reader at 450/630-700 nm or 450 nm.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Follow the instructions for use carefully. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.
- Wear protective clothing and disposable gloves when dealing with samples and reagents. Wash hands after operations.
- Do not use reagents beyond the labeled expiry date.
- Do not mix or use components from kits with different batch codes.
- It is important to calibrate all the equipment e.g. micropipettes, and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Ensure that the bottom of the plate is clean and dry and that no bubbles are present on the surface of the liquid before reading the plate.
- Never eat, drink, smoke, or apply cosmetics in the assay laboratory. Never pipette solutions by mouth.
- Avoid cross contamination between reagents to ensure valid test results.
- Follow the wash procedure to ensure optimum assay performance.
- Use Plate Sealer to cover microwell plate during incubation to minimize evaporation.
- Use a new pipette tip for each specimen assayed.
- Do not touch or splash the rim of the well. Do not blow out from micropipettes.
- Do not allow sodium hypochlorite fumes from chlorine bleach or other sources to contact the microwell plate during the assay as the color reaction may be inhibited.

HEALTH AND SAFETY INFORMATION

- Collect samples in accordance with correct medical practices.
- Some reagents may cause toxicity, irritation, burns or have a carcinogenic effect as raw materials. Contact with the skin and the mucosa should be avoided but not limited to the following reagents: Stop solution, the Conjugate, and the Wash buffer, Extraction solution, Substrate.
- The Stop solution 0.5M H₂SO₄ is an acid. Use it with appropriate care. Wipe up spills immediately and wash with water if it comes into contact with the skin or eyes.
- ProClin™ 300 0.1% is used as a preservative; it can cause irritation of the skin. Wipe up spills immediately or wash with water if it comes into contact with the skin or eyes.
- Pipette tips, vials, strips and specimen containers should be collected and autoclaved for not less than 2 hours at 121°C before any further steps of disposal. Solutions containing sodium hypochlorite should NEVER be autoclaved.
- All specimens and materials should be considered as potentially infectious. Strict adherence to GLP (Good Laboratory Practice) regulations can ensure personal safety.
- Chemicals should be handled and disposed of only in accordance with the current GLP (Good Laboratory Practices) and the local or national regulations.
- Neutralized acids and other liquids should be decontaminated by adding sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1.0% sodium hypochlorite may be necessary to ensure effective decontamination.

STORAGE AND STABILITY

- Components of the kit will remain stable through the expiration date indicated on the label and package when stored between 2-8°C, Once opened; all reagents are stable for up to 3 months after the first opening date if stored between 2-8°C. Return reagents to 2-8°C immediately after use.
- Place unused wells in the zip-lock aluminum foiled pouch and return to 2-8 °C, under which conditions the wells will remain stable for 3 months from the opening date.
- Concentrated Wash Buffer may be stored at room temperature to avoid crystallization. If crystals are present, warm up the solution at 37°C. Working Wash Buffer is stable for 2 weeks at room temperature.
- Do not expose reagents especially the Substrate to strong light or hypochlorite fumes during storage or incubation steps.
- Do not store Stop Solution in a shallow dish or return it to the original bottle after use.

SPECIMEN COLLECTION AND PREPARATION

- This *H. pylori* Antigen ELISA Test can be performed using only human stool.
- Stool samples should be collected in clean containers. Samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample should be totally thawed and brought to room temperature before testing.
- The patient has to be asked to collect the specimen avoiding any possible contact with urine or water.
- The patient submitted to the test should not be under antibiotic or anti-bacterial treatments as this pharmaceutical therapy is known to affect *H. pylori* up to a certain extent, depending on the antibiotic used, giving rise to false interpretation.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

MATERIALS

MATERIALS PROVIDED

1. ***H. Pylori* Antigen Microwell Plate:** Microwell plate coated with anti-*H. Pylori* antibodies. **(1 plate: 96 wells/plate).**
2. ***H. pylori* Antigen Conjugate:** One red cap vial containing antibodies to *H. pylori* bound to peroxidase; Preservative: 0.1% ProClin™ 300. **(1 x 8 mL).**
3. **Concentrated Wash Buffer (25x):** One white cap bottle containing Tris-HCl buffer containing 0.1% Tween 20; Preservative: 0.1% ProClin™ 300. **(1 x 40 mL).**
4. **Extraction Solution:** One white cap bottle containing 0.9% NaCl buffer containing EDTA; Preservative: 0.1% ProClin™ 300. **(1 x100 mL).**
5. **Substrate A:** One white cap vial containing Citrate-phosphate buffer containing hydrogen peroxide; Preservative: 0.1% ProClin™ 300. **(1 x 8 mL).**
6. **Substrate B:** One blue cap vial containing Buffer containing tetramethylbenzidine (TMB); Preservative: 0.1% ProClin™ 300. **(1 x 8 mL).**
7. **Stop Solution:** One yellow cap vial containing 0.5M Sulfuric acid. **(1 x 8 mL).**

8. **H.pylori Antigen Calibrator 1:** One white cap vial containing Buffer non-reactive for *H. pylori* Antigen; Preservative: 0.1% ProClin™ 300. (1 x 1 mL).
9. **H.pylori Antigen Calibrator 2:** One white cap vial containing Buffer containing 0.1 µg/mL *H. pylori* Antigen; Preservative: 0.1% ProClin™ 300. (1 x 1 mL).
10. **H.pylori Antigen Calibrator 3:** One white cap vial containing Buffer containing 0.5 µg/mL *H. pylori* Antigen; Preservative: 0.1% ProClin™ 300. (1 x 1 mL).
11. **H.pylori Antigen Calibrator 4:** One white cap vial containing Buffer containing 1.0 µg/mL *H. pylori* Antigen; Preservative: 0.1% ProClin™ 300. (1 x 1 mL).
12. **Plate Sealers (2 pieces).**
13. **Package Insert (1 copy).**

MATERIALS REQUIRED BUT NOT PROVIDED

- Freshly distilled or deionized water.
- Sodium hypochlorite solution for decontamination.
- Absorbent paper or paper towel.
- Water bath or incubator capable of maintaining 15°C to 30°C.
- Calibrated automatic or manual microwell plate washer capable of aspirating and dispensing 350 µL/well.
- Disposable gloves.
- Automated processor (optional).
- Calibrated micropipettes with disposable tips capable of dispensing 50 and 100 µL.
- Graduated cylinders for wash buffer dilution.
- Vortex mixer for specimen mixing (optional).
- Disposable reagent reservoirs.
- Calibrated microplate reader capable of reading at 450 nm with a 630-700 nm reference filter, or reading at 450 nm without a reference filter.
- Timer.

DIRECTIONS FOR USE

1. Remove unused strips from the microwell plate, and store in the original resealable pouch at 2-8°C.
2. Allow reagents and specimens to reach room temperature (15-30°C) prior to testing.

WASH PROCEDURE

- a. The wash procedure is critical. Insufficient washing will Result in a poor precision and falsely elevated absorbance readings.
- b. Prepare working wash buffer by adding content of wash buffer bottle provided with the kit to distilled or deionized water to reach a final volume of 1 liter. The working wash buffer is stable for 2 weeks at 15-30°C.
3. Dispense 1 mL of Extraction Solution into Specimen Extraction Tube.

For Solid Stool Specimens:

- i. Take out the cap of the **Specimen Extraction Tube**
- ii. Randomly stab the specimen collection stick into the stool specimen in at least 3 different sites to collect approximately 30 mg of specimen (equivalent to 1/4 of a pea). Do not scoop the stool specimen.
- iii. Transfer into Specimen Extraction Tube.

For Liquid Stool Specimens:

- i. Hold the Liquid Specimen Dropper vertically.

- ii. Aspirate stool specimens and then dispense 2 drops (approximately 50 µL) into the Specimen Extraction Tube containing the Extraction Solution.
 - iii. Screw on and tighten the cap onto the Specimen Extraction Tube.
 - iv. Shake the Specimen Extraction Tube vigorously to mix the specimen and the Extraction Solution.
4. Leave A1 as Blank well.
 5. Dispense 50 µL of Calibrator 1 in wells B1 and C1. (Light Yellow Reagent)
 6. Dispense 50 µL of Calibrator 2 in wells D1 and E1. (Green Blue Reagent)
 7. Dispense 50 µL of Calibrator 3 in wells F1 and G1. (Light Blue Reagent)
 8. Dispense 50 µL of Calibrator 4 in wells H1 and A2. (Dark Blue Reagent)
 9. Hold the Specimen Extraction Tube upright and break off the tip of the tube. Invert the Specimen Extraction Tube and dispense 2 drops of the specimen Extraction Solution (approx. 50 µL) to assigned wells starting at B2. (Yellow Reagent)
 10. Dispense 50 µL of Conjugate to each well except for the Blank well. (Red Reagent)
 11. Mix gently by swirling the microwell plate on a flat bench for 30 seconds.
 12. Cover the microwell plate with the Plate Sealer and incubate at room temperature (15-30°C) in a room, a water bath, or an incubator for 60 minutes ± 5 minutes.
 13. Remove the Plate Sealer.
 14. Wash each well 5 times with 350 µL of Working Wash Buffer per well, and then remove the liquid.
 15. Turn the microwell plate upside down on absorbent tissue for a few seconds. Ensure that all wells have been completely washed and dried.
Note: Improper washing may cause false positive results.
 16. Dispense 50 µL of Substrate A to each well. (Clear Reagent)
 17. Dispense 50 µL of Substrate B to each well. (Clear Reagent) Then a blue color should develop in wells containing Positive specimens.
 18. Mix gently then cover microwell plate with Plate Sealer and incubate at room temperature (15-30°C) in a room, a water bath, or an incubator for 10 minutes ± 1 minute.
 19. Remove the Plate Sealer.
 20. Dispense 50 µL of Stop Solution to each well. (Clear Reagent) Then a yellow color should develop in wells containing Positive specimens.
 21. Read at 450/630-700 nm within 30 minutes.
 22. Note: Microwell plate can also be read at 450 nm, but it is strongly recommended to read it at 450/630-700 nm for better results.

VALIDATION REQUIREMENT AND QUALITY CONTROL

1. Calculate the Mean Absorbance of Calibrators 1-4 by referring to the table below.

Example of Calibrator 2 Calculation	
Item	Absorbent
Calibrator 2: Well D1	0.469
Calibrator 2: Well E1	0.507
Total Absorbance of Calibrator 2	0.976
Mean Absorbance of Calibrator 2	0.488

2. Check the validation requirements below to determine if the test results are valid.

Item	Validation Requirements
Blank Well	Blank Absorbance should be < 0.050 if read at 450/630-700 nm Note: It should be < 0.100 if read at 450 nm
Calibrator 1	Mean Absorbance after subtraction of Blank Absorbance should be < 0.100
Calibrator 2	Mean Absorbance after subtraction of Blank Absorbance should be > 0.150
Calibrator 3	Mean Absorbance after subtraction of Blank Absorbance should be > 0.500
Calibrator 4	Mean Absorbance after subtraction of Blank Absorbance should be > 1.000

NOTE: The test results are considered invalid if the above validation requirements are not met. Repeat the test or contact your local distributor.

INTERPRETATION OF RESULTS

Qualitative

Calculate the Index Value to obtain qualitative specimen results.

1. If the test is valid, obtain Cut-Off Value by subtracting the Blank Absorbance from the Mean Absorbance of 1/2× (Calibrator 2+Calibrator 1). See an example of Cut-Off Value calculation below.


Item	Absorbent
Blank Absorbance: Well A1	0.011
Cut-Off Value: 1/2× (Mean Absorbance of Calibrator 2+ Mean Absorbance of Calibrator 1) – Blank Absorbance	1/2× (0.488+0.012)-0.011=0.239

2. Calculate the Index Value by dividing the Specimen Absorbance by the Cut-Off Value, and then read the results by referring to the Interpretation of Results table below.

Item	Absorbent
Specimen: Well F2	0.968
Blank Absorbance: Well A1	0.011
Cut-Off Value	0.239
Index Value: Specimen/Cut-Off Value	(0.968-0.011)/0.239=4.0

Quantitative

Draw the calibration curve and obtain quantitative specimen results.

1. Subtract the Blank Absorbance from the Mean Absorbance of each Calibrator, and then plot them on the Y-axis against their Log₁₀ of the corresponding concentration in µg/mL on the X-axis on a linear graph paper and draw the calibration curve. Draw the best fitted line through data points to obtain a standard curve. Refer to an example of the calibration curve at right.

NOTE: Do not use the calibration curve at right to make any calculation. A calibration curve must be performed for each run.

2. Obtain quantitative specimen results from their absorbance by using the calibration curve.

NOTE: Specimens that have absorbance above Calibrator 4 should be pre-diluted using Extraction Solution and retested. The concentration must be multiplied by the dilution factor. Automated reading and calculation may also be performed using linear regression function on suitable computer programs.

Interpretation of Results - Qualitative and Quantitative

Results	Qualitative	Quantitative
	Index Value	Concentration
Negative	< 0.9	< 0.045 µg/mL
Positive	> 1.1	> 0.055 µg /mL
Equivocal*	≥ 0.9 and ≤ 1.1	0.045 – 0.055 µg/mL

***NOTE:** For Equivocal results, the specimen should be retested. Specimens that are repeatedly Equivocal after retest should be confirmed using an alternate method. If the results remain Equivocal, collect a new specimen in two weeks. If the new specimen is Positive, the specimen is presumed to be Positive.

LIMITATIONS

1. The *H. pylori* Antigen ELISA Test Kit is used for the detection of *H. pylori* antigen in human stool. Diagnosis of an infectious disease should not be established based on a single test result. Further testing, including confirmatory testing, should be performed before a specimen is considered positive. A negative test result does not exclude the possibility of exposure. Specimens containing precipitate may give inconsistent test results.
2. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
3. As with other sensitive immunoassays, a false positive result may arise due to inadequate washing from the initial test. The results may be affected due to procedural or instrument error.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The *H. pylori* Antigen ELISA Test Kit has been compared to a leading commercial *H. Pylori* Antigen ELISA test using clinical specimens. The results show that the clinical sensitivity of the *H. pylori* Antigen ELISA Test Kit is 98.6%, and the clinical specificity is 95.4%.

H. pylori Antigen ELISA vs. Other ELISA

Method		Other ELISA		Total Results
<i>H. pylori</i> Antigen ELISA	Results	Positive	Negative	
	Positive	70	6	76
	Negative	1	125	126
Total Results		71	131	202

Clinical Sensitivity: 98.6% (92.4-100.0%) *

Clinical Specificity: 95.4% (90.3-98.3%)*

Overall Agreement: 96.5% (93.0-98.6%)*

*95% Confidence Interval

Reproducibility

Intra-Assay: Within-run precision has been determined by using 10 replicates of two specimens: a low positive and a high positive.

Inter-Assay: Between-run precision has been determined by using 10 replicates on the same two specimens: a low positive and a high positive. Three different lots of the *H. pylori* Antigen ELISA Test Kit have been tested using these specimens.

Specimen	Intra-Assay			Inter-Assay		
	Mean Absorbance/Cut-Off	Standard Deviation	Coefficient of Variation (%)	Mean Absorbance/Cut-Off	Standard Deviation	Coefficient of Variation (%)
1	1.741	0.156	8.96	1.723	0.133	7.72
2	4.726	0.252	5.33	4.861	0.252	5.18

REFERENCES

1. Marshall, BJ, McGechie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Ansorg, R, Von Recklinghausen, G, Pomarius, R and Schmid, EN. Evaluation of techniques for isolation, subcultivation and preservation of *Helicobacter pylori*. J. Clin. Micro. (1991), 29:51-53.
4. Pronovost, AP, Rose, SL, Pawlak, J, Robin, H and Schneider, R. Evaluation of a new immunodiagnostic assay for *Helicobacter pylori* antibody detection: Correlation with histopathological and microbiological results. J. Clin. Micro. (1994), 32: 46-50.
5. Megraud, F, Bassens-Rabbe, MP, Denis, F, Belbouri, A and Hoa, DQ. Seroepidemiology of *Campylobacter pylori* infection in various populations. J. Clin. Micro. (1989), 27: 1870-3.

ATLAS MEDICAL

Ludwig-Erhard Ring 3

15827 Blankenfelde-Mahlow

Germany

Tel: +49 - 33708 – 3550 30

Email: Info@atlas-medical.com

Website: www.atlas-medical.com

PPI1532A01

Rev B (21.10.2021)


	Catalogue Number		Temperature limit
	In Vitro diagnostic medical device		Caution
	Contains sufficient for <n> tests and Relative size		Consult instructions for use (IFU)
	Batch code		Manufacturer
	Fragile, handle with care		Use-by date
	Manufacturer fax number		Do not use if package is damaged
	Manufacturer telephone number		Date of Manufacture
	Keep away from sunlight		Keep dry

anti-Trichinella spiralis

ELISA kit for the qualitative detection of antibodies to
Trichinella spiralis

Instructions for use

IVD

REF
EI-605

Σ 96
tests

CE

EQUI anti-*Trichinella spiralis*

ELISA kit for the qualitative detection of antibodies to *Trichinella spiralis*

1. INTENDED USE

The «EQUI anti-*Trichinella spiralis*» is ELISA kit intended to qualitatively detect anti-*Trichinella spiralis* antibodies in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose trichinosis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: villagers, hunters and their families, people whose diet includes meat, particularly pork.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Among helminthiasis, which are often caused by meat consumption, trichinosis has a special place. The causative agent of this disease is the *Trichinella spiralis* roundworm, the larvae of which were found even in the mummies of Ancient Egypt. The greatest danger of this helminthiasis is the possibility of severe damage to the central nervous system, which leads to death.

Adult *Trichinella spiralis* have a filamentous body up to 0.5 mm in diameter and 2 mm in length. *Trichinella* is characterized by live birth – while in the human body, one female can give birth to about one and a half thousand rod-shaped larvae. The larvae actively migrate to other organs, but only settle massively in the striated muscles. Connective tissue capsule with feeding vessels gets formed around the larva after 1-1.5 months. About a year later, the capsule begins to calcify, but the larva inside remains alive and invasive for many years. Most mammals can become infected with trichinosis in the natural environment by eating fresh or decomposed meat that contains *Trichinella* larvae.

Trichinella enters the human body through the meat infested with encapsulated larvae (pork, horse meat, nutria's meat, wild boar, bear, etc.). Acute (intestinal) trichinosis may be asymptomatic or show symptoms 1-2 days after eating raw or undercooked meat. This stage of the disease is characterized by abdominal pain, nausea or diarrhea. However, the classic symptoms of trichinosis appear at the stage of larval migration about 2 weeks after the invasion and can last up to two months. These include fever, fatigue, muscle aches and pain, swelling of the face and eyes, skin rashes, and more. *Trichinella spiralis* affects the most active human muscles: diaphragm, intercostal and abdominal muscles, tongue, as well as facial and skeletal muscles. The overall extent of clinical manifestations is directly correlated with the quantity of consumed parasites.

During larval migration, the organism is hyper-sensitized by highly immunogenic products of *Trichinella* metabolism. Patients have eosinophilia and an increase in IgE antibodies. However, a strong immune response does not lead to the elimination of parasites, although it may be crucial for protection from re-infection.

Thus, recurring cases of trichinosis are much easier to tolerate than the primary ones.

To diagnose this helminthiasis, specific tests are conducted on meat that could be a source of invasion, to test for the presence of encapsulated larvae of *Trichinella spiralis*. During the early stages of infection, PCR is used to detect *Trichinella* DNA. A few weeks after the probable invasion, a biopsy of the affected patient's muscle tissue may be performed to detect the parasite. But from the second week of the invasion, serodiagnosis with the detection of specific antibodies to antigens of larvae of the pathogen trichinosis and an increase in the titer of these antibodies in paired sera is particularly useful. 4-5 months after the invasion, the number of specific antibodies decreases, but may remain at a level sufficient for detection for several years. Immunoblotting is used to confirm the positive result of the detection of specific IgG antibodies. However, for a comprehensive diagnosis of trichinosis we should use not only laboratory-generated data but also clinical and epidemiological studies.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Trichinella spiralis* antibodies in «EQUI anti-*Trichinella spiralis*» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Trichinella spiralis* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Trichinella spiralis* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species (anti-IgG and anti-IgA) monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	Microplate
		Each plate well is coated with <i>Trichinella spiralis</i> larvae antigens. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months
CONTROL +	1 x 0,25 ml	Positive control
		Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C

		Negative control
CONTROL -	1 x 0,6 ml	Negative human serum with a preservative (yellow). Store at 2-8 °C
		Serum dilution solution
DIL SAMPLE	1 x 13 ml	Buffer solution with a milk extract, a detergent and a preservative (purple). Store at 2-8 °C
		Conjugate solution (ready to use)
SOLN CONJ	1 x 13 ml	Buffer solution of monoclonal antibodies to human IgG and IgA, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
		TMB solution (ready to use)
SOLN TMB	1 x 13 ml	TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
		Washing solution TWEEN (20x concentrated)
TWEEN WASH 20x	1 x 50 ml	20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
		Stop Solution (ready to use)
SOLN STOP	1 x 13 ml	0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Besure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI anti-Trichinella spiralis» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;

- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI anti-Trichinella spiralis» kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

- 9.4. Add 90 µL of [DIL|SAMPLE] into each plate well.
- 9.5. Add 10 µL of controls and test samples into the wells:
 [CONTROL|+] – into well A1,
 [CONTROL|-] – into wells B1, C1 and D1,
 and test samples into the remaining wells.
 At the time of adding, the solution changes its colour from purple to blue.
 Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
 – aspirate the content of all wells into a liquid waste container;
 – add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 – aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 – repeat the washing step 4 more times;
 – after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control ($\overline{Nc}$), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}} / CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +	$OD \geq 1,0$
-----------	---------------

CONTROL -	$OD \leq 0,150$
-----------	-----------------

CONTROL -	$\bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0$	where Ncn is the OD for each Nc run
-----------	---	-------------------------------------

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and Nc is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{\text{sample}} > 1,1$	POSITIVE
$0,9 \leq IP_{\text{sample}} \leq 1,1$	BORDERLINE*
$IP_{\text{sample}} < 0,9$	NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
G12	0,554	1,59	4,1
E5	0,999	2,87	6,5

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
G12	0,561	1,56	4,7
E5	1,042	2,89	7,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 μ mol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

Studies of the characteristics of the method in comparison with a similar commercial ELISA kit were performed on a sample of characterized serum and a sample of donor blood serum. The relative sensitivity of «EQUI anti-Trichinella spiralis» ELISA kits was determined from 24 serum samples that were tested for antibodies to anti-*Trichinella spiralis* and characterized as positive in a commercial ELISA kit. According to the results of the analysis, the relative sensitivity of the sets «EQUI anti-Trichinella spiralis» is 96%. Studies of the relative specificity of ELISA kits were performed on a complex of 260 samples of donor blood serum. According to research, the relative specificity of ELISA kits «EQUI anti-Trichinella spiralis» was 98%. The coincidence percentage is 97.86%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI anti-Trichinella spiralis» ELISA kit is an indication that the patient has antibodies of IgG and / or IgA specific to *Trichinella spiralis*. The presence of IgG antibodies in newborns is not an evidence of *Trichinella spiralis* invasion.

Uncertain results may indicate a history of *Trichinella spiralis* invasion. A negative result in the «EQUI anti-Trichinella spiralis» ELISA kit indicates the absence of antibodies specific for *Trichinella spiralis* in the test sample or the concentration of specific antibodies below the sensitivity limit of the analysis. Specific antibodies may not be detected at the beginning of the clinical manifestations of the invasion. In this case, it is recommended to re-obtain and test serum samples from patients with clinical signs of trichinosis in one or two weeks.

The final diagnosis cannot be established solely on the basis of serological test results. When making a diagnosis the results of a set of laboratory and instrumental studies, should be taken into account as well as clinical manifestations of the disease. Cross-reactions with antibodies to antigens of other helminths cannot be completely ruled out. To exclude a false-positive result, it is recommended to conduct a verification study of positive samples by immunoblotting.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
High background in all wells	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$

Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL][SAMPLE] into the wells (purple)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL][+], B1, C1, D1 – [CONTROL][-],
other wells – examined samples
(change of colour from purple to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN][CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN][TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN][STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\overline{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \overline{Nc} + 0,3;$$

$$IP_{sample} = OD_{sample} / CO$$

$\overline{Nc}$ - the average value of OD 3-x [CONTROL][-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{sample} > 1,1$	POSITIVE
$0,9 \leq IP_{sample} \leq 1,1$	BORDERLINE
$IP_{sample} < 0,9$	NEGATIVE

Ascaris lumbricoides IgG

ELISA kit for the qualitative detection of IgG
antibodies to *Ascaris lumbricoides*

Instructions for use

IVD

REF
EI-601

Σ 96
tests

CE

EQUI *Ascaris lumbricoides* IgG

ELISA kit for the qualitative detection of
IgG antibodies to *Ascaris lumbricoides*

1. INTENDED USE

The «EQUI *Ascaris lumbricoides* IgG» is ELISA kit intended to qualitatively detect anti-*Ascaris lumbricoides* IgG in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose lumbricosis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: children, rural people, summer visitors.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Ascaris lumbricoides is a human parasite resulting in lumbricosis — one of the most common helminthiases in the world. By some estimates, over a milliard of people infested with acaricides are on earth.

Human ascaris belongs to *Nematoda* roundworms infesting the small intestine of a man who is its exclusive host. *Ascaris lumbricoides* eggs are excreted in the environment with faeces of the infested man. In a warm, wet soil, ascaris larvae develops in the eggs, therefore eggs become invasive only after a maturation period (2 to 3 weeks at 25–30 °C, lower temperatures require longer term). After infestation, larvae leave eggs in the human intestine, penetrates blood circulation and migrate to the liver and lungs with blood flow. The larvae move to the pharynx from the lungs, and here they are re-ingested and further enter the small intestine. In 2 to 3 months, adult ascaris able to propagate develops from larvae in the small intestine.

The helminths are transferred by faecal-oral route upon injection of mature eggs of *Ascaris lumbricoides* with soil-contaminated vegetables, fruits, water, as well as through dirty hands after contact with soil. Lumbricosis is conditionally divided into the early stage (migration of larvae) and late stage (parasitism of adults in the intestine). Invasion is asymptomatic in most cases. Primary feeling of being unwell occurs as early as several days after infestation and is accompanied by weakness, abdominal pain, nausea. Migration of larvae to the lungs may manifest as rales and cough. In some cases, intense invasion may result in pneumonia and liver damage. However, the most common symptom of early lumbricosis are allergic reactions due to hypersensitivity to metabolic products of larvae.

Late stage manifests as decreases appetite, abdominal pain, vomiting, diarrhoea, constipation. Massive ascaris invasion may result in the intestinal obstruction with a lump of helminths or rupture of the walls with peritonitis. When ascarides penetrate other organs, complications may develop such as

hepatitis, cholangitis, pancreatitis and even asphyxia. Cases of neurological disorders sometimes develop in lumbricosis, namely: headache, irritability, sleep impairment, inattention, etc. If no timely treatment is started for intense invasion, it may lead to death, especially in younger children.

Strong immune response to *Ascaris lumbricoides* invasion develops as early as at the early stage. It includes cellular and humoral immunity. Antigens of ascaris larvae stimulate secretion of all-class specific immunoglobulins, however, the level of specific and total IgE antibodies is the highest. The intensity of the immune response (including increased IgG titres) correlates with the massiveness of the invasion.

For diagnosis of lumbricosis, parasitologic stool test for presence of ascaris larvae and eggs is the most common. X-ray imaging of the lungs is additionally applied at the early stage of invasion. Complete blood count (eosinophilia develops in lumbricosis) and detection of serum anti-*Ascaris lumbricoides* antibodies also is included in the set of exams. The presence of specific anti-ascaris antibodies may suggest asymptomatic invasion, and allows initiation of treatment before complications develop in conjunction with other diagnostic instruments.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Ascaris lumbricoides* IgG in «EQUI Ascaris lumbricoides IgG» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Ascaris lumbricoides* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Ascaris lumbricoides* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgG monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS

1 x 96
wells

Microplate

Each plate well is coated with *Ascaris lumbricoides* antigen. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months

		Positive control
CONTROL +	1 x 0,25 ml	Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
		Negative control
CONTROL -	1 x 0,6 ml	Negative human serum with a preservative (yellow). Store at 2-8 °C
		Serum dilution solution
DIL SAMPLE	1 x 13 ml	Buffer solution with a milk extract, a detergent and a preservative (brown). Store at 2-8 °C
		Conjugate solution (ready to use)
SOLN CONJ	1 x 13 ml	Buffer solution of monoclonal antibodies to human IgG, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
		TMB solution (ready to use)
SOLN TMB	1 x 13 ml	TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
		Washing solution TWEEN (20x concentrated)
TWEEN WASH 20x	1 x 50 ml	20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
		Stop Solution (ready to use)
SOLN STOP	1 x 13 ml	0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620–695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI Ascaris lumbricoides IgG» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;

- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI *Ascaris lumbricoides* IgG» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;

- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

- 9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.
- 9.2. Fill in the sample application plan.
- 9.3. Prepare the detergent as per para. 8.2.
- 9.4. Add 90 µL of [DIL|SAMPLE] into each plate well.
- 9.5. Add 10 µL of controls and test samples into the wells:
[CONTROL|+] – into well A1,
[CONTROL|-] – into wells B1, C1 and D1,
and test samples into the remaining wells.
At the time of adding, the solution changes its colour from brown to blue. Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
– aspirate the content of all wells into a liquid waste container;
– add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
– aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
– repeat the washing step 4 more times;
– after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN|CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN|TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN|STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN|TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only [SOLN|TMB] and [SOLN|STOP] must be added

in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control ($\bar{Nc}$), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \bar{Nc} + 0,3$$

$$IP_{\text{sample}} = OD_{\text{sample}}/CO, \text{ where } OD_{\text{sample}} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

$$\boxed{\text{CONTROL} +} \quad OD \geq 1,0$$

$$\boxed{\text{CONTROL} -} \quad OD \leq 0,150$$

$$\boxed{\text{CONTROL} -} \quad \bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0 \quad \text{where } Ncn \text{ is the OD for each } Nc \text{ run}$$

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and $\bar{Nc}$ is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$$\begin{array}{ll} IP_{\text{sample}} > 1,1 & \text{POSITIVE} \\ 0,9 \leq IP_{\text{sample}} \leq 1,1 & \text{BORDERLINE*} \\ IP_{\text{sample}} < 0,9 & \text{NEGATIVE} \end{array}$$

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated in 24 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
547	0,504	1,43	2,9
671	0,753	2,13	3,6
413	1,165	3,30	3,1

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
547	0,534	1,55	5,0
671	0,750	2,17	4,6
413	1,159	3,36	3,6

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To evaluate clinical sensitivity and specificity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits, 55 serum samples from patients with clinical symptoms typical for lumbricosis and 60 serum samples from patients without clinical manifestations (seronegative in terms of *Ascaris lumbricoides*) were used. Clinical sensitivity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits was 94.55 % and clinical specificity — 93.3 %.

Method characteristics in comparison with equal commercial ELISA kit was studied in target paediatric population (160 samples) and population of donors (346 samples). For paediatric population serum, relative specificity of «EQUI *Ascaris lumbricoides* IgG» ELISA kits was established at the level of 97.92 % and percent agreement was 95.51 %. For donor population serum, relative specificity of was 89.74 %, relative specificity — 96.30 % and percent agreement was 95.47 %.

12. LIMITATIONS OF ASSAY

Positive result in «EQUI *Ascaris lumbricoides* IgG» ELISA kit supports presence of anti-*Ascaris lumbricoides* specific IgG antibodies. Presence of this class antibodies in newborns is not an evidence of *Ascaris lumbricoides* invasion.

Inconclusive results may suggest a history of *Ascaris lumbricoides* invasion.

Negative result of «EQUI *Ascaris lumbricoides* IgG» ELISA kit supports the absence of anti- *Ascaris lumbricoides* IgG specific antibodies in the test sample or concentration of specific antibodies is below the sensitivity limit of the assay.

The results of serological test only are not the basis for final diagnosis. When establishing the diagnosis, the results of complex laboratory and instrumental tests, as well as clinical manifestations should be considered. Cross-reactions with antibodies to antigens of other helminths cannot be fully ruled out.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

REFERENCES

1. CDC - Ascariasis // <https://www.cdc.gov/parasites/ascariasis/index.html>.
2. Cooper P. J., Chico M. E. et al. Human Infection with *Ascaris lumbricoides* Is Associated with a Polarized Cytokine Response // *The Journal of Infectious Diseases*. - 2003. - Vol. 182 (4). - P. 1207–1213.
3. Gupta S., Kumar S. et al. *Ascaris lumbricoides*: an unusual aetiology of gastric perforation // *Journal of Surgical Case Reports*. - 2012. - Vol. 2012. - rjs008.
4. Li Q., Zhao D. et al. Life-threatening complications of ascariasis in trauma patients: a review of the literature // *World Journal of Emergency Medicine*. - 2014. - Vol. 5 (3). - P. 165–170.
5. McSharry C., Xia Y. et al. Natural Immunity to *Ascaris lumbricoides* Associated with Immunoglobulin E Antibody to ABA-1 Allergen and Inflammation Indicators in Children // *Infection and Immunity*. - 1999. - Vol. 67(2). - P. 484–489.
6. Palmer L. J., Celedón J. C. et al. *Ascaris lumbricoides* Infection Is Associated with Increased Risk of Childhood Asthma and Atopy in Rural China // *American Journal of Respiratory and Critical Care Medicine*. - 2002. - Vol. 165, No. 11. - P. 1489–1493.
7. Shalaby N. Effect of *Ascaris lumbricoides* infection on T helper cell type 2 in rural Egyptian children // *Therapeutics and Clinical Risk Management*. - 2016. - Vol. 12. - P. 379–385.
8. WHO. Water related diseases: ascariasis. 2013 // http://www.who.int/water_sanitation_health/diseases/ascariasis/en/.
9. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU.
10. Закон України «Про відходи» // Відомості Верховної Ради України. - 1998. - №36-37.
11. Наказ МОЗ України №325 від 08.06.2015 «Про затвердження Державних санітарно-протиепідемічних правил і норм щодо поводження з медичними відходами».
12. Постанова КМУ від 02 жовтня 2013р. №754 «Про затвердження технічного регламенту щодо медичних виробів для діагностики in vitro».
13. Hanna Tolonen, Kari Kuulasmaa, Tiina Laatikainen, Hermann Wolf and the European Health Risk Monitoring Project. Recommendation for indicators, international collaboration, protocol and manual of operations for chronic disease risk factor surveys Part 4.Storage and transfer of serum/plasma samples// Finnish National Public Health Institute 2002// https://thl.fi/publications/ehrm/product2/part_iii4.htm
14. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Annex 3.Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results// Geneva: World Health Organization; 2012 Dec.

Manufacturer

Authorized Representative in the European Community

In vitro diagnostic medical device

Catalogue number

Date of manufacture

Use by date

Batch code

Temperature limit

Contains sufficient for <n> tests

Caution

Non-Sterile

Consult instructions for use

Keep away from sunlight

Keep dry

Compliance with EU safety requirements

Edition 8, 10.02.2022

For questions and suggestions regarding the ELISA kit contact:

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL][SAMPLE] into the wells (brown)

Add to 10 µl of controls and samples into the wells:

A1 – [CONTROL][+], B1, C1, D1 – [CONTROL][-],

other wells – examined samples

(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN][CONJ] into all wells (green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN (300 µl per well)

Add 100 µl of [SOLN][TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN][STOP] into all wells (change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at 450/620-695 nm

CALCULATION OF RESULTS

$$\overline{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \overline{Nc} + 0,3;$$

$$IP_{sample} = OD_{sample} / CO$$

$\overline{Nc}$ - the average value of OD 3-x [CONTROL][-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{sample} > 1,1$	POSITIVE
$0,9 \leq IP_{sample} \leq 1,1$	BORDERLINE
$IP_{sample} < 0,9$	NEGATIVE

HAV IgM

**ELISA kit for the qualitative detection of IgM
antibodies to hepatitis A virus**

Instruction for use

IVD

REF
EI-031

Σ 96
tests

CE

EQUI HAV IgM

ELISA kit for the qualitative detection
of IgM antibodies to hepatitis A virus

1. INTENDED USE

The «EQUI HAV IgM» is ELISA kit intended for the qualitative detection of IgM antibodies to hepatitis A virus in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) to diagnose acute hepatitis A. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: blood or organ donors; pregnant women and children born to infected mothers; patients with symptoms of liver disease.

Usage: ELISA kit is used in clinical diagnostic laboratories, blood transfusion stations, as well as in other institutions working in the field of *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

One of the most common foodborne infections is hepatitis A. The hepatitis A virus (HAV) causes acute liver disease, which can be mild or severe. Unlike hepatitis B and C, this hepatitis does not become chronic, but can lead to acute liver failure, which is characterized by high mortality.

Hepatitis A virus is a small shellless RNA virus from the *Picornaviridae* family. It is characterized by being highly stable in different environments and can be stored at + 4 ° C for several months, but becomes inactive after 5 minutes at 100 ° C. The virus replicates in liver cells, and then is released through the bile into the environment with the fecal masses of the infected person. The cellular immune response to HAV infection leads to the destruction of hepatocytes, liver dysfunction and the development of symptoms, typical for other types of hepatitis.

Acute hepatitis A, even with the clinical manifestations, does not differ from other viral hepatitis. Therefore, serological markers of infection are used for diagnosis, namely the detection of specific antibodies to HAV antigens. IgM antibodies are detected in the serum 1-2 weeks after infection, with the onset of symptoms or a few days before. In the maximum titer of anti-CAA IgM are detected in the jaundice period, after that their level gradually decreases. In most patients, specific IgM ceases to be detected after 6 months, and may occasionally circulate in the blood for more than a year. IgG antibodies to HAV antigens begin to be released shortly after IgM antibodies and stay in the blood throughout life. Also, specific IgG is produced after a vaccination. Detection of IgG antibodies to CAA indicates the formation of a stable immunity due to infection or immunization.

3. ANALYSIS PRINCIPLE

Detection of specific IgM antibodies to hepatitis A virus in the «EQUI HAV IgM» ELISA kit is based on the principle of «IgM capture» of solid-phase ELISA in a two-stage incubation. Monoclonal antibodies specific for human IgM immunoglobulins are adsorbed into the wells of the plate. During the first step of incubation of the

test samples in the wells of the ELISA plate, IgM immunoglobulins, if present in the samples, bind to monoclonal antibodies in the solid phase. The wells are washed to remove unbound components, leaving only specific antibody-antibody complexes. A mixture of HAV antigen and peroxidase conjugate of HAV-specific antibodies that bind to solid-phase immune complexes is then added. Unbound components are removed during washing. Immune complexes are detected by adding a solution of chromogen 3,3', 5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After a 30-minute incubation, the reaction is halted by adding a stop solution. Optical density (OG) in the wells is determined on a spectrophotometer at a wavelength of 450 / 620-695 nm. The intensity of the yellow color is proportional to the number of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

STRIPS	1 x 96 wells	Microplate Each plate well is coated with monoclonal antibodies specific for human IgM immunoglobulins. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months
CONTROL +	1 x 0,25 ml	Positive control The solution of human IgM immunoglobulins crosslinked with monoclonal antibodies specific for horseradish peroxidase, with a preservative (pink). Store at 2-8 °C
CONTROL -	1 x 0,6 ml	Negative control Negative human serum with a preservative (yellow). Store at 2-8 °C
DIL SAMPLE	1 x 13 ml	Serum dilution solution Buffer solution with monoclonal antibodies to human IgG, milk extract, detergent and preservative (brown). Store at 2-8 °C
CONJ 11x	1 x 1,3 ml	Conjugate (11x concentrated) 11-fold concentrate of conjugate of antibodies to hepatitis A virus with horseradish peroxidase in buffer solution with stabilizers (purple). Dilute the conjugate (11x) 1:11 with the conjugate dilution solution before use (eg 100 µl concentrate + 1 ml conjugate dilution solution, enough for 8 wells). Diluted solution should be stored at 2-8 °C for no more than 1 day
DIL CONJ	1 x 13 ml	Conjugate dilution solution Buffer solution of inactivated hepatitis A virus antigen with detergent and preservative (yellow). Store at 2-8 °C
SOLN TMB	1 x 13 ml	TMB solution (ready to use) TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C

TWEEN WASH 20x	1 x 50 ml	Washing solution TWEEN (20x concentrated) 20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
SOLN STOP	1 x 13 ml	Stop Solution (ready to use) 0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620-695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI HAV IgM» ELISA kit;
- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- **SOLN|TMB** solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of **SOLN|TMB** with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for conjugate solution and **SOLN|TMB**;
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI HAV IgM» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and conjugate solution with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;
- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by

centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

8.3. Conjugate solution preparation

Working dilution of the conjugate is prepared as follows: dilute **CONJ|11x** (purple) in a clean vial of solution **DIL|CONJ** (yellow) in the ratio 1:11 (ie, 1 + 10), the solution turns green. For example, for 8 well analysis add to 1 ml **DIL|CONJ** 100 µl **CONJ|11x**. The solution of the conjugate in the working dilution is stable during the day when stored at 2-8 °C.

9. ASSAY PROCEDURE

9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame.

Be sure to add control wells in every test run.

9.2. Fill in the sample application plan.

9.3. Prepare the detergent as per para. 8.2.

9.4. Add 90 µL of **DIL|SAMPLE** into each plate well.

9.5. Add 10 µL of controls and test samples into the wells:

CONTROL|+ – into well A1,

CONTROL|- – into wells B1, C1 and D1,

and test samples into the remaining wells.

At the time of adding, the solution changes its colour from brown to blue.

Pipette the mix in the wells carefully to avoid foaming.

9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.

- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
- aspirate the content of all wells into a liquid waste container;
 - add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
 - aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
 - repeat the washing step 4 more times;
 - after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Prepare conjugate solution as per para. 8.3.
- 9.9. Add 100 µL of conjugate solution into each well. Cover the strips with a new piece of adhesive film and incubate for **60 minutes at 37 °C**.
- 9.10. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.11. Add 100 µL of **SOLN|TMB** into the wells; do not touch the bottom and the walls of the plate wells.
- 9.12. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.13. Add 100 µL of **SOLN|STOP** into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding **SOLN|TMB**. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.14. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

*Measurement at the single wavelength of 450 nm is possible, in that case, it is needed to leave one well for blank (only **SOLN|TMB** and **SOLN|STOP** must be added in blank well).*

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD of the negative control ($\overline{Nc}$), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\overline{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \overline{Nc} + 0,3$$

$$IP_{sample} = OD_{sample} / CO, \text{ where } OD_{sample} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +

$$OD \geq 1,2$$

CONTROL -

$$OD \leq 0,150$$

CONTROL -

$$\bar{Nc} \times 0,5 \leq Ncn \leq \bar{Nc} \times 2,0 \quad \text{where } Ncn \text{ is the OD for each } Nc \text{ run}$$

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and $\bar{Nc}$ is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$$IP_{\text{sample}} > 1,1$$

POSITIVE

$$0,9 \leq IP_{\text{sample}} \leq 1,1$$

BORDERLINE*

$$IP_{\text{sample}} < 0,9$$

NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated in 32 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
37s/2	1,576	4,85	5,3
24s	2,462	7,57	4,3

Inter assay reproducibility

The coefficient of variation (CV) for two sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
37s/2	1,600	4,78	6,3
24s	2,463	7,36	7,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 μmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

The diagnostic characteristics of the ELISA kit were evaluated by examining a set of 30 samples containing IgM antibodies to hepatitis A virus, a set of donor serum samples (188 samples) and samples of PHT202 Anti-Hepatitis A Virus (HAVed) panel. Performance Panel (contains 8 positive and 13 negative samples) - a total of 239 samples - was compared to similar commercial kits. For this set (239 samples) the relative sensitivity of the «EQUI HAV IgM» ELISA kit is 100%, the relative specificity - 100%, the percentage of coincidence - 100%.

12. LIMITATIONS OF ASSAY

A positive result in the «EQUI HAV IgM» ELISA kit is the evidence that the patient has IgM antibodies specific for hepatitis A virus. Anti-CAA specific IgM antibodies are usually markers of active replication of hepatitis A virus.

In order to counteract the false-positive results caused by the presence of autoantibodies specific for class G immunoglobulins (rheumatoid factor) in human serum samples, the kit uses a special block component that prevents the formation of immune complexes with anti-human antibodies in the solid phase.

The final diagnosis cannot be established solely on the basis of serological test results. When making a diagnosis the results of a set of laboratory and instrumental studies, as well as clinical manifestations of the disease should all be taken into account.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer

<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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Manufacturer

Authorized Representative in the European Community

In vitro diagnostic medical device

Catalogue number

Date of manufacture

Use by date

Batch code

Temperature limit

Contains sufficient for <n> tests

Caution

Non-Sterile

Consult instructions for use

Keep away from sunlight

Keep dry

Compliance with EU safety requirements

Edition 6, 01.11.2021

For questions and suggestions regarding the ELISA kit contact:

Obelis s.a.
Bd Général Wahis 53
1030 Brussels
Belgium
Tel: +(32)2 732-59-54
Fax: +(32)2 732-60-03
mail@obelis.net

Ekvittestlab LLC
Velyka Vasylykivska St. 114, Kyiv, Ukraine, 03150
Tel: 0(800)31-89-87, +38 (044)334-89-87,
e-mail: info@equitest.com.ua, www.equitest.com.ua

ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL][SAMPLE] into the wells
(brown)

Add to 10 µl of controls and samples into the wells:

A1 – [CONTROL+], B1, C1, D1 – [CONTROL-],

other wells – examined samples

(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of prepared 1:11 (1+10) conjugate solution into all wells
(green)

Cover strips with an adhesive film, incubate for **60 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN][TMB] into all wells

Incubate for **30 min in the dark at 18-25°C**

Add 100 µl of [SOLN][STOP] into all wells

(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at
450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{sample} = OD_{sample} / CO$$

$\bar{Nc}$ - the average value of OD 3-x [CONTROL-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS


$IP_{sample} > 1,1$	POSITIVE
$0,9 \leq IP_{sample} \leq 1,1$	BORDERLINE
$IP_{sample} < 0,9$	NEGATIVE

Toxocara canis IgG

ELISA kit for the qualitative detection of IgG
antibodies to *Toxocara canis*

Instructions for use

IVD

REF
EI-603

Σ 96
tests

CE

EQUI *Toxocara canis* IgG

ELISA kit for the qualitative detection of
IgG antibodies to *Toxocara canis*

1. INTENDED USE

The «EQUI *Toxocara canis* IgG» is ELISA kit intended to qualitatively detect anti-*Toxocara canis* IgG in human serum or plasma by enzyme-linked immunosorbent assay (ELISA) in order to diagnose toxocariasis. The testing procedure is designed for both manual arrangement with automatic pipettes and standard equipment, and for automated «open» immunoassay analysers.

Target group: children, pet owners, rural people, summer visitors, forest guards, veterinarians.

Usage: ELISA kit is used in clinical diagnostic laboratories and other institutions engaged in *in vitro* diagnostics.

2. CLINICAL SIGNIFICANCE

Toxocariasis is a common disease induced by *Toxocara* helminth which is transmitted from animals to human. Toxocariasis is spread throughout the world, however, it is more common in depressed areas with poor hygienic conditions. In some regions, up to 90 % of puppies and up to 10 % of adult domesticated dogs are infested with toxocara. The risk of infestation is higher for owners of cats and dogs and for children due to playing in the sandpits and on the playgrounds contaminated with animal faeces.

Toxocara are threadworms belonging to *Nematoda*. Human conditions are mostly caused by *Toxocara canis*, which infested canids, rare - *Toxocara cati*, which is more common in felids. Adult toxocara in the body of infested animals reaches 5–15 cm in length; their propagation takes place here. Female helminths lay about 200 thous eggs daily, which are excreted in the environment with faeces. If conditions are favourable, following several weeks of maturation in the soil they become invasive — a larva is developed in the eggs. In the paratenic host (mice, poultry, cows, pigs, etc.). larva develops without propagation. If the conditions are unfavourable, larvae are encapsulated and may maintain viability for a long time (up to 10 years). They may also be the source of invasion.

People are infested through faecal-oral route when ingesting *Toxocara canis* mature eggs with soil-contaminated vegetables, fruits, berries, via dirty hands or when consuming meat of paratenic hosts. In the small intestine, larvae leave their cover and penetrates blood circulation through the intestinal walls. The larvae migrate to other organs and tissues with blood, namely: liver, lungs, muscles, eyes, CNS, etc. In the most of the infested, toxocariasis is asymptomatic. Clinical manifestations of this disease are associated with the site of larvae migration and depend on the intensity of invasion and age of the host. Visceral syndrome larva migrans is typical after infestation of the internal organs with *Toxocara canis* and ocular

toxocariasis, when eye and optic nerve are involved. Symptoms of visceral toxocariasis: fever, fatigue, abdominal pain, anorexia, hepatomegaly, cough and others. Heart and respiratory failure may develop in severe cases. Due to a strong immune response to larvae antigens, immediate and delayed hypersensitivity reactions develop. Granulomatosis in ocular toxocariasis may result in retinal detachment and loss of vision.

Diagnosis of toxocariasis is complicated due to the lack of specific manifestations of the disease, even upon intense invasion. Furthermore, a man is an intermediate host of *Toxocara canis* and does not excrete parasites in the environment, whereas it is difficult to localise larvae in certain organs via non-invasive methods. Eosinophilia may appear in blood tests, however, serological tests are more common to detect toxocariasis (immunofluorescence reaction, ELISA and immunoblotting). Detection of specific anti-*Toxocara canis* IgG to larvae antigens may suggest current or previous invasion. High titer of IgE antibodies is also typical for active invasion. However, the combination of clinical manifestations and laboratory findings are necessary for diagnosis.

3. ANALYSIS PRINCIPLE

The procedure of testing for anti-*Toxocara canis* IgG in «EQUI *Toxocara canis* IgG» ELISA kit is based on «indirect» solid-phase ELISA with a two-stage incubation. Antigens of *Toxocara canis* larvae are entrapped in the wells. During the first step of incubation of ELISA plate wells with test samples, specific anti-*Toxocara canis* antibodies (if present in the samples) bind to the solid-phase antigens. The wells are washed to remove unbound antibodies and have only specific antigen-antibody complexes left. Then, a conjugate of anti-species IgG monoclonal antibodies with horseradish peroxidase is added, which binds to solid-phase immune complexes. Unbound components are removed by washing. Antigen-antibody complexes are detected by adding a solution of chromogen 3,3',5,5'-tetramethylbenzidine (TMB) with hydrogen peroxide. After 30-minute incubation, the reaction is stopped by adding the stop solution. The optical density (OD) in the wells is determined using a spectrophotometer at 450/620-695 nm. The intensity of the yellow colour is proportional to the level of antibodies in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Contents of the ELISA kit

Microplate

Each plate well is coated with *Toxocara canis* larval antigens. The wells are detachable. After the first opening, store unused strips in the package at 2-8 °C for a maximum of 6 months

STRIPS

1 x 96 wells

		Positive control
CONTROL +	1 x 0,25 ml	Conjugated specific monoclonal antibody solution with preservative (pink). Store at 2-8 °C
		Negative control
CONTROL -	1 x 0,6 ml	Negative human serum with a preservative (yellow). Store at 2-8 °C
		Serum dilution solution
DIL SAMPLE	1 x 13 ml	Buffer solution with a milk extract, a detergent and a preservative (brown). Store at 2-8 °C
		Conjugate solution (ready to use)
SOLN CONJ	1 x 13 ml	Buffer solution of monoclonal antibodies to human IgG, conjugated with horseradish peroxidase, with stabilizers and preservative (green). Store at 2-8 °C
		TMB solution (ready to use)
SOLN TMB	1 x 13 ml	TMB solution, H ₂ O ₂ , a stabilizer, a preservative (colourless). Store at 2-8 °C
		Washing solution TWEEN (20x concentrated)
TWEEN WASH 20x	1 x 50 ml	20-fold phosphate buffer concentrate with Tween-20 (colourless). Dilute TWEEN detergent (20x) at 1:20 with distilled or deionized water (e. g., 5 mL of concentrate + 95 mL of water for 8 wells) before use. Store the diluted solution at 2-8 °C for a maximum of 7 days
		Stop Solution (ready to use)
SOLN STOP	1 x 13 ml	0.5 mol H ₂ SO ₄ solution (colourless). Store at 2-8 °C

The ELISA kit also includes adhesive films (2 items), sample application plan (1 item), checklist, and instruction for use.

4.2. Optional reagents, materials and equipment

Automatic single and multichannel pipettes 10–1000 µL, tips, volumetric laboratory glassware (10–1,000 mL), deionized or distilled water, thermostat at 37 °C, automatic or semi-automatic plate washer, spectrophotometer (reader) for microplates at 450/620–695 nm, appropriate containers for potentially contaminated waste, timer, filter paper, disposable powder-free gloves, disinfectants.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

Be sure to read the instructions for use carefully before the test. The validity of the test results depends on strict following of the test procedure.

- do not use the ELISA kit components after the expiry date;
- do not use for analysis or mix components of different batches, components of kits for different nosologies, or reagents from other manufacturers with the «EQUI Toxocara canis IgG» ELISA kit;

- do not freeze the ELISA kit or its contents;
- after using a reagent, close each vial with its cap;
- when washing, control filling and complete aspiration of solution from the wells;
- use a new pipette tip each time you add samples or reagents;
- prevent direct sunlight from reaching the reagents from the ELISA kit;
- [SOLN|TMB] solution must be colourless before use. Do not use the solution if its colour is blue or yellow. Avoid contact of [SOLN|TMB] with metals or metal ions. Use only clean glassware thoroughly rinsed with distilled water;
- do not use reagents with colour not in line with para. 4.1;
- under no circumstances should the same glassware be used for [SOLN|CONJ] and [SOLN|TMB];
- do not evaluate the test results visually (without a reader);
- any optional equipment that is in direct contact with biological material or kit components should be considered contaminated and requires cleaning and decontamination;
- the ELISA kit includes materials for 96 tests. Dispose of the used components as well as any remaining unused components.

5.2. Safety requirements

- all reagents in the ELISA kit are for laboratory professional use for *in vitro* diagnosis only and may only be used by qualified personnel;
- conduct the tests in disposable powder-free gloves and goggles only;
- do not eat, drink, smoke, or apply make-up in the test room;
- do not mouth-pipette the solutions;
- controls from the «EQUI *Toxocara canis* IgG» ELISA kit have been tested and found to be for anti-HIV1/2, anti-HCV and anti-*Treponema pallidum* antibodies and HBsAg negative; however, controls and test samples should be handled as potentially hazardous infectious materials;
- some of the kit components contain low concentrations of harmful substances and can damage skin or mucosa. In case of contact of [SOLN|TMB], [SOLN|STOP] and [SOLN|CONJ] with mucous membranes or skin, immediately wash the affected area with plenty of water;
- in case of spillage of acid-free solutions, e. g. sera, treat the surface with a disinfectant solution and then wipe dry with filter paper. Otherwise first neutralize acid with sodium bicarbonate solution and then wipe the surface dry as described above.

5.3. Waste inactivation and disposal

- the liquid waste must be inactivated, for example, with hydrogen peroxide solution at the final concentration of 6% for 3 hours at room temperature, or with sodium hypochlorite at the final concentration of 5% for 30 minutes, or with other approved disinfectants;

- the solid waste must be inactivated by autoclaving at a temperature not less than 132°C;
- do not autoclave the solutions that contain sodium azide or sodium hypochlorite;
- disposal of inactivated waste must be conducted due to national laws and regulations.

6. STORAGE AND STABILITY

ELISA kit is stable up to the expiry date stated on the label when stored at 2-8°C. The kit should be transported at 2-8°C. Single transportation at a temperature up to 23°C for two days is possible.

7. SAMPLE COLLECTION, TRANSPORTATION AND STORAGE GUIDELINES

Collect blood from the vein into the sterile test tube. Test tube must be marked with patient ID and date of sample collecting. Blood before serum separation can be stored at 2-8 °C for 24 hours, avoiding freezing.

Serum or plasma can be stored at 2-8 °C for maximum 3 days. Frozen serum can be stored for longer periods of time at -20 °C or -70 °C. Thaw frozen samples and keep them at room temperature for 30 minutes before use. After thawing, the stir samples to achieve homogeneity. Avoid repeated freezing-thawing cycles for test samples. If serum (or plasma) is turbid, remove insoluble inclusions by centrifugation at 3000 rpm for 10-15 minutes. Do not use serum samples with hyperlipidemia, hemolysis, and bacterial growth.

Transport serum samples in insulated containers. To do that, put closed labelled tubes in a plastic bag, tightly seal it and place in the centre of an insulated container. Put the frozen cold packs on the bottom, along the side walls of the insulated container and on top of the serum samples.

8. REAGENT PREPARATION

NOTE! It is very important to keep all ELISA kit components for at least 30 min at room temperature 18-25 °C before the assay!

8.1. Microplate preparation

To prevent water condensation in the wells, keep the **STRIPS** for 30 minutes at a room temperature before opening. Open the vacuum pack, detach the appropriate number of wells, and carefully pack the remaining wells with a desiccant and store tightly zip-locked at 2-8 °C. Storing the packed plate this way ensures its stability for 6 months.

8.2. Washing solution preparation

To prepare detergent, dilute **TWEEN|WASH|20x** at 1:20 (1+19) with distilled or deionized water and stir. E. g., 5 mL of concentrate + 95 mL of water, which is enough for 8 wells. If there are crystals present in the detergent concentrate, heat the vial at 37 °C until the crystals dissolve completely (15–20 minutes). Store the diluted solution at 2-8 °C for a maximum of 7 days.

9. ASSAY PROCEDURE

- 9.1. Prepare the necessary number of wells (four wells for controls and a necessary number of wells for test samples) and insert them into the ELISA plate frame. Be sure to add control wells in every test run.
- 9.2. Fill in the sample application plan.
- 9.3. Prepare the detergent as per para. 8.2.
- 9.4. Add 90 µL of [DIL][SAMPLE] into each plate well.
- 9.5. Add 10 µL of controls and test samples into the wells:
[CONTROL +] – into well A1,
[CONTROL -] – into wells B1, C1 and D1,
and test samples into the remaining wells.
At the time of adding, the solution changes its colour from brown to blue. Pipette the mix in the wells carefully to avoid foaming.
- 9.6. Cover the strips up with adhesive film and incubate for **30 minutes at 37 °C**.
- 9.7. Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
– aspirate the content of all wells into a liquid waste container;
– add a minimum of 300 µl of diluted washing solution to each well, soak each well for 30 seconds;
– aspirate the content of all wells again. The residual volume after every aspiration should be less than 5 µl;
– repeat the washing step 4 more times;
– after the final aspiration, eliminate extra moisture by tapping the plate against a piece of filter paper.
- 9.8. Add 100 µL of [SOLN][CONJ] into each well. Cover the strips with a new piece of adhesive film and incubate for **30 minutes at 37 °C**.
- 9.9. Following incubation, remove the film carefully and wash the wells five times as described in para. 9.7.
- 9.10. Add 100 µL of [SOLN][TMB] into the wells; do not touch the bottom and the walls of the plate wells.
- 9.11. Incubate the strips for **30 minutes** in a dark place at a room temperature of 18-25 °C. Do not use adhesive film at this stage.
- 9.12. Add 100 µL of [SOLN][STOP] into each strip well to stop the enzymatic reaction; adhere to the same sequence of actions as when adding [SOLN][TMB]. At the time of adding, the solution colour changes from blue to yellow, and clear solution slightly changes its shade.
- 9.13. Measure the optical density (OD) of the wells at 450/620-695 nm wavelength using an ELISA microplate reader within 5 minutes after stopping the reaction. Pay attention to the cleanness of the plate bottom and the absence of bubbles in the wells before reading.

Measurement at the single wavelength of 450 nm is impossible, in that case, it is needed to leave one well for blank (only [SOLN][TMB] and [SOLN][STOP] must be added

in blank well).

10. CALCULATION AND INTERPRETATION OF RESULTS

10.1. Calculation of results

Calculate the average OD for the negative control ($\overline{Nc}$), Cut off (CO) and a sample positivity index (IP_{sample}).

$$\overline{Nc} = (Nc1 + Nc2 + Nc3)/3; \quad CO = \overline{Nc} + 0,3$$

$$IP_{sample} = OD_{sample} / CO, \text{ where: } OD_{sample} \text{ is the OD sample.}$$

10.2. Quality control (assay validation)

The test results are considered valid if they meet the following requirements:

CONTROL +

 $OD \geq 1,0$

CONTROL -

 $OD \leq 0,150$

CONTROL -

 $\overline{Nc} \times 0,5 \leq Ncn \leq \overline{Nc} \times 2,0$ where Ncn is the OD for each Nc run

If any of the OD values for the negative control is beyond the above interval, it should be discarded, and $\overline{Nc}$ is calculated based on the remaining OD values for the negative control. If several OD values for the negative control fail to meet the above requirements, the test is considered invalid and requires a new run.

10.3. Interpretation of results

$IP_{sample} > 1,1$	POSITIVE
$0,9 \leq IP_{sample} \leq 1,1$	BORDERLINE*
$IP_{sample} < 0,9$	NEGATIVE

* Uncertain samples are recommended to be re-examined in two wells of the ELISA kit. If the results are again uncertain, a new sample should be selected and analyzed in 2-4 weeks. In case of repeated indeterminate results, such samples shall be considered negative.

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical performance characteristics

Precision of measurement

Intra assay repeatability

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated in 24 replicates on one series of ELISA kits.

Sample No.	OD _{av}	IP _{av}	CV, %
669	0,927	2,81	4,8
544	1,503	4,56	1,4
666	1,694	5,14	4,5

Inter assay reproducibility

The coefficient of variation (CV) for three sera with different levels of specific antibodies was evaluated for 4 days in 4 sets of analysis, 8 replicates in each analysis.

Sample No.	OD _{av}	IP _{av}	CV, %
669	1,016	3,04	4,7
544	1,516	4,54	1,9
666	1,683	5,04	4,1

Analytical specificity

The test results are not affected by bilirubin at up to 0.21 mg/mL (361.8 µmol/L), haemoglobin at up to 10 mg/mL and triglycerides at up to 10 mg/mL (11.3 mmol/l) present in the sample.

11.2. Diagnostic characteristics

To evaluate diagnostic characteristics of «EQUI *Toxocara canis* IgG» ELISA kits, 78 serum samples from patients with clinical symptoms typical for toxocariasis and 60 serum samples from patients without clinical manifestations (seronegative in terms of *Toxocara canis*) were used. Clinical sensitivity of «EQUI *Toxocara canis* IgG» ELISA kits was 98.7 %, clinical specificity — 96.7 %.

Method characteristics in comparison with equal commercial ELISA kit was studied in target paediatric population (160 samples) and population of donors (298 samples). For paediatric population serum, relative specificity of «EQUI *Toxocara canis* IgG» ELISA kits was established at the level of 99.28 % and percent agreement was 97.45 %. For donor population serum, relative specificity of was 89.19 %, relative specificity — 93.55 % and percent agreement was 91.73 %.

12. LIMITATIONS OF ASSAY

Positive result in «EQUI *Toxocara canis* IgG» ELISA kit supports presence of anti-*Toxocara canis* specific IgG antibodies. Presence of this class antibodies in newborns is not an evidence of *Toxocara canis* invasion.

Inconclusive results may suggest a history of *Toxocara canis* invasion.

Negative result of «EQUI *Toxocara canis* IgG» ELISA kit supports the absence of anti-*Toxocara canis* specific IgG antibodies in the test sample or concentration of specific antibodies is below the sensitivity limit of the assay.

The results of serological test only are not the basis for final diagnosis. When establishing the diagnosis, the results of complex laboratory and instrumental tests, as well as clinical manifestations should be considered. Cross-reactions with antibodies to antigens of other helminths cannot be fully ruled out.

13. DIFFICULTIES THAT CAN OCCUR DURING THE ASSAY PROCEDURE

Possible reasons	Solution
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head and rinse according to the instructions for use
Poor quality or contaminated water	Use purified water with specific resistance $\geq 10 \text{ M}\Omega \cdot \text{cm}$
Use of poorly washed glassware	Use chemically clean utensils
Use of chlorinated disinfectants	Do not use chlorine disinfectants
Use of contaminated tips	Use new tips
Increased incubation times or change in the temperature conditions	Adhere to the incubation regime according to the instructions for use
<i>High background in a row of wells</i>	
Repeat application of TMB solution	TMB solution should be applied once
Contamination of the automatic pipette nozzle with conjugate solution	Clean the pipette and dial carefully liquid
Contamination of one of the washer's channel	Clean the flush channel, rinse washer
<i>Received OD of the positive control is below the border value</i>	
One of the reagents (conjugate solution or TMB solution) was not prepared in a correct way or was not added	Re-conduct ELISA, pay attention to the correctness of the introduction of these reagents
Reduced incubation times at any stage	Incubate according to instructions for use
<i>The colour density of the wells fails to meet the obtained optical density value</i>	
This may suggest that the optical beam has been displaced	Check the correct operation of the reader

14. TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

In case of technical problems, you can obtain assistance by contacting the manufacturer.

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ASSAY PROCEDURE SCHEME

Keep all reagents for 30 min at temperature 18-25°C before use

Dispense 90 µl [DIL][SAMPLE] into the wells
(brown)

Add to 10 µl of controls and samples into the wells:
A1 – [CONTROL][+], B1, C1, D1 – [CONTROL][-],
other wells – examined samples
(change of colour from brown to blue)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN][CONJ] into all wells
(green)

Cover strips with an adhesive film, incubate for **30 min at 37°C**

Rinse the wells 5 times with prepared 1:20 (1+19) washing solution TWEEN
(300 µl per well)

Add 100 µl of [SOLN][TMB] into all wells

Incubate for **30 min** in the dark at **18-25°C**

Add 100 µl of [SOLN][STOP] into all wells
(change of colour from blue to yellow)

Measure the optical density (OD) with an ELISA microplate reader at
450/620-695 nm

CALCULATION OF RESULTS

$$\bar{Nc} = (Nc1 + Nc2 + Nc3)/3;$$

$$CO = \bar{Nc} + 0,3;$$

$$IP_{sample} = OD_{sample} / CO$$

$\bar{Nc}$ - the average value of OD 3-x [CONTROL][-]

CO - Cut off

IP_{sample} - sample positivity index

INTERPRETATION OF RESULTS

$IP_{sample} > 1,1$	POSITIVE
$0,9 \leq IP_{sample} \leq 1,1$	BORDERLINE
$IP_{sample} < 0,9$	NEGATIVE

Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)

INTENDED USE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is a rapid lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of Malaria P.falciparum specific histidine rich protein-2 (Pf HRP-II) and Malaria P.vivax specific lactate dehydrogenase (Pv-LDH) in human blood specimen as an aid in the diagnosis of Malaria infection. It is for *In-Vitro* Diagnostic use only.

INTRODUCTION

Malaria is a serious, sometimes fatal, parasitic disease characterized by fever, chills, and anemia and is caused by a parasite that is transmitted from one human to another by the bite of infected Anopheles mosquitoes. There are four kinds of malaria that can infect humans: Plasmodium falciparum, P. vivax, P. ovale, and P. malariae. In humans, the parasites (called sporozoites) migrate to the liver where they mature and release another form, the merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At the present, malaria is diagnosed by looking for the parasites in a drop of blood. Blood will be put onto a microscope slide and stained so that the parasites will be visible under a microscope.

PRINCIPLE

The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) contains a membrane, which is precoated with mouse monoclonal antibodies specific to HRP-II of P. falciparum on test line Pf region and with mouse monoclonal antibodies specific to lactate dehydrogenase of P. vivax species on test line Pv region respectively. Conjugate pad is dispensed with monoclonal antibodies conjugated to colloidal gold, which are specific to P. falciparum histidine rich protein-2 (Pf HRP-II) and specific to the lactate dehydrogenase of P. vivax.

During the assay, an adequate volume of the blood specimen is dispensed into the sample well (S) of the test cassette, a lysis buffer is added to the buffer well (B). The buffer contains a detergent that lyses the red blood cells and releases various antigens, which migrate by capillary action across the strip held in the cassette. Pv-LDH if presents in the specimen will bind to the Pv-LDH-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-Pv-LDH antibody, forming a burgundy colored Pv band, indicating a Pv positive test result.

Alternatively, pHRP-II if presents in the specimen will bind to the pHRP-II-gold conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-pHRP-II antibodies, forming a burgundy colored Pf band, indicating a Pf positive test result.

Absence of any T bands suggests a negative result. The test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex of goat anti- mouse IgG I mouse IgG (anti-Pv-LDH and anti-pHRP-II)-gold conjugates regardless of the color development on any of the T bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

MATERIALS SUPPLIED

25 Sealed pouches each containing a test cassette, a dropper and a desiccant
1 Buffer, 7.0 mL
1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

1. Clock or timer
2. Collection by venipuncture: collection tube (containing EDTA, citrate or heparin)
3. Collection using a lancet: sterile lancet

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened, preferably at 2°C-30°C. Do not expose the kit over 30°C. Do not freeze the kit. Ensure that the test device is brought to room temperature before opening. The test device is stable through the expiration date printed on the sealed pouch if it is stored at 2°C-30°C.

WARNINGS AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only. Do not use after expiration date.
2. The instruction must be followed exactly to get accurate results. Failure to follow the insert gives inaccurate test results.
3. Do not eat, drink or smoke in the area where the specimens or kits are handled.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.

5. Hemolized blood may be used for the testing, but do not take precipitants.
6. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
7. Humidity and temperature can adversely affect results.
8. Do not perform the test in a room with strong air flow, ie. an electric fan or strong airconditioning.

SPECIMEN COLLECTION

Collection by venipuncture:

- 1) Collect whole blood into a collection tube (containing EDTA, citrate or heparin) by venipuncture.
- 2) If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage periods greater than three days, freezing is recommended. They should be brought to room temperature prior to use. Using the specimen after long-term storage of more than three days can cause non-specific reaction.
- 3) When stored at 2-8°C, the whole blood sample should be used within three days.

Collection using a lancet:

- 1) Clean the area to be lanced with an alcohol swab.
- 2) Squeeze the end of the fingertip and pierce with a sterile lancet.
- 3) Wipe away the first drop of blood with sterile gauze or cotton.
- 4) Using the dropper provided, while gently squeezing the tube, immerse the open end in the blood drop and then gently release the pressure to draw blood into the dropper.

TEST PROCEDURE

Allow the test device, specimen, buffer, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.


1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface. Be sure to label the device with specimen's ID number.
3. With a 5 µL mini plastic dropper provided, draw whole blood specimen to exceed the specimen line as showed in the following image and then transfer drawn whole blood into the sample well (S). Then add 3 drops (about 120 µL) of Lysis Buffer to the buffer well (B) immediately.

Note: Practice a few times prior to testing if you are not familiar with the mini dropper. For better precision, transfer specimen by pipette capable to deliver 5 µL of volume.

4. Set up timer.

If preferred, after 5 minutes of adding specimen and buffer, you may add one more drop of Lysis Buffer to help the background become clearer.

5. Results can be read in 20 to 30 minutes. It may take more than 20 minutes to have the background become clearer. Don't read results after 30 minutes. To avoid confusion, discard the test cassette after interpreting the result.

Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood)

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE:

P.f Positive: One line appears in the control region, and one line appears in P.f. line region.

P.v Positive: One line appears in the control region and one line appears in Pv line region.

P.f and P.v Positive: One line appears in the control region, one line appears in Pv line region and one line appears in P.f. line region.

NEGATIVE: Only one colored line appears in the control region.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) is for *in vitro* diagnostic use only. This test should be used for the detection of P.f and P.v antigens in whole blood specimens only. Neither the quantitative value nor the rate of increase in P.f and P.v concentration can be determined by this qualitative test.
2. The Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) will only indicate the presence of antigens of P.f and / or P.v in the specimen and should not be used as the sole criterion for the diagnosis of malaria infection.
3. As known relevant interference, haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples can lead to impair the test results.
4. The test is limited to the detection of antigen to Malaria Plasmodium sp. Although the test is very accurate in detecting HRP-II specific to P.f or pLDH specific to P.v, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained.
5. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of malaria infection.
6. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Clinical Performance for P.f Ag test:

A total of 352 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	50	4	54
	Negative	0	298	298
Total Results		50	302	352

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 98.9%

2. Clinical Performance for P.v Ag test:

A total of 289 samples from susceptible subjects were tested by the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) and by thick blood smear test.

Method		Smear Test		Total Results
Malaria Pf/Pv Ag Rapid Test	Results	Positive	Negative	
	Positive	63	3	66
	Negative	0	223	223
Total Results		63	226	289

Relative Sensitivity: 100%

Relative Specificity: 98.7%

Overall Agreement: 99.0%

3. Precision: Within-run and between-run have been determined by the testing 10 replicates of four specimens: a negative, a low positive, a medium positive and a strong positive. All values were correctly identified 100% of the time.

4. Interference: To evaluate the interference of Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood) with known relevant interfering specimens, the haemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples were investigated. In these studies, those specimens did not interfere with the Malaria P.f./P.v. Ag Rapid Test Cassette (Whole Blood).

REFERENCE

1. Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras.: Plasmodium falciparum and Plasmodium vivax: Lactate Dehydrogenase Activity and its Application for *in vitro* Drug Susceptibility Assay. Experimental Parasitology 80, 260-271 (1995)
2. David L. Vander Jagt, Lucy A. Hunsaker and John E. Heidrich : Partial Purification and Characterization of Lactate Dehydrogenase from Plasmodium falciparum. Molecular and Biochemical Parasitology, 4 (1981) 255-264
3. David J. Bzik, Barbara A. Fox and Kenneth Gonyer : Expression of Plasmodium falciparum lactate dehydrogenase in Escherichia coli Molecular and Biochemical Parasitology, 59(1993) 155-166
4. Histidine-Rich Protein II: a Novel Approach to Malaria Drug Sensitivity Testing ANTIMICROBIAL AGENTS AND CHEMOTHERAPY, June 2002, p. 1658©1664 Vol. 46, No. 6

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co., Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China
Tel: +86-572-5226111 Fax: +86-572-5226222
Website: www.orientgene.com

EC REP

QARAD BV
Cipalstraat 3, 2440 Geel BELGIUM

REF

GCMAL(pf/pv)-402a

Revision Date: 2022-09-26
B20885-03

Fecal Occult Blood Rapid Test Cassette (Feces)

INTENDED USE

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid chromatographic immunoassay for the qualitative detection of human occult blood in feces by professional laboratories or physician's offices. It is useful to detect bleeding caused by a number of gastrointestinal disorders, e.g., diverticulitis, colitis, polyps, and colorectal cancer.

Fecal Occult Blood Rapid Test Cassette (Feces) is recommended for use in 1) routine physical examinations, 2) hospital monitoring for bleeding in patients, and 3) screening for colorectal cancer or gastrointestinal bleeding from any source.

INTRODUCTION

Most of diseases can cause hidden blood in the stool. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood. Traditional guaiac-based method lacks sensitivity and specificity, and has diet-restriction prior to the testing.

Fecal Occult Blood Rapid Test Cassette (Feces) is a rapid test to qualitatively detect low levels of fecal occult blood in feces. The test uses double antibody-sandwich assay to selectively detect as low as 50 ng/mL of hemoglobin or 6 µg hemoglobin/g feces. In addition, unlike the guaiac assays, the accuracy of the test is not affected by the diet of the patients.

PRINCIPLE

Fecal Occult Blood Rapid Test Cassette (Feces) is a lateral flow chromatographic immunoassay based on the principle of the double antibody-sandwich technique. The membrane is pre-coated with anti-hemoglobin antibodies on the test line region of the device. During testing, the specimen reacts with the colloidal gold coated with anti-hemoglobin antibodies. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-hemoglobin antibodies on the membrane and generate a colored line. The presence of this colored line in the test region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS PROVIDED

20 Test cassettes
20 Specimen collection tubes with buffer
1 Package insert

MATERIALS REQUIRED BUT NOT PROVIDED

1. Specimen collection containers 2. Clock or timer

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test is not stable out of the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. This package insert must be read completely before performing the test. Failure to follow the insert gives inaccurate test results.
3. Do not use it if the tube/pouch is damaged or broken.
4. Test is for single use only. Do not re-use under any circumstances.
5. **Do not use specimen with visible blood for the testing.**
6. Handle all specimens as if they contain infectious agents. Observe established standard procedure for proper disposal of specimens.
7. Specimen extraction buffer contains Sodium Azide (0.1%). Avoid contact with skin or eyes. Do not ingest.
8. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
9. Humidity and temperature can adversely affect results.
10. Do not perform the test in a room with strong air flow, i.e. electric fan or strong air conditioning.

PATIENT PREPARATION

1. A specimen should not be collected from a patient with following conditions that may interfere with the test results:

- Menstrual bleeding
 - Bleeding hemorrhoids
 - Constipating bleeding
 - Urinary bleeding.
2. Dietary restrictions are not necessary.
 3. Alcohol and certain medications such as aspirin, indomethacin, phenylbutazone, reserpine, cortocosteroids, and nonsteroidal anti-inflammatory drugs may cause gastrointestinal irritation and subsequent bleeding, thus gives positive reactions. On the advice of the physician, such substances should be discontinued at least 48 hours prior to testing.

SPECIMEN COLLECTION AND PREPARATION

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.


1. Collect a random sample of feces in a clean, dry receptacle.
2. Unscrew the top of the collection tube and remove the applicator stick.
3. Randomly pierce the fecal specimen in at least five (5) different sites.
4. Remove excess sample off the shaft and outer grooves. Be sure sample remains on inside grooves.
5. Replace the stick in the tube and tighten securely.
6. Shake the specimen collection bottle so that there is proper homogenisation of feces in buffer solution.

Note: Specimens prepared in the specimen collection tube may be stored at room temperature (15-30°C) for 3 days maximum, at 2-8°C for 7 days maximum or at -20°C for 3 months maximum if not tested within 1 hour after preparation.

TEST PROCEDURE

Allow the test cassette, specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean, flat surface.
3. Shake the specimen collection tube several times.
4. Hold the specimen collection tube upright and then unscrew and open the upper cap.
5. Squeeze 3 drops (~90 µL) of the sample solution in the sample well of the cassette and start the timer.
6. Wait for the colored line(s) to appear. Read results in 5 minutes. Do not interpret the result after 5 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

Positive: Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

Negative: One colored line appears in the control line region (C). No line appears in the test line region (T).

Invalid: Control line fails to appear. The test should be repeated using a new cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

NOTE:

1. The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered positive. Note that this is a qualitative test only, and

Fecal Occult Blood Rapid Test Cassette (Feces)

cannot determine the concentration of analytes in the specimen.

2. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

QUALITY CONTROL

An internal procedural control is included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. This test kit is to be used for the qualitative detection of human hemoglobin in fecal samples. A positive result suggests the presence of human hemoglobin in fecal samples. In addition to intestinal bleeding the presence of blood in stools may have other causes such as hemorrhoids, blood in urine etc.
2. Not all colorectal bleedings are due to precancerous or cancerous polyps. The information obtained by this test should be used in conjunction with other clinical findings and testing methods, such as colonoscopy gathered by the physician.
3. Negative results do not exclude bleeding since some polyps and colorectal region cancers can bleed intermittently or not at all. Additionally, blood may not be uniformly distributed in fecal samples. Colorectal polyps at an early stage may not bleed.
4. Urine and excessive dilution of sample with water from toilet bowl may cause erroneous test results. The use of a receptacle is recommended.
5. Feces specimens should not collect during the menstrual period and not three day before or afterwards, at bleeding due to constipation, bleeding haemorrhoids, or at taking rectally administered medication. It could cause false positive results.
6. This test may be less sensitive for detecting upper g.i. Bleeding because blood degrades as it passes through the g.i. Track.
7. The Fecal Occult Blood Rapid Test Cassette (Feces) is to aid diagnosis and is not intended to replace other diagnostic procedures such as G.I. fibroscope, endoscopy, colonoscopy, or X-ray analysis. Test results should not be deemed conclusive with respect to the presence or absence of gastrointestinal bleeding or pathology. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source for the occult blood in the feces.

PERFORMANCE CHARACTERISTICS

1. Sensitivity: 99.6%

Fecal Occult Blood Rapid Test Cassette (Feces) can detect the levels of human occult blood as low as 50 ng/mL hemoglobin or 6 µg hemoglobin/g feces.

2. Prozone Effect:

It is observed that this FOB test can detect 2 mg/mL hemoglobin.

3. Specificity: 99.9%

Fecal Occult Blood Rapid Test Cassette (Feces) is specific to human hemoglobin. Specimen containing the following substances at the standard concentration was tested on both positive and negative controls and showed no effects on test results at standards concentration.

Substances	Concentrations (Diluted with the extraction buffer)
Beef hemoglobin	2 mg/mL
Chicken hemoglobin	0.5 mg/mL
Pig hemoglobin	0.5 mg/mL
Goat hemoglobin	0.5 mg/mL
Horse hemoglobin	20 mg/mL
Rabbit hemoglobin	0.06 mg/mL

REFERENCES

1. Simon J.B. Occult Blood Screening for Colorectal Carcinoma: A Critical Review, Gastroenterology, Vol. 1985;88:820.
2. Blebea J. and Ncpherson RA. False-Positive Guaiac Testing With Iodine, Arch Pathol Lab Med, 1985;109:437-40.

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co.,Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China
Tel: +86-572-5226111 Fax: +86-572-5226222
Website: www.orientgene.com

Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany

GEFOB-602b

H. pylori Ag Rapid Test Cassette (Feces)

INTENDED USE

H. pylori Ag Rapid Test Cassette (Feces) is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of H. Pylori antigen in feces. It is for professional *in vitro* diagnostic use only.

INTRODUCTION

H. Pylori is associated with a variety of gastrointestinal diseases included non-ulcer dyspepsia, duodenal and gastric ulcer and active, chronic gastritis.^{1,2} The prevalence of H. pylori infection could exceed 90% in patients with signs and symptoms of gastrointestinal diseases. Recent studies indicate an association of H. Pylori infection with stomach cancer.³ H. Pylori colonizing in the gastrointestinal system elicits specific antibody responses^{4,5,6} which aids in the diagnosis of H. Pylori infection and in monitoring the prognosis of the treatment of H. Pylori related diseases. Antibiotics in combination with bismuth compounds have been shown to be effective in treating active H. Pylori infection. Successful eradication of H. pylori is associated with clinical improvement in patients with gastrointestinal diseases providing a further evidence.⁷

PRINCIPLE

H. pylori Ag Rapid Test Cassette (Feces) is a lateral flow chromatographic immunoassay based on the principle of the double antibody-sandwich technique. The test cassette consists of: 1) a burgundy colored conjugate pad containing H. Pylori antibodies conjugated with color particles (H. Pylori conjugates). 2) a nitrocellulose membrane strip containing a test band (T band) and a control band (C band). The T band is pre-coated with non-conjugated H. Pylori antibodies.

When an adequate volume of test specimen is dispensed into the sample well of the cassette, the specimen migrates by capillary action across the cassette. The antigen of H. Pylori if present in the specimen will bind to the H. Pylori antibodies conjugates. The immunocomplex is then captured on the membrane by the pre-coated H. Pylori antibodies, forming a burgundy colored T band, indicating a H. Pylori antigen positive test result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred. Otherwise, the test result is invalid and the specimen must be retested with another device.

PRODUCT CONTENTS

H. pylori Ag Rapid Test Cassette (Feces) containing anti- H. pylori antibodies particles and anti-H. pylori antibodies coated on the membrane.

MATERIALS SUPPLIED

20 Sealed pouches each containing a test cassette and a desiccant
20 Specimen collection tubes with extraction buffer, 2.0 mL
1 Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

1. Clock or timer
2. Specimen collection containers.

STORAGE AND STABILITY

All reagents are ready to use as supplied. Store unused test device unopened at 2°C-30°C. If stored at 2°C-8°C, ensure that the test device is brought to room temperature before opening. The test is not stable out off the expiration date printed on the sealed pouch. Do not freeze the kit or expose the kit over 30°C.

WARNINGS AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. Do not use it if the tube/pouch is damaged or broken.
3. Test is for single use only. Do not re- use under any circumstances.
4. Handle all specimens as if they contain infectious agents. Observe established standard procedure for proper disposal of specimens
5. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assay.
6. Humidity and temperature can adversely affect results

SPECIMEN COLLECTION

Collect sufficient quantity of feces (1-2 mL or 1-2 g) in a clean, dry specimen collection container to obtain maximum antigens (if present). Best results will be obtained if the assay is performed within 6 hours after collection. Specimen collected may be stored for 3 days at 2-8°C if not tested within 6 hours. For long term storage, specimens should be kept below -20°C.


To process fecal specimens:

• For Solid Specimens:

Unscrew the cap of the specimen collection tube, then randomly stab the specimen collection applicator into the fecal specimen in at least 3 different sites to collect approximately 50 mg of feces (equivalent to 1/4 of a pea). Do not scoop the fecal specimen.

• For Liquid Specimens:


Hold the dropper vertically, aspirate fecal specimens, and then transfer 2 drops (approximately 80 µL) into the specimen collection tube containing the dilution buffer. Screw on and tighten the cap onto the specimen collection tube, then shake the specimen collection tube vigorously to mix the specimen and the dilution buffer. Leave the tube alone for 2 minutes.

TEST PROCEDURE

1. Remove the test device from its foil pouch by tearing along the notch and use it as soon as possible.
2. Specimen collection. See also specimen collection.
3. Holding the sample collection device upright, carefully break off the tip of collection device.
4. Squeeze 2 drops (~80 µL) of the sample solution in the sample well of the cassette, as in the illustration.
5. Read the test results in 10 minutes. It is important that the background is clear before the result is read. Do not read results after 10 minutes. To avoid confusion, discard the test device after interpreting the result.

INTERPRETATION OF RESULTS

H. pylori Ag Rapid Test Cassette (Feces)

Positive: Two lines appear. One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

Negative: One colored line appears in the control line region(C). No line appears in the test line region (T).

Invalid: Control line fails to appear.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The Assay Procedure and the Assay Result Interpretation must be followed closely when testing the presence of H. Pylori antigen in feces from individual subjects. Failure to follow the procedure may give inaccurate results.

2. H. pylori Ag Rapid Test Cassette (Feces) is limited to the qualitative detection of H. Pylori antigen in feces. The intensity of the test band does not have linear correlation with the antigen titer in the specimen.

3. A negative result for an individual subject indicates absence of detectable H. Pylori antigen. However, a negative test result does not preclude the possibility of exposure to or infection with H. Pylori.

4. A negative result can occur if the quantity of the H. Pylori antigen present in the specimen is below the detection limits of the assay, or the antigen that are detected are not present during the stage of disease in which a sample is collected.

5. The results obtained with this test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

PERFORMANCE CHARACTERISTICS

A study was performed with 165 patient feces samples including both symptomatic gastrointestinal disorders and samples from non-symptomatic patients and 100 normal feces samples.Comparison for all subjects with H. pylori Ag Rapid Test Cassette (Feces) and reference ELISA kit is showed in the following table:

Method		EIA		Total Results
H.P Test Cassette	Results	Positive	Negative	
	Positive	163	0	163
	Negative	2	100	102
Total Results		165	100	265

Relative sensitivity: 98.8%

Relative specificity: 100%

Accuracy:98.9%

REFERENCE

1. Marshall,B.J.et.al. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia.149:439-44, 1985.

2. Marshall,B.J.et.al. Prospective double-blind trial of duodenal ulcer relapse after eradication of Campylobacter pylori. Lancet. Dec.1437-42,1988.

3. Megraud,F.et.al. Seroepidemiology of Campylobacter pylori infection in virious populations J.Clin.Microbiology. 27:1870-3,1989.

4. Soll,A.H. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med.322:909-916,1990.

5. Parsonnet,J.et.al. Helicobacter pylori infection and the risk of gastric carcinoma. New England J.Med. 325:1127-31,1991.

6. Ansong,R. et.al. Evaluation of techniques for isolation, subcultivation and preservation of Helicobacter pylori. J.Clin.Micro. 29:51-53,1991.

7. Pronovost,A.P.et.al. Evaluation of a new immunodiagnostic assay for Helicobacter pylori antibody detection: Correlation with histopathological and microbiological results. J.Clin.Microbiol.32:46-50,1994.

INDEX OF SYMBOLS

	Consult instructions for use		Tests per kit		Authorized Representative
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2~30°C		Lot Number		Catalog#

Zhejiang Orient Gene Biotech Co.,Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China.
TEL: +86-572-5226111 FAX: +86-572-5226222
Website: www.orientgene.com

Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany

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Vitrotest Helicobacter pylori IgA

ELISA test kit for the qualitative and semiquantitative determination of IgA class antibodies to *Helicobacter pylori*

REF ELO99-96

IVD

Σ 96

1. INTENDED USE

The test kit Vitrotest *Helicobacter pylori* IgA is an enzyme linked immunosorbent assay (ELISA) for the qualitative and semiquantitative determination of IgA antibodies to *Helicobacter pylori* in human serum or plasma.

The test kit might be applied for the ELISA using both automatic pipettes and standard equipment as well as open system automated ELISA analyzers.

2. CLINICAL VALUE

Helicobacter pylori is a widespread microorganism infecting half of the world's population. Its prevalence is very high in developing countries and is quite low in the developed world. According to the World Gastroenterology Organization, infection of adults in Eastern Europe and Asia is 70-80%.

Studies during recent decades have shown the key role of bacterium *Helicobacter pylori* in the pathogenesis of stomach and duodenal lesions. *H. pylori* is detected almost in 100 % of adult patients with duodenal ulcer, approximately in 80 % of patients with peptic ulcer, in 92 % of patients with gastric cancer and in 92 % of patients with active chronic gastritis. Research has demonstrated that elimination of helicobacter leads to the disappearance of gastritis and significant reduction in the incidence of duodenal ulcer recurrence.

Helicobacteriosis is a chronic infection with long, often asymptomatic course. Its symptoms do not differ from clinical manifestations of gastro-duodenitis since usual constant pain in the epigastrium occurs. *H. pylori* is often present in patients with no clinical manifestations of disease.

H. pylori strains are extremely heterogeneous and divided into two groups: strains that express both antigens VacA and CagA (type I) and strains that do not express these antigens (type II). Strains of the first group dominate in patients with peptic ulcer disease and gastric cancer. CagA protein penetrates the epithelial cells of the mucous membrane, causes abnormal mitosis and induces chromosome instability. Organism infected with expressing CagA *H. pylori* produces antibodies, which are specific for CagA. Antibodies to CagA protein are found in 80-100 % of patients with duodenal ulcer and in 94 % of patients with gastric cancer. Thus, detection of the antibodies to CagA protein is an informative marker in the diagnosis of duodenal ulcer and gastric cancer.

Strains of *H. pylori* II type which do not express antigens CagA and VacA are associated with severe gastric and duodenal ulcers, particularly peptic ulcer disease and cancer.

H. pylori infection of the patient can be identified with both invasive and non-invasive diagnostic methods. Invasive methods include bioptic studies of the gastrointestinal tract mucosa with histological or culture methods or a rapid urease test. However, heterogeneous spread of *H. pylori* in tissues often leads to false-negative results. Non-invasive diagnostic methods include serologic studies of patient sera for the presence of specific antibodies to *H. pylori* and respiratory urease test using radioactively labeled urea. ELISA for the detection of specific IgG/IgA/IgM antibodies is the least invasive, rapid, highly sensitive and informative method for the diagnosis of *H. pylori* infection.

3. PRINCIPLE OF THE TEST

Determination of IgA antibodies to *H. pylori* in the test kit Vitrotest *Helicobacter pylori* IgA is based on a solid phase, indirect ELISA in a two-step incubation procedure. Microwells are coated a mixture of recombinant and native *H. pylori* antigens. During the first incubation step, the specific antibodies to *H. pylori*, if present in the sample, will be bound to the solid phase precoated antigens. The wells are washed to remove unbound antibodies. A secondary antibody (anti-IgA), which is conjugated to horseradish peroxidase (HRP), is added next and binds to the immune complexes on the solid phase. Unbound components are removed by washing. Antigen-antibody complexes are revealed by addition of chromogen solution containing 3,3',5,5'- tetramethylbenzidine (TMB) and hydrogen peroxide. After 15 minutes the reaction has been stopped, the absorbance values are read using a spectrophotometer at 450/620-695 nm. The colour intensity is proportional to the amount of the antibodies present in the sample.

4. MATERIALS AND EQUIPMENT

4.1. Composition of the test kit

ELISA STRIPS	1x96 wells	Microplate (12 strips x 8 wells) Each well is coated with a mixture of recombinant and native <i>H. pylori</i> antigens. The wells can be separated.
CONTROL +	1x0.5 ml	Positive control Solution of specific monoclonal immunoglobulins with preservative (pink).

CONTROL –	1x0.5 ml	Negative control Buffer solution with detergent and preservative (yellow).
CONTROL CUT-OFF	1x0.5 ml	Cut-off control Solution of specific monoclonal immunoglobulins with preservative (orange).
SAMPLE DILUENT	1x12 ml	Sample diluent Buffer solution with detergent and preservative (brown-green).
CONJUGATE SOLUTION	1x12 ml	Conjugate solution Buffer solution of monoclonal antibodies to human IgA conjugated to HRP with stabilizers and preservative (orange), ready to use.
TMB SOLUTION	1x12 ml	TMB solution TMB, H ₂ O ₂ , stabilizers, preservative (colourless), ready to use.
WASH TWEEN 20X	1x50 ml	Washing solution Tw20 (20x concentrate) 20X concentrated of PBS buffer with Tween-20 and NaCl (colourless)
STOP SOLUTION	1x12 ml	Stop Solution 0.5 mol/l H ₂ SO ₄ (colourless), ready to use.

Adhesive films (2), sera identification plan (1), instruction for use and certificate of analysis.

4.2. Material required but not provided

- variable volume automatic pipettes (10 µl–1000 µl) and disposable pipette tips;
- plate reader (single wavelength 450 nm or dual wavelength 450/620–695 nm);
- volumetric laboratory glassware (10–1000 ml);
- distilled/DI water;
- incubator thermostatically controlled at 37 °C;
- automatic/semiautomatic plate washer;
- appropriate waste containers for potentially contaminated materials;
- timer;
- absorbent paper;
- disposable gloves;
- disinfectants;
- protective clothes.

5. PRECAUTIONS AND SAFETY

5.1. Precautions

The ELISA assays are time and temperature sensitive. Strictly follow the test procedure and do not modify it.

- do not use expired reagents;
- do not use for analyses and do not mix reagents from different lots or from test kits of different nosology as well as other manufacturer's reagents with Vitrotest kits;

*Note: it is possible to use **WASH TWEEN 20X**, **TMB SOLUTION** and **STOP SOLUTION** from other Vitrotest ELISA kits.*

- close reagents after use only with appropriate caps;
- control the filling and full aspiration of the solution in the wells;
- use a new tip for each sample and reagent;
- avoid exposure of kit reagents to direct sunlight;
- **TMB SOLUTION** must be colourless before use. If **TMB SOLUTION** is blue or yellow it cannot be used. Avoid any contact of **TMB SOLUTION** with metals or metal ions. Use glassware thoroughly washed and rinsed with distilled/DI water;
- never use the same glassware for **CONJUGATE SOLUTION** and **TMB SOLUTION**.

The manufacturer is not responsible or liable for any incorrect results and/or incidents taking place as a result of any violation of the instruction. The manufacturer is not responsible for visual readings of samples (without using a plate reader).

5.2. Safety

- all components of test kit are intended for *in vitro* diagnostic use only;
- all materials of human or animal origin should be regarded and handled as potentially infectious;
- the ELISA is only designed for qualified personnel;

- disposable gloves and safety glasses must be worn at all times while performing analysis;
- never eat, drink, smoke or apply cosmetics in the assay laboratory;
- never pipette solutions by mouth;
- controls do not contain of human origin components;
- avoid contact with [STOP SOLUTION] containing 0.5 mol/l H_2SO_4 . It may cause skin irritation and burns;
- some components of the test kit contain low concentrations of harmful compounds and could cause irritation of the skin and the mucosa. In the case of contact of [TMB SOLUTION], [STOP SOLUTION] or [CONJUGATE SOLUTION] with skin or mucosa, the place of contact should be immediately rinsed with large amounts of water;
- in case of spilling of solutions that do not contain acid, e.g. sera, rinse the surface with disinfectant, then dry it with absorbent paper. In other case acid first must be neutralized by sodium bicarbonate and then wiped out as described above;
- for information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request.

5.3. Waste treatment

Patient specimens, controls, and incubated microplate strips should be treated as infectious waste, residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

6. STORAGE AND STABILITY

Reagents are stable until stated expiration date on the label when stored refrigerated (2-8°C). Do not freeze. The kit should be shipped at 2-8°C. Single transportation at the temperature up to 23°C for two days is possible.

After the first opening of the packaging, the components of the ELISA kits are stable within 3 months, except for those specified in p. 8 of this Instruction.

7. SPECIMEN COLLECTION

The fresh serum or plasma (EDTA, lithium-heparin, sodium citrate, potassium fluoride) samples can be stored for 3 days at 2-8 °C, or frozen for longer periods at -20 – -70 °C. Frozen samples must be thawed and kept at room temperature for at least 30 min before use. Do not use preheated samples. Mix thawed samples thoroughly to homogeneity. Avoid repeated freezing/thawing. Samples containing aggregates must be clarified by centrifugation (3000 rpm for 10-15 min). Do not use hyperlipemic, hyperhaemolysed or contaminated by microorganisms serum specimens. The presence of bilirubin up to concentration of 0.21 mg/ml (361.8 µmol/l), haemoglobin up to concentration of 10 mg/ml and triglycerides up to concentration of 10 mg/ml (11.3 mmol/l) are allowed.

8. REAGENT PREPARATION

It is very important to keep all test components for at least 30 min at room temperature (18-25 °C) before the assay!

8.1. [ELISA STRIPS] preparation

Before opening the bag with the [ELISA STRIPS], keep it at room temperature for 30 minutes to avoid water condensation inside the wells. Open the vacuum bag and take out the necessary number of wells. Once opened, the bag with the remaining strips and desiccant must be **resealed with the zip-lock** immediately and kept refrigerated at 2-8°C for no more than 3 months.

8.2. Washing solution preparation

Check the [WASH TWEEN 20X] for the presence of salt crystals. If crystals have formed, re-solubilise them by warming the vial at 37°C, until crystals have been fully dissolved (15-20 min). Dilute the [WASH TWEEN 20X] 1:20 (1+19) with distilled/DI water before use. For example, 4ml concentrate + 76 ml water is sufficient for 8 wells. Once diluted it is stable at 2-8°C for 7 days.

9. ASSAY PROCEDURE

- Take out from the protective bag the support frame and the necessary number of the wells [ELISA STRIPS] (the number of specimens + 4 for controls). Place the wells into the frame. Wells with the controls must be included in every test.
- Complete the sera identification plan.
- Prepare washing solution (see 8.2.).
- Dispense 90 µl of [SAMPLE DILUENT] into each well.
- Dispense 10 µl of controls and patient samples into the wells in the following order: A1 – [CONTROL +], B1, C1 – [CONTROL CUT-OFF], and D1 – [CONTROL -], other wells – patient samples. Mix gently to avoid foaming. The colour of the sample diluent changes from brown-green to blue.
- Cover strips with an adhesive film and incubate for 30 min at 37°C.

- 9.7 Remove and discard the adhesive film and wash all wells 5 times with automatic washer or 8-channel pipette as follows:
- aspirate the contents of all wells into a liquid waste container and add immediately a minimum of 300 µl of diluted Washing solution to each well;
 - soak each well for 30 seconds between each wash cycle;
 - aspirate again. The residual volume must be lower than 5µl.
 - repeat the washing step 4 times;
 - after the final washing cycle, turn down the plate onto an absorbent paper and tap it to remove any residual buffer.
- 9.8 Dispense 100 µl of [CONJUGATE SOLUTION] per well. Cover strips with a new adhesive film, incubate for 30 min at 37°C.
- 9.9 Remove and discard the adhesive film and wash all wells five times as described above (see 9.7).
- 9.10 Dispense 100 µl [TMB SOLUTION] into all wells. Do not touch the walls and bottoms of the wells to avoid contamination.
- 9.11 Incubate the strips for 15 minutes at room temperature (18-25°C) in the dark. Do not use adhesive film in this step.
- 9.12 Add 100 µl [STOP SOLUTION] to each well in the same order and at the same rate as for [TMB SOLUTION].
- 9.13 Read the optical density (OD) of the wells at 450/620-695 nm using a microplate reader within 5 minutes after adding the [STOP SOLUTION]. Pay attention to the cleanness of the plate bottom and absence of bubbles in the wells before reading.
- Measurement in the single-wave procedure at 450 nm is possible. Reserve blank well to adjust spectrophotometer in such analysis. Only [TMB SOLUTION] and [STOP SOLUTION] must be added in blank well).*

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10.1. Validation of the test

The test run may be considered valid provided the following criteria are met:

[CONTROL +]	OD ≥ 1.2
[CONTROL CUT-OFF]	OD in a range 0.25-0.65
[CONTROL -]	OD ≤ 0.150

If one of the control cut-off absorbances does not match the above criteria, this value should be discarded and a mean value should be calculated using the remaining cut-off value. If both control cut-off absorbance do not meet the criteria, the test is invalid and must be re-tested.

10.2. Calculation of results

The cut-off (CO) is the mean optical density (OD) of the wells containing [CONTROL CUT-OFF]:

$$CO = (OD_{\text{CONTROL CUT-OFF } 1} + OD_{\text{CONTROL CUT-OFF } 2})/2;$$

The sample result is reported as a Ratio:

$$\text{Ratio}_{\text{sample}} = OD_{\text{sample}} / CO, \quad OD_{\text{sample}} - \text{optical density of the well containing sample}$$

10.3. Interpretation of results

Ratio _{sample} > 1.1	POSITIVE
0.9 ≤ Ratio _{sample} ≤ 1.1	DOUBTFUL*
Ratio _{sample} < 0.9	NEGATIVE

* If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

11. PERFORMANCE CHARACTERISTICS

11.1. Specificity and sensitivity

Sensitivity of the test kit Vitrotest Helicobacter pylori IgA was evaluated using 20 sera which contain IgG antibodies to *H. pylori*. In the test kit Vitrotest Helicobacter pylori IgA among IgG positive sera 17 samples were identified as positive. The data obtained were confirmed by an CE-marked enzyme immunoassay for the detection of IgA antibodies against *H. pylori*. In the test kit Vitrotest Helicobacter pylori IgA sensitivity was 85%.

In the study of 23 negative for antibodies to *H. pylori* sera, specificity rate of the test kit Vitrotest Helicobacter pylori IgA was above 91.3%.

11.2. Accuracy

Intra assay repeatability

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 32-replicate determinations using 1 lot of the test kit.

Serum No.	OD	Ratio	CV, %
H3	0.917	2.94	7.5
H4	0.526	1.69	5.1

Inter assay reproducibility

Coefficient of variation (CV) was calculated by measuring 2 samples with various specific antibody levels in 4 ELISA performances during 4 days, in 8-replicate determinations.

Serum No.	OD	Ratio	CV, %
H3	0.851	2.73	7.4
H4	0.529	1.70	5.7

12. LIMITATIONS OF THE PROCEDURE

A positive result in the test kit Vitrotest *Helicobacter pylori* IgA indicates the presence of specific antibodies IgA to *H. pylori*.

A negative result in the test kit Vitrotest *Helicobacter pylori* IgA indicates the absence in serum of IgA antibodies specific to *H. pylori*.

For a comprehensive serological diagnosis of *H. pylori*-infection it is recommended to test the samples for the presence IgG specific antibodies to *H. pylori* using, for example, test kit Vitrotest *Helicobacter pylori* IgG.

It should be noted that in case of early *H. pylori* infection the result may be negative due to the absence or low concentration of antibodies. In the presence of clinical manifestations of the disease it is recommended to repeat the test at least in two weeks.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis, in fact, should take into consideration as well as clinical history, symptomatology and serological data.

13. TROUBLESHOOTING

Possible causes	Solutions
<i>High background in all wells</i>	
Contaminated washer	Clean the washer head, then rinse it with 30% ethanol and distilled water
Low quality water or contaminated water	Use distilled/DI with resistivity $\geq 10 \text{ M}\Omega\cdot\text{cm}$.
Using contaminated glassware	Use clean glassware
Using chlorine based disinfectants	Use disinfectants without chlorine
Using contaminated tips	Use new tips
Increased time of incubation or temperature regimen was changed	Follow incubation regimen according to instruction for use
<i>High background in a few wells</i>	
TMB solution was added more than once	Add TMB solution once
Pipette shaft was contaminated with conjugate solution	Clean the pipette; pipette the liquids carefully
One the channels of the washer was contaminated	Clean the washer channel, clean the washer
<i>OD of the positive control below normal</i>	
Conjugate solution/tmb solution was prepared improperly or not added	Run ELISA repeatedly, prepare conjugate solution / TMB solution properly

Reduced incubation time in one of the stages

Follow incubation regimen according to the instruction for use

Visual colour intensity of the wells does not correspond to optical density

The optical beam or another component of the reader is misaligned or malfunctioning

Test the absorbance reader's performance

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Vitrotest Helicobacter pylori IgA

ASSAY PROCEDURE

Keep all reagents and specimens for at least 30 min at 18-25 °C before use

Dispense 90 µl of [SAMPLE DILUENT] into the wells
(brown-green colour)

Dispense 10 µl of controls and samples into the wells in the following order:

A1 – [CONTROL +],

B1, C1 – [CONTROL CUT-OFF],

D1 – [CONTROL -],

E1 and other wells – patient samples

(colour changes from brown-green to blue)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20
(300 µl per well)

Add 100µl of [CONJUGATE SOLUTION] into the wells
(orange colour)

Cover wells with an adhesive film, incubate for 30 min at 37 °C

Rinse the wells 5 times with diluted 1:20 (1+19) washing solution Tween-20
(300 µl per well)

Add 100 µl of [TMB SOLUTION] into the wells

Incubate for 15 min in the dark at 18-25 °C

Add 100 µl of [STOP SOLUTION]
(colour changes from blue to yellow)

Determine the optical density (OD) at 450/620-695 nm

CALCULATION

$$CO = (OD_{\text{CONTROL CUT-OFF 1}} + OD_{\text{CONTROL CUT-OFF 2}})/2;$$

$$\text{Ratio}_{\text{sample}} = OD_{\text{sample}}/CO$$

INTERPRETATION

Ratio _{sample} > 1.1	POSITIVE
0.9 ≤ Ratio _{sample} ≤ 1.1	DOUBTFUL
Ratio _{sample} < 0.9	NEGATIVE

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redkin@xema-medica.com

www.xema-medica.com

Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUALITATIVE DETERMINATION OF IgM ANTIBODIES TO *HELICOBACTER PYLORI* IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the qualitative determination of IgM antibodies to *Helicobacter pylori* in serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 46 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

H. pylori is detectable in nearly 100% of adult patients with duodenal ulcer, about 80% of patients with gastric ulcer, more than 90% of gastric cancer patients and more than 90% of patients with active chronic gastritis. An association between *H. pylori* and gastric cancer is confirmed. In developing countries, where prevalence of *H. pylori* infection is high (ca. 70-80% of population for Eastern Europe and Asia), most children become infected by the age of 10, and gastric cancer rates are very high. In the USA and other developed countries, standards of hygiene and the increasing socioeconomic status of the population have reduced the incidence of infection, and in parallel, the rates of peptic ulcers and gastric cancer have declined. There is an excellent correlation between the clinical presentation of gastritis, the presence of *H. pylori* in the stomach and elevated serum *H. pylori* IgG antibody. Moreover, eradication of *H. pylori* from a patient leads to a recovery from gastritis and a significant reduction in duodenal ulcer recurrence. *H. pylori*-specific IgG antibody also falls significantly after successful antibacterial therapy. *H. pylori* strains are classified into two broad groups - those that express both VacA (vacuolising toxin) and CagA (cytotoxin-associated gene) (type I) and those that produce neither (type II). Type I strains are predominate in patients with ulcers and cancer, while type II strains are mostly associated with asymptomatic infections. 80-100% of patients with duodenal ulcer disease produce CagA antibodies against a 128 kd antigen compared with 60-63% of *H. pylori*-infected persons with gastritis only, indicating that serologic responses to the 128 kd protein are more prevalent among *H. pylori*-infected persons with duodenal ulcers than infected persons without peptic ulceration. In *H. pylori*-infected patients who develop gastric cancer or duodenal ulcer, serum IgG against CagA are found in 94% and 80-100% of patients, resp., thus indicating that detection of antibodies to CagA is a useful marker for diagnosis of gastric cancer and duodenal ulcer. For gastritis, prevalence of IgG-antibodies to *H. pylori* is lower – ca. 60-63%.

The IgM test detects antibodies formed 2-4 weeks after infection with *H. pylori*. IgG and IgA antibodies require up to two months after infection before becoming detectable. If the immune system fails to eradicate the bacterium, it will begin to colonise the mucosal lining of the upper gastrointestinal tract. In the absence of *H. pylori* IgG, IgM positivity should be considered as a reason for treatment only in case of existing active symptoms of gastric ulcer. If treatment is not chosen in the presence of *H. pylori* IgM antibodies, it is advisable to carry out IgG testing after 2-3 months to confirm infection eradication by the immune system.

3. PRINCIPLE OF THE TEST

This test is based on indirect enzyme immunoassay principle. Tested specimen is placed into the microwells coated with CagA protein of *Helicobacter pylori*. Antibodies from the specimen bind coated CagA protein of *Helicobacter pylori* on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific IgG, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Positivity index (PI, %) is calculated by the formula (see Calculations). Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

- 4.1.** For professional use only.
- 4.2.** This kit is intended for in vitro diagnostic use only.
- 4.3. INFECTION HAZARD:** There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.
- 4.4.** Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.
- 4.5.** Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.
- 4.6.** Do not use the kit beyond the expiration date.
- 4.7.** All indicated volumes have to be performed according to the protocol. Optimal test results are obtained only when using calibrated pipettes and microplate readers.
- 4.8.** Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
- 4.9.** Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.
- 4.10.** Do not mix reagents from different lots.
- 4.11.** Replace caps on reagents immediately. Do not swap caps.
- 4.12.** Do not pipette reagents by mouth.
- 4.13.** Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.
- 4.14.** Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.
- 4.15.** The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5. KIT COMPONENTS

5.1. Contents of the Kit

	Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP	<i>Helicobacter pylori</i> IgM EIA strips, 8x12 wells	1	pcs		until exp. date
2	CONTROL - CONTROL +	Control sera (0.5 ml and 0.2 ml, resp.)	2	pcs	colourless; red	2 months
3	CONJ HRP	Conjugate, 14 ml	1	pcs	green	until exp. date
4	DIL	EIA Sample buffer 14 ml	1	pcs	purple	until exp. date
5	SUBS TMB	Substrate solution, 14 ml	1	pcs	colourless	until exp. date
6	BUF WASH 26X	Washing solution concentrate 26x, 22 ml	1	pcs	colourless	Concentrate - until exp. date
7	STOP	Stop solution, 14 ml	1	pcs	Diluted washing solution - 1 month at 2-8 °C or 5 days at RT	until exp. date
8	N003	Plate sealing tape	2	pcs	colourless	until exp. date
9	K119MI	Instruction <i>Helicobacter pylori</i> IgM EIA	1	pcs	-	N/A
10	K119MQ	QC data sheet <i>Helicobacter pylori</i> IgM EIA	1	pcs	-	N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 90-250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10-250 µl;
- Dry thermostat for 37 °C ±2 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0-3.0.

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE.

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2-8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the vials with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilution in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units:

7.5. Assay procedure:

1	Put the desired number of microstrips into the frame; allocate 4 wells for control samples CONTROL - CONTROL + (3 and 1 wells resp.) and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Pipet 90 µl of EIA Sample buffer into each well.
3	Pipet 10 µl of control samples CONTROL - CONTROL + and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 30 minutes at 37 °C.
5	Prepare washing solution by 26x dilution of washing solution concentrate BUF WASH 26X with distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
6	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
7	Incubate 30 minutes at 37 °C.
8	Wash the strips 5 times.
9	Dispense 100 µl of SUBS TMB into the wells.
10	Incubate 10-20 minutes at 18-25 °C.
11	Dispense 100 µl of STOP into the wells.
12	Measure OD (optical density) at 450 nm.
13	Set photometer blank on air.

For the assay to be valid, the following requirements should be met:

1. OD450 for CONTROL+ should be nlt 0.6 AU.
2. OD450 for CONTROL- should not be more than 0.15 AU for all replicates.

8. EXPECTED VALUES

CALCULATION OF RESULTS.

1. Calculate the mean absorbance values (OD450) for CONTROL- in triplicates and each pair of samples.
2. Calculate the cut-off value: (mean OD450 for CONTROL-) + 0.25
3. Calculate Positivity Index (PI) for each sample:

$$PI = \text{mean OD450}(\text{sample}) / \text{Cut-off}$$

If PI value is greater than 1.1, the result is POSITIVE.

If PI value is less than 0.9, the result is NEGATIVE.

If PI value is between 0.9 and 1.1, the result is EQUIVOCAL. Such samples should be retested. If the result is equivocal again, a new sample should be obtained 2-4 weeks later and tested again. If the result remains equivocal, the sample should be considered negative. Therapeutic consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutic measures. NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

11. PERFORMANCE CHARACTERISTICS

11.1 Analytical specificity.

Specificity of the test was evaluated on 93 negative serum specimens. Based on these data, specificity of the test is 98%.

11.2 Analytical sensitivity.

Specificity of kit was evaluated using internal panel "H.Pylori Positive Control".

11.3 Precision.

Intra-assay precision for two different lots (CV1, CV2) is shown below:

Serum, no	Replicates	Value, IP overage	CV1, %	CV2, %
1	32	3.4	2.9	2.8
2	32	9.1	4.8	5

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Russian Diagnostic Manufacturers Association

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RUSSIAN ASSOCIATION OF MEDICAL LABORATORY DIAGNOSTICS

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xemahelp

xemahelp@gmail.com

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma

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
ASSAY PROCEDURE

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	8
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma
Total IgE EIA

1. INTENDED USE

The Total IgE EIA kit is an enzyme immunoassay, intended for the quantitative determination of total IgE concentration in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Total immunoglobulin E (IgE) serum level is widely reported as the laboratory marker of atopic diseases such as atopic asthma, atopic dermatitis, and pollenosis. An atopic (IgE-dependent) mechanism can also underlie gastroenterocolitis, urticaria, other forms of vasculitis (including systemic), cholecystitis, vulvovaginitis, and cystitis. Part of the drug allergy (mainly to penicillin and protein drugs) also develops according to the IgE-dependent mechanism. In all of the conditions listed above, the production of high titers of specific IgE antibodies can lead to an increase in the level of total IgE in the serum. A particularly high level of total IgE is characteristic of atopic dermatitis. In addition to atopic diseases, total serum IgE is significantly increased in parasitic infestations and mycoses (especially systemic), rarely in systemic autoimmune diseases and immunodeficiency states (especially in hyper-IgE syndrome), as well as in mastocytosis (mast cell tumor) and extremely rare IgE-myeloma. A decrease in the level of total IgE in serum (below 15 IU/ml in adults) is a rare and little-studied phenomenon described in hypogammaglobulinemia, some autoimmune diseases, ulcerative colitis, and primary biliary cirrhosis.

3. TEST PRINCIPLE

The determination of the total IgE is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human IgE. Second antibodies – rabbit polyclonal antibodies to IgE conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes three stages of incubation:

- during the first stage the total IgE from the specimen is captured by the monoclonal antibodies coated onto the microwell surface;
- during the second stage horseradish peroxidase-conjugated with rabbit polyclonal antibodies bind to free epitopes of immobilized total IgE, fixed in the formed at the previous stage complexes;
- during the third stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured total IgE in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of total IgE in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P200Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to total IgE; ready to use
C200Z	CAL 1	Calibrator C1	0.8 mL	1	Solution based on phosphate buffer, free of total IgE, with preservative, ready to use (yellow liquid)
C200Z	CAL 1-5	Calibrators	0.8 mL	5	Solutions based on phosphate buffer, containing 50; 200; 500 and 1000 IU/mL of total IgE, ready to use (red liquids)
Q200Z	CONTROL	Control serum	0.8 mL	1	Solution based on human serum, containing of known total IgE content, with preservative, ready to use (colourless liquid)
T200Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of rabbit polyclonal antibodies to human total IgE conjugated to the horseradish peroxidase; ready to use (red liquid)
S011Z	DIL	EIA Buffer	14 mL	1	Buffer solution with detergent and preservative, ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	30 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (3 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- dry thermostat for 37 °C±2 °C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The Total IgE EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The Total IgE EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- EIA Buffer, Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing solution preparation

Add the contents of the 30 mL washing solution concentrate vial to 750 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the washing solution concentrate, mL	2.5	5	7.5	10	12.5	15	17.5	20	22.5	25	27.5	30
Volume of water, mL	62.5	125	187.5	250	312.5	375	437.5	500	562.5	625	687.5	750

10. ASSAY PROCEDURE

10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).

10.2 Dispense **50 µL of EIA Buffer** to all wells.

10.3 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5	SAMP13	SAMP13						
D	CAL4	CAL4	SAMP6	SAMP6	SAMP14	SAMP14						
E	CAL5	CAL5	SAMP7	SAMP7	SAMP15	SAMP15						
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 3 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Conjugate Solution** to all wells.
- 10.7 Cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.8 At the end of the incubation period, aspirate and wash each well 5 times as described in 10.5.
- 10.9 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 20 minutes**.
- 10.10 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.11 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the stop solution. Set photometer blank on CAL1.
- 10.12 Plot a calibration curve in linear coordinates: (x) is the concentration of total IgE in the Calibrators IU/mL, (y) – OD versus concentration of total IgE (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.13 Determine the corresponding concentration of total IgE in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on the results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgE. Based on data obtained by XEMA LLC, the following normal range is recommended (see below).

NOTE: values of total IgE concentrations in the tested samples that are below the LoD (3 IU/mL) and also exceed the value of the upper calibrator (1000 IU/mL) should be provided in the following form : «the total IgE concentration of tested sample X is «lower than 3 IU/mL» or «higher than 1000 IU/mL».

12.2. The calibrators concentration values of the Total IgE EIA kit are expressed in IU/mL. To calculate concentrations in ng/mL, the received concentration value in IU/mL shall be multiplied by 2.4.

$$1 \text{ IU/mL} = 2.4 \text{ ng/mL.}$$

Sex, age	Units, IU/mL		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
< 6 months	-	12	-	28.8
6-12 months	-	30	-	72.0
1-3 yrs	-	45	-	108.0
4-6 yrs	-	70	-	168.0
7-9 yrs	-	90	-	216.0
10-15 yrs	-	120	-	288.0
>15 yrs	-	130	-	312.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, IU/mL	CV, %
1	10.6	4.33
2	116.2	5.47

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, IU/mL	CV, %
1	12.5	8.36
2	113.4	1.47

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, IU/mL	Concentration2, IU/mL	Concentration3, IU/mL	CV, %
1	12.7	13.3	12.3	3.66
2	115.5	117.8	115.1	1.25

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known total IgE concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 50-1000 IU/mL $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest total IgE concentration in the serum or plasma sample that is detected by the Total IgE EIA kit is no lower than 3 IU/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for Total IgE EIA kit is 50IU/mL.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of total IgE with other analytes is shown in the table:

Analyte	Concentration, IU/mL	Cross-reactivity, %
IgA	1000	Not detected
IgM	1000	Not detected
IgG	1000	Not detected

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC

Akademika Yefremova St. 23

03179, Kyiv, Ukraine

tel.:+38 044 422-62-16

tel.:+38 044 294-69-78

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www.xema.in.ua