FOR INFORMATION ONLY.
WHEN PERFORMING
THE ASSAY ALWAYS REFER
TO PACKAGE INSERT
SUPPLIED
WITH THE KIT

RO

CanAg CA242 EIA

REF

101-10

IVD

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EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIFICAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Könlilik kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunaudoti iki/Izletot līdz/Houdbaar
tot/Brukes innen/Użyć przed/
Prazo de validade/Expiră la/
Использовать до/Använd före/
Použite ně do/ Uporabno do/
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Πατίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzszzám/Codice del lotto/
Partijos kodus/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Datuma na proizvodstvo/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de
fabricación/Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcji/Data de fabrico/Data
fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum
izdelave/Datum proizvodnje/Üretim tarihi

IVD

In Vitro Diagnostic Medical Device/
Μεδικινικις υρειδ ζα διαγνστικα
ιν ιντρο/Diagnostický zdravotnícký
prostředek in vitro/Medicinsk udstyr til
in vitro-diagnostik/In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para
diagnóstico in vitro/In vitro diagnostiline
meditsiiniseade/Dispositif médical de
diagnostic in vitro/Diagnostički medicinski
uređaj In Vitro/In vitro orvosdiagnostikai
eszköz/Dispositivo medico per test
diagnostici in vitro/In Vitro Diagnostinė
Medicinos Priemonė/Mediciniska ierice
in vitro diagnostikai/In vitro-diagnostisch
medisch instrument/In vitro diagnostisk
medisinsk utstyr/Wyrób medyczny do
diagnostyki in vitro/Dispositivo Médico
de Diagnóstico In Vitro/Dispozitiv medical
pentru diagnostic in vitro/Только для
диагностики In Vitro/Endast för in
vitro-diagnostik/ Zdravotnicka pomůcka na
diagnostiku in vitro/In vitro diagnostični
pripomoček/Diagnostički medicinski
uređaj In Vitro/<96> testleri için yeterlilik
içerir

REF

Catalogue number/Καταλογζην νομερ/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Katalogi number/
Número de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Nomer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalógové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Temperature limitation/
Τεμπερατυρνι γρانيци/
Teplotní omezení/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/
Temperatūriniai apribojimai/
Temperatūras ierobežojums/
Temperatuurbeperking/
Temperaturbegrensninger/
Temperatury graniczne/
Limite de temperatura/
Limite de temperatură/
Τεμπερατυρньй режим/
Temperaturbegrænsning/
Teplotné obmedzenie
Omejitve temperature/
Temperaturno ograničenje/
Sicaklik sinirlamasi/

Contains sufficient for <96> tests/Съдържа
достатъчно количество за тестове
<96>/Lze použít pro <96> testů/Ineholder
tilstrækkeligt/Inhalt ausreichend für <96>
Prüfungen/Περιεχόμενο επαρκές για
«96» εξετάσεις/Contenido suficiente para
<96> ensayos/Kogusest piisab <96> testi
läbiviimiseks/Contenu suffisant pour «96»
tests/Sadržaji dovoljno za <96> testova/A
doboz tartalma <96> vizsgálat elvégzéséhez
elegendő/Contenuto sufficiente per «96»
saggi/Turins skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/Inhoud
voldoende voor «96» testen/til «96» test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para «96» ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для
«96» определений/Innehåller tillräckligt
till «96» antal tester/Obsah postačuje na
tento počet testov: <96>/Vsebinsa zadostuje
za <96> testov/Sadržina dovoljna za <96>
testova/<96> testleri için yeterlilik içerir

Consult Instructions for Use/
Прочетете инструкцията за
употреба/Konzultujte s návodem
k použití/Se brugsanvisning/Siehe
Gebrauchsanweisung/Συμβουλευτείτε
τις Οδηγίες σχετικά με τη χρήση/
Consulte las instrucciones de uso/
Vt kasutusjuhendit/Consulter le mode
d'emploi/Pročitajte upute za uporabu/
Olvassa el a használati utasítást/
Consultare le istruzioni per l'uso/Dél
naudojimo žiūrėkite instrukcijas/Izlasiet
lietošanas instrukciju/Raadpleeg de
instructies voor gebruik/Les instruksene
før bruk/Sprawdzić w instrukcji użycia/
Consulte as Instruções de Utilização/
Consultați instrucțiunile de utilizare/
Обратитесь к инструкции по
применению/Se bruksanvisning/
Prečítajte si návod na používanie/
Pročitajte uputstvo za upotrebu/
Kullanım Talimatlarını Bakınız

CONT

Contents of kit/Съдържание на набора/
Obsah soupravy/Kittets indhold/Inhalt
des Kits/Περιεχόμενα του kit/Contenido
del kit/Komplekt sisaladab/Contenu du
kit/Sadržaj opreme/A készlet tartalma/
Contenuto del kit/Rinkinio turinys/
Komplekta saturs/Inhoud van de set/
Settets innhold/Zawartość zestawu/
Conteúdo do kit/Conținutul setului/
Компоненты набора/Kit innehåll/
Obsah súpravy/Vsebina kompleta/Sadržaj
opreme/Kitin içindekiler

Biological risks/Биологическа
опасност/Biologická rizika/Biologisk
fare/Biologische Gefahren/Biológicoi
kivõnnoi/Riesgos biológicos/
Bioloogilised ohud/Risques biologiques/
Biološki rizici/Biológiai kockázatok/Rischi
biologici/Biologinis pavojus/Biologiskais
risks/Biologische risico's/Biologische
risikoer/Zagroženie biologické/Riscos
biológicos/ Biologisk risk/Pericole
biologice/Биологическая опасность/
Biologický rizikové/Biologické riziká/
Biološki rizici/Bijolojik riskler

ORIG HUM

Human/C човешки произход/Lidské/
Human/Human/ἄνθρωπος ἀναφορᾶς/
Humano/Inimăpăritolu/Humaine/Ljudskog
porjekla/Humán/Origine Umana/
Žmogaus kilmės/Cilvēku izcelsmes/
Human/Menneske/Ludzka/Humano/
Origine umană/Человеческого
происхождения/Human/Ludské/
Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/C миши произход/Мыši/
Fra mus/Maus/από ποντίκι/de ratón/
Hiirtelt/De souris/Mišijeg porjekla/
Egérbdli/Murino/Pelės kilmės/No peles/
Van muizen/Fra mus/Mysia/Do rato/De
la șoareci/Мышиного происхождения/
Från mus/Myšie/Mišjega izvora/Mišijeg
porekla/Farenden

ORIG BOV

Bovine/C говежди произход/
Hovězí/Bovinn/Rind/από βοοειδή/
Bovino/Veistelt/Bovine/Rogate stoke/
Szarvasmarha/Bovino/Jaučio/No
liellopa/Bovien/Bovinn/Wolowy/Bovino/
Origine bovină/крупного рогатого
скота/Från ko/Hovädzie/Govejega
izvora/Rogate krupne stoke/Bovinn

Reconstitute with/Пазтваряне с/
Rozfeďte pomocí/Rekonstituere med/
Rekonstituieren mit/Ανασύσταση με/
Reconstituer con/Lahjendamine/
Reconstituer avec/Rekonstituiraite s/
Feloldáshoz/Ricostituire con/Atkurti,
ištirpdant su/Atšķaidīt ar/Reconstitutie
met/Rekonstituere med/Odtworzyć
za pomocą/Reconstituer com/A
se reconstitui cu/Растворить в/
Rekonstituera med/Rozriedzte pomocou/
Rekonstituiraite z/s/Ponovno formiranje
sa/Yeniden oluşturulur

Manufacturer/Производител/Výrobce/
Producent/Hersteller/Κατασκευαστής/
Fabricante/Tootja/Fabricant/Proizvođač/
Gyártó/Fabbrikante/Gamintojas/
Ražotājs/Fabrikant/Produsent/
Producent/Fabricante/Producător/
Производител/Tilverkare/ Výrobca/
Izdelovalec/Proizvođač/Üretici

INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE

Please visit our website **www.fdi.com/ifu** to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website **www.fdi.com**. If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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- V, Roselli M, Guadagni F. (2001) Comparative analysis of CA 242 and CA 19-9 serum tumor markers in colorectal cancer patients. A longitudinal evaluation. *Anticancer Res.* Mar-Apr;21(2B):1263-70.
8. Carpelan-Holmström, M., Haglund, K., Lundin, J., Järvinen, H and Roberts, P., (1996) Pre-operative serum levels of CA 242 and CEA predict outcome in colorectal cancer. *Eur. J. Cancer* 32(7), 1156-1161.
 9. Engarås B. (2003) Individual cutoff levels of carcinoembryonic antigen and CA 242 indicate recurrence of colorectal cancer with high sensitivity. *Dis Colon Rectum.* Mar;46(3):313-21.
 10. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
 11. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
 12. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
 13. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
 14. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

CanAg[®] is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB

Elof Lindälv's gata 13

SE-414 58 Göteborg

Sweden

Phone + 46 31 85 70 30

Fax + 46 31 85 70 40

info@fdab.com

www.fdab.com