Product List – CE Marked

Certified by

ISO 13485:2016

EC – Directive 98 / 79 EC
For In-Vitro-Diagnostics

2020-02-1

NovaLisa[®] Virology

Prod. No.	Name
ADVA0010	Adenovirus IgA
ADVG0010	Adenovirus IgG
ADVM0010	Adenovirus IgM
CHIG0590	Chikungunya Virus IgG capture
CHIM0590	Chikungunya Virus IgM μ -capture
CMVG0110	Cytomegalovirus (CMV) IgG
ACMV7110	Avidity Cytomegalovirus (CMV) IgG
CMVM0110	Cytomegalovirus (CMV) IgM
DENG0120	Dengue Virus IgG
DENM0120	Dengue Virus IgM
DVM0640	Dengue Virus IgM μ -capture
NS1D4020	Dengue Virus NS1 Antigen
EBVA0150	Epstein-Barr Virus (VCA) IgA
EBVG0150	Epstein-Barr Virus (VCA) IgG
AEBV7150	Avidity Epstein-Barr Virus (VCA) IgG
EBVM0150	Epstein-Barr Virus (VCA) IgM
EBVG0580	Epstein-Barr Virus (EBNA) IgG
HANG0670	Hantavirus IgG
HANM0670	Hantavirus IgM
HEVG0780	Hepatitis E Virus (HEV) IgG
HEVM0780	Hepatitis E Virus (HEV) IgM
HSVG0250	Herpes simplex Virus 1+2 (HSV) IgG
HSVM0250	Herpes simplex Virus 1+2 (HSV) IgM
HSV1G0500	Herpes simplex Virus 1 (HSV 1) IgG
HSV1M0500	Herpes simplex Virus 1 (HSV 1) IgM
HSV2G0540	Herpes simplex Virus 2 (HSV 2) IgG
HSV2M0540	Herpes simplex Virus 2 (HSV 2) IgM
INFA0290	Influenza Virus A IgA
INFG0290	Influenza Virus A IgG
INFM0290	Influenza Virus A IgM
INFA0300	Influenza Virus B IgA
INFG0300	Influenza Virus B IgG
INFM0300	Influenza Virus B IgM
MEAG0330	Measles Virus IgG
AMEA7330	Avidity Measles Virus IgG
MEAM0330	Measles Virus IgM
MUMG0340	Mumps Virus IgG
MUMM0340	Mumps Virus IgM
PAIA0360	Parainfluenza Virus 1,2,3 IgA
PAIG0360	Parainfluenza Virus 1,2,3 IgG
PARG0370	Parvovirus B 19 IgG
PARM0370	Parvovirus B 19 IgM
RSVA0380	Respiratory syncytial Virus IgA
RSVG0380	Respiratory syncytial Virus IgG
RSVM0380	Respiratory syncytial Virus IgM
RUBG0400	Rubella Virus IgG

ARUB7400	Avidity Rubella Virus IgG
RUBM0400	Rubella Virus IgM μ -capture
TICG0440	TBE / FSME IgG
TICM0440	TBE / FSME IgM
PTICG044	TBE / FSME IgG plus
VZVA0490	Varicella-Zoster Virus (VZV) IgA
VZVG0490	Varicella-Zoster Virus (VZV) IgG
VZVM0490	Varicella-Zoster Virus (VZV) IgM
ZVG0790	Zika Virus IgG capture
ZVM0790	Zika Virus IgM μ -capture

NovaLisa[®] Bacteriology

Prod. No.	Name
BAR0900	Bartonella
BOPA0030	Bordetella pertussis IgA
BOPG0030	Bordetella pertussis IgG
BOPM0030	Bordetella pertussis IgM
BPTA0610	Bordetella pertussis toxin (PT) IgA
BPTG0610	Bordetella pertussis toxin (PT) IgG
BORG0040	Borrelia burgdorferi IgG
BORM0040	Borrelia burgdorferi IgM
BRUG0050	Brucella IgG
BRUM0050	Brucella IgM
CHLA0070	Chlamydia trachomatis IgA
CHLG0070	Chlamydia trachomatis IgG
CHLM0070	Chlamydia trachomatis IgM
CHLA0510	Chlamydia pneumoniae IgA
CHLG0510	Chlamydia pneumoniae IgG
CHLM0510	Chlamydia pneumoniae IgM
CORG0090	Corynebacterium diphtheriae toxin IgG
CORG5009	Corynebacterium diphtheriae toxin 5S IgG
PCORG009	Corynebacterium diphtheriae toxin 5S IgG plus
COX1G0600	Coxiella burnetii (Q-Fever) Phase 1 IgG
COX2G0600	Coxiella burnetii (Q-Fever) Phase 2 IgG
COX2M0600	Coxiella burnetii (Q-Fever) Phase 2 IgM
HELA0220	Helicobacter pylori IgA
HELG0220	Helicobacter pylori IgG
PHELA022	Helicobacter pylori IgA plus
PELG022	Helicobacter pylori IgG plus
LEGG0650	Legionella Pneumophila IgG
LEGM0650	Legionella Pneumophila IgM
LEPG0660	Leptospira IgG
LEPM0660	Leptospira IgM

MYCA0350	Mycoplasma pneumoniae IgA
MYCG0350	Mycoplasma pneumoniae IgG
MYCM0350	Mycoplasma pneumoniae IgM
TETG0430	Clostridium tetani toxin IgG
TETG5043	Clostridium tetani toxin 5S IgG
PTETG043	Clostridium tetani toxin 5S IgG plus

NovaLisa® Parasites

Prod. No.	Name
CHAG0560	Chagas (Trypanosoma cruzi) IgG
TRYP0570	Chagas
ENTG0140	Entamoeba histolytica IgG
LEIG0310	Leishmania infantum IgG
MAL0620	Malaria
TOXA0460	Toxoplasma gondii IgA
TOXG0460	Toxoplasma gondii IgG
ATOX7460	Avidity Toxoplasma gondii IgG
TOXM0460	Toxoplasma gondii IgM μ -capture

NovaLisa® Worms

Prod. No.	Name
ASCG0020	Ascaris lumbricoides IgG
ECHG0130	Echinococcus IgG
FIL0760	Filariasis
SCHG0410	Schistosoma mansoni IgG
SCHM0410	Schistosoma mansoni IgM
STRO0690	Strongyloides
TAEG0420	Taenia solium IgG
TOCG0450	Toxocara canis IgG
TRIG0480	Trichinella spiralis IgG

NovaLisa® Fungi

Prod. No.	Name
ASPG0680	Aspergillus fumigatus IgG
ASPM0680	Aspergillus fumigatus IgM
CANA0060	Candida albicans IgA
CANG0060	Candida albicans IgG
CANM0060	Candida albicans IgM

NovaLisa® Hormones

THYROID HORMONES

(ELISAs for the determination of thyroid hormones and antibodies)

Prod. No.	Name
ATG1010	Anti-TG
ATPO1020	Anti-TPO
FT41050	Free T4
TSH1030	TSH

Hormones

STEROID HORMONES

(ELISAs for the determination of steroid hormones in plasma and serum)

Prod. No.	Name
DNOV001	Cortisol
DNOV002	Testosterone
DNOV003	17 beta-Estradiol
DNOV004	17-OH Progesterone
DNOV005	DHEA-S
DNOV006	Progesterone
DNOV008	Androstenedione
DNOV009	Free Testosterone
DNOV011	Total Estriol
DNOV012	Aldosterone

STEROID HORMONES IN URINE

(ELISAs for the determination of steroid hormones in urine)

Prod. No.	Name
DNOV010	Urinary Cortisol

STEROID HORMONES IN SALIVA

(ELISAs for the determination of steroid hormones in saliva)

Prod. No.	Name
DSNOV20	Cortisol Saliva
DSNOV21	Testosterone Saliva
DSNOV24	DHEA-S Saliva
DSNOV27	Androstenedione Saliva

PROTEIN HORMONES

(ELISAs for the determination of proteins in plasma and serum)

Prod. No. Name

DNOV030	LH
DNOV031	FSH
DNOV032	Prolactin
DNOV033	AFP
DNOV034	beta HCG

THYROID HORMONES

(ELISAs for the determination of thyroid hormones and antibodies)

Prod. No. Name

DNOV051	Free T3
DNOV053	Total T3
DNOV054	Total T4
DNOV057	Thyroglobulin

DIABETES MONITORING

(ELISAs for the determination of specific analytes in plasma and serum)

Prod. No. Name

DNOV111	Insulin
DNOV112	C-Peptide

CIRCULATING IMMUNO COMPLEXES

(ELISAs for the determination of specific analytes in plasma and serum)

Prod. No. Name

DNOV093	CIC-C1q
DNOV094	CIC-C3d
DNOV096	CH-50

TUMOR MARKERS

(ELISAs for the determination of specific analytes in plasma and serum)

Prod. No. Name

DNOV 060	CEA
DNOV061	CA 125
DNOV062	CA 15-3
DNOV063	CA 19-9

MISCELLANEOUS

(ELISAs for the determination of specific analytes in plasma and serum)

Prod. No.	Name
------------------	-------------

DNOV100	Ferritin
DNOV101	HGH
DNOV102	IgE

NovoLisa[®] Autoimmune

Autoimmune

(ELISAs for the determination of specific autoimmune antibodies)

Prod. No.	Name
------------------	-------------

ATG1010	Anti-TG
ATPO1020	Anti-TPO

Rheumatology

(ELISAs for the determination of specific analytes in plasma and serum)

Prod. No.	Name
------------------	-------------

RFM3010	Rheumatoid Factor IgM
---------	-----------------------

NovoLisa[®] Recombinant Antigens

Prod. No.	Name
------------------	-------------

BORG0040	Borrelia burgdorferi IgG
BORM0040	Borrelia burgdorferi IgM
CHAG0560	Chagas (Trypanosoma cruzi) IgG
TRYP0570	Chagas
HANG0670	Hantavirus IgG
HANM0670	Hantavirus IgM
HELA0220	Helicobacter pylori IgA
PHELA022	Helicobacter pylori IgA plus
HEVG0780	Hepatitis E Virus (HEV) IgG
HEVM0780	Hepatitis E Virus (HEV) IgM
HSV1G0500	Herpes simplex Virus 1 (HSV 1) IgG
HSV1M0500	Herpes simplex Virus 1 (HSV 1) IgM
HSV2G0540	Herpes simplex Virus 2 (HSV 2) IgG
HSV2M0540	Herpes simplex Virus 2 (HSV 2) IgM
MAL0620	Malaria
STRO0690	Strongyloides
ZVG0790	Zika Virus IgG capture
ZVM0790	Zika Virus IgM μ -capture

NovaLisa[®] Quantitative Assays (WHO standardized)

Prod. No.	Name
BPTA0610	Bordetella pertussis toxin (PT) IgA
BPTG0610	Bordetella pertussis toxin (PT) IgG
CORG0090	Corynebacterium diphtheriae toxin IgG
CORG5009	Corynebacterium diphtheriae toxin 5S IgG
PCORG009	Corynebacterium diphtheriae toxin 5S IgG plus
RFM3010	Rheumatoid Factor IgM
RUBG0400	Rubella Virus IgG
TETG0430	Clostridium tetani toxin IgG
TETG5043	Clostridium tetani toxin 5S IgG
PTETG043	Clostridium tetani toxin 5S IgG plus
TOXG0460	Toxoplasma gondii IgG
ATOX7460	Avidity Toxoplasma gondii IgG
TSH1030	TSH

NovaLisa[®] Quantitative Assays

Prod. No.	Name
ATG1010	Anti-TG
ATPO1020	Anti-TPO
BPTA0610	Bordetella pertussis toxin (PT) IgA
BPTG0610	Bordetella pertussis toxin (PT) IgG
CORG0090	Corynebacterium diphtheriae toxin IgG
CORG5009	Corynebacterium diphtheriae toxin 5S IgG
PCORG009	Corynebacterium diphtheriae toxin 5S IgG plus
FT41050	Free T4
HELA0220	Helicobacter pylori IgA
HELG0220	Helicobacter pylori IgG
PHELA022	Helicobacter pylori IgA plus
PHELG022	Helicobacter pylori IgG plus
RFM3010	Rheumatoid Factor IgM
RUBG0400	Rubella Virus IgG
ARUB7400	Avidity Rubella Virus IgG
TETG0430	Clostridium tetani toxin IgG
TETG5043	Clostridium tetani 5S toxin IgG
PTETG043	Clostridium tetani toxin 5S IgG plus
TICG0440	TBE / FSME IgG
PTICG044	TBE / FSME IgG plus
TOXG0460	Toxoplasma gondii IgG
ATOX7460	Avidity Toxoplasma gondii IgG
TSH1030	TSH

Antigen Assays

Prod. No.	Name
NS1D4020	Dengue Virus NS1 Antigen

NovaLisa[®] IgM μ -capture Assays

Prod. No.	Name
CHIM0590	Chikungunya Virus IgM μ -capture
DVM0640	Dengue Virus IgM μ -capture
RUBM0400	Rubella Virus IgM μ -capture
TOXM0460	Toxoplasma gondii IgM μ -capture
ZVM0790	Zika Virus IgM μ -capture

NovaLisa[®] Antibody Assays

Prod. No.	Name
ASCG0020	Ascaris lumbricoides IgG
CHAG0560	Chagas (Trypanosoma cruzi) IgG
TRYP0570	Chagas
ENTG0140	Entamoeba histolytica IgG
LEIG0310	Leishmania infantum IgG
MAL0620	Malaria
STRO0690	Strongyloides
TAEG0420	Taenia solium IgG
TOCG0450	Toxocara canis IgG
TRIG0480	Trichinella spiralis IgG

NovaLisa[®] Avidity Assays

Prod. No.	Name
ACMV7110	Avidity Cytomegalovirus (CMV) IgG
AEBV7150	Avidity Epstein-Barr Virus (VCA) IgG
AMEA7330	Avidity Measles Virus IgG
ARUB7400	Avidity Rubella Virus IgG
ATOX7460	Avidity Toxoplasma gondii IgG

NovaLisa[®] Liquor Diagnostic

Prod. No.	Name
BORG0040	Borrelia burgdorferi IgG
BORM0040	Borrelia burgdorferi IgM

NovaLisa®

Bordetella pertussis IgM

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português.....	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía/ Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: BOPM0030 (96 Determinations)

ENGLISH

1. INTRODUCTION

Bordetella pertussis is a respiratory pathogen that causes pertussis, commonly known as whooping cough, a localized infection of the ciliated epithelium of the bronchial tree. Pertussis is characterized by a prolonged paroxysmal cough often accompanied by an inspiratory whoop.

The disease affects mainly children, but adults have also been increasingly reported to be affected. The pathogen produces toxins which cause local damage to the cilia of epithelial cells, which leads to prolonged illness and pertussis. Disease presentation varies with age and history of previous exposure or vaccination. Severe disease is infrequent in healthy, vaccinated persons. Infants, particularly those who have not received the primary vaccination series against pertussis, are at risk for complications and mortality.

In addition to *B. pertussis*, three other *Bordetella* species can cause disease in humans: *B. parapertussis*, *B. holmesii*, and *B. bronchiseptica*. *B. parapertussis* causes a pertussis-like illness that is generally milder than pertussis because the bacteria do not produce pertussis toxin. Co-infection of *B. pertussis* and *B. parapertussis* is not unusual.

B. pertussis is of worldwide prevalence. Globally, 20-40 million cases of pertussis occur each year, 90 % of which are in developing countries, and there are up to 400,000 fatalities each year, mostly in young infants.

Transmission of *B. pertussis* occurs primarily via close direct contact with an infected person or inhalation of airborne droplets. Symptoms develop following inhalation of the airborne pathogen. The organism is highly contagious, with up to 90 % of household contacts developing the disease. Infected persons are most contagious in the catarrhal and the paroxysmal stages.

The incubation period is usually seven to 10 days, with a range of 4-21 days.

Typical pertussis symptoms occur in three different stages: catarrhal, paroxysmal, and convalescent.

The catarrhal stage lasts for about 1-2 weeks, and is characterized by non-specific symptoms such as rhinorrhea, sneezing, low-grade fever and cough. The second stage is the paroxysmal stage, lasting for about 4-6 weeks, and is characterized by various pathognomonic symptoms of pertussis such as episodes of paroxysmal cough with a characteristic whooping sound. The final stage is the convalescent stage. During this stage, the respiratory symptoms gradually decrease although cough can continue for months.

Many factors can alter the usual course of pertussis, causing an atypical presentation. Previously vaccinated adolescents and adults may have less severe paroxysmal symptoms.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Bordetella pertussis</i>	Pertussis whooping cough	<u>1. Stadium catarrhale</u> : symptoms of a cold with slight fever (1-2 weeks) <u>2. Stadium convulsivum</u> : severe, spasmodic coughing; after deep inspiration follows a coughing staccato (2-6 weeks) <u>3. Stadium decrementi</u> : Ease of disease with symptoms of a bronchitis (up to 6 weeks)	Highly contagious droplet infection

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- Microscopy
- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The *Bordetella pertussis* IgM ELISA is intended for the qualitative determination of IgM class antibodies against *Bordetella pertussis* in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Bordetella pertussis* antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Bordetella pertussis* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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8. ASSAY PROCEDURE

Please read the instruction for use carefully before performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 μ L to 350 μ L to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

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Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100 % (95% confidence interval: 96.19% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 89.19 % (95% confidence interval: 74.58% - 96.97%).

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Product Number: BOPG0030 (96 Determinations)

ENGLISH

1. INTRODUCTION

Bordetella pertussis is a respiratory pathogen that causes pertussis, commonly known as whooping cough, a localized infection of the ciliated epithelium of the bronchial tree. Pertussis is characterized by a prolonged paroxysmal cough often accompanied by an inspiratory whoop.

The disease affects mainly children, but adults have also been increasingly reported to be affected. The pathogen produces toxins which cause local damage to the cilia of epithelial cells, which leads to prolonged illness and pertussis. Disease presentation varies with age and history of previous exposure or vaccination. Severe disease is infrequent in healthy, vaccinated persons. Infants, particularly those who have not received the primary vaccination series against pertussis, are at risk for complications and mortality.

In addition to *B. pertussis*, three other *Bordetella* species can cause disease in humans: *B. parapertussis*, *B. holmesii*, and *B. bronchiseptica*. *B. parapertussis* causes a pertussis-like illness that is generally milder than pertussis because the bacteria do not produce pertussis toxin. Co-infection of *B. pertussis* and *B. parapertussis* is not unusual.

B. pertussis is of worldwide prevalence. Globally, 20-40 million cases of pertussis occur each year, 90 % of which are in developing countries, and there are up to 400,000 fatalities each year, mostly in young infants.

Transmission of *B. pertussis* occurs primarily via close direct contact with an infected person or inhalation of airborne droplets. Symptoms develop following inhalation of the airborne pathogen. The organism is highly contagious, with up to 90 % of household contacts developing the disease. Infected persons are most contagious in the catarrhal and the paroxysmal stages.

The incubation period is usually seven to 10 days, with a range of 4-21 days.

Typical pertussis symptoms occur in three different stages: catarrhal, paroxysmal, and convalescent.

The catarrhal stage lasts for about 1-2 weeks, and is characterized by non-specific symptoms such as rhinorrhea, sneezing, low-grade fever and cough. The second stage is the paroxysmal stage, lasting for about 4-6 weeks, and is characterized by various pathognomonic symptoms of pertussis such as episodes of paroxysmal cough with a characteristic whooping sound. The final stage is the convalescent stage. During this stage, the respiratory symptoms gradually decrease although cough can continue for months.

Many factors can alter the usual course of pertussis, causing an atypical presentation. Previously vaccinated adolescents and adults may have less severe paroxysmal symptoms.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Bordetella pertussis</i>	Pertussis whooping cough	<u>1. Stadium catarrhale</u> : symptoms of a cold with slight fever (1-2 weeks) <u>2. Stadium convulsivum</u> : severe, spasmodic coughing; after deep inspiration follows a coughing staccato (2-6 weeks) <u>3. Stadium decrementi</u> : Ease of disease with symptoms of a bronchitis (up to 6 weeks)	Highly contagious droplet infection

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- Microscopy
- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The *Bordetella pertussis* IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Bordetella pertussis* in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Bordetella pertussis* antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Bordetella pertussis* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

It is 93.02% (95% confidence interval: 80.94% - 98.54%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

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Product Number: CORG0090 (96 Determinations)

ENGLISH

1. INTRODUCTION

Corynebacteria are aerobic non spore-forming gram-positive rods of irregular shape (0.5 –1 µm thick and 2-6 µm long). They comprise skin commensals, opportunist pathogens and several major pathogens, including *Corynebacterium diphtheriae*. In general, they are isolated from throat swabs on selective media containing tellurite. The bacterial infection caused by *C. diphtheriae*, Diphtheria, has two forms. Respiratory diphtheria is typically caused by toxin-producing (toxigenic) strains; cutaneous disease can be caused by either toxigenic or nontoxigenic strains. In the respiratory form of the disease, a membrane is formed; this membrane is usually visible on the throat or tonsils. Persons may die from asphyxiation when the membrane obstructs breathing. Other complications are caused by remote effects of the diphtheria toxin (myocarditis, nerve paralysis) Cutaneous diphtheria is usually mild, typically consisting of non-distinctive sores or shallow ulcers and only rarely involving toxic complications (1-2% of infections with toxigenic strains). Diphtheria was one of the most common causes of death among children during the prevaccine era.

Since the introduction and widespread use of diphtheria toxoid vaccine (formalin-inactivated diphtheria toxin) in most industrialized countries the disease is now characterized by sporadic cases and intermittent outbreaks of low intensity. But recent large epidemics of diphtheria in several eastern European countries have again drawn attention to this „forgotten“ disease – and, the majority of these cases have occurred among adolescents and adults instead of children.

The only effective way to control diphtheria is by prophylactic immunization with diphtheria toxoid. Antibody to the toxoid protects against the action of the toxin; immunized persons can be infected by toxin-producing strains of diphtheria, but the systemic manifestations of diphtheria do not occur. The outcome of the disease improves with early, appropriate treatment. Prompt recognition of the disease is important to assure early, appropriate treatment with diphtheria anti-toxin

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Corynebacterium diphtheriae</i>	Diphtheria (respiratory)	with malaise, sore throat, anorexia, low-grade fever and swelling of the neck (“bull neck”) from inflammation. Complications: exotoxin-induced damage to other organs.	Transmission from person to person through close physical and respiratory contact Transmission is increased in overcrowded and poor socio-economic conditions

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- Microscopy
- Serology: e.g. by ELISA

2. INTENDED USE

The *Corynebacterium diphtheriae* toxin IgG ELISA is intended for the quantitative determination of IgG class antibodies against *Corynebacterium diphtheriae* toxin in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break apart 8-well snap-off strips coated with *Corynebacterium diphtheriae* toxin (toxoid) antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0,1 %; ready to use; yellow cap.
- **Standards:** 4 vials, each containing 2 mL standard; coloured yellow; ready to use; ≤ 0.02% (v/v) MIT.
Standard A: 0.000 IU/mL; blue cap
Standard B: 0.015 IU/mL; green cap
Standard C: 0.075 IU/mL; yellow cap
Standard D: 0.150 IU/mL; red cap
The standards are calibrated in accordance with the "1st International Standard for Diphtheria Antitoxin Human IgG (WHO, 2012).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Corynebacterium diphtheriae* toxin (toxoid) antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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8. ASSAY PROCEDURE

Please read the instruction for use carefully before performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature(20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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
9.2. Calculation of Results

In order to obtain **quantitative results in IU/mL** plot the (mean) absorbance values of the 4 Standards A, B, C and D on (linear/linear) graph paper in a system of coordinates against their corresponding concentrations (0.000, 0.015, 0.075, 0.150 IU/mL) and draw a standard calibration curve (absorbance values on the y-axis, concentrations on the x-axis).

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9.4. Interpretation of Results

according to: RKI 1999

< 0.01 IU/mL	No protective antibody level! Immediate full course of basic immunization is recommended!
0.01 - 0.09 IU/mL	No reliable protection! Immediate booster injection is recommended.
0.1 – 1.0 IU/mL	Reliable protection!
> 1.0 IU/mL	Reliable long term protection: After about 10 years after last booster control and booster injection is recommended. It is recommended that the basic immunisation or booster is checked 4-6 weeks after immunisation and to record the data on the certificate of vaccination.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 89.42% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 95.44% - 100%).

10.4. Analytical Sensitivity

The analytical sensitivity (according to CLSI EP17-A) is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator. It is 0.00092 IU/mL.

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Reagents may contain CMIT/MIT (3:1) or MIT (refer to 4.1)
Therefore, the following hazard and precautionary statements apply.

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ENGLISH

1. INTRODUCTION

Echinococci are microscopic cestodes (tapeworms) of 1-6 mm which are dependent on their genus found either in dogs or other canids (*E. granulosus*) or in foxes, coyotes and wolves (*E. multilocularis*). In their larval stage they are the causative agent of human echinococcosis (Hydatidosis, or hydatid disease). The adult tapeworms reside in the small bowel of the definitive hosts, and gravid proglottids release eggs that are passed in the feces. After ingestion of a suitable intermediate host, the egg hatches in the small bowel and releases an oncosphere that penetrates the intestinal wall and through the circulatory system into various organs, especially the liver and lungs, where it develops into a cyst. Echinococcus infections remain silent for years before the enlarging cysts cause symptoms in the affected organs (liver, lung, and less commonly other organs as brain, bone, heart). *E. granulosus* occurs practically worldwide, and more frequently in rural, grazing areas where dogs ingest organs from infected animals. *E. multilocularis* occurs in the northern hemisphere, including central Europe and the northern parts of Europe, Asia, and North America. Although human cases are rare, infection in humans causes parasitic tumors to form in the liver, the lungs, and less commonly, the brain, and other organs. If left untreated, infection can be fatal.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>E. granulosus</i>	Cystic echinococcosis (Cystic Hydatid Disease, CHD)	(Depends on localization size, and number of cysts) Liver: Upper abdominal pain, hepatomegaly, cholestasis, jaundice, etc.	"hand-to-mouth" transmission. Infection by oral uptake of eggs.
<i>E. multilocularis</i>	Alveolar Echinococcosis (Alveolar Hydatid Disease, AHD)	Lungs: Thoracic pains, cough, expectoration, dyspnea, etc. CNS: Neurological symptoms	e.g.: contaminated wild berries.

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- Microscopy
- Serology: e.g by ELISA

2. INTENDED USE

The Echinococcus IgG ELISA is intended for the qualitative determination of IgG class antibodies against Echinococcus in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with Echinococcus antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap;
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with Echinococcus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

It is 98.82% (95% confidence interval: 95.81% - 99.86%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

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	Warning	H317	May cause an allergic skin reaction.
		P261	Avoid breathing spray.
		P280	Wear protective gloves/ protective clothing.
		P302+P352	IF ON SKIN: Wash with plenty of soap and water.
		P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
		P362+P364	Take off contaminated and Wash it before reuse.

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Product Number: MEAM0330 (96 Determinations)

ENGLISH

1. INTRODUCTION

Measles or morbilli virus belongs to the RNA viruses of the family Paramyxoviridae. The virions are spherical particles of 150-250 nm in diameter consisting of the ribonucleoprotein with helical symmetry and an envelope with spikes containing the strain-specific and hemagglutinating antigens. Morbilli viruses have no neuraminidase activity. Measles is a classic childhood disease. The virus is endemic: at the age of 20 about 90% of the population has had immunological experience with it. Newborns are protected by maternal antibodies for the first 3-4 months of life; the active disease leaves lifelong immunity. The measles virus has a contagiousness index of about 96%, is worldwide distributed, and can be serious. Bacterial superinfection was a serious threat in the preantibiotic era, but the prognosis of uncomplicated measles is now good. CNS complications such as encephalomyelitis (0.1%) which may occur after the acute phase of measles infection subsides, however still have a high mortality (10%). Prognosis of recovery in these patients is poor. Between 10-30% of all cases are fatal; 20-50% develop significant damages. Subacute sclerosing panencephalitis (SSPE) is a rare (1:1000) degenerative disease of the CNS which is thought to be a slow virus infection.

Species	Disease	Symptoms (e.g.)	Transmission route
Measles Virus	Measles	Fever, malaise, productive cough and runny nose, headache, abdominal pain, typical lesions in the mouth (Koplik spots), characteristic exanthema	By air or contact with saliva or nasal secretions
	SSPE (Subacute sclerosing panencephalitis)	Complications: Encephalomyelitis, decrease in intellectual skills that progresses to an almost complete loss of brain function and death	

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- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The Measles Virus IgM ELISA is intended for the qualitative determination of IgM class antibodies against Measles Virus in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate (IgM):** 12 break-apart 8-well snap-off strips coated with Measles Virus antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with Measles Virus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 98.71% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 91.19% - 100%).

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13. ORDERING INFORMATION

Prod. No.: MEAM0330 Measles Virus IgM ELISA (96 Determinations)

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Product Number: MEAG0330 (96 Determinations)

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- PCR
- Serology: e.g. by ELISA

2. INTENDED USE

The Measles Virus IgG ELISA is intended for the qualitative determination of IgG class antibodies against Measles Virus in human serum or plasma (citrate, heparin).

Measles Virus IgG avidity can be determined with assay Avidity Measles Virus IgG (Product code: AMEA7330).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with Measles Virus antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

	NTU	mIU/mL (3 rd International Standard)	
Cut-off	10	-	-
Positive	> 11	> 220	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11	120 – 220	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9	< 120	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 90.0% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 97.01% (95% confidence interval: 93.93% - 98.79%).

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13. ORDERING INFORMATION

Prod. No.: MEAG0330 Measles Virus IgG ELISA (96 Determinations)

For avidity testing:

AMEA7330 Avidity Measles Virus IgG ELISA (48 Determinations)

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Product Number: AMEA7330 (48 Determinations)

ENGLISH

1. INTRODUCTION

The presence of IgG antibodies to Measles virus indicates the occurrence of the infection but does not distinguish between recent and past infection. Specific IgM antibodies are first detected approximately in ten days and peak at about four weeks post infection. They may persist for several months after acute infections. Based on the evidence that antibody avidity gradually increases after exposure to an immunogen, avidity of IgG antibodies can be used as a marker for distinguishing recent primary from long-term infections. Avidity describes the binding strength of a specific antibody to its antigen. Low-avidity IgG antibodies indicate a primary infection, whereas the presence of IgG antibodies with high avidity points to persistency or reactivation of infection.

2. INTENDED USE

The Avidity Measles Virus IgG ELISA is intended to indicate the Measles-specific IgG avidity in human serum or plasma (citrate, heparin) to differentiate between acute and past infection.

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Microplates are coated with specific antigens to bind corresponding antibodies of the sample (dual pipetting). After washing the wells to remove all unbound sample material, one well is incubated with avidity reagent and the corresponding well with washing buffer. The avidity reagent removes the low-avidity antibodies from the antigens whereas the high-avidity ones are still bound to the specific antigens. After second washing step to remove the rest of avidity reagent and low-avidity antibodies, a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a third washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

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4. MATERIALS

4.1. Reagents supplied

- **Avidity Reagent:** 1 bottle containing 15 mL of an urea solution, coloured blue, ready to use; black cap.
- **Control Low:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; $\leq 0.02\%$ (v/v) MIT.
- **Control High:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; $\leq 0.02\%$ (v/v) MIT.

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4.2. Materials supplied

- 1 instruction for use Avidity Measles Virus IgG ELISA (Product Number: AMEA7330)
- 1 Instruction for use Measles Virus IgG ELISA (Product Number: MEAG0330)
- 1 empty labelled bottle (white with white cap) for ready to use Washing Buffer

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6.1. Avidity Reagent

If crystals have formed in the reagent warm up to 37 °C e.g. in a water bath and mix gently until they disappear.

6.2. Washing Buffer

It is recommended to fill 15 mL ready to use Washing Buffer into supplied bottle (s. 4.2) to use it in step Step 5 of the test preparation.

Note: Ready to use Washing Buffer is stable for 5 days at room temperature (20...25 °C).

7. SPECIMEN COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) **Measles Virus IgG positive** samples with this assay.

Note: For samples with high absorbance values (OD > 2.000) appropriate higher dilutions should be used.

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Adjust the incubator to 37 ± 1 °C.

For avidity determination dual pipetting of standards/controls and diluted samples is needed.

1. Dispense 100 μ L standards/controls and diluted samples into their respective wells. Leave wells A1/A2 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 μ L of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 μ L of Avidity Reagent in wells B1, C1, D1, E1 etc, except for the Substrate Blank well A1.
Dispense 100 μ L of Washing Buffer in wells B2, C2, D2, E2 etc, except for the Substrate Blank well A2.
6. **Incubate for exactly 10 min at 37 ± 1 °C.**
7. Repeat step 4.
8. Dispense 100 μ L Conjugate into all wells except in the blank wells (A1/A2).
9. **Incubate for 30 min at room temperature ($20...25$ °C).** Do not expose to direct sunlight.
10. Repeat step 4.
11. Dispense 100 μ L TMB Substrate Solution into all wells.
12. **Incubate for exactly 15 min at room temperature ($20...25$ °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
13. Dispense 100 μ L Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate, thereby a colour change from blue to yellow occurs.
14. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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9.2. Calculation of Results

For each patient sample or control calculate the ratio between the absorbance of the well dispensed with Avidity Reagent and the absorbance of the well dispensed with Washing Buffer multiplied by 100:

$$\frac{\text{Absorbance (sample or control) Avidity Reagent}}{\text{Absorbance (sample or control) Washing Buffer (diluted 1+19)}} \times 100 = \text{Avidity (\%)}$$

Note: For samples with high absorbance values (OD > 2.000) appropriate higher dilutions should be used.

9.3. Interpretation of Results

Result	Avidity	Interpretation
Low-avidity IgG	< 45 %	An avidity index of less than 45 % indicates a primary infection acquired within the past 2 months.
Equivocal	45 – 55 %	No clinical interpretation can be deduced from an equivocal result. It is recommended to take a second sample within an appropriate period of time (e.g. 2 weeks) and repeat testing. If the result of the repeated test is still equivocal, precise statements regarding the time of infection cannot be made.
High-avidity IgG	> 55 %	The presence of high-avidity IgG indicates a past infection or reinfection.

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A result of high avidity can not exclude the possibility of a recent infection.

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IgG	IgM	IgG-Avidity	Probable result
+	-	low	Vague, further investigation necessary
+	-	high	Indicatives of a past infection
+	+	low	Suggests a current or very recent infection
+	+	high	Suggests a past infection with persisting IgM or reactivation of infection

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10.1. Diagnostic Performance

The evaluation of the diagnostic performance of the Avidity Measles Virus IgG test was performed in comparison to well defined samples. The resulting relative agreement was 100 %.

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13. ORDERING INFORMATION

Prod. No.: AMEA7330 Avidity Measles Virus IgG (48 Determinations)

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Product Number: MUMG0340 (96 Determinations)

ENGLISH

1. INTRODUCTION

Mumps viruses are RNA viruses of the family Paramyxoviridae. The virions are spherical particles of 150-250 nm in diameter consisting of a ribonucleoprotein with helical symmetry and enveloped by matrix protein and a lipid bilayer which contains two spike structures: viral hemagglutinin (H) and viral neuraminidase (N). Mumps virus involves primarily the parotid and related salivary glands; however infection can lead to CNS disease and accumulation of the virus in CSF. Mumps (Epidemic Parotitis) is an acute contagious viral disease mostly occurring in children. Nearly 50% of all infections are subclinical. The highest incidence of clinical manifestations is found in the age group of 4 to 15 years. Secondary infections are rare because of long-lasting immunity.

10 to 35 % of mumps cases develop orchitis which occurs nearly always after puberty. The process is mostly unilateral and the prognosis usually good. Mumps virus has been one of the most important causes of viral CNS disease (meningitis and encephalitis) in USA; vaccine administration has greatly reduced its incidence

Species	Disease	Symptoms (e.g.)	Transmission route
Mumps Virus	Mumps	Fever and unilateral or bilateral swelling of the parotid gland; the sublingual and submaxillary glands may also be involved Complications: Orchitis, Meningoencephalitis, Pancreatitis	Virus transmission occurs by droplet infection

Infection or presence of pathogen may be identified by:

- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The Mumps Virus IgG ELISA is intended for the qualitative determination of IgG class antibodies against Mumps Virus in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with Mumps Virus antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-ssoff Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with Mumps Virus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgG Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value **< 0.100**
- **Negative Control:** Absorbance value **< 0.200 and $< \text{Cut-off}$**
- **Cut-off Control:** Absorbance value **$0.150 - 1.300$**
- **Positive Control:** Absorbance value **$> \text{Cut-off}$**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.411	4.14
#2	24	1.173	4.32
#3	24	1.338	1.41
Interassay	n	Mean (NTU)	CV (%)
#1	12	30.07	4.58
#2	12	31.66	5.52
#3	12	2.87	9.27

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95.83% (95% confidence interval: 85.75% - 99.49%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 93.55% (95% confidence interval: 89.41% - 96.43%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: MUMG0340 Mumps Virus IgG ELISA (96 Determinations)

NovaLisa®

Mumps Virus IgM

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português.....	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: MUMM0340 (96 Determinations)

ENGLISH

1. INTRODUCTION

Mumps viruses are RNA viruses of the family Paramyxoviridae. The virions are spherical particles of 150-250 nm in diameter consisting of a ribonucleoprotein with helical symmetry and enveloped by matrix protein and a lipid bilayer which contains two spike structures: viral hemagglutinin (H) and viral neuraminidase (N). Mumps virus involves primarily the parotid and related salivary glands; however infection can lead to CNS disease and accumulation of the virus in CSF. Mumps (Epidemic Parotitis) is an acute contagious viral disease mostly occurring in children. Nearly 50% of all infections are subclinical. The highest incidence of clinical manifestations is found in the age group of 4 to 15 years. Secondary infections are rare because of long-lasting immunity.

10 to 35 % of mumps cases develop orchitis which occurs nearly always after puberty. The process is mostly unilateral and the prognosis usually good. Mumps virus has been one of the most important causes of viral CNS disease (meningitis and encephalitis) in USA; vaccine administration has greatly reduced its incidence.

Species	Disease	Symptoms (e.g.)	Transmission route
Mumps Virus	Mumps	Fever and unilateral or bilateral swelling of the parotid gland; the sublingual and submaxillary glands may also be involved Complications: Orchitis, Meningoencephalitis, Pancreatitis	Virus transmission occurs by droplet infection

Infection or presence of pathogen may be identified by:

- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The Mumps Virus IgM ELISA is intended for the qualitative determination of IgM class antibodies against Mumps Virus in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with Mumps Virus antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with Mumps Virus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgM Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgM Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 μ L to 350 μ L to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 μ L standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 μ L of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 μ L Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 μ L TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 μ L Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value **< 0.100**
- **Negative Control:** Absorbance value **< 0.200 and $< \text{Cut-off}$**
- **Cut-off Control:** Absorbance value **$0.150 - 1.300$**
- **Positive Control:** Absorbance value **$> \text{Cut-off}$**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.600	6.05
#2	24	1.039	4.81
#3	24	0.639	3.71

Interassay	n	Mean (NTU)	CV (%)
#1	12	21.38	7.50
#2	12	13.88	8.95
#3	12	1.39	10.35

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 97.68% (95% confidence interval: 95.03% - 99.15%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 94.44% (95% confidence interval: (72.71% - 99.86%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: MUMM0340 Mumps Virus IgM ELISA (96 Determinations)

Certificate

Quality Management System
EN ISO 13485:2016

Registration No.: SX 1624072-1

Organization: Seramun Diagnostica GmbH
Spreenhagener Str. 1
15754 Heidesee
Germany

Scope: Design and development, manufacture and distribution of in-vitro diagnostics and laboratory systems, consisting of laboratory equipment and associated analysis software, for the detection of autoimmune diseases and infectious diseases

Design and development, manufacture and distribution of substrates and protein stabilizers for in-vitro diagnostics, as well as antibody production and contract development

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Report No.: 3324901-90
Effective date: 2021-03-10
Expiry date: 2024-03-09
Issue date: 2021-03-03

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Serazym[®] Clostridium difficile Toxin A+B

Enzyme immunoassay for detection of *Clostridium difficile* Toxin A and B
in stool specimens and culture suspensions

REF E-040 ▽ 96 REF E-040-A2 ▽ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Clostridium difficile is a bacterium causing nosocomial diarrhea in adults during or after the treatment with antibiotics such as 3rd generation cephalosporines (1). Although 2 - 3% of healthy adults and 20 - 50% of healthy children are colonized with *Clostridium difficile*, the infection is usually of exogenous origin and results from the contact either to hospital staff or to *Clostridium difficile* spores which may contaminate toilets, bed clothes etc. Both exotoxins A and B of this spore-forming bacteria cause the depolymerisation of actin filaments due to the intracellular enzymatic modification of rho-proteins. Consequently, the permeability of cell membrane is raised and neutrophils may invade leading to expression of the clinical picture of the so-called *Clostridium difficile*-associated diarrhea and colitis or finally the pseudomembraneous colitis (PMC) (1). As the production of toxins and the outbreak of disease is correlated, diagnosis of *Clostridium difficile* infection is based mainly on a direct detection of the toxins in stool specimens. To date the cytotoxicity test has been considered as the gold standard for detection of *Clostridium difficile* toxins. Recently it has been replaced to a large extent by immunological tests such as enzyme immunoassay (2).

References:

1. Rambaud J-C., LaMont J-T. (Hrsg.): "Ökosystem Darm Special - Updates on Clostridium difficile" Springer Verlag 1995
2. Wilkins T.D. and Lyster D.M. (2003): „Clostridium difficile Testing: after 20 Years, Still Challenging“ Journal Of Clinical Microbiology, Vol. 41, No. 2, p. 531-534

Intended Use

The *Serazym*[®] Clostridium difficile Toxin A+B is an *in-vitro*-diagnostic medical device for direct detection of the toxins A and B of *Clostridium difficile* in stool specimens and culture suspensions.

Principle of The Test

Serazym[®] Clostridium difficile Toxin A + B is an indirect two-site-immunoassay for the qualitative determination of both *Clostridium difficile* toxins A and B based on polyclonal and monoclonal antibodies. *Clostridium difficile* toxins of stool specimens or culture suspensions and the positive control react with monoclonal anti-toxin A and polyclonal anti-toxin B antibodies coated on the solid phase of the microplate. After incubation non-bound material is removed by a washing step. Subsequently bound toxins specifically react with biotinylated polyclonal anti-toxin A and monoclonal anti-toxin B antibodies during a second incubation period. Non-bound material is separated from the solid-phase immune complexes by a subsequent washing step. During the next incubation period horseradish peroxidase (HRP) conjugated streptavidin reacts with the bound biotinylated antibodies. Unbound conjugate is removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-tetramethylbenzidine (TMB) into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of *Clostridium difficile* toxin A and B. After consideration of the cut-off value, results are interpreted as positive or negative.

Test Components

		For 96 Wells	For 2x 96 Wells	
1	WELLS	Microtitration plate coated with monoclonal anti-Toxin A- (mouse) and polyclonal anti-Toxin B-antibodies (rabbit)	12 single breakable 8-well strips colour coding red vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding red vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control <i>C. difficile</i> Toxin reactive sample	2.0 ml · ready to use coloured blue red cap	4.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>C. difficile</i> Toxin negative sample	2.0 ml · ready to use coloured blue green cap	4.0 ml · ready to use coloured blue green cap
6/1	CONJ BIOTIN	Biotin-conjugate Biotinylated, polyclonal anti-Toxin A- (rabbit) and monoclonal anti-Toxin B-antibodies (mouse)	15 ml · ready to use coloured green white cap	30 ml · ready to use coloured green white cap
6/2	CONJ STREPT	Streptavidin-HRP-conjugate	15 ml · ready to use coloured red brown cap	30 ml · ready to use coloured red brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation And Storage Of Samples

Toxin detection from stool specimens

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours or frozen at at least -20°C. Storage at -20°C as well as repeated freezing and thawing of samples should be avoided.

Formalin-preserved stool samples should not be used in this assay. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing in the ELISA.

Preparation

Warm samples to room temperature and mix thoroughly. Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. If necessary spin down floating particles in a micro centrifuge at maximum speed for 1 min.

Toxin detection from culture suspensions (toxigenic culture)

Colonies of *Clostridium difficile* grown on blood or CCFA agar for 48 hours can be tested directly in the *Serazym*[®] Clostridium difficile Toxin A+B. Prepare a bacterial suspension according to Mc Farland standard 1 (OD value at 600 nm: 0.20 - 0.25 after zero compensation to the yellow coloured *Serazym*[®] sample diluent):

Pipette 1000 µl of sample diluent into a clean tube.

Transfer 2 - 4 inoculating loops of a *C. difficile* culture into the sample diluent and suspend on a vortex mixer. Read OD value at 600 nm as described above where required.

Use 100 µl for ELISA testing. If selective culture media are used the detectable amount of toxins may be reduced due to inhibitory components of such media resulting in decreased OD values in the ELISA. Therefore using selective media for toxigenic culture requires the preparation of a bacterial suspension at least according to Mc Farland standard 4 (OD 600 nm > 1.0 after zero compensation to the yellow coloured *Serazym*[®] sample diluent). In this case the *Clostridium difficile* colonies of at least half of a densely grown agar plate have to be harvested. Where required the recommendations and instructions of the medium manufacturers are to be observed.

Materials Required But Not Provided

Micropipettes · multi-channel pipette or multi-pipette · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer or 8-channel pipette · microplate reader with optical filters for 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation · orbital shaker for performance of test variant 2

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

The *Serazym*[®] Clostridium difficile Toxin A+B can be performed in two ways:

1. Incubation without shaking; complete test duration 2 hours and 15 minutes
2. Incubation with shaking; complete test duration 1 hour and 15 minutes

Dilute samples with sample diluent (3) 1 : 6, e.g. 200 mg or 200 µl stool + 1.0 ml sample diluent (3) or transfer 2 - 4 inoculation loops of a *C. difficile* colony into a tube with 1.0 ml sample diluent (3) and mix thoroughly on a vortex.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle!

Avoid direct light exposure of the TMB substrate solution!

Working steps variant 1: without shaking

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted stool specimen** or **culture suspension**.
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 120 µl) **CONJ BIOTIN** biotin-conjugate (6/1) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 120 µl) **CONJ STREPT** streptavidin-HRP-conjugate (6/2) per well.
9. Cover plate and incubate for 30 min at RT
10. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
11. Dispense 3 drops (or 120 µl) **SUBSTR TMB** substrate (7) per well.
12. Incubate for 15 min at RT protected from light.
13. Dispense 3 drops (or 120 µl) **STOP** stop solution (8) per well, mix gently.
14. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Working steps variant 2: with shaking

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted stool specimen** or **culture suspension**.
3. Cover plate and incubate for 30 min at RT on an orbital shaker with a frequency of 500-700/min.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 120 µl) **CONJ BIOTIN** biotin-conjugate (6/1) per well.
6. Cover plate and incubate for 15 min at RT on an orbital shaker with a frequency of 500-700/min.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 120 µl) **CONJ STREPT** streptavidin-HRP-conjugate (6/2) per well.
9. Cover plate and incubate for 15 min at RT on an orbital shaker with a frequency of 500-700/min.
10. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
11. Dispense 3 drops (or 120 µl) **SUBSTR TMB** substrate (7) per well.
12. Incubate for 15 min at RT without shaking protected from light.
13. Dispense 3 drops (or 120 µl) **STOP** stop solution (8) per well, mix gently.
14. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.20

Samples with absorbances higher than the cut-off value are considered positive, samples with absorbances 10% below the cut-off value are considered negative for *Clostridium difficile* toxin A and B antigen. Samples within 10 % below the cut-off up to the cut-off value have to be considered borderline and should be repeatedly tested. In case of repeated borderline result a second sample of the corresponding patient should be investigated.

Reference Values

Serazym® C. difficile Toxin A+B	
Positive	> Cut-off
Borderline	0,9 x Cut-off – Cut-off
Negative	< 0,9 x Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.20 (manual performance)
 ≤ 0.30 (automatic performance)
- the mean OD of the positive control is ≥ 1.00

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. Formalin treated samples may cause false positive results. A negative test result in the *Serazym*[®] Clostridium difficile Toxin A+B does not exclude an infection: The overall interpretation of the ELISA results should always consider the microbiological examination as well as clinical findings.


Automatic Processing

Performing the *Serazym*[®] Clostridium difficile Toxin A+B on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x-8x.

Correlation: manual – automatic processing

A panel of 125 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.976$.

Serazym® Clostridium difficile Toxin A+B

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*® Clostridium difficile Toxin A+B calculated from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.386	0.042	3.0
2	0.506	0.017	3.3
3	0.332	0.028	8.5

Inter-assay coefficient of variation (CV) in the *Serazym*® Clostridium difficile Toxin A+B in 5 different test runs on 2 different days from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.321	0.102	7.7
2	0.485	0.034	6.9
3	0.345	0.037	10.8

Specificity and sensitivity

A total of 154 stool specimens were tested in parallel with the *Serazym*® Clostridium difficile Toxin A+B and another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym</i> ® ELISA positive	103	4
<i>Serazym</i> ® ELISA negative	2	45

Specificity: 91.8% · Sensitivity: 98.0%

Cross reactivity

Faecal samples positive for one of the following intestinal bacteria did not show any cross reaction in the *Serazym*® Clostridium difficile Toxin A+B:

Staphylococcus aureus, enterotoxin negative; *Staphylococcus aureus*, enterotoxin positive; EHEC; *Pseudomonas aeruginosa*; *Salmonella typhimurium*; *Salmonella enteritidis*; *Salmonella spec.* *Aeromonas hydrophila*; *Aeromonas caviae*; *Campylobacter spec.*; *Hafnia alvei*; *Yersinia enterocolitica* O:3.

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units of the following microorganisms and tested negative with the *Serazym*® ELISA (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)
<i>Bacillus cereus</i>	(ATCC 11778)
<i>Bacillus subtilis</i>	(ATCC 6633)
<i>Bacteroides fragilis</i>	(ATCC 25285)
<i>Candida albicans</i>	(ATCC 10231)
<i>Campylobacter coli</i>	(ATCC 33559)
<i>Campylobacter jejuni</i>	(ATCC 33291)
<i>Citrobacter freundii</i>	(ATCC 8090)
<i>Clostridium sordellii</i>	(ATCC 9714)
<i>Enterobacter aerogenes</i>	(ATCC 13048)
<i>Enterobacter cloacae</i>	(ATCC 13047)
<i>Enterococcus faecalis</i>	(ATCC 29212)

<i>Escherichia coli</i>	(ATCC 25922)
<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Salmonella enterica</i> Serovar <i>enteritidis</i>	(ATCC 13076)
<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	(ATCC 14028)
<i>Shigella flexneri</i>	(ATCC 12022)
<i>Shigella sonnei</i>	(ATCC 25931)
<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Vibrio parahaemolyticus</i>	(ATCC 17802)

The *C. sordellii* strain ATCC 9714 did not cross react in the *Serazym*[®] Clostridium difficile Toxin A+B although some publications describe cross reactivities of toxins of some *C. sordellii* strains with anti-*C. difficile* toxin antibodies.

Common Advices And Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*[®] stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!
Always use protective gloves!
Never pipette material by mouth!
Note safety precautions of the single test components!

History Of Changes

Version	Section	Modifications
2020-07-21	Intended Use	Correction
	Test Components	Correction
	Preparation And Storage Of Samples	Update
	Assay Procedure	Update
	Common Advices And Precautions	Update
	History of Changes	New section "History of Changes"

Incubation Scheme *Serazym*[®] Clostridium difficile Toxin A+B (E-040)

1.  pipette
 100 µl **CONTROL +** (4)
 100 µl **CONTROL -** (5)
 100 µl **diluted stool specimen** or culture suspension
 60 min incubation (RT) **alternatively 30 min while shaking**
 5 x wash with wash solution
2.  3 drops (or 120 µl) **CONJ BIOTIN** (6/1)
 30 min incubation (RT) **alternatively 15 min while shaking**
 5 x wash with wash solution
3.  3 drops (or 120 µl) **CONJ STREPT** (6/2)
 30 min incubation (RT) **alternatively 15 min while shaking**
 5 x wash with wash solution
4.  3 drops (or 120 µl) **SUBSTR TMB** (7)
 15 min incubation protected from light (RT) **without shaking**
5.  3 drops (or 120 µl) **STOP** (8)

Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

Serazym[®] Clostridium difficile GDH

Enzyme immunoassay for detection of *Clostridium difficile* GDH in faecal specimens

REF E-107 ▽ 96 REF E-107-A2 ▽ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

The spore forming anaerobic bacterium *Clostridium difficile* is the most common cause of nosocomial diarrhea predominantly developing in the course of antibiotic therapy. The severity of the disease in dependence of the health status of the patient but also of the pathogenicity of the *C. difficile* strain may range from bloody diarrhea to Colitis and pseudomembranous Colitis up to toxic megacolon. The two exotoxins A and B are responsible for the development of the symptoms. Most *Clostridium difficile* strains produce both, toxin A and B, but some are characterized by isolated toxin B production. Atoxigenic strains are generally considered as nonpathogenic and therefore diagnostically irrelevant. The *Clostridium difficile* specific enzyme Glutamatdehydrogenase, also described as “common antigen” is produced by both: toxigenic and atoxigenic strains.

The human intestine is colonized with *Clostridium difficile* in 2 - 3% of healthy adults and in up to 50% of children < 2 years of age. Therefore, the diagnosis of a *Clostridium difficile* associated diarrhea (CDAD) has to be confirmed by the detection of the toxins A and B.

Changes in the epidemiological situation are observed since several years: Increasingly also young adults and not hospitalized patients get sick developing the symptoms of CDAD. An increasing number of severe developments of the disease caused by highly virulent strains also occur. In this context testing of patient's status as carrier seems reasonable.

Different approaches have developed in laboratory diagnosis of *C. difficile* infections. Usually direct toxin A and B detection from stool or enrichment cultures is performed by enzyme immunoassay. In addition, a new two-step procedure has been established in several laboratories that comprise a first screening of *C. difficile* GDH and the subsequent testing for toxins A and B of GDH positive samples. The isolated detection of GDH is neither sufficient for diagnosing a CDAD nor suitable for a therapeutic decision, but always needs to be completed by the testing for toxins A and B.

References:

1. Wilkins TD and Lyerly DM (2003): „Clostridium difficile Testing: after 20 Years, Still Challenging“ Journal Of Clinical Microbiology, Vol. 41, No.2, p531-534
2. Von Eichel-Streiber C und Braun V (2008): “Das difficile Clostridium” Journal Of Laboratory Medicine, Vol. 32, No. 4, p219-234
3. Lyerly DM, Barroso LA, Wilkins TD (1991): “Identification of the latex test-reactive protein of Clostridium difficile as glutamate dehydrogenase.” Journal Of Clinical Microbiology, Vol. 29, p2639-2642.
4. Zhen L, Keller SF, Lyerly DM et al. (2004): „Multicenter Evaluation of a New Screening Test That Detects Clostridium difficile in Faecal Specimens“ Journal Of Clinical Microbiology, Vol. 42, No. 8, p3837-3840.
5. Peterson LR and Robicsek A (2009): “Does My Patient Have Clostridium difficile Infection?” Annals of Internal Medicine, Vol. 151, No. 3, p176-180.

Intended Use

Serazym® Clostridium difficile GDH is an *in-vitro* diagnostic medical device for direct detection of Glutamatdehydrogenase (GDH) of *Clostridium difficile* in faecal specimens.

Principle of the Test

Serazym® Clostridium difficile GDH is a one-step enzyme immunoassay on the basis of polyclonal and monoclonal antibodies to Glutamatdehydrogenase (GDH) of *Clostridium difficile*. Diluted stool specimens and monoclonal anti-GDH antibodies are dispensed simultaneously into the wells of a microtitration plate coated with polyclonal anti-GDH antibodies. After an incubation time of 60 min at room temperature (RT) unbound components are removed by a washing step. Alternatively, the first incubation can be shortened to 30 min by shaking. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-Tetramethylbenzidine (TMB) within a 10 min reaction time into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of GDH antigen in the respective sample. Considering the cut-off value results are interpreted as positive or negative.

Test Components

			For 96 Wells	For 2x 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti- <i>C.difficile</i> GDH antibodies (sheep)	12 single breakable 8-well strips colour coding green vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding green vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control Recombinant <i>C. difficile</i> GDH	2.0 ml · ready to use coloured blue red cap	4.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>C. difficile</i> GDH negative sample	2.0 ml · ready to use coloured blue green cap	4.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, monoclonal anti- <i>C.difficile</i> GDH antibodies (mouse)	15 ml · ready to use coloured green brown cap	25 ml · ready to use coloured green brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation and Storage of Samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at least at -20°C. Repeated (> 3x) freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8 °C before testing in the ELISA.

Preparation

Mix samples thoroughly. Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. Spin down floating particles in a micro centrifuge at maximum speed for 1 min if necessary.

Materials required but not provided

Micropipettes · multi-channel pipette or multi-pipette · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer or 8-channel pipette · microplate reader with optical filters for 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation · orbital shaker for performance of test variant 2

Preparation and Storage of Reagents

Kit size and expiry

One kit is designed for 1 x 96 or 2 x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

The *Serazym*[®] Clostridium difficile GDH can be performed in two ways:

1. Incubation without shaking; complete test duration 1 h and 10 minutes
2. Incubation with shaking; complete test duration 40 minutes

Dilute samples with sample diluent (3) 1 : 6 e.g. 200 mg or 200 µl stool + 1.0 ml sample diluent (3) and mix thoroughly.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that residual fluid is completely drained in every single wash cycle!

Avoid direct light exposure of the TMB substrate solution!

Working steps variant 1: without shaking

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Dispense 3 drops (or 100 µl) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted sample**, mix gently
4. Cover plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
6. Dispense 3 drops (or 100 µl) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 3 drops (or 100 µl) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Working steps variant 2: with shaking

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Dispense 3 drops (or 100 μ l) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 100 μ l **CONTROL +** positive control (4)
100 μ l **CONTROL -** negative control (5)
100 μ l **diluted sample**, mix gently
4. Cover plate and incubate for 30 min at RT on an orbital shaker with a frequency of 500 – 700 / min.
5. Decant, then wash each well 5x with 300 μ l wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
6. Dispense 3 drops (or 100 μ l) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 3 drops (or 100 μ l) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.10

Samples with OD values equal to or higher than the cut-off value are considered positive, samples with OD values below the cut-off value are considered negative in the *Serazym*[®] C. difficile GDH.

Reference Values

<i>Serazym</i>[®] Clostridium difficile GDH	
Positive	\geq Cut-off
Negative	$<$ Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is \leq 0.20 (manual performance)
 \leq 0.30 (automatic performance)
- the mean OD of the positive control is \geq 1.00

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

The qualitative determination of *Clostridium difficile* GDH in stool specimens by enzyme immunoassay is not equivalent to the diagnosis of a *C. difficile* associated disease (CDAD). CDADs are caused by the toxins A and B of pathogenic *C. difficile* strains. Therefore, a positive GDH test result has to be supplemented by testing for toxins A and B, in order to confirm or exclude a toxigenic strain. On the other hand a negative GDH test result does not necessarily exclude a *C. difficile* infection. After dilution in sample diluent the samples should be tested by ELISA as quickly as possible, at least within 72 hours, because storage-depending antigen degradation may cause false negative results. More than 3 freeze-thaw cycles may also cause false negative results due to antigen degradation. As a result of inhomogeneous antigen distribution in some stool samples insufficient homogenizing may also cause false negative results. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. The overall interpretation of the ELISA results should always consider the microbiological examination as well as clinical findings.

Automatic Processing

Performing the *Serazym*[®] *Clostridium difficile* GDH on fully automated microplate processors (e.g. DS2 or DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x - 8x.

Correlation: manual – automatic processing

A panel of 90 stool specimens (45 positive and 45 negative samples) was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.999$.

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) calculated from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	2.800	0.040	1.52
2	1.960	0.045	2.46
3	0.611	0.041	7.19
4	0.352	0.017	5.15

Inter-assay coefficient of variation (CV) in 5 different test runs calculated from twofold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.707	0.250	14.6
2	1.155	0.148	12.8
3	0.895	0.098	10.9
4	0.381	0.068	17.9

Lot-to-Lot reproducibility in 3 different production lots calculated from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.840	0.166	9.0
2	1.285	0.172	13.4
3	0.349	0.093	26.6
4	0.026	0.005	17.6

Lower detection limit

The lower detection limit of Glutamatdehydrogenase (GDH) in the *Serazym*[®] Clostridium difficile GDH was determined 10 ng/ml by titration of recombinant GDH antigen.

Sensitivity and specificity

Sensitivity in comparison to PCR

Ninety eight out of 102 stool specimens characterized *Clostridium difficile* positive by PCR were tested positive with the *Serazym*[®] ELISA corresponding to a sensitivity of 96.1%.

Comparative sensitivity and specificity

In two comparative studies 235 and 170 stool samples were tested in parallel in the *Serazym*[®] Clostridium difficile GDH and 2 other commercially available ELISAs respectively.

Study 1

n = 235	comparative ELISA 1 positive	comparative ELISA 1 negative
<i>Serazym</i> [®] ELISA positive	101	3**
<i>Serazym</i> [®] ELISA negative	12*	119

Sensitivity: 89.4% · Specificity: 97.5%

* 10 out of 12 *Serazym*[®] ELISA negative and comparative ELISA 1 positive samples were tested negative in 2 other commercial ELISAs. Sensitivity amended: 98.1%

** One sample was confirmed true positive by PCR. Specificity amended: 98.3%

Study 2

n = 170	comparative ELISA 2 positive	comparative ELISA 2 negative
<i>Serazym</i> [®] ELISA positive	69	1*
<i>Serazym</i> [®] ELISA negative	3	97

Sensitivity: 95.8% · Specificity: 98.9%

* This sample was confirmed true positive by PCR.

Specificity amended: 100%

Cross reactivity

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units of the following microorganisms and tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)
<i>Bacillus cereus</i>	(ATCC 11778)
<i>Bacillus subtilis</i>	(ATCC 6633)
<i>Bacteroides fragilis</i>	(ATCC 25285)
<i>Candida albicans</i>	(ATCC 10231)
<i>Campylobacter coli</i>	(ATCC 33559)
<i>Campylobacter jejuni</i>	(ATCC 33291)
<i>Citrobacter freundii</i>	(ATCC 8090)
<i>Clostridium sordellii</i>	(ATCC 9714)
<i>Enterobacter aerogenes</i>	(ATCC 13048)
<i>Enterobacter cloacae</i>	(ATCC 13047)
<i>Enterococcus faecalis</i>	(ATCC 29212)
<i>Escherichia coli</i>	(ATCC 25922)

<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Salmonella enterica</i> Serovar <i>enteritidis</i>	(ATCC 13076)
<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	(ATCC 14028)
<i>Shigella flexneri</i>	(ATCC 12022)
<i>Shigella sonnei</i>	(ATCC 25931)
<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Vibrio cholerae</i>	Clinical isolate
<i>Yersinia enterocolitica</i> Serotyp <i>O3, O9</i>	Clinical isolates

Interference

None of the following substances added to GDH positive and negative stool samples showed a significant impact on the test result:

barium sulphate (5%), Buscopan® (2 mg/ml), cyclamate (5%), Diclofenac (2 mg/ml), haemoglobine (5 mg/ml), Hylak® N (5%), Immodium® akut duo (0.2/12.5 mg/ml), Iberogast® (5%), loperamide (0.2 mg/ml), metronidazole (2 mg/ml), mucin (5 mg/ml), Nexium® (2 mg/ml), palmitic acid (20%), Pentofuryl® (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie® (8 mg/ml), Simagel® (2 mg/ml), stearic acid (20%), vancomycin hydrochloride (0.5%).

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*® stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2021-02-02	Preparation And Storage Of Samples, Assay Procedure	Update

Incubation Scheme *Serazym*[®] Clostridium difficile GDH (E-107)

1.  3 drops (or 100 µl) **CONJ HRP** (6)
 +
 pipette
 100 µl **CONTROL +** (4)
 100 µl **CONTROL -** (5)
 100 µl **stool sample**, mix gently
 60 min incubation (room temperature)
alternatively 30 min while shaking

 5x wash with wash solution
2.  3 drops (or 100 µl) **SUBSTR TMB** (7)
 10 min incubation protected from light (room temperature)
without shaking
3.  3 drops (or 100 µl) **STOP** (8)
Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

NovaLisa[®]

Clostridium tetani toxin IgG

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	8
Français	14
Italiano	20
Español	26
Bibliography / Literatur / Bibliographie / Bibliografía / Bibliografia	32
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	32
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: TETG0430 (96 Determinations)

ENGLISH

1. INTRODUCTION

Clostridia are spore-forming gram-positive bacteria. The round spores are built at the terminal end which results in the microscope in a "tennis racket" like shape.

Tetanus develops only when spores of *Clostridium tetani* germinate under strict anaerobic conditions after gaining access to wounds and small lacerations. The clinical manifestation of the disease is primarily not caused by the invasion of the exciter but by the secretion of a powerful neurotoxin (tetanospasmin). This toxin blocks the inhibition of the signal transduction and has a high affinity to the central nervous system. The consequence is hyper excitability of the muscles to external stimuli in combination with a principal increase of the muscle tonus without influence of consciousness. It starts with tonic spasm of muscles (trismus), mimic muscles and gullet muscles. Neck, back and abdominal musculature follow. At the same time the appearance of refractory spasm of whole muscle groups can hamper breathing. Hyper salivation and swallowing problems cause aspiration and pneumonia with the next breath.

Clostridium tetani is ubiquitous present in soil and intestine of humans and animals. Ingestion of bacteria or growth in the intestine of man or animal is without harm. The spores are extremely resistant towards heat and can stay infectious for a long period. The bacteria can get under the skin by even smallest wounds. In Europe tetanus mainly occurs after injuries and sometimes postoperative whereas in developing countries Tetanus is widely disseminated. The WHO assumes that one million people die because of tetanus worldwide per year.

Tetanus toxin is an excellent immunogen in man - only one antigenic type of toxin. The only effective way to control tetanus is by prophylactic active immunization.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Clostridium tetani</i>	Tetanus	Trismus, dysphagia, severe, painful spasms of whole muscle groups, hypersalivation, aspiration, asphyxia	Injury (Infection of the wound with <i>Clostridium tetani</i>)

The presence of pathogen or infection may be identified by

- Microscopy
- Serology: e.g. ELISA

2. INTENDED USE

The *Clostridium tetani* toxin IgG ELISA is intended for the quantitative determination of IgG class antibodies against *Clostridium tetani* toxin in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The quantitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Clostridium tetani toxin Coated Microplate (IgG):** 12 break apart 8-well snap-off strips coated with Clostridium tetani toxin (toxoid) antigens; in resealable aluminium foil.
- **IgG Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Clostridium tetani toxin anti-IgG Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap; < 5% NMP.
- **Clostridium tetani toxin IgG Standards:** 4 vials, each containing 2 ml standard (human serum or plasma); coloured yellow; ready to use.

Standard A:	0.0	IU/ml; blue cap
Standard B:	0.1	IU/ml; green cap
Standard C:	0.5	IU/ml; yellow cap
Standard D:	1.0	IU/ml; red cap

The standards are calibrated in accordance with the Who International Standard; "1st International Standard for Tetanus Immunoglobulin, Human"; NIBSC Code: TE-3.

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Coated Microplate

The break-apart snap-off strips are coated with Clostridium tetani toxin (toxoid) antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 ml Washing Buffer + 190 ml distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37°C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. For CSF please use the instruction for use ABVL0001. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 10 µl sample and 1 ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three to five and the volume of Washing Buffer from 300 µl to 350 µl to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µl standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µl Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25°C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µl TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.2. Measurement

Adjust the ELISA microwell plate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA microwell plate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrat-Blank:** Absorbance value < 0.100
- **Standard A:** Absorbance value < 0.200
- **Standard B:** Absorbance value > 0.150
- **Standard C:** Absorbance value > 0.500
- **Standard D:** Absorbance value > 1.000

Standard A < Standard B < Standard C < Standard D

If these criteria are not met, the test is not valid and must be repeated.


9.2. Calculation of Results

In order to obtain **quantitative results in IU/ml** blot the (mean) absorbance values of the 4 Standards A - D on (linear/linear) graph paper in a system of coordinates against their corresponding concentrations (0.0 / 0.1 / 0.5 and 1.0 IU/ml) and draw a standard curve (absorbance values on the y-axis, concentrations on the x-axis).

Read results from this standard curve employing the (mean) absorbance values of each patient sample.

For the calculation of the standard-curve mathematical Point to Point function should be used.

9.3. Typical standard Curve

9.4. Interpretation of Results and Recommendations [IU/ml]

< 0.1 IU/ml	No protective antibody level or no reliable protection! Immediate full course of basic immunization or booster injection and control of antibody concentration 4 to 6 weeks later is recommended.
0.11 - 0.5 IU/ml	Reliable protection! Booster injection and control of antibody concentration 4 to 6 weeks later is recommended.
0.51 - 1.0 IU/ml	Reliable protection; control of antibody concentration after about 2 years is recommended. Booster injection is not required. Note: In cases of antibody concentrations greater than 0.5 IU/ml vaccination can cause side effects!
1.1 - 5.0 IU/ml	Range of long term protection: Control after 5 to 10 years
> 5.0 IU/ml	Range of long term protection: Control after 10 years immunisation and to record the data on the certificate of vaccination.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.

In immunocompromised patients and newborns serological data only have restricted value.

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	Cv (%)
#1	24	1.306	3.60
#2	24	1.805	3.46
#3	24	1.591	5.34
Interassay	n	Mean (IU/ml)	Cv (%)
#1	12	0.060	9.62
#2	12	0.084	11.33
#3	12	0.658	13.99

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100.0% (95% confidence interval: 76.84% - 100.0%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 99.22% (95% confidence interval: 95.76% - 99.98%).

10.4. Analytical Sensitivity

The analytical sensitivity (according to CLSI EP17-A) is defined as the apparent concentration of the analyte that can be distinguished from the zero calibrator. It is 0.01 IU/ml.

10.5. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/ml hemoglobin, 5 mg/ml triglycerides and 0.5 mg/ml bilirubin.

10.6. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

10.7. Measurement range

The measurement range is 0.01 IU/ml – 1 IU/ml.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: TETG0430 Clostridium tetani toxin IgG ELISA (96 Determinations)

NovaLisa[®]

Coxiella burnetii (Q-Fever) Phase 1 IgG

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português.....	27
Bibliography / Literatur / Bibliographie / Bibliografía / Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: COX1G0600 (96 Determinations)

ENGLISH

1. INTRODUCTION

Q-Fever is a disease that results from infection with small, polymorph and gram-negative bacteria called *Coxiella burnetii*. After an outbreak in Brisbane, Australia, the responsible organism was isolated and named *Coxiella burnetii* in honour of Dr. Herald Rae Cox and Sir Frank Burnet. New molecular research demonstrated a close relationship to *Legionella*. The zoonosis Q-Fever is found everywhere except New Zealand (no data available). There is an extensive reservoir (mainly ticks) of *C. burnetii*. Ticks are an important vector of the pathogen in the transmission between domestic and wildlife animals. But the ticks are unimportant in the direct infection of humans. Cattles, sheep and goats are usually the source of transmission of this microorganism to humans. However cats, dogs and rabbits are also important in this regard. In most instances humans become infected with *Coxiella burnetii* following inhalation of contaminated aerosols (respiratory tract). The incubation period for Q-Fever in humans is about 2 weeks. The resulting illness can be divided into acute and chronic varieties. During the acute phase of illness antibodies to the phase 2-antigen are formed. Anti phase-1 antibodies in high titers are typical for a chronic disease. In areas where Q-Fever is endemic, 12% or more of the population have antibodies to *C. burnetii*. Most of the infections are subclinical or undiagnosed.

The acute infection shows symptoms of high fever, shivers, muscle pain and headache. Later on more severe diseases such as pneumonia or hepatitis can occur. Infections during pregnancy can lead to an abort or premature birth. Approximately 1% of all infections become chronic. The most frequent organ manifestation in Q-Fever is endocarditis.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Coxiella burnetii</i>	Q-Fever	Pneumonia white fever, headache, muscle pain; Hepatitis; Myocarditis and Endocarditis.	Inhalation of contaminated aerosols (respiratory tract). Transmission by ingestion of contaminated products such as milk or meat is likely.

Infection or presence of pathogen may be identified by:

- Cell culture
- PCR
- Serologia: e.g. ELISA, IFT

2. INTENDED USE

The *Coxiella burnetii* (Q-Fever) Phase 1 IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Coxiella burnetii* (Q-Fever) Phase 1 in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Coxiella burnetii* (Q-Fever) Phase 1 antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with *Coxiella burnetii* (Q-Fever) Phase 1 antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37°C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgG Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value **< 0.100**
- **Negative Control:** Absorbance value **< 0.200 and $< \text{Cut-off}$**
- **Cut-off Control:** Absorbance value **$0.150 - 1.300$**
- **Positive Control:** Absorbance value **$> \text{Cut-off}$**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1 Antibody Isotypes and State of Infection

Serology	Significance
Phase II IgM	Characteristic of the primary antibody response May persist for several months
Phase II IgG	Characteristic of the primary antibody response May persist for several years
Phase I IgG	Characteristic of the chronic infection Occasionally during convalescent phase

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.303	6.75
#2	24	0.601	7.03
#3	24	1.380	11.19
Interassay	n	Mean (NTU)	CV (%)
#1	12	21.41	4.84
#2	12	47.55	6.23
#3	12	2.13	8.84

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 86.77% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 59.04% - 100%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: COX1G0600 Coxiella burnetii (Q-Fever) Phase 1 IgG ELISA (96 Determinations)

NovaLisa[®]

Coxiella burnetii (Q-Fever) Phase 2 IgG

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português.....	27
Bibliography / Literatur / Bibliographie / Bibliografía / Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: COX2G0600 (96 Determinations)

ENGLISH

1. INTRODUCTION

Q-Fever is a disease that results from infection with small, polymorph and gram-negative bacteria called *Coxiella burnetii*. After an outbreak in Brisbane, Australia, the responsible organism was isolated and named *Coxiella burnetii* in honour of Dr. Herald Rae Cox and Sir Frank Burnet. New molecular research demonstrated a close relationship to *Legionella*. The zoonosis Q-Fever is found everywhere except New Zealand (no data available). There is an extensive reservoir (mainly ticks) of *C. burnetii*. Ticks are an important vector of the pathogen in the transmission between domestic and wildlife animals. But the ticks are unimportant in the direct infection of humans. Cattles, sheep and goats are usually the source of transmission of this microorganism to humans. However cats, dogs and rabbits are also important in this regard. In most instances humans become infected with *Coxiella burnetii* following inhalation of contaminated aerosols (respiratory tract). The incubation period for Q-Fever in humans is about 2 weeks. The resulting illness can be divided into acute and chronic varieties. During the acute phase of illness antibodies to the phase 2-antigen are formed. Anti phase-1 antibodies in high titers are typical for a chronic disease. In areas where Q-Fever is endemic, 12% or more of the population have antibodies to *C. burnetii*. Most of the infections are subclinical or undiagnosed.

The acute infection shows symptoms of high fever, shivers, muscle pain and headache. Later on more severe diseases such as pneumonia or hepatitis can occur. Infections during pregnancy can lead to an abort or premature birth. Approximately 1% of all infections become chronic. The most frequent organ manifestation in Q-Fever is endocarditis.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Coxiella burnetii</i>	Q-Fever	Pneumonia white fever, headache, muscle pain; Hepatitis; Myocarditis and Endocarditis.	Inhalation of contaminated aerosols (respiratory tract). Transmission by ingestion of contaminated products such as milk or meat is likely.

Infection or presence of pathogen may be identified:

- Cell culture
- PCR
- Serology: e.g. ELISA, IFT

2. INTENDED USE

The *Coxiella burnetii* (Q-Fever) Phase 2 IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Coxiella burnetii* (Q-Fever) Phase 2 in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Coxiella burnetii* (Q-Fever) Phase 2 antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37°C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with *Coxiella burnetii* (Q-Fever) Phase 2 antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgG Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value **< 0.100**
- **Negative Control:** Absorbance value **< 0.200 and $< \text{Cut-off}$**
- **Cut-off Control:** Absorbance value **$0.150 - 1.300$**
- **Positive Control:** Absorbance value **$> \text{Cut-off}$**

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
Phase II IgM	Characteristic of the primary antibody response May persist for several months
Phase II IgG	Characteristic of the primary antibody response May persist for several years
Phase I IgG	Characteristic of the chronic infection Occasionally during convalescent phase

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.428	9.69
#2	24	0.794	8.51
#3	24	0.645	6.27
Interassay	n	Mean (NTU)	CV (%)
#1	12	17.91	8.48
#2	12	14.30	12.09
#3	12	1.93	7.38

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 91.4% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 39.76% - 100%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet.

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: COX2G0600 Coxiella burnetii (Q-Fever) Phase 2 IgG ELISA (96 Determinations)

NovaLisa®

Coxiella burnetii (Q-Fever) Phase 2 IgM

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía/ Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procedure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: COX2M0600 (96 Determinations)

ENGLISH

1. INTRODUCTION

Q-Fever is a disease that results from infection with small, polymorph and gram-negative bacteria called *Coxiella burnetii*. After an outbreak in Brisbane, Australia, the responsible organism was isolated and named *Coxiella burnetii* in honour of Dr. Herald Rae Cox and Sir Frank Burnet. New molecular research demonstrated a close relationship to *Legionella*. The zoonosis Q-Fever is found everywhere except New Zealand (no data available). There is an extensive reservoir (mainly ticks) of *C. burnetii*. Ticks are an important vector of the pathogen in the transmission between domestic and wildlife animals. But the ticks are unimportant in the direct infection of humans. Cattles, sheep and goats are usually the source of transmission of this microorganism to humans. However cats, dogs and rabbits are also important in this regard. In most instances humans become infected with *Coxiella burnetii* following inhalation of contaminated aerosols (respiratory tract). The incubation period for Q-Fever in humans is about 2 weeks. The resulting illness can be divided into acute and chronic varieties. During the acute phase of illness antibodies to the phase II-antigen are formed. Anti phase-I antibodies in high titers are typical for a chronic disease. In areas where Q-Fever is endemic, 12% or more of the population have antibodies to *C. burnetii*. Most of the infections are subclinical or undiagnosed.

The acute infection shows symptoms of high fever, shivers, muscle pain and headache. Later on more severe diseases such as pneumonia or hepatitis can occur. Infections during pregnancy can lead to an abort or premature birth. Approximately 1% of all infections become chronic. The most frequent organ manifestation in Q-Fever is endocarditis.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Coxiella burnetii</i>	Q-Fever	Pneumonia white fever, headache, muscle pain; Hepatitis; Myocarditis and Endocarditis.	Inhalation of contaminated aerosols (respiratory tract). Transmission by ingestion of contaminated products such as milk or meat is likely.

Infection or presence of pathogen may be identified by:

- Cell culture
- PCR
- Serology: e.g. ELISA, IFT

2. INTENDED USE

The *Coxiella burnetii* (Q-Fever) Phase 2 IgM ELISA is intended for the qualitative determination of IgM class antibodies against *Coxiella burnetii* (Q-Fever) Phase 2 in the early stages of infection in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Coxiella burnetii* (Q-Fever) Phase 2 antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 12.1.

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Microtiterplate

The break-apart snap-off strips are coated with *Coxiella burnetii* (Q-Fever) Phase 2 antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37°C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgM Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgM Sample Dilution Buffer into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.1. Measurement

Adjust the ELISA Microtiterplate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay run to be considered valid, these Instructions for Use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < 0.100
- **Negative Control:** Absorbance value < 0.200 and $< \text{Cut-off}$
- **Cut-off Control:** Absorbance value $0.150 - 1.300$
- **Positive Control:** Absorbance value $> \text{Cut-off}$

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
Phase II IgM	Characteristic of the primary antibody response May persist for several months
Phase II IgG	Characteristic of the primary antibody response May persist for several years
Phase I IgG	Characteristic of the chronic infection Occasionally during convalescent phase

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24	0.485	3.11
#2	24	0.824	3.86
#3	24	0.794	3.23
Interassay	n	Mean (NTU)	CV (%)
#1	12	20.76	3.16
#2	12	14.60	8.05
#3	12	1.02	6.76

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95.95 % (95 % confidence interval: 88.61 % - 99.16 %).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 95.74 % (95 % confidence interval: 85.46 % - 99.48 %).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

10.5. Cross Reactivity

Cross reactivity with antibodies against CMV, Dengue virus, EBV and Mycoplasma cannot be excluded.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

12.1. Safety note for reagents containing hazardous substances

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Therefore, the following hazard and precautionary statements apply.

Warning

H317	May cause an allergic skin reaction.
P261	Avoid breathing spray.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P362+P364	Take off contaminated and Wash it before reuse.

Further information can be found in the safety data sheet

12.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

13. ORDERING INFORMATION

Prod. No.: COX2M0600 Coxiella burnetii (Q-Fever) Phase 2 IgM ELISA (96 Determinations)

Serazym[®] Adenovirus

Enzyme immunoassay for detection of *Adenovirus* in faecal samples

REF E-017 ∇ 96 REF E-017-A2 ∇ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Adenovirus is the causative agent of respiratory tract infections, conjunctivitis or enteritis. *Adenovirus* infections are spread via faecal-oral transmission or aerosols (1). Most infections are mild to moderate and do not last longer than one week. *Adenovirus* is responsible for 2 - 8% of respiratory tract infections and 7 - 17% of diarrheal diseases in children (2). Thirty percent of viral diarrheas in immunocompromised patients are caused by *Adenovirus* infections (1). The diagnosis of *Adenovirus* infections is preferably done by direct detection in faecal specimens or swabs. Due to the difficulty of *Adenovirus* culture on tissue or cell culture virus detection so far has been performed by electron microscopy. Meanwhile immunological methods like enzyme immunoassay (ELISA) have been developed for antigen detection (3). ELISA techniques are based on poly- and/or monoclonal antibodies to the protein hexon representing the major part of the virus capsid.

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Intended Use

Serazym® Adenovirus is an *in-vitro*-diagnostic medical device for direct detection of *Adenovirus* in faecal samples.

Principle Of The Test

Serazym® Adenovirus is a one-step enzyme immunoassay on the basis of monoclonal antibodies against an epitope of the capsid protein hexon, common to all human pathogenic *Adenovirus* serotypes. Diluted stool specimens and horseradish peroxidase (HRP) labelled monoclonal anti-Adenovirus-antibodies are dispensed simultaneously into the wells of a microtitration plate coated with monoclonal anti-Adenovirus-antibodies. After an incubation time of 60 min at room temperature (RT) unbound components are removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-Tetramethylbenzidine (TMB) within a 10 min reaction time at room temperature protected from light into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of *Adenovirus*. Considering the cut-off value results are interpreted as positive or negative.

Test Components

			For 96 Wells	For 2x 96 Wells
1	WELLS	Microtitration plate coated with monoclonal anti-Adenovirus-antibodies (mouse)	12 single breakable 8-well strips colour coding violet vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding violet vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control <i>Adenovirus</i> antigen Ad 41 (inactivated)	1.5 ml · ready to use coloured blue red cap	3.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>Adenovirus</i> negative sample	1.5 ml · ready to use coloured blue green cap	3.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, monoclonal anti-Adenovirus-antibodies (mouse)	12 ml · ready to use coloured green brown cap	24 ml · ready to use coloured green brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation And Storage Of Samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing in the ELISA.

Preparation

Quickly thaw frozen samples. Warm samples to room temperature and mix well.

The *Serazym*[®] Adenovirus can be performed with 1 : 6 or 1 : 11 diluted specimens. In case of additional testing of the same sample in the *Serazym*[®] Norovirus, the *Serazym*[®] Campylobacter or the *Serazym*[®] Clostridium difficile Toxin A+B the 1 : 6 dilution is recommended.

Preparation of a 1 : 11 sample dilution:

Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 100 mg (diameter about 2 - 3 mm) of faeces if solid or pipette 100 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Preparation of a 1 : 6 sample dilution:

Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Materials Required But Not Provided

Micropipettes · multi-channel pipette or multi-pipette · reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least 30 days when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute samples with sample diluent (3) 1 : 11 or 1 : 6, e.g. 100 mg or 100 µl stool + 1.0 ml (1 : 11) sample diluent (3) or 200 mg or 200 µl stool + 1.0 ml (1 : 6) sample diluent (3).
Avoid any time shift during dispensing of reagents and samples.
Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!
Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Dispense 2 drops (or 75 µl) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 75 µl **CONTROL +** positive control (4)
75 µl **CONTROL -** negative control (5)
50 µl **diluted sample**, mix gently.
4. Cover plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
6. Dispense 2 drops (or 75 µl) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 2 drops (or 75 µl) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction Stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.20

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for *Adenovirus* antigen.

Reference Values

Serazym® Adenovirus	
Positive	\geq Cut-off
Negative	$<$ Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.20 (manual performance)
 ≤ 0.30 (automatic performance)
- the mean OD of the positive control is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. A negative test result not necessarily excludes an *Adenovirus* infection. Inhomogeneous virus distribution in the sample can cause false negative results. The investigation of samples that were taken beyond the acute phase of the disease can cause false negative results, because the number of virus particles has decreased under the detection limit of the test. It is therefore recommended to take samples within the acute phase of the disease where a maximum number of excreted virus particles are to be expected. A final interpretation of the test results should consider clinical findings as well.

Automatic Processing

Performing the *Serazym*[®] Adenovirus on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x - 8x.

Correlation: manual – automatic processing

A panel of 110 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.974$.

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*[®] Adenovirus from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	2.792	0.160	5.7
2	2.059	0.167	8.1
3	1.368	0.094	6.9
4	0.718	0.068	9.4

Inter-assay coefficient of variation (CV) in the *Serazym*[®] Adenovirus in 6 different test runs from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.850	0.107	5.8
2	1.057	0.069	6.5
3	0.574	0.042	7.3
4	0.312	0.030	9.6

Lower detection limit

The lower detection limit of *Adenovirus* antigen in the *Serazym*[®] Adenovirus was determined by titration of purified *Adenovirus* antigen (hexon). Lower detection limit: 6 ng / ml

Specificity and sensitivity

A total of 330 stool samples were investigated in parallel in the *Serazym*[®] Adenovirus and in another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym</i> [®] ELISA positive	55	1
<i>Serazym</i> [®] ELISA negative	2	272

Specificity: 99.6% · Sensitivity: 96.5%

Cross reactivity

Stool samples positive for one of the subsequent pathogens have been tested with the *Serazym*[®] Adenovirus and showed no cross reactivity:

Rotavirus (n = 10), *Astrovirus* (n = 8), *Norovirus* (n = 31), *Clostridium difficile* (n = 11), *Campylobacter jejuni* (n = 7), *Campylobacter coli* (n = 1), *Salmonella enteritidis* (n = 18), *Giardia lamblia* (n = 1).

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units of the following microorganisms and tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)	<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Bacillus cereus</i>	(ATCC 11778)	<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Bacillus subtilis</i>	(ATCC 6633)	<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Bacteroides fragilis</i>	(ATCC 25285)	<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Candida albicans</i>	(ATCC 10231)	<i>Salmonella enterica Serovar enteritidis</i>	(ATCC 13076)
<i>Campylobacter coli</i>	(ATCC 33559)	<i>Salmonella enterica Serovar typhimurium</i>	(ATCC 14028)
<i>Campylobacter jejuni</i>	(ATCC 33291)	<i>Shigella flexneri</i>	(ATCC 12022)
<i>Citrobacter freundii</i>	(ATCC 8090)	<i>Shigella sonnei</i>	(ATCC 25931)
<i>Clostridium sordellii</i>	(ATCC 9714)	<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Enterobacter aerogenes</i>	(ATCC 13048)	<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Enterobacter cloacae</i>	(ATCC 13047)	<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Enterococcus faecalis</i>	(ATCC 29212)	<i>Vibrio cholerae</i>	clinical isolate
<i>Escherichia coli</i>	(ATCC 25922)	<i>Yersinia enterocolitica Serotyp 03, 09</i>	clinical isolates

Interference

None of the following substances added to positive and negative stool samples showed a significant impact on the test result:

Barium sulfate (5%), Buscopan® (2 mg/ml), Cyclamate (5%), Diclofenac (2 mg/ml), Hemoglobine human (5 mg/ml), Blood human (1,25%), Hylak® N (5%), Iberogast® (5%), Immodium® akut duo (0.2/12.5 mg/ml), Loperamide (0.2 mg/ml), Metronidazole (2 mg/ml), Mucin (5 mg/ml), Nexium® (2 mg/ml), Palmitic acid (20%), Pentofuryl® (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie® (8 mg/ml), Simigel® (2 mg/ml), Stearic acid (20%), Vancomycin (2 mg/ml).

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the Serazym® stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!
Always use protective gloves!
Never pipette material by mouth!
Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2020-11-04	Common Advices and Precautions	Update
	History of Changes	Adding section "History of Changes "
	Test Components	Update

Incubation Scheme *Serazym*[®] Adenovirus (E-017)

1.  2 drops (or 75 µl) **CONJ HRP** (6)
 +
 pipette
 75 µl **CONTROL +** (4)
 75 µl **CONTROL -** (5)
 50 µl **diluted stool sample**, mix gently
 60 min incubation (room temperature)

5 x wash with wash solution
2.  2 drops (or 75 µl) **SUBSTR TMB** (7)
 10 min incubation (room temperature) protected from light
3.  2 drops (or 75 µl) **STOP** (8)

Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

Serazym[®] Astrovirus

Enzyme immunoassay for detection of *Astrovirus* in faecal samples

REF E-045 ∇ 96 REF E-045-A2 ∇ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Astrovirus was firstly described in 1975 and named according to its star-shaped structure visible under the electron microscope. *Astrovirus* belongs to the family *Astroviridae*. Human *Astroviruses* are subdivided into 7 serotypes (1). Together with Rotavirus and Adenovirus *Astrovirus* is one of the most common causes of non-bacterial gastroenteritis in children under 5 years of age all over the world. Thus 80% of children between 5 and 10 years of age are anti-Astrovirus-antibody positive. *Astrovirus* caused gastroenteritis in adults and nosocomial infections are observed as well (2). The course of the disease is usually self-limiting and of short duration. After the incubation time of 1 - 2 days a 1 - 4 days lasting gastroenteritis develops accompanied by vomiting, diarrhea, fever and abdominal pain finally causing dehydration. Although occurring all over the year *Astrovirus* infections are mainly observed during the winter months (3, 4). *Astrovirus* infections are spread via faecal-oral transmission from person to person or via contaminated things or food. Infected persons excrete high amounts of *Astrovirus* particles with their faeces (1, 2). The detection of *Astrovirus* may be performed by electron microscopy or by molecular biology techniques such as polymerase chain reaction (PCR). Meanwhile immunological methods like enzyme immunoassay have established as preferential methods for routine laboratory diagnosis since these methods are fast, safe and automation is possible (1).

References:

1. Rohwedder, A. (2000): "Virale Gastroenteritiden, Erreger und Diagnostik", *Mikrobiologie*, 10. Jg. p.121-126.
2. Palombo, E. A. and Bishop, R. F. (1996): "Annual Incidence, Serotype Distribution and Genetic Diversity of Human Astrovirus Isolates from Hospitalized Children in Melbourne, Australia"; *Journal of Clinical Microbiology*, Vol. 34, No. 7, p. 1750-1753.
3. Cukor, G. and Blacklow, N. R. (1984): "Human Viral Gastroenteritis", *Microbiological Reviews*, June, Vol. 48 No. 2, p. 157-179.

4. Gaggero, A.; O’Ryan, M. et al. (1998): “Prevalence of Astrovirus Infection among Chilean Children with Acute Gastroenteritis”, Journal of Clinical Microbiology, Vol. 36 No. 12, p. 3691-3693.

Intended use

Serazym® Astrovirus is an *in-vitro*-diagnostic medical device for direct detection of *Astrovirus* in faecal samples.

Principle of the test

Serazym® Astrovirus is a one-step enzyme immunoassay on the basis of polyclonal and monoclonal antibodies against *Astrovirus* antigens. Diluted stool specimens and horseradish peroxidase (HRP) labelled monoclonal anti-Astrovirus-antibodies are dispensed simultaneously into the wells of a microtitration plate coated with polyclonal anti-Astrovirus-antibodies. After an incubation time of 60 min at room temperature (RT) unbound components are removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-Tetramethylbenzidine (TMB) within a 10 min reaction time at room temperature protected from light into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of *Astrovirus*.

Test components

		For 96 Wells	For 2x 96 Wells	
1	WELLS	Microtitration plate coated with polyclonal anti-Astrovirus-antibodies (sheep)	12 single breakable 8-well strips colour coding light blue vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding light blue vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer, 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control <i>Astrovirus</i> reactive sample	1.5 ml · ready to use coloured blue red cap	3.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>Astrovirus</i> negative sample	1.5 ml · ready to use coloured blue green cap	3.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, monoclonal anti-Astrovirus-antibodies	12 ml · ready to use coloured green brown cap	24 ml · ready to use coloured green brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation and storage of samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 48 h at 2...8°C before testing in the ELISA.

Preparation

Quickly thaw frozen samples. Warm samples to room temperature and mix well.

The *Serazym*[®] Astrovirus can be performed with 1 : 6 or 1 : 11 diluted specimens. In case of additional testing of the same sample in the *Serazym*[®] Norovirus, the *Serazym*[®] Campylobacter or the *Serazym*[®] Clostridium difficile Toxin A+B the 1 : 6 dilution is recommended.

Preparation of a 1 : 11 sample dilution: Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 100 mg (diameter about 2 - 3 mm) of faeces if solid or pipette 100 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Preparation of a 1 : 6 sample dilution: Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Materials required but not provided

Micropipettes · multi-channel pipette or multi-pipette · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least 30 days when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay procedure

Dilute samples with sample diluent (3) 1 : 11 or 1 : 6, e.g. 100 mg or 100 µl stool + 1.0 ml (1 : 11) sample diluent (3) or 200 mg or 200 µl stool + 1.0 ml (1 : 6) sample diluent (3).

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use.
Mix gently without causing foam.
2. Dispense 2 drops (or 75 µl) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 75 µl **CONTROL +** positive control (4)
75 µl **CONTROL -** negative control (5)
50 µl **diluted sample**, mix gently.
4. Cover plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
6. Dispense 2 drops (or 75 µl) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 2 drops (or 75 µl) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

Result interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.10

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for *Astrovirus* antigen.

Reference values

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.15 (manual test performance)
 ≤ 0.30 (automatic test performance)
- the mean OD of the positive control is ≥ 1.00

If the above-mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. A negative test result not necessarily excludes an *Astrovirus* infection. Inhomogeneous virus distribution in the sample can cause false negative results. The investigation of samples that were taken beyond the acute phase of the disease can cause false negative results, because the number of virus particles has decreased under the detection limit of the test. It is therefore recommended to take samples within the acute phase of the disease where a maximum number of excreted virus particles are to be expected. A final interpretation of the test results should consider clinical findings as well.


Automatic Processing

Performing the *Serazym*[®] *Astrovirus* on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary, the number of washing steps can be enhanced from 5x to 7x - 8x.

<i>Serazym</i>[®] <i>Astrovirus</i>	
Positive	\geq Cut-off
Negative	$<$ Cut-off

Correlation: manual – automatic processing

A panel of 96 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.993$.

Performance characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym®* Astrovirus from 8-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	1.667	0.148	8.9
2	0.994	0.063	6.4
3	0.443	0.027	6.1
4	0.185	0.018	9.8

Inter-assay coefficient of variation (CV) in the *Serazym®* Astrovirus in 6 different test runs from 8-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	1.853	0.071	3.8
2	1.019	0.059	5.8
3	0.583	0.069	11.9
4	0.350	0.034	9.7

Lower detection limit

The lower detection limit of *Astrovirus* antigen in the *Serazym®* Astrovirus was determined by titration of purified *Astrovirus*-antigen. Lower detection limit: 6 ng / ml.

Specificity and sensitivity

A total of 98 stool samples were investigated in parallel in the *Serazym®* Astrovirus and in another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym®</i> ELISA positive	49	0
<i>Serazym®</i> ELISA negative	2	47

Specificity: 100% Sensitivity: 96%

Cross reactivity

Stool samples positive for one of the subsequent pathogens have been tested with the *Serazym*[®] Astrovirus and showed no cross reactivity:

Rotavirus (n = 10), *Adenovirus* (n = 20), *Norovirus* (n = 31), *Clostridium difficile* (n = 11), *Campylobacter jejuni* (n = 7), *Campylobacter coli* (n = 1), *Salmonella enteritidis* (n = 18), *Giardia lamblia* (n = 1).

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units of the following microorganisms and tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)	<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Bacillus cereus</i>	(ATCC 11778)	<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Bacillus subtilis</i>	(ATCC 6633)	<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Bacteroides fragilis</i>	(ATCC 25285)	<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Candida albicans</i>	(ATCC 10231)	<i>Salmonella enterica Serovar enteritidis</i>	(ATCC 13076)
<i>Campylobacter coli</i>	(ATCC 33559)	<i>Salmonella enterica Serovar typhimurium</i>	(ATCC 14028)
<i>Campylobacter jejuni</i>	(ATCC 33291)	<i>Shigella flexneri</i>	(ATCC 12022)
<i>Citrobacter freundii</i>	(ATCC 8090)	<i>Shigella sonnei</i>	(ATCC 25931)
<i>Clostridium sordellii</i>	(ATCC 9714)	<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Enterobacter aerogenes</i>	(ATCC 13048)	<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Enterobacter cloacae</i>	(ATCC 13047)	<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Enterococcus faecalis</i>	(ATCC 29212)	<i>Vibrio cholerae</i>	Clinical isolate
<i>Escherichia coli</i>	(ATCC 25922)	<i>Yersinia enterocolitica Serotyp O3, O9</i>	Clinical isolates

Interference

None of the following substances added to positive and negative stool samples showed a significant impact on the test result:

barium sulphate (5%), Buscopan[®] (2 mg/ml), cyclamate (5%), Diclofenac (2 mg/ml), haemoglobine (5 mg/ml), Hylak[®] N (5%), Immodium[®] akut duo (0.2/12.5 mg/ml), Iberogast[®] (5%), loperamide (0.2 mg/ml), metronidazole (2 mg/ml), mucin (5 mg/ml), Nexium[®] (2 mg/ml), palmitic acid (20%), Pentofuryl[®] (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie[®] (8 mg/ml), Simagel[®] (2 mg/ml), stearic acid (20%).

* The addition of 5% (v/v) Hylak[®] N (lactic acid containing preparation against digestive complaints) to Astrovirus positive stool suspensions may decrease OD values.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution are universally applicable for the *Serazym*[®] stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!
Always use protective gloves!
Never pipette material by mouth!
Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2020-07-21	Intended use	Correction
	Test Components	Correction
	Assay Procedure	Update
	Common Advices and Precautions	Update
	History of Changes	New section "History of Changes"

Incubation scheme *Serazym*[®] Astrovirus (E-045)


1.  2 drops (or 75 µl) **CONJ HRP** (6)
 +
 pipette
 75 µl **CONTROL +** (4)
 75 µl **CONTROL -** (5)
 50 µl **stool sample**, mix gently

 60 min incubation (room temperature)

 5 x wash with wash solution
2.  2 drops (or 75 µl) **SUBSTR TMB** (7)

 10 min incubation (room temperature) protected from light
3.  2 drops (or 75 µl) **STOP** (8)

Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

Serazym[®] Norovirus

Enzyme immunoassay for detection of *Norovirus* in stool specimens

REF E-061 ▾ 96 REF E-061-A2 ▾ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Noroviruses belong to the family *Caliciviridae*, single stranded RNA viruses of 30 - 40 nm in size characterized by a typical cup-shaped capsid. Within the genus *Norovirus* the two human pathogenic genogroups GG1 and GGII have been identified. GG1 and GGII strains can be further subclassified into at least 15 and 18 genotypes resp. The genetic heterogeneity of Noroviruses causes distinct capsid protein divergences between different genogroups (about 60%) as well as between different genotypes within one genogroup (about 20 - 30%). Since 1994 genotype GGII.4 is predominantly circulating. Noroviruses are very resistant to environmental conditions and highly contagious. The infection is transmitted by direct contact to already infected people either by faecal-oral transmission or by ingestion of aerosols from vomit or by contaminated food, drinking water or objects. After a short, 10 - 50 hours lasting incubation time fulminant diarrhea and often vomiting develop as the characteristic symptoms. The infection is usually self-limiting and symptoms disappear after 2 - 3 days. Norovirus infections are characterized by seasonal fluctuations with a climax during the winter months. They are considered as the most common cause of non-bacterial gastroenteritis outbreaks worldwide, but may also be responsible for single cases of viral gastroenteritis. The high sequence variability of the capsid proteins circumvents the production of protective antibodies and hampers diagnostic detection. Methods like PCR (usually as "Real time Reverse Transcriptase PCR – Rt RT-PCR) and enzyme immunoassay are commonly used for laboratory diagnosis.

References:

1. Venkataram, B.V. et al. (1999): "X-ray Crystallographic Structure of the Norwalk Virus Capsid" *Science* Vol 286: 287-290.
2. Künkel, U. und Schreier, E. (2002): "Caliciviren-Virale Auslöser akuter Gastroenteritiden" *Bundesgesundheitsbl-Gesundheitsforsch-Gesundheitsschutz* 45: 534-542.
3. Hansman, G.S. et al. (2006): "Genetic and antigenic diversity among noroviruses" *Journal of General Virology* 87:909-919.
4. Marshall, J.A. and Bruggink, L.D. (2006): "Laboratory Diagnosis of Norovirus" *Clin. Lab.* 52: 571-581
5. Tomoyuki Shiota et al. (2007): "Characterization of a Broadly Reactive Monoclonal Antibody against Norovirus Genogroups I and II: Recognition of a Novel Conformational Epitope" *Journal of Virology*, Vol. 81, No. 22, p.12298-12306.
6. Lindesmith, L.C. et al. 2008: „Mechanisms of GII.4 Norovirus Persistence in Human Populations“ *PLoS Medicine* Vol 5 (2): 269-289

Intended Use

The *Serazym*[®] Norovirus is an *in-vitro*-diagnostic medical device for direct detection of *Norovirus* specific antigens in stool specimens.

Principle Of The Test

Norovirus antigens from stool specimens and the positive control react with polyclonal anti-Norovirus antibodies coated on the solid phase of the microplate. After incubation non-bound material is removed by a washing step. Subsequently bound antigens specifically react with horseradish peroxidase (HRP) labeled polyclonal anti-Norovirus antibodies during a second incubation period. Unbound conjugate is removed by a washing step. HRP converts the subsequently added colourless chromogenic substrate solution (TMB / H₂O₂) into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of *Norovirus* antigens. By means of the calculated cut-off value, results are interpreted either as positive or negative.

Test Components

			For 96 Wells	For 2x 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti-Norovirus antibodies (sheep)	12 single breakable 8-well strips colour coding silver vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding silver vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control Recombinant Norovirus capsid proteins	2.0 ml · ready to use coloured blue red cap	4.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control Norovirus negative sample	2.0 ml · ready to use coloured blue green cap	4.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled polyclonal anti-Norovirus antibodies (sheep)	15 ml · ready to use coloured green brown cap	25 ml · ready to use coloured green brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation And Storage Of Samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing in the ELISA.

Preparation

Warm samples to room temperature and mix thoroughly. Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. If necessary spin down floating particles in a micro centrifuge for one min at maximum speed.

Materials Required But Not Provided

Micropipettes · multi-channel pipette or multi-pipette · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with 450 nm filter for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute samples with sample diluent (3) 1 : 6, e.g. 200 mg or 200 µl stool + 1.0 ml sample diluent (3) and mix thoroughly on a vortex.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid direct light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted stool specimen**.
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 100 µl) **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 100 µl) **SUBSTR TMB** substrate (7) per well.
9. Incubate for 10 min at RT protected from light.
10. Dispense 3 drops (or 100 µl) **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.10

Samples with OD values higher than or equal to the cut-off value are considered positive, samples with OD values below the cut-off value are considered negative for *Norovirus* antigen.

Reference Values

<i>Serazym</i> [®] Norovirus	
Positive	≥ Cut-off
Negative	< Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.20 (manual test performance)
≤ 0.30 (automatic test performance)
- the mean OD of the positive control is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. A negative test result not necessarily excludes a *Norovirus* infection. Inhomogeneous virus distribution in the sample can cause false negative results. The investigation of samples that were taken beyond the acute phase of the disease can cause false negative results, because the number of virus particles has decreased under the detection limit of the test. It is therefore recommended to take samples within the acute phase of the disease where a maximum number of excreted virus particles are to be expected. Genetic recombination between different *Norovirus* strains may cause antigenic shift finally leading to the occurrence of virus variants that are not detected by ELISA. The overall interpretation of the ELISA results should always consider clinical findings.


Automatic Processing

Performing the *Serazym*[®] Norovirus on fully automated microplate processors (e.g. DS2 or DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control.

It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x-8x.

Correlation: manual – automatic processing

A panel of 90 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.999$.

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym®* Norovirus calculated from 8-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	2.132	0.043	2.00
2	0.902	0.035	3.93
3	0.534	0.017	3.28
4	0.217	0.013	5.92

Inter-assay coefficient of variation (CV) in the *Serazym®* Norovirus in 8 different test runs on 2 different days from 8-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	1.924	0.146	7.59
2	0.813	0.026	3.17
3	0.562	0.021	3.74
4	0.247	0.007	2.76

Lower detection limit

The lower detection limit of the *Serazym®* Norovirus was determined < 10 ng/ml capsid protein for genogroup I and II.

Specificity and sensitivity

One retrospective study with 159 stool specimens was performed to compare the *Serazym®* Norovirus with another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym®</i> ELISA positive	111	3
<i>Serazym®</i> ELISA negative	6	39

Specificity: 92.9% Sensitivity: 94.9%

Cross reactivity

Faecal samples positive for one of the following pathogens did not show any cross reaction in the *Serazym*[®] Norovirus:

Adenovirus (n = 6), *Astrovirus* (n = 8), *Rotavirus* (n = 6), *Clostridium difficile* (n = 8), *Campylobacter jejuni* (n = 6), *Helicobacter pylori* (n = 5), *Giardia lamblia* (n = 8), *Cryptosporidium parvum* (n = 7), *Entamoeba histolytica/dispar* (n = 6).

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units and virus particles per ml stool suspension resp. All microorganisms were tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-Off):

<i>Adenovirus</i>	typ e41	<i>Peptostreptococcus anaerobius</i>	ATCC 27337
<i>Aeromonas hydrophila</i>	ATCC 7966	<i>Proteus vulgaris</i>	ATCC 8427
<i>Astrovirus</i>	serotype 4	<i>Pseudomonas aeruginosa</i>	ATCC 10145
<i>Bacillus cereus</i>	ATCC 117788	<i>Rotavirus</i>	strain SA11
<i>Bacillus subtilis</i>	ATCC 6633	<i>Salmonella enterica</i> Serovar <i>thyphimurium</i>	ATCC 14028
<i>Bacteroides fragilis</i>	ATCC 25285	<i>Salmonella enterica</i> ssp. <i>galolyticus</i>	ATCC 13076
<i>Campylobacter jejuni</i>	ATCC 33291	<i>Shigella flexneri</i>	ATCC 12022
<i>Candida albicans</i>	clinical isolate	<i>Shigella sonnei</i>	ATCC 25931
<i>Citrobacter freundii</i>	ATCC 8090	<i>Staphylococcus aureus</i>	ATCC 25923
<i>Clostridium sordelli</i>	ATCC 9714	<i>Staphylococcus epidermidis</i>	ATCC 12228
<i>Enterobacter aerogenes</i>	ATCC 13048	<i>Vibrio cholerae</i>	RV 2011/ST5
<i>Enterobacter cloacae</i>	ATCC 13047	<i>Vibrio parahaemolyticus</i>	ATCC 17802
<i>Enterococcus faecalis</i>	ATCC 29212	<i>Yersinia enterocolitica</i> 0:9	clinical isolate
<i>Escherichia coli</i>	ATCC 25922	<i>Yersinia enterocolitica</i> 0:3	clinical isolate
<i>Klebsiella pneumonia</i>	ATCC 13883		

Interference

None of the following substances added to positive and negative stool samples showed a significant impact on the test result:

Barium sulfate (5%), Buscopan[®] (2 mg/ml), Cyclamate (5%), Diclofenac (2 mg/ml), Hemoglobine human (5 mg/ml), Blood human (5%), Hylak[®] N (5%), Iberogast[®] (5%), Immodium[®] akut duo (0.2/12.5 mg/ml), Loperamide (0.2 mg/ml), Metronidazole (2 mg/ml), Mucin (5 mg/ml), Nexium[®] (2 mg/ml), Palmitic acid (20%), Pentofuryl[®] (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie[®] (8 mg/ml), Simagel[®] (2 mg/ml), Stearic acid (20%), Vancomycin (2 mg/ml).

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*[®] stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!
Always use protective gloves!
Never pipette material by mouth!
Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2021-02-02	Common Advices and Precautions	Update
	History of Changes	Adding section "History of Changes"
	Test Components	Update
	Assay Procedure	Update

Incubation Scheme *Serazym*[®] Norovirus (E-061)

- | | | | |
|----|---|---|--|
| 1. |  | pipette
100 µl
100 µl
100 µl
60 min
5 x wash | CONTROL + (4)
CONTROL - (5)
diluted stool sample
incubation (room temperature)
with wash solution |
| 2. |  | 3 drops (or 100 µl)
30 min
5 x wash | CONJ HRP (6)
incubation (room temperature)
with wash solution |
| 3. |  | 3 drops (or 100 µl)
10 min | SUBSTR TMB (7)
incubation (room temperature) protected from light |
| 4. |  | 3 drops (or 100 µl) | STOP (8) |
- Read OD at 450 / ≥ 620 nm**

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

NovaLisa[®]

Brucella IgG

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: BRUG0050 (96 Determinations)

ENGLISH

1. INTRODUCTION

Brucella is a small Gram-negative bacterium (0.4-0.8 µm in diameter and 0.4-3.0 µm in length) which is non-flagellated, and non-spore-forming. They are named after the military doctor David Bruce, who on Malta in 1887 isolated the pathogens from the spleen of a soldier who died of undulating fever. Four species are pathogenic to human: Brucella abortus, Brucella melitensis, Brucella suis and Brucella canis. All four species are excitors of Brucellosis, a disease characterized by undulating fever. Depending on exciter the disease is also called Morbus Bang (B. abortus) or Malta fever (B. melitensis).

The pathogens are transmitted from animals, which are mainly affected. The infection is caused by contact with ill animals or their excrements as well as by non-pasteurized milk and milk products like fresh cheese from sheep or goat. Main entrances are skin wounds, conjunctives and digestive tract. The intact pathogens are transported by granulocytes into local lymph nodes, from where they spread haematogenous. All kind of organs can be infected. Symptoms depend up on the infected organ, where an inflammation takes place. Until now the pathogenic mechanism is not completely understood.

Brucellosis appears worldwide. In non-pasteurized milk and milk products Brucella is viably and infectious for weeks. Bovine brucellosis caused by Brucella abortus is still the most widespread form, although reported incidence and prevalence of the disease vary widely from country to country (from <0.01 to >200 per 100,000 population). Brucella melitensis is endemic in areas where keeping of sheeps and goats is frequent. It causes serious human infections. Working with these bacteria in laboratories necessitates highest carefulness because of high contagiousity. Risk groups include abattoir workers, meat inspectors, animal handlers, veterinarians, and laboratorians. Brucellosis is a nationally notifiable disease and reportable to the local health authority.

Species	Disease	Symptoms (e.g.)	Transmission route
B. abortus (cattle) B. melitensis (sheep, goats) B. suis (pigs) B. canis (dogs)	Brucella	Fever, chills (undulating fever), malaise, arthritis, hepatitis, endocarditis, hepatomegalie, osteomyelitis (OM)	Oral (non-pasteurized milk and milk products) Percutan (contact with ill animals or their excrements) In general no transmission from human to human

The presence of pathogen or infection may be identified by

- Histology
- Serology: e.g. ELISA

2. INTENDED USE

The Brucella IgG ELISA is intended for the qualitative determination of IgG class antibodies against Brucella in human serum or plasma (citrate, heparin).

3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA microwell plate reader.

4. MATERIALS

4.1. Reagents supplied

- **Brucella Coated Microplate (IgG):** 12 break-apart 8-well snap-off strips coated with Brucella antigens; in resealable aluminium foil.
- **IgG Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Brucella anti-IgG Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap; < 5 % NMP.
- **Brucella IgG Positive Control:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap.
- **Brucella IgG Cut-off Control:** 1 vial containing 3 ml control (human serum or plasma); coloured yellow; ready to use; green cap.
- **Brucella IgG Negative Control:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap.

For potential hazardous substances please check the safety data sheet.

4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

4.3. Materials and Equipment needed

- ELISA microwell plate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Pipettes to deliver volumes between 10 and 1000 µl
- Vortex tube mixer
- Distilled water
- Disposable tubes

5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

6.1. Coated Microplate

The break-apart snap-off strips are coated with Brucella antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 + 19; e. g. 10 ml Washing Buffer + 190 ml distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37°C e.g. in a water bath. Mix well before dilution.

6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. For CSF please use the instruction for use ABVL0001. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

Heat inactivation of samples is not recommended.

7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgG Sample Diluent. Dispense 10 µl sample and 1 ml IgG Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

8. ASSAY PROCEDURE

8.1. Test Preparation

Please read the instruction for use carefully before performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three to five and the volume of Washing Buffer from 300 µl to 350 µl to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µl standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µl Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µl TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

8.2. Measurement

Adjust the ELISA microwell plate reader **to zero** using the **Substrate Blank**.

If - due to technical reasons - the ELISA microwell plate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at **450 nm** and record the absorbance values for each standard/control and sample in the plate layout.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the **mean absorbance values** of all duplicates.

9. RESULTS

9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate Blank:** Absorbance value < 0.100
- **Negative Control:** Absorbance value < 0.200 and $< \text{Cut-off}$
- **Cut-off Control:** Absorbance value $0.150 - 1.300$
- **Positive Control:** Absorbance value $> \text{Cut-off}$

If these criteria are not met, the test is not valid and must be repeated.

9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43
Cut-off = 0.43

9.2.1. Results in Units [NTU]

$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$

Example: $\frac{1.591 \times 10}{0.43} = 37 \text{ NTU (Units)}$

9.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative .
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.
In immunocompromised patients and newborns serological data only have restricted value.

9.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

10.1. Precision

Intraassay	n	Mean (OD)	CV (%)
#1	24	0.577	4.14
#2	24	1.276	3.34
#3	24	1.200	2.75
Interassay	n	Mean (NTU)	CV (%)
#1	12	23.22	4.97
#2	12	20.13	6.05
#3	12	5.10	8.55

10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte.
It is 98.78% (95% confidence interval: 93.39% - 99.97%)

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte.
It is 100.0% (95% confidence interval: 66.37% - 100.0%).

10.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/ml hemoglobin, 5 mg/ml triglycerides and 0.5 mg/ml bilirubin.

10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal evidence of false-positive results due to cross-reactions.

11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

12. PRECAUTIONS AND WARNINGS

- In compliance with article 1 paragraph 2b European directive 98/79/EC the use of the in vitro diagnostic medical devices is intended by the manufacturer to secure suitability, performances and safety of the product. Therefore the test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or strips of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel who are familiar with good laboratory practice.

12.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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Product Number: BRUM0050 (96 Determinations)

ENGLISH

1. INTRODUCTION

Brucella is a small Gram-negative bacterium (0.4-0.8 µm in diameter and 0.4-3.0 µm in length) which is non-flagellated, and non-spore-forming. They are named after the military doctor David Bruce, who on Malta in 1887 isolated the pathogens from the spleen of a soldier who died of undulating fever. Four species are pathogenic to human: Brucella abortus, Brucella melitensis, Brucella suis and Brucella canis. All four species are exciter of Brucellosis, a disease characterized by undulating fever. Depending on exciter the disease is also called Morbus Bang (B. abortus) or Malta fever (B. melitensis).

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- Histology
- Serology: e.g. ELISA

2. INTENDED USE

The Brucella IgM ELISA is intended for the qualitative determination of IgM class antibodies against Brucella in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Brucella Coated Microplate (IgM):** 12 break-apart 8-well snap-off strips coated with Brucella antigens; in resealable aluminium foil.
- **IgM Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2 ; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2 , for washing the wells; white cap.
- **Brucella anti-IgM Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap; < 5 % NMP.
- **Brucella IgM Positive Control:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap.
- **Brucella IgM Cut-off Control:** 1 vial containing 3 ml control (human serum or plasma); coloured yellow; ready to use; green cap.
- **Brucella IgM Negative Control:** 1 vial containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap.

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7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgM Sample Diluent. Dispense 10 µl sample and 1 ml IgM Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100.0% (95% confidence interval: 97.49% - 100.0%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100.0% (95% confidence interval: 81.47% - 100.0%).

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13. ORDERING INFORMATION

Prod. No.: BRUM0050 Brucella IgM ELISA (96 Determinations)

NovaLisa®

Hantavirus IgG

ELISA

CE

Only for in-vitro diagnostic use

English	2
Deutsch	7
Français	12
Italiano	17
Español	22
Português	27
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía/ Bibliografia	34
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	34
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	35
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste	36

Product Number: HANG0670 (96 Determinations)

ENGLISH

1. INTRODUCTION

Hantaviruses are negative sense RNA viruses in the Bunyaviridae family. Humans may be infected with Hantaviruses through urine, saliva or contact with rodent waste products. Some Hantaviruses may lead to serious diseases in humans, such as hemorrhagic fever with renal syndrome (HFRS) and hantavirus pulmonary syndrome (HPS).

Human infections of Hantaviruses have almost entirely been linked to human contact with rodent excrement, but recent human-to-human transmission has been reported with the Andes virus in South America.

Hantavirus has an incubation time of two to four weeks in humans before symptoms of infection occur. The symptoms of HFRS can be split into five phases:

- Febrile phase: Symptoms include fever, chills, sweaty palms, diarrhea, malaise, headaches, nausea, abdominal and back pain, respiratory problems such as the ones common in influenza virus infection, as well as gastro-intestinal problems. These symptoms normally occur for three to seven days and arise about two to three weeks after exposure.
- Hypotensive phase: This occurs when the blood platelet levels drop and symptoms can lead to tachycardia and hypoxemia. This phase can last for 2 days.
- Oliguric phase: This phase lasts for three to seven days and is characterized by the onset of renal failure and proteinuria occurs.
- Diuretic phase: This is characterized by diuresis of three to six liters per day, which can last for a couple of days up to weeks.
- Convalescent phase: This is normally when recovery occurs and symptoms begin to improve.

Regions especially affected by HFRS include China, the Korean Peninsula, Russia (Hantaan, Puumala and Seoul viruses), and northern and western Europe (Puumala and Dobrava virus).

Species	Disease	Symptoms (e.g.)	Transmission route
Puumala virus Dobrava virus Hantaan virus Seoul virus	Hemorrhagic fever with renal syndrome (HFRS)	Initial: suddenly occurring symptoms like intense headache, back and abdominal pain, fever, chills, nausea, and blurred vision. Late: low blood pressure, acute shock, vascular leakage, and acute kidney failure	After exposure to aerosolized urine, droppings, or saliva of infected rodents or their nests (airborne transmission). Also by direct contact with these materials to broken skin or onto mucous membranes.
Andes virus Sin-Nombre-virus (New world strains)	Hantavirus pulmonary syndrome (HPS)	Initial: universal symptoms include fatigue, fever and muscle aches, especially in the large muscle groups - thighs, hips, back, and sometimes shoulders. There may also be headache, dizziness, chills, and abdominal problems, such as nausea, vomiting, diarrhea, and abdominal pain. Late: coughing and shortness of breath, lungs fill with fluid.	Bites by infected rodents. Human to human transmission can not be excluded (for New World strains).

The presence of pathogen or infection may be identified by:

- PCR
- Serology (e. g. ELISA)

2. INTENDED USE

The Hantavirus IgG ELISA is intended for the qualitative determination of IgG antibodies against Hantavirus in human serum or plasma (citrate or heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Hantavirus Coated Microplate (IgG):** 12 breakapart 8-well snap-off strips coated with recombinant Hantavirus antigens in resealable aluminium foil.
- **IgG Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M); pH 7.2 ± 0.2; for washing the wells; white cap.
- **Hantavirus anti-IgG Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgG; in phosphate buffer (10 mM); coloured blue, ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1%; ready to use; yellow cap; < 5% NMP.
- **Hantavirus IgG Positive Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap.
- **Hantavirus IgG Cut-off Control:** 1 bottle containing 3 ml control (human serum or plasma); coloured yellow; ready to use; green cap.
- **Hantavirus IgG Negative Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap.

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6.1. Coated Microplate

The break-apart snap-off strips are coated with recombinant Hantavirus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Use human serum or plasma (citrate or heparin) samples with this assay. For CSF please use the instruction for use ABVL0001. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 96.59% (95% confidence interval: 90.36% - 99.29%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 99.16% (95% confidence interval: 95.41% - 99.98%).

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10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal significant evidence of false-positive results due to cross-reactions.

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Product Number: HANM0670 (96 Determinations)

ENGLISH

1. INTRODUCTION

Hantaviruses are negative sense RNA viruses in the Bunyaviridae family. Humans may be infected with Hantaviruses through urine, saliva or contact with rodent waste products. Some Hantaviruses may cause serious diseases in humans, such as hemorrhagic fever with renal syndrome (HFRS) and hantavirus pulmonary syndrome (HPS).

Human infections of Hantaviruses have almost entirely been linked to human contact with rodent excrement, but recent human-to-human transmission has been reported with the Andes virus in South America.

Hantavirus has an incubation time of two to four weeks in humans before symptoms of infection occur. The symptoms of HFRS can be split into five phases:

- Febrile phase: Symptoms include fever, chills, sweaty palms, diarrhea, malaise, headaches, nausea, abdominal and back pain, respiratory problems such as the ones common in influenza virus infection, as well as gastro-intestinal problems. These symptoms normally occur for three to seven days and arise about two to three weeks after exposure.
- Hypotensive phase: This occurs when the blood platelet levels drop and symptoms can lead to tachycardia and hypoxemia. This phase can last for 2 days.
- Oliguric phase: This phase lasts for three to seven days and is characterized by the onset of renal failure and proteinuria occurs.
- Diuretic phase: This is characterized by diuresis of three to six liters per day, which can last for a couple of days up to weeks.
- Convalescent phase: This is normally when recovery occurs and symptoms begin to improve.

Regions especially affected by HFRS include China, the Korean Peninsula, Russia (Hantaan, Puumala and Seoul viruses), and northern and western Europe (Puumala and Dobrava virus).

Species	Disease	Symptoms (e.g.)	Transmission route
Puumala virus Dobrava virus Hantaan virus Seoul virus	Hemorrhagic fever with renal syndrome (HFRS)	Initial: suddenly occurring symptoms like intense headache, back and abdominal pain, fever, chills, nausea, and blurred vision. Late: low blood pressure, acute shock, vascular leakage, and acute kidney failure	After exposure to aerosolized urine, droppings, or saliva of infected rodents or their nests (airborne transmission). Also by direct contact with these materials to broken skin or onto mucous membranes.
Andes virus Sin-Nombre-virus (New world strains)	Hantavirus pulmonary syndrome (HPS)	Early: universal symptoms include fatigue, fever and muscle aches, especially in the large muscle groups - thighs, hips, back, and sometimes shoulders. There may also be headaches, dizziness, chills, and abdominal problems, such as nausea, vomiting, diarrhea, and abdominal pain. Late: coughing and shortness of breath, lungs fill with fluid.	Bites by infected rodents. Human to human transmission can not be excluded (for New World strains).

The presence of pathogen or infection may be identified by:

- PCR
- Serology (e. g. ELISA)

2. INTENDED USE

The Hantavirus IgM ELISA is intended for the qualitative determination of IgM antibodies against Hantavirus in human serum or plasma (citrate or heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Hantavirus Coated Microplate (IgM):** 12 breakapart 8-well snap-off strips coated with recombinant Hantavirus antigens in resealable aluminium foil.
- **IgM Sample Diluent:** 1 bottle containing 100 ml of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap.
- **Stop Solution:** 1 bottle containing 15 ml sulphuric acid, 0.2 mol/l; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 ml of a 20-fold concentrated phosphate buffer (0.2 M); pH 7.2 ± 0.2; for washing the wells; white cap.
- **Hantavirus anti-IgM Conjugate:** 1 bottle containing 20 ml of peroxidase labelled antibody to human IgM; in phosphate buffer (10 mM); coloured red, ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 ml 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1%; ready to use; yellow cap; < 5% NMP.
- **Hantavirus IgM Positive Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; red cap.
- **Hantavirus IgM Cut-off Control:** 1 bottle containing 3 ml control (human serum or plasma); coloured yellow; ready to use; green cap.
- **Hantavirus IgM Negative Control:** 1 bottle containing 2 ml control (human serum or plasma); coloured yellow; ready to use; blue cap.

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6.1. Coated Microplate

The break-apart snap-off strips are coated with recombinant Hantavirus antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Use human serum or plasma (citrate or heparin) samples with this assay. For CSF please use the instruction for use ABVL0001. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing.

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7.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with IgM Sample Diluent. Dispense 10 µl sample and 1 ml IgM Sample Diluent into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µl standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µl of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µl Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µl TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µl Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 99.03% (95% confidence interval: 94.71% - 99.98%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 96.3% (95% confidence interval: 89.56% - 99.23%).

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10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal significant evidence of false-positive results due to cross-reactions.

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13. ORDERING INFORMATION

Prod. No.: HANM0670 Hantavirus IgM ELISA (96 Determinations)

Serazym[®] Rotavirus

Enzyme immunoassay for detection of *Rotavirus* in faecal samples

REF E-020 ∇ 96 REF E-020-A2 ∇ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Group A Rotaviruses are the most common cause of non-bacterial gastroenteritis in children aged between 4 months and 3 years (1 - 5). Rotavirus is excreted into the intestine in large amounts (10^9 – 10^{11} virus particles per g faeces). Nosocomial infections therefore cause problems especially on baby wards and in childrens hospitals (3).

Rotavirus may also be responsible for travellers diarrhea in adults and have been detected in stool specimens of asymptomatic carriers as well (1). Rotavirus is spread by faecal-oral transmission from person to person or via contaminated staff. In temperate climates Rotavirus infections are mainly observed during the winter months (1).

Since virus culture on primary monkey kidney cells or permanent cell lines is time-consuming, these methods are of no diagnostic relevance. The golden standard is direct virus detection by electron microscopy (1, 2). Meanwhile antigen detection methods based on immunological techniques like agglutination tests or enzyme immunoassays using polyclonal or monoclonal antibodies against the group A specific antigen (VP-6) have been established (1 - 5).

References:

1. Böthig, B. und Diedrich, S. (1996): "Rotaviren" Diagnostische Bibliothek Band 1 Virusdiagnostik, Hrsg. Tomas Porstmann, Blackwell Wissenschafts-Verlag Berlin, Wien 1996, S. 441-451
2. Grauballe, B.F. et al. (1981): "Optimized Enzyme-Linked Immunosorbent Assay for Detection of Human and Bovine Rotavirus in Stools: Comparison with Electron-Microscopy, Immunoelectro-Osmophoresis and Fluorescent Antibody Techniques." *Journal of Medical Virology* 7: 29-40
3. Coulson, B.S. and I.H. Holmes (1984): "An Improved Enzyme-Linked Immunosorbent Assay For The Detection Of Rotavirus In Faeces Of Neonates." *Journal of Virological Methods*, 8: 165-179
4. Cukor G. et al. (1984): "Detection of Rotavirus in Human Stools by Using Monoclonal Antibody." *Journal of Clinical Microbiology* 19: 888-892
5. Cukor G. and N.R. Blacklow (1984): "Human Viral Gastroenteritis", *Microbiological Reviews* 48 No.2,p. 157-179.

Intended Use

The *Serazym*[®] Rotavirus is an *in-vitro*-diagnostic medical device for direct detection of Rotavirus in faecal samples.

Principle Of The Test

Serazym[®] Rotavirus is a one-step enzyme immunoassay on the basis of polyclonal antibodies to the group specific VP-6 antigen, the major protein of group A Rotaviruses. Diluted stool specimens and horseradish peroxidase (HRP) labelled polyclonal anti-Rotavirus-antibodies are dispensed simultaneously into the wells of a microtitration plate coated with polyclonal anti-Rotavirus antibodies. After an incubation time of 60 min at room temperature unbound components are removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-Tetramethylbenzidine (TMB) within a 10 min reaction time into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / 620 nm is directly proportional to the specifically bound amount of Rotavirus. Considering the cut-off value results are interpreted as positive or negative.

Test Components

			For 96 Wells	For 2x 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti-Rotavirus-antibodies (sheep)	12 single breakable 8-well strips colour coding dark blue vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding dark blue vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap	2x 100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control Rotavirus reactive sample	1.5 ml · ready to use coloured blue red cap	3.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control Rotavirus negative sample	1.5 ml · ready to use coloured blue green cap	3.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, polyclonal anti-Rotavirus-antibodies (rabbit)	12 ml · ready to use coloured green brown cap	24 ml · ready to use coloured green brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

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Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8 °C before testing in the ELISA.

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The *Serazym*[®] Rotavirus can be performed with 1 : 6 or 1 : 11 diluted specimens. In case of additional testing of the same sample in the *Serazym*[®] Norovirus, the *Serazym*[®] Campylobacter or the *Serazym*[®] Clostridium difficile Toxin A+B the 1 : 6 dilution is recommended.

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Preparation And Storage Of Reagents**Kit size and expiry**

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least 1 month when stored at 2...8°C.

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Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Dispense 2 drops (or 75 µl) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 75 µl **CONTROL +** positive control (4)
75 µl **CONTROL -** negative control (5)
50 µl **diluted sample**, mix gently.
4. Cover plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
6. Dispense 2 drops (or 75 µl) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 2 drops (or 75 µl) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

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Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for Rotavirus antigen.

Reference Values

Serazym® Rotavirus	
Positive	\geq Cut-off
Negative	$<$ Cut-off

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Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. A negative test result not necessarily excludes a *Rotavirus* infection. Inhomogeneous virus distribution in the sample can cause false negative results. The investigation of samples that were taken beyond the acute phase of the disease can cause false negative results, because the number of virus particles has decreased under the detection limit of the test. It is therefore recommended to take samples within the acute phase of the disease where a maximum number of excreted virus particles are to be expected. Faecal samples from vaccinated children may contain vaccine virus causing positive ELISA results. A final interpretation of the test results should consider clinical findings as well.

Automatic Processing

Performing the *Serazym*[®] Rotavirus on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary, the number of washing steps can be enhanced from 5x to 7x - 8x.

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A panel of 133 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.96$.

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*[®] Rotavirus from 12-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.841	0.137	7.45
2	1.208	0.078	6.50
3	0.620	0.040	6.38
4	0.463	0.024	5.26

Inter-assay coefficient of variation (CV) in the *Serazym*[®] Rotavirus in 10 different test runs from 3-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	2.720	0.128	4.71
2	1.647	0.122	7.38
3	0.968	0.074	7.69
4	0.409	0.019	4.58

Lower detection limit

The lower detection limit of *Rotavirus* antigen in the *Serazym*[®] Rotavirus was determined by titration of purified *Rotavirus* antigen SA-11: < 10 ng / ml corresponding to 10⁶ virus particles / g faeces.

Specificity and sensitivity

A total of 488 stool samples were investigated in parallel in the *Serazym*[®] Rotavirus and in another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym</i> [®] ELISA positive	246	0
<i>Serazym</i> [®] ELISA negative	4	238

Specificity: 100% · Sensitivity: 98.4%

Cross reactivity

Stool samples positive for one of the subsequent pathogens have been tested with the *Serazym*[®] Rotavirus and showed no cross reactivity:

Adenovirus (n = 20), Astrovirus (n = 8), Norovirus (n = 31), *Clostridium difficile* (n = 11), *Campylobacter jejuni* (n = 7), *Campylobacter coli* (n = 1), *Salmonella enteritidis* (n = 18), *Giardia lamblia* (n = 1) and stool samples (n = 93) with detectable levels (> 10 µg / g) of haemoglobin.

Negative stool specimens have been spiked with ≥ 10⁸ colony forming units of the following microorganisms and tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)	<i>Klebsiella pneumoniae</i>	(ATCC 13883)
<i>Bacillus cereus</i>	(ATCC 11778)	<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Bacillus subtilis</i>	(ATCC 6633)	<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Bacteroides fragilis</i>	(ATCC 25285)	<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Candida albicans</i>	(ATCC 10231)	<i>Salmonella enterica</i> Serovar <i>enteritidis</i>	(ATCC 13076)
<i>Campylobacter coli</i>	(ATCC 33559)	<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	(ATCC 14028)
<i>Campylobacter jejuni</i>	(ATCC 33291)	<i>Shigella flexneri</i>	(ATCC 12022)
<i>Citrobacter freundii</i>	(ATCC 8090)	<i>Shigella sonnei</i>	(ATCC 25931)
<i>Clostridium sordellii</i>	(ATCC 9714)	<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Enterobacter aerogenes</i>	(ATCC 13048)	<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Enterobacter cloacae</i>	(ATCC 13047)	<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Enterococcus faecalis</i>	(ATCC 29212)	<i>Vibrio cholerae</i>	clinical isolate
<i>Escherichia coli</i>	(ATCC 25922)	<i>Yersinia enterocolitica</i> Serotyp <i>O3, O9</i>	clinical isolates

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History of Changes

Version	Section	Modifications
2020-11-16	Common Advices and Precautions	Update
	History of Changes	Adding section "History of Changes "
	Test Components	Update

Incubation Scheme *Serazym*[®] Rotavirus (E-020)

1. 

2 drops (or 75 µl)

CONJ HRP (6)

+

pipette

75 µl **CONTROL +** (4)

75 µl **CONTROL -** (5)

50 µl **diluted stool sample**, mix gently

60 min incubation (room temperature)

 5 x Wash with wash solution
2. 

2 drops (or 75 µl)


SUBSTR TMB (7)

10 min incubation (room temperature) protected from light
3. 

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Product Number: LEGM0650 (96 Determinations)

ENGLISH

1. INTRODUCTION

Legionellae are aerobic gram-negative facultative intracellular parasites of certain protozoa. They are found in freshwater environments worldwide and can cause respiratory disease (legionellosis) in humans.

Legionella was first identified after an outbreak of pneumonia involving delegates of the 1976 American Legion Convention at a Philadelphia hotel.

The genus Legionella currently has at least 50 species comprising 70 distinct serogroups. One species of Legionella, *L. pneumophila*, is the aetiological agent of approximately 90 % of legionellosis cases, and serogroup 1 (Sg1) accounts for about 84 % of these cases.

L. pneumophila multiplies itself at temperatures between 25 and 42 °C, with an optimal growth temperature of 35 °C. Legionella thrives in warm, stagnant water in the environment and in artificial systems such as cooling towers, evaporative condensers, hot and cold water systems and spa pools that mimic the natural environment in which the organism thrives. These systems also provide the means by which aerosols/droplets are generated and the organism dispersed into the atmosphere.

Legionellosis can be acquired by the inhalation of aerosols containing Legionella bacteria or by micro-aspiration of ingested water contaminated with Legionella. Person-to-person transmission is not thought to be a risk.

The likelihood of contracting Legionnaires' disease depends on the level of contamination in the water source, the susceptibility of the person exposed, and the intensity of exposure. Legionnaires' disease is characterized as an "opportunistic" disease that attacks individuals who have an underlying illness or a weakened immune system. Predisposing risks include increasing age, being male, heavy smoking, alcohol abuse, chronic lung disease, immunosuppressive therapy, cancer chemotherapy, organ or bone marrow transplant, and corticosteroid therapy.

Legionellosis can appear in two distinct clinical presentations: Legionella pneumonia (Legionnaires' disease) with an incubation period of approx. 2-10 days (may extend up to 16-20 days) and Pontiac fever (incubation period: normally 12-48 hours).

Legionella pneumonia (Legionnaires' disease) is a serious form of pneumonia that carries with it a case-fatality ratio of 10-15 %. Legionnaires' disease patients initially present with cough, fever and nonspecific symptoms including malaise, myalgia and headache. Some patients develop shaking chills, chest pain, diarrhea, delirium or other neurologic symptoms. Extra pulmonary involvement is rare.

Pontiac fever is a milder form of the disease without manifestations of pneumonia and presents as an influenza-like illness. Symptoms may include headache, chills, muscle aches, a dry cough and fever. It is usually self-limiting and typically does not require treatment. The attack rate is much higher than for Legionnaires' disease (up to 95 % of those exposed).

Species	Disease	Symptoms (e.g.)	Transmission route
Legionella pneumophila	Legionella pneumonia (Legionnaires' disease)	Cough, fever and nonspecific symptoms (malaise, myalgia, headache). Some patients develop shaking chills, chest pain, diarrhea, delirium or other neurologic symptoms.	Inhalation of aerosols containing Legionella bacteria or micro-aspiration of ingested water contaminated with Legionella
	Pontiac fever	Influenza-like illness (headache, chills, muscle aches, a dry cough and fever) without manifestations of pneumonia	

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- Culture
- Urinary antigen detection
- PCR
- Serology: Detection of antibodies by IF, ELISA

2. INTENDED USE

The Legionella pneumophila IgM ELISA is intended for the qualitative determination of IgM class antibodies against Legionella pneumophila in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with Legionella pneumophila antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microplate

The break-apart snap-off strips are coated with Legionella pneumophila antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95.65% (95% confidence interval: 85.16% - 99.47%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 66.37% - 100%).

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10.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters did not reveal significant evidence of false-positive results due to cross-reactions.

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Product Number: LEGG0650 (96 Determinations)

ENGLISH

1. INTRODUCTION

Legionellae are aerobic gram-negative facultative intracellular parasites of certain protozoa. They are found in freshwater environments worldwide and can cause respiratory disease (legionellosis) in humans.

Legionella was first identified after an outbreak of pneumonia involving delegates of the 1976 American Legion Convention at a Philadelphia hotel.

The genus Legionella currently has at least 50 species comprising 70 distinct serogroups. One species of Legionella, *L. pneumophila*, is the aetiological agent of approximately 90 % of legionellosis cases, and serogroup 1 (Sg1) accounts for about 84 % of these cases.

L. pneumophila multiplies itself at temperatures between 25 and 42 °C, with an optimal growth temperature of 35 °C. Legionella thrives in warm, stagnant water in the environment and in artificial systems such as cooling towers, evaporative condensers, hot and cold water systems and spa pools that mimic the natural environment in which the organism thrives. These systems also provide the means by which aerosols/droplets are generated and the organism dispersed into the atmosphere.

Legionellosis can be acquired by the inhalation of aerosols containing Legionella bacteria or by micro-aspiration of ingested water contaminated with Legionella. Person-to-person transmission is not thought to be a risk.

The likelihood of contracting Legionnaires' disease depends on the level of contamination in the water source, the susceptibility of the person exposed, and the intensity of exposure. Legionnaires' disease is characterized as an "opportunistic" disease that attacks individuals who have an underlying illness or a weakened immune system. Predisposing risks include increasing age, being male, heavy smoking, alcohol abuse, chronic lung disease, immunosuppressive therapy, cancer chemotherapy, organ or bone marrow transplant, and corticosteroid therapy.

Legionellosis can appear in two distinct clinical presentations: Legionella pneumonia (Legionnaires' disease) with an incubation period of approx. 2-10 days (may extend up to 16-20 days) and Pontiac fever (incubation period: normally 12-48 hours).

Legionella pneumonia (Legionnaires' disease) is a serious form of pneumonia that carries with it a case-fatality ratio of 10-15 %. Legionnaires' disease patients initially present with cough, fever and nonspecific symptoms including malaise, myalgia and headache. Some patients develop shaking chills, chest pain, diarrhea, delirium or other neurologic symptoms. Extra pulmonary involvement is rare.

Pontiac fever is a milder form of the disease without manifestations of pneumonia and presents as an influenza-like illness. Symptoms may include headache, chills, muscle aches, a dry cough and fever. It is usually self-limiting and typically does not require treatment. The attack rate is much higher than for Legionnaires' disease (up to 95 % of those exposed).

Species	Disease	Symptoms (e.g.)	Transmission route
Legionella pneumophila	Legionella pneumonia (Legionnaires' disease)	Cough, fever and nonspecific symptoms (malaise, myalgia, headache). Some patients develop shaking chills, chest pain, diarrhea, delirium or other neurologic symptoms	Inhalation of aerosols containing Legionella bacteria or micro-aspiration of ingested water contaminated with Legionella
	Pontiac fever	Influenza-like illness (headache, chills, muscle aches, a dry cough and fever) without manifestations of pneumonia	

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- Culture
- Urinary antigen detection
- PCR
- Serology: Detection of antibodies by ELISA

2. INTENDED USE

The Legionella pneumophila IgG ELISA is intended for the qualitative determination of IgG class antibodies against Legionella pneumophila in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Legionella pneumophila* antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Legionella pneumophila* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 88.78% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 90.0% (95% confidence interval: 68.3% - 98.77%).

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Species	Disease	Symptoms (e.g.)	Transmission route
<i>Leptospira</i> spp.	Leptospirosis	Fever accompanied by chills, intense headache, severe myalgia (muscle ache), abdominal pain, conjunctival suffusion (red eye), and occasionally a skin rash	Direct or indirect contact with the urine of an infected animal (via cuts, skin abrasions or mucous membranes)
	Weil's disease	jaundice, renal failure, meningitis, haemorrhage and myocarditis with arrhythmias	

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- Microscopy
- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The *Leptospira* IgM ELISA is intended for the qualitative determination of IgM class antibodies against *Leptospira* spp. in human serum or plasma (heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Leptospira* antigens; in resealable aluminium foil.
- **IgM Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent); coloured green; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM); coloured red; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 96.0% (95% confidence interval: 79.65% - 99.9%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 93.51% - 100%).

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10.5. Cross Reactivity

It cannot be excluded that Cytomegalovirus, Treponema pallidum and Coxiella specimens may result in false-positive IgM antibody results. In addition, it should be noted that IgM class antibodies directed against Leptospira generally remain detectable for months or even years but at low titer.

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Product Number: LEPG0660 (96 Determinations)

ENGLISH

1. INTRODUCTION

Leptospirosis (also known as Weil's syndrome) is probably the most widespread zoonosis in the world. It is caused by infection with spirochete bacteria of the genus *Leptospira* and affects humans as well as a broad spectrum of animal hosts. The incidence is significantly higher in warm climate countries than in temperate regions. The disease is seasonal, with peak incidence occurring in summer or fall in temperate regions, where temperature is the limiting factor in survival of leptospires, and during rainy seasons in warm climate regions, where rapid desiccation would otherwise prevent survival.

Natural reservoirs for the pathogenic *Leptospira interrogans* include rodents as well as a large variety of domesticated mammals (e. g. pigs, cattle and dogs). Leptospires occupy the lumen of nephritic tubules in their natural host and are shed into the urine.

Transmission can occur when humans are directly or indirectly exposed to the urine of infected animals or a urine-polluted environment. Leptospires gain entry into the human blood stream via cuts, skin abrasions or mucous membranes through contact with moist soil, vegetation, and contaminated waters; handling infected animal tissues; and ingestion of food and water. Leptospires are rarely transmitted from human to human.

The incubation period is usually 5-14 days, with a range of 2-30 days.

The spectrum of clinical symptoms is extremely wide. The vast majority of leptospiral infections are either subclinical or result in very mild illness and recover without any complications. Clinical manifestations of leptospirosis range from mild influenza-like symptoms to severe life-threatening disease forms, characterized by jaundice, renal failure, bleeding and severe pulmonary hemorrhage.

The clinical presentation of leptospirosis is biphasic, with the acute or septicemic phase lasting about a week, followed by the immune phase, characterized by antibody production and excretion of leptospires in the urine. Most of the complications of leptospirosis are associated with localization of leptospires within the tissues during the immune phase and thus occur during the second week of the illness. The classical syndrome of Weil's disease represents only the most severe presentation. It is characterized by jaundice, renal failure, hemorrhage and myocarditis with arrhythmias.

Species	Disease	Symptoms (e.g.)	Transmission route
Leptospira spp.	Leptospirosis	Fever accompanied by chills, intense headache, severe myalgia (muscle ache), abdominal pain, conjunctival suffusion (red eye), and occasionally a skin rash	Direct or indirect contact with the urine of an infected animal (via cuts, skin abrasions or mucous membranes)
	Weil's disease	jaundice, renal failure, meningitis, haemorrhage and myocarditis with arrhythmias	

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- Microscopy
- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The *Leptospira* IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Leptospira* spp. in human serum or plasma (heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Leptospira* antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Leptospira* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Use human serum or plasma (heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 97.37% (95% confidence interval: 86.19% - 99.93%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 93,4% - 100%).

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10.5. Cross Reactivity

It cannot be excluded that polyclonal B-cell activation induced by Epstein-Barr virus (EBV) or the presence of Rheumatoid Factors may result in false-positive Leptospira IgG antibody results.

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Only for in-vitro diagnostic use

Instructions for use

English	2
Abbreviations	10
Symbols Key	11
Summary of Test Procedure	12

REF

YERG0990 (96 Determinations)

ENGLISH

1. INTENDED USE

The *Yersinia enterocolitica* IgG ELISA is intended for the qualitative determination of IgG class antibodies against antigens of the 70 kb virulence plasmid of *Yersinia enterocolitica* in human serum or plasma (citrate, heparin). The ELISA is intended for use as an aid in identifying individuals with an adaptive immune response to *Yersinia enterocolitica*. The determination of increased antibody levels contributes to diagnosis of *Yersinia*-induced reactive arthritis. The test is not intended for diagnosing acute enteric diseases.

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3. MATERIALS

3.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with specific antigen; in resealable aluminium foil.
- **DIL:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02 % (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02 % (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).

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7. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems, we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

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Example: $\frac{1.591 \times 10}{0.43} = 37$ NTU (Units)

8.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks.
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

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9.1.1. Single-Site Study

The precision study was performed at a single site. A negative, a high negative, a low positive, and a moderate positive sample were run in 4 replicates, two times per day for 12 days for a total of 96 results.

Sample	n	Mean (NTU)	Repeatability		Between Run		Within Day		Between Day		Within lab	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
moderate positive	96	21.50	0.8883	4.1	1.6690	7.8	1.8907	8.8	0.9469	4.4	2.1145	9.8
low positive	96	15.28	0.6923	4.5	1.0913	7.1	1.2924	8.5	0.0000	0.0	1.2924	8.5
high negative	96	8.56	0.5784	6.8	0.7467	8.7	0.9445	11.0	0.2395	2.8	0.9744	11.4
negative	96	4.60	0.4439	9.7	0.3461	7.5	0.5629	12.2	0.1438	3.1	0.5810	12.6

9.1.2. Multisite-Study

The precision study was performed at two different sites. A negative, a high negative, a low positive, and a moderate positive sample were run in 5 replicates, once a day for 5 days for a total of 50 results.

Sample	n	Mean (NTU)	Repeatability		Within Site		Reproducibility	
			SD	%CV	SD	%CV	SD	%CV
moderate positive	50	21.68	0.6338	2.9	1.4163	6.5	2.1946	10.1
low positive	50	14.58	0.5295	3.6	1.1098	7.6	1.2922	8.9
high negative	50	8.60	0.4094	4.8	0.9157	10.6	1.2857	14.9
negative	50	3.93	0.3859	9.8	0.7162	18.2	0.9696	24.7

9.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 98.82 % (95 % confidence interval: 95.81 % - 99.86 %).

9.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 95.65 % (95 % confidence interval: 85.16 % - 99.47 %).

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9.5. Cross Reactivity

A minimum of 5 samples with antibody activities to potentially cross-reacting parameters (Adenovirus, Borrelia burgdorferi, Brucella, Campylobacter jejuni, Chlamydia pneumoniae, Enterovirus, Epstein-Barr virus, Helicobacter pylori, Parvovirus B19) or samples positive for ANA or rheumatoid factors, and samples from pregnant women were tested to evaluate the cross reactivity of the assay. Positive findings were additionally analyzed with a CE-marked reference assay. The results are shown in the following table.

Pathogen/Condition	Samples tested	Number of positive samples	
		Yersinia enterocolitica IgG	CE-marked reference assay
Adenovirus	15	2	4
Antinuclear antibodies (ANA)	15	3	2
Borrelia burgdorferi	15	4	4
Brucella	10	9	9
Campylobacter jejuni	8	4	5
Chlamydia pneumoniae	12	4	2
Enterovirus	15	7	5
Epstein-Barr virus (EBV)	15	3	3
Helicobacter pylori	14	6	6
Parvovirus B19	14	6	1
Pregnancy samples	14	1	0
Rheumatoid factor (RF)	13	4	5

Cross-reactions with antibodies against Adenovirus, Borrelia burgdorferi, Brucella, Campylobacter jejuni, Chlamydia pneumoniae, Enterovirus, Epstein-Barr virus, Helicobacter pylori, Parvovirus B19 as well as with samples from pregnant women or samples positive for antinuclear antibodies (ANA) or rheumatoid factor (RF) cannot be excluded.

Based on the prevalence values of 4-35 % for anti-Yersinia enterocolitica IgG antibodies reported in the literature, it cannot be excluded that some of the tested samples are correctly positive^{1,2}.

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Only for in-vitro diagnostic use

Instructions for use

English	2
Abbreviations	10
Symbols Key	11
Summary of Test Procedure.....	12

REF

YERA0990 (96 Determinations)

ENGLISH

1. INTENDED USE

The *Yersinia enterocolitica* IgA ELISA is intended for the qualitative determination of IgA class antibodies against antigens of the 70 kb virulence plasmid of *Yersinia enterocolitica* in human serum or plasma (citrate, heparin). The ELISA is intended for use as an aid in identifying individuals with an adaptive immune response to *Yersinia enterocolitica*. The determination of increased antibody levels contributes to diagnosis of *Yersinia*-induced reactive arthritis. The test is not intended for diagnosing acute enteric diseases.

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3. MATERIALS

3.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with specific antigen; in resealable aluminium foil.
- **DIL:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap; 0.2 % (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgA in phosphate buffer (10 mM); coloured violet; ready to use; black cap.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02 % (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02 % (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015 % (v/v) CMIT/MIT (3:1).

For hazard and precautionary statements see 11.1

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5.1. Microtiterplate

The break-apart snap-off strips are coated with specific antigen. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

5.2. WASH BUF 20x

Dilute WASH BUF 20x 1 + 19; e.g. 10 mL WASH BUF 20x + 190 mL distilled water. The diluted buffer (WASH BUF 1x) is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g., in a water bath. Mix well before dilution.

5.3. SUB TMB

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If SUB TMB turns into blue, it may have become contaminated and should be thrown away.

6. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

6.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with DIL. Dispense 10 µL sample and 1 mL DIL into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

7. ASSAY PROCEDURE

Please read the instruction for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems, we recommend increasing the washing steps from three up to five and the volume of WASH BUF 1x from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour ± 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of WASH BUF 1x. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL SUB TMB into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL SOLN STOP into all wells in the same order and at the same rate as for the SUB TMB, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the SOLN STOP.

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8.2.1. Results in Units [NTU]

Sample (mean) absorbance value x 10 / Cut-off = [NovaTec Units = NTU]

Example: $\frac{1.591 \times 10}{0.43} = 37$ NTU (Units)

8.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks.
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.
Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data. In immunocompromised patients and newborns serological data only have restricted value.		

8.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection
IgA	Produced in mucosal linings throughout the body (⇒ protective barrier) Usually produced early in the course of the infection

9. SPECIFIC PERFORMANCE CHARACTERISTICS

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9.1. Precision

Evaluation of precision of the assay was performed according to "CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline - Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014".

9.1.1. Single-Site Study

The precision study was performed at a single site. A negative, a high negative, an equivocal, and a positive sample were run in 4 replicates, two times per day for 12 days for a total of 96 results.

Sample	n	Mean (NTU)	Repeatability		Between Run		Within Day		Between Day		Within lab	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
positive	96	25.03	1.3174	5.3	1.3537	5.4	1.8889	7.5	0.8140	3.3	2.0568	8.2
equivocal	96	10.01	0.7302	7.3	0.7926	7.9	1.0777	10.8	0.4165	4.2	1.1554	11.5
high negative	96	6.81	0.6429	9.4	0.0912	1.3	0.6494	9.5	0.3696	5.4	0.7472	11.0
negative	96	2.82	0.4468	15.8	0.1387	4.9	0.4678	16.6	0.2013	7.1	0.5093	18.0

9.1.2. Multisite-Study

The precision study was performed at three different sites. A negative, a high negative, an equivocal, and a moderate positive sample were run in 5 replicates, once a day for 5 days for a total of 75 results.

Sample	n	Mean (NTU)	Repeatability		Within Site		Reproducibility	
			SD	CV	SD	CV	SD	CV
moderate positive	75	23.10	1.0818	4.7	1.8064	7.8	2.0543	8.9
equivocal	75	9.98	0.6410	6.4	0.9572	9.6	1.3167	13.2
high negative	75	6.13	0.4330	7.1	0.6231	10.2	0.6688	10.9
negative	75	2.14	0.2894	13.5	0.4488	20.9	0.6611	30.8

9.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 97.32 % (95 % confidence interval: 92.37 % - 99.44 %).

9.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 93.94 % (95 % confidence interval: 79.77 % - 99.26 %).

9.4. Interferences

The assay was evaluated for interferences according to guideline EP07-A3 ("Interference Testing in Clinical Chemistry" from the Clinical and Laboratory Standards Institute). Three samples, covering the relevant measuring range, were spiked with high levels of interferents and were tested along with the unspiked sample. The following table shows the tested substances added to patient samples at the indicated concentrations. These correspond to the recommendations in the CLSI guideline to represent pathological elevated concentrations in patient samples.

Interferent	Concentration tested
Albumin	60 mg/mL
Bilirubin, unconjugated	0.4 mg/mL
Bilirubin, conjugated	0.4 mg/mL
Cholesterol	4 mg/mL
Hemoglobin	10 mg/mL
Triglycerides	15 mg/mL

No clinically significant interference effect was found for all tested substances.

9.5. Cross Reactivity

A minimum of 5 samples with antibody activities to potentially cross-reacting parameters (Adenovirus, Borrelia burgdorferi, Brucella, Campylobacter jejuni, Chlamydia pneumoniae, Enterovirus, Epstein-Barr virus, Helicobacter pylori, Parvovirus B19, Salmonella typhi) or samples positive for ANA or rheumatoid factors, and samples from pregnant women were tested to evaluate the cross reactivity of the assay. Positive findings were additionally analyzed with a CE-marked reference assay. The results are shown in the following table.

Pathogen/Condition	Samples tested	Number of positive samples	
		Yersinia enterocolitica IgA	CE-marked reference assay
Adenovirus	15	3	0
Antinuclear antibodies (ANA)	15	5	1
Borrelia burgdorferi	10	1	1
Brucella	10	6	5
Campylobacter jejuni	9	1	0
Chlamydia pneumoniae	12	4	2
Enterovirus	7	2	2
Epstein-Barr virus (EBV)	13	4	0
Helicobacter pylori	13	2	1
Pregnancy samples	13	2	0
Rheumatoid factor (RF)	11	3	0
Salmonella typhi	9	3	1

Cross-reactions with antibodies against Borrelia burgdorferi, Chlamydia pneumoniae, Epstein-Barr virus, Helicobacter pylori, Enterovirus, Adenovirus, Campylobacter jejuni, Salmonella typhi, Brucella as well as with samples from pregnant women or samples positive for rheumatoid factor (RF) or antinuclear antibodies (ANA) cannot be excluded.

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11.1. Safety note for reagents containing hazardous substances

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Reagents may contain 5-Bromo-5-nitro-1,3-dioxane (refer to 3.1)

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Warning

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.

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ENGLISH

1. INTRODUCTION

Taenia solium is a tapeworm of 2-7 m in length which resides in the small intestine of humans but also other animal species (monkeys, hamsters). The tapeworms produce proglottids (less than 1,000, and each with 50,000 eggs) which mature, become gravid, detach from the tapeworm, and migrate to the anus or are passed in the stool. The eggs contained in the gravid proglottids and passed with the faeces can survive for months to years in the environment. After ingestion of a suitable intermediate host (pigs and other animals, including humans) the eggs release the oncosphere, invade the intestinal wall and migrate to the striated muscles, into the brain, liver and other tissues of the host where they develop in cysticerci. In the human intestine, a cysticercus develops over 2 months into an adult tapeworm, which can survive for up to 25 years. The important parasitic infection caused by *Taenia solium* is cysticercosis which may involve the eye and the central nervous system. The swine tapeworm *Taenia solium* is worldwide in distribution. Prevalence is higher in poorer communities where humans live in close contact with pigs and eat undercooked pork, and is very rare in Muslim countries. The main symptom of Taeniasis (only mild) is often the passage (passive) of proglottids.

The most important feature of *Taeniasis solium* is the risk of development of Cysticercosis.

Species	Disease	Symptoms (e.g.)	Transmission route
<i>Taenia solium</i>	Taeniasis Cysticercosis (Neurocysticercosis)	Abdominal pain; Nausea; Weakness and fatigue; Weight loss; Flatulence (gases); Diarrhea or constipation; Appetite changes (too much hunger or loss of appetite) Cysticerci in the brain may cause increased cranial pressure, convulsions and altered mental states	Ingestion of undercooked pork meat containing cysticerci or ingestion of <i>Taenia solium</i> eggs via fecally contaminated food or water

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- Histology
- Microscopy
- PCR
- Serology: e.g. ELISA

2. INTENDED USE

The *Taenia solium* IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Taenia solium* in human serum or plasma (citrate, heparin).

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4. MATERIALS

4.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Taenia solium* antigens; in resealable aluminium foil.
- **IgG Sample Dilution Buffer:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **Stop Solution:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **Washing Buffer (20x conc.):** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap.
- **Conjugate:** 1 bottle containing 20 mL of peroxidased Protein A in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **TMB Substrate Solution:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

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4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplate
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

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6.1. Microtiterplate

The break-apart snap-off strips are coated with *Taenia solium* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

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Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 1 hour \pm 5 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Protein A conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

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10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100% (95% confidence interval: 96.67% - 100%).

10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 78.2% - 100%).

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10.5. Cross Reactivity

Cross reaction of the antigens with antibodies against Echinococcus and Entamoeba is possible.

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13. ORDERING INFORMATION

Prod. No.: TAEG0420 Taenia solium IgG ELISA (96 Determinations)

	ZAO "Vector-Best"	Rev. 01
	EC Declaration of conformity	Page 1 of 4

EC DECLARATION OF CONFORMITY

ZAO "Vector-Best" hereby ensures under own responsibility and declares that the products listed on pages 2-4 are in conformity with applicable provisions and fulfill the essential requirements of Annex I Directive 98/79/EC of 27 October 1998 regarding in vitro diagnostic medical devices.

Classification of products:

Other devices (all devices except Annex II and self-testing devices)

Conformity assessment procedure:

Annex III (not including section 6).

Manufacturer:

ZAO "Vector-Best"
 Address: AHC, Koltsovo,
 Novosibirsk Region, 630559, Russia,
 Tel. +7 (383) 363 20 60,
 Fax: +7 (383) 363 35 55

European authorized representative:

Bioron GmbH,
 Rheinhorststr. 18, D-67071
 Ludwigshafen, Germany.
 tel.: +49 (0) 621 5720 915,
 fax: +49 (0) 621 5720 916

Date: 2013/04/12

Murat Khusainov
 General Director ZAO «Vector-Best»

	ZAO "Vector-Best"	Rev. 01
	EC Declaration of conformity	Page 2 of 4

No.	Product name	Identification data	REF
1.	Vectohep A-IgM	ELISA kit for determination of IgM to hepatitis A virus	D-0352
2.	Vectohep A-IgG	ELISA kit for quantitative and qualitative determination of IgG to hepatitis A virus	D-0362
3.	Vectohep TTV-IgG	ELISA kit for determination of IgG to TT virus	D-0802
4.	Vectohep E-IgG	ELISA kit for determination of IgG to hepatitis E virus	D-1056
5.	Vectohep E-IgM	ELISA kit for determination of IgM to hepatitis E virus	D-1058
6.	Vectohep G-IgG	ELISA kit for determination of IgG to hepatitis G virus	D-1252
7.	LymeBest-IgG	ELISA kit for determination of IgG to infectious borreliosis agents	D-1452
8.	LymeBest-IgM	ELISA kit for determination of IgM to infectious borreliosis agents	D-1454
9.	RecombiBest antipallidum-IgG	ELISA kit for determination of IgG to Treponema pallidum	D-1852
10.	RecombiBest antipallidum-total antibodies	ELISA kit for determination of total antibodies to Treponema pallidum	D-1856
11.	RecombiBest antipallidum-IgM	ELISA kit for determination of IgM to Treponema pallidum	D-1858
12.	RecombiBest antipallidum-total antibodies	ELISA kit for determination of total antibodies to Treponema pallidum	D-1857
13.	VectoHSV-1,2 - IgG	ELISA kit for determination of IgG to herpes simplex virus types 1 and 2	D-2152
14.	VectoHSV - IgM	ELISA kit for determination of IgM to herpes simplex virus types 1 and 2	D-2154
15.	VectoHHV-8 - IgG	ELISA kit for determination of IgG to human herpes virus type 8	D-2160
16.	VectoHHV-6 - IgG	ELISA kit for determination of IgG to human herpes virus type 6	D-2166
17.	Ureaplasma urealyticum – IgG-EIA-BEST	ELISA kit for determination of IgG to Ureaplasma urealyticum antigens	D-2254
18.	Ureaplasma urealyticum – IgA-EIA-BEST	ELISA kit for determination of IgA to Ureaplasma urealyticum antigens	D-2258
19.	VectoParotitis-IgG	ELISA kit for determination of IgG to parotitis virus	D-2602
20.	VectoParotitis-IgM	ELISA kit for determination of IgM to parotitis virus	D-2604
21.	Toxocara-IgG-EIA-BEST	ELISA kit for determination of IgG to toxocara antigens	D-2752
22.	Opisthorchiasis – IgG-EIA-BEST	ELISA kit for determination of IgG to opisthorchiasis antigens	D-2952
23.	Echinococcus-IgG-EIA-BEST	ELISA kit for determination of IgG to Echinococcus	D-3356

		antigens	
24.	Ascarid-IgG-EIA-BEST	ELISA kit for determination of IgG to Ascaris lumbricoides	D-3452
25.	Lambliia-antibodies-EIA-BEST	ELISA kit for determination of IgG, IgM and IgA to Lambliia antibodies	D-3552
26.	Lambliia-IgM-EIA-BEST	ELISA kit for determination of IgM to Lambliia antibodies	D-3554
27.	Lambliia-antigen-EIA-BEST	ELISA kit for determination of Lambliia antigen	D-3556
28.	Helicobacter pylori-CagA-antigen-EIA-BEST	ELISA kit for determination of total antibodies to CagA Helicobacter pylori	D-3752
29.	TSH-EIA-BEST	ELISA kit for determination of concentration of thyroid-stimulating hormone	X-3952
30.	T3 total-EIA-BEST	ELISA kit for determination of concentration of total triiodothyronine	X-3954
31.	T4 total-EIA-BEST	ELISA kit for determination of concentration of total thyroxine	X-3956
32.	Anti-TPO-EIA-BEST	ELISA kit for determination of antibody concentration to thyroperoxidase	X-3968
33.	PAPP-A-EIA-BEST	ELISA kit for determination of concentration of pregnancy-associated plasma protein A	D-4160
34.	Mycoplasma hominis-IgG-EIA-BEST	ELISA kit for determination of IgG to Mycoplasma hominis	D-4352
35.	Mycoplasma hominis-IgA-EIA-BEST	ELISA kit for determination of IgA to Mycoplasma hominis	D-4358
36.	Mycoplasma pneumoniae-IgG-EIA-BEST	ELISA kit for determination of IgG to Mycoplasma pneumoniae	D-4362
37.	Mycoplasma pneumoniae-IgM-EIA-BEST	ELISA kit for determination of IgM to Mycoplasma pneumoniae	D-4366
38.	Vectocrimean – CHF – IgG	ELISA kit for determination of IgG to Crimean-Congo hemorrhagic fever virus	D-5052
39.	Vectocrimean – CHF – IgM	ELISA kit for determination of IgM to Crimean-Congo hemorrhagic fever virus	D-5054
40.	CEA-EIA-BEST	ELISA kit for determination of concentration of carcinoembryonic antigen	T-8454
41.	AFP-EIA-BEST	ELISA kit for determination of concentration of Alpha-Fetal Protein	T-8456
42.	CA-125-EIA-BEST	ELISA kit for determination of concentration of oncomarker CA-125	T-8466
43.	CA 19-9-EIA-BEST	ELISA kit for determination of concentration of CA 19-9	T-8470
44.	CA 15-3-EIA-BEST	ELISA kit for determination of concentration of oncomarker CA 15-3	T-8472
45.	NSE-EIA-BEST	ELISA kit for determination of concentration of neuron specific enolase	T-8476

	ZAO "Vector-Best"	Rev. 01
	EC Declaration of conformity	Page 4 of 4

46.	Ferritin-EIA-BEST	ELISA kit for determination of concentration of ferritin	T-8552
47.	IgE total-EIA-BEST	ELISA kit for determination of concentration of total IgE	A-8660
48.	IgG total-EIA-BEST	ELISA kit for determination of concentration of total IgG	A-8662
49.	IgM total-EIA-BEST	ELISA kit for determination of concentration of total IgM	A-8664
50.	IgA total-EIA-BEST	ELISA kit for determination of concentration of total IgA	A-8666
51.	Gamma-Interferon-EIA-BEST	ELISA kit for determination of concentration of gamma-interferon	A-8752
52.	Interleukine-4-EIA-BEST	ELISA kit for determination of concentration of Interleukine-4	A-8754
53.	Alpha-TNF-EIA-BEST	ELISA kit for determination of concentration of alpha-tumor necrosis factor	A-8756
54.	Alpha-Interferon-EIA-BEST	ELISA kit for determination of concentration of alpha-interferon	A-8758
55.	Interleukine-6-EIA-BEST	ELISA kit for determination of concentration of Interleukine-6	A-8768
56.	Interleukine-2-EIA-BEST	ELISA kit for determination of concentration of Interleukine-2	A-8772
57.	Procalcitonin-EIA-BEST	ELISA kit for determination of concentration of procalcitonin	A-9004
58.	NTproBNP-EIA-BEST	ELISA kit for determination of concentration of N-terminal prohormone of brain natriuretic peptide	A-9102
59.	Troponin I-EIA-BEST	ELISA kit for determination of concentration of troponin I	A-9106

Certificate

mdc medical device certification GmbH
certifies that

VECTOR

**AO Vector-Best
Research and Production Area
Building 36, Office 211, Koltsovo
630559 Novosibirsk region
Russian Federation**

with the locations listed in the attachment
for the scope

**Design and development, production and distribution of
medical devices for in vitro diagnostics (PCR, ELISA, Biochemistry)**

has introduced and applies a

Quality Management System

The mdc audit has proven that this quality management system
meets all requirements of the following standard

EN ISO 13485

Medical devices – Quality management systems –
Requirements for regulatory purposes

EN ISO 13485:2016 + AC:2016 - ISO 13485:2016

Valid from	2020-07-04
Valid until	2023-07-03
Registration no.	D1213100019
Report no.	P20-00568-173687
Stuttgart	2020-06-02

Head of Certification Body

Attachment of the certificate

No. D1213100019

date 2020-06-02

Page 1 of 1

Location	Scope
AO Vector-Best Arbuzova str. 1/1, 630117 Novosibirsk Russian Federation	design and development, production and distribution of medical devices for in vitro diagnostics
AO Vector-Best Research and Production area, building 36, Koltsovo, 630559 Novosibirsk region Russian Federation	design and development, production of medical devices for in vitro diagnostics
AO Vector-Best Pasechnaya str, 3, 630117 Novosibirsk Russian Federation	design and development, production of medical devices for in vitro diagnostics

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Scope: Design and development, manufacture, installation, service and distribution of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents, and instruments / software for in vitro diagnostics

TÜVRheinland

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the abovementioned standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.: 3313978-90

Effective date: 2020-05-19

Expiry date: 2023-05-18

Issue date: 2020-05-14

D. Wiedemuth

Dipl.-Ing. (FH) D. Wiedemuth
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany

Certificate

Quality Management System EN ISO 13485:2016

Registration No.: SX 1483000-1

Organization: EUROIMMUN
Medizinische Labordiagnostika AG
Seekamp 31
23560 Lübeck
Germany

No.	Facility	Scope
/01	EUROIMMUN Medizinische Labordiagnostika AG Seekamp 31 23560 Lübeck Germany	Design and development and manufacture of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents, and instruments / software for instruments for in vitro diagnostics
/02	EUROIMMUN Medizinische Labordiagnostika AG Werkstr. 1 23942 Dassow Germany	Design, development, manufacture and distribution of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents and instruments / software for in vitro diagnostics
/03	EUROIMMUN Medizinische Labordiagnostika AG An der Trave 1 23923 Selmsdorf Germany	Design, development, manufacture, service and distribution of immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents and instruments / software for in vitro diagnostics

Report No.: 3313978-90

Effective date: 2020-05-19

Expiry date: 2023-05-18

Issue date: 2020-05-14

D. Wiedemuth

Dipl.-Ing. (FH) D. Wiedemuth
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany

Certificate

Quality Management System EN ISO 13485:2016

Registration No.: SX 1483000-1

Organization: EUROIMMUN
Medizinische Labordiagnostika AG
Seekamp 31
23560 Lübeck
Germany

- | | | |
|-----|--|--|
| /04 | EUROIMMUN
Medizinische Labordiagnostika AG
Am Sonnenberg 9
23627 Groß Grönau
Germany | Design, development and manufacture of immunofluorescence test systems, molecular diagnostic / genetic test systems and, test systems for the determination of infectious agents for in vitro diagnostics |
| /05 | EUROIMMUN
Medizinische Labordiagnostika AG
Am Born 24
23627 Groß Grönau
Germany | Design and development of software for in vitro diagnostics |
| /06 | EUROIMMUN
Medizinische Labordiagnostika AG
Im Kreppel 1
02747 Herrnhut
Germany | Manufacture of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems and test systems for the determination of infectious agents for in vitro diagnostics |
| /07 | EUROIMMUN
Medizinische Labordiagnostika AG
Am Pließnitztal 1
02748 Bernstadt
Germany | Manufacture of immunobiochemical test systems, test systems for the determination of infectious agents and instruments for in vitro diagnostics |

Report No.: 3313978-90

Effective date: 2020-05-19

Expiry date: 2023-05-18

Issue date: 2020-05-14

D. Wiedemuth

Dipl.-Ing. (FH) D. Wiedemuth
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany

Certificate

Quality Management System EN ISO 13485:2016

Registration No.: SX 1483000-1

Organization: EUROIMMUN
Medizinische Labordiagnostika AG
Seekamp 31
23560 Lübeck
Germany

- | | | |
|-----|--|--|
| /08 | EUROIMMUN
Medizinische Labordiagnostika AG
Schloßstr. 11
91257 Pegnitz
Germany | Manufacture of immunofluorescence test systems, installation and service of instruments / software for in vitro diagnostic |
| /09 | EUROIMMUN
Medizinische Labordiagnostika AG
Am Flugplatz 4
23560 Lübeck
Germany | Installation, service and distribution of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents, and instruments / software for in vitro diagnostics |
| /10 | EUROIMMUN
Medizinische Labordiagnostika AG
Gewerbestr. 19
23942 Dassow
Germany | Manufacture of sheet metal and other components for instruments for in vitro diagnostics |

Report No.: 3313978-90

Effective date: 2020-05-19

Expiry date: 2023-05-18

Issue date: 2020-05-14

D. Wiedemuth

Dipl.-Ing. (FH) D. Wiedemuth
TÜV Rheinland LGA Products GmbH
Tillystraße 2 · 90431 Nürnberg · Germany

Certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810000**

Certificate Holder: **EUROIMMUN
Medizinische Labordiagnostika AG**
Seekamp 31
23560 Lübeck
Germany

including the locations according to annex

Scope: Design, development, manufacture, installation, service and sales of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents, and instruments / software for in vitro diagnostics in humans and animals; trainings

Proof has been furnished by means of an audit that the requirements of ISO 9001:2015 are met.

Validity: The certificate is valid from 2020-05-19 until 2023-05-18.
First certification 2018

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Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810000**

No.	Location	Scope
/01	c/o EUROIMMUN Medizinische Labordiagnostika AG Seekamp 31 23560 Lübeck Germany	Design, development, manufacture, installation, service and sales of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents, and instruments / software for in vitro diagnostics in humans and animals; trainings
/02	c/o EUROIMMUN Medizinische Labordiagnostika AG Werkstr. 1 23942 Dassow Germany	Design, development, manufacture and sales of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents and instruments / software for in vitro diagnostics in humans and animals
/03	c/o EUROIMMUN Medizinische Labordiagnostika AG An der Trave 1 23923 Selmsdorf Germany	Design, development, manufacture, service and sales of immunobiochemical test systems, immunofluorescence test systems and instruments / software for in vitro diagnostics for humans
/04	c/o EUROIMMUN Medizinische Labordiagnostika AG Am Sonnenberg 9 23627 Groß Grönau Germany	Design, development and manufacture of immunofluorescence test systems for in vitro diagnostics for humans and animals
/05	c/o EUROIMMUN Medizinische Labordiagnostika AG Am Born 24 23627 Groß Grönau Germany	Design and development of software for in vitro diagnostics for humans and animals

Annex to certificate

Standard **ISO 9001:2015**

Certificate Registr. No. **01 100 1810000**

- | | | |
|-----|--|---|
| /06 | c/o EUROIMMUN
Medizinische Labordiagnostika AG
Im Kreppel 1
02747 Herrnhut
Germany | Manufacture of immunobiochemical test systems and immunofluorescence test systems for in vitro diagnostics for humans |
| /07 | c/o EUROIMMUN
Medizinische Labordiagnostika AG
Am Pließnitztal 1
02748 Bernstadt
Germany | Manufacture of immunobiochemical test Systems, test systems for the determination of infectious agents and instruments for in vitro diagnostics for humans |
| /08 | c/o EUROIMMUN
Medizinische Labordiagnostika AG
Schloßstr. 11
91257 Pegnitz
Germany | Manufacture of immunofluorescence test systems, installation and service of instruments / software for in vitro diagnostics in humans, trainings |
| /09 | c/o EUROIMMUN
Medizinische Labordiagnostika AG
Am Flugplatz 4
23560 Lübeck
Germany | Design, development, installation, service and sales of immunobiochemical test systems, immunofluorescence test systems, molecular diagnostic / genetic test systems, test systems for the determination of infectious agents and instruments / software for in vitro diagnostics for humans and animals; trainings |
| /10 | c/o EUROIMMUN
Medizinische Labordiagnostika AG
Gewerbestr. 19
23942 Dassow
Germany | Manufacture of sheet metal and other components for instruments for in vitro diagnostics in humans and animals |

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Seekamp 31, D-23560 Lübeck, Germany

declare under our sole responsibility that the WESTERBLOT products

Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)

DY 2132-#### G

Anti-Borrelia burgdorferi-WESTERNBLOT (IgM)

DY 2132-#### M

(product name, order no)

meet the demands of

*Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998 and its transpositions in national laws which apply to it.*

Conformity assessment procedure: Annex III

Lübeck, 04.10.2016

(Place and date of issue)

Dr. Wolfgang Schlumberger
- Member of the Board -

Susanne Aleksandrowicz
- Member of the Board -

Anti-Borrelia burgdorferi-WESTERNBLOT (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DY 2132-3001 G DY 2132-24001 G	Borrelia burgdorferi (complete antigens)	IgG	Antigen coated membrane strips	30 x 01 (30) 240 x 01 (240)

Indication: The Westernblot test kit provides a qualitative in vitro assay for human autoantibodies of the immunoglobulin class IgG against Borrelia in serum or plasma for the diagnosis of Lyme borreliosis and associated diseases (Erythema chronicum migrans, lymphadenosis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculoneuritis and neuroborreliosis).

Principles of the test: The test kit contains test strips with electrophoretically separated antigen extracts of Borrelia burgdorferi. The blot strips will be blocked and incubated in the first reaction step with diluted patient samples. In the case of positive samples, specific antibodies of the class IgG (and IgA, IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Format	Format	Symbol
1. Test strips Single strips with electrophoretically separated Borrelia burgdorferi antigens	30 x 1	240 x 1	STRIPS
2. Evaluation matrix with control strip Test strip incubated with a positive control serum	1 pattern	1 pattern per test strip lot	---
3. Enzyme conjugate Alkaline phosphatase-labelled anti-human IgG (goat), 10x concentrate	2 x 3 ml	16 x 3 ml	CONJUGATE 10x
4. Universal buffer 10x concentrate	1 x 100 ml	8 x 100 ml	BUFFER 10x
5. Substrate solution Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NBT/BCIP), ready for use	1 x 50 ml	8 x 50 ml	SUBSTRATE
6. Adhesive foil	---	8 sheets	---
7. Test instruction	1 booklet	1 booklet	---

LOT Lot description

IVD In vitro diagnostic medical device

Storage temperature

Unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Undiluted patient samples and incubated blot strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

The following components are not provided in the test kits but can be ordered at EUROIMMUN under the respective order numbers. Performance of the test requires an **incubation tray**:

ZD 9895-0130 Incubation tray with 30 channels (black)

ZD 9898-0144 Incubation tray with 44 channels (black, for the EUROBlotOne and EUROBlotCamera system)

For the creation of work protocols and the evaluation of incubated test strips using **EUROLineScan** green paper and adhesive foil are required:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

If a **visual evaluation** is to be performed in individual cases, the required evaluation protocol can be ordered under: ZD 2132-0101 Evaluation protocol visual Anti-Borrelia burgdorferi-WESTERNBLOT.

Preparation and stability of the reagents

Note: The bag containing the blot strips is printed with a number in addition to the test kit lot number. This number refers to the strip batch and is also printed on the corresponding evaluation template. These two numbers must match to ensure correct evaluation of test results.

All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, if earlier, unless stated otherwise in the instructions. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- **Coated test strips:** Ready for use. Open the packing with the test strips only when the strips have reached room temperature to prevent condensation on the strips. After removal of the strips the packing should be sealed tightly and stored at +2°C to +8°C. To ensure correct evaluation of results, the lot number on the bag must match the lot number on the strips as well as on the evaluation matrix.
- **Enzyme conjugate:** The enzyme conjugate is supplied as a 10x concentrate. For the preparation of the ready for use enzyme conjugate the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with ready for use diluted universal buffer. For 1 test strip dilute 0.15 ml anti-human IgG concentrate with 1.35 ml ready for use diluted universal buffer. The ready for use diluted enzyme conjugate should be used at the same working day.
- **Universal buffer:** The universal buffer is supplied as a 10x concentrate. For the preparation of the ready for use universal buffer the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with deionised or distilled water. For the incubation of 1 test strip 1.5 ml buffer concentrate should be diluted with 13.5 ml deionised or distilled water. The ready for use diluted universal buffer should be used at the same working day.
- **Substrate solution:** Ready for use. Close bottle immediately after use, as the contents are sensitive to light ☀.

Warning: The control of human origin has tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the serum or plasma samples

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: The **patient samples** for analysis are diluted **1:51** in ready for use diluted universal buffer. For example, add 30 µl of sample to 1.5 ml ready for use diluted universal buffer and mix well by vortexing. Sample pipettes are not suitable for mixing.

Incubation

Blocking: According to the number of serum samples to be tested fill each channel of the incubation tray with 1.5 ml ready for use diluted universal buffer and a blot strip. Remove the required amount of blot strips from the packing using a pair of tweezers. The number on the test strip should be visible. Incubate for **15 minutes** at room temperature (+18°C to +25°C) on a rocking shaker. Afterwards aspirate off all the liquid.

Sample incubation: (1st step) Fill each channel with 1.5 ml of the diluted serum samples. Incubate at room temperature (+18°C to +25°C) for **30 minutes** on a rocking shaker.

Wash: Aspirate off the liquid from each channel and wash **3 x 5 minutes** each with 1.5 ml working strength universal buffer on a rocking shaker.

Conjugate incubation: (2nd step) Pipette 1.5 ml diluted enzyme conjugate (alkaline phosphatase-conjugated anti-human IgG) into each channel. Incubate for **30 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.


Wash: Aspirate off the liquid from each channel. Wash as described above.

Substrate incubation: (3rd step) Pipette 1.5 ml substrate solution into the channels of the incubation tray. Incubate for **10 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

Stopping: Aspirate off the liquid from each channel and wash each strip **3 x 1 minute** with deionised or distilled water.

For automated incubation with the EUROBlotMaster select the program **Euro02 Inf WB30**.

For automated incubation with the EUROBlotOne select the program **Euro 01/02**.

**Anti-Borrelia burgdorferi-WESTERNBLOT (IgG)****Incubation protocol**

Anti-Borrelia burgdorferi-WESTERNBLOT (IgG) Evaluation and interpretation

Handling: For evaluation of incubated test strips we generally recommend using the **EUROLineScan** software. After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. After they have dried, the test strips will be stuck to the adhesive foil. The dry test strips are then scanned using a flatbed scanner (EUROIMMUN AG) and evaluated with **EUROLineScan**. Alternatively, imaging and evaluation is possible directly from the incubation trays (EUROBlotCamera and EUROBlotOne). For general information about the EUROLineScan program please refer to the EUROLineScan user manual (EUROIMMUN AG). The code for entering the **Test** in EUROLineScan is **B.b._WB_IgG**.

If a visual evaluation must be performed in exceptional cases, hold the evaluation matrix next to the stuck-on blot strips and position it so that the black band above the number on the blot strips lines up with the alignment bar of the evaluation matrix. **The lot number on the evaluation matrix must match the lot number on the blot strips.** Clearly recognisable bands on the blot strips which concur with the labelled bands on the evaluation matrix are noted in the evaluation protocol.

Antigens: The antigen source for the EUROIMMUN Anti-Borrelia burgdorferi-WESTERNBLOT is provided by a particularly suitable Borrelia burgdorferi strain. The cultured Borrelia have been solubilised using sodium dodecyl sulphate followed by a separation of the solubilised protein using discontinuous polyacrylamide gel electrophoresis according to molecular mass and transfer of the separated proteins to nitrocellulose. From each test kit, 2 control test strips have been removed and incubated with a reference serum. One of these stained strips is included in the kit, the other remains with EUROIMMUN for documentation purposes. Diagnostically relevant antigens have been characterised with monoclonal reference antibodies from the German National Reference Laboratory for Borreliae.

Specificity of the antigens on the test strips:

Band	Antigen	Specificity
83 kDa	Membrane-vesical protein, p83	Degradation product of p100, high specificity
75 kDa	Heat shock protein, p75	Unspecific
62 kDa	Heat shock protein, p62	Unspecific
57/59 kDa	p57 and p59	Unspecific
50 kDa	p50	Unspecific
47 kDa	p47	Probably genus specific
43 kDa	p43	Unspecific
41 kDa	Flagellin, p41	Genus specific, cross reactivity to other spirochaetaceae and bacteria having flagella
39 kDa	Bmp A, p39	High specificity
34 kDa	Osp B, p34	Outer surface protein B, high specificity
32 kDa	p32	Unspecific
31 kDa	Osp A, p31	Outer surface protein A, high specificity
29 kDa	p29	Probably specific, poorly investigated
28 kDa	p28	Unspecific
25 kDa	Osp C, p25	Outer surface protein C, high specificity
21/22 kDa	p21/22	High specificity
18 kDa	p18	Probably specific
17 kDa	p17	Poorly investigated

In the lower part of the test strip there is a conjugate control membrane chip (IgA, IgG and IgM). Below the conjugate control, there is a membrane chip with a control band (Control).

Attention: A correctly performed determination of antibodies of class IgG against the antigens described above is indicated by a positive reaction of the control band and a positive reaction of the IgG band.

If one of these bands only shows a very weak reaction or none at all, the result is not valid.

Specificity of the antigens: *Borrelia burgdorferi* antigens can generally be divided into three categories.

Category	Antigens
1	Cross-reacting and undefined antigens with the molecular mass 17 kDa, 28 kDa, 32 kDa, 43 kDa, 47 kDa, 50 kDa, 57 kDa, 59 kDa, 62 kDa and 75 kDa.
2	A genus-specific antigen with a molecular mass of 41 kDa (flagellin).
3	Species-specific and highly specific antigens with the molecular mass 18 kDa, 21/22 kDa, 25 kDa, 31 kDa, 34 kDa, 39 kDa and 83 kDa.

IgG class antibodies against *Borrelia burgdorferi*

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration, as negative sera sometimes produce weak signals in individual bands. Based on experience, the results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

Result	Characteristics
Negative	No bands, or weak intensities of some antigens from categories 1 and 2.
Borderline	A distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!) and several distinctive signals from categories 1 and 2. It is recommended that a new sample be taken and the test repeated after a few weeks.
Positive	More than one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!). In addition, and particularly in the case of patients in the late stage of the disease, numerous bands from categories 1 and 2 can be observed.

IgM class antibodies against *Borrelia burgdorferi*

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration. The results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

In the early phase of a *Borrelia* infection IgM antibodies are typically directed against Osp C (p25). IgM antibodies against other specific *Borrelia* antigens are not considered definitive indicators of a fresh *Borrelia* infection.

IgM antibodies against flagellin (p41) can represent the initial response of the body to *Borrelia burgdorferi*. However, an unspecific reaction cannot be excluded, since it is known that antibodies directed against other microorganisms cross react with *Borrelia burgdorferi* flagellin (p41). For this reason, a single band at the position of flagellin (p41) in IgM detection should not be considered proof of a fresh infection with *Borrelia burgdorferi*. If only the flagellin band (p41) reacts positively, the test should be repeated several weeks later with a fresh blood sample.

In the serological investigation of a *Borrelia* infection the determination of antibodies of class IgM often yields unclear results. IgM antibodies can sometimes be found in serum years after an infection or following antibiotic treatment. Therefore, the detection of IgM antibodies does not necessarily indicate a fresh infection. A negative IgM results does not exclude a fresh infection. With a second infection, only antibodies of class IgG and not IgM can be formed.

In the late stage of borreliosis a positive IgM result does not provide any additional information, due to the antibody persistence mentioned above. The cause of such false-positive IgM results often remains unclear. They are observed, for example, in infectious mononucleosis, herpes virus infections and various autoimmune diseases.

For the diagnosis of a fresh *Borrelia* infection, a positive IgM result should be confirmed with a positive IgG result using a fresh blood sample 3 to 6 weeks later.

Result	Characteristics
Negative	No antigen bands recognisable or weak intensities of some bands of category 1
Borderline	One antigen band of category 2 (flagellin, p41) or a weak band of category 3. It is recommended that a fresh sample be taken and the test repeated after a few weeks. Antibodies against Osp C are characteristic for a fresh infection.
Positiv	At least one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!). Antibodies against Osp C are characteristic for a fresh infection.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Measurement range: The Westernblot is a qualitative method. No measurement range is provided.

Inter- and intra-assay variation: The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characterised samples on one day. In every case, the intensity of the bands was within the specified range. This Westernblot displays excellent inter- and intra-assay reproducibility.

Interference: Haemolytic, lipaemic and icteric sera showed no effect on the analytical results.

Prevalence: Sera from 156 clinically characterised patients and 517 healthy blood donors were investigated with the EUROIMMUN Anti-Borrelia burgdorferi-Westernblot.

Clinically characterised sera	n	Prevalence		
		IgG	IgM	IgM/IgG
Erythema migrans	108	64%	67%	84%
Neuroborreliosis	32	85%	44%	88%
Arthritis	10	100%	30%	100%
Acrodermatitis	6	83%	50%	100%

Serum samples	n	Prevalence	
		IgG	IgM
Healthy blood donors *	517	5%	4%

*Medical University of Luebeck

The prevalence of anti-Borrelia antibodies from the healthy blood donor samples agrees with the values from literature [14].

Cross reactivity: The quality of the antigen used (whole antigen, SDS extract) and the antigen source (*Borrelia burgdorferi sensu stricto*) ensure high specificity of the Westernblot. The determination of cross reactivity is not necessary with the Westernblot, since specific reactions and unspecific or cross reactions can be differentiated directly with this test system.

Clinical significance

The history of Lyme disease, a contagious condition caused by *Borrelia burgdorferi* and transmitted to humans by ticks, offers infectiologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Hauser noted the affection was transmitted by ticks. Independently, Afzelius, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemia of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemia, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "treponema-like spirochaetes", which were later authenticated as the causal agent: *Borrelia burgdorferi*.

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: *Borrelia*, *spirochetes*, *crisistospira*, *treponema* and *leptospira*. Among about 30 tick species feeding on humans, *Ixodes ricinus* is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi*, which causes Lyme disease, of *Anaplasma phagocytophilum* and the tick-borne encephalitis virus. *I. ricinus* ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. afzelii*, and *B. garinii*. The human seroprevalence rate of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries is about 8% (for IgG); in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-*Borrelia* antibodies in about 40%, hunters in more than 50% of cases.

A *Borrelia burgdorferi* infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi*. The main antigen VlsE, which is exclusively expressed *in vivo*, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears at one side and after several weeks on the opposite side. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiological manifestations such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% to 30%. VlsE displays the highest sensitivity of all antigens tested for IgG detection.

Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titers are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*: ELISA, indirect immunofluorescence (IIFT), passive haemagglutination and immunoblot. Generally, ELISA or IIFT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it is recommended performing ELISA/IIFT and immunoblot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. VlsE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, VlsE possesses the highest sensitivity for the detection of a *Borrelia* infection.

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Anti-Borrelia burgdorferi-WESTERNBLOT (IgM) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DY 2132-3001 M DY 2132-24001 M	Borrelia burgdorferi (complete antigens)	IgM	Antigen coated membrane strips	30 x 01 (30) 240 x 01 (240)

Indication: The Westernblot test kit provides a qualitative in vitro assay for human autoantibodies of the immunoglobulin class IgM against Borrelia in serum or plasma for the diagnosis of Lyme borreliosis. Associated diseases: Erythema chronicum migrans, lymphadenosis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculoneuritis and neuroborreliosis.

Principles of the test: The test kit contains test strips with electrophoretically separated antigen extracts of Borrelia burgdorferi. The blot strips will be blocked and incubated in the first reaction step with diluted patient samples. In the case of positive samples, specific antibodies of the class IgM (and IgA, IgG) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Format	Format	Symbol
1. Test strips Single strips with electrophoretically separated Borrelia burgdorferi antigens	30 x 1	240 x 1	STRIPS
2. Evaluation matrix with control strip Test strip incubated with a positive control serum	1 pattern	1 pattern per test strip lot	---
3. Enzyme conjugate Alkaline phosphatase-labelled anti-human IgM (goat), 10x concentrate	2 x 3 ml	16 x 3 ml	CONJUGATE 10x
4. Universal buffer 10x concentrate	1 x 100 ml	8 x 100 ml	BUFFER 10x
5. Substrate solution Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NBT/BCIP), ready for use	1 x 50 ml	8 x 50 ml	SUBSTRATE
6. Adhesive foil	---	8 sheets	---
7. Test instruction	1 booklet	1 booklet	---
LOT Lot description	CE		Storage temperature
IVD In vitro diagnostic medical device			Unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Undiluted patient samples and incubated blot strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

The following components are not provided in the test kits but can be ordered at EUROIMMUN under the respective order numbers.

Performance of the test requires an **incubation tray**:

ZD 9895-0130 Incubation tray with 30 channels (black)

ZD 9898-0144 Incubation tray with 44 channels (black, for the EUROBlotOne and EUROBlotCamera system)

For the creation of work protocols and the evaluation of incubated test strips using **EUROLineScan** green paper and adhesive foil are required:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

If a **visual evaluation** is to be performed in individual cases, the required evaluation protocol can be ordered under: ZD 2132-0101 Evaluation protocol visual Anti-Borrelia-burgdorferi-WESTERNBLOT

Preparation and stability of the reagents

Note: The bag containing the blot strips is printed with a number in addition to the test kit lot number. This number refers to the strip batch and is also printed on the corresponding evaluation template. These two numbers must match to ensure correct evaluation of test results.

All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, if earlier, unless stated otherwise in the instructions. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- **Coated test strips:** Ready for use. Open the packing with the test strips only when the strips have reached room temperature to prevent condensation on the strips. After removal of the strips the packing should be sealed tightly and stored at +2°C to +8°C. To ensure correct evaluation of results, the lot number on the bag must match the lot number on the strips as well as on the evaluation matrix.
- **Enzyme conjugate:** The enzyme conjugate is supplied as a 10x concentrate. For the preparation of the ready for use enzyme conjugate the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with ready for use diluted universal buffer. For 1 test strip dilute 0.15 ml anti-human IgG concentrate with 1.35 ml ready for use diluted universal buffer. The ready for use diluted enzyme conjugate should be used at the same working day.
- **Universal buffer:** The universal buffer is supplied as a 10x concentrate. For the preparation of the ready for use universal buffer the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with deionised or distilled water. For the incubation of 1 test strip 1.5 ml buffer concentrate should be diluted with 13.5 ml deionised or distilled water. The ready for use diluted universal buffer should be used at the same working day.
- **Substrate solution:** Ready for use. Close bottle immediately after use, as the contents are sensitive to light ☼.

Warning: The control of human origin has tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the serum or plasma samples

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: The **patient samples** for analysis are diluted **1:51** in ready for use diluted universal buffer. For example, add 30 µl of sample to 1.5 ml ready for use diluted universal buffer and mix well by vortexing. Sample pipettes are not suitable for mixing.

Incubation

Blocking: According to the number of serum samples to be tested fill each channel of the incubation tray with 1.5 ml ready for use diluted universal buffer and a blot strip. Remove the required amount of blot strips from the packing using a pair of tweezers. The number on the test strip should be visible. Incubate for **15 minutes** at room temperature (+18°C to +25°C) on a rocking shaker. Afterwards aspirate off all the liquid.

Sample incubation: (1st step) Fill each channel with 1.5 ml of the diluted serum samples and incubate at room temperature (+18°C to +25°C) for **30 minutes** on a rocking shaker.

Wash: Aspirate off the liquid from each channel and wash **3 x 5 minutes** each with 1.5 ml working strength universal buffer on a rocking shaker.

Conjugate incubation: (2nd step) Pipette 1.5 ml diluted enzyme conjugate (alkaline phosphatase-conjugated anti-human IgM) into each channel and incubate for **30 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.


Wash: Aspirate off the liquid from each channel. Wash as described above.

Substrate incubation: (3rd step) Pipette 1.5 ml substrate solution into the channels of the incubation tray. Incubate for **10 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

Stopping: Aspirate off the liquid from each channel and wash each strip **3 x 1 minute** with deionised or distilled water.

For automated incubation with the EUROBlotMaster select the program **Euro02 Inf WB30**.

For automated incubation with the EUROBlotOne select the program **Euro 01/02**.

**Anti-Borrelia-burgdorferi WESTERNBLOT (IgM)****Incubation protocol**

Evaluation and Interpretation of the results of the Anti-Borrelia-burgdorferi WESTERNBLOT (IgM)

Handling: For evaluation of incubated test strips we generally recommend using the **EUROLineScan** software. After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. After they have dried, the test strips will be stuck to the adhesive foil. The dry test strips are then scanned using a flatbed scanner (EUROIMMUN AG) and evaluated with **EUROLineScan**. Alternatively, imaging and evaluation is possible directly from the incubation trays (EUROBlotCamera and EUROBlotOne). For general information about the EUROLineScan program please refer to the EUROLineScan user manual (EUROIMMUN AG). The code for entering the **Test** in EUROLineScan is **B.b._WB_IgM**.

If a visual evaluation must be performed in exceptional cases, hold the evaluation matrix next to the stuck-on blot strips and position it so that the black band above the number on the blot strips lines up with the alignment bar of the evaluation matrix. **The lot number on the evaluation matrix must match the lot number on the blot strips.** Clearly recognisable bands on the blot strips which concur with the labelled bands on the evaluation matrix are noted in the evaluation protocol.

Antigens: The antigen source for the EUROIMMUN Anti-Borrelia-burgdorferi WESTERNBLOT is provided by a particularly suitable *Borrelia burgdorferi* strain [2, 3]. The cultured borrelia have been solubilised using sodium dodecyl sulphate followed by a separation of the solubilised protein using discontinuous polyacrylamide gel electrophoresis according to molecular mass and transfer of the separated proteins to nitrocellulose. From each test kit, 2 control test strips have been removed and incubated with a reference serum. One of these stained strips is included in the kit, the other remains with EUROIMMUN for documentation purposes. Diagnostically relevant antigens have been characterised with monoclonal reference antibodies from the German National Reference Laboratory for Borreliae [39].

Specificity of the antigens on the test strips: [41, 42, 43]

Band	Antigen	Specificity
83 kDa	Membrane-vesical protein, p 83	Degradation product of p 100, high specificity.
75 kDa	Heat shock protein, p 75	Unspecific.
62 kDa	Heat shock protein, p 62	Unspecific.
57/59 kDa	p 57 and p 59	Unspecific.
50 kDa	p 50	Unspecific.
47 kDa	p 47	Probably genus specific.
43 kDa	p 43	Unspecific.
41 kDa	Flagellin, p 41	Genus specific, cross reactivity to other spirochaetaceae and bacteria having flagella.
39 kDa	Bmp A, p 39	High specificity.
34 kDa	Osp B, p 34	Outer surface protein B, high specificity.
32 kDa	p 32	Unspecific.
31 kDa	Osp A, p 31	Outer surface protein A, high specificity.
29 kDa	p 29	Probably specific, poorly investigated.
28 kDa	p 28	Unspecific.
25 kDa	Osp C, p 25	Outer surface protein C, high specificity.
21/22 kDa	p 21/22	High specificity.
18 kDa	p 18	Probably specific.
17 kDa	p 17	Poorly investigated.

In the lower part of the test strip there is a conjugate control membrane chip (IgA, IgG and IgM). Below the conjugate control, there is a membrane chip with a control band (Control).

Attention: A correctly performed determination of antibodies of class IgM against the antigens described above is indicated by a positive reaction of the control band and a positive reaction of the IgM band.

If one of these bands only shows a very weak reaction or none at all, the result is not valid.

Specificity of the antigens: *Borrelia burgdorferi* antigens can generally be divided into three categories [40].

Category	Antigens
1	Cross-reacting and undefined antigens with the molecular mass 17 kDa, 28 kDa, 32 kDa, 43 kDa, 47 kDa, 50 kDa, 57 kDa, 59 kDa, 62 kDa and 75 kDa.
2	A genus-specific antigen with a molecular mass of 41 kDa (flagellin).
3	Species-specific and highly specific antigens with the molecular mass 18 kDa, 21/22 kDa, 25 kDa, 31 kDa, 34 kDa, 39 kDa and 83 kDa.

IgG class antibodies against *Borrelia burgdorferi*

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration, as negative sera sometimes produce weak signals in individual bands. Based on experience, the results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

Result	Characteristics
Negative	No bands, or weak intensities of some antigens from categories 1 and 2.
Borderline	A distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!) and several distinctive signals from categories 1 and 2. It is recommended that a new sample be taken and the test repeated after a few weeks.
Positive	More than one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!). In addition, and particularly in the case of patients in the late stage of the disease, numerous bands from categories 1 and 2 can be observed.

IgM class antibodies against *Borrelia burgdorferi*

Interpretation of results: In order to evaluate the signals, the band positions and intensity of staining must be taken into consideration. The results of the *Borrelia burgdorferi* WESTERNBLOT test can be divided into negative, borderline and positive results.

In the early phase of a *Borrelia* infection IgM antibodies are typically directed against Osp C (p 25). IgM antibodies against other specific *Borrelia* antigens are not considered definitive indicators of a fresh *Borrelia* infection.

IgM antibodies against flagellin (p 41) can represent the initial response of the body to *Borrelia burgdorferi*. However, an unspecific reaction cannot be excluded, since it is known that antibodies directed against other microorganisms cross react with *Borrelia burgdorferi* flagellin (p 41). For this reason, a single band at the position of flagellin (p 41) in IgM detection should not be considered proof of a fresh infection with *Borrelia burgdorferi*. If only the flagellin band (p 41) reacts positively, the test should be repeated several weeks later with a fresh blood sample.

In the serological investigation of a *Borrelia* infection the determination of antibodies of class IgM often yields unclear results. IgM antibodies can sometimes be found in serum years after an infection or following antibiotic treatment. Therefore, the detection of IgM antibodies does not necessarily indicate a fresh infection. A negative IgM result does not exclude a fresh infection. With a second infection, only antibodies of class IgG and not IgM can be formed.

In the late stage of borreliosis a positive IgM result does not provide any additional information, due to the antibody persistence mentioned above. The cause of such false-positive IgM results often remains unclear. They are observed, for example, in infectious mononucleosis, herpes virus infections and various autoimmune diseases [39].

For the diagnosis of a fresh *Borrelia* infection, a positive IgM result should be confirmed with a positive IgG result using a fresh blood sample 3-6 weeks later.

Result	Characteristics
Negative	No antigen bands recognisable or weak intensities of some bands of category 1
Borderline	One antigen band of category 2 (flagellin, p 41) or a weak band of category 3. It is recommended that a fresh sample be taken and the test repeated after a few weeks. Antibodies against Osp C are characteristic for a fresh infection.
Positive	At least one distinctive band from category 3 (antigens in category 3 are shaded grey in the evaluation protocol!). Antibodies against Osp C are characteristic for a fresh infection.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Measurement range: The Westernblot is a qualitative method. No measurement range is provided.

Inter- and intra-assay variation: The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characterised samples on one day. In every case, the intensity of the bands was within the specified range. This Westernblot displays excellent inter- and intra-assay reproducibility.

Interference: Haemolytic, lipaemic and icteric sera showed no effect on the analytical results.

Prevalence: Sera from 156 clinically characterised patients and 517 healthy blood donors were investigated with the EUROIMMUN Anti-Borrelia burgdorferi-WESTERNBLOT.

Clinically characterised sera	n	Prevalence		
		IgG	IgM	IgM/IgG
Erythema migrans	108	64%	67%	84%
Neuroborreliosis	32	85%	44%	88%
Arthritis	10	100%	30%	100%
Acrodermatitis	6	83%	50%	100%

Serum samples	n	Prevalence	
		IgG	IgM
Healthy blood donors *	517	5%	4%

*Medical University of Luebeck

The prevalence of anti-Borrelia antibodies from the healthy blood donor samples agrees with the values from literature [7].

Cross reactivity: The quality of the antigen used (whole antigen, SDS extract) and the antigen source (*Borrelia burgdorferi sensu stricto*) ensure high specificity of the Westernblot. The determination of cross reactivity is not necessary with the Westernblot, since specific reactions and unspecific or cross reactions can be differentiated directly with this test system.

Clinical significance

The history of Lyme disease, a contagious condition, caused by *Borrelia burgdorferi*, transmitted to humans by ticks, offers infectiologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis [1]. Fifty years later Hauser noted the affection was transmitted by ticks [1]. Independently, Afzelius, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemic of arthritis occurred in children in Lyme, Connecticut, USA [1]. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemic, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that tick bites transmitted a Spirochaeta which was later authenticated as the causal agent: *Borrelia burgdorferi* [1, 2, 3].

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: borrelia, spirochaetes, cristispira, treponema and leptospira. Among about 30 tick species feeding on humans, *Ixodes ricinus* is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi* which causes Lyme disease, and of the tick-borne encephalitis virus [4, 5, 6, 7]. *I. ricinus* ticks pass through three developmental stages: larvae, nymphs and adults (females and males) [8]. The density of this tick species may be very high, reaching in some places more than 300 ticks/100 m² [8]. Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere [9, 10, 11]. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. afzelii*, and *B. garinii* [4, 6, 10, 11, 12]. The human seroprevalence rates of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries ranges from 18% to 52% [13]. In East Asia, e.g. in China, there are 26% on the average [14]. Persons working in the area of forestry display anti-borrelia antibodies in about 40%, hunters in more than 50% of cases [7, 15].

A *Borrelia burgdorferi* infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine [7, 16]. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner [17]. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting [18]. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular, hepatitis, regional lymphadenitis, and mixed [19, 20, 21]. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients in Stage I. The prevalence of specific IgG antibodies is considerably lower [22, 23]. However, serological tests often provide negative results in this stage of the illness [3, 19]. In order to overcome this problem a new test was created including the newly identified VlsE (variable major protein-like sequence, expressed), which can be considered as the major antigen for *Borrelia* serology. Over 85% of IgG-positive sera can be identified at a glance by assessing the VlsE band in incubated Westernblot strips. VlsE allows detection of antibodies against all *Borrelia* species. The risk of a false negative reaction due to species differences is ten times lower [24]. Humans also produce highly specific borrelicidal antibodies against outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi sensu stricto*: immunoglobulin M (IgM) OspC and immunoglobulin G (IgG) OspC borrelicidal antibodies [25]. VlsE is the most sensitive antigen for IgG antibody detection, OspC for IgM antibody detection [3, 26, 27].

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis [20, 28, 29]. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears at one side and after several weeks on the opposite side [20]. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found [30]. Rarer are cardiological manifestations, such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II [24]. In the early phase of this stage, mainly IgM antibodies are found, but in the late phase, often only IgG antibodies occur. However, IgM antibodies can persist for a long time [3, 30]. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20%-30% [31]. VlsE displays the highest sensitivity of all antigens tested for IgG detection [26].

Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibodies are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable [3, 24].

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing [21, 29, 32]. An Antibody Index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in blood/CSF barrier function [33].

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi* [34, 35, 36, 37, 38, 39]: ELISA, indirect immunofluorescence, passive haemagglutination and immunoblot [26]. Many investigators employ an ELISA test in parallel with an indirect immunofluorescence test for preliminary characterisation of the serum sample [26, 34]. According to guidelines of the USA and Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as first step: if reactive, followed by Westernblot [3, 24, 40]. In fresh infections it is recommended performing ELISA/IIFT and Westernblot in parallel, since some weak reactions become detectable earlier in Westernblot than in the screening tests [12, 24]. VlsE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection, should be included in the test [3]. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% compared to whole extract Westernblots and by 30% compared to recombinant antigen Westernblots [24, 27, 31]. Of all recombinant antigens tested, VlsE possesses the highest sensitivity for the detection of a *Borrelia* infection [31, 32, 33].

Literature references

1. Belaich S. **Lyme disease**. Presse Med 14 (1995) 81-87.
2. Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Grunwaldt E, Davis JP. **Lyme disease - a tick-borne Spirochetosis?** Science 216 (1982) 1317-1319.
3. Johnson R, Schmid GP, Hyde FW, Steigerwalt AG, Brenner DJ. **Borrelia burgdorferi sp. nov.: Etiologic agent of Lyme disease**. Int J Syst Bacteriol 34 (1984) 496-497.
4. Wilske B. **Epidemiology and diagnosis of Lyme borreliosis**. Ann Med 37 (2005) 568-579.
5. Derdakova M, Lencakova D. **Association of genetic variability within the Borrelia burgdorferi sensu lato with the ecology, epidemiology of Lyme borreliosis in Europe**. Ann Agric Environ Med 12 (2005) 165-172.
6. Steffen I, Hirsch HH. **Diagnostic tests of Lyme borreliosis**. Ther Umsch 62 (2005) 737-744.
7. Robert-Koch-Institut. **Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg**. Epidemiologisches Bulletin 14 (1998) 93-97.
8. Gern L. **The biology of the Ixodes ricinus tick**. Ther Umsch 62 (2005) 707-712.
9. Holzer BR. **Tick borne diseases**. Ther Umsch 62 (2005) 757-763.
10. Bazovska S, Machacova E, Spalekova M, Kontrosova S. **Reported incidence of Lyme disease in Slovakia and antibodies to B. burgdorferi antigens detected in healthy population**. Bratisl Lek Listy 106 (2005) 270-273.
11. Egger M. **Lyme borreliosis - an overview**. Ther Umsch 62 (2005) 731-736.
12. Hauser U, Lehnert G, Lobentanzer R, Wilske B. **Interpretation Criteria for Standardized Western Blots for Three European Species of Borrelia burgdorferi Sensu Lato**. J Clin Microbiol 35 (1997) 1433-1444.
13. Oehme R, Hartelt K, Backe H, Brockmann S, Kimmig P. **Foci of tick-borne diseases in southwest Germany**. Int J Med Microbiol 291 (2002) 22-29.
14. Ai CX, Zhang WF, Zhao JH. **Sero-epidemiology of Lyme disease in an endemic area in China**. Microbiol Immunol 38 (1994) 505-509.
15. Cetin E, Sotoudeh M, Auer H, Stanek G. **Paradigm Burgenland: risk of Borrelia burgdorferi sensu lato infection indicated by variable seroprevalence rates in hunters**. Wien Klin Wochenschr 118 (2006) 677-681.

16. Escudero-Nieto R, Guerrero-Espejo A. **Diseases produced by Borrelia.** *Enferm Infecc Microbiol Clin* 23 (2005) 232-240.
17. Vorob'eva NN, Sumlivaia ON. **Clinical manifestations in acute Lyme diseases.** *Med Parazitol (Mosk)* 4 (2003) 3-7.
18. Roscoe C, Epperly T **Tick-borne relapsing fever.** *Am Fam Physician* 15 (2005) 2039-2044.
19. Coulter P, Lema C, Flayhart D, Linhardt AS, Aucott JN, Auwaerter PG, Dumler JS. **Two-year evaluation of Borrelia burgdorferi culture and supplemental tests for definitive diagnosis of Lyme disease.** *J Clin Microbiol* 43 (2005) 5080-5084.
20. Mlodzikowska-Albrecht J, Zarowski M, Steinborn B, Winczewska-Wiktor A, Gurda B, Wigowska-Sowinska J. **Bilateral facial nerve palsy in the course of neuroborreliosis in children-dynamics, laboratory tests and treatment.** *Rocz Akad Med Bialymst* 50 Suppl 1 (2005) 64-69.
21. Christen HJ, Hanefeld F, Eiffert H, Thomssen R. **Epidemiology and clinical manifestations of Lyme borreliosis in childhood. A prospective multicentre study with special regard to neuroborreliosis.** *Acta Paediatr Suppl* 386 (1993) 1-75.
22. Bryksin AV, Godfrey HP, Carbonaro CA, Wormser GP, Aguero-Rosenfeld ME, Cabello FC. **Borrelia burgdorferi BmpA, BmpB, and BmpD proteins are expressed in human infection and contribute to P39 immunoblot reactivity in patients with Lyme disease.** *Clin Diagn Lab Immunol* 12 (2005) 935-940.
23. Christova I. **Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis.** *Int J Immunopathol Pharmacol* 16 (2003) 261-268.
24. Wilske B, Fingerle V, Schulte-Spechtel U. **Microbiological and serological diagnosis of Lyme borreliosis.** *FEMS Immunol Med Microbiol* 49 (2007) 13-21.
25. Lovrich SD, Jobe DA, Schell RF, Callister SM. **Borrelia burgdorferi OspC antibodies specific for a highly conserved epitope are immunodominant in human Lyme disease and do not occur in mice or hamsters.** *Clin Diagn Lab Immunol* 12 (2005) 746-751.
26. Goettner G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. **Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of VlsE and DbpA homologues.** *J Clin Microbiol* 43 (2005) 3602-3609.
27. Zhang JR, Hardham JM, Barbour AG, Norris SJ. **Antigenic Variation in Lyme Disease Borreliae by Promiscuous Recombination of VMP-like Sequence Cassettes.** *Cell* 89 (1997) 275-285.
28. Brooks CS, Vuppala SR, Jett AM, Akins DR. **Identification of Borrelia burgdorferi outer surface proteins.** *Infect Immun* 74 (2006) 296-304.
29. Steenhoff AP, Smith MJ, Shah SS, Coffin SE. **Neuroborreliosis with progression from pseudotumor cerebri to aseptic meningitis.** *Pediatr Infect Dis J* 25 (2006) 91-92.
30. Shi C, Wolfe J, Russell JQ, Fortner K, Hardin N, Anguita J, Budd RC. **Fas ligand deficiency impairs host inflammatory response** *Infect Immun* 74 (2006) 1156-1160.
31. Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimerl C, Johnson BJ, Wilske B. **Significant improvement of the recombinant Borrelia-specific immunoglobulin G immunoblot test by addition of VlsE and a DbpA homologue derived from 2003** 1299-1303.
32. Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimerl C, Johnson B, Wilske B. **Significant improvement of the recombinant Borrelia IgG immunoblot for serodiagnosis of early neuroborreliosis.** *Int J Med Microbiol* 293 (2004) 158-160.
33. Reiber H, Lange P. **Virus-spezifische Antikörper in Liquor und Serum. ELISA-Analytik und Auswertung mittels Antikörper-Index und Quotientendiagramm.** *Lab Med* 15 (1991) 204-207.
34. Jovicic VLj, Grego EM, Lako BL, Ristic BM, Lepsanovic ZA, Stajkovic NT. **Improved serodiagnosis of early Lyme borreliosis: immunoblot with local Borrelia afzelii strain.** *APMIS* 111 (2003) 1053-1059.

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declare under our sole responsibility that the ELISA-products

Anti-Borrelia plus VlsE ELISA (IgG)

EI 2132-9601-2 G

CSF: Anti-Borrelia plus VlsE ELISA (IgG)

EI 2132-9601-L G

CSQ pair of controls Anti-Borrelia (IgG)

EI 2132-0208-8 L G

(product name, order no)

in combination with automated analyzer for ELISA

EUROIMMUN Analyzer I

meet the demands of

*Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998*

Lübeck, 26.05.2010

(Place and date of issue)

Wolfgang Schlumberger, PhD
- Member of the Board -

Susanne Aleksandrowicz
- Member of the Board -

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declare under our sole responsibility that the ELISA products

Anti-Borrelia ELISA (IgG)

EI 2132-9601 G

Anti-Borrelia ELISA (IgM)

EI 2132-9601 M

CSF: Anti-Borrelia ELISA (IgM)

EI 2132-9601-L M

(product name, order no)

meet the demands of

*Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998*

Lübeck, 25.04.2008

(Place and date of issue)

Wolfgang Schlumberger, PhD
- Member of the Board -

Susanne Aleksandrowicz
- Member of the Board -

Anti-Borrelia plus VlsE ELISA (IgG)

Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2132-9601-2 G	Borrelia burgdorferi, Borrelia afzelii, Borrelia garinii, Borrelia VlsE Antigen	IgG	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the IgG class against Borrelia antigens in serum or plasma for the diagnosis of infection with Borrelia and associated diseases: Erythema chronicum migrans, lymphadenosis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculitis and neuro-borreliosis.

Application: Clinical diagnosis of borreliosis can be achieved by determination of Borrelia-specific antibodies of classes IgG and IgM. The Anti-Borrelia plus VlsE ELISA (IgG) is based on an optimised lysate mixture of the most relevant human pathogenic Borrelia strains and also contains recombinant VlsE. Due to its complete antigen spectrum, the ELISA offers a very high sensitivity and is therefore ideally suited for antibody screening. Samples with positive or borderline ELISA results should always be further investigated using immunoblot.

Principle of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with a mix of whole antigen extracts of Borrelia burgdorferi sensu stricto, Borrelia afzelii, Borrelia garinii and recombinant VlsE of Borrelia burgdorferi. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	STRIPS
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	CAL 1
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	CAL 2
4. Calibrator 3 2 RU/ml (IgG, human), ready for use	light red	1 x 2.0 ml	CAL 3
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	POS CONTROL
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	NEG CONTROL
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	CONJUGATE
8. Sample buffer ready for use	light blue	1 x 100 ml	SAMPLE BUFFER
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	WASH BUFFER 10x
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	SUBSTRATE
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	STOP SOLUTION
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
LOT Lot description			 Storage temperature
IVD In vitro diagnostic medical device			 Unopened usable until

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bags).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light ☀. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples must be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** with sample buffer.

Example: Add 10 µl of sample to 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: The calibrators and controls are prediluted and ready for use, do not dilute them.

Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation:
(1st step) Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.
Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing:
Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation:
(2nd step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.
Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Empty the wells. Wash as described above.

Substrate incubation:
(3rd step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells.
Incubate for **15 minutes** at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping: Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement: **Photometric measurement** of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient sample (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient sample (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. Therefore, the number of tests performed can be matched to the number of samples, minimising reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of the calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive

Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use “point-to-point” plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.

If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as “>200 RU/ml”. It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4. The upper limit of the normal range of non-infected persons (cut-off value) recommended by EUROIMMUN is 20 relative units (RU)/ml. EUROIMMUN recommends interpreting results as follows:

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For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends to retest the samples.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7 to 10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against *Borrelia*, the calibration is performed in relative units (RU)/ml.

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is provided by particularly suitable *Borrelia* strains (*Borrelia burgdorferi sensu stricto*, *Borrelia afzelii*, *Borrelia garinii*) as well as recombinant VlsE (variable major protein-like sequence, expressed) of *Borrelia burgdorferi*. VlsE is a newly characterised surface protein of *Borrelia* which is expressed exclusively *in vivo* and which contains conserved and highly immunogenic epitopes. The cultured *borrelia* have been solubilised by using sodium dodecyl sulphate. The used antigen mixture contains all relevant proteins.

Linearity: The linearity of the Anti-*Borrelia* plus VlsE ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95 . The Anti-*Borrelia* plus VlsE ELISA (IgG) is linear at least in the tested concentration range (8 RU/ml to 124 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-*Borrelia* plus VlsE ELISA (IgG) is 0.7 RU/ml.

Cross reactivity: Cross reactivity of the Anti-*Borrelia* plus VlsE ELISA (IgG) was evaluated in a study performed on 256 patient sera with antibodies against *Treponema pallidum* and 16 patient sera with antibodies against *Leptospira*. 95 sera of the Anti-*Treponema pallidum*-positive sera and 1 serum of the Anti-*Leptospira*-positive sera presented reactivity. Therefore, the cross reactivity amounts to 37.1% and 6.3%, respectively. Cross reactions with antibodies against other spirochaetes cannot be excluded. The Anti-*Borrelia* plus VlsE ELISA is designed as a screening test which should be followed by a confirmatory test (Western- or line blot). Further sera from patients with different possible influencing factors were investigated. An overview of these results can be found in the following table.

Possible influencing factors	n	Anti- <i>Borrelia</i> plus VlsE ELISA (IgG) positive
Acute EBV infection	33	0%
Anti-HSV 1	12	0%
Anti-CMV	12	0%
Anti-VZV	12	0%
Anti-Adenovirus	12	0%
Anti-RSV	12	0%
Anti-Parainfluenza	12	0%
Anti-Influenza-A virus	12	0%
Anti-Influenza-B virus	12	0%
Anti- <i>Mycoplasma pneumoniae</i>	12	0%
Anti-Measles virus	12	0%
Anti-Mumps virus	12	0%
Anti-Rubella virus	12	0%
Anti- <i>Toxoplasma gondii</i>	11	0%
Anti- <i>Chlamydia pneumoniae</i>	11	0%
Anti- <i>Helicobacter pylori</i>	12	0%
ANA + DNS (AAb)	36	2.8%
Rheuma factor	37	0%
Neurological diseases	54	3.7%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	52	3.5
2	97	3.5
3	128	3.6

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	55	3.6
2	103	3.7
3	138	3.8

Specificity and sensitivity: Sera of 165 patients with suspected borreliosis were analysed using the EUROIMMUN Anti-Borrelia plus VisE ELISA (IgG) and the EUROIMMUN Anti-Borrelia EUROLINE-Westernblot (IgG) as a reference method. The test showed a specificity of 90.2% and a sensitivity of 100%.

n = 165		EUROIMMUN EUROLINE-Westernblot		
		positive	borderline	negative
ELISA	positive	60	11	9
	borderline	0	1	1
	negative	0	0	83

Borderline results are not included by the calculation of specificity and sensitivity.

Clinical study: 138 patients with clinically characterised borreliosis in different disease stages were screened with the Anti-Borrelia plus VisE-ELISA (IgG) and Anti-Borrelia-ELISA (IgM). Sensitivities between 89 to 96% were found.

n = 138		EUROIMMUN Anti-Borrelia plus VisE-ELISA (IgG), EUROIMMUN Anti-Borrelia-ELISA (IgM)		
		IgG positive	IgM positive	IgG and/or IgM positive
Clinic	erythema migrans, n = 97	60	68	86 (88.7%)
	mono- or poly- arthritis, n = 26	22	5	25 (96.2%)
	neuroborreliosis, n = 15	13	4	14 (93.3%)

All ELISA negative results were confirmed as negative by an Anti-Borrelia Westernblot.

Reference range: The levels of the anti-Borrelia antibodies (IgG) were analysed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 20 RU/ml, 5% of the blood donors were anti-Borrelia positive (IgG), which reflects the known percentage of infections in adults.

Clinical significance

The history of Lyme disease, a contagious condition caused by *Borrelia burgdorferi* and transmitted to man by ticks, offers infectiologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Hauser noted the affection was transmitted by ticks. Independently, Afzelius, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemic of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemic, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "treponema-like spirochaetes", which were later authenticated as the causal agent: *Borrelia burgdorferi*.

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: borrelia, spirochaetes, cristispira, treponema and leptospira. Among about 30 tick species feeding on humans, *Ixodes ricinus* is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi*, which causes Lyme disease, of *Anaplasma phagocytophilum* and the tick-borne encephalitis virus. *I. ricinus* ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. afzelii*, and *B. garinii* (in rare cases also by *B. spielmanii* and *B. bavariensis*). The human seroprevalence rate of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries is about 8% (for IgG); in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-*Borrelia* antibodies in about 40%, hunters in more than 50% of cases.

A *Borrelia burgdorferi* infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi*. The main antigen VlsE, which is exclusively expressed in vivo, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears on one side and after several weeks on the opposite side. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiological manifestations such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% to 30%. VlsE displays the highest sensitivity of all antigens tested for IgG detection.

Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titers are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*: ELISA, indirect immunofluorescence (IIFT), passive haemagglutination and immunoblot. Generally, ELISA or IIFT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it is recommended performing ELISA/IIFT and immuno-blot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. VlsE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, VlsE possesses the highest sensitivity for the detection of a *Borrelia* infection.

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10. Goettner G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. **Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of VlsE and DbpA homologues.** J Clin Microbiol 43 (2005) 3602-3609.
11. Hauser U, Lehnert G, Lobentanzer R, Wilske B. **Interpretation Criteria for Standardized WesternBlots for Three European Species of Borrelia burgdorferi Sensu Lato.** J Clin Microbiol 35 (1997) 1433-1444.
12. Johnson R, Schmid GP, Hyde FW, Steigerwalt AG, Brenner DJ. **Borrelia burgdorferi sp. Nov.: Etiologic agent of lyme disease.** Int J Syst Bacteriol 34 (1984) 496-497.
13. Oehme R, Hartelt K, Backe H, Brockmann S, Kimmig P. **Foci of tick-borne diseases in south-west Germany.** Int J Med Microbiol 291 (2002) 22-29.
14. Robert-Koch-Institut. **Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg.** Epidemiologisches Bulletin 14 (1998) 93-97.
15. Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimerl C, Johnson BJ, Wilske B. **Significant improvement of the recombinant Borrelia-specific immunoglobulin G immuno-blot test by addition of VlsE and a DbpA homologue derived from Borrelia garinii for diagnosis of early neuroborreliosis.** J Clin Microbiol 41 (2003) 1299-1303.
16. Tilton RC, Ryan RW. **The laboratory diagnosis of Lyme disease.** J Clin Immunoassay 16 (1993) 208-214.
17. Wilske B, Fingerle V, Schulte-Spechtel U. **Microbiological and serological diagnosis of Lyme borreliosis.** FEMS Immunol Med Microbiol 49 (2007) 13-21.
18. Wilske B. **Epidemiology and diagnosis of Lyme borreliosis.** Ann Med 37 (2005) 568-579.
19. Wilske B, Zöller L, Brade V, Eiffert H, Göbel UB, Stanek G, Pfister HW. **MiQ 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik: Lyme-Borreliose.** Urban & Fischer Verlag München Jena (2000).
20. Zhang JR, Hardham JM, Barbour AG, Norris SJ. **Antigenic Variation in Lyme Disease Borreliae by Promiscuous Recombination of VMP-like Sequence Cassettes.** Cell 89 (1997) 275-285.

Anti-Borrelia ELISA (IgM) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2132-9601 M	Borrelia burgdorferi, Borrelia afzelii, Borrelia garinii	IgM	Ag-coated microplate wells	96 x 01 (96)

Indications: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the IgM class against Borrelia antigens in serum or plasma for the diagnosis of Infection with Borrelia. Associated diseases: Erythema chronicum migrans, lymphadenosis cutis benigna, acrodermatitis chronica atrophicans, arthritis, carditis, lymphocytic meningoradiculitis, neuroborreliosis.

Application: Clinical diagnosis of borreliosis can be achieved by determination of Borrelia-specific antibodies of classes IgG and IgM. The Anti-Borrelia ELISA (IgM) is based on an optimised lysate mixture of the most relevant human pathogenic Borrelia strains. Specific cultivation methods help to ensure a large proportion of OspC in the test. Due to its complete antigen spectrum, the ELISA offers a very high sensitivity and is therefore ideally suited for antibody screening. Samples with positive or borderline ELISA results should always be further investigated using immunoblot.

Principle of the test: The test kit contains microtiter strips each with 8 break-off reagent wells coated with antigen extracts of Borrelia burgdorferi sensu stricto, Borrelia afzelii and Borrelia garinii. In the first reaction step, diluted patient samples are incubated with the wells. In the case of positive samples, specific IgM antibodies (also IgA and IgG) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 , 200 RU/ml (IgM, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 , 20 RU/ml (IgM, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 , 2 RU/ml (IgM, human), ready for use	light red	1 x 2.0 ml	
5. Positive control , (IgM, human), ready for use	blue	1 x 2.0 ml	
6. Negative control , (IgM, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgM (goat), ready for use	red	1 x 12 ml	
8. Sample buffer containing IgG/RF absorbent (anti-human IgG antibody preparation obtained from goat), ready for use	green	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
 Lot description			 Storage temperature
 In vitro diagnostic medical device			 Unopened usable until

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature (+18°C to +25°C) to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bags).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use. The green coloured sample buffer contains IgG/RF absorbent. Serum or plasma samples diluted with this sample buffer are only to be used for the determination of IgM antibodies.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionised or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light . The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Warning: The calibrators and controls of human origin have tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Introduction: Before the determination of specific antibodies of class IgM, antibodies of class IgG should be removed from the patient sample. This procedure must be carried out in order to prevent any rheumatoid factors from reacting with specifically bound IgG, which would lead to false positive IgM test results, and to prevent that specific IgG displace IgM from the antigen, which would lead to false IgM negative test results.

Functional principle: The sample buffer (green coloured!) contains an anti-human antibody preparation from goat. IgG of a serum or plasma sample is bound with high specificity by these antibodies and precipitated. If the sample also contains rheumatoid factors, these will be absorbed by the IgG/anti-human IgG complex.

Separation properties:

- All IgG subclasses are bound and precipitated by the anti-human IgG antibodies.
- Human serum IgG in concentrations of up to 15 mg per ml are removed (average serum IgG concentration in adults: 12 mg per ml).
- Rheumatoid factors are also removed.
- The recovery rate of the IgM fraction is almost 100%.

Performance: The **patient samples** for analysis are diluted **1:101** with green coloured sample buffer. For example, add 10 µl sample to 1.0 ml sample buffer and mix well by vortexing. Sample pipettes are not suitable for mixing. Incubate the mixture for at least **10 minutes** at room temperature (+18°C to +25°C). Subsequently, it can be pipetted into the microplate wells according to the pipetting protocol.

Notes:

- Antibodies of the class IgG should not be analysed with this mixture.
- It is possible to check the efficacy of the IgG/RF absorbent for an individual patient sample by performing an IgG test in parallel to the IgM test using the mixture. If the IgG test is negative, the IgM result can be considered as reliable.
- The calibrators and controls are ready for use, do not dilute them.

Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step) Transfer 100 µl of the calibrators, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Manual: Empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.
Automatic: Wash the reagent wells 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) remaining in the reagent wells after washing can interfere with the substrate and lead to false low extinction values.

Insufficient washing (e.g. less than 3 wash cycles, too small wash buffer volumes, or too short residence times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step) Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgM) into each of the microplate wells. Incubate for **30 minutes** at room temperature (+18°C to +25°C).

Washing: Empty the wells. Wash as described above.

Substrate incubation: (3rd step) Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes** at room temperature (+18°C to +25°C), protect from direct sunlight.

Stopping: Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement: **Photometric measurement** of the colour intensity should be made at a **wavelength of 450 nm** and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, carefully shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using an analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on enquiry.

Automated test performance using other fully automated, open system analysis devices is possible. However, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example of the **semiquantitative determination** of antibodies in 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example of the **quantitative determination** of antibodies in 24 patient samples (P 1 to P 24).

Calibrators (C 1 to C 3), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimises reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator 2. Calculate the ratio according to the following formula:


$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive

Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibrators against the corresponding units (linear/linear). Use “point-to-point” plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.

If the extinction for a patient sample lies above the value of calibrator 1 (200 RU/ml), the result should be reported as “>200 RU/ml”. It is recommended that the sample be retested at a dilution of e.g. 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4. The upper limit of the reference range of non-infected persons (cut-off value) recommended by EUROIMMUN is 20 relative units RU/ml. EUROIMMUN recommends interpreting results as follows:

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For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another, EUROIMMUN recommends to retest the samples.

A negative serological result does not exclude an infection. Particularly in the early phase of an infection, antibodies may not yet be present or are only present in such small quantities that they are not detectable. In case of a borderline result, a secure evaluation is not possible. If there is a clinical suspicion and a negative test result, we recommend clarification by means of other diagnostic methods and/or the serological investigation of a follow-up sample. A positive result indicates that there has been contact with the pathogen. In the determination of pathogen-specific IgM antibodies, polyclonal stimulation of the immune system or antibody persistence may affect the diagnostic relevance of positive findings. Significant titer increases (exceeding factor 2) and/or seroconversion in a follow-up sample taken after 7-10 days can indicate an acute infection. To investigate titer changes, sample and follow-up sample should be incubated in adjacent wells of the ELISA microplate within the same test run. For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Calibration: As no international reference serum exists for antibodies against *Borrelia*, the calibration is performed in relative units (RU)/ml.

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratio values determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The binding activity of the antibodies and the activity of the enzyme used are temperature-dependent. It is therefore recommended using a thermostat in all three incubation steps. The higher the room temperature (+18°C to +25°C) during the incubation steps, the greater will be the extinction values. Corresponding variations apply also to the incubation times.

However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is provided by particularly suitable *Borrelia* strains (*Borrelia burgdorferi sensu stricto*, *Borrelia afzelii* and *Borrelia garinii*). The cultured bacteria were solubilised using sodium dodecyl sulphate. The used antigen mixture contains all relevant proteins.

Linearity: The linearity of the Anti-*Borrelia* ELISA (IgM) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95. The Anti-*Borrelia* ELISA (IgM) is linear at least in the tested concentration range (3 RU/ml to 177 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-*Borrelia* ELISA (IgM) is 1.3 RU/ml.

Cross reactivity: Cross reactivity of the Anti-*Borrelia* ELISA (IgM) was evaluated in a study performed on 263 patient sera with antibodies against *Treponema pallidum* and 18 patient sera with antibodies against *Leptospira*. 26 sera of the Anti-*Treponema pallidum*-positive sera and 3 sera of the Anti-*Leptospira*-positive sera presented reactivity. Therefore, the cross reactivity amounts to 9.9% and 16.7%, respectively. Cross reactions with antibodies against other spirochaetes cannot be excluded. The Anti-*Borrelia* ELISA is designed as a screening test which should be followed by a confirmatory test (Western- or line blot). Further sera from patients with different possible influencing factors were investigated. An overview of these results can be found in the following table.

Possible influencing factors	n	Anti- <i>Borrelia</i> ELISA (IgM) positive
Acute EBV infection	27	0%
Anti-CMV	19	0%
Anti-Measles virus	13	0%
Anti-Mumps virus	12	0%
Anti- <i>Toxoplasma gondii</i>	14	0%
Anti-VZV	15	0%
Anti-Rubella virus	10	0%
Anti-HSV	5	0%
Anti-HAV	7	0%
Anti-HBV	8	0%
Anti-HAMA	4	0%
ANA + DNS (AAb)	45	0%
Rheuma factor	39	0%
Neurological diseases	54	1.9%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation (CV) using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	75	4.3
2	103	2.9
3	115	3.4

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	78	6.9
2	109	6.2
3	122	4.7

Specificity and sensitivity: Sera from 150 patients with suspected borreliosis were analysed using the EUROIMMUN Anti-Borrelia ELISA (IgM) and the EUROIMMUN Anti-Borrelia burgdorferi Westernblot (IgM) as a reference method. The test showed a specificity of 96.4% and a sensitivity of 100%.

n = 150		EUROIMMUN Westernblot		
		positive	borderline	negative
ELISA	positive	34	5	4
	negative	0	1	106

Clinical study: 138 patients with clinically characterised borreliosis in different disease stages were screened with the Anti-Borrelia plus VisE-ELISA (IgG) and Anti-Borrelia-ELISA (IgM). Sensitivities between 89-96% were found.

n = 138		EUROIMMUN Anti-Borrelia plus VisE-ELISA (IgG), EUROIMMUN Anti-Borrelia-ELISA (IgM)		
		IgG positive	IgM positive	IgG und/oder IgM positive
Clinic	erythema migrans, n = 97	60	68	86 (88.7%)
	mono- or poly- arthritis, n = 26	22	5	25 (96.2%)
	neuroborreliosis, n = 15	13	4	14 (93.3%)

All ELISA negative results were confirmed as negative by an Anti-Borrelia Westernblot.

Reference range: The levels of the anti-Borrelia antibodies (IgM) were analysed with this EUROIMMUN ELISA in 500 healthy blood donors. With a cut-off of 20 RU/ml, 1.6% of the blood donors were anti-Borrelia positive (IgM).

Clinical significance

The history of Lyme disease, a contagious condition caused by *Borrelia burgdorferi* and transmitted to man by ticks, offers infectiologists a formidable lesson on how medicine progresses. Clinical description started in Europe at the turn of the 19th/20th century with Pick's description of what was then labelled chronic atrophic acrodermatitis. Fifty years later, Hauser noted the affection was transmitted by ticks. Independently, Afzelius, then Lipschutz, described erythema migrans and its relationship with tick bites. Neurological involvement was also described with the skin signs. These early dermatological descriptions suddenly came into the limelight in 1975 when an epidemia of arthritis occurred in children in Lyme, Connecticut, USA. Many of the affected children had erythema migrans. Based on these observations and an epidemiological analysis of the epidemia, Steele and co-workers defined "Lyme disease" as a rheumatological disorder commonly associated with erythema migrans and sometimes with multiple organ involvement. In 1982 Burgdorfer suggested that ticks transmitted "Treponema-like spirochaetes", which were later authenticated as the causal agent: *Borrelia burgdorferi*.

Borrelia burgdorferi as a species of bacteria is belonging to the spirochaetaceae family. This family includes five genera: *borrelia*, *spirochaetes*, *cristispira*, *treponema* and *leptospira*. Among about 30 tick species feeding on humans, *Ixodes ricinus* is the most frequent tick species biting humans in Europe. It is the vector of *Borrelia burgdorferi*, which causes Lyme disease, of *Anaplasma phagocytophilum* and the tick-borne encephalitis virus. *I. ricinus* ticks pass through three developmental stages: larvae, nymphs and adults (females and males). The density of this tick species can be very high, reaching in some places more than 300 ticks/100 m². Lyme borreliosis is the most frequent tick-transmitted disease in the northern hemisphere. In Europe, especially in Central Europe and Scandinavia there are up to 155 cases per 100,000 individuals caused by the species *B. burgdorferi sensu stricto*, *B. afzelii*, and *B. garinii* (in rare cases also by *B. spielmanii* and *B. bavariensis*). The human seroprevalence rate of antibodies against *Borrelia burgdorferi* in the normal population of Germany and other Central European countries is about 8% (for IgG); in highly endemic areas even considerably higher with ranges from 18% to 52%. In East Asia, e.g. in an endemic area in China, the seroprevalence amounts to 26%. Persons working in the area of forestry display anti-*Borrelia* antibodies in about 40%, hunters in more than 50% of cases.

A *Borrelia burgdorferi* infection can manifest itself in the areas of dermatology, neurology, ophthalmology, rheumatology and internal medicine. The clinical expression of borreliosis can be divided into three stages:

Stage I: The typical primary manifestation of a *Borrelia burgdorferi* infection is erythema migrans (80%), a reddening of the skin which appears around the area of the tick bite and spreads in a circular manner. The erythema is accompanied by influenza-like general symptoms with fever, shaking chill, headaches and vomiting. Lymphadenopathies are observed in a few cases (lymphadenosis cutis benigna). The leading clinical types of the erythema-free borreliosis (20%) are neurological, arthromyalgic, influenza-like, cardiovascular borreliosis and borreliosis with hepatitis, with regional lymphadenitis and mixtures of these. Stage I can result in spontaneous healing or can develop into a generalised borreliosis. The transition phase is generally symptom-free. IgM antibodies against *Borrelia burgdorferi* can be detected serologically in 50% to 90% of patients during stage I. Humans produce highly specific antibodies against the outer surface protein C (OspC) shortly after infection with *Borrelia burgdorferi*. The main antigen VlsE, which is exclusively expressed *in vivo*, is the most sensitive antigen for IgG antibody detection, whereas OspC is the most sensitive antigen for IgM antibody detection.

Stage II: A variety of symptoms can develop several weeks to some months after receiving the tick bite. In the foreground of these are neurological manifestations: meningitis, encephalitis, asymmetric polyneuritis, cranial nerve pareses, Bannwarth's lymphocytic meningoradiculoneuritis. One of the most frequent neurological symptoms of borreliosis in children is the peripheral facial palsy. Usually the paresis of the facial nerve appears on one side and after several weeks on the opposite side. Also arthritis, particularly of the knee joints, and non-localised pains in the bones, joints and muscles are frequently found. Less frequently are cardiological manifestations such as myocarditis and pericarditis. Antibodies against *Borrelia burgdorferi* can be detected in 50% to 90% of patients in stage II. In the early phase of this stage mainly IgM antibodies are found but in the late phase often only IgG antibodies occur. However, IgM antibodies can persist for a long time. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% to 30%. VlsE displays the highest sensitivity of all antigens tested for IgG detection.

Stage III: The typical manifestations of a *Borrelia burgdorferi* infection in stage III are chronic-relapsing erosive arthritis, acrodermatitis chronica atrophicans and progressive encephalomyelitis, which can proceed in a similar way to multiple sclerosis. Without treatment, the tertiary stage can develop over a period of years to decades after the original infection. In this stage, IgG antibody titers are significantly increased in 90% to 100% of patients, while IgM antibodies are only rarely detectable.

For the diagnosis of neuroborreliosis, the determination of antibodies in cerebrospinal fluid is of decisive importance and far superior to serological testing. An antibody index (AI) discriminates between a blood-derived and a pathological, brain-derived specific antibody fraction in cerebrospinal fluid (CSF) and takes into account individual changes in the blood/CSF barrier function.

Various techniques come into question for the detection of antibodies against *Borrelia burgdorferi*: ELISA, indirect immunofluorescence (IIFT), passive haemagglutination and immunoblot. Generally, ELISA or IIFT are used as screening test for preliminary characterisation of the serum sample. According to guidelines in the USA and in Germany, serological diagnosis should follow the principle of a two-step procedure: ELISA or indirect immunofluorescence (IIFT) as a first step, which, if reactive, is followed by an immunoblot. In fresh infections it is recommended performing ELISA/IIFT and immuno-blot in parallel, since some weak reactions become detectable earlier in the blot than in the screening test. VlsE, the most sensitive antigen for IgG antibody detection, and OspC for IgM antibody detection should be included in the test. By additionally determining antibodies against VlsE the serological hit rate can be increased by 20% compared to whole extract Westernblots. Of all recombinant antigens tested, VlsE possesses the highest sensitivity for the detection of a *Borrelia* infection.

Literature references

1. Burgdorfer W, Barbour AG, Hayes SF, Benach JL, Grunwaldt E, Davis JP. **Lyme disease - a tick-borne Spirochetosis?** Science 216 (1982) 1317-1319.
2. Christova I. **Enzyme-linked immunosorbent assay, immunofluorescent assay, and recombinant immunoblotting in the serodiagnosis of early Lyme borreliosis.** Int J Immunopathol Pharmacol 16 (2003) 261-268.
3. EUROIMMUN AG. Komorowski L, Janssen A, Stöcker W, Probst C. **Dimerisation of recombinant OspC leads to an antigen with enhanced potential for active vaccination.** 12. International Conference on Lyme Borreliosis and other Tick-Borne Diseases, Ljubljana, Slowenien (2010).
4. EUROIMMUN AG. Komorowski L, Probst C, Janssen A, Stöcker W. **Polypeptide und Verfahren zur spezifischen Detektion von Antikörpern bei Patienten mit einer Borrelien-infektion.** Europäische Patentanmeldung EP2199303 (2008).
5. EUROIMMUN AG. Meyer W, Janssen A, Scheper T, Schlumberger W, Stöcker W. **Lyme borreliosis: Prevalence of antibodies against non-proteinic (lipid) antigens.** Int J Med Microbiol 298S2 (Suppl. 45) (2008) 39.
6. EUROIMMUN AG. Probst C, Ott A, Scheper T, Meyer W, Stöcker W, Komorowski L. **N-Terminal disulfide-bridging of Borrelia outer surface protein C increases its diagnostic and vaccine potentials.** Ticks and Tick-Borne Diseases 3 (2012) 1-7.
7. EUROIMMUN AG. Meyer W, Scheper T, Schlumberger W, Stöcker W. **Antibodies Against the Newly Identified Major Antigen VlsE: A milestone in the Serological Diagnosis of Lyme Borreliosis.** Poster zum EUREGIO. Congress of Clinical Chemistry and Laboratory Medicine, Aachen (2003).
8. EUROIMMUN AG. Steinhagen K, Schlumberger W, Stöcker W. **Nachweis einer spezifischen intrathekalen Antikörpersynthese mit modernen ELISA-Testsystemen: Hohe Trefferquote bei Multipler Sklerose und Neuroborreliose.** J Lab Med 25 (2001) 135-149.
9. EUROIMMUN AG. Meyer W, Janssen A, Kühn B, Scheper T, Suer W, Komorowski L, Probst C, Schlumberger W, Stöcker W. **EUROLINE Borrelia RN-AT: Marriage of authentic recombinant and native borrelia antigens in a line immunoassay format.** International Journal of Medical Microbiology 298S2 (Suppl. 45): 39-40 (2008).
10. Goettner G, Schulte-Spechtel U, Hillermann R, Liegl G, Wilske B, Fingerle V. **Improvement of Lyme borreliosis serodiagnosis by a newly developed recombinant immunoglobulin G (IgG) and IgM line immunoblot assay and addition of VlsE and DbpA homologues.** J Clin Microbiol 43 (2005) 3602-3609.
11. Hauser U, Lehnert G, Lobentanzer R, Wilske B. **Interpretation Criteria for Standardized WesternBlots for Three European Species of Borrelia burgdorferi Sensus Lato.** J Clin Microbiol 35 (1997) 1433-1444.

12. Johnson R, Schmid GP, Hyde FW, Steigerwalt AG, Brenner DJ. **Borrelia burgdorferi sp. Nov.: Etiologic agent of lyme disease.** Int J Syst Bacteriol 34 (1984) 496-497.
13. Oehme R, Hartelt K, Backe H, Brockmann S, Kimmig P. **Foci of tick-borne diseases in south-west Germany.** Int J Med Microbiol 291 (2002) 22-29.
14. Robert-Koch-Institut. **Surveillance der Lyme-Borreliose am Beispiel des Bundeslandes Brandenburg.** Epidemiologisches Bulletin 14 (1998) 93-97.
15. Schulte-Spechtel U, Lehnert G, Liegl G, Fingerle V, Heimerl C, Johnson BJ, Wilske B. **Significant improvement of the recombinant Borrelia-specific immunoglobulin G immuno-blot test by addition of VlsE and a DbpA homologue derived from Borrelia garinii for diagnosis of early neuroborreliosis.** J Clin Microbiol 41 (2003) 1299-1303.
16. Tilton RC, Ryan RW. **The laboratory diagnosis of Lyme disease.** J Clin Immunoassay 16 (1993) 208-214.
17. Wilske B, Fingerle V, Schulte-Spechtel U. **Microbiological and serological diagnosis of Lyme borreliosis.** FEMS Immunol Med Microbiol 49 (2007) 13-21.
18. Wilske B. **Epidemiology and diagnosis of Lyme borreliosis.** Ann Med 37 (2005) 568-579.
19. Wilske B, Zöller L, Brade V, Eiffert H, Göbel UB, Stanek G, Pfister HW. **MiQ 12/2000 Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik: Lyme-Borreliose.** Urban & Fischer Verlag München Jena (2000).
20. Zhang JR, Hardham JM, Barbour AG, Norris SJ. **Antigenic Variation in Lyme Disease Borreliae by Promiscuous Recombination of VMP-like Sequence Cassettes.** Cell 89 (1997) 275-285.

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declares under its sole responsibility as manufacturer that the Blot product:

Anti-Echinococcus EUROLINE WB (IgG)

DY 2321-####-1 G

(product name, order number)

meet the following demands of:

Directive 98/79/EC on in vitro diagnostic medical devices of 27 October 1998 and its transpositions in national laws which apply to it.

Conformity assessment procedure: Annex III

This Declaration of Conformity is valid based on the respective currently valid version of technical documentation.

Lübeck, 25.04.2022

(Place and date of issue)

Dr. Ewald Müller-Kunert
- Head of Quality Management
and Regulatory Affairs -

Susanne Aleksandrowicz
- Member of the Executive Board -

Anti-Echinococcus EUROLINE WB (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
DY 2321-1601-1 G	Echinococcus antigen extract (whole antigen) plus Em18, Em95 and EgAgB	IgG	Antigen coated membrane strips	16 x 01 (16)

Indication: The EUROLINE-WB test kit provides a qualitative in vitro assay for human antibodies of the immunoglobulin class IgG against Echinococcus granulosus and Echinococcus multilocularis antigens in serum or plasma for the diagnosis of echinococcosis.

Application: The Anti-Echinococcus EUROLINE-WB (IgG) is excellently suited for result confirmation of serological screening tests for the detection of Echinococcus-specific antibodies (IHA, IIFT, ELISA). The use of an Echinococcus whole antigen extract and biochemically produced single antigens of Echinococcus granulosus and Echinococcus multilocularis in the Anti-Echinococcus EUROLINE-WB (IgG) also enables differentiation between cystic and alveolar echinococcosis, caused by Echinococcus granulosus and multilocularis, respectively.

Test principle: The test kit contains blot strips with electrophoretically separated antigen extract of Echinococcus. Additionally, each blot has a membrane chip coated with the biochemically produced antigens Em18, Em95 and EgAgB. In the first reaction step the blot strips are incubated with diluted patient samples. In the case of positive samples, specific antibodies of class IgG (and IgA, IgM) will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Format	Symbol
1. Test strips Single strips with electrophoretically separated Echinococcus antigen extract plus Em18, Em95 and EgAgB	16 x 1	STRIPS
2. Evaluation matrix with control strip Test strip incubated with a positive control serum	1 pattern	---
3. Enzyme conjugate Alkaline phosphatase-labelled anti-human IgG (goat), 10x concentrate	1 x 3 ml	CONJUGATE 10x
4. Universal buffer 10x concentrate	1 x 50 ml	BUFFER 10x
5. Substrate solution Nitroblue tetrazolium chloride/5-Bromo-4-chloro-3-indolylphosphate (NBT/BCIP), ready for use	1 x 30 ml	SUBSTRATE
6. Incubation tray	2 x 8 channels	TRAY
7. Test instruction	1 booklet	---
LOT Lot description		 Storage temperature
IVD In vitro diagnostic medical device		 Unopened usable until

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C, do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Undiluted patient sera and incubated blot strips should be handled as infectious waste. Other reagents do not need to be collected separately, unless stated otherwise in official regulations.

The following components are not provided in the test kits but can be ordered at EUROIMMUN under the respective order numbers. Performance of the test requires an **incubation tray**:

ZD 9895-0130 Incubation tray with 30 channels (black)

ZD 9898-0144 Incubation tray with 44 channels (black, for the EUROBlotOne and EUROBlotCamera system)

For the creation of work protocols and the evaluation of incubated test strips using **EUROLineScan** green paper and adhesive foil are required:

ZD 9880-0101 Green paper (1 sheet)

ZD 9885-0116 Adhesive foil for approx. 16 test strips

ZD 9885-0130 Adhesive foil for approx. 30 test strips

If a **visual evaluation** is to be performed in individual cases, the required evaluation protocol can be ordered under: ZD 2321-0101-1 Evaluation protocol visual Anti-Echinococcus EUROLINE-WB.

Preparation and stability of the reagents

Note: The bag containing the blot strips is printed with a number in addition to the test kit lot number. This number refers to the strip batch and is also printed on the corresponding evaluation template. These two numbers must match to ensure correct evaluation of test results.

All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. Unopened, reagents are stable until the indicated expiry date when stored at +2°C to +8°C. After initial opening, reagents are stable for 12 months or until the expiry date, if earlier, unless stated otherwise in the instructions. Opened reagents must also be stored at +2°C to +8°C and protected from contamination.

- **Coated test strips:** Ready for use. Open the packing with the test strips only when the strips have reached room temperature to prevent condensation on the strips. After removal of the strips the packing should be sealed tightly and stored at +2°C to +8°C. To ensure correct evaluation of results, the lot number on the bag must match the lot number on the strips as well as on the evaluation matrix.
- **Enzyme conjugate:** The enzyme conjugate is supplied as a 10x concentrate. For the preparation of the ready for use enzyme conjugate the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with ready for use diluted universal buffer. For 1 test strip dilute 0.15 ml anti-human IgG concentrate with 1.35 ml ready for use diluted universal buffer. The ready for use diluted enzyme conjugate should be used at the same working day.
- **Universal buffer:** The universal buffer is supplied as a 10x concentrate. For the preparation of the ready for use universal buffer the amount required should be removed from the bottle using a clean pipette and diluted 1:10 with deionised or distilled water. For the incubation of 1 test strip 1.5 ml buffer concentrate should be diluted with 13.5 ml deionised or distilled water. The ready for use diluted universal buffer should be used at the same working day.
- **Substrate solution:** Ready for use. Close bottle immediately after use, as the contents are sensitive to light ☀.

Warning: The control of human origin has tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the agent sodium azide in a non-declarable concentration. Avoid skin contact.

Preparation and stability of the serum or plasma samples

Sample material: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: The **patient samples** for analysis are diluted **1:51** in the ready for use diluted universal buffer. For example, add 30 µl of sample to 1.5 ml ready for use diluted universal buffer and mix well by vortexing. Sample pipettes are not suitable for mixing.

Incubation

Blocking: According to the number of serum samples to be tested fill each channel of the incubation tray with 1.5 ml ready for use diluted universal buffer and a blot strip. Remove the required amount of blot strips from the packing using a pair of tweezers. The number on the test strip should be visible. Incubate for **15 minutes** at room temperature (+18°C to +25°C) on a rocking shaker. Afterwards aspirate off all the liquid.

Sample incubation: (1st step) Fill each channel with 1.5 ml of the diluted serum samples and incubate at room temperature (+18°C to +25°C) for **30 minutes** on a rocking shaker.

Wash: Aspirate off the liquid from each channel and wash **3 x 5 minutes** each with 1.5 ml working strength universal buffer on a rocking shaker.

Conjugate incubation: (2nd step) Pipette 1.5 ml ready for use diluted enzyme conjugate (alkaline phosphatase-conjugated anti-human IgG) into each channel and incubate for **30 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.


Wash: Aspirate off the liquid from each channel. Wash as described above.

Substrate incubation: (3rd step) Pipette 1.5 ml substrate solution into the channels of the incubation tray. Incubate for **10 minutes** at room temperature (+18°C to +25°C) on a rocking shaker.

Stopping: Aspirate off the liquid from each channel and wash each strip **3 x 1 minute** with deionised or distilled water.

For automated incubation with the EUROBlotMaster select the program **Euro02 Inf WB 30**.

For automated incubation with the EUROBlotOne select the program **Euro 01/02**.

**Anti-Echinococcus EUROLINE-WB (IgG)****Incubation protocol**

Evaluation and Interpretation of the results of the Anti-Echinococcus EUROLINE-WB (IgG)

Handling: For evaluation of incubated test strips we generally recommend using the **EUROLineScan** software. After stopping the reaction using deionised or distilled water, place the incubated test strips onto the adhesive foil of the green work protocol using a pair of tweezers. The position of the test strips can be corrected while they are wet. As soon as all test strips have been placed onto the protocol, they should be pressed hard using filter paper and left to air-dry. After they have dried, the test strips will be stuck to the adhesive foil. The dry test strips are then scanned with a flatbed scanner (EUROIMMUN AG) and evaluated with **EUROLineScan**. Alternatively, imaging and evaluation is possible directly from the incubation trays (EUROBlotCamera and EUROBlotOne). For general information about the EUROLineScan program please refer to the EUROLineScan user manual (EUROIMMUN AG). The code for entering the **Test** in EUROLineScan is **Echino_EL-WB_IgG**.


If a visual evaluation must be performed in exceptional cases, hold the evaluation matrix next to the stuck-on blot strips and position it so that the black band above the number on the blot strips lines up with the alignment bar of the evaluation matrix. **The lot number on the evaluation matrix must match the lot number on the blot strips.** Clearly recognisable bands on the blot strips which concur with the labelled bands on the evaluation matrix are noted in the evaluation protocol.

Antigens: The antigen source for the **EUROIMMUN Anti-Echinococcus EUROLINE-WB** is provided by a particularly suitable Echinococcus antigen preparation. This antigen has been solubilised using sodium dodecyl sulphate followed by a separation of the solubilised protein using discontinuous polyacrylamide gel electrophoresis according to molecular mass and transfer of the separated proteins to nitrocellulose. Additional membrane strips containing the antigens Em 95, Em 18 and EgAgB are applied to each nitrocellulose membrane.

Control blot strips were taken from each nitrocellulose membrane and incubated with a reference serum. One of these stained strips is included in the kit.

Specificity of the antigens on the test strips:

Band	Antigen	Specificity
24-26 kDa	p25/26	Unspecific
21 kDa	p21	Specific for Echinococcus and other parasites
16-18 kDa	p16/18	Specific for Echinococcus
7 kDa	p7	Specific for Echinococcus
Em 95	Em 95	Specific for Echinococcus multilocularis
Em 18	Em 18	Specific for Echinococcus multilocularis
EgAgB	EgAgB	Specific for Echinococcus

In the lower part of the test strip there is a conjugate control membrane chip (IgA, IgG and IgM). Below the conjugate control, there is a membrane chip with a control band (Control).

Attention: A correctly performed determination of antibodies of class IgG against the antigens described above is indicated by a positive reaction of the control band and a positive reaction of the IgG band.

If one of these bands only shows a very weak reaction or none at all, the result is not valid.

IgG class antibodies against Echinococcus

In order to evaluate the signals, the band position and the intensity of the staining must be taken into account, as negative sera can also in some cases produce weak signals at individual band positions.

The defined antigens on the Anti-Echinococcus EUROLINE-WB can be divided into 5 categories:

Category	Antigens
1	antigen: p25/26
2	genus-specific Echinococcus antigen: EgAgB
3	Echinococcus antigen: p21
4	Echinococcus antigens: p7 and p16/18
5	Echinococcus multilocularis antigens: Em18 and Em95

Result interpretation: The results obtained with the Anti-Echinococcus EUROLINE-WB can be classified into negative, borderline and positive.

Result	Characteristics
Negative	No bands, or one positive antigen band of category 1, or one borderline antigen band of category 2.
Borderline	Positive antigen band of category 2, or at least one borderline antigen band of category 3, 4 or 5. It is recommended that a fresh sample be taken and the test be repeated after a few weeks.
Positive	At least one positive antigen band of category 3, 4 or 5. If one antigen band of category 3 is positive, either alone or in combination with an antigen band of category 1, cross reactivity may have occurred due to an Ascaris or Anisakis infection.

The Anti-Echinococcus EUROLINE-WB allows a serological differentiation between an infection with Echinococcus granulosus and Echinococcus multilocularis in many cases. The differentiation is automatically made during evaluation with EUROLIneScan.

For differentiation the following applies:

Result	Characteristics
Echinococcus granulosus	Positive antigen band of category 2 and additionally at least one positive antigen band of category 3 or 4.
Echinococcus multilocularis	At least one positive antigen band of category 5. Further positive antigen bands from the other categories may also occur.

For diagnosis, the clinical picture of the patient always needs to be taken into account along with the serological findings.

Test characteristics

Measurement range: The EUROLINE-WB is a qualitative method. No measurement range is provided.

Inter- and intra-assay variation: The inter-assay variation was determined by multiple analyses of characterised samples over several days. The intra-assay variation was determined by multiple analyses of characterised samples on one day. In every case, the intensity of the bands was within the specified range. This EUROLINE-WB displays excellent inter- and intra-assay reproducibility.

Interference: Haemolytic, lipaemic and icteric sera showed no effect on the analytical results.

Cross reactions: The cross reactivity of the EUROLINE-WB was investigated using patient samples with the following parasite infections:

Species	Number of samples	Number of positive results
<i>Ascaris lumbricoides</i>	10	3
<i>Anisakis simplex</i>	16	1
Filarioidea	5	0
<i>Strongyloides stercoralis</i>	10	0
<i>Schistosoma</i> ssp.	9	0
Multi-Helminth infection	2	0
<i>Plasmodium</i> ssp.	7	0
<i>Toxocara canis</i>	10	0
<i>Taenia solium</i>	7	0
<i>Trichinella spiralis</i>	13	0
<i>Fasciola hepatica</i>	17	0
<i>Entamoeba histolytica</i>	11	0
<i>Leishmania</i> ssp.	5	0

Specificity and sensitivity: A panel of 107 defined patient samples with a positive Echinococcus result and a control panel with serum samples from 50 healthy blood donors and 50 tumour patients were analysed using the Anti-Echinococcus EUROLINE-WB (IgG). The samples had been provided by the Institute of Parasitology of the University of Bern, Switzerland.

n = 207		Characterisation: Institute of Parasitology, University of Bern	
		positive	negative
Anti-Echinococcus EUROLINE-WB (IgG)	positive	99	0
	negative	8	100

The Anti-Echinococcus EUROLINE-WB (IgG) has a specificity of 100% and a sensitivity of 93%.

Of the 99 Echinococcus samples that were positive in the Anti-Echinococcus EUROLINE-WB (IgG) 80 samples could be differentiated further:

Echinococcus IgG positive: n = 99	Number of samples
<i>Echinococcus granulosus</i>	47
<i>Echinococcus multilocularis</i>	33
No differentiation possible	19

The rate of differentiation between *Echinococcus granulosus* and *Echinococcus multilocularis* is 81%.

Clinical significance

Echinococcosis is an infectious disease caused by parasites of the genus *Echinococcus*. In Europe, the dog tapeworm (*E. granulosus*), causing cystic echinococcosis (CE), and the fox tapeworm (*E. multilocularis*), causing alveolar echinococcosis (AE), are most important from the medical point of view.

The development of all *Echinococcus* species includes an obligatory host change: the definitive hosts are carnivores, the intermediate hosts mostly herbivores. The mature worms of the fox tapeworm, *E. multilocularis*, live in the intestine of their definitive hosts (in Europe, mainly red foxes, seldom dogs and cats) who release the mature cestode eggs with their droppings. These eggs are very robust and may be infectious for several months if the conditions are good. The intermediate hosts (small mammals like field mice or European water voles) take the eggs in with their food. The larvae of the worms then develop in their inner organs (mostly in the liver).

Humans can be infected with *E. granulosus* by smear infection, through dealing with contaminated soil or consuming contaminated foods. Also here, the eggs are dispersed with the definitive hosts' droppings (mostly dogs, at times cats) and stay infectious for months. The larvae develop in liquid-filled blisters (hyatides) which are found in liver, lung, other organs and also in the skeletal system of the intermediate host (e.g. hooved animals such as cows or sheep).

Humans are incidental hosts for fox and dog tapeworms. Human infection usually takes place by intake of cestodes after contact with infected hosts (e.g. smear infections or via the animal's fur).

The clinical image of both echinococcoses (CE and AE) differs in the different developmental behaviours of both parasites in the human body. The clinical appearance of AE corresponds to that of a malignoma. After haematogenous transport of the cestodes into the liver, infection of the liver takes place (usually unnoticed for a long time) and an alveolar tumour develops.

In infections with *E. granulosus*, larvae are released in the human intestine. They are transported haematogenously via the portal vein first into the liver and then into other organs, e.g. the lungs. The clinical course may vary significantly and is characterised by the slowly growing cysts and their different localisation.

In humans, both diseases remain asymptomatic for many years, until after 10 to 15 years they show symptoms like cholestatic icterus, epigastric pains, fatigue, weight loss and hepatomegaly. By invasion and destruction of the healthy liver tissue, an untreated echinococcosis may lead to the patient's death. Differential diagnosis from cysts, malignant and benign tumours, abscesses and the distinction between AE and CE are important for a diagnosis.

Imaging techniques, such as sonography, CT and MRT are used for diagnosis. The use of serological test systems for the detection of parasite-specific antibodies in serum or plasma helps to confirm results from imaging techniques. If whole antigens of *Echinococcus* are used, an echinococcosis can be detected with good sensitivity by ELISA or IFT tests.

The use of species-specific antigens in Westernblot and ELISA tests enables in many cases the serological differentiation of *E. granulosus* and *E. multilocularis*.

A negative result in this serological investigation does not exclude an infection. In liver CE, serum antibodies can be detected in 80 to 94% of cases. In lung echinococcosis, the prevalence amounts to only 65 to 70%.

In practice, molecular detection methods (PCR for the detection of *Echinococcus* DNA and RNA) have not proven to deliver.

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Seekamp 31, D-23560 Lübeck, Germany

declare under our sole responsibility that the ELISA products

Anti-West Nile Virus ELISA (IgG)

EI 2662-9601 G

Anti-West Nile Virus ELISA (IgM)

EI 2662-9601 M

Avidity: Anti-West Nile Virus ELISA (IgG)

EI 2662-9601-1 G

(product name, order no)

meet the demands of

*Directive 98/79/EC on in vitro diagnostic medical devices
of 27 October 1998 and its transpositions in national laws which apply to it.*

Conformity assessment procedure: Annex III

Lübeck, 13.10.2016

(Place and date of issue)

Susanne Aleksandrowicz
- Member of the Board -




Dr. Wolfgang Schlumberger
- Member of the Board -

Anti-West Nile Virus ELISA (IgG) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2662-9601 G	West Nile virus	IgG	Ag-coated microplate wells	96 x 01 (96)

Principle of the test: The ELISA test kit provides a semiquantitative or quantitative in vitro assay for human antibodies of the IgG class against West Nile virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with West Nile virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgG (also IgA and IgM) antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgG (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	
2. Calibrator 1 200 RU/ml (IgG, human), ready for use	dark red	1 x 2.0 ml	
3. Calibrator 2 20 RU/ml (IgG, human), ready for use	red	1 x 2.0 ml	
4. Calibrator 3 2 RU/ml, (IgG, human), ready for use	light red	1 x 2.0 ml	
5. Positive control (IgG, human), ready for use	blue	1 x 2.0 ml	
6. Negative control (IgG, human), ready for use	green	1 x 2.0 ml	
7. Enzyme conjugate peroxidase-labelled anti-human IgG (rabbit), ready for use	green	1 x 12 ml	
8. Sample buffer ready for use	light blue	1 x 100 ml	
9. Wash buffer 10x concentrate	colourless	1 x 100 ml	
10. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	
11. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	
12. Test instruction	---	1 booklet	
13. Quality control certificate	---	1 protocol	
14. Protective foil	---	2 pieces	
 Lot description		 Storage temperature	
 In vitro diagnostics		 Unopened usable until	

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

The thermostat adjusted ELISA incubator must be set at 37°C +/- 1°C.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrators and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallisation occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: The controls used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Sample dilution: Patient samples are diluted **1:101** in sample buffer. For example: dilute 10 µl serum in 1.0 ml sample buffer and mix well by vortexing (sample pipettes are not suitable for mixing).

NOTE: Calibrators and controls are prediluted and ready for use, do not dilute them.

Incubation

For **semiquantitative analysis** incubate **calibrator 2** along with the positive and negative controls and patient samples. For **quantitative analysis** incubate **calibrators 1, 2 and 3** along with the positive and negative controls and patient samples.

(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrator, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.

For manual processing of microplate wells, cover the finished test plate with the protective foil. When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate **60 minutes at 37°C ± 1°C**.

Washing:

Manual: Remove the protective foil and empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Remove the protective foil and empty the wells and subsequently wash 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgG) into each of the microplate wells.

When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate **30 minutes at room temperature (+18°C to +25°C)**

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3rd step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes at room temperature (+18°C to +25°C)** (protect from direct sunlight).

Stopping the reaction:

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the EUROIMMUN Analyzer I, Analyzer I-2P or the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry. Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C 2	P 6	P 14	P 22			C 1	P 4	P 12	P 20		
B	pos.	P 7	P 15	P 23			C 2	P 5	P 13	P 21		
C	neg.	P 8	P 16	P 24			C 3	P 6	P 14	P 22		
D	P 1	P 9	P 17				pos.	P 7	P 15	P 23		
E	P 2	P 10	P 18				neg.	P 8	P 16	P 24		
F	P 3	P 11	P 19				P 1	P 9	P 17			
G	P 4	P 12	P 20				P 2	P 10	P 18			
H	P 5	P 13	P 21				P 3	P 11	P 19			

The pipetting protocol for microtiter strips 1-4 is an example for the **semiquantitative analysis** of 24 patient samples (P 1 to P 24).

The pipetting protocol for microtiter strips 7-10 is an example for the **quantitative analysis** of 24 patient samples (P 1 to P 24).

The calibrators (C 1 to C 3), the positive (pos.) and negative (neg.) controls, and the patient samples have each been incubated in one well. The reliability of the ELISA test can be improved by duplicate determinations for each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of the calibrator 2. Calculate the ratio according the following formula:

$$\frac{\text{Extinction of the control or patient sample}}{\text{Extinction of calibrator 2}} = \text{Ratio}$$


EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive

In cases of borderline test results, an additional patient sample should be taken 7 days later and re-tested in parallel with the first patient sample. The results of both samples allow proper evaluation of titer changes.

Quantitative: The standard curve from which the concentration of antibodies in the patient samples can be taken is obtained by point-to-point plotting of the extinction values measured for the 3 calibration sera against the corresponding units (linear/linear). Use "point-to-point" plotting for calculation of the standard curve by computer. The following plot is an example of a typical calibration curve. Please do not use this curve for the determination of antibody concentrations in patient samples.

If the extinction of a serum sample lies above the value of calibrator 1 (200 RU/ml). The result should be given as ">200 RU/ml". It is recommended that the sample be re-tested at a dilution of 1:400. The result in RU/ml read from the calibration curve for this sample must then be multiplied by a factor of 4.

The upper limit of the normal range of non-infected persons (**cut-off value**) recommended by EUROIMMUN is **20 relative units (RU)/ml**. EUROIMMUN recommends interpreting results as follows:

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Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Test characteristics

Calibration: As no quantificated international reference serum exists for antibodies against West Nile virus, the calibration is performed in relative units (RU).

For every group of tests performed, the extinction values of the calibrators and the relative units and/or ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the controls are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibration sera are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is a recombinant, detergent-extracted glycoprotein E of WNV from the membrane fraction of human cells.

Linearity: The linearity of the Anti-West Nile Virus ELISA (IgG) was determined by assaying 4 serial dilutions of different patient samples. The coefficient of determination R^2 for all sera was > 0.95 . The Anti-West Nile Virus ELISA (IgG) is linear at least in the tested concentration range (10 RU/ml to 160 RU/ml).

Detection limit: The lower detection limit is defined as the mean value of an analyte-free samples plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-West Nile Virus ELISA (IgG) is 0.4 RU/ml.

Cross reactivity: Cross reactivities to other flaviviruses cannot be excluded. They were recognized with anti-TBE positive and anti-Dengue virus positive samples.

Antibodies against	n	Anti-West Nile Virus ELISA (IgG)
Adenovirus	12	0%
Chlamydia pneumoniae	12	0%
CMV	12	0%
EBV-CA	12	0%
Helicobacter pylori	12	0%
HSV-1	12	0%
Influenza virus A	12	0%
Influenza virus B	12	0%
Measles virus	12	0%
Mumps virus	12	0%
Mycoplasma pneumoniae	12	0%
Parainfluenza virus Pool	12	0%
RSV	12	0%
Rubella virus	12	0%
Toxoplasma gondii	9	0%
VZV	12	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (RU/ml)	CV (%)
1	101	6.1
2	104	3.3
3	149	2.7

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (RU/ml)	CV (%)
1	103	5.1
2	120	4.9
3	171	4.2

Sensitivity and specificity: Samples from 295 patients (origin: Europe) were investigated using the EUROIMMUN Anti-West Nile Virus ELISA and a neutralization test (NT) (performed by RKI, Berlin) as a reference method. The specificity was 96.9%, with a sensitivity of 99.5%. Values for 4 of the samples were borderline and were not included in the calculation.

Reference range: The levels of anti-West Nile virus antibodies (IgG) were analyzed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off of 20 RU/ml, 1.0% of the blood donors were anti-West Nile virus positive (IgG).

Clinical significance

West Nile virus (WNV) is an enveloped single-stranded RNA virus of the Flaviviridae family [1]. This family comprises around 100 virus types that are presently categorized into the three known species Flavivirus, Pestivirus and Hepacivirus [1, 2, 3, 4, 5]. West Nile virus received its name in 1937 when it was first isolated from a blood sample of an elderly woman living in the West Nile district in Uganda, who had fever of unknown cause accompanied by neurological disorders [6]. Further isolates were achieved only in 1951 from the sera of children with weak, unspecific symptoms, namely in Egypt where the virus is endemically distributed. At that time mice and embryonated hen's eggs were used for virus detection [1].

WNV is present not only in tropical areas, but also in moderate climate regions [2, 3, 4, 5]. Significant epidemics were observed in 1951/52 and 1957 in Israel and 1974 and 1983/84 in South Africa [1]. In the mid 90's the virus changed its virulence causing an epidemic accumulation of WNV encephalitis in Algeria (1994), Rumania (1996/97), the Czech Republic (1997), the Democratic Republic of Congo (1998), Russia, North America (1999) and Israel (2000) [2, 3, 7, 8, 9, 10, 11, 12].

In the USA 149 infections with 18 cases of death were recorded from 1999 to 2001. In 2002 this number rose to 4156 infections and 284 deaths, in 2003 to 9858 infections and 262 deaths [7, 8, 10]. Currently the virus has been detected in seven Canadian provinces, in 48 USA states and in Mexico, as well as in Puerto Rico, the Dominican Republic, Jamaica, Guadeloupe and El Salvador [3, 9, 11].

Since 1958, when antibodies against the WNV were first detected in the sera of two Albanians, repeated outbreaks of West Nile fever have occurred in Southern and Eastern Europe and meanwhile also in Central and West Europe [8, 9, 10, 12, 14, 15]. Its emergence and rapid spread is credited to world climate change, long-distance travel and globalization of economic trade [9, 12, 13, 15, 16, 17]. Consequent monitoring of West Nile activity by controlling sera of exposed persons is essential [15, 17, 18, 19]. Seroprevalence studies in endemic regions have shown an infection spread of up to 40% [14]. West Nile virus is therefore the Flavivirus with the largest distribution area [2, 3, 4, 5, 9].

WNV is transmitted by a number of mosquitoes. In the Mediterranean region and in Africa mosquitoes of the *Culex univittatus* complex species are the main arthropod hosts, while in North America WNV could be detected in 37 mosquito species, with *Culex pipiens* being the main vector [1, 2, 3, 4, 7, 9]. In India *Culex vishnui* and in France *Culex modestus* were identified as the main vectors [1]. In total WNV could be found in more than 40 mosquito and in several tick species [3, 4, 5].

Birds represent the vertebrate reservoir [1, 20, 21]. Alone in the USA WNV has been found in more than 162 species of birds. Many of them showed clinical symptoms and thousands of birds died after contracting a natural infection [7, 8, 10, 20, 21]. Birds that survive develop lifelong immunity. Acting as coincidental hosts mammals can also become infected when bitten by an infected mosquito [10, 13, 14]. Transmission has additionally been documented via breast milk, bone marrow transplantations, liver and heart transplants, blood transfusions, lab accidents such as open wounds during handling of infected brain tissue as well as transplacental transmission [2, 3, 4, 5, 8, 22, 23, 24, 25, 26, 27, 28]. Other than humans, mostly only horses became ill after an infection [2, 13, 21]. As well as WNV infection via a mosquito bite, a second natural infection source is possible in animals, namely via feeding on infected prey [3, 5]. An experimental infection of cats was successfully achieved by feeding them infected mice [20, 21].

70% to 80% of the humans infected with WNV showed no symptoms [8]. In the remaining 20%-30%, signs of sudden flu-like symptoms appear after an incubation period of 2-6 days with fever ranging from 38.5 to 40°C lasting for 3-5 days, nausea, shivering, head and back aches, joint and muscle pain and other unspecific symptoms such as loss of appetite, dizziness, vomiting, diarrhoea, coughing and a sore throat [1, 2, 3, 4, 5].

Typical for epidemical occurring fever are exanthema on the breast, back and upper extremities and general lymph node swelling [1, 3, 26]. Severe clinical cases of WNV infections are characterized by myocarditis, pancreatitis and hepatitis and since 1996, also neurological disorders, as WNV is now capable of crossing the blood-brain barrier [3]. The neurological symptoms begin after a short febrile prodrome phase approximately 1-7 days after infection and become manifest in the form of encephalitis and meningoencephalitis accompanied by stiffness, spasms and shivering as the result of damage done to the basal ganglia [3, 8, 29, 30].

Another widespread symptom is general muscle weakness similar to the Guillain-Barré syndrome and also polio-like paralysis [29, 30]. Approximately 4%-14% of the hospitalized patient cases are fatal [29]. High risk factors are old age and a weak immune system [1, 3, 29].

An infection with WNV during pregnancy can cause miscarriages, congenital meningitis, birth defects in approx. 10% of the cases and in an additional 10% of newborns growth disturbances [3, 8, 27, 28].

The diagnosis of WNV can be performed by virus detection or by detection of specific antibodies [1, 4, 17, 18, 19]. As virus isolation from serum or cerebrospinal fluid or virus detection using Real-Time Reverse Transcription-Polymerase Chain Reaction (RT-PCR) is usually unsuccessful due to short viraemia and low virus titers, the detection of specific WNV antibodies using ELISA and IFA has gained importance [11, 17, 18, 19, 31, 32, 33, 34, 35, 36, 37, 38].

Specific IgM antibodies in serum can be determined using ELISA or IIFT [19]. Antibodies of class IgM are detectable in serum from the second day after initial symptoms of the illness occur. A four-fold increase in titer of the respective class of antibody is considered proof of a WNV infection.

If the IgM test is negative, even though the symptoms indicate a WNV infection, a second serum sample should be taken and tested for IgM antibodies a few days later. A combination of ELISA and IIFT provides close to 100% reliability [17, 18, 31, 32, 33, 34]. Anti-WNV IgM antibodies persist for 2 to 3 months, often for more than a year [4, 17, 18, 19, 22].

Antibodies of class IgG are detectable approx. 2 days after the appearance of IgM antibodies [11, 19, 32, 33, 34, 40]. Two to four weeks after a positive IgM result the infection can be confirmed and its severity and prognosis evaluated using a qualitative and quantitative test for the detection of specific WNV IgG antibodies in the patient serum [34].

For the reliable differentiation between acute and past infections the detection of low-avidity IgG antibodies gives evidence for a primary or an acute WNV infection, while high-avidity antibodies indicate a past or reactivated WNV infection [39, 40, 41]. EUROIMMUN offers additional test systems for determination of IgG avidity in both ELISA and IIFT formats. The detection of low-avidity antibodies using ELISA and IIFT in parallel is possible for WNV as it is for *Toxoplasma gondii*, rubella virus, EBV-EA, EBV-CA and Corona virus [39, 40, 41].

As the degree of similarity within the Flavivirus family is high antibody cross reactions can occur [19, 31, 42]. Therefore samples that are positive for specific IgM and/or IgG antibodies against WNV should be titrated and investigated on all relevant Flavivirus IIFT substrates for cross reactions. By comparing the titer strengths the initial result can be confirmed or disproved by the second detection and an infection with another Flavivirus identified as the source of illness [17, 32, 34].

To supplement and extend the current Anti-West Nile Virus ELISA and Anti-West Nile Virus IIFT (each IgG or IgM or avidity) BIOCHIP Mosaics and Profiles for the detection (IIFT) of infections with Flaviviruses and the BIOCHIP Mosaic Fever Profile 1: South-East Asia have been developed. With these tests specific antibodies (IgG and IgM) against several infectious agents can be investigated simultaneously [36, 37, 38, 43, 44, 45]. These supplementary tests allow similar or ambiguous disease symptoms and potential cross reactions to be clarified and differential diagnostic issues to be addressed [1].

A specific antiviral therapy for WNV encephalitis is not available at present [1, 2, 3, 4, 5, 14, 35]. Intensive medical care is the only possibility to positively influence the illness. Eradication of WNV is impossible due to the natural bird-mosquito cycle [1]. A vaccine with formalin inactivated WNV is only available for horses [2, 10]. Therefore public education, individual precautionary measures and protection against insect bites are essential contributions to preventing WNV infections [2, 3, 4, 5, 14, 16].

Literature references

1. Falke D. **Flavi-Viren**. In: Hahn H, Falke D, Kaufmann SHE, Ullmann U. **Medizinische Mikrobiologie und Infektiologie**. Springer-Verlag, 5. Auflage (2004) 523-529.
2. Campbell GL, Marfin AA, Lanciotti RS, Gubler DJ. **West Nile Virus**. *Lancet Infect Dis* 2 (2002) 519-529.
3. Nosal B, Pellizzari R. **West Nile virus**. *CMAJ* 168 (2003) 1443-1444.
4. Granwehr BP, Lillibridge KM, Higgs S, Mason PW, Aronson JF, Campbell GA, Barrett ADT. **West Nile virus: where are we now?** *The Lancet Infect Dis* 4 (2004) 547-551.
5. Stock I. **West nile virus. An unusual flavivirus with increased importance**. [Article in German] *Med Monatsschr Pharm* 27 (2004) 112-120.
6. Smithburn KC, Hughes TP, Burke AW et al. **A neurotropic virus isolated from the blood of a native of Uganda**. *Am J Trop Med Hyg* 20 (1940) 471-472.
7. Kilpatrick AM, Kramer LD, Jones MJ, Marra PP, Daszak P. **West Nile virus epidemics in North America are driven by shifts in mosquito feeding behavior**. *PLoS Bio* 4 (2006) 82.
8. Zeller HG, Schuffenecker I. **West Nile virus: an overview of its spreads in europe and the mediterranean basin in contrast to its spread in the americas**. *Eur J Clin Microbiol Infects Dis* 23 (2004) 147-156.
9. Rappole JH, Compton BW, Leimgruber P, Robertson J, King DI, Renner SC. **Modeling movement of West Nile virus in the Western hemisphere**. *Vector Borne Zoonotic Dis* 6 (2006) 128-139.
10. Hubalek Z. **European experience with the West Nile virus ecology and epidemiology: could it be relevant for the New World?** *Viral Immunol* 13 (2000) 415-426.
11. Tilley PA, Walle R, Chow A, Jayaraman GC, Fonseca K, Drebot MA, Preiksaitis J, Fox J. **Clinical utility of commercial enzyme immunoassays during the inaugural season of West Nile virus activity, Alberta, Canada**. *J Clin Microbiol* 43 (2005) 4691-4695.

12. Hubalek Z, Halouzka J. **West Nile Fever – a re-emerging mosquito-borne viral disease in europe.** Emerg Infect Dis 5 (1999) 643-650.
13. Pauli G. **West Nile virus. Prevalence and significance as a zoonotic pathogen.** [Article in German] Bundesgesundheitsblatt Gesundheitsforschung Gesundheitsschutz 47 (2004) 653-660.
14. Robert-Koch-Institut. **West-Nil-Fieber.** Epidemiol Bull 13 (2003) 96-99.
15. Epstein PR. **West Nile virus and the climate.** J Urban Health 78 (2001) 367-371.
16. Gubler DJ. **The global emergence/resurgence of arboviral diseases as public health problems.** Arch Med Res 33 (2002) 330-342.
17. Koraka P, Zeller H, Niedrig M, Osterhaus AD, Groen J. **Reactivity of serum from patients with flavivirus infection measured by immunofluorescence assay and ELISA.** Microbes Infect 12 (2002) 1209-1215.
18. Shi PY, Wong SJ. **Serologic diagnosis of West-Nile virus infection.** Expert Rev Mol Diagn 6 (2003) 733-741.
19. Tardei G, Ruta S, Chitu V, Rossi C, Tsai TF, Cernescu C. **Evaluation of immunoglobulin M (IgM) and IgG enzyme immunoassays in serologic diagnosis of West Nile virus Infection.** J. Clin. Microbiol 38 (2000) 2232-2239.
20. Rappole JH, Hubalek Z. **Migratory birds and West Nile virus.** J Appl Microbiol 94 (2003) 47-58.
21. Kilpatrick AM, Daszak P, Jones MJ, Marra PP, Kramer LD. **Host heterogeneity dominates West Nile virus transmission.** Proc Biol Sci 22 (2006) 2327-2333.
22. Prince HE, Tobler LH, Lape-Nixon M, Foster GA, Stramer SL, Busch MP. **Development and persistence of West Nile virus-specific immunoglobulin M (IgM), IgA, and IgG in viremic blood donors.** J Clin Microbiol 43 (2005) 4316-4320.
23. Iwamoto M, Jernigan DB, Guasch A, Trepka MJ, Blackmore CG, Hellinger WC, Pham SM, Zaki S, Lanciotti RS, Lance-Parker SE, DiazGranados CA, Winqvist AG, Perlino CA, Wiersma S, Hillyer KL, Goodman JL, Marfin AA, Chamberland ME, Petersen LR. **West Nile Virus in Transplant Recipients Investigation Team. Transmission of West Nile virus from an organ donor to four transplant recipients.** N Engl J Med 348 (2003) 2196-2203.
24. Montgomery SP, Brown JA, Kuehnert M, Smith TL, Crall N, Lanciotti RS, Macedo de Oliveira A, Boo T, Marfin AA, 2003 West Nile Virus Transfusion-Associated Transmission Investigation Team. **Transfusion-associated transmission of West Nile virus, United States 2003 through 2005.** Transfusion 46 (2006) 2038-2046.
25. Cameron C, Reeves J, Antonishyn N, Tilley P, Alport T, Towns D, Lane D, Saldanha J. **West Nile virus in Canadian blood donors.** Transfusion 45 (2005) 487-491.
26. Anderson RC, Horn KB, Hoang MP, Gottlieb E, Bennin B. **Punctate exanthem of West Nile Virus infection: report of 3 cases.** J Am Acad Dermatol 51 (2004) 820-823.
27. O'Leary DR, Kuhn S, Kniss KL, Hinckley AF, Rasmussen SA, Pape WJ, Kightlinger LK, Beecham BD, Miller TK, Neitzel DF, Michaels SR, Campbell GL, Lanciotti RS, Hayes EB. **Birth outcomes following West Nile Virus infection of pregnant women in the United States: 2003-2004.** Pediatrics 117 (2006) 537-545.
28. Paisley JE, Hinckley AF, O'Leary DR, Kramer WC, Lanciotti RS, Campbell GL, Hayes EB. **West Nile virus infection among pregnant women in a northern Colorado community, 2003 to 2004.** Pediatrics 117 (2006) 814-820.
29. Guarner J, Shieh WJ, Hunter S, Paddock CD, Morken T, Campbell GL, Marfin AA, Zaki SR. **Clinicopathologic study and laboratory diagnosis of 23 cases with West Nile virus encephalomyelitis.** Hum Pathol 35 (2004) 983-990.
30. Bouffard JP, Riudavets MA, Holman R, Rushing EJ. **Neuropathology of the brain and spinal cord in human West Nile virus infection.** Clin Neuropathol 23 (2004) 59-61.

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Anti-West Nile Virus ELISA (IgM) Test instruction

ORDER NO.	ANTIBODIES AGAINST	IG CLASS	SUBSTRATE	FORMAT
EI 2662-9601 M	West Nile Virus	IgM	Ag-coated microplate wells	96 x 01 (96)

Principle of the test: The ELISA test kit provides a semiquantitative in vitro assay for human antibodies of the IgM class against West Nile virus in serum or plasma. The test kit contains microtiter strips each with 8 break-off reagent wells coated with West Nile virus antigens. In the first reaction step, diluted patient samples are incubated in the wells. In the case of positive samples, specific IgM (also IgA and IgG) antibodies will bind to the antigens. To detect the bound antibodies, a second incubation is carried out using an enzyme-labelled anti-human IgM (enzyme conjugate) catalysing a colour reaction.

Contents of the test kit:

Component	Colour	Format	Symbol
1. Microplate wells coated with antigens: 12 microplate strips each containing 8 individual break-off wells in a frame, ready for use	---	12 x 8	STRIPS
2. Calibrator (IgM, human), ready for use	dark red	1 x 2.0 ml	CAL
3. Positive control (IgM, human), ready for use	blue	1 x 2.0 ml	POS CONTROL
4. Negative control (IgM, human), ready for use	green	1 x 2.0 ml	NEG CONTROL
5. Enzyme conjugate peroxidase-labelled anti-human IgM (goat), ready for use	red	1 x 12 ml	CONJUGATE
6. Sample buffer buffer containing IgG/RF-Absorbent (Anti-human IgG antibody preparation obtained from goat), ready for use	green	1 x 100 ml	SAMPLE BUFFER
7. Wash buffer 10x concentrate	colourless	1 x 100 ml	WASH BUFFER 10x
8. Chromogen/substrate solution TMB/H ₂ O ₂ , ready for use	colourless	1 x 12 ml	SUBSTRATE
9. Stop solution 0.5 M sulphuric acid, ready for use	colourless	1 x 12 ml	STOP SOLUTION
10. Test instruction	---	1 booklet	
11. Quality control certificate	---	1 protocol	
12. Protective foil	---	3 pieces	
LOT Lot description	CE	 Storage temperature	
IVD In vitro diagnostics		 Unopened usable until	

Storage and stability: The test kit has to be stored at a temperature between +2°C to +8°C. Do not freeze. Unopened, all test kit components are stable until the indicated expiry date.

Waste disposal: Patient samples, calibrators, controls and incubated microplate strips should be handled as infectious waste. All reagents must be disposed of in accordance with local disposal regulations.

Preparation and stability of the reagents

Note: All reagents must be brought to room temperature (+18°C to +25°C) approx. 30 minutes before use. After first use, the reagents are stable until the indicated expiry date if stored at +2°C to +8°C and protected from contamination, unless stated otherwise below.

The thermostat adjusted ELISA incubator must be set at 37°C +/- 1°C.

- **Coated wells:** Ready for use. Tear open the resealable protective wrapping of the microplate at the recesses above the grip seam. Do not open until the microplate has reached room temperature to prevent the individual strips from moistening. Immediately replace the remaining wells of a partly used microplate in the protective wrapping and tightly seal with the integrated grip seam (Do not remove the desiccant bag).
Once the protective wrapping has been opened for the first time, the wells coated with antigens can be stored in a dry place and at a temperature between +2°C and +8°C for 4 months.
- **Calibrator and controls:** Ready for use. The reagents must be mixed thoroughly before use.
- **Enzyme conjugate:** Ready for use. The enzyme conjugate must be mixed thoroughly before use.
- **Sample buffer:** Ready for use. The green coloured sample buffer contains IgG/RF absorbent. Serum or plasma samples diluted with this sample buffer are only to be used for the determination of IgM antibodies.
- **Wash buffer:** The wash buffer is a 10x concentrate. If crystallization occurs in the concentrated buffer, warm it to 37°C and mix well before diluting. The quantity required should be removed from the bottle using a clean pipette and diluted with deionized or distilled water (1 part reagent plus 9 parts distilled water).
For example: For 1 microplate strip, 5 ml concentrate plus 45 ml water.
The working strength diluted wash buffer is stable for 4 weeks when stored at +2°C to +8°C and handled properly.
- **Chromogen/substrate solution:** Ready for use. Close the bottle immediately after use, as the contents are sensitive to light. The chromogen/substrate solution must be clear on use. Do not use the solution if it is blue coloured.
- **Stop solution:** Ready for use.

Warning: The control sera used have been tested negative for HBsAg, anti-HCV, anti-HIV-1 and anti-HIV-2 using enzyme immunoassays and indirect immunofluorescence methods. Nonetheless, all materials should be treated as being a potential infection hazard and should be handled with care. Some of the reagents contain the toxic agent sodium azide. Avoid skin contact.

Preparation and stability of the patient samples

Samples: Human serum or EDTA, heparin or citrate plasma.

Stability: Patient samples to be investigated can generally be stored at +2°C to +8°C for up to 14 days. Diluted samples should be incubated within one working day.

Introduction: Before the determination of specific antibodies of class IgM, antibodies of class IgG should be removed from the patient sample. This procedure must be carried out in order to prevent any rheumatoid factors from reacting with specifically bound IgG, which would lead to false positive IgM test results, and to prevent specific IgG displacing IgM from the antigen, which would lead to false IgM negative test results.

Functional principle: The sample buffer (green coloured!) contains an anti-human antibody preparation from goat. IgG from a serum sample is bound with high specificity by these antibodies and precipitated. If the sample also contains rheumatoid factors, these will be absorbed by the IgG/anti-human IgG complex.

Separation properties:

- All IgG subclasses are bound and precipitated by the anti-human IgG antibodies.
- Human serum IgG in concentrations of up to 15 mg per ml are removed (average serum IgG concentration in adults: 12 mg per ml).
- Rheumatoid factors are also removed.
- The recovery rate of the IgM fraction is almost 100%.

Performance: The **patient samples** for analysis are diluted **1:101** with sample buffer. For example, add 10 µl serum to 1.0 ml sample buffer and mix well. Incubate the mixture for at least **10 minutes** at room temperature. Subsequently, it can be pipetted into the microplate wells according to the pipetting protocol.

Notes:

- Antibodies of the class IgG should not be analyzed with this mixture.
- It is possible to check the efficacy of the IgG/RF absorbent for an individual patient sample by performing an IgG test in parallel to the IgM test using the mixture. If the IgG test is negative, the IgM result can be considered as reliable.
- The calibrator and controls containing IgM antibodies are pre-diluted and ready for use, do not dilute them.

Incubation

(Partly) manual test performance

Sample incubation: (1st step)

Transfer 100 µl of the calibrator, positive and negative controls or diluted patient samples into the individual microplate wells according to the pipetting protocol.

For manual processing of microplate wells, cover the finished test plate with the protective foil. When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate **60 minutes at 37°C ± 1°C**.

Washing:

Manual: Remove the protective foil and empty the wells and subsequently wash 3 times using 300 µl of working strength wash buffer for each wash.

Automatic: Remove the protective foil and empty the wells and subsequently wash 3 times with 450 µl of working strength wash buffer (program setting: e.g. TECAN Columbus Washer "Overflow Modus").

Leave the wash buffer in each well for 30 to 60 seconds per washing cycle, then empty the wells. After washing (manual and automated tests), thoroughly dispose of all liquid from the microplate by tapping it on absorbent paper with the openings facing downwards to remove all residual wash buffer.

Note: Residual liquid (> 10 µl) in the reagent wells after washing can interfere with the substrate and lead to false low extinction values. Insufficient washing (e.g., less than 3 wash cycles, too small wash buffer volumes, or too short reaction times) can lead to false high extinction values.

Free positions on the microplate strip should be filled with blank wells of the same plate format as that of the parameter to be investigated.

Conjugate incubation: (2nd step)

Pipette 100 µl of enzyme conjugate (peroxidase-labelled anti-human IgM) into each of the microplate wells.

When using an automated microplate processor for incubation, follow the instrument manufacturer's recommendations with regard to microwell plate sealing.

Incubate **30 minutes at room temperature (+18°C to +25°C)**.

Washing:

Empty the wells. Wash as described above.

Substrate incubation: (3. step)

Pipette 100 µl of chromogen/substrate solution into each of the microplate wells. Incubate for **15 minutes at room temperature (+18°C to +25°C)** (protect from direct sunlight).

Stopping the reaction: (3rd step)

Pipette 100 µl of stop solution into each of the microplate wells in the same order and at the same speed as the chromogen/substrate solution was introduced.

Measurement:

Photometric measurement of the colour intensity should be made at a wavelength of 450 nm and a reference wavelength between 620 nm and 650 nm **within 30 minutes of adding the stop solution**. Prior to measuring, slightly shake the microplate to ensure a homogeneous distribution of the solution.

Test performance using fully automated analysis devices

Sample dilution and test performance are carried out fully automatically using the analysis device. The incubation conditions programmed in the respective software authorised by EUROIMMUN may deviate slightly from the specifications given in the ELISA test instruction. However, these conditions were validated in respect of the combination of the Analyzer I, Analyzer I-2P and the DSX from Dynex and this EUROIMMUN ELISA. Validation documents are available on inquiry. Automated test performance using other fully automated, open system analysis devices is possible, however, the combination should be validated by the user.

Pipetting protocol

	1	2	3	4	5	6	7	8	9	10	11	12
A	C	P 6	P 14	P 22								
B	pos.	P 7	P 15	P 23								
C	neg.	P 8	P 16	P 24								
D	P 1	P 9	P 17									
E	P 2	P 10	P 18									
F	P 3	P 11	P 19									
G	P 4	P 12	P 20									
H	P 5	P 13	P 21									

The above pipetting protocol is an example of the **semiquantitative analysis** of antibodies in 24 patient samples (P 1 to P 24).

Calibrator (C), positive (pos.) and negative (neg.) control as well as the patient samples have been incubated in one well each. The reliability of the ELISA test can be improved by duplicate determinations of each sample.

The wells can be broken off individually from the strips. This makes it possible to adjust the number of test substrates used to the number of samples to be examined and minimizes reagent wastage.

Both positive and negative controls serve as internal controls for the reliability of the test procedure. They should be assayed with each test run.

Calculation of results

The extinction value of the calibrator defines the upper limit of the reference range of non-infected persons (**cut-off**) recommended by EUROIMMUN. Values above the indicated cut-off are to be considered as positive, those below as negative.

Semiquantitative: Results can be evaluated semiquantitatively by calculating a ratio of the extinction value of the control or patient sample over the extinction value of calibrator. Use the following formula to calculate the ratio:

$$\frac{\text{Extinction of the control/patient sample}}{\text{Extinction of calibrators}} = \text{Ratio}$$

EUROIMMUN recommends interpreting results as follows:

Ratio <0.8:	negative
Ratio ≥0.8 to <1.1:	borderline
Ratio ≥1.1:	positive

Evaluation information: For duplicate determinations the mean of the two values should be taken. If the two values deviate substantially from one another the sample should be retested.

For the interpretation of borderline results an investigation using further tests (e.g. avidity determination of antibody class IgG) can be helpful. Diagnosis can be secured by the determination of the titer change in two serum samples taken at an interval of at least 7 days and analysed in parallel.

For diagnosis, the clinical symptoms of the patient should always be taken into account along with the serological results.

Test characteristics

Calibration: As no quantificated international reference serum exists for antibodies of the IgM class against West Nile virus, results are provided in the form of ratios which are a relative measure for the concentration of antibodies.

For every group of tests performed, the extinction values of the calibrators and the ratios determined for the positive and negative controls must lie within the limits stated for the relevant test kit lot. A quality control certificate containing these reference values is included. If the values specified for the control sera are not achieved, the test results may be inaccurate and the test should be repeated.

The activity of the enzyme used is temperature-dependent and the extinction values may vary if a thermostat is not used. The higher the room temperature during substrate incubation, the greater will be the extinction values. Corresponding variations apply also to the incubation times. However, the calibrators are subject to the same influences, with the result that such variations will be largely compensated in the calculation of the result.

Antigen: The antigen source is a recombinant, detergent-extracted glycoprotein E of WNV from the membrane fraction of human cells.

Detection limit: The lower detection limit is defined as the mean value of an analyte-free sample plus three times the standard deviation and is the smallest detectable antibody titer. The lower detection limit of the Anti-West Nile ELISA (IgM) is a ratio value of 0.03.

Cross reactivity: Cross reactivities to other flaviviruses cannot be excluded. They were recognized with anti-TBE positive and anti-Dengue virus positive samples.

Antibodies against	n	Anti-West Nile ELISA (IgM)
Borrelia burgdoferi	9	0%
CMV	8	0%
EBV-CA	9	0%
HSV-1/2	2	0%
Measles virus	10	0%
Mumps virus	9	0%
Parvovirus B19	8	0%
Rubella virus	10	0%
Toxoplasma gondii	10	0%
VZV	4	0%

Interference: Haemolytic, lipaemic and icteric samples showed no influence on the result up to a concentration of 10 mg/ml for haemoglobin, 20 mg/ml for triglycerides and 0.4 mg/ml for bilirubin in this ELISA.

Reproducibility: The reproducibility of the test was investigated by determining the intra- and inter-assay coefficients of variation using 3 sera. The intra-assay CVs are based on 20 determinations and the inter-assay CVs on 4 determinations performed in 6 different test runs.

<i>Intra-assay variation, n = 20</i>		
Serum	Mean value (Ratio)	CV (%)
1	1.9	5.9
2	2.0	2.2
3	2.4	4.0

<i>Inter-assay variation, n = 4 x 6</i>		
Serum	Mean value (Ratio)	CV (%)
1	1.9	3.9
2	2.1	9.4
3	2.4	6.8

Sensitivity and specificity: Study I: For the determination of the sensitivity 18 clinically and serologically precharacterised sera (Robert Koch Institute, Berlin, Germany) were tested with the EUROIMMUN Anti-West Nile Virus ELISA (IgM). The sensitivity amounted to 94.4%.

Study II: 99 patient sera, characterised as positive at the Saskatchewan Disease Control Laboratory (Canada) using several serological methods, were investigated with the EUROIMMUN Anti-West Nile Virus ELISA (IgM). The sensitivity was 93.7%. For the determination of the specificity 500 blood donor samples were investigated. The specificity was 99.8%.

Reference range: The levels of anti-West Nile virus antibodies (IgM) were analyzed with this EUROIMMUN ELISA in a panel of 500 healthy blood donors. With a cut-off ratio of 1.0, 0.2% of the blood donors were anti-West Nile virus positive (IgM).

Clinical significance

West Nile virus (WNV) is an enveloped single-stranded RNA virus of the Flaviviridae family [1]. This family comprises around 100 virus types that are presently categorized into the three known species Flavivirus, Pestivirus and Hepacivirus [1, 2, 3, 4, 5]. West Nile virus received its name in 1937 when it was first isolated from a blood sample of an elderly woman living in the West Nile district in Uganda, who had fever of unknown cause accompanied by neurological disorders [6]. Further isolates were achieved only in 1951 from the sera of children with weak, unspecific symptoms, namely in Egypt where the virus is endemically distributed. At that time mice and embryonated hen's eggs were used for virus detection [1].

WNV is present not only in tropical areas, but also in moderate climate regions [2, 3, 4, 5]. Significant epidemics were observed in 1951/52 and 1957 in Israel and 1974 and 1983/84 in South Africa [1]. In the mid 90's the virus changed its virulence causing an epidemic accumulation of WNV encephalitis in Algeria (1994), Rumania (1996/97), the Czech Republic (1997), the Democratic Republic of Congo (1998), Russia, North America (1999) and Israel (2000) [2, 3, 7, 8, 9, 10, 11, 12]. In the USA 149 infections with 18 cases of death were recorded from 1999 to 2001. In 2002 this number rose to 4156 infections and 284 deaths, in 2003 to 9858 infections and 262 deaths [7, 8, 10]. Currently the virus has been detected in seven Canadian provinces, in 48 USA states and in Mexico, as well as in Puerto Rico, the Dominican Republic, Jamaica, Guadeloupe and El Salvador [3, 9, 11].

Since 1958, when antibodies against the WNV were first detected in the sera of two Albanians, repeated outbreaks of West Nile fever have occurred in Southern and Eastern Europe and meanwhile also in Central and West Europe [8, 9, 10, 12, 14, 15]. Its emergence and rapid spread is credited to world climate change, long-distance travel and globalization of economic trade [9, 12, 13, 15, 16, 17]. Consequent monitoring of West Nile activity by controlling sera of exposed persons is essential [15, 17, 18, 19]. Seroprevalence studies in endemic regions have shown an infection spread of up to 40% [14]. West Nile virus is therefore the Flavivirus with the largest distribution area [2, 3, 4, 5, 9].

WNV is transmitted by a number of mosquitoes. In the Mediterranean region and in Africa mosquitoes of the *Culex univittatus* complex species are the main arthropod hosts, while in North America WNV could be detected in 37 mosquito species, with *Culex pipiens* being the main vector [1, 2, 3, 4, 7, 9]. In India *Culex vishnui* and in France *Culex modestus* were identified as the main vectors [1]. In total WNV could be found in more than 40 mosquito and in several tick species [3, 4, 5].

Birds represent the vertebrate reservoir [1, 20, 21]. Alone in the USA WNV has been found in more than 162 species of birds. Many of them showed clinical symptoms and thousands of birds died after contracting a natural infection [7, 8, 10, 20, 21]. Birds that survive develop lifelong immunity. Acting as coincidental hosts mammals can also become infected when bitten by an infected mosquito [10, 13, 14]. Transmission has additionally been documented via breast milk, bone marrow transplantations, liver and heart transplants, blood transfusions, lab accidents such as open wounds during handling of infected brain tissue as well as transplacental transmission [2, 3, 4, 5, 8, 22, 23, 24, 25, 26, 27, 28]. Other than humans, mostly only horses became ill after an infection [2, 13, 21]. As well as WNV infection via a mosquito bite, a second natural infection source is possible in animals, namely via feeding on infected prey [3, 5]. An experimental infection of cats was successfully achieved by feeding them infected mice [20, 21].

70% to 80% of the humans infected with WNV showed no symptoms [8]. In the remaining 20%-30%, signs of sudden flu-like symptoms appear after an incubation period of 2-6 days with fever ranging from 38.5 to 40°C lasting for 3-5 days, nausea, shivering, head and back aches, joint and muscle pain and other unspecific symptoms such as loss of appetite, dizziness, vomiting, diarrhoea, coughing and a sore throat [1, 2, 3, 4, 5].

Typical for epidemical occurring fever are exanthema on the breast, back and upper extremities and general lymph node swelling [1, 3, 26]. Severe clinical cases of WNV infections are characterized by myocarditis, pancreatitis and hepatitis and since 1996, also neurological disorders, as WNV is now capable of crossing the blood-brain barrier [3]. The neurological symptoms begin after a short febrile prodrome phase approximately 1-7 days after infection and become manifest in the form of encephalitis and meningoencephalitis accompanied by stiffness, spasms and shivering as the result of damage done to the basal ganglia [3, 8, 29, 30].

Another widespread symptom is general muscle weakness similar to the Guillain-Barré syndrome and also polio-like paralysis [29, 30]. Approximately 4%-14% of the hospitalized patient cases are fatal [29]. High risk factors are old age and a weak immune system [1, 3, 29].

An infection with WNV during pregnancy can cause miscarriages, congenital meningitis, birth defects in approx. 10% of the cases and in an additional 10% of newborns growth disturbances [3, 8, 27, 28].

The diagnosis of WNV can be performed by virus detection or by detection of specific antibodies [1, 4, 17, 18, 19]. As virus isolation from serum or cerebrospinal fluid or virus detection using Real-Time Reverse Transcription-Polymerase Chain Reaction (RT-PCR) is usually unsuccessful due to short viraemia and low virus titers, the detection of specific WNV antibodies using ELISA and IFA has gained importance [11, 17, 18, 19, 31, 32, 33, 34, 35, 36, 37, 38].

Specific IgM antibodies in serum can be determined using ELISA or IIFT [19]. Antibodies of class IgM are detectable in serum from the second day after initial symptoms of the illness occur. A four-fold increase in titer of the respective class of antibody is considered proof of a WNV infection.

If the IgM test is negative, even though the symptoms indicate a WNV infection, a second serum sample should be taken and tested for IgM antibodies a few days later. A combination of ELISA and IIFT provides close to 100% reliability [17, 18, 31, 32, 33, 34]. Anti-WNV IgM antibodies persist for 2 to 3 months, often for more than a year [4, 17, 18, 19, 22].

Antibodies of class IgG are detectable approx. 2 days after the appearance of IgM antibodies [11, 19, 32, 33, 34, 40]. Two to four weeks after a positive IgM result the infection can be confirmed and its severity and prognosis evaluated using a qualitative and quantitative test for the detection of specific WNV IgG antibodies in the patient serum [34].

For the reliable differentiation between acute and past infections the detection of low-avidity IgG antibodies gives evidence for a primary or an acute WNV infection, while high-avidity antibodies indicate a past or reactivated WNV infection [39, 40, 41].

EUROIMMUN offers additional test systems for determination of IgG avidity in both ELISA and IIFT formats. The detection of low-avidity antibodies using ELISA and IIFT in parallel is possible for WNV as it is for *Toxoplasma gondii*, rubella virus, EBV-EA, EBV-CA and Corona virus [39, 40, 41].

As the degree of similarity within the Flavivirus family is high antibody cross reactions can occur [19, 31, 42]. Therefore samples that are positive for specific IgM and/or IgG antibodies against WNV should be titrated and investigated on all relevant Flavivirus IIFT substrates for cross reactions. By comparing the titer strengths the initial result can be confirmed or disproved by the second detection and an infection with another Flavivirus identified as the source of illness [17, 32, 34].

To supplement and extend the current Anti-West Nile Virus ELISA and Anti-West Nile Virus IIFT (each IgG or IgM or avidity) BIOCHIP Mosaics and Profiles for the detection (IIFT) of infections with Flaviviruses and the BIOCHIP Mosaic Fever Profile 1: South-East Asia have been developed. With these tests specific antibodies (IgG and IgM) against several infectious agents can be investigated simultaneously [36, 37, 38, 43, 44, 45]. These supplementary tests allow similar or ambiguous disease symptoms and potential cross reactions to be clarified and differential diagnostic issues to be addressed [1].

A specific antiviral therapy for WNV encephalitis is not available at present [1, 2, 3, 4, 5, 14, 35]. Intensive medical care is the only possibility to positively influence the illness. Eradication of WNV is impossible due to the natural bird-mosquito cycle [1]. A vaccine with formalin inactivated WNV is only available for horses [2, 10]. Therefore public education, individual precautionary measures and protection against insect bites are essential contributions to preventing WNV infections [2, 3, 4, 5, 14, 16].

Literature references

1. Falke D. **Flavi-Viren**. In: Hahn H, Falke D, Kaufmann SHE, Ullmann U. **Medizinische Mikrobiologie und Infektiologie**. Springer-Verlag, 5. Auflage (2004) 523-529.
2. Campbell GL, Marfin AA, Lanciotti RS, Gubler DJ. **West Nile Virus**. *Lancet Infect Dis* 2 (2002) 519-529.
3. Nosal B, Pellizzari R. **West Nile virus**. *CMAJ* 168 (2003)1443-1444.
4. Granwehr BP, Lillibridge KM, Higgs S, Mason PW, Aronson JF, Campbell GA, Barrett ADT. West Nile virus: where are we now? *The Lancet Infect Dis* 4 (2004) 547-551.
5. Stock I. **West nile virus. An unusual flavivirus with increased importance**. [Article in German] *Med Monatsschr Pharm* 27 (2004) 112-120.
6. Smithburn KC, Hughes TP, Burke AW et al. **A neurotropic virus isolated from the blood of a native of Uganda**. *Am J Trop Med Hyg* 20 (1940) 471-472.
7. Kilpatrick AM, Kramer LD, Jones MJ, Marra PP, Daszak P. **West Nile virus epidemics in North America are driven by shifts in mosquito feeding behavior**. *PLoS Bio* 4 (2006) 82.
8. Zeller HG, Schuffenecker I. **West Nile virus: an overview of its spreads in europe and the mediterranean basin in contrast to its spread in the americas**. *Eur J Clin Microbiol Infects Dis* 23 (2004) 147-156.
9. Rappole JH, Compton BW, Leimgruber P, Robertson J, King DI, Renner SC. **Modeling movement of West Nile virus in the Western hemisphere**. *Vector Borne Zoonotic Dis* 6 (2006) 128-139.
10. Hubalek Z. **European experience with the West Nile virus ecology and epidemiology: could it be relevant for the New World?** *Viral Immunol* 13 (2000) 415-426.
11. Tilley PA, Walle R, Chow A, Jayaraman GC, Fonseca K, Drebot MA, Preiksaitis J, Fox J. **Clinical utility of commercial enzyme immunoassays during the inaugural season of West Nile virus activity, Alberta, Canada**. *J Clin Microbiol* 43 (2005) 4691-4695.
12. Hubalek Z, Halouzka J. **West Nile Fever – a re-emerging mosquito-borne viral disease in europe**. *Emerg Infect Dis* 5 (1999) 643-650.

13. Pauli G. **West Nile virus. Prevalence and significance as a zoonotic pathogen.** [Article in German] Bundesgesundheitsblatt Gesundheitsforschung Gesundheitsschutz 47 (2004) 653-660.
14. Robert-Koch-Institut. **West-Nil-Fieber.** Epidemiol Bull 13 (2003) 96-99.
15. Epstein PR. **West Nile virus and the climate.** J Urban Health 78 (2001) 367-371.
16. Gubler DJ. **The global emergence/resurgence of arboviral diseases as public health problems.** Arch Med Res 33 (2002) 330-342.
17. Koraka P, Zeller H, Niedrig M, Osterhaus AD, Groen J. **Reactivity of serum from patients with flavivirus infection measured by immunofluorescence assay and ELISA.** Microbes Infect 12 (2002) 1209-1215.
18. Shi PY, Wong SJ. **Serologic diagnosis of West-Nile virus infection.** Expert Rev Mol Diagn 6 (2003) 733-741.
19. Tardei G, Ruta S, Chitu V, Rossi C, Tsai TF, Cernescu C. **Evaluation of immunoglobulin M (IgM) and IgG enzyme immunoassays in serologic diagnosis of West Nile virus Infection.** J. Clin. Microbiol 38 (2000) 2232-2239.
20. Rappole JH, Hubalek Z. **Migratory birds and West Nile virus.** J Appl Microbiol 94 (2003) 47-58.
21. Kilpatrick AM, Daszak P, Jones MJ, Marra PP, Kramer LD. **Host heterogeneity dominates West Nile virus transmission.** Proc Biol Sci 22 (2006) 2327-2333.
22. Prince HE, Tobler LH, Lape-Nixon M, Foster GA, Stramer SL, Busch MP. **Development and persistence of West Nile virus-specific immunoglobulin M (IgM), IgA, and IgG in viremic blood donors.** J Clin Microbiol 43 (2005) 4316-4320.
23. Iwamoto M, Jernigan DB, Guasch A, Trepka MJ, Blackmore CG, Hellinger WC, Pham SM, Zaki S, Lanciotti RS, Lance-Parker SE, DiazGranados CA, Winquist AG, Perlino CA, Wiersma S, Hillyer KL, Goodman JL, Marfin AA, Chamberland ME, Petersen LR. **West Nile Virus in Transplant Recipients Investigation Team. Transmission of West Nile virus from an organ donor to four transplant recipients.** N Engl J Med 348 (2003) 2196-2203.
24. Montgomery SP, Brown JA, Kuehnert M, Smith TL, Crall N, Lanciotti RS, Macedo de Oliveira A, Boo T, Marfin AA, 2003 West Nile Virus Transfusion-Associated Transmission Investigation Team. **Transfusion-associated transmission of West Nile virus, United States 2003 through 2005.** Transfusion 46 (2006) 2038-2046.
25. Cameron C, Reeves J, Antonishyn N, Tilley P, Alport T, Towns D, Lane D, Saldanha J. **West Nile virus in Canadian blood donors.** Transfusion 45 (2005) 487-491.
26. Anderson RC, Horn KB, Hoang MP, Gottlieb E, Bennin B. **Punctate exanthem of West Nile Virus infection: report of 3 cases.** J Am Acad Dermatol 51 (2004) 820-823.
27. O'Leary DR, Kuhn S, Kniss KL, Hinckley AF, Rasmussen SA, Pape WJ, Kightlinger LK, Beecham BD, Miller TK, Neitzel DF, Michaels SR, Campbell GL, Lanciotti RS, Hayes EB. **Birth outcomes following West Nile Virus infection of pregnant women in the United States: 2003-2004.** Pediatrics 117 (2006) 537-545.
28. Paisley JE, Hinckley AF, O'Leary DR, Kramer WC, Lanciotti RS, Campbell GL, Hayes EB. **West Nile virus infection among pregnant women in a northern Colorado community, 2003 to 2004.** Pediatrics 117 (2006) 814-820.
29. Guarner J, Shieh WJ, Hunter S, Paddock CD, Morken T, Campbell GL, Marfin AA, Zaki SR. **Clinicopathologic study and laboratory diagnosis of 23 cases with West Nile virus encephalomyelitis.** Hum Pathol 35 (2004) 983-990.
30. Bouffard JP, Riudavets MA, Holman R, Rushing EJ. **Neuropathology of the brain and spinal cord in human West Nile virus infection.** Clin Neuropathol 23 (2004) 59-61.

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Design, development and manufacturing of reagents, reagent products, calibrators and control materials for immunochemistry, microbiology, infectious immunology and molecular biology techniques.

Design and development, management of production and technical servicing of instruments and software for "in vitro" diagnostic.

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ORGANISMO NOTIFICADO 0318

Ea IgG

**Enzyme ImmunoAssay (ELISA) for
the quantitative/qualitative
determination of IgG antibodies to
Epstein Barr Virus Early Antigen
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

Ea IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Epstein Barr Virus Early Antigen in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Epstein Barr Virus or EBV is the principal etiological agent of infectious mononucleosis, as well as a contributory factor in the etiology of Burkitt's lymphoma and nasopharyngeal carcinoma, or NPC. A member of the family Herpesviridae, it has a worldwide distribution, such that 80 to 90% of all adults have been infected. Primary infections usually occur during the first decade of life. While childhood infections are mostly asymptomatic, 50 to 70% of young adults undergoing primary EBV infections show mild to severe illness. EBV may cause a persistent, latent infection which can be reactivated under immunosuppression or in AIDS affected patients.

As humoral responses to primary EBV infections are quite rapid, the level and class of antibodies raised in most cases allow classification as to whether the patient is still susceptible, has a current or recent primary infection, had a past infection or may be having reactivated EBV infection. The detection of EBV-specific IgG, IgM and IgA antibodies to its major immunodominant antigens (Nuclear Antigen, Viral Capsid Antigen, Early Antigen) has become therefore an important and useful determination for the monitoring and follow-up of EBV infected patients.

C. PRINCIPLE OF THE TEST

Microplates are coated with EBV-specific affinity purified Early Antigen or EA.

In the 1st incubation, the solid phase is treated with diluted samples and anti-EA IgG are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti-EA IgG are detected by the addition of anti hlgG antibody, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti EA IgG antibodies present in the sample.

IgG in the sample may be quantitated by means of a standard curve calibrated in arbitrary units per milliliter (Uarb/ml) as no international standard is available.

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1. Microplate: MICROPLATE

12 strips x 8 microwells coated with **affinity purified EBV Ea**. Plates are sealed into a bag with desiccant.

Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: CAL N° ...

6 vials. Ready to use and color coded standard curve ranging:

4 ml CAL1 = 0 arbU/ml
4 ml CAL2 = 5 arbU/ml
2 ml CAL3 = 10 arbU/ml
2 ml CAL4 = 20 arbU/ml
2 ml CAL 5 = 50 arbU/ml
4 ml CAL6 = 100 arbU/ml.

Standards are calibrated against an internal Gold Standard or IGS as no international one is defined.

Contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

4. Enzyme conjugate : CONJ

1x16ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

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7. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

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1. Calibrated Micropipettes (1000, 100 and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.

2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.

5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.

6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.

7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.

8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.

9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 6 months.

11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.

13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

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4. Sera and plasma can be stored at +2°.8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples

should not be freezed/thawed more than once as this may generate particles that could affect the test result.

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6. Samples whose anti-Ea IgG antibody concentration is expected to be higher than 100 arbU/ml should be diluted before use, either 1:10 or 1:100 in the Calibrator 0 arbU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 ul of each specimen with 450 ul of Cal 0 (1:10). Then 50 ul of the 1:10 dilution are diluted with 450 ul of the Cal 0 (1:100). Mix tubes thoroughly on vortex and then proceed toward the dilution step reported in section M.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of storing.

In this case call Dia.Pro's customer service.

Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°.8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the

requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dilute all the content of the 20x concentrated Wash Solution as described above.
6. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
7. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
8. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
9. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
10. Check that the micropipettes are set to the required volume.
11. Check that all the other equipment is available and ready to use.
12. In case of problems, do not proceed further with the test and advise the supervisor.

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M1. QUANTITATIVE DETERMINATION:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of Microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
3. Dispense 100 µl of Calibrators in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

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- Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 100 µl of Calibrator 0 arbU/ml and Calibrator 5 arbU/ml in duplicate and Calibrator 100 arbU/ml in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **60 min at +37°C**.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from yellow to blue.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-

630nm (background subtraction, strongly recommended), blanking the instrument on A1.

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N. ASSAY SCHEME

Method	Operations
Calibrators	100 µl
Samples diluted 1:101	100 µl
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
CAL 1 0 arbU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
CAL 2 5 arbU/ml OD450nm < OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 100 arbU/ml < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead) ; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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
Then on the calibration curve calculate the concentration of anti Eα IgG antibody in samples.

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Example of calculation:

Note: The following data must not be used instead of real figures obtained by the user.

Calibrator 0 arbU/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted
 Calibrator 5 arbU/ml: 0.250 – 0.270 OD450nm
 Mean Value: 0.260 OD450nm
 Higher than Cal 0 + 0.100 – Accepted
 Calibrator 100 arbU/ml: 2.045 OD450nm
 Higher than 1.000 – Accepted

The OD450nm of the Calibrator 5 arbU/ml is considered the cut-off (or Co) of the system.

The ratio between the OD450nm value of the sample and the OD450nm of the Calibrator 5 arbU/ml (or S/Co) can provide a semi-quantitative estimation of the content of specific IgG in the sample.

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 5 arbU/ml are considered negative for anti-Ea IgG antibody.

Samples with a concentration higher than 5 arbU/ml are considered positive for anti-Ea IgG antibody.

Ea IgG results alone are not, anyway, enough to provide a clear diagnosis of EBV infection.

At least EBV VCA IgG and EBV VCA IgM results, possibly together with EBNA IgG, are necessary in combination.

A reference range of the minimum essential serological markers of Epstein-Barr infection, derived from Infectious Diseases Handbook, 3rd edition, published by Lexi-Comp Inc., USA, is reported schematically below:

VCA IgM	EBNA (or VCA) IgG	Interpretation
negative	negative	No history of EBV infection
positive	negative	Acute primary infection
negative	positive	History of previous infection
positive	positive	Reactivation

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R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in an external clinical center on negative and positive samples with reference to a FDA approved commercial kit.

1. Limit of detection

No international standard for Ea IgG Antibody detection has been defined so far by the European Community.

In its absence, an Internal Gold Standard (or IGS), derived from a patient with an history of past mononucleosis infection, has been defined in order to provide the device with a constant and excellent sensitivity.

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The diagnostic performances were evaluated in a performance evaluation study conducted in an external centre, with excellent experience in the diagnosis of infectious diseases.

The diagnostic sensitivity was studied on samples, pre-tested positive with a different reference kit of European origin in use at the laboratory. Positive samples were collected from patients that experienced mononucleosis infection.

The diagnostic specificity was determined on panels of negative samples from normal individuals and blood donors, classified negative with the reference kit, including potentially interfering specimens.

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Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.

The Performance Evaluation provided the following values :

Sensitivity	≥ 98 %
Specificity	≥ 98 %

3. Reproducibility:

Data obtained from a study conducted on three samples of different Ea IgG reactivity, examined in 16 replicates in three separate runs show CV% values ranging 3-16% depending on OD450nm readings.

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S. LIMITATIONS

False positivity has been assessed as less than 2-5% of the normal population depending on the reference kit used.

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All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with EN ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	Ea IgG CODE: EAG.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – JUNE 2010
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

EBNA IgG

**Enzyme ImmunoAssay (ELISA) for
the quantitative/qualitative
determination of IgG antibodies to
Epstein Barr Virus Nuclear Antigen
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

EBNA IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Epstein Barr Virus Nuclear Antigen in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Epstein Barr Virus or EBV is the principal etiological agent of infectious mononucleosis, as well as a contributory factor in the etiology of Burkitt's lymphoma and nasopharyngeal carcinoma, or NPC. A member of the family Herpesviridae, it has a worldwide distribution, such that 80 to 90% of all adults have been infected. Primary infections usually occur during the first decade of life. While childhood infections are mostly asymptomatic, 50 to 70% of young adults undergoing primary EBV infections show mild to severe illness. EBV may cause a persistent, latent infection which can be reactivated under immunosuppression or in AIDS affected patients. As humoral responses to primary EBV infections are quite rapid, the level and class of antibodies raised in most cases allow classification as to whether the patient is still susceptible, has a current or recent primary infection, had a past infection or may be having reactivated EBV infection. The detection of EBV-specific IgG, IgM and IgA antibodies to its major immunodominant antigens (mainly Nuclear Antigen or EBNA and Viral Capsidic Antigen or VCA) has become therefore an important and useful determination for the monitoring and follow-up of EBV infected patients.

C. PRINCIPLE OF THE TEST

In order to get rid of crossreactions with other viruses of the same family, microplates are coated with affinity purified native EBNA antigen, capable to provide the assay with the highest specificity.

In the 1st incubation, the solid phase is treated with diluted samples and anti-EBNA IgG are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti-EBNA IgG are detected by the addition of anti hIgG antibody, labeled with peroxidase (HRP). The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti EBNA IgG antibodies present in the sample.

IgG in the sample may therefore be quantitated by means of a standard curve calibrated in arbitrary units per milliliter (arbU/ml) as no international standard is available.

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1. Microplate: MICROPLATE

12 strips x 8 microwells coated with affinity purified native EBNA antigen. Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: CAL N° ..

Ready to use and color coded standard curve ranging:
4 ml CAL1 = 0 arbU/ml
4 ml CAL2 = 5 arbU/ml
2 ml CAL3 = 10 arbU/ml
2 ml CAL4 = 20 arbU/ml
2ml CAL 5 = 50 arbU/ml
4 ml CAL6 = 100 arbU/ml.

Standards are calibrated against an internal Gold Standard or IGS as no international one is defined.

Contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized.

It contains fetal bovine serum proteins, human IgG antibodies to EBNA at 20 arbU/ml±20%, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

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4. Enzyme conjugate : CONJ

1x16ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

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7. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

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1. Calibrated Micropipettes (1000, 100 and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National

Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.

5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.

6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.

7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.

8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.

9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.

11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.

13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

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3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results.

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6. Samples whose anti-EBNA IgG antibody concentration is expected to be higher than 100 arbU/ml should be diluted before use, either 1:10 or 1:100 in the Calibrator 0 arbU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 ul of each specimen with 450 ul of Cal 0 (1:10). Then 50 ul of the 1:10 dilution are diluted with 450 ul of the Cal 0 (1:100). Mix tubes thoroughly on vortex and then proceed toward the dilution step reported in section M.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of storing. In this case call Dia.Pro's customer service.

Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2..8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Note: Once diluted, the wash solution is stable for 1 week at +2..8°C.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the Control Serum as reported.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.
12. Check that all the other equipment is available and ready to use.
13. In case of problems, do not proceed further with the test and advise the supervisor.

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M1. QUANTITATIVE DETERMINATION:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use.

Mix carefully all the liquid components on vortex and then proceed as described below.

- Place the required number of Microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
- Dispense 100 µl of Calibrators and 100 µl Control Serum in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

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- Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 100 µl of Calibrator 0 arbU/ml and Calibrator 10 arbU/ml in duplicate and Calibrator 100 arbU/ml in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **60 min at +37°C**.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from yellow to blue.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

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N. ASSAY SCHEME

Method	Operations
Calibrators & Control(*)	100 µl
Samples diluted 1:101	100 µl
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not got contaminated during the assay
CAL 1 0 arbU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.

CAL 2 5 arbU/ml OD450nm < OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 3 10 arbU/ml OD450nm < OD450nm CAL1 + 0.200	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 100 arbU/ml < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Problem	Check
Control Serum Different from Expected value	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the control has occurred.

Anyway, if all other parameters (Blank, CAL1, CAL2, CAL 6), match the established requirements, the test may be considered valid.

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
Then on the calibration curve calculate the concentration of anti EBNA IgG antibody in samples.

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P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 10 arbu/ml and then check that the assay is valid.

Example of calculation:

Note: The following data must not be used instead or real figures obtained by the user.

Calibrator 0 arbu/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted

Calibrator 10 arbu/ml: 0.450 – 0.470 OD450nm
 Mean Value: 0.460 OD450nm
 Higher than Cal 0 + 0.200 – Accepted

Calibrator 100 arbu/ml: 2.045 OD450nm
 Higher than 1.000 – Accepted

The OD450nm of the Calibrator 10 arbu/ml is considered the cut-off (or Co) of the system.

The ratio between the OD450nm value of the sample and the OD450nm of the Calibrator 10 arbu/ml (or S/Co) can provide a semi-quantitative estimation of the content of specific IgG in the sample.

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 5 arbu/ml are considered negative for anti EBNA IgG antibody.

Samples with a concentration ranging 5-10 arbu/ml are considered in the gray-zone. Samples with a concentration higher than 10 arbu/ml are considered positive for anti EBNA IgG antibody.

EBNA IgG results alone are not, anyway, enough to provide a clear diagnosis of EBV infection. At least EBV VCA IgM results are necessary in combination.

A reference range of the minimum essential serological markers of Epstein-Barr infection, derived from Infectious Diseases

Handbook, 3rd edition, published by Lexi-Comp Inc., USA, is reported schematically below:

VCA IgM	EBNA IgG	Interpretation
negative	negative	No history of EBV infection
positive	negative	Acute primary infection
negative	positive	History of previous infection
positive	positive	Reactivation

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R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in an external clinical center on negative and positive samples with reference to a FDA approved commercial kit.

1. Limit of detection

No international standard for EBNA IgG Antibody detection has been defined so far by the European Community.

In its absence, an Internal Gold Standard (or IGS), derived from a patient with an history of past mononucleosis infection, has been defined in order to provide the device with a constant and excellent sensitivity.

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The method is based on the use of an affinity purified native EBNA antigen to provide the assay with the highest specificity to EBV.

The diagnostic performances were evaluated in a performance evaluation study conducted in an external centre, with excellent experience in the diagnosis of infectious diseases and in particular in EBV infection.

The Diagnostic Sensitivity was studied on more than 50 samples, pre-tested positive with two reference kits of European origin in use at the laboratory. Positive samples were collected from patients that experienced mononucleosis infection.

The diagnostic specificity was determined on panels of more than 50 negative samples from normal individuals and blood donors, classified negative with the reference kit, including potentially interfering specimens.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity.

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Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.

The Performance Evaluation provided the following values :

Sensitivity	≥ 98 %
Specificity	≥ 98 %

3. Reproducibility:

Data obtained from a study conducted on three samples of different EBNA IgG reactivity, examined in 16 replicates in three separate runs show CV% values ranging 5-20% depending on OD450nm readings.

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Depending on the reference kit in use, due to some heterogeneity among different devices, the presence of 2-5% false reactivity may be seen.

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PLACE & DATE OF FIRST ISSUE	MILANO – SEPTEMBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	VCA IgG CODE: VCAG.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – SEPTEMBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

VCA IgG

**Enzyme ImmunoAssay (ELISA) for
the quantitative/qualitative
determination of IgG antibodies to
Epstein Barr Virus Capsidic Antigen
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italia**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

VCA IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Epstein Barr Virus Capsidic Antigen in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Epstein Barr Virus or EBV is the principal etiological agent of infectious mononucleosis, as well as a contributory factor in the etiology of Burkitt's lymphoma and nasopharyngeal carcinoma, or NPC. A member of the family Herpesviridae, it has a worldwide distribution, such that 80 to 90% of all adults have been infected. Primary infections usually occur during the first decade of life. While childhood infections are mostly asymptomatic, 50 to 70% of young adults undergoing primary EBV infections show mild to severe illness. EBV may cause a persistent, latent infection which can be reactivated under immunosuppression or in AIDS affected patients. As humoral responses to primary EBV infections are quite rapid, the level and class of antibodies raised in most cases allow classification as to whether the patient is still susceptible, has a current or recent primary infection, had a past infection or may be having reactivated EBV infection. The detection of EBV-specific IgG, IgM and IgA antibodies to its major immunodominant antigens (mainly Nuclear Antigen or EBNA and Viral Capsidic Antigen or VCA) has become therefore an important and useful determination for the monitoring and follow-up of EBV infected patients.

C. PRINCIPLE OF THE TEST

In order to get rid of crossreactions with other viruses of the same family, microplates are coated with affinity purified native VCA antigen, to provide the assay with the highest specificity and sensitivity.

In the 1st incubation, the solid phase is treated with diluted samples and anti-VCA IgG are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti-VCA IgG are detected by the addition of anti hIgG antibody, labeled with peroxidase (HRP). The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti-VCA IgG antibodies present in the sample.

IgG in the sample may therefore be quantitated by means of a standard curve calibrated in arbitrary units per milliliter (arbU/ml) as no international standard is available.

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1. Microplate: MICROPLATE

12 strips x 8 microwells coated with affinity purified native VCA antigen. Plates are sealed into a bag with desiccant.

Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve ranging:

4 ml CAL1 = 0 arbU/ml
4 ml CAL2 = 5 arbU/ml
2 ml CAL3 = 10 arbU/ml
2 ml CAL4 = 20 arbU/ml
2 ml CAL 5 = 50 arbU/ml
4 ml CAL6 = 100 arbU/ml.

Standards are calibrated against an internal Gold Standard or IGS as no international one is defined.

Contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized. It contains bovine serum proteins, human IgG antibodies to VCA at 20 arbU/ml±20%, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

3. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

4. Enzyme conjugate : CONJ

1x16ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

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7. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

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1. Calibrated Micropipettes (1000, 100 and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.

2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 6 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
14. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.
3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

4. Sera and plasma can be stored at +2..8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be freezed/thawed more than once as this may generate particles that could affect the test result.
5. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.
6. Samples whose anti-VCA IgG antibody concentration is expected to be higher than 100 arbU/ml should be diluted before use, either 1:10 or 1:100 in the Calibrator 0 arbU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 ul of each specimen with 450 ul of Cal 0 (1:10). Then 50 ul of the 1:10 dilution are diluted with 450 ul of the Cal 0 (1:100). Mix tubes thoroughly on vortex and then proceed toward the dilution step reported in section M.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of storing. In this case call Dia.Pro's customer service.

Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2..8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Note: Once diluted, the wash solution is stable for 1 week at +2..8°C.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the

liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the Control Serum as reported.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.
12. Check that all the other equipment is available and ready to use.
13. In case of problems, do not proceed further with the test and advise the supervisor.

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M1. QUANTITATIVE DETERMINATION:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of Microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
3. Dispense 100 µl of Calibrators and 100 µl Control Serum in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

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- Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 100 µl of Calibrator 0 arbU/ml and Calibrator 5 arbU/ml in duplicate and Calibrator 100 arbU/ml in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **60 min at +37°C**.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from yellow to blue.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

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N. ASSAY SCHEME

Method	Operations
Calibrators	100 µl
Control Serum (*)	100 µl
Samples diluted 1:101	100 µl
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm

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Legenda: BLK = Blank
S = Sample
CAL = Calibrators

CAL 2 5 arbU/ml OD450nm < OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 100 arbU/ml < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
CAL 1 0 arbU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.

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Problem	Check
Control Serum Different from Expected value	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the control serum has occurred.

Anyway, if all other parameters (Blank, CAL1, CAL2, CAL 6), match the established requirements, the test may be considered valid.

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
Then on the calibration curve calculate the concentration of anti VCA IgG antibody in samples.

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Example of calculation:

Note: The following data must not be used instead of real figures obtained by the user.

Calibrator 0 arbU/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted

Calibrator 5 arbU/ml: 0.250 – 0.270 OD450nm
 Mean Value: 0.260 OD450nm
 Higher than Cal 0 + 0.100 – Accepted

Calibrator 100 arbU/ml: 2.045 OD450nm
 Higher than 1.000 – Accepted

The OD450nm of the Calibrator 5 arbU/ml is considered the cut-off (or Co) of the system.

The ratio between the OD450nm value of the sample and the OD450nm of the Calibrator 5 arbU/ml (or S/Co) can provide a semi-quantitative estimation of the content of specific IgG in the sample.

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 5 arbU/ml are considered negative for anti-VCA IgG antibody.

Samples with a concentration higher than 5 arbU/ml are considered positive for anti-VCA IgG antibody.

VCA IgG results alone are not, anyway, enough to provide a clear diagnosis of EBV infection. At least EBV VCA IgM results, possibly together with EBNA IgG, are necessary in combination. A reference range of the minimum essential serological markers of Epstein-Barr infection, derived from Infectious Diseases Handbook, 3rd edition, published by Lexi-Comp Inc., USA, is reported schematically below:

VCA IgM	EBNA (or VCA) IgG	Interpretation
negative	negative	No history of EBV infection
positive	negative	Acute primary infection
negative	positive	History of previous infection
positive	positive	Reactivation

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R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in an external clinical center on negative and positive samples with reference to a FDA approved commercial kit.

1. Limit of detection

No international standard for VCA IgG Antibody detection has been defined so far by the European Community. In its absence, an Internal Gold Standard (or IGS), derived from a patient with an history of past mononucleosis infection, has been defined in order to provide the device with a constant and excellent sensitivity.

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Microplates are coated with with affinity purified native VCA antigen capable to provide the assay with the highest specificity and sensitivity.

The diagnostic performances were evaluated in a performance evaluation study conducted in an external centre, with excellent experience in the diagnosis of infectious diseases.

The diagnostic sensitivity was studied on more than 50 samples, pre-tested positive with a different reference kit of European origin in use at the laboratory. Positive samples were collected from patients that experienced mononucleosis infection.

The diagnostic specificity was determined on panels of more than 50 negative samples from normal individuals and blood donors, classified negative with the reference kit, including potentially interfering specimens.

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Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.

The Performance Evaluation provided the following values :

Sensitivity	≥ 98 %
Specificity	≥ 98 %

3. Reproducibility:

Data obtained from a study conducted on three samples of different VCA IgG reactivity, examined in 16 replicates in three separate runs show CV% values ranging 3-16% depending on OD450nm readings.

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S. LIMITATIONS

False positivity has been assessed as less than 2-5% of the normal population depending on the reference kit used.

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VCA IgM

**“Capture” Enzyme ImmunoAssay
(ELISA) for the quantitative/qualitative
determination of IgM class antibodies to
Epstein Barr Virus Capsidic Antigen
in human plasma and sera**

- for “in vitro” diagnostic use only -

DIA.PRO

Diagnostic Bioprobes Srl

Via Carducci n° 27

20099 – Sesto San Giovanni

Milano – Italy

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

VCA IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative or qualitative determination of IgM class antibodies to Epstein Barr Virus (EBV) Capsidic Antigen in human plasma and sera with the "capture" system.

The kit is intended for the classification of the viral infective agent and the follow-up of EBV infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Epstein Barr Virus or EBV is the principal etiological agent of infectious mononucleosis, as well as a contributory factor in the etiology of Burkitt's lymphoma and nasopharyngeal carcinoma, or NPC.

A member of the family Herpesviridae, it has a worldwide distribution, such that 80 to 90% of all adults have been infected. Primary infections usually occur during the first decade of life. While childhood infections are mostly asymptomatic, 50 to 70% of young adults undergoing primary EBV infections show mild to severe illness.

EBV may cause a persistent, latent infection which can be reactivated under immunosuppression or in AIDS affected patients. As humoral responses to primary EBV infections are quite rapid, the level and class of antibodies raised in most cases allow classification as to whether the patient is still susceptible, has a current or recent primary infection, had a past infection or may be having reactivated EBV infection.

The detection of EBV-specific IgG, IgM and IgA antibodies to its major immunodominant antigens has become therefore an important and useful determination for the monitoring and follow-up of EBV infected patients.

C. PRINCIPLE OF THE TEST

The assay is based on the "IgM Capture" method and on affinity purified native VCA antigen.

Microplates are coated with a polyclonal anti-IgM antibody that in the 1st incubation "captures" specifically this class of antibodies.

After washing out all the other components of the sample, in the 2nd incubation bound anti EBV-VCA IgM are detected by the addition of a complex formed by biotinylated affinity purified native VCA antigen and Streptavidine, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of IgM antibodies present in the sample and can be detected by an ELISA reader.

Quantification of IgM is made possible by a standard curve calibrated in arbitrary units, in absence of an international standard to refer to.

D. COMPONENTS

Each kit contains sufficient reagents to carry out 96 tests.

1. Microplate: MICROPLATE

12 strips x 8 breakable wells coated with affinity-purified anti human IgM specific (u-chain) goat polyclonal antibody and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve ranging:

4 ml CAL1 = 0 arbU/ml

4 ml CAL2 = 10 arbU/ml

2 ml CAL3 = 20 arbU/ml

2 ml CAL4 = 50 arbU/ml

4 ml CAL5 = 100 arbU/ml.

Standards are calibrated against an internal Gold Standard or IGS as no international one is defined.

Contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized. Contains fetal bovine serum proteins, human anti EBV VCA IgM antibodies at $20 \pm 20\%$ arbU/ml, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label .

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Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

4. Enzyme conjugate: CONJ 20X

1x0.8 ml/vial. 20x concentrated solution. It contains peroxidase (HRP) labeled Streptavidine, dissolved into a buffered solution of 10 mM Tris buffer pH 6.8 +/-0.1, 5% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

5. Antigen Diluent : AG DIL

n° 1 vial of 16 ml. Protein buffer solution for the preparation of the working EBV VC antigen. The solution contains 10 mM Tris buffer pH 6.8 +/-0.1, 2% BSA, 0.1% Kathon GC and 0.2 mg/ml gentamicine sulphate as preservatives. The reagent is code coloured with 0.01% red alimentary dye

6. EBV VCA Antigen : Ag VCA

1x6 vials. Lyophilized reagent to be dissolved with 1.9 ml of Antigen Diluent as reported in the proper section. It contains biotinylated affinity purified native VCA antigen, 25 mM Tris buffer pH 7.8 +/-0.1 and 5% BSA as proteic carrier.

7. Specimen Diluent: DILSPE

2x60.0 ml/vial. Buffered solution for the dilution of samples. It contains 2% casein, 0.2 M Tris buffer pH 6.0 +/-0.1, 0.2% Tween 20, 0.1% Kathon GC and 0.09% sodium azide as preservatives. The component is blue color coded.

8. Chromogen/Substrate: SUBS TMB

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide of H₂O₂.

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1. Calibrated Micropipettes in the range 10-1000 ul and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.

4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and if with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB/H₂O₂) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at +2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic labware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
14. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2..8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

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Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2..8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Control Serum:

Lyophilized reagent to be dissolved with EIA grade water as reported in the label.

Note: In order to maintain its reactivity fully preserved, upon dissolution keep the excess frozen in aliquots at -20°C and use just once. Do not freeze again.

Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200 ml and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

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Antigen-Conjugate Complex:

Proceed carefully as follows:

1. Dissolve the content of a lyophilized vial with 1.9 ml of Antigen Diluent. Let fully dissolved the lyophilized content and then gently mix on vortex.
2. Gently mix the concentrated Enzyme Conjugate on vortex. Then add 0.1 ml of it to the vial of the dissolved EBV VC Ag and mix gently on vortex.

Important Notes:

1. Dissolve and prepare only the number of vials necessary to the test. The complex obtained is not stable. Store any residual solution frozen in aliquots at -20°C.
2. The preparation of the complex has to be done **right before** the dispensation of samples and controls into the plate. Mix again on vortex gently just before its use.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.
2. The ELISA incubator has to be set at +37°C (tolerance of ±0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water

baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.

3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container).
5. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Dissolve the Control Serum as described above and gently mix.

8. Prepare the Antigen/Conjugate complex as reported before.
9. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
10. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
11. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
12. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
13. Check that the micropipettes are set to the required volume.
14. Check that all the other equipment is available and ready to use.
15. In case of problems, do not proceed further with the test and advise the supervisor.

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M.1 Quantitative analysis

1. Place the required number of strips in the microplate holder. Leave A1 and B1 wells empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2..8°C, sealed.
2. Dilute samples 1:101 dispensing 1 ml Specimen Diluent into a disposable tube and then 10 ul sample; mix on vortex before use. Do not dilute the calibrators and the control serum as they are ready-to-use.
3. Prepare the Antigen/Conjugate complex as reported in Section H.
4. Pipette 100 µl of all the Calibrators and 100 µl of Control Serum in duplicate; then dispense 100 µl of samples. The Control Serum is used to verify that the whole analytical system works as expected. Check that Calibrators Control Serum and samples have been correctly added. Then incubate the microplate at **+37°C for 60 min**.

Important note: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

5. Wash the microplate as reported in section I.3.
6. In all the wells, except A1 and B1, pipette 100 µl Antigen/ Conjugate Complex. Check that the reagent has been correctly added. Incubate the microplate at **+37°C for 60 minutes**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Complex. Contamination might occur.

7. Wash the microplate as described in section I.3.
8. Pipette 100 µl TMB/H₂O₂ mixture in each well, the blank wells A1+B1 included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 minutes**.

Important note: Do not expose to strong direct light as a high background might be generated.

9. Stop the enzymatic reaction by pipette 100 µl Sulphuric Acid into each well and using the same pipetting sequence as in step 8. Then measure the color intensity with a microplate reader at 450nm (reading) and at 620-630nm (blinking, strongly recommended), blanking the instrument on A1, or B1 or both wells.

M.2 Qualitative analysis

1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2..8°C, sealed.
2. Dilute samples 1:101 dispensing 1 ml Specimen Diluent into a disposable tube and then 10 ul sample; mix on vortex before use. Do not dilute the calibrators as they are ready-to-use. Then prepare the Antigen/Conjugate complex as reported in Section H.
3. Pipette 100 µl CAL 1 in duplicate, 100 µl CAL 2 in duplicate, 100 µl CAL 5 in single. Then dispense 100 µl of samples. Check that Calibrators and samples have been correctly added. Then incubate the microplate at **+37°C for 60 min**.

Important note: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

4. Wash the microplate as reported in section I.3.
5. In all the wells, except A, pipette 100 µl Antigen/ Conjugate Complex. Check that the reagent has been correctly added. Incubate the microplate at **+37°C for 60 minutes**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Complex. Contamination might occur.

6. Wash the microplate as described in section I.3.
7. Pipette 100 µl TMB/H₂O₂ mixture in each well, the blank well A1 included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 minutes**.

Important note: Do not expose to strong direct light as a high background might be generated.

8. Stop the enzymatic reaction by pipette 100 µl Sulphuric Acid into each well and using the same pipetting sequence as in step 7. Then measure the color intensity with a microplate reader at 450nm (reading) and at 620-630nm (blinking, strongly recommended), blanking the instrument on A1.

Important general notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.

N. ASSAY SCHEME

Calibrators	100 ul
Control Serum (*)	100 ul
Samples diluted 1:101	100 ul
1st incubation	60 min
Temperature	+37°C
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
TMB/H ₂ O ₂ mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm & 620nm

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
CAL 1 OD450nm > 0.200 coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure when the dispensation of calibrators is carried out; 4. that no contamination of the Cal 1 or of the wells where it was dispensed has occurred due to spills of positive samples or Antigen/Conjugate complex; 5. that micropipettes have not become contaminated with positive samples or with the Antigen/Conjugate complex 6. that the washer needles are not blocked or partially obstructed.
CAL 2 OD450nm < Cal 1 + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of a wrong calibrator); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 5 OD450nm < 1.000	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibration has occurred.

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Problem	Check
Control Serum Different from Expected value	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (e.g.: dispensation of a wrong calibrator); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the control has occurred.


Anyway, if all other parameters (Blank, CAL1, CAL2, CAL 5), match the established requirements, the test may be considered valid.

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Then on the calibration curve calculate the concentration of anti EBV VCA IgM antibody in samples.
An example of Calibration curve is reported below.

Note: Do not use these data to calculate the real assay results. The figures above are reported only as an example.

P.2 Qualitative method

Check that the assay is valid.
An example is provided below:

Note: The following data must not be used instead of real figures obtained by the user.

Calibrator 0 arbU/ml: 0.020 – 0.024 OD450nm
Mean Value: 0.022 OD450nm
Lower than 0.200 – Accepted
Calibrator 10 arbU/ml: 0.250 – 0.270 OD450nm
Mean Value: 0.260 OD450nm
Higher than CAL 1 + 0.100 – Accepted
Calibrator 100 arbU/ml: 2.045 OD450nm
Higher than 1.000 – Accepted

The OD450nm of the Calibrator 10 arbU/ml is considered the cut-off (or Co) of the system.

The ratio between the OD450nm value of the sample and the OD450nm of the Calibrator 10 arbU/ml (or S/Co) can provide a semi-quantitative estimation of the content of specific IgM in the sample.

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 10 arbU/ml are considered negative for anti EBV VCA IgM antibody.
Samples with a concentration higher than 10 arbU/ml are considered positive for anti EBV VCA IgM antibody. The patient is likely to be in the acute phase of infection (mononucleosis).

VCA IgM results alone are not, anyway, enough to provide a clear diagnosis of EBV infection. At least EBNA IgG results are necessary in combination.

A reference range of the minimum essential serological markers of Epstein-Barr infection, derived from Infectious Diseases Handbook, 3rd edition, published by Lexi-Comp Inc., USA, is reported schematically below:

VCA IgM	EBNA IgG	Interpretation
negative	negative	No history of EBV infection
positive	negative	Acute primary infection
negative	positive	History of previous infection
positive	positive	Reactivation

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R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in an external clinical center on panels of negative and positive samples with reference to a commercial kit.

1. Limit of detection

No international standard for EBV VCA IgM Antibody detection has been defined so far by the European Community. In its absence, an Internal Gold Standard (or IGS), derived from a patient in the acute phase of mononucleosis infection, has been defined in order to provide the device with a constant and excellent sensitivity.

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The assay is based on the "IgM Capture" method and on affinity purified native VCA antigen in order to provide the highest specificity and sensitivity.

The diagnostic sensitivity was studied on more than 50 samples, pre-tested positive with the reference kit of European origin in use at the laboratory. Positive samples were collected from patients undergoing acute mononucleosis infection.

The diagnostic specificity was determined on panels of more than 250 negative samples from normal individuals and blood donors, classified negative with the reference kit, including potentially interfering specimens.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed. Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.

The Performance Evaluation provided the following values :

Sensitivity	> 98 %
Specificity	> 98 %

3. Reproducibility:

Data obtained from a study conducted on three samples of different VCA IgM reactivity, examined in 16 replicates in three separate runs showed CV% results ranging 2-8%, depending on the OD450nm readings.

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S. LIMITATIONS

False positivity has been assessed as less than 2 % of the normal population, mostly due to high titers of Rheumatoid Factor. IgM capture systems, even if acknowledged to be more specific than sandwich assays, may in fact be influenced by this kind of interfering substance..

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T. CONFIRMATION TEST

In order to provide the medical doctor with the best accuracy in testing for EBV infection, a confirmation assay is reported.

The confirmation test has to be carried out on any positive sample before a diagnosis of primary infection of EBV is released to the doctor.

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1. Prepare the Antigen/Conjugate Complex as described in the proper section.
2. The well A1 of the strip is left empty for blanking.
3. CAL 2 (10 arbU/ml) is dispensed in the strip in positions B1+C1.
4. The positive sample to be confirmed, diluted 1:101, is dispensed in the strip in position D1+E1.
5. The strip is incubated for 60 min at +37°C.
6. After washing, the blank well A1 is left empty.
7. 100 µl of Antigen/Conjugate Complex are dispensed in wells B1+C1+D1.
8. Then 100 µl of Enzyme Conjugate (**CONJ**) alone are added to well E1. **Note:** *This material does not contain any VCA antigen, only the conjugate*
9. The strip is incubated for 60 min at +37°C.
10. After washing, 100 µl Chromogen/Substrate are added to all the wells and the strip is incubated for 20 min at r.t.
11. 100 µl Sulphuric Acid are added to all the wells and then their color intensity is measured at 450nm (reading filter) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

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1. If the sample in position D1 shows an OD450nm lower than the one of CAL 2, a problem of dispensation or contamination in the first test is likely to be occurred. The Assay Procedure in Section M has to be repeated to double check the analysis.
2. If the sample in position D1 shows an OD450nm value higher than the one of CAL 2 and in position E1 shows an OD450nm value still higher than the one of CAL 2, the sample is considered a **false positive**. The reactivity of the sample is in fact not dependent on the specific presence of EBV VCA antigens and a crossreaction with the enzyme conjugate has occurred.
3. If the sample in position D1 shows an OD450nm value higher than the one of CAL 2 and in position E1 shows an OD450nm value lower the one of CAL 2, the sample is considered a **true positive**. The reactivity of the sample is in fact dependent on the specific presence of EBV VCA antigens and not due to any crossreaction with the conjugate alone.

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	VCA IgM CODE: VCAM.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – SEPTEMBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HP Ag CODE: HPAG.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – MARCH 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

HP Ag

**Enzyme Immunoassay for the
qualitative/quantitative determination
of Helicobacter pylori Antigen
in human stools**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF HPAG.CE.96
96 Tests

HP Ag

A. INTENDED USE

Enzyme Immunoassay (ELISA) for the one-step qualitative/quantitative determination of *Helicobacter pylori* Antigen (HP Ag) in human stools. The kit may be used for the follow-up of HP-infected patients and their pharmacological treatment. For "in vitro" diagnostic use only.

B. INTRODUCTION

Helicobacter pylori (Hp) is a Gram negative bacterium, firstly isolated in gastric mucosa by Marshall and Warren in 1983.

This bacterium is widely diffused in men, without limitations of sex and age; it has been found that infections can be transmitted directly by contact with contaminated biological fluids (saliva, stool, body secretions) and also from contaminated food and beverages.

H.pylori, and in particular some pathogenic strains (CagA +), is the etiological agent responsible of most of active infections and lesions of the gastric mucosa in man.

H.pylori infection moreover acts as cofactor in the development of tumoral pathologies of the gastric apparatus and it is suspected to be associated to some inflammatory pathologies of the genital female apparatus, evolving toward neoplastic transformation.

At the present time, the identification of *Helicobacter pylori* is mostly made with invasive histochemical techniques, with the determination of its urease activity on a isotopic substrate (breath test and mass analysis), with time-consuming bacteriological culture systems and with expensive molecular biology techniques (PCR).

ELISA for HP Ag have been only recently introduced as a specific, fast, non invasive (analysis of stools) and cheaper method of detection.

C. PRINCIPLE OF THE TEST

Stools from patients are used as a source of sample for the determination of HP antigen.

Microplates are coated with a cocktail of affinity purified mouse monoclonal antibodies directed to the most specific *Helicobacter pylori* antigens.

In the 1st incubation, the solid phase is treated with the sample, previously extracted from stools, and simultaneously with a mixture of monoclonal antibodies to Hp, conjugated with peroxidase (HRP).

After washing out all the other components of the sample, in the 2nd incubation the bound enzyme specifically present on the solid phase generates an optical signal that is proportional to the amount of *H.pylori* antigens present in the sample.

D. COMPONENTS

Code HPAG.CE.96 contains reagents to perform 96 tests.

1. Microplate MICROPLATE

8x12 microwell strips coated with anti HP Ag specific affinity purified mouse monoclonal antibodies and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Set: CAL ...

n° 1 set of 4 vials - Lyophilized calibrators. To be dissolved with EIA grade water. When dissolved, Calibrators have the following concentrations: **0 - 0.1 - 0.5 - 1.0 ug/ml HP Ag**

They contain fetal bovine serum, inactivated HP Ag, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Important Note: Calibrators when dissolved are not stable. Proceed as described in the proper section for storage.

3. Wash buffer concentrate WASHBUF 20X

1x60ml/bottle - 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2 and 0.05% Tween 20 and 0.045% ProClin 300

4. Enzyme Conjugate CONJ

1x16ml/vial - Ready to use component. It contains Horseradish Peroxidase (HRP) labeled mouse monoclonal antibodies to HP Ag, 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The Enzyme Conjugate is color coded red.

5. Chromogen/Substrate SUBS TMB

1x25ml/vial - It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine (TMB) and 0.02% hydrogen peroxide (H₂O₂).

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6. Specimen Diluent: DILSPE

2x60ml/vial - Buffered solution for the extraction of HP Ag from the specimen and preparation of the sample. It contains 10 mM Tris-HCl buffer pH 7.4+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The component is color coded blue.

7. Sulphuric Acid H₂SO₄ 0.3 M

1x15ml/vial - It contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

8. Plate sealing foils: n° 2

9. Package insert: n° 1

Upon request:

HP Ag Extraction kit n° 1

The kit contains all what is necessary to prepare n° 50 samples extracted from stools collected by patients.

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1. Calibrated variable volume Micropipettes ranging 1000 ul and 200 ul; disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.
9. Disposable plastic micro-spoon stools collection container (available upon request from Dia.Pro s.r.l.)

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for

Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.

4. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.

5. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.

6. Check that the liquid components of the kit are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.

7. Avoid cross-contamination between samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.

8. Avoid cross-contamination between kit components by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.

9. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.

10. Treat all specimens as potentially infective, according to national regulations and laws concerning biological sample handling and wasting.

11. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.

12. Wastes produced during the use of the kit have to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

13. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

14. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water

15. Other waste materials generated from the use of the kit (example: tips used for samples and controls, tools for the extraction of the sample from specimens, used microplates, etc.) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: COLLECTION, PREPARATION AND WARNINGS

1) It is recommended to collect fresh stools in the morning with the plastic collector provided on request together with the kit. Alternatively a conical bottomed disposable tube, provided by the laboratory to the patient, may be used.

2) The patient submitted to the test should not be under antibiotic or anti-bacterial treatments as this pharmaceutical therapy is known to affect H.pylori up to a certain extent, depending on the antibiotic used, giving origin to false interpretation.

3) The patient has to be asked to collect the specimen avoiding any possible contact with urine or water using the plastic spoon present in the stool collector and taking just the amount of specimen necessary to fill up the cavity of the spoon.

4) The patient is asked to deliver the specimen the same day to the laboratory. From the time of collection, the specimen can be stored in the laboratory up to 24 hr at 2..8°C or kept frozen at -20°C for longer time.

5) Specimens, and then samples derived from them, have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is recommended when the number of samples on testing is pretty high.

Important Note: Degradation of HP antigen heavily occurs in stools after 24 hrs generating false negative results, even if the specimen is stored at 2..8°C.

The next following operations are described and represented in figures in the Instructions for Use of the Stool Extraction Kit provided together with the kit.

Operate according to the following instructions:

1) Open the stool collection device and introduce the extraction brush deeply into the specimen. Rotate the brush 2-4 times in order to collect the right amount of biological material (about 0.2 gr).

2) Transfer the brush carefully into the test tube supplied in the kit and then add 1 ml Specimen Diluent. Keeping the brush inside the tube, mix vigorously on vortex for 1 min +/-10% in order to dissolve H.pylori into solution.

3) Discard the brush and insert the filtering piston, supplied with the kit, into the tube. Push gently the piston down into the tube in order to collect not more than 150-200 ul of the liquid phase of the suspension, volume enough to carry out the test.

Important Notes:

a) Be careful not to apply a too strong manual pressure on the piston. The piston could break the tube and spills could be generated. If this should happen, use a paper towel soaked with an hospital disinfectant to clean up the contaminated surfaces.

b) Avoid any addition of preservatives to samples, especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.

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Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C.

Important Note: After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Calibration Set:

Add the volume of ELISA grade water reported in the label to the lyophilized powder of each Calibrator. Let fully dissolve the content and then gently mix on vortex.

Important Note: When dissolved, Calibrators are not stable. Store Calibrators frozen in aliquots at -20°C, carefully labeled with the content of HP Ag present in each of them.

Wash buffer concentrate:

The concentrated solution has to be diluted 20x with ELISA grade water and mixed gently end-over-end before use. During

preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Important Note: *Once diluted, the wash solution is stable for 1 week at +2..8° C.*

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
2. The ELISA incubator has to be set at +37°C (tolerance of $\pm 0.5^\circ\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).

5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.

An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions. The washer has to be correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests.

Important Note: *Due to the nature of the sample used and the possible presence of particles in the sample, be careful to control that the needles of the washer do not get blocked by the presence of stool bodies.*

4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Internal Quality Control". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated workstation, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

Important Note: *Due to the nature of the sample used and the possible presence of particles in the sample, be careful to control that the needles of the workstation do not get blocked by the presence of stool bodies. We strongly suggest to use disposable sample tips in order to avoid any block or damage of fix probes.*

7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.
8. Upon request, Dia.Pro srl offers a sample preparation device able to produce a particle free sample showing excellent performances in the assay. Please inquire.

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1. Prepare the sample from stools as described in section G and represented in the Instructions for Use of the HP Ag Extraction Kit.
2. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
3. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.

4. Dilute all the content of the 20x concentrated Wash Solution as described above.
 5. Dissolve the Calibrator Set as described above.
 6. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
 7. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles reported in the specific-section.
 8. Check that the ELISA reader has been turned on at least 20 minutes before reading.
 9. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
 10. Check that the micropipettes are set to the required volume.
 11. Check that all the other equipment is available and ready to use.
 12. In case of problems, do not proceed further with the test and advise the supervisor.
9. Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

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Two procedures are available: a quantitative method able to provide a quantification of HP Ag in the specimen and a qualitative method.

A. Quantitative Assay

1. Place the required number of strips in the plastic holder and carefully identify the wells for calibrators and samples. Leave A1+B1 wells empty for blanking purposes.
2. Pipette 100 µl Calibrators in duplicate into the calibration wells (see the example of dispensation reported below).
3. With the Pasteur pipette supplied aspirate the liquid filtered up into the inner chamber of the piston and dispense 2 drops (about 100 µl) of sample into the sample well. Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm. (samples show OD values higher than 0.100).
4. Dispense then 100 µl Enzymatic Conjugate in all wells, except for A1+B1, used for blanking operations.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination might occur.

5. Following addition of the conjugate, check that the color of the samples has turned from brown to pale reddish and incubate the microplate for **120 min at +37°C**.

Important notes: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA workstations.

6. When the first incubation is over, wash the microwells as previously described (section I.3)
7. Pipette 200 µl Chromogen/Substrate into all the wells, A1+B1 included. Incubate the microplate protected from light at **room temperature (18-24°C) for 20 min**.

Important note: Do not expose to strong direct light as a high background might be generated.

8. Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction, using the same pipetting sequence as in step 6.

B. Qualitative Assay

1. Place the required number of strips in the plastic holder and carefully identify the wells for calibrators and samples. Leave A1well empty for blanking purposes.
2. Pipette 100 µl Calibrator 1 in duplicate, 100 µl Calibrator 2 in duplicate, 100 µl Calibrator 4 in single and then 100 µl samples. Check for the presence of samples in wells as reported before.
3. Dispense 100 µl Enzymatic Conjugate in all wells, except for A1, used for blanking operations.

Important note: Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination might occur.

4. Following addition of the conjugate, check that the color of the samples has turned from brown to pale reddish and then incubate the microplate for **120 min at +37°C**.

Important notes: Strips have to be sealed with the adhesive sealing foil only when the test is performed manually. Do not cover strips when using ELISA workstations.

5. When the first incubation is over, wash the microwells as previously described (section I.3)
6. Pipette 200 µl Chromogen/Substrate into all the wells, A1 included. Incubate the microplate protected from light at **room temperature (18-24°C) for 20 min**.

Important note: Do not expose to strong direct light as a high background might be generated.

7. Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction, using the same pipetting sequence as in step 6.
8. Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1.

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Legenda: BLK = Blank CAL = Calibrator S = Sample

Important notes:

1. Ensure that no fingerprints or dust are present on the external bottom of the microwell before reading. They could generate false positive results on reading.
2. Reading should ideally be performed immediately after the addition of the acid solution but definitely no longer than 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.

N. ASSAY SCHEME

Operations	Procedure
Calibrators&samples	100 ul
Enzyme Conjugate	100 ul
1st incubation	120 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20'' of soaking OR n° 6 cycles without soaking
Chromogen/Substrate	200ul
2nd incubation	20 min
Temperature	room
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

O. INTERNAL QUALITY CONTROL

A check is performed on the controls/calibrator any time the kit is used in order to verify whether the expected OD450nm/620-630nm or S/Co values have been matched in the analysis. Ensure that:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
CAL 0 ug/ml	< 0.200 mean OD450nm value after blanking
CAL 0.1 ug/ml	OD450nm > OD450nm CAL 0 ug/ml + 0.100
CAL 1 ug/ml	> 1.000 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
CAL 0 ug/ml > 0.200 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive calibrators instead of the negative one); 4. that no contamination of the calibrator or of the wells where the calibrator was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
CAL 0.1 ug/ml OD450nm < Cal 0 ug/ml + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative calibrator instead) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 1 ug/ml < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the calibrator (dispensation of negative calibrator instead). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

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Important note:

The analysis must be done proceeding as the reading step described in the section M, points A8 and B8.

P. CALCULATION OF RESULTS

Quantitative Assay:

Calculate the mean OD450nm/620-630nm value of the calibrators. Then draw a calibration curve possibly using a 4 parameters fitting curve system. Then calculate on the curve the concentration of HP antigen in the sample.

Qualitative Assay:

The test results are calculated by means of a cut-off value determined from the O450nm/ 620-630nm value of the CAL 0 ug/ml (CAL 0) and the OD450nm/620-630nm of the CAL 0.1 ug/ml (CAL 0.1) with the following formula:

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Q. INTERPRETATION OF RESULTS

In the quantitative assay, samples showing a concentration of H.pylori antigen higher than 0.05 ug/ml are considered positive. For the qualitative assay, test results are interpreted as a ratio of the sample OD450nm/620-630nm (S) and the Cut-Off value (Co), mathematically S/Co, according to the following table:

S/Co	Interpretation
< 1.0	Negative
1.0 – 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient is not infected by H.pylori.

Any patient showing an equivocal result should be retested on a second sample.

A positive result is indicative of HP infection and therefore the patient should be treated accordingly.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed first by repeating the test and then, if still positive, by an alternative method before a diagnosis of HP infection is confirmed.
3. When test results are transmitted from the laboratory to another department, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of HP infection has to be taken and released to the patient by a suitably qualified medical doctor. This should be done taking also into account other diagnostic evidences of infection.

An example of qualitative method is reported below. (data obtained proceeding as the the reading step described in the section M, point B8):


Note: The following data must not be used instead of real figures obtained by the user.

Cal 0 ug/ml: 0.040–0.060 OD450nm
 Mean Value: 0.050 OD450nm
 Lower than 0.200 – Accepted

Cal 0.1 ug/ml: 0.210-0.230 OD450nm
 Mean Value: 0.220 OD450nm
 Higher than Cal 0 ug/ml + 0.100 – Accepted
 Cut-Off = (CAL 0 + CAL 0.1) / 2 = 0.135
 Calibrator 1 ug/ml: 2.000 OD450nm
 OD450nm higher than 1.000 – Accepted

Sample 1: 0.028 OD450nm
 Sample 2: 1.690 OD450nm
 Sample 1 S/Co < 1.0 = negative
 Sample 2 S/Co > 1.1 = positive

An Example of Calibration curve is reported below:

R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances was conducted by testing negative and positive samples in an external clinical site and in the own laboratories as well.

1. Limit of detection

The limit of detection of the assay was calculated by examining serial dilution of HP antigen in Sample Diluent.

Results show that the analytical sensitivity of the assay is better than **0.05 ug/ml** when the limit of dilution is considered mean OD450nm/620-630nm CAL 0 ug/ml + 5 SD.

2. Diagnostic performances:

The diagnostic performances of the kit were studied in some external trials against (a) the "breath test", considered by the medical literature the Gold Standard for HP Ag determination and (b) against a commercial CE-marked ELISA, as well.

In the studies were the breath test was used as reference stools were collected from patients the same day of the breath test.

In the study carried out with a commercial ELISA as reference the same samples were extracted with the specific devise supplied in combination with kits.

In the studies conducted in reference to the **breath test** the following mean values were obtained from two centers, one in Italy and one in Spain:

- Sensitivity: 98%
- Specificity: 96%

No confirmation by ELISA was possible on discrepant samples.

In the studies conducted in reference to a CE-marked commercial **ELISA**, similarly based on monoclonal antibodies to H.pylori, the following mean values were obtained:

- Sensitivity: 100%
- Specificity: 93%

The discrepant samples (7% "false positive") turned to be positive in the breath test, revealing a better overall performance of our kit respect to the commercial reference one, when compared on the gold standard method.

Due to the highest specificity shown by the monoclonal antibodies used in the kit no cross-reaction was observed with Campilobacter species.

3. Precision:

The variability shown in the tables did not result in sample misclassification. CV values ranging 4-8%, depending on OD450nm/620-630nm values were observed.

S. LIMITATIONS

False negative results were obtained from samples extracted from specimens stored for more than 1 day at 2..8°C. This finding is explain by the fact that the devise detects "live" H.pylori as confirmed by the results of excellent comparison to the breath test.

False positive results were mostly obtained from samples of liquid stools, difficult to handle and extract.

ASSAY GRAPHICAL SCHEME


Add 100 µl Calibrators, samples and conjugate to the plate and then incubate for 120 min at +37°C


Wash as described in the proper section

Add 200 µl Chromogen/Substrate and incubate 20 min at r.t.

Add 100 µl Sulphuric Acid

Read the plate at 450nm (reading) and at 620-630nm (blanking)

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EC DESIGN-EXAMINATION CERTIFICATE
in accordance with Annex IV, Section 4, Directive 98/79/EC

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**WE HEREBY DECLARE THAT THE ABOVE MENTIONED
PRODUCT MEETS THE PROVISIONS OF THE COUNCIL
DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

NOTIFIED BODY	AEMPS – n° 0318
(EC) CERTIFICATE(S)	<ul style="list-style-type: none">• FULL QUALITY ASSURANCE SYSTEM N° 2003 12 0388 CT (in accordance with Annex IV – except Section IV) of the Directive 98/79/EC), RELEASED BY EC NOTIFIED BODY N° 0318• DESIGN CERTIFICATE N° 2008 02 0539 ED RELEASED BY EC NOTIFIED BODY N° 0318• UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY EC NOTIFIED BODY N° 0318

PLACE & DATE OF FIRST ISSUE	MILANO – FEBRUARY 2008
PLACE & DATE OF CURRENT EMISSION	SESTO SAN GIOVANNI (MI) – MARCH 2018
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 0318

HIV Ab&Ag

**Fourth generation Enzyme Immunoassay
for the determination of antibodies to
Human Immunodeficiency Virus or HIV
type 1&2&O and P24 HIV-1 Antigen
in human serum and plasma**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

REF IVCOMB.CE
96/192/480/960 Tests

HIV Ab&Ag

A. INTENDED USE

The kit is a solid phase enzyme immunoassay for the in-vitro diagnostic screening of antibodies to all subtypes of HIV-1 and HIV-2 and HIV-1 antigen (p24) in human serum or plasma.

This kit is intended exclusively for *In vitro* diagnostic use in an authorized clinical laboratory and the test has to be carried out by specifically trained health-care professional personnel.

B. INTRODUCTION

Epidemiological evidence indicates that an infectious agent transmitted through intimate contact, intravenous drug use or use of infected blood or blood products leads to Acquired Immunodeficiency Syndrome (AIDS).

This disease affects T-cell mediated immunity, resulting in severe lymphopenia and a reduced subpopulation of helper T-lymphocytes. Destruction of this T-lymphocyte population by the virus causes an immune deficiency, resulting in a reduced or deficient response to subsequent infections.

Consequently, infections become more severe and may cause death. At present, there is no successful treatment for AIDS.

The etiological agent has been identified as a retrovirus, human immunodeficiency virus type 1 (HIV-1).

A closely related, but distinct type of immunodeficiency virus, designated HIV-2, has also been isolated. This virus causes a disease that is indistinguishable from AIDS.

Serological cross-reactivity between HIV-1 and HIV-2 has been shown to be highly variable from sample to sample.

This variability requires the inclusion of antigens to both HIV-1 and HIV-2 for the screening of antibodies to HIV-1 and HIV-2.

The presence of anti-HIV-1 and/or anti-HIV-2 and/or HIV p24 antigen in the blood indicates potential infection with HIV-1 and/or HIV-2 and consequently this blood should not be used for transfusion or for manufacture of injectable products.

C. PRINCIPLE OF THE TEST

Synthetic peptides representing immunodominant epitopes of HIV-1 and HIV-2 together with a monoclonal antibody to p24 HIV-1 antigen are coated onto wells of a microplate.

The peptides and the antibody have been carefully selected to ensure the screening of antibody and p24 antigen to all HIV-1 subtypes, including subtype O and HIV-2. Serum or plasma samples are added to these wells and, if antibodies specific to HIV-1 and/or HIV-2 (IgG, IgM or IgA) are present in the sample, they will form stable complexes with the HIV peptide antigens in the well. In case HIV-1 p24 is present in the sample, the antigen will be captured by the specific monoclonal antibody.

Antigen-antibody complexes are then identified through the successive addition of: (1) biotinylated peptides, a biotinylated monoclonal antibody to HIV-1 p24, and; (2) horseradish peroxidase HRP Streptavidin conjugate.

The hydrolytic activity of horseradish peroxidase allows for the quantification of these antibody-antigen complexes.

Peroxidase substrate solution is then added.

During incubation, a blue color will develop in proportion to the amount of anti-HIV-1/2 antibodies or HIV-1 p24 antigen bound to the well, thus establishing their presence or absence in the sample.

Wells containing samples negative for anti-HIV antibody and/or p24 antigen remain colorless.

A stop solution is added to each well and the resulting yellow color is read on a microplate reader at 450 nm.

D. COMPONENTS

The standard format of the product code IVCOMB.CE contains reagents for 192 tests.

1. Microplate **MICROPLATE**

n° 2 microplates. 12 strips of 8 breakable wells coated with HIV specific gp36, gp41 and gp120 peptides and with a Monoclonal Antibody specific to the HIV-1 p24 Ag. Plates are sealed into a bag with desiccant.

2. Negative Control **CONTROL -**

1x4.0ml/vial. Ready to use control. It contains animal serum negative for HIV antibodies and for p24 antigen, and 0.1% Kathon GC as preservatives. The negative control is yellow-brown color coded.

3. Positive Control HIV-1 Ab **CONTROL 1+**

1x4.0ml/vial. Ready to use control. It contains inactivated HIV 1 antibody positive serum, filtered HIV Ab&Ag negative animal serum and 0.1% Kathon GC as preservative. The Positive Control is light green color coded.

Important Note: *The positive control has been inactivated using B-propionolactone BPL/UV. This does not fully ensure the absence of viable pathogens, and therefore, the control should be handled as potentially biohazardous, in accordance with good laboratory practices.*

4. Positive Control HIV-2 Ab **CONTROL 2+**

1x4.0ml/vial. Ready to use control. It contains inactivated HIV2 antibody positive serum, filtered HIV Ab&Ag negative animal serum and 0.1% Kathon GC as preservatives. The Positive Control is dark green color coded.

Important Note: *The positive control has been inactivated using B-propionolactone BPL/UV. This does not fully ensure the absence of viable pathogens, and therefore, the control should be handled as potentially biohazardous, in accordance with good laboratory practices.*

5. HIV-1 P24 Ag Calibrator **CAL Ag**

2 vials. Lyophilized. It contains not infectious recombinant p24 antigen in a 10 mM phosphate buffer pH 7.0+/-0.2 with 0.3 mg/ml Gentamicine Sulphate and 0.1% Kathon GC as stabilizers. This component is calibrated against the NIBSC 1st International HIV-1 p24 Ag reference sample 90/636 (diluted 1:256) as well as the EFS HIV Ag performance panel (3015-3022).

Important Notes:

- 1) *The Calibrator contains p24 recombinant Ag with a concentration of about 100 pg/ml, corresponding to about 4 IU/ml.*
- 2) *The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.*

6. Wash buffer concentrate **WASHBUF 20X**

2x60ml/bottle. 20x concentrated solution. It contains 0.1% Kathon GC. Once diluted, the wash solution contains 10 mM phosphate buffer saline pH 7.0+/-0.2 and 0.05% Tween 20.

7. Conjugate # 1 **CONJ 1**

8 vials. The vial contains lyophilized biotinylated HIV1&2&0 gp36, gp41 and gp120 peptides and a biotinylated monoclonal antibody specific for HIV 1 p24 antigen. Vials are to be resuspended with 6 ml of the Conjugate # 1 diluent.

8. Conjugate 1 Diluent **CONJ 1 DIL**

1x60ml/bottle. Used to dissolve the lyophilized powder of Conjugate # 1, it contains Tris saline Buffer supplemented with 0.05% Kathon GC, Tween 20 and BSA.

9. Conjugate # 2 **CONJ 2**

1x25ml/bottle The solution contains HRP conjugated with streptavidin in Tris saline Buffer supplemented with 0.05% Kathon GC, Tween 20 and BSA. This component is color coded in pink/red.

10. Chromogen/Substrate SUBS TMB

1x45ml/bottle Ready-to-use component. It contains 50 mM citrate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

Note: To be stored protected from light as sensitive to strong illumination.

11. Sulphuric Acid H₂SO₄ 0.3 M

1x25ml/bottle It contains 0.3 M H₂SO₄ solution.
Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363)

12. Sample Diluent: DILSPE

1x14ml/vial Contains Tris saline buffer supplemented with 0.05% Kathon GC, anti HAMA blocker, and Tween 20; used for specimen dilution. This component is color coded in light blue.

13. Plate sealing foils n° 4

14. Package insert n° 1

Important note: Upon request, Dia.Pro can supply reagents for 96, 480, 960 tests, as reported below :

1.Microplate	n°1	n°5	n°10
2.Negative Control	1x2.0ml/vial	1x10ml/vial	1x20ml/vial
3.Positive Control 1	1x2.0ml/vial	1x10ml/vial	1x20ml/vial
4.Positive Control 2	1x2.0ml/vial	1x10ml/vial	1x20ml/vial
5.HIV p24 Calibrator	n°1 vial	n° 5 vials	n° 10 vials
6.Wash buff conc	1x60ml/bottle	5x60ml/bottles	4x150ml/bottles
7.Conjugate # 1	n°4 vials	n°20 vials	n°40 vials
8.Conjugate 1 Diluent	1x30ml/vial	3x50ml/bottles	2x150ml/bottles
9.Conjugate # 2	1x15ml/vial	2x38ml/bottle	2x75ml/bottle
10.Chrom/Substrate	1x25ml/vial	3x42ml/bottle	2x125ml/bottles
11.Sulphuric Acid	1x15ml/vial	2x40ml/bottles	2x80ml/bottles
12.SampleDiluent	1x7ml/vial	1x35ml/bottle	1x70ml/bottle
Plate seal foils	n° 2	n° 10	n° 20
Pack. insert	n° 1	n°1	n° 1
Number of tests	96	480	960
Code	IVCOMB.CE.96	IVCOMB.CE.480	IVCOMB.CE.960

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (200ul and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.

6. Upon receipt, store the kit at 2.8°C into a temperature controlled refrigerator or cold room.

7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.

8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.

9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.

10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

11. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.

12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.

14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are

treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

15. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

16. The Sulphuric Acid is irritant. In case of spills, wash the surface with plenty of water.

17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMENDATIONS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.

3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.

4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results.

Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

5. Sera and plasma can be stored at +2°.8°C for up to seven days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

6. If particles are present filter using 0.2-0.8µ filters to clean up the sample for testing.

7. Do not use heat inactivated samples as they could give origin to false reactivity.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not shown any relevant loss of activity up to 2 months.

Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the pouch is not broken or that some defect is present indicating a problem of storage. In this case call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, in presence of desiccant supplied, firmly zipped and stored at +2°.8°C. When opened the first time, residual strips are stable up to two months.

Negative Control:

Ready to use. Mix well on vortex before use.

Positive Controls Ab:

Positive controls are ready to use. Handle Positive Controls Ab as potentially infective, even if HIV, if present in the control, has been chemically inactivated.

Calibrator Ag

The Lyophilized Calibrator Ag contains a non-infectious recombinant p24 antigen. The volume of EIA grade water to be used for its dissolution and to reach the appropriate p24 concentration (about 100 pg/ml) is written on the vial label. To help dissolve the lyophilized pellet, vortex a few times, at regular intervals. Complete dissolution should be achieved within 2-5 minutes.

Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.

Wash buffer concentrate:

The 20x concentrated solution has to be diluted with EIA grade water up to 1200 ml and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution.

In the preparation avoid foaming as the presence of bubbles could give origin to a bad washing efficiency.

Important Note: Once diluted, the wash solution is stable for 1 week at +2°.8°C.

Conjugate # 1:

The Conjugate # 1 mix solution must be prepared immediately before starting the test. Add 6 ml of Conjugate 1 diluent directly to one Conjugate # 1 vial to dissolve the lyophilized powder. This preparation (a total of 6 ml in one vial) is sufficient for 32 tests, or 4 complete vertical strips of the microplate. To help dissolve the lyophilized powder, vortex a few times, at regular intervals.

Important Note: Any unused portion of this reconstituted Conjugate # 1 Solution may be stored at +2°.8°C for no more than 12 hours. After this time it has to be discarded and the empty, used container has to be washed with EIA grade water and kept dry for any following re-use.

Conjugate # 2:

Ready to use reagent. Mix well end-over-end before use.

Chromogen/Substrate:

Ready to use. Mix well end-over-end before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container.

Sulphuric Acid:

Ready to use. Mix well end-over-end before use.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P364).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 - Take off contaminated clothing and wash it before reuse.

Sample Diluent:

Ready to use. Mix well end-over-end before use.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. **Micropipettes** have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The **ELISA incubator** has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the sections "Validation of Test" and "Assay Performances". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The **ELISA microplate reader** has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; (d) repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure

that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

6. When using an **ELISA automated work station**, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended for blood screening when the number of samples to be tested exceed 20-30 units per run.
7. When using automatic devices, in case the vial holder of the instrument does not fit with the vials supplied in the kit, transfer the solution into appropriate containers and label them with the same label peeled out from the original vial. This operation is important in order to avoid mismatching contents of vials, when transferring them. When the test is over, return the secondary labeled containers to 2..8°C, firmly capped.
8. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator Ag.
5. Dissolve the Conjugate # 1 vial containing lyophilized powder with the Conjugate 1 Diluent (1 lyophilized Conjugate # 1 + 6ml Conjugate # 1 Diluent) to obtain the Conjugate # 1 Mix as described in the proper section.
6. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
7. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
8. Check that the ELISA reader has been turned on at least 20 minutes before reading.
9. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
10. Check that the micropipettes are set to the required volume.
11. Check that all the other equipment is available and ready to use.
12. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

1. Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument dispense 50 ul Sample Diluent first and then 150 ul controls and samples.

Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples or tips have to be changed.

For the next operations follow the operative instructions reported below for the Manual Assay.

It is strongly recommended to check that the time lap between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

The correct number of lyophilized Conjugate # 1 must be dissolved each with 6 ml Conjugate # 1 Diluent. Once the lyophilized powders are dissolved and mixed well, they are to be mixed together into a plastic container and the assay may begin.

2. Manual assay:

1. Dissolve the right number of lyophilized Conjugate # 1 with Conjugate # 1 Diluent before starting to dispense the samples and controls of the test.
2. Place the required number of strips in the microwell holder. Leave the 1st well empty for the operation of blanking.
3. Dispense 50 ul Sample Diluent in all the wells, except A1 used for blanking.
4. Dispense 150 ul of Negative Control in triplicate, 150 ul HIV1 Positive Control, 150 ul HIV2 Positive Control and 150 ul of Calibrator Ag in duplicate in proper wells.
5. Dispense 150 ul of Sample in each properly identified well. Mix gently the plate on the work surface, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into the diluent.
6. Incubate the microplate for **60 min at +37°C**.

Important note: *Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.*

7. Wash the microplate with an automatic washer by delivering and aspirating 350ul/well of diluted washing solution as reported previously (section I.3).
8. Pipette 150 ul Conjugate # 1 mix, prepared as described before, into each well, except the 1st blanking well, and cover with the sealer.

Important note: *Be careful not to touch the plastic inner surface of the well with the tip filled with the Conjugate. Contamination might occur.*

9. Incubate the microplate for **30 min at +37°C**.
10. Pipette 100 ul of Conjugate # 2 in all the wells, except A1, and gently agitate the microplate to mix the two conjugates.

Important Note: *This solution must be added to the bottom of each well to ensure proper performance. Inadequate mixing of the two solutions (Conjugate 1 and Conjugate 2) may reduce the binding of streptavidin HRP (Conjugate 2) to the biotinylated reagents and consequently affect the performance of the assay. Be sure to provide an adequate mixing when the Conjugate # 2 is added, both in the manual and in the automated procedures.*

11. Incubate the microplate sealed for **30 min at +37°C**.
12. Wash as in section 7.
13. Dispense 200 ul of Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-25°C) for 30 minutes**. Start the timing immediately after addition of this component to the first well.

Important note: Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 ul Sulphuric Acid into all the wells using the same pipetting sequence as in step 13 to stop the enzymatic reaction. Addition of acid will turn the positive controls and positive samples from blue to yellow/brown.
- Measure the color intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

Important notes:

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 30 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.
- The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

Method	Operations
Sample Diluent	50 ul
Controls and calibrator	150 ul
Samples	150 ul
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Conjugate # 1	150 ul
2nd incubation	30 min
Temperature	+37°C
Conjugate # 2	100 ul
3rd incubation	30 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂	200 ul
4th incubation	30 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm

An example of dispensation scheme is reported below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL Ag											
B	NC	CAL Ag											
C	NC	S1											
D	NC	S2											
E	POS 1 Ab	S3											
F	POS 1 Ab	S4											
G	POS 2 Ab	S5											
H	POS 2 Ab	S6											

Legenda: BLK = Blank NC = Negative Control POS 1 Ab = HIV -1 Ab Positive Control, POS 2 Ab = HIV -2 Ab Positive, CAL Ag = HIV p24 Ag Calibrator, S = Sample

O. INTERNAL QUALITY CONTROL

A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450nm values are as expected and reported in the table below.

Check	Requirements
Blank well	≤ 0.100 OD450nm value
Negative Control (NC)	≤ 0.200 mean OD450nm value after blanking Absorbance of individual negative control values must be less than or equal to 0.200. If one value is outside this range, discard this value and recalculate mean. If two values are outside this range the run should be repeated.
HIV-1 Ab Positive Control	Mean OD450nm ≥ 0.700.
HIV-2 Ab Positive Control	Mean OD450nm ≥ 0.700.
HIV Ag Calibrator	S/Co > 1

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
Negative Control (NC) > 0.200 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of their wells has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Controls < 0.700 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.200, too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.
HIV Ag Calibrator S/Co < 1	1. that the procedure has been correctly executed; 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of Calibrator Ag. In this case, the negative control will have an OD450nm value > 0.200, too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred. 5.that the lyophilize powder was dissolved correctly with the correct volume of water written on the vial label.

Should these problems happen, after checking, report any residual problem to the supervisor for further actions.

P. CALCULATION OF THE CUT-OFF

The tests results are calculated by means of a cut-off value determined with the following formula on the mean OD450nm value of the Negative Control (NC):

$$NC + 0.125 = \text{Cut-Off (Co)}$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Important note: When the calculation of results is done by the operative system of an ELISA automated work station be sure that the proper formulation is used to calculate the cut-off value and generate the right interpretations of results.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 1	Negative
> 1	Positive

A negative result indicates that the patient has not been infected by HIV.

If the initial absorbance value is equal to or greater than the cut-off value, retest the sample in duplicate. If both retest values are less than the cut-off, the interpretation is not reactive for HIV antibody and/or antigen (negative).

If one or both retest values are equal to or greater than the cut-off the interpretation of the test results is repeatedly reactive. The sample should be considered reactive or positive for HIV antibody and/or antigen according to the criteria of this HIV ELISA test.

A positive result is indicative of HIV infection and therefore the patient should be treated accordingly.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
2. Repeatedly reactive specimens should be submitted to a Confirmation Assay before diagnosis of HIV infection is released.
3. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
4. Diagnosis of HIV infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead or real figures obtained by the user.

Negative Control: 0.110 – 0.120 - 0.115
 OD450nm Mean Value: 0.115 OD450nm
 Lower than 0.200 – Accepted

HIV 1 Ab Positive Control: 2.000 OD450nm mean value
 Higher than 0.700 – Accepted

HIV 2 Ab Positive Control: 2.100 OD450nm mean value
 Higher than 0.700 – Accepted

Calibrator Ag: 0.322 OD450nm mean value
 S/Co > 1 - Accepted

Cut-Off = 0.115 + 0.125 = 0.240

Sample 1: 0.070 OD450nm

Sample 2: 1.690 OD450nm

Sample 1 S/Co < 1 = negative

Sample 2 S/Co > 1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

The performance evaluation was carried out both in an external centre of excellence for HIV diagnosis, that examined the device on a population of antibody positive and negative samples against a CE-marked kit used in the laboratory as reference, and in DiaPro's laboratories as well to complete the study.

R.1 ANALYTICAL SENSITIVITY

The limit of detection (or analytical sensitivity) of the assay has been calculated by means of preparations specific for HIV-1 and HIV-2 antibody and HIV-1 p24 Ag detection, supplied by NIBSC Blanche Lane South Mimms Potters Bar Hertfordshire EN6 3QG, UK.

Samples were diluted in HIV Ab&Ag negative plasma to generate limiting dilution curves and examined in duplicate.

The tables below reports the mean OD450nm values and the S/Co index:

**NIBSC anti-HIV 2 monitor sample
code 99/674 – 005**

Sample	Lot #	0506	Lot #	0706	Lot #	0906
Dilution	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1x	over	>14.5	over	>15.0	over	>15.1
2x	3.838	14.1	3.765	14.5	3.774	14.4
4x	2.371	8.7	2.268	8.7	2.319	8.9
8x	1.253	4.6	1.097	4.2	1.140	4.4
16x	0.700	2.6	0.712	2.7	0.735	2.8
32x	0.462	1.7	0.439	1.7	0.483	1.8
64x	0.281	1.0	0.260	1.0	0.294	1.1
128x	0.189	0.7	0.171	0.7	0.174	0.7
diluent	0.140	0.5	0.122	0.5	0.131	0.5

The device shows a limiting dilution value at 64x.

**NIBSC British working standard for anti HIV 1
code 99/750 –007**

Sample	Lot #	0506	Lot #	0706	Lot #	0906
Dilution	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1x	2.206	8.1	2.086	8.0	2.142	8.2
2x	0.999	3.7	0.925	3.6	1.027	3.9
4x	0.475	1.7	0.486	1.9	0.486	1.9
8x	0.295	1.1	0.301	1.2	0.289	1.1
16x	0.212	0.8	0.206	0.8	0.220	0.8
diluent	0.140	0.5	0.122	0.5	0.131	0.5

The device shows a limiting dilution value at 8x.

**NIBSC British working standard for anti HIV 1
code 99/710 – 007**

Sample	Lot #	0506	Lot #	0706	Lot #	0906
Dilution	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1x	0.588	2.2	0.611	2.4	0.607	2.3
2x	0.301	1.1	0.309	1.2	0.312	1.2
4x	0.198	0.7	0.210	0.8	0.192	0.7
diluent	0.140	0.5	0.122	0.5	0.131	0.5

The device shows a limiting dilution value at about 2x.

**NIBSC HIV-1 p24 Antigen Monitor Sample
code 02/146-002**

Sample	Lot #	0506	Lot #	0706	Lot #	0906
Dilution	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1x	3.664	13.4	3.650	14.0	3.552	13.5
2x	2.151	7.9	2.133	8.2	2.086	8.0
4x	1.209	4.4	1.178	4.5	1.214	4.6
8x	0.734	2.7	0.729	2.8	0.780	3.0
16x	0.388	1.4	0.342	1.3	0.351	1.3
32x	0.259	0.9	0.236	0.9	0.229	0.9
diluent	0.140	0.5	0.122	0.5	0.131	0.5

The device shows a limiting dilution value at about 16x.

**NIBSC 1st International reference Reagent for HIV 1 Ag
code 90/636 – (Version 4, 12 May 2009)**

Sample	Lot # 1		Lot # 2		Lot # 3	
	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
16	2.053	14.8	2.048	14.8	2.053	14.5
8	1.118	8.0	1.121	8.1	1.124	8.0
4	0.571	4.1	0.574	4.1	0.576	4.1
2	0.290	2.1	0.291	2.1	0.290	2.0
1	0.160	1.2	0.162	1.2	0.160	1.1
0.5	0.104	0.7	0.105	0.8	0.103	0.7
diluent	0.014	//	0.014	//	0.014	//

Note: Lot 1 = C2E2T2/6 – Lot 2 = C2E2T2/5 – Lot 3 = C2E2T2/4

The devise shows a sensitivity ≤ 2 IU/ml as required by CTS:2009.

**BBI Anti-HIV 1 Low Titer
Performance Panel - PRB 107
(modified version)**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
2	5.6	3.7
3	2.1	4.5
5	0.5	0.1
6	4.4	7.3
7	1.4	1.0
8	6.3	7.1
10	1.8	3.5
12	5.9	5.0
13	7.1	3.3
15	3.2	2.7

R.2 DIAGNOSTIC SPECIFICITY AND SENSITIVITY

R.2.1 Diagnostic Specificity:

In addition to the first study, where a total of more than 5000 unselected blood donors, more than 200 hospitalized patients, (under examination for non HIV pathologies) and more than 100 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.) were tested, the diagnostic specificity was recently assessed by testing a total of 3268 negative samples on four different lots. A value of diagnostic specificity of 100% was observed.

No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity.

Frozen specimens have been tested, as well, to check for interferences due to collection and storage.

No interference was observed.

R.2.2 Diagnostic Sensitivity

The diagnostic sensitivity of the test was determined on a population of HIV positive specimens.

Results are reported in the tables below.

**BBI anti HIV-1 Low Titer
Performance Panel - PRB 106**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	4.5	7.5
2	2.6	1.8
3	5.2	2.0
4	6.1	8.9
5	11.4	17.6
6	0.5	0.1
7	3.9	7.9
8	13.2	15.7
9	>14.5	16.6
10	6.1	12.7
11	4.8	6.8
12	3.0	2.4
13	5.9	10.5
14	7.7	6.0
15	>14.5	5.5

**Etablissement Francais du Sang
Mixed titer Performance
Panel HIV Ab (1-6) Lot # 11**

ID	Composition	Lot # 0506	Lot # 0706	Lot # 0906	S/Co
		S/Co	S/Co	S/Co	mean
1	HIV2(1/200)	>14.5	>15.0	>15.1	>14.9
2	HIV2(1/800)	10.0	10.2	10.2	10.1
3	negative	0.4	0.4	0.4	0.4
4	HIV1(1/700)	8.9	9.0	9.1	9.0
5	HIV1(1/160)	>14.5	>15.0	>15.1	>14.9
6	HIV1(1/200)	1.9	1.9	1.8	1.9

**BBI Anti-HIV 1 Low Titer
Performance Panel - PRB 108**

Member ID #	IVCOMB.CE S/Co	IVAB.CE S/Co
1	3.5	3.1
2	0.5	0.2
3	2.3	2.5
4	11.4	10.6
5	3.8	2.9
6	6.9	5.3
7	4.0	3.7
8	7.5	6.3
9	6.3	3.2
10	2.4	1.5
11	8.8	7.2
12	4.1	2.9
13	3.7	2.3
14	7.2	6.1
15	10.9	6.3

**BBI Anti-HIV 1 Mixed Titer
Performance Panel - PRB 204**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	1.2	1.0
2	>14.5	65.2
3	0.3	0.2
4	>14.5	65.1
5	>14.5	66.1
6	>14.5	67.5
7	>14.5	66.0
8	>14.5	64.6
9	1.3	2.4
10	>14.5	65.1
11	>14.5	65.4
12	>14.5	63.8
13	4.7	67.0
14	>14.5	63.8
15	>14.5	63.6
16	>14.5	66.6
17	>14.5	68.8
18	>14.5	71.2
19	>14.5	71.5
20	>14.5	68.5
21	>14.5	71.7
22	>14.5	68.5
23	0.3	0.2
24	2.3	2.3
25	0.8	0.9

**BBI anti HIV 1/2 Combo
Performance Panel - PRZ 206**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	10.2	>15.7
2	0.4	0.2
3	>14.5	>15.7
4	>14.5	>15.7
5	6.6	9.2
6	>14.5	9.1
7	0.4	0.2
8	>14.5	>15.7
9	11.4	15.5
10	>14.5	>15.7
11	>14.5	>15.7
12	>14.5	12.8
13	>14.5	15.7

**BBI HIV p24 Antigen Mixed Titer
Performance Panel - PRA 203**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	2.3	4.3
2	4.1	12.2
3	1.4	6.0
4	1.9	3.0
5	Nd	Nd
6	1.0	1.3
7	2.9	5.7
8	3.0	1.5
9	0.4	0.5
10	3.6	12.2
11	2.4	3.1
12	1.6	3.5
13	1.2	2.1
14	0.4	0.9
15	3.3	2.5
16	1.1	1.8
17	2.5	6.3
18	1.7	2.4
19	0.4	0.5
20	4.6	17.0

**BBI HIV-1 Incidence/Prevalence
Performance Panel - PRB 601**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	>14.5	18.5
2	>14.5	17.4
3	>14.5	17.4
4	>14.5	17.4
5	>14.5	15.8
6	>14.5	17.4
7	>14.5	18.5
8	>14.5	17.4
9	>14.5	17.4
10	>14.5	17.4
11	>14.5	18.5
12	>14.5	17.4
13	>14.5	17.4
14	>14.5	17.4
15	>14.5	17.4

Moreover, in the external Performance Evaluation a total of 651 positive samples, including HIV type 1, HIV type 2, HIV type 1 mixed subtypes (including O), HIV 1 Antigen, more than 40 early seroconversion HIV samples and cell culture supernatants were evaluated and a value of 100% was found.

Finally, more than 30 panels of seroconversion containing samples of HIV 1/2/0 Antibodies and/or HIV-1 p24 Antigen positive, obtained from BBI, USA, were evaluated using IVCOMB.CE lot # 0506. In the table below results are reported.

**BBI Worldwide HIV
Performance Panel - WWRB 303**

Member ID #	IVCOMB.CE S/Co	Ref. Kit S/Co
1	>14.5	>15.9
2	>14.5	>15.9
3	>14.5	>15.9
4	>14.5	>15.9
5	4.1	>15.9
6	>14.5	>15.9
7	>14.5	>15.9
8	>14.5	>15.9
9	1.1	2.4
10	>14.5	>15.9
11	11.3	0.3
12	0.3	0.4
13	>14.5	>15.9
14	>14.5	>15.9
15	>14.5	>15.9

Seroconversion Panel	IVCOMB.CE 4 th Generation HIV Ab&Ag	IVAB.CE 3 rd Generation HIV Ab
ID	First specimen detected positive in the panel	
PRB 910 (J)	2	3
PRB 916 (P)	4	5
PRB 917 (Q)	2	5
PRB 919 (S)	1	2
PRB 922 (V)	1	1
PRB 924 (X)	5	6
PRB 926 (Z)	3	5
PRB 927 (AB)	2	3
PRB 928 (AC)	2	3
PRB 929 (AD)	4	6
PRB 930 (AE)	2	3
PRB 931 (AF)	6	6
PRB 933 (AH)	2	2
PRB 935 (AJ)	6	7
PRB 937 (AL)	4	not detected
PRB 938 (AM)	1	3
PRB 939 (AN)	7	9
PRB 940 (AP)	2	4
PRB 941 (AQ)	4	5
PRB 942 (AR)	4	not detected
PRB 944 (AT)	3	5
PRB 946 (AV)	3	not detected
PRB 947 (AW)	2	4
PRB 948 (AX)	4	not detected
PRB 949 (AY)	4	5
PRB 950 (AZ)	3	4
PRB 952 (BB)	3	5
PRB 953 (BC)	3	4
PRB 955 (BE)	3	4
PRB 956 (BF)	5	not detected

**Etablissement Francais du Sang
Performance Panel
HIV Ag (3015-3022) lot 2004**

Sample	IVCOMB.CE Lot # 0506 S/Co	Concentration HIV 1 p24 Ag [pg/ml]
3015	8.7	500
3016	4.2	250
3017	1.8	100
3018	1.2	50
3019	0.9	25
3020	0.6	10
3021	0.6	5
3022(diluent)	0.5	diluent

The device shows a better sensitivity than the previous generation as it is able to detect the p24 antigen.

The results of the Performance Evaluation, correlate perfectly with what stated by EU CTS and show an overall value of diagnostic sensitivity of 100%

R.3 PRECISION

The precision of the device was assessed by determining its values in a within and between runs. In the table below results are reported for a negative sample and a low positive sample.

Average values N = 48	Negative Sample	Low Positive
OD450nm	0.136	0.916
Std.Deviation	0.011	0.022
CV %	7.6	4.0

S. LIMITATIONS

The user of this kit is advised to carefully read and understand this package insert. Strict adherence to the protocol is necessary in order to obtain reliable test results. In particular, accurate sample and reagent pipetting, along with careful washing and timing of incubation steps is essential for accurate and reproducible detection of HIV-1 and HIV-2 antibodies and HIV-1 p24 antigen.

After the EIA test is performed, repeatedly reactive specimens should be submitted for additional testing using Western Blot (WB), Immunofluorescence Assay (IFA), Radioimmunoprecipitation Assay (RIPA) tests and PCR for HIV nucleic acid.

The determination that a person's serum contains antibodies or p24 antigen to HIV has extensive medical, social, psychological and economic implications.

It is recommended that confidentiality, appropriate counselling and medical evaluation be considered an essential aspect of the testing sequence. AIDS and AIDS-related conditions are clinical diseases and their diagnosis can only be established clinically.

EIA testing alone cannot be used to diagnose AIDS.

A non-reactive test result at any point in the testing sequence does not preclude the possibility of exposure to or infection with HIV. The risk of an asymptomatic person, who is repeatedly reactive, developing AIDS and/or AIDS-related conditions is not known.

Falsely reactive test results can be observed with a test kit of this nature. The proportion of reactive samples will depend on the sensitivity and specificity of the test kit and on the prevalence of HIV-1 and HIV-2 antibodies in the population to be screened.

Antibodies to HIV may occur due to voluntary participation in an HIV vaccine study.

Interpretation of this diagnostic test will depend on the type of vaccine given. Correlation with the medical history and additional testing may be necessary to accurately diagnose HIV in vaccine volunteers.

BIBLIOGRAPHY

- Alizon, M., Sonigo, P., Barré-Sinoussi, F., Chermann, J.-C., Tiollais, P., Montagnier, L. and Wain-Hobson, S., 1984. Molecular Cloning of Lymphadenopathy-Associated Virus. *Nature* 312: 757-760.
- Hahn, B.N., Shaw, G.M., Arya, S.K., Popovic, M., Gallo, R.C. and Wong-Staal, F., 1984. Molecular Cloning and Characterization of the HTLV-III Virus Associated with AIDS. *Nature* 312: 166-169.
- Luciw, P.A., Potter, S.J., Steimer, K., Dina, D. and Levy, J.A., 1984. Molecular Cloning of AIDS-Associated Retrovirus. *Nature* 312: 760-763.
- Popovic, M., Sarngadharan, M.G., Read, E. and Gallo, R.C., 1984. Detection, Isolation, and Continuous Production of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and Pre-AIDS. *Science* 224: 497-500.
- Sarngadharan, M.G., Popovic, M., Bruch, L., Schüpbach, J. and Gallo, R.C., 1984. Antibodies Reactive with Human T-Lymphotropic Retroviruses (HTLV-III) in the Serum of Patients with AIDS. *Science* 224: 506-508.
- Vézinet-Brun, F., Barré-Sinoussi, F., Saimot, A.G., Christol, D., Montagnier, L., Rouzioux, C., Klatzmann, D., Rozenbaum, W., Gluckmann, J.C. and Chermann, J.-C., 1984. Detection of IgG Antibodies to Lymphadenopathy-Associated Virus in Patients with AIDS or Lymphadenopathy Syndrome. *Lancet*: 1253-1256, June 9.
- Spire, B., Montagnier, L., Barré-Sinoussi, F. and Chermann, J.-C., 1984. Inactivation of Lymphadenopathy Associated Virus by Chemical Disinfectants. *Lancet*: 899-901, Oct. 20.
- Sehulster, L.M., Hollinger, F.B., Dreesman, G.R. and Melnick, J.L., 1981. Immunological and Biophysical Alteration of Hepatitis B Virus

Antigens by Sodium Hypochlorite Disinfection. *App. and Envir. Microbiol.* 42: 762-767.

- Barre-Sinoussi, F., Chermann, J.C., Rey, F., Nugeyre, M.T., Chamaret, S., Gruest, J., Dauguet, C., Axler-Blin, C., Vézinet-Brun, F., Rouzioux, C., Rozenbrum, W. et Montagnier, L. 1983. Isolation of a T-lymphotropic retrovirus from a patient at risk for acquired immune deficiency syndrome (AIDS). *Science* 220: 868-871.
- Gallo, R.C., Salahuddin, S.Z., Popovic, M., Shearer, G.M., Kaplan, M., Haynes, B.F., Palker, T.J., Redfield, R., Oleske, J., Safai, B., White, G., Foster, P. and Markham, P.D. 1984. Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 224:500-503.
- Gold, J., Dwyer, J. 1994. A short history of AIDS. *Med. J. Aust.* 160:251-252.
- Saville R. D., Constantine N. T., Cleghorn F. R. and Al. Fourth-Generation Enzyme-Linked Immunosorbent Assay for the simultaneous Detection of Human Immunodeficiency Virus Antigen and Antibody. *J. of Clin Microbiology*, July 2001, p.2518-2524.
- Bernard Weber, El Hadji Mbargane Fall, Annemarie Berger, Hans Wilhelm Doerr. Reduction of diagnostic Window by New Fourth-Generation Human Immunodeficiency Virus Screening Assays. *J. of Clin. Microb.* Aug. 1998, p. 2235-2239.
- Clark J., Coates T. J., Lescano A. G., et Al. Different Positive Predictive Values of commercially available Human Immunodeficiency Virus Enzyme-Linked Immunosorbent Assays. *Clin. And Vacc. Immunology*. Feb 2006 p.302-303.
- Novack L., Galai N., Yaari A., Orgel M., Shinar E., Sarov B. Use of Seroconversion Panels to estimate Delay in the detection of Anti-Human Immunodeficiency Virus Antibodies by Enzyme-Linked Immunosorbent Assay of pooled Compared to Singleton Serum Samples. *J. of Clin Microbiol.* Aug. 2006 p.2909-2913.
- Apetrei C, Buzdugan I, Mitroi I, Duca M. The clinical and immunological correlations between the p24 antigenemia levels and those of anti-p24 antibodies in HIV-seropositive children. *Bacteriol Virusol Parazitol Epidemiol.* 1995 Apr-Jun;40(2):141-4. Romanian.
- Shumai EP, Vorob'ev SM, Makarova NE, Tugizov ShM, Zverev VV, Kushch AA. The immunoenzyme detection of the HIV-1 antigen by using monoclonal antibodies to protein p24Vopr *Virusol.* 1992 Sep-Dec;37(5-6):229-32. Russian.
- d'Arminio Monforte A, Novati R, Marchisio P, Zanchetta N, Uberti-Foppa C, Tornaghi R, Massironi E, Lazzarin A, Principi N. Early diagnosis of HIV infection in infants. *AIDS.* 1989 Jun;3(6):391-5.
- Borghj V, De Rienzo B, Pietrosomoli P, Pecorari M, Mongiardo N, Pellegrino F, Zanchetta GP, Lami G, Squadrini F. Detection of serum HIV-Ag related to the major core protein (p24) in persons at risk for AIDS. *Microbiologica.* 1989 Jan;12(1):81-3.
- Goudsmit J, Lange JM, Krone WJ, Teunissen MB, Epstein LG, Danner SA, van den Berg H, Breederveld C, Smit L, Bakker M, et al. Pathogenesis of HIV and its implications for serodiagnosis and monitoring of antiviral therapy. *J Virol Methods.* 1987 Aug;17(1-2):19-34.
- Lyamuya E, Bredberg-Raden U, Massawe A, Urassa E, Kawo G, Msemo G, Kazimoto T, Ostborn A, Karlsson K, Mhalu F, Biberfeld G. Performance of a modified HIV-1 p24 antigen assay for early diagnosis of HIV-1 infection in infants and prediction of mother-to-infant transmission of HIV-1 in Dar es Salaam, Tanzania. *J Acquir Immune Defic Syndr Hum Retrovirol.* 1996 Aug 1;12(4):421-6.
- Clavel, F., Mansinho, K., Chamaret, S. et al. 1987. Human immunodeficiency virus type 2 infection associated with AIDS in West Africa. *N. Engl. J. Med.* 316:1180-1185.
- Sinicco, A., Fora, R., Sciadra, M., Lucchini, A., Caramello, P. and Giannini, P. 1993. Risk of developing AIDS after primary acute HIV-1 infection. *J. of Acquired Immune Deficiency Syndromes* 6:575-581.
- Ju Lin, Hsiang. 1995. Laboratory tests for human immunodeficiency viruses. *J. Intl. Fed. Clin. Chem.* 7(2):61-66.
- IVD Directive 98/79/CE, Common Technical Specifications (CTS) – Annex II, List A.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

HSV1 IgG

**Enzyme ImmunoAssay (ELISA) for the
quantitative/qualitative determination
of IgG antibodies to
Herpes Simplex Virus type 1
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF HSV1G.CE
96 Tests

HSV1 IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Herpes Simplex Virus type 1 in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Herpes Simplex Virus type 1 (HSV1) and type 2 (HSV2) are large complex DNA-containing viruses which have been shown to induce the synthesis of several proteins during infection, possessing an high number of crossreactive determinants and just a few of type-specific sequences.

The majority of primary and recurrent genital herpetic infections are caused by HSV2; while non genital infections, such as common cold sores, are caused primarily by HSV1.

The detection of virus specific IgG and IgM antibodies are important in the diagnosis of acute/primary virus infections or reactivations of a latent one, in the absence of evident clinical symptoms.

Asymptomatic infections may happen for HSV in apparently healthy individuals and during pregnancy. Severe herpetic infections may happen in immunocompromised and suppressed patients in which the disease may evolve toward critical pathologies.

The determination of HSV specific antibodies has then become important in the monitoring of "risk" patients and in the follow up of acute and severe infections.

C. PRINCIPLE OF THE TEST

Microplates are coated with native inactivated HSV1.

The solid phase is first treated with the diluted sample and IgG to HSV are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti HSV1 IgG are detected by the addition of polyclonal specific anti hlgG antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti HSV1 IgG antibodies present in the sample. A Calibration Curve, calibrated against an internal Gold Standard, makes possible a quantitative determination of the IgG antibody in the patient.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate: MICROPLATE

n° 1. 12 strips x 8 microwells coated with native UV inactivated HSV1 in presence of bovine proteins.

Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve derived from human plasma positive for HSV1 IgG ranging:

4ml CAL1 = 0 arbU/ml
4ml CAL2 = 5 arbU/ml
2ml CAL3 = 10 arbU/ml
2ml CAL4 = 20 arbU/ml
2mlCAL5 = 50 arbU/ml
4ml CAL6 = 100 arbU/ml.

Standards are calibrated in arbitrary units against an internal Gold Standard (or IGS).

It contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized. It contains fetal bovine serum proteins, human IgG antibodies to HSV1 at about 20 arbU/ml ±20%, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

4. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

5. Enzyme conjugate : CONJ

2x8ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300, 0.02% gentamicine sulphate as preservatives and 0.01% red alimentary dye.

6. Chromogen/Substrate: SUBS TMB

1x16ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H₂O₂).

Note: To be stored protected from light as sensitive to strong illumination.

7. Sulphuric Acid: H₂SO₄ 0.3 M

1x15ml/vial. It contains 0.3 M H₂SO₄ solution.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363)

8. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide 0.1% and 0.045% ProClin 300 as preservatives. The reagent is blue colour coded.

9. Plate sealing foils n°2

10. Package insert n°1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (1000 ul, 100 ul and 10 ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-0.5°C tolerance)..
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

- The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
- Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
- Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
- Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
- Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
- Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
- Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
- Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
- Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
- Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
- The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water
- Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND WARNINGS

- Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.
- Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storing.

In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C.

After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Calibration Curve

Ready to use component. Mix carefully on vortex before use.

Control Serum

Add the volume of ELISA grade water, reported on the label, to the lyophilised powder; let fully dissolve and then gently mix on vortex.

Note: *The control after dissolution is not stable. Store frozen in aliquots at -20°C.*

Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: *Once diluted, the wash solution is stable for 1 week at +2..8°C.*

Enzyme conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

If this component has to be transferred use only plastic, possibly sterile disposable containers.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container

Sample Diluent

Ready to use component. Mix carefully on vortex before use.

Sulphuric Acid:

Ready to use. Mix well on vortex before use.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 - Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.
An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the lyophilised Control Serum as reported in the proper section.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.
12. Check that all the other equipment is available and ready to use.
13. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

The kit may be used for quantitative and qualitative determinations as well.

M1. QUANTITATIVE DETERMINATION:

Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 1000 µl Sample Diluent and then 10 µl sample (1:101 dilution factor). The whole content is then dispensed into a properly defined dilution tube. Before the next sample is aspirated, needles have

to be duly washed to avoid any cross-contamination among samples. When all the samples have been diluted make the instrument dispense 100 µl samples into the proper wells of the microplate.

This procedure may be carried out also in two steps of dilutions of 1:10 each (90 µl Sample Diluent + 10 µl sample) into a second dilution platform. Make then the instrument aspirate first 100 µl Sample Diluent, then 10 µl liquid from the first dilution in the platform and finally dispense the whole content in the proper well of the assay microplate.

Do not dilute Calibrators and the dissolved Control Serum as they are ready to use.

Dispense 100 µl calibrators/control in the appropriate calibration/control wells.

For the next operations follow the operative instructions reported below for the Manual Assay.

It is strongly recommended to check that the time lap between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

Manual assay:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
3. Dispense 100 µl of Calibrators and 100 µl Control Serum in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

5. Wash the microplate with an automatic washer as reported previously (section I.3).
6. Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

7. Incubate the microplate for **60 min at +37°C**.
8. Wash microwells as in step 5.
9. Pipette 100 µl Chromogen/Substrate mixture into each well, the blank wells A1 and B1 included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

10. Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
11. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

M2. QUALITATIVE DETERMINATION

If only a qualitative determination is required, proceed as described below:

Automated assay:

Proceed as described in section M1.

Manual assay:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
3. Dispense 100 µl of Calibrator 0 arbU/ml and Calibrator 5 arbU/ml in duplicate and Calibrator 100 arbU/ml in single. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

5. Wash the microplate with an automatic washer as reported previously (section I.3).
6. Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

7. Incubate the microplate for **60 min at +37°C**.
8. Wash microwells as in step 5.
9. Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

10. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
11. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1.

General Important notes:

1. Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
2. Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

N. ASSAY SCHEME

Method	Operations
Calibrators & Control (*)	100 µl
Samples diluted 1:101	100 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

(*) Important Notes:

- The Control Serum (CS) it does not affect the test's results calculation.
- The Control Serum (CS) used only if a laboratory internal quality control is required by the Management.

An example of dispensation scheme for Quantitative Analysis is reported below:

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL4	S 1									
B	BLK	CAL4	S 2									
C	CAL1	CAL5	S 3									
D	CAL1	CAL5	S 4									
E	CAL2	CAL6	S 5									
F	CAL2	CAL6	S 6									
G	CAL3	CS(*)	S 7									
H	CAL3	CS(*)	S 8									

Legenda: BLK = Blank CAL = Calibrator
CS(*) = Control Serum - Not mandatory S = Sample

An example of dispensation scheme in qualitative assays is reported below:

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S 3	S 11									
B	CAL1	S 4	S 12									
C	CAL1	S 5	S 13									
D	CAL2	S 6	S 14									
E	CAL2	S 7	S 15									
F	CAL6	S 8	S 16									
G	S 1	S 9	S 17									
H	S 2	S 10	S 18									

Legenda: BLK = Blank CAL = Calibrators
S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the calibrators any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

Check	Requirements
Blank well	< 0.050 OD450nm value
CAL 1 0 arbU/ml	< 0.150 mean OD450nm value after blanking coefficient of variation < 30%
CAL 2 5 arbU/ml	OD450nm > OD450nm CAL1 + 0.100
CAL 6 100 arbU/ml	OD450nm > 1.000

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.050 OD450nm	1. that the Chromogen/Substrate solution has not got contaminated during the assay
CAL 1 0 arbU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
CAL 2 5 arbU/ml OD450nm < OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 100 arbU/ml < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

Should one of these problems have happened, after checking, report to the supervisor for further actions.

** Note:

If Control Serum has used, verify the following data:

Check	Requirements
Control Serum	Mean OD450nm CAL 4 ± 20%

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Control Serum Different from expected value	<ol style="list-style-type: none"> 1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the control serum has occurred.

Anyway, if all other parameters (Blank, CAL1, CAL2, CAL6), match the established requirements, the test may be considered valid.

P. RESULTS


P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm (4-parameters interpolation is suggested).

Then on the calibration curve calculate the concentration of anti Herpes Simplex Virus type 1 IgG antibody in samples.

An example of Calibration curve is reported in the next page.

Example of Calibration Curve :

Important Note:

Do not use the calibration curve above to make calculations.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 5 arbU/ml and then check that the assay is valid.

Example of calculation:

The following data must not be used instead of real figures obtained by the user.

Calibrator 0 arbU/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted

Calibrator 5 arbU/ml: 0.350 – 0.370 OD450nm
 Mean Value: 0.360 OD450nm
 Higher than Cal 0 + 0.100 – Accepted

Calibrator 100 arbU/ml: 2.245 OD450nm
 Higher than 1.000 – Accepted

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 5 arbU/ml are considered negative for anti HSV1 IgG antibody.

Samples with a concentration higher than 5 arbU/ml are considered positive for anti HSV1 IgG antibody.

Particular attention in the interpretation of results has to be used in the follow-up of pregnancy for a primary infection of HSV due to the risk of neonatal malformations.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. In the follow-up of pregnancy for HSV infection a positive result (presence of IgG antibody > 5 arbU/ml) should be confirmed to ruled out the risk of a false positive result and a false definition of protection.

R. PERFORMANCES

1. Limit of detection

The limit of detection of the assay has been calculated by means of an internal Gold Standard in absence of an international preparation to refer to.

The limit of detection has been calculated as mean OD450nm Calibrator 0 arbU/ml + 5 SD.

The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined in the assay for three lots.

Mean OD450nm values (n = 2)

IgG arbU/ml	HSV1G.PU Lot # 0703	HSV1G.PU Lot # 1203	HSV1G.PU Lot # 0204/2
0	0.077	0.034	0.043
5	0.355	0.404	0.318
10	0.742	0.713	0.516
20	1.254	1.216	0.944
50	1.952	1.928	1.728
100	2.623	2.261	2.072

The assay shows a limit of detection far better than 5 arbU/ml.

In addition the preparation code Accurun n° 150, produced by Boston Biomedica Inc., BBI, USA, was tested in dilutions to determine the limit of its detection and provide a further value of analytical sensitivity

Mean OD450nm values (n = 2)

Dilution	HSV1G.CE Lot # 1004	HSV1G.PU Lot # 1203	HSV1G.PU Lot # 0204/2
1 X	1.248	1.218	1.300
2 X	0.860	0.848	0.876
4 X	0.545	0.526	0.583
8 X	0.315	0.300	0.329
16 X	0.164	0.152	0.148
32 X	0.082	0.064	0.072
0 arbU/ml	0.057	0.050	0.047
5 arbU/ml	0.288	0.355	0.318

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in a performance evaluation study on panels of samples classified positive by a kit US FDA approved. Positive samples from different stage of HSV infection were tested. The value, obtained from the analysis of more than 300 specimens, has been > 98%.

3. Diagnostic specificity:

The diagnostic specificity has been determined on panels of negative samples from not infected individuals, classified negative with a kit US FDA approved.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity.

Frozen specimens have been tested, as well, to check for interferences due to collection and storage.

No interference was observed.

Potentially interfering samples derived from patients with different pathologies (mostly ANA, AMA and RF positive) and from pregnant women were tested.

No crossreaction was observed.

An overall value > 98% of specificity was found when examined on more than 100 specimens.

3. Precision:

It has been calculated on the Calibrator 5 arbU/ml, considered the cut-off of the assay, examined in 16 replicates in three separate runs for three lots.

Results are reported as follows:

HSV1G.CE Lot # 1004

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.292	0.290	0.285	0.289
Std.Deviation	0.024	0.024	0.027	0.025
CV %	8.24	8.28	9.42	8.65

HSV1G.PU: lot 1203

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.365	0.382	0.378	0.375
Std.Deviation	0.022	0.029	0.018	0.023
CV %	6.02	7.59	4.76	6.12

HSV1G.PU: Lot 0204/2

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.322	0.298	0.304	0.308
Std.Deviation	0.018	0.019	0.016	0.018
CV %	5.59	6.38	5.26	5.74

The variability shown in the tables above did not result in sample misclassification.

S. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

Frozen samples containing fibrin particles or aggregates after thawing may generate some false results.

This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

REFERENCES

- Engvall E. and Perlmann P.. J.Immunochemistry 8: 871-874, 1971
- Engvall E. and Perlmann P.. J.Immunol.. 109: 129-135, 1971
- Remington J.S. and Klein J.O.. (1996) In "Infectious diseases of fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
- Volk W.A. (1982) In "essential of Medical Microbiology". 2nd ed., pp 729, G.B. Lippincott Co. Philadelphia, New York, S.Josè, Toronto.
- Leinikki P.O. et al.. J.Clin.Microbiol.. 8:418, 1978
- Piroid E. et al.. Révue Méd.Vet.. 131:25, 1980.
- Vaheri A. et al.. J.Med.Virol.. 5:171, 1980.
- Vejtorp M. et al.. Acta Path.Microbiol.Scand.. 88:349, 1980.
- Voller A. et al.. Brit.J.Exp.Pathol.. 56:338, 1975

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HSV1 IgG CODE: HSV1G.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED PRODUCT MEETS
THE PROVISIONS OF THE COUNCIL DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

ISO CERTIFICATE	UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY AEMPS (AGENCIA ESPAÑOLA DE MEDICAMENTOS Y PRODUCTOS SANITARIOS)
------------------------	---

PLACE & DATE OF FIRST ISSUE	MILANO – OCTOBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

HSV1 IgM

**“Capture” Enzyme Immuno Assay
(ELISA) for the determination
of IgM antibodies to
Herpes Simplex Virus type 1
in human plasma and sera**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

REF HSV1M.CE
96 tests

HSV1 IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of IgM antibodies to Herpes Simplex Virus types 1 in human plasma and sera with the "capture" system. The device is intended for the follow-up of HSV1 infected patients and for the monitoring of risk of neonatal defects due to HSV infection during pregnancy.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Herpes Simplex Virus type 1 (HSV1) and type 2 (HSV2) are large complex DNA-containing viruses which have been shown to induce the synthesis of several proteins during infection, possessing an high number of cross-reactive determinants and just a few of type-specific sequences.

The majority of primary and recurrent genital herpetic infections are caused by HSV2; while non genital infections, such as common cold sores, are caused primarily by HSV1.

The detection of virus specific IgG and IgM antibodies are important in the diagnosis of acute/primary virus infections or reactivations of a latent one, in the absence of evident clinical symptoms.

A-symptomatic infections may happen for HSV in apparently healthy individuals and during pregnancy. Severe herpetic infections may happen in immuno-compromised and suppressed patients in which the disease may evolve toward critical pathologies.

The determination of HSV specific antibodies has then become important in the monitoring of "risk" patients and in the follow up of acute and severe infections.

C. PRINCIPLE OF THE TEST

The assay is based on the principle of "IgM capture" where IgM class antibodies in the sample are first captured by the solid phase coated with anti hIgM antibody.

After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a preparation of inactivated HSV1, labeled with a HSV1 specific antibody conjugated with peroxidase (HRP).

After incubation, microwells are washed to remove unbound conjugate and then the chromogen/substrate is added.

In the presence of bound conjugate the colorless substrate is hydrolyzed to a colored end-product, whose optical density may be detected and is proportional to the amount of IgM antibodies to HSV1 present in the sample.

A system is described how to control whether the positivity shown by a sample is true or not (Confirmation Test), helpful for the clinician to make a correct interpretation of results.

D. COMPONENTS

The kit contains reagents for 96 tests.

1. Microplate: MICROPLATE

12 strips x 8 microwells coated with anti human IgM affinity purified goat antibody, in presence of bovine proteins. Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Negative Control: CONTROL -

1x4.0 ml/vial. Ready to use control. It contains 1% human serum proteins, 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1%

Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives.

The negative control is cpale yellow color coded..

3. Positive Control: CONTROL +

1x4.0 ml/vial. Ready to use control. It contains 1% human serum positive for HSV1 IgM, 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives.

The positive control is green colour coded.

4. Calibrator: CAL ...ml

N° 1 lyophilized vial. To be dissolved with EIA grade water as reported in the label. It contains anti HSV1 IgM, fetal bovine serum, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

5. Lyophilized HSV1 Ag: AG HSV1

N° 6 lyophilized vials. The vials contain gamma-ray inactivated HSV1 in protein buffer. The solution contains 2% bovine proteins, 10 mM Tris HCl buffer pH 6.8+/-0.1, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC. To be dissolved with 1.9 ml of Antigen Diluent as reported in the specific section.

6. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle. 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

7. Enzyme conjugate: CONJ 20X

1x0.8 ml/vial. 20x concentrated solution of a HSV1-specific antibody, labeled with HRP and diluted in a protein buffer containing 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.1% Kathon GC and 0.2 mg/ml gentamicine sulphate as preservatives.

8. Antigen Diluent : AG DIL

n° 1 vial of 16 ml. Protein buffer solution for the preparation of the Immunocomplex. The solution contains 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.1% Kathon GC and 0.2 mg/ml gentamicine sulphate as preservatives. The reagent is code coloured with 0.01% red alimentary dye

9. Specimen Diluent : DILSPE

2x60.0 ml/vial. Proteic buffered solution for the dilution of samples. It contains 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives.

The reagent is color coded with 0.01% blue alimentary dye.

10. Chromogen/Substrate : SUBS TMB

1x16ml/vial. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine (TMB), 0.02% hydrogen peroxide (H₂O₂) and 4% dimethylsulphoxide.

Note: To be stored protected from light as sensitive to strong illumination.

11. Sulphuric Acid: H₂SO₄ 0.3M

1x15ml/vial. It contains 0.3 M H₂SO₄ solution.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363)

12. Plate sealing foils n° 2

13. Package insert n° 1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (1000 ul, 100 ul and 10 ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-0.5°C tolerance)..
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

14. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water

16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND WARNINGS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.
3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
4. Sera and plasma can be stored at +2°..8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
5. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storing. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Negative Control:

Ready to use. Mix well on vortex before use.

Positive Control:

Ready to use. Mix well on vortex before use.

Calibrator:

Add the volume of ELISA grade water reported on the label to the lyophilized powder. Let fully dissolve and then gently mix on vortex.

Important Note: *The solution is not stable. Store the Calibrator frozen in aliquots at -20°C.*

Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: *Once diluted, the wash solution is stable for 1 week at +2..8° C.*

Ag/Ab Immunocomplex:

Proceed carefully as follows:

1. Dissolve the content of a lyophilized vial with 1.9 ml of Conjugate/Antigen Diluent. Let fully dissolved the lyophilized content and then gently mix on vortex.
2. Gently mix the concentrated Enzyme Conjugate on vortex. Then add 0.1 ml of it to the vial of the dissolved HSV1 Ag and mix gently on vortex.

Important Notes:

1. Dissolve and prepare only the number of vials necessary to the test. The Immunocomplex obtained is not stable. Store any residual solution frozen in aliquots at -20°C .
2. The preparation of the Immunocomplex has to be done **right before** the dispensation of samples and controls into the plate. Mix again on vortex gently just before its use.

Specimen Diluent:

Ready to use. Mix well on vortex before use

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container

Sulphuric Acid:

Ready to use. Mix well on vortex before use.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 - Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation

of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label (primary container). Do not use the device if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator as described above and gently mix.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.

- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

M.1 Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 1000 µl Specimen Diluent and then 10 µl sample (1:101 dilution factor). The whole content is then dispensed into a properly defined dilution tube. Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples. When all the samples have been diluted make the instrument dispense 100 µl diluted samples into the proper wells of the microplate.

This procedure may be carried out also in two steps of dilutions of 1:10 each (90 µl Specimen Diluent + 10 µl sample) into a second dilution platform. Make then the instrument aspirate first 100 µl Specimen Diluent, then 10 µl liquid from the first dilution in the platform and finally dispense the whole content in the proper well of the assay microplate.

Do not dilute controls/calibrator as they are ready to use.

Dispense 100 µl calibrators/controls in the appropriate calibration/control wells.

For the next operations follow the operative instructions reported below for the Manual Assay.

It is strongly recommended to check that the time lap between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

M. 2 Manual assay:

- Dilute samples 1:101 by dispensing first 10 µl sample and then 1 ml Specimen Diluent into a dilution tube; mix gently on vortex.
- Place the required number of Microwells in the microwell holder. Leave the well in position A1 empty for the operation of blanking.
- Dispense 100 µl of Negative Control and 100 µl of Calibrator in the proper wells in duplicate. Dispense 100 µl of Positive Control in single into the proper well. Do not dilute controls and the calibrator as they are ready to use!
- Dispense 100 µl diluted samples in the proper sample wells and then check that all the samples wells are blue colored and that controls and calibrator have been dispensed.
- Incubate the microplate for **60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 µl of the **Ag/Ab Immunocomplex** into each well, except the blanking well A1, and cover with the sealer. Check that all wells are red colored, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the **Ag/Ab Immunocomplex**. Contamination might occur.

- Incubate the microplate for **60 min at +37°C**.
- Wash microwells as in step 6.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10. Addition of acid will turn the positive control and positive samples from blue to yellow.
- Measure the color intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

Important notes:

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

N. ASSAY SCHEME

Controls&calibrator(*)	100 ul
Samples diluted 1:101	100 ul
1st incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
Immunocomplex	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
TMB/H2O2 mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm

(*) Important Notes:

- The Calibrator (CAL) does not affect the Cut Off calculation, therefore it does not affect the test's results calculation.
- The Calibrator (CAL) used only if a laboratory internal quality control is required by the Management.

An example of dispensation scheme is reported below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3											
B	NC	S4											
C	NC	S5											
D	CAL(*)	S6											
E	CAL(*)	S7											
F	PC	S8											
G	S1	S9											
H	S2	S10											

Legenda: BLK = Blank NC = Negative Control
 CAL(*) = Calibrator-Not mandatory PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

Parameter	Requirements
Blank well	< 0.05 OD450nm value
Negative Control mean value (NC)	< 0.200 OD450nm value after blanking coefficient of variation < 30%
Positive Control	> 1.000 OD450nm

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.05 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.200 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

** Important Note:

If the Calibrator has used, verify the following data:

Check	Requirements
Calibrator	S/Co > 1.2

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Calibrator S/Co < 1.2	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (e.g.: dispensation of negative control instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

Anyway, if all other parameters (Blank, Negative Control, Positive Control), match the established requirements, the test may be considered valid.

P. CALCULATION OF THE CUT-OFF

The test results are calculated by means of the mean OD450nm value of the Negative Control (NC) and a mathematical calculation, in order to define the following cut-off formulation:

$$\text{Cut-Off} = \text{NC} + 0.250$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as a ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 1.0	Negative
1.0 - 1.2	Equivocal
> 1.2	Positive

A negative result indicates that the patient is not undergoing an acute infection of Herpes Simplex Virus type 1.

Any patient showing an equivocal result, should be re-tested by examining a second sample taken from the patient after 1-2 weeks from first testing.

A positive result is indicative of a Herpes Simplex Virus type 1 infection.

An example of calculation is reported below:

Important Note: The following data must not be used instead or real figures obtained by the user.

Negative Control: 0.100 – 0.120 – 0.080 OD450nm

Mean Value: 0.100 OD450nm

Lower than 0.150 – Accepted

Positive Control: 1.850 OD450nm

Higher than 1.000 – Accepted

$$\text{Cut-Off} = 0.110 + 0.250 = 0.360$$

Calibrator: 1.000 - 0.900 OD450nm

Mean value: 0.950 OD450nm

S/Co = 2.6

S/Co higher than 1.2 – Accepted

Sample 1: 0.075 OD450nm

Sample 2: 1.580 OD450nm

Sample 1 S/Co < 1 = negative

Sample 2 S/Co > 1.2 = positive

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Particular attention in the interpretation of results has to be used in the follow-up of pregnancy for an infection of HSV due to the risk of severe neonatal malformations.
3. In pregnancy monitoring, it is strongly recommended that any positive result is confirmed first with the procedure described below and secondly with a different device for HSV IgM detection, before taking any preventive medical action.
4. Any positive sample should be submitted to the Confirmation Test reported in section T before giving a result of positivity. By carrying out this test, false reactions, leading to a misinterpretation of the analytical result, can be revealed and then ruled out.
5. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
6. Diagnosis of infection has to be taken and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCE CHARACTERISTICS

1. Limit of detection

No international standard for HSV1&2 IgM Antibody detection has been defined so far by the European Community.

In its absence, an Internal Gold Standard (or IGS), calibrated on the preparation named "Accurun – Anti HSV2 IgM plasma" produced by Boston Biomedica Inc., USA, code 9106072, has been defined in order to provide the device with a constant and excellent sensitivity.

The limit of detection of the assay has been therefore calculated on the IGS. A limiting dilution curve was prepared in the Negative Control (NC).

Results of Quality Control are given in the following table:

OD450nm values

IGS	HSV1M.CE Lot # RD1	HSV1M.CE Lot # RD2	HSV1M.CE Lot # RD3
1X	0.450	0.460	0.455
2X	0.277	0.300	0.288
4X	0.216	0.198	0.185
NC	0.115	0.085	0.086

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in a performance evaluation study on panels of 40 samples classified positive by a CE marked kit. The value obtained from the analysis was > 98%.

3. Diagnostic specificity:

The diagnostic specificity has been determined in the performance evaluation on panels of more than 300 specimens, negative with the reference kit, derived from normal individuals of European origin.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples.

A study conducted on more than 60 potentially cross-reactive samples has not revealed any interference in the system.

No cross reaction were observed.

The Performance Evaluation has provided a value > 98%.

False positive reactions may be anyway pointed out and then ruled out in the interpretation of results with the procedure reported in section T, able to verify whether or not a positive result is real.

4. Precision:

Results are reported as follows:

HSV1M.CE: lot # RD1

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.083	0.107	0.116	0.102
Std.Deviation	0.004	0.017	0.013	0.011
CV %	5.12	15.82	11.59	10.84

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.393	0.436	0.421	0.417
Std.Deviation	0.031	0.019	0.007	0.019
CV %	7.93	4.38	1.68	4.66

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.469	1.530	1.541	1.513
Std.Deviation	0.034	0.055	0.037	0.042
CV %	2.31	3.60	2.39	2.77

HSV1M.CE: lot # RD2

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.101	0.099	0.097	0.099
Std.Deviation	0.009	0.011	0.013	0.011
CV %	8.91	11.11	13.40	11.14

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.412	0.395	0.420	0.409
Std.Deviation	0.015	0.009	0.012	0.012
CV %	3.64	2.27	2.86	2.92

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.512	1.498	1.534	1.515
Std.Deviation	0.042	0.035	0.028	0.035
CV %	2.78	2.34	1.83	2.31

HSV1M.CE: lot # RD3

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.095	0.112	0.092	0.100
Std.Deviation	0.012	0.009	0.010	0.011
CV %	12.6	8.04	10.86	10.50

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.405	0.398	0.412	0.405
Std.Deviation	0.012	0.015	0.014	0.014
CV %	2.96	3.77	3.40	3.37

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.489	1.475	1.518	1.494
Std.Deviation	0.025	0.032	0.028	0.028
CV %	1.68	2.17	1.84	1.90

S. LIMITATIONS

Frozen samples containing fibrin particles or aggregates may generate false positive results.

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

T. CONFIRMATION TEST

In order to provide the medical doctor with the best accuracy in the follow-up of pregnancy, where a false positive result could lead to an operation of abortion, a confirmation test is reported. The confirmation test has to be carried out on any positive sample before a diagnosis of primary infection of HSV is released to the doctor.

Proceed for confirmation as follows:

1. Prepare the Antigen/Conjugate Complex as described in the proper section. This reagent is called Solution A.
2. Then 25 ul concentrated Enzymatic Conjugate are diluted in 500 ul Antigen Diluent and mixed gently on vortex. Do not use any lyophilized antigen vial for this procedure ! This solution is called Solution B.
3. The well A1 of the strip is left empty for blanking.
4. The Negative Control is dispensed in the strip in positions B1+C1. This is used for the calculation of the cut-off and S/Co values.
5. The positive sample to be confirmed, diluted 1:101, is dispensed in the strip in position D1+E1.
6. The strip is incubated for 60 min at +37°C.
7. After washing, the blank well A1 is left empty.
8. 100 µl of Solution A are dispensed in wells B1+C1+D1.
9. Then 100 µl of Solution B are added to well E1.
10. The strip is incubated for 60 min at +37°C.
11. After washing, 100 µl Chromogen/Substrate are added to all the wells and the strip is incubated for 20 min at r.t.
12. 100 µl Sulphuric Acid are added to all the wells and then their color intensity is measured at 450nm (reading filter) and at 620-630nm (background subtraction,), blanking the instrument on A1.

Interpretation of results is carried out as follows:

1. If the sample in position D1 shows a S/Co value lower than 1.0 a problem of dispensation or contamination in the first test is likely to be occurred. The Assay Procedure in Section M has to be repeated to double check the analysis.
2. If the sample in position D1 shows a S/Co value higher than 1.2 and in position E1 shows a S/Co value still higher than 1.2 the sample is considered a **false positive**. The reactivity of the sample is in fact not dependent on the specific presence of HSV1 and a crossreaction with enzymatic conjugate has occurred.
3. If the sample in position D1 shows a S/Co value higher than 1.2 and in position E1 shows a S/Co value lower than 1.0 the sample is considered a **true positive**. The reactivity of the sample is in fact dependent on the specific presence of HSV and not due to any crossreaction.

The following table is reported for the interpretation of results

Well	S/Co		
	< 1.0	> 1.2	> 1.2
D1	< 1.0	> 1.2	> 1.2
E1	< 1.0	> 1.2	< 1.0
Interpretation	Problem of contam.	False positive	True positive

REFERENCES

1. Engvall E. and Perlmann P.. J.Immunochemistry 8: 871-874, 1971
2. Engvall E. and Perlmann P.. J.Immunol.. 109: 129-135, 1971
3. Remington J.S. and Klein J.O.. (1996) In "Infectious diseases of fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
4. Volk W.A. (1982) In "essential of Medical Microbiology". 2nd ed., pp 729, G.B. Lippincott Co. Philadelphia, New York, S.Josè, Toronto.
5. Leinikki P.O. et al.. J.Clin.Microbiol.. 8:418, 1978
6. Piroid E. et al.. Rêvue Méd.Vet.. 131:25, 1980.
7. Vaheeri A. et al.. J.Med.Virol.. 5:171, 1980.
8. Vejtorp M. et al.. Acta Path.Microbiol.Scand.. 88:349, 1980.
1. Voller A. et al.. Brit.J.Exp.Pathol.. 56:338, 1975

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HSV1 IgM CODE: HSV1M.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED PRODUCT MEETS
THE PROVISIONS OF THE COUNCIL DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

ISO CERTIFICATE	UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY AEMPS (AGENCIA ESPAÑOLA DE MEDICAMENTOS Y PRODUCTOS SANITARIOS)
------------------------	---

PLACE & DATE OF FIRST ISSUE	MILANO – OCTOBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

HSV2 IgG

**Enzyme ImmunoAssay (ELISA) for the
quantitative/qualitative determination
of IgG antibodies to
Herpes Simplex Virus type 2
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF HSV2G.CE
96 Tests

HSV2 IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Herpes Simplex Virus type 2 in human plasma and sera.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Herpes Simplex Virus type 1 (HSV1) and type 2 (HSV2) are large complex DNA-containing viruses which have been shown to induce the synthesis of several proteins during infection, possessing an high number of crossreactive determinants and just a few of type-specific sequences.

The majority of primary and recurrent genital herpetic infections are caused by HSV2; while non genital infections, such as common cold sores, are caused primarily by HSV1.

The detection of virus specific IgG and IgM antibodies are important in the diagnosis of acute/primary virus infections or reactivations of a latent one, in the absence of evident clinical symptoms.

Asymptomatic infections may happen for HSV in apparently healthy individuals and during pregnancy. Severe herpetic infections may happen in immunocompromised and suppressed patients in which the disease may evolve toward critical pathologies.

The determination of HSV specific antibodies has then become important in the monitoring of "risk" patients and in the follow up of acute and severe infections.

C. PRINCIPLE OF THE TEST

Microplates are coated with synthetic HSV2 specific glycoprotein G or gG.

The solid phase is first treated with the diluted sample and IgG to HSV2 are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti HSV2 IgG are detected by the addition of polyclonal specific anti hlgG antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti HSV2 IgG antibodies present in the sample. A Calibration Curve, calibrated against an internal Gold Standard, makes possible a quantitative determination of the IgG antibody in the patient.

D. COMPONENTS

Each kit contains sufficient reagents to perform 96 tests.

1. Microplate: MICROPLATE

n° 1. 12 strips x 8 microwells coated with synthetic HSV2-specific gG in presence of bovine proteins.

Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve derived from human plasma positive for HSV2 IgG ranging:

4ml CAL1 = 0 arbU/ml
4ml CAL2 = 5 arbU/ml
2ml CAL3 = 10 arbU/ml
2ml CAL4 = 20 arbU/ml
2ml CAL5 = 50 arbU/ml
4ml CAL6 = 100 arbU/ml.

Standards are calibrated in arbitrary units against an internal Gold Standard (or IGS).

It contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and ProClin 300 0.045% as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized. It contains fetal bovine serum proteins, human IgG antibodies to HSV2 at about 20 arbU/ml ± 20%, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

4. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

5. Enzyme conjugate : CONJ

2x8ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300, 0.02% gentamicine sulphate as preservatives and 0.01% red alimentary dye.

6. Chromogen/Substrate: SUBS TMB

1x16ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H₂O₂).

Note: To be stored protected from light as sensitive to strong illumination.

7. Sulphuric Acid: H₂SO₄ 0.3 M

1x15ml/vial. It contains 0.3 M H₂SO₄ solution.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363)

8. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide 0.1% and 0.045% ProClin 300 as preservatives. The reagent is blue colour coded.

9. Plate sealing foils n°2

10. Package insert n°1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (1000 ul, 100 ul and 10 ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-0.5°C tolerance)..
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses.

The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
14. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND WARNINGS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.

3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

4. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection.

Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

5. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storage.

In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C.

After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Calibration Curve

Ready to use component. Mix carefully on vortex before use.

Control Serum

Add the volume of ELISA grade water, reported on the label, to the lyophilised powder; let fully dissolve and then gently mix on vortex.

Note: *The control after dissolution is not stable. Store frozen in aliquots at -20°C.*

Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: *Once diluted, the wash solution is stable for 1 week at +2..8° C.*

Enzyme conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

If this component has to be transferred use only plastic, possibly sterile disposable containers.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container

Sample Diluent

Ready to use component. Mix carefully on vortex before use.

Sulphuric Acid:

Ready to use. Mix well on vortex before use.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 - Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.
An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, mandatory) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and

validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the lyophilised Control Serum as reported in the proper section.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.
12. Check that all the other equipment is available and ready to use.
13. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

The kit may be used for quantitative and qualitative determinations as well.

M1. QUANTITATIVE DETERMINATION:

Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 1000 µl Sample Diluent and then 10 µl sample (1:101 dilution factor). The whole content is then dispensed into a properly defined dilution tube. Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples. When all the samples have been diluted make the instrument dispense 100 µl samples into the proper wells of the microplate.

This procedure may be carried out also in two steps of dilutions of 1:10 each (90 µl Sample Diluent + 10 µl sample) into a second dilution platform. Make then the instrument aspirate first

100 µl Sample Diluent, then 10 µl liquid from the first dilution in the platform and finally dispense the whole content in the proper well of the assay microplate.

Do not dilute Calibrators and the dissolved Control Serum as they are ready to use.

Dispense 100 µl calibrators/control in the appropriate calibration/control wells.

For the next operations follow the operative instructions reported below for the Manual Assay.

It is strongly recommended to check that the time lap between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

Manual assay:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
3. Dispense 100 µl of Calibrators and 100 µl Control Serum in duplicate. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

Important note: *Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.*

5. Wash the microplate with an automatic washer as reported previously (section I.3).
6. Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

Important note: *Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.*

7. Incubate the microplate for **60 min at +37°C**.
8. Wash microwells as in step 5.
9. Pipette 100 µl Chromogen/Substrate mixture into each well, the blank wells A1 and B1 included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: *Do not expose to strong direct illumination. High background might be generated.*

10. Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
11. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

M2. QUALITATIVE DETERMINATION

If only a qualitative determination is required, proceed as described below:

Automated assay:

Proceed as described in section M1.

Manual assay:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
3. Dispense 100 µl of Calibrator 0 arbU/ml and Calibrator 5 arbU/ml in duplicate and Calibrator 100 arbU/ml in single. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

Important note: *Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.*

5. Wash the microplate with an automatic washer as reported previously (section I.3).
6. Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: *Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.*

7. Incubate the microplate for **60 min at +37°C**.
8. Wash microwells as in step 5.
9. Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: *Do not expose to strong direct illumination. High background might be generated.*

10. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
11. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1.

General Important notes:

1. *Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.*
2. *Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.*

N. ASSAY SCHEME

Method	Operations
Calibrators & Control (*)	100 µl
Samples diluted 1:101	100 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H ₂ O ₂	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

(* Important Notes:

- The Control Serum (CS) it does not affect the test's results calculation.
- The Control Serum (CS) used only if a laboratory internal quality control is required by the Management.

An example of dispensation scheme for Quantitative Analysis is reported below:

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL4	S 1									
B	BLK	CAL4	S 2									
C	CAL1	CAL5	S 3									
D	CAL1	CAL5	S 4									
E	CAL2	CAL6	S 5									
F	CAL2	CAL6	S 6									
G	CAL3	CS(*)	S 7									
H	CAL3	CS(*)	S 8									

Legenda: BLK = Blank CAL = Calibrator
CS(*) = Control Serum - Not mandatory S = Sample

An example of dispensation scheme in qualitative assays is reported below:

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S 2	S 10									
B	CAL1	S 3	S 11									
C	CAL1	S 4	S 12									
D	CAL2	S 5	S 13									
E	CAL2	S 6	S 14									
F	CAL6	S 7	S 15									
G	S 1	S 8	S 16									
H	S 2	S 9	S 17									

Legenda: BLK = Blank CAL = Calibrators
S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the calibrators any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

Check	Requirements
Blank well	< 0.050 OD450nm value
CAL 1 0 arbU/ml	< 0.150 mean OD450nm value after blanking coefficient of variation < 30%
CAL 2 5 arbU/ml	OD450nm ≥ OD450nm CAL1 + 0.100
CAL 6 100 arbU/ml	OD450nm ≥ 1.000

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.050 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
CAL 1 0 arbU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
CAL 2 5 arbU/ml OD450nm ≤ OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 100 arbU/ml ≤ 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

Should one of these problems have happened, after checking, report to the supervisor for further actions.

**** Note:**

If Control Serum has used, verify the following data:

Check	Requirements
Control Serum	Mean OD450nm CAL 4 ± 20%

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Control Serum Different from expected value	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the control serum has occurred.

Anyway, if all other parameters (Blank, CAL1, CAL2, CAL6), match the established requirements, the test may be considered valid.

P. RESULTS


P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm (4-parameters interpolation is suggested).

Then on the calibration curve calculate the concentration of anti Herpes Simplex Virus type 2 IgG antibody in samples.

An example of Calibration curve is reported in the next page.

Example of Calibration Curve :

Important Note:

Do not use the calibration curve above to make calculations.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 5 arbU/ml and then check that the assay is valid.

Example of calculation:

The following data must not be used instead of real figures obtained by the user.

Calibrator 0 arbU/ml: 0.020 – 0.024 OD450nm
 Mean Value: 0.022 OD450nm
 Lower than 0.150 – Accepted

Calibrator 5 arbU/ml: 0.350 – 0.370 OD450nm
 Mean Value: 0.360 OD450nm
 Higher than Cal 0 + 0.100 – Accepted

Calibrator 100 arbU/ml: 2.245 OD450nm
 Higher than 1.000 – Accepted

Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 5 arbU/ml are considered negative for anti HSV2 IgG antibody.

Samples with a concentration higher than 5 arbU/ml are considered positive for anti HSV2 IgG antibody.

Particular attention in the interpretation of results has to be used in the follow-up of pregnancy for a primary infection of HSV due to the risk of neonatal malformations.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. In the follow-up of pregnancy for HSV infection a positive result (presence of IgG antibody > 5 arbU/ml) should be confirmed to ruled out the risk of a false positive result and a false definition of protection.

R. PERFORMANCES

1. Limit of detection

The limit of detection of the assay has been calculated by means of an internal Gold Standard in absence of an international preparation to refer to.

The limit of detection has been calculated as mean OD450nm Calibrator 0 arbU/ml + 5 SD.

The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined in the assay for three lots.

Mean OD450nm values (n = 2)

IgG arbU/ml	HSV2G.PU Lot # 1203	HSV2G.PU Lot # 1103	HSV2G Lot # 0304/2
0	0.022	0.030	0.014
5	0.353	0.384	0.269
10	0.596	0.606	0.557
20	1.169	1.471	0.895
50	2.030	2.276	1.776
100	3.102	3.353	2.893

The assay shows a limit of detection far better than 5 arbU/ml.

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in a performance evaluation study on panels of samples classified positive by a kit US FDA approved. Positive samples from different stage of HSV infection were tested. The value, obtained from the analysis of more than 300 specimens, has been $\geq 98\%$.

3. Diagnostic specificity:

The diagnostic specificity has been determined on panels of negative samples from not infected individuals, classified negative with a kit US FDA approved.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity.

Frozen specimens have been tested, as well, to check for interferences due to collection and storage.

No interference was observed.

Potentially interfering samples derived from patients with different pathologies (mostly ANA, AMA and RF positive) and from pregnant women were tested.

No crossreaction was observed.

An overall value $> 98\%$ of specificity was found when examined on more than 100 specimens.

3. Precision:

It has been calculated on the Calibrator 5 arbU/ml, considered the cut-off of the assay, examined in 16 replicates in three separate runs for three lots.

Results are reported as follows:

HSV2G.CE: lot 1004

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.286	0.303	0.256	0.282
Std.Deviation	0.022	0.037	0.020	0.026
CV %	7.7	12.4	7.74	9.28

HSV2G.PU: lot 1103

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.375	0.384	0.394	0.384
Std.Deviation	0.019	0.022	0.015	0.019
CV %	5.07	5.73	3.81	4.87

HSV2G.PU: lot 1203

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.352	0.345	0.332	0.343
Std.Deviation	0.017	0.020	0.024	0.020
CV %	4.83	5.78	7.23	5.95

The variability shown in the tables above did not result in sample misclassification.

S. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

Frozen samples containing fibrin particles or aggregates after thawing may generate some false results.

This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history,

symptomatology, as well as other diagnostic data should be considered.

REFERENCES

- Engvall E. and Perlmann P.. J.Immunochemistry 8: 871-874, 1971
- Engvall E. and Perlmann P.. J.Immunol.. 109: 129-135, 1971
- Remington J.S. and Klein J.O.. (1996) In "Infectious diseases of fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
- Volk W.A. (1982) In "essential of Medical Microbiology". 2nd ed., pp 729, G.B. Lippincott Co. Philadelphia, New York, S.Josè, Toronto.
- Leinikki P.O. et al.. J.Clin.Microbiol.. 8:418, 1978
- Piroid E. et al.. Révue Méd.Vet.. 131:25, 1980.
- Vaheri A. et al.. J.Med.Virol.. 5:171, 1980.
- Vejtorp M. et al.. Acta Path.Microbiol.Scand.. 88:349, 1980.
- Voller A. et al.. Brit.J.Exp.Pathol.. 56:338, 1975.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia.Pro Diagnostic Bioprobes S.r.l.
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HSV2 IgG CODE: HSV2G.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

**WE HEREBY DECLARE THAT THE ABOVE MENTIONED PRODUCT MEETS
THE PROVISIONS OF THE COUNCIL DIRECTIVE 98/79/EC
FOR IN VITRO DIAGNOSTIC DEVICES.**

ISO CERTIFICATE	UNE EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY AEMPS (AGENCIA ESPAÑOLA DE MEDICAMENTOS Y PRODUCTOS SANITARIOS)
------------------------	---

PLACE & DATE OF FIRST ISSUE	MILANO – OCTOBER 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

HSV2 IgM

**“Capture” Enzyme Immuno Assay
(ELISA) for the determination
of IgM antibodies to
Herpes Simplex Virus type 2
in human plasma and sera**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF HSV2M.CE
96 tests

HSV2 IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of IgM antibodies to Herpes Simplex Virus types 2 in human plasma and sera with the "capture" system. The device is intended for the follow-up of HSV2 infected patients and for the monitoring of risk of neonatal defects due to HSV infection during pregnancy.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Herpes Simplex Virus type 1 (HSV1) and type 2 (HSV2) are large complex DNA-containing viruses which have been shown to induce the synthesis of several proteins during infection, possessing an high number of cross-reactive determinants and just a few of type-specific sequences.

The majority of primary and recurrent genital herpetic infections are caused by HSV2; while non genital infections, such as common cold sores, are caused primarily by HSV1.

The detection of virus specific IgG and IgM antibodies are important in the diagnosis of acute/primary virus infections or reactivations of a latent one, in the absence of evident clinical symptoms.

A-symptomatic infections may happen for HSV in apparently healthy individuals and during pregnancy. Severe herpetic infections may happen in immuno-compromised and suppressed patients in which the disease may evolve toward critical pathologies.

The determination of HSV specific antibodies has then become important in the monitoring of "risk" patients and in the follow up of acute and severe infections.

C. PRINCIPLE OF THE TEST

The assay is based on the principle of "IgM capture" where IgM class antibodies in the sample are first captured by the solid phase coated with anti hIgM antibody.

After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a preparation of inactivated HSV2, labeled with a HSV2 specific antibody conjugated with peroxidase (HRP).

After incubation, microwells are washed to remove unbound conjugate and then the chromogen/substrate is added.

In the presence of bound conjugate the colorless substrate is hydrolyzed to a colored end-product, whose optical density may be detected and is proportional to the amount of IgM antibodies to HSV2 present in the sample.

A system is described how to control whether the positivity shown by a sample is true or not (Confirmation Test), helpful for the clinician to make a correct interpretation of results.

D. COMPONENTS

The kit contains reagents for 96 tests.

1. Microplate: MICROPLATE

12 strips x 8 microwells coated with anti human IgM affinity purified goat antibody, in presence of bovine proteins. Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Negative Control: CONTROL -

1x4.0 ml/vial. Ready to use control. It contains 1% human serum proteins, 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1%

Tween 20, 0.09% sodium azide and 0.045% ProClin 300 as preservatives.

The negative control is pale yellow color coded.

3. Positive Control: CONTROL +

1x4.0 ml/vial. Ready to use control. It contains 1% human serum positive for HSV2 IgM, 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.045% ProClin 300 as preservatives.

The positive control is green colour coded.

4. Calibrator: CAL ...ml

N° 1 lyophilized vial. To be dissolved with EIA grade water as reported in the label. It contains anti HSV2 IgM, fetal bovine serum, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

5. Lyophilized HSV2 Ag: AG HSV2

N° 6 lyophilized vials. The vials contain lyophilized gamma-ray inactivated HSV2 in protein buffer. The solution contains 2% bovine proteins, 10 mM Tris HCl buffer pH 6.8+/-0.1, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300. To be dissolved with 1.9 ml of Antigen Diluent as reported in the specific section.

6. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle. 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

7. Enzyme conjugate: CONJ 20X

1x0.8 ml/vial. 20x concentrated solution of a HSV2-specific antibody, labeled with HRP and diluted in a protein buffer containing 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.2 mg/ml gentamicine sulphate as preservatives.

8. Antigen Diluent : AG DIL

n° 1 vial of 16 ml. Protein buffer solution for the preparation of the Immunocomplex. The solution contains 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.2 mg/ml gentamicine sulphate as preservatives. The reagent is code coloured with 0.01% red alimentary dye

9. Specimen Diluent : DILSPE

2x60.0 ml/vial. Proteic buffered solution for the dilution of samples. It contains 2% casein, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.045% ProClin 300 as preservatives.

The reagent is color coded with 0.01% blue alimentary dye.

10. Chromogen/Substrate : SUBS TMB

1x16ml/vial. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine (TMB), 0.02% hydrogen peroxide (H₂O₂) and 4% dimethylsulphoxide.

Note: To be stored protected from light as sensitive to strong illumination.

11. Sulphuric Acid: H₂SO₄ 0.3 M

1x15ml/vial. It contains 0.3 M H₂SO₄ solution.
Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363)

12. Plate sealing foils n° 2

13. Package insert n° 1

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (1000 ul, 100 ul and 10 ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are

treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

14. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water

16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND WARNINGS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.

3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

4. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

5. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 3 months.

Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storing. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°.8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

Negative Control:

Ready to use. Mix well on vortex before use.

Positive Control:

Ready to use. Mix well on vortex before use.

Calibrator:

Add the volume of ELISA grade water reported on the label to the lyophilized powder. Let fully dissolve and then gently mix on vortex.

Important Note: *The solution is not stable. Store the Calibrator frozen in aliquots at -20°C.*

Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before

use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at +2..8° C.

Ag/Ab Immunocomplex:

Proceed carefully as follows:

1. Dissolve the content of a lyophilized vial with 1.9 ml of Conjugate/Antigen Diluent. Let fully dissolved the lyophilized content and then gently mix on vortex.
2. Gently mix the concentrated Enzyme Conjugate on vortex. Then add 0.1 ml of it to the vial of the dissolved HSV2 Ag and mix gently on vortex.

Important Notes:

1. Dissolve and prepare only the number of vials necessary to the test. The Immunocomplex obtained is not stable. Store any residual solution frozen in aliquots at -20°C.
2. The preparation of the Immunocomplex has to be done **right before** the dispensation of samples and controls into the plate. Mix again on vortex gently just before its use.

Specimen Diluent:

Ready to use. Mix well on vortex before use

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container

Sulphuric Acid:

Ready to use. Mix well on vortex before use.

Attention: Irritant (H315, H319; P280, P302+P352, 332+P313, P305+ P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.

P332 + P313 – If skin irritation occurs: Get medical advice/attention.

P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 – If eye irritation persists: Get medical advice/attention.

P362 + P363 - Take off contaminated clothing and wash it before reuse.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination

of spills or residues of kit components should also be carried out regularly.

2. The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested). 5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing. An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter de 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label (primary container). Do not use the device if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.

- Dissolve the Calibrator as described above and gently mix.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

M.1 Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 1000 µl Specimen Diluent and then 10 µl sample (1:101 dilution factor). The whole content is then dispensed into a properly defined dilution tube. Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples. When all the samples have been diluted make the instrument dispense 100 µl diluted samples into the proper wells of the microplate.

This procedure may be carried out also in two steps of dilutions of 1:10 each (90 µl Specimen Diluent + 10 µl sample) into a second dilution platform. Make then the instrument aspirate first 100 µl Specimen Diluent, then 10 µl liquid from the first dilution in the platform and finally dispense the whole content in the proper well of the assay microplate.

Do not dilute controls/calibrator as they are ready to use.

Dispense 100 µl calibrators/control in the appropriate calibration/control wells.

For the next operations follow the operative instructions reported below for the Manual Assay.

It is strongly recommended to check that the time lap between the dispensation of the first and the last sample will be calculated by the instrument and taken into consideration by delaying the first washing operation accordingly.

M. 2 Manual assay:

- Dilute samples 1:101 by dispensing first 10 µl sample and then 1 ml Specimen Diluent into a dilution tube; mix gently on vortex.
- Place the required number of Microwells in the microwell holder. Leave the well in position A1 empty for the operation of blanking.
- Dispense 100 µl of Negative Control and 100 µl of Calibrator in the proper wells in duplicate. Dispense 100 µl of Positive Control in single into the proper well. Do not dilute controls and the calibrator as they are ready to use !
- Dispense 100 µl diluted samples in the proper sample wells and then check that all the samples wells are blue colored and that controls and calibrator have been dispensed.
- Incubate the microplate for **60 min at +37°C** .

Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.

- Wash the microplate with an automatic as reported previously (section I.3).
- Pipette 100 µl of the **Ag/Ab Immunocomplex** into each well, except the blanking well A1, and cover with the sealer. Check that all wells are red colored, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the **Ag/Ab Immunocomplex**. Contamination might occur.

- Incubate the microplate for **60 min at +37°C** .
- Wash microwells as in step 6.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10. Addition of acid will turn the positive control and positive samples from blue to yellow .
- Measure the color intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1.

Important notes:

- Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

N. ASSAY SCHEME

Controls&calibrator (*)	100 ul
Samples diluted 1:101	100 ul
1st incubation	60 min
Temperature	+37°C
Washing	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Immunocomplex	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2 mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

(*) Important Notes:

- The Calibrator (CAL) does not affect the Cut Off calculation, therefore it does not affect the test's results calculation.
- The Calibrator (CAL) used only if a laboratory internal quality control is required by the Management.

An example of dispensation scheme is reported below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3											
B	NC	S4											
C	NC	S5											
D	CAL(*)	S6											
E	CAL(*)	S7											
F	PC	S8											
G	S1	S9											
H	S2	S10											

Legenda: BLK = Blank NC = Negative Control
 CAL(*) = Calibrator-Not Mandatory PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified.

Control that the following data are matched:

Parameter	Requirements
Blank well	< 0.05 OD450nm value
Negative Control mean value (NC)	< 0.200 OD450nm value after blanking coefficient of variation < 30%
Positive Control	> 1.000 OD450nm

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.05 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.200 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

**** Important Notes:**

The analysis must be done proceeding as the reading step described in the section M, point 12.

If the Calibrator has used, verify the following data:

Check	Requirements
Calibrator	S/Co > 1.2

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Calibrator S/Co < 1.2	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (e.g.: dispensation of negative control instead) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

Anyway, if all other parameters (Blank, Negative Control, Positive Control), match the established requirements, the test may be considered valid.

P. CALCULATION OF THE CUT-OFF

The test results are calculated by means of the mean OD450nm/620-630nm value of the Negative Control (NC) and a mathematical calculation, in order to define the following cut-off formulation:

$$\text{Cut-Off} = \text{NC} + 0.250$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as a ratio of the sample OD450nm/620-630nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 1.0	Negative
1.0 - 1.2	Equivocal
> 1.2	Positive

A negative result indicates that the patient is not undergoing an acute infection of Herpes Simplex Virus type 2.

Any patient showing an equivocal result, should be re-tested by examining a second sample taken from the patient after 1-2 weeks from first testing.

A positive result is indicative of a Herpes Simplex Virus type 2 infection.

An example of calculation is reported below (data obtained proceeding as the reading step described in the section M, point 12).

Important Note: The following data must not be used instead of real figures obtained by the user.

Negative Control: 0.090 – 0.110 – 0.070 OD450nm

Mean Value: 0.100 OD450nm

Lower than 0.200 – Accepted

Positive Control: 1.850 OD450nm

Higher than 1.000 – Accepted

Cut-Off = 0.100+0.250 = 0.350

Calibrator: 0.900 – 1.100 OD450nm

Mean value: 1.000 OD450nm S/Co = 2.8

S/Co higher than 1.2 – Accepted

Sample 1: 0.070 OD450nm

Sample 2: 1.690 OD450nm

Sample 1 S/Co < 1 = negative

Sample 2 S/Co > 1.2 = positive

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Particular attention in the interpretation of results has to be used in the follow-up of pregnancy for an infection of HSV due to the risk of severe neonatal malformations.
3. In pregnancy monitoring, it is strongly recommended that any positive result is confirmed first with the procedure described below and secondly with a different device for HSV IgM detection, before taking any preventive medical action.
4. Any positive sample should be submitted to the Confirmation Test reported in section T before giving a result of positivity. By carrying out this test, false reactions, leading to a misinterpretation of the analytical result, can be revealed and then ruled out.
5. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
6. Diagnosis of infection has to be taken and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCE CHARACTERISTICS

1. Limit of detection

No international standard for HSV1&2 IgM Antibody detection has been defined so far by the European Community.

In its absence, an Internal Gold Standard (or IGS), calibrated on the preparation named "Accurun – Anti HSV2 IgM plasma" produced by Boston Biomedica Inc., USA, code 9106072, has been defined in order to provide the device with a constant and excellent sensitivity..

The limit of detection of the assay has been therefore calculated on the IGS. A limiting dilution curve was prepared in Negative Control (NC).

Results of Quality Control are given in the following table:

OD450nm values

IGS	HSV2M.CE Lot # RD1	HSV2M.CE Lot # RD2	HSV2M.CE Lot # RD3
1X	0.560	0.572	0.590
2X	0.343	0.324	0.348
4X	0.239	0.218	0.225
NC	0.145	0.132	0.139

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in a clinical trial on panels of 40 samples classified positive by a kit US FDA approved. The value obtained from the analysis was > 98%.

3. Diagnostic specificity:

The diagnostic specificity has been determined in a performance evaluation study on panels of more than 300 specimens, negative with the reference kit, derived from normal individuals of European origin.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples.

A study conducted on more than 60 potentially cross-reactive samples has not revealed any interference in the system.

No cross reaction were observed.

The Performance Evaluation has provided a value > 98%.

False positive reactions may be anyway pointed out and then ruled out in the interpretation of results with the procedure reported in section T, able to verify whether or not a positive result is real.

4. Precision:

Results are reported as follows:

HSV2M.CE: lot # RD1

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.092	0.113	0.097	0.101
Std.Deviation	0.011	0.019	0.010	0.013
CV %	12.25	16.83	10.24	13.11

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.451	0.471	0.435	0.452
Std.Deviation	0.018	0.000	0.033	0.017
CV %	3.92	0.00	7.48	3.8

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.530	1.574	1.527	1.543
Std.Deviation	0.023	0.052	0.006	0.027
CV %	1.48	3.33	0.37	1.73

HSV2M.CE: lot # RD2

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.095	0.101	0.097	0.098
Std.Deviation	0.006	0.008	0.005	0.006
CV %	6.30	7.92	5.15	6.45

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.431	0.428	0.453	0.437
Std.Deviation	0.023	0.018	0.023	0.021
CV %	5.3	4.2	5.10	4.9

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.558	1.552	1.541	1.550
Std.Deviation	0.031	0.025	0.039	0.032
CV %	1.98	1.61	2.53	2.04

HSV2M.CE: lot # RD3

Negative (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.104	0.108	0.099	0.104
Std.Deviation	0.015	0.010	0.011	0.012
CV %	14.4	9.2	11.11	11.57

Low reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.425	0.436	0.440	0.434
Std.Deviation	0.008	0.006	0.009	0.008
CV %	1.8	1.4	2.0	1.7

High reactive (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.571	1.562	1.558	1.564
Std.Deviation	0.040	0.034	0.024	0.033
CV %	2.54	2.17	1.54	2.08

Important note:

The performance data have been obtained proceeding as the reading step described in the section M, point 12.

S. LIMITATIONS

Frozen samples containing fibrin particles or aggregates may generate false positive results.

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.

This test is suitable only for testing single samples and not pooled ones.

Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.

T. CONFIRMATION TEST

In order to provide the medical doctor with the best accuracy in the follow-up of pregnancy, where a false positive result could lead to an operation of abortion, a confirmation test is reported. The confirmation test has to be carried out on any positive sample before a diagnosis of primary infection of HSV is released to the doctor.

Proceed for confirmation as follows:

1. Prepare the Antigen/Conjugate Complex as described in the proper section. This reagent is called Solution A.
2. Then 25 ul concentrated Enzymatic Conjugate are diluted in 500 ul Antigen Diluent and mixed gently on vortex. Do not use any lyophilized antigen vial for this procedure ! This solution is called Solution B.
3. The well A1 of the strip is left empty for blanking.
4. The Negative Control is dispensed in the strip in positions B1+C1. This is used for the calculation of the cut-off and S/Co values.
5. The positive sample to be confirmed, diluted 1:101, is dispensed in the strip in position D1+E1.
6. The strip is incubated for 60 min at +37°C.
7. After washing, the blank well A1 is left empty.
8. 100 µl of Solution A are dispensed in wells B1+C1+D1.
9. Then 100 µl of Solution B are added to well E1.
10. The strip is incubated for 60 min at +37°C.
11. After washing, 100 µl Chromogen/Substrate are added to all the wells and the strip is incubated for 20 min at r.t.
12. 100 µl Sulphuric Acid are added to all the wells and then their color intensity is measured at 450nm (reading filter) and at 620-630nm (background subtraction), blanking the instrument on A1.

Interpretation of results is carried out as follows:

1. If the sample in position D1 shows a S/Co value lower than 1.0 a problem of dispensation or contamination in the first test is likely to be occurred. The Assay Procedure in Section M has to be repeated to double check the analysis.
2. If the sample in position D1 shows a S/Co value higher than 1.2 and in position E1 shows a S/Co value still higher than 1.2 the sample is considered a **false positive**. The reactivity of the sample is in fact not dependent on the specific presence of HSV2 and a crossreaction with enzymatic conjugate has occurred.
3. If the sample in position D1 shows a S/Co value higher than 1.2 and in position E1 shows a S/Co value lower than 1.0 the sample is considered a **true positive**. The reactivity of the sample is in fact dependent on the specific presence of HSV and not due to any crossreaction.

The following table is reported for the interpretation of results

Well	S/Co		
D1	< 1.0	> 1.2	> 1.2
E1	< 1.0	> 1.2	< 1.0
Interpretation	Problem of contam.	False positive	True positive

REFERENCES

1. Engvall E. and Perlmann P.. J.Immunochemistry 8: 871-874, 1971
2. Engvall E. and Perlmann P.. J.Immunol.. 109: 129-135, 1971
3. Remington J.S. and Klein J.O.. (1996) In "Infectious diseases of fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
4. Volk W.A. (1982) In "essential of Medical Microbiology". 2nd ed., pp 729, G.B. Lippincott Co. Philadelphia, New York, S.Josè, Toronto.
5. Leinikki P.O. et al.. J.Clin.Microbiol.. 8:418, 1978
6. Piroid E. et al.. Révue Méd.Vet.. 131:25, 1980.
7. Vaheri A. et al.. J.Med.Virol.. 5:171, 1980.
8. Vejtorp M. et al.. Acta Path.Microbiol.Scand.. 88:349, 1980.
9. Voller A. et al.. Brit.J.Exp.Pathol.. 56:338, 1975

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:
Dia.Pro Diagnostic Bioprobes S.r.l.
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HSV2 IgM CODE: HSV2M.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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Rev: 05/2018

RUB IgG

**Enzyme ImmunoAssay (ELISA) for the
quantitative/qualitative determination of
IgG antibodies to Rubella Virus
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

Rub IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Rubella Virus in human plasma and sera.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Rubella is a small spherical enveloped virus, 55-60nm in diameter, and is the only member of the genus Rubivirus of the family Togaviridae.

The virus contains a single positive stranded 42s RNA molecule and only one serotype is known. The virus encodes for at least three envelope glycoproteins, E1, E2a, E2b; a nucleocapsid-associated protein, C; and two nonstructural proteins.

The detection of Rubella-specific IgG and IgM antibodies is very important for the serological diagnosis of both congenital and primary postnatal rubella infections as they can lead to severe birth defects.

The absence of Rubella-specific IgG antibodies in sera, characteristically of long-term duration after primary infections, in presence of virus-specific IgM, is indicative for the risk of defects in newborn infants.

Highly specific Rubella IgG assays provide the clinician with a helpful and reliable test for the monitoring of these risks in pregnancy and for the monitoring of the immunological response upon vaccination.

C. PRINCIPLE OF THE TEST

Microplates are coated with native Rubella Virus, highly purified by sucrose gradient centrifugation and inactivated.

The solid phase is first treated with the diluted sample and IgG to Rubella Virus are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti Rubella IgG are detected by the addition of polyclonal specific anti hlgG antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti Rubella Virus IgG antibodies present in the sample. A Calibration Curve, calibrated against the 1st W.H.O international standard for anti-Rubella immunoglobulin code RUBI-1-94, makes possible a quantitative determination of the IgG antibody in the patient.

D. COMPONENTS

The standard kit contains reagents to perform 96 tests.

1. Microplate: MICROPLATE

12 strips x 8 microwells coated with highly purified and UV inactivated Rubella Virus in presence of bovine proteins.

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2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve derived from human plasma positive for Rubella IgG and titrated on WHO standard ranging:

4 ml CAL1 = 0 WHO IU/ml
 4 ml CAL2 = 10 WHO IU/ml
 2ml CAL3 = 20 WHO IU/ml
 2 ml CAL4 = 50 WHO IU/ml
 2 ml CAL 5 = 100 WHO IU/ml
 4 ml CAL6 = 250 WHO IU/ml.

Standards are calibrated against the 1st W.H.O international standard for anti-Rubella immunoglobulin code RUBI-1-94.

It contains human serum proteins, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. Standards are blue colored.

3. Control Serum: CONTROL ...ml

1 vial. Lyophilized.

It contains fetal bovine serum proteins, human IgG antibodies to Rubella Virus calibrated at 20 WHO IU/ml \pm 10%, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

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4. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

5. Enzyme conjugate : CONJ

2x8ml/vial..Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.1% Kathon GC, 0.02 mg/ml gentamicine sulphate as preservatives and 0.01% red alimentary dye.

6. Chromogen/Substrate: SUBS TMB

1x16ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H₂O₂) and 4% dimethylsulphoxide.

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8. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. The reagent is blue colour coded.

9. Plate sealing foils n°2

10. Package insert n°1

Important note: Only upon specific request, Dia.Pro can supply reagents for 192 and 480 tests, as reported below :

1. Microplate	n°2	n°5
2. Calibration curve	6x4.0ml/vial	6x10ml/vial
3. Control Serum	n° 2 vials	n° 5 vials
4. Wash buff conc	2x60ml/bottle	5x60ml/bottles
5. Enz. Conjugate	2x16ml/vial	2x40ml/bottles
6. Chromog/Subs	2x16ml/vial	2x40ml/bottles
7. Sulphuric Acid	1x30ml/vial	2x40ml/bottles
8. Specimen Diluent	2x120ml/vial	10x60ml/bottles
9. Plate seal foils	n° 4	n° 10
10. Pack. insert	n° 1	n° 1
Number of tests	192	480
Code	RUBG.CE.192	RUBG.CE.480

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.
3. Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
4. Sera and plasma can be stored at +2°..8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
5. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.
6. Samples whose anti-Rubella IgG antibody concentration is expected to be higher than 250 IU/ml should be diluted before use, either 1:10 or 1:100 in the Calibrator 0 IU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 ul of each specimen with 450 ul of Cal 0 (1:10). Then 50 ul of the 1:10 dilution are diluted with 450 ul of the Cal 0 (1:100). Mix tubes thoroughly on vortex and then proceed toward the dilution step reported in section M.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing.

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Wash buffer concentrate:

The 20x concentrated solution has to be diluted with EIA grade water up to 1200 ml and mixed gently end-over-end before use.

As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. In the preparation avoid foaming as the presence of bubbles could give origin to a bad washing efficiency.

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P362 + P363 – Take off contaminated clothing and wash it before reuse.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory

tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the sections "Validation of Test" and "Assay Performances". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. Incubation times have a tolerance of ±5%.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. When using automatic devices, in case the vial holder of the instrument does not fit with the vials supplied in the kit, transfer the solution into appropriate containers and label them with the same label peeled out from the original vial. This operation is important in order to avoid mismatching contents of vials, when transferring them. When the test is over, return the secondary labeled containers to 2..8°C, firmly capped.
8. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the lyophilised Control Serum as reported in the proper section.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.

8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.
12. Check that all the other equipment is available and ready to use.
13. In case of problems, do not proceed further with the test and advise the supervisor.

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The kit may be used for quantitative and qualitative determinations as well

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5. Wash the microplate with an automatic washer by delivering and aspirating 350 µl/well of diluted washing solution as reported previously (section I.3).

6. Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

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10. Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
11. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

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Proceed as described in section M1.

Manual assay:

1. Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
2. Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
3. Dispense 100 µl of Calibrator 1 (0 IU/ml) and 100 µl Calibrator 2 (10 IU/ml) in duplicate, and 100 µl Calibrator 6 (250 IU/ml) in single. Then dispense 100 µl of diluted samples in each properly identified well.
4. Incubate the microplate for **60 min at +37°C**.

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5. Wash the microplate with an automatic washer by delivering and aspirating 350 µl/well of diluted washing solution as reported previously (section I.3).
6. Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

General Important notes:

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.
- The Control Serum (CS) does not affect the test results calculation. The Control Serum may be used only when a laboratory internal quality control a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

Method	Operations
Calibrators & Control	100 µl
Samples diluted 1:101	100 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H2O2	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm

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Legenda: BLK = Blank CAL = Calibrator CS = Control Serum S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified. Control that the following data are matched:

Check	Requirements
Blank well	< 0.050 OD450nm value
CAL 1 0 IU/ml	< 0.150 mean OD450nm value after blanking coefficient of variation < 30%
CAL 2 10 IU/ml	OD450nm > OD450nm CAL1 + 0.100
CAL 6 250 IU/ml	OD450nm > 1.000
Control Serum	20 IU/ml ±10%

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.050 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
CAL 1 0 IU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.

CAL 2 10 IU/ml OD450nm < OD450nm CAL1 + 0.100	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 6 250 IU/ml < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.
Control Serum Different from expected value	First verify that: 1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (ex.: dispensation of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the standard has occurred. 5. the Control Serum has been dissolved with the right volume reported on the label. If a mistake has been pointed out, the assay has to be repeated after eliminating the reason of this error. If no mistake has been found, proceed as follows: a) a value up to +/-20% is obtained: the overall Precision of the laboratory might not enable the test to match the expected value +/-10%. Report the problem to the Supervisor for acceptance or refusal of this result. b) a value higher than +/-20% is obtained: in this case the test is invalid and the DiaPro's customer service has to be called.

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
P. RESULTS

P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm (4-parameters interpolation is suggested). Then on the calibration curve calculate the concentration of anti Rubella Virus IgG antibody in samples. An example of Calibration curve is reported in the next page.

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P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 10 IU/ml and then check that the assay is valid.

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Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 10 WHO IU/ml are considered negative for anti Rubella Virus IgG antibody by most of the international medical literature.

Samples with a concentration higher than 10 WHO IU/ml are considered positive for anti Rubella Virus IgG antibody.

This titer is considered the lowest concentration of IgG to provide an effective immunological protection against a second infection of Rubella Virus by NCCLS, USA.

Particular attention in the interpretation of results has to be used in the follow-up of pregnancy for an infection of Rubella Virus due to the risk of severe neonatal malformations.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. In the follow-up of pregnancy for Rubella Virus infection a positive result (presence of IgG antibody > 10 IU/ml) should be confirmed to ruled out the risk of a false positive result and a false definition of protection.

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what suggested in NCCLS's approved guideline for Rubella IgG testing (I/LA6-A)

1. Limit of detection

The limit of detection of the assay has been calculated by means of the 1st W.H.O international standard for anti-Rubella immunoglobulin code RUBI-1-94. The limit of detection has been calculated as mean OD450nm Calibrator 0 IU/ml + 5 SD. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined in the assa for three lots.

WHO IU/ml	RUBG.CE Lot. 0303	RUBG.CE Lot. 0503	RUBG.CE Lot. 0603
50	1.292	1.301	1.354
20	0.701	0.742	0.724
10	0.402	0.451	0.425
5	0.211	0.241	0.231
Std 0	0.024	0.032	0.038

The assay shows a limit of detection far better than 10 IU/ml.

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in an external study of performance evaluation (University Hospital, Microbiology Department, Salamanca, Spain) on panels of samples classified positive by a kit US FDA approved. Positive samples from different stage of Rubella Virus infection were tested. The value, obtained from the analysis of more than 300 specimens, has been > 98%.

3. Diagnostic specificity:

The diagnostic specificity has been determined in the same centre on panels of negative samples from not infected individuals, classified negative with a kit US FDA approved. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage. No interference was observed. Potentially interfering samples derived from patients with different pathologies (mostly ANA, AMA and RF positive) and from pregnant women were tested. No crossreaction was observed. An overall value > 98% of specificity was found when examined on more than 100 specimens.

4. Precision:

It has been calculated on three samples, a negative, a low positive and a positive, examined in 16 replicates in three separate runs for three lots. Results are reported as follows:

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5. Accuracy

The assay accuracy has been checked by the dilution and recovery tests. Any "hook effect", underestimation likely to happen at high doses of analyte, was ruled out up to 1.000 IU/ml.

S. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte. Frozen samples containing fibrin particles or aggregates after thawing may generate some false results.

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NOTIFIED BODY	AEMPS – n° 0318
(EC) CERTIFICATE(S)	<ul style="list-style-type: none">• FULL QUALITY ASSURANCE SYSTEM N° 2004 05 0442 CT (in accordance with Annex IV – except Section IV) of the Directive 98/79/EC), RELEASED BY EC NOTIFIED BODY N° 0318• UNI EN ISO 13485 N° 2013 11 0039 EN, RELEASED BY EC NOTIFIED BODY N° 0318

PLACE & DATE OF FIRST ISSUE	MILANO – MAY 2004
PLACE & DATE OF CURRENT EMISSION	SESTO SAN GIOVANNI (MI) – MAY 2018
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

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PLACE & DATE OF FIRST ISSUE	SESTO SAN GIOVANNI (MI) – MAY 2013
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Parvovirus B19 IgG

**Enzyme ImmunoAssay (ELISA) for the
quantitative/qualitative determination of
IgG antibodies to Parvovirus B19
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
Milano - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

Parvovirus B19 IgG

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgG antibodies to Parvovirus B19 in human plasma and sera.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The **B19 virus**, generally referred to as **parvovirus B19** was the first (and until 2005 the only) known human virus in the family of parvoviruses, genus erythrovirus. Parvovirus B19 is a non-enveloped, icosahedral virus that contains a single-stranded linear DNA genome. It is classified as erythrovirus because of its capability to invade red blood cell precursors in the bone marrow. Three genotypes (with subtypes) have been recognized. The viral capsid is composed of two structural proteins, namely VP1 (83kD) and VP2 (53 kD). Infection by Parvovirus B19 spreads through respiratory secretions but also through blood or blood products. The infection causes a mild illness characterized by an erythematous maculopapular facial rash called fifth disease or erythema infectiosum. It is typical in childhood and is also seen in adults. A person usually gets sick within 4 to 14 days after getting infected with parvovirus B19 but about 20% of children and adults who get infected with this virus will not have any symptoms. Infection during pregnancy presents the risk of transmission to the fetus that may result in hydrops fetalis. People with weakened immune systems caused by leukemia, cancer, organ transplants, or HIV infection are at risk for serious complications from fifth disease. It can cause chronic anemia that requires medical treatment. Therefore the detection of Parvovirus B19-specific antibodies becomes very important.

C. PRINCIPLE OF THE TEST

Microplates are coated with Parvovirus B19 antigen.

The solid phase is first treated with the diluted sample and IgG to Parvovirus B19 are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti Parvovirus IgG are detected by the addition of polyclonal specific anti hlgG antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti Parvovirus IgG antibodies present in the sample. A Calibration Curve, calibrated against the 2nd W.H.O international standard for Anti-Parvovirus B19 code 01/602, makes possible a quantitative determination of the IgG antibody in the patient.

D. COMPONENTS

The kit contains reagents to perform 96 tests.

1. Microplate: MICROPLATE

12 strips x 8 microwells coated with Parvovirus B19 antigens. Plates are sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Calibration Curve: CAL N° ...

Ready to use and color coded standard curve derived from human plasma positive for Parvovirus B19 IgG and titrated on WHO standard ranging:

4ml/vial CAL1 = 0 WHO IU/ml

4ml/vial CAL2 = 3 WHO IU/ml

2ml/vial CAL3 = 6 WHO IU/ml

2ml/vial CAL4 = 12 WHO IU/ml

2ml/vial CAL5 = 20 WHO IU/ml

4ml/vial CAL6 = 40 WHO IU/ml.

Standards are calibrated against the 2nd W.H.O international standard for Anti-Parvovirus B19 code 01/602.

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3. Control Serum: CONTROL ...ml

1 vial. Lyophilized.

It contains bovine serum proteins, human plasma positive to Parvovirus B19 calibrated at 12 WHO IU/ml ± 10%, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

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4. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

5. Enzyme conjugate : CONJ

2x8ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300, 0.02 mg/ml gentamicine sulphate as preservatives and 0.01% red alimentary dye.

6. Chromogen/Substrate: SUBS TMB

1x16ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H₂O₂) and 4% dimethylsulphoxide.

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8. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. The reagent is blue colour coded.

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- Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. Bar code labeling and electronic reading is strongly recommended.
- Haemolysed ("red") and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.
- Samples whose anti-Parvovirus B19 IgG antibody concentration is expected to be higher than 40 IU/ml should be diluted before use, either 1:10 in the Specimen Diluent. Dilutions have to be done in clean disposable tubes by diluting 50 µl of each specimen with 450 µl of Specimen Diluent (1:10). Mix tubes thoroughly on vortex and then proceed toward the dilution step reported in section M.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call Dia.Pro's customer service. Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Note: The control after dissolution is not stable. Store frozen in aliquots at -20°C. Do not thaw and freeze the aliquot again.

Wash buffer concentrate:

The 20x concentrated solution has to be diluted with EIA grade water up to 1200 ml and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. In the preparation avoid foaming as the presence of bubbles could give origin to a bad washing efficiency.

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Enzyme conjugate:

Ready to use. Mix well on vortex before use.
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.
If this component has to be transferred use only plastic, possibly sterile disposable containers.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.
Do not expose to strong illumination, oxidizing agents and metallic surfaces.
If this component has to be transferred use only plastic, possible sterile disposable container

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Sulphuric Acid:

Ready to use. Mix well on vortex before use.
Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.
H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.
P332 + P313 – If skin irritation occurs: Get medical advice/attention.
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313 – If eye irritation persists: Get medical advice/attention.
P362 + P363 – Take off contaminated clothing and wash it before reuse.

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- Micropipettes** have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
- The **ELISA incubator** has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests and the right temperature of +37°C is assured to the microplate.
- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution.

The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).

- 5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.
An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of $\pm 5\%$.
5. The **ELISA microplate reader** has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an **ELISA automated work station**, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the lyophilised Control Serum as reported in the proper section.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturer instructions. Set the right number of washing cycles as reported in the specific section.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.

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M1. QUANTITATIVE DETERMINATION:

- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibration Set as calibrators are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of microwells in the microwell holder. Leave the A1 and B1 empty for the operation of blanking.
- Dispense 100 µl of Calibrators in duplicate. Then dispense 100 µl of diluted samples in each properly identified well. The Control Serum doesn't have to be used in every single analysis; it may be used whenever an internal quality control is required by the management to check the overall performances of the laboratory itself. In case, dispense 100 µl of the Control Serum, prepared according to instructions, in duplicate into a proper well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer by delivering and aspirating 350 µl/well of diluted washing solution as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except A1+B1 blanking wells, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1 and B1.

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- Incubate the microplate for **60 min at +37°C**.
- Wash microwells as in step 5.
- Pipette 100 µl Chromogen/Substrate mixture into each well, the blank wells A1 and B1 included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.

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- Pipette 100 µl Sulphuric Acid to stop the enzymatic reaction into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction), blanking the instrument on A1 or B1 or both (mandatory).

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Calibrator 1 (0IU/ml) and Calibrator 5 (20IU/ml) as they are ready to use. The Control Serum doesn't have to be used in every single analysis; it may be used whenever an internal quality control is required by the management to check the overall performances of the laboratory itself. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 100 µl of Calibrator 1 (0 IU/ml) and 100 µl Calibrator 5 (20 IU/ml) in duplicate, and 100 µl Control Serum, prepared according to instructions, in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer by delivering and aspirating 350 µl/well of diluted washing solution as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **60 min at +37°C**.
 - Wash microwells as in step 5.
 - Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 20 minutes**.
- Important note:** Do not expose to strong direct illumination. High background might be generated.

- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.

- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction), blanking the instrument on A1 or B1 or both (mandatory).

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N. ASSAY SCHEME

Method	Operations
Calibrators & Control(*)	100 µl
Samples diluted 1:101	100 µl
1 st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 µl
2 nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H ₂ O ₂	100 µl
3 rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

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O. ASSAY QUALITY CONTROL

A validation check is carried out on the Calibrators any time the kit is used in order to verify whether the performances of the assay are as expected and required by the IVDD directive 98/79/EC.

The Control Serum is used only when required by the management for an internal verification of the performances of the laboratory itself.

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Problem	Check
Blank well > 0.050 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
CAL 1 0 IU/ml > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive calibrator instead of the negative one); 4. that no contamination of the negative calibrator or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.

CAL 2 3 IU/ml OD450nm < OD450nm CAL1 + 0.100	<ol style="list-style-type: none"> 1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
CAL 5 20 IU/ml < 0.750 OD450nm	<ol style="list-style-type: none"> 1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.
CAL 6 40 IU/ml < 1.000 OD450nm	<ol style="list-style-type: none"> 1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong calibrator instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Check	Requirements
Control Serum	Mean OD450nm CAL4 ± 20%

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Problem	Check
Control Serum Different from expected value	First verify that: <ol style="list-style-type: none"> 1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (e.g.: dispensation of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the control has occurred. 5. the Control Serum has been dissolved with the right volume reported on the label. If a mistake has been pointed out, the assay has to be repeated after eliminating the reason of this error.


Anyway, if all other parameters (Blank, CAL1, CAL2, CAL5, CAL6), match the established requirements, the test may be considered valid.

P. RESULTS

P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm (4-parameters interpolation is suggested). Then on the calibration curve calculate the concentration of anti Parvovirus B19 IgG antibody in samples.

Example of Calibration Curve :

Important Note:
Do not use the calibration curve above to make calculations.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for Calibrator 1(0 IU/ml), and for Calibrator 5(20IU/ml) and then check that the assay is valid. In this case the results are calculated by means of a cut-off value determined with the following formula:

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Q. INTERPRETATION OF RESULTS

Q.1 Quantitative method

Samples with a concentration lower than 3 WHO IU/ml are considered negative for anti Parvovirus B19 IgG antibody by most of the international medical literature. Samples with a concentration between 3 and 5 WHO IU/ml are considered equivocal for anti Parvovirus B19 IgG antibody. Samples with a concentration higher than 5 WHO IU/ml are considered positive for anti Parvovirus B19 IgG antibody. This titer is considered the lowest concentration of IgG to provide an effective immunological protection.

Q.2 Qualitative method

Results are interpreted as ratio between the sample OD450nm and the cut-off value or S/Co.

Results are interpreted according to the following table:

S/Co	Interpretation
< 0.8	Negative
0.8 – 1.2	Equivocal
> 1.2	Positive

Any patient showing an equivocal result should be retested on a second sample taken 1-2 weeks after the initial sample.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. Diagnosis has to be done and released to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

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Cut-Off = 1.517 / 5 = 0.303

Sample 1: 0.028 OD450nm
 Sample 2: 1.890 OD450nm
 Sample 1 S/Co < 0.9 negative
 Sample 2 S/Co > 1.0 positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what suggested in NCCLS's approved guideline C24-A2.

1. Limit of detection

The limit of detection of the assay has been calculated by means of the 2nd W.H.O international standard for Anti-Parvovirus B19 code 01/602. The limit of detection has been calculated as mean OD450nm Calibrator 0 IU/ml + 5 SD. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined in the assay for two lots.

WHO IU/ml	PARVOG.CE Lot P1	PARVOG.CE Lot P2
40	3.041	3.110
20	1.686	1.570
12	0.954	0.925
6	0.473	0.549
3	0.266	0.233
1.5	0.112	0.087
Std 0	0.030	0.055

The assay shows a limit of detection far better than 3 IU/ml.

2. Diagnostic sensitivity & Specificity :

The Diagnostic **Sensitivity** was calculated on a panel of 50 samples classified positive for the IgG anti parvovirus B19 by a reference kit CE marked.

A value of $\geq 98\%$ was observed when referring to the reference device.

The Diagnostic **Specificity** was calculated on a panel of more than 100 samples classified negative with the reference device. A value $\geq 98\%$ was observed.

These findings are summarized in the following table.

Sensitivity	$\geq 98\%$
Specificity	$\geq 98\%$

4. Precision:

It has been calculated on three samples, a negative, a low positive and a positive, examined in 16 replicates in three separate runs for two lots. Results are reported as follows:

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5. Accuracy

The assay accuracy has been checked by the dilution and recovery tests. Any "hook effect", underestimation likely to happen at high doses of analyte, was ruled out up to 77 IU/ml.

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	Parvovirus B19 IgM CODE: PARVOM.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	SESTO SAN GIOVANNI (MI) – MAY 2013
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Parvovirus B19 IgM

**Enzyme ImmunoAssay (ELISA) for the
qualitative determination of
IgM antibodies to Parvovirus B19
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

Parvovirus B19 IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the qualitative determination of IgM antibodies to Parvovirus B19 in human plasma and sera.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The **B19 virus**, generally referred to as **parvovirus B19** was the first (and until 2005 the only) known human virus in the family of parvoviruses, genus erythrovirus. Parvovirus B19 is a non-enveloped, icosahedral virus that contains a single-stranded linear DNA genome. It is classified as erythrovirus because of its capability to invade red blood cell precursors in the bone marrow. Three genotypes (with subtypes) have been recognized. The viral capsid is composed of two structural proteins, namely VP1 (83kD) and VP2 (53 kD). Infection by Parvovirus B19 spreads through respiratory secretions but also through blood or blood products. The infection causes a mild illness characterized by an erythematous maculopapular facial rash called fifth disease or erythema infectiosum. It is typical in childhood and is also seen in adults. A person usually gets sick within 4 to 14 days after getting infected with parvovirus B19 but about 20% of children and adults who get infected with this virus will not have any symptoms. Infection during pregnancy presents the risk of transmission to the fetus that may result in hydrops fetalis. In particular the presence of IgM antibodies is reported to be correlated to the acute phase of illness, while IgG antibodies become present at different titers shortly after primary infections and last in blood for many years.

People with weakened immune systems caused by leukemia, cancer, organ transplants, or HIV infection are at risk for serious complications from fifth disease. It can cause chronic anemia that requires medical treatment.

Therefore the detection of Parvovirus B19-specific antibodies becomes very important.

C. PRINCIPLE OF THE TEST

Microplates are coated with Parvovirus B19 antigens.

The solid phase is first treated with the diluted sample and IgM to Parvovirus B19 are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti Parvovirus IgM are detected by the addition of polyclonal specific anti hIgM antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti Parvovirus IgM antibodies present in the sample. The presence of IgM in the sample may therefore be determined by means of a cut-off value able to discriminate between negative and positive samples.

Neutralization of IgG anti-Parvovirus, carried out directly in the well, is performed in the assay in order to block interferences due to this class of antibodies in the determination of IgM.

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1. Microplate: MICROPLATE

12 strips x 8 microwells coated with Parvovirus B19 antigens.

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2. Negative Control: CONTROL -

1x4.0 ml/vial. Ready to use. It contains, human plasma negative to Parvovirus B19, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives.

The Negative Control is pale yellow color coded.

3. Positive Control: CONTROL +

1x4.0 ml/vial. Ready to use. It contains human plasma positive to Parvovirus B19, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives.

The Positive Control is green yellow color coded.

4. Calibrator: CAL ...

n° 1 vial. Lyophilized reagent to be dissolved with EIA grade water as reported in the label. It contains bovine serum proteins, human plasma positive to Parvovirus, 0.2 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

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5. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.

6. Enzyme conjugate : CONJ

1x16ml/vial. Ready to use and red colour coded. It contains Horseradish peroxidase conjugated polyclonal antibodies to human IgM, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

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9. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.5% NP40, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

10. Neutralizing Reagent: SOLN NEUT

1x8ml/vial. It contains goat anti hIgG, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives.

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1. Calibrated Micropipettes (1000 ul, 100 ul and 10 ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-0.5°C tolerance)..
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 3 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are

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4. Sera and plasma can be stored at +2..8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of manufacturing. In this case call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminium pouch, in presence of desiccant supplied, firmly zipped and stored at +2..8°C. When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.

Negative Control

Ready to use components. Mix carefully on vortex before use.

Positive Control

Ready to use components. Mix carefully on vortex before use.

Calibrator

Add the volume of ELISA grade water, reported on the label, to the lyophilized powder; let fully dissolve and then gently mix on vortex.

Note: The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.

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Enzyme conjugate:

Ready to use. Mix well on vortex before use.
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.
If this component has to be transferred use only plastic, possibly sterile disposable containers.

Chromogen/Substrate:

Ready to use. Mix well on vortex before use.
Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.
Do not expose to strong illumination, oxidizing agents and metallic surfaces.
If this component has to be transferred use only plastic, possible sterile disposable container

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Sulphuric Acid:

Ready to use. Mix well on vortex before use.
Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

Legenda:

Warning H statements:

H315 – Causes skin irritation.
H319 – Causes serious eye irritation.

Precautionary P statements:

P280 – Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.
P332 + P313 – If skin irritation occurs: Get medical advice/attention.
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313 – If eye irritation persists: Get medical advice/attention.
P362 + P363 – Take off contaminated clothing and wash it before reuse.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to

set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal quality Control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates.
3. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
4. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
5. Dissolve the content of the Calibrator as reported.
6. Dilute all the content of the 20x concentrated Wash Solution as described above.
7. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
8. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
9. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
10. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
11. Check that the micropipettes are set to the required volume.

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the Controls/Calibrator as they are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Not dispense Neutralizing Reagent in A1 used for blanking operations and in the wells used for the Controls and the Calibrator.
- Dispense 50 µl Neutralizing Reagent in all the samples wells.
- Dispense 100 µl of Negative Control in duplicate, 100 µl of Positive Control in single, 100 µl of Calibrator in duplicate and 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

General Important notes:

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before

reading at 450nm. Finger prints could generate false positive results on reading.

- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

N. ASSAY SCHEME

Method	Operations
Controls & Calibrator (*)	100 µl
Neutralizing Reagent (only for samples)	50 µl
Samples diluted 1:101	100 µl
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

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An example of dispensation scheme is reported in the table below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S3											
B	NC	S4											
C	NC	S5											
D	CAL(*)	S6											
E	CAL(*)	S7											
F	PC	S8											
G	S1	S9											
H	S2	S10											

Legenda: BLK = Blank NC = Negative Control
 CAL = Calibrator - Not mandatory PC = Positive Control
 S = Sample

O. INTERNAL QUALITY CONTROL

A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as expected and required by the IVDD directive 98/79/EC. Control that the following data are matched:

Check	Requirements
Blank well	< 0.100 OD450nm value
Negative Control	< 0.150 mean OD450nm value after blanking coefficient of variation < 30%
Positive Control	OD450nm > 1.000

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not got contaminated during the assay
Negative Control > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of a positive control instead of the negative one); 4. that no contamination of the negative control or of their wells has occurred due spills of positive samples or the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong control) ; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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**** Note:**

If Calibrator has used, verify the following data:

Check	Requirements
Calibrator	S/Co > 1.0

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Problem	Check
Calibrator S/Co < 1.0	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (dispensation of a wrong control instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

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P. RESULTS

If the test turns out to be valid, results are calculated from the mean OD450nm value of the Negative Control (NC) by means of a cut-off value (Co) determined with the following formula:

$$\text{Cut-Off} = \text{NC} + 0.250$$

Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as a ratio of the sample OD450nm value (S) and the cut-off value (Co), or S/Co, according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not developed IgM antibodies to Parvovirus.

Any patient showing an equivocal result should be retested on a second sample taken 1-2 weeks after the initial sample.

A positive result is indicative of an ongoing Parvovirus infection and therefore the patient should be treated accordingly.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. Diagnosis has to be done and released to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

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Positive Control: 1.500 OD450nm
Higher than 0.1000 – Accepted

$$\text{Cut-Off} = 0.100 + 0.250 = 0.350$$

Calibrator: 0.500 – 0.540 OD450nm
Mean value: 0.520 OD450nm
S/Co higher than 1.0 – Accepted

Sample 1: 0.080 OD450nm
Sample 2: 1.800 OD450nm
Sample 1 S/Co < 1.0 = negative
Sample 2 S/Co > 1.2 = positive

R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in accordance to what suggested in NCCLS's approved guideline C24-A2.

1. Limit of detection

No international standard for Parvovirus B19 IgM Antibody detection has been defined so far by the European Community. In its absence, an Internal Gold Standard (or IGS), has been defined in order to provide the device with a constant and excellent sensitivity.

I.G.S. Dilution	PARVOM.CE Lot P1	PARVOM.CE Lot P2
1x	1.418	1.219
2x	0.804	0.665
4x	0.407	0.383
8x	0.225	0.212
Negative Control	0.065	0.070

Positive Control (N = 16)				
Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	3.350	2.937	3.170	3.152
Std.Deviation	0.174	0.199	0.197	0.190
CV %	5.2	7.0	7.1	6.4

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2. Diagnostic Sensitivity and Specificity:

The Diagnostic **Sensitivity** was calculated on a panel of 50 samples classified positive for the IgM anti Parvovirus B19 by a reference kit CE marked.

A value of $\geq 98\%$ was observed when referring to the reference device.

The Diagnostic **Specificity** was calculated on a panel of more than 100 samples classified negative with the reference device.

A value $\geq 98\%$ was observed.

These findings are summarized in the following table.

Sensitivity	$\geq 98\%$
Specificity	$\geq 98\%$

4. Precision:

It has been calculated on three samples, a negative, a low positive and a positive, examined in 16 replicates in three separate runs for two lots. Results are reported as follows:

PARVOM.CE: lot P1

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.149	0.136	0.136	0.140
Std.Deviation	0.013	0.018	0.018	0.016
CV %	8.5	13.0	13.0	11.5

Low Positive sample (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.996	0.977	0.950	0.974
Std.Deviation	0.032	0.057	0.056	0.048
CV %	3.2	5.8	5.9	5.0

Positive Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	3.032	2.671	3.362	3.022
Std.Deviation	0.221	0.221	0.184	0.209
CV %	7.3	8.3	5.5	7.0

PARVOM.CE: lot P2

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.101	0.091	0.094	0.095
Std.Deviation	0.013	0.010	0.013	0.012
CV %	12.5	11.2	13.7	12.5

Low Positive sample (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.227	1.261	0.970	1.153
Std.Deviation	0.085	0.085	0.090	0.087
CV %	6.9	6.7	9.2	7.6

5. Accuracy

The assay accuracy has been checked by the dilution and recovery tests.

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Produced by
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

HBcAb

**Competitive Enzyme Immunoassay for
the determination of antibodies
to Hepatitis B core Antigen
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

HBcAb

A. INTENDED USE

Competitive Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis B core Antigen in human plasma and sera.

The kit is intended for the screening of blood units and the follow-up of HBV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) defines Hepatitis B as follows:

"Hepatitis B is one of the major diseases of mankind and is a serious global public health problem. Hepatitis means inflammation of the liver, and the most common cause is infection with one of 5 viruses, called hepatitis A,B,C,D, and E. All of these viruses can cause an acute disease with symptoms lasting several weeks including yellowing of the skin and eyes (jaundice); dark urine; extreme fatigue; nausea; vomiting and abdominal pain. It can take several months to a year to feel fit again. Hepatitis B virus can cause chronic infection in which the patient never gets rid of the virus and many years later develops cirrhosis of the liver or liver cancer.

HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. The main ways of getting infected with HBV are: (a) perinatal (from mother to baby at the birth); (b) child-to-child transmission; (c) unsafe injections and transfusions; (d) sexual contact.

Worldwide, most infections occur from infected mother to child, from child to child contact in household settings, and from reuse of un-sterilized needles and syringes. In many developing countries, almost all children become infected with the virus. In many industrialized countries (e.g. Western Europe and North America), the pattern of transmission is different. In these countries, mother-to-infant and child-to-child transmission accounted for up to one third of chronic infections before childhood hepatitis B vaccination programmes were implemented. However, the majority of infections in these countries are acquired during young adulthood by sexual activity, and injecting drug use. In addition, hepatitis B virus is the major infectious occupational hazard of health workers, and most health care workers have received hepatitis B vaccine.

Hepatitis B virus is not spread by contaminated food or water, and cannot be spread casually in the workplace. High rates of chronic HBV infection are also found in the southern parts of Eastern and Central Europe. In the Middle East and Indian sub-continent, about 5% are chronically infected. Infection is less common in Western Europe and North America, where less than 1% are chronically infected.

Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 30% to 50% of children infected between 1 to 4 years of age develop chronic

infection. The risk of death from HBV-related liver cancer or cirrhosis is approximately 25% for persons who become chronically infected during childhood.

Chronic hepatitis B in some patients is treated with drugs called *interferon or lamivudine*, which can help some patients. Patients with cirrhosis are sometimes given liver transplants, with varying success. It is preferable to prevent this disease with vaccine than to try and cure it.

Hepatitis B vaccine has an outstanding record of safety and effectiveness. Since 1982, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 95% effective in preventing children and adults from developing chronic infection if they have not yet been infected. In many countries where 8% to 15% of children used to become chronically infected with HBV, the rate of chronic infection has been reduced to less than 1% in immunized groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programmes."

Hepatitis B core Antigen (or HBcAg) is the major component of the core particles of HBV.

HBcAg is composed of a single polypeptide of about 17 kD that is released upon disaggregating the core particles; the antigen contains at least one immunological determinant.

Upon primary infection, anti HBcAg antibodies are one of the first markers of HBV hepatitis appearing in the serum of the patient, slightly later than HBsAg, the viral surface antigen.

Anti HBcAg antibodies are produced usually at high titers and their presence is detectable even years after infection. Isolated HBcAb, in absence of other HBV markers, have been observed in infected blood units, suggesting the use of this test for screening HBV, in addition of HBsAg.

The determination of HBcAb has become important for the classification of the viral agent, together with the detection of the other markers of HBV infection, in sera and plasma.

C. PRINCIPLE OF THE TEST

The assay is based on the principle of competition where the antibodies in the sample compete with a monoclonal antibody for a fixed amount of antigen on the solid phase.

A purified recombinant HBcAg is coated to the microwells.

The patient's serum/plasma is added to the microwell together with an additive able to block interferences present in the sample.

In the second incubation after washing, a monoclonal antibody, conjugated with Horseradish Peroxidase (HRP) and specific for HBcAg is added and binds to the free rec-HBcAg coated on the plastic.

After incubation, microwells are washed to remove any unbound conjugate and then the chromogen/substrate is added. In the presence of peroxidase enzyme the colorless substrate is hydrolyzed to a colored end-product.

The color intensity is inversely proportional to the amount of antibodies to HBcAg present in the sample.

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1. Microplate **MICROPLATE**

8x12 microwell strips coated with recombinant HBcAg and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Negative Control CONTROL -

1x1.0ml/vial. Ready to use. Contains 5% bovine serum albumin, 10 mM phosphate buffer pH 7.4 +/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The negative control is pale yellow color coded.

3. Positive Control CONTROL +

1x1.0ml/vial. Ready to use. Contains 5% bovine serum albumin, anti HBcAg antibodies at a concentration of about 10 PEI U/ml, (calibrated on PEI HBc Reference Material 82), 10 mM phosphate buffer pH 7.4 +/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The positive control is green color coded.

4. Calibrator CAL ...

n° 1 vial. Lyophilised. To be dissolved with EIA grade water as reported in the label. Contains fetal bovine serum, human antibodies to HBcAg at a concentration of 2 PEI U/ml +/-10% (calibrated on PEI HBc Reference Material 82) and 0.1% Kathon GC as preservative.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label .

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6. Enzyme Conjugate CONJ

1x16ml/vial. Ready-to-use solution. Contains 5% bovine serum albumine, 10 mM tris buffer pH 6.8 +/-0.1, Horseradish peroxidase conjugated mouse monoclonal antibody to HBcAg in presence of 0.3 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives. The component is red colour coded .

7. Chromogen/Substrate SUBS TMB

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.6 +/-0.1, 0.03% tetra-methyl-benzidine (TMB), 0.02% hydrogen peroxide (H₂O₂) and 4% dimethylsulphoxide

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8. Specimen Diluent DILSPE

4x3ml/vial. 10 mM tris buffered solution pH 8.0 +/-0.1 containing 0.1% Kathon GC for the pre-treatment of samples and controls in the plate, blocking interference.

Note: Use all the content of one vial before opening a second one. The reagent is sensitive to oxidation.

9. Sulphuric Acid H₂SO₄ 0.3 M

1x15ml/vial. Contains 0.3 M H₂SO₄ solution. Attention: Irritant (Xi R36/38; S2/26/30)

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1. Calibrated Micropipettes (100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
6. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
11. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
13. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
15. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and

disposed according to national directives and laws concerning laboratory wastes.

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- Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
- Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2°.8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

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A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-uses of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in storage. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°.8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Note: *The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C.*

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Note: *Once diluted, the wash solution is stable for 1 week at +2°.8°C.*

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Avoid contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

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Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, and if possible, sterile disposable container.

8. Specimen Diluent

Ready to use solution. Mix gently on vortex before use. Use all the content of one vial before opening a second one. The reagent is sensitive to oxidation.

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Attention: Irritant (Xi R36/38; S2/26/30)

Legenda: R 36/38 = Irritating to eyes and skin.

S 2/26/30 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

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- Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.
- The ELISA incubator has to be set at +37°C (tolerance of ±0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350 µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well characterized negative and positive reference samples, and check to match the values reported below in the sections "Validation of Test" and "Assay Performances". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
- Incubation times have a tolerance of ±5%.
- The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.
- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Dissolve the Calibrator as described above and gently mix.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be performed according to the procedure given below, taking care to maintain the same incubation time for all the samples being tested.

- Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- Leave the A1 well empty for blanking purposes.
- Dispense 50 µl Specimen Diluent into all the control and sample wells.
- Pipette 50 µl of the Negative Control in triplicate, 50 µl of the Calibrator in duplicate and then 50 µl of the Positive Control in single. Then dispense 50 µl of each of the samples.
- Incubate the microplate for **60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Enzyme Conjugate in all the wells, except A1; incubate the microplate for **60 min at +37°C**.

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- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

Important note: Do not expose to strong direct light. as a high background might be generated.

- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with negative control and negative samples will turn from clear to blue (competitive method).
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to stop the enzymatic reaction. Addition of the stop solution will turn the negative control and negative samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1.

Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.

N. ASSAY SCHEME

Specimen Diluent	50 µl
Controls&calibrator and samples	50 µl
1st incubation	60 min
Temperature	+37°C
Wash	n° 4-5
Enzyme Conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash	n° 4-5
TMB/H2O2 mix	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm

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Legenda: BLK = Blank NC = Negative Control
CAL = Calibrator PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A check is performed on the controls/calibrator any time the kit is used in order to verify whether the expected OD450nm or Co/S values have been matched in the analysis.

Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.050 OD450nm value
Negative Control (NC)	> 1.000 OD450nm after blanking coefficient of variation < 20%
Calibrator (about 2 PEI U/ml)	Co/S > 1
Positive Control	< 0.200 OD450nm

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Problem	Check
Blank well > 0.050 OD450nm	that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) < 1.000 OD450nm after blanking coefficient of variation > 20%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator Co/S < 1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control > 0.200 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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P. RESULTS

The results are calculated by means of a cut-off value determined with the following formula:

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Q. INTERPRETATION OF RESULTS

Results are interpreted as ratio between the cut-off value and the sample OD450nm or Co/S.

Results are interpreted according to the following table:

Co/S	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HBV.

Any patient showing an equivocal result should be re-tested on a second sample taken 1-2 weeks after the initial sample. The blood unit should not be transfused.

A positive result is indicative of HBV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

An example of calculation is reported below.

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Positive Control: 0.100 OD450nm
Lower than 0.200 – Accepted

$$\text{Cut-Off} = (2.100 + 0.100) / 5 = 0.440$$

Calibrator: 0.400-0.360 OD450nm
Mean value: 0.380 OD450nm
Co/S > 1 – Accepted

Sample 1: 0.028 OD450nm
Sample 2: 1.890 OD450nm
Sample 1 Co/S > 1.1 positive
Sample 2 Co/S < 0.9 negative

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1. LIMIT OF DETECTION:

The sensitivity of the assay has been calculated by means of the reference preparation for HBcAb supplied by Paul Erlich Institute (PEI HBc Reference Material 82). The assay shows a sensitivity of about 1.25 PEI U/ml.

The table below reports the Co/S values shown by the PEI standard diluted as suggested by the manufacturer to prepare a limiting dilution curve in Fetal Calf Serum (FCS).

PEI U/ml	Lot 1001	Lot 0702	Lot 0702/2	Lot 1202
5	22.6	18.0	19.0	17.7
2.5	8.0	5.5	5.4	5.0
1.25	1.1	1.3	1.0	1.0
0.625	0.4	0.4	0.4	0.4

In addition Accurun 1 – series 3000 – supplied by Boston Biomedica Inc., USA, was tested to determine its Co/S value. Results are reported in the table below:

Accurun 1 – series 3000

Value	Lot 1001	Lot 0702	Lot 1202
Co/S	2.9	2.3	2.2

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The Performance Evaluation of the device was carried out in a trial conducted on more than total 6000 samples.

2.1 Diagnostic Specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. A total of more 5000 unselected donors, including 1st time donors, were examined.

In a first study 2023 samples were tested against a US company as reference. A specificity of 99.5% was found. In a second study 1588 samples were examined against a European company. A specificity of 99.7% was found. In the last study 1565 samples were assayed against the same US company; a value of 99.8% was found.

In addition to the above population, 206 samples from hospitalized patients were tested against the European company. A value of 99.3% specificity was found.

Moreover, diagnostic specificity was assessed by testing 164 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolyzed, lipemic, etc.) against the European company. A value of specificity of 100% was assessed.

Finally, both human plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and human sera have been used to determine the specificity.

No false reactivity due to the method of specimen preparation has been observed.

2.2 Diagnostic Sensitivity

It defined as the probability of the assay of scoring positive in the presence of specific analyte.

373 positive specimens were tested against the European company; a diagnostic sensitivity of 99.7 was found.

3. PRECISION

The mean values obtained from a study conducted on three lots and on two samples of different anti-HBcAg reactivity, examined in 16 replicates in three separate runs is reported below:

BCAB.CE: lot # 1202

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	1.943	1.939	1.924	1.935
Std.Deviation	0.081	0.078	0.103	0.087
CV %	4.2	4.0	5.3	4.5

Calibrator (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.143	0.147	0.148	0.146
Std.Deviation	0.014	0.017	0.018	0.016
CV %	9.8	11.4	12.1	11.1
Co/S	2.8	2.7	2.6	2.7

BCAB.CE: lot # 0702

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.163	2.110	2.106	2.126
Std.Deviation	0.105	0.088	0.139	0.111
CV %	4.9	4.2	6.6	5.2

Calibrator (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.182	0.193	0.195	0.190
Std.Deviation	0.018	0.023	0.019	0.020
CV %	10.0	12.0	9.9	10.6
Co/S	2.5	2.2	2.3	2.3

BCAB.CE: lot # 0702/2

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	2.278	2.098	2.130	2.169
Std.Deviation	0.135	0.126	0.159	0.140
CV %	5.9	6.0	7.5	6.5

Calibrator (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.193	0.190	0.199	0.134
Std.Deviation	0.023	0.023	0.027	0.025
CV %	12.1	12.3	13.5	12.6
Co/S	2.4	2.2	2.2	2.3

The variability shown in the tables did not result in sample misclassification.

S. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte. This test is suitable only for testing single samples and not pooled ones.

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Produced by
Dia.Pro Diagnostic Bioprobes Srl
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

CE
0318

HBc IgM

**“Capture” Enzyme ImmunoAssay (ELISA)
for the quantitative/qualitative
determination of IgM class antibody to
Hepatitis B Virus core Antigen
in human plasma and sera**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

HBc IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IgM class antibodies to Hepatitis B Virus core Antigen in human plasma and sera with the "capture" system.

The kit is intended for the classification of the viral agent and for the follow-up of chronic patients under therapy.
For "in vitro" diagnostic use only.

B. INTRODUCTION

Hepatitis B core Antigen (or HBcAg) is the major component of the core particles of Hepatitis B virus (or HBV).

Particles have a size of 27nm and contain a circular double-stranded DNA molecule, a specific DNA-polymerase and HBcAg. HBcAg is composed of a single polypeptide of about 17 kD that is released upon disaggregation of the core particles ; the antigen contains at least one immunological determinant.

Upon primary infection, anti HBcAg IgM antibodies are one of the first markers of HBV hepatitis appearing in the serum of the patient, together or slightly later than HBsAg, the viral surface antigen.

Anti HBcAg IgM titers, very high during the acute phase, decrease along the illness, as IgG antibodies appear, down to undetectable levels in convalescent patients.

In chronic hepatitis, however, spikes of anti HBcAg IgM synthesis are present, confirming reactivation of HBV in hepatocytes and giving origin to permanent IgM low titers.

The determination of anti HBcAg IgM antibodies has become very important for the fast classification of the virus, of the phase of the illness and for the monitoring of patients under treatment with interferon.

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After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a purified preparation of recombinant HBcAg, labelled with a monoclonal antibody conjugated with peroxidase (HRP).

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In the presence of peroxidase the colourless substrate is hydrolysed to a coloured end-product, whose optical density may be detected and is proportional to the amount of IgM antibodies to HBcAg present in the sample.

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1. Microplate: **MICROPLATE**

8x12 microwell strips coated with purified anti human IgM specific mouse monoclonal antibody, post-coated with bovine serum proteins and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Calibration Curve: **CAL N° ...**

6x2.0 ml/vial. Ready to use and color coded standard curve calibrated on the HBcIgM reference preparation supplied by Paul Ehrlich Institute (HBc-Referenzserum-IgM 84), ranging: CAL1 = 0 U/ml // CAL2 = 5 U/ml // CAL3 = 10 U/ml // CAL4 = 20 U/ml // CAL 5 = 50 U/ml // CAL 6 = 100 U/ml.

It contains chemical inactivated HBcIgM positive human plasma, 100 mM Tris buffer pH 7.4+/-0.1, 0.5% Tween 20, 0.09% sodium azide and 0.045% ProClin 300 as preservatives. The Calibration Curve is coded with blue alimentary dye.

Important Note: Even if plasma has been chemically inactivated, handle this component as potentially infectious.

3. Wash buffer concentrate: **WASHBUF 20X**

1x60ml/bottle. 20x concentrated solution.

Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

4. Enzyme Conjugate (Immunocomplex) : **CONJ**

1x16.0 ml/vial. Ready-to-use solution. Contains an immunocomplex formed by a specific mouse monoclonal antibody, labelled with HRP, and a purified recombinant HBcAg. The reagent is dissolved into a buffer solution 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The component is red colour coded.

5. Specimen Diluent : **DILSPE**

2x60.0 ml/vial. Buffered solution for the dilution of samples; it contains 100 mM Tris buffer pH 7.4+/-0.1, 0.5% Tween 20, 2% Casein, 0.045% ProClin 300 and 0.09% sodium azide as preservatives. The component is blue color coded.

6. Control Serum : **CONTROL ...ml**

1 vial. Lyophilized. Contains fetal bovine serum, human HBcIgM positive human plasma calibrated at 20 ± 10% PEI U/ml. 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Important Notes

1. The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label .

2. Important Note: Even if plasma has been chemically inactivated, handle this component as potentially infectious.

7. Chromogen/Substrate : **SUBS TMB**

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

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1. Calibrated Micropipettes (150ul, 100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

- The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
- Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
- Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
- Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
- Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
- Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
- Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
- Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
- Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
- Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
- The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water
- Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

G. SPECIMEN: PREPARATION AND RECOMMANDATIONS

- Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been

observed in the preparation of the sample with citrate, EDTA and heparin.

- Avoid any addition of preservatives; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results.
- Haemolysed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

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Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call Dia.Pro's customer service. Unused strips have to be placed back into the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200ml and mixed gently end-over-end before use.

During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

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Avoid contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

Specimen Diluent

Ready to use. Mix on vortex before use.

Control Serum

Dissolve the content of the vial with EIA grade water as reported in the label. Mix well on vortex before use. The dissolved control serum is ready to use.

Note: The control after dissolution is not stable. Store frozen in aliquots at -20°C.

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Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong light, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, and if possible, sterile disposable container

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P362 + P363 – Take off contaminated clothing and wash it before reuse.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested). 5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing. An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Blanking is carried out on the well identified in the section "Assay Procedure". The optical

system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB+H₂O₂) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Control Serum as described above and gently mix.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.

In case of problems, do not proceed further with the test and advise the supervisor.

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Two procedures can be carried out with the device according to the request of the clinician.

M.1 Quantitative analysis

1. Place the required number of strips in the plastic holder and carefully identify the wells for standards and samples.

- Dilute samples **1:101** dispensing 1 ml Sample Diluent into a disposable tube and then 10 µl sample; mix on vortex before use. Do not dilute the Calibrators and the dissolved Control Serum as they are ready-to-use.
- Leave the A1+B1 wells empty for blanking purposes.
- Pipette 100 µl of the Calibrators in duplicate, 100 µl dissolved Control Serum in duplicate followed by 100 µl of diluted samples. The Control Serum is used to verify that the whole analytical system works as expected. Check that Calibrators, Control Serum and samples have been correctly added.
- Incubate the microplate **for 60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- In all the wells except A1+B1, pipette 100 µl Enzyme Conjugate. Incubate the microplate **for 60 min at +37°C**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.

- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1+B1 included.

Important note: Do not expose to strong direct light. as a high background might be generated.

- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive samples, the control serum and the positive calibrators, as well, will turn from clear to blue.
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to block the enzymatic reaction.. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

M.2 Qualitative analysis

- Place the required number of strips in the plastic holder and carefully identify the wells for standards and samples.
- Dilute samples **1:101** dispensing 1 ml Sample Diluent into a disposable tube and then 10 µl sample; mix on vortex before use. Do not dilute the Calibrators as they are ready-to-use.
- Leave the A1 well empty for blanking purposes.
- Pipette 100 µl Calibrator 0 U/ml in duplicate, 100 µl Calibrator 10 U/ml in duplicate and 100 µl Calibrator 100 U/ml in single. Then dispense 100 µl diluted samples in proper sample wells. Check that Calibrators and samples have been correctly added.
- Incubate the microplate **for 60 min at +37°C**.

Important note: Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

- When the first incubation is finished, wash the microwells as previously described (section I.3)
- In all the wells except A1, pipette 100 µl Enzyme Conjugate. Incubate the microplate **for 60 min at +37°C**.

Important note: Be careful not to touch the inner surface of the well with the pipette tip and not to immerse the top of it into samples or controls. Contamination might occur.

- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

Important note: Do not expose to strong direct light. as a high background might be generated.

- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive samples, the control serum and the positive calibrators, as well, will turn from clear to blue.
- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 9 to block the enzymatic reaction. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1 or B1 or both.

Important notes:

- Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
- Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
- The Control Serum (CS) does not affect the cut-off calculation and therefore the test results calculation. The Control Serum may be used only when a laboratory internal quality control is required by the management

N. ASSAY SCHEME

The assay protocol can be summarized in the table below:

Calibrators & diluted samples & dissolved Control Serum	100 ul
1st incubation	60 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing steps	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Chromogen/Substrate	100ul
3rd incubation	20 min
Temperature	room
Sulphuric Acid	100 ul
Reading OD	450nm /620-630nm

An example of dispensation scheme in quantitative assays is reported below:

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Calibrator 0 U/ml > 0.150 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study;
coefficient of variation > 30%	2. that the proper washing solution has been used and the washer has been primed with it before use;
	3. that no mistake has been done in the assay procedure (dispensation of positive calibrators instead of Cal 0);
	4. that no contamination of the Cal 0, or of the wells where this was dispensed, has occurred due to positive samples, to spills or to the enzyme conjugate;
	5. that micropipettes have not become contaminated

	with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator 5 U/ml < CAL 0 + 5SD or < CAL 0 + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Calibrator 10 U/ml < CAL 0 + 0.200	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Calibrator 100 U/ml < 1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the calibrator; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Control Serum Different from expected value	First verify that: 1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (ex.: dispensation of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the standard has occurred. 5. the Control Serum has been dissolved with the right volume reported on the label. If a mistake has been pointed out, the assay has to be repeated after eliminating the reason of this error. If no mistake has been found, proceed as follows: a) a value up to +/-20% is obtained: the overall Precision of the laboratory might not enable the test to match the expected value +/-10%. Report the problem to the Supervisor for acceptance or refusal of this result. b) a value higher than +/-20% is obtained: in this case the test is invalid and the DiaPro's customer service has to be called.

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
P. RESULTS

P.1 Quantitative method

If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm/620-630nm (4-parameters interpolation is suggested).

Then on the calibration curve calculate the concentration of anti HBc IgM antibody in samples.

An example of Calibration curve is reported below.

Important Note: Do not use this example to make real calculations on samples.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm/620-630nm values for the Calibrators 0 and 10 U/ml and then check that the assay is valid.

Example of calculation (data obtained proceeding as the the reading step described in the section M, point 12).

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Q. INTERPRETATION OF RESULTS

Q.1 Qualitative results

For qualitative interpretations, the medical literature generally considers positive samples showing a concentration of HBc IgM ≥ 10 PEI U/ml.

Test results are therefore interpreted as a ratio of the sample OD450nm and the OD450nm/620-630nm of the Cal 10 PEI U/ml (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

Q.2 Quantitative results

The calibration curve is used to determine the concentration of IgM antibodies to HBcAg in samples.

Samples with a concentration lower than 5 PEI U/ml are considered negative for HBcIgM.

Samples with a concentration between 5 and 10 PEI U/ml are considered in a gray-zone.

In the follow up of chronic hepatitis, however, values higher of 5 PEI U/ml may be considered positive for HBcIgM, when in presence of other clinical signs.

Samples with a concentration higher than 10 PEI U/ml are considered positive for HBcIgM.

Important general notes:

- When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to produce the calibration curve, calculate sample concentration and generate the correct interpretation of results.
- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
- A positive result is indicative of HBV infection and therefore the patient should be treated accordingly.
- When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
- Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

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1. Limit of detection

The limit of detection of the assay has been calculated by means of :

- the HBcIgM reference preparation supplied by Paul Erlich Institute, Germany (HBc-Referenzserum-IgM 84), on which the Standard Curve has been calibrated.
- Accurun 113 (cat. N° A113-5001) supplied by Boston Biomedica Inc., USA

Results of Quality Control for three lots are given in the following tables:

BCM.CE	Lot #	0103	Lot #	0103/2	Lot #	0303
PEI U/ml	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
100	2.752	8.9	2.883	9.7	2.911	9.1
50	1.917	6.2	1.972	6.7	2.053	6.4
20	0.980	3.2	0.914	3.1	1.095	3.4
10	0.544	1.8	0.513	1.7	0.592	1.8
5	0.310	1.0	0.296	1.0	0.321	1.0
2.5	0.155	0.5	0.149	0.5	0.161	0.5
1.25	0.084	0.3	0.084	0.3	0.093	0.3
negative	0.040		0.035		0.044	

BBI Accurun # 113 lot # 48-9999-0621

BCM.CE	Lot #	0103	Lot #	0103/2	Lot #	0303
BBI 113	OD450nm	S/Co	OD450nm	S/Co	OD450nm	S/Co
1 x	3.336	10.8	3.195	10.4	3.269	10.3
2 x	2.472	8.0	2.385	7.8	2.385	7.5
4 x	1.467	4.7	1.413	4.6	1.429	4.5
8 x	0.865	2.8	0.807	2.6	0.856	2.7
16 x	0.430	1.4	0.427	1.4	0.410	1.3
32 x	0.234	0.8	0.234	0.8	0.248	0.8
64 x	0.129	0.4	0.133	0.4	0.122	0.4
128 x	0.086	0.3	0.082	0.3	0.089	0.3
negative	0.040		0.040		0.052	

Moreover the BBI's panel # PHE 102 was also examined in three lots of product; data are reported below with reference to a European kit (BBI's results).

BBI – Panel code PHE 102

	Lot # 0103	Lot # 0103/2	Lot # 0303	Sorin EIA
Member	S/Co	S/Co	S/Co	S/Co
01	6.7	6.3	6.5	2.0
02	11.3	10.0	10.7	6.1
03	9.5	7.2	8.4	3.0
04	5.8	3.4	4.1	2.1
05	11.3	11.4	11.2	3.1
06	12.1	11.6	11.8	4.1
07	0.1	0.1	0.1	0.2
08	9.2	8.5	8.8	2.3
09	12.2	11.7	11.9	4.2
10	11.7	10.2	10.8	2.8
11	5.9	5.8	5.8	2.1
12	12.7	11.4	11.7	5.2
13	11.6	11.0	11.3	3.6
14	7.0	6.3	6.6	2.3
15	12.4	11.5	11.8	4.5

2. Diagnostic Sensitivity:

It is defined as the probability of the assay of scoring positive in the presence of the specific analyte.

The diagnostic sensitivity has been tested internally and externally in a qualified Clinical Laboratory on panels of samples classified positive by a US FDA approved kit.

Positive samples were collected from different patients and from different HBV pathologies (acute and chronic hepatitis).

An overall value > 98% has been found in the study conducted on a total number of more than 200 samples.

A Seroconversion panel produced by BBI, USA, code # PHM 935A, has also been studied; results are reported below with reference to two commercial kits (BBI's results).

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3. Diagnostic Specificity:

It is defined as the probability of the assay of scoring negative in the absence of the specific analyte.

The diagnostic specificity has been determined internally and externally in a qualified Clinical Laboratory on panels of negative samples from normal individuals and blood donors, classified negative with a US FDA approved kit.

A total number of more than 400 negative specimens were tested. A diagnostic specificity > 98% has been found.

Moreover, the diagnostic specificity was assessed by testing more than 50 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolyzed, lipemic, etc.). No interference was observed in the study.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been

used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

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4. Precision:

It has been calculated on three samples examined in 16 replicate in three different runs, carried out on three different lots. The values found were as follows:

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0318

HBe Ag&Ab

**Enzyme Immunoassay (ELISA) for the
determination of Hepatitis B Virus
"e" Antigen and Antibody
in human plasma and sera.**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

HBe Ag&Ab

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of Hepatitis B Virus "e" Antigen and Antibody in human plasma and sera.

The kit is intended for the follow-up of acute infection and of chronic patients under therapy.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Hepatitis B "e" Antigen or HBeAg is known to be intimately associated with Hepatitis B Virus or HBV replication and the presence of infectious Dane particles in the blood.

Recently, it has been found that HBeAg is a product of proteolytic degradation of Hepatitis B core Antigen or HBcAg, occurring in hepatocytes, whose expression is under the control of the precore region of HBV genome.

If HBeAg is considered a specific marker of infectivity, the presence of anti HBeAg antibodies in blood is recognised to be a clinical sign of recovery from infection to convalescence.

The determination of these two analytes in samples from HBV patients has become important for the classification of the phase of illness and as a prognostic value in the follow up of infected patients.

C. PRINCIPLE OF THE TEST

HBeAg:

HBeAg, if present in the sample, is captured by a specific monoclonal antibody, in the 1st incubation.

In the 2nd incubation, after washing, a tracer, composed of a mix of two specific anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP), is added to the microplate and binds to the captured HBeAg.

The concentration of the bound enzyme on the solid phase is proportional to the amount of HBeAg in the sample and its activity is detected by adding the chromogen/substrate in the 3rd incubation.

The presence of HBeAg in the sample is determined by means of a cut-off value that allows for the semiquantitative detection of the antigen.

HBeAb

Anti HBeAg antibodies, if present in the sample, compete with a recombinant HBeAg preparation for a fixed amount of an anti HBeAg antibody, coated on the microplate wells.

The competitive assay is carried out in two incubations, the first with the sample and rechBeAg, and the second with a tracer, composed of two anti HBeAg monoclonal antibodies, labeled with peroxidase (HRP).

The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti HBeAg antibodies in the sample and its activity is detected by adding the chromogen/substrate in the third incubation.

The concentration of HBeAg specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti HBeAg antibodies.

D. COMPONENTS

The kit contains reagents for total 96 tests.

1. Microplate: MICROPLATE

n° 1 coated microplate

12 strips of 8 breakable wells coated with anti HBeAg specific monoclonal antibody, postcoated with bovine serum proteins and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 2..8°C.

2. Negative Control: CONTROL

1x2.0ml/vial. Ready to use control. It contains bovine serum, 0.09% sodium azide and 0.045% ProClin 300 as preservatives. The negative control is colorless.

3. Antigen Positive Control: CONTROL + Ag

1x1.0ml/vial. Ready to use control. It contains 2% bovine serum albumin, non infectious recombinant HBeAg, 100 mM tris buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.045% ProClin 300 as preservatives.

The positive control is green color coded.

4. Antibody Positive Control: CONTROL + Ab

1x1.0ml/vial. Ready to use control. It contains 2% bovine serum albumin, human anti HbeAg positive plasma at about 10 PEI U/ml, 100 mM tris buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.045% ProClin 300 as preservatives. The label is red colored.

The positive control is yellow color coded.

5. Antigen Calibrator: CALAG ...ml

n° 1 vial. Lyophilized calibrator for HBeAg. To be dissolved with EIA grade water as reported in the label. It contains fetal bovine serum, non infectious recombinant HBeAg at 1 PEI U/ml +/-10%, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

6. Antibody Calibrator: CALAB ...ml

n° 1 vial. Lyophilized calibrator for anti HBeAg antibody. To be dissolved with EIA grade water as reported in the label. It contains fetal bovine serum, positive plasma at 0.25 PEI U/ml +/-10%, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives. The label is red colored.

Important Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

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8. Enzyme conjugate: CONJ

1x16ml/vial. Ready to use conjugate. It contains Horseradish peroxidase conjugated with a mix of monoclonal antibodies to HBeAg, 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The reagent is red color coded.

9. HBe Antigen: Ag-HBe

1x10ml/vial. Ready to use reagent. It contains recombinant HBeAg, fetal bovine serum, buffered solution pH 8.0+/-0.1, 0.045% ProClin 300 and 0.09% sodium azide as preservatives. The reagent is blue color coded.

10. Chromogen/Substrate: SUBS TMB

1x16ml/vial. Ready-to-use component. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

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F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are

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14. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

15. The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results.
4. Haemolysed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

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1. Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in manufacturing. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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7. Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200 ml and mixed gently end-over-end before use.

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Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

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Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

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- Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
- The ELISA incubator has to be set at +37°C (tolerance of $\pm 0.5^\circ\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing. An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
- Incubation times have a tolerance of $\pm 5\%$.
- The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB+H₂O₂) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation

and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.

- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Dissolve the Calibrator as described above and gently mix.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.

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A) HBe Antigen:

- Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- Leave the A1 well empty for blanking purposes.
- Pipette 100 µl of the Negative Control in triplicate, 100 µl of the Antigen Calibrator in duplicate and then 100 µl of the Antigen Positive Control in single.
- Then dispense 100 µl of samples in the proper wells.
- Check for the presence of samples in wells by naked eye (there is a marked colour difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0.100).
- Incubate the microplate for **60 min at +37°C**.

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- When the first incubation is finished, wash the microwells as previously described (section I.3)
- Dispense 100 µl Enzyme Conjugate in all wells, except for A1, used for blanking operations.

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- Check that the reagent has been dispensed properly and then incubate the microplate for **60 min at +37°C**.
- When the second incubation is finished, wash the microwells as previously described (section I.3)
- Pipette 100 µl Chromogen/Substrate into all the wells, A1 included.

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- Incubate the microplate protected from light at **room temperature (18-24°C) for 20 minutes**. Wells dispensed with positive control and positive samples will turn from clear to blue.

- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 11. Addition of the stop solution will turn the positive control and positive samples from blue to yellow.

- Measure the color intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1.

B) HBe Antibody:

- Place the required number of strips in the plastic holder and carefully identify the wells for controls, calibrator and samples.
- Leave the A1 well empty for blanking purposes.
- Pipette 50 µl of the Negative Control in triplicate, 50 µl of the Antibody Calibrator in duplicate and then 50 µl of the Antibody Positive Control in single.
- Then dispense 50 µl of samples in the proper wells.
- Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm (samples show OD values higher than 0.100).
- Dispense then 50 µl of HBe Antigen in all the wells, except for A1.
- Incubate the microplate for **60 min at +37°C**.

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- When the first incubation is finished, wash the microwells as previously described (section I.3)
- Finally proceed as described for the HBeAg assay from point 8 to the last one.

Important notes:

- Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
- Reading should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
- The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

HBe antigen test

Controls and calibrator	100 ul
Samples	100 ul
1st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2 mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm/620-630nm

HBe antibody test

Controls and calibrator	50 ul
Samples	50 ul
Neutralising antigen	50 ul
1st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzymatic conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2 mixture	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm/620-630nm

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HBe Antibody

Check	OD450nm
Blank well	< 0.100 OD450nm
Negative Control (NC)	> 1.000 OD450nm after blanking coefficient of variation < 10%
Antibody Calibrator	OD450nm < NC/1.5
Positive Control (PC)	OD450nm < NC/10

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If they do not, don't proceed any further and perform the following checks:

HBeAg

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co < 2	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control < 1.500 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

HBe antibody

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) < 1.000 OD450nm after blanking coefficient of variation > 10%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (e.g.: dispensation of positive control instead of negative control; no dispensation of the Neutralizing Antigen; no dispensation of the Enzyme Conjugate); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred; 5. that micropipettes have not become contaminated with positive samples; 6. that the washer needles are not blocked or partially obstructed.

Calibrator OD450nm > NC/1.5	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead; no dispensation of the Neutralizing Antigen; no dispensation of the Enzyme Conjugate); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control OD450nm > NC/10	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Important note:

The analysis must be done proceeding as the reading step described in the section M, point 14.

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HBeAb:

$$(NC + PC) / 3 = \text{Cut-Off (Co)}$$

Important note: *When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.*

Q. INTERPRETATION OF RESULTS

Results are interpreted as follows:

HBeAg:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

HBeAb:

Co/S	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

Note:

$$S = OD450nm/620-630nm \text{ of the sample}$$

$$Co = \text{cut-off value}$$

An example of calculation for HBeAg assay is reported below (data obtained proceeding as the the reading step described in the section M, point 14):

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Sample 1 S/Co < 0.9 = negative

Sample 2 S/Co > 1.1 = positive

An example of calculation for HBeAb is reported below (data obtained proceeding as the the reading step described in the section M, point 14):

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Sample 2 Co/S < 0.9 negative

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory director to reduce the risk of judgment errors and misinterpretations.
2. The Identification of the clinical status of a HBV patient (acute, chronic, asymptomatic hepatitis) has to be done on the basis also of the other markers of HBV infection (HBsAg, HBsAb, HBcAb, HBcIgM);
3. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCE CHARACTERISTICS

A) HBeAg

1. Limit of detection

The limit of detection of the assay has been calculated by means of the International Standard for HBeAg, supplied by Paul Erlich Institute (PEI).

The data obtained by examining the limit of detection on three lots is reported in the table below.

HBE.CE Lot ID	PEI U/ml HBeAg
0103	0.25
0103/2	0.25
0303	0.25

In addition the preparation Accurun # 51, produced by Boston Biomedica Inc., USA, has been tested, upon dilution in FCS. Results are reported for three lots of products.

BBI's Accurun 51 (S/Co)

HBE.CE Lot ID	1 x	2 x	4 x	8 x	16x
0103	4.1	1.6	0.9	0.6	0.4
0103/2	4.1	1.7	0.9	0.6	0.4
0303	4.0	1.6	0.9	0.5	0.4

2. Diagnostic Sensitivity:

The diagnostic sensitivity has been tested on panels of samples classified positive by a US FDA approved kit.

Positive samples were collected from different HBV pathologies (acute, chronic) bearing HBeAg reactivity.

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Moreover the Panel of Seroconversion code PHM 935B, produced by BBI, was examined.

Data are reported below and compared with those reported by BBI for two other commercial products.

Sample ID	HBE.CE S/Co	Abbott EIA S/Co	Sorin EIA S/Co
21	5.4	4.5	6.3
22	3.7	4.3	5.4
23	1.9	3.2	3.1
24	1.1	2.4	1.5
25	1.0	2.1	1.2
26	0.6	1.7	0.7
27	0.2	0.8	0.3
28	0.2	0.6	0.2
29	0.2	0.4	0.2
30	0.2	0.3	0.2
31	0.1	0.3	0.2
32	0.1	0.3	0.2

Finally the Performance Panel code PHJ 201, produced by BBI, was tested. Data are reported below and compared with those reported by BBI for an other commercial product.

Member	PEI U/ml	HBE.CE	Sorin EIA
1	3	3.3	7.0
2	6	17.5	21.9
3	26	30.1	37.1
4	31	29.4	23.5
5	1	1.1	2.2
6	2	2.3	6.9
7	35	30.1	24.6
8	38	29.2	31.9
9	4	16.6	10.8
10	-	0.3	0.2

11	1	3.4	3.6
12	< 1	0.2	1.2
13	< 1	0.9	1.4
14	-	0.2	0.2
15	-	0.4	0.1
16	-	0.5	0.1
17	-	0.3	0.2
18	-	0.2	0.2
19	-	0.2	0.1
20	-	0.2	0.1
21	-	0.3	1.0
22	-	0.3	0.1
23	-	0.4	0.1
24	-	0.2	0.2
25	-	0.3	0.2

3. Diagnostic Specificity:

The diagnostic specificity has been determined on panels of negative samples from normal individuals and blood donors, classified negative with a FDA approved kit.

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Samples derived from patients with different viral (HCV and HAV) and non viral pathologies of the liver that may interfere with the test were examined.

No cross reaction were observed.

The Performance Evaluation study conducted in a qualified external reference center on more than 500 samples has provided a value > 98% .

4. Precision

It has been calculated on two samples examined in 16 replicate in three different runs on three lots.

The values found were as follows:

HBE.CE: lot # 0103

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.030	0.027	0.032	0.029
Std.Deviation	0.002	0.002	0.003	0.002
CV %	7.4	8.2	7.9	7.8

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.569	0.559	0.575	0.568
Std.Deviation	0.027	0.029	0.028	0.028
CV %	4.7	5.3	4.9	4.9
S/Co	4.4	4.4	4.4	4.4

HBE.CE: lot # 0103/2

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.033	0.031	0.030	0.032
Std.Deviation	0.003	0.003	0.002	0.003
CV %	7.9	8.5	7.4	8.0

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.565	0.573	0.568	0.569
Std.Deviation	0.026	0.025	0.024	0.025
CV %	4.7	4.3	4.2	4.4
S/Co	4.2	4.4	4.4	4.3

HBE.CE: lot # 0303

Negative Control (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.029	0.034	0.038	0.034
Std.Deviation	0.003	0.003	0.004	0.003
CV %	9.7	9.8	9.2	9.6

PEI 1 U/ml (N = 16)

Mean values	1st run	2nd run	3 rd run	Average value
OD 450nm	0.579	0.573	0.564	0.572
Std.Deviation	0.023	0.028	0.025	0.025
CV %	4.1	4.8	4.5	4.5
S/Co	4.5	4.3	4.1	4.3

B) HBe Antibody

1. Limit of detection

The limit of detection of the assay has been calculated by means of the International Standard for HBeAb, supplied by Paul Erlich Institute (PEI).

The data obtained by examining the limit of detection on three lots is reported in the table below.

HBE.CE Lot ID	PEI U/ml HBeAb
0103	0.25
0103/2	0.25
0303	0.25

In addition the preparation Accurun # 52, produced by Boston Biomedica Inc., USA, has been tested, upon dilution in FCS. Results are reported for three lots of products.

Accurun 52 (Co/S)

HBE.CE Lot ID	1 x	2 x	4 x	8 x	16x
0103	1.0	0.8	0.6	0.4	0.4
0103/2	1.0	0.8	0.6	0.5	0.4
0303	1.0	0.8	0.6	0.4	0.4

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested on panels of samples classified positive for HBeAb by a US FDA approved kit. Positive samples were collected from different HBV pathologies bearing anti HBeAg antibody reactivity.

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Moreover the Panel of Seroconversion code PHM 935B, produced by BBI, was examined.

Data are reported below and compared with those reported by BBI for two other commercial products.

Sample ID	HBE.CE Co/S	Abbott EIA Co/S	Sorin EIA Co/S
21	0.4	0.4	0.5
22	0.4	0.5	0.6
23	0.4	0.6	0.5
24	0.4	0.5	0.6
25	0.4	0.6	0.5
26	0.5	0.6	0.6
27	0.6	0.8	0.7
28	0.7	0.9	0.7
29	0.6	0.9	0.7
30	0.8	1.0	0.9
31	1.0	1.3	1.1
32	1.0	1.2	1.0

Finally the Performance Panel code PHJ 201, produced by BBI, was tested. Data are reported below and compared with those reported by BBI for another commercial product.

Member	PEI U/ml	HBE.CE	Sorin EIA
1	-	0.3	0.5
2	-	0.2	0.5
3	-	0.2	0.5
4	-	0.2	0.5
5	-	0.3	0.6
6	-	0.3	0.6
7	-	0.2	0.4
8	-	0.2	0.4
9	-	0.2	0.5
10	-	1.9	0.6
11	-	0.3	0.5
12	-	0.4	0.9
13	2	4.4	9.1
14	1	3.8	2.9
15	< 1	1.0	1.5
16	> 50	4.3	120.9
17	< 1	1.0	1.0
18	5	5.6	21.8
19	1	2.7	6.4
20	11	5.0	47.3
21	2	1.9	10.0
22	26	28.1	90.7
23	-	0.3	0.5
24	< 1	0.8	1.3
25	50	28.1	167.4

3. Diagnostic specificity:

The clinical specificity has been determined as described before for HBeAg.

The Performance Evaluation study conducted in a qualified external reference center on more than 500 samples has provided a value > 98%.

4. Precision:

It has been calculated on two samples examined in 16 replicate in three different runs on three lots.

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0318

Important note:

The performance data have been obtained proceeding as the reading step described in the section M, point 14.

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HBsAb

**Enzyme Immunoassay for
qualitative/quantitative determination of
antibodies to Hepatitis B surface Antigen
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

REF SAB.CE
96 Tests

HBs Ab

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for both the quantitative and qualitative determination of antibodies to the Surface Antigen of Hepatitis B Virus in human plasma and sera.
For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) defines Hepatitis B Virus infection as follows:

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Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 30% to 50% of children infected between 1 to 4 years of age develop chronic infection. The risk of death from HBV-related liver cancer or

cirrhosis is approximately 25% for persons who become chronically infected during childhood.

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Hepatitis B surface Antigen (HBsAg) is the major structural polypeptide of the envelope of the Hepatitis B Virus (HBV).

This antigen is composed mainly of the type common determinant "a" and the type specific determinants "d" and "y", present only on the specific serotypes.

Upon infection, a strong immunological response develops firstly against the type specific determinants and in a second time against the "a" determinant.

Anti "a" antibodies are however recognised to be most effective in the neutralisation of the virus, protecting the patient from other infections and leading it to convalescence.

The detection of HBsAb has become important for the follow up of patients infected by HBV and the monitoring of recipients upon vaccination with synthetic and natural HBsAg.

C. PRINCIPLE OF THE TEST

Microplates are coated with a preparation of highly purified HBsAg that in the first incubation with sample specifically captures anti HBsAg antibodies to the solid phase.

After washing, captured antibodies are detected by an HBsAg, labelled with peroxidase (HRP), that specifically binds the second available binding site of these antibodies.

The enzyme specifically bound to wells, by acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of HBsAb in the sample and can be detected by an ELISA reader.

The amount of antibodies may be quantitated by means of a standard curve calibrated against the W.H.O reference preparation.

Samples are pre treated in the well with an specimen diluent able to block interference present in vaccinated individuals.

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1. Microplate: MICROPLATE

8x12 microwell strips coated with purified heat-inactivated HBsAg of both subtypes (ad and ay) from human origin and sealed into a bag with desiccant.

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2. Calibration Curve: **CAL N° ...**

5x2.0 ml/vial. Ready to use and colour coded standard curve, derived from HBsAg positive plasma titrated on WHO standard for anti HBsAg (1st reference preparation 1977, lot 17-2-77), ranging: CAL1 = 0 mIU/ml // CAL2 = 10 mIU/ml // CAL3 = 50 mIU/ml // CAL4 = 100 mIU/ml // CAL5 = 250 mIU/ml. Contains human serum proteins, 5% BSA, 10 mM phosphate buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. Standards are blue coloured.

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4. Enzyme conjugate : **CONJ**

1x16.0 ml/vial. Ready-to-use solution and red color coded. It contains inactivated purified HBsAg of both subtypes ad and ay, labelled with HRP, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.3 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

5. Chromogen/Substrate: **SUBS TMB**

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetramethyl-benzidine (TMB) and 0.02% hydrogen peroxide (H₂O₂).
Note: To be stored protected from light as sensitive to strong illumination.

6. Sulphuric Acid: **H2SO4 0.3 M**

1x15ml/vial. Contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

7. Specimen Diluent: **DILSPE**

1x8ml. 10 mM Tris Buffered solution pH 7.4 +/-0.1, suggested to be used in the follow up of vaccination. It contains 0.09% sodium azide as preservatives.

8. Control Serum: **CONTROL ...ml**

1 vial. Lyophilized. Contains fetal bovine serum proteins, human anti HBsAg antibodies calibrated at 50 ± 10% WHO mIU/ml. 0.3 mg/ml gentamicine sulphate and 0.1% Kathon GC as preservatives.

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1. Calibrated Micropipettes (100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), set at +37°C (+/-1°C tolerance)..
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking, strongly recommended) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

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1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

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4. Sera and plasma can be stored at +2°.8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be freeze/thawed more than once as this may generate particles that could affect the test result.

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6. Samples whose anti-HBsAg antibody concentration is expected to be higher than 250 mIU/ml should be diluted before use either 1:10 or 1:100 in the Calibrator 0 mIU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 µl of each specimen with 450 µl of Cal 0 (1:10). Then 50 µl of the 1:10 dilution are diluted with 450 µl of the Cal 0 (1:100). Mix tubes thoroughly on vortex when preparing the diluted samples.

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1. Microplate:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned green, indicating a defect in conservation.

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Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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Avoid contamination of the liquid with oxidising chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.

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Avoid contamination of the liquid with oxidising chemicals, air-driven dust or microbes. Do not expose to strong light, oxidising agents and metallic surfaces.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.

2. The ELISA incubator has to be set at +37°C (tolerance of ±1°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.

3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well-characterized negative and positive reference samples, and check to match the values reported below in the sections "Validation of Test" and "Assay Performances". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. Incubation times have a tolerance of ±5%.

5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

6. When using an ELISA automated workstation, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the

values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.

7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Control Serum as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
6. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipments are available and ready to use.

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M.1 Quantitative analysis

1. Place the required number of strips in the microplate holder. Leave A1 and B1 wells empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2..8°C, sealed. Then Dispense in all the wells to be used for the test, except for A1 and B1, 50µl of the Specimen Diluent.

Important note: This additive is added before distributing samples and controls into specific wells and is particularly

intended for blocking some substances present in people undergoing vaccination and capable to mask antibodies.

2. Pipette 100µl of all the Calibrators, 100µl of Control Serum in duplicate and then 100ul of samples. The Control Serum is used to verify that the whole analytical system works as expected. Check that Calibrators, Control Serum and samples have been correctly added. Then incubate the microplate at **+37°C for 60 min.**

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3. Wash the microplate as reported in section I.3.
4. In all the wells except A1 and B1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at **+37°C for 60 minutes.**

Important note:

- 1) *Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.*
- 2) *Mix thoroughly the Enzyme Conjugate on vortex before use.*

5. Wash the microplate as described.
6. Pipette 100µl TMB/H₂O₂ mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 minutes.**

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7. Stop the enzymatic reaction by pipette 100µl Sulphuric Acid into each well and using the same pipetting sequence as in step 6. Then measure the colour intensity with a microplate reader at 450nm (reading) and at 620-630nm (blanking, strongly recommended), blanking the instrument on A1 and B1 wells.

M.2 Qualitative analysis

1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2..8°C, sealed.
2. Dispense 50 ul Specimen Diluent in all the wells, except for the blank A1. Then pipette 100µl of the Calibrator 0 mIU/ml in duplicate, 100µl of the Calibrator 10 mIU/ml in duplicate, 100µl of the Calibrator 250 mIU/ml in single, and then 100ul of samples. Check that Calibrators and samples have been correctly added. Then incubate the microplate at **+37°C for 60 min.**
3. Wash the microplate as reported in section I.3.
4. In all the wells except A1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at **+37°C for 60 minutes.**

Important note:

- 3) *Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.*
- 4) *Mix thoroughly the Enzyme Conjugate on vortex before use.*

5. Wash the microplate as described.
6. Pipette 100µl TMB/H₂O₂ mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 minutes.**

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7. Stop the enzymatic reaction by pipette 100µl Sulphuric Acid into each well and using the same pipetting sequence as in step 6. Then measure the colour intensity with a microplate reader at 450nm (reading) and at 620-630nm (blinking, strongly recommended), blanking the instrument on A1 and B1 wells.

Important general notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
3. The Control Serum (CS) does not affect the cut-off calculation and therefore the test results calculation. The Control Serum may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME (standard procedure)

Specimen Diluent	50 ul
Calibrators	100 ul
Control Serum	100 ul
Samples	100 ul
1st incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H ₂ O ₂ mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm & 620nm

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
Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Calibrator 0 mIU/ml > 0.200	1. that the washing procedure and the washer settings are as validated in the pre qualification study;
coefficient of variation > 30%	2. that the proper washing solution has been used and the washer has been primed with it before use;
	3. that no mistake has been done in the assay procedure when the dispensation of standards is carried out;
	4. that no contamination of the Cal 0 mIU/ml or of the wells where it was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate;
	5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate
	6. that the washer needles are not blocked or partially obstructed.

Calibrator 10 mIU/ml OD450nm < Cal 0 + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of a wrong calibrator); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the standard has occurred.
Calibrator 250 mIU/ml OD450nm < 1.500 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the standard has occurred.
Control Serum Different from expected value	First verify that: 1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (ex.: dispensation of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the standard has occurred. 5. the Control Serum has been dissolved with the right volume reported on the label. If a mistake has been pointed out, the assay has to be repeated after eliminating the reason of this error. If no mistake has been found, proceed as follows: a) a value up to +/-20% is obtained: the overall Precision of the laboratory might not enable the test to match the expected value +/-10%. Report the problem to the Supervisor for acceptance or refusal of this result. b) a value higher than +/-20% is obtained: in this case the test is invalid and the DiaPro's customer service has to be called.

Example of Calibration Curve :

Important Note:
 Do not use the calibration curve above to make calculations.

P.2 Qualitative method

In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 10 mIU/ml and then check that the assay is valid.

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Calibrator 250 mIU/ml: 2.845 OD450nm
 Higher than 1.500 – Accepted

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Then on the calibration curve calculate the concentration of anti HBsAg antibody in samples.

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Q. INTERPRETATION OF RESULTS

Samples with a concentration lower than 10 WHO mIU/ml are considered negative for anti HBsAg antibody by most of the international medical literature.

Samples with a concentration higher than 10 WHO mIU/ml are considered positive for anti HBsAg antibody.

In the follow up of vaccination recipients, however, the value of 20 WHO mIU/ml is usually accepted by the medical literature as the minimum concentration at which the patient is considered clinically protected against HBV infection.

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. Diagnosis has to be done and released to the patient by a suitably qualified medical doctor.

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1. LIMIT OF DETECTION:

The limit of detection of the assay has been calculated by means of the HBsAb international preparation supplied by CLB on behalf of WHO (1st reference preparation 1977, lot 17-2-77), on which Calibration Curve has been calibrated. HBV negative serum was used as diluent, as recommended by the supplier. Results of Quality Control are given in the following table:

WHO mIU/ml	SAB.CE Lot # 1002	SAB.CE Lot # 1001	SAB.CE Lot # 1002/2
50	0.933	0.812	0.846
10	0.219	0.192	0.194
5	0.110	0.096	0.104
2.5	0.057	0.058	0.067
Std 0	0.021	0.015	0.023

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A Performance Evaluation has been conducted on a total number of more than 700 samples.

2.1 Diagnostic Specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte.

More than 500 negative specimens were tested, internally and externally, against a European company.

A diagnostic specificity of 98.8% was assessed. .

Moreover, diagnostic specificity was assessed by testing 113 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.) against the European company. A value of specificity of 100% was assessed.

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106 vaccinated patients were evaluated providing a diagnostic sensitivity of 100%.

More than 100 HBV naturally infected patients were tested, internally and externally, against the European company; a diagnostic sensitivity of 100% was found.

3. PRECISION:

The mean values obtained from a study conducted on three samples of different anti-HBsAg reactivity, examined in 16 replicates in three separate runs is reported below:

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4. ACCURACY

The assay accuracy has been checked by the dilution and recovery tests. Any "hook effect", underestimation likely to happen at high doses of analyte, was ruled out up to 10.000 mIU/ml.

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HBsAg_{one}

Version ULTRA

**Fourth generation Enzyme
Immunoassay (ELISA)
for the determination of
Hepatitis B surface Antigen or HBsAg
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

HBsAg One version ULTRA

A. INTENDED USE

Fourth generation Enzyme Immunoassay (ELISA) for the one-step determination of Hepatitis B surface Antigen or HBsAg in human plasma and sera.

The kit is intended for the screening of blood units, is able to detect HBsAg mutants and finds application in the follow-up of HBV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) defines Hepatitis B Virus infection as follows:

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HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. The main ways of getting infected with HBV are: (a) perinatal (from mother to baby at the birth); (b) child- to-child transmission; (c) unsafe injections and transfusions; (d) sexual contact.

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Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 30% to 50% of children infected between 1 to 4 years of age develop chronic infection. The risk of death from HBV-related liver cancer or cirrhosis is approximately 25% for persons who become chronically infected during childhood. Chronic hepatitis B in some patients is treated with drugs called *interferon* or *lamivudine*, which can help some patients. Patients with cirrhosis are sometimes given liver transplants, with varying success. It is preferable to prevent this disease with vaccine than to try and cure it.

Hepatitis B vaccine has an outstanding record of safety and effectiveness. Since 1982, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 95% effective in preventing children and adults from developing chronic infection if they have not yet been infected. In many countries where 8% to 15% of children used to become chronically infected with HBV, the rate of chronic infection has been reduced to less than 1% in immunized groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programs."

Hepatitis B surface Antigen or HBsAg is the most important protein of the envelope of Hepatitis B Virus, responsible for acute and chronic viral hepatitis.

The surface antigen contains the determinant "a", common to all the known viral subtypes, immunologically distinguished by two distinct subgroups (ay and ad).

The ability to detect HBsAg with high sensitive immunoassays in the last years has led to an understanding of its distribution and epidemiology worldwide and to radically decrease the risk of infection in transfusion.

C. PRINCIPLE OF THE TEST

A mix of mouse monoclonal antibodies specific to the determinants "a", "d" and "y" of HBsAg is fixed to the surface of microwells. Patient's serum/plasma is added to the microwell together with a second mix of mouse monoclonal antibodies, conjugated with Horseradish Peroxidase (HRP) and directed against a different epitope of the determinant "a" and against "preS".

The specific immunocomplex, formed in the presence of HBsAg in the sample, is captured by the solid phase.

At the end of the one-step incubation, microwells are washed to remove unbound serum proteins and HRP conjugate.

The chromogen/substrate is then added and, in the presence of captured HBsAg immunocomplex, the colorless substrate is hydrolyzed by the bound HRP conjugate to a colored end-product. After blocking the enzymatic reaction, its optical density is measured by an ELISA reader.

The color intensity is proportional to the amount of HBsAg present in the sample.

The version ULTRA is particularly suitable for automated screenings and is able to detect "s" mutants.

D. COMPONENTS

The standard configuration contains reagents to perform 192 tests and is made of the following components:

1. Microplate MICROPLATE

n° 2. 12 strips of 8 breakable wells coated with anti HBsAg, affinity purified mouse monoclonal antibodies, specific to "a", "y" and "d" determinants, and sealed into a bag with desiccant.

2. Negative Control CONTROL -

1x4.0ml/vial. Ready to use control. It contains goat serum, 10 mM phosphate buffer pH 7.4+/-0.1, 0.09% Na-azide and 0.1% Kathon GC as preservatives. The negative control is pale yellow color coded.

3. Positive Control CONTROL +

1x4.0ml/vial. Ready to use control. It contains goat serum, non infectious recombinant HBsAg, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.1% Kathon GC as preservatives. The positive control is color coded green.

4. Calibrator CAL ...

n° 2 vials. Lyophilized calibrator. To be dissolved with EIA grade water as reported in the label. Contains fetal bovine serum, non infectious recombinant HBsAg at 0.5 IU/ml (2nd WHO international standard for HBsAg, NIBSC code 00/588), 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.1% Kathon GC as preservatives.

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5. Wash buffer concentrate WASHBUF 20X

2x60ml/bottle. 20X concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

6. Enzyme Conjugate Diluent CONJ DIL

2x16ml/vial. Ready to use and pink/red color coded reagent. It contains 10 mM Tris buffer pH 6.8+/-0.1, 1% normal mouse serum, 5% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives. The solution is normally opalescent.

7. Enzyme Conjugate CONJ 20X

2x1ml/vial. 20X concentrated reagent. It contains Horseradish Peroxidase (HRP) labeled mouse monoclonal antibodies to HBsAg, determinant "a" and "preS", 10 mM Tris buffer pH 6.8+/-0.1, 5% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

8. Chromogen/Substrate SUBS TMB

2x25ml/bottle. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine (TMB) and 0.02% hydrogen peroxide (H₂O₂).

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9. Sulphuric Acid H₂SO₄ 0.3 M

1x25ml/bottle. It contains 0.3 M H₂SO₄ solution.

Note: Attention: Irritant (H315; H319; P280; P302+P352; P332+P313; P305+P351+P338; P337+P313; P362+P363)

10. Plate sealing foils n° 4

11. Package insert

Important note:

Only upon specific request, Dia.Pro can supply reagents for 96, 480, 960 tests, as reported below:

	N°1	N°5	N°10
Microplates	N°1	N°5	N°10
Negative Control	1x2ml/vial	1x10ml/vial	1x20ml/vial
Positive Control	1x2ml/vial	1x10ml/vial	1x20ml/vial
Calibrator	N° 1 vial	N° 5 vials	N° 10 vials
Wash buffer concentrate	1x60ml/vial	5x60ml/vial	4x150ml/vial
Enzyme conjugate	1x0.8ml/vial	1x4ml/vial	2x4ml/vial
Conjugate Diluent	1x16ml/vial	2x40ml/vial	2x80ml/vial
Chromogen/Substrate	1x25ml/vial	3x42ml/vial	2x125ml/vial
Sulphuric Acid	1x15ml/vial	2x40ml/vial	2x80ml/vial
Plate sealing foils	N° 2	N° 10	N° 20
Package insert	N° 1	N° 1	N° 1
Number of tests	96	480	960
Code SAGIULTRA.CE	96	480	960

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1. Calibrated Micropipettes (150ul, 100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet), capable to provide shaking at 1300 rpm+/-150, set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
6. Upon receipt, store the kit at 2.8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
11. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-use of the device and up to 6 months.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
15. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water
17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and lipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as well as they could give rise to false positive results. Specimens with an altered pathway of coagulation, presenting particles after blood collection and preparation of serum/plasma as those coming from hemodialized patients, could give origin to false positive results.
5. Sera and plasma can be stored at +2°..8°C for up to seven days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen sample should not be frozen/thawed more than once as this may generate particles that could affect the test result. If some turbidity is present or presence of microparticles is suspected after thawing, filter the sample on a disposable 0.2-0.8µ filter to clean it up for testing or use the two-steps alternative method.

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1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned green, indicating a defect in conservation. In this case, call Dia.Pro's customer service. Unused strips have to be placed back inside the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°..8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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3. Positive Control:

Ready to use. Mix well on vortex before use. The positive control does not contain any infective HBV as it is composed of recombinant synthetic HBsAg.

4. Calibrator:

Add the volume of ELISA grade water, reported on the label, to the lyophilized powder; let fully dissolve and then gently mix on vortex. The solution is not stable. Store the Calibrator frozen in aliquots at -20°C.

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Note: Once diluted, the wash solution is stable for 1 week at +2°..8°C.

6. Enzyme conjugate:

The working solution is prepared by diluting the 20X concentrated reagent into the Conjugate Mix well on vortex before use.

Avoid any contamination of the liquid with oxidizing chemicals, dust or microbes. If this component has to be transferred, use only plastic sterile disposable containers.

Important note: The working solution is not stable. Prepare only the volume necessary for the work of the day. As an example when the kit is used in combination with other instruments or manually, dilute 0.1 ml 20X Conjugate with 1.9 ml Conjugate Diluent into a disposable plastic vial and mix carefully before use.

7. Chromogen/Substrate:

Ready to use. Mix well by end-over-end mixing.

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8. Sulphuric Acid:

Ready to use. Mix well by end-over-end mixing.

Attention: Irritant (H315; H319; P280; P302+P352; P332+P313; P305+P351+P338; P337+P313; P362+P363).

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1. **Micropipettes** have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.
2. The **ELISA incubator** has to be set at +37°C (tolerance of ±1°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. In case of **shaking** during incubations, the instrument has to ensure 350 rpm ±150. Amplitude of shaking is very important as a wrong one could give origin to splashes and therefore to some false positive result.
4. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls/calibrator and reference panels, before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350µl/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it

is recommended to run an assay with the kit controls/calibrator and well-characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

5. **Incubation times** have a tolerance of $\pm 5\%$.
6. The **microplate reader** has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; (d) repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
7. When using **ELISA automated workstations**, all critical steps (dispensation, incubation, washing, reading, shaking, data handling, etc.) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated by checking full matching the declared performances of the kit. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set paying particular attention to avoid carry over by the needles used for dispensing samples and for washing. The carry over effect must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.
8. When using automatic devices, in case the vial holder of the instrument does not fit with the vials supplied in the kit, transfer the solution into appropriate containers and label them with the same label peeled out from the original vial. This operation is important in order to avoid mismatching contents of vials, when transferring them. When the test is over, return the secondary labeled containers to 2..8°C, firmly capped.
9. **Dia.Pro's customer service** offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the essential requirements of the assay. Support is also provided for the installation of new instruments to be used in combination with the kit.

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1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dilute the 20X concentrated Enzyme Conjugate with its Diluent as reported.
5. Dissolve the Calibrator as described above.
6. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
7. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution,

according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.

8. Check that the ELISA reader has been turned on at least 20 minutes before reading.
9. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
10. Check that the micropipettes are set to the required volume.
11. Check that all the other equipment is available and ready to use.
12. In case of problems, do not proceed further with the test and advise the supervisor.

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Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument dispense first 150 ul controls & calibrator, then all the samples and finally 100 ul diluted Enzyme Conjugate.

For the pre-washing step (point 1 of the assay procedure) and all the next operations follow the operative instructions reported below for the Manual Assay.

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Manual Assay:

1. Place the required number of strips in the plastic holder and wash them once to hydrate wells. Carefully identify the wells for controls, calibrator and samples.

Important note: *Pre washing (1 cycle: dispensation of 350ul/well of washing solution+ aspiration) is fundamental to obtain reliable and specific results both in the manual and in the automatic procedures. Do not omit it !*

2. Leave the A1 well empty for blanking purposes.
3. Pipette 150ul of the Negative Control in triplicate, 150ul of the Calibrator in duplicate and then 150ul of the Positive Control in single followed by 150ul of each of the samples.
4. Check for the presence of samples in wells by naked eye (there is a marked color difference between empty and full wells) or by reading at 450/620nm. (samples show OD values higher than 0.100).
5. Dispense 100ul diluted Enzymatic Conjugate in all wells, except for A1, used for blanking operations.

Important note: *Be careful not to touch the inner surface of the well with the pipette tip when the conjugate is dispensed. Contamination might occur.*

6. Following addition of the conjugate, check that the color of the samples have changed from yellowish to pink/red and then incubate the microplate for **120 min at +37°C**.

Important notes:

- a. *Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.*
- b. *If the procedure is carried out on shaking, be sure to deliver the rpm reported for in Section I.3 as otherwise intra-well contamination could occur.*
7. When the first incubation is over, wash the microwells as previously described (section I.4)

8. Pipette 200 µl Chromogen/Substrate into all the wells, A1 included.

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9. Incubate the microplate protected from light at **18-24°C for 30 min**. Wells dispensed with the positive control, the calibrator and positive samples will turn from clear to blue.
10. Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction, using the same pipetting sequence as in step 8. Addition of the acid solution will turn the positive control, the calibrator and positive samples from blue to yellow/brown.
11. Measure the color intensity of the solution in each well, as described in section I.6 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1.

Important general notes:

- If the second filter is not available, ensure that no fingerprints or dust are present on the external bottom of the microwell before reading at 450nm. They could generate false positive results on reading.
- Reading should ideally be performed immediately after the addition of the acid solution but definitely no longer than 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.
- When samples to be tested are not surely clean or have been stored frozen, the assay procedure reported below is recommended as long as it is far less sensitive to interferences due to hemolysis, hyperlipaemia, bacterial contamination and fibrin microparticles. The assay is carried out in two-steps at +37°C on shaking at 350 rpm ±150 as follows:
 - dispense 100 µl of controls, calibrator and samples
 - incubate 60 min at +37°C on shaking
 - wash according to instructions (section I.4)
 - dispense 100 µl diluted enzyme tracer
 - incubate 30 min at +37°C on shaking
 - wash
 - dispense 100 µl TMB&H₂O₂ mix
 - incubate 30 min at r.t. on shaking
 - stop and read

In this procedure the pre-wash can be omitted. This method shows performances similar to the standard one and therefore can be used in alternative.
- The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

Operations	Procedure
Pre-Washing step	n° 1 cycle
Controls&Calibrator&samples	150 µl
Diluted Enzyme Conjugate	100 µl
1st incubation	120 min
Temperature	+37°C
Washing steps	n° 4-5
Chromogen/Substrate	200µl
2nd incubation	30 min
Temperature	room
Sulphuric Acid	100 µl
Reading OD	450nm

An example of dispensation scheme is reported in the following section:

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O. INTERNAL QUALITY CONTROL

A check is performed on the controls/calibrator any time the kit is used in order to verify whether the expected OD450nm or S/Co values have been matched in the analysis.

Ensure that the following results are met:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	< 0.050 mean OD450nm value after blanking
Calibrator 0.5 IU/ml	S/Co ≥ 2
Positive Control	> 1.000 OD450nm value

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.

Calibrator S/Co < 2	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control < 1.000 OD450nm	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.050). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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P. CALCULATION OF THE CUT-OFF

The test results are calculated by means of a cut-off value determined on the mean OD450nm value of the negative control (NC) with the following formula:

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Q. INTERPRETATION OF RESULTS

Test results are interpreted as a ratio of the sample OD450nm (S) and the Cut-Off value (Co), mathematically S/Co, according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 – 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient is not infected by HBV and that the blood unit may be transfused.

Any patient showing an equivocal result should be retested on a second sample taken 1-2 weeks after the initial sample; the blood unit should not be transfused.

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Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Any positive result must be confirmed first by repeating the test on the sample, after having filtered it on 0.2-0.8 µ filter to remove any microparticles interference. Then, if still

positive, the sample has to be submitted to a confirmation test before a diagnosis of viral hepatitis is released.

3. When test results are transmitted from the laboratory to another department, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis infection has to be taken and released to the patient by a suitably qualified medical doctor.

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Negative Control: 0.012 – 0.008 – 0.010 OD450nm
Mean Value: 0.010 OD450nm
Lower than 0.050 – Accepted
Positive Control: 2.489 OD450nm
Higher than 1.000 – Accepted
Cut-Off = 0.010+0.050 = 0.060
Calibrator: 0.350 - 0.370 OD450nm
Mean value: 0.360 OD450nm S/Co = 6.0
S/Co higher than 2.0 – Accepted
Sample 1: 0.028 OD450nm
Sample 2: 1.690 OD450nm
Sample 1 S/Co < 0.9 = negative
Sample 2 S/Co > 1.1 = positive

R. PERFORMANCE CHARACTERISTICS

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC). Version ULTRA proved to be at least equivalent to the original design in a study conducted for the validation of the new version.

1. Analytical Sensitivity

The limit of detection of the assay has been calculated on the 2nd WHO international standard, NIBSC code 00/588.

In the following table, results are given for three lots (P1, P2 and P3) of the version ULTRA in comparison with the reference device (Ref.):

WHO IU/ml	Lot # P1 S/Co	Lot # P2 S/Co	Lot # P3 S/Co	Ref. S/Co
0.4	4.6	4.8	4.6	4.6
0.2	2.3	2.4	2.4	2.4
0.1	1.4	1.4	1.5	1.2
0.05	0.8	0.8	1.0	0.7
0.025	0.6	0.6	0.6	0.4
FCS (NC)	0.3	0.2	0.3	0.1

The assay shows an Analytical Sensitivity better than 0.1 WHO IU/ml of HBsAg.

In addition two panels of sensitivity supplied by EFS, France, and by SFTS, France, were tested and gave in the best conditions the following results:

Panel EFS Ag HBs HB1-HB6 lot n° 04

Sample ID	Characteristics	ng/ml	S/Co
HB1	diluent	/	0,2
HB2	adw2+ayw3	0,05	0,6
HB3	adw2+ayw3	0,1	1,0
HB4	adw2+ayw3	0,2	1,8
HB5	adw2+ayw3	0,3	2,4
HB6	adw2+ayw3	0,5	4,2

Sensitivity panel SFTS, France, Ag HBs 2005

Sample ID	Characteristics	ng/ml	S/Co
171	Adw2 + ayw3	2.21 ± 0.15	15,4
172	Adw2 + ayw3	1.18 ± 0.10	8,7
173	Adw2 + ayw3	1.02 ± 0.05	6,1
174	Adw2 + ayw3	0.64 ± 0.04	4,0
175	Adw2 + ayw3	0.49 ± 0.03	3,4
176	Adw2 + ayw3	0.39 ± 0.02	2,6
177	Adw2 + ayw3	0.25 ± 0.02	2,0
178	Adw2 + ayw3	0.11 ± 0.02	1,3
179	Adw2 + ayw3	0.06 ± 0.01	0,9
180	Adw2 + ayw3	0.03 ± 0.01	0,8
181	Adw2	0.5 – 1.0	4,7
182	Adw4	0.5 – 1.0	3,6
183	Adr	0.5 – 1.0	4,5
184	Ayw1	0.5 – 1.0	5,1
185	Ayw2	0.5 – 1.0	6,4
186	Ayw3	0.5 – 1.0	7,3
187	Ayw3	0.5 – 1.0	5,8
188	Ayw4	0.5 – 1.0	6,9
189	Ayr	0.5 – 1.0	6,1
190	diluent	/	0,6

The panel # 808, supplied by Boston Biomedical Inc., USA, was also tested to define the limit of sensitivity. Results in the best conditions are as follows :

BBI panel PHA 808

Sample ID	Characteristics	ng/ml	S/Co
01	ad	2,49	10,2
02	ad	1,17	4,8
03	ad	1,02	4,3
04	ad	0,96	3,8
05	ad	0,69	2,9
06	ad	0,50	2,2
07	ad	0,41	1,5
08	ad	0,37	1,3
09	ad	0,30	1,2
10	ad	0,23	1,0
11	ay	2,51	11,2
12	ay	1,26	5,9
13	ay	0,97	4,1
14	ay	0,77	3,7
15	ay	0,63	2,0
16	ay	0,48	2,4
17	ay	0,42	2,0
18	ay	0,33	1,8
19	ay	0,23	1,6
20	ay	0,13	1,1
21	negative	/	0,6

2. Diagnostic Sensitivity:

The diagnostic sensitivity was tested according to what required by Common Technical Specifications (CTS) of the directive 98/79/EC on IVD for HBsAg testing.

Positive samples, including HBsAg subtypes and a panel of "s" mutants from most frequent mutations, were collected from different HBV pathologies (acute, a-symptomatic and chronic hepatitis B) or produced synthetically, and were detected positive in the assay.

All the HBsAg known subtypes, "ay" and "ad", and isoforms "w" and "r", supplied by CNTS, France, were tested in the assay and determined positive by the kit as expected.

An overall value of 100% has been found in a study conducted on a total number of more than 400 samples positive with the original reference IVD code SAG1.CE, CE marked.

A total of 30 sero-conversions were studied, most of them produced by Boston Biomedica Inc., USA.

Results obtained by examining eight panels supplied by Boston Biomedica Inc., USA, are reported below for the version ULTRA in comparison with the reference device code SAG1.CE.

Panel ID	1 st sample positive	HBsAg subtype	HBsAg ng/ml	Version ULTRA S/Co	Ref. device S/Co
PHM 906	02	ad	0.5	3.7	1.4
PHM 907 (M)	06	ay	1.0	4.4	2.9
PHM 909	04	ad	0.3	1.2	0.8
PHM 914	04	ad	0.5	1.1	1.1
PHM 918	02	ad	0.1	1.8	0.5
PHM 923	03	ay	< 0.2	2.2	1.2
PHM 925	03	Ind.	n.d.	1.4	0.9
PHM 934	01	ad	n.d.	1.0	0.8

3. Diagnostic Specificity:

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. In addition to the first study, where more than 5000 negative samples from blood donors (two blood centers), classified negative with a CE marked device in use at the laboratory of collection were examined, the diagnostic specificity was recently assessed by testing a total of 2288 negative blood donors on seven different lots. A value of specificity of 100% was found.

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Samples derived from patients with different viral (HCV, HAV) and non viral pathologies of the liver that may interfere with the test were examined. No cross reaction were observed.

4. Precision:

It has been calculated for the version ULTRA on two samples examined in 16 replicates in 3 different runs for three lots.

Results are reported in the following tables:

Average values	Negative Sample	Calibrator 0.5 IU/ml
Total n = 144	0.026	0.332
Std.Deviation	0.004	0.027
CV %	16%	8%

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S. LIMITATIONS

Repeatable false positive results were assessed on freshly collected specimens in less than 0.1% of the normal population, mostly due to high titers Heterophilic Anti Mouse Antibodies (HAMA).

Interferences in fresh samples were also observed when they were not particles-free or were badly collected (see chapter G).

Old or frozen samples, presenting fibrin clots, cryoglobulins, lipid-containing micelles or microparticles after storage or thawing, can generate false positive results.

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0318

HEV IgG

**Third generation Enzyme Immunoassay
for the determination of IgG antibodies
to Hepatitis E Virus
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF EVG.CE
96 Tests

HEV IgG

A. INTENDED USE

Third generation Enzyme ImmunoAssay (ELISA) for the qualitative determination of IgG antibodies to Hepatitis E Virus in human plasma and sera.

The kit is intended for the follow-up of HEV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

Hepatitis E Virus or HEV is a recently discovered agent of enterically transmitted viral hepatitis. HEV is an un-enveloped single-strand RNA virus, after being provisionally assigned to the Caliciviridae family, HEV was re-classified as the sole member of the genus Hepevirus, family Hepeviridae, in 2004. HEV is found in the stool of infected patients and present in 4 strains (1, 2, 3 and 4) differently spread geographically and virulent.

HEV is a serious problem in many developing countries since its first outbreak was reported in 1955 in New Delhi, India.

A high case-fatality rate has been found among pregnant women and chronic hepatitis carriers.

The cloning and sequencing of HEV genome have led to the development of serological tests for the detection of anti HEV antibodies based on recombinant immunodominant antigens derived from conservative regions of the four virus strains.

C. PRINCIPLE OF THE TEST

Microplates are coated with HEV-specific recombinant antigens encoding for conservative and immunodominant determinants of all the 4 subtypes.

The solid phase is first treated with the diluted sample and anti HEV IgG are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti-HEV IgG are detected by the addition of polyclonal specific anti hlgG antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti HEV IgG present in the sample. A cut-off value let optical densities be interpreted into anti-HEV IgG negative and positive results.

D. COMPONENTS

Code EVG.CE contains reagents for 96 tests.

1. Microplate MICROPLATE

n° 1 microplate. 12 strips of 8 microwells coated with HEV specific recombinant antigens. Plates are sealed into a bag with desiccant.

2. Negative Control CONTROL -

1x4.0ml/vial. Ready to use control. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. The negative control is olive green colour coded.

3. Positive Control CONTROL +

1x4.0ml/vial. Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HEV, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives.

The Positive Control is blue colour coded.

4. Calibrator CAL ...

n° 1 vial. Lyophilized calibrator. To be dissolved with the volume of EIA grade water reported on the label.

It contains foetal bovine serum proteins, human antibodies to HEV whose content is calibrated on 1st WHO reference reagent for HEV antibody, NIBSC code 95/584 at 1 IU/ml +/-20%, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.3 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

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6. Enzyme Conjugate CONJ

1x16ml/vial. Ready to use and red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives.

7. Chromogen/Substrate SUBS TMB

1x16ml/vial. Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

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8. Assay Diluent DILAS

1x8ml/vial. 10 mM tris buffered solution pH 8.0 +/-0.1 containing 0.045% ProClin 300 for the pre-treatment of samples and controls in the plate, blocking interference.

9. Sulphuric Acid H₂SO₄ 0.3 M

1x15ml/vial. It contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315, H319; P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P363).

10. Sample Diluent: DILSPE

1x50ml/vial. It contains 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. To be used to dilute the sample.

Note: The diluent changes colour from olive green to dark bluish green in the presence of sample.

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F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and

qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.

3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

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6. Upon receipt, store the kit at +2..8°C into a temperature controlled refrigerator or cold room.

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9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.

10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.

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14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

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16. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been

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4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

5. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection.

Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.

6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

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A study conducted on an opened kit has not pointed out any relevant loss of activity up to 6 re-use of the device and up to 6 months.

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Unused strips have to be placed back into the aluminium pouch, in presence of desiccant supplied, firmly zipped and stored at +2°..8°C.

When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.

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3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infective, even if HCV, eventually present in the control, has been chemically inactivated.

4. Calibrator:

Dissolve carefully the content of the lyophilized vial with the volume of EIA grade water reported on its label.

Mix well on vortex before use.

Handle this component as potentially infective, even if HEV, eventually present in the control, has been chemically inactivated.

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5. Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bidistilled water up to 1200 ml and mixed gently end-over-end before use.

Once diluted, the wash solution is stable for 1 week at +2..8° C.

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Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.
If this component has to be transferred use only plastic, possibly sterile disposable containers.

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- Micropipettes** have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
- The **ELISA incubator** has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution.

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An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.

- Incubation times** have a tolerance of $\pm 5\%$.
- The **ELISA reader** has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, mandatory) for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an **ELISA automated work station**, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended for blood screening when the number of samples to be tested exceed 20-30 units per run.
- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
- Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Dissolve the Calibrator as described above.
- Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
- Check that the ELISA reader has been turned on at least 20 minutes before reading.
- If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

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Automated assay:

In case the test is carried out automatically with an ELISA workstation, we suggest to make the instrument aspirate 200 ul Sample Diluent and then 10 ul sample.

All the mixture is then carefully dispensed directly into the appropriate sample well of the microplate. Before the next sample is aspirated, needles have to be duly washed to avoid any cross-contamination among samples.

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Dispense 200 ul controls/calibrator in the appropriate control/calibration wells.

Important Note: *Visually monitor that samples have been diluted and dispensed into appropriate wells. This is simply achieved by checking that the color of dispensed samples has turned to dark bluish-green while the color of the negative control has remained olive green.*

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Manual assay:

- Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
- Dispense 200 ul of Negative Control in triplicate, 200 ul Calibrator in duplicate and 200 ul Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use !
- Add 200 ul of Sample Diluent (DILSPE) to all the sample wells; then dispense 10 ul sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: *Check that the color of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.*

- Dispense 50 ul Assay Diluent (DILAS) into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
- Incubate the microplate for **45 min at +37°C**.

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- Wash the microplate with an automatic washer by delivering and aspirating 300ul/well of diluted washing solution as reported previously (section I.3).
- Pipette 100µl Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red colored component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **45 min at +37°C**.
- Wash microwells as in step 6.

- Pipette 100µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 15 minutes**.

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- Pipette 100µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of acid will turn the positive control and positive samples from blue to yellow.
- Measure the color intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1.

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N. ASSAY SCHEME

Method	Operations
Controls & Calibrator(*)	200 ul
Samples	200ul dil.+10ul
Assay Diluent (DILAS)	50 ul
1st incubation	45 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 ul
2nd incubation	45 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2	100 ul
3rd incubation	15 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm/620-630nm

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O. INTERNAL QUALITY CONTROL

A check is carried out on the controls any time the kit is used in order to verify whether their OD450nm/620-630nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	< 0.050 mean OD450nm value after blanking
Positive Control	> 1.000 OD450nm value

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not got contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of their wells has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150, too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Check	Requirements
Calibrator	S/Co > 1.1

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Problem	Check
Calibrator S/Co < 1.1	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (e.g.: dispensation of negative control instead of Calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

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P. CALCULATION OF THE CUT-OFF

The tests results are calculated by means of a cut-off value determined with the following formula:

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Important note: *When the calculation of results is done by the operative system of an ELISA automated work station be sure that the proper formulation is used to calculate the cut-off value and generate the right interpretations of results.*

Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD450nm/620-630nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HEV or that the blood unit may be transfused.

Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.

A positive result is indicative of HEV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

- 1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.*
- 2. By definition, any positive result should be confirmed by an alternative method before a diagnosis of viral hepatitis is formulated.*
- 3. When test results are transmitted from the laboratory to an informatics center, attention has to be done to avoid erroneous data transfer.*
- 4. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.*

An example of calculation is reported below (data obtained proceeding as the reading step described in the section M, point 12):

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*Negative Control: 0.019 – 0.020 – 0.021 OD450nm
Mean Value: 0.020 OD450nm
Lower than 0.050 – Accepted
Positive Control: 2.189 OD450nm
Higher than 1.000 – Accepted
Cut-Off = 0.020+0.350 = 0.370
Calibrator: 0.600 - 0.640 OD450nm
Mean value: 0.620 OD450nm S/Co = 1.7
S/Co higher than 1.1 – Accepted*

Sample 1: 0.070 OD450nm
Sample 2: 1.690 OD450nm
Sample 1 S/Co < 0.9 = negative
Sample 2 S/Co > 1.1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted on negative and positive samples in an external clinical center with reference to a FDA approved kit.

1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated by means of 1st WHO reference reagent for HEV antibody, NIBSC code 95/584. The assay shows an analytical sensitivity of about 0.1 WHO IU/ml.

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They were checked with a sensitivity set at 0.25 WHO IU/ml on more than total 700 samples.

2.1 Diagnostic Specificity:

It is defined as the probability of the assay of scoring negative in the absence of specific analyte.

A total number of more than one hundred 1-5 years old children, by definition negative for HEV antibodies as they never chanced to eat uncooked swine meat and get therefore infected, were tested; a value of 100% specificity (negativity) was assured at a sensitivity set at a cut-off of 0.25 WHO IU/ml.

In addition, a total of more 500 unselected donors, including 1st time donors, and HEV negative hospitalized patients, coming from Italy, were examined maintaining the sensitivity of 0.25 WHO IU/ml. About 5% of such population turned out to be repeatedly positive; confirmation was not carried out in absence of a commercial Confirmation kit. However, these samples were detected positive with a commercial CE-marked ELISA. From this study a diagnostic specificity of 100% was observed.

Moreover, the Diagnostic Specificity was also assessed by testing more than 100 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, haemolized, lipemic, etc.). A value of specificity of 100% was assessed.

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No interference was observed.

2.2 Diagnostic Sensitivity

It is defined as the probability of the assay of scoring positive in the presence of specific analyte.

The diagnostic sensitivity has been assessed externally (Institute of Virology, University of Milan) and internally on a total number of 200 positive specimens (maintaining a sensitivity set at a cut-off of 0.25 WHO IU/ml); a diagnostic sensitivity (or correlation with a commercial reference kit manufactured in Europe and CE marked) of 100% was found.

3. PRECISION:

It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs.

CV values ranging between 5-10%, depending on OD450nm/620-630nm values, were found. The variability seen did not result in sample misclassification.

S. LIMITATIONS

Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HEV IgG CODE: EVG.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – JULY 2004
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HEV IgM CODE: EVM.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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Rev: 05/2018

HEV IgM

**Enzyme Immunoassay (ELISA)
for the determination of IgM antibodies
to Hepatitis E Virus
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

HEV IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of IgM antibodies to Hepatitis E Virus in human plasma and sera. The kit may be used for the determination of the acute phase of infection where IgM antibodies are generated before the other immunoglobulins and for the follow-up of HEV-infected patients. For "in vitro" diagnostic use only.

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The solid phase is first treated with the diluted sample and anti HEV IgM are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound anti-HEV IgM are detected by the addition of polyclonal specific anti hIgM antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti HEV IgM present in the sample. A cut-off value let optical densities be interpreted into anti-HEV IgM negative and positive results.

Neutralization of IgG anti-HEV and Rheumatoid Factor, carried out directly in the well, is performed in the assay in order to block such kind of interferences.

D. COMPONENTS

Code EVM.CE contains reagents for 96 tests.

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2. Negative Control CONTROL -

1x4.0ml/vial. Ready to use control. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. The negative control is yellow colour coded.

3. Positive Control CONTROL +

1x4.0ml/vial. Ready to use control. It contains 1% goat serum proteins, human anti HEV IgM, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives.

The Positive Control is dark green colour coded.

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5. Enzyme Conjugate CONJ

1x16ml/vial. Ready to use and red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgM, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives.

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8. Specimen Diluent: DILSPE

2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. To be used to dilute the sample.

9. Neutralizing Reagent: SOLN NTR

1x8ml/vial. Ready-to-use Reagent. It contains goat anti hIgG, 2% casein, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.1% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives.

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1. Calibrated Micropipettes (100ul and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the

test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.

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8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.

9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.

10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.

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2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.

3. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

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Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be

frozen/thawed more than once as this may generate particles that could affect the test result.

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Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of conservation.

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- Micropipettes** have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
- The **ELISA incubator** has to be set at +37°C (tolerance of +/-0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing.
An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
- Incubation times have a tolerance of ±5%.
- The **ELISA reader** has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an **ELISA automated work station**, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the

possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates.
- Check that the Chromogen (TMB) is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
- Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
- Dissolve the content of the Control Serum as reported.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as as reported in the specific section.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

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- Dilute samples 1:101 into a properly defined dilution tube (example: 1000 µl Sample Diluent + 10 µl sample). Do not dilute the negative Control and the Positive Control as they are ready to use. Mix carefully all the liquid components on vortex and then proceed as described below.
- Place the required number of Microwells in the microwell holder. Leave A1 well empty for the operation of blanking.
- Dispense 50 µl of the Neutralizing Reagent (SOLN NTR) in all the wells of the samples. Do not add it in the wells used for Controls and in A1!

Important note: *The Neutralizing Reagent is able to block false positive reactions due to RF. Positive samples in internal QC panels might be detected negative if such samples were tested positive with an IVD that does not carry out any RF blocking reaction.*

- Dispense 100 µl of Negative Control in triplicate and 100 µl of Positive control in single. Then dispense 100 µl of diluted samples in each properly identified well.
- Incubate the microplate for **60 min at +37°C**.

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- Wash the microplate with an automatic washer by delivering and aspirating 300 µl/well of diluted washing solution as reported previously (section I.3).
- Pipette 100 µl Enzyme Conjugate into each well, except the A1 well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

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- Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10. Addition of acid will turn the positive calibrators, the control serum and the positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, mandatory), blanking the instrument on A1.

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N. ASSAY SCHEME

Method	Operations
Neutralizing Reagent (in sample wells only !)	50 µl
Negative and Positive Controls	100 µl
Samples diluted 1:101	100 µl
1st incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 µl
2nd incubation	60 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H2O2	100 µl
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 µl
Reading OD	450nm/620-630nm

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Legenda: BLK = Blank NC = Negative Control
PC = Positive Control S = Sample

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
Negative Control (NC) > 0.100 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of their wells has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 0.500 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150, too). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD_{450nm/620-630nm} and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 1.0	Negative
1.0 – 1.2	Equivocal
> 1.2	Positive

A negative result indicates that the patient has no detectable anti HEV IgM reactivity.

Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined.

A positive result is indicative of HEV infection and therefore the patient should be treated accordingly.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method before a diagnosis of viral hepatitis is formulated.
3. When test results are transmitted from the laboratory to an informatics center, attention has to be done to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis E infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below (data obtained proceeding as the the reading step described in the section M, point 12):

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1. LIMIT OF DETECTION

The limit of detection of the product has been checked on the international reference reagent for HEV antibody supplied by NIBSC/WHO with code n° 95/584. This material was assessed to be positive also for anti HEV IgM, low titer. The observed value is about 1 IU/ml.

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They were checked on about 700 sample derived from acute infections, HEV Ab positive patients, random individuals under diagnostic examination and healthy individuals with a sensitivity set at about 5 IU/ml in order to assure the highest sensitivity and be able to detect primary infection at the earliest phase.

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A total of more 500 unselected donors and HEV negative hospitalized patients, including 1st time donors, were examined. The diagnostic specificity was assessed against a kit US FDA approved. A diagnostic specificity of ≥ 95% was observed.

Moreover, diagnostic specificity was assessed by testing more than 100 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, haemolized, lipemic, etc.). A value of specificity of 100% was assessed.

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No interference was observed.

High reactive RF positive samples were observed to give origin to false positive results in not more than 5% of HEV Ab negative individuals.

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The diagnostic sensitivity has been assessed externally and internally on a total number of 100 positive specimens coming from Germany, Mexico and from Burma.

A diagnostic sensitivity of 100% was found.

3. PRECISION:

It has been calculated on two samples, one negative and one positive, examined in 16 replicates in three separate runs.

CV values ranging between 5-15%, depending on OD_{450nm} values, were found. The variability seen did not result in sample misclassification.

S. LIMITATIONS

Repeatable false positive results were assessed for high titer RF positive samples, escaping the effect of the Neutralizing Reagent.

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HDV Ab

**Competitive Enzyme Immunoassay
for the qualitative determination of
antibodies to Hepatitis Delta Virus
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF DAB.CE
96 Tests

HDV Ab

A. INTENDED USE

Competitive Enzyme ImmunoAssay (ELISA) for the qualitative determination of antibodies to Hepatitis Delta Virus or HDV in human plasma and sera with a "two-steps" methodology.

The kit is used for the follow-up of patients infected by HDV.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The Hepatitis Delta Virus or HDV is a RNA defective virus composed of a core presenting the delta-specific antigen, encapsulated by HBsAg, that requires the helper function of HBV to support its replication.

Infection by HDV occurs in the presence of acute or chronic HBV infection. When acute delta and acute HBV simultaneously occur, the illness becomes severe and clinical and biochemical features may be indistinguishable from those of HBV infection alone. In contrast, a patient with chronic HBV infection can support HDV replication indefinitely, usually with a less severe illness appearing as a clinical exacerbation.

The determination of HDV specific serological markers (HDV Ag, HDV Ab, HDV IgM and HDV IgG) represents in these cases an important tool to the clinician for the classification of the etiological agent, for the follow up of infected patients and their treatment. The detection of HDV total antibodies allows the classification of the illness and the monitoring of the seroconversion event.

C. PRINCIPLE OF THE TEST

Anti-HDV antibodies, if present in the sample, compete with a virus-specific polyclonal IgG, labeled with peroxidase (HRP), for a fixed amount of rec-HDV coated on the microplate. The test is carried out with a two steps incubation competitive system. First the sample is added to the plate and specific anti HDV antibodies bind to the adsorbed antigen. After washing, an enzyme conjugated antibody to HDV is added and binds to the free portion of the antigen coated. After washing a chromogen/substrate mixture is dispensed. The concentration of the bound enzyme on the solid phase becomes inversely proportional to the amount of anti-HDV antibodies in the sample and its activity is detected by the added chromogen/substrate. The concentration of HDV-specific antibodies in the sample is determined by means of a cut-off value that allows for the semi quantitative detection of anti-HDV antibodies.

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1. Microplate: **MICROPLATE**

8x12 microwell strips coated with recombinant HDV-specific antigen and sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening; reseal unused strips in the bag with desiccant and store at 4°C.

2. Negative Control: **CONTROL -**

1x2.0ml/vial. Ready to use. Contains goat serum proteins, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.09% Sodium Azide and 0.045% ProClin 300 as preservatives. The negative control is colour coded pale yellow.

3. Positive Control: **CONTROL +**

1x2.0ml/vial. Ready to use. Contains goat serum proteins, high titer anti HDV antibodies, 100 mM Tris-HCl buffer pH 7.4 +/-0.1, 0.09% Sodium Azide and 0.045% ProClin 300 as preservatives. The positive control is colour coded green.

4. Calibrator: **CAL**

n° 1 vial. Lyophilised. To be dissolved with EIA grade water as reported in the label. Contains bovine serum proteins, low titer human antibodies to HDV, 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

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6. Enzyme conjugate: **CONJ**

1x16ml/vial. Ready-to-use solution. Contains 5% bovine serum albumine, 10 mM tris buffer pH 6.8 +/-0.1, Horseradish peroxidase conjugated antibody to HDV in presence of 0.2 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives. The component is colour coded red.

7. Chromogen/Substrate: **SUBS TMB**

1x16ml/vial. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 4% DMSO, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide of H₂O₂.

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1. Calibrated Micropipettes in the range 10-1000 ul and disposable plastic tips.
2. EIA grade water (double distilled or deionized, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blinking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.

4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB/H₂O₂) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at +2..8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
12. The use of disposable plastic labware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
14. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

5. Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for at least 12 months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

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4. Calibrator:

Low positive control. Add precisely the volume of EIA grade water, reported on its label, to the lyophilized powder; let fully dissolve and then gently mix on vortex.

Note: *The dissolved calibrator is not stable. Store it frozen in aliquots at -20°C. When thawed do not freeze again; discard it.*

5. Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with EIA grade water up to 1200 ml and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

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1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing. An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to 4; (c) linearity to 4; repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the sections "Validation of Test" and "Assay Performances". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for

dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator as described above and gently mix.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

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1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at $+2.8^{\circ}\text{C}$, sealed.
2. Pipette 100 μl of Negative Control in triplicate, 100 μl Positive Control in single and then 100 μl of samples. Check that controls and samples have been correctly added. Then incubate the microplate at **$+37^{\circ}\text{C}$ for 60 min.**
3. Wash the microplate as reported in section I.3.
4. In all the wells except A1, pipette 100 μl Enzyme Conjugate. Check that the reagent has been correctly added. Then incubate the microplate at **$+37^{\circ}\text{C}$ for 60 min.**

Important note: Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.

5. Wash the microplate as described.

6. Pipette 100 µl TMB/H₂O₂ mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 min**.

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7. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step n° 6 to stop the enzymatic reaction. Addition of the stop solution will turn the negative control and negative samples from blue to yellow.

8. Measure the colour intensity of the solution in each well, as described in section I.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, mandatory), blanking the instrument on A1.

Important notes:

1. Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.
3. The use of the Calibrator, a low positive control, is not mandatory for the assay as the CAL does not enter into the cut-off calculation. The CAL may be used as a low titer positive control when a laboratory internal quality verification is required by the management. When used for such purpose, dispense 100 ul of it, possibly in duplicate.

N. ASSAY SCHEME

Controls/Calibrator	100 ul
Samples	100 ul
1st incubation	60 min
Temperature	+37°C
Washing step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H ₂ O ₂ mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

An example of dispensation scheme (including CAL) is reported in the table below:

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O. INTERNAL QUALITY CONTROL

A check is performed on the negative and positive controls any time, and on the Calibrator in addition when the kit is used for the first time, in order to verify whether the expected OD_{450nm} / 620-630nm or Co/S values have been matched in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.100 OD _{450nm} value
Negative Control (NC)	> 1.000 OD _{450nm} after blanking If lower carefully control the washing procedure and decrease the number of cycles or the soaking time coefficient of variation < 30%
Positive Control (PC)	OD _{450 nm} < NC/10
Calibrator (CAL)	PC ≤ OD _{450nm} < (NC+PC)/5

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If they don't, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD _{450nm}	that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) < 1.000 OD _{450nm} after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate; 6. that the washer needles are not blocked or partially obstructed.

Calibrator OD450nm Outside the range	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of Calibrator); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control OD450nm > NC/10	<ol style="list-style-type: none"> 1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Important note:

The analysis must be done proceeding as the reading step described in the section M, point 8.

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Q. INTERPRETATION OF RESULTS

Results are interpreted as ratio between the cut-off value and the sample OD450nm / 620-630nm or Co/S. Results are interpreted according to the following table:

Co/S	Interpretation
< 0.9	Negative
0.9 – 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HDV.

Any patient showing an equivocal result should be re-tested on a second sample taken 1-2 weeks after the initial sample.

A positive result is indicative of HDV infection and therefore the patient should be treated accordingly.

Important notes:

1. *Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.*
2. *When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.*
3. *Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.*

An example of calculation is reported below (data obtained proceeding as the the reading step described in the section M, point 8).

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Calibrator: 0.300-0.260 OD450nm

Mean value: 0.280 OD450nm

Within the range PC ≤ OD450nm < (NC+PC)/5 – Accepted

Sample 1: 0.020 OD450nm

Sample 2: 1.900 OD450nm

Sample 1 Co/S > 1.1 positive

Sample 2 Co/S < 0.9 negative

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC)

1. LIMIT OF DETECTION:

In absence of an international standard, the sensitivity of the assay has been calculated by means of the product named Accurun n° 127 supplied by Boston Biomedica Inc. – USA .

The table below reports the OD450nm shown by this preparation when diluted in Fetal Calf Serum to prepare a limiting dilution curve, in three different lots.

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The diagnostic performances were evaluated in a clinical trial conducted by the Department of Gastro-Hepatology, Prof. M.Rizzetto, S.Giovanni Battista hospital, Torino, Italy, on more than 400 samples against a reference kit.

Negative, positive and potentially interfering samples were examined in the trial.

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Results are briefly reported in the tables below:

Sensitivity	> 98 %
Specificity	> 98 %

3. PRECISION

The mean values obtained from a study conducted on two samples of different anti-HDV antibody reactivity, examined in 16 replicates in three separate runs for three lots of product, is reported below:

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Important note:

The performance data have been obtained proceeding as the reading step described in the section M, point 8.

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0318

HCV Ab

**Version 4.0 Enzyme Immunoassay
for the determination of
anti Hepatitis C Virus antibody
in human serum and plasma**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

REF CVAB.CE
96,192,480,960 Tests

HCV Ab

A. INTENDED USE

Version 4.0 Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis C Virus in human plasma and sera. The kit is intended for the screening of blood units and the follow-up of HCV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) define Hepatitis C infection as follows:

"Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post-transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. HCV is spread primarily by direct contact with human blood. The major causes of HCV infection worldwide are use of unscreened blood transfusions, and re-use of needles and syringes that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is too costly for most persons in developing countries to afford. Thus, from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviours (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the *flaviviridae* family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease.

An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high propensity (80%) of inducing chronic infection. HCV is clustered into several distinct genotypes which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Cirrhosis develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 5% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer who do not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with

alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, through the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug-users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioural practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils, or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used, the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or re-use unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (when prevalence data are available) is high compared to some countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme immunosorbent assays (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of active infection and a potential for spread of the infection and or/the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin, can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with ribavirin is effective in about 30% to 50% of patients. Ribavirin does not appear to be effective when used alone.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any protective immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HCV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation of plasma derived products; (c) implementation and maintenance of infection control practices in health care settings, including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce overuse of injections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices. “

The genome encodes for structural components, a nucleocapsid protein and two envelope glycoproteins, and functional constituents involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

C. PRINCIPLE OF THE TEST

Microplates are coated with HCV-specific antigens derived from “core” and “ns” regions encoding for conservative and immunodominant antigenic determinants (Core peptide, recombinant NS3, NS4 and NS5 peptides).

The solid phase is first treated with the diluted sample and HCV Ab are captured, if present, by the antigens.

After washing out all the other components of the sample, in the 2nd incubation bound HCV antibodies, IgG and IgM as well, are detected by the addition of polyclonal specific anti hlgG&M antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti HCV antibodies present in the sample. A cut-off value let optical densities be interpreted into HCV antibody negative and positive results.

D. COMPONENTS

Code CVAB.CE contains reagents for 192 tests.

1. Microplate **MICROPLATE**

n° 2 microplates

12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed into a bag with desiccant.

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3. Positive Control **CONTROL +**

1x4.0ml/vial. Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HCV, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. The Positive Control is blue colour coded.

4. Calibrator **CAL ...**

n° 2 vials. Lyophilized calibrator. To be dissolved with the volume of EIA grade water reported on the label. It contains foetal bovine serum proteins, human antibodies to HCV whose content is calibrated on the NIBSC Working Standard code 99/588-003-WI, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.3 mg/ml gentamicine sulphate and 0.045% ProClin 300 as preservatives.

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5. Wash buffer concentrate **WASHBUF 20X**

2x60ml/bottle. 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

6. Enzyme Conjugate **CONJ**

2x16ml/vial. Ready to use and pink/red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives.

7. Chromogen/Substrate **SUBS TMB**

2x16ml/vial. Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂.

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8. Assay Diluent **DILAS**

1x15ml/vial. 10 mM tris buffered solution pH 8.0 +/-0.1 containing 0.045% ProClin 300 for the pre-treatment of samples and controls in the plate, blocking interference.

9. Sulphuric Acid **H₂SO₄ 0.3 M**

1x32ml/bottle. It contains 0.3 M H₂SO₄ solution.

Attention: Irritant (H315; H319; P280; P302+P352; P332+P313; P305+P351+P338; P337+P313; P362+P363)

10. Sample Diluent: **DILSPE**

2x50ml/bottle. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 6.0 +/-0.1, 0.5% Tween 20, 0.09% Na-azide and 0.045% ProClin 300 as preservatives. To be used to dilute the sample.

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11. Plate sealing foils n° 4

12. Package insert n° 1

Important note: Only upon specific request , Dia.Pro can supply reagents for 96, 480, 960 tests , as reported below:

1. Microplate	n°1	n°5	n°10
2. Negative Control	1x2.0ml/vial	1x10ml/vial	1x20.ml/vial
3. Positive Control	1x2.0ml/vial	1x10ml/vial	1x20.ml/vial
4. Calibrator	n° 1 vial	n° 5 vials	n° 10 vials
5. Wash buff conc	1x60ml/bottle	5x60ml/bottles	4x150ml/bottles
6. Enz. Conjugate	1x16ml/vial	2x40ml/bottles	4x40ml/bottles
7. Chromog/Subs	1x16ml/vial	2x40ml/bottles	4x40ml/bottles
8. Assay Diluent	1x8ml/vial	1x40ml/bottle	1x80ml/bottle
9. Sulphuric Acid	1x15ml/vial	2x40ml/bottle	2x80ml/bottles
10. Sample Diluent	1x50ml/vial	5x50ml/bottles	4x125ml/bottles
11. Plate seal foils	n° 2	n° 10	n° 20
12. Pack. insert	n° 1	n° 1	n° 1
Number of tests	96	480	960
Code	CVAB.CE.96	CVAB.CE.480	CVAB.CE.960

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1. Calibrated Micropipettes (200ul and 10ul) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate from strong light and avoid vibration of the bench surface where the test is undertaken.
6. Upon receipt, store the kit at 2..8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
11. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated

before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

15. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water
17. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.
4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2°..+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.

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1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of manufacturing. In this case call Dia.Pro's customer service. Unused strips have to be placed back into the aluminium pouch, in presence of desiccant supplied, firmly zipped and stored at +2°..8°C. When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.

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4. Calibrator:

Dissolve carefully the content of the lyophilised vial with the volume of EIA grade water reported on its label.

Mix well on vortex before use.

Handle this component as potentially infective, even if HCV, eventually present in the control, has been chemically inactivated.

Note: *When dissolved the Calibrator is not stable. Store in aliquots at -20°C.*

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- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution. The instrument weekly has to be submitted to decontamination according to its manual (NaOH 0.1 M decontamination suggested).
5 washing cycles (aspiration + dispensation of 350ul/well of washing solution + 20 sec soaking = 1 cycle) are sufficient to ensure the assay with the declared performances. If soaking is not possible add one more cycle of washing. An incorrect washing cycle or salt-blocked needles are the major cause of false positive reactions.
- Incubation times have a tolerance of ±5%.
- The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to ≥ 2.0; (c) linearity to ≥ 2.0; (d) repeatability ≥ 1%. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended for blood screening when the number of samples to be tested exceed 20-30 units per run.
- When using automatic devices, in case the vial holder of the instrument does not fit with the vials supplied in the kit, transfer the solution into appropriate containers and label them with the same label peeled out from the original vial. This operation is important in order to avoid mismatching contents of vials, when transferring them. When the test is over, return the secondary labeled containers to 2..8°C, firmly capped.
- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a trueness of +/-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
- The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water

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- Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
- Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the

- aluminum pouch, containing the microplate, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
 - Dissolve the Calibrator as described above.
 - Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
 - Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.
 - Check that the ELISA reader has been turned on at least 20 minutes before reading.
 - If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
 - Check that the micropipettes are set to the required volume.
 - Check that all the other equipment is available and ready to use.
 - In case of problems, do not proceed further with the test and advise the supervisor.

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Automated assay:

In case the test is carried out automatically with an ELISA system, we suggest to make the instrument aspirate 200 ul Sample Diluent and then 10 ul sample.

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Do not dilute controls/calibrator as they are ready to use. Dispense 200 ul controls/calibrator in the appropriate control/calibration wells.

Important Note: *Visually monitor that samples have been diluted and dispensed into appropriate wells. This is simply achieved by checking that the colour of dispensed samples has turned to dark bluish-green while the colour of the negative control has remained olive green.*

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Important note: *Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.*

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- Wash the microplate with an automatic washer by delivering and aspirating 350ul/well of diluted washing solution as reported previously (section I.3).
- Pipette 100µl Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this pink/red coloured component has been dispensed in all the wells, except A1.

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- Incubate the microplate for **45 min at +37°C**.
- Wash microwells as in step 6.
- Pipette 100µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature (18-24°C) for 15 minutes**.

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- Pipette 100µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of acid will turn the positive control and positive samples from blue to yellow/brown.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction), blanking the instrument on A1 (mandatory).

Important notes:

- Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.
- Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.
- The Calibrator (CAL) does not affect the cut-off calculation and therefore the test results calculation. The Calibrator may be used only when a laboratory internal quality control is required by the management.

N. ASSAY SCHEME

Method	Operations
Controls & Calibrator	200 ul
Samples	200ul dil.+10ul
Assay Diluent (DILAS)	50 ul
1st incubation	45 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme conjugate	100 ul
2nd incubation	45 min
Temperature	+37°C
Wash step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H ₂ O ₂	100 ul
3rd incubation	15 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630nm

	samples, to spills or to the enzyme conjugate; 5. that micropipettes haven't got contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co < 1.1	1. that the procedure has been correctly executed; 2. that no mistake has been done in its distribution (ex.: dispensation of negative control instead of control serum) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly executed; 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150, too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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If the results of the test match the requirements stated above, proceed to the next section.
If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Sustrate solution has not got contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of their wells has occurred due to positive

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A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused.
Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.
A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies **earlier** than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result ! The sample has to be anyway submitted to a confirmation test (supplied upon request by DiaPro srl, code CCONF).
4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV antibodies - lacking anti hIgM conjugate in the formulation of the enzyme tracer and therefore missing IgM reactivity - may be present. The real positivity of the sample for antibodies to HCV should be then confirmed by examining also IgM reactivity, important for the diagnosis of HCV infection.
5. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
6. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.

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Negative Control: 0.019 – 0.020 – 0.021 OD450nm
 Mean Value: 0.020 OD450nm
 Lower than 0.050 – Accepted
 Positive Control: 2.189 OD450nm
 Higher than 1.000 – Accepted
 Cut-Off = 0.020+0.350 = 0.370
 Calibrator: 0.550 - 0.530 OD450nm
 Mean value: 0.540 OD450nm S/Co = 1.4
 S/Co higher than 1.1 – Accepted
 Sample 1: 0.070 OD450nm
 Sample 2: 1.690 OD450nm
 Sample 1 S/Co < 0.9 = negative
 Sample 2 S/Co > 1.1 = positive

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1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated by means of the British Working Standard for anti-HCV, NIBSC code 99/588-003-WI. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot # 1	Lot # 2
Factor	S/Co	S/Co
1 X	2.0	2.0
2 X	1.1	1.2
4 X	0.7	0.8
8 X	0.5	0.5
Negative plasma	0.3	0.3

In addition the sample coded Accurun 1 – series 3000 - supplied by Boston Biomedica Inc., USA, has been evaluated "in toto" showing the results below:

CVAB.CE Lot ID	Accurun 1 Series	S/Co
1201	3000	1.5
0602	3000	1.5
1202	3000	1.9

In addition, n° 7 samples, tested positive for HCV Ab with Ortho HCV 3.0 SAVE, code 930820, lot. # EXE065-1, were diluted in HCV Ab negative plasma in order to generate limiting dilutions and then tested again on CVAB.CE, lot. # 1202, and Ortho. The following table reports the data obtained.

Sample n°	Limit Dilution	CVAB.CE S/Co	Ortho 3.0 S/Co
1	256 X	1.9	1.3
2	256 X	1.9	0.7
3	256 X	2.4	1.0
4	128 X	2.5	3.2
5	85 X	3.3	1.4
6	128 X	2.2	0.8
7	135 X	3.2	2.2

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The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

2.1 Diagnostic specificity:

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. In addition to the first study, where a total of 5043 unselected blood donors, (including 1st time donors), 210 hospitalized patients and 162 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.) were examined, the diagnostic specificity was recently assessed by testing a total of 2876 negative blood donors on six different lots. A value of specificity of 100% was found.

No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage. No interference was observed.

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The diagnostic sensitivity has been assessed externally on a total number of 359 specimens; a diagnostic sensitivity of 100% was found. Internally more than other 50 positive samples were tested, providing a value of diagnostic sensitivity of again 100%. Positive samples from infections carried out by different genotypes of HCV were tested as well.

Furthermore, most of seroconversion panels available from Boston Biomedica Inc., USA, (PHV) and Zeptomatrix, USA, (HCV) have been studied.

Results are reported below for some of them.

Panel	N° samples	DiaPro*	Ortho* **
PHV 901	11	9	9
PHV 904	7	2	4
PHV 905	9	3	4
PHV 906	7	7	7
PHV 907	7	3	2
PHV 908	13	10	8
PHV 909	3	2	2
PHV 910	5	3	3
PHV 911	5	3	3
PHV 912	3	1	1
PHV 913	4	2	2
PHV 914	9	5	5
PHV 915	4	3	0
PHV 916	8	4	3
PHV 917	10	6	6
PHV 918	8	2	0
PHV 919	7	3	3
PHV 920	10	6	6
HCV 10039	5	2	0
HCV 6212	9	6	7
HCV 10165	9	5	4

Note: * Positive samples detected

** HCV v.3.0

Finally the Product has been tested on the panel EFS Ac HCV, lot n° 01/08.03.22C/01A, supplied by the Etablissement Francais Du Sang (EFS), France, with the following results:

EFS Panel Ac HCV

Sample	Lot # 1	Lot # 2	Lot # 2	Results expected
	S/Co	S/Co	S/Co	
HCV 1	2.2	2.4	2.6	positive
HCV 2	1.6	2.0	2.1	positive
HCV 3	1.5	1.7	1.6	positive
HCV 4	5.2	6.5	5.5	positive
HCV 5	1.6	1.8	1.6	positive
HCV 6	0.4	0.4	0.4	negative

3. PRECISION:

It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs. Results are reported as follows:

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S. LIMITATIONS

Repeatable false positive results, not confirmed by RIBA or similar confirmation techniques, were assessed as less than 0.1% of the normal population. Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

REFERENCES

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. MMWR 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. Hepatology 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. Viral Hepatitis and Liver Disease, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. Digestive diseases in the United States: epidemiology and impact. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. JAMA 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. Ann Intern Med 1972;77:691-9.

7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.

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CE
0318

HAV Ab

**Competitive Enzyme ImmunoAssay
(ELISA) for the determination of
antibodies to Hepatitis A Virus
in human plasma and sera**

- for "in vitro" diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via Carducci n° 27
20099 Sesto San Giovanni
(Milano) - Italy**

Phone +39 02 27007161

Fax +39 02 44386771

e-mail: info@diapro.it

HAV Ab

A. INTENDED USE

Competitive Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis A Virus in human plasma and sera. The kit is used for the follow-up of patients infected by HAV. For "in vitro" diagnostic use only.

B. INTRODUCTION

The Center for Disease Control or CDC, Atlanta, USA, defines Hepatitis A Virus as follows:

Hepatitis A continues to be one of the most frequently reported vaccine-preventable diseases in the world, despite the licensure of hepatitis A vaccine in 1995. Widespread vaccination of appropriate susceptible populations would substantially lower disease incidence and potentially eliminate indigenous transmission of hepatitis A virus (HAV) infection.

HAV, a 27-nm RNA agent classified as a picornavirus, can produce either asymptomatic or symptomatic infection in humans after an average incubation period of 28 days (range, 15-50 days). The illness caused by HAV infection typically has an abrupt onset of symptoms that can include fever, malaise, anorexia, nausea, abdominal discomfort, dark urine, and jaundice. The likelihood of having symptoms with HAV infection is related to the person's age. In children less than 6 years of age, most (70%) infections are asymptomatic; if illness does occur, it is not usually accompanied by jaundice. Among older children and adults, infection is usually symptomatic, with jaundice occurring in greater than 70% of patients. Signs and symptoms usually last less than 2 months, although 10%-15% of symptomatic persons have prolonged or relapsing disease lasting up to 6 months.

In infected persons, HAV replicates in the liver, is excreted in bile, and is shed in the stool. Peak infectivity of infected persons occurs during the 2-week period before onset of jaundice or elevation of liver enzymes, when the concentration of virus in stool is highest. The concentration of virus in stool declines after jaundice appears. Children and infants can shed HAV for longer periods than adults, up to several months after the onset of clinical illness. Chronic shedding of HAV in feces does not occur; however, shedding can occur in persons who have relapsing illness. Viremia occurs soon after infection and persists through the period of liver enzyme elevation.


Hepatitis A cannot be differentiated from other types of viral hepatitis on the basis of clinical or epidemiologic features alone. Serologic testing to detect immunoglobulin M (IgM) antibody to the capsid proteins of HAV (IgM anti-HAV) is required to confirm a diagnosis of acute HAV infection. In most persons, IgM anti-HAV becomes detectable 5-10 days before the onset of symptoms and can persist for up to 6 months after infection. Immunoglobulin G (IgG) anti-HAV, which appears early in the course of infection, remains detectable for the person's lifetime and confers lifelong protection against the disease. Commercial diagnostic tests are available for the detection of IgM and total (IgM and IgG) anti-HAV in serum.

HAV RNA can be detected in the blood and stool of most persons during the acute phase of infection by using nucleic acid amplification methods, and nucleic acid sequencing has been used to determine the relatedness of HAV isolates.

HAV infection is acquired primarily by the fecal-oral route by either person-to-person contact or ingestion of contaminated food or water. On rare occasions, HAV infection has been transmitted by transfusion of blood or blood products collected from donors during the viremic phase of their infection. In experimentally infected nonhuman primates, HAV has been detected in saliva during the incubation period; however, transmission by saliva has not been demonstrated.

Depending on conditions, HAV can be stable in the environment for months. Heating foods at temperatures greater than 185 F (85 C) for 1 minute or disinfecting surfaces with a 1:100 dilution of sodium hypochlorite (i.e., household bleach) in tap water is necessary to inactivate HAV.

Because most children have asymptomatic or unrecognized infections, they play an important role in HAV transmission and serve as a source of infection for others. In one study of adults without an identified source of infection, 52% of their households included a child less than 6 years old, and the presence of a young child was associated with HAV transmission within the household. In studies where serologic testing of the household contacts of adults without an identified source of infection was performed, 25%-40% of the contacts less than 6 years old had serologic evidence of acute HAV infection (IgM anti-HAV).

C. PRINCIPLE OF THE TEST

The assay is based on the principle of competition where the antibodies in the sample compete with an anti-HAV specific antibody, labeled with HRP, for a fixed amount of antigen on the solid phase.

A purified and inactivated HAV is coated to the microwells.

The patient's serum/plasma is added to the microwell and antibodies to HAV are captured by the solid phase.

After washing, the enzyme conjugate is added and binds to the free HAV antigen, if still present.

The plate is washed to remove unbound conjugate and then the chromogen/substrate is added.

In the presence of peroxidase the colorless substrate is hydrolysed to a coloured end-product, whose optical density may be detected and is inversely proportional to the amount of antibodies to HAV present in the sample.

An additive is added to the sample directly into the well to block interferences able to mask the presence of antibodies, mostly appearing in the follow up of vaccination.

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1. Microplate **MICROPLATE**

8x12 microwell strips coated with purified and inactivated HAV, sealed into a bag with desiccant. Allow the microplate to reach room temperature before opening. Reseal unused strips in the bag with desiccant and store at 2..8°C°.

2. Negative Control: **CONTROL -**

1x4.0ml/vial. Ready to use. Contains bovine serum proteins, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives. The negative control is color coded pale yellow.

3. Positive Control: **CONTROL +**

1x4.0ml/vial. Ready to use. Contains bovine serum proteins, anti HAV antibodies at a concentration higher than 100 WHO mIU/ml, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives. The positive control is colour coded green.

4. Calibrator: **CAL ...**

n° 1 vial. Lyophilized. To be dissolved with EIA grade water as reported in the label. Contains bovine serum proteins, anti HAV antibodies at a concentration of about 10 WHO mIU/ml, 10 mM phosphate buffer pH 7.4+/-0.1, 0.02% gentamicine sulphate and 0.045% ProClin 300 as preservatives.

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5. Wash buffer concentrate: **WASHBUF 20X**

1x60ml/bottle. 20x concentrated solution. to be diluted up to 1200ml with distilled water before use. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.045% ProClin 300.

6. Enzyme conjugate: **CONJ**

1x16ml/vial. Ready-to-use solution. Contains Horseradish peroxidase conjugated antibody, specific to HAV, in presence of 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.045% ProClin 300 and 0.02% gentamicine sulphate as preservatives. The reagent is colored with a red dye.

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8. Specimen Diluent: **DILSPE**

1x8ml. Buffered solution suggested to be used in the follow up of vaccination. It contains 0.09% sodium azide and 0.045% ProClin 300 as preservatives. The reagent is color coded dark green.

9. Sulphuric Acid: **H₂SO₄ 0.3 M**

1x15ml/vial. Contains 0.3 M H₂SO₄ solution. Attention: Irritant (H315; H319; P280; P302+P352; P332+P313; P305+P351+P338; P337+P313; P362+P363)

10. Plate sealing foils n° 2

11. Package insert n° 1

Upon request:

Calibration Curve: **CAL N° ...**

5x2.0 ml/vial. Ready to use and colour coded standard curve ranging: 0-5-10-50-100 WHO mIU/ml.

(CAL1=0mIU/ml, CAL2=5mIU/ML, CAL3=10mIU/ml, CAL4=50mIU/ml, CAL5=100mIU/ml).

Contains serum proteins, 0.3 mg/ml gentamicine sulphate and 0.045% ProClin 300 GC as preservatives.

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1. Calibrated Micropipettes (150ul, 100ul and 50ul) and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.1°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample. Do not reuse disposable tips.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one. Do not reuse disposable tips.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
11. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.

12. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..
13. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
14. The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water
15. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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- Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- Avoid any addition of preservatives; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
- Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and reading is strongly recommended.
- Haemolysed and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- Sera and plasma can be stored at +2°...+8°C in primary collection tubes for up to five days after collection. Do not freeze primary tubes of collection. For longer storage periods, sera and plasma samples, carefully removed from the primary collection tube, can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

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A study conducted on opened kit has pointed out no relevant loss of performances up to 3 months from first opening.

1. Antigen coated microwells:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in conservation. In this case, call Dia.Pro's customer service. Unused strips have to be placed back into the aluminum pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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4. Calibrator:

Add the volume of ELISA grade water, reported on the label, to the lyophilized powder; let fully dissolve and then gently mix on vortex. The dissolved calibrator is not stable; store it frozen in aliquots at -20°C.

5. Wash buffer concentrate:

The whole content of the 20x concentrated solution has to be diluted with bi-distilled water up to 1200ml and mixed gently end-over-end before use.

Once diluted, the wash solution is stable for 1 week at 2-8° C.

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- Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.
- The ELISA incubator has to be set at +37°C (tolerance of +/- 0.5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The **ELISA washer** is extremely important to the overall performances of the assay. The washer must be carefully validated in advance, checked for the delivery of the right

dispensation volume and regularly submitted to maintenance according to the manufacturer's instructions for use. In particular the washer, at the end of the daily workload, has to be extensively cleaned out of salts with deionized water. Before use, the washer has to be extensively primed with the diluted Washing Solution.

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4. Incubation times have a tolerance of $\pm 5\%$.
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter of 620-630nm, mandatory for blanking purposes. Its standard performances should be (a) bandwidth ≤ 10 nm; (b) absorbance range from 0 to 4; (c) linearity to 4; (d) repeatability $\geq 1\%$. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB+H₂O₂) is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Calibrator as described above and gently mix.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
6. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as reported in the specific section.

7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

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1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at +2..8°C, sealed.

2. Dispense 50 ul Specimen Diluent in all the wells identified for samples and controls/calibrator, except for A1. Then pipette 100 ul of Negative Control in triplicate, 100 ul of Calibrator in duplicate, 100 ul Positive Control in single and then 100 ul of samples.

Check that controls/calibrator and samples have been correctly added. Incubate the microplate at **+37°C for 60 min.**

3. Wash the microplate as reported in section I.3.

4. In all the wells except A1, pipette 100 ul Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at **+37°C for 60 minutes.**

Important note: *Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.*

5. Wash the microplate as described.

6. Pipette 100 ul TMB/H₂O₂ mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at **room temperature for 20 minutes.**

Important note: *Do not expose to strong direct light as a high background might be generated.*

7. Pipette 100 ul Sulphuric Acid into each well to stop the enzymatic reaction using the same pipetting sequence as in step 6. Then measure the color intensity with a microplate reader at 450nm (reading) and at 620-630nm (background subtraction), blanking the instrument on A1 well (mandatory).

Important notes:

1. *Ensure that no finger prints are present on the bottom of the microwell before reading. Finger prints could generate false positive results on reading.*
2. *Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.*

N. ASSAY SCHEME

Specimen Diluent	50 ul
Controls/Calibrator(*)	100 ul
Samples	100 ul
1st incubation	60 min
Temperature	+37°C
Washing step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
Enzyme Conjugate	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing step	n° 5 cycles with 20" of soaking OR n° 6 cycles without soaking
TMB/H ₂ O ₂ mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm / 620-630

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O. INTERNAL QUALITY CONTROL

A check is performed on the controls any time the kit is used in order to verify whether the expected OD450nm or Co/S values have been matched in the analysis.

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC)	1. that the washing procedure and the washer settings are as validated in the

< 0.750 OD450nm after blanking coefficient of variation > 30%	pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control > 0.300 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (ex.: dispensation of negative control instead of the positive one); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Problem	Check
Calibrator Co/S < 1.0	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

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P. CALCULATION OF THE CUT-OFF

The test results are calculated by means of a cut-off value determined with the following formula:

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Important note: When the calculation of results is performed by the operating system of an ELISA automated work station, make sure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as a ratio of the Cut-Off value and the OD450nm of the sample (or Co/S) according to the following table:

Co/S	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HAV.

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A positive result is indicative of a past or recent HAV infection and therefore the patient should be treated accordingly.

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Positive Control: 0.100 OD450nm
 Lower than 0.300 – Accepted

Cut-Off = (2.000 + 0.100) / 3 = 0.700

Calibrator: 0.400-0.360 OD450nm
 Mean value: 0.380 OD450nm
 Co/S > 1 – Accepted

Sample 1: 0.050 OD450nm
 Sample 2: 1.900 OD450nm
 Sample 1 Co/S > 1.1 positive
 Sample 2 Co/ < 0.9 negative

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R. PERFORMANCE CHARACTERISTICS

1. Limit of detection

The limit of detection of the assay has been calculated by means of the 2nd International Standard supplied by WHO.


Two control samples, supplied by Boston Biomedica Inc., USA, with code Accurun 52 and 120, were also examined.

The sensitivity shown by the assay is < 10 WHO mIU/ml or < 5 PEI mU/ml.

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Curves are reported below:

2. Diagnostic sensitivity:

The diagnostic sensitivity has been tested in a clinical trial on panels of samples classified positive by a US FDA approved kit. An overall value of 100% has been found in the study conducted on a total number of more than 200 samples.

Seroconversion and performance panels have also been studied. Results obtained by examining two panels supplied by Boston Biomedica Inc., USA, are reported below.

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3. Diagnostic specificity:

The diagnostic specificity has been determined on panels of negative samples from normal individuals and blood donors, classified negative with a US FDA approved kit.

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Samples derived from patients with different viral (HCV, HDV, HBV, HIV) and non viral pathologies of the liver that may interfere with the test were examined.

No cross reaction were observed.

The Performance Evaluation study conducted in the external reference center on more than 1000 samples has provided a value > 98% .

4. Precision

The mean values obtained from a study conducted on two samples of different anti-HAV reactivity, examined in 16 replicates in three separate runs is reported below:

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REFERENCES

- CDC. Summary of notifiable diseases, United States, 1997. *MMWR* 1998;46:1-87.
- CDC. Prevention of hepatitis A through active or passive immunization. Recommendations of the Advisory Committee on Immunization Practices (ACIP). *MMWR* 1996;45(RR-15).
- Krugman S, Giles JP. Viral hepatitis: new light on an old disease. *JAMA* 1970;212:1019-29.
- Hadler SC, Webster HM, Erben JJ, Swanson JE, Maynard JE. Hepatitis A in day-care centers: a communitywide assessment. *N Engl J Med* 1980;302:1222-7.
- Lednar WM, Lemon SM, Kirkpatrick JW, Redfield RR, Fields ML, Kelley PW. Frequency of illness associated with epidemic hepatitis A virus infection in adults. *Am J Epidemiol* 1985;122:226-33.
- Glikson M, Galun E, Oren R, Tur-Kaspa R, Shouval D. Relapsing hepatitis A. Review of 14 cases and literature survey. *Medicine* 1992;71:14-23.
- Skinh j P, Mathiesen LR, Kryger P, M Iler AM. Faecal excretion of hepatitis A virus in patients with symptomatic hepatitis A infection. *Scand J Gastroenterol* 1981;16:1057-9.
- Tassopoulos NC, Papaevangelou GJ, Ticehurst JR, Purcell RH. Faecal excretion of Greek strains of hepatitis A virus in patients with hepatitis A and in experimentally infected chimpanzees. *J Infect Dis* 1986;154:231-7.
- Rosenblum LS, Villarino ME, Nainan OV, et al. Hepatitis A outbreak in a neonatal intensive care unit: risk factors for transmission and evidence of prolonged viral excretion among preterm infants. *J Infect Dis* 1991;164:476-82.
- Sjogren MH, Tanno H, Fay O, et al. Hepatitis A virus in stool during clinical relapse. *Ann Intern Med* 1987;106:221-6.
- Lemon SM. The natural history of hepatitis A: the potential for transmission by transfusion of blood or blood products. *Vox Sang* 1994;67(suppl 4):19-23.
- Bower WA, Nainan OV, Margolis HS. Duration of viremia in naturally-acquired hepatitis A viral infections. [Abstract 103] In: Abstracts of the Infectious Diseases Society of America 35th Annual Meeting. Alexandria, VA: Infectious Diseases Society of America, 1997.
- Liaw YF, Yang CY, Chu CM, Huang MJ. Appearance and persistence of hepatitis A IgM antibody in acute clinical hepatitis A observed in an outbreak. *Infection* 1986;14:156-8.
- Stapleton JT. Host immune response to hepatitis A virus. *J Infect Dis* 1995;171(suppl 1):S9-14.
- Hutin YJF, Pool V, Cramer EH, et al. A multistate, foodborne outbreak of hepatitis A. *N Engl J Med* 1999;340:595-602.
- Soucie JM, Robertson BH, Bell BP, McCaustland KA, Evatt BL. Hepatitis A virus infections associated with clotting factor concentrate in the United States. *Transfusion* 1998;38:573-9.
- Cohen JI, Feinstone S, Purcell RH. Hepatitis A virus infection in a chimpanzee: duration of viremia and detection of virus in saliva and throat swabs. *J Infect Dis* 1989;160:887-90.
- McCaustland KA, Bond WW, Bradley DW, Ebert JW, Maynard JE. Survival of hepatitis A virus in feces after drying and storage for 1 month. *J Clin Microbiol* 1982;16:957-8.
- Favero MS, Bond WW. Disinfection and sterilization. In: Zuckerman AJ, Thomas HC, eds. *Viral hepatitis, scientific basis and clinical management*. New York, NY: Churchill Livingstone, 1993:565-75.
- Staes C, Schlenker T, Risk I, et. al. Source of infection among persons with acute hepatitis A and no identified risk factors, Salt Lake County, Utah, 1996 [Abstract 302]. *Clin Infect Dis* 1997;25:411.

21. Smith PF, Grabau JC, Werzberger A, et al. The role of young children in a community-wide outbreak of hepatitis A. *Epidemiol Infect* 1997;118:243-52.
22. Williams I, Bell B, Kaluba J, Shapiro C. Association between chronic liver disease and death from hepatitis A, United States, 1989-92 [Abstract A39]. IX Triennial International Symposium on Viral Hepatitis and Liver Disease. Rome, Italy, April 1996.
23. Akriadias EA, Redeker AG. Fulminant hepatitis A in intravenous drug users with chronic liver disease. *Ann Intern Med* 1989;110:838-9.
24. Willner IR, Uhl MD, Howard SC, Williams EQ, Riely CA, Waters B. Serious hepatitis A: an analysis of patients hospitalized during an urban epidemic in the United States. *Ann Intern Med* 1998;128:111-4.

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Manufactured by
Dia.Pro Diagnostic Bioprobes S.r.l.
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy

Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HAV Ab CODE: AVAB.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – SEPTEMBER 2003
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

HAV IgM

**“Capture” Enzyme Immuno Assay (ELISA)
for the determination of IgM class
antibodies to Hepatitis A Virus
in human plasma and sera**

- for “in vitro” diagnostic use only -

DIA.PRO

**Diagnostic Bioprobes Srl
Via G. Carducci n° 27
20099 Sesto San Giovanni
(Milano) – Italy**

Phone +39 02 27007161

Fax +39 02 26007726

e-mail: info@diapro.it

HAV IgM

A. INTENDED USE

Enzyme ImmunoAssay (ELISA) for the determination of IgM class antibodies to Hepatitis A Virus in human plasma and sera with the "capture" system. The kit may be used for the identification of the viral agent causing hepatitis in the patient and the follow up of the acute phase of the infection.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The Center for Disease Control or CDC, Atlanta, USA, defines Hepatitis A Virus as follows:

Hepatitis A continues to be one of the most frequently reported vaccine-preventable diseases in the world, despite the licensure of hepatitis A vaccine in 1995. Widespread vaccination of appropriate susceptible populations would substantially lower disease incidence and potentially eliminate indigenous transmission of hepatitis A virus (HAV) infection.

HAV, a 27-nm RNA agent classified as a picornavirus, can produce either asymptomatic or symptomatic infection in humans after an average incubation period of 28 days (range, 15-50 days). The illness caused by HAV infection typically has an abrupt onset of symptoms that can include fever, malaise, anorexia, nausea, abdominal discomfort, dark urine, and jaundice. The likelihood of having symptoms with HAV infection is related to the person's age. In children less than 6 years of age, most (70%) infections are asymptomatic; if illness does occur, it is not usually accompanied by jaundice. Among older children and adults, infection is usually symptomatic, with jaundice occurring in greater than 70% of patients. Signs and symptoms usually last less than 2 months, although 10%-15% of symptomatic persons have prolonged or relapsing disease lasting up to 6 months. In infected persons, HAV replicates in the liver, is excreted in bile, and is shed in the stool. Peak infectivity of infected persons occurs during the 2-week period before onset of jaundice or elevation of liver enzymes, when the concentration of virus in stool is highest. The concentration of virus in stool declines after jaundice appears. Children and infants can shed HAV for longer periods than adults, up to several months after the onset of clinical illness. Chronic shedding of HAV in feces does not occur; however, shedding can occur in persons who have relapsing illness. Viremia occurs soon after infection and persists through the period of liver enzyme elevation.


Hepatitis A cannot be differentiated from other types of viral hepatitis on the basis of clinical or epidemiologic features alone. Serologic testing to detect immunoglobulin M (IgM) antibody to the capsid proteins of HAV (IgM anti-HAV) is required to confirm a diagnosis of acute HAV infection. In most persons, IgM anti-HAV becomes detectable 5-10 days before the onset of symptoms and can persist for up to 6 months after infection. Immunoglobulin G (IgG) anti-HAV, which appears early in the course of infection, remains detectable for the person's lifetime and confers lifelong protection against the disease. Commercial diagnostic tests are available for the detection of IgM and total (IgM and IgG) anti-HAV in serum.

HAV RNA can be detected in the blood and stool of most persons during the acute phase of infection by using nucleic acid amplification methods, and nucleic acid sequencing has been used to determine the relatedness of HAV isolates.

HAV infection is acquired primarily by the fecal-oral route by either person-to-person contact or ingestion of contaminated food or water. On rare occasions, HAV infection has been transmitted by transfusion of blood or blood products collected from donors during the viremic phase of their infection. In experimentally infected nonhuman primates, HAV has been detected in saliva during the incubation period; however, transmission by saliva has not been demonstrated.

Depending on conditions, HAV can be stable in the environment for months. Heating foods at temperatures greater than 185 F (85 C) for 1 minute or disinfecting surfaces with a 1:100 dilution of sodium hypochlorite (i.e., household bleach) in tap water is necessary to inactivate HAV.

Because most children have asymptomatic or unrecognized infections, they play an important role in HAV transmission and serve as a source of infection for others. In one study of adults without an identified source of infection, 52% of their households included a child less than 6 years old, and the presence of a young child was associated with HAV transmission within the household. In studies where serologic testing of the household contacts of adults without an identified source of infection was performed, 25%-40% of the contacts less than 6 years old had serologic evidence of acute HAV infection (IgM anti-HAV).

C. PRINCIPLE OF THE TEST

The assay is based on the principle of "IgM capture" where IgM class antibodies in the sample are first captured by the solid phase coated with anti-IgM antibody.

After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a purified preparation of inactivated HAV, labelled with an antibody conjugated with peroxidase (HRP).

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In the presence of peroxidase the colorless substrate is hydrolysed to a colored end-product, whose optical density may be detected and is proportional to the amount of antibodies to HAV present in the sample.

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1. Microplate: MICROPLATE

12 strips of 8 breakable wells coated with anti human IgM antibody, affinity purified, and sealed into a bag with desiccant. Bring the microplate to room temperature before opening the bag. Unused strips have to be returned into the bag and the bag has to be sealed and stored back to 2..8°C, in presence of the desiccant.

2. Negative Control: CONTROL -

1x4.0 ml/vial. Ready to use control. It contains goat serum proteins, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The negative control is colourless.

3. Positive Control: CONTROL +

1x4.0 ml/vial. Ready to use control. It contains anti HAV IgM, goat serum proteins, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The positive control is green colour coded.

4. Calibrator: CAL ...

N° 1 lyophilized vial. To be dissolved with EIA grade water as reported in the label. It contains anti HAV IgM, 2% BSA, 10 mM tris buffer pH 6.0+/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives.

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5. Wash buffer concentrate: WASHBUF 20X

1x60ml/bottle. 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.

6. Enzyme conjugate 20X: CONJ

1x0.8 ml/vial. 20X concentrated solution. It contains Horseradish peroxidase conjugated antibody specific to HAV in presence of 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives.

7. HAV Antigen: Ag HAV

1x16 ml/vial. Ready-to-use solution. It contains inactivated and stabilised HAV in presence of 10 mM Tris buffer pH 6.8+/-0.1, 2% BSA, 0.1% Kathon GC and 0.02% gentamicine sulphate as preservatives. The reagent is red colour coded.

8. Specimen Diluent: DILSPE

2x60.0 ml/vial. Proteic buffered solution for the dilution of samples. It contains goat serum proteins, 10 mM tris buffer pH 6.0+/-0.1, 0.1% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives. The reagent is blue colour coded.

9. Chromogen/Substrate: SUBS TMB

1x16ml/vial. It contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.8, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide of H₂O₂.

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1. Calibrated Micropipettes of 10ul, 100ul and 1000ul and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.1°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB & H₂O₂) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
11. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
12. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and

inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min..

13. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.

14. The Stop Solution is an irritant. In case of spills, wash the surface with plenty of water

15. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

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2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and reading is strongly recommended.

3. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.

4. Sera and plasma can be stored at +2°...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be freeze/thawed more than once as IgM antibodies may get damaged and as this procedure may generate particles that could affect the test result.

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A study conducted on an opened kit has not pointed out any relevant loss of activity up to 3 months.

1. Antibody coated microwells:

Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in conservation. In this case, call Dia.Pro's customer service.

Unused strips have to be placed back into the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2°-8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

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3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious, even if HAV, eventually present in the control, has been chemically inactivated.

4. Calibrator:

Add the volume of ELISA grade water, reported on the label, to the lyophilised powder; let fully dissolve and then gently mix on vortex. The solution is not stable. Store the Calibrator frozen in aliquots at -20°C.

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5. Wash buffer concentrate:

The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before use.

Once diluted, the wash solution is stable for 1 week at 2-8° C. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

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6. Enzyme conjugate:

20X preparation. Mix well on vortex.

Avoid contamination of the liquid with oxidizing chemicals, dust or microbes when the reagent is aspirated to be used.

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Handle this component as potentially infectious, even if HAV has been chemically inactivated.

6+7. HAV Antigen/Antibody complex:

About 5-10 min before its use, dilute the 20X concentrated Enzyme Conjugate in the proper volume of HAV Antigen, necessary for the assay. Then mix on vortex carefully.

Example: To run 2 strips, dilute 100 µl Enzyme Conjugate 20X into 2 ml of HAV Antigen.

Note: This immunocomplex is not stable; discard the exceeding volume.

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Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong light, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

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1. Micropipettes have to be calibrated to deliver the correct volume (tolerance +/-5%) required by the assay and must be submitted to regular decontamination (household

alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$.

- The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
- The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests. 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section 'O'. Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
- Incubation times have a tolerance of $\pm 5\%$.
- The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm, strongly recommended) for blanking purposes. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section 'O' "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.
- Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

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- Check the expiration date of the kit printed on the external label (primary container). Do not use the device if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.

- Dissolve the Calibrator as described above and gently mix.
- Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at $+37^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the micropipettes are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

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- Dilute samples 1:101 by dispensing first 10 μl sample and then 1 ml Sample Diluent into a dilution tube; mix gently on vortex.
- Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
- Dispense 100 μl Negative Control in triplicate, 100 μl Positive Control in single and 100 μl Calibrator in duplicate in proper wells. Do not dilute controls and the calibrator as they are ready to use !
- Dispense 100 μl diluted samples in the proper sample wells and then check that all the samples wells are blue coloured and that controls and calibrator have been dispensed.
- Incubate the microplate for **60 min at $+37^{\circ}\text{C}$** .

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- About 5-10 minutes before use, prepare the HAV Antigen/Antibody immunocomplex as described previously.
- Wash the microplate with an automatic washer as reported previously (section I.3).
- Pipette 100 μl HAV Antigen/Antibody complex into each well, except the 1st blanking well, and cover with the sealer. Check that all wells are red coloured, except A1.

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- Incubate the microplate for **60 min at $+37^{\circ}\text{C}$** .
- Wash microwells as in step 7.
- Pipette 100 μl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at **room temperature ($18-24^{\circ}\text{C}$) for 20 minutes**.

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- Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction using the same pipetting sequence as in step 10. Addition of acid will turn the positive control and positive samples from blue to yellow.
- Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and at 620-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

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N. ASSAY SCHEME

Controls&Calibrator (*) samples diluted 1:101	100 ul
1st incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
HAV & Tracer	100 ul
2nd incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
TMB/H ₂ O ₂ mix	100 ul
3rd incubation	20 min
Temperature	r.t.
Sulphuric Acid	100 ul
Reading OD	450nm

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O. INTERNAL QUALITY CONTROL

A check is performed on the controls any time the kit is used in order to verify whether the expected OD450nm or S/Co values have been matched in the analysis.

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Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.150 OD450nm after blanking coefficient of variation > 30%	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Control < 0.500 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

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Problem	Check
Calibrator S/Co < 1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (e.g.: dispensation of negative control instead) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.

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A negative result indicates that the patient is not undergoing an acute infection by HAV.

Any patient showing an equivocal result, should be re-tested by examining a second sample after 1-2 weeks from first testing.

A positive result is indicative of an HAV infection event and therefore the patient should be treated accordingly.

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Sample 2 S/Co > 1.2 = positive

Important notes:

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method (confirmation test) before a diagnosis of viral hepatitis is confirmed.
3. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

R. PERFORMANCE CHARACTERISTICS

1. Limit of detection

In absence of a defined international standard for HAV IgM, the limit of detection of the assay has been calculated by means of the following preparations:

1. Accurun # 121 supplied by Boston Biomedica Inc. – USA
2. Accurun # 51 supplied by Boston Biomedica Inc., USA

These preparation were prepared according to the manufacturer's instructions, diluted in Sample Diluent (1:100) and then further diluted in Sample Diluent to generate a limiting curve (accurun # 121).

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Positive samples were collected from patients carrying HAV acute infection, confirmed by clinical symptoms and analysis.

An overall value of 100% has been found in the study conducted on a total number of more than 100 samples.

A seroconversion panel has also been studied.

Results obtained by examining a preparation supplied by Boston Biomedica Inc., USA, are reported below.

Seroconversion Panel : PHT 902

Sample	OD450nm	S/Co	DiaSorin Refer.	
			S/Co	Score
CTRL (-)	0,048	0,2		
CTRL (+)	1,736	5,8		
PHT902				
1	0,037	0,1	0,3	neg
2	0,042	0,1	0,3	neg
3	1,956	6,6	6,8	pos
4	1,988	6,7	6,7	pos
5	0,669	2,2	1,5	pos

3. Diagnostic Specificity:

The diagnostic specificity has been determined on panels of specimens, negative with the reference kit, derived from normal individuals and blood donors of European origin.

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Samples derived from patients with different viral (HCV, HDV, HBV, HEV) and non viral pathologies of the liver that may interfere with the test were examined.

No cross reaction were observed.

The Performance Evaluation study conducted in a qualified external reference centre on more than 500 samples has provided a value > 98%.

3. Precision:

It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs. Results are reported as follows:

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S. LIMITATIONS

False positivity has been assessed as less than 2% of the normal population, mostly due to high titers of RF.

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Dia.Pro
Diagnostic
Bio**Probes**

EC DECLARATION OF CONFORMITY

MANUFACTURER	DIA.PRO DIAGNOSTIC BIOPROBES S.R.L. VIA G. CARDUCCI N° 27 – 20099 SESTO SAN GIOVANNI (MILANO) – ITALY
PRODUCT	HAV IgM CODE: AVM.CE (96 tests)
CLASSIFICATION	GENERAL IVD
CONFORMITY ASSESSMENT ROUTE	SELF CERTIFICATION

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PLACE & DATE OF FIRST ISSUE	MILANO – SEPTEMBER 2003
PLACE & DATE OF CURRENT ISSUE	SESTO SAN GIOVANNI (MI) – MARCH 2019
SIGNATURE Legal Representative Dr.ssa Fiorenza Scozzesi	

Rev: 05/2018

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

EC Certificate - Full Quality Assurance System No. 11 0040 QS/NB

The quality system of manufacturer

Federal Budget Institute of Science “Central Research Institute for Epidemiology”

3a Novogireevskaya Street, Moscow 111123, Russia

has been certified as meeting the requirements of

Directive 98/79/EC

on in vitro diagnostic medical devices, Annex IV excluding (4, 6)

for the following product category(ies):

AmpliSens® PCR kits

The Notified Body No. 1023 declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with IVDD Annex IV. This quality assurance system conforms to the requirements of this Directive and is subjected to periodical surveillance. For placing on the market of List A devices covered by this certificate, an EC Design-Examination Certificate according to Annex IV (Section 4) is required.

Valid from: 2022-04-28
Valid until: 2025-05-26
First Issued: 2011-01-24
Revision: k

Date: 2022-04-28

A handwritten signature in blue ink, appearing to read 'Jiri Heš'.

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”**
3a Novogireevskaya Street, Moscow 111123, Russia

Product(s):

Name: **AmpliSens® Rubella virus-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-50 F

Classification: List B

GMDN: 30793

Name: **AmpliSens® Toxoplasma gondii-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-50 F

Classification: List B

GMDN: 52428

Name: **AmpliSens® CMV-FEP PCR kit**

Trade name(s): -

Model(s): variant FEP (0.2-ml tubes)

Classification: List B

GMDN: 30798

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB
issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”**
3a Novogireevskaya Street, Moscow 111123, Russia

Name: **AmpliSens® CMV-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 30798

Name: **AmpliSens® HSV / CMV-MULTIPRIME-FRT
PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 61348

Name: **AmpliSens® CMV-screen/monitor-FRT PCR
kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 30798

Date: 2022-04-28

Revision: k

Mgr. Jiří Heš

Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Name: **AmpliSens® EBV / CMV / HHV6-screen-FRT
PCR kit**

Trade name(s): -
Model(s): variant FRT-100 F
Classification: List B
GMDN: 61348

Name: **AmpliSens® Chlamydia trachomatis-FEP
PCR kit**

Trade name(s): -
Model(s): variant FEP (0.2-ml tubes)
Classification: List B
GMDN: 30677

Name: **AmpliSens® Chlamydia trachomatis-FRT
PCR kit**

Trade name(s): -
Model(s): variant FRT, variant FRT-100 F
Classification: List B
GMDN: 30677

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš

Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Name: **AmpliSens® *C.trachomatis* / *Ureaplasma* /
M.genitalium-MULTIPRIME-FEP PCR kit**

Trade name(s): -

Model(s): variant FEP (0.2 ml tubes)

Classification: List B

GMDN: 50409

Name: **AmpliSens® *C.trachomatis* / *Ureaplasma* /
M.genitalium-MULTIPRIME-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 50409

Name: **AmpliSens® *C.trachomatis* / *Ureaplasma* /
M.hominis-MULTIPRIME-FEP PCR kit**

Trade name(s): -

Model(s): variant FEP (0.2 ml tubes)

Classification: List B

GMDN: 50409

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Name: **AmpliSens® *C.trachomatis* / *Ureaplasma* /
M.hominis-MULTIPRIME-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 50409

Name: **AmpliSens® *C.trachomatis* / *Ureaplasma* /
M.genitalium / *M.hominis*-MULTIPRIME-FRT
PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 50409

Name: **AmpliSens® *N.gonorrhoeae* / *C.trachomatis* /
M.genitalium / *T.vaginalis*-MULTIPRIME-FRT
PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 50409

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Name: **AmpliSens® *N.gonorrhoeae / C.trachomatis /
M.genitalium*-MULTIPRIME-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 50409

Name: **AmpliSens® Genoscreen HLA B*5701-FRT
PCR kit**

Trade name(s): -

Model(s): variant FRT

Classification: List B

GMDN: 56403

Name: **AmpliSens® *Mycoplasma pneumoniae /
Chlamydomphila pneumoniae*-FEP PCR kit**

Trade name(s): -

Model(s): variant FEP (0.2 ml tubes)

Classification: List B

GMDN: 58957

Date: 2022-04-28
Revision: k

A handwritten signature in blue ink, appearing to read 'Mgr. Jiří Heš'.

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Name: **AmpliSens® *Mycoplasma pneumoniae* /
Chlamydomydia pneumoniae-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 58957

Name: **AmpliSens® *T.vaginalis* / *N.gonorrhoeae* /
C.trachomatis-MULTIPRIME-FRT PCR kit**

Trade name(s): -

Model(s): variant FRT-100 F

Classification: List B

GMDN: 61144

Facility(ies):

Federal Budget Institute of Science “Central Research Institute for Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”
3a Novogireevskaya Street, Moscow 111123, Russia**

Certificate History:

Revision	Date	Reference Number	Action
	2011-01-24	813600111	Certification
a	2011-07-21	813600161	Change of manufacturer name
b	2012-02-13	343601304	Product scope extension
c	2014-05-13	343602568	Product scope extension
d	2016-01-15	813600504a	Prolongation of certificate validity
e	2016-06-17	813600504	Re-certification process
f	2016-08-29	343603690	Change of manufacturer facility address
g	2017-11-30	343603888	Changes of product compositions, packaging and quality system documentation
h	2018-10-31	813600754	Change of product labelling, shelf life extension and quality system documentation
i	2019-05-09	813600859	Product shelf life extension
j	2021-04-27	813601045	Re-certification process

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

Annex to EC Certificate No. 11 0040 QS/NB

issued for manufacturer:

**Federal Budget Institute of Science “Central Research Institute for
Epidemiology”**
3a Novogireevskaya Street, Moscow 111123, Russia

Revision	Date	Reference Number	Action
k	2022-04-28	813601141	Extension of the certificate validity regarding to REGULATION (EU) 2022/112 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL, dated 25 th January 2022

Date: 2022-04-28
Revision: k

Mgr. Jiří Heš
Representative of the Notified Body No. 1023

FEDERAL SERVICE FOR SUPERVISION OF CONSUMER RIGHTS PROTECTION AND HUMAN WELFARE

FEDERAL BUDGET INSTITUTE OF SCIENCE
«CENTRAL RESEARCH INSTITUTE FOR EPIDEMIOLOGY»

111123, Moscow, 3A Novogireevskaya street, Tel.: +7 495 974 96 42, Fax: +7 495 305 54 23,
e-mail: obtk@pcr.ru

EC DECLARATION OF CONFORMITY

Directive 98/79/EC of the European Parliament and of the Council of 27th of October 1998 on
In Vitro Diagnostic Medical Devices

Federal Budget Institute of Science "Central Research Institute for Epidemiology" hereby under own responsibility declares that the products covered by the declaration conform with Essential Requirements listed in Annex I of EC Directive 98/79/EC (IVD Directive). Supporting documentation is retained under the premises of the manufacturer.

The quality management system meets the requirements of the standard EN ISO 13485 "Medical devices – Quality management systems – Requirements for regulatory purposes" and is certified by Institute for testing and certification, Inc. (certificate No. 21 0023 SJ, valid until 26.04.2024).

Manufacturer:	Federal Budget Institute of Science "Central Research Institute for Epidemiology"
Authorized Representative:	Ecoli Dx, s.r.o. Purkyňova 74/2 Praha 1, 110 00 Czech Republic Tel: +420 325 209 912 Cell: +420 739 802 523 E-mail: ecolix@ecoli.sk
Product Name:	Annex for this Declaration
Description:	Reagent kits for qualitative detection and quantification of DNA (RNA) of different infectious agents or HLA B*5701 DNA in human specimens
Classification:	Article 9, paragraph 3 of EC Council Directive 98/79/EC on <i>in Vitro</i> Diagnostic Devices Annex II List B IVDs (According to EC Declaration of Conformity List)
Conformity Assessment Route:	Annex IV (IVDD) Full QA System
Notified Body:	Institute for testing and certification, Inc. třída Tomáše Bati 299 Louky, 763 02 Zlín, Czech Republic E-mail: itc@itczlin.cz Notified Body No. 1023
EC Certificate:	No. 11 0040 QS/NB revision k, valid until 2025-05-26
Place, Date of Issue:	Zlín, Czech Republic, 2022-04-28

Signed

Full name Vasily G. Akimkin
Title Director

Valid from 2022-04-28

Valid until 2025-05-26

№№	Description	Model(s)
1.	AmpliSens® <i>Rubella virus</i> -FRT PCR kit	variant FRT-50 F
2.	AmpliSens® <i>Toxoplasma gondii</i> -FRT PCR kit	variant FRT-50 F
3.	AmpliSens® <i>CMV</i> -FEP PCR kit	variant FEP (0.2-ml tubes)
4.	AmpliSens® <i>CMV</i> -FRT PCR kit	variant FRT-100 F
5.	AmpliSens® <i>HSV / CMV</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
6.	AmpliSens® <i>CMV</i> -screen/monitor-FRT PCR kit	variant FRT-100 F
7.	AmpliSens® <i>EBV / CMV / HHV6</i> -screen-FRT PCR kit	variant FRT-100 F
8.	AmpliSens® <i>Chlamydia trachomatis</i> -FEP PCR kit	variant FEP (0.2-ml tubes)
9.	AmpliSens® <i>Chlamydia trachomatis</i> -FRT PCR kit	variant FRT variant FRT-100 F
10.	AmpliSens® <i>C.trachomatis / Ureaplasma / M.genitalium</i> -MULTIPRIME-FEP PCR kit	variant FEP (0.2-ml tubes)
11.	AmpliSens® <i>C.trachomatis / Ureaplasma / M.genitalium</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
12.	AmpliSens® <i>C.trachomatis / Ureaplasma / M.hominis</i> -MULTIPRIME-FEP PCR kit	variant FEP (0.2-ml tubes)
13.	AmpliSens® <i>C.trachomatis / Ureaplasma / M.hominis</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
14.	AmpliSens® <i>C.trachomatis / Ureaplasma / M.genitalium / M.hominis</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
15.	AmpliSens® <i>N.gonorrhoeae / C.trachomatis / M.genitalium / T.vaginalis</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
16.	AmpliSens® <i>N.gonorrhoeae / C.trachomatis / M.genitalium</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F
17.	AmpliSens® Genoscreen HLA B*5701-FRT PCR kit	variant FRT
18.	AmpliSens® <i>Mycoplasma pneumoniae / Chlamydophila pneumoniae</i> -FEP PCR kit	variant FEP (0.2-ml tubes)
19.	AmpliSens® <i>Mycoplasma pneumoniae / Chlamydophila pneumoniae</i> -FRT PCR kit	variant FRT-100 F
20.	AmpliSens® <i>T.vaginalis / N.gonorrhoeae / C.trachomatis</i> -MULTIPRIME-FRT PCR kit	variant FRT-100 F

EC DECLARATION OF CONFORMITY
Directive 98/79/EC of the European Parliament and of the Council of 27th of October 1998 on
In Vitro Diagnostic Medical Devices

Federal Budget Institute of Science "Central Research Institute for Epidemiology" hereby under own responsibility declares that the products covered by the declaration conform with Essential Requirements listed in Annex I of EC Directive 98/79/EC (IVD Directive). Supporting documentation is retained under the premises of the manufacturer.

The quality management system meets the requirements of the standard EN ISO 13485 "Medical devices – Quality management systems – Requirements for regulatory purposes" and is certified by Institute for testing and certification, Inc. (certificate No. 21 0023 SJ, valid until 26.04.2024).

Manufacturer:	Federal Budget Institute of Science "Central Research Institute for Epidemiology"
Authorised Representative:	Ecoli Dx, s.r.o Purkyňova 74/2 Praha 1, 110 00 Czech Republic Tel: +420 325 209 912 Cell: +420 739 802 523 Email: ecoli@ecoli.sk
Product Name:	Annex for this Declaration
Description:	Reagent kits for qualitative detection and quantification of DNA (RNA) of different infectious agents
Classification:	Article 9, paragraph 1 of EC Council Directive 98/79/EC on <i>in Vitro</i> Diagnostic Devices
Conformity Assessment Route:	Annex III (IVDD)

Signed _____

Full name: Vasily G. Akimkin
Title: Director

Valid from 25.05.2022

№№	Description	Product Code (for reference only)
1.	AmpliSens® All bacto-screen-FRT PCR kit	H-2631-1-CE H-2632-1-4-CE
2.	AmpliSens® All-screen-FRT PCR kit	R-B45(RG,iQ)-CE
3.	AmpliSens® All viro-screen-FRT PCR kit	H-2761-1-CE
4.	AmpliSens® ARVI-screen-FRT PCR kit	R-V57-100-F(RG,iQ,Dt)-CE
5.	AmpliSens® Ascariidosis-FRT PCR kit	H-1971-1-CE
6.	AmpliSens® <i>Bacillus anthracis</i> -FRT PCR kit	R-B41(RG)-CE
7.	AmpliSens® <i>Bordetella</i> multi-FRT PCR kit	R-B84-100-F(RG,iQ,Dt)-CE
8.	AmpliSens® <i>Borrelia burgdorferi sensu lato</i> -FRT PCR kit	R-B37(RG)-CE
9.	AmpliSens® <i>Borrelia miyamotoi</i> -FRT PCR kit	H-2791-1-CE H-2792-1-4-CE
10.	AmpliSens® BRCA1-FRT PCR kit	S-3901-1-CE
11.	AmpliSens® <i>Brucella</i> spp.-FRT PCR kit	R-B10(RG)-CE
12.	AmpliSens® <i>C.albicans</i> / <i>C.glabrata</i> / <i>C.krusei</i> -MULTIPRIME-FRT PCR kit	R-F3-F(RG,iQ)-CE
13.	AmpliSens® <i>Candida albicans</i> -FEP PCR kit	F1-100-R0,2-FEP-CE
14.	AmpliSens® <i>Candida albicans</i> -FRT PCR kit	R-F1-F(RG,iQ)-CE
15.	AmpliSens® CCHFV-FRT PCR Kit	R-V22-50-F(RG,iQ,Mx,Dt)-CE
16.	AmpliSens® <i>Corynebacterium diphtheriae</i> / tox-genes-FRT PCR kit	H-2842-1-CE H-2843-1-4-CE
17.	AmpliSens® Cov-Bat-FRT PCR kit	H-2242-1-CE
18.	AmpliSens® COVID-19-FL PCR kit	H-4094-1-1-CE
19.	AmpliSens® <i>Coxiella burnetii</i> -FRT PCR kit	R-B85-50-F(RG,iQ,Mx,Dt)-CE
20.	AmpliSens® <i>Cryptococcus neoformans</i> -FRT PCR kit	R-F4-F(RG,iQ)-CE
21.	AmpliSens® Dengue virus type-FRT PCR kit	R-V63(RG,CFX)-CE
22.	AmpliSens® Dengue virus-FRT PCR kit	H-2391-1-CE H-2392-1-4-CE
23.	AmpliSens® EBOV Zaire-FRT PCR kit	R-V69-50-F-CE
24.	AmpliSens® EBV-screen/monitor-FRT PCR kit	R-V9-100-S(RG,iQ,Mx)-CE
25.	AmpliSens® Enterovirus 71-FRT PCR kit	R-V64-F-CE
26.	AmpliSens® Enterovirus-FRT PCR kit	R-V16(RG)-CE
27.	AmpliSens® Enterovirus / Parechovirus-FRT PCR kit	H-3751-1-2-CE
28.	AmpliSens® ESBL CTX-M-FRT PCR kit	HN-3571-1-CE
29.	AmpliSens® Escherichioses-FRT PCR kit	R-B62(RG,iQ)-CE
30.	AmpliSens® F2/F5-SNP-FRT PCR kit	S-3451-1-CE S-3452-1-4-CE
31.	AmpliSens® FiloA-screen-FRT PCR kit	H-2781-1-4-CE
32.	AmpliSens® Florocenosis / Aerobes-FRT PCR kit	R-B88-100-FT-CE
33.	AmpliSens® Florocenosis / Bacterial vaginosis-FRT PCR kit	R-B74-100-FT(RG)-CE
34.	AmpliSens® Florocenosis / <i>Candida</i> -FRT PCR kit	R-F5-100-FT(RG,CFX)-CE
35.	AmpliSens® Florocenosis / <i>Mycoplasma</i> -FRT PCR kit	R-B75-100-FT(RG,iQ,Mx)-CE
36.	AmpliSens® <i>Giardia lamblia</i> -FRT PCR kit	H-2821-1-CE H-2822-1-4-CE

№№	Description	Product Code (for reference only)
37.	AmpliSens® <i>Gardnerella vaginalis</i> -FEP PCR kit	B7-100-R0,2-FEP-CE
38.	AmpliSens® <i>Gardnerella vaginalis</i> -FRT PCR kit	R-B7-F(RG,iQ)-CE
39.	AmpliSens® Genoscreen- <i>IL28B</i> -FRT PCR kit	R-O5-100-F(RG,iQ,Dt,CFX)-CE
40.	AmpliSens® <i>HAV</i> -FRT PCR kit	R-V4(RG,iQ)-CE
41.	AmpliSens® <i>Helicobacter pylori</i> -FRT PCR kit	R-B9(RG,iQ)-CE H-2831-1-CE H-2832-1-4-CE
42.	AmpliSens® Hemochromatosis-FRT PCR kit	S-2451-1-CE S-2452-1-4-CE
43.	AmpliSens® <i>HGV</i> -FRT PCR kit	R-V2-50-F(RG,iQ,Mx,Dt)-CE
44.	AmpliSens® <i>HHV6</i> -screen-titre-FRT PCR kit	R-V10-T(RG,iQ,Mx)-CE
45.	AmpliSens® <i>HHV7</i> -screen/monitor-FRT PCR kit	H-2431-1-1-CE
46.	AmpliSens® <i>HHV8</i> -screen/monitor-FRT PCR kit	H-3581-1-1-CE H-3582-1-14-CE
47.	AmpliSens® <i>HPV 16/18</i> -FRT PCR kit	R-V12-100-CE R-V12-F-CE R-V12(RG,iQ,Mx)-CE
48.	AmpliSens® <i>HPV 6/11</i> -FRT PCR kit	R-V11-100-CE R-V11-Mod(RG,iQ,Mx)-CE R-V11(RG,iQ,Mx)-CE
49.	AmpliSens® <i>HPV</i> HCR genotype-FRT PCR kit	R-V25(RG,iQ,Mx)-CE
50.	AmpliSens® <i>HPV</i> HCR genotype-titre-FRT PCR kit	R-V67-F-CE H-2261-1-13-CE
51.	AmpliSens® <i>HPV</i> HCR screen-EPh PCR kit	V31-100F-CE
52.	AmpliSens® <i>HPV</i> HCR screen-FEP PCR kit	V31-3x-FEP-CE V31-FEP-CE
53.	AmpliSens® <i>HPV</i> HCR screen-titre-14-FRT PCR kit	H-2311-1-13-CE
54.	AmpliSens® <i>HPV</i> HCR screen-titre-FRT PCR kit	R-V31-T-2x(RG,iQ,SC)-CE R-V31-T-4x(RG,iQ,Mx)-CE R-V31-F-CE
55.	AmpliSens® <i>HSV I, II</i> -FRT PCR kit	R-V8-F(RG,iQ)-CE
56.	AmpliSens® <i>HSV</i> -typing-FEP PCR kit	V38-100-R0,2-FEP-CE
57.	AmpliSens® <i>HSV</i> -typing-FRT PCR kit	R-V38-F(RG,iQ)-CE
58.	AmpliSens® <i>Human enterovirus</i> -FEP PCR kit	H-2771-2-2-CE H-2772-2-CE
59.	AmpliSens® <i>Human enterovirus</i> -FRT PCR kit	H-2771-1-2-CE H-2773-1-CE H-2773-1-4-CE
60.	AmpliSens® <i>Influenza virus A/H1-swine</i> -FEP PCR kit	V55-50-R0,2-FEP-CE
61.	AmpliSens® <i>Influenza virus A/H1-swine</i> -FRT PCR kit	R-V55(RG)-CE R-V55-F(SC)-CE
62.	AmpliSens® <i>Influenza virus A-type</i> -FRT PCR kit	R-V54-100-F(RG,iQ,Dt,SC)-CE
63.	AmpliSens® <i>Influenza virus A-type-H5, H7, H9</i> -FRT PCR kit	R-V66-F-CE
64.	AmpliSens® <i>Influenza virus A H5N1</i> -FRT PCR kit	R-V33(SC)-CE
65.	AmpliSens® <i>Influenza virus A/B</i> -FRT PCR kit	R-V36-100-F-Mod(RG,iQ,Dt,CFX,SC)-CE
66.	AmpliSens® <i>Influenza virus B-type</i> -FRT PCR kit	H-3991-1-23-CE H-3992-1-3-CE
67.	AmpliSens® <i>JCV-BKV</i> screen-monitor-FRT PCR kit	H-2441-1-1-CE

№№	Description	Product Code (for reference only)
68.	AmpliSens® <i>Legionella pneumophila</i> -FEP PCR kit	B50-R0,2-FEP-CE
69.	AmpliSens® <i>Legionella pneumophila</i> -FRT PCR kit	R-B50(RG)-CE
70.	AmpliSens® <i>Leptospira</i> -FRT PCR kit	R-B49(RG,iQ)-CE
71.	AmpliSens® Leucosis Quantum <i>M-bcr</i> -FRT PCR kit	TR-O1(RG,iQ,Mx,A)-CE
72.	AmpliSens® <i>Listeria monocytogenes</i> -screen/monitor-FRT PCR kit	H-2161-1-1-CE
73.	AmpliSens® MDR A.b.-OXA-FRT PCR kit	HN-3871-1-CE HN-3872-1-4-CE
74.	AmpliSens® MDR KPC/OXA-48-FRT PCR kit	R-C2(RG,CFX)-CE
75.	AmpliSens® MDR MBL-FRT PCR kit	R-C1(RG,CFX)-CE
76.	AmpliSens® MDR MCR-1-FRT PCR kit	HN-4171-1-CE HN-4172-1-4-CE
77.	AmpliSens® MDR VRE-FRT PCR kit	HN-3891-1-CE HN-3892-1-4-CE
78.	AmpliSens® <i>MRSA</i> -screen-titre-FRT PCR kit	R-B78-100-FT(RG,iQ)-CE
79.	AmpliSens® <i>MTC</i> -diff-FRT PCR kit	R-B80(RG,iQ,Dt,SC)-CE
80.	AmpliSens® <i>MTC</i> -MDR-FRT PCR kit	H-3611-1-CE H-3612-1-4-CE
81.	AmpliSens® <i>MTC</i> -FEP PCR kit	B57-FEP-CE
82.	AmpliSens® <i>MTC</i> -FRT PCR kit	R-B57(RG,iQ,SC,Dt)-CE
83.	AmpliSens® <i>MTHFR</i> -SNP-FRT PCR kit	S-3721-1-CE S-3722-1-4-CE
84.	AmpliSens® <i>Mycoplasma genitalium</i> -FEP PCR kit	B4-100-R0,2-FEP-CE
85.	AmpliSens® <i>Mycoplasma genitalium</i> -FRT PCR kit	R-B4(RG)-CE R-B4-F(RG,iQ)-CE
86.	AmpliSens® <i>Mycoplasma hominis</i> -FEP PCR kit	B3-100-R0,2-FEP-CE
87.	AmpliSens® <i>Mycoplasma hominis</i> -FRT PCR kit	R-B3(RG)-CE R-B3-F(RG,iQ)-CE
88.	AmpliSens® <i>M.genitalium</i> -ML/FQ-Resist-FRT PCR kit	H-3971-1-CE
89.	AmpliSens® <i>N.meningitidis</i> / <i>H.influenzae</i> / <i>S.pneumonia</i> -FRT PCR kit	R-B25(RG,iQ)-CE
90.	AmpliSens® <i>Neisseria gonorrhoeae</i> -screen-FEP PCR kit	B51-100-R0,2-FEP-CE
91.	AmpliSens® <i>Neisseria gonorrhoeae</i> -screen-FRT PCR kit	R-B51(RG)-CE R-B51-F(RG,iQ)-CE
92.	AmpliSens® NmABCW-FRT PCR kit	H-3861-1-3-CE
93.	AmpliSens® <i>Norovirus</i> GI / GII-FRT PCR kit	H-2751-1-3-CE
94.	AmpliSens® <i>Parvovirus</i> B19-FRT PCR kit	R-V49(RG,iQ,Mx)-CE
95.	AmpliSens <i>Plasmodium</i> spp. / <i>P.falciparum</i> / <i>P.vivax</i> -FRT PCR kit	H-3981-1-CE H-3982-1-4-CE
96.	AmpliSens® <i>Pneumocystis jirovecii (carinii)</i> -FRT PCR kit	R-F2-Mod(RG,iQ,Mx)-CE
97.	AmpliSens® Pneumo-quantum-FRT PCR kit	H-2811-1-1-CE H-2812-1-14-CE
98.	AmpliSens® <i>Poliovirus</i> -FRT PCR kit	R-V58(RG,iQ)-CE
99.	BRCA-screen kit	S-1619-6-CE
100.	PEERO-prep reagent kit for sample preparation	K15-1611-40-CE
101.	AmpliSens® Pyroscreen PHARMA-screen-Imatinib kit	S-16121-6-CE

№№	Description	Product Code (for reference only)
102.	AmpliSens® Pyroscreen PHARMA-screen-transport kit	S-16119-6-CE
103.	AmpliSens® Pyroscreen PHARMA-screen-Warfarin kit	S-16120-6-CE
104.	AmpliSens® <i>Rickettsia conorii</i> -FRT PCR kit	H-2741-1-CE H-2742-1-4-CE
105.	AmpliSens® <i>Rickettsia</i> spp. SFG-FRT PCR kit	H-3741-1-CE H-3742-1-4-CE
106.	AmpliSens® Rotavirus / Norovirus / Astrovirus-FRT PCR kit	R-V40(RG,iQ)-CE
107.	AmpliSens® SARS-CoV-2-IT reagent kit	H-4121-10-CE
108.	AmpliSens® SARS-CoV-2-N501Y-IT reagent kit	H-4161-10-CE
109.	AmpliSens® <i>Shigella</i> spp. and <i>EIEC</i> / <i>Salmonella</i> spp. / <i>Campylobacter</i> spp.-FRT PCR kit	R-B44(RG,iQ)-CE
110.	AmpliSens® <i>Streptococcus agalactiae</i> -screen-titre-FRT PCR kit	R-B77-100-FT(RG,iQ)-CE
111.	AmpliSens® <i>Streptococcus pyogenes</i> -screen/monitor-FRT PCR kit	H-2171-1-1-CE H-2172-1-14-CE
112.	AmpliSens® <i>T.vaginalis</i> / <i>N.gonorrhoeae</i> -MULTIPRIME-FRT PCR kit	R-B65-F(RG,iQ)-CE
113.	AmpliSens® TBE-FRT PCR kit	R-V52(RG)-CE
114.	AmpliSens® TBEV, <i>B.burgdorferi</i> sl, <i>A.phagocytophilum</i> , <i>E.chaffeensis</i> / <i>E.muris</i> -FRT PCR kit	R-V59(RG,iQ,Mx,Dt)-CE
115.	AmpliSens® <i>Treponema pallidum</i> -FRT PCR kit	R-B20-F(RG,iQ)-CE
116.	AmpliSens® <i>Trichomonas vaginalis</i> -EPH PCR kit	B6-100-R0,2-CE
117.	AmpliSens® <i>Trichomonas vaginalis</i> -FEP PCR kit	B6-100-R0,2-FEP-CE
118.	AmpliSens® <i>Trichomonas vaginalis</i> -FRT PCR kit	R-B6-F(RG,iQ)-CE
119.	AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i> -FEP PCR kit	B19-100-R0,2-FEP-CE
120.	AmpliSens® <i>U.parvum</i> / <i>U.urealyticum</i> -FRT PCR kit	R-B19(RG)-CE R-B19-F(RG,iQ)-CE
121.	AmpliSens® <i>Ureaplasma</i> spp.-FRT PCR kit	R-B2(RG)-CE R-B2-F(RG,iQ)-CE
122.	AmpliSens® <i>Ureaplasma</i> spp.-screen-titre-FRT PCR kit	R-B2-100-FT(RG,iQ,Mx)-CE
123.	AmpliSens® <i>Vibrio cholerae</i> -FRT PCR kit	R-B53(RG)-CE
124.	AmpliSens® VZV-FRT PCR kit	R-V61-50-F(RG)-CE
125.	AmpliSens® WNV-FRT PCR kit	R-V53(RG,iQ,Mx)-CE
126.	AmpliSens® Yellow fever virus-FRT PCR kit	H-2461-1-CE H-2462-1-4-CE
127.	AmpliSens® <i>Yersinia enterocolitica</i> / <i>Y.pseudotuberculosis</i> -FRT PCR kit	R-B64(RG,iQ)-CE
128.	AmpliSens® <i>Yersinia pestis</i> -FRT PCR kit	R-B79(RG,iQ,Dt)-CE
129.	AmpliSens® Zika virus-FRT PCR kit	H-2411-1-CE

Signed _____

Full name: Vasily Akimkin
Title: Director

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CERTIFICATE

No. 21 0023 SJ

We confirm on the basis of a performed audit that company

Federal Budget Institute of Science “Central Research Institute for Epidemiology“

3a, Novogireevskaya str., 111123 Moscow, Russian Federation
Company VAT No.: 7720024671

has implemented and documented a functional quality management system
in compliance with the requirements of the standard

EN ISO 13485:2016

Covering the following activities:

Design and development, manufacturing and final control of *in vitro* diagnostic
medical devices

The Certificate is issued on the basis of the results mentioned in Audit Report No. 233404948/2021.
The Certificate validity is conditioned by positive results of surveillance audits, which the certified
company committed to undergo.

During use of the Certificate the Certificate Holder undertakes to follow the Rules of Use of the Certificate. This
document is publicly available on www.itczlin.cz

Date of Issue: 27. 04. 2021
Valid until: 26. 04. 2024

Date of the first certification awarding: 20. 05. 2015

Ing. Pavel Vaněk
Head of Certification Body

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Instruction Manual

KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Sufficient for
	<i>In vitro</i> diagnostic medical device		Use-by Date
	Version		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Manufacturer	NCA	Negative control of amplification
	Date of manufacture	C-	Negative control of extraction
	Authorized representative in the European Community	C+	Positive control of amplification
PCE	Positive control of extraction	IC	Internal control

1. INTENDED USE

AmpliSens® *Leptospira*-FRT PCR kit is an *in vitro* nucleic acid amplification test for qualitative detection of 16S RNA of pathogenic *Leptospira* genospecies in the clinical material (blood and cerebrospinal fluid), autopsy material (brain, kidney, liver, lung tissue, and mesenteric lymph nodes) and biological material (material obtained from died animals (lung, brain, and kidney tissue) and animals suffering from acute leptospirosis (blood) or *Leptospira* persisting in kidneys (urine)) using real-time hybridization-fluorescence detection of amplified products.

NOTE: The results of PCR analysis are taken into account in complex diagnostics of disease.

2. PRINCIPLE OF PCR DETECTION

Detection of 16S RNA of pathogenic *Leptospira genospecies* by the polymerase chain reaction (PCR) is based on the amplification of the pathogen genome specific region using specific primers. In the real-time PCR, the amplified product is detected with the use of fluorescent dyes. These dyes are linked to oligonucleotide probes, which bind specifically to the amplified product during thermocycling. The real-time monitoring of the fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

AmpliSens® *Leptospira*-FRT PCR kit is a qualitative test that contains the Internal Control (Internal Control STI-87-rec (IC)). It must be used in the isolation procedure in order to control the extraction process of each individual sample and to identify possible reaction inhibition.

AmpliSens® *Leptospira*-FRT PCR kit uses "hot-start", which greatly reduces the frequency of nonspecifically primed reactions. "Hot-start" is guaranteed by separation of nucleotides and Taq-polymerase using chemically modified polymerase (TaqF). Chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The results of amplification are registered in the following fluorescence channels:

Table 1

Channel for fluorophore	FAM	JOE
DNA-target	IC STI-87-rec cDNA	<i>Leptospira</i> cDNA
Target gene	Artificially synthesized sequence	16S RNA

3. CONTENT

AmpliSens® *Leptospira*-FRT PCR kit is produced in 1 form:

variant FRT R-B49(RG,iQ)-CE.

Variant FRT includes:

Reagent	Description	Volume, ml	Quantity
RT-G-mix-2	colorless clear liquid	0.01	2 tubes
RT-PCR-mix-1-FRT <i>Leptospira</i>	colorless clear liquid	0.6	1 tube
RT-PCR-mix-2-FEP/FRT	colorless clear liquid	0.3	1 tube
Polymerase (TaqF)	colorless clear liquid	0.03	1 tube
TM-Revertase (MMIv)	colorless clear liquid	0.015	1 tube
Positive Control cDNA <i>Leptospira</i> (C+ <i>Leptospira</i>)	colorless clear liquid	0.1	1 tube
RNA-eluent	colorless clear liquid	0.07	2 tubes
Negative Control (C-)*	colorless clear liquid	1.6	2 tubes
Positive Control <i>Leptospira</i> -rec	colorless clear liquid	0.03	5 tubes
Internal Control STI-87-rec (IC)**	colorless clear liquid	0.12	5 tubes

* must be used in the extraction procedure as Negative Control of Extraction

** add 10 µl of Internal Control during the DNA extraction procedure directly to the sample/lysis mixture (see **RIBO-sorb** K2-1-Et-50-CE, **RIBO-zol-C** K2-13-50-CE, **RIBO-prep** K2-9-Et-50-CE protocols).

Variant FRT PCR kit is intended for 60 reactions, including controls.

4. ADDITIONAL REQUIREMENTS

- RNA extraction kit.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile RNase-free pipette tips with aerosol filters (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with a rotor for 2 ml tubes.
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia); iCycler iQ5 (Bio-Rad, USA); or equivalent).
- Disposable polypropylene PCR tubes (0.2-ml)
- Refrigerator for 2 to 8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol barriers and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distinctly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use the PCR kit if the internal packaging was damaged or its appearance was changed.
- Do not use the PCR kit if the transportation and storage conditions according to the Instruction Manual were not observed.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid samples and reagents contact with the skin, eyes, and mucous membranes. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice immediately.
- Safety Data Sheets (SDS) are available on request.
- The PCR kit is intended for single use for PCR analysis of specified number of samples (see the section "Content").
- The PCR kit is ready for use in accordance with the Instruction Manual. Use the PCR kit strictly for intended purpose.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.

Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING

Obtaining samples of biological materials for PCR-analysis, transportation and storage is described in manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

AmpliSens® Leptospira-FRT PCR kit is intended for analysis of RNA extracted with RNA extraction kits from the following material.

Human material

- clinical material (blood, cerebrospinal fluid);
- autopsy material (brain, kidney, liver, and lung tissue and mesenteric lymph nodes).

Animal material (biological material)

- urine;
- blood;
- brain, kidney, and lung tissue.

The material can be stored at 2-8 °C for 1 day. The autopsy material can be stored at the temperature not more than minus 16 °C for 1 week, at the temperature not more than minus 68 °C for a long time.

Sampling and pretreatment

6.1. **Blood and cerebrospinal fluid.** Whole blood is taken in the morning after overnight fasting into the tube with 6 % EDTA solution in the ratio 1 : 20. The closed tube with whole peripheral blood should be rotated several times. The tube with blood should be centrifuged at 1,000 g for 10 min to obtain blood plasma (if the blood was stored at 2-8 °C more than 1 hour, it should be mixed carefully by inverting the tube). Transfer 1 ml of plasma into two tubes. Two tubes with 1.0 ml of plasma should be centrifuged at 13,000 rpm for 10 min to concentrate bacterial cells. Then, 900 µl of the supernatant plasma should be removed with a filter tip into the container with disinfectant. The pellet and 100 µl of the supernatant are tested for the presence of *Leptospira* 16S RNA. The second pellet prepared in the same way should be stored at the temperature not more than minus 16 °C for repeated extraction (if any technological procedure is performed incorrectly). The pellet obtained from blood plasma can be stored at the temperature not more than minus 16 °C for 1 week or at the temperature not more than minus 68 °C for a long time.

When cerebrospinal fluid is analyzed, the pellet is obtained by the same procedure by centrifugation at 13,000 rpm for 10 min. The pellet and 100 µl of the supernatant are analyzed.

6.2. **Urine** for analyses is taken into a sterile container. If there is no chance to test material within 24 h after sampling, urine is transferred to a centrifuge tube or an Eppendorf tube. The contents of the tube is mixed with glycerol (~10 % v/v) and frozen. It can be stored at the temperature not more than minus 16 °C for 1 week or at the temperature not more than minus 68 °C for a long time.

If a cooling centrifuge (4 °C) with a speed of 9,000-10,000 g intended for 30-ml tubes is available, the following sample preparation procedure is used. The sample is centrifuged at 9,000-10,000 g for 10 min, the supernatant is transferred to a container with disinfectant. Leave ~1 ml of the supernatant over the pellet in the tube and resuspend it. Transfer the suspension to a new tube and concentrate it by centrifugation at 13,000 rpm for 10 min. Then, 900 µl of the supernatant is transferred to the container with disinfectant, and the pellet and 100 µl of the supernatant is used for RNA isolation. In case of large quantities of salts and mucus, 100 µl of the supernatant and the upper layer of cells should be carefully taken from the salt pellet and transferred into a new tube for RNA isolation.

If you have no centrifuge for 30-ml tubes and a speed of 9,000-10,000 g, bacteria are concentrated from 1 ml of urine as described above using 1.5-ml tubes and a microcentrifuge for Eppendorf tubes. The remaining urine should be decontaminated in a disinfectant.

6.3. **Animal internal organs and autopsy material** is to be homogenized in sterile porcelain mortars with pestles. Then, 10 % suspension in sterile saline or phosphate buffer is prepared; 30 µl of the suspension is taken for RNA extraction.

7. WORKING CONDITIONS

AmpliSens® Leptospira-FRT PCR kit should be used at 18-25 °C.

8. PROTOCOL

8.1. RNA Extraction

It is recommended to use the following nucleic acid extraction kits:

- RIBO-prep, **REF** K2-9-Et-50-CE,
- RIBO-zol-C, **REF** K2-13-50-CE,
- RIBO-sorb, **REF** K2-1-Et-50-CE.

NOTE: Extract the RNA according to the manufacturer's protocol.

NOTE: If using the **RIBO-prep kit** extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- If RNA is extracted from the tissues homogenates. Add 10 µl of **Internal Control STI-87-rec (IC)** to each tube and then add 300 µl of **Solution for Lysis**. Label test tubes. Then, add 50 µl of test suspensions to the tubes with **Solution for Lysis** and **Internal Control STI-87-rec (IC)**.
- If RNA is extracted from the cell pellets after *Leptospira* concentration, 300 µl of **Solution for Lysis** and 10 µl of **Internal Control STI-87-rec (IC)** should be added directly into the tubes with the pellets.
- Add only 10 µl of **Internal Control STI-87-rec (IC)** and 300 µl of **Solution for Lysis** into the tube for the Negative Control of extraction (C-).
- Add 300 µl of **Solution for Lysis**, and then add 90 µl of **Negative Control (C-)**, 10 µl of **Internal Control STI-87-rec (IC)** and 10 µl of **Positive Control Leptospira-rec** into the tube for the Positive Control of extraction (PCE).
- 40 µl of **RNA-buffer** is to be used.

NOTE: If using the **RIBO-zol-C** and **RIBO-sorb kit** extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- Add 300 µl of **Solution D**, 10 µl of **Internal Control STI-87-rec (IC)** and 100 µl of **Negative control (C-)** into the tube for the Negative Control of extraction (C-).
- Add 300 µl of **Solution D**, and then add 90 µl of **Negative Control (C-)**, 10 µl of **Internal Control STI-87-rec (IC)** and 10 µl of **Positive Control Leptospira-rec** into the tube for the Positive Control of extraction (PCE).
- After the first stage of RNA extraction with **RIBO-zol-C** kit the top phase of obtained samples (400 µl when extracting RNA from cerebrospinal fluid, urine pellets and tissues homogenates, 200 µl - from blood pellets, 450 µl - from C- and PCE samples) should be transferred into new 1.5-ml tubes with 400 µl of **Lysis Solution** (from RIBO-sorb kit).
- 40 µl of **RNA-buffer** is to be used.

8.2. Preparing reverse transcription and PCR

The total reaction volume is 25 µl, the volume of DNA sample is 10 µl.

NOTE: Only RNase-free, DNase-free disposable plastic consumables must be used when working with RNA.

8.2.1 Preparing tubes for reverse transcription and PCR

- 1 Prepare the required number of tubes for amplification of cDNA obtained from clinical and control samples. The type of tubes depends on the PCR instrument used for analysis. For carrying out N reactions with 2 controls, N+2 tubes are required.
- 2 Prepare the reaction mixture, calculating per one reaction:
 - 10 µl of **RT-PCR-mix-1-FRT Leptospira**
 - 5 µl of **RT-PCR-mix-2-FEP/FRT**
 - 0.5 µl of **polymerase (TaqF)**
 - 0.25 µl of **TM-Revertase (MMV)**
 - 0.25 µl of **RT-G-mix-2**
- 3 Transfer 15 µl of the prepared mixture to each tube. Discard the unused mixture.
- 4 Using filter tips add 10 µl of **RNA** samples obtained at the RNA extraction stage into prepared tubes.
- 5 Carry out the control amplification reactions:
 - NCA** - Add 10 µl of **RNA-eluent** to the tube labeled NCA (Negative Control of Amplification).
 - C+** - Add 10 µl of **Positive Control cDNA Leptospira (C+Leptospira)** to the tube labeled C+ (Positive Control of Amplification).
 - C-** - Add 10 µl of the **sample extracted from the Negative Control of Extraction sample** to the tube labeled C- (Negative control of Extraction).
 - PCE** - Add 10 µl of the **sample extracted from the Positive control of Extraction sample** to the tube labeled PCE (Positive control of Extraction).

NOTE: Amplification is to be carried immediately after mixing the reaction mixture, RNA-sample and controls. The time period between addition of RNA-samples into the reaction mixture and amplification starting is to be not more than 10-15 min.

8.2.2. Amplification

- 1 Create a temperature profile on your instrument as follows:

Table 2

Step	Rotor-type instruments ¹			Plate-type instruments ²		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	50	30 min	1	50	30 min	1
2	95	15 min	1	95	15 min	1
3	95	20 s	10	95	20 s	10
	65	50 s		65	50 s	
	72	20 s				
4	95	20 s	38	95	20 s	40
	61	Fluorescence acquiring		61	Fluorescence acquiring	
	72	20 s		65	20 s	

Fluorescent signal is detected in the channels for the FAM and JOE fluorophores.

- 2 Adjust the fluorescence channel sensitivity according to the Guidelines [2].
- 3 Insert tubes into the reaction module of the device.
- 4 Run the amplification program with fluorescence detection.
- 5 Analyze results after the amplification program is completed

9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in two channels:

- The signal of the IC cDNA amplification product is detected in the channel for the FAM fluorophore.
- The signal of the *Leptospira* cDNA amplification product is detected in the channel for the JOE fluorophore.

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a Ct value of the cDNA sample in the corresponding column of the results grid.

Principle of interpretation is the following:

- *Leptospira* cDNA is **detected** if the Ct value determined in the results grid in the channel for the JOE fluorophore is less than the boundary Ct value specified in the Guidelines.
- *Leptospira* cDNA is **not detected** in a sample if the Ct value is not determined (absent) in the channel for the JOE fluorophore, whereas the Ct value determined in the channel for the FAM fluorophore is less than the boundary Ct value specified in the Guidelines.
- The result is **invalid** if the Ct value is not determined (absent) in the channel for the JOE fluorophores whereas the Ct value in the channel for the FAM fluorophore is greater than the specified boundary Ct value. In such cases, the PCR analysis of this sample should be repeated starting from the RNA extraction stage.
- The result is **equivocal** if the Ct value determined in the channel for the JOE fluorophore is greater than the boundary Ct value specified in the Guidelines, whereas the Ct value determined in the channel for the FAM fluorophore is less than the boundary Ct value specified in the Guidelines. In such cases, the PCR analysis of this sample should be repeated two times starting from the RNA extraction stage.

NOTE: Boundary Ct values are specified in the Guidelines [2]

The result of the analysis is considered reliable only if the results obtained for the Positive and Negative Controls of amplification and extraction are correct (see Table 3).

Table 3

Control	Stage for control	Results for controls	
		Ct value in the channel for the fluorophore	
		FAM	JOE
PCE	DNA extraction	Present	< boundary value
C-	DNA extraction	Present	Absent
NCA	PCR	Absent	Absent
C+	PCR	Absent	< boundary value

¹ For example, Rotor-Gene 3000/Rotor-Gene 6000 (Corbett Research, Australia).

² For example, iCycler iQ5 (Bio-Rad, USA).

10. TROUBLESHOOTING

Results of analysis are not being registered in the following cases:

1. If the Ct value is determined for the Negative Control of extraction (C-) in the channel for the JOE fluorophore and/or for the Negative Control of amplification (NCA) in the channels for the FAM and JOE fluorophores in the results grid, it indicates contamination of reagents or samples. In such cases, the results of analysis are considered to be irrelevant. Analysis should be repeated and measures to detect and eliminate the source of contamination should be taken.
2. If no signal is detected for the Negative Control of extraction (C-) in the channel for the FAM fluorophore and/or for the Positive Control of extraction (PCE) in the channels for the FAM and JOE fluorophores, the results of analysis are considered invalid. Analysis of all samples should be repeated starting from the extraction stage.
3. If no signal is detected for Positive Control of amplification (C+) in the channel for the JOE fluorophore, the results of analysis are considered invalid. Analysis of all samples should be repeated starting from the RT-PCR stage.

If you have any further questions or if encounter problems, please contact our Authorized representative in the European Community.

11. TRANSPORTATION

AmpliSens® Leptospira-FRT PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the **AmpliSens® Leptospira-FRT** PCR kit are to be stored at 2–8 °C when not in use (except for RT-G-mix-2, RT-PCR-mix-1-FRT *Leptospira*, RT-PCR-mix-2-FEP/FRT, Polymerase (TaqF), and TM-Revertase (MMIv)). All components of the **AmpliSens® Leptospira-FRT** PCR kit are stable until the expiry date stated on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

NOTE: RT-G-mix-2, RT-PCR-mix-1-FRT *Leptospira*, RT-PCR-mix-2-FEP/FRT, Polymerase (TaqF), and TM-Revertase (MMIv) are to be stored at the temperature from minus 24 to minus 16 °C.

NOTE: RT-PCR-mix-1-FRT *Leptospira* is to be stored away from light

13. SPECIFICATIONS

13.1. Sensitivity

Analytical Sensitivity of **AmpliSens® Leptospira-FRT** PCR kit is not less than 5×10^3 copies per 1 ml of sample (copies/ml).

The claimed analytical features of **AmpliSens® Leptospira-FRT** PCR kit are guaranteed only when additional reagents kits **RIBO-sorb** and **RIBO-zol-C** or **RIBO-prep** (manufactured by Federal Budget Institute of Science "Central Research Institute for Epidemiology") are used.

13.2. Specificity

The analytical specificity of **AmpliSens® Leptospira-FRT** PCR kit is ensured by the selection of specific primers and probes as well as stringent reaction conditions. The primers and probes have been checked for possible homologies to all sequences published in gene banks by sequence comparison analysis.

The clinical specificity of **AmpliSens® Leptospira-FRT** PCR kit was confirmed in laboratory clinical trials.

14. REFERENCES

1. Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR diagnostics", developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being.
2. Guidelines to the **AmpliSens® Leptospira-FRT** PCR kit for qualitative detection of 16S RNA of pathogenic *Leptospira* genospecies in the clinical material, autopsy material and biological material by the polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology".

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of the **AmpliSens® Leptospira-FRT** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
13.12.10	Cover page	The phrase "For Professional Use Only" was added
	Intended use	The phrase "The results of PCR analysis are taken into account in complex diagnostics of disease" was added
	Content	New sections "Working Conditions" and "Transportation" were added The "Explanation of Symbols" section was renamed to "Key to Symbols Used"
	Stability and Storage	The information about the shelf life of reagents before and after the first use was added Information that RT-PCR-mix-1-FRT <i>Leptospira</i> is kept away from light was added
	Key to Symbols Used	The explanation of symbols was corrected
	Footer	Reference number was changed from R-B49(RG)-CE to R-B49(RG,iQ)-CE
	Sampling and Handling	Duration of the urine sample storage was changed Information about storage and disposal of urine samples was added
01.07.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"
16.10.15 ChA	3. Content	Quantity of Negative Control (C-) tubes was changed from 1 to 2
	8.1. RNA Isolation	For RIBO-prep the procedure of positive control of extraction preparation was changed: 90 µl of Negative Control (C-) is to be added additionally
17.01.17 ME	Text	Corrections according to the template
	1. Intended use	Types of biological material was specified
	8.1. RNA extraction	The sections were rewritten
	8.2.2. Amplification	
	9. Data analysis	
14. References	The reference to Guidelines was added	
23.04.20 KK	Through the text	The text formatting was changed
	2. Principle of PCR detection	The table with targets was added.
	Footer	The phrase "Not for use in the Russian Federation" was added

AmpliSens®

Ecolis s.r.o., Studenohorska 12
841 03 Bratislava 47
Slovak Republic
Tel.: +421 2 6478 9336
Fax: +421 2 6478 9040

Federal Budget Institute of Science "Central Research Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia

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Este certificado ampara todas las marcas de estos productos incluidas por el fabricante en su declaración de conformidad. / This certificate covers all trademarks of these products included by the manufacturer in his declaration of conformity.