

Application: Viruton® Pulver is a modern generation preparation for cleaning and high-level disinfection of invasive and non-invasive medical devices. It is recommended for surgical and dental instruments and all kinds of endoscopes, surgical probes, ventilator tubes and other anesthetic equipment. It can also be used for the disinfection of surfaces of medical equipment, such as incubators. It is recommended for manual disinfection, in ultrasonic cleaners, automatic washers and endoscope cleaning baths. The product does not require an activator. Uncontaminated solution is stable up to 30 hours. The innovative formula contains four enzymes (protease, lipase, amylase and mannanase) responsible for removing organic contaminants, such as blood, protein, fats and secretions. The composition of surfactants prevents the powder from dusting and guarantees excellent washing properties.

Effectiveness	EN 14885	10 min.	15 min.	30 min.	60 min.
Bacteria (S. aureus, P. aeruginosa, E. hirae, MRSA)	EN 13727		0,5%		
Bacteria (S. aureus, P. aeruginosa, E. hirae, MRSA)	EN 14561	1%		0,5%	
Bacteria (S. aureus, P. aeruginosa, E. hirae, MRSA)	EN 13697			0,5%	
Fungi (C. albicans)	EN 13624			0,5%	
Fungi (C. albicans)	EN 14562	1%		0,5%	
Fungi (A. Brasiliensis)	EN 13624	1%			0,5%
Fungi (A. Brasiliensis)	EN 14562	2%	1,5%	1%	
Mycobacteria (M. terrae, M. avium)	EN 14348		0,5%		
Mycobacteria (M. terrae)	EN 14563	1%			0,5%
Mycobacteria (M.avium)	EN 14563	1%			
Viruses (Polio, Adeno, Noro, Rota, HIV, HBV, HCV, wirus grypy A, B, C, Vaccina, BVDV, SARS-CoV-2, Herpes simplex, Ebola)	EN 14476	2%		1%	
Spores (Cl. difficile)	EN 13704	1%*,2%		0,5%	
Spores (B. subtilis)	EN 13704	1%*		0,5%	
Spores (Cl. sporogenes)	EN 13704	2%			
Spores (Cl. difficile r027)	EN 17126	1%*, 2%			