Anti-HBc

**For the Qualitative Detection
of Total Antibodies to Hepatitis B Core Antigen
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Anti-HBc

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the qualitative detection of total antibodies against hepatitis B core antigen (HBcAg) in human serum or plasma (heparinized, sodium citrate or EDTA), as an aid in the determination of a prior immune status to the hepatitis B virus.

Catalog Number: **L2KHC2** (200 tests)

Test Code: **aBC** Color: **Light Gray**

Summary and Explanation

Hepatitis B virus (HBV) is the sole human pathogen in the family of hepatitis-associated DNA viruses, and is found world-wide. Distribution of HBV infection will vary among geographical areas and population groups. Transmission of the virus is due to parenteral contact, through the exchange of blood or blood products, sexual contact, and perinatal spread from mother to newborn.^{1,2} Clinical manifestations range from mild asymptomatic infections to severe fulminant hepatitis.^{1,2,3} Over 90% of infected adults will have an acute self-limiting infection,¹ with jaundice and abnormal liver function. Recovery occurs without any chronic sequelae.^{1,2}

Chronic liver disease, a condition in which infection persists for more than six months, a known sequela of a hepatitis B infection, is usually progressive.^{1,2} The risk of developing the chronic carrier state is more likely to follow infection acquired in childhood than as an adult.^{4,5} In chronic HBV carriers, there is no evidence of continued hepatic damage,^{1,2} however, the infection persists and the carrier maintains the ability to transmit the virus.²

Availability of recombinant HBV vaccines, and the recommendation of universal immunization for infants and other high-risk persons has aided in the prevention of HBV infections. In addition, treatment with alpha-interferon to relieve symptoms is available. Results have shown positive response to treatment in 40–50% of

selected individuals with chronic active hepatitis B.^{4,5}

Classification of a hepatitis B infection requires the identification of several serological markers expressed during three phases (incubation, acute and convalescent) of the infection. The first marker to appear during the incubation phase is HBsAg, and indicates an ongoing infection with HBV.^{1,2,4} Anti-HBc appears shortly after the appearance of HBsAg and peaks during the acute phase prior to the appearance of anti-HBs. IgM antibody to the core antigen will decline in uncomplicated acute infection, whereas IgG antibody will persist for years.^{4,5} Anti-HBc is also elevated in chronic HBV infections.⁴

Presence of anti-HBc indicates an ongoing or previous HBV infection. When used in conjunction with tests for other HBV serological markers, diagnosis of viral hepatitis can be achieved.

Principle of the Procedure

IMMULITE 2000 Anti-HBc is a solid-phase, chemiluminescent two-step immunoassay.

Incubation Cycles: 2 × 30 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot

activators and/or anticoagulants. IMMULITE 2000 Anti-HBc has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 µL serum or plasma (heparinized, sodium citrate or EDTA).

Storage: 3 days at 2–8°C.⁶
For longer storage: at –20°C.⁷

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

The Anti-HBc Adjustor, Anti-HBc Low Positive Control and Anti-HBc Positive Control contain HBcAg which has been inactivated by proven, documented methods. However, always handle all controls as if capable of transmitting infectious agents.

The concentration of anti-HBc in a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Anti-HBc EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Anti-HBc Bead Pack (L2HC12)

With barcode. 200 beads coated with purified recombinant HBcAg. Stable at 2–8°C until expiration date.

L2KHC2: 1 pack.

Anti-HBc Reagent Wedge (L2HCA2)

With barcode. 11.5 mL of a protein-based buffer, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to purified monoclonal murine anti-HBc in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KHC2: 1 wedge.

Before use, tear off the top of the label at the perforations without damaging the barcode on the main label. Remove the foil seal from the top of the Reagent Wedge, and snap the sliding cover down into the ramps on the reagent lid.

Anti-HBc Adjustor (LHCR)

2 mL human serum reactive to HBcAg in a buffer, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHC2: 1 vial.

Anti-HBc Controls (LHCC1, LHCC2, LHCC3)

Three vials containing 4 mL each. **LHCC1 (Negative Control)**: human serum nonreactive to HBcAg, with preservative.

LHCC2, LHCC3 (Low Positive Control, Positive Control): human serum reactive to HBcAg, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KHC2: 1 set.

For the current control ratio ranges, please refer to the Control insert.

Aliquot Labels with barcodes are supplied with the kit, for use with the Adjustors and Controls. Before use, place the appropriate Aliquot Labels on test tubes, so the barcodes can be read by the barcode reader on the IMMULITE 2000.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, adjustment, assay and quality control procedures.

Adjustment Interval: 4 weeks.

Quality Control Samples: The controls supplied with the kit should be used as quality control material to monitor assay performance.

For the current control ratio ranges, please refer to the Control insert.

Calculation of Cutoff and CO/S Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry : Kit Entry.)

Calculation of a cutoff/signal (co/s) ratio is done by using the following formula:

$$\text{CO/S Ratio} = \frac{\text{Mean Adjustor cps} \times \text{P1}}{\text{Sample or Control cps}}$$

Calculation and reporting of qualitative (reactive/nonreactive/indeterminate) and co/s ratio results are handled automatically by the IMMULITE 2000.

The result for a sample is "indeterminate" if the counts per second for that sample fall within $\pm 15\%$ of the cutoff. The result is "reactive" if the sample's counts are *below* the indeterminate range, and "nonreactive" if *above* this range.

Interpretation of Results

- A result of "**Reactive**" (co/s ratio of ≥ 1.15) indicates that the patient sample is reactive and that anti-HBc antibodies were detected in the sample, which is indicative of either ongoing or previous infection.
- A result of "**Nonreactive**" (co/s ratio of < 0.85) indicates that the patient sample is nonreactive and that anti-HBc antibodies were not detected in the sample.
- Any result of "**Indeterminate**" (co/s ratio between 0.85 and < 1.15) should be retested. Samples which still test as "Indeterminate" should be tested by an alternate method, or a second sample should be taken — if possible — within a reasonable period of time (e.g., one week).

The results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Anti-HBc EIA. Values obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

Individuals acutely infected with the hepatitis B virus will exhibit anti-HBc between two weeks and four months after exposure, usually through the course of clinical illness.⁴ Antibody levels persist for life.

Limitations

The measurement of Anti-HBc may be affected by the presence of elevated levels of human albumin.

The results of the test must be taken within the context of the patient's clinical and vaccination history, symptomology and other laboratory findings.

A nonreactive result does not indicate that the patient was not infected with HBV. The patient sample should be tested for the presence of other serological markers, such as HBsAg or anti-HBc IgM.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed as a cutoff-to-signal ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Bilirubin: Presence of unconjugated bilirubin in concentrations up to 400 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 5 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.

Results from patients taking biotin supplements or receiving high-dose biotin therapy should be interpreted with caution due to possible interference with this test.

Hemolysis: Presence of hemoglobin in concentrations up to 504 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Blood was collected from 57 laboratory volunteers into plain, heparinized, sodium citrate and EDTA vacutainer tubes, spiked with anti-HBc and then assayed by the IMMULITE 2000 Anti-HBc procedure. Results are expressed as a cutoff-to-signal ratio.

(Heparin) = 0.98 (Serum) + 0.018
 $r = 0.99$

(NaCitrate) = 1.07 (Serum) + 0.027
 $r = 0.94$

(EDTA) = 0.97 (Serum) + 0.031
 $r = 0.99$

Means:
0.255 (Serum)
0.269 (Heparin)
0.301 (NaCitrate)
0.280 (EDTA)

Method Comparison 1: The assay was compared to IMMULITE Anti-HBc on 139 patient samples.

IMMULITE 2000				
IML	Reac	Non-Reac	Relative Sensitivity	Relative Specificity
Reactive	60	0	100%	98.7%
Non-reactive	1	78		

Agreement: 99.3%

Method Comparison 2: The IMMULITE 2000 Anti-HBc procedure was compared to the Elecsys Anti-HBc immunoassay on 2,219 samples, consisting of seroconversion panels, potentially cross-reactive samples, samples from patients with acute, chronic and resolved HBV infection, from blood donors, hospitalized patients, and routine samples. The table below presents the results of this study.

		<u>No. of samples</u>		Sensitivity	Specificity	PPV	NPV
		anti-HBc reactive	anti-HBc non-reactive				
IMMULITE 2000 Anti-HBc							
Reactive	263	1		98.9%	99.9%	99.6%	99.8%
Non-reactive	3	1,953					
Elecsys Anti-HBc							
Reactive	266	9		100%	99.5%	96.7%	100%
Non-reactive	0	1,944					

References

1) Locarnini SA, Gust ID. Hepadnaviridae: hepatitis B virus and the delta virus. In: Balows A, et al, editors. Laboratory diagnosis of infectious diseases: principles and practices. New York: Springer-Verlag, 1988: 750-96. 2) Follett EAC. Diagnosis of hepatitis B infection. In: Young H, McMillan A, editors. Immunological diagnosis of sexually transmitted diseases. New York: Marcel Dekker, 1988: 433-49. 3) Hollinger FB, Dienstag JL. Hepatitis B and D viruses. In: Lennette EH, et al, editors. Manual of clinical microbiology. 6th ed. Washington, D.C.: American Society for Microbiology, 1995:1033-49. 4) Nowicki MJ, Balistreri WF. Hepatitis A to E: building up the alphabet. Contemporary Peds 1992: 118-28. 5) Zuckerman AJ, et al. Hepatitis B virus and hepatitis D virus. In: Principles and practice of clinical virology. 2nd ed. New York: John Wiley & Sons, 1992: 153-72. 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:354-60. 7) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995:322-4. 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.
©2008-2020 Siemens Healthcare Diagnostics Inc. All rights reserved.