



en

DIRECT BILIRUBIN

REF 8G63-22

G95979R03

B8G6D0

ARCHITECT

DIRECT BILIRUBIN

This package insert contains information to run the Direct Bilirubin assay on the ARCHITECT cSystems.

Revised March 2022.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Customer Service: Contact your local representative or find country-specific contact information on www.corelaboratory.abbott

Key to Symbols

ISO 15223 Symbols	Other Symbols
Consult instructions for use	DISTRIBUTED IN THE USA BY Distributed in the USA by
Manufacturer	FOR USE WITH Identifies products to be used together
Sufficient for	INFORMATION FOR USA ONLY Information needed for United States of America only
Temperature limitation	MANUFACTURED FOR Manufactured for
Use by/Expiration date	PRODUCT OF CANADA Product of Canada
IVD <i>In Vitro</i> Diagnostic Medical Device	R1 Reagent 1
LOT Lot Number	R2 Reagent 2
REF List Number	Rx ONLY For use by or on the order of a physician only (applicable to USA classification only).
SN Serial number	

NAME

DIRECT BILIRUBIN

INTENDED USE

The Direct Bilirubin assay is used for the quantitative analysis of direct bilirubin in human serum or plasma.

SUMMARY AND EXPLANATION OF TEST

Red blood cells at the end of their circulating life are broken down in the reticuloendothelial system, mainly the spleen. The resulting heme, once the iron is removed, is then converted to bilirubin. This process accounts for about 80% of the 500 µmol (300 mg) of bilirubin formed daily. Other sources of bilirubin include the breakdown of myoglobin and cytochromes and the catabolism of immature red blood cells in the bone marrow.

Once formed, bilirubin is transported to the liver bound to albumin. This fraction of bilirubin is referred to as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid (mono- and diglucuronides) to form conjugated bilirubin by the enzyme uridyl diphosphate glucuronyl transferase. Conjugated bilirubin or direct bilirubin is excreted via the biliary system into the intestine, where it is metabolized by bacteria to a group of products known collectively as stercobilinogen. Elimination is almost complete and serum levels are normally negligible.

Direct bilirubin is the sum of the conjugated fractions. Direct bilirubin is elevated in conditions causing hepatic obstruction, hepatitis, cirrhosis, several inherited enzyme deficiencies, and inherited defects in canalicular excretion.

PRINCIPLES OF PROCEDURE

Bilirubin determination is generally based on the reaction of bilirubin with a diazotized sulfanilic acid, described by Ehrlich.¹ In this method, direct (conjugated fractions) bilirubin couples with a diazonium salt in the presence of sulfamic acid to form the colored compound azobilirubin. The increase in absorbance at 548 nm due to azobilirubin is proportional to the direct bilirubin concentration.

Methodology: Diazo Reaction

REAGENTS

Reagent Kit

REF 8G63-22 Direct Bilirubin is supplied as a liquid, ready-to-use, two-reagent kit which contains:

R1 10 x 39 mL

R2 10 x 13 mL

Estimated tests per kit: 2000

Calculation is based on the minimum reagent fill volume per kit.

Reactive Ingredients	Concentration
R1 Sulfamic acid	9.7 g/L
R2 2, 4-dichloroaniline	< 0.1 g/L
Sodium nitrite	< 0.1 g/L
HCl	33.54 g/L

REAGENT HANDLING AND STORAGE

Reagent Handling

Remove air bubbles, if present in the reagent cartridge, with a new applicator stick. Alternatively, allow the reagent to sit at the appropriate storage temperature to allow the bubbles to dissipate. To minimize volume depletion, do not use a transfer pipette to remove the bubbles.

CAUTION: Reagent bubbles may interfere with proper detection of reagent level in the cartridge, causing insufficient reagent aspiration which could impact results.

Reagent Storage

Unopened reagents are stable until the expiration date when stored at 2 to 8°C. The reagents should be clear.

Reagent stability is 28 days if the reagent is uncapped and onboard.

Indications of Deterioration

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, extreme turbidity, microbial growth, if calibration does not meet the appropriate package insert and/or **ARCHITECT System Operations Manual** criteria, or if controls do not meet the appropriate criteria.

WARNINGS AND PRECAUTIONS

Precautions for Users

- **IVD**
- For *In Vitro* Diagnostic Use.
- **Rx ONLY**
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens,² Biosafety Level 2³ or other appropriate biosafety practices^{4,5} should be used for materials that contain or are suspected of containing infectious agents.
- The following warnings and precautions apply to: **R1**



DANGER: Contains sulfamic acid.
H314 Causes severe skin burns and eye damage.

Prevention

P260 Do not breathe mist/vapors/spray.
P264 Wash hands thoroughly after handling.
P280 Wear protective gloves/protective clothing/eye protection.

Response

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P310 Immediately call a POISON CENTER or doctor/physician.

Disposal

P501 Dispose of contents/container in accordance with local regulations.

- The following warnings and precautions apply to: **R2**



DANGER: Contains hydrochloric acid.
H314 Causes severe skin burns and eye damage.
H332 Harmful if inhaled.
H290 May be corrosive to metals.



Prevention

P260 Do not breathe mist/vapors/spray.
P264 Wash hands thoroughly after handling.
P280 Wear protective gloves/protective clothing/eye protection.
P234 Keep only in original container.

Response

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P310 Immediately call a POISON CENTER or doctor/physician.
P390 Absorb spillage to prevent material damage.

Disposal

P501 Dispose of contents/container in accordance with local regulations.

- Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.
- For the most current hazard information, see the product Safety Data Sheet.
- Safety Data Sheets are available at www.corelaboratory.abbott.com or contact your local representative.
- For a detailed discussion of safety precautions during system operation, refer to the **ARCHITECT System Operations Manual**, Section 8.

SPECIMEN COLLECTION AND HANDLING

Suitable Specimens

Use serum or plasma specimens without visible hemolysis or lipemia. Refer to the SPECIFIC PERFORMANCE CHARACTERISTICS section of this package insert.

Abbott Laboratories has not verified the assay performance characteristics with neonatal specimens.

NOTE: Abbott Laboratories recommends the use of sample interference indices in the semi-quantitative mode to assist in the determination of sample integrity for all specimens. Refer to the Sample Interference Indices (HIL) application sheets.

- **Serum:** Use serum collected by standard venipuncture techniques into glass or plastic tubes with or without gel barriers. Ensure complete clot formation has taken place prior to centrifugation. Centrifuge according to tube manufacturer's instructions to ensure proper separation of serum from blood cells.
Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may take longer to complete their clotting processes. Fibrin clots may subsequently form in these sera and the clots could cause erroneous test results.
- **Plasma:** Use plasma collected by standard venipuncture techniques into glass or plastic tubes. Acceptable anticoagulants are lithium heparin (with or without gel barrier), sodium heparin, and EDTA. The use of tubes containing sodium fluoride/potassium oxalate is not recommended due to the potential of hemolysis formation with this anticoagulant. Ensure centrifugation is adequate to remove platelets. Centrifuge according to tube manufacturer's instructions to ensure proper separation of plasma from blood cells.

For total sample volume requirements, refer to the ASSAY PARAMETERS section of this package insert and *Section 5* of the **ARCHITECT System Operations Manual**.

Specimen Storage

Serum and Plasma: Specimens should be protected from bright light as bilirubin is photolabile.⁶ Direct bilirubin is stable in serum and plasma as follows:

Temperature	Maximum Storage	Bibliographic Reference
20 to 25°C	2 days	7
2 to 8°C	7 days	7, 8
-20°C	3 months	9
-80°C	3 months	9

Limitations of laboratory equipment make it necessary in practice for clinical laboratories to establish a range around -20°C and/or -80°C for specimen storage. These temperature ranges may be established from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

NOTE: Stored specimens must be inspected for particulates. If present, mix and centrifuge the specimen to remove particulates prior to testing.

PROCEDURE

Materials Provided

[REF] 8G63-22 Direct Bilirubin Reagent Kit

Materials Required but not Provided

- **[REF]** 1E66-05 Bilirubin Calibrator
- Control Material
- Saline (0.85% to 0.90% NaCl) for specimens that require dilution

Assay Procedure

For a detailed description of how to run an assay, refer to *Section 5* of the **ARCHITECT System Operations Manual**.

Specimen Dilution Procedures

The ARCHITECT cSystems have an automatic dilution feature; refer to *Section 2* of the **ARCHITECT System Operations Manual** for additional information.

Serum and Plasma: Specimens with direct bilirubin value exceeding 15.0 mg/dL (256.5 µmol/L) are flagged and may be diluted by following the Manual Dilution Procedure.

PROCEDURE (Continued)

Automated Dilution Protocol

If using the Automated Dilution Protocol, the system performs a dilution of the specimen and automatically corrects the concentration by multiplying the result by the appropriate dilution factor. To set up the automatic dilution feature, refer to *Section 2* of the **ARCHITECT System Operations Manual** for additional information.

Manual Dilution Procedure

Manual dilutions should be performed as follows:

- Use saline (0.85% to 0.90% NaCl) to dilute the sample.
- The operator must enter the dilution factor in the patient or control order screen. The system uses this dilution factor to automatically correct the concentration by multiplying the result by the entered factor.
- If the operator does not enter the dilution factor, the result must be multiplied by the appropriate dilution factor before reporting the result.

NOTE: If a diluted sample result is flagged indicating it is less than the linear low limit, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to *Section 5* of the **ARCHITECT System Operations Manual**.

CALIBRATION

Calibration is stable for approximately 14 days (336 hours) and is required with each change in reagent lot number. Verify calibration with at least two levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

For a detailed description of how to calibrate an assay, refer to *Section 6* of the **ARCHITECT System Operations Manual**.

For information on calibrator standardization, refer to the Bilirubin Calibrator package insert.

QUALITY CONTROL

The following is the recommendation of Abbott Laboratories for quality control. As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- Two levels of controls (normal and abnormal) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, patient values may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

RESULTS

Refer to *Appendix C* of the **ARCHITECT System Operations Manual** for information on results calculations.

Representative performance data are given in the EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert. Results obtained in individual laboratories may vary.

LIMITATIONS OF THE PROCEDURE

Refer to the SPECIMEN COLLECTION AND HANDLING and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

Some specimens may give a direct bilirubin result slightly greater than the total bilirubin result. During internal testing at Abbott Laboratories, specimens with total bilirubin concentrations of 0.2 mg/dL (3.4 µmol/L) or less occasionally gave a direct bilirubin result that slightly exceeded their respective total bilirubin result. This may be observed when nearly all reacting bilirubin is direct bilirubin.

For patients undergoing evaluations involving the administration of indocyanine green (ICG), it is recommended that samples are drawn after ICG has been eliminated. See the Interfering Substances section for additional information.^{10,11}

Abbott Laboratories has not verified the assay performance characteristics with neonatal specimens.

EXPECTED VALUES

Reference Range

Serum¹²

	Range (mg/dL)	Range (μmol/L)
Adult	0.0 to 0.5	0.0 to 8.6

To convert results from mg/dL to μmol/L, multiply mg/dL by 17.1.

A study was conducted using 135 serum samples from volunteers ranging in age from 25 to 66 years. Data were analyzed as described by Clinical and Laboratory Standards Institute (CLSI) protocol NCCLS C28-A.¹³ From this study, 95% of all specimens fell within 0.0 to 0.5 mg/dL, with samples ranging from 0.1 to 0.6 mg/dL.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

SPECIFIC PERFORMANCE CHARACTERISTICS

Linearity

Linearity for Direct Bilirubin is 0.1 to 15.0 mg/dL (1.7 to 256.5 μmol/L), with recovery within 10% or within the 95% confidence level of the predicted value. Linearity was verified using a modified CLSI protocol NCCLS EP6-P.¹⁴ A study performed on an ARCHITECT cSystem produced linear results for Direct Bilirubin up to 16.9 mg/dL (289.0 μmol/L).

Limit of Detection (LOD)

The LOD is the Limit of Absence (LOA*) + 1.645 SD, where SD = the pooled, within-run standard deviation of a low concentration sample. A study performed on an ARCHITECT cSystem produced an LOD for Direct Bilirubin of 0.04 mg/dL (0.69 μmol/L).

* LOA = mean concentration of analyte-free sample + 1.645 SD, where SD = pooled SD of analyte-free sample.

Limit of Quantitation (LOQ)

The LOQ for Direct Bilirubin is ≤ 0.10 mg/dL (1.71 μmol/L). The LOQ is the analyte concentration at which the CV = 20%. Performance studies produced an LOQ of 0.05 mg/dL (0.86 μmol/L).

Interfering Substances

Potential interference in the Direct Bilirubin assay from 62 mg/dL (0.62 g/L) hemoglobin, 125 mg/dL (1.25 g/L) Intralipid, 0.50 mmol/L Indican (indoxyl sulfate), or 6.3 mg/L (8.1 μmol/L) indocyanine green is ≤ 10% or ± 0.1 mg/dL, whichever is greater, at the medical decision level of the analyte.

Interference studies were conducted using a modified CLSI protocol NCCLS EP7-P.¹⁵ Interference effects were assessed by Dose Response and Paired Difference methods, at the medical decision level of the analyte.

Interfering Substance	Interferent	Concentration	N	Target (mg/dL)	Difference from Target (mg/dL)
Hemoglobin	31 mg/dL (0.31 g/L)		5	0.4	-0.1
	62 mg/dL (0.62 g/L)		5	0.4	-0.1
	125 mg/dL (1.25 g/L)		5	0.4	-0.2
	250 mg/dL (2.50 g/L)		5	0.4	-0.2
Human triglyceride	500 mg/dL (5.00 g/L)		5	0.4	-0.2
	519 mg/dL (5.86 mmol/L)		5	0.4	-0.1
Intralipid	1,034 mg/dL (11.68 mmol/L)		5	0.4	+0.3
	125 mg/dL (1.25 g/L)		5	0.4	-0.1
	250 mg/dL (2.50 g/L)		5	0.4	+0.1
Indocyanine Green	500 mg/dL (5.00 g/L)		5	0.4	+0.4
	6.3 mg/L (8.1 μmol/L)		3	0.3	+0.1
	12.5 mg/L (16.1 μmol/L)		3	0.3	+0.3
	18.8 mg/L (24.2 μmol/L)		3	5.1	+0.4
	25.0 mg/L (32.3 μmol/L)		3	5.1	+0.5

Interfering Substances (Continued)

Hemoglobin solutions at the above concentrations were prepared by addition of hemolysate to human serum pools. Human triglyceride solutions at the above concentrations were prepared by mixing an elevated triglyceride human serum pool with a normal triglyceride human serum pool. Intralipid solutions at the above concentrations were prepared by addition of Intralipid to human serum pools.

Taki et al. reported indoxyl sulfate concentrations up to 8.62 mg/dL (0.40 mmol/L), with an average of 3.52 mg/dL (0.17 mmol/L), in 224 hemodialysis (HD) patients.¹⁶ Indoxyl sulfate does not cause significant interference with this direct bilirubin method. Testing at Abbott Laboratories demonstrated that addition of 12.57 mg/dL (0.50 mmol/L) 3-indoxyl sulfate potassium salt to specimens increased the direct bilirubin concentration by a maximum of 0.1 mg/dL.

Indocyanine green solutions at the concentrations listed in the table were prepared by the individual addition of indocyanine green to two pools of plasma, one with a high concentration of bilirubin and one with a low concentration of bilirubin.

Interferences from medications or endogenous substances may affect results.¹⁷

Precision

The imprecision of the Direct Bilirubin assay is ≤ 5% Total CV. Representative data from studies using CLSI protocol NCCLS EP5-A¹⁸ are summarized below.

Control		Level 1	Level 2	Level 3	Level 4
N		80	80	80	80
Mean (mg/dL)		0.41	2.11	3.50	8.59
Within Run	SD	0.01	0.01	0.02	0.05
	%CV	2.1	0.6	0.6	0.5
Between Run	SD	0.00	0.04	0.06	0.14
	%CV	0.0	1.7	1.6	1.6
Between Day	SD	0.01	0.04	0.07	0.16
	%CV	3.6	1.8	2.0	1.9
Total	SD	0.02	0.05	0.09	0.22
	%CV	4.1	2.6	2.7	2.6

Method Comparison

Correlation studies were performed using a modified CLSI protocol NCCLS EP9-A.¹⁹

Serum results from the Direct Bilirubin assay on an AEROSSET System were compared with those from a commercially available methodology based on Jendrassik-Gróf procedure; factored application.

Serum results from the Direct Bilirubin assay on an ARCHITECT cSystem were compared with those from a commercially available methodology based on Jendrassik-Gróf procedure; factored application.

	AEROSSET vs. Comparative Method	ARCHITECT vs. Comparative Method
N	129	107
Y - Intercept	0.21	0.21
Correlation Coefficient	0.995	0.996
Slope	1.08	1.08
Range (mg/dL)	0.1 to 9.9	0.1 to 9.9

Serum results from the Direct Bilirubin assay on an AEROSSET System were compared with the manual Jendrassik-Gróf method.

Serum results from the Direct Bilirubin assay on an ARCHITECT cSystem were compared with the Direct Bilirubin assay ([REF] 8G63-22) on an AEROSSET System.

	AEROSSET vs. Comparative Method	ARCHITECT vs. AEROSSET
N	49	107
Y - Intercept	0.31	-0.01
Correlation Coefficient	0.998	0.999
Slope	0.983	1.01
Range (mg/dL)	-0.1 to 15.6	0.1 to 11.1

BIBLIOGRAPHY

1. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1136.
2. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
4. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
6. Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA: WB Saunders; 1999:1169.
7. Guder WG, Narayanan S, Wisser H, et al. List of Analytes—preanalytical variables. Annex In: *Samples: From the Patient to the Laboratory*. Darmstadt: GIT Verlag; 1996:Annex 8–9.
8. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*, 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
9. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 2nd ed. Washington, DC: AACC Press; 1997:3-88.
10. Donnachie EM, Seccombe DW, Urquhart NI. Indocyanine Green Interference in the Kodak Ektachem Determination of Total Bilirubin. *Clin Chem* 1989; 35(5): 899-900.
11. Meijer DKF, Weert B, Vermeer GA. Pharmacokinetics of Biliary Excretion in Man. VI. Indocyanine Green. *Eur J Clin Pharmacol* 1988; 35: 295-303.
12. Data on file at Abbott Laboratories.
13. Sasse EA, Aziz KJ, Harris EK, et al. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline (C28-A)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1995.
14. Passey RB, Bee DE, Caffo A, et al. *Evaluation of the Linearity of Quantitative Analytical Methods; Proposed Guideline (EP6-P)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1986.
15. Powers DM, Boyd JC, Glick MR, et al. *Interference Testing in Clinical Chemistry; Proposed Guideline (EP7-P)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1986.
16. Taki K, Tsuruta Y, Niwa T. Indoxyl sulfate and atherosclerotic risk factors in hemodialysis patients. *Am J Nephrol* 2007; 27:30–5.
17. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 4th ed. Washington, DC: AACC Press; 1995:3-105–3-108.
18. Kennedy JW, Carey RN, Coolen RB, et al. *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline (EP5-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 1999.
19. Kennedy JW, Carey RN, Coolen RB, et al. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline (EP9-A)*. Wayne, PA: The National Committee for Clinical Laboratory Standards, 1995.

TRADEMARKS

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

ARCHITECT and related brand marks are trademarks of Abbott. Other trademarks are the property of their respective owners



Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

MANUFACTURED FOR

Abbott Laboratories

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Revised March 2022.

©2017, 2022 Abbott Laboratories



ARCHITECT

Direct Bilirubin Serum/Plasma—Conventional and SI Units

Configure assay parameters — General			
<input checked="" type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation			
Assay: BilD	Type: Photometric	Version: †	
Number: 1065			
Run controls for onboard reagents by: Lot			
<input checked="" type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input type="radio"/> Validity checks			
Reaction mode: End up			
Primary		Secondary	
Wavelength: 548 / 660		Read times Main: 31 – 33	
Last required read: 33			
Absorbance range: ___ – ___		Color correction: ___ – ___	
Sample blank type: Self		Blank: 14 – 16	

Configure assay parameters — Results			
<input type="radio"/> Reaction definition <input checked="" type="radio"/> Reagent / Sample <input type="radio"/> Validity checks			
Reagent: BILD0		Reagent volume: 160 40	
Diluent: Saline		Water volume: ___	
Diluent dispense mode: Type 0		Dispense mode: Type 0 Type 0	
Dilution name	Sample	Diluted sample	Dilution factor
STANDARD	5.0	___	1:1.00
___	___	___	___
___	___	___	___
___	___	___	___

Configure assay parameters — Calibration			
<input type="radio"/> Reaction definition <input type="radio"/> Reagent / Sample <input checked="" type="radio"/> Validity checks			
Reaction check: None			
Maximum absorbance variation: ___			

Configure assay parameters — Calibration			
<input type="radio"/> General <input checked="" type="radio"/> Calibration <input type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation			
Assay: BilD		Calibration method: Linear	
<input checked="" type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks			
Calibrator set: Bil		Calibrator level: 0.0	
Blank: Water		Concentration: §	
Cal 1: Bil1		§	
Cal 2: Bil2		§	
Replicates: 3 [Range 1 – 3]			

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators <input checked="" type="radio"/> Volumes <input type="radio"/> Intervals <input type="radio"/> Validity checks			
Calibrator: Bil			
Blank:	Calibrator level	Sample	Diluted sample
Water	5.0	___	___
Cal 1: Bil1	5.0	___	___
Cal 2: Bil2	5.0	___	___

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators <input type="radio"/> Volumes <input checked="" type="radio"/> Intervals <input type="radio"/> Validity checks			
Calibration intervals:			
Full interval: 336		(hours)	
Calibration type:			
Adjust type: None			

Configure assay parameters — Calibration			
<input type="radio"/> Calibrators <input type="radio"/> Volumes <input type="radio"/> Intervals <input checked="" type="radio"/> Validity checks			
Blank absorbance range: ___ – ___			
Span: Blank – Blank			
Span absorbance range: ___ – ___			
Expected cal factor: 0.00			
Expected cal factor tolerance %: 0			

Configure assay parameters — SmartWash				
<input type="radio"/> General <input type="radio"/> Calibration <input checked="" type="radio"/> SmartWash <input type="radio"/> Results <input type="radio"/> Interpretation				
Assay: BilD				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates
R1	BILD0	Detergent A	345	1
R2	BILD0	Detergent A	345	1

Direct Bilirubin Serum/Plasma—Conventional Units

Configure assay parameters — Results			
<input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation			
Assay: BilD		Assay number: 1065	
Dilution default range:		Result units: mg/dL	
Low-Linearity: 0.1^{††}			
High-Linearity: 15.0			
Gender and age specific ranges:			
GENDER	AGE (UNITS)	NORMAL	EXTREME
Either	0 – 130 (Y)	0.0 – 0.5	

Configure result units	
Assay:	BilD
Version:	†
Result units:	mg/dL
Decimal places:	1 [Range 0 – 4]
Correlation factor:	1.0000
Intercept:	0.0000

Direct Bilirubin Serum/Plasma—SI Units

Configure assay parameters — Results			
<input type="radio"/> General <input type="radio"/> Calibration <input type="radio"/> SmartWash <input checked="" type="radio"/> Results <input type="radio"/> Interpretation			
Assay: BilD		Assay number: 1065	
Dilution default range:		Result units: µmol/L	
Low-Linearity: 1.8^{††}			
High-Linearity: 256.5			
Gender and age specific ranges:			
GENDER	AGE (UNITS)	NORMAL	EXTREME
Either	0 – 130 (Y)	0.0 – 8.6	

Configure result units	
Assay:	BilD
Version:	†
Result units:	µmol/L
Decimal places:	1 [Range 0 – 4]
Correlation factor:	1.0000
Intercept:	0.0000

† Due to differences in instrument systems and unit configurations, version numbers may vary.
 § Refer to the concentration specified on calibrator labeling or value sheet. These values are defined on the Configure calibrator set screen.
 †† The linear low value (Low-Linearity) is LOQ rounded up to the number of decimal places defined in the decimal places parameter field.