<i>Ocusil Silikon Yağı</i>	TR
-----------------------------------	----

Açıklama

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Description

Ocusil Silicone Oil is a colourless, homogeneous liquid, comprising ultrapure polydimethyl siloxane, designed for the retinal endo-tamponad after vitrectomy. Ocusil Silicone Oil used particulary in severe cases of retinal detachment for which conventional treatments would fail.

Composition

Contains highly pure sterile silicone oil. Ocusil Silicone Oil is available in six viscosities.

Physical Properties

Viscosity {25 °C} 1000 cSt 1300 cSt 2200 cSt 5000 cSt 5500 cSt 5700 cSt
Refraktive index {25 °C} 1,401-1,405
Appearance Viscous, Transparent, Odorless, liquid
Formula POLIDIMETHYL-SILOXANE (CH₃)₂SiO-(Si(CH₃)₂O)n-Si(CH₃)₃

Instructions For Use

Ocusil Silicone Oil is designed for use by ophthalmologists only. Before use, the pack preserving the product's sterility should be checked to make sure it is intact.

Vials

- Transfer in aseptic conditions the content of the vial into a sterile syringe.
- Connect the Ocusil Silicone Oil syringe to the infusion cannula.
- Inject slowly.

Syringe

- Ocusil Silicone Oil could be inject directly

Ocusil Silicone Oil acts as a mechanical internal tamponade that reattaches the retina to choroid. There is no diffusion or absorption of the silicone out of the eye. The implantation duration of the Ocusil Silicone Oil should be controlled and adapted to the retinal state of the patient. When removing Ocusil Silicone Oil from the eye it must be ensured that the silicone oil is completely removed from the eye.

Indications

The Ocusil Silicone Oil is indicated in prolonged tamponade after surgical treatment for severe retinal detachment, especially:

- Retinal detachment with proliferative vitreoretinopathy,
- Retinal detachment thet complicate a diabetic retinopathy,
- Retinal detachment with giant tears,
- Traumatic retinal detachment,
- Secondary retinal detachment to a viral retinist

Contraindications

The internal use of silicone oil is contraindicated for the following patients:

Are hypersensitive to silicone oil

Silicone oil must not be used in patients who wear an IOL made from silicone

Side effects

As in any surgical operation, surgery concerning the posterior segment using silicone oil has risks that are best evaluated by the surgery. The use of silicone oils can lead to the formation of cataracts and secondary glaucoma. As in any type of surgery, there are risks linked to the material or developments of the initial pathology. Potential complications associated with the implantation of Ocusil Silicone Oil include but are not limited to the following:

Intraoperative Complications

- Ocular hyportony
- Ocular Hypotony

Postoperative Complications

- Choroidal detachment
- Corneal epithelium edema
- Secondary retinal detachment
- Cystoid macular edema
- Infection, endophthalmis
- Closed angel glaucoma
- Keratitis
- Cataract
- Late emulsification

Interactions with other agents are not known

Presentation And Sterilisation

Ocusil Silicone Oil is available in 10/15 ml units in vials and syringes. It is packaged in a pouch, to make it easier to handle in aseptic conditions. It is steam sterilized in its secondary packaging.

Storage Conditions

Ocusil Silicone Oil must be stored at room temperature.

Warnings

Ocusil Silicone Oil is a single use liquid and must not be resterilized. The sterility and purity of the product are only guaranteed for a single use. It should be stored at room temperature, in a dry place and must not be used after the expiration date shown on the package. Glass syringes are potentially fragile due to their material structure, it is recommended to use appropriate application systems. We recommend removing the Ocusil Silicone Oil once the desired effects have been observed. Ocusil Silicone Oil may reduce the absorption of light which could cause damage to the retina in the case of photocoagulation. Ocusil Silicone Oil's refractive index is rather higher than that of intraocular fluid and causes refractive error, notably when using standart calculations for the power of intraocular lenses.

Adverse Reactions

Unexpected adverse reactions and complications related to Ocusil Silicone Oil must be reported to Atakan Dede-Miray Medikal.

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