AMPLATZER™ GUIDEWIRES

Ordering Information

Contents: 1 guidewire

Model/ Reorder No	Diam (in)	Body	Floppy Tip Lgth (cms)	Tip Description	Usable Lgth (cm)
9-GW-001	0.035	Super Stiff	6	7.5 mm, Modified J-tip	260
9-GW-002	0.035	Super Stiff	5	1.5 mm, Modified J-tip	260
9-GW-003	0.035	Super Stiff	20	6 mm, J-tip	300
9-GW-004*	0.035	Soft Tip, Fixed Core	NA	6 mm, J-tip	300

*9-GW-004, also referred to as “Noodlewire,” is a soft tip, flexible guidewire recommended for establishing an arterial-venous loop, facilitating closure of ventricular septal defects.

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
STERILE EO

LATEX

Does not contain
natural rubber
latex components

2022-01-24

CE
2797

Abbott Medical
5050 Nathan Lane North
Plymouth, MN
55442 USA
+1 855 478 5833
+1 651 756 5833

EC REP

Abbott Medical
The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem
Belgium
+32 2 774 68 11

Manufactured By

Lake Region Medical Limited,
Butlersland, New Ross, Co.
Wexford, Ireland.

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SJM Declaration of Conformity

St. Jude Medical (SJM) hereby declares that the following SJM facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/ECC. All supporting documentation is retained under the premises of SJM. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

Manufacturer Address:	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442, USA
European Representative:	St. Jude Medical Coordination Center BVBA The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium
Product Type:	Catheter Guide Wire
Product Name(s):	AMPLATZER Guidewires
Model Number(s):	9-GW-001, 9-GW-002, 9-GW-003, 9-GW-004
Classification:	Class III (Rule 7) Annex II, Section 4 GHTF Class D
GMDN Code(s):	35094
Original CE Mark Date:	22 March 2001 (1-3), 22 Sept 2003 (4, Noodlewire)
EC Certificate No and expiration date:	Certificate No: CE 594293 Expiration Date: 23 Feb 2023
Annex II, Clause 3:	Certificate No: CE 590631 Expiration Date: 23 Feb 2023
Applicable Quality System Standards:	ISO 13485
Notified Body:	BSI Kitemark Court Davy Avenue Knowlhill Milton Keynes MK5 8PP UK
Notified Body Number:	0086

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No.**CE 694955****Issued To:****Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA****In respect of:****AMPLATZER Guidewires**

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...making excellence a habit.™

Page 1 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

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Part Number	Description
9-GW-001	JFC-SS Modified J, Fixed, PTFE Coated, Super Stiff
9-GW-002	J1.5FC-SS, Modified J, Fixed, PTFE Coated, Super Stiff
9-GW-003	J9FC-FS-LLLT Fixed Core, Long (20 cm) PTFE Coated, Finger-Straightenable
9-GW-004	Noodlewire Guidewire

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Certificate History

Date	Reference Number	Action
03 September 2018	8957249	First Issue. Mirror certificate to CE 594293.
Current	8243107	Traceable to NB 0086.

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Eliberat pentru:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
SUA

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A handwritten signature in blue ink, appearing to be "Cristina", written next to the official stamp.

AMPLATZER™ PFO OCCLUDER

**A LANDMARK DEVICE.
AND A TURNING POINT
FOR PFO CLOSURE.**

FIRST TO MARKET, STILL FIRST WORLDWIDE

As the device that created the category, AMPLATZER™ PFO occluder has sustained leadership over decades of use by pursuing clinical evidence—even beyond an initial study end date—to become the first device supported by positive PFO trial results¹. Today, we continue to innovate around advancing patient safety and reducing risk for patients around the world.

A LANDMARK DEVICE. AND A TURNING POINT FOR PFO CLOSURE.

- Industry-leading device, designed for ease of use and effective closure
- Backed by the largest trial with the most extensive patient follow-up
- Demonstrated excellent safety and efficacy¹

INDICATIONS AND USAGE The AMPLATZER™ PFO Occluder is indicated for percutaneous transcatheter closure of a patent foramen ovale (PFO) to reduce the risk of recurrent ischemic stroke in patients, predominantly between the ages of 18 and 60 years, who have had a cryptogenic stroke due to a presumed paradoxical embolism, as determined by a neurologist and cardiologist following an evaluation to exclude known causes of ischemic stroke. See important safety information referenced within.

WHY IS AMPLATZER PFO OCCLUDER RELIED UPON BY THOUSANDS OF PHYSICIANS AROUND THE WORLD?

GLOBALY OVER
100K
IMPLANTS

WE SET THE STANDARD

- Pioneered treatment with a **PFO-specific** device
- **Over 100,000** devices implanted globally

5,810
PATIENT YEARS
OF DATA

WE RAISE THE BAR: THE LANDMARK RESPECT TRIAL¹

- Most extensive patient **follow-up**, **>2X more** than other PFO trials
- Only trial to include patients on anticoagulation therapy, a **real-world cross-section** of patients

990
PATIENTS IMPLANTED
WITH DEVICE
IN RCTS*

WE DEMONSTRATE EXCELLENCE

- **ZERO** device erosions, thrombus, or embolization events in **SIX**** published trials with **990** patients¹⁻⁶
- **94.2% effective** closure rates at 6 months¹

* RCTs=Randomized Clinical Trials

**Patients in device group of each trial implanted with AMPLATZER PFO occluder device: RESPECT = 465, PREMIUM = 119, PC = 191, CLOSE = 121, DEFENSE = 53, PRIMA = 41. RCTs=Randomized Clinical Trials

OFTEN IMITATED, NEVER MATCHED

Industry-leading device. Developed specifically for the treatment of PFO closure.


DESIGN ADVANTAGES THAT MAKE THE DIFFERENCE

DURABLE NITINOL WIRE MESH WITH POLYESTER FABRIC THREAD

Excellent visibility under fluoro

ASYMMETRIC DOUBLE DISC DESIGN

Minimizes material in the Left Atrium*


94%
CLOSURE RATE†
at 6 months in
RESPECT trial¹

INTAGLIO™ WIRE TREATMENT

Designed to reduce nickel leaching⁷

ONE STEP AHEAD: MINIMIZING COMPLEXITY IN CLOSURE

UNIQUE SELF-EXPANDING DISCS

Linked by a short-connecting waist, the discs align to the PFO without an additional “locking” step.

LOW PROFILE DELIVERY

8 F & 9 F introducer sheaths for multiple sizes (18, 25, 35 mm) enables treatment of patients with smaller vasculature.

FULLY RECAPTURABLE AND REPOSITIONABLE DESIGN

Allows confirmation of device placement prior to final release of the device.

*On the most commonly used sizes (25mm and 35mm devices)

†Effective Closure

See important safety information referenced within.

THE TURNING POINT FOR PFO CLOSURE

Three trials published concurrently in the NEJM provide conclusive evidence of the superiority of PFO closure versus medical management in reducing risk of recurrent stroke.

	RESPECT ¹	REDUCE ⁸	CLOSE ⁴
Devices Used	100% AMPLATZER™ PFO Occluder	39% GORE [‡] HELEX, 61% GORE Cardioform	51% AMPLATZER PFO Occluder; 49% other approved PFO devices
Patients	980	664	473
Follow-Up-Patient Years	5,810 (median 5.9 yrs)	2,232 (median 3.2 yrs)	NR* (mean 5.4 yrs)
Anticoagulant Allowed in Control Group?	Yes	No	No
Relative Risk Reduction	62% (Recurrent ischemic stroke of unknown mechanism)	77% (Recurrent ischemic stroke)	97% (Recurrent ischemic stroke)
Effective Closure	94.2% Freedom from >9 bubbles (Evaluated after 6 months)	94.5% Freedom from >25 bubbles (Evaluated after 12 months)	NR*

*Not reported
See important safety information referenced within.

PUBLISHED DATA HIGHLIGHTS EXCELLENT SAFETY

LOW RISK OF DEVICE-RELATED EVENTS

	RESPECT ¹	REDUCE ⁸
Device Embolization/ Dislocation	0	2
Aortic Erosion/ Dissection	0	1
Device Thrombus	0	2

CLOSE Trial data not included due to device- and procedure-related events reported in combination.

LOW RISK OF ATRIAL FIBRILLATION (AF)

RATE (PER 100 PT YRS)	RESPECT ¹	REDUCE ^{8*}
Serious AF**/ Flutter	0.22	0.65
Any AF/ Flutter	0.76	1.90

CLOSE Trial data not included as follow-up patient-years was not reported.

*Rates calculated based on data in final publication.
**In RESPECT, serious AF was adjudicated by an independent board of physicians. In REDUCE, it was determined by the local investigator.

AMPLATZER™ PFO OCCLUDER

Device Specifications

SIZING AND DEVICE SELECTION

AMPLATZER™ PFO Occluder

Model/Reorder Number	Right Atrial Disc Diameter [A]	Left Atrial Disc Diameter [B]	Minimum Recommend Sheath Size
9-PFO-018	18 mm	18 mm	8 F; 45° curve
9-PFO-025	25 mm	18 mm	8 F; 45° curve
9-PFO-035	35 mm	25 mm	9 F; 45° curve

DELIVERY SYSTEM

AMPLATZER™ TorqVue® 45° Delivery System

Model/Reorder Number	Sheath Size	Tip Angle	Sheath Inner Diameter	Sheath Outer Diameter	Usable Length
9-ITV08F45/80	8 F	45°	2.69 mm/0.11 inch	3.45 mm/0.14 inch	80 cm
9-ITV09F45/80	9 F	45°	3.00 mm/0.12 inch	3.81 mm/0.15 inch	80 cm

ANCILLARY PRODUCTS

AMPLATZER™ Guidewire

Model/Reorder Number	Diameter	Body	Tip Description	Usable Length
9-GW-002	0.035 inch	Super Stiff	1.5 mm, Modified J-tip	260 cm

AMPLATZER™ TorqVue® Exchange Systems

Ordering Information – Contents: 1 each delivery sheath, dilator, exchange wire, hemostasis valve, loader, and plastic vise

Model/Reorder Number	Sheath Size	Tip Angle	Sheath Inner Diameter	Sheath Outer Diameter	Usable Length
9-EITV09F45/80	9 F	45°	3.00 mm/0.12 inch	3.81 mm/0.15 inch	80 cm
9-EITV12F45/80	12 F	45°	3.99 mm/0.16 inch	4.80 mm/0.19 inch	80 cm

REFERENCES

1. Saver JL, Carroll JD, Thaler DE, et al. Long-term outcomes of patent foramen ovale closure or medical therapy after stroke. *N Engl J Med* 2017; 377: 1022-32. 2. Tobis J, Charles A, Silbertson D, et al. Prospective, randomized investigation to evaluate incidence of headache reduction in subjects with frequent migraine and PFO using the AMPLATZER PFO occluder to medical management. *J Am Coll Cardiol* 2017; 70:2766-74. 3. Meier B, Kalesan B, Mattle HP, et al. Percutaneous closure of patent foramen ovale in cryptogenic embolism. *N Engl J Med* 2013; 368: 1083 - 91. 4. Mas J-L, Derumeaux G, Guillon B, et al. Patent foramen ovale closure or anticoagulation vs. antiplatelets after stroke. *N Engl J Med* 2017;377:1011-21 and supplementary appendix. 5. Lee PH, Song J-K, Kim JS, et al. Cryptogenic Stroke and High-Risk Patent Foramen Ovale: The DEFENSE-PFO Trial, *Journal of the American College of Cardiology* (2018), doi: 10.1016/j.jacc.2018.02.046. 6. Heinrich P. Mattle, Stefan Evers, David Hildick-Smith, et al. Percutaneous closure of patent foramen ovale in migraine with aura, a randomized controlled trial, *European Heart Journal*, Volume 37, Issue 26, 7 July 2016, Pages 2029–2036. 7. Based on internal lab testing. Data on file at Abbott. 8. Søndergaard L, Kasner SE, Rhodes JF, et al. Patent foramen ovale closure or antiplatelet therapy for cryptogenic stroke. *N Engl J Med* 2017; 377: 1033-42.

For more information about the AMPLATZER™ PFO Occluder or the RESPECT trial, contact your Abbott sales representative, or visit CRYPTOGENICSTROKE.COM.

INDICATIONS AND IMPORTANT SAFETY INFORMATION

INDICATIONS AND USAGE The AMPLATZER™ PFO Occluder is indicated for percutaneous transcatheter closure of a patent foramen ovale (PFO) to reduce the risk of recurrent ischemic stroke in patients, predominantly between the ages of 18 and 60 years, who have had a cryptogenic stroke due to a presumed paradoxical embolism, as determined by a neurologist and cardiologist following an evaluation to exclude known causes of ischemic stroke. **CONTRAINDICATIONS** Patients with intra-cardiac mass, vegetation, tumor or thrombus at the intended site of implant, or documented evidence of venous thrombus in the vessels through which access to the PFO is gained. • Patients whose vasculature, through which access to the PFO is gained, is inadequate to accommodate the appropriate sheath size. • Patients with anatomy in which the AMPLATZER™ PFO device size required would interfere with other intracardiac or intravascular structures, such as valves or pulmonary veins. • Patients with other source of right-to-left shunts, including an atrial septal defect and/or fenestrated septum. • Patients with active endocarditis or other untreated infections. **WARNINGS** Patients who are at increased risk for venous thromboembolic events should be managed with thromboembolic risk reduction regimen after the PFO Closure following standard of care. • Do not use this device if the sterile package is open or damaged. • Prepare for situations that require percutaneous or surgical removal of this device. This includes availability of a surgeon. • Embolized devices must be removed as they may disrupt critical cardiac functions. Do not remove an embolized occluder through intracardiac structures unless the occluder is fully recaptured inside a catheter or sheath. • Patients who are allergic to nickel can have an allergic reaction to this device. • This device should be used only by physicians who are trained in standard transcatheter techniques. • Transient hemodynamic compromise may be encountered during device placement, which may require fluid replacement or other medications as determined by the physician. • Do not release the device from the delivery cable if the device does not conform to its original configuration, or if the device position is unstable or if the device interferes with any adjacent cardiac structure (such as Superior Vena Cava (SVC), Pulmonary Vein (PV), Mitral Valve (MV), Coronary Sinus (CS), aorta (AO)). If the device interferes with an adjacent cardiac structure, recapture the device and redeploy. If still unsatisfactory, recapture the device and either replace with a new device or refer the patient for alternative treatment. • Ensure there is sufficient distance from the PFO to the aortic root or SVC (typically defined as 9 mm or greater as measured. **PRECAUTIONS** The safety and effectiveness of the AMPLATZER™ PFO Occluder has not been established in patients (with): • Age less than 18 years or greater than 60 years because enrollment in the pivotal study (the RESPECT trial) was limited to patients 18 to 60 years old • A hypercoagulable state including those with a positive test for an anticardiolipin antibody (IgG or IgM), Lupus anticoagulant, beta-2 glycoprotein-1 antibodies, or persistently elevated fasting plasma homocysteine despite medical therapy • Unable to take antiplatelet therapy • Atherosclerosis or other arteriopathy of the intracranial and extracranial vessels associated with a ≥50% luminal stenosis • Acute or recent (within 6 months) myocardial infarction or unstable angina • Left ventricular aneurysm or akinesis • Mitral valve stenosis or severe mitral regurgitation irrespective of etiology • Aortic valve stenosis (mean gradient greater than 40 mmHg) or severe aortic valve regurgitation • Mitral or aortic valve vegetation or prosthesis • Aortic arch plaques protruding greater than 4 mm into the aortic lumen • Left ventricular dilated cardiomyopathy with left ventricular ejection fraction (LVEF) less than 35% • Chronic, persistent, or paroxysmal atrial fibrillation or atrial flutter • Uncontrolled hypertension or uncontrolled diabetes mellitus • Diagnosis of lacunar infarct probably due to intrinsic small vessel as qualifying stroke event • Arterial dissection as cause of stroke n Index stroke of poor outcome (modified Rankin score greater than 3) n Pregnancy at the time of implant • Multi-organ failure • Use on or before the last day of the expiration month that is printed on the product packaging label. • This device was sterilized with ethylene oxide and is for single use only. Do not reuse or re-sterilize this device. Attempting to re-sterilize this device can cause a malfunction, insufficient sterilization, or harm to the patient. • The AMPLATZER™ PFO Occluder device consists of a nickel-titanium alloy, which is generally considered safe. However, in vitro testing has demonstrated that nickel is released from this device for a minimum of 60 days. Patients who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may also develop an allergy to nickel if this device is implanted. • Store in a dry place. • Pregnancy – Minimize radiation exposure to the fetus and the mother. • Nursing mothers – There has been no quantitative assessment for the presence of leachables in breast milk. **ADVERSE EVENTS** Potential adverse events that may occur during or after a procedure using this device may include, but are not limited to: Air embolus; Allergic drug reaction; Allergic dye reaction; Allergic metal reaction: Nitinol (nickel, titanium), platinum/iridium, stainless steel (chromium, iron, manganese, molybdenum, nickel); Anesthesia reactions; Apnea; Arrhythmia; Bacterial endocarditis; Bleeding; Brachial plexus injury; Cardiac perforation; Cardiac tamponade; Cardiac thrombus; Chest pain; Device embolization; Device erosion; Deep vein thrombosis; Death; Endocarditis; Esophagus injury; Fever; Headache/migraine; Hypertension/hypotension; Myocardial infarction; Pacemaker placement secondary to PFO device closure; Palpitations; Pericardial effusion; Pericardial tamponade; Pericarditis; Peripheral embolism; Pleural effusion; Pulmonary embolism; Reintervention for residual shunt/device removal; Sepsis; Stroke; Transient ischemic attack; Thrombus; Valvular regurgitation; Vascular access site injury; Vessel perforation.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at manuals.sjm.com for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Abbott Structural Heart

3200 Lakeside Dr., Santa Clara, CA. 95054 USA, Tel: 1.800.227.9902

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

©2018 Abbott. All rights reserved. SJM-AMPLP-0418-0090 | Item approved for US use only.

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A handwritten signature in blue ink, appearing to be "Cristina", written in a cursive style.

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SJM Declaration of Conformity

St. Jude Medical (SJM) hereby declares that the following SJM facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/EEC. All supporting documentation is retained under the premises of SJM. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

Manufacturer Address:

AGA Medical Corporation
5050 Nathan Lane North
Plymouth, Minnesota 55442, USA

European Representative:

St. Jude Medical Coordination Center BVBA
The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem, Belgium

Product Type:

Cardiac Occluder

Product Name(s):

AMPLATZER Septal Occluder
AMPLATZER Multi-Fenestrated Septal Occluder -
"Cribiform"
AMPLATZER PFO Occluder

Model Number(s):

Product	Model Numbers			Original CE Mark Date
AMPLATZER Septal Occluder	9-ASD-004	9-ASD-013	9-ASD-024	24 February 1998
	9-ASD-005	9-ASD-014	9-ASD-026	
	9-ASD-006	9-ASD-015	9-ASD-028	
	9-ASD-007	9-ASD-016	9-ASD-030	
	9-ASD-008	9-ASD-017	9-ASD-032	
	9-ASD-009	9-ASD-018	9-ASD-034	
	9-ASD-010	9-ASD-019	9-ASD-036	
	9-ASD-011	9-ASD-020	9-ASD-038	
	9-ASD-012	9-ASD-022	9-ASD-040	
AMPLATZER Multi-Fenestrated Septal Occluder – "Cribiform"	9-ASD-MF-018	9-ASD-MF-035		2 September 2002
	9-ASD-MF-025	9-ASD-MF-040		
	9-ASD-MF-030			
AMPLATZER PFO Occluder	9-PFO-018	9-PFO-025		24 February 1998
	9-PFO-030	9-PFO-035		



Signature:

Lisa Becker
Senior Director, Regulatory Affairs

Issue Date

SJM Declaration of Conformity

Classification:	Class III per Annex II, Rule 8
GMDN Code(s):	45418
EC Design Certificate No and Expiration Date:	Certificate No: CE 594291 Expiration Date: 23 February 2023
Annex II, Clause 3 Certificate No and Expiration Date:	Certificate No: CE 590631 Expiration Date: 23 February 2023
Applicable Quality System Standards:	ISO 13485
Notified Body:	BSI Kitemark Court Davy Avenue Knowlhill Milton Keynes MK5 8PP UK
Notified Body Number:	0086
AMPLATZER Septal Occluder Manufacturing Facilities:	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA St. Jude Medical, Costa Rica Ltda. Edificio #44, Calle 0, Avenida 2 Zona Franca Coyol, El Coyol, Alajuela, Costa Rica
AMPLATZER Multi-Fenestrated Septal Occluder - "Cribiform" Manufacturing Facility:	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA
AMPLATZER PFO Occluder Manufacturing Facility:	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA
Signature:  Lisa Becker Senior Director, Regulatory Affairs	 Issue Date

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No.**CE 594291****Issued To:**

AGA Medical Corporation
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

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For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0086):

Stewart Brain, Head of Compliance & Risk -
Medical Devices

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Supplementary Information to CE 594291

Issued To:

AGA Medical Corporation
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

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Supplementary Information to CE 594291

Issued To:

AGA Medical Corporation
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Device Size (mm)	Model Number	Wire Diameter (mm)/Number of Wires	Distal Disc Diameter (mm)	Proximal Disc Diameter (mm)	Waist Diameter	Waist Length
17	9-ASD-017	.005/72	31	27	17	4
18	9-ASD-018	.006/72	32	28	18	4
19	9-ASD-019	.006/72	33	29	19	4
20	9-ASD-020	.006/72	34	30	20	4
22	9-ASD-022	.006/72	36	32	22	4
24	9-ASD-024	.006/72	38	34	24	4
26	9-ASD-026	.007/72	40	36	26	4
28	9-ASD-028	.007/72	42	38	28	4
30	9-ASD-030	.007/72	44	40	30	4
32	9-ASD-032	.007/72	46	42	32	4
34	9-ASD-034	.008/72	50	44	34	4
36	9-ASD-036	.008/72	52	46	36	4
38	9-ASD-038	.008/72	54	48	38	4
40	9-ASD-040	.008/72	56	50	40	4

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Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000
BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.
A member of BSI Group of Companies.

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Supplementary Information to CE 594291

Issued To:

AGA Medical Corporation
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer Multi-Fenestrated Septal Occluder – “Cribriform”:

Device Size (mm)	Model Number	Wire Diameter (mm)/Number of Wires	Disc Diameter (mm)	Diameter of Raised Disc (mm)
18	9-ASD-MF-018	.005/72	18	6
25	9-ASD-MF-025	.005/72	25	13
30	9-ASD-MF-030	.005/72	30	30
35	9-ASD-MF-035	.006/72	35	23
40	9-ASD-MF-040	.006/72	40	40

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Supplementary Information to CE 594291

Issued To:

AGA Medical Corporation
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Certificate History

Date	Reference Number	Action
24 February 2013	10139223	First Issue – Transfer from another Notified Body.
28 May 2013	10141409	Change the surface finish of the nitinol wire from black oxide to chemically etched.
03 February 2014	10144456	Update to Amplatzer Septal Occluder IFU. There are no changes to the intended use.
06 February 2015	10152724	Addition of St. Jude Medical Costa Rica Ltda. as an alternate manufacturing site for the Amplatzer Septal Occluder. Addition of Synergy Health AST, SRL as an alternate sterilization site for the Amplatzer Septal Occluder.
25 August 2015	10156926	DuPont Tyvek Medical Transition Project update.
01 February 2016	10160623	Addition of Sterigenics Willowbrook, IL as a sterilizer.
Current	8872315	Certificate renewal.

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A handwritten signature in blue ink, appearing to be "Dr.", located to the right of the official stamp.

AMPLATZER™ TORQVUE™ 45° AND 180° DELIVERY SYSTEMS

Ordering Information

Contents: 1 delivery cable, delivery sheath, dilator, hemostasis valve, loader and plastic vise

Model/ Reorder No	Sheath Sz	Tip Angle	Inner Diam of Sheath (mm/in)	Outer Diam of Sheath (mm/in)	Usable Lgth (cm)	Delivery Cable Lgth (cm)
9-ITV06F45/60	6 F	45°	2.11/0.08	2.79/0.11	60	110
9-ITV07F45/60	7 F	45°	2.44/0.10	3.18/0.13	60	110
9-ITV07F45/80	7 F	45°	2.44/0.10	3.18/0.13	80	120
9-ITV08F45/60	8 F	45°	2.69/0.11	3.45/0.14	60	110
9-ITV08F45/80	8 F	45°	2.69/0.11	3.45/0.14	80	120
9-ITV09F45/80	9 F	45°	3.00/0.12	3.81/0.15	80	120
9-ITV10F45/80	10 F	45°	3.30/0.13	4.14/0.16	80	120
9-ITV12F45/80	12 F	45°	3.99/0.16	4.80/0.19	80	120
9-ITV13F45/80	13 F	45°	4.32/0.17	5.13/0.20	80	120
9-ITV05F180/60	5 F	180°	1.83/0.07	2.51/0.10	60	110
9-ITV06F180/60	6 F	180°	2.11/0.08	2.79/0.11	60	110
9-ITV06F180/80	6 F	180°	2.11/0.08	2.79/0.11	80	120
9-ITV07F180/80	7 F	180°	2.44/0.10	3.18/0.13	80	120
9-ITV08F180/80	8 F	180°	2.69/0.11	3.45/0.14	80	120
9-ITV09F180/80	9 F	180°	3.00/0.12	3.81/0.15	80	120

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A handwritten signature in blue ink, appearing to be "Cristina", written in a cursive style to the right of the official stamp.

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QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

St. Jude Medical
177 County Road B East
St Paul
Minnesota
55117
USA

Holds Certificate No:

FM 558476

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design, development and Manufacture of and finished Mechanical Heart Valves, Tissue Heart Valves, Annuloplasty Rings, Valve and Annuloplasty Ring Sizer Sets, Mechanical Valve Leaflet Testers, Holder Rotators, Transcatheter Heart Valve Delivery and Loading Systems and related accessories along with manufacturing of intermediate components used in other medical devices.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2009-12-24

Latest Revision Date: 2020-01-13

Effective Date: 2020-02-29

Expiry Date: 2023-02-28

Page: 1 of 2

...making excellence a habit.™

Certificate No: **FM 558476**

Location	Registered Activities
St. Jude Medical 177 County Road B East St Paul Minnesota 55117 USA	Design, development and Manufacture of and finished Mechanical Heart Valves, Tissue Heart Valves, Annuloplasty Rings, Valve and Annuloplasty Ring Sizer Sets, Mechanical Valve Leaflet Testers, Holder Rotators, Transcatheter Heart Valve Delivery and Loading Systems and related accessories along with manufacturing of intermediate components used in other medical devices.
St. Jude Medical Brasil Ltda. Rua Professor Jose Vieira de Mendonca 1301 Bairro Engenho Nogueira Pampulha, Belo Horizonte Minas Gerais 31.310-026 Brasil	The Manufacture and final inspection of tissue made heart valves and tissue vascular prostheses.

Original Registration Date: 2009-12-24

Latest Revision Date: 2020-01-13

Effective Date: 2020-02-29

Expiry Date: 2023-02-28

Page: 2 of 2

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Abbott Medical Declaration of Conformity

Abbott Medical (Abbott) hereby declares that the following Abbott facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/EEC. All supporting documentation is retained under the premises of Abbott. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

Manufacturer Address:

Abbott Medical
5050 Nathan Lane North
Plymouth, Minnesota 55442, USA

European Representative:


Abbott Medical
The Corporate Village
Da Vincilaan 11 Box F1
1935 Zaventem, Belgium

Product Type:

Cardiac Occluder Delivery Kit

Product Name(s):

Amplatzer TorqVue Delivery System
Amplatzer TorqVue Exchange System
Amplatzer TorqVue 2 Delivery Sheath
Amplatzer TorqVue LP Delivery System
Amplatzer TorqVue LP Catheter
Amplatzer TorqVue Delivery System with Pusher Catheter
Amplatzer TorqVue 45°x45° Delivery Sheath
Amplatzer Amulet Delivery Sheath
Amplatzer Trevisio Intravascular Delivery System

Signature:

Michelle Grossman
Director, Regulatory Affairs


April 21, 2020
Issue Date

Abbott Medical Declaration of Conformity

Model Number(s):

Product	Model Numbers	Original CE Mark Date
Amplatzer TorqVue Delivery System	9-ITV06F45/60 9-ITV13F45/80 9-ITV07F45/60 9-ITV05F180/60 9-ITV07F45/80 9-ITV06F180/60 9-ITV08F45/60 9-ITV06F180/80 9-ITV08F45/80 9-ITV07F180/80 9-ITV09F45/80 9-ITV08F180/80 9-ITV10F45/80 9-ITV09F180/80 9-ITV12F45/80	10 October 2005
Amplatzer TorqVue Exchange System	9-EITV09F45/80 9-EITV12F45/80 9-EITV06F180/80 9-EITV08F180/80 9-EITV10F180/80	10 October 2005
Amplatzer TorqVue 2 Delivery Sheath	9-TV2-05F120 9-TV2-06F120 9-TV2-07F120	19 February 2010
Amplatzer TorqVue LP Delivery System	9-TVLP4F90/060 9-TVLP4F90/080 9-TVLP5F90/060 9-TVLP5F90/080	07 February 2008
Amplatzer TorqVue LP Catheter	9-TVLPC4F90/080	28 April 2011
Amplatzer TorqVue Delivery System with Pusher Catheter	9-ITVP07F-180/80 9-ITVP08F-180/80 9-ITVP09F-180/80	21 June 2011
Amplatzer TorqVue 45°x45° Delivery Sheath	9-TV45X45-09F-100 9-TV45X45-10F-100 9-TV45X45-12F-100 9-TV45X45-13F-100 9-TV45X45-14F-100	03 December 2008 (9-13 Fr) 24 February 2012 (14 Fr)
Amplatzer Amulet Delivery Sheath	DS-TV45X45-12F-080 DS-TV45X45-14F-080	08 February 2017
Amplatzer Trevisio Intravascular Delivery System	9-ATV06F45/60 9-ATV09F45/80 9-ATV07F45/60 9-ATV10F45/80 9-ATV07F45/80 9-ATV12F45/80 9-ATV08F45/60 9-ATV13F45/80 9-ATV08F45/80	20 April 2020

Signature:

 Michelle Grossman
 Director, Regulatory Affairs

April 21 2020

 Issue Date

Abbott Medical Declaration of Conformity

Classification:	Class III (Rule 7) Annex II, Section 4 GHTF Class D
GMDN Code(s):	45419
EC Design Certificate No and Expiration Date:	Certificate No: CE 694956 Expiration Date: 23 February 2023
Annex II, Clause 3 Certificate No and Expiration Date:	Certificate No: CE 694788 Expiration Date: 23 February 2023
Applicable Quality System Standards:	ISO 13485
Notified Body:	BSI Group The Netherlands B.V. Say Building John M. Keynesplein 9 1066 EP Amsterdam The Netherlands
Notified Body Number:	2727

Signature:

Michelle Grossman
Director, Regulatory AffairsApril 21, 2020
Issue Date

EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

No.**CE 694956****Issued To:**

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

In respect of:**Amplatzer Delivery Systems**

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):

Gary E Slack, Senior Vice President Medical Devices

First Issued: 2018-09-03**Date: 2020-04-20****Expiry Date: 2023-02-23**

...making excellence a habit.™

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ TorqVue™ Delivery System			
Intended purpose per IFU:			
The Amplatzer™ TorqVue™ Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices.			
Classification: Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ITV06F45/60	6	45°	60
9-ITV07F45/60	7	45°	60
9-ITV07F45/80	7	45°	80
9-ITV08F45/60	8	45°	60
9-ITV08F45/80	8	45°	80
9-ITV09F45/80	9	45°	80
9-ITV10F45/80	10	45°	80
9-ITV12F45/80	12	45°	80
9-ITV13F45/80	13	45°	80
9-ITV05F180/60	5	180°	60
9-ITV06F180/60	6	180°	60
9-ITV06F180/80	6	180°	80
9-ITV07F180/80	7	180°	80
9-ITV08F180/80	8	180°	80
9-ITV09F180/80	9	180°	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 2 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ TorqVue™ Exchange System			
Intended purpose per IFU: The Amplatzer™ TorqVue™ Exchange System is intended for removal of an Amplatzer™ Delivery Sheath and subsequent exchange for an Amplatzer™ Delivery Sheath of equal or larger diameter.			
Classification: Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-EITV09F45/80	9	45°	80
9-EITV12F45/80	12	45°	80
9-EITV06F180/80	6	180°	80
9-EITV08F180/80	8	180°	80
9-EITV10F180/80	10	180°	80

Amplatzer™ TorqVue™ 2 Delivery Sheath			
Intended purpose per IFU: The Amplatzer™ TorqVue™ 2 Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers and coronary vasculature of the heart or in the peripheral vasculature.			
Classification: Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-TV2-05F120	5	none	120
9-TV2-06F120	6	none	120
9-TV2-07F120	7	none	120

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 3 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
 A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ TorqVue™ Delivery System with Pusher Catheter			
Intended purpose per IFU: The Amplatzer™ TorqVue™ Delivery System with Pusher Catheter is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Membranous VSD Occluder device.			
Classification: Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ITVP07F180/80	7	180°	80
9-ITVP08F180/80	8	180°	80
9-ITVP09F180/80	9	180°	80

Amplatzer™ TorqVue™ LP Delivery System				
Intended purpose per IFU: The Amplatzer™ TorqVue™ LP Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of the Amplatzer™ devices.				
Classification: Class III				
Catalogue Number	Model, Type			
	Device Size (Fr)	Curve Dimension	Length (cm)	Delivery Wire Length (cm)
9-TVLP4F90/060	4	90°	60	160
9-TVLP4F90/080	4	90°	80	195
9-TVLP5F90/060	5	90°	60	160
9-TVLP5F90/080	5	90°	80	195

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 4 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ Trevisio™ Intravascular Delivery System

Intended purpose per IFU:

The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices.

Classification: Class III

Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ATV06F45/60	6	45°	60
9-ATV07F45/60	7	45°	60
9-ATV07F45/80	7	45°	80
9-ATV08F45/60	8	45°	60
9-ATV08F45/80	8	45°	80
9-ATV09F45/80	9	45°	80
9-ATV10F45/80	10	45°	80
9-ATV12F45/80	12	45°	80
9-ATV13F45/80	13	45°	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 5 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
 A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ TorqVue™ LP Catheter				
Intended purpose per IFU: The TorqVue™ LP Catheter is intended to facilitate the loading, delivery, and deployment of Amplatzer™ devices.				
Classification: Class III				
Catalogue Number	Model, Type			
	Device Size (Fr)	Usable Length (cm)	Tip Outer Diameter mm (in)	Tip Inner Diameter mm (in)
9-TVLPC4F90/080	4	80	1.40 (.055)	1.17 (0.046)

Amplatzer™ TorqVue™ 45x45 Delivery Sheath		
Intended purpose per IFU: The Amplatzer™ TorqVue™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.		
Classification: Class III		
Catalogue Number	Model, Type	
	Sheath Size (Fr)	Length (cm)
9-TV45X45-09F-100	9	100
9-TV45X45-10F-100	10	100
9-TV45X45-12F-100	12	100
9-TV45X45-13F-100	13	100
9-TV45X45-14F-100	14	100

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 6 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Amplatzer™ Amulet™ Delivery Sheath		
Intended purpose per IFU: The Amplatzer™ Amulet™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.		
Classification: Class III		
Catalogue Number	Model, Type	
	Sheath Size (Fr)	Length (cm)
DS-TV45X45-12F-080	12	80
DS-TV45X45-14F-080	14	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 7 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.
A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 694956

Issued To:

Abbott Medical
5050 Nathan Lane North
Plymouth
Minnesota
55442
USA

Certificate History

Date	Reference Number	Action
03 September 2018	8957249	First Issue. Mirror certificate to CE 594294.
20 February 2019	8243107	Traceable to NB 0086.
20 March 2019	9738457	Addition of Sterigenics US LLC, Salt Lake City, Utah for ETO Sterilization.
16 December 2019	3053900	Addition of Midwest Sterilization Corporation, Jackson, Missouri USA for ETO Sterilization in chambers 1, 2, 3, 6, and 13.
Current	9784335	Extension to scope to include the Amplatzer™ Trevisio™ Intravascular Delivery System. Revision to scope statement to remove "TorqVue" brand name. Reformat product tables. Correction of delivery wire length for Amplatzer™ TorqVue™ LP Delivery System. Correction of tip outer diameter for Amplatzer™ TorqVue™ LP Catheter.

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 8 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.
This certificate was issued electronically and is bound by the conditions of the contract.

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Abbott Medical
5050 Nathan Lane North
Plymouth, Minnesota 55442, SUA

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BSI Group The Netherlands B.V. înregistrată în Olanda cu numărul 33264284.