Product data sheet

GENERAL INFORMATION

Product: Vacuum manifold for Solid Phase Extraction SPE

Description: Compatible with SPHE and standard luer-end SPE cartridges. Made of transparent glass, high resistance to corrosion. Withstands more than 96 kPa (720 mm Hg) vacuum pressure with no deformation, long-term use. Individual valves allow vacuum control for each cartridge improving reproducibility. Joint highly resistant to corrosion and PTFE internal tube rack. Includes gauge, lid, collection rack, solvent guide needles and silicone tube

SPECIFICATIONS

Reference	pcs/ box	description
MNFL-S12-001	1	12 positions, square

PACKING

Type: Carboard box

Label:

BRANCHIA

MNFL-S12-001

Vacuum manifold for Solid Phase Extraction SPE, 1 pcs

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Batch n°.: xxxx

PROCESS

The manifold for the extraction in solid phase is used to carry out the process of the same name in order to concentrate a substance that is very diluted within the sample to be analyzed. Next we explain the Solid Phase Extraction:

ACTORS:

Absorbent: It is the material inside the SPE cartridge, solid in nature, that retains the analyte of the sample by passing it through it

Sample: Matrix from which we want to analyze a substance

Analyte: Substance to be analyzed within the matrix

Interferences: Substances that interfere in some way in the analyte, either interfering directly with it, or in any of the processes in which the analyte is involved


STEPS:

Conditioning: In some cases, you need to pass some type of solution through the solid absorber to activate it. Without this step the absorbent is not able to perform the following steps:

Loading: The sample is passed through the solid absorber to retain a few substances, including the analyte.

Washing: A solvent is passed through the absorber to remove all previously absorbed substances that interfere in some way in the analytical, the interference.

Elution: A solvent is passed through the absorbent to extract the analyte from the absorber and can be dissolved in a sample with a lower volume than the original sample.

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