## Pipeline™ Vantage Embolization Device with Shield Technology™

<b>fr</b>	Dispositif d'embolisation Pipeline™ Vantage avec Shield Technology™
<b>de</b>	Pipeline™ Vantage Embolisationsimplantat mit Shield Technology™
<b>it</b>	Dispositivo di embolizzazione Pipeline™ Vantage con Shield Technology™
<b>es</b>	Dispositivo de embolización Pipeline™ Vantage con Shield Technology™
<b>sv</b>	Pipeline™ Vantage emboliseringsinstrument med Shield Technology™
<b>nl</b>	Pipeline™ Vantage-embolisatiehulpmiddel met Shield Technology™
<b>pt</b>	Dispositivo de embolização Pipeline™ Vantage com Shield Technology™
<b>fi</b>	Shield Technology™ -käsitelty Pipeline™ Vantage -embolisaatiolaite
<b>da</b>	Pipeline™ Vantage-emboliseringsanordning med Shield Technology™
<b>el</b>	Συσκευή εμβολισμού Pipeline™ Vantage με Shield Technology™
<b>cs</b>	Zařízení pro embolizaci Pipeline™ Vantage s Shield Technology™
<b>hu</b>	Shield Technology™ alkalmazásával készült Pipeline™ Vantage embolizációs eszköz
<b>ru</b>	Устройство для эмболизации Pipeline™ Vantage с технологией Shield Technology™
<b>pl</b>	Urządzenie do embolizacji Pipeline™ Vantage wyposażone w technologię Shield Technology™
<b>tr</b>	Shield Technology™ özellikli Pipeline™ Vantage Embolizasyon Cihazı
<b>no</b>	Pipeline™ Vantage-emboliseringsenhet med Shield Technology™
<b>sk</b>	Embolizačné zariadenie Pipeline™ Vantage s technológiou Shield Technology™
<b>ro</b>	Dispozitiv de embolizare Pipeline™ Vantage cu Shield Technology™
<b>bg</b>	Устройство за емболизация Pipeline™ Vantage с технология Shield Technology™
<b>et</b>	Shield Technology™ tehnoloogiaga embolisatsiooniseade Pipeline™ Vantage
<b>hr</b>	Medicinski proizvod za embolizaciju Pipeline™ Vantage s tehnologijom Shield Technology™
<b>lt</b>	„Pipeline™ Vantage“ embolizavimo prietaisas su „Shield Technology™“
<b>lv</b>	Embolizācijas ierīce Pipeline™ Vantage ar tehnoloģiju Shield Technology™
<b>sl</b>	Pripomoček za embolizacijo Pipeline™ Vantage s tehnologijo Shield Technology™
<b>sr</b>	Pipeline™ Vantage sredstvo za embolizaciju sa Shield Technology™
<b>pt-br</b>	Dispositivo de embolização Pipeline™ Vantage com Shield Technology™
<b>ko</b>	Shield Technology™를 이용한 Pipeline™ Vantage 색전술 기구
<b>vi</b>	Thiết bị thuyên tắc mạch máu Pipeline™ Vantage với Shield Technology™
<b>kk</b>	Shield Technology™ технологиясы бар Pipeline™ Vantage эмболизациялау құрылғысы
<b>ar</b>	Shield Technology™ جهاز الانصمام Pipeline™ Vantage المزود بـ

## TABLE OF CONTENTS

# *Pipeline™ Vantage Embolization Device with Shield Technology™*

---

English.....	3
Français.....	6
Deutsch.....	9
Italiano.....	12
Español.....	15
Svenska.....	18
Nederlands.....	21
Português.....	24
Suomi.....	27
Dansk.....	30
Ελληνικά.....	33
Česky.....	36
Magyar.....	39
Русский.....	42
Polski.....	45
Türkçe.....	48
Norsk.....	51
Slovenčina.....	54
Română.....	57
Български.....	60
Eesti.....	63
hrvatski.....	65
Lietuvių.....	68
Latviski.....	71
Slovenščina.....	74
Srpski.....	77
Português (Brasil).....	80
한국어.....	83
Tiếng Việt.....	87
Қазақша.....	90
العربية.....	96
<b>Symbol Glossary.....</b>	<b>97</b>

## Instructions for Use

## Pipeline™ Vantage Embolization Device with Shield Technology™

### CAUTION

- This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

### DESCRIPTION

The Pipeline™ Vantage Embolization Device with Shield Technology™ consists of a permanent implant combined with a guidewire-based delivery system. The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is a braided, multi-alloy, mesh cylinder woven with cobalt-chromium-nickel and platinum wires. An image of the Pipeline™ Vantage Embolization Device with Shield Technology™ implant is shown in Figure 1 and the design of the device is shown in Figure 2. The woven wires of the device provide approximately 30% metal coverage of the arterial wall surface area. The implant is designed for placement in a parent vessel across the neck of an intracranial aneurysm (IA). The expanded or unconstrained diameter is 0.25 mm larger than the labeled diameter. Shield Technology™ is a surface-modification that is not derived from any animal or human sources.

The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is assembled on a guide-wire based delivery system that consists of a 304-stainless steel core wire and a 304L stainless steel hypotube. The implant is assembled over 304 stainless steel resheathing components. A Platinum-Iridium Restraint is distal to the resheathing components and is termed the Resheathing Marker. Refer to Figure 3 for the Resheathing Marker position.

The tip coil is made of platinum-tungsten alloy. The tip, distal, and proximal solder joints are a tin-silver. The ePTFE protective sleeves cover and protect the distal portion of the braid while the Pipeline™ Vantage Embolization Device with Shield Technology™ implant is advanced through the micro catheter. The Resheathing components allow the user to resheath the implant back into the micro catheter. The Resheathing Marker provides the user fluoroscopic visualization for the limit of resheathing the implant.

The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is compressed inside an introducer sheath. The Pipeline™ Vantage Embolization Device with Shield Technology™ implant is designed to be delivered through a compatible micro catheter of either 0.021 inch (0.53 mm) or 0.027 inch (0.69 mm) inside diameter and minimum 135 cm in length. Refer to Table 1 for micro catheter compatibility for each device size.


Figure 1. The Pipeline™ Vantage Embolization Device with Shield Technology™


Figure 2. The Pipeline™ Vantage Embolization Device with Shield Technology™ delivery system and implant (not to scale)

- |                      |                      |
|----------------------|----------------------|
| 1. Tip Coil          | 5. Fluorosafe Marker |
| 2. Proximal Bumper   | 6. Delivery Wire     |
| 3. Introducer Sheath | 7. Distal Marker     |
| 4. Braid             | 8. ePTFE Sleeves     |

Labeled Diameter (mm)	Compatible catheter inner diameter	Labeled Lengths (mm)
2.50	0.021 inch (0.53 mm)	10, 12, 14, 16, 18, 20
2.75		10, 12, 14, 16, 18, 20
3.00		10, 12, 14, 16, 18, 20, 25
3.25		10, 12, 14, 16, 18, 20, 25
3.50	0.027 inch (0.69 mm)	10, 12, 14, 16, 18, 20, 25, 30, 35
4.00		10, 12, 14, 16, 18, 20, 25, 30, 35, 40
4.50		10, 12, 14, 16, 18, 20, 25, 30, 35, 40
5.00		10, 12, 14, 16, 18, 20, 25, 30, 35, 40
5.50		10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50
6.00		10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50

### DEVICE COMPATIBILITY

Micro catheter compatibility is defined on the product label:

The Pipeline™ Vantage 021 system is designed to be delivered through a compatible microcatheter of 0.021 inch (0.53 mm) inside diameter at least 135 cm in length. Compatibility testing has been performed with the Phenom 21 Catheter.

The Pipeline™ Vantage 027 system is designed to be delivered through a compatible micro catheter of 0.027 inch (0.69 mm) inside diameter at least 135 cm in length. Compatibility testing has been performed with the Phenom 27 Catheter.

### INTENDED PURPOSE / INDICATIONS FOR USE

The Pipeline™ Vantage Embolization Device with Shield Technology™ is intended for endovascular embolization of cerebral aneurysms.

### CONTRAINDICATIONS

- Patients with active bacterial infection.
- Patients in whom antiplatelet therapy (i.e. aspirin and clopidogrel) is contraindicated.
- Patients who have not received antiplatelet agents prior to the procedure.
- The Pipeline™ Vantage Embolization device with Shield Technology™ should not be used alone as sole therapy for acutely ruptured aneurysms.

### PREPARATION FOR USE

- Choose a Pipeline™ Vantage device with a labeled diameter that is approximately equivalent to the target vessel diameter.
  - Select an appropriately sized Pipeline™ Vantage device such that its fully expanded diameter is equivalent to that of the largest target vessel diameter. An incorrectly sized Pipeline™ Vantage device may result in inadequate device placement, incomplete opening or migration.
- Choose a Pipeline™ Vantage device with labeled length that is at least 6 mm longer than the aneurysm neck.
  - Take device foreshortening into account when deploying the Pipeline™ Vantage device. The Pipeline™ Vantage device foreshortens 47 - 58% during deployment.
- Remove packaging hoop from the pouch and pull the distal end of the introducer sheath from the blue clip on the packaging hoop.
- Carefully remove device from the packaging hoop until the delivery wire is exposed.

### WARNING

- Pre-deploying the distal end of the device prior to introduction into the micro catheter may cause damage to the distal end of the braid
- Partially insert introducer sheath into the rotating hemostatic valve (RHV) at the micro catheter hub and close the RHV. Use a minimum flush pressure of 250 mmHg and confirm back flush of the saline at the proximal end of the introducer sheath prior to advancing the Pipeline™ Vantage device into the micro catheter.
  - Advance introducer sheath into the RHV; visually confirm the tip of the sheath is seated deeply in the hub of the micro catheter.

## DIRECTIONS FOR USE

- Using standard interventional radiographic technique, place the micro catheter tip at least 20 mm past the distal edge of the aneurysm. Gently retract the micro catheter to reduce slack in the micro catheter prior to inserting the Pipeline™ Vantage device.

**NOTE:** It is recommended to use a heparinized saline drip to continuously flush micro catheter during Pipeline™ Vantage device use.

- Secure introducer sheath to the hub by locking down the RHV tightly.

**CAUTION:** Avoid deploying the device prior to introduction into the micro catheter.

- Advance the proximal end of the delivery wire until it aligns with the proximal end of the introducer sheath.
- Remove the introducer sheath.

**NOTE:** The delivery wire has a fluorosafe marker no further than 125 cm from the distal end.

**CAUTION:** The fluorosafe marker is only compatible with microcatheters with a minimum length of 135 cm.

- Advance the Pipeline™ Vantage device into the micro catheter by pushing the delivery wire until the tip of the delivery wire aligns with the tip of the micro catheter.

**CAUTION:** If high forces or excessive friction are encountered during delivery, discontinue delivery of the device and identify the cause of the resistance, remove device and micro catheter simultaneously. Advancement of the Pipeline™ Vantage device against resistance may result in device damage or patient injury.

**CAUTION:** The presence of other indwelling endovascular stents may interfere with proper deployment and function of the Pipeline™ Vantage device.

- Once the tip of delivery wire and micro catheter are aligned, verify that the Pipeline™ Vantage implant is in the desired location. The distal end of Pipeline™ Vantage implant should be placed at least 3 mm past the distal edge of the aneurysm neck.
- Begin to deliver the Pipeline™ Vantage implant using a combination of unsheathing the Pipeline™ Vantage implant and pushing the delivery wire simultaneously.

## WARNING

- Pushing delivery wire without retracting the micro catheter at the same time will cause the open-end of the braid to move distally in the vessel. This may cause damage to the braid or vessel.
- Use in anatomy with severe tortuosity, stenosis or parent vessel narrowing may result in difficulty or inability to deploy the Pipeline™ Vantage device and can lead to damage to the Pipeline™ Vantage device and microcatheter. To mitigate potential problems as a result of increased delivery forces, reduce the load in the system by:
  - Unloading the microcatheter to the inner curves of vessel by pulling back on the system (i.e., the microcatheter and delivery wire together).
  - Continue unloading the system until advancement of the device (inside of microcatheter) is observed, while minimizing the distal tip movement to prevent loss of position.
  - Begin to re-advance the delivery wire while maintaining reduced load in the microcatheter. This process should be repeated until the device passes through tortuous area and the delivery force is decreased.

### 8. Resheathing Instructions:


During deployment of the Pipeline™ Vantage device resheathing can be performed by either:

- Advancing the micro catheter while pinning the delivery wire
  - Advancing the micro catheter while applying tension on the delivery wire
  - Advancing the micro catheter while gently pulling the delivery wire proximally
- During deployment, the point of no return/Resheathing limit is reached when the Resheathing marker aligns with the Distal marker of the micro catheter (see Figure 3). The Resheathing limit is the maximum length of the implant that can be deployed while maintaining the ability to fully resheath the device.
- The Pipeline™ Vantage device implant is fully resheathed when the distal marker is retracted completely inside the micro catheter. The system is designed to allow for a 2 full cycles of resheathing of the Pipeline™ Vantage device.
- After the distal end of the implant has successfully expanded, deploy the remainder of the implant by pushing and/or unsheathing the delivery wire. Manipulation of the micro catheter by locking down the delivery wire and moving both as a system may facilitate expansion of the implant.
- CAUTION:** Under fluoroscopy, carefully monitor the tip coil position during deployment of the Pipeline™ Vantage device.

## WARNING

- Resheathing the Pipeline™ Vantage device more than 2 full cycles may cause damage to the distal or proximal ends of the braid.
- Resheathing the Pipeline™ Vantage device past the distal marker of the delivery system may cause damage to the distal end of the braid.

**Figure 3. Pipeline™ Vantage Embolization Device with Shield Technology™ (Resheathing schematic as seen under fluoroscopy, image not to scale).**

- |                                 |                      |
|---------------------------------|----------------------|
| 1. Proximal End of device       | 5. Proximal Bumper   |
| 2. Micro Catheter               | 6. Delivery Wire     |
| 3. Micro Catheter Distal Marker | 7. Resheathing Limit |
| 4. Resheathing Marker           | 8. Device Detached   |

- After the entire implant is deployed, advance the micro catheter through the implant making sure not to dislodge the braid. When the micro catheter tip is distal to the implant, retract the delivery wire into the micro catheter tip.

**CAUTION:** Avoid advancing or retracting the Resheathing Marker within the implant without coverage of the microcatheter.

**CAUTION:** If the catheter cannot be advanced through the Pipeline™ Vantage implant, carefully withdraw the delivery wire through the implant.

**CAUTION:** If the delivery wire cannot be retracted into the micro catheter, carefully remove the delivery wire and micro catheter simultaneously as a system.

- Carefully inspect the deployed implant under fluoroscopy to confirm it is completely apposed to the vessel wall and not kinked/twisted. If the device is not fully apposed or is kinked/twisted, consider using a balloon catheter, micro catheter, or guidewire to fully open it.

**CAUTION:** It is not recommended to use the Pipeline™ Vantage delivery wire to influence apposition of the implant.

## POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to the following:

- Adverse reaction to antiplatelet/anticoagulation agents, anesthesia reactions such as pain, nausea, aspiration, or to contrast media such as burn sensation and organ damage or failure or due to radiation exposure such as alopecia, burns, skin reddening, ulcers, skin discoloration, cataracts, delayed neoplasia
- Access site complications such as edema, abscess, bleeding including retroperitoneal hemorrhage, tissues damage, hematoma, hemorrhage, and nerve damage
- Vascular complications such as vasospasm, stenosis, dissection, perforation, rupture, AV fistula formation, pseudo aneurysm, occlusion, thromboembolic complications including ischemia, occlusion, embolism (to unintended territory)
- Device malfunctions such as kink, stretching, friction, fracture, breakage, foreign body, misplacement, migration, inadequate deployment, premature deployment, non-detachment, reaction to device materials (such as hypersensitivity, hemolysis, fever, mutagenic effects, inflammation, granuloma, toxicity)
- Systemic complications such as infection, discomfort, pain, fever, shock, allergic reactions, organ damage, organ failure, hypertension, hypotension, arrhythmia, angina, myocardial infarction.
- Neurological deficits or dysfunctions including stroke, infarction, visual deficits, loss of vision, seizures, motor function, transient ischemic attack, headache, cranial neuropathy, confusion, emotional changes, coma
- Bleeding/ hemorrhagic complications.

- Visual complications include but are not limited to Amaurosis fugax/transient blindness, Blindness, Diplopia, Reduced visual acuity/field, Retinal artery occlusion, Retinal ischemia, Retinal infarction, Vision impairment including scintillations, blurred vision, eye floaters
- Decreased therapeutic response including need for target aneurysm retreatment
- Intra-cranial hemorrhage (including from aneurysm rupture), mass effect, brain edema, hydrocephalus
- Death

\* Consult instructions for use for other therapy devices and medications for additional potential complication information.

## WARNING

- Person with known allergy to cobalt/chromium alloy (including major elements cobalt, chromium, nickel, molybdenum) or platinum may suffer an allergic reaction to the Pipeline™ Vantage device implant.
- Person with known allergy to platinum alloy (including major elements platinum, tungsten, iridium), tin, silver, stainless steel or silicone elastomer may suffer an allergic reaction to the Pipeline™ Vantage device delivery system.
- Do not reprocess or resterilize. Reprocessing and resterilization increase the risk of patient infection and compromised device performance.
- Placement of multiple Pipeline™ Vantage devices may increase the risk of ischemic complications.
- Do not attempt to reposition the device after full deployment.

## PRECAUTIONS

- Physicians should undergo appropriate training prior to using the Pipeline™ Vantage device in patients.
- The Pipeline™ Vantage device is intended for single use only. Carefully inspect the sterile package and device components prior to use. Do not use if sterile package or device components are damaged.
- Use the Pipeline™ Vantage device system prior to the "Use-By-date" printed on the package.
- Do not use the Pipeline™ Vantage device in patients in whom angiography demonstrates inappropriate anatomy for endovascular treatment, such as severe pre- or post-aneurysmal narrowing or severe intracranial vessel tortuosity.
- The appropriate anti-platelet and anti-coagulation therapy should be administered in accordance with standard medical practice.
- A thrombosing aneurysm may aggravate pre-existing, or cause new, symptoms of mass effect and may require medical therapy.
- Use of implants with labeled diameter larger than the parent vessel diameter may result in decreased effectiveness and additional safety risk due to incomplete foreshortening resulting in an implant longer than anticipated.
- Take all necessary precautions with patients in whom a pre-existing stent is in place in the parent artery at the target aneurysm location.
- Take all necessary precautions to limit X-ray radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

## HOW SUPPLIED

This device is supplied STERILE using ethylene oxide. This device is Non-pyrogenic.

## STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

## DIAGNOSTIC MAGNETIC RESONANCE (MR) IMAGING

Non-clinical testing has demonstrated that the Pipeline™ Vantage device is MR Conditional for single and overlapping stents up to 70 mm in length. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla or 3-Tesla, only.
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less.
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2-W/kg (Normal Operating Mode of Operation for the MR system).
- Maximum head SAR of 3.2 W/kg.

After 15-minutes of continuous scanning the Pipeline™ Vantage device is expected to produce a maximum temperature rise up to 4.15°C.

### Artifact Information

In non-clinical testing, the image artifact caused by the Pipeline™ Vantage device extends approximately 20.2 mm from this implant when imaged using a T1-weighted spin echo pulse sequence and a 3-Tesla MR system.

Multilayer implant configuration of the Pipeline™ Vantage device does not affect its MRI compatibility, including temperature rise, torque, displacement, and artifact.

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Micro Therapeutics, Inc.
Medtronic B.V.