
# Medtronic

## Admiral™ Xtreme

PTA Balloon Catheter

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Instructions for Use • Mode d'emploi

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## PTA Balloon Catheter

### 1. Description

The Admiral Xtreme PTA balloon catheter is an over-the-wire (OTW) peripheral balloon catheter, specifically designed for percutaneous transluminal angioplasty (PTA) in stenosed vessel segments. The catheter has a dual-lumen shaft that is branched at the proximal end. One lumen forms the entrance to the central lumen for the guidewire, and the other lumen is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution. The catheter and balloon are designed to reach targeted inflation diameters, depending on the balloon size and defined pressure.

The Admiral Xtreme PTA balloon catheter is available in a variety of balloon sizes. The Compliance Chart included in the catheter packaging lists the balloon diameters at given pressures.

The balloon catheter is provided with hydrophilic coating.

### 2. Indications

The Admiral Xtreme PTA balloon dilatation catheter is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

### 3. Contraindications

The Admiral Xtreme PTA balloon catheter is contraindicated for use in coronary arteries and in the neurovasculature. It is also contraindicated for cases where the guidewire cannot cross the target lesion.

### 4. Warnings

- This device is designed and intended for single use only. DO NOT RESTERILIZE AND/OR REUSE. Reuse or resterilization may create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Reuse or resterilization may compromise the structural integrity of the device or lead to device failure, which, in turn, may result in patient injury, illness, or death. Medtronic is not responsible for any direct, incidental, or consequential damages resulting from resterilization or reuse.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not manipulate the balloon in its inflated state. The position of the balloon catheter may only be changed with the guidewire in place.
- If resistance occurs during manipulation, the cause must first be ascertained by fluoroscopy, road mapping, or DSA before moving the balloon backward or forward.
- Do not move the guidewire during inflation of the balloon.
- Completely deflate the balloon before removing it from the vascular system.
- Do not exceed the rated burst pressure. The rated burst pressure is based on the results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure. Use of a pressure-monitoring device is recommended to prevent over pressurization.
- Use only a mixture of contrast medium and saline solution to fill the balloon (1:1). Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the Use-by date specified on the package.

### 5. Precautions

- Only interventionalists who have sufficient experience in percutaneous transluminal angioplasty should use this PTA balloon catheter. A thorough understanding of the technical principles, clinical applications, and risks associated with percutaneous transluminal angioplasty is necessary before using this product.
- Administer appropriate drug therapy (anticoagulant, vasodilator, etc.) to the patient according to standard protocols for PTA before inserting the dilatation catheter.
- Exercise care during handling in order to avoid possible damage to the catheter. Do not use a catheter that has been damaged.
- To minimize the introduction of air, aspirate and flush the system and keep a tight catheter connection throughout the procedure.
- Use the device with caution for procedures involving calcified lesions or synthetic vascular grafts due to the abrasive nature of these lesions.

- Identify allergic reactions to contrast media before PTA.
- Observe the general technical requirements for catheter insertion at all times, including flushing the components with sterile, isotonic saline solution prior to use and the usual prophylactic, systemic heparinization.
- Catheter applications vary. Select the technique on the basis of the patient's condition and the experience of the interventionalist.
- The maximum diameter of the guidewire must not exceed over 0.035 in (0.89 mm). Certain coatings on guidewires may cause increased friction between the guidewire and the catheter.

**Caution:** Larger models of Admiral Xtreme PTA balloon catheters may exhibit slower deflation times particularly on long catheter shafts.

## 6. Potential complications / adverse effects

The complications that may result from a balloon dilatation procedure include:

Puncture related:

- Local hematoma
- Local hemorrhage
- Local or distal thromboembolic episodes
- Thrombosis
- Arteriovenous fistula
- Pseudoaneurysm
- Local infections

Dilatation related

- Dissection in the dilated artery wall
- Perforation of the artery wall
- Prolonged spasms
- Acute reocclusion necessitating surgical intervention
- Restenosis of the dilated artery
- Total occlusion of the peripheral artery

Angiography related

- Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection
- Systemic embolization
- Endocarditis
- Short-term hemodynamic deterioration
- Death
- Drug reactions
- Allergic reaction to contrast medium
- Pyrogenic reaction

## 7. Preparation technique

Before using the device, carefully examine the unit to verify that the catheter and sterile package have not been damaged in shipment. After moistening the balloon's protective sheath with physiological saline solution, carefully remove the sheath so as not to damage the balloon. In case of resistance, twist the sheath with one hand while holding the shaft with the other.

**Caution:** Do not use with Lipiodol™ or Ethiodol™ contrast media, or other such contrast media that incorporate the components of these agents.

**To displace air:**

1. Flush the guidewire lumen by connecting a syringe filled with saline solution. Flush until liquid flows out of the guidewire lumen at the distal tip of the balloon catheter.
2. Connect a syringe containing contrast medium to the Luer connector located at the proximal hub of the balloon catheter. Apply negative pressure for about 15 seconds until no bubbles appear in the contrast medium solution. Return to atmospheric pressure allowing contrast medium to flow into the catheter lumen. Remove the syringe, leaving a meniscus of contrast in the hub of the balloon lumen.
3. Remove air bubbles from the inflation device following the manufacturer's instructions.

4. Using a stopcock, attach the inflation device to the PTA balloon catheter. Do not allow air to enter the system.

**Caution:** Do not manipulate the balloon catheter during removal from packaging, nor while flushing the guidewire lumen.

The Admiral Xtreme PTA balloon catheter is provided with a hydrophilic coating. Before inserting the catheter, activate the coating by immersing the catheter in normal saline for approximately 30 to 60 seconds, or by wiping the catheter shaft with a saturated gauze sponge.

**Caution:** Do not wipe down the catheter surface with dry gauze.

## 8. Insertion technique

Access to the vessel must be sufficiently patent and sufficiently recanalized to proceed with further intervention.

If the guidewire has not already been positioned across the lesion, insert a 0.035 in (0.89 mm) guidewire across the lesion. Use fluoroscopy and follow standard PTA techniques.

Insert the proximal end of the guidewire into the distal tip of the catheter.

Carefully insert the PTA balloon catheter through the introducer sheath and, under fluoroscopic guidance, carefully advance it through the peripheral vasculature following the guidewire toward the lesion.

**Caution:** If resistance is met, do not advance the guidewire or the balloon dilatation catheter without first determining the cause of resistance and taking remedial action.

Position the balloon within the lesion using the 2 radiopaque markers located under the balloon as reference points for precise placement across the target lesion.

To avoid kinking, advance the dilatation catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.

## 9. Balloon inflation

**Caution:** Use only a mixture of contrast medium and saline solution to fill the balloon (1:1). Do not use any gaseous substances as a balloon inflation medium.

Purge all air bubbles from the catheter and the inflation device before balloon inflation.

Open the stopcock on the inflation device. Inflate the dilation balloon gradually to dilate the lesion to the targeted diameter in accordance with the Compliance Chart. Apply a constant pressure for about 30 seconds.

The standard inflation medium is a 1:1 mixture of contrast medium and normal saline.

**Note:** A 1:3 ratio of contrast to saline has yielded faster balloon inflation and deflation times.

**Caution:** Do not exceed the rated burst pressure (RBP).

## 10. Balloon deflation and removal

Deflate the balloon catheter in accordance with standard PTA procedures. Apply negative pressures to the balloon for about 60 to 120 seconds. Ensure that the balloon is fully deflated before carefully removing the catheter from the vessel.

After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable hospital, administrative, and government regulations.

## 11. How supplied

The Admiral Xtreme PTA balloon catheter is supplied sterile and intended for single use only. The Admiral Xtreme PTA balloon catheter is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. Use product before the labeled Use-by date.

The Admiral Xtreme PTA balloon catheter is not made with natural rubber latex. However, during the manufacturing, assembly, or packaging process, it may have incidental contact with materials, tools, or equipment containing or coming in contact with latex.

**Caution:** Do not use if the inner package is open or damaged.

## 12. Storage

Store at controlled room temperature in a dry place. Do not expose to organic solvents (such as alcohol), ionizing radiation, or ultraviolet light. Rotate the inventory so that catheters are used before the Use-by date on the package label.

## 13. Disclaimer of warranty

Although this product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic therefore disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Medtronic to any representation or warranty with

respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular term held to be invalid.

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© 2017 Medtronic
Trademarks may be registered and are the property of their respective owners.
Medtronic, Inc.