ALINITY H-SERIES

Diluent

INTENDED USE

The Alinity h-series Diluent is to be used in conjunction with Alinity hq analyzer as a primary diluent for the red blood cells (RBC), platelets (PLT) and reticulocyte (RETIC) measurement.

METHOD PRINCIPLE

- Prepares each cell for counting and sizing.
- Provides enough wetting action to prevent the accumulation of bubbles in the flow system.
- Primes pumps and flow path and rinses the tubing and other components at the end of each cycle.

REAGENT TYPE PRESENTATION

Liquid 1 cubitainer x 9.6 L

REAGENT STORAGE AND STABILITY

- Store reagent between 15 °C – 30 °C
- The shelf life of the reagent is printed on the bottle label
- On-board reagent stability is 90 days.

When storing reagents, calibrators, and controls, comply with the following requirements:

- Store reagents, calibrators, and controls as instructed on their labels or in the product documentation.
- Protect reagents from extreme heat and freezing during storage. Temperatures below 0 °C can cause layering which changes the tonicity and conductivity of the reagent. Do not use reagents that have been frozen.

If reagents, calibrators, or controls are received in a condition that is contrary to the label recommendation or the product documentation, or are damaged, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REAGENT PREPARATION

The Alinity h-series Diluent is ready to use; no special preparation is necessary.

WASTE HANDLING AND DISPOSAL

Please refer to the MSDS associated with the reagent.² All waste produced by the Alinity h-series must be disposed of in accordance with local, state, and national regulations that govern the treatment and disposal of medical waste. All waste containers must be labeled as biohazardous waste.

CORELABORATORY. ABBOTT/HEMATOLOGY

© 2023 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. ADD-125991 EMEA-EN.

REAGENT REPLACEMENT

The Alinity h-series Diluent is a liquid ready to use reagent to be used with Alinity hq hematology analyzer. Replace a Diluent cubitainer in the lower reagent drawer by using the LUI. Pull open the upper reagent drawer.

- On the LUI menu bar, tap Reagents
- On the Reagent Status screen, tap Load
- Pull open the lower reagent drawer
- Place the new cubitainer on the floor next to the empty cubitainer in the lower reagent drawer
- Unscrew the cap from the new cubitainer
- Unscrew the cap and tubing assembly from the empty cubitainer and insert the assembly into the new cubitainer.
- Remove the empty cubitainer.
- Remove the reagent bar code label from the new cubitainer.
- Scan the reagent bar code with the LUI bar code reader to enter the reagent information.
- Tap Save
- Place the bar code label on top of the new cubitainer
- When a message is displayed that provides instructions to replace the cubitainer, for Alinity hq, tap Diluent A or Diluent B waste disposal procedures for the laboratory.

NOTES

- Reagent status indicators are not present in the lower drawer.
- Ensure that the end of the tube does not touch anything during the transfer process
- The scanner illuminates as a blue box when the bar code being scanned is correctly positioned for reading.
- Reagent packaging may be labeled as recyclable. Noncontaminated packaging may be used for recycling. Refer to applicable local regulations and laboratory waste disposal procedures. For disposal of contaminated packaging, refer to applicable local regulations and laboratory waste disposal procedure.

ALARMS

- Please refer to Alinity h-series Operations Manual
- Section 3: Principles of operation - System messages and flags and Section 5: Reagent Inventory notifications

REFERENCES

1. Alinity h-series Operations Manual
2. Technical Library: www.corelaboratory.abbott

ALINITY H-SERIES

WBC Reagent

INTENDED USE

The Alinity h-series WBC Reagent is to be used in conjunction with Alinity hq analyzer for counting and differentiation of white blood cells (WBC) in human peripheral blood.

METHOD PRINCIPLE

The Alinity h-series WBC Reagent is a ready to use reagent and uses Abbott's proprietary fluorescent nuclear dye membrane-permeable. WBC reagent lyses red blood cells and maintains the WBC integrity and morphology during the measurement. The dye facilitates the separation of nucleated cells (WBC and NRBC-nucleated red blood cell) from nonnucleated events, such as lyse-resistant red blood cells (RBC) and cellular debris. Lyse-resistant RBC do not interfere with WBC counts. Fluorescent emission from cell nuclei is measured with a detector in the fluorescence channel.

Alinity hq provides a six-part WBC differential. WBC subpopulations, including immature granulocytes (IG), are reported both in percent and as absolute number. NRBC are counted using fluorescence and light scatter characteristics and are expressed as both concentration and Nucleated red blood cells per 100 white blood cells (NR/W).

REAGENT TYPE PRESENTATION

Liquid, concentrated 4 bottles x 975 mL

COMPOSITION

- Potassium Bicarbonate..... <0.3%
- Ammonium Chloride..... <0.3%
- Formaldehyde..... <0.09%
- Sodium Acetate..... <0.2%
- Nuclear stain..... <0.0004%
- Saponin..... <0.02%
- Preservatives..... <0.07%

REAGENT STORAGE AND STABILITY

- Store reagent between 15°C – 30°C
- The shelf life of the reagent is printed on the bottle label
- On-board reagent stability is 90 days

When storing reagents, calibrators, and controls, comply with the following requirements:

- Store reagents, calibrators, and controls as instructed on their labels or in the product documentation.

- Protect reagents from extreme heat and freezing during storage. Temperatures below 0°C can cause layering which changes the tonicity and conductivity of the reagent. Do not use reagents that have been frozen.

If reagents, calibrators, or controls are received in a condition that is contrary to the label recommendation or the product documentation, or are damaged, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REAGENT PREPARATION

The Alinity h-series WBC Reagent is ready to use; no special preparation is necessary.

WASTE HANDLING AND DISPOSAL

Please refer to the MSDS associated with the reagent. (5)

All waste produced by the Alinity h-series must be disposed of in accordance with local, state, and national regulations that govern the treatment and disposal of medical waste. All waste containers must be labeled as biohazardous waste.

SPECIMEN COLLECTION AND PREPARATION

SPECIMEN TYPE

The Alinity h-series currently accepts venous and capillary whole blood specimens.

RECOMMENDED ANTICOAGULANTS

All performance claimed in the Alinity h-series Operations Manual are based on data from specimens that are collected in dipotassium EDTA (K2EDTA) or tripotassium EDTA (K3EDTA). Specimens that are collected in other anticoagulants are not recommended.

REQUIREMENTS FOR COLLECTION

When collecting specimens, comply with the following requirements:

- Collect all whole blood specimens according to laboratory procedures. Follow the recommendations of the collection tube manufacturer.
- Comply with all precautions for blood collection by venipuncture to prevent clotting or specimen hemolysis.

SPECIMEN COLLECTION

All blood samples must be collected using the proper technique. Considered as potentially infectious all samples, reagents, controls, etc. that contain extracts human samples and follow the instruction of biological safety. Sample collection tube must be filled with the exact amount of blood indicated on the tube itself to avoid variations in results.¹⁻³

SPECIMEN REQUIREMENTS FOR PREPARATION AND STORAGE

When preparing and storing whole blood specimens, comply with the following requirements:

- Use fresh whole blood specimens to achieve the most reliable result data. The International Council for Standardization in Hematology (ICSH) defines a fresh blood specimen as one that is processed within 4 hours after collection.
- For specimens that are processed within 8 hours after collection, room-temperature storage is recommended.
- For specimens that are processed more than 8 hours after collection, refrigerated (2°C to 8°C) storage is recommended.
- Before refrigerator-stored specimens are mixed and processed, bring the specimens to room temperature.

SPECIMEN STABILITY

- 24 hours when stored at 15°C - 30°C
- 48 hours when stored at 2°C - 10°C

REAGENT REPLACEMENT

The Alinity h-series WBC Reagent is a liquid ready to use reagent to be used with Alinity hq hematology analyzer. Replace a reagent bottle in the upper reagent drawer by using the LUI:

- Pull open the upper reagent drawer.
- View the reagent status indicators on the reagent bottle holder release buttons.
- Press the bottle holder release button to release the empty reagent bottle.
- Gently remove the empty reagent bottle from the bottle holder.
- On the menu bar, tap Reagents.
- On the Reagent Status screen, tap Load.
- On the New Reagent screen, scan the reagent bar code to enter the reagent information.
- Remove the lid from the cap on the new reagent bottle to expose the septum.
- Invert the new reagent bottle so that the cap faces downward.
- Firmly insert the new reagent bottle into the bottle holder until the bottle locks into position.
- When a confirmation message is displayed, tap OK.
- Close the upper reagent drawer.
- Discard the empty reagent bottle and the lid according to the waste disposal procedures for the laboratory.

NOTES

- Do not pierce the reagent bottle septum more than three times.
- A red indicator shows which reagent bottle requires replacement.

- Forceful removal of the empty reagent bottle may damage the septum piercers.
- The scanner illuminates as a blue box when the bar code being scanned is correctly positioned for reading.
- Reagent packaging may be labeled as recyclable. Noncontaminated packaging may be used for recycling.

PERFORMANCES CHARACTERISTICS

Please refer to Alinity h-series Operations Manual Section 4: Performance characteristics and specifications.

RESULTS CALCULATION

Please refer to Alinity h-series Operations Manual Section 3: Principles of operation - Data analysis and presentation (Alinity hq).

ALARMS

- Please refer to Alinity h-series Operations Manual
- Section 3: Principles of operation - System messages and flags and Section 5: Reagent Inventory notifications.

LIMITATIONS ON INTERPRETATION

Please refer to Alinity h-series Operations Manual Section 7 Operational precautions and limitations - Limitations of result interpretation.

QUALITY CONTROL PROCEDURES

Quality control on the Alinity h-series involves the monitoring of control results and patient data.

The quality control programs on the system help to assess precision and accuracy, identify shifts and trends, and determine the nature and cause of errors on the Alinity hq.

CONTROL MATERIAL NEEDED, BUT NOT INCLUDED

- 04U72-06 Alinity h-series Control 29P 3 levels, 6x3mL
- 04U72-16 Alinity h-series Control 29P 3 levels, 12x3mL

External quality control programs use resources outside the laboratory to assess system performance. For more information about external quality control programs, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REFERENCES

1. CLSI – OCT 2007: GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Sixth Edition
2. CLSI – Sept 2008: GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Sixth Edition
3. CLSI – May 2014: M29-A4 Protection of laboratory Workers from Occupationally Acquired Infections; Guideline – Forth Edition
4. Alinity h-series Operations Manual
5. Technical Library: www.corelaboratory.abbott

ALINITY H-SERIES

HGB Reagent

INTENDED USE

The Alinity h-series HGB Reagent is to be used in conjunction with Alinity hq analyzer to measure the hemoglobin in human peripheral blood.

METHOD PRINCIPLE

The Alinity h-series HGB Reagent is a ready to use reagent.

Absorption photometry is used to measure hemoglobin (HGB). This method is based on the linear relationship between the amount of light absorbed in a well-mixed, stationary sample and the concentration of the analyte in a sample (as defined by Beer's Law).

In the HGB analysis, red blood cells (RBC) are lysed using a cyanide-free HGB reagent. The HGB reagent destroys red blood cell membranes and white blood cells and converts the hemoglobin to a stable chromogenic complex that is measurable at 540 nm.

Alinity's HGB method is linear up to 250 g/L and is not affected by white blood cells (WBC) counts up to $250 \times 10^9/L$ (Kim, Y 1997).

REAGENT TYPE PRESENTATION

Liquid, concentrated 4 bottles x 975 mL

COMPOSITION

- Imidazole..... <2.8%
- Sodium Hydroxide..... <0.5%
- Surfactant..... <1.0%

REAGENT STORAGE AND STABILITY

- Store reagent between 15°C – 30°C
- The shelf life of the reagent is printed on the bottle label
- On-board reagent stability is 90 days.

When storing reagents, calibrators, and controls, comply with the following requirements:

Store reagents, calibrators, and controls as instructed on their labels or in the product documentation.

- Protect reagents from extreme heat and freezing during storage. Temperatures below 0°C can cause layering which changes the tonicity and conductivity of the reagent. Do not use reagents that have been frozen.

If reagents, calibrators, or controls are received in a condition that is contrary to the label recommendation or the product documentation, or are damaged, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REAGENT PREPARATION

The Alinity h-series HGB Reagent is ready to use; no special preparation is necessary.

WASTE HANDLING AND DISPOSAL

Please refer to the MSDS associated with the reagent.⁵ All waste produced by the Alinity h-series must be disposed of in accordance with local, state, and national regulations that govern the treatment and disposal of medical waste. All waste containers must be labeled as biohazardous waste.

SPECIMEN COLLECTION AND PREPARATION

SPECIMEN TYPE

The Alinity h-series currently accepts venous and capillary whole blood specimens.

RECOMMENDED ANTICOAGULANTS

All performance claimed in the Alinity h-series Operations Manual are based on data from specimens that are collected in dipotassium EDTA (K2EDTA) or tripotassium EDTA (K3EDTA). Specimens that are collected in other anticoagulants are not recommended.

REQUIREMENTS FOR COLLECTION

When collecting specimens, comply with the following requirements:

- Collect all whole blood specimens according to laboratory procedures. Follow the recommendations of the collection tube manufacturer.
- Comply with all precautions for blood collection by venipuncture to prevent clotting or specimen hemolysis.

SPECIMEN COLLECTION

All blood samples must be collected using the proper technique.

Considered as potentially infectious all samples, reagents, controls, etc. that contain extracts human samples and follow the instruction of biological safety. Sample collection tube must be filled with the exact amount of blood indicated on the tube itself to avoid variations in results.¹⁻³

SPECIMEN REQUIREMENTS FOR PREPARATION AND STORAGE

When preparing and storing whole blood specimens, comply with the following requirements:

- Use fresh whole blood specimens to achieve the most reliable result data. The International Council for Standardization in Hematology (ICSH) defines a fresh blood specimen as one that is processed within 4 hours after collection.
- For specimens that are processed within 8 hours after collection, room-temperature storage is recommended.
- For specimens that are processed more than 8 hours after collection, refrigerated (2°C to 8°C) storage is recommended.
- Before refrigerator-stored specimens are mixed and processed, bring the specimens to room temperature.

SPECIMEN STABILITY

- 24 hours when stored at 15°C - 30°C
- 48 hours when stored at 2°C - 10°C

REAGENT REPLACEMENT

The Alinity h-series HGB Reagent is a liquid ready to use reagent to be used with Alinity hq hematology analyzer. Replace a reagent bottle in the upper reagent drawer by using the LUI:

- Pull open the upper reagent drawer.
- View the reagent status indicators on the reagent bottle holder release buttons.
- Press the bottle holder release button to release the empty reagent bottle.
- Gently remove the empty reagent bottle from the bottle holder.
- On the menu bar, tap Reagents.
- On the Reagent Status screen, tap Load.
- On the New Reagent screen, scan the reagent bar code to enter the reagent information.
- Remove the lid from the cap on the new reagent bottle to expose the septum.
- Invert the new reagent bottle so that the cap faces downward.
- Firmly insert the new reagent bottle into the bottle holder until the bottle locks into position.
- When a confirmation message is displayed, tap OK.
- Close the upper reagent drawer.
- Discard the empty reagent bottle and the lid according to the waste disposal procedures for the laboratory.

NOTES

- Do not pierce the reagent bottle septum more than three times.
- A red indicator shows which reagent bottle requires replacement.
- Forceful removal of the empty reagent bottle may damage the septum piercers.
- The scanner illuminates as a blue box when the bar code being scanned is correctly positioned for reading.
- Reagent packaging may be labeled as recyclable. Noncontaminated packaging may be used for recycling.

CORELABORATORY. ABBOTT/HEMATOLOGY

PERFORMANCES CHARACTERISTICS

Please refer to Alinity h-series Operations Manual Section 4: Performance characteristics and specifications.

RESULTS CALCULATION

Please refer to Alinity h-series Operations Manual Section 3: Principles of operation - Data analysis and presentation (Alinity hq).

ALARMS

- Please refer to Alinity h-series Operations Manual
- Section 3: Principles of operation - System messages and flags and Section 5: Reagent Inventory notifications.

LIMITATIONS ON INTERPRETATION

Please refer to Alinity h-series Operations Manual Section 7: Operational precautions and limitations - Limitations of result interpretation.

QUALITY CONTROL PROCEDURES

Quality control on the Alinity h-series involves the monitoring of control results and patient data.

The quality control programs on the system help to assess precision and accuracy, identify shifts and trends, and determine the nature and cause of errors on the Alinity hq.

CONTROL MATERIAL NEEDED, BUT NOT INCLUDED

- 04U72-06 Alinity h-series Control 29P 3 levels, 6x3mL
- 04U72-16 Alinity h-series Control 29P 3 levels, 12x3mL

External quality control programs use resources outside the laboratory to assess system performance. For more information about external quality control programs, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REFERENCES

1. CLSI – OCT 2007: GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Sixth Edition
2. CLSI – Sept 2008: GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Sixth Edition
3. CLSI – May 2014: M29-A4 Protection of Laboratory Workers from Occupationally Acquired Infections; Guideline – Forth Edition
4. Alinity h-series Operations Manual
5. Technical Library: www.corelaboratory.abbott

ALINITY H-SERIES

Autoclean Solution

INTENDED USE

The Alinity h-series AutoClean Solution is to be used in conjunction with Alinity hq analyzer as cleaning solution.

METHOD PRINCIPLE

The major function is to clean the internal pathway of the Alinity hq module. The Alinity h-series AutoClean Solution is combined with CLRW compliant water from the external water purification system to provide the final 0.525% sodium hypochlorite concentration.

REAGENT TYPE

Liquid, concentrated

PRESENTATION

4 bottles x 1000 mL

COMPOSITION

Sodium Hypochlorite $\leq 6.3\%$ Sodium Hydroxide $<0.06 - 0.60\%$

REAGENT STORAGE AND STABILITY

- Store reagent between 15 – 30 °C
- The shelf life of the reagent is printed on the bottle label
- On-board reagent stability is 90 days.

When storing reagents, calibrators, and controls, comply with the following requirements:

- Store reagents, calibrators, and controls as instructed on their labels or in the product documentation
- Protect reagents from extreme heat and freezing during storage. Temperatures below 0°C can cause layering which changes the tonicity and conductivity of the reagent. Do not use reagents that have been frozen. If reagents, calibrators, or controls are received in a condition that is contrary to the label recommendation or the product documentation, or are damaged, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REAGENT PREPARATION

The Alinity h-series AutoClean Solution is ready to use; no special preparation is necessary.

WASTE HANDLING AND DISPOSAL

Please refer to the MSDS associated with the reagent. (2) All waste produced by the Alinity h-series must be disposed of in accordance with local, state, and national regulations that govern the treatment and disposal of medical waste. All waste containers must be labeled as biohazardous waste.

CORELABORATORY. ABBOTT/HEMATOLOGY

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. ADD-124253 EMEA/EN.

REAGENT REPLACEMENT

The Alinity h-series AutoClean Solution is a liquid ready to use reagent to be used with Alinity hq hematology analyzer. Replace a reagent bottle in the upper reagent drawer by using the LUI:

- Pull open the upper reagent drawer.
- View the reagent status indicators on the reagent bottle holder release buttons.
- Press the bottle holder release button to release the empty reagent bottle.
- Gently remove the empty reagent bottle from the bottle holder.
- On the menu bar, tap Reagents.
- On the Reagent Status screen, tap Load.
- On the New Reagent screen, scan the reagent bar code to enter the reagent information.
- Remove the lid from the cap on the new reagent bottle to expose the septum.
- Invert the new reagent bottle so that the cap faces downward.
- Firmly insert the new reagent bottle into the bottle holder until the bottle locks into position.
- When a confirmation message is displayed, tap OK.
- Close the upper reagent drawer.
- Discard the empty reagent bottle and the lid according to the waste disposal procedures for the laboratory.

NOTES

- Do not pierce the reagent bottle septum more than three times.
- A red indicator shows which reagent bottle requires replacement.
- Forceful removal of the empty reagent bottle may damage the septum piercers.
- The scanner illuminates as a blue box when the bar code being scanned is correctly positioned for reading.
- Reagent packaging may be labeled as recyclable. Noncontaminated packaging may be used for recycling.

ALARMS

- Please refer to Alinity h-series Operations Manual
- Section 3: Principles of operation - System messages and flags and Section 5: Reagent Inventory notifications

REFERENCES

1. Alinity h-series Operations Manual
2. Technical Library: www.corelaboratory.abbott

ALINITY H-SERIES

Retic Reagent

INTENDED USE

The Alinity h-series Retic Reagent is to be used in conjunction with Alinity hq analyzer to determine the concentration and percentage of reticulocytes in human peripheral blood.

METHOD PRINCIPLE

The Alinity h-series Retic Reagent is a ready to use reagent and uses Abbott's proprietary fluorescent reticulocyte dye that is excited at 488 nm and emits fluorescence at 530nm when binding to DNA or RNA. Reticulocytes are non-nucleated, immature RBC recently released from the bone marrow. They contain cytoplasmic ribosomal RNA, indicating active hemoglobin production. Mature RBC, on the other hand, do not contain RNA and no longer make hemoglobin. This principle is used in the Alinity hq reticulocyte method. The reticulocyte reagent uses a fluorescent dye which binds to the RNA found in reticulocytes. When a reticulocyte count is ordered, an aliquot is taken from the RBC/PLT dilution. The algorithm first separates RBC from platelets and WBC based on size and fluorescence properties. It then differentiates mature RBC, with very low fluorescence, from higher fluorescing reticulocytes. A second dynamic threshold is used for separating IRF (with highest fluorescent intensity) from more mature reticulocytes. The same dye is used to enumerate reticulated platelets.

REAGENT TYPE

Liquid, concentrated

PRESENTATION

4 bottles x 370 mL

COMPOSITION

- Imidazole..... <0.4%
- Sodium Chloride..... <0.7%
- Nuclear stain..... <0.00005%
- Surfactant..... <0.006%
- Preservatives..... <0.04%

REAGENT STORAGE AND STABILITY

- Store reagent between 2°C – 8°C
- The shelf life of the reagent is printed on the bottle label
- On-board reagent stability is 30 days

When storing reagents, calibrators, and controls, comply with the following requirements:

- Store reagents, calibrators, and controls as instructed on their labels or in the product documentation.
- Protect reagents from extreme heat and freezing during storage. Temperatures below 0°C can cause layering which changes the tonicity and conductivity of the reagent. Do not use reagents that have been frozen.

If reagents, calibrators, or controls are received in a condition that is contrary to the label recommendation or the product documentation, or are damaged, contact the local representative. Country-specific contact information is available at abbottdiagnostics.com

REAGENT PREPARATION

The Alinity h-series Retic Reagent is ready to use; no special preparation is necessary.

WASTE HANDLING AND DISPOSAL

Please refer to the MSDS associated with the reagent.⁵ All waste produced by the Alinity h-series must be disposed of in accordance with local, state, and national regulations that govern the treatment and disposal of medical waste. All waste containers must be labeled as biohazardous waste.

SPECIMEN COLLECTION AND PREPARATION

SPECIMEN TYPE

The Alinity h-series currently accepts venous and capillary whole blood specimens.

RECOMMENDED ANTICOAGULANTS

All performance claimed in the Alinity h-series Operations Manual are based on data from specimens that are collected in dipotassium EDTA (K2EDTA) or tripotassium EDTA (K3EDTA). Specimens that are collected in other anticoagulants are not recommended.

REQUIREMENTS FOR COLLECTION

When collecting specimens, comply with the following requirements:

- Collect all whole blood specimens according to laboratory procedures. Follow the recommendations of the collection tube manufacturer.
- Comply with all precautions for blood collection by venipuncture to prevent clotting or specimen hemolysis.

SPECIMEN COLLECTION

All blood samples must be collected using the proper technique. Considered as potentially infectious all samples, reagents, controls, etc. that contain extracts human samples and follow the instruction of biological safety. Sample collection tube must be filled with the exact amount of blood indicated on the tube itself to avoid variations in results.¹⁻³

SPECIMEN REQUIREMENTS FOR PREPARATION AND STORAGE

When preparing and storing whole blood specimens, comply with the following requirements:

- Use fresh whole blood specimens to achieve the most reliable result data. The International Council for Standardization in Hematology (ICSH) defines a fresh blood specimen as one that is processed within 4 hours after collection.
- For specimens that are processed within 8 hours after collection, room-temperature storage is recommended.
- For specimens that are processed more than 8 hours after collection, refrigerated (2°C to 8°C) storage is recommended.
- Before refrigerator-stored specimens are mixed and processed, bring the specimens to room temperature.

SPECIMEN STABILITY

- 24 hours when stored at 15°C - 30°C
- 48 hours when stored at 2°C - 10°C

REAGENT REPLACEMENT

The Alinity h-series Retic Reagent is a liquid ready to use reagent to be used with Alinity hq hematology analyzer. Replace a reagent bottle in the upper reagent drawer by using the LUI:

- Pull open the upper reagent drawer.
- View the reagent status indicators on the reagent bottle holder release buttons.
- Press the bottle holder release button to release the empty reagent bottle.
- Gently remove the empty reagent bottle from the bottle holder.
- On the menu bar, tap Reagents.
- On the Reagent Status screen, tap Load.
- On the New Reagent screen, scan the reagent bar code to enter the reagent information.
- Remove the lid from the cap on the new reagent bottle to expose the septum.
- Invert the new reagent bottle so that the cap faces downward.
- Firmly insert the new reagent bottle into the bottle holder until the bottle locks into position.
- When a confirmation message is displayed, tap OK.
- Close the upper reagent drawer.
- Discard the empty reagent bottle and the lid according to the waste disposal procedures for the laboratory.

NOTES

- Do not pierce the reagent bottle septum more than three times.
- A red indicator shows which reagent bottle requires replacement.
- Forceful removal of the empty reagent bottle may damage the septum piercers.
- The scanner illuminates as a blue box when the bar code being scanned is correctly positioned for reading.
- Reagent packaging may be labeled as recyclable. Noncontaminated packaging may be used for recycling.

PERFORMANCES CHARACTERISTICS

Please refer to Alinity h-series Operations Manual Section 4: Performance characteristics and specifications.

RESULTS CALCULATION

Please refer to Alinity h-series Operations Manual Section 3: Principles of operation - Data analysis and presentation (Alinity hq).

ALARMS

- Please refer to Alinity h-series Operations Manual
- Section 3: Principles of operation - System messages and flags and Section 5: Reagent Inventory notifications.

LIMITATIONS ON INTERPRETATION

Please refer to Alinity h-series Operations Manual Section 7 Operational precautions and limitations - Limitations of result interpretation.

QUALITY CONTROL PROCEDURES

Quality control on the Alinity h-series involves the monitoring of control results and patient data.

The quality control programs on the system help to assess precision and accuracy, identify shifts and trends, and determine the nature and cause of errors on the Alinity hq.

CONTROL MATERIAL NEEDED, BUT NOT INCLUDED

- 04U72-06 Alinity h-series Control 29P 3 levels, 6x3mL
- 04U72-16 Alinity h-series Control 29P 3 levels, 12x3mL

External quality control programs use resources outside the laboratory to assess system performance. For more information about external quality control programs, contact the local representative. Country-specific contact information is available at abbottddiagnostics.com

REFERENCES

1. CLSI – OCT 2007: GP41-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Sixth Edition
2. CLSI – Sept 2008: GP42-A6 Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Sixth Edition
3. CLSI – May 2014: M29-A4 Protection of laboratory Workers from Occupationally Acquired Infections: Guideline – Fourth Edition
4. Alinity h-series Operations Manual
5. Technical Library: www.corelaboratory.abbott

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BIBLIOGRAFIE

1. Clinical and Laboratory Standards Institute (CLSI). *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI Document H26-A2. Wayne, PA: CLSI; 2010.
2. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
3. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 6th ed. Washington, DC: US Government Printing Office; June 2020.
4. World Health Organization. *Laboratory Biosafety Manual*. 4th ed. Geneva: World Health Organization; 2020.
5. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
6. International Committee for Standardization in Haematology Expert Panel on Cytometry. The assignment of values to fresh blood used for calibrating automated blood cell counters. *Clin Lab Haematol* 1988;10(2):203-212.
7. The Expert Panel in Cytometry on the ICSH, J.M. England, et. al. Recommended Methods for the Visual Determination of White Cell and Platelet Counts. World Health Organization, WHO/LAB/88.3, 1988.
8. International Council for Standardization in Haematology Expert Panel on Cytometry. Reference method for the enumeration of erythrocytes and leucocytes. *Clin Lab Haematol* 1994;16(2):131-138.
9. Clinical and Laboratory Standards Institute (CLSI). *Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition*. CLSI Document H15-A3. Wayne, PA: CLSI; 2000.
10. Clinical and Laboratory Standards Institute (CLSI). *Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard—Third Edition*. CLSI Document H07-A3. Wayne, PA: CLSI; 2000.
11. International Council for Standardization in Haematology Expert Panel on Cytometry; International Society of Laboratory Hematology Task Force on Platelet Counting. Platelet counting by the RBC/platelet ratio method. *Am J Clin Pathol* 2001;115(3):460-464.

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 Abbott GmbH
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

MANUFACTURED FOR
Abbott Laboratories

DISTRIBUTED IN THE USA BY
Abbott Laboratories
Abbott Park, IL 60064 USA

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■ BIBLIOGRAFIE

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
5. The Expert Panel in Cytometry on the ICSH, J.M. England, et.al.: Recommended Methods for the Visual Determination of White Cell and Platelet Counts. World Health Organization, WHO/LAB/88.3, 1988.
6. International Council for Standardization in Haematology Expert Panel on Cytometry. Reference method for the enumeration of erythrocytes and leucocytes. *Clin Lab Haematol* 1994;16(2):131-138.
7. Clinical and Laboratory Standards Institute (CLSI). *Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard—Third Edition*. CLSI Document H15-A3. Wayne, PA: CLSI; 2000.
8. Clinical and Laboratory Standards Institute (CLSI). *Procedure for Determining Packed Cell Volume by the Microhematocrit Method; Approved Standard—Third Edition*. CLSI Document H07-A3. Wayne, PA: CLSI; 2000.
9. International Committee for Standardization in Haematology Expert Panel on Cytometry. The assignment of values to fresh blood used for calibrating automated blood cell counters. *Clin Lab Haematol* 1988;10(2):203-212.
10. Clinical and Laboratory Standards Institute (CLSI). *Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard—Second Edition*. CLSI Document H26-A2. Wayne, PA: CLSI; 2010.

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Abbott GmbH & Co. KG