Anti-HBs

**For the Quantitative Measurement
of Antibodies to Hepatitis B Surface Antigen
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Anti-HBs

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of *antibodies against hepatitis B surface antigen (Anti-HBs)* in human serum or plasma (heparinized or EDTA), as an aid in the determination of immune status to the hepatitis B virus.

Catalog Number: **L2KAH2** (200 tests)

Test Code: **aHB** Color: **Dark Pink**

Summary and Explanation

Hepatitis B virus (HBV) is the sole human pathogen in the family of hepatitis-associated DNA viruses, and is found world-wide. Distribution of HBV infection will vary among geographical areas and population groups. Transmission of the virus is due to parenteral contact, through the exchange of blood or blood products, sexual contact, and perinatal spread from mother to newborn.^{1,2} Clinical manifestations range from mild asymptomatic infections to severe fulminant hepatitis.^{1,2,3} Over 90% of infected adults will have an acute self-limiting infection,¹ with jaundice and abnormal liver function. Recovery occurs without any chronic sequelae.^{1,2}

Chronic liver disease, a condition in which infection persists for more than six months, a known sequela of a hepatitis B infection, is usually progressive.^{1,2} The risk of developing the chronic carrier state is more likely to follow infection acquired in childhood than as an adult.^{4,5} In chronic HBV carriers, there is no evidence of continued hepatic damage,^{1,2} however, the infection persists and the carrier maintains the ability to transmit the virus.²

Availability of recombinant HBV vaccines, and the recommendation of universal immunization for infants and other high-risk persons has aided in the prevention of HBV infections. In addition, treatment with alpha-interferon to relieve symptoms is available. Results have shown positive response to treatment in 40–50% of

selected individuals with chronic active hepatitis B.^{4,5}

Classification of a hepatitis B infection requires the identification of several serological markers expressed during three phases (incubation, acute and convalescent) of the infection. The first marker to appear during the incubation phase is HBsAg, and indicates an ongoing infection with HBV.^{1,2,4} Antibodies to HBsAg generally appear after HBsAg has been cleared from the blood stream, usually 6 months after infection, and its presence represents recovery and immunity. However, in a few patients known to have antibodies to HBsAg, subclinical infections have developed.⁵ The presence of HBsAg antibodies should not be used as the sole marker in determining a prior hepatitis B infection.

Principle of the Procedure

IMMULITE 2000 Anti-HBs is a solid-phase, two-step chemiluminescent enzyme immunoassay.

The solid phase, a polystyrene bead, is coated with purified HBsAg subtypes ad and ay. An alkaline phosphatase-labeled HBsAg, supplied in the reagent, is added to the polystyrene bead. Unbound enzyme conjugate is then removed by a centrifugal wash. Finally, chemiluminescent substrate is added and undergoes hydrolysis in the presence of alkaline phosphatase; the photon output, as measured by the luminometer, is related to the presence of antibodies to HBsAg in the sample.

Incubation Cycles: 2 × 30 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples should not be used in this assay.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving

anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Anti-HBs has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 50 μ L serum or plasma (heparinized or EDTA). (See Alternate Sample Types section.)

Storage: 3 days at 2–8°C.⁶
For longer storage: at –20°C.⁷

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Certain components contain HBsAg which has been inactivated by proven, documented methods. However, always handle them as if capable of transmitting infectious agents.

The concentration of anti-hepatitis B surface antigen in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Anti-HBs EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Anti-HBs Bead Pack (L2AH12)

With barcode. 200 beads coated with purified inactivated human HBsAg subtypes ad and ay. Stable at 2–8°C until expiration date.

L2KAH2: 1 Pack

Anti-HBs Reagent Wedge (L2AHA2)

With barcode. 11.5 mL of a protein-based buffer, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to purified inactivated human HBsAg subtypes ad and ay in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KAH2: 1 Wedge

Before use, tear off the top of the label at the perforations without damaging the barcode on the main label. Remove the foil seal from the top of the Reagent Wedge, and snap the sliding cover down into the ramps on the reagent lid.

Anti-HBs Adjustors (LAHL, LAHH)

Two vials (Low and High), 2 mL each, of human serum reactive to HBsAg in a buffer, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KAH2: 1 set

Anti-HBs Controls (LAHC1, LAHC2, LAHC3)

Three vials, 4 mL each, with barcodes.

LAHC1 (Negative Control): One vial containing human serum nonreactive to HBsAg, with preservative. **LAHC2, LAHC3 (Low Postive Control, Positive Control):** Two vials containing human serum reactive to HBsAg, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KAH2: 1 set

Refer to the control insert for concentration levels.

Aliquot Labels with barcodes are supplied with the kit, for use with Adjustors and Controls. Before use, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water, test tubes

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Adjustment Interval: 4 weeks

Quality Control Samples: The Controls supplied with the kit should be used as quality control material to monitor assay performance.

Interpretation of Results

The IMMULITE 2000 Anti-HBs calibration employs a stored master curve, generated by the four-parameter logistic mathematical model based on the dose-CPS (counts per second) relationship during the calibration process.

Reactive: A result of greater than or equal to 10 mIU/mL (WHO 1st IRP 26-1-77) indicates that antibodies to HBsAg are present and were detected in the patient sample. In vaccinated individuals this usually indicates protection against infection.

Nonreactive: A result of less than 10 mIU/mL indicates that antibodies to HBsAg were not detected in the patient sample, or are below the protective level for immunity.

The results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Anti-HBs EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

Individuals acutely infected with the hepatitis B virus will exhibit anti-HBs approximately two weeks after the disappearance of HBsAg. This antibody response will reach peak levels after several months and gradually decline over a period of years. Persons who have been vaccinated against HBV will also have detectable levels of anti-HBs.

Limitations

The measurement of Anti-HBs may be affected by the presence of elevated levels of human albumin.

The results of the test must be taken within the context of the patient's clinical and vaccination history, symptomology and other laboratory findings.

Hemolyzed or grossly contaminated samples may give erroneous results.

A nonreactive result does not indicate that the patient was not infected with HBV or was not vaccinated against HBV. The patient sample should be tested for the presence of other serological markers.

Testing for anti-HBs antibodies alone is not sufficient in determining previous infections.

Sodium citrate will affect the measurement of antibodies against hepatitis B surface antigen. Use only serum, EDTA plasma or heparinized plasma.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: 3–2000 mIU/mL (WHO 1st IRP 26-1-77)

Analytical Sensitivity: ≤ 3.0 mIU/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Bilirubin: Presence of unconjugated bilirubin in concentrations up to 400 mg/L has no effect on results, within the precision of the assay.

Hemolysis: May interfere with the assay causing inconsistent or erroneous results. (See "Hemolysis" table.)

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: Samples (n = 46) were collected into plain, heparinized and EDTA vacutainer tubes.

(Heparin) = 1.03 (Serum) + 14 mIU/mL
r = 0.97

(EDTA) = 1.05 (Serum) + 16 mIU/mL
r = 0.88

Means:

350 mIU/mL (Serum)
376 mIU/mL (Heparin)
384 mIU/mL (EDTA)

Method Comparison 1: The IMMULITE 2000 Anti-HBs procedure was compared to IMMULITE Anti-HBs Kit on 59 patient samples. (Concentration Range: approximately 3–2000 mIU/mL. See "Method Comparison 1" graph.)
By linear regression:

(IML 2000) = 1.00 (IML) – 3.4 mIU/mL
r = 0.993

Means:

343 mIU/mL (IMMULITE 2000)
346 mIU/mL (IMMULITE)

Method Comparison 2: The IMMULITE 2000 Anti-HBs procedure was compared to the Elecsys Anti-HBs immunoassay on 775 samples, consisting of seroconversion panels, potentially crossreactive samples, samples from patients with resolved HBV infection, from vaccinated individuals and hospitalized patients, and routine samples. The table below presents the results of this study.

	No. of samples		Sensitivity	Specificity	PPV	NPV
	anti-HBs reactive	anti-HBs non-reactive				
IMMULITE 2000 Anti-HBs						
Reactive	386	1	99.5%	99.7%	99.7%	99.5%
Non-reactive	2	386				
Elecsys Anti-HBs						
Reactive	383	5	99.0%	98.7%	98.7%	99.0%
Non-reactive	4	383				

References

- 1) Locarnini SA, Gust ID. Hepadnaviridae: hepatitis B virus and the delta virus. In: Balows A, et al, editors. Laboratory diagnosis of infectious diseases: principles and practices. New York: Springer-Verlag, 1988: 750-96.
- 2) Follett EAC. Diagnosis of hepatitis B infection. In: Young H, McMillan A, editors. Immunological

diagnosis of sexually transmitted diseases. New York: Marcel Dekker, 1988: 433-49. 3) Hollinger FB, Dienstag JL. Hepatitis B and D viruses. In: Lennette EH, et al, editors. Manual of clinical microbiology. 6th ed. Washington, D.C.: American Society for Microbiology, 1995:1033-49. 4) Nowicki MJ, Balistreri WF. Hepatitis A to E: building up the alphabet. Contemporary Peds 1992: 118-28. 5) Zuckerman AJ, et al. Hepatitis B virus and hepatitis D virus. In: Principles and practice of clinical virology. 2nd ed. New York: John Wiley & Sons, 1992: 153-72. 6) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: WB Saunders, 1995:354-60. 7) Tietz NW, editor. Clinical guide to laboratory tests. 3rd ed. Philadelphia: W.B. Saunders, 1995:322-4. 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8. 9) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

siemens.com/healthcare


The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2008 Siemens Healthcare Diagnostics. All rights reserved.