Guideline[®]

V3 catheter

English/Instructions for Use.....	1
Česky/Návod k použití.....	2
Dansk/Brugsanvisning.....	3
Nederlands/Gebruiksaanwijzing.....	3
Eesti/Kasutusjuhised.....	4
Suomi/Käyttöohjeet.....	5
Français/Mode d'emploi.....	6
Deutsch/Gebrauchsanweisung.....	7
Ελληνικά/Οδηγίες χρήσης.....	8
Magyar/Használati utasítás.....	9
Italiano/Istruzioni per l'uso.....	10
Latviski/Lietošanas instrukcija.....	11
Lietuvių/Naudojimo instrukcijos.....	12
Norsk/Bruksanvisning.....	12
Polski/Instrukcja obsługi.....	13
Português/Instruções de utilização.....	14
Русский/Инструкция по применению.....	15
Español/Instrucciones de uso.....	16
Svenska/Bruksanvisning.....	17
Türkçe/Kullanım Talimatları.....	18
Slovenčina/Návod na použitie.....	19
Tiếng Việt/Hướng dẫn sử dụng.....	20
Hrvatski/Upute za uporabu.....	21
Srpski/Uputstvo za upotrebu.....	21
Română/Instrucțiuni de utilizare.....	22
Українська/Інструкції для застосування.....	23
ภาษาไทย/คำแนะนำสำหรับการใช้งาน.....	24

GuideLiner[®] V3 catheter

Instructions For Use

USA CAUTION

Federal law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION

The GuideLiner catheter is a single lumen rapid exchange catheter offered in sizes compatible with 5F, 6F, 7F, and 8F guide catheters and may be placed over a standard length guidewire. The larger sizes of GuideLiner catheters are intended to be used within the proximal portions of the coronary vasculature to provide support and/or facilitate use of multiple interventional devices. The 150cm device has a 125cm stainless steel shaft section covered on the distal 17cm with a semi-circular polymer. The steel shaft is followed distally by a 25cm lumen section wiped with silicone.

The GuideLiner catheter has two platinum-iridium marker bands, which enable visibility while using standard fluoroscopic methods. The distal marker band is located on the distal tip. The proximal marker band is located near the collar. The device has two positioning marks located at 95cm (single mark) and 105cm (double mark) from the distal tip.


The GuideLiner catheter is delivered through a guide catheter resulting in an inner diameter that is approximately 1 French smaller than the guide catheter. The GuideLiner catheter has a proximal tab which indicates guide catheter compatibility and the resulting GuideLiner catheter inner diameter.

The GuideLiner catheter has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Model	Compatible Guide Catheter	GuideLiner Minimum I.D.	GuideLiner Tip O.D. (X)
5569 5F	≥ 5F (≥ 0.056" / 1.42mm I.D.)	0.046" (1.17mm)	0.053" (1.35mm)
5570 5.5F	≥ 6F (≥ 0.066" / 1.68mm I.D.)	0.051" (1.30mm)	0.063" (1.60mm)
5571 6F	≥ 6F (≥ 0.070" / 1.78mm I.D.)	0.056" (1.42mm)	0.067" (1.70mm)
5572 7F	≥ 7F (≥ 0.078" / 1.98mm I.D.)	0.062" (1.57mm)	0.075" (1.90mm)
5573 8F	≥ 8F (≥ 0.088" / 2.24mm I.D.)	0.071" (1.80mm)	0.085" (2.16mm)

INDICATIONS

GuideLiner catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.

CONTRAINDICATIONS

The GuideLiner catheters are contraindicated in vessels less than 2.5mm in diameter and in vessels in the neurovasculature or the venous system.

WARNINGS

The GuideLiner catheter is provided sterile for single use only. Reuse of single-use devices creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Never advance the GuideLiner catheter into a vessel without a leading guidewire as vessel damage may result.

Never advance the GuideLiner catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury and/or occlusion may result. If pressure in a vessel dampens after inserting the GuideLiner catheter, withdraw the catheter until the pressure returns to normal.

Due to the size and non-tapered tip of the GuideLiner catheter, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel damage.

PRECAUTIONS

Do not use the GuideLiner catheter if the packaging has been damaged. A damaged package could result in a breach of sterility or device damage.

Inspect the GuideLiner catheter prior to use for any bends or kinks. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumen should be flushed with sterile, heparinized saline prior to use to prevent clot formation and to ensure the catheter is free from debris which could be introduced into the body resulting in an embolism and/or occlusion.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system.

Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking. Do not apply torque to the catheter during delivery, as catheter damage may result.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response, as catheter damage or vessel injury may occur.

Never advance the GuideLiner catheter more than 10cm beyond the tip of the guide catheter as the GuideLiner catheter may become lodged in the guide catheter making it difficult to remove.

Do not withdraw an undeployed stent back into the GuideLiner catheter when the catheter is in the body, as it may result in dislodging the stent. Instead, simultaneously pull both the GuideLiner catheter and undeployed stent back into the guide and remove together.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the GuideLiner catheter include, but are not limited to, the following:

- Cardiac arrest
- Embolism
- Infection
- Myocardial infarction
- Slow-flow/Occlusion
- Stent dislodgement
- Thrombosis
- Vessel dissection
- Vessel perforation

CLINICAL PROCEDURE

The GuideLiner catheter should be used by physicians trained on the procedures for which it is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

PACKAGE CONTAINS:

1x Catheter

OTHER ITEMS REQUIRED BUT NOT PROVIDED:

- Guide catheter with an inner diameter large enough to accommodate the specific model of GuideLiner catheter in use
- Y-adaptor with hemostasis valve (Tuohy-Borst type)
- Guidewire with diameter ≤ 0.014" / 0.36mm
- Sterile syringe (for flushing)
- Sterile heparinized saline (for flushing)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the GuideLiner packaging and components for damage.
2. Using sterile technique, transfer the dispenser coil with the GuideLiner catheter into the sterile field.
3. Remove the GuideLiner catheter from the dispenser coil and thoroughly flush the lumen of the GuideLiner catheter from the distal tip with sterile, heparinized saline solution.

DEPLOYMENT PROCEDURE

Deploy the GuideLiner catheter according to the following steps:

1. Secure the previously inserted guidewire and backload the distal tip of the GuideLiner catheter onto the guidewire and advance until the catheter is just proximal to the hemostasis valve.
2. Open the hemostasis valve and advance the GuideLiner catheter through the hemostasis valve and into the guide catheter.
3. Under fluoroscopy, advance the GuideLiner catheter beyond the distal tip of the guide catheter and into the desired location within the vessel.

WARNING: Never advance the GuideLiner catheter into a vessel with an effective diameter less than 2.5mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the GuideLiner catheter, withdraw the catheter until the pressure returns to normal.

WARNING: Due to the size and non-tapered tip of the GuideLiner catheter, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Teleflex[®]

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis, MN 55369 USA
Phone: (866) 246-6990
Fax: (866) 804-9881

EC REP

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- Using fluoroscopy, confirm the desired position of the GuideLiner catheter in the vessel.
- If performing an interventional procedure, backload the interventional device over the existing guidewire and advance the device through the guide catheter and GuideLiner catheter into the desired vascular space.
NOTE: If a second wire is used during the intervention and encounters resistance within the guide catheter, pull the wire back several centimeters and slowly readvance.
- Tighten the Y-adaptor hemostasis valve securely on the proximal shaft of the GuideLiner catheter to prevent back-bleeding.
- Perform the catheterization procedure. After completing the procedure, remove the GuideLiner catheter prior to removing the guide catheter from the vessel.
- Dispose of the GuideLiner catheter following standard hospital procedures.

STORAGE & HANDLING

No special storage or handling conditions.

LIMITED WARRANTY

Vascular Solutions LLC warrants that the GuideLiner catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special or consequential damages arising from the use of the GuideLiner catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions LLC has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions LLC.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS LLC.

PATENTS AND TRADEMARKS

May be covered by one or more U.S. or international patents.

See: www.teleflex.com/patents-intv

GuideLiner is a registered trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

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DEVICE DESCRIPTION

The Twin-Pass dual access catheter is a dual lumen catheter designed for use in the coronary and/or peripheral vasculature. The Twin-Pass catheter consists of an over-the-wire (OTW) lumen that runs the length of the catheter and a rapid exchange (RX) delivery lumen on the distal segment. The Twin-Pass catheter has a working length of 135cm and is compatible with various guidewire diameters (see specification table). The Twin-Pass catheter also has hydrophilic coating on the distal 18cm of the catheter.

The Twin-Pass catheter has white positioning marks located at 95cm (single mark) and 105cm (double marks) from the distal tip, respectively. The Twin-Pass catheter has a marker band located 1mm from the distal tip and a second marker band located 11mm from the distal tip to identify the distal end of the OTW lumen.


The Twin-Pass catheter has been sterilized with ethylene oxide.

STERILE EO

SPECIFICATIONS

Model Number	5200
Max. guidewire diameter - RX lumen	0.014" / 0.36mm
RX lumen inner diameter	0.017" / 0.43mm
Max. guidewire diameter - OTW lumen	0.014" / 0.36mm
OTW lumen inner diameter	0.018" / 0.46mm
Proximal shaft outer diameter	0.038" / 0.97mm (2.9F)
Distal shaft outer diameter	0.047" / 1.19mm (3.5F)
Distal tip outer diameter	0.026" / 0.66mm (2F)
Minimum guide catheter inner diameter	0.066" / 1.68mm

Model 5200

INDICATIONS

The Twin-Pass catheter is intended to be used in conjunction with steerable guidewires in order to access discrete regions of the coronary and peripheral arterial vasculature, to facilitate placement and exchange of guidewires and other interventional devices, and for use during two guidewire procedures.

The Twin-Pass catheter is also used to sub selectively infuse/deliver diagnostic or therapeutic agents.

CONTRAINDICATIONS

The Twin-Pass catheter is contraindicated for high pressure injections and for use in the cerebral vasculature.

WARNINGS

The Twin-Pass catheter is provided sterile for single use only. Reuse of a single-use device creates a potential risk of patient or user infections and may compromise device functionality, which may lead to illness or serious patient injury.

Do not advance the Twin-Pass catheter without a guidewire in place through an RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in vessel injury.

Do not exceed the maximum recommended infusion pressure of 300 psi (2068 kPa) when injecting through the Twin-Pass catheter as catheter damage or vessel injury could result.

If a guidewire of appropriate diameter cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter damage or vessel injury could result. Identify and resolve the cause of the blockage, or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in catheter damage or vessel injury.

PRECAUTIONS

Do not use the Twin-Pass catheter if the packaging has been damaged. A damaged package could indicate a breach of sterility or device damage.

Inspect the catheter prior to use. Do not use a damaged catheter. Vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumens must be flushed with sterile, heparinized saline prior to use to ensure catheter is free from debris which could be introduced into the body.

Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the guidewire lumen, difficulty while inserting the catheter or guidewires, or pressure damping.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response as catheter damage or vessel injury may occur.

Check that all fittings and attachments are secure and observe that the syringe is free of air bubbles, prior to injection, so that air is not introduced into the catheter during use. Failure to do so may result in air being introduced into the vasculature.

ADVERSE EFFECTS

Potential adverse effects that may be associated with the use of the Twin-Pass catheter include, but are not limited to, the following:

- Effusion/tamponade
- Embolism
- Hypertension
- Hypotension
- Infection
- Inflammatory response
- Myocardial infarction
- Tachycardia
- Thrombosis
- Vessel dissection
- Vessel perforation
- Vessel rupture
- Vessel spasm

CLINICAL PROCEDURE

The Twin-Pass catheter should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

- Single-use Twin-Pass catheter

Other items required but not provided:

- Guide catheter with an inner diameter large enough to accommodate the specific model of Twin-Pass catheter in use
- Appropriate guidewires
- 10ml syringe (for system flushing)
- Sterile heparinized saline (for system flushing)
- Transducer (if measuring pressure)

PREPARATIONS FOR USE

1. Prior to use, carefully inspect the Twin-Pass catheter packaging and components for damage.
2. Using sterile technique, transfer the dispenser coil with the Twin-Pass catheter into the sterile field.
3. Attach a 10ml syringe filled with sterile, heparinized saline to the luer-lock guidewire entry port of the Twin-Pass catheter and thoroughly flush the catheter.
4. Flush the Twin-Pass catheter with sterile, heparinized saline to activate the hydrophilic coating.
5. Remove the Twin-Pass catheter from the dispenser coil and inspect it for any bends or kinks.

DEPLOYMENT PROCEDURE

Deploy the Twin-Pass catheter according to the following steps:

1. Backload the rapid exchange segment of the Twin-Pass catheter onto the guidewire that is already in place in the distal vasculature.

WARNING: Do not advance the Twin-Pass catheter without a guidewire in place through an RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in vessel injury.
2. Under fluoroscopy, carefully advance the Twin-Pass catheter into the desired distal vascular space.

WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in catheter damage or vessel injury.
3. To deliver an exchange-length guidewire through the OTW lumen:
 - Load the guidewire into the proximal lumen of the Twin-Pass catheter.
 - Advance the guidewire until it exits the OTW lumen into the distal vascular space.
 - Fix both guidewires using standard guidewire exchange techniques and carefully withdraw the Twin-Pass catheter until

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the distal tip exits the hemostatic valve and both wires can be secured.

- To deliver fluid through the OTW lumen:
 - Create an air-free, fluid-to-fluid connection between the OTW lumen and the injection syringe.
 - Aspirate prior to injection, and observe the syringe to remove any air bubbles prior to injecting into the catheter.
 - Using the injection syringe, manually infuse the diagnostic or therapeutic agent.

WARNING: Do not exceed the maximum recommended infusion pressure of 300 psi (2068 kPa) when injecting through the Twin-Pass catheter. Catheter damage or vessel injury could result.

Precaution: Check that all fittings and attachments are secure so that air is not introduced into the catheter during the procedure.

- Dispose of Twin-Pass catheter following standard hospital procedure.

STORAGE AND HANDLING

No special or handling conditions.

LIMITED WARRANTY

Vascular Solutions LLC warrants that the Twin-Pass catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions LLC to be defective in workmanship or materials. Vascular Solutions LLC shall not be liable for any incidental, special or consequential damages arising from the use of the Twin-Pass catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

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See: www.teleflex.com/patents-intv

Twin-Pass is a registered trademark of Teleflex Innovations S.à.r.l., Vascular Solutions LLC, or Teleflex Medical, each a part of Teleflex Incorporated.

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
EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

No.**CE 647900**

Issued To:

**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**


In respect of:

Design, Development, and Manufacturing of Guidewires, Hemostatic and Embolic Agents for Peripheral Use, Fluid Delivery Kits, Surgical Laser Procedure Kits, Intravascular Catheters Used for Diagnostic, Angiographic, Embolic Delivery and/or Embolectomy, Support and Exchange Applications, Vascular Access Devices including ultrasound guided devices, and Accessories for use in Catheterization Procedures, and Retrieval Devices

Those aspects of Annex II related to maintaining the sterility of accessory Torque Devices, SmartNeedle Monitor Covers, and Auto-Fill Syringe Kits. Those aspects of Annex II related to maintaining the sterility and metrology function of inflation devices

on the basis of our examination of the quality assurance system under the requirements of Council Directive 93/42/EEC, Annex II excluding section 4. The quality assurance system meets the requirements of the directive. For the placing on the market of class III products an Annex II section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):

Albert Roossien, Regulatory Lead

First Issued: **2016-06-10**Date: **2019-03-11**Expiry Date: **2022-07-01**

...making excellence a habit.™

Page 1 of 1

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

This certificate was issued electronically and is bound by the conditions of the contract.

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
Biomerics LLC 12560 Fletcher Lane, Suite 300 Rogers MN 55374 USA	Manufacture
Cartika Medical, Inc. 6550 Wedgwood Road North Suite 300 Maple Grove Minnesota 55311 USA	Manufacture
E-Beam Services, Inc. 2775 Henkle Drive Lebanon Ohio 45036 USA	E Beam Sterilization

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
E-Beam Services, Inc. 3400 Union Street Lafayette Indiana 47905 USA	E Beam Sterilization
Galt Medical Corp 2220 Merritt Drive Garland, TX 75041 USA	Manufacture
Harland Medical Systems 7418 Washington Ave S Eden Prairie MN 55344 USA	Manufacture

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
Harmac Medical Products IDA Business Park Castlerea Roscommon Ireland	Manufacture
Isomedix Operations, Inc. 380 90th Avenue NW Minneapolis MN 55433 USA	ETO Sterilization
Lake Region Medical 340 Lake Hazeltine Drive Chaska Minnesota 55318 USA	Manufacture

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
Merit Medical Systems, Inc. 1600 West Merit Parkway South Jordan UT 84095 USA	Manufacture
Merit Medical Systems, Inc. 14646 Kirby Drive Houston Texas 77047 USA	Manufacture
Phillips-Medisize Ireland High Road Letterkenny Co. Donegal Ireland	Manufacture

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
Polymicro Technologies a Subsidiary of Molex 18019 N. 25th Ave. Phoenix AZ 85023 USA	Manufacture
Precision Wire Components 10230 SW Spokane Court Tualatin Oregon 97062 USA	Manufacture
Sterigenics US, LLC 7775 South Quincy Street Willowbrook Illinois 60527 USA	ETO Sterilization

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:	Service(s) supplied
Synergy Health Westport Ltd Lodge Road Westport County Mayo Ireland	Gamma Sterilization
Teleflex Medical IDA Business and Technology Park Dublin Road Athlone Co. Westmeath Ireland	EU Representative
Vascular Solutions LLC 14005 13th Avenue Plymouth Minnesota 55441 USA	Manufacture

...making excellence a habit.™

EC Certificate - Full Quality Assurance System

Directive 93/42/EEC on Medical Devices, Annex II excluding Section 4

List of Significant Subcontractors

Recognised as being involved in services relating to the product covered by:

Certificate No: **CE 647900**
Date: **2019-03-11**
Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Subcontractor:

Service(s) supplied

Vascular Solutions LLC
6401 Sycamore Court North
Maple Grove
Minnesota
55369
USA

Regulatory Compliance

Vascular Solutions LLC
6420 Sycamore Lane North
Minneapolis
Minnesota
55369
USA

Manufacture

...making excellence a habit.™

EC Certificate - Full Quality Assurance System Certificate History

Certificate No: **CE 647900**
 Date: **2019-03-11**
 Issued To: **Vascular Solutions LLC**
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Date	Reference Number	Action
10 June 2016	8466178	First Issue – Transfer from another Notified Body.
15 December 2016	8649395	Addition of Sterigenics Willowbrook, IL as a significant subcontractor for EO sterilization.
27 January 2017	8588049	Amend the scope to include guidewires. Addition of Precision Wire Components to the certificate as significant subcontractors for manufacturing.
30 June 2017	8760886	Certificate renewal. Subcontractor name change from Steris Isomedix Services to Isomedix Operations, Inc.
26 July 2018	8940684	Change of legal manufacturer name from "Vascular Solutions, Inc." to "Vascular Solutions LLC." Change of multiple subcontractor names from "Vascular Solutions, Inc." to "Vascular Solutions LLC." Change EU Representative to Teleflex Medical (Ireland). Removal of Vascular Solutions Zerusa Limited. Clarifications to sterilization services supplied for several subcontractors, and corrections to several subcontractor addresses.
Current	8558535	Traceable to NB 0086.

EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

No.**CE 678704**

Issued To:


**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

In respect of:

GuideLiner V3 catheter

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):

Albert Roossien, Regulatory Lead

First Issued: **2018-10-11**Date: **2019-03-11**Expiry Date: **2023-10-10**

...making excellence a habit.™

Page 1 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 678704

Issued To:

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Model Number	Device Name	Size	Intended Purpose per IFU	Classification
5569	GuideLiner V3 Catheter	5F	The GuideLiner catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices	Class III
5570		5.5F		
5571		6F		
5572		7F		
5573		8F		

First Issued: **2018-10-11**

Date: **2019-03-11**

Expiry Date: **2023-10-10**

...making excellence a habit.™

Page 2 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 678704

Issued To:

**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

Certificate History

Date	Reference Number	Action
11 October 2018	8780387	Initial release as per up-classification review.
01 February 2019	8958037	Addition of alternate sterilizer Sterigenics US, LLC.
Current	8558535	Traceable to NB 0086.

First Issued: **2018-10-11**

Date: **2019-03-11**

Expiry Date: **2023-10-10**

...making excellence a habit.™

Page 3 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

No.**CE 678705**

Issued To:


**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

In respect of:

Turnpike Catheters

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):

Albert Roossien, Regulatory Lead

First Issued: **2018-09-26**Date: **2019-03-11**Expiry Date: **2023-09-25**

...making excellence a habit.™

Page 1 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 678705

Issued To:

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Turnpike Catheters

Catalog Number	Device Name	Model, Type	Intended purpose per IFU	Classification
5640 and 5641	Turnpike Spiral	Soft Tip, Spiral Shaft	The Turnpike Spiral catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.	Class III
5642 and 5643	Turnpike	Soft Tip, Smooth Shaft	The Turnpike catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.	Class III

First Issued: **2018-09-26**

Date: **2019-03-11**

Expiry Date: **2023-09-25**

...making excellence a habit.™

EC Design-Examination Certificate

Supplementary Information to CE 678705

Issued To:

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Catalog Number	Device Name	Model, Type	Intended purpose per IFU	Classification
5638 and 5639	Turnpike LP	Soft Tip, Smooth Shaft	The Turnpike LP catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.	Class III
5621 and 5622	Turnpike Gold	Threaded Tip, Spiral Shaft	The Turnpike Gold catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.	Class III

First Issued: **2018-09-26**

Date: **2019-03-11**

Expiry Date: **2023-09-25**

...making excellence a habit.™

Page 3 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Supplementary Information to CE 678705

Issued To:

**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

Certificate History

Date	Reference Number	Action
26 September 2018	8780390	First issue.
Current	8558535	Traceable to NB 0086.

First Issued: **2018-09-26**

Date: **2019-03-11**

Expiry Date: **2023-09-25**

...making excellence a habit.™

Page 4 of 4

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.

EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

No.

CE 678707

Issued To:


**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

In respect of:

Twin-Pass Dual Access Catheters

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0086):

Stewart Brain, Head of Compliance & Risk -
Medical Devices

First Issued: **2018-09-26**

Date: **2018-09-26**

Expiry Date: **2023-09-25**

...making excellence a habit.™

Page 1 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

EC Design-Examination Certificate

Supplementary Information to CE 678707

Issued To:

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

Twin-Pass Dual Access Catheters

Catalog Number	Device Name	Model, Type	Intended purpose per IFU	Classification
5200	Twin-Pass Dual Access Catheter Twin-Pass Torque	Twin-Pass	<p>The Twin-Pass catheters are intended to be used in conjunction with steerable guidewires in order to access discrete regions of the coronary and peripheral arterial vasculature, to facilitate placement and exchange of guidewires and other interventional devices, and for use during two guidewire procedures.</p> <p>The Twin-Pass is also used to subselectively infuse/deliver diagnostic or therapeutic agents.</p>	Class III
5201		Twin-Pass Torque	<p>The Twin-Pass catheter is intended to access discrete regions of the coronary and/or peripheral vasculature. It may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents.</p>	Class III

First Issued: **2018-09-26**

Date: **2018-09-26**

Expiry Date: **2023-09-25**

...making excellence a habit.™

EC Design-Examination Certificate

Supplementary Information to CE 678707

Issued To:

**Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA**

Certificate History

Date	Reference Number	Action
Current	8780393	First issue.

First Issued: **2018-09-26**

Date: **2018-09-26**

Expiry Date: **2023-09-25**

...making excellence a habit.™

Page 3 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that:

Vascular Solutions LLC
6464 Sycamore Court North
Minneapolis
Minnesota
55369
USA

DUNS Number: F000466

Holds Certificate No:

MDSAP 690601

Statement of Conformity: The company listed on this certificate has been audited to and found to conform with the following criteria: ISO 13485:2016 and Australia - Therapeutic Goods (Medical Devices) Regulations, 2002, Schedule 3 Part 1 (excluding Part 1.6) - Full Quality Assurance Procedure; Brasil - RDC ANVISA n. 16/2013, RDC ANVISA n. 23/2012, RDC ANVISA n. 67/2009; Canada - Medical Devices Regulations - Part 1 - SOR 98/282; Japan - MHLW Ministerial Ordinance 169, Article 4 to Article 68, PMD Act; USA - 21 CFR 820, 21 CFR 803, 21 CFR 806, 21 CFR 807 - Subparts A to D

Design, Development, Manufacturing and Distribution of Guidewires, Vascular and Topical Hemostasis Devices, Embolectomy Catheters, Surgical Lasers, Intravascular Catheters, Vascular Access Devices, Vascular Introducer Kits. Fluid Delivery Kits Embolization Devices and Retrieval Devices.

For and on behalf of BSI:

Stewart Brain, Head of Compliance & Risk - Medical Devices

Original Registration Date: 2019-05-07

Effective Date: 2019-10-03

Expiry Date: 2022-02-18

BSI Group America Inc. is an MDSAP authorized auditing organization

Page: 1 of 2

...making excellence a habit.™

Certificate No: **MDSAP 690601**

Location	Registered Activities
Vascular Solutions LLC 6464 Sycamore Court North Minneapolis Minnesota 55369 USA DUNS Number: F000466	Design, Development, Manufacturing and Distribution of Guidewires, Vascular and Topical Hemostasis Devices, Embolectomy Catheters, Surgical Lasers, Intravascular Catheters, Vascular Access Devices, Vascular Introducer Kits. Fluid Delivery Kits Embolization Devices and Retrieval Devices.
Vascular Solutions LLC 14005 13th Avenue Plymouth Minnesota 55441 USA DUNS Number: F000466	Receiving, Receiving Inspection, Non Sterile Kit Assembly, Final Inspection and Release, Warehousing / Distribution.
Vascular Solutions LLC 6401 Sycamore Court North Maple Grove Minnesota 55369 USA	Research and Design, HR, Marketing, Finance.
Vascular Solutions LLC 6420 Sycamore Lane North Minneapolis Minnesota 55369 USA	Catheter Component Manufacturing including extrusion, coils, snare loops, braiding, and laser processing.

Original Registration Date: 2019-05-07

Effective Date: 2019-10-03

Expiry Date: 2022-02-18

Page: 2 of 2

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](https://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Declaration of Conformity

Manufacturer:	Vascular Solutions LLC 6464 Sycamore Court North Minneapolis Minnesota 55369 USA
European Representative:	Teleflex Medical IDA Business and Technology Park Dublin Road Athlone Co. Westmeath Ireland
Product Name:	Turnpike Catheter Turnpike Gold Catheter Turnpike LP Catheter Turnpike Spiral Catheter
Reference/Catalog Number):	5642, 5643 5621, 5622 5638, 5639 5640, 5641
Generic Indication(s):	The Turnpike catheters are intended to be used to access discrete regions of the coronary and/or peripheral vasculature. They may be used to facilitate placement and exchange of guidewires and to subselectively infuse/deliver diagnostic and therapeutic agents
Classification:	Class III, Rule 6
GMDN Code	10691
Conformity Assessment Route:	Annex II Section 4 of the EC-Directive (Council Directive 93/42/EEC of 14 June 1993)
Declaration of Conformity:	We hereby declare that the distributed CE marked products, specified on this Declaration of Conformity, conform to the required documentation in accordance with Annex II Section 4 of the EC – Directive, the Council Directive 93/42/EEC of 14 June 1993, concerning medical devices.
Notified Body:	BSI Group the Netherlands BV Say Building John M. Keynesplein 9 1066 EP Amsterdam Notified Body Number: 2797
EC Certificate(s)	MD648016: Quality Management System – ISO 13485:2016 & EN ISO 13485:2016 CE647900: Full Quality Assurance System, 93/42/EEC Annex II, excluding Section 4 CE678705: Design Examination Certificate, 93/42/EEC Annex II, Section 4 MDSAP 690601: ISO 13485:2016
Date CE Mark Affixed:	Turnpike Catheter: December 12, 2014 Turnpike Gold Catheter: February 9, 2015 Turnpike LP Catheter: December 11, 2015 Turnpike Spiral Catheter: December 12, 2014

Place of Issuance:	Vascular Solutions LLC 6464 Sycamore Court North Minneapolis Minnesota 55369 USA		
Signature of Authorized Individual:		Date of Issue:	January 23, 2020
Name (printed)	Brenda Johnson, Director of Regulatory		

REVISION HISTORY

Rev	Reason
A	Initial release of Class III DoC
B	Update to reflect new Notified Body address and number. Removed list of standards as they will be maintained with the Essential Requirements Checklist. Updated GMDN code to reflect recent revisions to the codes.

Declaration of Conformity

Manufacturer:	Vascular Solutions LLC 6464 Sycamore Court North Minneapolis Minnesota 55369 USA		
European Representative:	Teleflex Medical IDA Business and Technology Park Dublin Road Athlone Co. Westmeath Ireland		
Product Name:	GuideLiner V3 catheter		
Reference/Catalog Number:	5569, 5570, 5571, 5572, 5573		
Generic Indication(s):	GuideLiner catheters are intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.		
Classification:	Class III, Rule 6 in Annex IX of 93/42/EEC (MDD)		
GMDN Code	17846		
Conformity Assessment Route:	Annex II, Section 4 of 93/42/EEC (MDD)		
Declaration of Conformity:	We hereby declare that the distributed CE marked products, specified on this Declaration of Conformity, conform to the required documentation in accordance with Annex II of the EC – Directive, the Council Directive 93/42/EEC of 14 June 1993, concerning medical devices.		
Notified Body:	BSI Group the Netherlands BV Say Building John M. Keynesplein 9 1066 EP Amsterdam The Netherlands Notified Body Number: 2797		
EC Certificate(s)	MDSAP690601: QMS ISO 13485:2016 CE647900: Full Quality Assurance System, 93/42/EEC Annex II, excluding Section 4 CE678704: Design Examination Certificate, 93/42/EEC Annex II, Section 4		
Date CE Mark Affixed:	April 8, 2013 (models 5570 through 5573) June 9, 2015 (model 5569)		
Place of Issuance:	Vascular Solutions LLC 6464 Sycamore Court North Minneapolis Minnesota 55369 USA		
Signature of Authorized Individual:		Date of Issue:	March 9, 2020
Name (printed)	Brenda Johnson, Director of Regulatory		

REVISION HISTORY

Rev	Reason
A	New release for Class III certification
B	Add MDSAP Certificate. Remove standards applied table.