Medtronic

ClosureFast™

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! USA Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

CE
2797
UK
CA
0086

Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.

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English

Instructions for use

Sterile, single use only. For use with 0.64 mm (0.025 in) guidewire and ClosureRFG generator.

Note: Thoroughly read all instructions, including the ClosureRFG radiofrequency generator operator's manual, before using the ClosureFast system. Observe all warnings, precautions, and cautions noted throughout these instructions. Failure to do so may result in patient complications.

ClosureFast catheters

The three models of ClosureFast catheters are the ClosureFast catheters with a 7 cm heating element (model numbers: CF7-7-60 and CF7-7-100) and the ClosureFast catheter with a 3 cm heating element (model number: CF7-3-60) devices. See *Table 1* for product configuration details.

Device description

The ClosureFast system consists of 2 main components: a ClosureFast catheter and a ClosureRFG radiofrequency generator. The catheter and the integrated connection cable are provided sterile, and are single-use, disposable devices. The catheter's function is to provide thermal energy to the desired treatment site via radiofrequency heating of the catheter heating element and to relay temperature back to the radiofrequency generator. The radiofrequency generator remains out of the sterile field during use and is provided nonsterile. The catheter is connected to the radiofrequency generator via the connection cable.

Intended purpose

The ClosureFast catheter is intended to be used in conjunction with the ClosureRFG generator in order to occlude the treated vessel via a fibrotic seal.

Intended patient population

This device is intended for use in adult patients with chronic venous insufficiency.

Indications for use

The ClosureFast catheter is intended for endovascular coagulation of blood vessels in patients with superficial vein reflux.

Contraindications

The ClosureFast catheter is contraindicated for use in patients with thrombus in the target vein segment.

Note: No data exists regarding the use of this catheter in patients with documented peripheral arterial disease. Take the same care in the treatment of patients with significant peripheral arterial disease as for a traditional vein ligation and stripping procedure.

Clinical benefits

The clinical benefits of vein closure with the ClosureFast catheter in the treatment of superficial vein reflux include the following:

- Resolution of the symptoms of chronic venous insufficiency
- Venous leg ulcer healing
- Improvement in quality of life
- Increased mobility

Intended users

The ClosureFast catheter is intended to be used by physicians who are familiar with duplex ultrasound, vascular interventions, and chronic venous insufficiency.

Performance characteristics

The ClosureFast system consists of 2 main components: a ClosureFast catheter and a ClosureRFG radiofrequency generator. The catheter provides thermal energy to deliver segmental ablation to the desired treatment site through radiofrequency heating of the catheter heating element. The catheter relays the temperature back to the radiofrequency generator in order for the radiofrequency generator to automatically control RF power and maintain the preset temperature. The catheter and the integrated connection cable are provided sterile and are single-use, disposable devices. The catheter is connected to the radiofrequency generator via the connection cable.

Warnings

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Treatment of a vein located near the skin surface can result in skin burn if the skin is not protected with fluid infiltration.
- Nerve injury can occur from thermal damage to adjacent nerves. The risk of nerve injury may be higher with treatment at or below the knee, or without perivenous fluid infiltration.

Precautions

- Do not bend the catheter into a tight radius. Kinking the shaft can damage the catheter.
- When inserting the catheter into the vein, ensure that the guidewire does not protrude from the catheter tip.
- The vein wall may be thinner in an aneurysmal segment. To effectively occlude a vein with an aneurysmal segment, additional tumescent infiltration may be needed over the aneurysmal segment. The treatment of the vein should include segments proximal and distal to the aneurysmal segment.

Potential adverse events

The potential complications include, but are not limited to, the following list:

- Adjacent nerve injury or nerve damage
- Hematoma or seroma
- Infection
- Inflammation
- Pain
- Phlebitis
- Pulmonary embolism
- Skin burn or discoloration
- Thrombosis
- Vessel perforation

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

How supplied

The device is supplied sterile and is sterilized using ethylene oxide.

Storage

Store the device in a dry area away from sunlight.

Recommended supplies and equipment

- Tilt table
- Duplex ultrasound scanner
- Sterile ultrasound gel
- Sterile ultrasound transducer cover
- Percutaneous access needle
- 7 Fr introducer sheath (7 cm or 11 cm length)
- 0.64 mm (0.025 in) guidewire

Generator set-up

Note: Refer to the ClosureRFG generator operator's manual.

- Plug the power supply into the generator.
- Turn on the power by using the power switch.
- Confirm the software version on screen. Refer to *Table 1* for software version compatibility.
- The default target temperature setting is 120 °C for the ClosureFast catheter.

Note: When using the RFG2 model, the user can select from a range of 95 °C to 120 °C. Note that the default settings do not appear until a catheter is connected to the radiofrequency generator. Treatment settings can be adjusted according to physician preference (RFG2 only). See the radiofrequency generator operator's manual for instructions on how to change the settings.

Table 1. ClosureFast system specifications

Catheter model	CF7-7-60	CF7-7-100	CF7-3-60
Introducer sheath (minimum ID size)	7 Fr (2.3 mm)	7 Fr (2.3 mm)	7 Fr (2.3 mm)
Insertable length	60 cm	100 cm	60 cm
Heating element diameter	2.3 mm	2.3 mm	2.3 mm
Heating element length	7 cm	7 cm	3 cm
Maximum power setting	40 W	40 W	18 W
Default target temperature setting	120 °C	120 °C	120 °C
Software version	RFG2 model	4.0.0 or higher	4.0.0 or higher
	RFG3 model	1.11.0 or higher	1.11.0 or higher

Directions for use

Note: Use aseptic technique.

Patient preparation

- Flush disposable accessories with sterile, physiologic saline (0.9% sodium chloride).
- If local anesthetic is used, administer local anesthetic at the vein access site. Mild sedation can also be given.
Note: Venospasm may hinder the ability to access the target vein. Avoid any factors that can induce venospasm such as certain drugs, a cold environment, or patient anxiety.
- Position the patient for vein access. Lowering the patient's legs below the level of the heart increases vein diameter, which may facilitate vein access.
- Access the vein to be treated via a percutaneous stick, using an access needle or via a small cut-down.
- Prepare and place an introducer sheath per manufacturer instructions for use.

Catheter inspection and preparation

- Inspect the outer box for any signs of visible damage. Remove the pouch from the box and inspect the pouch for damage (such as tears or punctures).

Caution: Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.

- Open the pouch and remove the tray.
- Inspect the catheter.
- Using aseptic technique, pass the cable connector out of the tray and connect the cable to the radiofrequency generator.
Caution: Do not allow cable connector to come in contact with any fluid.
- Remove the catheter from the tray and place the catheter into the sterile field.
- Using sterile, physiologic saline, flush and fill the catheter lumen, cap the lumen at the handle, and wipe the outer surface of the catheter shaft.
Caution: Use of a flush through the catheter while the heating element is active will heat the fluid exiting the end of the catheter. Avoid fluid delivery through the catheter when the tip of the catheter is near an area that should not be thermally coagulated.
- Insert the catheter into the introducer sheath and advance the catheter tip to the most proximal point of treatment. Navigate the catheter to the treatment site using ultrasound guidance, palpation, or a guidewire.
- If using a center lumen guidewire to assist catheter advancement, refer to the manufacturer's instructions for use of the guidewire. After removal of the guidewire, re-flush the catheter lumen with sterile, physiologic saline and cap the lumen at the end of the catheter.
Caution: Do not advance the catheter or guidewire against resistance or vein perforation can occur.

Tumescent infiltration and catheter tip position

- Use tumescent infiltration of dilute local anesthetic or saline into the perivascular space to create a circumferential fluid layer around the vessel to be treated. To achieve contact between the catheter heating element and the vein wall, an approximate volume of 10 cm³ per cm of vein to be treated is recommended. Infiltrate up to approximately 5 cm distal to the saphenofemoral junction (SFJ) or saphenopopliteal junction (SPJ); infiltration beyond the junction will be performed after confirmation of final tip position.
Note: When the vein is located near the skin surface, create a subcutaneous distance of >1 cm between the anterior vein wall and the skin by using tumescent infiltration of saline or dilute local anesthetic solution.
- Verify the catheter tip position using the measurement calipers of the ultrasound machine. When treating either the great saphenous vein (GSV) or small saphenous vein (SSV), place the tip 2 cm inferior to the junction.
- Infiltrate tumescent fluid beyond the junction using ultrasound guidance.

Treatment

- Place the patient's legs above the level of the heart to facilitate vein collapse, apposition, and exsanguination (Trendelenburg position).
- While maintaining catheter tip position, partially withdraw the introducer sheath until the sheath hub is aligned with the first visible shaft marker or draw a mark on the skin at the level of the first visible shaft marker. Secure sheath to skin (optional).
Caution: Do not treat with the heating element in the deep venous system.
- Create a near-bloodless field by applying external compression along the full length of the heating element using the ultrasound transducer longitudinally aligned with the heating element, plus 2-3 fingertips (to approximate the length of the heating element) of compression distal of the transducer.
Caution: Failure to compress the vein over the full length of the heating element may result in inconsistent effectiveness or possible catheter damage.
- Enable radiofrequency energy delivery by pressing the "RF Power" button on the radiofrequency generator. This action will cause the "RF Power" button to start blinking. If the "RF Power" button does not light or start blinking, observe any displayed message and respond. Refer to the radiofrequency generator operator's manual for further detail (RFG2).
- Initiate radiofrequency energy delivery by pressing the button on the catheter handle. During treatment, energy delivery can be turned off by pressing the button on the catheter handle again. (RFG2 only: radiofrequency energy can be initiated using the "START RF" button. Radiofrequency energy can be turned off using the "STOP RF" button or the "RF Power" button.)
Note: For the ClosureFast catheter with the 7 cm heating element, power typically begins at 40 W and drops below 20 W within 10 seconds if compression is located correctly and the vein segment being treated has been properly exsanguinated. For the ClosureFast catheter with the 3 cm heating element, power typically begins at 18 W and drops below 10 W within 10 seconds if compression is located correctly and the vein segment being treated has been properly exsanguinated.
Note: If the set temperature is not reached within 5 seconds after radiofrequency energy delivery initiation, or if the power level stays high, there may be flow within the vein that is cooling the treatment segment. Terminate radiofrequency energy delivery, verify effectiveness of exsanguination methods and proper tip position, correct as necessary, and reinitiate treatment of the segment.
Note: Continuous temperature readings below the set temperature can result in incomplete treatment. If this occurs, stop the treatment and reconfirm vessel apposition to the catheter heating element and absence of blood flow in the vessel segment to be treated. Apply more firm external compression, if needed, and retreat the segment.
Caution: If treatment is halted due to non-uniform temperature, remove the catheter and inspect the heating element for damage. If the catheter is damaged, replace it.
Caution: Failure to respond to alerts can result in severe damage to the catheter.
- After the treatment time interval, radiofrequency energy delivery will terminate automatically. When using a ClosureFast catheter with the 7 cm heating element, deliver a second energy cycle to the one segment closest to the SFJ. When using a ClosureFast catheter with the 3 cm heating element, a second energy cycle may be given at the physician's discretion.
- Radiofrequency energy delivery may be repeated in a given vein segment at the physician's discretion.
Caution: Do not administer more than 3 energy delivery cycles at any given vein segment because the amount of thermal ablation depth will not significantly increase.
Caution: Do not re-advance catheter through an acutely treated vein segment.
- Quickly withdraw the catheter until the next visible shaft marker is aligned with the hub of the sheath.
Note: Some friction between the vein wall and the catheter after a heating cycle is normal and may be noticed while withdrawing the catheter.
- Treat the next vein segment according to *Step 3* through *Step 8*, repeating the compression, treatment, and indexing sequence until all segments are treated. The diagonal lines and the printed numbers on the outside of the catheter shaft, which correlate to the introducer sheath length, indicate the last full treatment segment when they are fully visible.
Note: Use the triple shaft mark located 3 cm from the heating element to determine the minimum distance from the heating element to the puncture site.
Caution: Treatment with the heating element inside the sheath or outside the body can result in skin burn or catheter damage.
- Remove the catheter and introducer sheath from the vein and power off the radiofrequency generator. Evaluate treated vein segments with ultrasound to determine treatment outcome.
Caution: There is no re-treatment algorithm with the ClosureFast catheter; do not re-advance the catheter through an acutely treated vein segment.

11. Obtain hemostasis at the access site.
12. Apply a multilayer compression wrap from foot to groin.

Follow-up care

1. Instruct patient to ambulate frequently and refrain from strenuous activities or heavy lifting for several days.
2. Postoperative compression for at least 1 week is recommended.
3. A follow-up examination within 72 hours should include an assessment to ensure that there is no thrombus extension into deep veins.

Disposal

Caution: Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic, Inc.
Medtronic B.V.
Medtronic Limited