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 06-02806-11

IDEXX

Avian Influenza Virus Antibody Test Kit

For veterinary use only.

Name and Intended Use

IDEXX AI is IDEXX's enzyme immunoassay for the detection of antibody to Avian Influenza Virus (AI) in chicken serum.

General Information

Domestic and wild avian species are affected by avian influenza viruses. The disease is characterized by a wide range of responses that include virtually no clinical signs to high mortality. Respiratory signs are common, along with drop in egg production, greenish diarrhea, bloodstained nasal and oral discharges and cyanosis and edema of the head, comb and wattle. Due to the variation and severity of clinical symptoms, serological testing produces significant advantages to detection of infected birds. Monitoring for exposure of a flock to influenza is facilitated by the measurement of antibody to avian influenza in serum.

Descriptions and Principles

This assay is designed to measure the relative level of antibody to AI in chicken serum. Viral antigen is coated on 96-well plates. Upon incubation of the test sample in the coated well, antibody specific to AI forms a complex with the coated viral antigens. After washing away unbound material from the wells, a conjugate is added which binds to any attached chicken antibody in the wells. Unbound conjugate is washed away and enzyme substrate is added. Subsequent color development is directly related to the amount of antibody to AI present in the test sample.

Reagents

Volume

	Reagents	Volume
1	AI Antigen Coated Plate	5
2	Positive Control — diluted chicken anti-AI serum; preserved with ProClin™ 150	1 x 1.9 mL
3	Negative Control — diluted chicken serum non-reactive to AI, preserved with ProClin™ 150	1 x 1.9 mL
4	Conjugate — (goat) anti-chicken/(goat) anti-turkey: HRPO Conjugate; preserved with ProClin™ 150	1 x 50 mL
5	Sample Diluent — preserved with ProClin™ 150	1 x 235 mL
A	TMB Substrate	1 x 60 mL
B	Stop Solution	1 x 60 mL

Note: See table at the end of the insert for a description of symbols used on the insert and labels of this kit.

Storage

Store the reagents at 2–8°C. Reagents are stable until expiration date, provided they have been stored properly.

Materials Required but Not Provided

- Precision micropipettes and multi-dispensing micropipettes
- Disposable pipette tips
- 96-well microplate reader (equipped with 650 nm filter)
- Microplate washer (manual, semi-automatic or automatic system)
- Use only distilled or deionized water for preparation of the reagents used in the test
- Vortex or equivalent

Precautions and Warnings

- Handle all biological material as potentially infectious. The antigen used in the kit reagents may not be completely inactivated.
- Wear protective gloves / protective clothing / eye or face protection when handling samples and reagents.
- Refer to the product Material Safety Data Sheet for additional information.
- See the end of this insert for reagent hazard and precaution warnings.

Laboratory Practices

- Optimal results will be obtained by strict adherence to this protocol. Careful pipetting, timing, and washing throughout this procedure are necessary to maintain precision and accuracy. Use a separate pipette tip for each sample and control.
- Do not expose TMB solution to strong light or any oxidizing agents. Handle TMB solution with clean glass or plastic ware.
- All wastes should be properly decontaminated prior to disposal. Dispose of contents in accordance with local, regional, and national regulations.
- Care should be taken to prevent contamination of kit components. Do not pour unused reagents back into containers.
- Do not use kit past expiration date.

Preparation of Samples

Dilute test samples five hundred fold (1:500) with sample diluent prior to being assayed (e.g., by diluting 1 μL of sample with 500 μL of Sample Diluent). **NOTE: DO NOT DILUTE CONTROLS.** Be sure to change tips for each sample. Samples must be thoroughly mixed prior to dispensing into the coated plate.

Test Procedure

All reagents must be allowed to come to 18–26°C before use. Mix reagents by gentle inverting or swirling.

- 1 Obtain antigen-coated plate(s) and record the sample position.

- 2 Dispense 100 μL of UNDILUTED Negative Control (NC) into duplicate wells.

- 3 Dispense 100 μL of UNDILUTED Positive Control (PC) into duplicate wells.

- 4 Dispense 100 μL of DILUTED sample into appropriate wells. Samples may be tested in duplicate, but a single well is acceptable.

- 5 Incubate for 30 minutes (± 2 minutes) at 18–26°C.

- 6 Remove the solution and wash each well with approximately 350 μL of distilled or deionized water 3–5 times. Avoid plate drying between plate washings and prior to the addition of the next reagent. Tap each plate onto absorbent material after the final wash to remove any residual wash fluid.

- 7 Dispense 100 μL of Conjugate into each well.

- 8 Incubate for 30 minutes (± 2 minutes) at 18–26°C.

- 9 Repeat step 6

- 10 Dispense 100 μL of TMB Substrate into each well.

- 11 Incubate for 15 minutes (± 1 minute) at 18–26°C.

- 12 Dispense 100 μL of Stop Solution into each well.

- 13 Measure and record absorbance values at 650nm, A(650).

14 Calculation:

Controls

$$NC\bar{x} = \frac{NC1 A(650) + NC2 A(650)}{2}$$

$$PC\bar{x} = \frac{PC1 A(650) + PC2 A(650)}{2}$$

Validity criteria

$$PC\bar{x} - NC\bar{x} > 0.075$$

$$NC\bar{x} \leq 0.150$$

For invalid assays, technique may be suspect and the assay should be repeated following a thorough review of the package insert.

Samples

$$S/P = \frac{\text{Sample Mean} - NC\bar{x}}{PC\bar{x} - NC\bar{x}}$$

The presence or absence of antibody to AI is determined by relating the A(650) value of the unknown to the Positive Control mean. The Positive Control is standardized and represents significant antibody levels to AI in chicken serum. The relative level of antibody in the sample is determined by calculating the sample to positive (S/P) ratio.

15 Interpretation:

Negative

Positive

$$S/P \leq 0.50$$

$$S/P > 0.50$$

ELISA positive samples should be confirmed by additional serological testing such as Agar Gel Precipitation Test (AGP).

Note: U.S. Customers must submit all positive samples and/or results to the National Veterinary Services Laboratories for H and N titration.

Note: IDEXX has instrument and software systems available which calculate results and provide data summaries.

For technical assistance:

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: idexx.com/contacttpd

VLN/PCN: 313/5004.00

IDEXX and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries

© 2020 IDEXX Laboratories, Inc. All rights reserved.

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WARNING / ATTENTION / ATENCIÓN / ATENÇÃO / ACHTUNG

H317 / H412 / P280 / P302 + P352 / P333 + P313

Positive Control / Conjugate – May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

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H317 / P280 / P302 + P352 / P333 + P313

Negative Control / Sample Diluent – May cause an allergic skin reaction. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

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IDEXX Europe B.V.