Medtronic

Export™

Aspiration catheter
Катетър за аспирация
Aspirační katétr
Aspirationskateter
Aspirationskatheter
Καθετήρας αναρρόφησης
Catéter de aspiración
Aspiratsioonikateeter
Aspirointikatetri
Cathéter d'aspiration
Aspiracijski kateter
Aspirációs katéter
Kateter aspirasi
Catetere per aspirazione
Íšsiurbimo kateteris
Aspirācijas katetrs
Aspiratiekatheter
Aspirasjonskateter
Cewnik do aspiracji
Cateter de aspiração
Cateter de aspirație
Аспирационный катетер
Aspiračný katéter
Aspiracijski kateter
Aspirationskateter
Aspirasyon kateteri
Аспіраційний катетер

BRA

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
English

1 Description

Carefully read all instructions before use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

The Export aspiration catheter is a single operator design, dual lumen catheter, compatible with 0.36 mm (0.014 in) guidewires and minimum 6 Fr guide catheter (1.78 mm [0.070 in] minimum ID). It has a distal radiopaque tip marker, proximal luer-lock port, and preloaded stiffening stylet (as applicable). The proximal luer-lock port is used to connect the aspiration line (supplied) and the aspiration syringe (supplied). Also, an appropriately sized syringe filled with an infusate may be attached to the aspiration line for infusing fluids.

2 Intended use of the device

2.1 Intended purpose

The Export aspiration catheter's intended purpose is to perform aspiration and containment of soft thrombo-embolic material during nonroutine percutaneous transluminal angioplasty (PTA) or stenting procedures in coronary arteries, with the effect of improving perfusion. The aspiration catheter is also intended to improve perfusion by selectively infusing and delivering diagnostic or therapeutic agents to the desired vasculature site during PTA or stenting procedures.

2.2 Intended patient population

The Export aspiration catheter is intended for nonroutine use in adult patients presenting with cardiac symptoms resulting from the reduction of perfusion caused by thrombus formation.

2.3 Indications for use

Export aspiration catheters are indicated in clinical conditions in which a vascular occlusion has occurred in vessels ≥ 2.0 mm diameter, due to soft thrombo-embolic formation located in the coronary arteries.

2.4 Contraindications

There are no known contraindications.

2.5 Clinical benefits

The Export aspiration catheter aspirates and contains embolic material such as thrombus or debris during percutaneous transluminal angioplasty or stenting procedures in coronary arteries. The aspiration catheter can also selectively infuse and deliver diagnostic or therapeutic fluids during these procedures. The catheter removes thrombi to restore perfusion, thereby re-establishing coronary artery flow. Therefore, the primary clinical benefit resulting from reperfusion is the improvement of symptoms related to coronary occlusion (chest pain, shortness of breath, fatigue, nausea, headache, or crushing pressure in the shoulders, neck, or jaw).

The performance of Export aspiration catheters has been evaluated in multiple clinical trials, which included a variety of patient demographics, conditions, and pathologies. The average procedural success rate specific to Export aspiration catheters to aid in the restoration of blood flow in coronary procedures is 90.5 % as shown in *Table 1*.

Table 1. Clinical benefits

Patient-relevant outcomes	Export aspiration catheters; N=6584 patients
Performance measure	
Procedural success rate ^a	90.5%

^a Defined as delivery of the aspiration catheter to the target lesion; containing and aspirating embolic material; selective infusion or delivery of diagnostic or therapeutic fluids.

2.6 Intended users

Interventional procedures using the Export aspiration catheter should only be performed by physicians who have experience with interventional techniques in the vascular system.

3 Performance characteristics

The Export aspiration catheter extracts embolic material during percutaneous vascular interventional treatments. The aspiration syringe creates suction to pull embolic material through the lumen of the aspiration catheter, for containment and removal from the treatment site. The Export catheter can also deliver diagnostic or therapeutic fluids to the treatment site.

4 Summary of safety and clinical performance

The Summary of Safety and Clinical Performance (SSCP) for this device can be found at <https://ec.europa.eu/tools/eudamed>. Search for the SSCP using the manufacturer, device name, and the following elements as applicable: device model, reference number, and catalog number.

5 Warnings

- If the Export aspiration catheter is to be used with an 0.36 mm (0.014 in) guidewire, consult the guidewire instructions for use (IFU) before use.
- If using the Export aspiration catheter with a minimum 6 Fr guide catheter (1.78 mm [0.070 in] minimum ID), consult the guide catheter IFU before use.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Visually inspect all sterile-barrier packaging before use. If the device is damaged or the integrity of the sterilization barrier has been compromised, do not use the product. Contact your Medtronic representative for return information.
- Handle the Export aspiration catheter and components carefully. Before use, inspect the Export aspiration catheter carefully for bends, kinks, or other damage. Do not use a damaged Export aspiration catheter or damaged components.
- Do not attempt to advance or retract the Export aspiration catheter against resistance until the cause of resistance has been determined by fluoroscopy or other means. Maintain visualization of the distal end of the guidewire while using the Export aspiration catheter, specifically in the region distal to the guide catheter and proximal to the Export aspiration catheter marker band. Should a loop or excessive slack in the guidewire be visualized, do not attempt to retract the Export aspiration catheter into the guide catheter until the loop or excessive slack is eliminated. Retracting the Export aspiration catheter without straightening the guidewire may

potentially tear the guidewire lumen, which could cause dislodgement of device components into the vasculature.

- Due to the presence of DEHP (Di[2-ethylhexyl] phthalate) in the aspiration line of this device, it is not recommended for use in minors and pregnant or nursing women.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.
- Only physicians who have received appropriate training should use the Export aspiration catheter.
- **If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.**

6 Precautions

- The Export aspiration catheter and components should be used in conjunction with fluoroscopic guidance and proper anticoagulation agents.
- When using the Export aspiration catheter for fluid delivery, do not exceed the maximum flow rate specified for the device (see *Table 2*).
- This device is supplied sterile, and is designed and intended for single use only. **Do not** resterilize or reuse the device.
- As in any elective coronary intervention, it is recommended that the patient have a mean systolic blood pressure greater than or equal to 90 mmHg whether or not IV pressors or intra-aortic balloon pump augmentation are used.

7 Adverse events

As with most percutaneous interventions, other potential adverse events associated with the use of this device include the following conditions: death, myocardial infarction, coronary or bypass graft thrombosis or occlusion, myocardial ischemia, stroke/CVA, emergent or nonemergent surgery, emboli (air, tissue or thrombotic), dissection, arterial perforation, arterial rupture, ventricular fibrillation, hemorrhage, hypotension, pseudoaneurysm at the access site, arteriovenous fistula, infection, allergic reaction, and bleeding complications at the access site.

8 Device materials

Some patients are allergic to, or can become sensitive to, nickel.

Patients can be exposed to the following materials in the system, which are carcinogenic, mutagenic, and toxic for reproduction (CMR):

Material	CAS number
Benzyl butyl phthalate (BBP)	85-68-7
Diethylphthalate (DOP) / Di(2-ethylhexyl) phthalate (DEHP)	117-81-7
Cobalt	7440-48-4

9 How supplied

The Export aspiration catheter is supplied with an aspiration line, two locking aspiration syringes, a preloaded stylet (for Export Advance only), and a filter cup. The device and components are supplied sterile and nonpyrogenic in the sterile pouch. The Export aspiration catheter, including all of its components, is intended for single use only. Standard storage conditions are sufficient to safeguard the device. Store the device in the original packaging at room temperature in a dry place.

10 Required items (not supplied)

- Minimum 6 Fr arterial or venous introducer sheath and minimum 6 Fr guide catheter (1.78 mm [0.070 in] minimum ID) in the appropriate configuration to cannulate the vessel (preferably with side holes if ostial narrowing is present or the guide catheter is occlusive).
- 0.89 mm (0.035 in) guidewire may be used to support the aspiration lumen during introduction of the aspiration catheter into a tortuous target vessel. If the Export catheter includes a preloaded stylet, remove the stylet before inserting the 0.89 mm (0.035 in) guidewire into the aspiration lumen.
- Rotating hemostatic valve (RHV)
- Heparinized normal saline
- 10 mL syringe
- Flush needle
- Appropriately-sized syringe for fluid delivery (optional)

11 Instructions for use

11.1 Preparing the Export aspiration catheter

The Export aspiration catheter is supplied with an aspiration line, 2 locking aspiration syringes, a preloaded stylet (for Export Advance only), and a filter cup.

1. Remove the components and the hoop containing the Export aspiration catheter from the packaging.
2. Attach an aspiration syringe filled with approximately 10 mL of heparinized saline to the flushing luer port on the hoop and completely flush the hoop to activate the hydrophilic coating.
3. Remove the Export aspiration catheter from the hoop. Fill the aspiration syringe with approximately 5 to 10 mL of heparinized saline and attach the aspiration line and syringe to the hub of the Export aspiration catheter (or to the hub of the preloaded stylet, if using the stylet to deliver the Export aspiration catheter).
4. Open the stopcock on the aspiration line and flush the entire length of the Export aspiration catheter using all of the heparinized saline contained in the aspiration syringe. Close the stopcock.
5. Verify that the stopcock on the aspiration line is in the closed position. Retract the plunger of the aspiration syringe and pull until it locks at the fully extended position.

Note: When using the Export aspiration catheter for fluid delivery, perform steps 1 through 4, and then connect an infusion syringe of preferred volume or an intravenous line to the aspiration line stopcock. Do not use the locking aspiration syringe for fluid infusion.

6. Thoroughly flush the guidewire lumen with heparinized normal saline, using a 10 mL syringe and flush needle.

11.2 Using the Export aspiration catheter during an interventional procedure

1. Perform aspiration using the Export aspiration catheter:
 - a. Load the proximal end of an in situ 0.36 mm (0.014 in) guidewire into the distal end of the guidewire lumen of the prepared Export aspiration catheter. Optionally, disconnect the aspiration line from the catheter hub and load the proximal end of an in situ 0.89 mm (0.035 in) guidewire into the distal end of the aspiration lumen of the aspiration catheter. Advance the distal tip of the aspiration catheter to the distal opening of the in situ guide catheter.
 - b. Using fluoroscopy, advance the Export aspiration catheter beyond the guide catheter and position the distal tip marker proximal to the desired aspiration site. Stop advancing the Export aspiration catheter if any resistance is encountered.
 - c. Begin aspiration by opening the stopcock on the aspiration line. Blood will enter the aspiration syringe until the entire vacuum is gone (or the aspiration syringe is filled).

Notes: When a stiffening stylet is used, disconnect the aspiration line from the stylet hub, remove the stylet from the catheter, and connect the aspiration line to the catheter hub before beginning aspiration. When using a 0.89 mm (0.035 in) guidewire to support the aspiration lumen, remove the 0.89 mm (0.035 in) guidewire from the aspiration catheter and connect the aspiration line to the catheter hub before beginning aspiration.

If blood does not begin filling the syringe within 5 s, check placement of the guide catheter tip. Unseat the guide catheter, if necessary, to resume flow.

If no blood is aspirated as a result of unseating the guide catheter, turn the stopcock off and remove the Export aspiration catheter. Outside of the patient, either flush the aspiration lumen or use a new Export aspiration catheter. **Do not** flush the system while the catheter is still inside the patient vasculature.

During performance testing, the Export aspiration catheter demonstrated the ability to evacuate fluid and debris at a minimum rate of 0.5 mL/s.

Caution: When a stiffening stylet is used, remove it before aspiration.

- d. After completing the aspiration process, turn the aspiration line stopcock to close off the aspiration syringe.
2. Remove the Export aspiration catheter.

Slowly retract and remove the Export aspiration catheter. If necessary, loosen the Tuohy Borst adapter of the RHV (rotating hemostasis valve) to allow easy withdrawal of the distal shaft.

Note: If the aspiration catheter will be reinserted in the patient after its removal, remove the aspiration syringe and reflush the aspiration lumen and the guidewire lumen of the Export aspiration catheter with heparinized saline. If a stylet is included with the device, reinsert the stiffening stylet before reinserting the catheter into the patient. Empty the aspiration syringe, reattach it to the aspiration line, close the stopcock, and retract the plunger to the fully extended lock position.
 3. Remove the catheters and follow standard hospital practice for management of the insertion site.

11.3 Optional steps: Using the Export aspiration catheter for fluid delivery

1. After removing diagnostic or interventional devices, load and advance the prepared Export aspiration catheter to the tip of the guide catheter.

Note: The aspiration syringe is replaced by an appropriately sized syringe filled with the fluid to be delivered or by an intravenous line.

2. Advance the Export aspiration catheter, position the distal tip marker at the vessel site, and begin fluid delivery.

Caution: When delivering fluid through the Export aspiration catheter, do not exceed infusion pressures of 6.8 atm (100 psi) or the maximum recommended flow rate per the following chart:

Table 2. Export aspiration catheter flow rate

Maximum recommended Export aspiration catheter flow rate	
Saline	4.0 mL/s
60% Ionic Contrast Media	1.8 mL/s
76% Ionic Contrast Media	0.5 mL/s

12 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2022 Medtronic
Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
Medtronic, Inc.
Medtronic B.V.