

**P.2 Quantitative method**  
In the qualitative method, calculate the mean OD450nm values for the Calibrators 0 and 10 U/ml and then check that the assay is valid.

**Example of calculation:**  
The following data must not be used instead of real figures obtained by the user:

Calibrator 0 U/ml:	0.020 - 0.024 OD450nm
Mean Value:	0.022 OD450nm
Lower than 0.150 - Accepted	0.350 - 0.330 OD450nm
Calibrator 10 U/ml:	0.340 OD450nm
Mean Value:	Higher than Cal 0 + 0.200 - Accepted
Calibrator 100 U/ml:	2.845 OD450nm
Higher than 1.000 - Accepted	

**Q. INTERPRETATION OF RESULTS**

**Q.1 Quantitative results**  
For qualitative interpretations, the medical literature generally considers positive samples showing a concentration of HBeIgt  $\geq 10$  PEI U/ml.  
Test results are therefore interpreted as a ratio of the sample OD450nm and the OD450nm of the Cal 10 PEI U/ml (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

**Q.2 Quantitative results**

The calibration curve is used to determine the concentration of IgM antibodies to HBeIgt in samples.  
Samples with a concentration lower than 5 PEI U/ml are considered negative for HBeIgt.  
Samples with a concentration between 5 and 10 PEI U/ml are considered in a gray-zone.  
In the follow up of chronic hepatitis, however, values higher of 5 PEI U/ml may be considered positive for HBeIgt, when in presence of other clinical signs.  
Samples with a concentration higher than 10 PEI U/ml are considered positive for HBeIgt.

**Important general notes:**

1. When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to produce the calibration curve, calculate sample concentration and generate the correct interpretation of results.
2. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
3. A positive result is indicative of HBV infection and therefore the patient should be treated accordingly.
4. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
5. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

**R. PEI-CRM/ARTICES**  
Evaluation of Performance has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of MD Directive 98/79/EC).

**1. Limit of detection**  
The limit of detection of the assay has been calculated by means of:

- 1.1 the HBeIgt reference preparation supplied by Paul Ehrlich Institute Germany (HBe-Reference-serum-1GM 84), on which the Standard Curve has been calculated.
- 1.2 Accurin 113 (USA, N° A113-3001) supplied by Boston Biomedica Inc., USA

Results of Quality Control for three lots are given in the following tables:

Member	Lot #	Lot #	Lot #	Lot #	Lot #
PEI U/ml	OD450nm	S/Co	OD450nm	S/Co	OD450nm
100	2.752	8.9	2.883	9.7	2.911
50	1.517	6.2	1.972	6.7	2.053
20	0.980	3.2	0.914	3.1	1.005
10	0.544	1.6	0.513	1.7	0.592
5	0.310	1.0	0.295	1.0	0.321
2.5	0.185	0.5	0.149	0.5	0.181
1.25	0.084	0.3	0.084	0.3	0.093
negative	0.040	0.035	0.044	0.044	0.044

Member	Lot #	Lot #	Lot #	Lot #	Lot #
PEI U/ml	OD450nm	S/Co	OD450nm	S/Co	OD450nm
100	2.752	8.9	2.883	9.7	2.911
50	1.517	6.2	1.972	6.7	2.053
20	0.980	3.2	0.914	3.1	1.005
10	0.544	1.6	0.513	1.7	0.592
5	0.310	1.0	0.295	1.0	0.321
2.5	0.185	0.5	0.149	0.5	0.181
1.25	0.084	0.3	0.084	0.3	0.093
negative	0.040	0.035	0.044	0.044	0.044

Moreover the Bbis panel # PHE 102 was also examined in three lots of product; data are reported below with reference to a European kit (BBI's results).

Member	Lot #	Lot #	Lot #	Lot #	Lot #
S/Co	S/Co	S/Co	S/Co	S/Co	S/Co
01	11.3	10.0	8.7	6.1	2.0
02	11.3	10.0	8.7	6.1	2.0
03	9.5	7.2	4.4	3.0	1.0
04	5.8	3.4	4.4	3.0	1.0
05	11.3	11.4	11.2	5.1	1.1
06	12.1	11.8	11.6	4.1	1.1
07	0.1	0.1	0.1	0.2	0.2
08	9.2	8.5	8.8	2.3	2.3
09	11.9	10.2	11.9	4.2	4.2
10	11.7	10.7	10.8	2.8	2.8
11	5.9	5.8	5.8	2.1	2.1
12	12.1	11.4	11.7	3.2	3.2
13	7.0	6.0	6.0	3.6	3.6
14	7.0	6.0	6.0	3.6	3.6
15	12.4	11.5	11.8	4.5	4.5

**2. Diagnostic Sensitivity:**  
It is defined as the probability of the assay of scoring positive in the presence of the specific analyte.  
The diagnostic sensitivity has been tested internally and externally in a qualified Clinical Laboratory on panels of samples classified positive by a US FDA approved kit.  
Positive samples were collected from different patients and from different HBV pathologies (acute and chronic hepatitis).

An overall value  $\geq 98\%$  has been found in the study conducted on a total number of more than 200 samples.  
A Seroconversion panel produced by BBI, code # PHE 555A, has also been studied; results are reported below with reference to two commercial kits (BBI's results).

Member #	S/Co	Abbott EIA	Diagnost	S/Co	EIA
01	0.2	0.1	0.1	0.1	0.1
02	0.2	0.1	0.1	0.1	0.1
03	0.1	0.1	0.1	0.1	0.1
04	0.2	0.1	0.1	0.1	0.1
05	0.2	0.1	0.1	0.1	0.1
06	0.2	0.1	0.1	0.1	0.1
07	0.2	0.1	0.1	0.1	0.1
08	0.1	0.1	0.1	0.1	0.1
09	0.1	0.1	0.1	0.1	0.1
10	0.1	0.1	0.1	0.1	0.1
11	0.2	0.1	0.1	0.1	0.1
12	0.2	0.1	0.1	0.1	0.1
13	2.8	3.7	0.7	0.7	0.7
14	5.0	6.4	0.9	0.9	0.9
15	1.2	6.2	4.5	4.5	4.5
16	> 12	5.6	4.5	4.5	4.5
17	> 12	4.8	4.8	4.8	4.8
18	> 12	4.8	4.8	4.8	4.8
19	> 12	5.6	4.4	4.4	4.4
20	> 12	5.6	5.2	5.2	5.2

**3. Diagnostic Specificity:**

It is defined as the probability of the assay of scoring negative in the absence of the specific analyte.  
The diagnostic specificity has been determined internally and externally in a qualified Clinical Laboratory on panels of negative samples from normal individuals and blood donors, classified negative with a US FDA approved kit.  
A total number of more than 400 negative specimens were tested. A diagnostic specificity  $> 96\%$  has been found.  
Moreover, the diagnostic specificity was assessed by testing more than 50 potentially interfering specimens (other infectious diseases, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolyzed (hemic, etc.)).  
No interference was observed in the study.  
Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.  
Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples.

**4. Precision:**

It has been calculated on three samples examined in 16 replicate in three different runs, carried out on three different kits. The values found were as follows:

Member	1st run	2nd run	3rd run	Average Value
Cal 5 U/ml (N = 16)	0.324	0.308	0.321	0.318
Std.Deviation	0.022	0.018	0.024	0.022
CV %	6.8	5.7	7.5	6.7

Member	1st run	2nd run	3rd run	Average Value
Cal 50 U/ml (N = 16)	2.109	2.058	2.052	2.073
Std.Deviation	0.101	0.088	0.136	0.109
CV %	4.8	4.3	6.7	5.2

Member	1st run	2nd run	3rd run	Average Value
Cal 5 U/ml (N = 16)	0.257	0.253	0.244	0.252
Std.Deviation	0.025	0.024	0.024	0.024
CV %	8.3	8.0	7.3	8.0

Member	1st run	2nd run	3rd run	Average Value
Cal 50 U/ml (N = 16)	0.332	0.331	0.322	0.328
Std.Deviation	0.017	0.018	0.016	0.017
CV %	5.0	5.5	4.9	5.1

Member	1st run	2nd run	3rd run	Average Value
Cal 50 U/ml (N = 16)	2.311	2.208	2.212	2.245
Std.Deviation	0.110	0.090	0.095	0.098
CV %	4.7	4.1	4.3	4.4

Member	1st run	2nd run	3rd run	Average Value
Cal 5 U/ml (N = 16)	0.043	0.042	0.046	0.042
Std.Deviation	0.004	0.005	0.004	0.004
CV %	10.3	11.1	10.9	10.3

Member	1st run	2nd run	3rd run	Average Value
Cal 5 U/ml (N = 16)	0.320	0.326	0.314	0.320
Std.Deviation	0.023	0.024	0.026	0.024
CV %	7.1	7.4	8.2	7.6

Member	1st run	2nd run	3rd run	Average Value
Cal 50 U/ml (N = 16)	2.150	2.153	2.092	2.135
Std.Deviation	0.057	0.067	0.076	0.067
CV %	2.6	3.1	3.6	3.1

**5. LIMITATIONS**

Frozen samples containing fibrin particles or aggregates may generate false positive results.  
Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent distortion of the level of the analyte.  
This test is suitable only for testing single samples and not pooled ones.  
Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, seroprevalence, as well as other diagnostic data should be considered.

**REFERENCES**

1. Engvall E. and Perlman P., J. Immunol., 87: 874-875, 1971.
2. Engvall E. and Perlman P., J. Immunol. 109: 129-135, 1972.
3. Rerington J.S. and Klein J.O., in Infectious diseases of the fetus and newborn infant, Sanders, Philadelphia, London, Toronto.
4. Wik W.A., in "Essentials of Medical Microbiology", 2<sup>nd</sup> ed., pp. 723 G.B. Lippincott Company, Philadelphia, New York, St. José, Toronto.

**HBs Ab**

**A. INTENDED USE**  
Enzyme immunoassay (ELISA) for both the qualitative and quantitative determination of antibodies to the Surface Antigen of Hepatitis B Virus in human plasma and sera.  
For "in vitro" diagnostic use only.

**B. INTRODUCTION**  
The World Health Organization (WHO) defines Hepatitis B Virus infection as follows:

"Hepatitis B is one of the major diseases of mankind and is a serious global public health problem. Hepatitis means inflammation of the liver, and the most common cause is infection with one of 5 viruses, called hepatitis A, B, C, D, and E. All of these viruses can cause an acute disease with symptoms lasting several weeks including yellowing of the skin and eyes (jaundice), dark urine, extreme fatigue, nausea, vomiting and abdominal pain. It can take several months to a year to feel fit again. Hepatitis B virus can cause chronic infection in which the patient never gets rid of the virus and many years later develops cirrhosis of the liver or liver cancer."

HBV is the most serious type of viral hepatitis and the only type causing chronic hepatitis for which a vaccine is available. Hepatitis B virus is transmitted by contact with blood or body fluids of an infected person in the same way as human immunodeficiency virus (HIV), the virus that causes AIDS. However, HBV is 50 to 100 times more infectious than HIV. The main ways of getting infected with HBV are: (a) perinatal, from mother to baby at the birth; (b) child-to-child transmission; (c) unsafe injections and transfusions; (d) sexual contact.

Worldwide, most infections occur from infected mother to child, from child to child contact in household settings, and from reuse of unsterilized needles and syringes. In many developing countries, almost all children become infected with the virus. In many industrialized countries (e.g. Western Europe and North America), the pattern of transmission is different. In these countries, mother-to-infant and child-to-child transmission accounted for up to one third of chronic infections before childhood hepatitis B vaccination programmes were implemented. However, the majority of infections in these countries are acquired during young adulthood by sexual activity, and injecting drug use. In addition, hepatitis B virus is the major infectious occupational hazard of health workers, and most health care workers have received hepatitis B vaccine.

Hepatitis B virus is not spread by contaminated food or water, and cannot be spread casually in the workplace. High rates of chronic HBV infection are also found in the southern parts of Eastern and Central Europe, in the Middle East and Indian subcontinent, about 5% are chronically infected. Infection is less common in Western Europe and North America, where less than 1% are chronically infected.

Young children who become infected with HBV are the most likely to develop chronic infection. About 90% of infants infected during the first year of life and 50% of young children infected between 1 to 4 years of age develop chronic infection. The risk of death from HBV-related liver cancer or

cirrhosis is approximately 25% for persons who become chronically infected during childhood.

Chronic hepatitis B in some patients is treated with drugs called *interferon* or *lamivudine*, which can help some patients. Patients with cirrhosis are sometimes given liver transplants, with varying success. It is preferable to prevent this disease with vaccine than to try and cure it.

Hepatitis B vaccine has an outstanding record of safety and effectiveness. Since 1982, over one billion doses of hepatitis B vaccine have been used worldwide. The vaccine is given as a series of three intramuscular doses. Studies have shown that the vaccine is 95% effective in preventing children and adults from developing chronic infection if they have not yet been infected. In many countries where 8% to 15% of children used to become chronically infected with HBV, the rate of chronic infection has been reduced to less than 1% in immunized groups of children. Since 1991, WHO has called for all countries to add hepatitis B vaccine into their national immunization programmes."

Hepatitis B surface Antigen (HBsAg) is the major structural polypeptide of the envelope of the Hepatitis B Virus (HBV). This antigen is composed mainly of the type common determinant, "s" and the type specific determinants "r" and "y", present only on the specific serotypes.

Upon infection, a strong immunological response develops firstly against the type specific determinants and in a second time against the "s" determinant. Anti "s" antibodies are however recognised to be most effective in the neutralisation of the virus, protecting the patient from other infections and leading it to convalescence. The detection of HBsAb has become important for the follow up of patients infected by HBV and the monitoring of recipients upon vaccination with synthetic and natural HBsAg.

**C. PRINCIPLE OF THE TEST**

Microplates are coated with a preparation of highly purified HBsAg that in the first incubation with sample specifically captures anti HBsAg antibodies to the solid phase. After washing, captured antibodies are detected by an HBsAg labelled with peroxidase (HRP), that specifically binds the second available binding site of these antibodies. The enzyme specifically bound to wells, by acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of HBsAb in the sample and can be detected by an ELISA reader.

The amount of antibodies may be quantified by means of a standard curve calibrated against the WHO reference preparation. Samples are pre-treated in the well with an specimen diluent able to block interference present in vaccinated individuals.

**D. COMPONENTS**

Each kit contains sufficient reagents to perform 96 tests.


**1. Microplate: MICROPLATE**

96x12 microwell strips coated with purified heat-inactivated HBsAg of both subtypes (ad and ay) from human origin and sealed into a bag with desiccant.  
Allow the microplate to reach room temperature before opening; resal unused strips in the bag with desiccant and store at 4°C.

# HBsAb

## Enzyme Immunoassay for qualitative/quantitative determination of antibodies to Hepatitis B surface Antigen in human serum and plasma

For "in vitro" diagnostic use only



**DIAPRO**  
Diagnostici Bioprobes Srl  
Via G. Carducci n° 27  
20099 Sesto San Giovanni  
(Milano) - Italy  
Phone +39 02 27007161  
Fax +39 02 26007726  
e-mail: info@diapro.it

REF.SAB.CE  
96 Tests



2. Calibration Curve: **CAL-N<sup>o</sup>**

5x2.0 optical density to use and deliver coded standard curve derived from HBsAg positive plasma stored on WHO standard for anti-HBsAg<sup>14</sup> reference preparation 1877, lot 17-2-77) ranging: CAL-1 = 0.0 U/ml // CAL-2 = 10 mIU/ml // CAL-3 = 50 mIU/ml // CAL-4 = 100 mIU/ml // CAL-5 = 250 mIU/ml. Contains human serum proteins, 5% BSA, 10 mM phosphate buffer pH 7.4+/-0.1, 0.09% sodium azide and 0.1% Kathon GC as preservatives. Standards are blue coloured.

3. Wash buffer concentrate: **WASHBUF-20X**

1x60ml/kit, 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.1% Kathon GC.

4. Enzyme conjugate: **CONJ**

1x15.0 ml/vial Ready-to-use solution and red color coded. It contains inactivated purified HBsAg of both subtypes ad and ay, labelled with HR2, 5% BSA, 10 mM Tris buffer pH 6.8+/-0.1, 0.3 mg/ml gentamicin sulphate and 0.1% Kathon GC as preservatives.

5. Chromogen/Substrate: **SUBS-TMB**  
1x160ml/kit. Contains a 50 mM citrate-phosphate buffered solution at pH 3.5-3.3, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzene (TMB) and 0.02% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>). Note: To be stored protected from light as sensitive to strong illumination.

6. Sulphuric Acid: **H2SO4 0.3M**  
1x150ml/kit. Contains 0.3 M H2SO4 solution. Attention: Irritant (H314), H319, P280, P302+P352, P332+P313, P505+P531+P533, P537+P531, P662+P563).

7. Specimen Diluent: **DILSPF**  
1x60ml, 10 mM Tris buffered solution pH 7.4 +/-0.1, suggested to be used in the follow up of vaccination. It contains 0.09% sodium azide as preservative.

8. Control Serum: **CONTROL...m**

1 vial, Lyophilized.  
Contains fetal bovine serum proteins, human anti HBsAg antibodies, calibrated at 50 ± 10% WHO mIU/ml, 0.3 mg/ml gentamicin sulphate and 0.1% Kathon GC as preservatives.

9. Plate sealing foil n° 2

10. Package insert n° 1

E MATERIALS REQUIRED BUT NOT PROVIDED

- 1. Calibrated Micropipettes (100ul and 50ul) and disposable plastic-tips.
2. ELA grade water (double distilled or deionised, chemical disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA micropipette.
6. Well set at +37°C (+/-1°C tolerance).
7. Calibrated ELISA microplate reader with 450nm (reading) and with 620-630nm (blanking, strongly recommended) filters.
8. Calibrated ELISA microplate washer.
7. Vortex or similar mixing tools.

F WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.

2. All the personnel involved in performing the assay have to wear protective laboratory clothes, including gloves and shoes. The use of any sharp instrument (cutting blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents. When opening Kit vials and microtiter plates and when performing the test, protect the Chromogen (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cool room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

10. Do not use the kit after the expiration date stated on the external container and internal (vial) labels. A study conducted on an opened kit did not pointed out any relevant loss of activity up to six 6 uses of the device and up to 6 months.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.

12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste management during the use of the kit has to be discarded in accordance with national directives and laws concerning laboratory waste of chemical and biological substances. In general, liquid waste generated from the washing procedure from test kits of controls and from samples has to be treated as potentially infective material and inactivated before waste. Stipulated procedures of inactivation are treatment with a 10% final concentration of household bleach for 15-18 hrs or heat-inactivation by autoclave at 121°C for 20 min.

14. Accidental spills from samples and operators have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
16. Other waste material generated from the use of the kit (example: tips used for samples and controls, used micropipettes) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.

C SPECIMEN: PREPARATION AND WARNINGS

1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Samples have to be clearly identified with codes or numbers in order to avoid misinterpretation of results. Bar code labelling and electronic reading is strongly recommended.
3. Haemolysed (red) and visibly hyperlipemic (milky) samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bubbles should be discarded as they could give rise to false results.
4. Sera and plasma can be stored at +2...+8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be freeze/thawed more than once as this may generate particles that could affect the test result.
5. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.8u filters to clean up the sample for testing.
6. Samples whose anti-HBsAg antibody concentration is expected to be higher than 250 mIU/ml should be diluted before use either 1:10 or 1:100 in the Calibrator 0 mIU/ml. Dilutions have to be done in clean disposable tubes by diluting 50 ul of each specimen with 450 ul of Cal 0 (1:10). Then 50 ul of the 1:10 dilution are diluted with 450 ul of the Cal 0 (1:100). Mix tubes thoroughly on vortex when preparing the diluted samples.

H PREPARATION OF COMPONENTS AND WARNINGS

1. Microplate: Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned green, indicating a defect in conservation. In this case, call Dia.Pro's customer service.
2. Wash buffer concentrate: The wash buffer concentrate is supplied in a plastic bag and stored at +2...+8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

I Calibration Curve

3. Control Serum: Add the volume of ELISA grade water reported on the label to the lyophilised powder (it fully dissolve and then gently mix on vortex. The control after dissolution is not stable. Store frozen in aliquots at -20°C.
4. Wash buffer concentrate: The whole content of the concentrate solution has to be diluted 20X with distilled water and then used for the assay. During preparation, avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles. Note: Once diluted, the wash solution is stable for 7 week at +2...+8°C.

J Enzyme conjugate

5. Enzyme conjugate: Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidising chemicals, dust or microbes. If this component has to be transferred, use only plastic, and if possible, sterile disposable containers.
6. Specimen Diluent: Ready to use. Mix well on vortex before use.

K Chromogen/Substrate

7. Chromogen/Substrate: Ready to use. Mix well on vortex before use. Avoid contamination of the liquid with oxidising chemicals, agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

L Sulphuric Acid

8. Sulphuric Acid: Ready to use. Mix well on vortex before use.

Attention: Irritant (H314), H319, P280, P302, P352, P332+P313, P505+P531+P533, P537+P531, P662+P563).
Legend:
Warning H statements:
H315 - Causes skin irritation.
H319 - Causes serious eye irritation.

Precautionary P statements:
P280 - Wear protective gloves/protective clothing/eye protection/face protection.
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water.
P332 + P313 - If skin irritation occurs: Get medical advice/attention.
P505 + P531 + P533 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337 + P313 - If eye irritation persists: Get medical advice/attention.
P362 + P363 - Take off contaminated clothing and wash it before reuse.

I INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (70% ethanol, 10% solution of bleach, isopropyl alcohol) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a biasness of <2%.
2. The ELISA incubator has to be set at +37°C (tolerance of +/-1°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests. The ELISA washer is extremely important to the overall performance of the assay. The washer must be carefully validated and correctly optimized before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit control/calibrator and well-characterized negative and positive reference samples, and check to make sure they are reported below in the sections "Validation of Test Kits Assay Performance". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
3. The ELISA microplate reader has to be equipped with a reading filter of 450nm and with a second filter (620-630nm strongly recommended) for blanking purposes. Its standard performance should be (a) bandwidth <= 10 nm; (b) repeatability <= 1%; (c) linearity from 0 to >= 2.0; (d) absorbance range from 0 to >= 2.0. Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
4. When using an ELISA automated workstation, all critical steps (dispensation, incubation, washing, reading, soaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the



values reported in the sections 'Validation of Test' and 'Assay Performance'. The assay protocol has to be installed in the operating system of the protocol that is for the washer and the reader. In addition, the liquid handling part of the station (dispensation and mixing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated workstations is recommended for blood screening and when the number of samples to be tested exceeds 20-30 units per run.

7. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure full compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

**L. PRE ASSAY CONTROLS AND OPERATIONS**

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Chromogen/Substrate is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Solution as described above.
4. Dissolve the Control Serum as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
6. Set the ELISA incubator at +37°C and prepare the ELISA washer by printing with the diluted washing solution, according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipments are available and ready to use.

**In case of problems, do not proceed further with the test and advise the supervisor.**

**M. ASSAY PROCEDURE**

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

Two procedures can be carried out with the device according to the request of the clinician.

**M.1 Quantitative analysis**

1. Place the required number of strips in the microplate holder. Leave A1 and B1 wells empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2,8°C, sealed.
2. Then Dispense in all the wells to be used for the test, except for A1 and B1, 50µl of the Specimen Diluent.

**Important note:** This additive is added before distributing samples and control into specific wells and is particularly

intended for blocking some substances present in people undergoing vaccination and carbohydrate-antibodies.

2. Pipette 100µl in all the Calibrators, 100µl of Control Serum in duplicate and then 100µl of samples. The Control Serum is used to verify that the whole analytical system works as expected. Check that Calibrators, Control Serum and samples have been correctly added. Then incubate the microplate at +37°C for 60 min.

**Important note:** Strips have to be sealed with the adhesive sealing roll only when the test is performed manually. Do not cover stripes when using ELISA automatic instruments.

3. Wash the microplate as reported in section I.3.
4. In all the wells except A1 and B1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at +37°C for 60 minutes.

**Important note:**

- 1) Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.
- 2) Mix thoroughly the Enzyme Conjugate on vortex before use.
5. Wash the microplate as described.
5. Pipette 100µl TMB/H<sub>2</sub>O<sub>2</sub> mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at room temperature for 20 minutes.

**Important note:** Do not expose to strong direct light as a light background might be generated.

7. Stop the enzymatic reaction by pipette 100µl Sulphuric Acid into each well and using the same pipetting sequence as in step 6. Then measure the colour intensity with a microplate reader at 450nm (reading) and at 620-630nm (blanking, blanking strongly recommended), blanking the instrument on A1 and B1 wells.

**M.2 Qualitative analysis**

1. Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of blanking. Store the other strips into the bag in presence of the desiccant at 2,8°C, sealed.
2. Dispense 50 µl Specimen Diluent in all the wells, except for the blank A1. Then pipette 100µl of the Calibrator 0 mixture in duplicate, 100µl of the Calibrator 10 mixture in duplicate, 100µl of the Calibrator 250 mixture in duplicate, 100µl of the Control Serum in duplicate and samples have been correctly added. Then incubate the microplate at +37°C for 60 min.
3. Wash the microplate as reported in section I.3.
4. In all the wells except A1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at +37°C for 60 minutes.

**Important note:**

- 3) Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.
- 4) Mix thoroughly the Enzyme Conjugate on vortex before use.
5. Wash the microplate as described.
5. Pipette 100µl TMB/H<sub>2</sub>O<sub>2</sub> mixture in each well, the blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at room temperature for 20 minutes.

**Important note:** Do not expose to strong direct light as a light background might be generated.

7. Stop the enzymatic reaction by pipette 100µl Sulphuric Acid into each well and using the same pipetting sequence as in step 6. Then measure the colour intensity with a microplate reader at 450nm (reading) and at 620-630nm (blanking, blanking strongly recommended), blanking the instrument on A1 and B1 wells.

**Important general notes:**

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has to be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self-oxidation of the chromogen can occur leading to a higher background.
3. The Control Serum (CS) does not affect the cut-off calculation and therefore the test results calculation. The Control Serum may be used only when a laboratory internal quality control is required by the management.

**N. ASSAY SCHEME (standard procedure)**

Specimen Diluent	50 µl
Calibrators	100 µl
Control Serum	100 µl
Samples	100 µl
1. Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
Enzyme Conjugate	100 µl
2. Incubation	60 min
Temperature	+37°C
Wash step	4-5 cycles
TMB/H <sub>2</sub> O <sub>2</sub> mix	100 µl
3. Incubation	20 min
Temperature	RT
Sulphuric Acid	100 µl
Reading OD	450nm & 620nm

An example of dispensation scheme in qualitative assays is reported below:

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS
B	BLK	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS
C	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS	CS
D	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS	CS
E	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS	CS
F	CAL 0	CAL 10	CAL 250	CS	CS	CS	CS	CS	CS	CS	CS	CS
G	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
H	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
I	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
J	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
K	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
L	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS

An example of dispensation scheme in qualitative assays is reported below:

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
B	CAL 1	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
C	CAL 1	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
D	CAL 2	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
E	CAL 2	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
F	CAL 5	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
G	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
H	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
I	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
J	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
K	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS
L	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS	CS

Legend: BLK = Blank // CAL = Calibrators // CS = Sample

**O. INTERNAL QUALITY CONTROL**  
A validation check is carried out on the controls any time the kit is used in order to verify whether the performances of the assay are as qualified.  
Control that the following data are matched:


Parameters	Requirements
Blank well	< 0.100 OD450nm
Calibrator	< 0.200 OD450nm after blanking
Calibrator	OD450nm higher than the OD450nm of the Calibrator 0 mixture + 0.105
Control Serum	> 1.500 OD450nm
Control Serum	OD450nm = OD450nm CAL 50 mixture ± 10%
Coefficient of variation	< 30% for the Calibrator 0 mixture

If the results of the test match the requirements stated above, proceed to the next section.  
If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay
Calibrator 0 > 0.200	1. that the washing procedure and the washer settings are as validated in the pre-qualification study.
Coefficient of variation > 30%	2. that the proper washing solution has been used before use.
	3. that no mistle has been done in the assay procedure when the dispensation of standards is carried out.
	4. that no contamination of the CAL 0 mixture or of the wells where it was dispensed has occurred due to positive samples, 10 spills or 10 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate
	6. that the washer needles are not blocked or partially obstructed

Calibrator 16 mIU/ml OD450nm < Cal 0 + 0.100	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the calibration (ex.: dispersion of a wrong calibrator); 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the standard has occurred.
Calibrator 250 mIU/ml < 1500 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the standard has occurred.
Control Serum Diluted from expected value	1. the procedure has been correctly performed; 2. no mistake has occurred during its distribution (ex.: dispersion of a wrong sample); 3. the washing procedure and the washer settings are correct; 4. no external contamination of the standard has occurred; 5. the Control Serum has been dissolved with the right volume reported on the label; 6. the control serum, poured out, the assay has to be repeated after stimulating the reason of this error; If no mistake has been found, proceed as follows: a) a value up to +/20% is obtained: the overall precision of the laboratory might not enable the report of the mean expected value +/10% acceptance or rebal of this result; b) a value higher than +/20% is obtained: in this case the test is invalid and the Diabro customer service has to be called.

Example of Calibration Curve :



**Important Note:**  
Do not use the calibration curve above to make calculations.

**P 2 Qualitative method**  
In the qualitative method, calculate the mean OD450nm values for the Calibrators D and 10 mIU/ml and then check that the assay is valid.

**Example of calculation:**  
The following data must not be used instead of real figures obtained by the user:  
Calibrator 0 mIU/ml: 0.020 - 0.024 OD450nm  
Mean Value: 0.022 OD450nm  
Lower than 0.200 - Accepted  
Calibrator 10 mIU/ml: 0.250 - 0.270 OD450nm  
Mean Value: 0.260 OD450nm  
Higher than Cal 0 + 0.100 - Accepted  
Calibrator 250 mIU/ml: 2.845 OD450nm  
Higher than 1.500 - Accepted

**P. RESULTS**  
**P 1 Quantitative method**  
If the test turns out to be valid, use for the quantitative method an approved curve fitting program to draw the calibration curve from the values obtained by reading at 450nm (4-parameters interpolation is suggested).  
Then on the calibration curve calculate the concentration of anti HbsAg antibody in samples.  
An example of Calibration curve is reported in the next page.

**Q. INTERPRETATION OF RESULTS**  
Samples with a concentration lower than 10 WHO mIU/ml are considered negative for anti HbsAg antibody by most of the international medical literature.  
Samples with a concentration higher than 10 WHO mIU/ml are considered positive for anti HbsAg antibody.  
In the follow up of vaccination recipients, however, the value of 20 WHO mIU/ml is usually accepted by the medical literature as the minimum concentration at which the patient is considered clinically protected against HBV infection.

- Important notes:**
1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
  2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
  3. Diagnosis has to be done and released to the patient by a suitably qualified medical doctor.

**R. PERFORMANCE**

Evaluation of Performance has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

**1. LIMIT OF DETECTION:**  
The limit of detection of the assay has been calculated by means of the HBSAB International preparation supplied by CLB on behalf of WHO (1<sup>st</sup> reference preparation 1977, lot 17-2-77), on which Calibration Curve has been calibrated. HBV negative serum was used as diluent, as recommended by the supplier. Results of Quality Control are given in the following table:

WHO mIU/ml	SAB CE Lot # 1002	SAB CE Lot # 1001	SAB CE Lot # 1002/2
50	0.933	0.812	0.846
10	0.219	0.192	0.194
5	0.110	0.096	0.104
2.5	0.057	0.058	0.057
Std 0	0.021	0.015	0.023

**2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY**  
A Performance Evaluation has been conducted on a total number of more than 700 samples.

**2.1 Diagnostic Specificity**  
It is defined as the probability of scoring negative in the absence of specific analyte.  
More than 500 negative specimens were tested, internally and externally, against a European company.  
A diagnostic specificity of 98.8% was assessed.  
Moreover, diagnostic specificity was assessed by testing 113 potentially interfering specimens (other infectious diseases, pregnant women, hemodialyzed, etc.) against the European company. A value of specificity of 100% was assessed.  
Finally, both human plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and human sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

**2.2 Diagnostic Sensitivity**  
It is defined as the probability of the assay of scoring positive in the presence of specific analyte.  
106 vaccinated patients were evaluated providing a diagnostic sensitivity of 100%.  
More than 100 HBV naturally infected patients were tested, internally and externally, against the European company: a diagnostic sensitivity of 100% was found.

**3. PRECISION:**  
The mean values obtained from a study conducted on three samples of different anti-HbsAg reactivity, examined in 16 replicates in three separate runs is reported below:  
SAB CE lot # 1202

Calibrator 250 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	2.998	3.000	3.259	3.085
Std Deviation	0.152	0.151	0.156	0.153
CV %	5.1	5.0	4.8	5.2

SAB CE lot # 1002

Calibrator 6 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	0.048	0.148	0.050	0.149
Std Deviation	0.005	0.004	0.006	0.006
CV %	9.4	8.4	11.5	9.8

Calibrator 16 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	0.249	0.252	0.242	0.248
Std Deviation	0.021	0.020	0.023	0.021
CV %	8.3	7.9	9.8	8.6

Calibrator 250 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	3.544	3.153	3.812	3.093
Std Deviation	0.160	0.176	0.198	0.159
CV %	4.3	4.8	3.9	4.3

SAB CE lot # 1002/2

Calibrator 6 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	0.050	0.051	0.050	0.050
Std Deviation	0.005	0.006	0.006	0.005
CV %	10.0	10.9	11.9	10.9

Calibrator 10 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	0.226	0.236	0.239	0.234
Std Deviation	0.026	0.017	0.018	0.018
CV %	0.9	7.0	7.5	7.0

Calibrator 250 mIU/ml (N = 16)

Mean values	1 <sup>st</sup> run	2 <sup>nd</sup> run	3 <sup>rd</sup> run	Average value
OD 450nm	3.526	3.457	3.489	3.494
Std Deviation	0.137	0.143	0.162	0.147
CV %	3.9	4.1	4.6	4.2

The variability shown in the tables did not result in sample misclassification.

**4. ACCURACY**  
The assay accuracy has been checked by the dilution and recovery tests. Any "hook effect", underestimation (likely to happen at high doses of analyte, was ruled out up to 10,000 mIU/ml).

**S. LIMITATIONS OF THE PROCEDURE**  
Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.  
This test is suitable only for testing single samples and not pooled ones.  
Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, symptomatology, as well as other diagnostic data should be considered.



REFERENCES

1. Engvall E, et al., *Immunochemistry*, 8: 671-674, 1971.
2. Engvall E, et al., *Immunochem*, 109, 129-135, 1971.
3. Remington J.S. and Klein J.O. In "Infectious diseases of the fetus and newborn infant". Sanders, Philadelphia, London, Toronto.
4. Volkov A. In "Essential of Medical Microbiology", 2<sup>nd</sup> ed, p. 729. G.B. Lippincott Company, Philadelphia, New York, S. Jose, Toronto.
5. Sjögreen D.F. et al., *Am.J.Pathol.*, 83: 838, 1975.
6. Bisher L.F., Dodd R.J., Sandler S.G., In "Viral Hepatitis: Lab. and Clinical Science". F. Denhardt, J. Denhardt eds. M. Decker Inc., New York, 215-230, 1983.
7. Costantini V., *Br.Med.J.*, 28: 156, 1972.
8. Laner J.J. et al., *Immunochem*, 108: 1066, 1971.
9. Muschler K. et al., *Am.J.Clin.Pathol.*, 76: 773, 1981.
10. Howard C.R., *Immuno Today*, 5: 185, 1984.
11. Ajayi R.D., *Lancet* 7874, 190-193, 1974.
12. Jilg W., et al., *J.Hepatol*, 9: 201-207, 1988.
13. P. Crovari et al., *Soll. Ist. Sieroter. Milan*, 63: 14-18, 1984.
14. M. Davidson et al., *J.Natl.Cancer Inst.*, 59: 1451-1467, 1977.
15. F. Gyorkov et al., *J.Natl.Cancer Inst.*, 59: 1451-1467, 1977.
16. S. Hadler et al., *N.E.J.Med.*, 313: 209-214, 1986.
17. J.H. Hooper et al., *Hepatology*, 7: 758-763, 1987.
18. C.L. Howard, *J.Gen.Viro.*, 67: 1215-1225.
19. W. Jilg et al., *J.Hepatol*, 6: 201-207, 1988.
20. P. Michel et al., *Hepatology*, 7: 1144-117, 1986.
21. W. Sarnaus et al., *N.E.J.Med.*, 303: 833-836, 1980.
22. P. Toulats et al., *Nature*, 317: 468-469, 1985.
23. A.J. Zuckerman et al., In "Hepatitis Viruses of Man". Academic Press, London, 1979.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:  
Dia, Pro Diagnostico Bioprobes Srl  
Via G. Carducci n° 27 - Sesto San Giovanni (MI) - Italy



0318

**HCV IGM**

*Important Note: Even if plasma has been chemically inactivated, handle this component as potentially infectious.*

# HCV IGM

**Enzyme ImmunoAssay (ELISA) for the quantitative/qualitative determination of IGM antibodies to Hepatitis C Virus in human serum and plasma**

- for "in vitro" diagnostic use only -



**DIA PRO**  
**Diagnostica Bioprobes S.p.A.**  
 Via G. Carducci n° 27  
 20099 Sesto San Giovanni  
 (Milano) - Italy  
 Phone +39 02 27007161  
 Fax +39 02 26007226  
 e-mail: [info@diapro.it](mailto:info@diapro.it)

REF. CVM/CE  
 96 Tests

**A. INTENDED USE**  
 Enzyme Immunoassay (ELISA) for the quantitative/qualitative determination of IGM antibodies to Hepatitis C Virus in human plasma and sera. The kit is mainly intended for the follow-up of HCV chronic patients submitted to anti-viral pharmaceutical treatment. For "in vitro" diagnostic use only.

**B. INTRODUCTION**  
 Antiviral drugs, such as Interferon taken alone or in combination with Ribavirin, can be used for the treatment of persons with chronic viral hepatitis C.  
 Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with Ribavirin is effective in about 30% to 50% of patients; Ribavirin does not appear to be effective when used alone.  
 Active production of HCV antigens in the liver of chronic patients generates spikes of IGM antibodies production and release of liver specific enzymes, similar to what happens in HBV chronic patients. The presence of anti HCV IGM is usually correlated to a phase of suffering and cellular damage of the liver.  
 During the pharmaceutical treatment HCV IGM may represent a marker for the follow-up of the efficacy of the drug itself, monitoring the balance between its effectiveness and the side effects, that often may be heavy for the patient.

**C. PRINCIPLE OF THE TEST**  
 Microplates are coated with HCV immunodominant synthetic peptides (core peptide, recombinant NS3, NS4 and NS5 peptides).  
 In the 1<sup>st</sup> incubation, the solid phase is treated with diluted samples and anti HCV IGM are captured, if present, by the antigens. After washing out all the other components of the sample, in the 2<sup>nd</sup> incubation bound anti-HCV IGM are detected by the addition of anti IGM antibody, labeled with peroxidase (HRP). The enzyme captured on the solid phase, acting on the substrate/chromogen mixture, generates an optical signal that is proportional to the amount of anti-HCV IGM antibodies present in the sample.  
 The presence of IGM in the sample may therefore be quantitated by means of a calibration curve able to determine the content of the antibody in anti-U/ml.  
 Neutralization of IgG anti-HCV, carried out directly in the well is performed in the assay in order to block interferences due to this class of antibodies in the determination of IGM.

**D. COMPONENTS**  
 Each kit contains sufficient reagents to perform 96 tests.

**1. Microplate: MICROELATE**  
 12 strips x 8 microtiter wells with HCV-specific synthetic antigens (core, NS4 and NS5 peptides and recombinant NS3). Plates are sealed into a bag with desiccant.

**2. Calibration Curve: CAL N° 1**  
 6x2.0 ml/vial. Ready to use end color coded standard curve calibrated on an Internal Gold Standard (in absence of a defined international one) for NS3 ranging 10-800 IU/ml.  
 CAL 1 = 0 antiU/ml CAL 4 = 80 antiU/ml  
 CAL 3 = 25 antiU/ml CAL 5 = 250 antiU/ml  
 CAL 5 = 100 antiU/ml CAL 6 = 250 antiU/ml  
 It contains chemical inactivated HCV IGM positive human plasma, 100 mM Tris buffer pH 7.4+/-0.1, 0.2% Tween 20, 0.09% sodium azide and 0.1% Kathon GC as preservatives.  
 The Calibration Curve is coded with blue alphanumeric type.

**3. Wash buffer concentrate: WASHBUFF 20X**  
 1x60ml/vial 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0+/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.

**4. Enzyme conjugate: ENZ**  
 1x15ml/vial. Ready to use and red color coded. It contains horseradish peroxidase conjugated polyclonal antibodies to human IGM, 5% BSA, 10 mM Tris buffer pH 5.8+/-0.1, 0.1% Kathon GC and 0.02% gentamicin sulphate as preservatives.

**5. Chromogen/substrate: SUBS IMB**  
 1x15ml/vial. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine (or TMB) and 0.02% hydrogen peroxide (or H<sub>2</sub>O<sub>2</sub>).  
**Note: To be stored protected from light as sensitive to strong illumination.**

**6. Sulphuric Acid: H2SO4 0.3 M**  
 1x15ml/vial. It contains 0.3 M H2SO4 solution.  
 Attention: Irritant (H315, H319, P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P362+P353).

**7. Specimen Diluent: DILSPE**  
 2x60ml/vial. It contains 2% casein, 10 mM Na-citrate buffer pH 6.0+/-0.1, 0.2% Tween 20, 0.09% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

**8. Neutralizing Reagent: SOLINEUT**  
 1x60ml/vial. It contains 90% anti-HCG 2% casein, 10 mM Na-citrate buffer pH 6.0+/-0.1, 0.09% Na-azide and 0.1% Kathon GC as preservatives.

**9. Plate sealing foils n°2**

**10. Package insert n°1**

**E. MATERIALS REQUIRED BUT NOT PROVIDED**

1. Calibrated Micropipettes (1000, 100 and 10ul) and disposable plastic tips.
2. EIA grade water (distilled or deionised charcoal treated) to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (+/-0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-650nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

**F. WARNINGS AND PRECAUTIONS**

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
2. All the personnel involved in performing the assay have to wear protective laboratory clothes, lab-coat, gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in bioassay procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Bioassay in Microbiological and Biomedical Laboratories", ed. 1984.
3. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.



4. The laboratory environment should be controlled so as to avoid contamination such as dust or air-borne microbial agents when opening kit vials and microplates and when performing the test. Protect the ChromogenSubstrate (or TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
6. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not also be interchanged.
7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each component.
10. Do not use the kit after the expiration date stated on the external container and internal (vials) labels. A study conducted on one kit did not pointed out any relevant loss of activity.
11. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. In compliance with what is reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed 1989.
12. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
13. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residues of concentrated materials and from samples that to be treated as potentially infectious material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 15-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
14. Accidental spills from samples and operations have to be absorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
15. The Sulphuric Acid is an irritant. In case of spills, wash the surface with plenty of water.
16. Other waste materials generated from the use of the kit (example: tips used for samples and controls/calibrators, used microplates) should be handled as potentially infectious and disposed according to national directives and laws concerning laboratory wastes.

**H. PREPARATION OF COMPONENTS AND WARNINGS**

- Microplate:**  
Allow the microplate to reach room temperature (about 1 h) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of drying. In this case call Dia.Pro's customer service. Unused strips have to be placed back into the aluminium pouch. In presence of desiccant supplied, firmly zippled and stored at -27,-8°C. When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.
- Calibration Curve**  
Ready to use components. Mix carefully on vortex before use.
- Wash buffer concentrate:**  
The whole content of the concentrated solution has to be diluted 20x with bidistilled water and mixed gently end-over-end before use. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles. **Note: Once diluted, the wash solution is stable for 1 week at +2,-8°C.**
- Enzyme conjugate:**  
Ready to use. Mix well on vortex before use. Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component has to be transferred use only plastic, possibly sterile disposable containers.
- ChromogenSubstrate:**  
Ready to use. Mix well on vortex before use. Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and other reagents. If this component has to be transferred use only plastic, possibly sterile disposable container.
- Sample Diluent**  
Ready to use component. Mix carefully on vortex before use.
- Neutralizing Reagent**  
Ready to use component. Mix carefully on vortex before use.
- Sulphuric Acid:**  
Ready to use. Mix well on vortex before use. Attention: Irritant (H314, H319), 2001, P02+P302, P332+P313, P305+P351+P338, P307+P311, P308+P313.
- Legends:**  
Warning H statements:  
H315 – Causes skin irritation.  
H319 – Causes serious eye irritation.  
H335 – Causes respiratory irritation.
- Precautionary P statements:**  
P201 – Wear protective gloves/protective clothing/eye protection/face protection.  
P302 + P332 – IF ON SKIN: Wash with plenty of soap and water.  
P332 + P313 – If skin irritation occurs: Get medical advice/attention.

- P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P331 + P313 – If eye irritation persists: Get medical advice/attention.  
P332 + P313 – Take off contaminated clothing and wash it before reuse.
- 1 INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT**
1. Microplates have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained in order to show a precision of 1% and a thickness of +1-2%. Decontamination of spills or residues of kit components should also be carried out regularly.
  2. The ELISA incubator has to be set at +37°C (tolerance of +4,0-5°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
  3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimised using the kit controls/calibrators and reference panels before using the kit for routine laboratory tests. Usually 4-5 washing cycles (aspiration + dispersion of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested in order to set correctly their number. It is recommended to run an assay with the kit controls/calibrators and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (detergents and cleaning of needles) of the washer have to be carried out according to the instructions of the manufacturer.
  4. The ELISA microplate reader has to be equipped with a reading tier of 450nm and with a second filter (620-650nm, slightly recommended) for blanking purposes. Its standard absorbance should be (a) bandwidth ≤ 10 nm, (b) repeatability 5 ng/ml for 10, 2, 20, (c) linearity to ≥ 2,0, (d) linearity to 100%. Blanking is carried out on the optical system of the section "Assay Procedure". The optical density of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly calibrated according to the manufacturer's instructions.
  5. When using an ELISA automated workstation, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, checked, controlled and regularly serviced in order to avoid any error in the section "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the station in addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. The attention must be paid to avoid carry over by the slides used for dispensing and for washing. This must be the slides and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceeds 20-30 units per run.
  6. Dia.Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

**1. PRE ASSAY CONTROLS AND OPERATIONS**

1. Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
  2. Check that the liquid components are not contaminated by visible particles or aggregates.
  3. Check that the ChromogenSubstrate is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette.
  4. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
  5. Divide all the content of the 20x concentrated Wash Solution as described above.
  6. Allow all the other components to reach room temperature (about 1 h) and then mix gently on vortex all liquid reagents.
  7. Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution, according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for the use with the kit.
  8. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
  9. If using an automated work station, turn on, check settings and set up to use the right assay protocol.
  10. Check that the microplates are set to the required volume.
  11. Check that all the other equipment is available and ready to use.
  12. In case of problems, do not proceed further with the test and advise the supervisor.
- ASSAY PROCEDURE**  
The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.
- Two methods of analysis are possible, as described below:
- M1 QUANTITATIVE ASSAY**
1. Place the required number of strips in the plastic holder and carefully identify the wells for calibrators and samples.
  2. Dilute samples 1:101 dispensing 1 ml Sample Diluent into a disposable tube and then 10 µl sample, mix on vortex before use. Do not dilute the Calibrators as they are ready-to-use.
  3. Leave the A1+B1 wells empty for blanking purposes.
  4. Dispense 50 µl Neutralizing Reagent in all the wells, except A1+B1 wells, used for blanking operations and the wells used for the Calibration Curve.
  5. In the identified positions pipette 100 µl of the Calibrators in duplicate followed by 100 µl of diluted samples. Check that Calibrators and samples have been correctly added.
  6. Incubate the microplate for 60 min at +37°C.
- Important note:** Strips have to be sealed with the adhesive sealing film, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.
7. When the first incubation is finished wash the microplate as previously described (section 13)
  8. In all the wells, except A1+B1, pipette 100 µl ENZYME Conjugate. Incubate the microplate for 60 min at +37°C.
- Important note:** Be careful not to touch the plastic inner surface of the well with the tip held with the Enzyme Conjugate. Contamination might occur.



9. When the second incubation is finished, wash the microwells as previously described (section 1.3).  
 10. Pipette 100 µl Chromogen/Substrate into all the wells. A1-S1 included.

**Important note:** Do not expose to strong direct light, as a high background might be generated.

11. Incubate the microplate protected from light at room temperature (18-24°C) for 20 minutes. Wells dispensed with positive samples and with positive calibrators will turn from clear to blue.  
 12. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10 to block the positive reaction. Addition of the stop solution will turn the yellow calibrators and the positive samples from blue to yellow.  
 13. Measure the color intensity of the solution in each well as described in section 1.5 using a 450nm filter (reading) and a 620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1 or B1 or both.

**M.2 QUALITATIVE ASSAY**

- Place the received number of strips in the plastic holder and carefully identify the wells for calibrators and samples.
- Divide samples 1:101 diluting 1 ml Sample Diluent into a disposable tube and then 10 µl sample mix on vortex before use. Do not dilute the Calibrators as they are ready-to-use.
- Leave the A1 well empty for blanking purposes.
- Dispense 50 µl Neutralizing Reagent in all the wells except the Calibrators.
- Then pipette 100 µl of Calibrator 0 and/or 1 in duplicate, 100-µl of Calibrator 10 and/or 100 in duplicate and finally 100 µl of diluted samples. Check that Calibrators and samples have been correctly added.
- Incubate the microplate for 60 min at +37°C.

**Important note:** Strips have to be sealed with the adhesive sealing foil, only when the test is performed manually. Do not cover strips when using ELISA automatic instruments.

7. When the first incubation is finished, wash the microwells as previously described (section 1.3).  
 8. In all the wells, except A1, pipette 100 µl Enzyme Conjugate. Incubate the microplate for 60 min at +37°C.

**Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

9. When the second incubation is finished, wash the microwells as previously described (section 1.3).  
 10. Pipette then 100 µl Chromogen/Substrate into all the wells. A1 included.

**Important note:** Do not expose to strong direct light, as a high background might be generated.

11. Incubate the microplate protected from light at room temperature (18-24°C) for 20 minutes. Wells dispensed with positive samples and with positive calibrators will turn from clear to blue.  
 12. Pipette 100 µl Sulphuric Acid into all the wells using the same pipetting sequence as in step 10 to block the enzymatic reaction. Addition of the stop solution will turn the yellow calibrators and the positive samples from blue to yellow.  
 13. Measure the color intensity of the solution in each well as described in section 1.5 using a 450nm filter (reading) and a

620-630nm filter (background subtraction, strongly recommended), blanking the instrument on A1.

**General important notes:**

- If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self-oxidation of the TMB Chromogen can occur leading to high background.

**N. ASSAY SCHEME**

Method	Operations
Neutralizing Reagent Calibrators (no SOLVENT!)	50 µl 100 µl
Samples diluted 1:101	100 µl
1 <sup>st</sup> incubation	60 min +37°C
Wash step	4-5 cycles
Enzyme conjugate	100 µl
2 <sup>nd</sup> incubation	60 min +37°C
Wash step	4-5 cycles
TMB/H <sub>2</sub> O <sub>2</sub>	100 µl
3 <sup>rd</sup> incubation	20 min +37°C
Temperature	18-24°C
Sulphuric Acid	100 µl
Reading OD	450nm

An example of dispensation scheme in qualitative assays is reported below:

**Microplate**

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	CAL4	S3									
B	BLK	CAL4	S4									
C	CAL1	CAL5	S5									
D	CAL1	CAL5	S6									
E	CAL2	CAL6	S7									
F	CAL2	CAL6	S8									
G	CAL3	S1	S9									
H	CAL3	S2	S10									

Legend: BLK = Blank // CAL = Calibrators // S = Sample

An example of dispensation scheme in qualitative assays is reported below:

**Microplate**

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S4										
B	CAL1	S5										
C	CAL1	S6										
D	CAL2	S7										
E	CAL2	S8										
F	S1	S9										
G	S2	S10										
H	S3	S11										

Legend: BLK = Blank // CAL = Calibrators // CS = Control Serum // S = Sample

**Q. INTERNAL QUALITY CONTROL**

A validation check is carried out any time the kit is used in order to verify whether the performance of the assay are as qualified. Control that the following data are met:

Parameter	Requirements
Blank well	< 0.100 OD450nm
Calibrator 0 and/or 1	> 0.200 OD450nm after blanking
Calibrator 10 and/or 100	OD450nm > OD450nm CAL 0 and/or 1 + 0.100
Calibrator 250 and/or 1000	3.500 > OD450nm > 2.000

If the results of the test match the requirements stated above, proceed to the next section.  
 If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay.
Calibrator 0 and/or 1 < 0.200 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study. 2. that the proper washing solution has been before use. 3. that no mistake has been made when pipetting the reagents. 4. that no contamination of the Cal 0 and/or 1 of the wash, where this was dispensed, has occurred due to positive samples, to spots or to the strip. 5. that micropipettes have not become contaminated with positive samples or with the reagents. 6. that the washer needles are not blocked or partially obstructed.
Calibrator 10 and/or 100 < CAL 0 + 0.100	1. that the procedure has been correctly performed. 2. that no mistake has occurred during its distribution. 3. that the washing procedure and the washer settings are as validated in the pre qualification study. 4. that no external contamination of the calibrator has occurred. 5. that the procedure has been correctly performed. 6. that no mistake has occurred during the distribution of the calibrator. 7. that the washing procedure and the washer settings are as validated in the pre qualification study.
Calibrator 250 and/or 1000 > 3.500 OD450nm after blanking	1. that the external contamination of the calibrator has occurred during the assay. 2. that the washing procedure and the washer settings are as validated in the pre qualification study. 3. that the proper washing solution has been before use. 4. that no mistake has been done in the assay. 5. that no contamination of the Cal 250 and/or 1000 of the wells where this was dispensed, has occurred due to positive samples, to spots or to the enzymatic conjugate. 6. that micropipettes have not become contaminated with positive samples or with the reagents. 7. that the washer needles are not blocked or partially obstructed.

If any of the above problems has occurred, report the problem to the supervisor for further actions.

**P. RESULTS**

If the test runs out to be valid, interpretation of results is carried out in the quantitative assay from the mean OD450nm value of the Calibration Curve elaborated with an appropriate curve fitting system (suggested 4 parameters).

In the qualitative assay interpretation of results is done on the mean OD450nm value of the Calibrator 10 and/or 100 (or CAL 2), by means of the following formulation:

$$\text{Mean OD450nm CAL 2} = \text{cut-off} (Co)$$

**Important note:** When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to generate the correct interpretation of results.

**Q. INTERPRETATION OF RESULTS**

**Q.1 QUANTITATIVE ASSAY**

Concentrations in and/or 1 are obtained elaborating OD450nm of samples on the fitted calibration curve.  
 The concentration of IgG is from Literature correlated proportionally with the live damage produced by antibodies to HCV upon virus replication in hepatocytes.  
 A decrease in IgG concentration upon pharmacological treatment is usually difficultly acknowledged as a sign of recovery and therapeutic efficacy.

**Q.2 QUALITATIVE ASSAY**

Test results are interpreted as a ratio of the sample OD450nm value (S) and the cut-off value (Co), or S/Co, according to the following table:

S/Co	Interpretation
< 1.0	Negative
> 1.0	Positive

A negative result indicates that the patient has not developed IgM antibodies to HCV.  
 A positive result is indicative of an ongoing HCV active infection.

**Important notes:**

- Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
- When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
- Diagnosis has to be done and released to the patient by a suitably qualified medical doctor.
- The results of this ELISA assay should be anyway implemented with other diagnostic and clinical tests.

An example of calculation is reported below.

The following data must not be used instead or real figures obtained by the user:

- CAL 1: 0.060 - 0.080 OD450nm
- Mean Value: 0.070 OD450nm
- Lower than 0.200 = Accepted
- CAL 2: 0.200 - 0.220 OD450nm
- Mean Value: 0.210 OD450nm
- Higher than CAL1+0.100 = accepted
- Cut-Off or Co = 0.210

Sample 1:	0.080 OD450nm
Sample 2:	1.860 OD450nm
Sample 1	S/Cc < 1.0 = negative
Sample 2	S/Cc > 1.0 = positive

## R PERFORMANCE CHARACTERISTICS

Evaluation of Performance has been conducted on selected panels carried out in a clinical external center and internally.

### 1. Limit of detection

No international standard for HCV IgM Antibody detection has been defined so far by the European Community. In its absence an Internal Std Standard (or ICS), derived from a patient with an history of chronic HCV infection, has been defined in order to provide the device with a constant and excellent sensitivity.

### 2. Diagnostic Sensitivity and Specificity:

The diagnostic performances were evaluated in a study conducted in an external clinical center, with excellent experience in the diagnosis of infectious diseases and HCV. The Diagnostic Sensitivity was studied on about 200 samples, pre-tested positive with an analytical system developed in house by the clinical laboratory where the study was conducted. Positive samples were collected from patients with a clinical history of HCV infection (acute and chronic). In addition some Serono assays (Panab, purchased from Boston Biomedical Inc, USA, were examined). The diagnostic specificity was determined on panels of more than 300 negative samples from individuals and blood donors, classified negative for anti HCV antibodies with the reference kit in use in the laboratory, including potentially interfering specimens.

A panel of potentially interfering samples (RF+, hemolysed, lipemic, etc.) was also examined. No interference was observed on the samples examined.

Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed. Frozen specimens have also been tested to check whether samples freezing interferes with the performance of the test. No interference was observed on clean and particle free samples.

The Performance Evaluation provided the following values :

<b>Sensitivity</b>	<b>&gt; 98 %</b>
<b>Specificity</b>	<b>&gt; 98 %</b>

### 3. Reproducibility:

It has been calculated on two samples examined in replicates in different runs. Results are reported below summarized in a table.

Average values	Calibrator 2	Calibrators 5
N = 48	10 arbU/ml	100 arbU/ml
OD450nm	0.241	1.632
Std Deviation	0.027	0.113
CV %	11.3	6.9

## S. LIMITATIONS

False positivity has been assessed on less than 2% of the normal population, mostly due to high titers of RF.  
Frozen samples containing fibrin particles or aggregates may generate false positive results.

## BIBLIOGRAPHY

- Kisavtsev-El, Zhanovnok-SV, Mitcra VM, Demiche AP, Zh Mikrobiol Epidemiol Immunobiol. 2006 Mar-Apr;(2):57-61
- Papathodoridis GV, Daliadisima JK, Katsoulidou A, Syras V, Albrecht M, Michel G, Halzaks A, Tassopoulos NC. J Hepatol. 1997 Jul;(27):136-41.
- Martiniell AL, Brown D, Braun HB, Michel G, Dushenko GM. J Hepatol. 1996 Jan;(24):121-6.
- Pawlotsky JM, Darbyth F, Remire J, Pellet C, Udén L, Stuyver L, Roudot-Thoraval F, Duvoix C, Douvin C, Mallat A, et al. J Med Virol. 1995 Nov;(47):285-91.
- Stransky J, Honzakova E, Vardasova J, Horšova M, Krcal J, Nerneck V, Horak J. Acta Virol. 1995 Apr;(40):2151-5.
- Nikolaeva LI, Blokhina NP, Tsukova NN, Voronkova NV, Mironovskii MI, Braginskij DM, Yasrebova ON, Boyutskaya OB, Isaeva OV, Michailov MI, Archakov AI. Gid. 2000 Nov;(47):698-702.
- Bizalion T, Ahmed SN, Guichard S, Chevalier P, Adham M, Ducer C, Bauleux J, Trepo C. J Med Virol. 1998 Nov;(56):224-9.
- Tan A, Yang G, Dreyfus G, Rouquie P, Durant J, Rappal A, Rampal P, Benzaken S, Am J Gastroenterol. 1997 Oct;(92):1835-8.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:  
Dia Pro Diagnostic Bioprodes Srl  
Via G. Carducci n° 27 – Sesto San Giovanni (MI) – Italy



0318



## HCV Ab

**A. INTENDED USE**  
Version 4.0 Enzyme immunoassay (EUSA) for the determination of antibodies to Hepatitis C Virus in human plasma and sera. The kit is intended for the screening of blood units and the follow-up of HCV-infected patients.  
For in vitro diagnostic use only.

# HCV Ab

## Version 4.0 Enzyme Immunoassay for the determination of anti Hepatitis C Virus antibody in human serum and plasma

- for "in vitro" diagnostic use only -



**DIA.PRO**  
**Diagnostic-Bioprobes Srl**  
Via G. Carducci n° 27  
20099 Sesto San Giovanni  
(Milano) - Italy  
Phone +39 02 27007161  
Fax +39 02 26007726  
e-mail: [info@dialpro.it](mailto:info@dialpro.it)

REF: CVAB/CE  
96.192.480.960 TS85

**B. INTRODUCTION**  
The World Health Organization (WHO) define Hepatitis C

Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post-transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. The major causes of HCV infection worldwide are use of unsterilized blood transfusions and reuse of needles or syringes that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is not costly for most persons in developing countries to afford. Thus, from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviors (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the *Flaviviridae* family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease. An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high prognostic (80%) of inducing chronic infection. HCV is classified into several distinct genotypes, which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Cirrhosis develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 2% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer who do not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with

alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, though the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioral practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils, or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or reuse unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (where prevalence data are available) is high compared to other countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme immunoassay assays (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of active infection and a potential for spread of the infection and/or the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% in 24hrs of patients. Interferon combined with ribavirin is effective in about 50% to 50% of patients. Ribavirin doses and regimen in the future will be studied.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any productive immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HIV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation and testing of blood products; (c) implementation and maintenance of infection control practices in health care settings, including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce reverse of infections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices.

The genome encodes for structural components: a nucleocapsid protein and two envelope glycoproteins, and functional components involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

**C. PRINCIPLE OF THE TEST**

Microplates are coated with HCV-specific antigens derived from core and the regions encoding for conservative and immunodominant antigenic determinants (Core peptide, recombinant NS3, NS4 and NS5 peptides). (Core peptide, NS3, NS4 and NS5 peptides). Plates are sealed into a bag with desiccant.

12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed into a bag with desiccant.

1. Microplate **MICROPLATE**

Code CVAB CE contains reagents for 192 tests.

**D. COMPONENTS**

- 1. Microplate **MICROPLATE**
- 2. Positive Control **CONTROL**
- 3. Negative Control **CONTROL**
- 4. Calibrator **CAL**

Code	Number of tests	CVAB CE 56	CVAB CE 480	CVAB CE 960
1. Microplate	n° 1	96	480	960
2. Positive Control	n° 1	96	480	960
3. Negative Control	n° 1	96	480	960
4. Calibrator	n° 1	96	480	960

1. Microplate: Ready to use control. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 8.0 +/-0.1, 0.5% Tween 20, 0.05% Na-azide and 0.1% Kathon GC as preservatives. The negative control is olive green colour coded.

2. Positive Control **CONTROL**: It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 8.0 +/-0.1, 0.5% Tween 20, 0.05% Na-azide and 0.1% Kathon GC as preservatives. The Positive Control is blue colour coded.

3. Negative Control **CONTROL**: It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 8.0 +/-0.1, 0.5% Tween 20, 0.05% Na-azide and 0.1% Kathon GC as preservatives. The Negative Control is blue colour coded.

4. Calibrator **CAL**: Ready to use. It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 8.0 +/-0.1, 0.5% Tween 20, 0.05% Na-azide and 0.1% Kathon GC as preservatives.

Note: The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.

- 5. Wash buffer concentrate **WASH-BUF 20X**
- 6. Enzyme Conjugate **CONJ**
- 7. Chromogen/substrate **SUBS/TMB**
- 8. Assay Diluent **DILAS**
- 9. Sulphuric Acid **H2SO4 0.3 M**
- 10. Sample Diluent **DILSEP**
- 11. Plate sealing foils
- 12. Package insert

5. Wash buffer concentrate **WASH-BUF 20X**: 20x concentrated solution. Once diluted, the wash solution contains 10 mM phosphate buffer pH 7.0 +/-0.2, 0.05% Tween 20 and 0.05% Kathon GC.

6. Enzyme Conjugate **CONJ**: Ready to use and pink/red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Tris buffer, pH 6.8 +/-0.1, 0.1% Kathon GC and 0.02% gentamicin sulphate as preservatives.

7. Chromogen/substrate **SUBS/TMB**: Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 5.5-3.8, 4% dimethylsulphoxide, 0.05% tetra-methyl-ortho-quinone or TMB and 0.02% hydrogen peroxide or H2O2.

8. Assay Diluent **DILAS**: Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 5.5-3.8, 4% dimethylsulphoxide, 0.05% tetra-methyl-ortho-quinone or TMB and 0.02% hydrogen peroxide or H2O2.

9. Sulphuric Acid **H2SO4 0.3 M**: It contains 0.3 M H2SO4 solution.

10. Sample Diluent **DILSEP**: It contains 1% goat serum proteins, 10 mM Na-citrate buffer pH 8.0 +/-0.1, 0.5% Tween 20, 0.05% Na-azide and 0.1% Kathon GC as preservatives. To be used to dilute the sample.

Note: The diluent changes colour from olive green to dark bluish green in the presence of sample.

Important note: Only upon specific request, DiaPro can supply reagents for 96, 480, 960 tests, as reported below:

- 1. Microplate
- 2. Positive Control
- 3. Negative Control
- 4. Calibrator
- 5. Wash buff conc
- 6. Enz Conjugate
- 7. Chromog/Subs
- 8. Assay Diluent
- 9. Sulphuric Acid
- 10. Sample Diluent
- 11. Plate sealing foils
- 12. Package insert

**E. MATERIALS REQUIRED BUT NOT PROVIDED**

- 1. Calibrator-Microplate: (200µl and 100µl and 50µl)
- 2. EIA grade water (distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
- 3. Timer with 60 minute range or higher.
- 4. Absorbent paper tissues.
- 5. Calibrated ELISA microplate thermolabile incubator capable to provide a temperature of +37°C.
- 6. Calibrated ELISA microplate reader with 450nm (reading) and with 620-630nm (blanking) filters.
- 7. Calibrated ELISA microplate washer.
- 8. Vortex or similar mixing tools.

**F. WARNINGS AND PRECAUTIONS**

- 1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- 2. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
- 3. All the personnel involved in performing the assay have to wear protective laboratory clothes, face-face gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported by the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- 4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
- 5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/substrate from strong light and avoid vibration of the bench surface where the test is undertaken.
- 6. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
- 7. Do not interchange components between different lots of the kits. It is recommended that components between two lots of the same lot should not be interchanged.
- 8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
- 9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
- 10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
- 11. Do not use the kit after the expiration date stated on the external container and internal (vial) labels.
- 12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- 13. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
- 14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste.

**G. SPECIMEN PREPARATION AND RECOMMENDATIONS**

- 1. Blood is drawn aseptically by venipuncture and plasma or serum is prepared using standard techniques or preparation of samples (clinical laboratory assays). No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
- 2. Avoid any addition of preservatives to samples, especially sodium azide as this chemical would affect the enzymatic activity of the conjugate generating false negative results.
- 3. Samples have to be clearly identified with code or names in order to avoid misidentification of results. The kit is used for the screening of blood units, bar code labelling and electronic reading is strongly recommended.
- 4. Haemolyzed (red) and visibly hyperlipidemic (milky) samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
- 5. Sera and plasma can be stored at +2-8°C for up to seven days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
- 6. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.8µm filters to clean up the sample for testing.

**H. PREPARATION OF COMPONENTS AND WARNINGS**

- 1. Microplates: Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant is not turned to dark green, indicating a defect of manufacturing. In this case call DiaPro's customer service. Unused strips have to be placed back into the aluminium pouch, in presence of desiccant supplied, firmly zipped and stored at +2-8°C. When opened the first time, residual strips are stable till the indicator of humidity inside the desiccant bag turns from yellow to green.
- 2. Negative Control: Ready to use. Mix well on vortex before use.
- 3. Positive Control: Ready to use. Mix well on vortex before use. Handle this component as potentially infective, even if HCV, eventually present in the control, has been chemically inactivated.
- 4. Calibrator: Dissolve carefully the content of the lyophilized vial with the volume of EIA grade water reported on its label. Mix well on vortex before use.





Microplate														
A	BLK	S2	1	2	3	4	5	6	7	8	9	10	11	12
B	NC	S3												
C	NC	S4												
D	NC	S5												
E	CAL	S6												
F	CAL	S7												
G	PC	S8												
H	S1	S9												

BLK = Blank  
 NC = Negative Control  
 S = Sample  
 CAL = Calibrator  
 PC = Positive Control

**O. INTERNAL QUALITY CONTROL**

A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	< 0.050 mean OD450nm value after blanking
Calibrator	S/Co > 1.1
Positive Control	> 1.000 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section.  
 If they do not, do not proceed any further and operate as follows:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromopropylamine solution has not got contaminated during the assay.
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study. 2. that the proper washing solution has been used and the washer has been primed with it. 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control). 4. that no contamination of the negative control or of their wells has occurred due to positive samples, to spills or to the enzyme conjugate. 5. that the coplains haven't got contaminated with positive samples or with the enzyme conjugate. 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co < 1.1	1. that the procedure has been correctly executed. 2. that no mistake has been done in its control instead of dispensation of negative control. 3. that the washing procedure and the washer settings are as validated in the pre qualification study. 4. that no external contamination of the calibrator has occurred.
Positive Control < 1.000 OD450nm	1. that the procedure has been correctly executed. 2. that no mistake has been done in the distribution of controls (dispensation of negative control instead of positive control, in this case, the negative control will have an OD450nm value > 0.150, too). 3. that the washing-procedure and the washer settings are as validated in the pre qualification study. 4. that no external contamination of the positive control has occurred.

Should these problems happen, after checking, report any residual problem to the supervisor for further actions.

**P. CALCULATION OF THE CUT-OFF**

The tests results are calculated by means of a cut-off value determined with the following formula on the mean OD450nm value of the Negative Control (NC):

$$NC + 0.350 = \text{Cut-Off (CO)}$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

**Important note:** When the calculation of results is done by the operative system of an ELISA automated work station, be sure that the proper formulation is used to calculate the cut-off value and generate the right interpretations of results.

**Q. INTERPRETATION OF RESULTS**

Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused.  
 Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.  
 A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

**Important notes:**

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies earlier than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result. The sample has to be anyway submitted to a confirmation test (supplied upon request by DiaPro srl, code CO00N).
4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV/ antibodies - carrying anti IgM conjugate in the formulation of the enzyme tracer and therefore missing IgM reactivity - may be present. The real possibility of the sample (or antibodies to detect) should be confirmed by serology (anti IgM reactivity, IgM detection) or by HCV infection.
5. When test results are transferred from the laboratory to an information system, attention has to be done to avoid erroneous data transfer. Attention has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user:  
 Negative Control: 0.019 - 0.020 - 0.021 OD450nm  
 Mean Value: 0.020 OD450nm  
 Lower than 0.050 - Accepted  
 Positive Control: 2.189 OD450nm  
 Higher than 1.000 - Accepted  
 Cut-Off = 0.020 + 0.350 = 0.370  
 Calibrator: 0.550 - 0.539 OD450nm S/Co = 1.4  
 Mean value: 0.540 OD450nm  
 S/Co higher than 1.1 - Accepted  
 Sample 1: 0.070 OD450nm  
 Sample 2: 1.690 OD450nm  
 Sample 1 S/Co < 0.9 = negative  
 Sample 2 S/Co > 1.1 = positive

**R. PERFORMANCES**

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC).

**1. LIMIT OF DETECTION**

The limit of detection of the assay has been calculated by means of the British Working Standard for anti-HCV, NIBSC code 991368-005-SV1. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot # 1		Lot # 2	
	S/Co	S/Co	S/Co	S/Co
Factor	2.0	2.0	2.0	2.0
1 X	1.1	1.1	1.2	1.2
2 X	0.7	0.7	0.8	0.8
4 X	0.5	0.5	0.5	0.5
8 X	0.3	0.3	0.3	0.3
Negative plasma	0.3	0.3	0.3	0.3

In addition, the sample coded Accurun 1 - series 3000 - supplied by Boston Biomedica Inc., USA, has been evaluated in lot# showing the results below:

CVA/CE	Accurun 1	S/Co
Lot10	3000	1.5
Series	3000	1.5
Lot10	3000	1.9

In addition, n° 7 samples, tested positive for HCV Ab with Ortho HCV 3.0 S/AVE, code 930820, lot. # EXE065-1, were diluted in HCV Ab negative plasma in order to generate limiting dilutions and then tested again on CVA/CE lot. # 1202, and Ortho. The following table reports the data obtained.

Sample n°	Dilution	CVA/CE		Ortho 3.0	
		S/Co	S/Co	S/Co	S/Co
1	256 X	1.9	1.9	1.3	1.3
2	256 X	1.9	1.9	0.7	0.7
3	256 X	2.4	2.4	1.0	1.0
4	128 X	2.5	2.5	3.2	3.2
5	64 X	3.3	3.3	1.4	1.4
6	128 X	2.2	2.2	0.8	0.8
7	135 X	3.2	3.2	2.2	2.2

**2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY**

The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

**2.1 Diagnostic specificity:**

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. In addition to the first study, where a total of 5043 unrelated blood donors (including 1111 time donors), 210 hospitalized patients and 162 potentially reacting specimens (after infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemodialyzed, hepatic, etc.) were examined, the diagnostic specificity was recently assessed by testing a total of 2876 negative blood donors on six different lots. A value of specificity of 100% was found.  
 No false reactivity due to the method of specimen preparation has been observed. Both plasma derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well to check for interferences due to collection and storage.  
 No interference was observed.

**2.2 Diagnostic Sensitivity**

It is defined as the probability of the assay of scoring positive in the presence of specific analyte.  
 The diagnostic sensitivity has been assessed externally on a total number of 338 specimens, a diagnostic sensitivity of 100% was found. Internally, more than 600 positive samples were tested, providing a value of diagnostic sensitivity of 99.910%.  
 Positive samples from infections carried out by different genotypes of HCV were tested as well.  
 Furthermore, most of seroconversion panels available from Boston Biomedica Inc., USA, (PHV) and Zepionem, USA, (NCO) have been studied.  
 Results are reported below for some of them.

Panel	N° samples	DiaPro*	Ortho**
PHV 901	1	6	6
PHV 904	7	2	4
PHV 905	9	3	4
PHV 906	7	7	7
PHV 907	7	7	7
PHV 908	13	10	6
PHV 909	3	2	2
PHV 910	5	3	3
PHV 911	5	3	3
PHV 912	3	1	1
PHV 913	2	2	2
PHV 914	4	5	5
PHV 915	4	3	0
PHV 916	8	4	4
PHV 917	10	6	6
PHV 918	8	2	2
PHV 919	7	3	3
PHV 920	10	6	6
HCV 10039	5	2	2
HCV 6212	6	6	7
HCV 10165	9	5	4

\* HCV v.3.0  
 \*\* HCV v.3.0

Finally the Product, has been tested on the panel EFS AC HCV, lot n° 01108.03.22C017A, supplied by the Establishment Francis D.Serg (EFS), France, with the following results:



ESFS Panel Ac HCV

Sample	SiCo	Lot #1	Lot #2	Lot #3	Lot #4	Lot #5	Lot #6	Lot #7	Lot #8	Lot #9	Lot #10	Results
HCV 1	2.2	2.4	2.4	2.6	2.6	2.6	2.6	2.6	2.6	2.6	2.6	positive
HCV 2	1.6	2.0	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	2.1	positive
HCV 3	1.5	1.7	1.7	1.6	1.6	1.6	1.6	1.6	1.6	1.6	1.6	positive
HCV 4	5.2	6.5	5.5	5.5	5.5	5.5	5.5	5.5	5.5	5.5	5.5	positive
HCV 5	1.6	1.8	1.8	1.6	1.6	1.6	1.6	1.6	1.6	1.6	1.6	positive
HCV 6	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	0.4	negative

The variability shown in the tables above did not result in sample misclassification.

**3. PRECISION:**  
It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs. Results are reported as follows:

Lot # 1202

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.094	0.099	0.096
Std.Deviation	0.008	0.007	0.008
CV %	8.7	6.6	7.9
Average value	0.097		

Cal # 2 - 7K (N = 16)

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.366	0.403	0.418
Std.Deviation	0.023	0.029	0.027
CV %	5.9	7.1	6.4
Average value	0.402		

Lot # 0602

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.097	0.096	0.094
Std.Deviation	0.009	0.010	0.008
CV %	8.9	10.1	8.4
Average value	0.099		

Cal # 2 - 7K (N = 16)

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.400	0.395	0.393
Std.Deviation	0.021	0.025	0.026
CV %	5.4	6.2	6.6
Average value	0.402		

Lot # 0602/2

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.087	0.091	0.088
Std.Deviation	0.009	0.007	0.008
CV %	10.0	8.2	8.6
Average value	0.089		


Cal # 2 - 7K (N = 16)

Negative Sample (N = 16)			
Mean values	1st run	2nd run	3rd run
OD 450nm	0.386	0.390	0.391
Std.Deviation	0.023	0.021	0.023
CV %	6.0	5.3	5.8
Average value	0.389		

**REFERENCES**

- CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
- Aller MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:625-55.
- McCullian GM, Aller MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. *JAMA* 1996;275:2338-41.
- Dufour MC. Chronic liver disease and cirrhosis. In: *Eschbart JE, ed. Digestive diseases in the United States: epidemiology and impact.* US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases, Washington, DC: US Government Printing Office, 1994. NIH publication no. 94-1447, 615-45.
- Aller MJ, Hadler SC, Jenson FN, et al. Risk factors for acute non-A non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
- Aller HU, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis B antigen-hepatitis donors. *Ann Intern Med* 1972;77:917-19.
- Aller HU, Purcell RH, Holland PV, Feinstone SM, Wornow AG, Krawczynski K. Serologic analysis of transfusion-associated hepatitis. *Lancet* 1975;2:938-41.
- Saefi LB, Wright EC, Zimmerman HJ, McCollum RW. VA Cooperative Studies Group. VA cooperative study of posttransfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
- Feinstone SM, Kapikian AZ, Purcell RH, Aller HU, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
- Chang CL, Kuo G, Wheeler RL, Chen DJ, Hwang DJY. Isolation of a cDNA clone derived from a bloodborne non-A, non-B viral hepatitis genome. *Science* 1989;244:559-62.
- Kuo G, Chao QL, Alter HU, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
- Aller HU, Purcell RH, Shin JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1594-950.
- Aach RD, Sturgeon CE, Halperin EB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1823-9.
- Aller MJ, Margolis HS, Krawczynski K, Jenson FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1898-1905.
- Aller MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
- Donahue JG, Nelson KE, Madao A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:626-31.
- Zeldin JB, Jahn S, Krawczynski K, et al. Seroprevalence of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:50-5.
- Fingeroth MD, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
- Garten RS, Vlahov D, Galai N, Donnelly MC, Nelson KE. Viral infections in street-drug injection drug users: the prevalence of the hepatitis B virus, hepatitis C virus, human immunodeficiency virus, and human T-lymphotropic virus. *Am J Pub Health* 1993;83:557-61.
- Breider DB, Aller HU, Denning JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
- Troisi CL, Hollinger FB, Hoos WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:142-6.
- Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis in patients with hemophilia. *J Virol Methods* 1993;41:205-9.
- Tekers JL, Miller ER, Aller MJ, Ardoin M, National surveillance of dialysis associated diseases in the United States, 1985. *ASAIO Journal* 1988;44:98-107.
- Ormond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1992;165:755-71.
- Weinstock HS, Bolen G, Rainald AL, Pollan LB. Hepatitis C virus infection in a community clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
- Thomas DL, Cannon RO, Shingo CN, Hook EW, III, Aller MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
- Burchincher SP, Kasz WH, Hessel NA, Liu J, O'Malley PM, Aller MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect Dis* 1994;229:253-9.
- Thomas DL, Zankaman JM, Alter HU, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted disease clinics in Baltimore—an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
- Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E, Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
- Cooper BW, Kuzell A, Tilton RC, Goodman R, Lenz R. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;17:519-23.

All the IVD Products manufactured by the company are under the control of a certified Quality Management System approved by an EC-Notified Body. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

  
 0318  
 Manufacturer:  
 Dia Pro Diagnostico Bioprobes Srl  
 Via G. Carducci n° 27 - Sesto San Giovanni (MI) - Italy

### HAV IGM

**A. INTENDED USE**  
 Enzyme Immunoassay (ELISA) for the determination of IgM class antibodies to Hepatitis A Virus in human plasma and sera with the "capture" system. The kit may be used for the identification of the viral agent causing hepatitis in the patient and the follow up of the acute phase of the infection.  
 For "in vitro" diagnostic use only.

**B. INTRODUCTION**  
 The Centre for Disease Control or CDC, Atlanta, USA, defines Hepatitis A Virus as follows:  
 Reported vaccine-preventable diseases in the world despite the licensure of hepatitis A vaccine in 1985. Widespread vaccination of appropriate susceptible populations would substantially lower disease incidence and potentially eliminate indigenous transmission of hepatitis A virus (HAV) infection.

# HAV IGM

**"Capture" Enzyme Immuno Assay (ELISA)**  
**for the determination of IgM class**  
**antibodies to Hepatitis A Virus**  
**in human plasma and sera**

- for "in vitro" diagnostic use only -




**DI.A.PRO**  
**Diagnostic Bioprobes Srl**  
 Via G. Carducci n° 27  
 20099 Sesto San Giovanni  
 (Milano) – Italy  
 Phone +39 02 27007161  
 Fax +39 02 26007226  
 e-mail: info@diapro.it

REF\_AVM/CE  
 06/15/15

HAV infection is acquired primarily by the fecal-oral route by either person-to-person contact or ingestion of contaminated food or water. On rare occasions, HAV infection has been transmitted by transfusion of blood or blood products collected from donors during the viraemic phase of their infection. In experimentally infected nonhuman primates, HAV has been detected in saliva during the incubation period; however, transmission by saliva has not been demonstrated.  
 Depending on conditions, HAV can be stable in the environment for months. Heating foods at temperatures greater than 185 °F (65 °C) for 1 minute or disinfecting surfaces with a 1:100 dilution of sodium hypochlorite (i.e., household bleach) in tap water is necessary to inactivate HAV.

Because most children have asymptomatic or unrecognized infections, they play an important role in HAV transmission and serve as a source of infection for others. In one study of adults without an identified source of infection, 32% of their households included a child less than 6 years old, and the presence of a young child was associated with HAV transmission within the household. In studies where serologic testing of the household contacts of adults without an identified source of infection was performed, 25%–40% of the contacts less than 6 years old had serologic evidence of acute HAV infection (IgM and HAV).



**C. PRINCIPLE OF THE TEST**  
 The assay is based on the principle of "IgM capture" where IgM class antibodies in the sample are first captured by the solid phase coated with anti-IgM antibody.  
 After washing out all the other components of the sample and in particular IgG antibodies, the specific IgM captured on the solid phase are detected by the addition of a purified preparation of inactivated HAV, labelled with an antibody conjugated with peroxidase (HRP).  
 After incubation, microwells are washed to remove unbound conjugate and then the chromogen/substrate is added.  
 In the presence of peroxidase the colorless substrate is hydrolysed to a colored end-product, whose optical density may be detected and is proportional to the amount of antibodies to HAV present in the sample.



**B. COMPONENTS**  
The kit contains reagents for 96 tests.

- 11. Plate sealing tools n° 2
- 12. Package insert n° 1
- E. MATERIALS REQUIRED BUT NOT PROVIDED**
  - 1. Calibrated Micropipettes of 10µl, 100µl and 1000µl and disposable plastic tips.
  - 2. EIA grade water (double distilled or deionised, chemical treated to remove oxidizing chemicals used as disinfectants).
  - 3. Timer with 50 minute range or higher.
  - 4. Absorbent paper tissues.
  - 5. Calibrated ELISA microplate thermostat incubator (dry or wet) set at +37°C (+/-0.1°C tolerance).
  - 6. Calibrated ELISA microwell reader with 450nm (reading) and with 620-630nm (blanking) filters.
  - 7. Calibrated ELISA microplate washer.
  - 8. Vortex or similar mixing tools.

**F. WARNINGS AND PRECAUTIONS**

- 1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory.
- 2. All the personnel involved in performing the assay have to wear protective laboratory clothes, lab-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported by the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- 3. All the personnel involved in sample handling should be vaccinated for HEV and HAV, for which vaccines are available, safe and effective.
- 4. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-borne microbial agents, when opening kit vials and microplates and when performing the test. Protect the Chromogen/Substrate (TMB & H2O2) from strong light and avoid vibration of the bench surface where the test is undertaken.
- 5. Upon receipt, store the kit at 2-8°C into a temperature controlled refrigerator or cold room.
- 6. Do not interchange components between different lots of the kits. It is recommended that components from two kits of the same lot should not be interchanged.
- 7. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures.
- 8. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
- 9. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
- 10. Do not use the kit after the expiration date stated on external (primary container) and internal (vials) labels. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
- 11. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of autoclaved workstations in order to avoid contamination.
- 12. Waste produced during the use of the kit has to be considered in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure from residues of controls and from samples has to be treated as potentially infective material and

**G. SPECIMEN: PREPARATION AND RECOMMENDATIONS**

- 1. Blood is drawn aseptically by venopuncture and plasma or serum is prepared using standard techniques. No preservatives have been observed in the preparation of the sample with citrate, EDTA and heparin.
  - 2. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labelling and reading is strongly recommended.
  - 3. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
  - 4. Sera and plasma can be stored at +2° -8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be refreezed more than once as IgM antibodies may get denatured and as this procedure may generate particles that could affect the test result.
  - 5. If particles are present, centrifuge at 2,000 rpm for 20 min or filter using 0.2-0.8µm filters to clean up the sample for testing.
- H. PREPARATION OF COMPONENTS AND WARNINGS**  
A study conducted on an opened kit has not pointed out any relevant loss of activity up to 3 months.
- 1. **Antibody coated microwells:**  
Allow the microplate to reach room temperature (about 1 hr) before opening the container. Check that the desiccant has not turned dark green, indicating a defect in conservation. In this case, call Dia.Pro's customer service.  
Unused strips have to be placed back into the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2° -8°C. When opened the first time, unused strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.
  - 2. **Negative Control:**  
Ready to use. Mix well on vortex before use.
  - 3. **Positive Control:**  
Ready to use. Mix well on vortex before use. Handle this component as potentially infectious, even if HAV, eventually present in the control, has been chemically inactivated.
  - 4. **Calibrator:**  
Add the volume of ELISA-grade water reported on the label, to the lyophilized powder, let it fully dissolve and then gently mix on vortex at -20°C.  
The solution is not stable. Store the Calibrator frozen in aliquots at -20°C.  
*Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.*

**5. Wash buffer concentrate**

- The reagent concentrate for concentrated solution has to be diluted 20x with distilled water and mixed gently end-over-end before use.
  - Once diluted, the wash solution is stable for 1 week at 2-8°C. During preparation, avoid frothing as the presence of bubbles could impact on the efficiency of the washing cycles.  
*Note: Once diluted, the wash solution is stable for 1 week at +2-8°C.*
  - 6. **Enzyme conjugate**  
20X preparation: Mix well on vortex.  
Avoid contamination of the liquid with oxidizing chemicals, dust or microbes when the reagent is aspirated to be used.
  - 7. **HAV Antigen:**  
Ready to use. Mix well on vortex before use.  
Handle this component as potentially infectious, even if HAV has been chemically inactivated.
  - 8. **Sample Diluent:**  
Ready to use. Mix well on vortex before use.
  - 9. **Chromogen/Substrate:**  
Ready to use. Mix well on vortex before use.  
Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes.  
Do not expose to strong light, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic and if possible, sterile disposable container.
  - 10. **Suphonic Acid:**  
Ready to use. Mix well on vortex before use.  
Attention: (H315, H319, P280, P302+P352, P332+P313, P305+P351+P338, P337+P313, P602+P603).
- Legenda:**  
Warning H statements:  
H315 – Causes skin irritation.  
H319 – Causes serious eye irritation.  
Precautionary P statements:  
P280 – Wear protective gloves/protective clothing/eye protection/face protection.  
P302 + P352 – IF ON SKIN: Wash with plenty of soap and water.  
P332 + P313 – If skin irritation occurs: Get medical advice/attention.  
P305 + P351 + P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337 + P313 – If eye irritation persists: Get medical advice/attention.  
P362 + P363 – Take off contaminated clothing and wash it before reuse.

**1. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT**

- 1. Micropipettes have to be calibrated to deliver the correct volume (tolerance +0.5%) required by the assay and must be submitted to regular decontamination (household

- 10% solution of bleach, hospital grade disinfectant of these parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of ±2%.
2. The ELISA incubator has to be set at +37°C (tolerance of ±4.05°C) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performance of the assay. The washer must be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests. 4-5 washing cycles (aspiration & dispensation of 350ul/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested in order to set correctly their number. It is recommended to run an assay with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section O. Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of ±5%.
5. The ELISA reader has to be equipped with a reading filter of 450nm and with a second filter (520-630nm, strongly recommended) for blanking purposes. Blanking is carried out on the well identified in the section Assay Procedure. The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in the section O. Internal Quality Control. The assay protocol has to be installed in the operating system. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry-over. This has to be validated and controlled to minimize the possibility of contamination of adjacent wells. The use of ELISA automated work stations is recommended when the number of samples to be tested exceed 20-30 units per run.

4. Dissolve the Calibrator as described above and gently mix.
5. Allow all the other components to reach room temperature (allow 1 hr) and then mix gently on vortex.
6. Set the ELISA incubator at +37°C ±4.0-1°C and prepare the ELISA washer by priming with the dilution washing solution, according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the microplates are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

**ASSAY PROCEDURE**

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

1. Dilute samples 1:101 by dispensing first 10 µl sample and then 1 ml Sample Diluent into a dilution tube, mix gently on vortex.
  2. Place the required number of Microwells in the microwell holder. Leave the 1<sup>st</sup> well empty for the operation of blanking.
  3. Dispense 100 µl Negative Control in triplicate, 100 µl Positive Control in single and 100 µl Calibrator in duplicate in proper wells. Do not dilute controls and the calibrator as they are ready to use!
  4. Dispense 100 µl diluted samples in the proper sample wells and then check that all the samples wells are blue coloured and that controls and calibrator have been dispensed.
  5. Incubate the microplate for 60 min at +37°C.
- Important note:** Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. About 5-10 minutes before use, prepare the HAV Antigen/Antibody immunocomplex as described previously.
  7. Wash the microplate with an automatic washer as reported previously (section 1.3).
  8. Pipette 100 µl HAV Antigen/Antibody complex into each well, except the 1<sup>st</sup> blanking well, and cover with the sealer. Check that all wells are red coloured, except A1.
- Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the EnZyme Conjugate. Contamination might occur.
9. Incubate the microplate for 60 min at +37°C.
  10. Wash microwells as in step 7.
  11. Pipette 100 µl Chromogen/Substrate mixture into each well, the blank well included. Then incubate the microplate at room temperature (18-24°C) for 20 minutes.

- L. PRE ASSAY CONTROLS AND OPERATIONS**
1. Check the expiration date of the kit printed on the external label (primary container). Do not use the device if expired.
  2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
  3. Dilute all the content of the 20x concentrated Wash Solution as described above.

**Important note:** Do not expose to strong direct illumination. High background might be generated.

12. Pipette 100 µl Sulphuric Acid into all the wells to stop the enzymatic reaction using the same pipetting sequence as in step 10. Addition of acid will turn the positive control and positive samples from blue to yellow.
13. Measure the colour intensity of the solution in each well as described in section 1.5, at 450nm filter (reading) and at 630-630nm (background subtraction, strongly recommended), blanking the instrument on A1.

**Important notes:**

1. If the second filter is not available ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results or reading.
2. Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 20 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

**N. ASSAY SCHEME**

Controls/Calibrator (*)	100 µl samples diluted 1:101
1 <sup>st</sup> incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
HAU & Tracer	100 µl
2 <sup>nd</sup> incubation	60 min
Temperature	+37°C
Washing	4-5 cycles
TUBERCULOZ mix	100 µl
3 <sup>rd</sup> incubation	20 min
Temperature	RT
Sulphuric Acid	100 µl
Reading OD	450nm

**(\*) Important Notes:**

- The Calibrator (CAL) does not affect the Cut Off calculation, therefore it does not affect the test results calculation.
- The Calibrator (CAL) used is required by the Internal Quality Control as required by the Management.

An example of dispensation scheme is reported in the table below.

	Microplate											
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL(1)	S6										
F	CAL(1)	S7										
G	PC	S8										
H	SI	S9										

Legend: BLK = Blank, NC = Negative Control, PC = Positive Control, SI = Sample

**O. INTERNAL QUALITY CONTROL**

A check is performed on the controls any time the kit is used in order to verify whether the expected OD450nm or S/Co values have been matched in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
Negative Control	< 0.150 OD450nm value after Blanking
Mean value (NC)	Coefficient of variation < 50%
Positive Control	> 0.500 OD450nm

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/substrate solution has not become contaminated during the assay.
Negative Control > 0.150 OD450nm	1. that the washing procedure and the washer settings are as validated in the pre-qualification study.
coefficient of variation > 50%	2. that the proper washing solution has been used and the washer has been primed with it before use.
Positive Control < 0.500	1. that the procedure has been correctly performed.
OD450nm	2. that no mistake has occurred during the dispensation and reading (dispensation of sample and reading of positive control).
	3. that the washing procedure and the washer settings are as validated in the pre-qualification study.
	4. that no external contamination of the positive control has occurred.
	6. that the washer needles are not blocked or partially obstructed.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

If Calibrator has used, verify the following data:

Parameter	Requirements
Calibrator	S/Co > 1

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Calibrator S/Co < 1	1. that the procedure has been correctly performed.
	2. that no mistake has occurred during its dispensation (e.g.: dispensation of negative control).
	3. that the washing procedure and the washer settings are as validated in the pre-qualification study.
	4. that no external contamination of the calibrator has occurred.

Anyway, if all other parameters (Blank, Negative Control, Positive Control), match the established requirements, the test may be considered valid.

**P. CALCULATION OF THE CUT-OFF**

The last results are calculated by means of the mean OD450nm value of the Negative Control (NC) and a mathematical calculation, in order to define the following cut-off formulation.



**Cut-Off = NC + 0.250**

The value found for the test is used for the interpretation of results as described in the next paragraph.

**Important note:** When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

**Q. INTERPRETATION OF RESULTS**

Test results are interpreted as a ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.8	Negative
0.8 – 1.2	Equivocal
> 1.2	Positive

A negative result indicates that the patient is not undergoing an acute infection by HAV.  
Any patient showing an equivocal result, should be re-tested by examining a second sample after 1-2 weeks from first testing.  
A positive result is indicative of an HAV infection event and therefore the patient should be treated accordingly.  
An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user.

**Negative Control:** 0.050 – 0.060 – 0.070 OD450nm  
**Mean Value:** 0.060 OD450nm  
**Lower than 0.150 – Accepted**  
**Positive Control:** 2.193 OD450nm  
**Higher than 0.500 – Accepted**

Cut-Off = 0.060+0.250 = 0.310

**Calibrator:** 0.550 – 0.530 OD450nm  
**Mean value:** 0.540 OD450nm  
**S/Co higher than 1.0 – Accepted**

**Sample 1:** 0.070 OD450nm  
**Sample 2:** 1.690 OD450nm  
**Sample 1 S/Co < 0.8 = negative**  
**Sample 2 S/Co > 1.2 = positive**

**Important notes:**

1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgment errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method (confirmation test) before a diagnosis of viral hepatitis is confirmed.
3. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
4. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

**R. PERFORMANCE CHARACTERISTICS**

**1. Limit of detection**  
In absence of a defined international standard for HAV IgM, the limit of detection of the assay has been calculated by means of the following preparations:

1. Accurn # 121 supplied by Boston Biomedica Inc. – USA
2. Accurn # S1 supplied by Boston Biomedica Inc. USA

These preparations were prepared according to the manufacturer's instructions, diluted in Sample Diluent (1:100) and then further diluted in Sample Diluent to generate a limiting curve (accurn # 121).  
Results of Quality Control are given in the following table:

Preparation	Dilutions	S/Co
Accurn # 121	1:100	5.4
	1:200	4.1
	1:400	2.8
	1:800	1.9
Accurn # S1	1:1600	1.0
	1:100	4.2

**2. Diagnostic Sensitivity:**

The diagnostic sensitivity has been tested on panels of samples classified positive by a US FDA approved kit.  
Positive samples were collected from patients carrying HAV acute infection, confirmed by clinical symptoms and analysis:

An overall value of 100% has been found in the study conducted on a total number of more than 100 samples.  
A seroconversion panel has also been studied.

Results obtained by examining a preparation supplied by Boston Biomedica Inc., USA, are reported below.

**Seroconversion Panel : PH1 902**

Sample	OD450nm	S/Co	Diasorin Refer.	S/Co	Score
CTRL (-)	0.048	0.2			
CTRL (+)	1.736	5.8			
PH1902	1	0.037	0.1	0.3	neg
	2	0.042	0.1	0.3	neg
	3	1.956	6.6	6.8	pos
	4	1.988	6.7	6.7	pos
	5	0.669	2.2	1.5	pos

**3. Diagnostic Specificity:**

The diagnostic specificity has been determined on panels of specimens, negative with the reference kit, derived from normal individuals and blood donors of European origin.  
Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the specificity. No false reactivity due to the method of specimen preparation has been observed.

Frozen specimens have also been tested to check whether this interferes with the performance of the test. No interference was observed on clean and particle free samples.  
Samples derived from patients with different viral (HCV, HDV, HBV, HEV) and non viral pathologies of the liver that may interfere with the test were examined.  
No cross reaction were observed.  
The Performance Evaluation study conducted in a qualified external reference centre on more than 500 samples has provided a value > 98%.

**3. Precision:**

It has been calculated on two samples, one negative and one low positive, examined in 16 replicates in three separate runs. Results are reported as follows:

Sample	Negative	Low Pos.
OD450nm	0.058	0.719
Std. Deviation	0.008	0.052
CV %	14.3	7.2

**Test # 1**

Sample	Negative	Low Pos.
OD450nm	0.048	0.709
Std. Deviation	0.007	0.053
CV %	13.9	8.9

**Test # 3**

Sample	Negative	Low Pos.
OD450nm	0.050	0.713
Std. Deviation	0.007	0.055
CV %	13.4	7.7

**5. LIMITATIONS**

False positivity has been assessed as less than 2% of the normal population, mostly due to high titers of RF.  
Frozen samples containing fibrin particles or aggregates may generate false positive results.  
Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analyte.  
This test is suitable only for testing single samples and not pooled ones.  
Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, serology, as well as other diagnostic data should be considered.

**REFERENCES**

1. Dienstag J.L., "Hepatitis A Virus: Identification, characterization and epidemiologic investigations", Progress in liver disease VI, Poppe E., Schaffner F., (eds), pp. 343-370, New York, Grune and Stratton, 1979.
2. Duemeyer W., Van der Veen J., Koster B., "ELISA in Hepatitis A", Lancet, I, 823-824, 1978
3. Parry J.V., (1981) "Hepatitis A infection guidelines for the development of satisfactory assays for laboratory diagnosis", The Institute of Medical Laboratory Sciences, 38, 303-311.
4. Lindberg J., Foster G., Hansson B.G. et al., "Seroologic markers of hepatitis A and B in chronic active hepatitis", Scandinavian Journal of Gastroenterology, 13:525-527, 1978.
5. Barbara J.A., Howell D.R., Biggs M., Parry J.V., "Post transmission hepatitis A", Lancet (1982), 1-738.
6. Zachoval R., Dienstag J.L., Purcell R.H., "Tests for hepatitis A virus antigen and antibody", in "Hepatitis A", Gevery R.J. (Ed), pp 33-46, Orlando, Academic Press, Inc, 1984

All the IVD Products manufactured by the company are under the control of a certified Quality Management System in compliance with ISO 13485 rule. Each lot is submitted to a quality control and released into the market only if conforming with the EC technical specifications and acceptance criteria.

Manufacturer:  
Dia-Pho Diagnostic Bioprobes Srl  
Via G. Carducci n° 27 - Sesto San Giovanni (MI) - Italy



## HAV Ab

# HAV Ab

### Competitive Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis A Virus in human plasma and sera

- for "in vitro" diagnostic use only -



**DIA.PRO**  
**Diagnostica Bioprobes Srl**  
Via Carducci n° 27  
20099 Sesto San Giovanni  
(Milano) - Italy  
Phone +39 02 27007161  
Fax +39 02 26007726  
e-mail: [info@diapro.it](mailto:info@diapro.it)

REF AVAB CE  
96 Tests

**A. INTENDED USE**  
Competitive Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis A Virus in human plasma and sera. The kit is used for the follow-up of patients infected by HAV. For "in vitro" diagnostic use only.

**B. INTRODUCTION**  
The Center for Disease Control or CDC, Atlanta, USA defines Hepatitis A Virus as follows:

Hepatitis A continues to be one of the most frequently reported vaccine-preventable diseases in the world, despite the licensure of hepatitis A vaccine in 1985. Widespread vaccination of appropriate susceptible populations would substantially lower disease incidence and potentially eliminate indigenous transmission of hepatitis A virus (HAV) infection.

HAV, a 27-nm RNA virus, classified as a picornavirus, can produce either asymptomatic or symptomatic infection in humans after an average incubation period of 28 days (range, 15-50 days). The illness caused by HAV infection typically has an abrupt onset of symptoms that can include fever, malaise, anorexia, nausea, abdominal discomfort, dark urine, and jaundice. The likelihood of having symptoms with HAV infection is related to the person's age; in children less than 6 years of age, most (70%) infections are asymptomatic; if illness does occur, it is not usually accompanied by jaundice. Among older children and adults, infection is usually symptomatic, with jaundice occurring in greater than 70% of patients. Signs and symptoms usually last less than 2 months, although 10%-15% of symptomatic persons have prolonged or relapsing disease lasting up to 6 months. In infected persons, HAV replicates in the liver; is excreted in bile, and is shed in the stool. Peak infectivity of infected persons occurs during the 2-week period before onset of jaundice or elevation of liver enzymes, when the concentration of virus in stool is highest. The concentration of virus in stool declines after jaundice appears. Children and infants can shed HAV for longer periods than adults; up to several months after the onset of clinical illness. Chronic shedding of HAV in feces does not occur; however, shedding can occur in persons who have relapsing illness. Viremia occurs soon after infection and persists through the period of liver enzyme elevation.

Hepatitis A cannot be differentiated from other types of viral hepatitis on the basis of clinical or epidemiologic features alone. Serologic testing to detect immunoglobulin M (IgM) antibody to the capsid proteins of HAV (IgM anti-HAV) is required to confirm a diagnosis of acute HAV infection. In most persons, IgM anti-HAV becomes detectable 5-10 days before the onset of symptoms and can persist for up to 6 months after infection. Immunoglobulin G (IgG) anti-HAV, which appears early in the course of infection, remains detectable for the patient's lifetime and confers lifelong protection against the disease. Commercial diagnostic tests are available for the detection of IgM anti-HAV and IgG anti-HAV in serum.


HAV RNA can be detected in the blood and stool of most persons during the acute phase of infection by using nucleic acid amplification methods, and nucleic acid sequencing has been used to determine the relatedness of HAV isolates.

HAV infection is acquired primarily by the fecal-oral route by either person-to-person contact or ingestion of contaminated food or water. On rare occasions, HAV infection has been transmitted by transfusion of blood or blood products collected from donors during the viremic phase of their infection. In experimentally infected nonhuman primates, HAV has been detected in saliva during the incubation period; however, transmission by saliva has not been demonstrated.

Depending on conditions, HAV can be stable in the environment for months. Heating foods at temperatures greater than 185 F (85 C) for 1 minute or disinfecting surfaces with a 1:100 dilution of sodium hypochlorite (i.e., household bleach) in tap water is necessary to inactivate HAV.

Because most children have asymptomatic or unrecognized infections, they play an important role in HAV transmission and serve as a source of infection for others. In one study of adults without an identified source of infection, 52% of their households included a child less than 6 years old, and the presence of a young child was associated with HAV transmission within the household. In studies where serologic testing of the household contacts of adults without an identified source of infection was performed, 23%-40% of the contacts less than 6 years old had serologic evidence of acute HAV infection (IgM anti-HAV).

**Hepatitis A**  
**Typical sequence of events following infection**



### C. PRINCIPLE OF THE TEST

The assay is based on the principle of competition, where the antibodies in the sample compete with an anti-HAV specific antibody labeled with HRP, for a fixed amount of antigen on the solid phase.

A purified and inactivated HAV is added to the microwell. The patient's serum/plasma is added to the microwell and antibodies to HAV are captured by the solid phase.

After washing, the enzyme conjugate is added and binds to the free HAV antigen, if still present.

The plate is washed to remove unbound conjugate and then the chromophore substrate is added.

In the presence of peroxidase, the colorless substrate is hydrolyzed to a colored end-product, whose optical density may be detected and is inversely proportional to the amount of antibodies to HAV present in the sample.

An additive is added to the sample directly into the well to block interferences able to mask the presence of antibodies, mostly appearing in the follow up of vaccination.







assay, with the kit controls and well characterized negative and positive reference samples, and check to match the values reported below in the section O Internal Quality Control. Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of the needles) of the washer has to be carried out according to the instructions of the manufacturer.

- The ELISA reader has a tolerance of 15%.
- The ELISA reader has to be equipped with a reading filter of 450nm (added) in a second filter (620-630nm, coming out recommended) for blanking purposes. Blanking is carried out on the well(s) defined in the section Assay Procedure. The optical system of the reader has to be calibrated regularly to ensure the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.
- When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to obtain the values reported in the section O Internal Quality Control. The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and aspiration) must be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing and for washing. This must be validated and controlled to minimize the possibility of contain holder adjacent wells. The use of ELISA automated work station is recommended when the number of samples to be tested exceeds 20-30 units per run.
- Dia-Pro's customer service offers support to the user in the setting and checking of instruments used in combination with the kit, in order to assure compliance with the requirements described. Support is also provided for the installation of new instruments to be used with the kit.

**L. PRE ASSAY CONTROLS AND OPERATIONS**

- Check the expiration date of the kit printed on the external label (primary container). Do not use if expired.
- Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Chromogen/Substrate (TMB-H<sub>2</sub>O<sub>2</sub>) is colorless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no leakage occurred in transportation and no spillage of liquid is present inside the box (primary container). Check that the aluminum pouch, containing the Chromogen, is not punctured or damaged.
- Dilute all the content of the 20x concentrated Wash Solution as described above.
- Disassemble the calorator as described above and gently mix (about 10s) and then mix gently on vortex all liquid reagents.
- Set the ELISA incubator at +37°C and prepare the ELISA washer by priming with the diluted washing solution according to the manufacturer's instructions. Set the right number of washing cycles as found in the validation of the instrument for use with the kit.
- Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
- If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
- Check that the microplates are set to the required volume.
- Check that all the other equipment is available and ready to use.
- In case of problems, do not proceed further with the test and advise the supervisor.

**M. ASSAY PROCEDURE**

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

- Place the required number of strips in the microplate holder. Leave A1 well empty for the operation of Blanking. Store the other strips into the bag in presence of the desiccant at +2,8°C, sealed.
- Dispense 50 µl Specimen Diluent in all the wells identified for samples and control/calibrator, except for A1. Then pipette 100 µl of Negative Control in triplicate, 100 µl of Calibrator in duplicate, 100 µl Positive Control in single and then 100 µl of samples. Check that control/calibrator and samples have been correctly added. Incubate the microplate at +37°C for 60 min.
- Wash the microplate as reported in section I.3.
- In all the wells except A1, pipette 100 µl Enzyme Conjugate. Check that the reagent has been correctly added. Incubate the microplate at +37°C for 60 minutes.

**Important note:** Be careful not to touch the inner surface of the well with the pipette tip when dispensing the Enzyme Conjugate. Contamination might occur.

- Wash the microplate as described.
- Pipette 100 µl TMB-H<sub>2</sub>O<sub>2</sub> mixture in each well. The blank wells included. Check that the reagent has been correctly added. Then incubate the microplate at room temperature for 20 minutes. **Important note:** Do not expose to strong direct light as a high background might be generated.
- Pipette 100 µl Sulphuric Acid into each well to stop the enzymatic reaction using the same pipetting sequence as in step 6. Then measure the color intensity with a microplate reader at 450nm (reading area at 620-630nm background subtraction, strongly recommended), blanking the instrument on A1 well.

**N. ASSAY SCHEME**

**Important notes:**

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading has to be ideally performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background.

Specimen Diluent	50 µl
Controls/Calibrator (*)	100 µl
1 <sup>st</sup> incubation	60 min
Temperature	+37°C
Washing step	4-5 cycles
Enzyme Conjugate	100 µl
2 <sup>nd</sup> incubation	60 min
Temperature	+37°C
Washing step	4-5 cycles
TMB/H <sub>2</sub> O <sub>2</sub> mix	100 µl
3 <sup>rd</sup> incubation	20 min
Temperature	rt
Sulphuric Acid	100 µl
Reading OD	450nm & 620nm

**(\*) Important Notes:**

- The Calibrator (CAL) does not affect the Cut-Off calculation, therefore it does not alter the test results calculation.
  - The Calibrator (CAL) used only if a laboratory internal quality control is required by the Management.
- An example of dispensation scheme is reported in the table below.

	Microplate											
	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL (*)	S6										
F	CAL (*)	S7										
G	PC	S8										
H	S1	S9										

Legend: BLK = Blank, NC = Negative Control, CAL (\*) = Calibrator - Not mandatory, PC = Positive Control, S = Sample

**O. INTERNAL QUALITY CONTROL**

A check is performed on the controls any time the kit is used in order to verify whether the expected OD450nm or CoS values have been reached in the analysis. Ensure that the following parameters are met:

Parameter	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	> 0.750 mean OD450nm value after blanking
Positive Control	< 0.300 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Chromogen/Substrate solution has not become contaminated during the assay.
Negative Control (NC) < 0.750 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre-qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to the enzyme conjugate; 5. that microplates have not become contaminated with positive samples or with the enzyme conjugate; 6. that the washer needles are not blocked or partially obstructed.
coefficient of variation > 30%	

**Positive Control > 0.300 OD450nm**

- that the procedure has been correctly performed;
- that no mistake has occurred during the distribution of the control (ex: dispensation of negative control instead of the positive one);
- that the washing procedure and the washer settings are as validated in the pre-qualification study;
- that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

If Calibrator has used, verify the following data:

Parameter	Requirements
Calibrator 10 mIU/ml (WHO)	CoS > 1.0

If the results of the test doesn't match the requirements stated above, operate as follows:

Problem	Check
Calibrator	1. that the procedure has been correctly performed;
CoS < 1.0	2. that no mistake has occurred during its distribution; 3. that the washing procedure and the washer settings are as validated in the pre-qualification study; 4. that no external contamination of the calibrator has occurred.

Anyway, if all other parameters (Blank, Negative Control, Positive Control), match the established requirements, the test may be considered valid.

**P. CALCULATION OF THE CUT-OFF**

The test results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off} = (\text{NC} + \text{PC}) / 3$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

**Important note:** When the calculation of results is performed by the operating system of an ELISA automated work station, make sure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

**Q. INTERPRETATION OF RESULTS**

Test results are interpreted as a ratio of the Cut-Off value and the OD450nm of the sample (or CoS) according to the following table:

CoS	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by rAVV.



Any pattern showing an equivocal result should be retested on a second sample taken 1-2 weeks after the initial sample. A positive result is indicative of a past or recent HAV infection and therefore the patient should be treated accordingly.

An example of calculation is reported below.  
The following data must not be used instead of real figures obtained by the user:

Negative Control: 1.900 - 2.000 - 2.100 OD450nm  
Mean Value: 2.000 OD450nm  
Higher than 0.750 - Accepted  
Positive Control: 0.100 OD450nm  
Lower than 0.300 - Accepted

Cut-Off = (2.000 + 0.100) / 3 = 0.700  
Calibrator: 0.400-0.360 OD450nm  
Mean value: 0.380 OD450nm  
CoS > 1 - Accepted

Sample 1: 0.050 OD450nm  
Sample 2: 1.900 OD450nm  
Sample 1 CoS > 1.1 positive  
Sample 2 CoS < 0.9 negative

**Important notes:**


1. Interpretation of results should be done under the supervision of the laboratory supervisor to reduce the risk of judgement errors and misinterpretations.
2. When test results are transmitted from the laboratory to another facility, attention must be paid to avoid erroneous data transfer.
3. Diagnosis of viral hepatitis infection has to be taken by and released to the patient by a suitably qualified medical doctor.

**R. PERFORMANCE CHARACTERISTICS**

**1. Limit of detection**  
The limit of detection of the assay has been calculated by means of the 2<sup>nd</sup> International Standard supplied by WHO. Two control samples, supplied by Boston Biomedica Inc., USA, with code Accuam 52 and 120, were also examined. The sensitivity shown by the assay is < 10 WHO mIU/ml or < 5 PEI mIU/ml.  
Results of Quality Control are given in the following table:

WHO mIU/ml	OD450 CoS	PEI mIU/ml	OD450 CoS
50	0.099	7.7	0.093
25	0.197	3.9	0.137
10	0.543	1.4	0.304
5	0.943	0.8	0.587
2.5	1.015	0.7	0.949
Neg. Control	2.217	Neg. Control	2.217
Accuam 52	0.060	12.7	Accuam 120
			0.115
			6.6

Curves are reported below:



**2. Diagnostic sensitivity:**  
The diagnostic sensitivity has been tested in a clinical trial on panels of samples classified positive by a US FDA approved kit. An overall value of 100% has been found in the study conducted on a total number of more than 200 samples. Seroreconversion and performance panels have also been studied. Results obtained by examining two panels supplied by Boston Biomedica Inc., USA, are reported below.

**Seroreconversion Panel: PHT 902**

Sample	OD450nm	CoS	Diagnosis
CTRL (-)	1.968	0.3	neg
CTRL (+)	0.084	8.1	pos
Calibrator	0.470	1.5	pos
PHT902			
1	1.878	0.4	neg
2	1.501	0.5	neg
3	0.090	7.6	pos
4	0.123	5.6	pos
5	0.120	5.7	pos

**Performance Panel: PHT 201**

Sample	OD450nm	CoS	Diagnosis	Sample	OD450nm	CoS	Diagnosis
1	0.189	4.0	pos	14	0.139	4.9	pos
2	0.132	5.2	pos	15	0.115	5.9	pos
3	0.143	4.8	pos	16	0.167	4.1	pos
4	0.104	6.6	pos	17	0.086	8.0	pos
5	0.438	1.6	pos	18	0.160	4.3	pos
6	0.121	5.7	pos	19	0.175	3.9	pos
7	0.137	5.4	pos	20	1.772	0.4	neg
8	0.150	4.6	pos	21	0.090	7.6	pos
9	0.115	5.9	pos	22	0.281	3.4	pos
10	0.094	7.3	pos	23	0.221	2.4	pos
11	0.070	9.8	pos	24	0.134	5.1	pos
12	1.814	0.4	neg	25	0.142	4.6	pos
13	0.097	7.1	pos	26	1.860	0.4	neg

**4. Precision**

The mean values obtained from a study conducted on two samples of different anti-HAV reactivity, examined in 16 replicates in three separate runs is reported below:

- Test # 1**
- | Sample         | Negative | Low Pos. |
|----------------|----------|----------|
| OD450nm        | 2.425    | 0.698    |
| Std. Deviation | 0.095    | 0.023    |
| CV %           | 2.7      | 3.9      |
- Test # 2**
- | Sample         | Negative | Low Pos. |
|----------------|----------|----------|
| OD450nm        | 2.373    | 0.573    |
| Std. Deviation | 0.107    | 0.034    |
| CV %           | 4.5      | 6.0      |
- Test # 3**
- | Sample         | Negative | Low Pos. |
|----------------|----------|----------|
| OD450nm        | 2.478    | 0.554    |
| Std. Deviation | 0.108    | 0.023    |
| CV %           | 4.4      | 4.2      |

The variability shown in the tables did not result in sample misclassification.

**S. LIMITATIONS OF THE PROCEDURE**

Bacterial contamination or heat inactivation of the specimen may affect the absorbance values of the samples with consequent alteration of the level of the analysis. This test is suitable only for testing single samples and not pooled ones. Diagnosis of an infectious disease should not be established on the basis of a single test result. The patient's clinical history, serology, as well as other diagnostic data should be considered.

**REFERENCES**

1. CDC. Summary of notifiable diseases - United States - 1997. *MMWR* 1999;46:1-47.
2. CDC. Prevention of hepatitis A through active or passive immunization. Recommendations of the Advisory Committee on Immunization Practices (ACIP). *MMWR* 1996;45:1R-15.
3. Krogan S, Giler JP. Viral hepatitis: new light on an old disease. *JAMA* 1970;222:1919-29.
4. Heller SC, Webster RM, Eron JJ, Swenson JE, Manard JE, Eng J, Wee 1980;302:1222-7.
5. Lethner WM, Lemon SM, Kirkpatrick JW, Reichard RR, Fields ML, Kelley PW. Frequency of illness associated with epidemic hepatitis A virus infection in adults. *Am J Epidemiol* 1985;122:285-33.
6. Gleason M, Galun E, Oren R, Tur-Kaspa R, Strouval D. Relapsing hepatitis A. Review of 14 cases and literature survey. *Medicine* 1992;71:16-25.
7. Shih J, P. Matheson LR, Kruger P, M. Iyer AM. Focal excretion of hepatitis A virus in patients with symptomatic hepatitis A infection. *Scand J Gastroenterol* 1981;16:1057-9.
8. Tsiapras NC, Papapanagiotou GI, Tzechuraj JR, Purcell RH. Focal excretion of Greek strains of hepatitis A virus in patients with hepatitis A and in experimentally infected chimpanzees. *J Infect Dis* 1985;154:2317.
9. Rosenblum LS, Villanar ME, Nelson CV, et al. Hepatitis A outbreak in a community in the United States: evidence of prolonged viral excretion among patients. *J Infect Dis* 1991;164:875-82.
10. Slogren MH, Tammo H, Fay O, et al. Hepatitis A virus in stool during clinical relapse. *Ann Intern Med* 1987;106:221-6.
11. Lemon SM. The natural history of hepatitis A: the potential for transmission by transfusion of blood or blood products. *Vox Sang* 1984;67(suppl 4):19-23.
12. Boer VA, Nelson CV, Menges HS. Duration of viremia in naturally acquired hepatitis A virus infection. *Proceedings of the 45th Annual Meeting of the Infectious Diseases Society of America*, 55th Annual Meeting, Alexandria, VA: Infectious Diseases Society of America, 1997.
13. Lau JF, Yang CY, Chu CM, Hung MJ. Appearance and persistence of hepatitis A IgM antibody in acute clinical hepatitis A observed in an outbreak. *Infection* 1986;14:155-9.
14. Sagarin JT. Host immune response to hepatitis A virus. *J Infect Dis* 1985;177(suppl 1):S9-14.
15. Huh JY, Paul V, Connor EH, et al. A multistate, foodborne outbreak of hepatitis A. *N Engl J Med* 1989;340:283-92.
16. Squire JM, Robertson BH, Bell BP, McCausland KA, East BJ. Hepatitis A virus infections associated with drinking water concentrate in the United States. *Transfusion* 1989;38:573-9.
17. Cohen JL, Faneone S, Purcell RH. Hepatitis A virus infection in a chimpanzee: duration of viremia and detection of virus in saliva and throat swabs. *J Infect Dis* 1989;160:867-90.
18. McCausland KA, Bond WW, Bradley DW, Ebert JW, Maynard JE. Survival of hepatitis A virus in tissues after drying and storage for 7 months. *J Clin Microbiol* 1982;16:557-9.
19. Favero MS, Bond WW. Disinfection and sterilization. In: Zuckerman AJ, Thomas HC, eds. *Viral hepatitis: scientific basis and clinical management*. New York, NY: Churchill Livingstone, 1993:565-75.
20. Sikes C, Scherker T, Rick L, et al. Source of infection among persons with acute hepatitis A and no identified risk factors. Salt Lake County, Utah, 1996 [Abstract 302]. *Clin Infect Dis* 1997;25:411.