


MENDEC[®] Spine

MV SYSTEM

MENDEC[®] Spine

MV System

IT

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
To the medical staff attention

CAUTION

Users should be trained in use of percutaneous vertebral body injection procedures using bone cement for the treatment of vertebral pathologies and in the management of complications resulting from procedures.

Read this document carefully. Failure to follow instructions may cause malfunction and/or cause patient injury.

DESCRIPTION

Mendec Spine MV System is radiopaque acrylic resin (PMMA based) pre-dosed in a sterile and disposable closed vacuum system. The package includes also a tube with luer-lock connections for the connection of the system body to a bone injection needle (not provided).

The acrylic resin consists of 9,2g sterile liquid and 20g sterile powder.

The liquid is sterilised by means of filtration and the powder, the system and the tube by means of ethylene oxide treatment.

The package consists of a Tyvek-sealed blister contained in an aluminium sachet.

Mendec Spine MV System allows:

- the mixing/delivery of the resin;
- the injection of the resin through a bone injection needle;
- preparation of the resin in a system totally closed allowing:
 - the limitation of the risk of microbiological contamination;
 - the reduction of the inhalation of monomer fumes by the operator, as well as the reduction of the possible mixing mistakes.
- preparation of the resin in a vacuum mixing chamber.

COMPOSITION of the acrylic resin

Liquid component	
Methylmethacrylate	99,10% w/w
N,N-dimethyl-p-toluidine	0,90% w/w
Hydroquinone	75 ppm
Powder component	
Polymethylmethacrylate	69,00% w/w
Barium sulphate	29,00% w/w
Benzoyl peroxide	2,00% w/w

INDICATIONS

Mendec Spine MV System is indicated for the filling of pathologic vertebral body in relation to osteoporotic compression fractures, metastasis, myeloma (Vertebroplasty).

CONTRAINDICATIONS

The use of Mendec Spine MV System is contraindicated in the following conditions:

- presence of active infections;
- hemorrhagic diathesis;
- extended vertebral destruction, with collapse superior to 2/3 of the standard thickness, destruction of the posterior wall with epidural extension of the pathologic tissue and clinical signs of medullar compression
- hypersensitivity to monomer or to any of the other components of the resin.

PRECAUTIONS FOR USE

The manufacturer does not provide any information on the use of the device in pregnancy or in children since no clinical data are available on this patient population; the use of the device in these cases requires that the potential benefits be weighed against the possible hazards by the surgeon.

For safe and effective use of Mendec Spine MV System, the surgeon should have specific training and experience to be thoroughly familiar with the properties, its preparation, handling characteristics, application limitations and correct application technique of the device.

Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events.

The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures. The product should not be used after the expiration date. As the liquid monomer is highly volatile and flammable, the operating room should be provided with adequate ventilation to eliminate the maximum amount of monomer vapour. The liquid component is a powerful lipid solvent. Avoid monomer contact with the skin and mucous membranes. It has caused contact dermatitis in susceptible individuals. Wearing of a second pair of surgical gloves and strict adherence to the mixing instructions may diminish the possibility of hypersensitive reactions. The liquid component is absorbed by rubber, that is why it should not be allowed to come in direct contact with rubber including rubber gloves during mixing. The vapours of the liquid component may affect contact lenses. Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient. Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the device's morphology and the mechanical features, causing a malfunction of the same with serious risks for the patient's health. All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation. Follow handling, mixing and preparation instructions with care. It is important to maintain strictly aseptic surgical procedures. It is important to note that deep wound infection will present a serious risk for the successful outcome of the surgical operation.

INSTRUCTIONS FOR USE

The acrylic resin is prepared by mixing the liquid and the powder included in the system.

Acrylic resins are temperature sensitive. Any increase in temperature of the working environment, the product and the mixing materials over 23°C reduces the waiting and hardening times of the acrylic resin. Equally, lower temperatures increase such times.

Before using Mendec Spine MV System it is strongly advised to make sure that the package is stored at a temperature of 23°C ± 1°C for the previous 24 hours.

NEVER add other substances or foreign bodies to the acrylic resin.

Ensure that the inner package is not damaged and that the components are integer, i.e. the powder should not present yellow or brown discolouring and the liquid should not be syrupy.

Preparation

For a correct preparation of Mendec Spine MV System the following steps have to be followed:

1

Grip the device on the vial holder (B) side and loosen the powder by tapping the mixing chamber (C) against the palm of the hand.

2

With the vial holder (B) facing upwards, push the thrust knob (A) smartly until the vial containing the monomer breaks.


3

Pull the thrust knob (A) upwards and then fully push it down. Repeat the cycle until the monomer has completely transferred into the mixing chamber (C).

4

When the monomer transfer is completed, lock the device by pushing the thrust knob (A) fully down and screw it in a clockwise direction for two complete turns.


5

Upturn the device and hold it in a vertical position with the vial holder (B) facing downwards. Start mixing by repeatedly pushing the mixing knob (E) down and then pulling it up, for about 1 minute.

Warning: do not rotate the mixing knob (E) while the mixing rod (D) is in the maximum upwards extension position, in order to prevent the mixing rod (D) from being removed before the correct time.

6

When mixing is complete, bring the mixing rod (D) into the maximum extension upwards and rotate the mixing knob (E) in the anticlockwise direction until the mixing rod (D) is unscrewed.

7

While keeping the device in the vertical position, remove the air from the mixing chamber (C) by screwing and slowly pushing the thrust knob (A) in the clockwise direction until the resin reaches the mouth of the device.

8

Connect the device to a needle via the luer-lock connector using the provided connection tube (F). The extrusion of the acrylic resin must occur within **5 minutes** from the end of the mixing procedure. The necessary amount of acrylic resin for the particular clinical application is defined by the surgeon once the components have been mixed. A graduated scale is present on the syringe in order to visualize the amount of the injected acrylic resin.

Application

Once the acrylic resin has reached the desired viscosity according to the needle used (check the viscosity by extruding the cement from the connection tube, and eliminate the first few ml of extruded resin - purging), the injection may start.

It is not advisable to use the product with needles which are thinner than 13G.

The injection has to be executed under continuous radioscopic guidance and must be done slowly (max injection speed 2ml/ min), checking the exit of the resin and its distribution inside the vertebral body, and any possible extra-vertebral diffusion (leakage). If any resin leakage is detected, immediately interrupt the injection and unscrew the push-knob.

Caution

Do not use Mendec Spine MV System as a delivery device to inject the resin directly in the vertebral bodies. Mendec Spine MV System must be connected with a proper bone injection needle with a luer-lock connection, using the provided connection tube.

Polymerization of the acrylic resin is an exothermic reaction and reach temperatures higher than physiological temperatures. Make sure that the preparation accessories and bone needles are specifically compatible with the product.

ADVERSE EFFECTS

Listed below are the most serious and frequent adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty. The surgeon should be aware of these reactions and be prepared to treat such reactions if they are encountered.

Associated to acrylic resin application	Procedural adverse events
Transient hypotension	Infection
BCIS	Local trauma (injury to pleura, kidneys, nerve roots, spinal cord, lamina and pedicle fracture)
Allergic reaction	Rib fracture
Vertebral fracture	Bleeding
Radiculopathy	Hemorrhage, Hematoma
Spinal cord compression	Short term cardiac irregularity
Cardiac perforation	Myocardial infarction
Pulmonary and cardiac embolism	Sudden death
Paraplegia, tetraplegia	
Cerebrovascular accident	
Death	

Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.

Other reported adverse events for acrylic bone cements intended for vertebroplasty include:

- Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.

IMPORTANT PHYSICIAN'S INFORMATION

The use of bone cement requires extremely close cooperation between the surgeon and the anaesthetist. The anaesthetist must be told when the bone cement is implanted during the operation. In some cases, events defined as bone cement implantation syndrome (BCIS) may occur, which are characterised by a series of clinical characteristics that include hypoxia, hypotension, cardiac arrhythmia, increased pulmonary vascular resistance and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology. Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the perivertebral venous system and these may migrate to the right side of the heart and into pulmonary circulation (pulmonary embolism) resulting in BCIS. BCIS has a very wide spectrum of severity and in many cases patients with pulmonary embolism detected during a chest X-ray remain asymptomatic.

It is recommended to verify the proper filling of the treated vertebra in order to avoid a further vertebral fracture. In osteoporotic patients the risk of further vertebral fractures is a possible event.


Storage

Store below 25°C in a dry environment and protect from light.

Sterility is guaranteed only if the unit container is undamaged.

TEMPERATURE EFFECTS ON PREPARATION AND APPLICATION TIMING MENDEC SPINE MV SYSTEM

The temperature of the storage and of the operative area affects the preparation and application of the resin. Temperature's effect on cement setting time was evaluated with a laboratory test. A graph on setting time according to temperature is shown below for ease of reference. The data provided was obtained in controlled environmental and storage conditions subjected to standard deviation.

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TECRES S.p.A.