

BACTICID AF WIPES

ALCOHOL BASED RAPID DISINFECTION WIPES FOR MEDICAL DEVICES

- ready-to-use disinfection wipes pre-soaked with Bactacid AF
- aldehyde, colourant and fragrance free
- rapid disinfection with broad biocidal activity
- effective against Noro, Adeno and Polio viruses
- residue free



Bactacid AF wipes are ethanol-based rapid disinfection wipes for non-invasive medical devices and alcohol resistant surfaces. Disposable wipes are pre-soaked with Bactacid AF solution and come in a convenient and ready to use packaging. Perfect product in conditions with higher hygiene requirements, where rapid and regular surface disinfection is significant. It does not contain aldehydes, fragrance and colorants and is suitable for all alcohol-resistant surfaces.

Bactacid AF wipes have a broad spectrum of bactericidal (incl. MRSA, VRE), mycobactericidal (incl. Mycobacterium avium and Mycobacterium terrae), fully virucidal (incl. all enveloped viruses, Rota, Adeno, Noro and Polio), fungicidal and yeasticidal activity.

Bactacid AF have been dermatologically tested.

APPLICATION AREAS:

Alcohol-based rapid disinfection wipes for medical devices and other alcohol resistant surfaces, where short contact times are required, for example in the following areas:

- Hospitals and other medical institutions
- Operating theatres and procedure rooms
- Intensive care unit
- In accident and emergency departments
- Dialysis departments
- Family doctor offices, first-aid stations
- Ambulances
- Laboratories and diagnostics.

- Dentist's offices
- For home care of patient
- Food industry and food processing areas
- Pharmaceutical and cosmetics industry

USAGE:

BACTICID AF wipes are ready for use. Open the cover of the bucket and pull out the impregnated wipe up to the perforation. To tear off the wipe, pull it suddenly and at a small angle. After tearing off the wipe, close the cover immediately to avoid drying of the wipes.

Thoroughly wipe the surface with the impregnated wipe, ensuring complete coverage. The surfaces must remain visibly moistened for 30 seconds to take effect.

Rinsing of the surfaces is not required. Mark the date of the first use on the product packaging.

Bactacid AF wipes is a perfect product for rapid disinfection of surfaces like patient nearby surfaces, operating tables, examination couches, stretchers, patient chairs, wheelchairs, medical equipment, trays, trolleys, hospital and maternity beds.

If the surface tolerance is unknown, try it out on an inconspicuous area first. This product is not suitable for acrylic glass, plexiglas and alcohol soluble varnish.

MICROBIOLOGY:

Efficacy	Concentration	Contact time
Escherichia coli	Ready to use	30 sec
Pseudomonas aeruginosa	Ready to use	30 sec
Staphylococcus aureus	Ready to use	30 sec
Enterococcus hirae	Ready to use	30 sec
Staphylococcus aureus MRSA	Ready to use	30 sec
Vancomycin-Resistant Enterococci (VRE)	Ready to use	30 sec
Candida albicans	Ready to use	30 sec
Mycobacterium terrae	Ready to use	30 sec
Mycobacterium avium	Ready to use	30 sec
All enveloped viruses (incl. HIV, HBV, HCV, Vaccinia)	Ready to use	15 sec
Rotavirus	Ready to use	15 sec
Norovirus	Ready to use	15 sec
Adenovirus	Ready to use	30 sec
Poliovirus	Ready to use	120 sec

BACTICID AF wipes are tested according to the following European standards: EN 13624, EN 13697, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, EN16615.



INGREDIENTS:

100 g of product contains active ingredients: ethanol 57 g, isopropanol 6 g.

PRECAUTIONS:

Highly flammable liquid and vapour


Keep out of reach of children. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with regulations on hazardous waste or packaging and packaging waste respectively.

pH 6,5 – 7,5

PACKAGING:



100 pcs bucket
16x20cm wipes



150 pcs bucket
20x20cm wipes



200 pcs bucket
20x20cm wipes



600 pcs bucket
18x20cm wipes