# Medtronic

## Reliant™

Stent Graft Balloon Catheter  
Балонен катетър за стент-графт  
Balónkový katétr pro stentgraft  
Stentgraftballonkateter  
Gefäßprothesen-Ballonkatheter  
Καθετήρας μπαλονιού μοσχεύματος stent  
Catéter balón para endoprótesis  
Stentsiiriku balloonkateeter  
Stentisiirrepallokatetri  
Cathéter à ballonnet d'endoprothèse  
Balonski kateter za stent-graftove  
Sztentgraft-ballonkatéter  
Kateter Balon Cangkok Stent  
Catetere a palloncino per endoprotesi  
Стент-графтқа арналған баллондық катетер  
스텐트 그래프트 풍선 카테터  
Stento transplantato balioninis kateteris  
Stentprotēzes balonkatetrs  
Балон катетер за стент-графтови  
Stentgraftballonkatheter  
Ballongkateter for stentgraft  
Cewnik balonowy do modelowania stentgraftu  
Cateter de balão para endoprótese  
Cateter cu balon pentru grefe de tip stent  
Баллонный катетер для стент-графта  
Balónový katéter pre stentgraft  
Balonski kateter pokrite žilne opornice (stent-graft)  
Balon-kateter za stent graftove  
Ballongkateter för stentgraft  
Stent Graft Balonlu Kateteri  
Балонний катетер для стент-графта  
Ống thông bóng nong giá đỡ nội mạch

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
## 1 Device description

The Reliant stent graft balloon catheter is designed for temporary occlusion of large vessels during an endovascular procedure and to assist in the expansion of self-expanding stent grafts used for the treatment of abdominal aortic aneurysms (AAA) and thoracic aortic aneurysms (TAA). Ballooning of a stent graft may improve modeling of the graft material and fixation of the stent graft to the vessel wall. Suboptimal expansion of self-expanding stent grafts may also be improved by inflating the balloon at the site of the stent.

**Note:** Read the entire Reliant stent graft balloon catheter Instructions for Use prior to performing a ballooning procedure.

**Warning:** Only physicians and teams trained in vascular surgery, interventional radiology, or cardiology, and who have completed training or have experience with stent grafts and balloon catheters should consider using this device.

**Figure 1.** Reliant Stent Graft Balloon Catheter

- 1 Extension with 3-Way Stopcock
- 2 Y-Connector
- 3 Strain Relief
- 4 Guidewire Lumen

- 5 Balloon Marker Bands
- 6 Catheter
- 7 Balloon

The Reliant stent graft balloon catheter consists of the following components:

- Multi-lumen placement catheter
- Compliant polyurethane balloon
- Radiopaque marker bands
- Strain relief
- Extension with 3-way stopcock
- Y-Connector (backend hub)

The Reliant stent graft balloon catheter has a compliant polyurethane balloon with a maximum diameter of 46 mm. The usable length for this device is 100 cm. This device is designed to accommodate a 0.97 mm (0.038 in) diameter or smaller guidewire. Two radiopaque marker bands are placed within the balloon to facilitate balloon placement, prior to inflation.

## 2 Intended purpose

The Reliant stent graft balloon catheter is intended to temporarily occlude large vessels to control hemorrhage and to improve endovascular stent graft function by aiding in graft expansion and enhancing graft-to-vessel and graft-to-graft interfaces.

### 2.1 Intended users

The Reliant stent graft balloon catheter should only be used by physicians and teams trained in vascular interventional techniques. The device should only be used in a clinical setting in a sterile environment.

### 2.2 Indications for use

The Reliant stent graft balloon catheter is indicated for use in patients suffering from aortic injury or rupture in order to control hemorrhage and in patients undergoing endovascular repair of aortic or aortoiliac injuries or lesions where self-expanding stent grafts are used.

The device is intended to assist in the expansion of self-expanding stent grafts.

### 2.3 Contraindications

The Reliant stent graft balloon catheter is contraindicated in patients who have the following traits:

- Are contraindicated to contrast media or anticoagulants
- Have an arterial entry site that cannot accommodate a 12 Fr (4 mm) introducer
- Are less than 18 years old
- Are pregnant or nursing

## 2.4 Intended patient population

The Reliant stent graft balloon catheter is intended for any patient undergoing a vascular procedure in which there is a need to temporarily occlude large vessels or to expand a vascular prosthetic device in endovascular procedures.

This device has not been tested in a pediatric population.

## 2.5 Clinical benefits

The Reliant stent graft balloon catheter can temporarily occlude large vessels to control hemorrhage from an injury or ruptured aneurysm; it can also be used to improve endovascular stent graft function by aiding in graft expansion and enhancing graft-to-vessel and graft-to-graft interfaces.

**Note:** After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00008788T.

## 2.6 Device performance characteristics

The Reliant stent graft balloon catheter is a single-piece device intended to temporarily occlude large vessels during an endovascular procedure and to improve the sealing performance of aortic stent grafts. This is achieved by tracking the device over a guidewire with a diameter of 0.97 mm (0.038 in) or less to the desired location using radiopaque markers located in the proximal and distal ends of the balloon. The balloon can then be inflated by injecting heparinized saline into the device with a syringe until the balloon has expanded to the required size (see *Section 8.4* for the balloon compliance chart). By inflating the stent graft at the proximal, distal, or overlapped areas, it may improve modelling of the graft material and fixation of the stent graft to the vessel wall.

## 3 Warnings and precautions

**Caution:** Federal (USA) law restricts this device to sale by or on order of a physician.

- Preparations should be made and a trained vascular surgical team available in the event that conversion to open surgery is required.
- Carefully monitor the patient's blood pressure throughout the procedure.
- If an obstruction in the vessel (eg, a tortuous bend, stenosis, calcification) prevents advancement of the catheter, use standard techniques to dilate or straighten the vessel before continuing to advance the catheter.
- Balloon rupture may occur under certain anatomical, procedural, or clinical circumstances so it is recommended that backup Reliant stent graft balloon catheters are available.
- Do not overinflate the balloon when modeling graft in vessels. Operator should visualize the stent graft at all times during balloon inflation to detect any stent graft movement.
- The Reliant stent graft balloon catheter is not intended for use as an angioplasty balloon.
- Do not use Reliant stent graft balloon catheter for more than 20 inflation or deflation cycles.
- **If the proximal and distal radiopaque markers of the Reliant device are not completely within the covered stent graft, there is an increased risk of vessel injury, rupture, or possible patient death.**
- The device is not for use in endovascular prostheses contraindicated for use with compliant balloon catheters.
- Do not use the Reliant stent graft balloon catheter in treatment of aortic dissection disease.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.
- Increased manipulation or duration of the procedure may lead to microembolization.
- Do not exceed the maximum inflation diameter (46 mm) because the balloon could rupture.
- **Overinflation of a balloon can cause damage to the vessel wall or stent graft, which can lead to vessel rupture.** Care should be taken when inflating a balloon in vessel, particularly when inflating in the distal-most spring, calcified, stenotic, or otherwise diseased vessels.
- **Not for use in the heart or coronary arteries.**
- Do not use a high-pressure injector with the Reliant stent graft balloon catheter. Hand injections are recommended.

- Do not use a power injector for injection of contrast medium through the distal catheter lumen because rupture may occur.
- Medtronic is not responsible for mis-sizing, misuse, or misplacement of the device.

## 4 Potential Complications and Adverse Events

Complications may occur with any stent graft balloon catheterization procedure. Therefore, only physicians trained in vascular surgery, interventional radiology, or cardiology, and who have completed training or have experience with stent grafts and balloon catheters should consider using this device. Possible complications associated with this type of procedure include, but are not limited to:

### Clinically Related:

- Vessel perforation or dissection
- Paresthesia
- Entry site infection
- Entry site hematoma
- Stroke
- Aneurysm rupture
- Cardiac events
- Respiratory failure
- General malaise
- Emboli
- Hemorrhage
- Renal complications
- Death

### Device Related:

- Balloon catching inside of graft during catheter movement
- Stent graft migration caused by balloon catheter device
- Balloon rupture
- Inability to inflate/deflate balloon
- Inability to insert guidewire

**Note:** If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority.

## 5 Patient Counseling

There are risks involved with any medical procedure. Both physician and patient should fully understand the risks and benefits associated with this endovascular device and procedure.

## 6 How supplied

The device is sterilized using ethylene oxide.

## 7 Additional equipment recommended

- Diluted contrast (75% sodium chloride/25% Renografin) for balloon inflation or deflation
- Heparinized saline solution
- Backup Reliant stent graft balloon catheters
- Freely-angled c-arm with:
  - High resolution fluoroscopy
  - High quality angiography
  - Digital subtraction angiography (DSA)
- Surgical suite standby, in the event that emergency surgery is necessary
- Inventory of introducer sheaths/dilators, syringes, and guidewires

## 8 Procedural Instructions

In accordance with standard technique, it is recommended that patients be anticoagulated during the procedure.

### 8.1 Preparation

**Note:** Entry site preparation to be performed prior to device preparation.

1. Attach syringe to guidewire lumen.
2. Use the syringe to flush the guidewire lumen with heparinized saline solution.
3. Remove the syringe from the guidewire lumen and connect it to the opened stopcock on balloon lumen. Draw a vacuum in the balloon and close stopcock.
4. Gently remove the protective sleeve from the balloon.
5. Fill syringe with heparinized saline solution and open stopcock.

6. Hold catheter with distal tip and balloon pointing down.
7. Partially inflate balloon.
8. Draw back on the syringe to deflate the balloon.
9. Repeat step 7 and step 8. Each time this is repeated, more volume of air is displaced with liquid. Attempt to remove all of the air in the balloon. Some changes in the catheter orientation may be necessary to vent all the air.
10. Draw a vacuum in the balloon and close stopcock.
11. Insert guidewire and advance it beyond the stent graft.
12. Place 12 Fr (4 mm) introducer/dilator in the femoral artery, over the guidewire.
13. Remove dilator from introducer.

## 8.2 Introduction

1. Verify that the balloon is deflated.
2. Slowly insert the Reliant stent graft balloon catheter over the guidewire.

**Caution:** When aligning the position of the catheter so that the balloon is in the proper position for expansion within the vessel, pay careful attention to the fluoroscope location to avoid parallax or other sources of visualization error.

## 8.3 Confirm position

Use fluoroscopy and angiography to ensure that the balloon is at the target expansion site. Radiopaque markers located near the bonds of the balloon should be placed just proximal and distal of the target location.

## 8.4 Balloon Inflation/Deflation

1. When the Reliant stent graft balloon catheter is positioned properly, it may be appropriate to decrease the patient's blood pressure to approximately 80 mmHg to facilitate ballooning.  
**Caution:** Failure to lower the blood pressure during thoracic cases may cause balloon movement and may result in stent graft misplacement.
2. Inflate the balloon to expand the target location. Continuously monitor fluoroscope screen, watching for stent movement. Proper modeling should show very slight outward expansion of stent with balloon inflation. Do not overinflate. Stop after stent graft expansion is evident. The table below is a guideline for determining the volume of solution (75% sodium chloride/25% Renografin) required to obtain a given balloon expansion diameter:

46 mm Balloon	
Diameter	mL (cc)
10 mm	3
20 mm	9
30 mm	19
40 mm	41
46 mm	60

**Caution:** This chart is only a guide. Balloon expansion should be carefully monitored with the use of fluoroscopy.

**Caution:** Do not exceed the maximum inflation diameter (46 mm) because the balloon could rupture. Overinflation can also cause damage to the vessel wall or stent graft or lead to vessel rupture.

**Note:** The balloon has been tested to temporarily occlude vessels up to 37 mm in diameter.

3. Fully deflate balloon and then verify it is deflated.
4. Move the balloon distally; inflate and deflate it to smooth wrinkles in the stent graft material and to enhance contact between the stent graft and vessel wall.  
**Warning:** Ensure the balloon is fully deflated before moving the Reliant stent graft balloon catheter.
5. Repeat until all target areas of the stent graft have been modeled. Inflate the balloon in the distal spring area with sufficient pressure to firmly embed the spring against the vessel.  
**Caution:** Balloon is highly compliant. Inflate slowly and do not overexpand. Use special care in areas of diseased vessels to avoid rupture or vessel trauma

## 8.5 Remove Catheter and Seal Entry Site

1. Draw a vacuum in the balloon. Verify that the balloon is fully deflated.
2. Withdraw the Reliant stent graft balloon catheter back through the introducer sheath.
3. Gently remove the Reliant stent graft balloon catheter and guidewire. Use fluoroscopy to ensure that the stent graft does not move and to track the movement of the Reliant stent graft balloon catheter throughout the withdrawal.

**Note:** Close entry site with standard surgical closure technique.

## **9 Disclaimer of warranty**

**Note: This disclaimer of warranty does not apply in Australia.**

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
**Medtronic, Inc.**
Medtronic B.V.
Medtronic Limited