# Progesterone

For use on IMMULITE® 2000 systems

**SIEMENS**

# IMMULITE® 2000 Progesterone

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of progesterone in serum, as an aid in the diagnosis and treatment of disorders of the ovaries or placenta.

Catalog Numbers: **L2KPW2** (200 tests), **L2KPW6** (600 tests)

Test Code: **PRG**

Color Code: **Light Green**

### Summary and Explanation

Progesterone is a steroid hormone which plays an important role in the preparation for and maintenance of pregnancy. It is synthesized from cholesterol via pregnenolone — then rapidly metabolized to pregnanediol, for the most part, in the liver.<sup>1,8,12</sup> The ovary and placenta are the major production sites; but a small amount is also synthesized by the adrenal cortex in both men and women.

Circulating progesterone levels, which are characteristically low during the follicular phase, increase sharply during the luteal phase of menstrual cycles, reaching a maximum some 5 to 10 days after the midcycle LH peak.<sup>11</sup> Unless pregnancy occurs, a steep decline to follicular levels sets in about 4 days before the next menstrual period. This pattern constitutes the rationale behind the well-established use of serum progesterone measurements as a simple and reliable method for ovulation detection.<sup>2,3,15</sup>

There is growing literature on luteal phase defects.<sup>4,7,10,16,17</sup> Daily progesterone levels are considered the most accurate means for documenting a defective luteal phase.<sup>4,7</sup> However, some investigators have found that three samples or even a single sample<sup>13,17</sup> (if well timed) can provide valuable information on the adequacy of the luteal phase.

Measurements of serum progesterone have also been used to check the effectiveness of ovulation induction,<sup>15</sup> to monitor progesterone replacement

therapy<sup>15</sup> and to detect and evaluate patients at risk for abortion during the early weeks of pregnancy.<sup>5,6,14</sup> On the other hand, although progesterone levels increase throughout pregnancy, they are not considered a suitable means for monitoring fetal well-being during the third trimester.<sup>8</sup>

### Principle of the Procedure

IMMULITE 2000 Progesterone is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 30 minutes

### Specimen Collection

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Progesterone has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**EDTA:** Because EDTA would have a significant effect on results, it should not be used as an anticoagulant.

**Gel Barrier Tubes:** Time-dependent decreases in progesterone levels have been reported when serum samples are collected and stored in gel barrier tubes.<sup>19,20,21</sup>

**Lipemia:** Gross lipemia has been shown to cause a decrease in the apparent progesterone concentration. The use of an ultracentrifuge is recommended to clear lipemic samples.

**Volume Required:** 25 µL serum

**Storage:** 7 days at 2–8°C, or 3 months at –20°C.<sup>18</sup>

**Dilution of High Samples:** All samples expected to have levels above the assay's calibration range should be diluted on-board.

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>26-28</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Progesterone Bead Pack (L2PW12)

With barcode. 200 beads, coated with polyclonal rabbit anti-progesterone. Stable at 2–8°C until expiration date.

**L2KPW2:** 1 pack    **L2KPW6:** 3 packs

### Progesterone Reagent Wedge (L2PWA2)

With barcode. 21 mL alkaline phosphatase (bovine calf intestine) conjugated to progesterone, in buffer, dispensed equally into chambers A and B. Stable at 2–8°C until expiration date.

**L2KPW2:** 1 wedge    **L2KPW6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Progesterone Adjustors (LPWL, LPWH)

Two vials (Low and High), 3 mL each, of progesterone in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KPW2:** 1 set    **L2KPW6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Multi-Diluent 1 (L2M1Z)

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) processed, normal human serum, containing undetectable to low levels of progesterone, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2M1Z:** 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M1Z:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

### Recommended Adjustment Interval:

2 weeks

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of progesterone.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

A multi-national study involving women in apparent good health (age: 16–44 years), who volunteered to have blood samples drawn, on a daily basis, throughout one complete ovulatory cycle,<sup>23</sup> yielded the following results.

Progesterone, ng/mL

Ovulatory Cycles	n*	Median	Central 95%
Follicular Phase	27 (382)	0.47	ND – 1.13
Midfollicular Days 5 to 11	27 (186)	0.43	ND – 0.98
Midcycle	27 (27)	1.06	0.48 – 1.72
Luteal Phase	27 (323)	8.9	0.95 – 21
Midluteal, Days 7 to 8 of Luteal Phase	27 (54)	13.1	6.0 – 24

\*Number of subjects (total number of results)

Progesterone, nmol/L

Ovulatory Cycles	n*	Median	Central 95%
Follicular Phase	27 (382)	1.5	ND – 3.6
Midfollicular Days 5 to 11	27 (186)	1.4	ND – 3.1
Midcycle	27 (27)	3.4	1.5 – 5.5
Luteal Phase	27 (323)	28	3.0 – 68
Midluteal, Days 7 to 8 of Luteal Phase	27 (54)	42	19 – 76

\*Number of subjects (total number of results)

See Menstrual Cycle Graph (see Tables and Graphs).

Another study yielded the following results.

Mass Units (ng/mL)	Absolute		n
	Median	Range	
<b>Males:</b>	0.52	0.27 – 0.90	63
<b>Females:</b>			
Follicular Phase	0.67	0.33 – 1.2	29
Luteal Phase	4.8	0.72 – 17.8	29
Post-menopausal	0.36	ND – 1.0	34
Oral Contraceptives	0.70	0.34 – 0.92	19
<b>Pregnant Females:</b>			
First Trimester	22.2	9.3 – 33.2	28
Second Trimester	35.4	29.5 – 50.0	10
Third Trimester	102	83.1 – 160	8

ND: not detectable

S.I. Units (nmol/L)	Absolute		n
	Median	Range	
<b>Males:</b>	1.7	0.86 – 2.9	63
<b>Females:</b>			
Follicular Phase	2.1	1.0 – 3.8	29
Luteal Phase	15.3	2.3 – 56.6	29
Post-menopausal	1.1	ND – 3.2	34
Oral Contraceptives	2.2	1.1 – 2.9	19
<b>Pregnant Females:</b>			
First Trimester	70.6	29.6 – 106	28
Second Trimester	113	93.8 – 159	10
Third Trimester	324	264 – 509	8

ND: not detectable

In pregnancy, the general tendency is for values to increase. There is considerable interpersonal variation in progesterone values, particularly in groups associated with elevated levels. (Note that measurement of progesterone levels is generally considered unsuitable for monitoring fetal well-being in the later weeks of pregnancy.<sup>8)</sup>)

A cross-sectional study of pediatric fertility values at a “wellness” clinic in the southwestern United States yielded the following results.

Group	Age (yr)	n	Progesterone, ng/mL	
			Median	Central 95%
Females	Cord	27	570	465 – 755
	0.1 – 0.4	24	1.2	0.25 – 17
	0.5 – 1	19	0.8	0.2 – 1.6
	1.1 – 9	38	0.4	ND – 1.4
Males	Cord	27	520	345 – 650
	0.1 – 0.4	33	1.5	0.3 – 14
	0.5 – 1	14	0.8	ND – 2
	1.1 – 9	42	0.4	ND – 1.3
Combined	Cord	54	550	350 – 750
	0.1 – 0.4	57	1.5	0.25 – 17
	0.5 – 1	33	0.8	ND – 2
	1.1 – 9	80	0.4	ND – 1.3

ND: nondetectable

Group	Age (yr)	n	Progesterone, nmol/L	
			Median	Central 95%
Females	Cord	27	1813	1479 – 2401
	0.1 – 0.4	24	3.8	0.8 – 54
	0.5 – 1	19	2.5	0.6 – 5.1
Males	1.1 – 9	38	1.3	ND – 4.5
	Cord	27	1654	1097 – 2067
	0.1 – 0.4	33	4.8	1.0 – 45
Combined	0.5 – 1	14	2.5	ND – 6.4
	1.1 – 9	42	1.3	ND – 4.1
	Cord	54	1749	1113 – 2385
Males	0.1 – 0.4	57	4.8	0.8 – 54
	0.5 – 1	33	2.5	ND – 6.4
	1.1 – 9	80	1.3	ND – 4.1

ND: nondetectable

**Pediatric:** Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Progesterone assay in accordance with CLSI guideline EP28-A3C.<sup>24</sup> Samples were collected prospectively from apparently healthy pediatric subjects, using predefined inclusion criteria. Reference values were generated for subpopulations based on age and Tanner stage subgroups based on physiological development. The study was designed to establish reference values for females within each age or Tanner stage subgroup. The subject's Tanner stage was assessed based on pubic hair and genitalia/breast development. The scale proposed by Neinstein and Kaufman was used for the determination of the Tanner stages.<sup>25</sup>

The reference intervals and Tanner values are based on the central 90% (5th and 95th percentiles). Where sample sizes were insufficient to calculate the 5th or 95th percentile, the minimum or maximum observed values are presented in the Reference Intervals tables.

### IMMULITE 2000/2000 XPI Progesterone Pediatric Reference Intervals

Age (Years)	n	Median	Range
		ng/mL	
12	38	< 0.2	< 0.2 – 1.4
13 – 21	127	0.4	< 0.2 – 10.0

Age (Years)	n	Median	Range
		nmol/L	
12	38	< 0.6	< 0.6 – 4.3
13 – 21	127	1.2	< 0.6 – 31.7

### IMMULITE 2000/2000 XPI Progesterone Pediatric Reference Intervals by Tanner Stage

Tanner Stage	n	Median	Range
		ng/mL	
1	1	Not Applicable	
2	18	< 0.2	< 0.2* – > 6.3 <sup>†</sup>
3	49	0.2	< 0.2 – 8.3
4	45	0.3	< 0.2 – 7.0
5	52	0.6	< 0.2 – 11.3

Tanner Stage	n	Median	Range
		nmol/L	
1	1	Not Applicable	
2	18	< 0.6	< 0.6* – > 20.1 <sup>†</sup>
3	49	0.6	< 0.6 – 26.3
4	45	0.9	< 0.6 – 22.2
5	52	2.0	< 0.6 – 35.8

\* Value presented is the minimum reportable value observed; insufficient sample size to calculate a 5th percentile limit.

<sup>†</sup> Value presented is the maximum value observed; insufficient sample size to calculate a 95th percentile limit.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

### Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of

interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

The supplement Dehydroepiandrosterone (DHEA) may cause falsely elevated progesterone results in immunoassays. At an initial concentration of Progesterone of 1.16 ng/mL, an 80.2% change in concentration was observed at the supraphysiologic level of DHEAS (DHEA metabolite) of 20,000 ng/mL. For patients being treated with DHEA, an alternate method that is not expected to show cross reactivity to DHEAS (DHEA metabolite), such as Liquid Chromatography-Mass Spectrometry (LC-MS), should be used.

With the advent of new steroid based medications (analogues) with similar chemical structures to progesterone, there is the possibility of cross-reactivity and falsely elevated results. For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. If the progesterone results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

### Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as ng/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

#### Conversion Factor:

ng/mL × 3.18 → nmol/L

**Calibration Range:** 0.2 to 40 ng/mL (0.6 to 127 nmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 0.1 ng/mL (0.3 nmol/L)

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See “Precision” table.)

**Linearity:** Samples were assayed under various dilutions (See “Linearity” table for representative data.)

**Recovery:** Samples spiked 1-to-19 with three progesterone solutions (25, 100 and 200 ng/mL) were assayed. (See “Recovery” table for representative data.)

**Specificity:** The assay is highly specific for progesterone (See “Specificity” table.)

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides interferes with the assay, causing a depression of values. (See “Lipemia” table.)

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 19 volunteers into plain glass and plastic serum tubes.

(Serum Plastic) = 1.03 (Serum Glass) – 0.06 ng/mL  
r = 0.997

Means:  
11.2 ng/mL (Serum Glass)  
11.5 ng/mL (Serum Plastic)

In another study, blood was collected from 20 volunteers into plain, heparinized, EDTA and Becton Dickinson SST<sup>®</sup> vacutainer glass tubes. Equal volumes of the matched samples were spiked with various concentrations of progesterone, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Progesterone procedure.

(Heparin) = 0.84 (Serum) + 0.4 ng/mL  
r = 0.974

(EDTA) = 2.21 (Serum) + 4.8 ng/mL  
r = 0.893

(SST) = 1.04 (Plain Tubes) + 0.02 ng/mL  
r = 0.987

Means:  
9.6 ng/mL (Serum)  
8.5 ng/mL (Heparin)

21.2 ng/mL (EDTA)  
10.2 ng/mL (SST)

**Method Comparison 1:** The assay was compared to Coat-A-Count<sup>®</sup> Progesterone on 162 patient samples. (Concentration range: approximately 0.2 to 40 ng/mL. See Method Comparison 1 graph.)  
By linear regression:

(IML 2000) = 0.71 (CAC) + 0.33 ng/mL  
r = 0.979

Means:  
8.2 ng/mL (IMMULITE 2000)  
11.1 ng/mL (Coat-A-Count)

**Method Comparison 2:** IMMULITE 2000 Progesterone (L2KPW) was also compared to IMMULITE/IMMULITE 1000 Progesterone (LKPG) on 182 patient samples (Concentration range: approximately 0.2 to 20 ng/mL. See “Method Comparison 2” graph.)  
By linear regression:

(IML 2000 – L2KPW) = 0.79 (IML – LKPG)  
+ 0.82 ng/mL  
r = 0.969

Means:  
5.9 ng/mL (IMMULITE 2000 – L2KPW)  
5.5 ng/mL (IMMULITE – LKPG)

## References

- 1) Aufrere MB, Benson H. Progesterone: an overview and recent advances. J Pharm Sci 1976; 65:783-800.
- 2) Bauman J. Basal body temperature: unreliable method of ovulation detection. Fertil Steril 1981; 36:729-33.
- 3) Brown JB. Timing of ovulation. Med J Austral 1977; ii: 780-83.
- 4) Gautray JP, et al. Clinical investigation of the menstrual cycle: clinical, endometrial and endocrine aspects of luteal defect. Fertil Steril 1981; 35:296-303.
- 5) Hensleigh PA, Fainstat T. Corpus luteum dysfunction: serum progesterone levels in diagnosis and assessment of therapy for recurrent and threatened abortion. Fertil Steril 1979; 32:396-9.
- 6) Hernandez Horta JL, et al. Direct evidence of luteal insufficiency in women with habitual abortion. Obstet Gynecol 1977; 49:705-8.
- 7) Jones G. Luteal phase defects. In: Behrman SJ, Kistner RW, editors. Progress in infertility. Boston: Little & Brown, 2nd Edition, 1975: 299-324.
- 8) Klopper A, Fuchs F. Progestagens. In: Fuchs F, Klopper A, editors. Endocrinology of pregnancy. Hagerstown: Harper & Row, 1977: 99-122.
- 9) Lehmann F, Bettendorf G. The endocrine shift from a normal cycle to anovulation. In: Insler V, Bettendorf G, editors. Advances in diagnosis and treatment of infertility. Amsterdam: Elsevier/North Holland, 1981: 105-13.
- 10) March CM. Luteal phase defects. In: Mishell DR, Davajan V, editors. Reproductive endocrinology, infertility and contraception. Philadelphia: F.A. Davis, 1979: 469-76.
- 11) March CM, Goebelsmann U,

Nakamura RM and Mishell DR. Roles of estradiol and progesterone in eliciting the midcycle luteinizing hormone and follicle-stimulating hormone surges. *J Clin Endocrinol Metab* 1979; 49:507-13. 12) Progesterone (Rochester: Bioeducational Publications, 1981). A BIO-ED slide/seminar educational program.

13) Radwanska E, et al. Plasma progesterone and oestradiol estimations in the diagnosis and treatment of luteal insufficiency in menstruating infertile women. *Acta Eur Fertil* 1976; 7:39-47.

14) Radwanska E, et al. Plasma progesterone levels in normal and abnormal early human pregnancy. *Fertil Steril* 1978; 30:398-402.

15) Radwanska E, et al. Single midluteal progesterone assay in the management of ovulatory infertility. *J Reprod Med* 1981; 26:85-89.

16) Sheehan KL, Casper RF, Yen SSC. Luteal phase defects induced by an agonist of luteinizing hormone-releasing factor: a model for fertility control. *Science* 1982; 215:170-72. 17) Wentz A. Pathophysiology of luteal phase inadequacy. In: Tozzini RI, Reeves G and Pineda RL, editors. *Endocrine physiopathology of the ovary*. Amsterdam: Elsevier/North Holland, 1980, 257-74.

18) Burtis CA, Ashwood ER, editors. *Tietz textbook of clinical chemistry*. 2nd ed. Philadelphia: W.B. Saunders, 1994. 19) Hilborn S, Krahn J. Effect of time of exposure to gel-barrier tubes on results for progesterone and some other endocrine tests. *Clin Chem* 1987; 33:204. 20) Reimers TJ, et al. Effect of storage times and temperature on T3, T4, LH, prolactin, insulin, cortisol, and progesterone concentrations in blood samples from cows. *J Anim Sci* 1983; 57:683-691.

21) Smith RL. Effect of serum-separating gels on progesterone assays. *Clin Chem* 1985; 31:1239.

22) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard*. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

23) Vankrieken L. *IMMULITE reproductive hormone assays: multicenter reference range data*. Los Angeles: Diagnostic Products Corporation, 2000. Document No. ZB157-D.

24) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 25) Neinstein LS and Kaufman FR, Chapter 1: Normal Physical Growth and Development in Neinstein L.S. *Adolescent Health Care: A Practical Guide*, 4<sup>th</sup> ed. 26) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8. 27) Clinical and Laboratory Standards Institute (formerly NCCLS). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 28) Federal

Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

## Technical Assistance

Contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

---

## Tables and Graphs

---

### Precision (ng/mL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	0.46	0.08	17.4%	0.10	21.7%
2	1.36	0.12	8.8%	0.17	12.5%
3	2.54	0.26	10.2%	0.27	10.6%
4	3.92	0.38	9.7%	0.41	10.5%
5	11.6	0.92	7.9%	1.23	10.6%
6	14.8	1.04	7.0%	1.49	10.1%
7	21.3	1.49	7.0%	2.03	9.5%

**Linearity (ng/mL)**

	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8 <sup>5</sup>	4.56	—	—
	4 in 8	2.35	2.28	103%
	2 in 8	1.28	1.14	112%
	1 in 8	0.70	0.57	123%
2	8 in 8	9.74	—	—
	4 in 8	4.86	4.87	100%
	2 in 8	2.54	2.44	104%
	1 in 8	1.36	1.22	111%
3	8 in 8	11.5	—	—
	4 in 8	6.48	5.75	113%
	2 in 8	3.38	2.88	117%
	1 in 8	1.71	1.44	119%
4	8 in 8	17.3	—	—
	4 in 8	9.78	8.66	113%
	2 in 8	4.60	4.33	106%
	1 in 8	2.43	2.16	113%
5	8 in 8	19.3	—	—
	4 in 8	11.3	9.65	117%
	2 in 8	5.83	4.83	121%
	1 in 8	3.24	2.41	134%
6	8 in 8	30.6	—	—
	4 in 8	13.2	15.3	86%
	2 in 8	7.02	7.65	92%
	1 in 8	3.43	3.83	90%

**Recovery (ng/mL)**

	Solution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	—	0.26	—	—
	A	1.09	1.50	73%
	B	3.86	5.25	74%
	C	8.53	10.3	83%
2	—	0.50	—	—
	A	1.66	1.73	96%
	B	4.65	5.48	85%
	C	9.96	10.5	95%
3	—	2.30	—	—
	A	3.13	3.44	91%
	B	5.79	7.19	81%
	C	10.8	12.2	89%
4	—	4.20	—	—
	A	4.86	5.24	93%
	B	8.38	8.99	93%
	C	12.5	14.0	89%
5	—	8.21	—	—
	A	8.75	9.05	97%
	B	12.3	12.8	96 <sup>^</sup>
	C	15.8	17.8	89%
6	—	12.7	—	—
	A	13.0	13.3	98%
	B	16.5	17.1	96%
	C	17.9	22.1	81%
7	—	14.8	—	—
	A	15.5	15.3	101%
	B	21.1	19.0	111%
	C	24.4	24.0	102%

## Specificity




Compound <sup>1</sup>	Amount Added <sup>2</sup> (ng/mL)	Apparent <sup>3</sup> ng/mL	% Cross reactivity <sup>4</sup>
Androstenedione	3000	2.28	0.076%
Corticosterone	3000	12.5	0.417%
Cortisol	50,000	1.7	0.003%
Danazol	10,000	ND	ND
DHEA sulfate*	5000	0.31	0.01%
11-Deoxy-corticosterone	250	4.55	1.82%
11-Deoxycortisol	250	ND	ND
Estradiol	10,000	ND	ND
17 $\alpha$ -Hydroxy-progesterone	500	2.22	0.444%
Medroxy-progesterone	10,000	2.9	0.029%
Pregnenolone	1000	0.47	0.047%
Testosterone	1000	1.19	0.119%

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.