# Total Testosterone

For use on IMMULITE® 2000 systems

**SIEMENS**

# IMMULITE® 2000 Total Testosterone

---

## English

---

**Intended Use:** For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of total testosterone in serum and heparinized plasma, as an aid in the diagnosis and management of conditions involving excess or deficiency of this androgen.

Catalog Number: **L2KTW2** (200 tests), **L2KTW6** (600 tests)

Test Code: **TES** Color: **Dark Green**

### Summary and Explanation

In the male, testosterone is mainly synthesized in the interstitial Leydig cells of the testis, and is regulated by the interstitial cell stimulating hormone (ICSH), or luteinizing hormone (LH) of the anterior pituitary (the female equivalent of ICSH). Testosterone is responsible for the development of secondary sex characteristics, such as the accessory sex organs, the prostate, seminal vesicles and the growth of facial, pubic and axillary hair. Testosterone measurements have been very helpful in evaluating hypogonadal states. Increased testosterone levels in males can be found in complete androgen resistance (testicular feminization).<sup>9</sup>

Common causes of decreased testosterone levels in males include: hypogonadism, orchidectomy, estrogen therapy, Klinefelter's syndrome, hypopituitarism, and hepatic cirrhosis.

In the female, testosterone levels are normally found to be much lower than those encountered in the healthy male. Testosterone in the female comes from three sources. It is secreted in small quantities by both the adrenal glands and the ovaries, and in healthy women 50–60% of the daily testosterone production arises from peripheral metabolism of prehormones, chiefly androstenedione. Common causes of increased serum testosterone levels in females include polycystic ovaries (Stein-Leventhal syndrome), ovarian tumors, adrenal tumors and adrenal

hyperplasia. Virilization in women is associated with the administration of androgens and endogenous overproduction of testosterone. There appears to be a correlation between serum testosterone levels and the degree of virilization in women, although approximately 25% of women with varying degrees of virilism have serum testosterone levels that fall within the female reference range.

### Principle of the Procedure

IMMULITE 2000 Total Testosterone is a solid-phase, competitive chemiluminescent enzyme immunoassay.

**Incubation Cycles:** 1 × 60 minutes

**Time to First Result:** 65 minutes

### Specimen Collection

**EDTA tubes** should not be used in the IMMULITE 2000 Total Testosterone procedure.

Lipemic or grossly contaminated samples may give erroneous results. The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence, the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Total Testosterone has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 20  $\mu$ L serum or plasma

**Storage:** 7 days at 2–8°C, or 2 months at –20°C.<sup>8</sup>

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>10-12</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

## Total Testosterone Bead Pack (L2TW12)

With barcode. 200 beads, coated with polyclonal rabbit anti-testosterone antibody. Stable at 2–8°C until expiration date.

**L2KTW2:** 1 pack **L2KTW6:** 3 packs

## Total Testosterone Reagent Wedge (L2TWA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to testosterone in buffer. Stable at 2–8°C until expiration date.

**L2KTW2:** 1 wedge **L2KTW6:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

## Total Testosterone Adjustors (LTWL, LTWH)

Two vials (Low and High), 4.0 mL each, of testosterone in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2KTW2:** 1 set **L2KTW6:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

### Multi-Diluent 1 (L2M1Z)

For the on-board dilution of high patient samples. One vial of concentrated (ready-to-use) processed, normal human serum, with preservative, containing undetectable to low levels of testosterone. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

**L2M1Z:** 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

**L2M1Z:** 3 labels

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**L2ZT:** 250 Sample Diluent Test Tubes (16 × 100 mm)

**L2ZC:** 250 Sample Diluent Tube Caps

Also required

Distilled or deionized water; test tubes; controls.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

**Recommended Adjustment Interval:**  
1 week

**Quality Control Samples:** Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of total testosterone.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

## Expected Values

To determine testosterone reference ranges for men and nonpregnant women, a total of 186 serum samples were collected from apparently healthy adult males and females. All samples were assayed by the IMMULITE 2000 Total Testosterone procedure, with the results tabulated below.

	Median N	Median (ng/dL)	Central 95% Range	Absolute Range
<b>Females</b>				
Ovulating	40	25	ND–73	ND–73
Post-menopausal	28	26	ND–43	ND–43
<b>Males</b>				
20–49 yrs	95	323	160–726	72–853
≥ 50 years	23	317	129–767	129–767

ND: not detectable

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitation

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

### Conversion Factor:

ng/dL × 0.03467 → nmol/L

**Calibration Range:** 20–1600 ng/dL (0.7–55 nmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

**Analytical Sensitivity:** 15 ng/dL (0.5 nmol/L)

**Precision:** Samples were processed over the course of 5 days, two runs per day, for a total of 10 runs and 80 replicates. (See "Precision" table.)

**Linearity:** Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

**Recovery:** Samples spiked 1 to 19 with three testosterone solutions (3000, 6000 and 12,000 ng/dL) were assayed. (See "Recovery" table for representative data.)

**Specificity:** The antibody is highly specific for testosterone. (See "Specificity" table.)

**Bilirubin:** May cause elevation of values. (See "Bilirubin" table.)

**Hemolysis:** Presence of hemoglobin in concentrations up to 270 mg/dL has no apparent effect on results, within the precision of the assay. (See "Hemolysis" table.)

**Lipemia:** Presence of triglycerides in concentrations above 500 mg/dL may interfere with the assay. (See "Lipemia" table.)

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from 17 volunteers into plain plastic serum tubes, plastic heparinized and EDTA tubes, and plastic gel barrier tubes (SST®). All tubes were from Becton Dickinson. All samples were assayed by the IMMULITE 2000 Total Testosterone procedure.

(EDTA Plastic) = 1.78 (Serum Plastic) + 47 ng/dL  
r = 0.981

(Na Heparin Plastic) = 0.89 (Serum Plastic) + 2.3 ng/dL  
r = 0.977

(SST Plastic) = 0.96 (Serum Plastic) + 4.0 ng/dL  
r = 0.976

Means:

296 ng/dL (Serum Plastic)  
266 ng/dL (Na Heparin Plastic)  
575 ng/dL (EDTA Plastic)  
288 ng/dL (SST Plastic)

EDTA plasma is not suitable for use.  
Sodium heparin causes a depression in values.

**Method Comparison:** The IMMULITE 2000 Total Testosterone assay was compared to Coat-A-Count® Total Testosterone on 197 samples.  
(Concentration range: approximately

20 to 1200 ng/dL. See graph.) By linear regression:

(IML 2000) = 1.00 (CAC) + 13.3 ng/dL  
r = 0.976

Means:

326 ng/dL (IMMULITE 2000)  
312 ng/dL (Coat-A-Count)

## References

- 1) Abraham GE, editor. Handbook of radioimmunoassay. Marcel Dekker, 1977.
- 2) DeGroot LJ, editor. Endocrinology. Grune and Stratton, 1980.
- 3) Jaffe BM, Behrman NR, editors. Methods of hormone radioimmunoassay. Academic Press, 1974.
- 4) Yen SSC, Jaffe RB, editors. Reproductive endocrinology. W.B. Saunders, 1978.
- 5) Testosterone. Slide/Seminar Program – Bioeducational Publications, 1980.
- 6) Testosterone. Mini-Seminar. Diagnostic Products Corporation, 1982;5(2).
- 7) Ismail AAA, et al. The role of testosterone measurements in the investigation of androgen disorders. Ann Clin Biochem 1986;23:113-34.
- 8) Demetriou JA. Testosterone. In: Pesce AJ, Kaplan LA, editors. Methods in clinical chemistry. St. Louis: The C. V. Mosby Company, 1987: 268.
- 9) Wilson JD, Foster DW, editors. Williams textbook of endocrinology. Philadelphia: Saunders, 1992: 923-6.
- 10) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 11) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 12) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

## Technical Assistance

In the United States, contact  
Siemens Healthcare Diagnostics  
Technical Services department.  
Tel: 877.229.3711. Outside the United States, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

# Tables and Graphs

## Precision (ng/dL)

	Mean <sup>3</sup>	Within-Run <sup>1</sup>		Total <sup>2</sup>	
		SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	27.1	4.41	16.3%	6.58	24.3%
2	86.1	10.1	11.7%	11.2	13.0%
3	152	15.2	10.0%	15.6	10.3%
4	280	23.3	8.3%	25.6	9.1%
5	414	29.8	7.2%	33.8	8.2%
6	991	50.2	5.1%	71.5	7.2%

## Linearity (ng/dL)


	Dilution <sup>1</sup>	Observed <sup>2</sup>	Expected <sup>3</sup>	%O/E <sup>4</sup>
1	8 in 8	449	—	—
	4 in 8	230	225	102%
	2 in 8	120	112	107%
	1 in 8	63.8	56.2	114%
2	8 in 8	533	—	—
	4 in 8	302	266	114%
	2 in 8	143	133	108%
	1 in 8	75.6	66.6	114%
3	8 in 8	636	—	—
	4 in 8	304	318	96%
	2 in 8	156	159	98%
	1 in 8	80.2	79.5	101%
4	8 in 8	654	—	—
	4 in 8	319	327	98%
	2 in 8	171	164	104%
	1 in 8	81.9	81.8	100%
5	8 in 8	1137	—	—
	4 in 8	569	569	100%
	2 in 8	291	284	102%
	1 in 8	152	142	107%

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.