**CALIBRADOR BIOQUÍMICA
(Humano)**

COD 18044

USO PREVISTO

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COD 18044

INTENDED USE

The Biochemistry Calibrator is intended for calibration of the BioSystems measurement methods specified in the value sheets. For in vitro diagnostic professional use in the clinical laboratory only.

CONTENTS AND COMPOSITION

Biochemistry Calibrator. For 5 x 5 mL. Lyophilized human serum containing component concentrations suitable for the calibration of measurement procedures.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

PREPARATION AND USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 5.00 mL of distilled water into the vial. The constituent values depend on the accuracy of this reconstitution step.
3. Close vial with the stopper and let stand for 20 minutes at room temperature. For alkaline phosphatase measurement, it is recommended to let the reconstituted material stand for 1 hour at room temperature before use.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. If the material is to be used for analysis of trace elements, avoid contact of the reconstituted material with the stopper to prevent a possible contamination.
6. The reconstituted control serum is to be treated like the patient serum.

STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Calibrator serum is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C and 30 days at -20°C (when frozen once), excepting:

- AST and alkaline phosphatase are stable 5 hours at 2-8°C and 30 days at -20°C.
- Bilirubin (total) and bilirubin (direct) are stable 4 days at 2-8°C and 30 days at -20°C.
- NEFAs are stable 2 days at 2-8°C and 30 days at -20°C.

CK and bilirubin are sensitive to light. Store the vials protected from light.

Discard the vial if there are signs of microbial contamination or excessive turbidity in the reconstituted product.

WARNINGS AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

SPECIFIED VALUES

The assigned concentration values for components are batch-dependent and their traceability is shown in the enclosed value sheets. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

NOTICE TO USER

The Summary of Safety and Performance is available from EUDAMED database: <https://ec.europa.eu/tools/eudamed>.

Language versions for IFU and SDS are available in the BioSystems website <https://einfo.bio>.

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