REF

12 mL × 4: 105-012298-00
48 mL × 1: 105-012297-00

Instructions for Use

Hematology Analyzer Reagent (M-6FD DYE)

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Find IFUs in more languages on our website:
<https://www.mindray.com/en/resources-center/ivd-edocs>

COMPANY CONTACT / KONTAKT MIT DEM UNTERNEHMEN/CONTACTO CON LA EMPRESA / COORDONNEES DE LA SOCIETE/INFORMAZIONI RELATIVE ALLA SOCIETÀ / ŞİRKET İLETİŞİM / CONTATO DA EMPRESA / КОНТАКТНАЯ ИНФОРМАЦИЯ

	Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
	Mindray Building, Keji 12 th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China
	www.mindray.com
	service@mindray.com
Tel /Tél/Тел:	+86 755 81888998
Fax/Télécopie/ Faks/Факс:	+86 755 26582680

	Shanghai International Holding Corp. GmbH(Europe)
	Eiffestraße 80, Hamburg 20537, Germany
Tel /Tél/Тел:	0049-40-2513175
Fax/Télécopie/ Faks/Факс:	0049-40-255726

Revision bar indicates update to previous version / Das Revisionsbarsymbol weist auf eine Aktualisierung der vorigen Version / La barra de revisión indica una actualización a la versión anterior / La barre de révision indique la mise à jour vers la version précédente / La barra di revisione indica l'aggiornamento alla versione precedente / Değişiklik çubuğu, önceki sürüme yapılan güncellemeyi gösterir / A barra de revisão indica a atualização para a versão anterior / Строка обновления указывает на обновление предыдущей версии

M-6FD DYE

PRODUCT

M-6FD DYE

MODEL

M-6FD

SPECIFICATIONS

48mLx1,12mLx4

INTENDED USE

The M-6FD DYE participates in WBC differentiation in the DIFF channel together with M-6LD LYSE.

PRINCIPLE

The M-6FD DYE is formulated to dye the white blood cells with the assist of M-6LD LYSE. It participates in the measurement of WBC-related parameters in the DIFF channel which is accomplished by using SF Cube cell analysis technology (3D analysis using information from scatter of laser light at two angles and fluorescence signals).

ACTIVE INGREDIENTS

Fluorochrome..... ≤0.01%
Ethylene Glycol.....94-99 %

STORAGE AND EXPIRATION DATE

This product shall be stored in ventilated room free from corrosive gas and light at 2 °C~30 °C with the humidity lower than 90%. The shelf life is 1 year. The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

APPLICABLE INSTRUMENT

This product applies to BC-6000 / BC-6000Plus / BC-6100 / BC-6100Plus / BC-6200 / BC-6600Plus / BC-6700Plus / BC-6800Plus / BC-7500[B] CRP / BC-7500[N] CRP / BC-7500[R] CRP / BC-7500[NR] CRP / BC-7300[B] CRP / BC-7300[N] CRP / BC-7300[R] CRP / BC-7300[NR] CRP / BC-7500[B] CS / BC-7500[N] CS / BC-7500[R] CS / BC-7500[NR] CS / BC-700[B] / BC-700[R] / BC-720[R] / BC-760[B] / BC-760[R] / BC-780[R] / BC-700[B] CS / BC-700[R] CS / BC-760[B] CS / BC-760[R] CS / BC-780[B] CS / BC-780[R] CS / BC-770[B] CS
EN

/ BC-770[R] CS / BC-7600 / BC-7800 / BC-7900 Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

SAMPLE REQUIREMENT

Anticoagulated fresh human whole blood samples, body fluid samples. Do not use contaminated samples.

MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

INSTRUCTIONS FOR USE

1. Restore the product to usage temperature;
2. Open the external packing and insert the pickup tube into the reagent pack based on the matching between the reagent pack label and the inlet label of the analyzer cap assembly connector;
3. Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
4. Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis See the operator's manual of the analyzer for details.

CUT-OFF VALUE/REFERENCE INTERVALS

Not applicable.

RESULT ELABORATION

Not applicable.

LIMITATIONS

Not applicable.

PRODUCT SPECIFICATIONS

1. Appearance: transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results of the reagent tested on the Mindray Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1. Blank Count Requirements

Parameter	Blank Count Requirements
WBC-D	$\leq 0.2 \times 10^9/L$

PRECAUTIONS

1. For professional use in vitro diagnosis only.
2. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
3. Do not use the reagent if it is frozen.
4. If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
5. If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
6. Dispose of the waste, residual and contaminated packing based on local regulations.
7. The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
8. Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
9. Disposal of waste liquid and materials should be in accordance with local guidelines.
10. The Material Safety Data Sheet (SDS) is available upon request.
11. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
12. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
13. This product contains components classified as follows in accordance

with the Regulation (EC) No. 1272/2008:

 Warning	
Statement(s)	
H302	Causes severe skin burns and eye damage.
Prevention	
P264	Wash all exposed external body areas thoroughly after handling.
P270	Do not eat, drink or smoke when using this product.
Response	
P301+P312	IF SWALLOWED: Call a POISON CENTER/doctor/physician/first aider/if you feel unwell.
P330	Rinse mouth.
Disposal	
P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.

GRAPHICAL SYMBOLS

		
Batch code	Use-by date	Temperature limit

Humidity Limitation

Manufacturer

In vitro diagnostic
medical device

Catalogue Number

Authorized
representative
in the European
Community

Consult
instructions for use

European Conformity

Keep away from
sunlight.
Indicates a medical
device that needs
protection from light
sources.

Unique Device
Identifier

REFERENCES

Not applicable.

Approval Date of the Operator's Manual

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